Shellee Kim Lamb March 9, 2017

Gloria Calero
Regulatory Affairs Specialist
Cardiac Rhythm Heart Failure

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K.v.K. Zuid Limburg : 090.30.178

V.A.T. : NL 0013.07.915.8.01

**Dimensions
(mm)
[Height ×
Length ×
Thickness]**

44.7 × 47.9 × 7.5

Mass (g)

27.1

Volume

12.1 cm³

Surface area

30.7 cm²

Connector

IS-1 dual chamber

Case material

Titanium

Connector materials

Polyurethane, silicone rubber

Power source

Lithium iodine

**Average available
capacity (Ah)**

1.18

Initial voltage (V)

2.8

Batteries

Battery and power consumption

Parameter	Capability	Notes
Model	Sigma 263	
Type	Lithium-Iodine	
Manufacturer	Medtronic Energy and Component Center	
Voltage	2.8 V	
Usable capacity	1.3 Ah	Capacity available from BOS to EOS
Minimum capacity at RRT/ERI	0.08 Ah	Minimum capacity is measured at nominal settings, 50 Ω load in unipo- lar or bipolar polarity, 100% pacing
Device input impedance	> 150 K Ω	Device input impedance is measured to CENELEC standard (EN 45502-2-1)

Recommended Replacement Time (RRT/ERI)

Vitatron Holding B.V.

Mailing address : P.O. Box 1220

6201 MP Maastricht

Visiting Address : Endepolsdomein 5

6229 GW Maastricht

Tel : +31 43 356 6551

Fax : +31 43 356 6542

K.v.K. Zuid Limburg : 090.30.178

V.A.T. : NL 0013.07.915.B.01

The Recommended Replacement Time (RRT/ERI) warns when the device battery is nearing depletion. When the battery voltage has dropped below a defined limit, the device sets the RRT/ERI status and reverts to VVI operation at a rate of 65 min⁻¹. At most programmed settings, the majority of devices will function for a minimum of three months after RRT/ERI. The devices indicators of RRT/ERI status are listed in the following table:

Operation	Indicator
Nonmagnet mode	VVI mode at 65 min ⁻¹
Magnet mode	VOO mode at 65 min ⁻¹
Telemetry	Replacement message on the programmer
Battery/Lead information	Replacement message and displayed battery voltage on programmer

Projected service life Vitatron A60DR

Pacing	A Amplitude, V Amplitude	Rate, Pulse Width	Lead impedance	
			500 Ω	1000 Ω
			Longevity (years)	
DDDR or DDD, 0%	1.5 V, 2.0 V ^a	60 min ⁻¹ , 0.4 ms	12.8	12.8
	2.5 V, 2.5 V		11.9	11.9
	3.5 V, 3.5 V		12.5	12.5
DDDR or DDD, 50%	1.5 V, 2.0 V ^a	60 min ⁻¹ , 0.4 ms	11.4	12.0
	2.5 V, 2.5 V		10.2	11.0
	3.5 V, 3.5 V		9.1	10.4

Vitatron Holding B.V.

Mailing address : P.O. Box 1220
6201 MP Maastricht
Visiting Address : Endepolsdomein 5
6229 GW Maastricht

Tel : +31 43 356 6551
Fax : +31 43 356 6542
K.v.K. Zuid Limburg : 090.30.178
V.A.T. : NL 0013.07.915.B.01

Pacing	A Amplitude, V Amplitude	Rate, Pulse Width	Lead impedance	
			500 Ω	1000 Ω
			Longevity (years)	
DDDR or DDD, 100%	1.5 V, 2.0 V*	60 min ⁻¹ , 0.4 ms	10.2	11.3
	2.5 V, 2.5 V		9.0	10.2
	3.5 V, 3.5 V		7.1	8.9
DDDR or DDD, 0%	2.5 V, 2.5 V	70 min ⁻¹ , 0.5 ms	11.8	----
	5.0 V, 5.0 V		11.5	----
DDDR or DDD, 100%	2.5 V, 2.5 V	70 min ⁻¹ , 0.5 ms	7.9	----
	5.0 V, 5.0 V		3.8	----
DDDR or DDD, 100%	5.0 V, 5.0 V	70 min ⁻¹ , 1.0 ms	2.5	----
DDDR or DDD, 100%	5.0 V, 5.0 V	100 min ⁻¹ , 1.0 ms	1.8	----

*For Atrial Output Management, the Minimum Adapted Amplitude is 1.5 V (nominal). For Ventricular Output Management, the Minimum Adapted Amplitude is 2.0 V (nominal)

Radiopaque dimensions:

Vitatron Holding B.V.

Mailing address : P.O. Box 1220

6201 MP Maastricht

Visiting Address : Endepolsdomein 5


6229 GW Maastricht

Tel : +31 43 356 6551

Fax : +31 43 356 6542

K.v.K. Zuid Limburg : 090.30.178

V.A.T. : NL 0013.07.915.8.01

Connector

Warning: Verify lead and connector compatibility before using a lead with this device. Using an incompatible lead may damage the connector, result in electrical current leakage, or result in an intermittent electrical connection.

Note: If you are using a lead that requires an adaptor for this device, please contact your Vitatron representative for information about compatible lead adaptors. Select a compatible lead. Refer to Table 1.

Table 1. Lead and connector compatibility

Vitatron Holding B.V.

Mailing address : P.O. Box 1220

6201 MP Maastricht

Visiting Address : Endepolsdomein 5
6229 GW Maastricht

Tel : +31 43 356 6551

Fax : +31 43 356 6542

K.v.K. Zuid Limburg : 090.30.178

V.A.T. : NL 0013.07.915.8.01

Connector	Polarity	Compatibility
Ventricular	Bipolar/Unipolar	Compatible with the IS-1 standard for unipolar and bipolar pacing and sensing ^a
Atrial	Bipolar	Compatible with the mechanical requirements of the IS-1 standard, but is only designed for sensing ^a

a) IS-1 refers to the International Connector Standard (see Document No. ISO 5841-3) whereby pulse generators and leads so designated are assured of meeting the electrical and mechanical parameters specified in the IS-1 International Standard.

Connector block geometry

model	Lead tip bore diameter	Overall bore depth
A60DR	1.93 mm	25.2 mm

All requirements for connector sensing/pacing configuration (bipolar, unipolar, other); connector geometry (bore depths and diameters in millimeters), or reference to published connector standards, an explanation of any markings used to identify the connector on the IPG were met in accordance with ISO 5841-3:2000

Requirement ISO 5841	pass/ not pass
1. Connection lead:	
A) location and number of insulating rings	pass
B) size	pass
B) continuity of the electrical circuit	continuous - pass
2. strength of the introduction / extraction flow gauge	<14H - pass
3. electrical impedance between the conductive parts	≥50 kOhm - pass
4. Deformation and its influence on the connection terminal on	<14H - pass
5. Impact of bipolar clamping screw connector sockets on Unipolar connection terminal	pass
6. Socket outlet	
A) size	pass
B) max force introduction	≤9 H - pass

Vitatron Holding B.V.

Mailing address : P.O. Box 1220

6201 MP Maastricht

Visiting Address : Endepolsdomein 5

6229 GW Maastricht

Tel : +31 43 356 6551

Fax : +31 43 356 6542

K.v.K. Zuid Limburg : 090.30.178

V.A.T. : NL 0013.07.915.8.01

7. Labeling sign («IS-1» clear markin)

pass

Vitatron Leads for use with Vitatron A-Series IPGs

Device Name	Device Model Number(s)
Impulse II	IHP09B, IHP09JB
Crystalline	ICL08, ICL08B, ICL08JB, ICM09, ICM09B, ICM09JB
Crystalline ActFix	ICF09B, ICQ09B
Excellence PS+	IMK49B, IMK49JB
Brilliant S+	IMW14Q, IMW15Q, IMW16Q, IMW17Q, IMW18Q

Wrench –Torque



included in the sterile package. The wrench is designed to prevent damage to the device from overtightening a setscrew.

Materials (not in contact with the human body):

1)Haft - Polyetherimide (G.E. Ultem 1100F-8104)

2)Hex Shaft - 455 Corrosion resistant steel

Dimensional (inches):

All A...

Mass: 8g

The torque is limited to 12 to17 in-oz. during first clockwise revolution; it will click when torque value in this range is achieved.

Vitatron Holding B.V.

Mailing address : P.O. Box 1220

6201 MP Maastricht

Visiting Address : Endepolsdomein 5

6229 GW Maastricht

Tel : +31 43 356 6551

Fax : +31 43 356 6542

K.v.K. Zuid Limburg : 090.30.178

V.A.T. : NL 0013.07.915.8.01

5. Diagnostic and emergency parameters

Automatic and clinician-selectable diagnostics:

The automatic diagnostics listed in Table are automatically collected by the pacemaker. Note that some of the diagnostics have programmable parameters that are accessed by selecting Data Collection Setup under the Data icon.

Automatic diagnostics

Parameter	Settings
Heart Rate Histograms (Short and Long Term, Atrial and Ventricular)	
Include Refractory Senses	Include; Exclude
AV Conduction Histograms (Short and Long Term)	
Reduced VP+ Histogram	
Atrial High Rate Episodes (Mode Switch On)	
Collection Delay after Mode Switch	0; 1; 2 ... 20; 25; 30 ... 60 s
Collection Method ^b	Frozen; Rolling
Atrial High Rate Episodes (Mode Switch Off)	
Detection Rate	80; 85; 90 ... 180; 200; 220; 240 ... 320; 330; 350; 370; 380; 400 min ⁻¹
Detection Duration	1; 2; 3 ... 20; 25; 30 ... 50; 55; 60 s
Termination Beats	5; 6; 7 ... 20 beats
Collection Method	Frozen; Rolling
Ventricular High Rate Episodes	
Detection Rate	80; 85; 90 ... 180; 200; 220; 240 ... 320; 330; 350; 370; 380; 400 min ⁻¹

Vitatron Holding B.V.

Mailing address : P.O. Box 1220

6201 MP Maastricht

Visiting Address : Endepolsdomein 5

6229 GW Maastricht

Tel : +31 43 356 6551

Fax : +31 43 356 6542

K.v.K. Zuid Limburg : 090.30.178

V.A.T. : NL 0013.07.915.B.01

Detection Duration	2; 3; 4 ... 198; 199; 200 beats
Termination Beats	5; 6; 7 ... 20 beats
SVT Filter	Off; On
Collection Method	Frozen; Rolling
AT/AF Trend	
AT/AF Durations	
Ventricular Rate During AT/AF	
Chronic Lead Trends	
Lead Monitor Counters	
Output management trend	Based on use of Output Management
Ventricular Output Management	
Key Parameter History	

^a Heart Rate Histograms can be programmed to include or exclude refractory sensed events.

^b Collection Method applies to Atrial High Rate Episodes and Ventricular High Rate Episodes.

Common clinician-selectable diagnostics

Diagnostic and parameters	Parameter settings
Custom Rate Trend	
Duration	Beat-to-Beat; 1 Hour; 24 hours
Collection Method	Frozen; Rolling
Include Refractory Senses	Include; Exclude
Ventricular Output Management Detail	
EGM Collection	Off; Atrial EGM; Ventricular EGM; Summed EGM ^a

a- For the G20 SR, E10 S, A20 SR, and A10 S models, EGM Collection options are limited to OFF and EGM.

High Rate Detail for A-series pacemakers ^a

High Rate Detaila	
High Rate Typec	AHR; VHR
EGM Type	Off; Atrial EGM; Ventricular EGM; Summed

Vitatron Holding B.V.

Mailing address : P.O. Box 1220

6201 MP Maastricht

Visiting Address : Endepolsdomein 5

6229 GW Maastricht

Tel : +31 43 356 6551

Fax : +31 43 356 6542

K.v.K. Zuid Limburg : 090.30.178

V.A.T. : NL 0013.07.915.8.01

	EGM ^d
Allocation (Collection Methode ^b = Frozen, High Rate Type = AHR only or VHR only)	1 for 0/16; 1 for 16/0; 1 for 8/8; 2 for 0/8; 2 for 8/0; 2 for 4/4; 4 for 0/4; 4 for 4/0; 4 for 2/2 (number of episodes for pre-onset seconds/ post onset seconds collected)
Allocation (Collection Method ^b = Rolling, High Rate Type = AHR only or VHR only)	1 for 8/0; 1 for 4/4; 1 for 0/16; 2 for 4/0; 2 for 2/2; 2 for 0/8; 4 for 2/0; 4 for 1/1; 4 for 0/4 (number of episodes for pre-onset seconds/ post onset seconds collected)
Pre-detection Timeout	1; 2; 3 ... 12; 14; 16 ... 24 weeks

^a High rate detection rate, detection duration, and termination criteria are set by parameters for the automatic diagnostic.

^b Collection Method is set in the High Rate automatic diagnostic.

Emergency pacing:
Emergency pacing provides VVI pacing at high output settings in emergency situations for pacemaker-dependent patients. Table lists the emergency settings. The settings apply to all pacemakers.

Parameter	Setting
Mode	VVI
Pacing Rate	70 min ⁻¹

Vitatron Holding B.V.

Mailing address : P.O. Box 1220

6201 MP Maastricht

Visiting Address : Endepolsdomein 5

6229 GW Maastricht

Tel : +31 43 356 6551

Fax : +31 43 356 6542


K.v.K. Zuid Limburg : 090.30.178

V.A.T. : NL 0013.07.915.B.01

	Ventricular Amplitude 7.5 V Pulse Width 1.5 ms Sensitivity 2.8 mV Pace Polarity Unipolar Sense Polarity Unipolar Lead Monitor Monitor Only Ventricular Refractory Period 330 ms Single Chamber Hysteresis Off Ventricular Output Management Off	
6. Materials	Material (longterm contact with internal environment of the body)	Device Component
	Titanium, grade 1	Left-hand shield
	Fe to 0,12	Right-hand shield
	C to 0,05	
	Si to 0,08	
	N to 0,04	
	Ti 99,5 - 99,99	
	O to 0,1	
	H to 0,008	
	Silicone Rubber Adhesive MED 2000	Adhesive
	Polyurethane, 75 durometer	Connector module
	Elasthane	
	ETR silicone rubber, 35 durometer	Split grommet
7. Labeling	IPG marking	

Vitatron Holding B.V.
Mailing address : P.O. Box 1220
6201 MP Maastricht
Visiting Address : Endepolsdomein 5
6229 GW Maastricht

Tel : +31 43 356 6551
Fax : +31 43 356 6542
K.v.K. Zuid Limburg : 090.30.178
V.A.T. : NL 0013.07.915.B.01

Shelf box label

Vitatron Holding B.V.

Mailing address : P.O. Box 1220

6201 MP Maastricht

Visiting Address : Endepolsdomein 5

6229 GW Maastricht

Tel : +31 43 356 6551

Fax : +31 43 356 6542


K.v.K. Zuid Limburg : 090.30.178

V.A.T. : NL 0013.07.915.B.01

A60 DR

vitatron

Reduced VP™ +

IPG device
Model A60A1
NBG DDR
Adaptive

SN 350001196G

2017-11-28
YYYY-MM-DD

2016-06
YYYY-MM

CE0344

STERILE EO
-18 °C / 0 °F
55 °C / 131 °F

DDDR

130 min¹
60 min¹
Upper tracking rate / lower rate

150 ms
120 ms
A-V interval (paced / sensed)

230 ms
Ventricular refractory period

0.5 mV
Atrial sensitivity

2.8 mV
Ventricular sensitivity

3.5 V
0.4 ms
Atrial amplitude / pulse width

2.5 V
0.4 ms
RV amplitude / pulse width

AUTO
Sensing polarity

AUTO
Pacing polarity

80 min¹
Lower rate

AUTO
PWP

Vitatron Holding B.V.
Endepolsdomein 5
6229 GW Maastricht
The Netherlands

Manufactured in
Singapore
Product of Singapore

NR54387A001 A - NR54377A001 A

Vitatron Holding B.V.

Mailing address : P.O. Box 1220

6201 MP Maastricht

Visiting Address : Endepolsdomein 5


6229 GW Maastricht

Tel : +31 43 356 6551

Fax : +31 43 356 6542

K.v.K. Zuid Limburg : 090.30.178

V.A.T. : NL 0013.07.915.8.01

Sterile package label

Vitatron Holding B.V.
Mailing address : P.O. Box 1220
6201 MP Maastricht
Visiting Address : Endepolsdomein 5
6229 GW Maastricht

Tel : +31 43 356 6551
Fax : +31 43 356 6542
K.v.K. Zuid Limburg : 090.30.178
V.A.T. : NL 0013.07.915.B.01

Packaging description:
Table 1: Primary Packaging – Sterile Package

Package part	Material
Inner Tray	Coextruded copolyester material consisting of a cap/core/cap construction. The core is virgin polyethylene terephthalate glycol (PETG) and the cap layers are virgin PETG with an antiblocking additive.
Outer Tray	Thermoformed polyethylene terephthalate glycol (PETG)
Inner Lid	Spunbonded polyolefin
Outer Lid	Dot-patterned hot melt adhesive coating

Table 2: Secondary Packaging –Shelf Box

Package part	Material
Shelf Box	White vat non-bending chip Box wrap – white litho Adhesive – Davis NW300C or Wisdom 6567, 6578 or equivalent Shrink-wrap of a polyethylene (PE) shrink film

Sticker for final package will be added locally in Russia after approval with information below (information can be changed in accordance with Russian legislation):

Vitatron Holding B.V.

Mailing address : P.O. Box 1220

6201 MP Maastricht

Visiting Address : Endepolsdomein 5

6229 GW Maastricht

Tel : +31 43 356 6551

Fax : +31 43 356 6542

K.v.K. Zuid Limburg : 090.30.178

V.A.T. : NL 0013.07.915.B.01

		Product name (in Russian)	
		Manufacturer (name)	UPN (unique product number)
		Country of origin	CFN n(customer facing number)
		See IFU	
		RC #	date of issue
8. Methods and conditions of disinfection and sterilization	Sterile. Single use only. Products cannot be resterilized		
9. Operating conditions	<u>Checking and opening the package</u> – Before opening the sterile package tray, visually check for any signs of damage that might invalidate the sterility of the package contents. <u>If the package is damaged</u> – The device packaging consists of an outer tray and inner tray. Do not use the device or accessories if the outer packaging tray is wet, punctured, opened, or damaged. Return the device to Vitatron because the integrity of the sterile packaging or the device functionality may be compromised. This device is not intended to be resterilized. <u>Dropped device</u> – Do not implant the device if it has been dropped on a hard surface from a height of 30 cm or more after it is removed from its packaging. <u>Use by date</u> – Do not implant the device after the “Use by” date on the package label. Battery longevity may be reduced. <u>For single use only</u> – Do not resterilize and reimplant an explanted device that has been contaminated by contact with body fluids. <u>Sterilization</u> – Vitatron has sterilized the package contents with ethylene oxide before shipment. This device is for single use only and is not intended to be resterilized. <u>Avoid magnets</u> – To avoid damaging the device, store the device in a clean area away from magnets, kits containing magnets, and any sources of electromagnetic interference. <u>Temperature limits</u> – Store and transport the package between –18 °C and +55 °C. Electrical reset may occur at temperatures below –18 °C. Device longevity may decrease and performance may be affected at temperatures above +55 °C.		

Vitatron Holding B.V.

Mailing address : P.O. Box 1220
6201 MP Maastricht
Visiting Address : Endepolsdomein 5
6229 GW Maastricht

Tel : +31 43 356 6551
Fax : +31 43 356 6542
K.v.K. Zuid Limburg : 090.30.178
V.A.T. : NL 0013.07.915.8.01

	Device temperature – Allow the device to reach room temperature before it is programmed or implanted. Device temperature above or below room temperature may affect initial device function. Protection of the IPG from damage caused by atmospheric pressure changes: Implantable parts of the IPG withstand changes of pressure. All devices passed the requirements in accordance with ISO 14708 (for detailed information please see <u>nEw3 Mechanical DAU Design Verification Test Report</u> (section 9.6) During the procedure the product is inside of the human body and is operated at a body temperature from 32 to 42 °C, pressure 50-150 kPa; Humidity up to 100%.
10. Storage and Transportation conditions	Temperature limits – Store and transport the package between –18 °C and +55 °C (ISO 14708-1). Electrical reset may occur at temperatures below –18 °C. Device longevity may decrease and performance may be affected at temperatures above +55 °C. Device temperature – Allow the device to reach room temperature before it is programmed Basic rate, test pulse rate, pulse duration, and pulse amplitude remain within expected tolerances when the device temperature is between 20 C to 43 C. Protection of the IPG from damage caused by atmospheric pressure changes: Implantable parts of the IPG withstand changes of pressure. all devices passed the requirements in accordance with ISO 14708 Leak Rate Each device is tested in manufacturing to assure an actual leak rate less than or equal to 1x 10-8 standard cm3/sec of helium. Dust protection The IPGs and ICDs are hermetically sealed, so dust will not affect them internally. Dust protection classification Basically does not apply to IPGs. Packed by the manufacturer products are transported by all modes of transport in closed vehicles in accordance with the rules of carriage of goods, operating in each mode of transport. Do not store and transport this product above + 55 °C. Transport in a manner that protects the integrity of the package and the sterile barrier. Keep in a dry location away from heat and sunlight.
11. Software	Software Programmers and software – Use only Medtronic programmers and Vitatron application software to communicate with the device. Programmers and software from other manufacturers are not compatible with Vitatron devices. Version - VSF20 Date of issue – 2013 Safety class – C (regarding IEC 62304)

Vitatron Holding B.V.

Mailing address : P.O. Box 1220

6201 MP Maastricht

Visiting Address : Endepolsdomein 5

6229 GW Maastricht

Tel : +31 43 356 6551

Fax : +31 43 356 6542

K.v.K. Zuid Limburg : 090.30.178

V.A.T. : NL 0013.07.915.8.01

User Characteristics - The primary audience for the programmer software is intended to be clinicians and nurses for the purpose of programming and interrogating the IPG. A secondary audience are the in-house Medtronic software and hardware development engineers and verification test engineers. The environment that the IPG application software and the 2090 programmer will be used in is a follow-up setting.

Session Description

This feature provides the capability of defining a patient session and controlling activities that occur at the start of a patient session. A patient session is generally defined as the period of use of the programmer for one patient involving one implanted device during

which data collection and/or programming of the device is performed. This is termed an “active” or “programming” session.

Common to all NG IPG sessions is the need to identify the device with which the programmer is communicating either through an auto-id or manual model selection. The device may also be identified during the process in which demographics information is

provided to fill in registration form information. The information gathered during a session will be stored in a patient session record. This

data will include the patient demographic information. At the start of a NG IPG session, after the device has been identified and an initial interrogate has occurred, the state of the device (i.e. the state of the device parameters and its components relative to patient history) will be analyzed through the quick scan diagnostic procedure (TACHY). The state of the device also will be saved at this time within a diagnostics. The patient session record will become the repository of other data collected during the

patient session (i.e. test results, episodes, etc.) If the session is communicating with an implantable device, model identification of the

device will occur prior to the start of a session. The model information will be used by the session to make available only the model applicable aspects of things like protocol and nominal selection lists.

Connection type Programmer - Device:

Device uses the Telemetry (175kHz) for connection with programmer through the programming head. Range of coverage – 6cm.

Vitatron Holding B.V.

Mailing address : P.O. Box 1220

6201 MP Maastricht

Visiting Address : Endepolsdomein 5

6229 GW Maastricht

Tel : +31 43 356 6551

Fax : +31 43 356 6542

K.v.K. Zuid Limburg : 090.30.178


V.A.T. : NL 0013.07.915.8.01

Terms of Limited Warranty

- A. This Limited Warranty provides the following assurance to a patient who receives this Vitatron pulse generator:
- (1) Should the pulse generator fail to function within normal tolerances at any time due to a defect in materials or workmanship, Vitatron will issue a credit as defined in Subsection A.3. Battery depletion will occur with time and, as used in Subsection A.2., is not considered a defect in material or workmanship.
 - (2) A credit as defined in Subsection A.3. will be issued if the pulse generator is explanted within 60 months due to battery depletion, in accordance with the pulse generator's Recommended Replacement Time Indicator. Each warranty period commences on the implant date. The warranty on battery depletion is made void if the pulse generator has been programmed to output values higher than the standard factory settings for an accumulated period of longer than 3 months.
 - (3) The credit is equal to the purchase price and is issued against the purchase of any Vitatron pulse generator requested as a replacement, or, at the discretion of Vitatron, is a functionally comparable replacement Vitatron pulse generator at no charge. Normally the credit will be issued only to the hospital or physician purchasing the replacement pulse generator.
- B. To qualify for the Limited Warranty, the following conditions must be met:
- (1) The pulse generator must be implanted prior to its "Use By" date.
 - (2) The pulse generator for which credit is sought must be returned to Vitatron in a shockproof package within 30 days of its removal and shall then become the property of Vitatron. Furthermore, a "Request for analysis" form should be completed and returned with the pulse generator.
- C. This Limited Warranty is limited to its express terms. In particular:
- (1) Vitatron shall not be liable to any person for any medical expenses or any direct or consequential damages resulting from implant or removal of any pulse generator or caused by any defect, failure or malfunction of any pulse generator, whether a claim for such damage is based upon warranty, contract, tort or otherwise. The exclusive remedy for any defect, failure or malfunction of any pulse generator shall be as stated.
 - (2) This Limited Warranty is made only to the patient in whom the pulse generator was originally implanted and is expressly in lieu of any implied warranty of merchantability or fitness for a particular purpose, whether arising from statute, common law, custom or otherwise. The remedies set forth in this Limited Warranty shall be the exclusive remedies available to any person. No person has any authority to bind Vitatron to any representation, condition or warranty except as contained in this Limited Warranty.
 - (3) This Limited Warranty is only valid as to the pulse generator and is not applicable, in whole or part, to any other Vitatron product including, but not limited to, Vitatron pacing leads.

Vitatron Holding B.V.
Mailing address : P.O. Box 1220
6201 MP Maastricht
Visiting Address : Endepolsdomein 5
6229 GW Maastricht

Tel : +31 43 356 6551
Fax : +31 43 356 6542
K.v.K. Zuid Limburg : 090.30.178
V.A.T. : NL 0013.07.915.8.01

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V.A.T. : NL 0013.07.915.B.01

Tel. +31 43 356 6551

Representative company of Vitatron B.V. in Moscow Russia

123317 Moscow,

Presnenskaya Naberezhnaya , 10

Tel.: +7 495 580 7377; +7 495 580 7377 FREE

Fax: + 7 495 580 7378

E-mail: info.russia@medtronic.ru

Internet: www.medtronic.eu, www.medtronic.ru

Vitatron Holding B.V.

Mailing address : P.O. Box 1220

6201 MP Maastricht

Visiting Address : Endepolsdomein 5

6229 GW Maastricht

Tel : +31 43 356 6551

Fax : +31 43 356 6542

K.v.K. Zuid Limburg : 090.30.178

V.A.T. : NL 0013.07.915.B.01

- (4) Vitatron shall inspect the returned pulse generator and determine whether credit is due under the terms of this warranty. No credit will be issued if Vitatron finds evidence of improper handling or of material alteration to the pulse generator.
- (5) This Limited Warranty is only applicable to pulse generators sold through a sales organization authorized by Vitatron. Vitatron as used in this warranty shall mean Vitatron B.V. in The Netherlands.
- (6) The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. Users may benefit from statutory warranty rights under the legislation governing the sales of consumer goods. If any part or term of this Limited Warranty is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the Limited Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Limited Warranty did not contain the particular part or term held to be invalid.

Measure unit of the device technical parameters (mentioned in IFU/Technical manual):

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