Instruction for Use

INDICATIONS

Bone tacks are indicated for stabilization during tissue regeneration and bone repair in the maxillofacial region.

USAGE

The MC Bio OSTEOSYNT system is designed to stabilize barrier membranes onto cortical plate bone. This may be used in craniofacial, maxillofacial or mandibular bone. General patient health, bone type and quality, and functional loads exerted should be considered and carefully evaluated prior to use.

TERMS OF TRANSPORTATION

The medical devices of the OSTEOSYNT system shall be conveyed using all types of covered vehicles in accordance with the requirements and rules of cargo transportation applicable to the corresponding mode of transport.

Type of shipping : by parcel, motor vehicles, planes, containers.

Transportation conditions: temperature; -10°C to +30°C, relative humidity not more than 75%.

Product boxes should not be exposed to excessive shocks and atmospheric precipitation during handling and transportation, in accordance with handling marks bearing the respective signs "Fragile, handle with care", "Top side", "Protect from moisture".

TEMPERATURE OF USE

Room temperature .

DISPOSAL

If not used, the MC Bio Supertack devices and their packaging , are produced in recyclable materials.

If the devices are used and contaminated, for a correct disposal, they should be previously cleaned by a disinfectant solution and then sterilized by a steamclave .

Products should be disposed of in accordance with the rules and regulations approved by the Russian legislation

WARRANTY

In accordance with the Directive 1999/44/EC.

The devices are covered by a warranty for a time of 1 year.

The warranty is considered only for production mistakes that don't allow the correct functioning of the device.

The compliance of one device to the standard of production is a pertinence of MC Bio company, so only the devices that will be send to the MC Bio company , might be investigated for the application of the warranty.

The loosing of the original properties in function of the use are not covered by warranty.

The warranty include only the replacing or the repairing of the device.

For the tacks, no warranty can be considered against the possible bending during the use.

The Court of Como (Italy) is exclusively competent for any dispute .

The manufacturer guarantees the conformity of the medical device Osteosynt System, to all regulatory requirements subject to the conditions of operation, transportation and storage.

The Manufacturer's warranty applies to the product prior to the expiration of its shelf life indicated on the packaging and relates to the unopened, undamaged and correctly stored product. The guaranteed shelf life subject to observance of storage conditions is :

DESCRIPTION

The MC Bio Osteosynt system is comprised of tacks (pins) , supplied STERILE, with three different lengths: 3mm, 4mm and 5mm, and instruments for the pins installation

The different tacks are:

- Pin 3.0
- Pin 4.0
- Pin 5.0
- Pin 3.0 Extra Strong
- Pin 4.0 Extra Strong
- Pin 5.0 Extra Strong

The drawing of the tacks include a 2.5mm head and a shaft with a barb for a perfect fixing in the bone.

The instruments of the MC Bio Osteosynt system includes:

- Lingual applicator Specialist
- Applicator Tip - Long
- Applicator Tip - Short
- Stop - Hard
- Stop - Medium
- Container for pins (Box)
- Set for the installation of the Supertack pins (Set SUPERTACK) , consisting of: surgical mallet, pin applicator, tool box, pin storage container.
- Specialist set for the installation of the Supertack pins (Set SPECIALIST) , consisting of: surgical mallet, pin applicator, stops (hard and medium), tool box, pin storage container and long applicator tip.

The pin applicator and the tips have a retaining head for the 2.5mm pins; the grip of larger diameter pin heads might damage the instrument.

The mallet weight and sizes are balanced for a correct fitting of the pins; its head has a diameter of 2cm and its weight is about 140gr.

The tool boxes (Supertack and Specialist) can be easily put in the standard steam claves (18-22 l of volume); their sizes are: length - 20cm, width - 9cm, depth - 2,5cm.

All the instruments can be easily removed from the tool box, during the use, but also for an easy cleaning and disinfection .


All the instruments and boxes are supplied NON-STERILE, and can be sterilized by steam autoclaves.

MATERIALS

- All the tacks are made of TiGr5
- The instruments are made of stainless steel
- The box is made of plastics medical grade.

SURGICAL PROCEDURE

1. Remove the tacks from their STERILE packaging and place in the organizer .
2. Firmly press the applicator onto the head of the desired tack in the organizer. This will grip the tack to be carried to the placement site.
3. Seat the tack through the barrier material only with the applicator perpendicular to the bone surface or the tack may be damaged. Placement should be accomplished with firm taps on the base of the applicator with the hammer. To disengage the tack, roll the applicator to one side.
4. Repeat this process until the barrier is firmly anchored. Always store the applicator in the surgical box; otherwise it may become distorted and not able to grip and pick up the tacks..
5. Every effort should be made to assume primary closure of all areas. Suture the soft tissue in the usual manner when the surgical procedure is complete.
6. The Patient's Record should include radiographs and written documentation indicating number and position of the tacks.

HOW SUPPLIED

The tacks of the SUPERTACK system are provided STERILE . Tacks are inside a glass vial, enclosed in a transparent blister sealed by a Tyvek cover.

EXPIRY DATE

It's recommended to respect the expiry date written on the label.

STORAGE

The storage of the tacks need temperature and conditions of Relative Humidity in accordance with the European Pharmacopeia. The finished product shall be stored in closed rooms at temperatures from +10°C to +30 °C and relative humidity not more than 75%. Store the Osteosynt devices in consumer packaging and remove from consumer packaging only immediately prior to use. Preservation of the package integrity is an indispensable prerequisite. The maximum time of storage is in function of the shelf life. On each label is present the expiration date of the STERILE condition.

MAX DURATION OF IMPLANTATION

The tacks (Pin MC Bio OSTEOSYNT) should be removed upon completion of the healing process (on average about 6-10 months) . In instances where the tack removal is required, after normal surgical exposure, the tacks may be dislodged by using a periosteal elevator, scalpel or other suitable thin flat instrument. The tacks should be removed, inventoried and discarded. The surgical site is now re-sutured.

REUSE

The tacks (Pin MC Bio OSTEOSYNT) are for single use only therefore reuse of the tacks is NOT recommended.

CONTRAINDICATIONS

General patient evaluation is critical prior to any procedure. Contraindications include, but are not limited to, local or systemic infection, clotting disorders, vascular impairment, radiation, steroid or anticoagulant therapy, diabetes or other systemic or metabolic limitations which would compromise healing.

COMPLICATIONS

Complications include those associated with any osseous surgical procedure and an esthetic usage. Severity and type of complications may indicate tack removal at the discretion of the clinician

WARNING


It's the clinician's responsibility to become familiar with proper technique for use this device. The clinician should have training and experience with bone fixation techniques.

The applicator must be held perpendicular to the bone during placement to avoid bending the tack. Complete seating should be attempted with a single strike of the hammer. The tacks are not designed to endure excessive or abnormal functional forces..

SURGICAL PROCEDURE

The instrumentation are provided NON STERILE and must be sterilized prior to use. An aseptic technique is required for predictable success.

1. Place the tacks in the organizer.
2. Firmly press the applicator onto the head of the desired pin in the organizer. This will grip the tack to be carried to the placement site.
3. Seat the tack through the barrier material only with the applicator perpendicular to the bone surface or the tack may be damaged. Placement should be accomplished with firm taps on the base of the applicator with the hammer. To disengage the tack, roll the applicator to one side.
4. Repeat this process until the barrier is firmly anchored. Always store the applicator in the surgical box; otherwise it may become distorted and not able to grip and pick up the tacks..

USE

-TOOL BOX

The Supertack and Specialist tool boxes allow the sterilization of the tools by steam claves.

For the sterilization, introduce the tools in the box, close the cover, and fit the box in aclave. It's recommended the use of steam sterilizer with Dynamic Air Removal (class B Clave)

-APPLICATOR

The applicator is made of two parts, the tip and the grip. When the tip doesn't allow the correct retaining of the pin, it can be replaced only disassembling it by the thread that connect the two parts.

To use correctly the applicator, firmly press it onto the head of the pin. In this way the pin can be carried the placement site. The pin will be fitted in the bone by firm taps on the base of the applicator. To disengage the tack, roll the applicator to one side.

-STOPS

The stops change the pin retaining of the applicator. Using the green stop on the Specialist tip, the retaining is increased to a medium level; an harder retaining is obtained using the blue stop.

-PIN ORGANIZER (BOX)

After the sterilization, the pin organizer can be removed from the box or eventually leaved inside. Removed the cover, the pins can be placed in the organizer, ready for the use.

-HAMMER

The hammer has a head diameter and a weight calibrated for transfer to the applicator a force that allow the correct insertion of the pins Supertack.

REUSE

The MC Bio SUPERTACK instrumentation are made for several uses.

The instrumentation must be sterilized prior the use.

WARNING

It's the clinician's responsibility to become familiar with proper technique for use these devices. The clinician should have training and experience with bone fixation techniques.

The applicator must be held perpendicular to the bone during placement to avoid bending the tack. Complete seating should be attempted with a single strike of the hammer. The tacks are not designed to endure excessive or abnormal functional forces.

HOW SUPPLIED

All components are provided NON- STERILE.

The surgical box is provided NON- STERILE.

The surgical box is studied for a sterilization by a steam sterilizer with Dynamic Air Removal (class B Clave)

STERILIZING

MC Bio recommends sterilization of the instruments, by a steam autoclave method and in a steam sterilizer with Dynamic Air Removal. Parameters are in accordance at the European Pharmacopeia and the standard EN ISO 17665-1.

STEAM AUTOCLAVE METHOD

- ☐ Recommended a steam class B Clave
- ☐ Sterilization in accordance with the EN ISO 17665 or the European Pharmacopeia

STORAGE

The finished product shall be stored in closed rooms at temperatures from +10°C to + 30 °C and relative humidity not more than 75%.. Store the Osteosynt devices in consumer packaging and remove from consumer packaging only immediately prior to use. Preservation of the package integrity is an indispensable prerequisite.

SERVICE

If the tools loose their original mechanical conditions , due to the use or fall down, promptly replace the part or the tool.

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