MEDPOR®

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MEDPOR

Surgical Implants including TITAN[™]
and BARRIER[™] Implant

CE 0086
STERILEO

IMPLANT, TITAN, BARRIER, FOR MEDICAL USE ONLY

POREX

SURGICAL PRODUCTS GROUP

15 Dart Road, Newnan, GA 30265-1017 USA www.porexmedical.com
1-678-479-1610 • 1-800-521-7321 • FAX 1-678-423-1435

Sterility guaranteed unless inner package is opened, damaged or wet.
Caution: Federal law (U.S.A.) restricts this device to sale by or on
the order of a physician.

Product Information

Authorized European Representative
Representante Autorizado en Europa
Representant agréé pour l'Europe

Rappresentante autorizzato per l'Europa
Representante autorizada na Europa
Autorisierte Vertretung für Europa

POREX SURGICAL GmbH

Lindwurmstrasse 101
80337 München, Germany
tel: +49 (0)89 232415-0
fax: +49 (0)89 232415-15

STERILEO

Product Information - English

Description

MEDPOR[®] Surgical Implants, MEDPOR[®] TITAN[™] Surgical Implants, and MEDPOR[®] BARRIER[™] Surgical Implants are manufactured from a linear, high-density polyethylene biomaterial. MEDPOR TITAN Surgical Implants have titanium mesh embedded in the polyethylene biomaterial. Polyethylene and titanium have long histories of use in surgical implantable products. MEDPOR Surgical Implants, MEDPOR TITAN Surgical Implants and MEDPOR BARRIER Surgical Implants (hereafter together referred to as MEDPOR Implants unless otherwise specifically noted) in block, sheet and pre-formed shapes are intended for non-weight bearing applications of craniofacial reconstruction/cosmetic surgery and repair of craniofacial trauma. MEDPOR Surgical Implants and MEDPOR BARRIER Implants are also intended for the augmentation or restoration of contour in the craniofacial skeleton.

MEDPOR Implants are manufactured in a wide variety of shapes and sizes to meet many applications. The interconnecting open pore structure of the MEDPOR Biomaterial allows for tissue ingrowth except in the case of the MEDPOR BARRIER Implants. The porosity of the porous MEDPOR Biomaterial is maintained large, with average pore sizes larger than 100 micrometers and pore volume in the 50 percent range (measured by mercury intrusion porosimetry). The MEDPOR BARRIER Implant is a non-porous high-density polyethylene material. The BARRIER material acts as a barrier to prevent tissue ingrowth and may be heat bonded to one side of a porous implant. The porous side of the implant becomes vascularized in the same manner as a regular MEDPOR Implant. The BARRIER Implant designs are intended to be modified by the surgeon to position the barrier area in the location where tissue ingrowth is not desired.

Animal data has demonstrated that the porous MEDPOR Biomaterial permits tissue ingrowth. Histologic analyses of biopsies from human implants have also demonstrated tissue ingrowth. The clinical significance of tissue ingrowth may vary with the application and implant site. The long term clinical effect of this characteristic has not been established. Because of the white color of the biomaterial, it will not show through overlying tissue. The firm nature of the biomaterial allows the implant to be readily carved and shaped without fear of collapsing the pores. In vitro and in vivo biocompatibility studies have shown the MEDPOR Biomaterial to be free from any observable systemic or cytotoxic effects.

Instructions for Use

MEDPOR Implants are provided STERILE and NON-PYROGENIC in double peel pouch packaging. The implants should remain in their protective packaging until the surgical site has been prepared by the surgeon. Implants should remain in the sterile field throughout handling,

modification and implantation. MEDPOR Implants are intended as basic shapes which can be modified by the surgeon for appropriate use for an individual patient's situation.

The success of any implant is dependent upon careful handling and good surgical technique.

Surgical Procedure

This document is not intended to be a comprehensive document on surgical techniques. Surgeons should utilize proper surgical techniques for which they were trained, and their clinical experience, to determine appropriate surgical procedures. Successful implantations are technique sensitive. Sound surgical judgment should be used in the selection and the implantation of all MEDPOR Implants.

When planning the implant site, the incision, whenever possible, should be made so as to avoid closure directly over the implant.

The implant may be shaped, cut or carved with the techniques outlined in the Implant Preparation section.

Implant Preparation

MEDPOR Implants are sold STERILE and SHOULD NOT BE RESTERILIZED. Porous materials are particularly susceptible to contamination either by micro-organisms or foreign material. Prior to handling the implants, the operative personnel should thoroughly rinse their gloves to remove glove powder and other foreign material. Following removal of the implant from its protective package, it should be placed in a basin of clean, sterile physiologic saline. This greatly reduces the chance of contaminating the implant with lint and other materials.

MEDPOR Surgical Implants and MEDPOR BARRIER Surgical Implants can be readily carved to shape using a sharp scalpel without fear of collapsing the pores. Large blades can be used for the rough cuts while final trimming may be more readily accomplished with a smaller blade. If considerable carving is required, an occasional blade change may be necessary. MEDPOR Implants can be readily cut to shape using scissors without collapsing the pores. When shaping MEDPOR Implants, sharp edges should be avoided to prevent trauma to the surrounding tissues. Before implantation, examine the tapered edges of the implant to be sure that there is sufficient material to ensure edge integrity. Trim any poorly attached material from the edge of the implant. Following carving or cutting of any MEDPOR Biomaterial, thoroughly wash the implant in sterile saline to remove any loose particles generated by trimming.

As an aid to contouring and shaping, the implant can be submerged in a hot, sterile saline bath (82° - 100° C, 180° - 212°F) for several minutes. The hot saline should relax the memory of the material, enabling modification of the shape. Test the implant for flexibility, then gently bend to modify the shape. The material can be returned to the hot saline if too much resistance is encountered. Ensure the saline bath remains hot. Repeat these steps if further modification is required. The material should be held in the modified shape and allowed to cool. The use of a cold sterile bath can be used to accelerate the cooling process.

Prior to implantation, the MEDPOR Surgical Implants and MEDPOR BARRIER Implants may be rinsed in a sterile antibiotic solution.

When fixation of the implant is desired, stabilization may be accomplished with suture, surgical wire, or rigid fixation screws. In the case of screw fixation, tightening the screw will compress the implant to the bone and will enable the surgeon to sink the screw head flush with the implant surface. One advantage of stabilizing the implant is the ability to delicately shape and feather the edges of the implant in situ after fixation. Care should be taken to remove all carved debris from the surgical site.

Multiple pieces of implant material may be stacked and sutured together. A cutting needle can easily penetrate the surface of the implant enabling the surgeon to re-suspend tissue or muscle by attaching it to the implant. Proper rigid fixation techniques allow for stabilization of the implant.

The porous nature of MEDPOR Surgical Implants and MEDPOR TITAN Implants allows for soft tissue ingrowth and vascularization of the implant. In patients that may require later surgical revision, the surgeon should be aware of this vascular and soft tissue ingrowth. In the event revision or removal of the implant is required after ingrowth has occurred, the surrounding soft tissue should be raised off the implant surface.

Indications

MEDPOR Implants in block, sheet and pre-formed shapes are intended for non-weight bearing applications of craniofacial reconstruction/cosmetic surgery and repair of craniofacial trauma. MEDPOR Surgical Implants and MEDPOR BARRIER Implants are intended for the augmentation or restoration of bony contour in the craniofacial skeleton. A reference list of many articles, presentations and publications on MEDPOR Implants and porous polyethylene is available upon request from Porex Surgical, Inc.

Indications for the MEDPOR BARRIER Surgical Implant include situations in which a biomaterial may be used for augmentation or restoration of bony contour in the craniofacial skeleton, when the surgeon wishes to prevent tissue ingrowth into one area of the implant.

Contraindications

1. MEDPOR Implants should not be used in areas where there is inadequate coverage of healthy, well vascularized tissue.
2. MEDPOR Implants should not be used in areas of recent or active infection until the surgeon has determined that the infection has been completely eradicated.
3. MEDPOR Implants should not be used in patients with systemic disorders which cause poor wound healing or will lead to soft tissue deterioration over the implant.
4. Since tissue irradiated in cancer therapy has been compromised, using MEDPOR Implants in these areas may be problematic.
5. MEDPOR Implants should not be used in areas exposed to the outside environment.
6. MEDPOR Implants should not be used in areas where there are prolonged, significant, compressive stresses.
7. MEDPOR Implants should not be used in weight bearing areas such as joints like the TMJ (temporal mandibular joint).
8. There is no known clinical data to suggest that MEDPOR Implants are safe and efficacious in the following areas:
 - a. non obliterated frontal sinuses
 - b. alveolar ridge
 - c. nasopharyngeal applications

Therefore, it is not recommended for use in these applications.

Precaution

When trimming or otherwise modifying MEDPOR TITAN Surgical Implants, sharp metal edges can be exposed. Handle MEDPOR TITAN Surgical Implants carefully to avoid cutting or puncturing sterile gloves or hands.

Do not place or carve MEDPOR Implants on cloth or any other surface which may contaminate the implant with lint and other particulate matter.

Do not use electrosurgery devices, such as battery or electric cautery units, lasers and other devices to modify the implant material.

When deemed necessary, initial stabilization of the implant can be accomplished by suture or craniofacial screw anchorage.

If the implants are to be placed using the intraoral route, precautions should be taken to reduce, as far as possible, the risk of contaminating the implant. The use of impervious plastic film or a rubber dam to isolate the incision is recommended.

MEDPOR Surgical Implants and MEDPOR BARRIER Implants may be impregnated with an appropriate antibiotic. This can readily be accomplished in the sterile field by adding the antibiotic to the saline in the implant storage basin.

Proper surgical procedures and techniques are necessarily the responsibility of the medical professional. Each surgeon must evaluate the appropriateness of the procedure based on his or her own medical training and experience.

PRECAUTION: As with any surgery, infections and other complications not related to the implant may result. Because of the high surface area of porous implants, there is an increased opportunity of contamination with this type of implant. It is most important, therefore, that the highest level of aseptic care be used when handling and preparing MEDPOR Implants.

Porex Surgical, Inc.
15 Dart Road, Newnan, GA 30265-1017 USA
1-678-479-1610 • 1-800-521-7321 • FAX 1-678-423-1435 • www.porexurgical.com

©2007, Porex Surgical, Inc. ECO 076-85-M
MEDPOR TITAN and MEDPOR BARRIER are trademarks and MEDPOR is a registered trademark of Porex Surgical, Inc. The POREX SURGICAL service mark is owned by or used under the license and authority of Porex Corporation. MEDPOR TITAN Patent Pending USA Patent Nos 6,187,041; 6,530,953; 5,454,226; 5,433,748; 5,466,258 and 5,876,435 Canadian Patent No. 2536940 Australia Patent No. D428,992

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