To whom it may concern:

This letter serves to declare that the attached document is a true original document, "Instruction for Use for a Medical Device," produced for Mentor.

For and on behalf of Mentor

Martin G. Sprunck
Associate Director, Regulatory Affairs

11/17/2016

Date

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document, to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.



State of California
County of Orange

On November 17, 2016, before me, Jamie L. Byun, Notary Public, personally appeared Martin G. Sprunck, who proved to me on the basis of satisfactory evidence to be the person whose name is subscribed to the within instrument and acknowledged to me that he executed the same in his authorized capacity, and that by his signature on the instrument the person, or the entity upon behalf of which the person acted, executed the instrument.

I certify under Penalty of Perjury under the laws of the State of California that the foregoing paragraph is true and correct.

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**INSTRUCTION FOR USE FOR A
MEDICAL DEVICE**

**CPX4 Breast Tissue Expander Contour
Profile, Textured, with Centerscope
Magnetic Detection Device and Winged
Infusion Set**

manufactured by
Ментор / Mentor
3041, Skyway Circle North Irving, Texas,
USA, 75038 / США

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1. NAME OF A MEDICAL DEVICE

CPX4 Breast Tissue Expander Contour Profile, Textured, with Centerscope Magnetic Detection Device and Winged Infusion Set

Configurations:

1. CPX4 Breast Tissue Expander Contour Profile, Textured, Low Height (Style 8100)
2. CPX4 Breast Tissue Expander Contour Profile, Textured, Med Height (Style 8200)
3. CPX4 Breast Tissue Expander Contour Profile, Textured, Tall Height (Style 8300)
4. CPX4 Breast Tissue Expander Contour Profile, Textured, with Suture Tabs, Low Height (Style 9100)
5. CPX4 Breast Tissue Expander Contour Profile, Textured, with Suture Tabs, Med Height (Style 9200)
6. CPX4 Breast Tissue Expander Contour Profile, Textured, with Suture Tabs, Tall Height (Style 9300)

материалом, высокий профиль (CPX4 Breast Tissue Expander Contour Profile, Textured, with Suture Tabs, Tall Height) Style 9300.

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Accessories:

1. Centerscope Magnetic Detection Device
2. Instructions for Use
3. Stickers - not more than 8

2. INFORMATION ON MEDICAL DEVICE MANUFACTURER

Mentor, 3041, Skyway Circle North, Irving, Texas 75038, USA

3. PURPOSE OF A MEDICAL DEVICE ESTABLISHED BY THE MANUFACTURER, INDICATIONS, CONTRAINDICATIONS

CPX4 Breast Tissue Expander Contour Profile, Textured, with Centerscope Magnetic Detection Device and Winged Infusion Set (“the expander”) used for breast reconstruction following mastectomy, and are intended for temporary subcutaneous or submuscular implantation.

Indications:

It is the responsibility of the surgeon to advise the prospective patients or their representatives, prior to surgery, of the indications associated with the use of these products.

Breast reconstruction after mastectomy, correction of an underdeveloped breast, scar revision and tissue defect procedures. The prosthesis are intended for temporary subcutaneous or submuscular implantation; they should be used within a time frame determined by the physician to achieve the

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clinically desired degree of tissue expansion in term four to six months.

Contraindications:

It is the responsibility of the surgeon to advise the prospective patients or their representatives, prior to surgery, of the contraindications associated with the use of these products.

Patient Groups in which the product is contraindicated:

The use of expanders is contraindicated in patients who have any of the following conditions:

- Implanted devices such as pacemakers, drug infusion devices, artificial sensing devices, etc., that would be affected by a magnetic field.
- Active infection anywhere in the body.
- Existing malignant or pre-malignant breast cancer without adequate treatment.

Surgical practices in which expander use is contraindicated due to compromise of expander integrity:

- Do not alter the expander shell or valve.
- Do not place inside the implant any drugs or substances other than sterile saline for injection.
- Do not allow the device to come into contact with Betadine.

NOTE: The satisfactory use of either tissue expander for tissue replacement following mastectomy or trauma may require special reconstructive procedures.

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Method of Application: surgical practice.

Conditions of the Medical Device Application: at health care organizations.

4. DESCRIPTION OF PRINCIPLES DEFINING FUNCTIONING OF A MEDICAL DEVICE AND THEIR FEATURES

Prosthesis is used for breast reconstruction following mastectomy, and is intended for temporary subcutaneous or submuscular implantation. Tissue expansion in breast reconstruction typically requires four to six months. For this purpose, the prosthesis is placed subcutaneously in a patient and filled with pyrogen-free sodium chloride solution for injection.

5. DESCRIPTION

In order to provide these tissue expanders with elasticity and integrity, the shells are made with successive cross-linked layers of silicone elastomer. Superior and anterior reinforcement allows for directional expansion in the lower pole of the devices; self-sealing Patch area with self-sealing technology (containing silicone gel) to the front patch of the device to minimize and/or prevent leakage in the event of an accidental needle puncture. The textured shell provides a disruptive surface for collagen interface.

Identification of the injection dome site can be accomplished by use of the CENTERSCOPE® Magnetic Injection Port Locator provided with the Tissue Expander. To use the

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CENTERSCOPE® Magnetic Injection Port Locator, follow the instructions provided with the specific dome detection device.

Injections must be made using sterile, pyrogen-free Sodium Chloride U.S.P. Solution and into the injection dome area. If injections are made on or outside the injection dome, leakage can occur.

6. APPLICATION

Expander selection

Some of the important surgical and implant sizing variables that have been identified include the following:

- The tissue expander should not be too small or too large in comparison to the patient's chest wall dimensions.
- Available tissue must provide adequate coverage of the device.
- Submuscular placement of the expander may be preferable in patients with thin or poor quality tissue.
- A well-defined, dry pocket of adequate size and symmetry must be created to allow the implant to be placed flat on a smooth surface.
- Avoid too small of an incision.

Testing procedures for expanders

The device should be tested for patency and shell integrity immediately prior to use. This can be accomplished by the following steps:

1. Using a 21 gauge standard needle, partially inflate the device with air through the injection dome.
2. Submerge the air-filled prosthesis in sterile, pyrogen-free testing fluid (water or saline).
3. Apply mild pressure and check for possible punctures or leaks.

Maintaining hemostasis/avoiding fluid accumulation

Careful hemostasis is important to prevent

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postoperative hematoma formation. Should excessive bleeding persist, the implantation of the device should be delayed until bleeding is controlled.

Postoperative evacuation of hematoma or seroma must be conducted with care to avoid expander contamination, or damage from sharp instruments.

Magnetic detection device Centerscope (hereinafter – detection device)

Instructions For Use

Using sterile technique, remove the detection device from the sterile pouch. Grip locator with either hand, ensuring that the free-swinging magnetic arm is pointing away from the hand holding it. Place the base of the unit (the flat side) on the patient. Move the unit in a circular motion until the magnetic arm detects (points towards) the location of the injection dome. Follow the direction that the arm points towards until the arm is pointing straight towards the hole in the base of locator (the target). When the arm is centered perfectly in the target, the injection site has been located.

To mark the injection, follow one of these options:

Option 1: With the magnetic arm centered in the target, make a mark with a surgical sterile marker in each of the three notches around the anterior perimeter of the base. Next, make a fourth mark in the hole located behind the magnetic arm that runs through the center base of the locator. After all four marks have been made, lift the device from the patient.

Using the same marker, carefully connect the opposing dots with a line creating an injection 'crosshair'. Then, make a clear mark at the point where the two perpendicular lines intersect. This point is the injection site.

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
Option 2: With the magnetic arm centered in the target, gently but firmly press the device against the patient. Hold for several seconds. Lift the device from the patient. The raised 'X' on the bottom of the locator will leave a clearly imprinted 'crosshair' mark on the patient's skin. Mark the center of the crosshair with a surgical sterile marker. This point is the injection site.

Option 3: Utilize a combination of options 1 and 2 for further confirmation of the best point for injection.


Expander filling procedure

To inflate either tissue expander:

1. Identify the injection dome site using the detection device Centerscope included with the tissue expander.

2. Once the center of the dome has been identified, a skin marker can be used to identify the area of injection.

3. Inflation is accomplished by premarking the skin, inserting a 21 gauge standard needle into the top of the

injection site, perpendicular $\pm 30^\circ$ to the base and filling the device using sterile, pyrogen-free Sodium Chloride U.S.P. Solution for injection.

4. Injections must be made into the injection dome. If injections are made on or outside the injection dome leakage can occur. Although the device has a selfsealing around the area of the injection dome, DO

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NOT ATTEMPT TO INJECT OUTSIDE THE DOME AS LEAKAGE MAY OCCUR.

Postoperative care

Mentor recommends that the patient be wrapped superiorly with an elastic bandage, taped laterally, and wear a surgical bra 24 hours a day to help prevent shifting of the device.

Recording procedure for expander

Each expander is supplied with a patient record label showing the catalog number and lot number for that unit. One of these pressure sensitive labels should be attached directly to the patient's chart. The date of placement, expansion data (date and volume), and date of explant should be indicated on the label.

7. SPECIFICATIONS OF A MEDICAL DEVICE

The expander is based on the principle of positioning the injection dome beneath the muscular tissue of a patient, and its subsequent filling with saline solution. To locate the injection dome, a special magnetic detection device, Centerscope, is used. It helps make crosshair marks on the patient's skin and thus locate the injection site.

текстурированной оболочкой (CPX4 Breast Tissue Expander Contour Profile, Textured)/ Figure 1. CPX4 Breast Tissue Expander Contour Profile, Textured

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After choosing an injection site and verifying it with the Centerscope magnetic detection device, fill the expander injection dome of perpendicularly (at right angle) $\pm 30^\circ$ to the base, using pyrogen-free sodium chloride solution for injection

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8. BASIC PARAMETERS AND CHARACTERISTICS OF A MEDICAL DEVICE

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** The product has been tested to meet the requirements of ASTM F1441 for over-expansion to 200% of maximum recommended volume and checked for leakage.

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
(Injection Port Competence): Apply 120-mm Hg intraluminal pressure to the port using water or test media with demonstrated equivalence. Using the prescribed gauge hypodermic needle, puncture the port 5 consecutive times within 1 mm² at a site near the center of the port. The port is considered leaking and fails the test if beads of fluid on the port surface are not static after 30 s.

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SET OF DELIVERING

- 1 Expander
- 1 Centerscope Magnetic Detection Device
- 1 Winged Infusion Set: 1 piece
- Stickers – not more than 8
- 1 Instruction for Use

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9. SAFETY REQUIREMENTS

1. Magnetic Fields

- DO NOT use these expanders in patients that have a previously implanted device that could be affected by a magnetic field. MRT should not be performed on a patient implanted with these expanders because the risk of movement may exist. It may cause pains or lead to displacement of the expander that may require a revision surgery.

2. Radiation Therapy

- Mentor has not evaluated *in vivo* effects of radiation therapy during the use of expanders - therefore, the safety of such use is not assured. The decision regarding the use of these expanders in patients, who will undergo the radiation therapy, should be made by a surgeon and a radiation oncologist.

3. Extrusion of the expander

- The incidence of extrusion of any expanders has been shown to increase when the expander is placed in injured areas: scarred, heavily irradiated or burned tissue, crushed bone areas or where severe surgical reduction of the area has previously been performed.

4. These expanders are temporary devices.

- They should be used within a time frame determined by the physician to achieve the clinically desired degree of tissue expansion.
- Expander use in breast reconstruction typically requires four to six months

5. Reuse

- Expanders are for single use only. Do not

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resterilize.

6. Avoid Damage During a Surgery

- Care should be taken to prevent damaging the expander with surgical instruments.
- Do not insert or repair a damaged expander.
- Take care during such procedures as tissue expansion, open capsulotomy, breast pocket revision, hematoma or seroma aspiration, and biopsy/lumpectomy to avoid damage to the expander shell, injection dome or valve.
- Do not contact the expander with disposable, capacitor-type cautery devices.

7. Proper Filling

- Surgeons should ensure themselves of the position of the fill dome prior to adding or withdrawing fluid. Needle punctures on or outside the palpation ring of the expander may penetrate the shell causing deflation or compromise the fill dome and necessitate replacement of the expander. Although the expander has a self-sealing BUFFERZONE area around the injection dome, **DO NOT ATTEMPT TO INJECT INTO THE AREA AROUND THE DOME**, as damage to the expander may occur.
- Mentor relies on the surgeon to select the optimum incision and pocket size for the chosen tissue expander configuration, style and projected volume.

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9. Instructions to Patient

- The patient should be advised that a vigorous body movement (e.g., physical exercise) or excessive manipulation or trauma in the region of the expander may cause the stress to the expander and result in a subsequent deflation.

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10. PRECAUTIONS

- It is the responsibility of the surgeon to advise the prospective patients or their representatives, prior to surgery, of the

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possible complications associated with the use of these products.

- Preexisting infection should be treated and resolved before implantation of either tissue expander.
- Any surgeon performing reconstructive mammaplasty with tissue expanders should be familiar with the currently available techniques for measuring the patient, determining the tissue expander size and performing the surgery.
- Lint, dust, talc, surgical glove powder, drape and sponge lint, fingerprints, skin oils and other surface contaminants deposited on a tissue expander by improper handling may cause foreign body reactions. Strict adherence to clean, aseptic techniques should be maintained to prevent tissue expander contamination and possible complications. Surgical instruments and gloves should be rinsed clean of any impurities before handling a tissue expander.
- The silicone elastomer shell may easily be cut by a scalpel or ruptured by excessive stress, manipulation with blunt instruments or penetration by a needle. Subsequent deflation and/or rupture will result. All prostheses should be carefully inspected for structural integrity prior to and during implantation.
- Any subsequent surgical procedures in the tissue expander area should be undertaken with extreme caution as damage could occur. In the event that a tissue expander is damaged, it must be removed.
- Each device should be checked for patency prior to surgery and continuously monitored throughout the surgical procedure to ensure the structural integrity of the device is not compromised in any way. A standby tissue expander should be available at the time of surgery.

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- Potential for contamination exists when fluid is added or removed from the device. Use the aseptic technique in the introduction of saline into a tissue expander, a single-use, sterile saline container is recommended.

11. ENVIRONMENTAL CONTROL REQUIREMENTS DURING THE USE OF A MEDICAL DEVICE

This medical device does not make a negative impact on humans and environment while used, transported and stored.

12. MATERIALS USED IN PRODUCT MANUFACTURE

Expander: Shell – MED 4750 (NuSil) (polydimethylsiloxane) silicone elastomer; base palpation ring, base, domelike spring washer – MED 4750 (NuSil) (polydimethylsiloxane) silicone elastomer. Needle protector – AISI 304 stainless steel (low carbon austenitic steel). Injection port magnet – neodymium rare-earth magnet.

Self-sealing technology, front patch of the device – dacron-reinforced MED 4750 (NuSil) (polydimethylsiloxane) silicone elastomer coating. Textured coating, dome, filled disk patch, bubble positioning loop – MED 4750 (NuSil) (polydimethylsiloxane) silicone elastomer. Placement loop – MED 4750 (NuSil) (polydimethylsiloxane) silicone elastomer, 6118 polyester mesh. Needle damper – AMOCO UDEL P-1700 polysulfone. Silicone adhesive – RTV cold-setting adhesive.

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Centerscope Magnetic Detection Device: base – Zytel 101 NC010 nylon (general purpose polyamide 66, for die casting and extrusion), magnet – neodymium (a rare-earth metal).

Winged Infusion Set: protective sheath – Bormed HE7541-PH polyethylene, L.L. connector, tube, wings – DAIREN GREEN PVC 9154, needle – AISI (The American Iron and Steel Institute) 304 stainless steel, safety-lock – Bormed HE7541-PH polyethylene.

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13. METHODS AND CONDITIONS FOR STERILIZATION

The products are provided sterile. The expander is dry heat sterilized and is for single use only. Do not resterilize.

Centerscope Magnet Detection Device: sterilized by irradiation.

Winged Infusion Set: sterilized by ethylene oxide.

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14. PACKAGING, LABELING, STORAGE AND TRANSPORTATION

The expander is placed in a polyethylene terephthalate glycol cellular tray sealed with a Tyvek film. The Centerscope Magnetic Detection Device is packed in two sterile Tyvek pouches. The Winged Infusion Set is packed in one Tyvek pouch.

Individual packaging ensures that the device is sterile.

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Detection Device and Winged Infusion Set

REF: xxxxxxxx

LOT see package

Intended purpose: surgical practice

Delivery set:

CPX4 Breast Tissue Expander (Style 8300): 1 unit;
Centerscope Magnetic Detection Device: 1 unit; Winged Infusion Set: 1 unit

Product description: The expander consists of the following elements: Silicone Elastomer Expander and Silicone Elastomer Integral Injection Dome.

LOT see package

Shelf life: 4 years

MENTOR, 3041 Skyway Circle North Irving, TX 75038-3540, USA

Importer: Johnson & Johnson LLC, 121614, Moscow, Krylatskaya str., 17/2, phone: +7(495) 580 77 77

Marketing Authorization No. -

Sterile unless the package has been opened or damaged. Do not resterilize

3548121-xxx

All device components are packed in individual, sterile packages which are placed into laminated cardboard boxes with instructions for use.

The medical devices may be transported by all means of transport in covered transportation vehicles as required by shipping rules and regulations applicable for each type of transport.

When transported and stored in its original packaging, the medical devices are resistant to climatic factors within the temperature range from +5 to +43°C, at relative humidity not exceeding 80%, without condensation, at atmospheric pressure from 200 to 1,500 Mpa.

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15. USE

The expanders may be used within the temperature range from +32 to + 42°C.

The Centerscope Magnetic Detection Device and the Winged Infusion Set may be used within the temperature range from +10 to + 40°C, at a relative humidity of 80%.

16. WARRANTY

The Manufacturer undertakes a warranty obligation to end users in relation to a medical device specified in the Transfer Protocol, and ensures that it functions faultlessly, has no existing defects and that such defects (if any) will be surely handled, within 12 months from the sales date, provided that all recommendations are observed.

Merchandise returned must have all manufacturer's seals intact and be returned within 60 days from the date of invoice to be eligible for credit or replacement. Please contact the Mentor Customer Service Department for details. Returned products may be subject to restocking charges. Authorization for return of merchandise should be obtained from your local Mentor representative. Other conditions noted above also apply.

17. DURABILITY (SHELF LIFE) AND PROHIBITION FOR USE OF A MEDICAL DEVICE UPON EXPIRATION OF DURABILITY PERIOD (SHELF LIFE)

The extender shall maintain sterility in its original package during 4 years.

The implantation period shall be determined by the physician to achieve the clinically desired

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degree of tissue expansion. Prosthesis use in breast reconstruction typically requires four to six months

18. PROCEDURE FOR UTILIZATION OR DISPOSAL OF A MEDICAL DEVICE

Utilization and disposal of the medical device shall be performed in compliance with the requirements of SanPiN 2.1.7. 2790-10 “Sanitary and Epidemiological Requirements for Managing Medical Waste” and the approved regulations of the healthcare institution as class B medical waste.

19. POTENTIAL ADVERSE REACTIONS ASSOCIATED WITH USE OF A MEDICAL DEVICE

The complications which may result from the use of this expander include the risks associated with the medication and methods used in the surgical procedures as well as the patient’s degree of intolerance to any foreign object placed in the body. The complications may include, but are not limited to, the following:

Additional Surgeries

- Additional surgery will be required to either replace a deflated tissue expander and/or to complete the breast reconstruction procedure.

Cancer

- Published studies indicate that the breast cancer is no more common in women with implants than in those without implants.
- Capsule Formation and Contracture

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- Postoperative formation of a fibrous tissue capsule around an implanted expander is a normal physiologic response to the implantation of a foreign object. Capsule formation occurs in all patients in varying degrees. Capsules range from thin to heavily-thickened.

- Contracture of the fibrous capsule may occur, independent of its thickness.

- Discomfort, pain, excessive tissue firmness and misshapen expanded tissue, deflation, increased palpability and wrinkling and/or displacement of the expander may occur and may require surgical intervention. In some patients, tissue firmness may recur subsequent to corrective surgical procedures.

Complications of Tissue Expansion

- Tissue thinning or necrosis.
- Sloughing of poorly vascularized tissue.
- Gross postoperative hematoma, manifested by enlargement, tenderness and discoloration leading, if untreated, to extrusion of the expander.
- Undue pressure on the tissue located over the expander or trauma to surrounding tissues which may lead to venous thrombosis, the breakdown of skin over the expander and subsequent extrusion. Deflation or removal of the expander may be necessary for tissue repair.

Connective Tissue Disease

- Concern over the association of expanders or any other implants to the development of autoimmune or connective tissue diseases,

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such as lupus, scleroderma, or rheumatoid arthritis, was raised because of cases reported in the literature with small numbers of women with implants. A review of several large epidemiological studies of women with and without implants indicates that these diseases are no more common in women with implants than those in women without implants.

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Deflation/Rupture/Leakage

- Saline-filled expanders deflate when saline solution leaks through an unsealed or damaged valve or through a break in the expander shell. Deflation can occur immediately or progressively over a period of days and is noticed by loss of size or shape of the expander. Additional surgery is needed to remove deflated expanders.

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Dissatisfaction with Cosmetic Results

- Incorrect expander size, inappropriate scar location or appearance and misplacement or migration of expanders may interfere with a satisfactory cosmetic result. These complications are generally associated with the surgical procedure and technique.

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Extrusion of Expander/Interruption of Wound Healing

- Skin necrosis and/or sloughing may result from undue tension on the skin overlying the expander, trauma to the skin during surgical procedures or inadequate tissue thickness inhibiting circulation. Subsequent

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exposure and/or extrusion of the expander may occur.

- Displacement, twisting, fracture or extrusion may occur from improper expander sizing and/or placement, e.g., when the expander is too large or the pocket is too small or when there has been inadequate preoperative assessment of stresses causing movement of the expander.
- The incidence of extrusion of the expander has been shown to increase when the expander has been placed in injured areas: scarred, heavily irradiated or burned tissue or crushed bone areas; where severe surgical reduction of the area has been performed; and where steroids are used in the surgical pocket.

Fluid Accumulation

- Excessive postoperative fluid accumulation and transient reaccumulation of fluid around the expander as a result of trauma and after vigorous exercise have been reported.

Hematoma

- Careful hemostasis is important to prevent postoperative hematoma formation. Should excessive bleeding persist, it is recommended that the expander not be used until bleeding is controlled.
- Gross postoperative hematoma, manifested by enlargement, tenderness and discoloration of tissue may, if untreated, lead to extrusion of the expander.

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Infection

- Infection, manifested by swelling, tenderness, pain and fever, may appear in the immediate postoperative period or at any time after insertion of the expander. In the absence of classic symptoms, subacute or chronic infections may be difficult to diagnose. If infection does not subside promptly with the appropriate treatment, removal of the expander is indicated.
- Toxic Shock Syndrome has been reported as a complication of both augmentation and reconstructive mammoplasty.

Magnetic Field

- The potential for adverse events could occur in patients who have any of the following conditions: Implanted devices such as pacemakers, drug infusion devices, artificial sensing devices, etc., that would be affected by a magnetic field. The tissue expander is specifically contraindicated in such patients to prevent any potential adverse events related to interference with these types of devices.
- The potential for adverse events such as pain or expander displacement requiring revision surgery could occur in patients undergoing Magnetic Resonance Imaging (MRI). The labeling includes warnings and symbols stating that the tissue expander is NOT MRI Safe and that MRI is NOT to be used on a patient implanted with these devices.

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Pain

- Pain may be felt of varying severity (degrees) and duration (length of time) during the tissue expansion process.

Wrinkling of the Expander

- Surgeons have reported that in some patients, visible or palpable wrinkling of the envelope, usually associated with textured prostheses, has occurred. Folds in the envelope can be visible beneath the overlying skin. This is reported to occur more frequently with: thin-skinned patients with little or no subcutaneous fat; subglandular rather than submuscular placement of the expander; a tissue expander that is too large relative to the pocket size or frame of the patient; overlying tissue that is minimal or of poor quality; or where there is contracture and/or insufficient fill volume.

20. AUTHORISED REPRESENTATIVE OF MANUFACTURER ON THE TERRITORY OF RF

Should any claims occur, please do not hesitate to contact Johnson & Johnson, Limited Liability Company; Johnson & Johnson, LLC, Russia, 121614, Moscow, Krylatskaya str. 17, building 2, phone (495) 580 77 77, fax (495) 580 78 78.