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INSTRUCTIONS

- For the Use of Medical Product
- Zirconia Blanks for Denture Frameworks
- Manufactured by BETTINI S.p.a., Italy

2. Moer

1. **Production.Designation**

The Product is intended for the fabrication of denture frameworks.

The application area of this medical product is dentistry.

NOTE: Dentures may be fitted in the patient's mouth only by medical personnel (dentists) working in accordance with the applicable laws.

It is applied in patient care and prophylactic institutions.

2. **Product Design**

Zirconia (zirconium dioxide) blanks are shaped like discs, cylinders or blocks (Fig. 1-6). The shape of billet depends on the machine used to mill off the denture framework.



Fig.1 Disc

d — Diameter
 s — Thickness

Fig. 2 Disc Drawing



Fig.3 Cylinder

// - Diameter
 - Thickness

Fig. 4 Cylinder Drawing



Fig.5 Block

b – width;
 h – height;
 l – length.

Fig. 6 Block Drawing

3. Basic Technical Information

Parameters	Discs	Cylinders	Blocks
Dimensions ($\pm 2\%$)	Diameter 98 mm; Thickness: 10 mm, 12 mm, 14 mm, 16 mm, 18 mm, 20 mm, 22 mm, 25 mm	Diameter: 16 mm, 20 mm; Height: 16 mm, 20 mm	Length: 40mm, 42mm, 43mm, 62mm, Height: 14 mm, 16 mm, 20 mm, Width: 15mm, 20mm, 25 mm
Weight	DB+: 238 g, 284 g, 326 g, 368 g, 410 g, 452 g, 494 g, 536 g DB+LUX: 238 g, 290 g, 334 g, 378 g, 424 g, 468 g, 512 g, 556 g	10 g, 20 g	26 g, 40 g, 42 g, 52 g, 54 g, 70 g, 96 g, 112 g, 142 g, 176 g
Shape	Disc	Cylinder	Rectangular parallelepiped
Color	White (DB+); Matt white (opal) (DB+LUX)	White (DB+)	White (DB+)
Material	95 % ZrO ₂ + 5 % Y ₂ O ₃ ; 94% ZrO ₂ + 5% Y ₂ O ₃ + 1% Al ₂ O ₃	95 % ZrO ₂ + 5 % Y ₂ O ₃	95 % ZrO ₂ + 5 % Y ₂ O ₃

Sterility	Non-sterile	Non-sterile	Non-sterile
Sintering temperature	1450° C	1450° C	1450° C
Flexural strength	Over 1400 MPa	Over 1400 MPa	Over 1400 MPa
Compressive strength	Over 2800 MPa	Over 2800 Mpa	Over 2800 Mpa
Firing density	6,05 g/cm ³	6,05 g/cm ³	6,05 g/cm ³
Fabrication technology	CAD/CAM technology		
Storage temperature	+5° C - +50° C		
Resistance to staining. Number of spots (after storage in staining solution)	0		
Porosity (number of pores)	≤ 16 (with diameter more than 30 μm) on the surface 1 mm		
Thermal expansion coefficient, 1/ 10 ⁶	25-400° C: 10,0 (DB+); 10,3 (DB+LUX) 25-800° C: 10,5 (DB+); 11,0 (DB+LUX)	25-400° C: 10,0 25-800° C: 10,5	25-400° C: 10,0 25-800° C: 10,5
Linear firing shrinkage	≤ 16 %		
Chemical solubility	≤ 0,5 %		

Blanks of various standard sizes are made of zirconia (zirconium dioxide) which is a high-tech ceramic material of the following composition: DB +: 95% ZrO₂ + 5% Y₂O₃, DB + LUX: 94% ZrO₂ + 5% Y₂O₃ + 1% Al₂O₃.

The shape of billet to be used depends on the machine used to mill off the denture framework.

Zirconia is the most up-to-date and advanced material which is currently used in dentistry. Materials based on zirconia have high fractural strength and crack resistance, as well as low thermal conductivity.

In general, zirconia is white, but it has the ability to scatter light which makes it looks like natural teeth.

Indications and Contraindications

Indications:

Lack of defect of a dental crown

Contraindications:

The product can be used only to fabricate frameworks for dentures intended to be worn in the mouth. They may not be inserted into the gum tissue or jaw bone, as well as to be fixed with screws or other similar components. DB+ and DB+LUX products may not be used for the manufacture of screws or other similar components because typology and classification of their application differ from typology and classification of denture frameworks. Even if so far there have been no reports of any severe allergies, it is the dentist's responsibility to check whether potentially allergic patient can wear DB+ or DB+LUX product in his/her mouth.

Possible Side Effects / Undesirable Reactions:

No possible side effects and undesirable reactions were observed in any patient.

Preparation to Work and Procedure of Work

Application Conditions

The Product may be applied at a temperature of +5°C - +50°C and humidity lower than 30-40%. Never expose it to moisture and frost (intensely low temperature).

Product Handling:

DB+ and DB+LUX products are very fragile; you must handle them with extreme care.

Do not let your patient swallow DB+ or DB+LUX product.

DB+ and DB+LUX products may be fitted in the mouth only if sintered as described in paragraph 4 herein.

- 1) Fabricate a master cast based on impression made by the dentist (Fig. 7, 8)

Fig. 7


Fig. 8

- 2) Scan the master cast and make modeling of the future framework on the computer (Fig. 9, 10)


Fig. 9

Fig. 10

- 3) Fix the billet in the holder of milling machine (Fig. 11)


Fig. 11

- 4) Program the milling process in accordance with the desired dimensions of final product. Then carve the framework of future denture from the billet on an automated milling machine (Fig. 12, 13)


Fig. 12

Fig. 13

- 5) Work on the framework with a diamond or tungsten carbide cutter.
6) While working on the framework, clean it with an air jet.
7) Dry the framework, using an infra-red lamp or drying kiln.
8) Prior to sintering, make a visual inspection of the framework to make sure that it is intact (Fig. 14)


Fig. 14

- 9) Place the framework in a high-temperature sintering furnace, and start the sintering cycle in accordance with the manufacturer's instructions.
- 10) On completion of the sintering, remove the framework from the furnace, compare it with the reference model to check its correspondence to the reference model, and work it with the diamond or tungsten carbide cutter, when needed (Fig. 15)


Fig. 15

- 11) Sand blast the sintered framework. If necessary, use silica alumina with the grain size of 60/100 μm at a pressure of 1-2 bar.
- 12) Make hygienic treatment of the framework with water.
- 13) Then apply dental ceramics to the framework and place the framework into the furnace for firing. Polish the framework, when needed (Fig. 16).


Fig. 16

RECOMMENDATIONS ON THE TREATMENT OF DENTAL ARCH AND SINGLE TOOTH

Should you need to fabricate a denture which is a set of several artificial teeth in a row or an entire dental arch (dental curve), it is essential to ensure homogeneity of all surrounding materials. This can be done by leaving small "junctions" or "connections" made of the same material as the basic workpiece, and cut at the completion of the sintering. This will also help to prevent the strain of the workpiece machined during the sintering, as well as to place the workpiece in the furnace so that to avoid undesirable contacts.

To fabricate several teeth frameworks or single tooth framework from the same billet, it is necessary to allow some time for filling cavities with materials that would dampen vibrations during their/its milling (such as wax, gypsum, silicone, etc.). Vibrations can adversely affect the homogeneity of ceramic material (dental ceramics) not yet exposed to sintering.

SHRINKAGE CALCULATION

Ceramic materials experienced pre-sintering are subject to shrinkage during sintering, which results in reduction of their initial size. Prior to milling, it is necessary to know the final size (final dimensions) of denture:

(a) if the method of multiplication factor is used, the final dimension of denture must be multiplied by this factor. The factor is approximately $1,250 \pm 0,007$, and each lot of blanks has its own factor indicated on the package;

(b) if the method of shrinkage percentage calculated with the formula $[(\text{pre-sintering} - \text{sintering}) / (\text{pre-sintering})] * 100$ is used, each final dimension of denture is multiplied by the following value: $100 / (100 - \text{shrinkage percentage})$, indicated on the package of each lot of material.

The values obtained in this manner (for each size) by either method are the dimensions of denture before its sintering in the furnace.

Shrinkage factor values were obtained by Bettini S.p.a. through a number of tests for machinability of materials and their density before and after the sintering.

The shrinkage factor is defined (refer to flow chart) as a difference in shrinkage between pre-sintered and sintered samples which are subsequently compared in appearance, weight and dimensions.

• WARNINGS:

Machining requires the installation of appropriate exhaust ventilation to prevent inhalation of dust. Continuous use of air jet for cleaning the frameworks in the course of machining is not recommended so that to avoid dust blowing.

Wear a face mask with fine filter to protect the respiratory apparatus.

Wear goggles during machining.

Do not ingest the product.

DB+ and DB+LUX products may be fitted in the mouth only if sintered as described in the Operating Instructions.

The products labeled DB+ or DB+LUX are supplied unsterile.

For sterilization:

Zirconia (yttrium-stabilized zirconium) is polymorphic material. To avoid deterioration of the qualitative characteristics during transformation, do not subject the material to sterilization with hot steam (autoclaving) or chemicals (except ethylene oxide). The application of gamma radiation may cause change in its color (discoloration).

Zirconia components are sterilized during the sintering cycle which takes place at a temperature of 1450 °C for two hours (the sintering medium is air). This process ensures complete decomposition and elimination of any bacteria or virus (which may be present in the ceramic parts) under high temperature, as well as complete burning of any residues of pyrolyzed carbon when exposed to air in the furnace.

Absolute sterility of personnel (dental technicians), sintering devices, denture packaging plants, including packaging material itself, furnaces and workplaces must comply with applicable health standards.

Note: zirconia components must be sterilized just once.

NOTE: Final machining of the sintered material can change biological, chemical, physical, and mechanical properties of the product: temperature that can be reached during machining (milling) of pre-sintered DB + SINT and DB + LUXSINT; and DB + DB + LUX materials (if they need minor correction) may not exceed 70 °C.

When machining use a jet of water to cool the ceramic product.

Sintering Cycle Proposed for DB+ and DB+LUX Components (Fig. 17).

Ambience: air

Exposure Time at 1450 °C: 2 hours

- Mechanical properties specified in the technical data were determined for the specific cycle; when different cycle is used, it should be checked by authorized personnel.


Fig. 17

6. Maintenance and Repair

The Product must not be repaired or reused. In case of damage (scratches, splits, etc.), the Product must not be further machined and fitted in the mouth.

7. Packaging and Transport. Delivery set

Each disc is packed in individual box that ensures its safety and prevents its environmental contamination (Fig. 18).


Fig. 18

Blocks are packed in boxes, 2, 11 or 12 blocks each (Fig. 19).


Fig. 19

Cylinders are packed in boxes, 28 blocks each (Fig. 20).


Fig. 20

Each piece is packed to guarantee safe storage.

The piece in inner package (consumer container) is put in transportation packing.

Transportation is carried out by all kinds of covered vehicles in compliance with the current regulations applied to this type of transportation means at a temperatures between $+ 5^{\circ} \text{C}$ and $+ 50^{\circ} \text{C}$.

8. Marking


Fig. 21

The inner package (consumer container) must contain the following information:

- name and address of manufacturer;
- name of Product;
- brand of material (DB+ and DB+ LUX);
- number of lot (LOT);
- number of pieces in a box (package);
- shrinkage percentage;
- product dimensions;
- bar-code;
- and symbols:

CE marking

Keep away from water and cold

Store at the indicated temperature

Read the instructions

Non-sterile

9. **Storage and Application**

Do not expose the Product to moisture and cold (low temperature).

Store the blanks at a temperature of +5 °C - +50 °C, and humidity lower than 30-40%.

The shelf-life of the Product is unlimited.

Do not store the Product together with aggressive media (corrosive fluids).

The Product must not be repaired or reused. In case of damage (scratches, splits, etc.), the

Product must not be further machined and fitted in the mouth.

Information on Inadequate Quality

To maintain a high level of customer satisfaction, we ask dentists and dental technicians to report all cases of inadequate quality of DB+ and DB+LUX products specifying the lot number, code number and essence of the claim. The sample should be sent to the manufacturer which shall, at its own expense, make the relevant examination.

Disposal

The waste must be collected in a yellow or yellow-marked disposable container, and disposed by specialized organizations in compliance with the applicable regulations and current laws.

Class of hazard medical wastes in accordance with the characteristics of morphological composition (SanPiN 2.1.7.2790-10): class B.

10. **Quality Guarantee**

The guarantee period is 90 days from the date of sale. The manufacturer -BETTINI S.p.A., Italy / BETTINI S.p.A, MONTEMARENZO (LC), Via Industriale, 11 CAP 23804, Italy – is liable for quality.

11. **Manufacturer**

BETTINI S.p.A., Italy / BETTINI S.p.A, MONTEMARENZO (LC), Via Industriale, 11 CAP 23804, Italy

Mr Edoardo Ferraris
General Direction

Annex to the Instructions for the Use of Medical Product
Zirconium Dioxide Blanks for Denture Frameworks

Discs of zirconia (zirconium dioxide), DB+:

- disc, diameter 98 mm, thickness 10 mm;
- disc, diameter 98 mm, thickness 12 mm;
- disc, diameter 98 mm, thickness 14 mm;
- disc, diameter 98 mm, thickness 16 mm;
- disc, diameter 98 mm, thickness 18 mm;
- disc, diameter 98 mm, thickness 20 mm;
- disc, diameter 98 mm, thickness 22 mm;
- disc, diameter 98 mm, thickness 25 mm.

II. Cylinders of zirconia (zirconium dioxide), DB+:

- cylinder, diameter 16 mm, height 16 mm;
- cylinder, diameter 20 mm, height 20 mm.

III. Blocks of zirconia (zirconium dioxide), DB+:


- block, length 40 mm, width 15 mm, height 14 mm;
- block, length 40 mm, width 20 mm, height 16 mm;
- block, length 42 mm, width 20 mm, height 16 mm;
- block, length 42 mm, width 20 mm, height 20 mm;
- block, length 43 mm, width 25 mm, height 16 mm;
- block, length 43 mm, width 25 mm, height 20 mm;
- block, length 62 mm, width 25 mm, height 20 mm;
- block, length 62 mm, width 25 mm, height 22 mm;
- block, length 75 mm, width 36 mm, height 16 mm;
- block, length 75 mm, width 36 mm, height 20 mm.

IV. Discs of zirconia (zirconium dioxide), DB+ LUX:

- disc, diameter 98 mm, thickness 10 mm;
- disc, diameter 98 mm, thickness 12 mm;
- disc, diameter 98 mm, thickness 14 mm;
- disc, diameter 98 mm, thickness 16 mm;
- disc, diameter 98 mm, thickness 18 mm;
- disc, diameter 98 mm, thickness 20 mm;
- disc, diameter 98 mm, thickness 22 mm;
- disc, diameter 98 mm, thickness 25 mm.

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