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
NOTARY PUBLIC

I, the undersigned Gabriel Rudbeck, Notary Public of the City of Stockholm, Sweden, hereby certify,

that this document is a true and correct copy of the document produced to me.

Stockholm, Sweden 2018-06-25

Ex officio:

Gabriel Rudbeck

Sveavägen 24-26, P O Box 3098, S-103 61 Stockholm, +46 (0) 8 20 59 90

Dermal implant **Sculptra™** (poly-L-lactic acid) manufactured by Q-Med AB, Sweden

DEVICE DESCRIPTION

SCULPTRA™ is a poly-L-lactic acid implant in the form of a sterile non-pyrogenic suspension, which is reconstituted from a sterile dry powder by the addition of sterile water for injections. This suspension contains microparticles of poly-L-lactic acid, the crystalline form of polylactic acid. Poly- L-lactic acid is a biocompatible, biodegradable, synthetic polymer from the alpha-hydroxy-acid family.

COMPOSITION OF SCULPTRA

Each vial of dry powder contains:
150 mg of poly-L-lactic acid
90 mg of sodium carboxymethylcellulose
127.5 mg of non-pyrogenic mannitol

PRIMARY PACKAGING

Glass Vial:	Clear colourless Type I Ph.Eur./USP glass vial
Rubber Stopper:	Red Ph.Eur./USP Bromobutyl rubber stopper
Aluminium cap provided by plastic flip-off:	Grey Polypropylene and aluminum flip-off cap

HOW SUPPLIED


The **SCULPTRA** dry powder is supplied after aseptic filtration sterilization in a sterile elongated clear glass vial with an aluminium ring at one end, which is hermetically sealed by a rubber bung, covered by a flip-off cap. **SCULPTRA** powder is reconstituted with 5 ml of sterile water for injections to form a sterile non-pyrogenic suspension.

CONFIGURATION

Vial containing **Sculptra™** – 1 piece
Instructions for use – 1 piece
Labeled carton – 1 piece

Guaranteed critical parameters /characteristics of the Device

Critical Parameter / Characteristic	Specification limit
Description	White to off-white lyophilized powder, in a clear glass vial, with a red rubber stopper and with an aluminium crimp with a grey flip-off cap
Reconstituted Product Evaluation	Formation of a well-dispensed suspension with some foam
pH (reconstituted solution)	5.9 - 7.2
Average Weight of content (on 20 vials) (mg)	273.0 – 406.3

 GALDERMA	Title Instructions for Use of Medical Device Sculptra	Doc id
		MA-35583

Critical Parameter / Characteristic	Specification limit
Uniformity of mass	Complies
Water content (%)	≤ 1.3 at release ≤ 2.0 at end of shelf life
Sterility	Sterile

INDICATIONS FOR USE

SCULPTRA is suitable for increasing the volume of depressed areas, particularly to correct skin depressions, such as in skin creases, wrinkles, folds, scars, and for skin aging.

SCULPTRA is also suitable for large volume corrections of the signs of facial fat loss (lipoatrophy).

Injection technique: the depth of injection and quantity of **SCULPTRA** used depend on the area to be treated and the result expected. Over-corrections should be avoided, but if they occur, the area concerned should be thoroughly massaged to ensure proper distribution of the product. See section "INSTRUCTIONS FOR USE" for additional information.

MODE OF ACTION

SCULPTRA is injected into the deep dermis or subcutaneous layer with a 26 G needle. The tight granulometric distribution of the microparticles of poly-L-lactic acid, its slow degradation kinetics, and a viscosity which is suitable for both deep intradermal or subcutaneous injections, gives **SCULPTRA** its mechanical properties and prolonged resorbability, which make this implant suitable for filling areas of depressed skin.

BIODEGRADATION

Sculptra consists of PLLA particles that slowly will degrade after injection. *In vitro* studies have shown degradation rate between 2% up to 10% after 24 months. The main metabolite lactate will be removed from the tissue by the cell glycolytic pathway.

The clinical experience is that effect of the treatment starts to slowly subside after 24 months. Biopsies from patients shows the presence of PLLA particles 12 months after injection, while no particles were observed in biopsy samples obtained 30 months after last injection; that might be an indication of complete biodegradation of PLLA particles (Vlegaar D & Bauer U 2004. J. Drugs Dermatol. 3(5):542-547).

Once injected, the implant is not possible to remove or replace.

CONTRAINDICATIONS

Do not use in the case of acute or chronic skin disease (infection or inflammation) in or near the area to be treated. **SCULPTRA** should not be used by any person who has an allergy to any of the constituents of the product.

Effective date: 2018-06-17

Effective

Version: 2.0

 GALDERMA	Title Instructions for Use of Medical Device Sculptra	Doc id MA-35583
--	--	--------------------

WARNINGS

SCULPTRA should only be used in the deep dermis or subcutaneous layer.

Avoid superficial injections in order to avoid the appearance of early papules or nodules at the injection site, which could be suggestive of improper injection techniques (superficial placement, excessive amount of product, incorrect reconstitution). In addition, massaging the treatment area to ensure proper distribution of the product may also minimize the appearance of papules or nodules.

SCULPTRA vials are for single patient and single session use only in order to avoid contamination. Do not re-use the vial or do not re-sterilize the vial.

Discard immediately after use. Do not use if package or vial is opened or damaged.

Caution must be taken to avoid injection into the blood vessels. An introduction into the vasculature may occlude the vessels and could cause skin infarction or embolism.

Do not overcorrect (overfill), a contour deficiency because the depression should gradually improve within several weeks as the treatment effect of **SCULPTRA** occurs.

Do not inject into the red area of the lip (vermillion).

Always mix the powder with sterile water for injections.

For injection of **SCULPTRA**, use the 26 G needles with single-use sterile syringes.

PRECAUTIONS FOR USE

The following precautions should be observed:

SCULPTRA should only be used by health care providers with expertise in the correction of volume defects after fully familiarizing themselves with the product, the product educational materials and its complete instruction leaflet.

The injection site should be cleaned with an antiseptic and free from inflammation or infection. As with all injections, patients treated with anti-coagulants may run the risk of a haematoma or localized bleeding at the injection site.

The safety of **SCULPTRA** for use during pregnancy, in breastfeeding females or in patients under 18 years has not been established.


No studies of interactions of **SCULPTRA** with drugs or other substances or implants have been made.

Caution must be taken when injecting **SCULPTRA** in areas of thin skin, such as the periorbital area.

An increased risk of papules and nodules in the periorbital area has been reported (see section “SIDE EFFECTS”). Refer to the product educational materials for instructions regarding injection techniques.

The safety of using **SCULPTRA** in patients with susceptibility to keloid formation and hypertrophic scarring has not been established. **SCULPTRA** should not be used in patients with known history of or susceptibility to keloid formation or hypertrophic scarring.

If laser treatment, chemical peeling or any other procedure based on active dermal response is considered after treatment with **SCULPTRA**, there is a possible risk of eliciting an inflammatory

reaction at the implant site. This also applies if **SCULPTRA** is administered before the skin has healed completely after such a procedure.

The patient should be informed that he or she should minimize exposure of the treatment area to sun and avoid UV lamp exposure until any initial swelling and redness has resolved

SIDE EFFECTS OF THE TREATMENT

The side effects usually resulting from the injections are transient bleeding from an area the size of the point of the needle or transient pain, localised redness at the injection site, ecchymosis, hematoma, or even slight oedema, which generally resolve within 2-6 days.

Based on data obtained through post-marketing surveillance and clinical studies, nodules have also been reported.

Subcutaneous papules, invisible but palpable, or visible nodules including periorbital nodules, or areas of induration have been noted in the injection area and may be due to over-correction.

Nodules are occasionally associated with inflammation or discoloration.

The early occurrence of subcutaneous nodules at the injection site (within 3 to 6 weeks after treatment) may be minimized by adhering to proper dilution and injection techniques (see section "WARNINGS").

Delayed occurrences of subcutaneous nodules at the injection site (within 1 to 14 months post injection) have been reported with sometimes a prolonged duration of up to 2 years.

For nodular areas or late granuloma formation, in some cases, they resolved spontaneously or following treatment with multiple intralesional injections of corticosteroids or other such agents.

Surgical excision of the nodules was sometimes required when they were larger in size, occurring in difficult anatomical regions (e.g. lower eyelid) or persisting after other treatments.

Other rarely reported adverse events include injection site abscess, injection site infection including cellulitis (facial), staphylococcal infection, late granuloma formation, allergic reaction, injection site urticaria, skin hypertrophy and injection site atrophy, angioedema, telangiectasis, skin sarcoidosis, scarring, skin discoloration.

ANY SIDE EFFECTS OR PRODUCT COMPLAINTS SHOULD BE NOTIFIED TO YOUR LOCAL GALDERMA REPRESENTATIVE OR DISTRIBUTOR.


INSTRUCTIONS FOR USE

Refer to the product educational materials for instructions regarding injection techniques.

Reconstitution prior to use

SCULPTRA is reconstituted extemporaneously in the following way:

1. Remove the flip-off cap from the vial and clean the penetrable stopper of the vial with an antiseptic. If the vial seal, or flip-off cap are damaged, do not use, and contact your local Galderma representative or distributor.

Version: 2.0

Effective

Effective date: 2018-06-19 09:40

	Title Instructions for Use of Medical Device Sculptra	Doc id MA-35583
--	---	---------------------------

2. Take a sterile single-use 5 ml syringe (graduated in 1/10 ml) and an 18 G sterile needle and attach the needle to the syringe.
3. Draw 5 ml of sterile water for injections with the syringe.
4. Slowly add the 5 ml of water to the dry powder, by introducing the 18 G needle into the rubber bung of the vial by a single puncture.

5. Let it stand for at least 2 hours to ensure complete hydration; do not shake during this period. **SCULPTRA** can be stored at refrigeration (2-8°C) or room temperature up to 30°C during hydration.

6. When ready to inject the product, product should be gently agitated immediately prior to use. Agitate the vial until a homogeneous translucent suspension is obtained. A single vial swirling agitator may be used. It is ready for use.

The reconstituted product must be injected within 72 hours following reconstitution.

If not used within 72 hours, it must be discarded (See "SPECIAL STORAGE CONDITIONS").

7. Immediately prior to use, clean the rubber bung of the vial with an antiseptic, and use a second 18 G needle to withdraw the appropriate amount of the suspension into a single-use sterile syringe. Do not store reconstituted product in the syringe.

8. Before injecting the product, change the 18 G needle to the 26 G needle, to inject the product into the deep dermis or subcutaneous layer.

9. To withdraw remaining contents of the vial, repeat steps 6 through 8.

10. Discard immediately after single session/patient use.

Patient treatment


1. A complete medical history should be taken to determine if the treatment is appropriate. Before treatment with **SCULPTRA**, the patient should be informed completely of the indications, contraindications, warnings, precautions for use, possible side effects and mode of administration of **SCULPTRA**. Each patient should be informed that the amount of **SCULPTRA** and the number of injection sessions will depend on the patient's need and the severity of the depressed area.

2. Respect the rules of asepsis and hygiene. Clean the injection site with an antiseptic.

3. **SCULPTRA** powder should be reconstituted extemporaneously with 5 ml of sterile water for injections, prior to subcutaneous or intradermal injection. **SCULPTRA** should be injected into the deep dermis or subcutaneous layer.

Do not inject with needles smaller than 26 G. Do not bend the needle.

To maintain a uniform suspension throughout the procedure, intermittently agitate the product in the syringe. Before initial injection, expel a few drops of the product through the attached 26 G needle to eliminate air and to check for needle blockage. If the 26 G needle becomes occluded or dull during an injection session, replacement may be necessary. If clogging occurs, remove the needle, expel a small amount of product, attach a new sterile 26 G needle, then expel a few drops of **SCULPTRA** to eliminate the air and re-check for needle blockage.

Effective date: 2018-06-19 09:40

Effective

	Title Instructions for Use of Medical Device Sculptra	Doc id MA-35583
--	--	--------------------

In order to control the injection depth of **SCULPTRA**, stretch/pull the skin opposite to the direction of the injection to create a firm injection surface. The 26 G sterile needle, bevel up, should be introduced into the skin at an angle of approximately 30-40 degrees, until the desired skin depth is reached.

A change in tissue resistance is felt when the needle crosses from the dermis into subcutaneous layer. If the needle is inserted at too shallow (small) an angle or if the needle tip is not sufficiently advanced, then the needle tip may be in the mid or superficial papillary dermis, the needle bevel may be visible through the skin. If product is injected too superficially the injected area will blanch immediately or shortly after injection. If this occurs, the needle should be removed and the treatment area gently massaged. In the event that the blanching does not disappear, the patient should not be re-injected.

4. During the first treatment session with **SCULPTRA**, only a limited correction should be made. The patient should then be evaluated at no sooner than four weeks posttreatment to determine if additional correction is needed. The patient should be informed of the potential need for additional treatments at the first consultation.

5. After the injection session, ice (in a suitable cloth, avoiding any direct contact with the skin) should be applied to the treatment area in order to reduce swelling and/or bruising.

6. It is important to thoroughly massage the treated area(s) to evenly distribute the product (use of an appropriate cream may help to reduce the friction on the skin surface during massaging).

7. The patient should periodically massage the treated areas for five minutes, five times per day for five days after the treatment to ensure a natural-looking correction.

8. **SCULPTRA** may be visualized with ultrasound imaging and MRI. It is not observed with CT scans and radiography.

SPECIAL STORAGE CONDITIONS

SCULPTRA powder should be stored at room temperature protected from humidity (store below 30°C). Upon reconstitution, **SCULPTRA** can be stored up to 72 hours at refrigeration (2-8°C) or room temperature up to 30°C. Do not freeze.

STERILIZATION

The sterility of the finished product is obtained through the aseptic manufacturing process.

IF THE VIAL, SEAL, OR FLIP-OFF CAP ARE DAMAGED, DO NOT USE, AND CONTACT YOUR LOCAL GALDERMA REPRESENTATIVE OR DISTRIBUTOR.

The shelf-life of Sculptra™ is 36 months.

	Title Instructions for Use of Medical Device Sculptra	Doc id MA-35583
--	---	---------------------------

DISPOSAL

After use, treatment syringes and needles may be potential biohazards. Handle accordingly and dispose of in accordance with accepted medical practice and applicable local requirements.

Dispose as special industrial waste at an authorized facility.

Disposal in the territory of the Russian Federation is carried out in accordance with SanPiN 2.1.7.2790-10 "Collection rules; storage and disposal of waste treatment and prevention facilities".

The class of danger of medical waste – B.

LEGAL MANUFACTURER:

Q-Med AB

Seminariegatan 21

SE-752 28 Uppsala

Sweden

For any information about this product, please contact your local representative:

Limited Liability Company GALDERMA, 31 A/1 Leningradskiy Avenue, Moscow

125284 Russian Federation

PACKAGING

The vial together with instructions for use is placed in a labeled carton.

Symbols on packaging


Caution

Consult instructions for use

Single patient use

Storage conditions

Sterile using aseptic processing techniques

 GALDERMA	Title Instructions for Use of Medical Device Sculptra	Doc id MA-35583
--	--	--------------------

LIST OF RELEVANT STANDARDS

Standard	Title
ISO 14644-1	Cleanrooms and associated controlled environments
ISO 10993-1	Biological evaluation of medical devices - Part 1: Evaluation and testing of within a risk management process
ISO 13408	Aseptic process of health care products: - Part 1: General requirements - Part 2: Filtration - Part 3: Lyophilization - Part 4: Clean-in-place technologies - Part 5: Sterilization in place - Part 6: Isolators
ISO 14155	Clinical investigation of medical devices for human subjects - Good clinical practice
ISO 14630	Non-active surgical implants - General requirements
ISO 15223-1	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
SS-EN 556-2	Sterilization of medical devices - Requirements for medical devices to be designated "STERILE" - Part 2: Requirements for aseptically processed medical devices
SS-EN 1041 +A1	Information supplied by the manufacturer of medical devices
SS-EN ISO 13485	Medical devices - Quality management systems - Requirements for regulatory purposes
SS-EN ISO 14971	Medical devices - Application of risk management to medical devices

WARRANTY LIABILITIES

Technical advice of Q-Med AB, both orally and in writing, is designed to help physicians use products of the company. Such advice does not extend limited warranty and does not eliminate the need for a physician to check products of Q-Med AB in order to ensure that they are suitable for the required procedures. A physician assumes all risks and liability for damage that may occur as a result of misuse of products of Q-Med AB. In the event of a defect in material or production, the liability of Q-Med AB is limited, at the discretion of Q-Med AB, to a replacement of a defective product or a refund of the actual cost of the defective product. To implement this limited warranty, the defective product must be returned to Q-Med AB. Q-Med AB is not liable under any circumstances for any indirect, incidental or consequential damages.

Q-Med AB declares compliance of medical device with the listed above international and national standards and guarantees that the medical device shall conform to the specifications during the whole shelf life on a condition that the transportation, storage and use are according to the conditions established by the manufacturer and described in instruction for use.

Q-Med AB makes no other warranties, either explicit or implied, except for the ones specified above.

SIGNATURES PAGE

Date	Signed by
2018-06-19 08:37	LUNDGREN Björn
Justification	Approved by Technical Expert

2018-06-19 08:45	KOPTINA Anna (EXT)
Justification	Compiled by

2018-06-19 09:40	NILSSON Lilian
Justification	Approved by Owner

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