INDIVIDUAL ACKNOWLEDGMENT

State/Commonwealth of Minnesota }
County of Hennepin } ss.

On this the 24th day of May, 2017, before me,
Deborah L. Shaffer, the undersigned Notary Public,
Name of Notary Public
personally appeared Julie Rivard,
Name(s) of Signer(s)

☒ personally known to me – OR –

☐ proved to me on the basis of satisfactory evidence

to be the person(s) whose name(s) is/are subscribed to the within instrument, and acknowledged to me that he/she/they executed the same for the purposes therein stated.

WITNESS my hand and official seal.

Place Notary Seal and/or Stamp Above

[Signature]
Signature of Notary Public
Plymouth, MN
Other Required Information (Printed Name of Notary, Residence, etc.)

OPTIONAL

Although the information in this section is not required by law, it may prove valuable to persons relying on the document and could prevent fraudulent removal and reattachment of this form to another document.

Description of Any Attached Document

Title or Type of Document: Operating Document-IFU for EverFlex with Entrust

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4600 Nathan Lane North
Plymouth, MN 55442, USA
Tel: +1 763 398 7000
www.ev3.net

М.П. / Stamp

"Operating document – Instruction for use" for
submission in the Russian Federation

**EverFlex Self-expanding Peripheral Stent with
Entrust Delivery System**

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EverFlex™

Self-expanding Peripheral Stent with
Entrust™ Delivery System

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INSTRUCTIONS FOR USE

DEVICE DESCRIPTION

The EverFlex Self-expanding Peripheral Stent with Entrust Delivery System is a self-expanding Nitinol stent system intended for permanent implantation. The self-expanding stent is made of a nickel titanium alloy (Nitinol) and comes pre-mounted on a 5 F, 0.035" over-the-wire delivery system. The stent is cut from a Nitinol tube that has an open lattice design, and has tantalum radiopaque markers at the proximal and distal ends of the stent. Upon deployment, the stent achieves its predetermined diameter and exerts a constant, gentle outward force to establish patency.

The Entrust delivery system is shown in Figure 1. It has a triaxial shaft design comprised of an inner shaft (1), retractable outer sheath (2), and an isolation sheath (3) for ease of use during deployment. The ergonomic deployment handle (8) has a removable safety lock (9), deployment wheel (10), and luer hub (6).

The distal portion of the delivery system (enlarged section) contains two radiopaque markers; one marker distal (4) and one marker proximal (5) to the constrained stent. The distal marker is mounted on the retractable outer sheath and the proximal marker is mounted on the inner shaft. The EverFlex self-expanding stent is constrained within the outer sheath. The space is flushed prior to the procedure through the luer hub. The outer sheath terminates distally in a flexible catheter tip (7).

Stent positioning at the target lesion is achieved prior to deployment utilizing the two radiopaque markers and the radiopaque markers on the stent. For stent deployment, the safety lock is removed from the handle to unlock the deployment wheel and outer sheath. The outer sheath retraction is achieved by slowly rotating the deployment wheel in the direction indicated by the arrow on the deployment handle. Complete deployment of the stent is achieved when the radiopaque marker on the outer sheath passes the proximal radiopaque marker on the inner shaft.


Figure 1: EverFlex Entrust™ Delivery System

INDICATION FOR USE

Peripheral

The stent is indicated for use in occlusions, lesions at high risk for abrupt closure or threatened closure following percutaneous transluminal angioplasty (PTA), or lesions believed to be at high risk for restenosis following PTA in the common iliac, external iliac, superficial femoral, proximal popliteal or subclavian arteries. Stenting is intended to improve and maintain artery lumen diameter.

CONTRAINDICATIONS

- Patients in whom anticoagulant and/or antiplatelet therapy is contraindicated.
- Patients who exhibit persistent acute intraluminal thrombus of the proposed lesion site, post thrombolytic therapy.
- Perforation at the angioplasty site evidenced by extravasation of contrast medium.
- Aneurysm of the artery to be treated.
- Patients with known hypersensitivity to nickel titanium.
- Patients with uncorrected bleeding disorders.
- All of the customary contraindications for PTA.
- The EverFlex Self-expanding Peripheral Stent with Entrust Delivery System is contraindicated for use in the Carotid artery.

WARNINGS

- The device is provided STERILE for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risk of patient infection and risk of compromised device performance.
- If abnormally high resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may cause damage to the stent or vessel, pre-deployment, or result in deployment complications. Carefully withdraw the stent system without deploying the stent.
- If high resistance is felt when rotating the thumbwheel, do not force deployment. Carefully withdraw the stent system without deploying the stent.

PRECAUTIONS

- Carefully inspect the sterile package and device prior to use to verify that no damage occurred during shipment.
- Do not exceed 300 psi / 20 ATM while flushing the delivery system.
- Do not use if the stent is partially deployed upon removal from the package, or before starting the deployment procedure.
- Always use a sheath during the implant procedure to protect both the vessel and puncture site. Support from a sheath is necessary to minimize lengthening or shortening during stent deployment.
- Using a 0.014" or 0.018" guidewire in a highly tortuous vessel, highly stenosed or calcified lesion could result in crossing, tracking, or deployment complications.
- Crossing and/or deployment complications may result from a lesion lumen diameter less than 2 mm.
- Torsing the device (multiple rotations) may lead to difficulty or inability to deploy and/or move device over the guidewire.
- Any slack in the system, inside or outside of the patient, could result in deploying the stent beyond the target site.
- Failure to completely remove the locking pin and tube assembly can result in the inability to deploy, increased deployment force, or partial deployment.
- The stent system is not designed for recapturing or repositioning after establishing vessel apposition.
- Failure to hold the handle in a fixed position may result in foreshortening or lengthening of the stent.
- The stent is not designed to be lengthened or shortened past its nominal length. Excessive stent lengthening or shortening may increase the risk of stent fracture.
- Use caution when crossing a deployed stent with any adjunct device.
- Stent should not be expanded past its nominal diameter.

POTENTIAL ADVERSE EVENTS

The potential adverse effects (e.g., complications) that may occur and/or require intervention with the use of this device include, but are not limited to:

- Abrupt or sub-acute closure
- Allergic reaction to device materials or procedure medications
- Allergic reaction to Nitinol
- Amputation
- Aneurysm
- Angina
- Arrhythmia
- Arterio-venous fistula
- Artery perforation or rupture
- Bleeding requiring transfusion
- Bruising
- Contrast medium reaction/renal failure
- Death
- Device breakage
- Dissection or intimal flap
- Edema
- Embolism
- Failure to deploy stent
- Fever
- Gastrointestinal bleeding due to anticoagulation
- Hematoma
- Hypertension/hypotension
- Infection
- Inflammation
- Intraluminal thrombus
- Myocardial infarction
- Pain
- Partial stent deployment
- Pseudoaneurysm
- Renal failure requiring dialysis
- Renal insufficiency (new or worsening)
- Restenosis
- Sepsis
- Shock
- Stent collapse or fracture
- Stent migration
- Stent misplacement
- Stroke
- Surgical or endovascular intervention
- Thrombosis/occlusion of the stent
- Transient ischemic attack
- Venous thromboembolism
- Vessel spasm

PROCEDURE

PREPARATION PROCEDURES

WARNING: The device is provided STERILE for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risk of patient infection and risk of compromised device performance.

1. Required Items for Implantation Procedure

- 5-10 cc syringe filled with heparinized saline
- 0.035" guidewire
- Hemostatic sheath
- PTA balloon

2. Select Stent Size

Measure the diameter of the reference vessel (proximal and distal to lesion). Refer to the table below for stent diameter sizing. Measure the length of the target lesion. Choose a stent length that will extend proximal and distal to the target lesion.

Stent Diameter and Length Sizing

Stent Diameter (mm)	Reference Vessel Diameter (mm)	Stent Length (mm)
5	3.5 – 4.5	20, 40, 60, 80, 100, 120, 150
6	4.5 – 5.5	20, 40, 60, 80, 100, 120, 150
7	5.5 – 6.5	20, 40, 60, 80, 100, 120, 150
8	6.5 – 7.5	20, 40, 60, 80, 100, 120, 150

3. Preparation of Stent Delivery System

- Open the shelf box to reveal the pouch containing the stent and delivery catheter.
- After careful inspection of the pouch, looking for damage to the sterile barrier, carefully peel open the outer pouch and extract the tray with contents.
- Set the tray on a flat surface. Carefully remove the stent delivery system.

CAUTION: Carefully inspect the sterile package and device prior to use to verify that no damage occurred during shipment. Verify that the red locking pin is still engaged with the handle.

CAUTION: Do not exceed 300 psi / 20 ATM while flushing the delivery system.

- Attach a 5-10 cc syringe filled with heparinized saline to the proximal luer lock injection hub. Inject the saline solution through the guidewire lumen until it comes out the catheter tip.

- Examine the distal end of the catheter to ensure the stent is fully contained within the tip of the catheter.

CAUTION: Do not use if the stent is partially deployed upon removal from the package, or before starting the deployment procedure.

STENT DEPLOYMENT PROCEDURE

1. Insertion of Sheath and Guidewire

- Gain access at the appropriate site using a sheath with a hemostatic valve that is compatible with a 5 F delivery system. If the iliac artery is treated, the access site is the common femoral artery. If the subclavian artery is treated, the access site is the brachial or axillary artery.

CAUTION: Always use a sheath during the implant procedure to protect both the vessel and puncture site. Support from a sheath is necessary to minimize lengthening or shortening during stent deployment.

- Insert a 0.035" guidewire of appropriate length across the target lesion via the sheath.

CAUTION: Using a 0.014" or 0.018" guidewire in a highly tortuous vessel, highly stenosed or calcified lesion could result in crossing, tracking, or deployment complications.

2. Dilate of Lesion

Pre-dilate the lesion using standard PTA techniques. Remove the PTA balloon from the patient while maintaining lesion access with the guidewire.

CAUTION: Crossing and/or deployment complications may result from a lesion lumen diameter less than 2 mm.

3. Introduction of Stent Delivery System

Advance the device over the guidewire through the hemostatic valve and sheath.

WARNING: If abnormally high resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may cause damage to the stent or vessel, pre-deployment, or result in deployment complications. Carefully withdraw the stent system without deploying the stent.


CAUTION: Torquing the device (multiple rotations) may lead to difficulty or inability to deploy and/or move device over the guidewire.

4. Stent Deployment


- Advance the delivery system until the distal (loading) radiopaque marker is distal to the target lesion.

- Pull back on the delivery system until there is no slack in the delivery system and the radiopaque markers extend distal and proximal to the target lesion.

CAUTION: Any slack in the system, inside or outside of the patient, could result in deploying the stent beyond the target site.

- Completely remove the locking pin and tube assembly by gently pulling in the direction indicated by the arrow on the pin. Discard the locking pin and attached tube assembly.

CAUTION: Failure to completely remove the locking pin and tube assembly can result in the inability to deploy, increased deployment force, or partial deployment.

- Visually confirm that the radiopaque markers still extend distal and proximal to the target lesion. Reposition if necessary.
- Initiate stent deployment by pinning down (holding) the handle in a fixed position and slowly rotating the thumbwheel in the direction indicated by the arrow on the handle (Figure 1).

NOTE: If a second hand is used to hold the delivery system, gently support the delivery system at the gold isolation sheath. Do not constrict the silver outer sheath during deployment.

CAUTION: The stent system is not designed for recapturing or repositioning after establishing vessel apposition.

- Once the stent distal radiopaque markers are approximately 1 mm to 3 mm distal to the outer sheath marker, and prior to achieving vessel apposition, reposition stent as needed using stent radiopaque markers.
- During release of the stent, the whole length of the flexible deployment system should be kept as straight as possible. In order to ensure that no slack is introduced to the delivery system, hold the handle stationary and fixed. Deployment is complete when the outer sheath marker reaches the inner shaft marker and the stent is released.

WARNING: If high resistance is felt when initially rotating the thumbwheel, do not force deployment. Carefully withdraw the stent system without deploying the stent.

CAUTION: Failure to hold the handle in a fixed position may result in foreshortening or lengthening of the stent.

CAUTION: The stent is not designed to be lengthened or shortened past its nominal length. Excessive stent lengthening or shortening may increase the risk of stent fracture.

NOTE: If more than one stent is needed, place the distal stent first. If overlap of sequential stents is necessary, the amount of overlap should be kept to a minimum.

5. Post Stent Deployment

- While using fluoroscopy following stent deployment, withdraw the entire delivery system as one unit, over the guidewire, into the catheter sheath and out of the body. Remove the delivery system from the guidewire.
- Using fluoroscopy, visualize the stent to verify full deployment.
- If incomplete expansion exists within the stent at any point along the lesion, post deployment balloon dilation may be performed.
- CAUTION: Use caution when crossing a deployed stent with any adjunct device.**
- CAUTION: Stent should not be expanded past its nominal diameter.**
- To dilate the stent, select an appropriate size PTA balloon catheter and dilate with conventional technique. The inflation diameter of the PTA balloon should approximate the diameter of the reference vessel.
- Confirm full stent expansion is complete, then remove the PTA balloon from the patient.
- Remove the guidewire and sheath from the body.
- Close entry wound as appropriate.
- Discard the delivery system, guidewire and sheath.

MRI INFORMATION

MR CONDITIONAL

Non-clinical testing demonstrated that the EverFlex Self-expanding Peripheral Stent with Entrust Delivery System stent is MR Conditional. A patient may be scanned safely, immediately after stent placement under the following conditions:

- Static magnetic field of 3-Tesla or 1.5-Tesla
- Maximum spatial gradient magnetic field of 720-Gauss/cm or less
- Normal operating mode (maximum WBA-SAR of 2.0 W/kg) for 15 minutes of scanning

MRI-RELATED TEMPERATURE RISE

In non-clinical testing, the EverFlex Self-expanding Peripheral Stent with Entrust Delivery System stent produced the following temperature rises during MRI performed for 15-minutes of scanning (i.e., per pulse sequence) in 1.5-Tesla/64-MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numana4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner) and 3-Tesla/128-MHz (Excite, HDx, Software 14X M5, General Electric Healthcare, Milwaukee, WI) MR systems:

Single, 8mm x 80mm	1.5-T	3.0-T
MR system reported, whole body averaged SAR (W/kg)	2.9	2.9
Calorimetry measured values, whole body averaged SAR (W/kg)	2.1	2.7
Highest temperature rise (°C)	2.0	2.4

Single, 8mm x 120mm	1.5-T	3.0-T
MR system reported, whole body averaged SAR (W/kg)	2.9	2.9
Calorimetry measured values, whole body averaged SAR (W/kg)	2.1	2.7
Highest temperature rise (°C)	2.9	3.7

These temperature changes will not pose a hazard to a patient under the conditions indicated above. It is recommended that patients register conditions under which the implant may be scanned safely with the MedAlert Foundation (www.medicalert.org) or equivalent organization.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the stent. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

WARRANTY DISCLAIMER

Although this product has been manufactured under carefully controlled conditions, ev3 Inc. has no control over the conditions under which this product is used. ev3 Inc. therefore disclaims all warranties, both express and implied, with respect to the product including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. ev3 Inc. shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind ev3 Inc. to any representation or warranty with respect to the product. The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

SYMBOL LEGEND

	Manufacturer
	Authorized representative in the European community
	Consult instructions for use
	Sterilized using ethylene oxide
	Catalogue number
	Batch code
	Keep dry
	Keep away from sunlight
	Use by
	Do not reuse
	Do not use if package is damaged
	Telephone
	Facsimile
	For prescription use only

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©2008 Национальная ассоциация нотариусов • 9350 Де Сото Авеню, п/я 2402 • Чатсуорт, Калифорния