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
CERTIFIED

5. at Indianapolis, Indiana
6. this *Thirtieth* of *May*, 2017
7. by the Secretary of State of Indiana
8. No. A2017-0530104106
9. Seal/Stamp:
10. Signature:

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CERTIFICATION STATEMENT

As a responsible representative of Cook Incorporated, I confirm this is an accurate version of the Instructions for Use for Fascial Incising Needles.

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

Fascial Incising Needle (090070, 090070-S5, 090071)


Fascial Incising Needles are used over a wire guide for cutting the renal fascia to facilitate dilation of a musculofascial tract. The Fascial Incising Needle is manufactured from stainless steel with a size of 18G/1.25mm and a length of 5 or 10 centimeters. The needle is designed to accept instrumentations up to 0.038 inch/0.97mm diameters, and is designed with a cut width of 4.5 or 10 millimeters.

Product Catalogue Number	Needle		Cut	
	Gage (Ga)/mm	Length (cm)	Length (cm)	Width (mm)
Tolerance	See Technical Characterists above	±10%	-	-
090070	18/1,25	5 cm	1	4,5
090070-S5	18/1,25	10 cm	1	4,5
090071	18/1,25	5 cm	1	10

Parameters of needle hub:

- “A” – length – 0.608 inch/ 15.443 (±0.005 inch/0.127mm)
- “B” – outer diameter – 0.220 inch /5.588mm (±0.005 inch/0.127mm)
- “C” – inner diameter – 0.151 inch/ 3.836mm (±0.001 inch/0.0254mm)
- “D” – diameter – 0.041 inch/1.041 mm (±0.005 inch/0.127mm)

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
PACKAGE

These devices is supplied, comprising:
Inserts – 2 pc;

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INTENDED USE/INDICATION FOR USE

Fascial Incising Needle (090070, 090070-S5, 090071) -- Used over a wire guide for cutting the renal fascia to facilitate dilation of a musculofascial tract.

Amplatz Needle Holder with silicone inserts (090000), Amplatz Needle Holder Accessory Inserts (091000) -- Used to direct a standard disposable needle during needle access procedure. It prevents the stylet from backing out of the cannula during puncture and also keeps the hands from radiation exposure during the puncture procedure.

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INSTRUCTIONS FOR USE

Fascial Incising Needle (090070, 090070-S5, 090071):

1. Make an incision using the fascial incision needle in the muscle or fascia tracking a 0.038 inch guidewire.

Warning: Do not use a Teflon-coated guidewire for tracking, as the needle blade will detach Teflon from the surface of the guidewire. To properly preserve the fascial incision needle,

rotate the fascial incision needle by 90 degrees after the initial incision and make a further incision.

- 2. Advance the incision using the guidewire, maintaining a fixed angle.
 - 3. Remove the incision needle and begin expansion of the muscle or fascial tissue.
- Warning:** Expansion of the muscle or fascial tissue may be difficult due to resistance from fascial tissue which is tough due to cross fibers.

Amplatz Needle Holder with silicone inserts (090000), Amplatz Needle Holder Accessory Inserts (091000)

- 1. Secure the Cook trocar needle in the silicone insert within the clear plastic needle handle.
- 2. Position your hands to the outside of the device to ensure minimal radiation exposure during the procedure.

TABLE OF MATERIALS

Fascial Incising Needle (090070, 090070-S5, 090071)

Device	Component	Material	Material Specification/ Compound Number	Patient Contacting
090071, 090070 and 090070-S5				
Needle	Needle Hub	Polycarbonate	Compound # 80020 CI Solvent Red 207	Indirect Contact
	Needle Cannula	Stainless Steel	AISI 304	Direct Contact
	Blade	Stainless Steel	AISI 304	Direct Contact
Product Protector	Product Protector	Polypropylene Non-Radiopaque Tubing	#30070	Non-Contact

Amplatz Needle Holder with silicone inserts (090000), Amplatz Needle Holder Accessory Inserts (091000)

RPNs Included	Device Components	Material	Material Specification #/ Compound #	Patient Contacting
090000	Needle Holder	Polycarbonate	Polycarbonate CH-10	Non Contact
	Cap	Polycarbonate Rod	Polycarbonate CH-10	Non Contact
	Button	Polycarbonate Rod	Polycarbonate CH-10	Non Contact
091000	Silicone Insert	Silicone	Silicone SILASTIC Q7-4735	Non Contact
	Silicone Insert	Silicone	Silicone, Blue SILASTIC Q7-4735 with blue pigment	Non Contact

		A35001BLU1	
--	--	------------	--

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- 1. Disassemble the Amplatz Needle Holder by removing the silicone inserts from the device.
- 2. Activate CIDEX® Solution per the manufacturer's instructions. Test the activated solution prior to use with CIDEX® Solution Test Strips.
- 3. Immerse the disassembled device completely, filling all hard-to-reach areas and eliminating air pockets, in CIDEX® Activated Dialdehyde Solution for a minimum of 45 minutes at 25°C.
- 4. Rinse thoroughly in a room temperature Sterile Water for Injection (SWFI) water bath. Rinse by immersing the device completely for a minimum of one minute. Manually flush all hard-to-reach areas with SWFI water.
- 5. Agitate the device under SWFI water, bring above water level, then re-immersed.
- 6. Repeat steps 4 and 5 two additional times using fresh SWFI for each of the three rinses.

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

UPN REF G14473

2123456

Get registered
ASCE member
Member ID: 1742-133
NSD, USA

COOK MEDICAL
Rushmore Medical Park
Birmingham, UK
www.cookmedical.com

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For Russian Federation – in accordance with national requirements. Risk class of the medical waste: class B

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Handwritten signature of T.A. Pestretsova in blue ink.

Handwritten signature of the notary in blue ink.

SECRETARY OF STATE
STATE OF INDIANA

Connie Lawson
Secretary of State


APOSTILLE

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1. Country: United States of America
2. This Public Document has been signed *Jamie L Ridner*
3. Acting in the capacity of notary public in & for *Monroe County*
4. and bears the seal/stamp of certification clerk in & for Sec. of State Office

CERTIFIED

5. at Indianapolis, Indiana
6. this *Thirtieth* of *May*, 2017
7. by the Secretary of State of Indiana
8. No. A2017-0530104103
9. Seal/Stamp:
10. Signature:

Connie Lawson

Connie Lawson
Indiana Secretary of State

*Effective May 1, 2003 all apostilles from the Indiana Secretary of State will have an electronically printed seal.
This document was revised March 17, 2012.*

COOK INCORPORATED
750 DANIELS WAY
BLOOMINGTON, IN 47404 USA
PHONE: 812.339.2235 TOLL FREE: 800.457.4500
WWW.COOKMEDICAL.COM

CERTIFICATION STATEMENT

As a responsible representative of Cook Incorporated, I confirm this is an accurate version of the Instructions for Use for Disposable Two-Part Trocar Needle.

Kris Weathers
Regulatory Affairs Manager
Cook Incorporated


State of Indiana)

ss

County of Monroe)

Subscribed and sworn to before me the 26 day of May month 2017 year.

My Commission expires August 6, 2021

Jamie L. Ridner

INSTRUCTIONS FOR USE

Puncture needles in separate packages, with accessories

INTENDED USE

Used for the initial percutaneous incision to allow access to the kidney and guide wire access
(Please see Intended Use for more information)

DEVICE DESCRIPTION



Conversion Table

Parameter	Equivalent
1 Fr	1/3 mm
1 inch	25,4 mm


TECHNICAL CHARACTERISTICS:

Tolerance for needles: 10 G /3.4 mm (± 0.1 mm); 11 G /3.05 mm (± 0.1 mm); 14 G /2.05 mm (± 0.1 mm); 15 G /1.83 mm (± 0.1 mm); 16 G /1.65 mm (± 0.4 mm); 18 G /1.25 mm (± 0.05 mm); 19 G /1.07 mm (± 0.03 mm); 20 G /0.89 mm (± 0.03 mm); 21 G /0.82 mm (± 0.01 mm); 22 G /0.71 mm (± 0.02 mm); 23 G /0.64 mm (± 0.03 mm); 24 G /0.57 mm (± 0.01 mm); 25 G / 0.52 mm (± 0.01 mm)

Disposable Two-Part Trocar Needle


The Two-Part Trocar Needle is used for percutaneous nephrostomy. The Two-Part Trocar Needle is comprised of a stainless steel beveled wire which can be inserted through a hubbed stainless steel cannula. Some needles are Echotip needles. The Echotip feature consists of a series of very small dimples in the surface, which causes that area of the needle to be easily visible under ultrasound.

Order Number	Cannula		Trocar Stylet		Comments
	Diameter (G/mm)	Length (cm)	Diameter (mm)	Length (cm)	
Tolerance	Refer to technical characteristics	$\pm 10\%$	$\pm 10\%$	$\pm 10\%$	
090020	18/1.25	20	1.016	21.5	-
090020-ET	18/1.25	20	1.016	21.5	EchoTip
090020-S1	18/1.25	15	1.016	16.6	-


Parameters of needle hub:

- “A” – length – 0.608 inch/ 15.443mm (± 0.005 inch/0.127mm)
- “B” – outer diameter – 0.220 inch /5.588mm (± 0.005 inch/0.127mm)
- “C” – inner diameter – 0.151 inch/ 3.836mm (± 0.001 inch/0.0254mm)
- “D” – diameter – 0.041 inch/ 1.041 mm (± 0.005 inch/0.127mm)

Parameters of stylet hub:

- “A” – length – 0.437 inch/ 11.0998 mm (±0.005 inch/0.127mm)
- “B” – outer diameter – 0.166 inch /4.216 mm (±0.003 inch/0.076mm)
- “C” – inner diameter – 0.150 inch/ 3.81mm (±0.001 inch/0.0254mm)


PACKAGE

These devices is supplied, comprising:
needle – 1 pc.;
stylet – 1 pc;
instruction manual – 1 pc.

RPN	Bevel angle degrees	Type of bevel	Blunt-tip radius, mm
Tolerance	±10%	N/A	no more than
090020	22	Trocar	0,03
090020-ET	22	Trocar	0,03
090020-S1	22	Trocar	0,03

Length of cutting edge for each design option- not applicable
Width of cutting edge for each design option- not applicable
Rockwell hardness of needle material – 44...48 HRC
Bending resistance: the needle tube shall not show a deviation of more than 0.45 mm in the range of 25.0 mm and a bending force is 20.0 N
The force of the puncture plastic film with a thickness of 0,150 mm – not more than 2,60 N
Fracture strength: the tube of the needle should not break if the distance between rigid support and point of application of force of bending – 31.5 mm
Product weight 4,0 g (±10%)

Skinny Needle with Chiba Tip (090010, 090010-ET)

The Skinny Needle with Chiba Tip needle is used for positioning of the renal collecting system throughout a percutaneous nephrostomy. The Skinny Needle with Chiba Tip is comprised of a stainless steel stylet wire which can be inserted through a hubbed stainless steel cannula. The cannula comes in 22G/0.71mm sizes with length of 20 centimeters. The distal end of the cannula is a Chiba tip, beveled at a 30° angle with Echotip detail for some needles. The Echotip feature consists of a series of very small dimples in the surface, which causes that area of the needle to be easily visible under ultrasound.

Order Number	Needle		Stylet		Comments
	Diameter G/mm	Length cm	Diameter mm	Length cm	
Tolerance		±10%	±10%	±10%	
090010	22/0,71	20	0.46	21.5	-
090010-ET	22/0,71	20	0.46	21.5	Echotip


Parameters of needle hub:

“A” – length – 0.975 inch/ 24.765 mm (±0.005 inch/0.127mm)

“B” – outer diameter – 0.214 inch /5.436 mm (±0.005 inch/0.127mm)

“C” – inner diameter – 0.141 inch/ 3.581mm (±0.001 inch/0.0254mm)

“D” – diameter – 0.207 inch/ 5.258 mm (±0.005 inch/0.127mm)


Parameters of stylet hub:

“A” – length – 0.439 inch/ 11.151mm (±0.001 inch/0.0254mm)

“B” – outer diameter – 0.332 inch /8.433mm (±0.005 inch/0.127mm)

“C” – inner diameter – 0.144 inch/ 3.658mm (± 0.005 inch/0.127mm)

PACKAGE

These devices is supplied, comprising:

needle – 1 pc.;

stylet – 1 pc;

instruction manual – 1 pc.

RPN	Bevel angle, degrees	Type of bevel	Blunt-tip radius, mm
Tolerance	$\pm 10\%$	N/A	no more than
090010	30	Chiba	0.03
090010-ET	30	Chiba	0.03

Length of cutting edge for each design option- not applicable

Width of cutting edge for each design option- not applicable

Rockwell hardness of needle material – 44...48 HRC

Bending resistance: the needle tube shall not show a deviation of more than 0.45 mm in the range of 25.0 mm and a bending force is 20.0 N


The force of the puncture plastic film with a thickness of 0,150 mm – not more than 2,60 N

Fracture strength: the tube of the needle should not break if the distance between rigid support and point of application of force of bending – 31.5 mm

Product weight 3.1 g ($\pm 10\%$)

Accessories for Puncture needles in separate packages

Amplatz Needle Holder with silicone inserts (090000)

The Amplatz Needle Holder is designed to prevent the end users hands from being directly exposed to radiation during x-ray. The Amplatz Needle Holder is comprised of a needle holder and silicone inserts. The needle holder is manufactured in a length of 30.2 centimeters and a width of 3.8 centimeters. A cap is situated in the center of the needle holder with a button/plunger that should be able to move freely in and out of the cap. There are two thru holes, located to the left and right of the cap, 10.2 centimeters distal from the cap when measuring from the centers of the hole and cap. One clear and one blue silicone insert is provided with the Amplatz Needle Holder. The silicone inserts hold the needle in place while preventing the stylet from backing out of the cannula during introductions.

Product Catalogue Number	Length of Needle Holder (cm)	Width of Needle Holder (cm)	Remarks	Specifications
090000	30.2	3.8	Needle holder with silicone inserts	Blue insert designed for use with Cook Disposable Chiba Needle. Clear insert designed for use with standard disposable Trocar needle.

PACKAGE

These devices is supplied, comprising:

Needle Holder – 1 pc.;

Inserts – 2 pc;

instruction manual – 1 pc.


Amplatz Needle Holder Accessory Inserts (091000)

The Amplatz Needle Holder Accessory Inserts are indicated for use with the Amplatz Needle Holder and Disposable Chiba and Standard Disposable Needles. The inserts are manufactured from translucent and blue silicone. The silicone insert is designed to fit into the preformed hole of the Amplatz Needle Holder. The inserts are round in shape with the center of the insert formed into an eight sided well. The insert is dimensionally tapered over a ½ inch space.

Parameters of silicone insert:

A – 0.500 inch/12.7 mm (± 0.005 inch /0.127 mm); B – 0.343 inch/8.71 mm; C – 0.327 inch/8.31 mm.

Dimensions for the translucent insert: Top; 0.093 inches/2.36mm, Mid; 0.270 inches/6.86mm, Bottom: 0.150 inches/3.81mm

Dimensions for the blue insert: Top; 0.093 inches/2.36mm, Mid; 0.270 inches/6.86mm, Bottom: 0.165 inches/4.19mm

PACKAGE

These devices is supplied, comprising:

Inserts – 2 pc;

The blue insert is designed for use with disposable Chiba needles as the design of the well fits snugly against the stylet knob/needle hub. The translucent insert is designed for use with standard two-part disposable needles as the design of the well fits snugly against the stylet knob/needle hub.

INTENDED USE/INDICATION FOR USE

Disposable Two-Part Trocar Needle (090020, 090020-ET, 090020-S1)- Used for percutaneous nephrostomy.

Skinny Needle with Chiba Tip (090010, 090010-ET) -- Used under fluoroscopy for initial location and positioning of the renal collecting system throughout a percutaneous nephrostomy.

Amplatz Needle Holder with silicone inserts (090000), Amplatz Needle Holder Accessory Inserts (091000) – Used to direct a standard disposable needle during needle access procedure . It prevents the stylet from backing out of the cannula during puncture and also keeps the hands from radiation exposure during the puncture procedure.

CONTRAINDICATIONS:

1. Coagulopathy
2. Acute inflammatory processes

POTENTIAL ADVERSE EVENTS/ COMPLICATIONS:

1. Fever
2. Urinary tract infection
3. Haematuria
4. Post-operative pain
5. Dysuria
6. Ureteral obstruction caused by stone fragments
7. Obstructive pyelonephritis
8. Urinary urgency
9. Haemorrhage requiring transfusion
10. Loss of kidney
11. Colon injury
12. Perforation of the collecting system
13. Hydrothorax
14. Perforation
15. Pneumothorax,
16. Sepsis
17. Prolonged urinary leak

INSTRUCTIONS FOR USE

Disposable Two-Part Trocar Needle (090020, 090020-ET, 090020-S1); Skinny Needle with Chiba Tip (090010, 090010-ET):

1. Insert internal needle into the external needle to combine them, and puncture any site.
2. Remove the internally held needle, and prepare for guidewire insertion.

Amplatz Needle Holder with silicone inserts (090000), Amplatz Needle Holder Accessory Inserts (091000)

1. Secure the Cook trocar needle in the silicone insert within the clear plastic needle handle.
2. Position your hands to the outside of the device to ensure minimal radiation exposure during the procedure.

TABLE OF MATERIALS

Disposable Two-Part Trocar Needle (090020, 090020-ET, 090020-S1)

Device	Device Components	Material	Material Spec #/ Compound #	Patient Contacting
090020, 090020-ET, 090020-S1	Needle			
	Stylet Knob	Polypropylene	N/A	Indirect
	Stylet	Stainless Steel	AISI304, MS#205	Non-Contact
	Stylet Locking Knob	Polypropylene	S-5827 GPPP TRANSL PINK	Non-Contact
	Needle Hub	Polypropylene	S-5827 GPPP TRANSL PINK	Indirect
	Cannula	Stainless Steel	AISI304, MS#206	Direct
	Product Protector	Polypropylene Non-Radiopaque Tubing	MS#137, Compound 30070	Non-Contact

Skinny Needle with Chiba Tip (090010, 090010-ET)

Device	Device Components	Material	Material Spec #/ Compound #	Patient Contacting
090010, 090010-ET	Needle			
	Cannula	Stainless Steel	AISI304, MS#206	Direct
	Stylet Wire	Stainless Steel	AISI304, MS#205	Non-Contact
	Needle Hub	Polycarbonate	Compound # 80020, CI Solvent Red 207	Indirect
	Stylet Knob	Polycarbonate	Compound # 80020, CI Solvent Red 207	Indirect
	Product Protector	Polypropylene Non-Radiopaque Tubing	MS#137, Compound 30070	Non-Contact

**Amplatz Needle Holder with silicone inserts (090000), Amplatz Needle Holder
Accessory Inserts (091000)**

RPNs Included	Device Components	Material	Material Specification #/ Compound #	Patient Contacting
090000	Needle Holder	Polycarbonate	Polycarbonate CH-10	Non Contact
	Cap	Polycarbonate Rod	Polycarbonate CH-10	Non Contact
	Button	Polycarbonate Rod	Polycarbonate CH-10	Non Contact
091000	Silicone Insert	Silicone	Silicone SILASTIC Q7-4735	Non Contact
	Silicone Insert	Silicone	Silicone, Blue SILASTIC Q7-4735 with blue pigment A35001BLU1	Non Contact

Material of Packaging

Components	Material Specification #/ Compound #	Patient Contacting
Packaging	Part of the paper- medical paper TYVEK: 1073B/109P Film part – Polyethylene 2505SG	Non contact

INFORMATION ABOUT PACKAGING.

The surface density of 50 to 80 (± 2) g/m²

Tensile strength - the tensile Strength is not less than 0.4 kN/m

Weld strength the weld Strength not less than 0.12 MPa

The width of the weld Width of the adhesive layer is not less than 8 mm

Primary packaging is the only.

PERFORMANCE CHARACTERISTICS

Technical specifications for needles described in the TDS:

The connection end of the needle and tubing should not break apart when applying minimal force by the end of the tube and in the direction of their separation, which is 40 N. The resilience: the application of force 10.0 N maximum deviation - 0.45 mm. Straightness: The deviation from straightness shall not exceed 0.2 mm. The maximum permissible deviation from the end of the axis is 3 °. In the event of the end of the needle at a distance of 17.5 mm - needles should not break. The radius of the working parts of the tools should not exceed 0.03 mm. Surface roughness - less than 0.16 microns.

Technical characteristics of the wire guides described in this TDS:

Flexural strength of the conductor is not less than 22 N - if deviation not more than 0.45 mm. Bond strength and coating the conductor wire is not less than 35 N. The strength of the core compound and coating the conductor is not less than 30 N.

PERFORMANCE CHARACTERISTICS

Puncture needles in separate packages, with accessories in packaging are stable under vibration loads in the range 10-55 Hz, amplitude displacement of 0.35 mm and shock loads at peak acceleration of 100 m / s² (10g), and duration of shock acceleration of 16 ms.

Puncture needles in separate packages, with accessories, in packaging are stable at temperatures in the range of +10 to + 42 ° C, relative humidity - 80% at + 25 ° C.

Puncture needles in separate packages, with accessories, in packaging are stable at a transportation and storage temperature in the in the range of -50 to + 50 ° C and relative humidity - 100% (non-condensing).

Puncture needles in separate packages, with accessories are stable at temperatures in the range of +10 to + 42 ° C, relative humidity - 80% at + 25 ° C.

Puncture needles in separate packages, with accessories are stable at a transportation and storage temperature in the in the range of -50 to + 50 ° C and relative humidity - 100% (non-condensing).

BIOCOMPATIBILITY TESTING

Device Name/Description	Biocompatibility Test
Puncture needles in separate packages, with accessories	Cytotoxicity, Sensitization, Intracutaneous reactivity, Hemolysis, Acute systemic toxicity

COMPLIANCE WITH RELEVANT REGULATIONS AND TECHNICAL STANDARDS

Cook Incorporated manufactures all medical products in accordance with Medical Devices Directive 93/42/EEC. Cook Incorporated manufacture all medical products in accordance with the United States Food and Drug Administration (FDA) Quality System Regulations, 21 CFR Part 820, *Quality System Regulations*.

All of Cook facilities are also certified to ISO 13485 (2003) *Quality Management System*.

SHELF LIFE

The “use by” date as denoted by symbols on the labelling. The shelf life is shown by indication the time of manufacture (factory symbol) to the stated expiration date (hourglass symbol). The difference between both dates is the shelf life.

Disposable Two-Part Trocar Needle (090020, 090020-ET, 090020-S1) – 3 year shelf-life

Skinny Needle with Chiba Tip (090010, 090010-ET) -- 3 year shelf-life

Amplatz Needle Holder with silicone inserts (090000), Amplatz Needle Holder Accessory Inserts (091000) – 5 year shelf life

STERILIZATION DETAILS

The sterile devices included in this technical data summary are supplied sterile in peel open packages utilizing an EtO (Ethylene Oxide) validated method to a sterility assurance level (SAL) of 10^{-6} . These devices are intended for one-time use, sterile if unopened or undamaged. Do not use the sterile product if there is doubt as to whether the product is sterile.

Amplatz Needle Holder with silicone inserts (090000) – Reusable device.

Cleaning, Disinfection and Sterilization

NOTE: All cleaning, disinfection and sterilization methods listed herein have been validated.

To Clean

1. Disassemble the Amplatz Needle Holder by removing the silicone inserts from the device.
2. Remove all residual organic matter, blood and irrigation solution to allow the sterilization or disinfection medium to contact the surface of the device. **NOTE:** Ethylene oxide gas and disinfection solutions do not effectively penetrate residual organic matter.
3. Prepare a mild, neutral-pH enzymatic detergent (i.e., Enzol using 1 ounce per gallon) and warm (30-43°C) tap or Deionized (DI) water. Immerse disassembled device and soak for 20 minutes. Use a soft brush to remove all residual organic matter, blood and irrigation solution from the device.
4. Prepare a fresh, neutral-pH enzymatic detergent (i.e., Enzol using 1 ounce per gallon) and warm (30-43°C) tap or DI water. Immerse disassembled device and perform ultrasonic cleaning for 20 minutes.
5. Rinse thoroughly with clean running tap or DI water.
6. Dry all exterior surfaces with a soft cloth.
7. Prior to disinfection/sterilization and reuse, this product should be inspected and reassembled to verify that it functions properly. All surfaces should be smooth and free of abrasions. The device must be disassembled again prior to disinfection and/or sterilization.

NOTE: Automated cleaning methods are not compatible with the Amplatz Needle Holder.

To Disinfect

NOTE: CIDEX® Solution is a high-level disinfection solution that has been found to be compatible with the Amplatz Needle Holder.

- 1. Disassemble the Amplatz Needle Holder by removing the silicone inserts from the device.
- 2. Activate CIDEX® Solution per the manufacturer’s instructions. Test the activated solution prior to use with CIDEX® Solution Test Strips.
- 3. Immerse the disassembled device completely, filling all hard-to-reach areas and eliminating air pockets, in CIDEX® Activated Dialdehyde Solution for a minimum of 45 minutes at 25°C.
- 4. Rinse thoroughly in a room temperature Sterile Water for Injection (SWFI) water bath. Rinse by immersing the device completely for a minimum of one minute. Manually flush all hard-to-reach areas with SWFI water.
- 5. Agitate the device under SWFI water, bring above water level, then re-immersed.
- 6. Repeat steps 4 and 5 two additional times using fresh SWFI for each of the three rinses.

NOTE: Automated disinfection methods are not compatible with the Amplatz Needle Holder.

To Sterilize

NOTE: Autoclave sterilization is not recommended for the Amplatz Needle Holder.

This device can be sterilized by ethylene oxide (EO) gas. EO sterilization should be performed in accordance with industry standards with a validated sterilization cycle. The device must be disassembled prior to sterilization. The recommended cycle for sterilizing this device using 100% EO in a double-pouched configuration is:

- Temperature = 54.4°C
- Humidity = 50-80%RH
- EO concentration = 735 mg/L EO
- Gas dwell time = minimum 60 minutes

STORAGE RULES

These devices should be stored in a dark, dry, cool place; avoid extended exposure to light.

USAGE CONDITIONS

These devices should be used in the patient care medical institutions.

PACKAGING PARAMETERS

Puncture needles in separate packages, with accessories are supplied in individual packaging with instructions for use enclosed. Individual packages are placed into corrugated cardboard boxes and sealed with tape.

TRANSPORTATION RULES

Transportation shall account for the type of transportation employed, the destination, the vendor used for transportation and an estimation of the frequency and harshness of any handling.

PACKAGING PARAMETERS

RPN	Weigth (g)	Height (cm)	Width (cm)
Tolerance	±10%	±10%	±10%
090020	11.2	7.6	40.5

090020-ET	12	7.7	45.7
090020-S1	12	7.7	45.7
090010	9.6	7.7	40.4
090010-ET	9.3	7.6	40.5
090000	109.9	15.2	53.1
091000	20	20.5	38.1


LABELING

The label of the device shows the following device parameters:

The label of the device shows the following device parameters:

- company-manufacturer and manufacturing address;
- device name (trade mark);
- serial number;
- marking about CE certification;
- address of the EU representative;
- trade mark of the manufacturer;
- manufacturing date;
- sizes;
- graphic picture;
- lot number;
- marking for "Sterile";
- marking for "Keep dry";
- marking for "Keep away from sunlight";
- shelf life

Amplatz
Needle
Holder
Accessory
Inserts
(091000)


DISPOSAL RULES

Upon completion of procedure, the devices are to be disposed of as per institutional guidelines for biohazardous medical waste.

For EU countries -in accordance with the World Health Organization guidance document " Safe management of waste from Health Care Activities"

For US countries – in accordance with the US Department of Labor of Occupational Safety and Health Administration (OSHA)

For Russian Federation – in accordance with national requirements. Risk class of the medical waste: class B.

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
CERTIFIED

5. at Indianapolis, Indiana
6. this *Thirtieth* of *May*, 2017
7. by the Secretary of State of Indiana
8. No. A2017-0530104104
9. Seal/Stamp:
10. Signature:

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
Effective May 1, 2003 all apostilles from the Indiana Secretary of State will have an electronically printed seal.
This document was revised March 17, 2012.

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CERTIFICATION STATEMENT

As a responsible representative of Cook Incorporated, I confirm this is an accurate version of the Instructions for Use for Mitty Pollack Needle.

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

Mitty Pollack Needle (090030)

The Mitty-Pollack needle is designed as a one-part stainless needle to allow guide wire access (supplied separately). The needle is manufactured from 18G/1.25mm stainless steel in a length of 14 centimeters. The needle bevel is sloped at a 45 degree angle and is rounded to a smooth finish.

Order Number	Outer Needle		Echotip(Y/N)
	Diameter (Ga)/mm	Length (cm)	
Tolerance	See Technical characteristics	±10%	
090030	18/1.25	14	No

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
Length of cutting edge - not applicable

Width of cutting edge - not applicable

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INTENDED USE/INDICATION FOR USE

Mitty Pollack Needle (090030) -- Used as a percutaneous entry needle to introduce wire guides into the kidney.

Amplatz Needle Holder with silicone inserts (090000), Amplatz Needle Holder

Accessory Inserts (091000) – Used to direct a standard disposable needle during needle access procedure . It prevents the stylet from backing out of the cannula during puncture and also keeps the hands from radiation exposure during the puncture procedure.

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1. Fever
2. Urinary tract infection
3. Haematuria
4. Post-operative pain

5. Dysuria
6. Ureteral obstruction caused by stone fragments
7. Obstructive pyelonephritis
8. Urinary urgency
9. Haemorrhage requiring transfusion
10. Loss of kidney
11. Colon injury
12. Perforation of the collecting system
13. Hydrothorax
14. Perforation
15. Pneumothorax,
16. Sepsis
17. Prolonged urinary leak

INSTRUCTIONS FOR USE

Mitty Pollack Needle (090030):

1. Cover the skinny needle and perform puncturing.
2. After puncturing, remove the skinny needle, and prepare for guidewire insertion.

Amplatz Needle Holder with silicone inserts (090000), Amplatz Needle Holder Accessory Inserts (091000):

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TABLE OF MATERIALS

Mitty Pollack Needle (090030)

Device	Device Components	Material	Material Spec #/ Compound #	Patient Contacting
090030,	Needle			
	Hub	Polycarbonate	Compound # 80020 CI Solvent Red 207	Indirect
	Cannula	Stainless Steel	AISI304, MS#206	Direct
	Silicone Coating	Aliphatic and isopropanol solvents	Dow Corning MDX4-4159 Fluid	Direct Contact
	Product Protector	Polypropylene Non-Radiopaque Tubing	#30070	Non-Contact

**Amplatz Needle Holder with silicone inserts (090000), Amplatz Needle Holder
Accessory Inserts (091000)**

RPNs Included	Device Components	Material	Material Specification #/ Compound #	Patient Contacting
090000	Needle Holder	Polycarbonate	Polycarbonate CH-10	Non Contact
	Cap	Polycarbonate Rod	Polycarbonate CH-10	Non Contact
	Button	Polycarbonate Rod	Polycarbonate CH-10	Non Contact
091000	Silicone Insert	Silicone	Silicone SILASTIC Q7-4735	Non Contact
	Silicone Insert	Silicone	Silicone, Blue SILASTIC Q7-4735 with blue pigment A35001BLU1	Non Contact

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Mitty Pollack Needle (090030) -- 3 year shelf-life

Amplatz Needle Holder with silicone inserts (090000), Amplatz Needle Holder Accessory Inserts (091000) – 5 year shelf life

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- 1. Disassemble the Amplatz Needle Holder by removing the silicone inserts from the device.
- 2. Activate CIDEX® Solution per the manufacturer’s instructions. Test the activated solution prior to use with CIDEX® Solution Test Strips.
- 3. Immerse the disassembled device completely, filling all hard-to-reach areas and eliminating air pockets, in CIDEX® Activated Dialdehyde Solution for a minimum of 45 minutes at 25°C.
- 4. Rinse thoroughly in a room temperature Sterile Water for Injection (SWFI) water bath. Rinse by immersing the device completely for a minimum of one minute. Manually flush all hard-to-reach areas with SWFI water.
- 5. Agitate the device under SWFI water, bring above water level, then re-immerses.
- 6. Repeat steps 4 and 5 two additional times using fresh SWFI for each of the three rinses.

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

- company-manufacturer and manufacturing address;
- device name (trade mark);
- serial number;
- marking about CE certification;
- address of the EU representative;
- trade mark of the manufacturer;
- manufacturing date;
- sizes;
- graphic picture;
- lot number;
- marking for “Sterile”;
- marking for “Keep dry”;
- marking for “Keep away from sunlight”;
- shelf life


Product	Sample Label
---------	--------------

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