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I, PHILLIP JONES, of the Town of Windsor, Berkshire, Notary Public duly admitted and sworn, practising in the said Town,

DO HEREBY CERTIFY AND ATTEST:

That I have no reason to doubt the authenticity of the document hereunto annexed namely Directions of Use signed for and on behalf of Allergan Limited by Denise Wheeler, Senior Regulatory Affairs Analyst, Medical Devices, EMEA Region.

IN TESTIMONY WHEREOF I have hereunto set my hand and affixed my Seal of Office, in the Town of Windsor aforesaid, this day Wednesday, 10th February 2016.

MR PHILLIP JONES

Notary Public
Charsley Harrison LLP
Windsor House
Victoria Street
Windsor
Berkshire
SL4 1EN

08 February 2016

To Whom it may Concern:

I certify that the attached copy Directions for Use for the Natrelle 133 Plus Tissue Expanders L3833 Rev. 2, is an accurate and true copy of the original.

Yours faithfully,

Denise Wheeler
Senior Regulatory Affairs Analyst,
Medical Devices, EMEA Region

Allergan Limited
Marlow International,
Parkway,
Marlow,
Buckinghamshire,
SL7 1YL,
UK

T +44 (0) 1628 494 444
F +44 (0) 1628 494 449

DFU, Style 133Plus TE (INTL)

Allergan
Marlow International, Parkway
Marlow, Bucks, SL7 1YL, UK

Natrelle

ALLERGAN

Marlow International, Parkway, Marlow,
Bucks, SL7 1YL, United Kingdom.

T. +44 (0)1628 494456

F. +44 (0)1628 494956


E. Productsupport_eame@allergan.com

www.allergan.com

NATRELLE® and ALLERGAN® and the Representation
of an Eye Logo are registered trademarks of Allergan,
Inc. All other trademarks mentioned herein are
property of their respective owners.

© 2015 Allergan. All rights reserved.

3833-01 (L3833 Rev.02) 10/2015

ALLERGAN

L3833

Release Date: 12 No

- (EN) 133 PLUS TISSUE EXPANDER
- (FR) EXPANSEUR TISSULAIRE 133 PLUS
- (DE) 133 PLUS-GEWEBEEXPANDER
- (IT) ESPANSORE TISSUTALE 133 PLUS
- (ES) EXPANSOR TISULAR 133 PLUS
- (PT) EXPANSOR TECIDULAR 133 PLUS
- (NL) 133 PLUS-WEEFSELEXPANDER
- (SV) 133 PLUS,VÄVNADSEXPANDER
- (NO) 133 PLUS VEVSEKSPANDERER
- (DA) 133 PLUS VÆVSEKSPANDER
- (FI) 133 PLUS-KUDOSVENYTIM
- (EL) ΔΙΑΤΑΗΡΑΣ ΙΣΤΩΝ 133 PLUS
- (TR) 133 PLUS DOKU GENİŞLETİCİ
- (RU) ТКАНЕВОЙ ЭКСПАНДЕР 133 PLUS
- (UK) ТКАНИННИЙ ЕКСПАНДЕР 133 PLUS

Natrelle

ALLERGAN

PROPRIETARY: This document and information contained may not be produced, used, or disclosed without the written permission of Allergan.
Page 1 of 117

DESCRIPTION

Natrelle® 133 Plus tissue expanders are designed for use in 2-stage reconstruction mammoplasty. All Natrelle® expanders are constructed of a silicone elastomer shell and are latex free.


FIGURE 1

- Natrelle® 133 Plus tissue expanders are intended for temporary subcutaneous or submuscular implantation to develop surgical flaps and additional tissue coverage.

TISSUE EXPANDERS DESIGN FEATURES

- The stable base in the Natrelle® 133 Plus tissue expanders provides greater control over expansion direction.
- Natrelle® 133 Plus tissue expanders are constructed from silicone elastomer and consist of an expansion envelope with a BIOCELL™ textured surface, designed to promote tissue adherence and a MAGNA-SITE™ injection site.
- Natrelle® 133 Plus tissue expanders include an orientation line on the lower pole of the device to aid in properly orienting the device intraoperatively.
- Natrelle® 133 Plus tissue expanders are supplied with the Natrelle® MAGNA-FINDER™ Xact external locating device.
- The MAGNA-SITE™ and MAGNA-FINDER™ Xact contain rare-earth, permanent magnets for accurate location of the MAGNA-SITE™ injection site. In vitro tests show that the MAGNA-SITE™ is detectable through 60mm of phantom tissue.
- All injection sites are self-sealing and contain a titanium needle guard to prevent inadvertent puncture through the base of the injection site (see Figure 2).


FIGURE 2

MAGNA-SITE™ & MAGNA-FINDER™ Xact Locating System

EACH PATIENT MUST BE INDIVIDUALLY EVALUATED FOR IMPLANT SURGERY BASED ON THE CLINICAL JUDGEMENT OF A QUALIFIED SURGEON.

INDICATIONS

- Breast reconstruction in patients with adequate tissue covering following mastectomy or trauma.
- Congenital deformity of the breast.
- Treatment of soft tissue deformities.

CONTRAINDICATIONS

- Tissue covering at expansion site determined inadequate or unsuitable by the surgeon.
- Active infection, that is either local or systemic.
- Existing carcinoma of the breast without mastectomy and residual gross local tumour of the breast after mastectomy.
- Advanced fibrocystic disease considered to be pre-malignant without mastectomy.
- Use of drugs that might result in high surgical risk and/or significant postoperative complications, including drugs that would interfere with blood clotting.
- A patient that demonstrates or shows signs of psychological instability (i.e., an inappropriate attitude or motivation).
- Natrelle® 133 Plus tissue expanders contain a MAGNA-SITE™ and should not be used in patients who already have implanted devices that would be affected by a magnetic field (e.g., pacemakers, drug infusion devices).
- Diagnostic testing with MRI is contraindicated in patients with any Natrelle® tissue expanders in place. The MRI equipment could cause movement of the MAGNA-SITE™ and tissue expander and result in not only patient discomfort but also tissue expander displacement, requiring revision surgery. Also, the MAGNA-SITE™ magnet can interfere with MRI and x-ray detection capabilities. All other Natrelle® implants may be safely imaged using MRI.

RELATIVE CONTRAINDICATIONS

- To varying degrees, radiation damage, ulceration, compromised vascularity, or history of compromised wound healing which may affect tissue covering suitability.
- Previous repeated contour correction failures.
- Patients about to undergo radiation therapy as this may make the use of tissue expanders and breast implants more difficult and increase the risk of complications.
- Physiological condition determined by the surgeon to pose unduly high risk of surgical and/or postoperative complications. To varying degrees, obesity, smoking, diabetes, autoimmune disease, coagulopathy, chronic lung or severe cardiovascular disease may affect patient suitability for surgical implantation.

INFORMATION THAT SHOULD BE PROVIDED TO THE PATIENT

All patients should be well informed on all the potential benefits and risks (see WARNINGS below) associated with the procedure prior to surgery.

Expected benefits include facilitating emotional healing after cancer, eliminating external prostheses, regaining body symmetry, allowing freedom in clothing and physical activities and improving sexual or interpersonal relationships.

Patients should be informed about the the tissue expander/implant reconstructive options available, including the surgical procedure, as well as non tissue expander/implant surgical reconstructive procedure options. As this surgery will most likely be carried out under a general anaesthetic patients should be made aware of the risks associated with anaesthesia. Patients should discuss with their surgeon any history which may indicate a contraindication (relative or absolute) to surgery. Post surgery care should be discussed, including time for wound healing, the need for any drainage tubes, recuperation duration and the need for implant integrity to be evaluated on a regular basis after primary surgery. Patients should be advised to consult a physician or pharmacist before using topical medicines (e.g. steroids) in the breast area and if any clinical examination or surgery in the breast

FU, Style 133Plus TE (INTL)

area is planned the patient should inform the doctor or nurse of the presence of an implant. The surgeon should advise the patient to consult a physician should she suspect any complications. All patients should receive a patient information booklet provided by Allergan.

WARNINGS

Tissue expander patients should be advised that tissue expanders are only to be used for a short term until the tissue has expanded sufficiently.

1. Deflation

Patients should be advised that the breast tissue expanders can deflate at any time and require replacement or revision surgery.

Causes of deflation include:

- Damage by surgical instruments
- Other trauma during surgery, such as improper handling or manipulation
- Underfilling below the recommended fill volume range for the individual implant may result in folds, fold abrasion and potentially, crease-fold failure; overfilling above the recommended fill volume range for the individual implant, which may compromise shell integrity
- Capsular contracture, or abrasive calcifications in the fibrous capsule
- Closed or external capsulotomy
- Stressors such as trauma, intense physical activity, vigorous massage and/or manipulation
- Excessive compression during mammographic imaging
- Leakage through an unsealed or damaged valve

2. Infection

Infection around a tissue expander may occur within days, weeks or even months after surgery. Signs of acute infection reported in association with tissue expanders include erythema, tenderness, fluid accumulation, pain and fever. Erythema may also occur as a normal response to expansion. Infection that is unresponsive to treatment may require tissue expander removal.

Very rarely, Toxic Shock Syndrome has been reported as a possible complication of breast implant surgery and may also be associated with other types of implant surgery.

3. Necrosis

Necrosis may inhibit wound healing and require surgical correction and/or explantation. Permanent scar deformity may occur as a result of necrosis.

4. Haematoma/Seroma

Haematoma/seroma may occur in the postoperative period inhibiting wound healing, or have delayed onset, either of which may require surgical correction and/or explantation.

5. Inflammatory Reaction

Studies evaluating the capsules around textured tissue expanders have reported possible silicone particles within giant cells, indicative of a local (and non-specific) foreign body reaction and silicone granuloma formation. Another study suggests that certain types of capsule cells, including some perceived as giant cells, may actually be secretory cells that form in response to the frictional forces of the tissue expander, providing lubrication at the capsule-expander interface. In case of an inflammatory reaction, the surgeon is advised to remove the device from the patient's body and to secure any evidence on the possible cause of the inflammatory reaction and treat the patient correspondingly. It is advised not to replace the tissue expander until the inflammatory reaction has passed completely and its cause has been eliminated.

Release Date: 12 Nov 2015 00:09:56 GMT -07:00

Revision: 2.0

L3833

6. Extrusion

Unstable or compromised tissue covering and/or interruption of wound healing may result in extrusion of the implant. In case of an extrusion, the device should be regarded as contaminated and should be removed. It may be replaced with another device after the wound has sufficiently healed.

7. Pain

As expected following any invasive surgical procedure, pain of varying intensity and duration may occur following implantation. In addition, improper size, placement, surgical technique, or capsular contracture may result in pain associated with nerve entrapment or interference with muscle motion. Unexplained pain must be promptly investigated.

The expansion process may cause some discomfort, but should not cause excessive pain. Pain may indicate expansion beyond tissue tolerance, which could result in ischemia and necrosis. Further expansion should be discontinued until the pain is resolved.

8. Breast Feeding

Tissue expanders may impact the ability to breast feed, though there is no conclusive clinical study data to support this. The periareolar incision may be associated with a higher likelihood of breast feeding difficulties than other incision sites. The risk of temporary or permanent changes in breast sensation resulting from breast surgery could interfere with the patient's ability to breast feed. Nerve traction and compression have been reported in rare cases in association with tissue expansion.

9. Dissatisfaction with Cosmetic Results

Scar deformity, hypertrophic scarring, capsular contracture, asymmetry, displacement, incorrect size, unanticipated contour, palpability and sloshing may occur. In some cases, cosmetic concerns may also lead to medical concerns. Careful surgical planning and technique can minimise, but not preclude, the risk of such results. Pre-existing asymmetry may not be entirely correctable. Revision surgery may be indicated to maintain patient satisfaction but carries additional considerations and risks.

10. Rotation

Rotation of an anatomical tissue expander may occur. Proper placement and pocket dissection reduces the risk of occurrence. Revision surgery may be necessary to correct rotation. In case of rotation, it is advised to rotate the device back into its correct position in an open surgical procedure. Reshaping of the implant pocket may be necessary to avoid any further rotation in the future.

11. Distortion

Tissue expansion is a time and labour intensive process that may cause temporary discomfort and distortion. If distortion occurs, the cause should be analysed and eliminated.

12. Tissue Atrophy/Chest Wall Deformity

Pressure of a tissue expander may cause tissue atrophy. In rare cases chest wall deformity has also been reported in association with the use of breast implants and tissue expanders.

13. Adulterated Fill

Do not use adulterated fill. *Natrelle® 133 Plus* tissue expanders are to be filled only with sterile saline for injection and only as described in Instructions for Use.

14. Inadequate Tissue Flap

Inadequate tissue flap following expansion may occur and may require additional surgery and expansion.

INSTRUCTIONS FOR USE

SURGICAL PROCEDURE

Allergan relies on the surgeon to know and follow the proper surgical procedures with *Natrelle*® tissue expanders. The surgeon must carefully evaluate implant size and contour, incision placement, pocket dissection and tissue expander placement criteria with respect to the patient's anatomy and desired physical outcome. Planning should include clear delineation of aesthetic goals to ensure mutual understanding between surgeon and patient. The surgeon should observe current and accepted techniques to minimise the risk of adverse and potentially disfiguring reactions.

Natrelle® products are designed and tested for compatibility with sterile water and saline solution. Other substances, such as alcohol or other chemical agents have not been tested in combination with the *Natrelle*® products.

Natrelle® products should not be exposed to extreme heat, cold or pressure. No excessive force should be used while implanting or removing an expander, and, accordingly, the skin incision should be planned for a sufficient size. No sharp objects such as knives or needles should be used in direct vicinity of any *Natrelle*® product, with the exception of the instruments used to fill the expander via the specially designed valve. Back-up tissue expanders must be available during the procedure.

RISK OF REUSE OF SINGLE USE DEVICES

These products are intended for single use only. DO NOT reuse explanted products. Reuse of the device poses a high risk of infection to the patient and damage to the tissue expander resulting in deflation. Reuse of the 21 Ga. Winged Infusion Set may also lead to contamination and/or impairment of functional capacity. Contamination or limited functionality may lead to injury, illness or death of the patient.

PRODUCT IDENTIFICATION

Product labels are supplied within the internal product packaging of each *Natrelle*® implant. The product labels provide specific information which allows product identification.

Important: These labels must be attached to the patient and hospital/doctors records to ensure product identification and device traceability.

STERILE PRODUCT

Each implant is sterilised by dry heat sterilisation and is supplied in a sealed, double primary package.

STORAGE CONDITIONS

Avoid prolonged exposure to extreme storage conditions. Store these devices at ambient room temperatures and at atmospheric pressure and in dry conditions away from direct sunlight.

HOW TO OPEN STERILE PRODUCT PACKAGE

Remove the implant and accessories (where applicable) from their packages in an aseptic environment and using talc-free gloved hands.

DO NOT expose the implant to lint, talc, sponge, towel, skin oils or other surface contaminants.

1. A non-sterile team member peels open the outer package.
2. The surgeon/scrub nurse removes the inner package and places it into the sterile field.
3. Peel open the inner package.
4. Gently retrieve the implant.

L3833

Release Date: 12 Nov 2015 00:09:56 GMT -07:00

Revision: 2.0

Prior to use, keep the implant covered in the inner package, to prevent contact with airborne and surgical field particulate contaminants.

PRELIMINARY PRODUCT EXAMINATION

Prior to use, examine the implant for any evidence of damage or particulate contamination.

TECHNIQUE FOR USING *Natrelle*® 133 PLUS TISSUE EXPANDERS WITH MAGNA-SITE™ INTEGRATED INJECTION SITES

Natrelle® 133 Plus tissue expanders are temporary devices and are not intended to be used for permanent implantation. Tissue expanders should be removed once adequate tissue has been developed. Tissue expansion in breast reconstruction typically requires four to six months. The total expansion period will vary depending on patient tolerance and desired flap size.

1. *Natrelle*® Tissue Expander Placement

Ensure incision is sufficiently large to facilitate insertion and avoid damage to the device. Inadequate pocket dissection increases the risk of deflation and tissue expander malposition. If the pocket is too small, the tissue expander may not have adequate room to unfold, increasing the risk of tissue erosion. Place the tissue expander flat end in the exact desired location, as the BIOCELL™ textured surface promotes tissue adherence to help immobilise the device. The MAGNA-SITE™ should be situated anteriorly, adjacent to the skin surface. The expander is oriented properly when the blue orientation line on the lower pole of the anterior of the device is vertically positioned in the surgical pocket.

2. Location of MAGNA-SITE™ using the MAGNA-FINDER™ Xact

The *Natrelle*® MAGNA-FINDER™ Xact is designed to assist with the location of the MAGNA-SITE™ injection site of the *Natrelle*® 133 Plus tissue expanders. While the injection site can be generally identified by palpation, always verify the location and orientation of the injection site with the MAGNA-FINDER™ Xact, as described below, before each filling.

- a) DO NOT store or use the *Natrelle*® MAGNA-FINDER™ Xact near any loose metal particles, as they may attach themselves to the magnet.
- b) Be sure the magnet inside the *Natrelle*® MAGNA-FINDER™ Xact moves freely without obstruction.
- c) Place the *Natrelle*® MAGNA-FINDER™ Xact on the tissue overlaying the implanted tissue expander. Slowly move the base of the *Natrelle*® MAGNA-FINDER™ Xact in a circular motion (as shown in Figure 3) until the magnet of the locating device detects (points toward) the location of the injection site in one plane.


FIGURE 3

- d) When the magnet on the *Natrelle*® MAGNA-FINDER™ Xact locating device is pointing straight towards the hole in the base of the device, then the MAGNA-SITE™ injection site has been located.
- e) With the magnet centered in the MAGNA-SITE™ injection site, gently press the *Natrelle*® MAGNA-FINDER™ Xact into the skin covering the injection site. The protrusions on the underside of the *Natrelle*® MAGNA-FINDER™ Xact will leave temporary indentations in the skin when the device is lifted. The injection site is located in the center of the circle created by the indentations (Figure 4).


FIGURE 4

As an alternative to the procedure described above, a mark may be made with a marker in each of the four notches around the perimeter of the base of the *Natrelle*® MAGNA-FINDER™ Xact (Figure 5). After all four marks have been made, lift the device from the patient, using the same pen, carefully connect the opposing dots/marks with a line, creating an injection "cross-hair". Then, make a clear mark at the point where the two perpendicular lines intersect. This is the location of the MAGNA-SITE™ injection site.


FIGURE 5

- f) Prepare the injection site for filling using an appropriate antiseptic swab.
3. ***Natrelle*® 133 Plus Tissue Expander Filling**
If the incision site is remote from the site of expansion, the tissue expander may be filled to tissue tolerance at the time of surgery. If the incision site is not remote from the site of expansion, the wound should be stable before tissue expansion begins. However, a slight amount of inflation to fill the pocket space without tension to the tissue may be possible initially.
- a) Insert a new, sterile 21-gauge (or smaller) standard 12° bevel hypodermic needle into injection site. Ideally, the needle should enter perpendicular to the top of the injection site as shown in Figure 6.
- DO NOT penetrate the injection site at an angle.

- b) Penetrate the injection site until the needle is stopped by the needle guard.


FIGURE 6

- c) Fill the tissue expander only with sterile saline for injection and only through the injection site after precise location of the MAGNA-SITE™ with the MAGNA-FINDER™ Xact. Fill carefully to a volume within the recommended fill range specified on the device package labelling and only to patient and tissue tolerance. The suggested fill volume is located on the patch of the tissue expander and is also listed on the device data sheet.

Filling is typically performed at weekly intervals. A *Natrelle*® Expander Fill Volume Record is provided with each device for recording fill volumes and monitoring the expansion process (see back of this book).

DIRECTIONS FOR USE:

Winged Infusion Set Ref: 7B3050

21 Ga. X 1 1/4 in (0.8mm x 44mm)

Instructions for USE: Use Aseptic Technique

Caution: Use care not to puncture the expansion envelope

- Replace per CDC guidelines and/or institutional protocol

1. Remove end of protector cap from the female luer of Winged Infusion Set.
2. Connect the female luer of Winged Infusion Set to male luer of delivery system.
3. Remove the protective sheath from needle of the Winged Infusion Set and insert needle into injection reservoir according to the instructions for the tissue expander.

Contents are sterile until package is opened or damaged.

INSTRUCTIONS AND PRECAUTIONS FOR REMOVAL

In case that removal of the product is necessary, the device should be removed from the patient without harming or breaking the outer shell of the product, if possible. If breakage of the outer shell occurs, be sure to remove all pieces of shell from the pocket. All devices removed from a patient's body should be treated as potentially contaminated and either safely disposed of or processed and returned to Allergan as per the instructions below.

RETURNED GOODS POLICY

Product returns and exchanges must be authorised through your Allergan representative. Exchange value is based on time limitations. All package seals must be intact to be eligible for return or exchange. Returned products may be subject to a restocking charge. Certain products are non-returnable. For more information, please contact your Allergan representative.

REPORTING AND RETURN OF EXPLANTED DEVICES

Explanted devices associated with a complaint or injury must be reported and returned to your local Allergan representative with a Product Field Note (PFN). In order for the explanted product to be returned to the manufacturer it must first be decontaminated and a Decontamination Certificate completed and returned with the explanted device within an explant return kit. Do not return device if the patient has HIV or hepatitis or is known or suspected to carry another infectious agent. Please notify your Allergan representative in these circumstances.

PRE-DISINFECTING INSTRUCTIONS:

Remove any remaining saline solution from the device and vent the device by puncturing the shell with a needle, or nicking the shell with a scalpel.

- Mark the puncture by circling the site with methylene blue, skin marking pencil, or permanent marker.

DISINFECTING METHODS:

Autoclave:

- DO NOT use a prevacuum autoclave or ethylene oxide sterilizer.
- Set autoclave on the "slow exhaust" or "liquid" setting.
- Autoclave by the following gravity displacement cycle: minimum of 70 minutes at 121 °C, 1 kg/cm² (250 °F, 15 psi).
- Open door slowly once the cycle is completed and allow device to cool to room temperature before preparing for shipment.

Bleach (use only if autoclave is not available):

Note: Do not use alcohol, Cidex, formaldehyde or other solvents for disinfecting.*

- Mix one part of household bleach (10% sodium hypochlorite) with 9 parts of water.
- Completely submerge the explanted device in the solution for 60-120 minutes.
- Rinse thoroughly with water and dry the device.

A Product Field Note, Returned Goods Authorisation (RGA) number, shipping instructions, a decontamination certificate within an Allergan return kit is required prior to dispatch of each explanted device. These can be obtained by contacting your local Allergan office/local distributor.

PRODUCT REPLACEMENT POLICY

Tissue expanders are covered under Allergan's Product Replacement Policy. To receive a replacement device, a Product Field Note (PFN) must be completed and submitted to your Allergan representative. Explanted devices must be returned in accordance with section above "Reporting and Return of explanted devices". The product replacement policy does not cover surgical or other expenses related

to rupture, deflation, cosmetic revision, capsular contracture, or other adverse events. Contact your local Allergan representation for specific warranty details for your region.

LIMITED WARRANTY, LIMITATION OF LIABILITY, AND DISCLAIMER OF OTHER WARRANTIES

Allergan warrants that reasonable care was used in the manufacture and production of this product. Allergan has no control over the conditions of use, patient selection, surgical procedure, post-surgical stresses, or handling of the device after it leaves our possession. Allergan does not warrant either a good effect or against an ill effect following its use. Allergan shall not be responsible for any incidental or consequential loss, damage or expenses directly or indirectly arising from use of this product. Allergan's sole responsibility in the event that Allergan determines the product was defective when shipped by Allergan, shall be replacement of the product. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law, or otherwise, including, however not limited to, any implied warranties of merchantability or fitness for use.

INFORMED CONSENT

An Informed Consent Form is provided (see back of book):

Please ensure that the patient receives the information from the section "INFORMATION THAT SHOULD BE PROVIDED TO THE PATIENT" and understands the information provided. The patient must realise that the surgical and post surgical risks associated with tissue expanders cannot be completely predicted, even with the best medical manufacturing, technology and surgical care and accept these conditions and limitations. Patients must fully inform their physician of their medical history, including any and all conditions that would contraindicate tissue expansion or implant surgery. Failure to inform their physician could result in significant surgical and post-surgical complications. The patient must decide herself on whether the expected benefits outweigh the said risks. If the patient decides that the expected benefits of the proposed tissue expander surgery outweigh the risks, then she must take full responsibility for her choice to proceed with tissue expander surgery. The two-part form located at the back of this book should be completed and signed once the patient has decided to proceed with tissue expander surgery. This form allows both the patient and the surgeon to retain copies for their records.

ID CARD

As part of device tracking, Allergan is providing the patient with an ID Card. The information on the inside of this card (when folded) is specific to the device(s) the patient received (patients should keep this card for their records and carry it at all times to facilitate medical care in case of emergency). When the device is replaced, another card will be provided for the new device. The ID Card is located at the back of this book.

Instructions for the Surgeon:

Place one label from each product in the appropriate space on the ID Card (L or R). These labels can be found attached to the bottom of the main inner label on top of the inner product packaging. In the absence of a Patient Record Label, copy the Product code, REF, SN and LOT from the package label to the space provided. Fill in all the remaining sections. Give this entire document to the patient for her records. See "Graphical Symbols" for explanation of symbols.

U, Style 133Plus TE (INTL)

GRAPHICAL SYMBOLS

SFR
SFV

SN
LOT
REF
ID CARD

21 Ga. WINGED INFUSION SET ONLY

SALINE FILL RANGE
SALINE FILL VOLUME

CHECK FOR LEFT BREAST IMPLANTATION

CHECK FOR RIGHT BREAST IMPLANTATION

SERIAL NUMBER

LOT NUMBER

CATALOGUE NUMBER

IDENTITY CARD

STERILE, DRY HEAT STERILISED,
DATE OF STERILISATION, YEAR & MONTH

SINGLE USE ONLY — DO NOT REUSE

USE BY, YEAR & MONTH

ATTENTION, SEE INSTRUCTIONS FOR USE

ANATOMICAL IMPLANT / TISSUE EXPANDER HEIGHT

ANATOMICAL IMPLANT / TISSUE EXPANDER WIDTH

ANATOMICAL IMPLANT / TISSUE EXPANDER PROJECTION

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CONFIDENTIAL: This document and information contained may not be reproduced, used, or disclosed without the written permission of Allergan.
Page 111 of 117 Effective

Style 133Plus TE (INTL)

EN Natrelle® Informed Consent Form

- I confirm that I have read the Natrelle® implant/tissue expander patient information brochure.
- I confirm that I have had ample opportunity to discuss my concerns and options with my surgeon.
- I confirm that my surgeon has given me time to consider Natrelle® implant/tissue expander surgery.
- I fully understand the effects and risks associated with having Natrelle® implant/tissue expander surgery.
- I consent to undergoing Natrelle® implant/tissue expander surgery as discussed with my surgeon.

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PROPRIETARY: This document and information contained may not be reproduced, used, or disclosed without the written permission of Allergan. Effective

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PROPRIETARY: This document and information contained may not be produced, used, or disclosed without the written permission of Allergan. Effective
Page 113 of 117

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PROPRIETARY: This document and information contained may not be produced, used, or disclosed without the written permission of Allergan.
Page 114 of 117 Effective

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PROPRIETARY: This document and information contained may not be produced, used, or disclosed without the written permission of Allergan. Effective

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Name/ Nom/ Name/ Nome/ Nombre/ Nome/ Naam/ Namn:

Address/ Adresse/ Anschrift/ Indirizzo/ Dirección/ Morada/ Adres/ Adress:

ALLERGAN, Marlow International, Parkway, Marlow,
Bucks, SL7 1YL, United Kingdom

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