

Instructions for Use for registration of surgical lights ADMECO Lux LED

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ADMECO LUX Camera Manual

Please read the operating and installation instructions carefully before using the camera. This prevents accidents and allows for optimal use of the camera.

ADMECO

ADMECO AG

Turbistrasse 17
CH-6281 Hochdorf
Switzerland


phone 0041 41 914 70 30
fax 0041 41 910 61 19
email info@admeco.ch

web www.admeco.ch

11/2015 StS

1	Used Symbols	4
2	Important information for safe use	5
	2.1 Identify the components	5
	2.2 Environment conditions	5
	2.3 Combination with other medical devices	5
	2.4 Proper use	6
	2.5 Improper use	6
	2.6 Maintenance	6
	2.7 Cleaning and disinfection	6
	Sterilize the Sterile Camera Handle	7
	2.8 Changing the device fuses	7
	2.9 Disposal	8
3	Safety Instructions	8
	3.1 Location requirements	8
	3.2 Requirements for the risk management	9
4	Description of the Vidi OP system	10
	4.1 The control unit	10
	4.1.1 Operational functions	10
	4.1.2 Start-up behaviour	13
	4.1.3 Connections of the control unit	13
	4.1.4 Cable length	14
	4.2 Mounting the camera module	15
	4.2.1 Insert the camera module	15
	4.2.2 Verify	15
	4.2.3 Lock the module	15
	4.3 Replacing the sterile camera handle	16
	4.3.1 Insert the sterile camera handle	16
	4.3.2 Remove the sterile camera handle	16
5	Installing the Vidi OP System	17
	5.1 The Control Unit	17
	5.2 The lighthouse camera	17
	5.3 The Adaption Unit	17
	5.4 Electric scheme	18
	5.5 Option separat arm	19
6	Technical Data	20
7	Appendix	21
	7.1 Questions	21
	7.2 Replacement Equipment	21
	7.3 EN 60601-1-2 Compliance	21
	7.3.1 Tested cable lengths	23
	7.3.2 Recommended Safety Distances	23
	7.4 Pin assignment	24
8	Declaration of conformity	25

1 Used symbols

	Gas explosion: warns of the danger of the explosive ignition of gas mixtures.
	Electric shock: warns of an electric shock which may result in severe or even fatal injury.
	Warning: indicates a potential hazard. If not avoided, this may cause death or serious injury.
	Damage to surface: warns of the risk of damage to the surface through the use of incorrect agents or disinfectants.
	Disposal: indicates, that this device is not allowed to be disposed in domestic waste.
	Refer to operation instructions: advises users / operators to refer to these operation manual.
	CE-Conformity mark: confirms conformity of the device to the guidelines of the European Medical Device Directive (MDD).
	Pay attention to the operating instructions: This symbol is found on the labels of the devices. It indicates that the operators have to consult the operating instructions.

2 Important informations for safe use

2.1 Identify the components

Each component of the Vidi OP Systems has a unique identification label. For the control unit, this label can be found at the bottom of the device. This label also consists of the serial number and the MAC address of the Ethernet interface.

Product designation: „Vidi OP System“

Model designation control unit: HD-VVD-ST-001

Model designation adaption unit: HD-VVD-AA-001

Model designation camera module: HD-VVD-CA-001

2.2 Environment conditions

For normal operation:

Temperature: 10°C to 40°C

Relative humidity: 5% to 95%

Air pressure: 700 hPa to 1060 hPa

For interim storage or transport:

Temperature: -15°C to 60°C

Relative humidity: 5% to 95%

Air pressure: 500 hPa to 1060 hPa

2.3 Combination with other medical devices

Read the operating instructions for combined medical products

- The Vidi OP system can be combined with medical devices from other manufacturers. For the operation of these devices, please refer to the respective operating instructions.
- Only medical devices approved in accordance with IEC 60601-1 or UL 60601-1 may be connected to the system. If a medical device is installed afterwards, the installation must be performed as specified in IEC 60601-1. Compliance with this standard must be ensured by the service technician responsible.
- The usage of invalid or unadmitted devices could result in a reduction of the immunity of the Vidi OP System!

2.4 Proper use

- The Vidi OP system is used to record the progress of surgery in the surgery field in HD quality.
- The video image is available not only for the purpose of documentation, but also for teaching and training purposes in the form of live transmissions.
- The Vidi OP system may only be operated, cleaned and disinfected by qualified personnel / authorised technicians.

2.5 Improper use

- The Vidi OP system is not suitable for use as an exclusive visualization medium for surgery. It must be ensured that the surgery staff have a direct view of the area of surgery at all times.
- The Vidi OP system is not suitable for diagnostic purposes!
- The Vidi OP system is not suitable for use in rooms or environments with atmospheres containing combustible mixtures of gases like anaesthetics with air, oxygen or N2O (laughing gas).

2.6 Maintenance

There is no need for any maintenance.

Cleaning and disinfection

Only use disinfectants approved by the manufacturer for use on the following materials: Polycarbonate (PC), polyamide (PA), polyvinyl chloride (PVC), acrylonitrile-butadiene-styrene copolymer (ABS), Polyetherimide (PEI) and on silicones.

Risk of contamination and infection of patients caused by excess disinfecting / cleaning agent

If cleaning agent or disinfectant is allowed to penetrate the device, surplus cleaning agent / disinfectant may drip into open wounds:

- Use cleaning agents / disinfectants sparingly so that no liquid enters the device.
- The cloth used to clean / disinfect the device should only be damp, not wet.
- Only use a wipe-down disinfectant to disinfect the device.

Damage to surfaces / lacquered surfaces

- Only use agents that do not contain chlorides or halides.
- Do not use benzine, paint thinner, alkaline / acidic agents or alcohol (e.g., ethanol, propanol, etc.) or aldehyde-based cleaning agents.
- Do not use abrasive cleaners.

Electric shock hazard

Touching live device parts may result in an electric shock! Always disconnect the device from the main power supply prior to cleaning and disinfecting:

- Switch off the mains switch in the OR. The camera system / lighting system is disconnected from the mains.
- Ensure the main switch is safeguarded against being inadvertently switched back on.

2.7 Sterilize the Sterile Camera Handle

Sterilization must be carried out in accordance with valid steam sterilization procedures, e.g. in a sterilizer in accordance with EN 285:2008 and validated in accordance with ISO 17665-1:2006.

Risk of damage to property from improper steam sterilization.

The sterile camera handles can last approx. 100 steam sterilization cycles with proper steam sterilization. Improper steam sterilization can damage the sterile handles and make the surface porous or crack it.

- Never exceed the maximum temperature of 134 °C.
- Sterilize the sterile camera handles individually in packaging suitable for steam sterilization.

2.8 Changing the device fuses

The fuse holder can be found on the back of the device. It is part of the mains socket. To be able to change the fuses, unplug the mains cord.

1. To remove the fuse holder, carefully pry the clamping bracket off using a screwdriver

Fuse specification

size	5x20
specifications	T3.15A / 250V H

- and pull the fuse holder out of the fuse compartment.
2. Remove both fuses from the fuse holder and insert new fuses.
3. Insert the fuse holder back into the fuse compartment and push it in until the clamping bracket engages.

Important:

- Use only fuses with high breaking capacity (H) according to IEC 60127-2!
- Always replace both fuses simultaneously!
- The fuses may only be changed by authorised and qualified technicians!

2.9 Disposal

- The Vidi OP system should be disposed of in accordance with the pertinent national regulations and at a suitable waste disposal point for the recycling of electrical and electronic devices.
- The Vidi OP system meets the requirements of Directive 2002/95/EC RoHS (restriction of the use of certain hazardous substances in electrical and electronic equipment).

3 Safety Instructions

3.1 Location requirements

- The Vidi OP system is not suitable for operation in potentially explosive atmospheres.
- The ambient temperature must be within the range 10°C to 40°C!
- The device must be set up or installed in such a way that sufficient air circulation is assured to cool it. Without such an sufficient air circulation, the device may be overheated or damaged.
- The device is not suitable for use in rooms or environments with atmospheres containing combustible mixtures of anaesthetics with air, oxygen or N2O (laughing gas).
- The device must be placed in a way that it can be easily disconnected from the mains at any time.

3.2 Requirements for the risk management

- Ensure, that under no circumstances, the patient and one of the components of the Vidi OP System is touched at the same time. This includes wires to the monitor or the camera too.
- The power switch of the Vidi OP System is not designed to disrupt the power! To completely disconnect the device from mains, unplug the mains cord.
- The Vidi OP system may only be operated in a electricity supply network with protective ground connections.
- Ensure, that under no circumstances, the camera module is changed while a surgery is in progress.
- Ensure, that under no circumstances, the system is used without an annual inspection in accordance with EN62353.
- Ensure, that under no circumstances, a camera module with integration into the light head is used without sterile camera handle.
- Ensure that all components and sterile handles are sterilized and disinfected. For details see chapter "Cleaning and disinfection" on page 6.
- Ensure that worn or damaged components are replaced and not used.
- It is not permitted to perform any modifications on the system or on individual components of the system.

4 Description of the Vidi OP system

The Vidi OP system is designed to “convert” the progress of surgery into high quality video images. These video images are available not only for the purpose of documentation, but also for teaching and training purposes in the form of live transmissions. Therefore the Vidi OP System provides plenty of standard interfaces.

A Vidi OP system consists of the following components.

- the HD camera module offering high resolution images with 1920x1080 pixels.
- the control unit for controlling the camera module.
- several connections for networking and output media.
- the adaption unit for adapting different cabling types.

4.1 The control unit

The control unit combines the operational functions to remote control the camera module and the connections of peripherals like monitors, recording units and connections to the local Ethernet. All of the video interfaces are active at the same time, so there is no need to toggle between them. The serial port of the control unit can be used to remote control the device from an external device or computer.

Important note:


To support further Ethernet based components, the control unit integrates an Ethernet switch!

4.1.1 Operational functions

The front panel consists of the power switch and the keypad. Figure 1 shows the whole front panel. The power switch physically disrupt the power! This is not a stand-by mode.

figure 1: the front panel of the Vidi OP system.

Position: This function is not provided by using the camera in the light head.

Rotation

Rotation: the rotation functions can be used to rotate the camera module clockwise or anticlockwise. The range of rotation is approximately 370°.

The camera keeps on rotating as long as the key is pressed!

Zoom

Zoom: use the zoom keys to enlarge (+) or reduce (-) the image section. The camera module allows an enlargement of up to 10x (optical zoom ratio).

The camera keeps on zooming as long as the key is pressed!

Focus

Focus: the Vidi OP system supports automatic focussing (AF) on distinctive image areas with higher contrast. To focus on parts of image outside these areas, the focus can be adjusted manually.

Press the **AF** key to toggle the automatic focus feature on or off. An active auto focus is indicated by a lighted LED.

Use **Near** and **Far** to adjust the focus as needed. The camera keeps on adjusting as long as the key is pressed!

Exposure: the Vidi OP system allows controlling the visual impression of the brightness of the image. Based on the current image content, the automatic exposure metering calculates an optimal brightness value. For the very few cases, where the automatic exposure scores the wrong areas, the exposure settings can be adjusted manually.

Press the AE key to toggle the automatic exposure on or off. An active auto exposure is indicated by a lighted LED. Use the "**Sun**" and "**Moon**" keys to adjust the exposure gradually as needed.

White Balance

White Balance: this feature is used to adapt the colour temperature setting of the camera to the given environment.

Info

Colour temperature is a characteristic of visible light. It is measured in Kelvin [K]. Higher values are called cool colours (blueish white), and lower values are called warm colours (yellowish, reddish white).

Use the followings steps to adjust the colour temperature:

1. adjust the surgery light as needed.
2. Place a white piece of paper in the lighted field.
3. Align the camera. Zoom into the paper to make sure, that the whole camera image shows white.
4. Press the "White Balance" key. The adjustment is performed automatically. This takes about 5 seconds. The new white balance settings are stored by the camera module.

Freeze

Freeze: the Vidi OP system is able to show a still image of the camera.

Pressing the "Freeze" key toggles between the live video and a still image. The still image is recorded each time the key is pressed. An active still image is indicated by a lighted LED.

4.1.2 Start-up behaviour

After turning the power on, all LEDs in the front panel start flashing. As soon as the device becomes ready, the LEDs stop flashing. They are now displaying the state of the corresponding functionality.

The defaults are:

- auto exposure "AE" is active
- auto focus "AF" is active
- zoom is set to a middle range
- image rotation is unchanged

To achieve the commands described above, the control unit communicates with the camera module. If this communication fails, the LEDs "AE" and "AF" start flashing at the same time. As soon as the communication is established again, the previous state of the LEDs is restored.

4.1.3 Connections of the control unit

The Vidi OP system control units has several connections at the rear of the device. These connections covers networking and output media for several SD and HD standart

figure 2: available connections on the back plane of the device.

Connections at the rear (from left to right):

1. FBAS / CVBS: composite video output (SD). Only connect MDD compliant devices!
2. Y/C: S-Video output (SD). Only connect MDD compliant devices!
3. Two DVI-D connections for two HD-compatible widescreen full HD monitors. Only connect MDD compliant devices!
4. Two HD-SDI connections for two HD-SDI-compatible widescreen full HD monitors. Only connect MDD compliant devices!
5. YUV / YCbCr: analogue component signal (HD). Only connect MDD compliant devices!

6. Two Ethernet LAN connections. These RJ45 sockets are connected to the internal Ethernet switch. Use one socket for the uplink.
7. Camera: current-carrying Binder connector for the camera module of the Vidi OP System.
8. Serial: optional RS-232 connection to remote control the Vidi OP system. Only connect MDD compliant devices!
9. 50Hz / 60Hz: switch the frame rate of the complete video paths between 50Hz and 60Hz. In the case of the FBAS/CVBS and Y/C interfaces, this means PAL (50Hz) and NTSC (60Hz).
10. Mains socket for non-heating appliances for the connection of a power supply cable. Includes a fuse holder for two device fuse, type: 5x20 T3,15A H 250 V
11. Ground stud

4.1.4 Cable length


The length of all cables for video or data interfaces are defined by the corresponding standards. For additional informations, see section "Tested cable lengths" on page 19.

4.2 Mounting the camera module

The camera module of the Vidi OP System can be mounted and unmounted with ease. The steps below describes the process of mounting the camera module into the light head unit.

4.2.1 Insert the camera module

To insert the camera module into the mount ring of the light head, ensure that the two latches are open. If the latches are open, the red marks are visible. The camera base is coded to ensure the correct orientation. Insert the camera module into the mount.

4.2.2 Verify

Press the camera module into the mount ring. Ensure that the module is inserted completely.

4.2.3 Lock the module

Press both latches down until the red marks disappear!

Important note:

The module is fully mounted, if the red marks are no longer visible! This is necessary to ensure a secure operation of the module.

4.3 Replacing the sterile camera handle

If the camera module is mounted into the light head, a sterile camera handle must be used.

4.3.1 Insert the sterile camera handle

Ensure that the sterile handle was properly disinfected and sterilized!

Push the sterile camera handle onto the mounted camera module until notice an perceptible resistance. Turn the handle until the securing button audibly snaps into place. After installing the handle make sure that it is secure by pulling it towards you without pressing the securing button.

4.3.2 Remove the sterile camera handle


Press the securing button and pull the sterile handle from the camera module. Check the sterile handle for wear and/or damage. Dispose and replace defective components!

5 Installing the Vidi OP System

The “Vidi OP System” consists of the following three units.

5.1 The Control Unit

The control unit is the main component of the Vidi OP System. The video signal of the camera is received and processed here. All common video signal are generated by this component. The user is able control all features of the separate pendant unit from one single place.

5.2 The light head camera

The light head holds the camera module. The light head is mounted to a spring arm. The sterile handle (3) used to manually adjust the point of view is a standard part too.

5.3 The Adaption Unit

The adaption unit is used to adapt the different physical characteristics of the two types of cables. The adaption unit is powered by the light voltage. It is important, that the metal case of the adaption unit is grounded, e.g. to the ceiling mount. To ensure a RF compliant grounding, use at least a flexible grounding strap ($\geq 10\text{mm}^2$).

PIN	colour	Description
1		
2		
3	white	Video line 2 (Video -)
4	brown	Video line 1 (Video +)

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


connection on building site
Installation bauseits

Pro Dimmer eine geschirmte Leitung nutzen!
Use 1 shielded cable per dimmer!
Veuillez installer un câble blindé par variateur!

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
6 Technical Data

Camera module HD-VVD-CA-001

Resolution	Full HD, 1920 x 1080 pixel
Display refresh rate	PAL/NTSC / 1080i/50 / 1080i/60
Zoom	10x optical
Lens Wide	3,4 mm / Tele 33,9 mm
Signal-to-noise ratio	50 dB
Electric shutter	1/2 to 1/10.000s
Min. illuminance	12,0 lux
White balance	automatic/manual
Motorized Turn	> 360°
Still image	digital
Focusing	automatic/manual
Amplifier	automatic/manual (-3 bis 28 dB, 16 Steps)
Power supply	24 V from the light head
Power Consumption	< 8W
Type of protection	IP41

Control unit HD-VVD-ST-001

Voltage	110 V - 240 V
Fusing	2 Stück T3.15 A / 250 V H / 5x20
HD-SDI	1080i/50 / 1080i/60
DVI	1080i/50 / 1080i/60
CVBS	PAL/NTSC
Y/C	PAL/NTSC
YUV	1080i/50 / 1080i/60
Power Consumption	< 25 W
Type of protection	IP40

data interfaces (control unit)

External remote control	RS-232 / V.24
Ethernet port	2x 10/100 MBit/s

7 Appendix

7.1 Contact Address

For questions about the Video OP System or about this documentation, feel free to contact

Admeco AG
Turbistrasse 17
CH-6280 Hochdorf
Tel: +41 41 914 70 30
Mail: info@admeco.ch

7.2 Replacement Equipment

For the Vidi OP System, there is no replacement equipment available. Faulty systems must be sent to the manufacturer as a whole.

7.3 Zulassung EN 60601-1

4.2.2 Emission

(according to EN 60601-1-2, medical electrical equipment, class B)

Standard Table	Type, line	Test	Normative reference Limits (required)	Executed	Limits (met)	Notes
EN 60601-1-2 6.1.1	AC power supply	Conducted emissions 0,15-30MHz	CISPR 11/22 class B	yes	class B	
EN 60601-1-2 6.1.3.1		Harmonic current (input current up to 16A / line) (for input power 75W...1KW)	IEC 61000-3-2 class A	no		4)
EN 60601-1-2 6.1.3.2		Voltage fluctuations and flicker (input current up to 16A / line)	IEC 61000-3-3	no		5)
EN 60601-1-2 6.1.1	Enclosure, lines	Radiated emissions 30-1000MHz	CISPR 11/22 class B	yes	class B	6)
EN 60601-1-2 5.2.1.1a	Manual	Examination of specification of failure criteria	-	no		
In making these tests only the EMC properties were examined without making a statement concerning the safety of the EUT.						

- 4) Not tested because the power consumption of the EUT is <75W.
 5) Not tested because there's no significant disturbance or influence estimated..
 6) The measurement was executed on an alternative test site according to prEN 50147-3.
 All other tests were executed in the accredited area.

4.2.1 Immunity

(according to EN 60601-1-2, medical electrical equipment)

Standard Table	Type, line	Test	Norm. ref. Criterion *) (required)	Executed	Criterion *) (achieved)	Notes
EN 60601-1-2 6.2.2	Enclosure, lines	Electrostatic discharge (ESD) $\pm 2/4/6$ kV contact-discharge (conduct. parts and indirect) $\pm 2/4/8$ kV air-discharge	IEC 61000-4-2 acc. to 6.2.1.10 ³⁾ acc. to 6.2.1.10 ³⁾	yes yes	yes yes	
EN 60601-1-2 6.2.3		Radiated RF 80–2500 MHz, 3V/m, 80% AM / 1 kHz, dwell time 1s (10V/m if the EUT is a life-sustaining device) (2Hz if physiological data processing)	IEC 61000-4-3 acc. to 6.2.1.10 ³⁾	yes	yes	1)
EN 60601-1-2 6.2.8		Power frequency magnetic field 50 Hz, 3A/m	IEC 61000-4-8 acc. to 6.2.1.10 ³⁾	yes	yes	2)
EN 60601-1-2 6.2.4		AC/DC power supply, signal/data lines	Fast transients (Burst) ± 2 kV, asym., 5/50 ns, 5 kHz, power supply lines, 60s ± 1 kV, asym., 5/50 ns, 5 kHz, signal lines, 60s (if line > 3m and is not coupled to patient)	IEC 61000-4-4 acc. to 6.2.1.10 ³⁾ acc. to 6.2.1.10 ³⁾	yes yes	yes yes
EN 60601-1-2 6.2.5	High energy impulses (Surge) (for mains supply) $\pm 0.5/1$ kV, symmetrical, 1.2/50 μ s (8/20), 5x 0°/90°/ 80°/270° $\pm 0.5/1/2$ kV, unsymmetrical, 1.2/50 μ s (8/20), 5x 0°/90°/180°/270°		IEC 61000-4-5 acc. to 6.2.1.10 ³⁾ acc. to 6.2.1.10 ³⁾	yes yes	yes yes	
EN 60601-1-2 6.2.6	Conducted RF 0.15–80 MHz, 3V, 80% AM / 1 kHz, dwell time 1s (10V at ISM-frequency if life-maintaining device) (2Hz if physiological data processing)		IEC 61000-4-6 acc. to 6.2.1.10 ³⁾	yes	yes	
EN 60601-1-2 6.2.7	Power variations/dropouts (deviation allowed if not life-maintaining device with >1kVA and max. 16A per phase) >95% for 0.5 periods, at zero point, 3x /10s pause 60% for 5 periods, at zero point, 3x/10s pause 30% for 25 periods, at zero point, 3x/10s pause		IEC 61000-4-11 acc. to 6.2.1.10 ³⁾ acc. to 6.2.1.10 ³⁾ acc. to 6.2.1.10 ³⁾	yes yes yes	yes yes yes	
	Power variations/dropouts (deviation allowed, but alarm required if life-maintaining device) >95% for 5s, at zero point, 3x/10s Pause		IEC 61000-4-11 acc. to 6.2.1.10 ³⁾ (temp. deviation allowed)	yes	yes	

*) The compliance criteria are defined on page 6.

1) The EUT was tested in the frequency range from 80-2700 MHz with 3 V/m.

2) The EUT was tested with a value of >110 A/m, but only in the X-axis in fact of the too large physical dimensions.

3) See the match of criterion item 6.2.1.10 of the used standard EN 60601-1-2 (see page 6).

7.3.1 Cable lengths

The Vidi OP System accomplishes the EN 60601-1-2. The table below shows the cable properties used for all tests.

Warning: the usage of different cable types could reduce the immunity of the whole system! The use of longer cable length may cause an increased emission or a reduced interference resistance.

Connector	shielded	Length (m)
230V AC main cord	no	3.0
Ethernet	yes	5.0
DVI-D	yes	4.0
Y/C	yes	4.0
Camera	yes	1.5
HD-SDI	yes	4.0
YUV	yes	4.0
FBAS/CVBS	yes	4.0

7.3.2 Recommended Safety Distances Between Portable and Mobile HF Telecommunication Devices and the Vidi OP System

The Vidi OP System is intended for use in an electromagnetic environment with controlled RF disturbances. The user of the Vidi OP System can help to avoid electromagnetic disturbances by keeping the minimum distance (safety distance) between portable and mobile telecommunication devices (transmitters) and the device - depending on the output power of the telecommunication devices as described below.

7 Appendix




Nominal power transmitter (W)	safety distance (m)		
	150 kHz – 80 Mhz	80 MHz – 800 Mhz	800 MHz – 2,5 Ghz
0.01	0.1	0.1	0.2
0.1	0.4	0.4	0.7
2	1.6	1.6	3.3
10	3.7	3.7	7.4
100	12	12	23

Typical cell phones use 2W transmission power. The safety distance for this devices is 3.3m!

7.4 Pin Assignment

Camera mount SUB-D 9pin (female)	description	adaption unit terminal block	RJ45 - Pin	signal
1	Twisted-Pair Video line 1	4	8	Video+
2				
3	shield / GND-supply voltage	1.2	3.6	GND
4				RS232-TxD
5				
6	Twisted-Pair Video line 2	3	7	Video-
7				
8	supply voltage +15...35V		1.2	+U _b (24V)
9				RS232-RxD

topview of the sockets

SUB-D 9pin (female)	terminal block	RJ45-socket
		

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is corresponding to the Essential Requirements (ERs) of the European Community
medical device directives (MDD) 93/42/EEC and 2007/47/EG and to the special requirements of the standard
EN60601.

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
Seite 1 von 1

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Medizintechnik
Medical Engineering

ADMECO AG
Turbistrasse 17
CH-6281 Hochdorf
Switzerland

Telefon +41 41 914 70 30
Telefax +41 41 910 61 19
E-Mail info@admeco.ch
www.admeco.ch

The document contains 27 pages bound, numbered and sealed pages.

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ADMECO LUX LED EL

Manual

Please read the operating and installation instructions carefully before using the surgical light. This prevents accidents and allows for optimal use of the surgical light.

1	Used Symbols	3
2	Description of the surgical light	4
3	Instructions for installing the surgical light to the „Multi“ ceiling system	5
3.1	Example: Arrangement and attachment to „Multi“ system	5
3.2	Installing the ceiling anchor ring	6
3.3	Installing the ceiling flange and stand tube	6
3.4	Installing the adapter with bom and spring arm	6
3.5	Installing the surgical light on the spring arm	7
3.6	Removing the spring arm	8
3.7	Checking the operating voltage	8
3.8	Connections on the power supply unit	9
4	Operating instructions	10
4.1	Operating conditions	10
4.2	Positioning before starting an operation	10
4.3	Sterilisable handle	10
4.4	Switching on and off	11
4.5	Adjusting the brightness	11
4.6	Endo mode	11
4.7	Centering the light beam	11
4.8	Repositioning during the operation	11
4.9	Moving the arms	11
5	Light care and disinfection	12
5.1	Cleaning and disinfecting the surgical light	12
6	Maintenance	13
6.1	Maintenance of the surgical light	13
6.2	Replacing the LED module	13
6.3	Adjusting the bracket brake	14
6.4	Ajusting the spring arm brake	15
7	Safety	16
8	Troubleshooting	17
8.1	Electrical causes	17
8.2	Mechanical causes	17
9	Miscellaneous	18
9.1	Serial number and type plate	18
9.2	General information on the LED modules	18
9.3	Storage and transportation	18
9.4	Disposal and environmental protection	18
10	Technical data	19
10.1	power supply unit	19
10.2	Light data	20
10.3	Electromagnetical Compability	20
10.4	Electronical scheme	21
10.5	Dimensions	22
11	Spare parts	23
11.1	Projector (focusable)	23
11.2	Focusable light	24
11.3	Ceiling system „Multi“	25
11.4	Wall bearing	26
11.5	Ceiling anchor ring	27

1 Used symbols

Warning: Indicates a potential hazard. If not avoided, this can cause a death or serious injury.

Electric shock: warns of an electric shock which may result in severe or even fatal injury.

Damage to surface: warns of the risk of damage to the surface through the use of incorrect agents or disinfectants.

Disposal: indicates, that this device is not allowed to be disposed in domestic waste.

CE-Conformity mark: confirms conformity of the device to the guidelines of the European Medical Device Directive (MDD).

Pay attention to the operating instructions: It indicates that the operators have to consult the operating instructions.

2 Description of the surgical light

The ADMECO LUX surgical light is an airstream-optimised lighting system for use in operating theatres, and is based on individual high-performance projectors. LEDs are used in these projectors. Thanks to this concept, ADMECO LUX surgical lights can be combined in a range of sizes for all types of operation according to modular design principles.

Originally developed specially for use under laminar supply-air systems in operating theatres, the system is now highly valued by surgeons worldwide thanks to its excellent lighting properties and exceptionally easy handling. Despite the compact projector design, extremely high light yields can be achieved. The surgical light is easy to operate and focus thanks to the sterilisable handle.

The ADMECO LUX surgical light has an extremely low air-resistant surface. As a result, the build-up of upward-moving airstreams and aerodynamic stagnation points underneath the surgical light are minimised.

Sterile surgical lights are positioned in the operating theatre so that the sterilisable handle cannot come into contact with the head of the surgeon at any time. This means that ADMECO LUX lights can be used at greater distances to the surgical site compared to lights where sterile manipulation is not possible.

Practical experience has shown that ADMECO LUX surgical lights are used approximately 1.2 to 1.4 metres away from the wound area. The projector optics are designed so that a light beam diameter of 22 cm is achieved at a distance of 1.2 metres.

Intended use

The surgical light is used in the vicinity of patients for the local lighting of surgical sites on the body and generates the appropriate illumination intensity. It is intended for use in operating theatres and as an aid in patient treatment and diagnosis.

The surgical lights are intended for use under ambient conditions according to DIN EN 60601, section 10.2.1 in medical rooms, which are equipped according to the guidelines set out in VDE 0100, part 710. The surgical lights are installed by the manufacturer or an authorised representative of ADMECO AG. Circuit diagrams and national regulations must be observed during installation and maintenance.

The surgical lights must be operated by trained specialist personnel, and not by patients. The operator must read the operating manual and ensure that the device is working properly before using the surgical light.

The surgical light is intended for continuous use under the following conditions:


temperature: 0...40°C

humidity: 10...75%

air pressure: 500...1060hPa

3 Instructions for installing the surgical light to the “Multi” ceiling system

3.1 Example: Arrangement and attachment to “Multi” system

3.2 Installing the ceiling anchor ring

Due to the strict safety requirements inside the operating theatre, the equipotential bonding between the ceiling / building reinforcements (or the steel building structure) and the surgical light must be isolated correctly. For this reason, the surgical light is insulated when suspended.

A ceiling with high static load resistance is required when suspending from the ceiling. The ceiling anchor ring must be attached to the ceiling with the anchor plugs according to the installation instructions from the plug manufacturer (UPAT, HILTI etc.), together with the washers and nuts. The pull-out strength of each anchor screw (plug) must be tested to a load of 12.5 KN.

3.3 Installing the ceiling flange and stand tube

The ceiling flange and stand tube are delivered as a complete, pre-assembled unit together with the connection terminals.



Installing the ceiling flange:

Unscrew the top nut on the bolt of the ceiling anchor ring. Insert the insulation sleeve and washer into the flange hole. Insert the ceiling flange into the bolt of the ceiling anchor ring, align it perpendicularly at a distance of approx. 40 mm using the washers and nuts, then fasten into place.

3.4 Installing the adapter with boom and spring arm


The adapter is pre-assembled with the booms and spring arms on delivery. Before installing the adapter, the cylindrical surface where the stand tube is inserted must be coated with graphitic grease.

The plug pair is then joined together. Ensure that the plug for the main light and the plug for the satellite lights are both assigned correctly. Insert the adapter into the stand tube and fasten it with the six grub screws. Align the adapter so that the long connecting port for the main light points towards the head of the patient (anaesthetic) and the short connecting port (satellite) points towards the feet, and are flush with the axis of the operating table. Tighten all of the grub screws. Different adapter arrangements are possible according to customer demands.


3.5 Installing the surgical light on the spring arm

The sliding part on the spring arm must be opened first, as shown below (1/2). Coat the pins of the surgical light with graphitic grease (3), insert into the spring arm (4) and attach the nose segment (5). Push the sliding part over the nose segment to secure it in place, then fasten with the screw and lock washer (6).

1.

2.


3.

4.

5.

6.

3.6 Removing the spring arm

If the spring arm ever needs to be replaced, then proceed as follows: Remove the surgical light as described in chapter 2.5, but in reverse order.

Loosen the screw on the boom of the spring arm and pull out the boom cap with plug. Remove the circlip and both narrow washers. Pull the spring arm downwards and remove the wide washer.

Installation of the new spring arm is made in reverse order.

3.7 Checking the operating voltage

In order to prevent a possible reduction in the service life of the LED modules, it is important to check the operating voltage.

The voltage can be measured directly on the terminals in the central joint (1).

Nominal V = 22...28 VDC

1.

3.8 Connections on the power supply unit

The wiring diagram for connecting the input can be found in chapter 9.

Terminal X1: INPUT

INPUT - Main power supply	
Voltage	95...250Vac
Output	max. 120 W
Frequency	50/60 Hz
Terminal cross-section	4mm ²

Terminal X2: AUXIN

OUTPUT - to the surgical lights	
Volatage	24Vdc
Terminal cross-section	4mm ²

VOUT - surgical light output			
Voltage	24Vdc		
Power	max. 110W		
diameter wiring	<10m	2.5 mm ²	
	10-20m	4 mm ²	
	20-30m	6 mm ²	!!! > connection diameter !!! - 4 mm ²
	30-50m	10 mm ²	!!! > connection diameter !!! - 4 mm ²

4 Operating instructions

4.1 Operating conditions

See also chapter 1

The surgical light is intended for continuous use under the following conditions:

temperature: 0...40°C

humidity: 10...75%

air pressure: 500...1060hPa

4.2 Positioning before starting an operation

We recommend that the surgical light is positioned by the staff beforehand so that any adjustments during the operation are avoided as far as possible. The surgeon then only has to make fine adjustments on the sterilisable handle during the operation, which can be made very easily thanks to the simple handling of the surgical light.

4.3 Sterilisable handle

The sterilisable handle is inserted and locked into place as shown, with the groove pointing towards the release lever on the surgical light. To remove the handle, push the release lever on the surgical light upwards and pull out the handle. The sterilisable handle is delivered non-sterile.

Spindle extension:

Installation of the spindle extension ensures optimal handling of the ADMECO sterilisable handle. When required, the cap can be removed and the extension screwed on. The cap prevents scratches to the sterilisable handles.

4.4 Switching on and off

The surgical light is switched on and off by shortly pressing the head button.

4.5 Adjusting the brightness

Depending on the depth of the wound, the light needed can increase during an operation. In order to prevent the surgeons from glare during the incision and suturing phases, the illumination intensity of the surgical lights is reduced before the operation. The brightness can be increased or decreased accordingly during the operation.

You can dim the lights by pressing the head button for a longer time.

4.6 Endo mode

To set the Endo mode the light has to be switched off.

Press the button for a longer time. The light will start up with the Endo mode.

4.7 Centering the light beam

It is recommended to position (centre) the surgical light on the wound area before the incision using the sterilisable handle. This is made by rotating the handle, which adjusts the pitch angle of the projectors. It is then possible to either concentrate the light beams on top of one another (minimum field size) or to increase the light spot so that the beam is opened sufficiently enough so that no visible dark points appear in the centre of the lit area (maximum field size). An adjustment of the surgical light during the operation is unnecessary, unless the distance between the wound area and the surgical light changes significantly.

4.8 Repositioning during the operation

The compact design and easy movement of the surgical lights mean that angle adjustments can be made without any problems during operations. This means that obstruction of the surgical light is largely prevented, depending on the position of the surgeon.

4.9 Moving the arms

The rotary resistance between the boom and the adapter is greater than that between the boom and spring arm. This then makes it possible to rotate the spring arm without the boom following.

5 Light care and disinfection

Important: The light is not waterproof. In order to prevent the penetration of liquids, the light should not be handled with excessive amounts of cleaning agent. Do not use organic solutions during cleaning.

Notice: Exposure to intensive, continuous ultraviolet light causes decay to plastics and coatings. The light should be protected from direct exposure to UV light (caused by germicidal lights, direct sunlight etc.).

The following care recommendations should be adhered to:

- use a cleaning agent in a water solution
- Prevent moisture from penetrating the housing interior
- Dry the cleaned surfaces with a soft cloth

To sterilise the handle, proceed as follows:

Remove the handle and wash the handle with water containing cleaning solution.

Sterilization must be carried out in accordance with valid steam sterilization procedures, e.g. in a sterilizer in accordance with EN 285:2016 and validated in accordance with ISO 17665-1:2006.

The temperature during sterilization should not exceed 134°C. The grip can last approximately 100 sterilization cycles.

Regular sterilisation and use can lead to natural material fatigue (cracks, scratches).

Replace the handle if you detect any cracks, discolouration or breakages.

Thanks to the design, the surgical light can be cleaned using a non-staining disinfectant without any problems.

6 Maintenance

6.1 Maintenance of the surgical light

Annual maintenance can be made by the in-house technical service team or ADMECO service team. We suggest making safety inspections annually. Standard EN 60601 applies for this inspection.

Maintenance:

- Greasing (lubrication) and any setting / adjustment of the articulated arms.
- Greasing (lubrication) of the boom / adapter connection.
- Greasing (lubrication) of the boom / spring arm connection.
- Greasing (lubrication) of the spring arm / light body connection.
- Check of the illumination intensity (in lux) on the surgical light. If there is a difference of more than -20% between the actual value and nominal value, then this may be due to clouding of the protective glass on the projectors.
- Check of the operating voltage. When at maximum output voltage on the dimmer, the voltage on the terminals in the central joint must not deviate strongly from the rated value of 24 VDC. A higher or lower voltage can lead to malfunctions on some of the projectors.

6.2 Replacing the LED module

Top of the projector:


Unscrew the cap (1) from the projector. Pull out the plug (3) used for the power supply using a suitable gripping tool (**do not pull the end of the cable**).

Bottom of the projector:

Remove the clamping ring (3) and pull out the protective glass using a vacuum lifter (4, included with delivery). Using the new LED module, the old module can be released and removed by rotating to the left several times (5/6/7).

Important: The lens array must only be touched around the edges!


Installation of the new LED module is made in reverse order.

1.

2.

3.

4.

5.

6.

7.

6.3 Adjusting the bracket brake

Adjustment of the bracket brake is necessary when rotation of the light around the centre of gravity on the bracket becomes too easy or too difficult.

Procedure:

Remove the cover lid and adjust the grub screw (hexagon socket, M8 x 8) on the bracket base using a 4 mm Allen key.

6.4 Adjusting the spring arm break

Adjustment of the bracket brake is necessary when rotation of the light around the centre of gravity on the bracket becomes too easy or too difficult.


increase the break force

Adjust the break force with a screwdriver in + direction on both sides with the same power

decrease the break force

Adjust the break force with a screwdriver in - direction on both sides with the same power

7 Safety

The following safety instructions should be read and applied carefully in order to ensure safe installation and handling.

The electrical system must correspond to DIN VDE 0100 and local regulations. This means that no additional measures are required relating to the safety of the ADMECO LUX surgical light.

Important: In order to avoid the risk of electric shocks, this device must only be connected to a power supply network with protective conductors.

The ADMECO surgical light is not suitable for use in explosive atmospheres.

If a single LED module fails, then the ADMECO LUX surgical light can still be used.

Light intensity in watts per m²

The standard states that a maximum value of 1000 watts per m² should not be exceeded in order to prevent burns to the patient. This value is exceeded when more than ten LED projectors on an ADMECO LUX LED combination are used to illuminate the wound area. A working distance of 1 metre is used as a basis. When this working distance is reduced (e.g. 80 cm), then the light intensity increases by 1.56 times! An increase in working distance to 1.2 metres leads to a 30% reduction in light intensity.

As absorbed radiation is transformed into heat, this can result in burns.

LEDs are extremely concentrated light sources. Avoid looking directly into the light.

ADMECO accepts no liability for the safety characteristics as a result of misuse, insufficient care or maintenance. All repairs or modifications must be made exclusively with original ADMECO parts. Modification of the device is not permitted.

The light has been designed for a service life of 25 years.

8 Troubleshooting

Malfunctions on the surgical light may be due to the following causes:

8.1 Electrical causes

- a) Defective power supply unit
- b) Defective dimmer control
- c) Operating voltage too high or too low

8.2 Mechanical causes



- a) **The spring arm is unstable**, meaning the surgical light does not remain at the set height.

- **Surgical light rises:**

The screw is located behind the little cover. Turn the allen wrench to right as it is shown in the picture.


- **Surgical light sinks:**

Turn the allen wrench to left as it is shown in the picture.

- b) **The surgical light can be moved too far or not far enough upwards.** The height stop on the surgical light can be adjusted by removing the covers.

- Adjusting the height stop upwards:
Nut must be turned higher (+, downwards)
- Adjusting the height stop downwards:
Nut must be turned lower (-, upwards)

9 Miscellaneous

9.1 Serial number and type plate

The serial number and type plate are located on the cover of the surgical light.

Note: Always give the serial number and the location of the surgical light in the building when contacting our customer service team.

Meaning of the serial number:

EL1000001 -> 00 01 0000
 Year Generation Sequential number

9.2 General information on the LED modules

LEDs have a very long service life (50,000 hours). However, with increased usage the illumination intensity is reduced when under the same loads (approximately 75% brightness after 50,000 hours). An assessment must then be made as to whether further use is possible, based on a measurement of the illumination intensity and subjective evaluation.

9.3 Storage and transportation

The following storage conditions apply for periods of up to 15 weeks:

- Temperature: -25 to +70 °C
- Rel. humidity: 10 to 75%
- Air pressure: 500 to 1060 hPa

Storage is made in closed or covered rooms only. After this, the values specified for the operating conditions in the instructions apply.

9.4 Disposal and environmental protection

The light is returned to the manufacturer for disposal. A special recycling company can also be appointed for disposal. Any defective or worn parts detected during maintenance must be disposed of correctly.

10 Technical data

10.1 EL power supply unit

Rated values	
Type:	EL
Item no.:	40.20.0519
Power:	max. 120 W
Frequency:	50/60 Hz
Primary voltage:	95... 250VAC
Secondary voltage:	24 VDC, max. 4A
Connections:	Clamps (see schema) 4mm ²
Protection class:	IP 20 / Dry installation
Ambient temperature	0°C...40°C

The installation plan is found in Chapter „Dimensions“.

10.2 Light data

criteria	2EL no focus	3EL	4EL
illuminance in 1m distance (lx)	80000	120000	160000
irradiance (W/m²)	250	375	500
diameter light field d50 (cm)	17		
diameter light field d10 (cm)	26		
remaining illuminance with 1 mask (%)	1	1	1
remaining illuminance with 2 masks (%)	55	34	36
remaining illuminance with tube (%)	99	99	99
remaining illuminance with tube and 1 mask (%)	0.5	0.5	0.5
remaining illuminance with tube and 2 masks (%)	54	34	35
Depth of illuminance L1+L2 (m)	1.10		
Ra	86		
R9	26		

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







In order to prevent incorrect deliveries, please specify the serial number of the surgical light.

11.2 Focusable light

Projektor

Abdeckung

Zentralgelenk

Spindel kurz

13

1) O-ring for pins, 21.95 x 1.78 40.20.0330

2) Pin, processed, disassembled 40.20.0090

3) Grub screw, hexagon socket, M6 x 6 91.12.3006

4) Cover 40.20.0006

5) Flat-head slider for bracket 40.20.0062

6) Sleeve, white 40.20.0005

7) Release lever (short), for spindle (short) 40.20.0012

8) Dowel pin for release lever, 3 x 20 91.55.9001

9) Grub screw, M6 x 8, unscrewed 40.21.0039

10) Pressure spring for release lever 91.55.9001

11) Pressure spring for bracket brake 91.55.9002

12) Grub screw, M8 x 8, for bracket brake 91.11.4008

13) Sterilisable handle 40.20.0048

14) Plug, three-pole 40.21.0034

15) Cap 40.20.0354


16) Spindle extension 40.20.0333

Vacuum lifter, Ø 60 mm, for protective glass 40.20.0099

Gleitmo 810 grease for plugs, sockets etc. 40.20.0284

In order to prevent incorrect deliveries, please specify the serial number of the surgical light.



11.3 Ceiling system „Multi“

1) Casing on boom, grey/white	40.20.0403
2) Washer, Ø 48/32 x 0.5	40.20.0224
3) Rondelle circlips ø39/32 x 0.5	40.20.0226
4) Washer with external tab, Ø 39/32 x 1	40.20.0226
5) V-circlip, d = 32	94.61.9005
6) Cover for angle adjustment, grey/white	40.20.0227
7) Cover for spring adjustment, grey/white	40.21.0167
8) Raised countersunk head screw, M4 x 8	90.90.9019
9) Sliding flange, grey/white	40.20.0229
10) Oval-head screw with Phillips head, M4 x 10	90.90.9017
11) Lock washer, M4	91.61.9001
12) Deflection leaf spring, white	40.20.0230
13) Nose segment	40.20.0231
14) Round plate, Ø 35 mm, matt silver	40.20.0089
15) Stand tube, Ø 95/80 mm, grey/white	40.21.0048
16) "Multi 2" adapter, assembled, grey/white	40.21.0221
17) "Multi 3" adapter, assembled, grey/white	40.21.0214
18) "Multi 4" adapter, assembled, grey/white	40.21.0215
19) Spring arm, 4.5 – 7.5 kg, grey/white	40.20.0396


21) Boom, 2000 V, grey/white, L=750 m, NL	40.20.0464
22) Boom, 2000 S, grey/white, L=750 m, NL	40.20.0397
23) Ceiling flange, Ø 320 mm, six-hole, grey/white	40.20.0043
24) Ceiling anchor ring, Ø 330 mm, assembled, M12, g/w	40.21.0095
25) Hexagon screws, M8 x 20	90.10.4020
26) Grub screws, M8 x 12, with base	91.13.4012
27) Plug	40.20.0118
28) Socket	40.20.0117
29) Cap, white	40.20.0067
30) Socket, three-pole	40.21.0161
31) Socket, three-pole	40.90.0029
32) Plug, three-pole	40.90.0015
33) Holder, W	40.20.0425
Power supply unit, EL	40.20.0519
Gleitmo 810 grease for plugs, sockets etc.	40.20.0284

11.4 Wall bearing

1) Casing on boom, grey/white	40.20.0403	17) Boom, 2000 S, grey/white, L=750 m, NL	40.20.0397
2) Washer, Ø 48/32 x 0.5	40.20.0224	18) Boom, 2000 V, grey/white, L=750 m, NL	40.20.0464
3) Rondelle circlips ø39/32 x 0.5	40.20.0225	19) Wall bearing, assembled, grey/white	40.21.0089
4) Washer with external tab, Ø 39/32 x 1	40.20.0226	20) Wall bearing cover, grey/white	40.20.0114
5) V-circlip, d = 32	94.61.9005	21) Wall bearing panel, grey/white	40.20.0115
6) Cover for angle adjustment, grey/white	40.20.0227	22) Cap for wall bearing cover, grey/white	40.20.0116
7) Cover for spring adjustment, grey/white	40.21.0167	23) Socket, three-pole	40.90.0004
8) Raised countersunk head screw, M4 x 8	90.90.9019	24) Socket, three-pole	40.90.0029
9) Sliding flange, grey/white	40.20.0229	25) Plug, three-pole	40.90.0015
10) Oval-head screw with Phillips head, M4 x 10	90.90.9017	26) Holder, W	40.20.0425
11) Lock washer, M4	91.61.9001	Power supply unit, EL	40.20.0519
12) Deflection leaf spring, white	40.20.0230		
13) Nose segment	40.20.0231		
14) Round plate, Ø 35 mm, matt silver	40.20.0089		
15) Spring arm, 4.5 – 7.5 kg, grey/white	40.20.0396		
		Gleitmo 810 grease for plugs, sockets etc.	40.20.0284

11.5 Ceiling anchor ring

Order number

Ceiling anchor ring, disassembled

40.20.0042

Individual assembly on threaded bolt, M12 (6 pcs. required per ceiling anchor ring)

40.21.0037

ADMECO AG

Turbistrasse 17
CH-6280 Hochdorf
Switzerland

phone 0041 41 914 70 30
fax 0041 41 910 61 19
email info@admeco.ch

web www.admeco.ch

The document contains 31 pages bound, numbered and sealed pages.

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1	Used symbols	3
2	Description of the surgical light	4
3	Instructions for installing the surgical light to the „Multi“ ceiling system	5
3.1	Example: Arrangement and attachment to „Multi“ system	5
3.2	Installing the ceiling anchor ring	6
3.3	Installing the ceiling flange and stand tube	6
3.4	Installing the adapter with boom- and spring arm	6
3.5	Installing the light head on the spring arm	7
3.6	Removing the spring arm	8
3.7	Checking the operating voltage	8
3.8	Connections on the power supply unit	9
4	Operating instructions	10
4.1	Operating conditions	10
4.2	Positioning before starting an operation	10
4.3	Sterilisable handle	10
4.4	Overview Controller	11
4.5	Switching on and off	11
4.6	Adjusting the brightness	11
4.7	Endoscopic light	11
4.8	Change the color temperatur	12
4.9	Centering the light beam	12
4.10	Repositioning during the operation	12
4.11	Moving the arms	12
5	Light care and disinfection	13
6	Maintenance	14
6.1	Maintenance of the surgical light	14
6.2	Replacing the LED module	14
5.3	Adjusting the bracket brake	15
5.4	Adjusting the spring arm break	16
7	Safety	17
8	Troubleshooting	18
8.1	Electrical causes	18
8.2	Mechanical causes	18
9	Miscellaneous	19
9.1	Serial number and type plate	19
9.2	General information on the LED modules	19
9.3	Storage and transportation	19
9.4	Disposal and environmental protection	19
10	Technical data	20
10.1	Dimmer	20
10.2	APS power supply unit	20
10.3	Light data	22
10.4	Electromagnetic Compability	22
10.5	Electrical scheme	23
10.6	Dimensions	25
11	Spare parts	26
11.1	Projector	26
11.2	Light head	27
11.3	Ceiling system „Multi“	28
11.4	Wall bearing	29

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
Surgical lights are positioned in the operating theatre so that the sterilisable handle cannot come into contact with the head of the surgeon at any time. This means that ADMECO LUX lights can be used at greater distances to the surgical site compared to lights where sterile manipulation is not possible.

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
humidity: 10...75%

air pressure: 500...1060hPa

3 Instructions for installing the surgical light to the "Multi" ceiling system

3.1 Example: Arrangement and attachment to "Multi" system

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


3.5 Installing the light head on the spring arm

- 1 + 2 remove controller / cover
- 3 + 4 remove the security- ring and plate
- 5 + 6 pull out the protection cover and und carefully plug in light head
- 7 + 8 mount the security plate and cover; adjust the force of the break with the screw
- 9 + 10 connect the controller and mount the covers

3.6 Removing the spring arm

If the spring arm ever needs to be replaced, then proceed as follows: Remove the surgical light as described in chapter „Installing the light head on the spring arm”, but in reverse order.
Loosen the screw on the boom of the spring arm and pull out the boom cap with plug. Remove the circlip and both narrow washers. Pull the spring arm downwards and remove the wide washer.

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
3.7 Checking the operating voltage


In order to prevent a possible reduction in the service life of the LED modules, it is important to check the polarity of the operating voltage.

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
Nominal V = 24...28 VDC (see photo in chapter “Replacing the LED module”).

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Terminal X1: INPUT

INPUT - Main power supply	
Voltage	100...240Vac
Power	max. 115 W
Frequency	50/60 Hz
Terminal cross-section	4mm ²

If the main power supply is not used, then "AUX" is automatically active.

Terminal X2: AUX



AUXIN - 24V Input / Emergency power supply	
Voltage	24V ac/dc
Power	max. 115 W
Terminal cross-section	4mm ²

The in-house 24 V emergency power supply system can be connected to the emergency power input. If there is no input voltage present on the INPUT, then the power supply unit switches automatically to AUX. **The dc-voltage has to be connected in the way it is indicated!**

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The light output is connected with the terminals on the stand leg of the surgical light.

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

reducing the brightness:

* push the button „darker“

the Led-beam indicates the position

increasing the brightness:

* push the button „brighter“

the Led-beam indicates the position

4.7 Endoscopic light

wall controller:

* switching on: push the buttons „brighter“ and „darker“ at the same time for 2 sec.

* switching off: push the buttons „brighter“ and „darker“ at the same time again for 2 sec.

the LED-beam will change his indication

spring arm controller:

* switching on: push the button „Endo“ for 2 sec.

* switching off: push the button „Endo“ for 2 sec.

4.8 Change the color temperatur

4.8.1 Wall controller

change the color temperatur

- * push the button „color temp“
- The LED beam will change his indication.
- * push the button „darker“ to reduce the color temperatur (warmer light)
- * push the button „brighter“ to increase the color temperatur (colder light)
- * leave the color temperature mode by pushing the button „color temp“. The mode will automatically stop after 10sec.

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4.10 Repositioning during the operation

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Procedure:

Remove the cover lid and adjust the screw.

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
Light intensity in watts per m²

The standard states that a maximum value of 1000 watts per m² should not be exceeded in order to prevent burns to the patient. This value is exceeded when more than ten LED projectors on an ADMECO LUX LED combination are used to illuminate the wound area. A working distance of 1 metre is used as a basis. When this working distance is reduced (e.g. 80 cm), then the light intensity increases by 1.56 times! An increase in working distance to 1.2 metres leads to a 30% reduction in light intensity. As absorbed radiation is transformed into heat, this can result in burns.

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Meaning of the serial number:

LED14000001 -> 00 01 0000
 Year Generation Sequential number

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

10 Technical data

10.1 Dimmer

This dimmer is especially suited for controlling the ADMECO APS power supply units.
Installation: See wiring diagram
Function: See chapter 3


Technical data		
Type:	Controller ARC	
Item no.:	40.20.0558	
Voltage:	max. 24VDC	
Connections:	Screw clamp 1,5mm²	
Dimensions: Installation depth:	88x88x49mm 35mm	for SEV und DIN device sockets
Protection class:	IP20	Dry installation
Ambient temperature:	0°C bis 45°C	

The ADP-1310 is mainly used as a controller for the APS power supply unit:

- Sends control commands (buttons) to the APS.
- Regularly queries the current brightness value in the APS.
- Displays the current brightness value on the LED bar.

10.2 APS power supply unit

Nenndaten	
Item No.:	40.20.0557
Type:	APS
Power:	max. 125VA
Frequency:	50/60 Hz
Primary voltage:	100... 240VAC 50/60 Hz
Secondary voltage:	24 VDC, max. 5A
incl. temperatur switch:	120°C
Connections:	clamps (see scheme)
Protection class:	IP 20 dry installation
Ambient temperature:	0°C....40°C
primary fuse:	230VAC 3A
secondary fuse:	24VDC 5A

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10.3 Light data

criteria	3SL	4SL	5SL
illuminance in 1m distance (lx)	96000	125000	160000
irradiance (W/m²)	312	416	520
diameter light field d50 (cm)	17		
diameter light field d10 (cm)	26		
remaining illuminance with 1 mask (%)	1	1	1
remaining illuminance with 2 masks (%)	34	35	37
remaining illuminance with tube (%)	99	99	99
remaining illuminance with tube and 1 mask (%)	0.5	0.5	0.5
remaining illuminance with tube and 2 masks (%)	34	34	37
Depth of illuminance L1+L2 (m)	1.10		
Ra	89		
R9	30		

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

Projector	
1) Projector, assembled	40.20.0477
3) LED module	40.20.0394
4) Protective glass, assembled	40.21.0187
5) Clamping ring, punched	40.20.0020
6) Projector cap, assembled	40.21.0232
7) Sleeve	40.20.0005


11.2 Light head

In order to prevent incorrect deliveries,
please specify the serial number of the
surgical light.


1) O-ring for pins, 21.95 x 1.78	40.20.0330
2) Pin, processed, disassembled	40.20.0090
3) Grub screw, hexagon socket, M6 x 6	91.12.3006
4) Cover	40.21.0206
5) Flat-head slider for bracket	40.20.0062
6) Sleeve, white	40.20.0005
7) Release lever (short), for spindle (short)	40.20.0012
8) Dowel pin for release lever, 3 x 20	91.65.9001
9) Grub screw, M6 x 8, unscrewed	40.21.0039
10) Pressure spring for release lever	91.55.9001
11) Pressure spring for bracket brake	91.55.9002
12) Grub screw, M8 x 8, for bracket brake	91.11.4008
13) Sterilizable handle	40.20.0048
14) Plug, three-pole	40.21.0034
15) Cap	40.20.0354
16) Spindle extension	40.20.0333
Vacuum lifter, Ø 60 mm, for protective glass	40.20.0099
Gleitmo 810 grease for plugs, sockets etc.	40.20.0284

11.3 Ceiling system „Multi“

1) Casing on boom, grey/white	40.20.0403	21) Boom, 2000 V, grey/white, L=750 m, NL	40.20.0464
2) Washer, Ø 48/32 x 0.5	40.20.0224	22) Boom, 2000 S, grey/white, L=750 m, NL	40.20.0397
3) Rondelle circlips ø39/32 x 0.5	40.20.0225	23) Ceiling flange, Ø 320 mm, six-hole, grey/white	40.20.0043
4) Washer with external tab, Ø 39/32 x 1	40.20.0226	24) Ceiling anchor ring, Ø 330 mm, assembled, M12, g/w	40.21.0095
5) V-circlip, d = 32	94.61.9005	25) Hexagon screws, M8 x 20	90.10.4020
6) Cover for angle adjustment, grey/white	40.20.0227	26) Grub screws, M8 x 12, with base	91.13.4012
7) Cover for spring adjustment, grey/white	40.21.0167	27) Plug	40.20.0118
8) Raised countersunk head screw, M4 x 8	90.90.9019	28) Socket	40.20.0117
9) Sliding flange, grey/white	40.20.0229	29) Cap, white	40.20.0067
10) Oval-head screw with Phillips head, M4 x 10	90.90.9017	30) Socket, three-pole	40.21.0161
11) Lock washer, M4	91.61.9001	31) Socket, three-pole	40.90.0029
12) Deflection leaf spring, white	40.20.0230	32) Plug, three-pole	40.90.0015
13) Nose segment	40.20.0231	33) Holder, W	40.20.0425
14) Round plate, Ø 35 mm, matt silver	40.20.0089	34) Baldachin	40.20.0428
15) Stand tube, Ø 95/80 mm, grey/white	40.21.0048	Power Supply Unit, APS	40.20.0557
16) "Multi 2" adapter, assembled, grey/white	40.21.0221	Dimmer APS	40.20.0558
17) "Multi 3" adapter, assembled, grey/white	40.21.0214	Gleitmo 810 grease for plugs, sockets etc.	40.20.0284
18) "Multi 4" adapter, assembled, grey/white	40.21.0215		
19) Spring arm, 4.5 – 7.5 kg, grey/white	40.20.0396		

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phone 0041 41 914 70 30
fax 0041 41 910 61 19
email info@admeco.ch

web www.admeco.ch
30

The document contains 32 pages bound, numbered and sealed pages.

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ADMECO LUX LED EL mobile Manual

Please read these instructions carefully before fitting or using the surgical light. Following these instructions will prevent unnecessary damage and optimize use.

1	Used Symbols	3
2	The ADMECO LUX Surgical Lights	4
3	Assembly	5
3.1	Assembly of the flexible OR light stand	5
3.2	Assembling the OR lights	6
3.3	Sterilisable handle	6
4	Advice on use	7
4.1	Operating conditions	7
4.1	Positioning prior to surgery	7
4.2	Switching on and off	7
4.3	Adjusting the brightness	7
4.4	Endo mode	7
4.5	Centering the beam	7
4.6	Repositioning during surgery	8
5	Care and disinfecting	8
6	Maintenance	9
6.1	Maintenance of surgical lights	9
6.2	Regulating the spring power	9
6.3	Changing the LED module	9
6.4	Ajusting the bracket brake	10
7	Safety	11
8	Troubleshooting	12
9	Miscellaneous	12
9.1	Labels	12
9.2	Storage and transportation	12
9.3	Disposal and environmental protection	13
10	Appendix	13
10.1	Model type	13
10.2	Spare parts	13
10.3	Technical data	14
10.3.1	EL power supply unit	14
10.3.2	Statics	14
10.3.3	Light data	14
10.3.4	Electromagnetic compability	15

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The ADMECO LUX Surgical Lights

The ADMECO LUX surgical lights use single high-power LED projector light heads to provide an air-stream optimized surgical lights system. ADMECO's single projector system is realized in a sophisticated modular construction allowing an entire range of surgical lights in various sizes for all types of surgery and operations.

Originally developed for use with laminar flow systems in operating rooms, the ADMECO LUX surgical lights combine outstanding lighting quality with incomparable simplicity in use – two reasons why they enjoy such a good reputation among so many operating room staff throughout Europe and in many other countries worldwide.

The lights generate extremely high levels of light thanks to the compact projector design. The surgical lights are easy to service and focus using the sterilizable grip.

ADMECO LUX surgical lights offer very little surface for air resistance. Consequently, they avoid problems related both to reverse air-streaming or a stagnant air-flow point under the lamp itself. Surgical lights for sterile use are positioned in the operating room so as to rule out any possible contact between the surgeon's head and the sterilized grip.

As a result, these lights are placed at a greater distance to the operating field than lights that allow non-sterile handling. Experience has shown that ADMECO LUX surgical lights perform best at approx. 1.20 - 1.40 m from the wound area. The projectors are optically designed to generate a light beam of 22 cm at a distance of 1.20 m.

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
temperature: 0...40°C
humidity: 10...75%
air pressure: 500...1060hPa

3 Assembly

Key

1. Transformer with spring arm and main plug
2. 5-arm base with rollers, weight, cylinder head screw and washer
3. Light head
4. Support tube
5. Main tube
Optional: UPS module
6. Tool set
7. Keys

3.1 Assembly of the flexible OR lights stand


First assemble the base with rollers and the weight (4), with the adjustable tube (3) using the screw provided. Finally, slide the support tube (2) and adjustable tube together and fix them using the set screw.

Option (if ordered): attach the UPS module to the main support tube before the adaptor and spring arm (1) are assembled.

Afterwards, the transformer with the spring arm (1) is attached to the main support tube.


1.

2.



3.


4.

3.2 Assembling the OR lights

First open the adjustable blue flange (1) on the spring arm, then lightly oil the trunnion of the surgical light with graphitic oil, introduce it into the spring arm and fix the nose-piece (2).

Then adjust the blue flange (1) so that the brake screw (3) can be inserted from below. This also secures the nose-piece(2).

3.3 Sterilisable handle

The sterilizable handle is fitted as shown, with the groove facing the lever of the surgical lights, and then clicked on and locked into place. To remove the handle, push up the lever of the surgical light and slide the handle out. The sterilisable handle is delivered non-sterile.

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air pressure: 500...1060hPa

4.2 Positioning prior to surgery

We recommend that the surgery personnel pre-position the surgical lights prior to the operation to avoid having to move the lamps during actual surgery. In this way, the surgeon only has to finely adjust the lamp using the sterilisable grip, a simple procedure thanks to the lamp's easy handling.

4.3 Switching on and off

The light is switched on or off by pressing the blue push-button on the dimmer switch. This switch regulates the strength of the light.

4.4 Adjusting the brightness

The need for more light during an operation depends on the wound depth. But to avoid blinding the surgeon during the initial process of covering and wound opening, it is advisable to use the dimmer to reduce the intensity of light before the operation. During surgery, the amount of light can then be regulated as needed. You can dim the lights by pressing the head button for a longer time.

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4.6 Centering the beam

It is advisable to use the sterilisable grip to set the lamp on the wound area before incision (centering the light). This can be simply done by turning the grip so that the angle of the projector heads is altered. This makes it possible either to concentrate the beams of lights completely together (giving the minimal light field size) or else to increase the light spot to the point where the beams of the projector heads cover an area as large as possible but there is still no visible dark spot in the center (maximum light field size). It is not normally necessary to re-position the lamps during an operation unless the distance between the wound area and the OR lights alters substantially.

4.7 Repositioning during surgery

The compact form of the surgical lights and their easy handling allow changes of angle to be made simply and as needed during an operation, largely avoiding the problem of the surgeon's working position blocking light access.

5 Care and disinfecting

Please note: The light is not watertight. To prevent liquid entering the lights, they should not be treated with an excessive amount of cleaning agent. Do not use any organic cleaning solutions.

Please note: Plastic and paint can suffer from an intensive and sustained exposure to ultraviolet light. The lights need to be protected from direct UV radiation (emitted by, for example, anti-bacterial lights or from direct sunlight, etc.).

To ensure the OR lights are correctly cleaned and disinfected, please follow these recommendations:

- use a cleaning agent in a water solution
- ensure no liquid enters the casing
- after the surface has been cleaned, wipe it dry with a soft cloth

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6 Maintenance

6.1 Maintenance of the surgical lights

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
Maintenance tasks:


- Lubrication (oiling) and adjustment as needed / adjustment of the articulated arm
- Lubrication (oiling) the bracket / adapter connection
- Lubrication (oiling) the bracket / spring arm connection
- Lubrication (oiling) the spring arm / lamp body connection
- Checking the luminous power (in lux) of the lamp. The lux value has to be assessed in relation to a light temperature of 4500 Kelvin, the size of the light field and the distribution of light within it.
- Checking the operating voltage

6.2 Regulating the spring power

- remove the cover

- regulate the spring power

6.3 Changing the LED module

Top of the light projector


Remove the cover (1) from the projector and draw the plug upwards. Using a suitable gripping device, remove the power supply plug (3) (do not pull on the cable).

Lower surface of the light projector

Remove the tension ring (4) and then use a suction lifter (5, included in product contents) to remove the protective glass (6/7/8). Dismantle the lens array by loosening the requisite screws (do not touch the lens surfaces) and remove the LED module by turning it several times to the left. (9) Take the LED module out through the top, carefully avoiding contact with the LED diodes.

Please note: The lens array can only be handled on the outside edge!


To fit a replacement LED module, simply reverse the steps described above.

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
7 Safety

Please read the following informations carefully for a safety installation and operation.

Electrical installation has to take place according to the German Standards Authority DIN VDE 0100 and local on-site regulations; for that reason, the ADMECO LUX surgical lights do not require any additional safety measures.

The ADMECO LUX surgical lights are not designed for hazardous locations where there is a danger of explosion.

If one LED-Module ceases to function, the ADMECO LUX surgical lights can safely continue to be used.

ADMECO is not liable for the safety specifications and conditions where there has been misuse of materials, or a lack of maintenance and careful handling. Any repairs or alterations needed are only to be carried out using original ADMECO parts.

Absorbed light converts into heat
„Irradiation energy“ in watts per m²

The standard states that the safety limit of a maximum value of 1000 watts per m² should not be exceeded since otherwise there is a danger of patients suffering from burns. This value will be exceeded if more than ten LED projector heads are used in any ADMECO LUX LED combination to illuminate the wound area. This calculation assumes a working distance of 1 meter. If the working distance is closer, for example, 80 centimetres, the irradiance strength is multiplied 1.56 times! A working distance set further away at 1 meter 20 cm produces an irradiance reduction of around 30%.

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8 Troubleshooting


If a problem arises in surgical lights use, the following are the most likely causes:

- a) Fault in the dimmer regulator (pulse switch)
- b) Fault in the power supply
- c) Faulty LED-Module

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9.3 Disposal and environmental protection

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10 Appendix

10.1 Model type

The model type is noted on the central joint of the OR lights.

10.2 Spare parts

The following spare parts can be ordered directly from ADMECO:

Item	Part Number
LED-module	40.20.0518
ADMECO LUX sterile grip For SL lights focusable (sterilisable)	40.20.0048
Vacuum lift Veribor 60mm: for plano-convex lens assembly	40.20.0099
Impulse switch (4 part)	40.21.0055
Spring arm nose-piece	40.20.0145
Spring arm brake screw	91.90.5008
Roller for base (without brake)	40.20.0204
Roller for base (with brake)	40.20.0203

10.3 Technical data

10.3.1 EL power supply unit

Rated values

Type:	EL
Item no.:	40.20.0519
Power:	max. 120 W
Frequency:	50/60 Hz
Primary voltage:	95... 250VAC
Secondary voltage:	24 VDC

10.3.2 Statics

Model	Description	Weight
mobile lights		
2 EL-F	LED surgical light, mobile, 2-reflector, unfocusable, 80 000 Lux (24 LEDs)	51.29kg
3 EL-F	LED surgical light, mobile,, 3-reflector, 120 000 Lux (36 LEDs)	54.10kg

10.3.3 Light data

criteria	2EL -F (no focus)	3EL -F
illuminance in 1m distance (lx)	80000	120000
irradiance	250	375
diameter light field d50 (cm)	17	
diameter light field d10 (cm)	26	
remaining illuminance with 1 mask (%)	1	1
remaining illuminance with 2 masks (%)	55	34
remaining illuminance with tube (%)	99	99
remaining illuminance with tube and 1 mask (%)	0.5	0.5
remaining illuminance with tube and 2 masks (%)	54	34
Depth of illuminance L1+L2 (m)	1.10	
Ra	86	
R9	26	

10.3.4 Electromagnetic Compability

Title	Standard	Executed	Verdict	Notes
RF disturbance	CISPR 11 / EN55011	yes	passed	Class B
RF conducted disturbances	IEC 61000-4-6 / EN 61000-4-6	yes	passed	
Radiated electromagnetic field	IEC 61000-4-3 / EN 61000-4-3	yes	passed	
Power frequency magnetic field	IEC 61000-4-8 / EN 61000-4-8	yes	passed	
Electrical fast transient burst	IEC 61000-4-4 / EN 61000-4-4	yes	passed	
Electrostatic discharge	IEC 61000-4-2 / EN 61000-4-2	yes	passed	
Surge	IEC 61000-4-5 / EN 61000-4-5	yes	passed	
Voltage dips - interruptions	IEC 61000-4-11 / EN 61000-4-11	yes	passed	
Harmonic current emissions	IEC 61000-3-2 / EN 61000-3-2	yes	passed	
Voltage fluctuations and flicker	IEC 61000-3-3 / EN 61000-3-3	yes	passed	

ADMECO AG

Turbistrasse 17
CH-6281 Hochdorf
Schweiz

phone 0041 41 914 70 30
fax 0041 41 910 61 19
email info@admeco.ch

web www.admeco.ch

The document contains 17 pages bound, numbered and sealed pages.

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