SECRETARY OF STATE
STATE OF INDIANA

Connie Lawson
Secretary of State

APOSTILLE

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1. Country: United States of America
2. This Public Document has been signed *Jamie L Ridner*
3. Acting in the capacity of notary public in & *Monroe County*
4. and bears the seal/stamp of certification clerk in & for Sec. of State Office

CERTIFIED

5. at Indianapolis, Indiana
6. this *Seventeenth* of *August, 2017*
7. by the Secretary of State of Indiana
8. No. *A2017-0817102536*
9. Seal/Stamp:
10. Signature:

Connie Lawson

Connie Lawson
Indiana Secretary of State

*Effective May 1, 2003 all apostilles from the Indiana Secretary of State will have an electronically printed seal.
This document was revised March 17, 2012.*

COOK INCORPORATED
750 DANIELS WAY
BLOOMINGTON, IN 47404 USA
PHONE: 812.339.2235 TOLL FREE: 800.457.4500
WWW.COOKMEDICAL.COM

CERTIFICATION STATEMENT

As a responsible representative of Cook Incorporated, I confirm this is an accurate version of the Instructions for Use for "Balloon Dilation Catheters for the Urinoexcretory System in Sets" Submission.

Kris Weathers
Regulatory Affairs Manager
Cook Incorporated


State of Indiana)

SS

County of Monroe)

Subscribed and sworn to before me the 15th day of August month 2017 year.

My Commission expires August 6, 2021

Instructions for Use

Balloon Dilation Catheters for the urinoexcretory system in sets

INTENDED USE


Used for urinoexcretory system dilation (see Instructions for Use).


CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner.)

DEVICE DESCRIPTION

Parameter	Conversion
1 French	1 mm = 1Fr / 3
1 Inch	25,4 mm

Balloon Catheter in set with endoscopic can and 10ml syringe (010004; 010005; 010010)

The balloon catheter is manufactured in 5.0 Fr/ 1.67 mm or 7.0 Fr/ 2.33 mm and in a length of 65 centimeters. The balloon catheter is manufactured in inflatable diameters of 5 and 6 millimeters and in lengths of 4 and 10 centimeters. The maximum balloon inflation pressures are 2.7 or 4.0 atm depending on the size of the balloon.


The proximal end of the catheter is characterized by a manifold with two lumens. One lumen serves as the inflation lumen, and the second lumen serves as the distal lumen. The inflation lumen is designed with a side angle adapter with a stopcock. The distal lumen is designed to accept a 0.028 /0.71 mm (± 0.005 inch/ 0.127 mm) or 0.038 inch /0.96mm (± 0.005 inch/ 0.127 mm). diameter wire guide. Two radiopaque markers are placed 0.061 inch/1.55 mm (± 0.009 inch/ 0.229 mm) from the proximal and distal ends of the balloon to provide fluoroscopic visibility.

Parameters of the two lumen manifold:

Length: A – 0.890 inch/22.606 mm (± 0.005 inch/0.127mm); B – 1.350 inch/ 34.29 mm (± 0.010 inch/0.254mm);

Inner diameters: C – 0.134 inch/3.404 mm (± 0.003 inch/0.076mm); D – 0.112 inch/2.84 mm (± 0.005 inch/0.127mm);

Outer diameters: E – 0.250 inch/6.35 mm (± 0.005 inch/0.127mm); F – 0.250 inch/6.35 mm (± 0.005 inch/0.127mm); G – 0.194 inch/4.93 mm (± 0.005 inch/0.127mm); H – 0.210 inch/5.33 mm (± 0.005 inch/0.127mm).

Stopcock parameters:

Length «A» - 1.594 inch/ 40.49mm (± 0.005 inch/0.127mm);

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
Inner diameter «B» - 0.098 inch/ 2.49 mm (± 0.005 inch/0.127mm);

Outer diameter: «D» - 0.430 inch/10.92mm (± 0.005 inch/0.127mm);

Luer-lock «Female»:

Inner diameter «C» - 0.169 inch/ 4.3 mm (± 0.005 inch/0.127mm);

Outer diameter «E» - 0.300 inch/7.62mm (± 0.005 inch/0.127mm).

Syringe parameters:

Syringe volume - 10 ml ($\pm 5\%$);

Graduation parameters - 1 ml (major graduations) 1/5 ml (minor graduations);

Tip length - 0.3 - 0.35 inch / 7.62 - 8.9 mm (± 0.01 inch/0.254mm);

Tip diameter - 0.071 inch / 1.8 mm (± 0.05 inch/0.127mm);


Dead space volume - 0.10 ml.

Endoscopic cap parameters:

Length: A - 0.552 inch/14.021mm (± 0.006 inch/ 0,152 mm); B - 0.039 inch/0.991mm (± 0.005 inch/ 0,127 mm); C - 0.118 inch/ 2.99mm (± 0.005 inch/ 0,127 mm); D - 0.079 inch/2.001mm (± 0.005 inch/ 0,127 mm).


Diameter: E - 0.394 inch/10.001mm (± 0.004 inch/ 0,101 mm); F - 0.197 inch/5.003mm (± 0.005 inch/ 0,127 mm); G - 0.355 inch/9.12mm (± 0.004 inch/ 0,101 mm).

The endoscopic cap is designed with a 45° ($\pm 5^\circ$) chamfer and an opening with a diameter of 0,075 inch/1,905 mm (± 0.005 inch/ 0,127 mm).

Parameters of FLLA:

Diameters: A – 0.035 inch/ 0.889 mm (± 0.002 inch/0.051mm); D – 0.210 inch/5.334mm (± 0.005 inch/0.127mm); E – 0.344 inch/ 8.738mm (± 0.005 inch/0.127mm); Length: B – 0.800 inch/ 20.32 mm (± 0.01 inch/0.254mm); C – 0.62 inch/15.748mm (± 0.005 inch/0.127mm); F – 0.172 inch/ 4.369mm (± 0.005 inch/0.127mm); G – 0.281inch/7.137mm (± 0.005 inch/0.127mm); H – 0.312 inch/7.925mm (± 0.005 inch/0.127mm).


Parameters of folding sleeve:

Order Number	A	B	C
010004	7 cm	0.133 inch/ 3.378 mm	0.113 inch/ 2.870 mm
010005	7 cm	0.107 inch/ 2.718 mm	0.087 inch/ 2.210mm
010010	14 cm	0.133 inch/ 3.378 mm	0.113 inch/ 2.870 mm

Order Number	Balloon Catheter				
	French Size Fr/mm	Catheter Length (cm)	Inflated Balloon Diameter (mm)	Balloon Length (cm)	Maximum Inflation Pressure (atm)
Tolerance	±10%	±10%	±10%	±10%	-
010005	5.0/ 1.67	65	5.0	4.0	2.7
010004	7.0/ 2.33	65	6.0	4.0	4.0
010010	7.0/ 2.33	65	6.0	10.0	4.0

Balloon Ureteral Dilator Set (010014; 010014-S1; 010014-S6; 010015)

The Balloon Ureteral Dilator Sets are comprised of a balloon catheter, wire guide, syringe and endoscopic cap.

The balloon catheter is manufactured in 5.0 Fr/ 1.67 mm or 7.0 Fr/ 2.33 mm and in a length of 65 centimeters. The balloon catheter is manufactured in inflatable diameters of 5 and 6 millimeters and in lengths of 4 and 10 centimeters. The maximum balloon inflation pressures are 2.7 or 4.0 atm depending on the size of the balloon.

The proximal end of the catheter is characterized by a manifold with two lumens. One lumen serves as the inflation lumen, and the second lumen serves as the distal lumen. The inflation lumen is designed with a side angle adapter with a stopcock. The distal lumen is designed to accept a 0.028 /0.71mm (± 0.005 inch/ 0.127 mm) or 0.038 inch /0.96mm (± 0.005 inch/ 0.127 mm). diameter wire guide. Two radiopaque markers are placed 0.061 inch/1.55 mm (± 0.009 inch/ 0.229 mm) from the proximal and distal ends of the balloon to provide fluoroscopic visibility.

The wire guide component is designed as a straight safety wire guide with a 0.028 /0.71 mm or 0.038 inch /0.96mm and measures 145 centimeters. A pin-vise handle is provided to stabilize the wire guide during introduction.

The syringe provided in this set is a 10 milliliter Luer lock syringe.

Syringe parameters:

Syringe volume – 10 ml ($\pm 5\%$);

Graduation parameters – 1 ml (major graduations) 1/5 ml (minor graduations);

Tip length – 0.3 - 0.35 inch / 7.62 – 8.9 mm (± 0.01 inch/0.254mm);

Tip diameter - 0.071 inch / 1.8 mm (± 0.05 inch/0.127mm);


Dead space volume – 0.10 ml.

Parameters of the two lumen manifold:

Length: A – 0.890 inch/22.606 mm (± 0.005 inch/0.127mm); B – 1.350 inch/ 34.29 mm (± 0.010 inch/0.254mm);

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




The endoscopic cap is designed with a 45° ($\pm 5^\circ$) chamfer and an opening with a diameter of 0,075 inch/1,905 mm ($\pm 0,005$ inch/ 0,127 mm).


Parameters of FLLA:

Diameters: A – 0.035 inch/ 0.889 mm (± 0.002 inch/0.051mm); D – 0.210 inch/5.334mm (± 0.005 inch/0.127mm); E – 0.344 inch/ 8.738mm (± 0.005 inch/0.127mm);
Length: B – 0.800 inch/ 20.32 mm (± 0.01 inch/0.254mm); C – 0.62 inch/15.748mm (± 0.005 inch/0.127mm); F – 0.172 inch/ 4.369mm (± 0.005 inch/0.127mm); G – 0.281 inch/7.137mm (± 0.005 inch/0.127mm); H – 0.312 inch/7.925mm (± 0.005 inch/0.127mm).

Parameters of folding sleeve:

Order Number	A	B	C
010014	7cm	0.133 inch/ 3.378 mm	0.113 inch/ 2.870 mm
010014-S1	14cm	0.133 inch/ 3.378 mm	0.113 inch/ 2.870 mm
010014-S6	7cm	0.133 inch/ 3.378 mm	0.113 inch/ 2.870 mm
010015	7cm	0.107 inch/ 2.718 mm	0.087 inch/ 2.210mm

Order Number	Balloon Catheter					Wire Guide	
	French Size Fr/mm	Catheter Length (cm)	Inflated Balloon Diameter (mm)	Balloon Length (cm)	Maximum Inflation Pressure (atm)	Diameter	Length (cm)
Tolerance	±10%	±10%	±10%	±10%	±10%	±10%	±10%
010015	5.0/ 1.67	65	5.0	4.0	2.7	0.028 inch (0.71 mm)	145
010014	7.0/ 2.33	65	6.0	4.0	4.0	0.038 inch (0.97 mm)	145
010014-S1	7.0/ 2.33	65	6.0	10.0	4.0	0.038 inch (0.97 mm)	145
010014-S6	7.0/ 2.33	65	6.0	4.0	4.0	0.038 inch (0.97 mm)	145

MATERIALS

Balloon Catheter in set with 10ml syringe and endoscopic cap (010004; 010005; 010010)

Device	Material	Patient Contacting
Balloon Catheter	Catheter Tubing : Polyethylene Radiopaque #30413 Markers: 5Fr -- Stainless Steel 7Fr -- Rhenium, Tungsten Re (6% or more), W (balance) Protective sheath : dark grey	Direct Contact
		Non-Contact
		Non-Contact
Luer Lock Adapter	Polyamide 6 #40340 Red pigments (for 010005): CI Pigment Brown 24 CI Solvent Red 225 CI Solvent Orange 63 Titanium Dioxide Blue pigments (for 010004 and 010010): Compound # 40340 Titanium Dioxide CI Pigment Blue 29 Zinc Stearate CI Pigment Green 7	Non-Contact

Balloon Catheter in set with 10ml syringe and endoscopic cap (010004; 010005; 010010)

(cont'd)

Device	Material	Patient Contacting
Balloon	Polyethylene Gruntzig Tubing #30040	Non-Contact
Stopcock	Body & Collar: Polycarbonate, Clear Handle: High Density Polyethylene, Blue	Non-Contact
Endoscopic cap	Silicone 8381-50-Blue	Indirect contact
Syringe	Syringle barrel: Polycarbonate #800200 Stopper – Chlorobutyl BD STOPPER Plunger - HYPAK PR3ML PSTYR	Indirect contact

Balloon Ureteral Dilator Sets (010014; 010014-S1; 010014-S6; 010015)

Component	Material	Patient Contacting
Catheter	Catheter shaft: Radiopaque polyethylene #30413 Markers: 5Fr -- Stainless Steel 7Fr -- Rhenium, Tungsten Re (6% or more), W (balance)	Direct contact
Adapter	Polyamide 6 #40340	Direct contact
Female Luer Lock Adapter	Polyamide 6 #40340 Red pigments (for 010015): CI Pigment Brown 24 CI Solvent Red 225 CI Solvent Orange 63 Titanium Dioxide Blue pigments (for 010014 (-S1/-S6): Compound # 40340 Titanium Dioxide CI Pigment Blue 29 Zinc Stearate CI Pigment Green 7	Indirect Contact
Balloon	Polyethylene #30040	Direct contact
Stopcock	Body & Collar: Polycarbonate, Clear Handle: High Density Polyethylene, Blue	Indirect contact
Wire guide	Stainless steel AISI 304 Solder: Silver All-State № 430 Wire guide coating: Teflon Dupont 850G-204	Direct contact
Pin wise handle	Polycarbonate #80110	Non-Contact
Endoscopic cap	Silicone 8381-50-Blue	Indirect contact
Syringe	Syringle barrel: Polycarbonate #800200 Stopper – Chlorobutyl BD STOPPER Plunger - HYPAK PR3ML PSTYR	Indirect contact

TECHNICAL CHARACTERISTICS

Balloon Dilation Catheters for the urinoexcretory system in sets are stable under vibration loads in the range 10-55 Hz, amplitude displacement of 0.15 mm.

Balloon Dilation Catheters for the urinoexcretory system in sets in packaging are stable under vibration loads in the range 10-55 Hz, amplitude displacement of 0.35 mm and shock loads at peak acceleration of 100 m / s² (10g), and duration of shock acceleration of 16 ms.

Balloon Dilation Catheters for the urinoexcretory system in sets are stable at temperatures in the range of +32 to +42°C.

Balloon Dilation Catheters for the urinoexcretory system in sets in packaging are stable at a transportation and storage temperatures in the in the range of -50°C to +50°C.

Bond tensile strenght (for diameters, mm): – $\geq 1,15 < 1,85$ - no less than 10 N;
 $\geq 1,85$ - no less than 15 N.

Balloon Dilation Catheters for the urinoexcretory system in sets are resistant to corrosion.

Balloon Dilation Catheters for the urinoexcretory system in sets are resistant to drastic change in the ambient temperature.

Products are supplied sterile, for one-time use and should not be re-sterilized.

INTENDED USE/INDICATION FOR USE

Balloon Catheter in set with 10ml syringe and endoscopic cap (010004; 010005; 010010) - used for ureteral dilation prior to ureteral stone manipulation or ureteroscopy, and dilation of the intramural ureter. Intended for one-time use.

Balloon Ureteral Dilator Sets (010014; 010014-S1; 010014-S6; 010015) -- used for ureteral dilation prior to ureteral stone manipulation or ureteroscopy, and dilation of the intramural ureter. Intended for one-time use.

CONTRAINDICATIONS

1. Acute urethritis
2. Acute prostatitis
3. Acute cystitis
4. Urethral abnormalities (diverticula, doubling, etc) because of the risk of perforation
5. Uncorrected coagulopathy
6. Untreated urinary infection

POTENTIAL ADVERSE EVENTS/COMPLICATIONS:

1. Urinary tract infections
2. Pyelonephritis
3. Hemorrhage

4. Trauma to the surrounding tissue (due to excessive dilation of the balloon)
5. Bleeding
6. Dysuria
7. Fever
8. Postoperative pains

INSTRUCTIONS FOR USE

1. Remove the protective sleeve from the balloon and discard. **CAUTION: Do not pre-inflate the balloon.**
 2. Lock the syringe onto the side-arm stopcock of the balloon catheter and aspirate until all air is removed. Close the stopcock.
 3. Lubricate the deflated balloon with a water-soluble lubricant.
 4. Attach the endoscopic cap to the accessory channel of the cystoscope.
 5. Under X-ray or fluoroscopic control, cystoscopically pass the wire guide the desired distance into the ureter.
 6. Pass the balloon catheter over the wire guide. Position the tip of the balloon within the distal ureter, being sure to leave a small portion of the balloon in the ureteral orifice and extending into the bladder.
- NOTE:** Alternatively, the balloon and catheter can be passed through the cystoscope as a unit, following the instructions below:
- a. Attach the pin vise handle to the female Luer lock fitting of the balloon catheter.
 - b. Place the wire guide, flexible tip leading, into the central lumen and pass it the length of the catheter until only the 3-cm flexible tip extends beyond the catheter tip.
 - c. Maintaining position of the wire guide, twist the pin vise handle clockwise to lock the wire guide in place.
 - d. Cystoscopically pass the wire guide/catheter assembly so that the tip of the balloon is within the distal ureter, being sure to leave a small portion of the balloon in the ureteral orifice and extending into the bladder.
7. Fill the syringe with dilute (30% solution) contrast medium and eliminate air from the syringe. **CAUTION: Always inflate the balloon with a sterile liquid. Never inflate with air, carbon dioxide or any other gas.**
 8. Open the stopcock on the inflation device and inflate the balloon. **NOTE:** To ensure proper regulation of balloon pressure, use of a balloon inflation device and pressure gauge is recommended. **CAUTION: Refer to product label for maximum inflation pressure. Do not exceed the maximum inflation pressure for this balloon device.**
 9. Close the stopcock for at least 30 seconds to maintain pressure and allow adequate dilation. **NOTE:** Time, rather than excessive pressure, is the key factor.
 10. Deflate the balloon by syringe aspiration. Reposition the balloon at a higher level and repeat the process until the desired location is reached.
 11. Aspirate the balloon completely, then gently remove the device. **CAUTION: Removal of the catheter is facilitated by rotating the shaft counterclockwise during withdrawal. Using excessive force to withdraw the balloon can inflict trauma to tissue and/or damage the cystoscope.**

COMPLIANCE WITH RELEVANT REGULATIONS AND TECHNICAL STANDARDS

Cook Incorporated manufactures all medical products in accordance with Medical Devices Directive 93/42/EEC. Cook Incorporated manufacture all medical products in accordance

with the United States Food and Drug Administration (FDA) Quality System Regulations, 21 CFR Part 820, *Quality System Regulations*.

All of Cook facilities are also certified to ISO 13485 (2003) *Quality Management System*.

WARRANTY INFORMATION

Manufacturer Cook Incorporated., Cook Incorporated, 750 Daniels Way, Bloomington, IN 47404, USA, phone: +812 339 2235, fax: +800 554 8335.

Authorized representative on the territory of the Russian Federation -- Company "SCHAG" (ZAO "SCHAG"). Russia, 119002, Moscow, Karmanitzki per. 9. Office 501A

Tel. (495) 956-13-09;
Fax (495) 956-1310,
e-mail: info@schag.ru

Manufacturer's warranty extends to the quality of products.

SHELF LIFE

The "use by" date as denoted by symbols on the labeling. The shelf life is shown by indication the time of manufacture (factory symbol) to the stated expiration date (hourglass symbol). The difference between both dates is the shelf life.

The shelf-life of the products falling under the scope of this technical data summary - is 3 (three) years.

STERILIZATION DETAILS

These devices are supplied sterile in peel open packages utilizing an EtO (Ethylene Oxide) validated method to a sterility assurance level (SAL) of 10^{-6} . These devices are intended for one-time use, sterile if unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile.

STORAGE RULES

These devices should be stored in a dark, dry, cool place; avoid extended exposure to light.

RULES FOR USE

Balloon Dilation Catheters for the urinoexcretory system in sets are intended for use in hospitals and health care facilities.

TRANSPORTATION RULES

Transportation shall account for the type of transportation employed, the destination, the vendor used for transportation and an estimation of the frequency and harshness of any handling.

PACKAGING

Balloon Dilation Catheters for the urinoexcretory system in sets are packaged in individual sterile peel-away pouch with instructions for use enclosed.

PACKAGING PARAMETERS


Order number	Height, cm	Width, cm	Weight, g	Tolerance
010004	50.2	20.3	61.2	±10%
010005	50.4	20.3	53.6	±10%
010010	50.2	20.3	61.2	±10%
010014	52.5	25.4	199.4	±10%
010014-S1	53.3	25.4	199	±10%
010014-S6	38.1	20.3	33	±10%
010015	52.2	25.2	198.2	±10%

LABELING

Balloon Dilation Catheters for the urinoexcretory system in sets are supplied in the individual packaging, which should show no signs of damage or opening. Upon detecting the bends, breaks and other defects the device cannot be used.

The label of the device shows the following device parameters:

- company-manufacturer and manufacturing address;
- device name (trade mark);
- serial number;
- marking about CE certification;
- marking about sterilization;
- address of the EU representative;
- trade mark of the manufacturer;
- manufacturing date;
- sizes;
- schematic picture;
- lot number;
- marking for "Single use only"
- marking for "Attention. See Instructions for Use"
- marking for consult Instructions for Use
- shelf life
- marking for phthalate

Product	Sample Label
Balloon Catheter in set with 10ml syringe and endoscopic cap (010004)	
Balloon Ureteral Dilator Sets (010014)	

DISPOSAL RULES

Upon completion of procedure, the devices are to be disposed of as per institutional guidelines for biohazardous medical waste.

For EU countries -in accordance with the World Health Organization guidance document " Safe management of waste from Health Care Activities"

For US countries – in accordance with the US Department of Labor of Occupational Safety and Health Administration (OSHA)

For Russian Federation – in accordance with national requirements. Risk class of the medical waste: class B.

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CERTIFIED

5. at Indianapolis, Indiana
6. this *Seventeenth* of *August, 2017*
7. by the Secretary of State of Indiana
8. No. *A2017-0817102537*
9. Seal/Stamp:
10. Signature:

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


Ascend AQ Ureteral Dilation Balloon Catheter in set with one-way stopcock (AUB)

The Ascend AQ Ureteral Dilation Balloon Catheter is comprised of a balloon catheter with batwing manifold. A one-way stopcock is supplied with the balloon catheter, supplied assembled. The balloon catheter is manufactured in diameter 5.0 French/ 1.67 mm with a working length of 65 centimeters. The balloon is manufactured in lengths of 4.0, 6.0 and 10 centimeters and with diameters of 4.0, 5.0, 6.0, 7.0 and 10.0 millimeter diameters. The Ascend AQ Ureteral Dilation Balloon Catheter has a maximum inflation pressure of 20 atmospheres. Nominal balloon pressure 14atm.

The distal tip of the balloon catheter is characterized by a 0.040 inch/ 1.02 mm (± 0.005 inch/ 0.127 mm) diameter endhole rounded to a smooth tip. Two radiopaque markers are placed on the catheter shaft beneath the balloon. Depending upon the specific balloon catheter, the radiopaque marker bands are spaced 4.0, 6.0 or 10 centimeters apart. Hydrophilic coating is applied to the surface of the balloon catheter shaft device for ease of access.

The proximal end of the balloon catheter is characterized by a high pressure one-way stopcock which is attached to the inflation lumen of the batwing manifold.

Stopcock parameters:

Length «A» - 1,594 inch/ 40,49mm (± 0.005 inch/0.127mm);

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
Inner diameter «B» - 0,098 inch/ 2,49 mm (± 0.005 inch/0.127mm);

Outer diameter: «D» - 0.430inch/10.92mm (± 0.005 inch/0.127mm);

Luer-lock «Female»:

Inner diameter «C» - 0,169 inch/ 4,3 mm (± 0.005 inch/0.127mm);

Outer diameter «E» - 0.300 inch/7.62mm (± 0.005 inch/0.127mm).


Parameters of batwing manifold:

Lengths: "A" – 2.275 inch/ 57.785mm (± 0.005 inch/0.127mm); "B" – 1.675 inch/ 42.545 mm (± 0.005 inch/0.127mm); "C" – 0.600 inch/15.24mm (± 0.005 inch/0.127mm); "D" – 1.575 inch/40.005 mm (± 0.005 inch/0.127mm).

Diameter: "E" – 0.136 inch/3.454 mm (± 0.005 inch/0.127mm);

Angle "F" -- 40°

Order number	Catheter		Balloon	
	Diameter Fr / mm	Length, cm	Diameter mm	Length, cm
Tolerance	$\pm 10\%$	$\pm 10\%$	$\pm 10\%$	$\pm 10\%$
AUB-4-4	5,0/1,67	65	4,0	4,0
AUB-4-6	5,0/1,67	65	4,0	6,0
AUB-4-10	5,0/1,67	65	4,0	10,0
AUB-5-4	5,0/1,67	65	5,0	4,0
AUB-5-6	5,0/1,67	65	5,0	6,0
AUB-5-10	5,0/1,67	65	5,0	10,0
AUB-6-4	5,0/1,67	65	6,0	4,0
AUB-6-6	5,0/1,67	65	6,0	6,0
AUB-6-10	5,0/1,67	65	6,0	10,0
AUB-7-4	5,0/1,67	65	7,0	4,0
AUB-7-6	5,0/1,67	65	7,0	6,0
AUB-7-10	5,0/1,67	65	7,0	10,0
AUB-10-4	5,0/1,67	65	10,0	4,0

MATERIALS

Ascend AQ Ureteral Dilation Balloon Catheter in set with one-way stopcock (AUB)

Component	Material	Patient Contacting
Catheter Tubing	Nylon № 40063	Direct Contact
Ink	PolyKing 3112 ink White Ink	Non-Contact
AQ coating	Polyvinylpyrrolidon PV05-0010	Direct contact
Balloon Tubing	Nylon #41143	Direct Contact
Marker	90% Pt, 10% Ir	Direct Contact
Batwing manifold	Polyamide 12-Compound # 40960 Red CI Pigment Red 149 CI Pigment Yellow 138 Calcium Carbonate Carbon Black	
One way stop cock	Body: Polycarbonate, Clear #80020 Rotating Connector: Polycarbonate, Clear Stem: Celcon, Grey	Non-Contact

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Bond tensile strenght (for diameters, mm): $\geq 1,15 < 1,85$ - no less than 10 N;
 $\geq 1,85$ - no less than 15 N.

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INTENDED USE/INDICATION FOR USE

Ascend AQ Ureteral Dilation Balloon Catheter in set with one-way stopcock (AUB) - used for transluminal dilation of ureteral strictures or ureteral dilation prior to ureteroscopy or stone manipulation. Intended for one-time use.

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1. Urinary tract infections
2. Pyelonephritis
3. Hemorrhage
4. Trauma to the surrounding tissue (due to excessive dilation of the balloon)
5. Bleeding
6. Dysuria
7. Fever
8. Postoperative pains

INSTRUCTIONS FOR USE

Ascend AQ Ureteral Dilation Balloon Catheter in set with one-way stopcock (AUB-)

BALLOON PREPARATION

1. Upon removal from the package, inspect the balloon catheter to ensure no damage has occurred during shipment.
2. Remove the protective sleeve from the balloon and discard. Lock a 10 ml syringe filled with 5 ml of dilute contrast medium onto the balloon lumen and apply negative pressure for 20 to 30 seconds. Release negative pressure on syringe allowing medium to draw into balloon. Detach the syringe leaving medium in the balloon lumen.

BALLOON CATHETER INTRODUCTION AND INFLATION

1. Under fluoroscopic control, cystoscopically pass a wire guide (up to .038 inch in diameter) the desired distance into the ureter beyond the area of planned dilation.
2. Activate the AQ hydrophilic coating by immersing the balloon and catheter in a sterile physiological solution for approximately 30 seconds.
3. Carefully advance the balloon catheter over the previously placed wire guide under fluoroscopic guidance utilizing the radiopaque markers to ensure proper positioning. Advance catheter so the entire balloon extends beyond the end of the cystoscope before inflating the balloon.

4. Prepare inflation device in standard manner and purge all air from syringe and tubing. To ensure proper regulation of balloon pressure, use of a balloon inflation device and pressure gauge is recommended.
5. Attach the prepared inflation device to the balloon lumen, open the stopcock on the inflation device and inflate the balloon.
6. Close the stopcock on the inflation device for at least 30 seconds to maintain pressure and allow adequate dilation. Time, rather than excessive pressure, is the key factor in transluminal dilation.

BALLOON DEFLATION AND WITHDRAWAL

1. Aspirate the balloon completely.
2. When the balloon is completely deflated, remove the device by withdrawing the catheter slowly. Removal of the catheter is facilitated by rotating the shaft counterclockwise during withdrawal. Using excessive force to withdraw the balloon can inflict trauma to tissue and/or damage the cystoscope.

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Cook Incorporated manufactures all medical products in accordance with Medical Devices Directive 93/42/EEC. Cook Incorporated manufacture all medical products in accordance with the United States Food and Drug Administration (FDA) Quality System Regulations, 21 CFR Part 820, *Quality System Regulations*.

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
- company-manufacturer and manufacturing address;

- device name (trade mark);
- serial number;
- marking about CE certification;
- marking about sterilization;
- address of the EU representative;
- trade mark of the manufacturer;
- manufacturing date;
- sizes;
- schematic picture;
- lot number;
- marking for "Single use only"
- marking for "Attention. See Instructions for Use"
- marking for consult Instructions for Use
- shelf life
- marking for phthalate

Product	Sample Label
Ascend AQ Ureteral Dilation Balloon Catheter in set with one- way stopcock (AUB-)	 

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For EU countries -in accordance with the World Health Organization guidance document "Safe management of waste from Health Care Activities"

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CERTIFIED

5. at Indianapolis, Indiana
6. this *Seventeenth* of *August, 2017*
7. by the Secretary of State of Indiana
8. No. *A2017-0817102538*
9. Seal/Stamp:
10. Signature:

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CERTIFICATION STATEMENT

As a responsible representative of Cook Incorporated, I confirm this is an accurate version of the Instructions for Use for Balloon Dilation Catheters for the Urinoexcretory System in Sets.

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


The Urethral Dilator Balloon Catheters are comprised of a balloon catheter and a syringe. Based upon the actual product, (UDBS-070029-OW) the set may also include a fixed core wire guide.

The balloon catheter has a diameter of 7.0 French/ 2.33 mm size and a working length of 29 centimeters. The distal lumen of the balloon catheter is characterized by a 0.040 inch/ 1.02 mm (± 0.005 inch/ 0.127 mm) diameter endhole, a coudé tip or a filiform tip. Because of the size of the coudé tip and filiform tip the device is labeled as a 9.0 French/ 3 mm device.

The Urethral Dilator Balloon Catheter designed for use over the wire, is manufactured with four evenly spaced ports placed on the catheter shaft inside the balloon at 4 centimeter increments.

The proximal end of the balloon catheter designed for an over the wire (UDBS-070029-OW) has a length of 18 centimeters and is designed with an inflatable diameter of 8.7 millimeters.

The Urethral Dilator Balloon Catheter designed with a coudé tip (UDBS-090029) or filiform tip (UDBS-090029-FF) manufactured in a length of 18 centimeters and is designed with an inflatable diameter of 8.7 millimeters.

Within the balloon segment of the catheter are 5 to 6 sideports placed 3 centimeters apart, then rotated at 180° there is a second set of 5 to 6 sideports. The proximal end of the balloon catheter designed with a coudé tip or a filiform tip are designed with connector cap and Luer lock adapter. A stopcock is attached to the proximal fitting assembly for inflation of the balloon.

The filiform fitting extends beyond the balloon catheter 0.243 inches/ 6.17 mm (± 0.005 inch/ 0.127 mm) and measures 0.105 inches/ 2.67 mm (± 0.005 inch/ 0.127 mm) in diameter.

The coudé tip curves upward. The distal end of the coudé tip forms a dome.

The wire guide measures 80 centimeters and measures 0.038 inches/ 0.97 mm in diameter.

The syringe is a 10 milliliter Luer lock syringe.

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Tip diameter - 0.071 inch / 1.8 mm (± 0.05 inch/0.127mm);

Dead space volume – 0.10 ml.

Order Number	Catheter		Balloon		Tip Configurations	Comments
	French Size Fr/ mm	Length (cm)	Length (cm)	Inflated Balloon Diameter (mm)		
Tolerance	$\pm 10\%$	$\pm 10\%$	$\pm 10\%$	$\pm 10\%$		
Urethral Dilation Balloon Catheter in set with wire guide and 10 ml syringe UDBS-070029-OW	7.0/ 2.33	29.0	18.0	8.7	Straight with Endhole	Includes a wire guide
Urethral Dilation Balloon Catheter in set with 10 ml syringe UDBS-090029	9.0/ 3	29.0	18.0	8.7	Coudé Tip	w/o wire guide
Urethral Dilation Balloon Catheter in set with 10 ml syringe UDBS-090029-FF	9.0/ 3	29.0	18.0	8.7	Filiform Tip	w/o wire guide

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INTENDED USE/INDICATION FOR USE

Urethral Dilation Balloon Catheters in set (UDBS-070029-OW; UDBS-090029; UDBS-090029-FF) -used to dilate the male or female urethra to treat stricture disease

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INSTRUCTIONS FOR USE

Urethral Dilation Balloon Catheters in sets (UDBS)

1. Remove protective sleeve from balloon and discard.

NOTE: when using filiform configuration, attach filiform using clockwise rotations prior to insertion.

2. Insert balloon catheter tip into the urethra and gently advance balloon segment into the urethra.

NOTE: lubrication of the balloon is optional. Use of lubricant may ease placement. Use of K-Y® Jelly or another water soluble lubricant is recommended.

NOTE: to ensure proper regulation of balloon pressure, use of a balloon inflation device and pressure gauge is recommended.*

3. Using an inflation gauge, inflate the balloon to a pressure of 90 psi, or less.

4. Inflate for 20 to 40 seconds, or until desired dilation is achieved.

5. Deflate balloon by syringe aspiration; remove and discard.

NOTE: aspirate the balloon completely before gently removing the device. Using excessive force to withdraw the balloon can inflict trauma to tissue.

INSTRUCTIONS FOR PLACING URETHRAL DILATION BALLOON CATHETER OVER-THE-WIRE

1. Remove protective sleeve from balloon and discard.

2. Insert wire guide into bladder.

3. Load catheter onto wire guide. Gently advance catheter over wire into the urethra.

NOTE: lubrication of the balloon is optional. Use of lubricant may ease placement. Use of K-Y® Jelly or another water soluble lubricant is recommended.

4. Remove wire guide.

5. Using an inflation gauge, inflate the balloon to a pressure of 90 psi, or less.

6. Inflate for 20 to 40 seconds, or until desired dilation is achieved.

7. Deflate balloon by syringe aspiration; remove and discard.

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




Urethral
Dilation
Balloon
Catheter with
10 ml Syringe
(UDBS-
090029)

Urethral Dilation Balloon Catheter With Coude' Tip

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