COPY CERTIFICATION BY DOCUMENT CUSTODIAN

I, PRIYA NAVKAR, of legal age, being first duly sworn, hereby
Name of document custodian

Swear (or affirm) that the attached reproduction of

IFU Flex
Description of documents

Is a true and exact copy of the correct and complete original document.

Priya Navkar
Signature of document custodian (affiant)

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.

State of California


County of Orange

Subscribed and sworn to (or affirmed) before me on this 9th day of Feb.

20 17 by PRIYA NAVKAR proved

to me on the basis of satisfactory evidence to be the person(s) who appeared before me.

Steve Oberlander
Signature (Seal)

INSTRUCTIONS FOR USE

PipelineTM Flex Embolization Device

Prinya Narkar

TABLE OF CONTENTS

Pipeline™ Flex Embolization Device

English	4
Français	5
Deutsch	6
Italiano	8
Español	9
Svenska	10
Nederlands	12
Português	13
Suomi	14
Dansk	16
Ελληνικά	17
Čeština	18
Magyar	20
Русский	21
Polski	23
Türkçe	24
Norsk	25
Slovenčina	27
Română	28
Български	29
한국어	31
عربي	34
Symbol Glossary	35

English Instructions for Use

EN

Pipeline™ Flex Embolization Device

CAUTION

- This device should be used only by physicians trained in percutaneous, intravascular techniques and procedures at medical facilities with the appropriate fluoroscopic equipment.
- The Pipeline™ Flex embolization device should be used by physicians who have received appropriate training for this device.
- Carefully inspect the sterile package and the Pipeline™ Flex embolization device prior to use to verify that neither has been damaged during shipment. Do not use kinked or damaged components.
- The Pipeline™ Flex embolization device is not to be used after the expiration date imprinted on the product label.
- Do not reprocess or resterilize. Reprocessing and resterilization increase the risks of patient infection and may compromise device performance.

DESCRIPTION

The Pipeline™ Flex embolization device consists of a permanent implant combined with a guide-wire-based delivery system. The Pipeline™ Flex embolization device implant is a braided, multi-alloy, mesh cylinder woven from platinum/tungsten and cobalt-chromium-nickel alloy wires. A photograph of Pipeline™ Flex embolization device is shown in Figure 1a and the design of the distal delivery system is shown in Figure 1b. The woven wires of the device provide approximately 30% metal coverage of the arterial wall surface area. The implant is designed for placement in a parent vessel across the neck of an intracranial aneurysm (IA). The expanded or un-constrained diameter is 0.25 mm larger than the labeled diameter.


The tip coil is made of platinum-tungsten alloy, the proximal bumper is platinum-iridium alloy, and the tip, distal and proximal solder joints are composed of tin-silver. The protective sleeves are designed to protect the distal portion of the braid while the Pipeline™ Flex embolization device is advanced through the micro catheter. The proximal bumper and resheathing pad allow the user to push the Pipeline™ Flex embolization device out of the micro catheter when the delivery system is advanced. The resheathing pad also allows the user to resheath the Pipeline™ Flex embolization device back into the micro catheter. The resheathing marker provides the user fluoroscopic visualization for the limit of resheathing the Pipeline™ Flex embolization device.

The Pipeline™ Flex embolization device implant is mounted on a 304 stainless steel micro guide-wire approximately 200 cm long and compressed inside an introducer sheath. The Pipeline™ Flex embolization device is designed to be delivered only through a compatible micro catheter of 0.027 inch (0.69 mm) inside diameter at least 135 cm in length.

Figure 1a: The Pipeline™ Flex embolization device

Figure 1b: The Pipeline™ Flex embolization device

1. Tip Coil
2. Proximal Bumper
3. Introducer Sheath
4. Braid
5. Fluorosafe Marker
6. Delivery Wire
7. Distal Marker
8. PTFE Sleeves
9. Resheathing Pad
10. Resheathing Marker

INDICATIONS FOR USE

The Pipeline™ Flex embolization device is intended for endovascular embolization of cerebral aneurysms.

CONTRAINDICATIONS

- Patients with active bacterial infection.
- Patients in whom antiplatelet and/or anticoagulation therapy is contraindicated.
- The Pipeline™ Flex embolization device should not be used alone as sole therapy for acutely ruptured aneurysms.
- Patients who have not received antiplatelet agents prior to the procedure.

POTENTIAL COMPLICATIONS

Potential complications include, but are not limited to, the following:

- Adverse reaction to antiplatelet/anticoagulation agents or contrast media
- Ischemia
- Blindness
- Mass effect
- Coma
- Neurological deficits
- Death
- Parent Artery Stenosis
- Device fracture
- Perforation
- Device migration or misplacement
- Perforator occlusion
- Dissection of the parent artery
- Post-procedure bleeding
- Distal embolization including to a previously uninvolved territory
- Ruptured or perforated aneurysm
- Embolism
- Seizure
- Groin injury
- Stroke
- Headache
- Thromboembolism
- Hemorrhage
- Transient Ischemic Attack (TIA)
- Hematoma or hemorrhage at the puncture site
- Vasospasm
- Hydrocephalus
- Vessel occlusion
- Infection
- Vessel perforation
- Intracerebral bleeding
- Vision impairment
- Chronic pain
- AV fistula
- Pseudoaneurysm

PRECAUTIONS

- Do not use in patients in whom the angiography demonstrates the anatomy is not appropriate for endovascular treatment, due to conditions such as severe intracranial vessel tortuosity or stenosis.
- Do not attempt to re-position the Pipeline™ Flex embolization device after deployment past the resheathing marker.
- Placement of multiple Pipeline™ Flex embolization devices may increase the risk of ischemic complications.
- The appropriate anti-platelet and anti-coagulation therapy should be administered in accordance with standard medical practice.
- A thrombosed aneurysm may aggravate pre-existing, or cause new, symptoms of mass effect and may require medical therapy.
- Do not place the Pipeline™ Flex embolization device in patients whom a pre-existing stent is in place in the parent artery at the target aneurysm location.
- Use of implants with labeled diameter larger than the parent vessel diameter may result in decreased effectiveness and additional safety risks due to incomplete foreshortening resulting in an implant longer than anticipated.
- Person with known allergy to cobalt/chromium alloy (including major elements of cobalt, chromium, nickel, molybdenum or tungsten) may suffer an allergic reaction to the Pipeline™ Flex embolization device implant.
- Person with known allergy to tin, silver, stainless steel or silicone elastomer may suffer an allergic reaction to the Pipeline™ Flex embolization device delivery system.

COMPATIBILITY

The Pipeline™ Flex embolization device is designed to be delivered through a compatible 0.027 (0.69 mm) inside diameter micro catheter. Unconstrained diameter of the Pipeline™ Flex embolization device is 0.25 mm greater than the labeled diameter (on the packaging). Do not use the Pipeline™ Flex embolization device in vessel diameters that are larger than the labeled diameter.

MAGNETIC RESONANCE IMAGING

Non-clinical testing has demonstrated that the Pipeline™ Flex embolization device is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 3 Tesla or less.
- Spatial gradient field of 720 Gauss/cm or less.
- Maximum whole-body-averaged specific absorption rate (SAR) of 4.0 W/kg for 15 minutes of scanning.

In non-clinical testing, the Pipeline™ Flex embolization device produced a temperature rise of less than 0.6°C at a maximum whole body averaged specific absorption rate (SAR) of 4.0 W/kg assessed by calorimetry for 15 minutes of MR scanning in a 3 Tesla MR 750 GE Signa 20 S System MR Scanner.

The Pipeline™ Flex embolization device may create local field inhomogeneity and susceptibility artifacts which may degrade the diagnostic quality of the MRI images. Based on the non-clinical testing of the 5.0 mm device using standard views, the worst case maximum artifact was < 4 mm when subjected to 3.0 Tesla. Local field artifact from the Pipeline™ Flex embolization device may decrease the accuracy of MR angiogram in assessing vessel luminal patency.

MR image quality may be compromised if the area is in the exact same area or relatively close to the position of the Pipeline™ Flex embolization device. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic implant.

DIRECTIONS FOR USE

- Using standard interventional radiographic technique, place the micro catheter tip at least 20 mm past the distal edge of the aneurysm. Gently retract the micro catheter to reduce slack in the micro catheter prior to inserting Pipeline™ Flex embolization device.
Note: It is recommended to use a heparinized saline drip to continuously flush micro catheter during Pipeline™ Flex embolization device use.
- Choose a Pipeline™ Flex embolization device with labeled diameter that approximates the target vessel diameter.
 - Select an appropriately sized Pipeline™ Flex embolization device such that its fully expanded diameter is equivalent to that of the largest target vessel. An incorrectly sized Pipeline™ Flex embolization device may result in inadequate device placement, incomplete opening, or migration.
 - The Pipeline™ Flex embolization device foreshortens substantially (50-60%) during deployment. Take device foreshortening into account when deploying the Pipeline™ Flex embolization device.

- Choose a Pipeline™ Flex Embolization Device with labeled length that is at least 6 mm longer than the aneurysm neck.
- Remove packaging hoop from the pouch and pull the distal end of the introducer sheath from the blue clip on the packaging hoop.
- Carefully remove system from packaging hoop until the delivery wire is exposed.
- Partially insert introducer sheath into the rotating hemostatic valve (RHV) at the catheter hub and close the RHV. Confirm back flush of the saline at the proximal end of the introducer sheath prior to advancing the Pipeline™ Flex embolization device into the micro catheter.
- Advance introducer sheath into the RHV; visually confirm the tip of the sheath is seated deeply in the hub of the micro catheter.
- Secure introducer sheath to the hub by locking down the RHV tightly.
- Advance the proximal end of delivery wire until it aligns with proximal end of introducer sheath.
- Remove the introducer sheath.

Note: The delivery wire has a fluorosafe marker no further than 125 cm from the distal end.

Caution: The fluorosafe marker is only compatible with micro catheters with a minimum length of 135 cm.

- Advance the Pipeline™ Flex embolization device into the micro catheter by pushing the delivery wire until the tip of the delivery wire aligns with the tip of the micro catheter.
Caution: If high forces or excessive friction is encountered during delivery, discontinue delivery of the device and identify the cause of the resistance, remove device and micro catheter simultaneously. Advancement of the Pipeline™ Flex embolization device against resistance may result in damage or patient injury.
Caution: The presence of other indwelling endovascular stents may interfere with proper deployment and function of the Pipeline™ Flex embolization device.
- Once the tip of delivery system and micro catheter are aligned, verify that the Pipeline™ Flex embolization device is in the desired location. Distal end of Pipeline™ Flex embolization device should be placed at least 3 mm past the distal edge of the aneurysm.
- Begin to deliver the Pipeline™ Flex embolization device using a combination of unsheathing the Pipeline™ Flex embolization device and pushing the delivery wire simultaneously.

WARNING

Pushing delivery wire without retracting the micro catheter at the same time will cause the open end braid to move distally in the vessel. This may cause damage to the braid or vessel.

- After the distal end of the Pipeline™ Flex embolization device has successfully expanded, deploy the remainder of Pipeline™ Flex embolization device by pushing the delivery wire and/or unsheathing the Pipeline™ Flex embolization device. Resheathing and/or manipulation of the micro catheter, by locking down the delivery wire and moving both as a system, may facilitate expansion of the Pipeline™ Flex embolization device.

Caution: Under fluoroscopy, carefully monitor the tip coil during Pipeline™ Flex embolization device deployment.

Caution: If the delivery wire cannot be retracted into the micro catheter, carefully remove the delivery wire and micro catheter simultaneously.


- Resheathing Instructions: During deployment of the Pipeline™ Flex embolization device resheathing can be performed by advancing the micro catheter while pulling delivery wire.

- The Pipeline™ Flex embolization device can be resheathed until the resheathing marker has reached the distal marker of the micro catheter (see Figure 2 below).
- The Pipeline™ Flex embolization device is fully resheathed when the distal marker is retracted completely inside the micro catheter. The system is designed to allow for 2 full cycles of resheathing of the Pipeline™ Flex embolization device.

WARNING

Resheathing the Pipeline™ Flex embolization device more than 2 full cycles may cause damage to the distal or proximal ends of the braid.

Figure 2. Pipeline™ Flex embolization device
(Resheathing schematic as seen under fluoroscopy,
image not to scale).

16. After the entire Pipeline™ Flex embolization device is deployed, advance the micro catheter through the device making sure not to dislodge the braid. When the micro catheter tip is distal to the Pipeline™ Flex embolization device, retract the delivery wire into the micro catheter tip.

Caution: If the catheter cannot be advanced through the Pipeline™ Flex embolization device, carefully remove the delivery wire through the Pipeline™ Flex embolization device construct.

Carefully inspect the deployed Pipeline™ Flex embolization device under fluoroscopy to confirm that it is completely apposed to the vessel wall and not kinked. If the device is not fully apposed or is kinked, consider using a balloon catheter, micro catheter, or guide-wire to fully open it.

DISPOSAL: The implant and/or delivery system should be disposed of or returned to the manufacturer per institutional guidelines.

QUESTIONS AND ANSWERS

Q If excessive friction is experienced during the insertion of delivery system at any time during the delivery of Pipeline™ Flex embolization device, what should I do?

A Carefully remove the entire system simultaneously (micro catheter and delivery system).

Q Can I retrieve the Pipeline™ Flex embolization device if the distal end of the Pipeline™ Flex embolization device has expanded at an undesirable location?

A Yes. A partially deployed Pipeline™ Flex embolization device can be resheathed per resheathing instructions, step 15 in the Directions for Use.

Q Can I retrieve a fully deployed Pipeline™ Flex embolization device?

A Once fully deployed, the Pipeline™ Flex embolization device cannot be removed. A second Pipeline™ Flex embolization device can be deployed if needed.

Q Can I place a second Pipeline™ Flex embolization device inside another Pipeline™ Flex embolization device?

A Yes. A second Pipeline™ Flex embolization device can be placed inside another Pipeline™ Flex embolization device. After placing the first Pipeline™ Flex embolization device, advance the micro catheter over the delivery wire while keeping the delivery core wire across the Pipeline™ Flex embolization device. Position the micro catheter at the desired location and retrieve the delivery wire. Select a new appropriate Pipeline™ Flex embolization device and deploy it as normal.

Caution: Placement of multiple Pipeline™ Flex embolization devices may increase the risk of ischemic complications.

Q If there is a difference between the proximal and distal diameter, which Pipeline™ Flex embolization device diameter do I choose?

A Choose a Pipeline™ Flex embolization device that matches larger (typically proximal) vessel diameter to ensure proper anchoring.

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The figure consists of three horizontal diagrams illustrating the repair process of a concrete beam.
 - The top diagram shows a beam with a crack (1) and a repair area (2). Labels 3, 4, 5, and 6 point to different sections of the beam.
 - The middle diagram shows the application of repair material (7) and reinforcement (8) to the beam.
 - The bottom diagram shows the final repaired beam with reinforcement (8) and repair material (7).

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Priya Navkar

Signature of document custodian (affiant)

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State of California

County of Orange


Subscribed and sworn to (or affirmed) before me on this 9th day of Feb.
20 17 by PRIYA NAVKAR proved

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
- Select an appropriately sized Pipeline Flex embolization device such that its fully expanded diameter is equivalent to that of the largest target vessel. An incorrectly sized Pipeline Flex embolization device may result in inadequate device placement, incomplete opening, or migration.

- The Pipeline Flex embolization device foreshortens substantially (50-60%) during deployment. Take device foreshortening into account when deploying the Pipeline Flex embolization device.

3. Choose a Pipeline Flex Embolization Device with labeled length that is at least 6 mm longer than the aneurysm neck.

4. Remove packaging hoop from the pouch and pull the distal end of the introducer sheath from the blue clip on the packaging hoop.

5. Carefully remove system from packaging hoop until the delivery wire is exposed.

6. Partially insert introducer sheath into the rotating hemostatic valve (RHV) at the catheter hub and close the RHV. Confirm back flush of the saline at the proximal end of the introducer sheath prior to advancing the Pipeline Flex embolization device into the micro catheter.

7. Advance introducer sheath into the RHV; visually confirm the tip of the sheath is seated deeply in the hub of the micro catheter.

8. Secure introducer sheath to the hub by locking down the RHV tightly.

9. Advance the proximal end of delivery wire until it aligns with proximal end of introducer sheath.

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11. Advance the Pipeline Flex embolization device into the micro catheter by pushing the delivery wire until

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
15. Resheathing Instructions: During deployment of the Pipeline Flex embolization device resheathing can be performed by advancing the micro catheter while pulling delivery wire.

- The Pipeline Flex embolization device can be resheathed until the resheathing marker has reached the distal marker of the micro catheter.
- The Pipeline Flex embolization device is fully resheathed when the distal marker is retracted completely inside the micro catheter. The system is designed to allow for 2 full cycles of resheathing of the Pipeline Flex embolization device.

Warning

Resheathing the Pipeline Flex embolization device more than 2 full cycles may cause damage to the distal or proximal ends of the braid.

Resheathing schematic as seen under fluoroscopy, image not to scale:

1. Proximal End of the device
2. Micro Catheter
3. Micro Catheter Distal Marker

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4. Resheathing Marker
5. Proximal Bumper
6. Delivery Wire
7. Resheathing Limit
8. Device Detached

16. After the entire Pipeline Flex embolization device is deployed, advance the micro catheter through the device making sure not to dislodge the braid. When the micro catheter tip is distal to the Pipeline Flex embolization device, retract the delivery wire into the micro catheter tip.


Caution: If the catheter cannot be advanced through the Pipeline Flex embolization device, carefully remove the delivery wire through the Pipeline Flex embolization device construct.

17. Carefully inspect the deployed Pipeline Flex embolization device under fluoroscopy to confirm that it is completely apposed to the vessel wall and not kinked. If the device is not fully apposed or is kinked, consider using a balloon catheter, micro catheter, or guide-wire to fully open it.

Disposal: The implant and or delivery system should be disposed of or returned to the manufacturer per institutional guidelines.

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minutes of MR scanning in a 3 Tesla MR 750 GE Signa 20.0 System MR Scanner.

The Pipeline Flex embolization device may create local field inhomogeneity and susceptibility artifacts which may degrade the diagnostic quality of the MRI images. Based on the non-clinical testing of the 5.0 mm device using standard views, the worst case maximum artifact was < 4 mm when subjected to 3.0 Tesla. Local field artifact from the Pipeline Flex embolization device may decrease the accuracy of MR angiogram in assessing vessel luminal patency.

MR image quality may be compromised if the area is in the exact same area or relatively close to the position of the Pipeline Flex embolization device. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic implant.

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Figure 1. The Pipeline Flex embolization device


Figure 2. The Pipeline Flex embolization device

- | | | |
|--------------------|----------------------|--------------------|
| 1. Tip Coil | 4. Braid | 8. PTFE Sleeves |
| 2. Proximal Bumper | 5. Fluorosafe Marker | 9. Resheathing Pad |

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3. Introducer Sheath 6. Delivery Wire 10. Resheathing Marker


Figure 3. Pipeline Flex Embolization Device

8	Wire, distal bump
9	Coil, distal tip
10	Restraints
11	Restraints
12	Subassembly, distal protective
13	PED-48 wire braid subassembly
14	Subassembly, PTFE covered proximal delivery system
15	Restraint, proximal pad
16	Introducer sheath, tapered
17	Subassembly, pad

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Each outer package contains the device in sterile package, instructions for use, MRI-card and sticker (identification label), on which are written the product type and lot number. This sticker should be pasted in the patient's medical card to document the exact details of the implant. The shipping box is embedded multilingual instructions for use in paper form. Membrane cover of sterile primary package is made of Tyvek material.

The device is provided sterile and non-pyrogenic. The device and its packaging should be carefully examined prior to use to verify that no damage has occurred during shipment. The device should not be used if the pouch is open or damaged. The device is intended for single use only and should be discarded after use.

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The expanded or unconstrained diameter is larger than the labeled diameter on	0.25	mm
Number of braid wires	48	pcs
Wire length	1980±28	mm
Wire diameter	0.508 ^{+0.013} _{-0.018}	mm
Tip coil length	15	mm
Resheathing Pad length	2.10 – 2.50	mm
Distal marker length	0.40±0.01	mm
Proximal bumper length	1±0.1	mm
Resheathing Marker length	0.05±0.01	mm
Outer introducer diameter	1.143±0.025	mm
Inner introducer diameter	0.686±0.013	mm
Introducer length	1060±13	mm
Packaging Hoop length	1562±13	mm
Coverage of the arterial wall surface area, no less	30	%
Wire break force	>3.1	N
Resistance to bends	Radius 9 mm	
The radial force, not more than	0.07	N/mm
Compatible microcatheter inner diameter	0.69	mm
Compatible microcatheter length	1350	mm
Protective sleeves width	9±3	mm
Maximum force, required for catheter delivery	≤9.9	N
Radial pressure over contact surface when device implanted in the largest recommended vessel size	>33.3	Pa
Coil bending angle	55±10	°

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ISO 11607-1:2009
ISO 11607-2:2006

Sterilization cycle № 70 parameters (CS-012)

Characteristic	Specification
Process Location	Sterigenics 4900 S. Gifford Avenue, Los Angeles CA 90058
Vessel ID No.	B
Cycle control program	Check Value 4607/Safe
Gas inject by pressure rise (Ethylene Oxide)	Inject gas to a pressure of at least 36.23 kPa and a maximum of 39.62 kPa.
Maintain pressure with inert during EO gas dwell	Maintain pressure at a setpoint of 78.56 kPa to a maximum of 83.64 kPa. Minimum allowed pressure is 76.87 kPa
Gas Dwell Time	Gas phase is a Minimum of 1 hour to 1 hour 30 minutes
Gas Dwell Temperature	Maintain temperature during gas dwell a minimum of 54.4°C to 60°C

Sterilization cycle № 570 parameters (CS-007)

Characteristic	Specification
Process Location	Sterigenics 4900 S. Gifford Avenue, Los Angeles CA 90058
Vessel ID No.	#7

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11. Advance the Pipeline Flex embolization device into the micro catheter by pushing the delivery wire until the tip of the delivery wire aligns with the tip of the micro catheter.

Caution: If high forces or excessive friction is encountered during delivery, discontinue delivery of the device and identify the cause of the resistance, remove device and micro catheter simultaneously. Advancement of the Pipeline Flex embolization device against resistance may result in damage or patient injury.

Caution: The presence of other indwelling endovascular stents may interfere with proper deployment and function of the Pipeline Flex embolization device.

12. Once the tip of delivery system and micro catheter are aligned, verify that the Pipeline Flex embolization device is in the desired location. Distal end of Pipeline Flex embolization device should be placed at least 3 mm past the distal edge of the aneurysm.

13. Begin to deliver the Pipeline Flex embolization device using a combination of unsheathing the Pipeline Flex embolization device and pushing the delivery wire simultaneously.

Warning

Pushing delivery wire without retracting the micro catheter at the same time will cause the open end braid to move distally in the vessel. This may cause damage to the braid or vessel

14. After the distal end of the Pipeline Flex embolization device has successfully expanded, deploy the remainder of Pipeline Flex embolization device by pushing the delivery wire and/or unsheathing the Pipeline Flex embolization device. Resheathing and/or manipulation of the micro catheter, by locking down the delivery wire and moving both as a system, may facilitate expansion of the Pipeline Flex embolization device.

Caution: Under fluoroscopy, carefully monitor the tip coil during Pipeline Flex embolization device deployment.

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
1. Proximal End of the device
2. Micro Catheter
3. Micro Catheter Distal Marker
4. Resheathing Marker
5. Proximal Bumper

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6. Delivery Wire
7. Resheathing Limit
8. Device Detached

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
may degrade the diagnostic quality of the MRI images. Based on the non-clinical testing of the 5.0 mm device using standard views, the worst case maximum artifact was < 4 mm when subjected to 3.0 Tesla. Local field artifact from the Pipeline Flex embolization device may decrease the accuracy of MR angiogram in assessing vessel luminal patency.

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- | | | |
|----------------------|----------------------|------------------------|
| 1. Tip Coil | 4. Braid | 8. PTFE Sleeves |
| 2. Proximal Bumper | 5. Fluorosafe Marker | 9. Resheathing Pad |
| 3. Introducer Sheath | 6. Delivery Wire | 10. Resheathing Marker |
| 7. Distal Marker | | |

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Figure 4. (Left) PFED proximal bumper and resheathing pad.
(Right) Protective sleeves.


Figure 5. Pipeline Flex delivery system


Figure 6. Pipeline Flex delivery system with braid

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Vessel ID No.	B
Cycle control program	Check Value 4607/Safe
Gas inject by pressure rise (Ethylene Oxide)	Inject gas to a pressure of at least 36.23 kPa and a maximum of 39.62 kPa.
Maintain pressure with inert during EO gas dwell	Maintain pressure at a setpoint of 78.56 kPa to a maximum of 83.64 kPa. Minimum allowed pressure is 76.87 kPa
Gas Dwell Time	Gas phase is a Minimum of 1 hour to 1 hour 30 minutes
Gas Dwell Temperature	Maintain temperature during gas dwell a minimum of 54.4°C to 60°C

Sterilization cycle № 570 parameters (CS-007)

Characteristic	Specification
Process Location	Sterigenics 4900 S. Gifford Avenue, Los Angeles CA 90058
Vessel ID No.	#7
Cycle control program	Check Value 34433/Safe
Gas inject by pressure rise (Ethylene Oxide)	Inject gas to a pressure of at least 37.93 kPa and a maximum of 39.62 kPa.
Maintain pressure with inert during EO gas dwell	Maintain pressure at a setpoint of 78.56 kPa to a maximum of 81.95 kPa. Minimum allowed pressure is 36.23 kPa

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