

ATTESTATION.

STATE OF TEXAS
TARRANT COUNTY

On this 13th day of February 2017, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following document(s), Integrated Power Console (IPC®) System User's Guide, issued for Medtronic Powered Surgical Solutions, presented to me by the document's custodian, Jenna Groves, and, to the best of my knowledge, that the photocopied document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

Official Notary Signature

United States of America
State of Texas
Tarrant County

Subscribed and sworn to before me
This 13th day of February, 2017.

Integrated Power Console (IPC®) System

MODEL: 1898001

User's Guide

RxOnly

CUSTOMER SERVICE

For further information regarding the use of this product or to report any problems, please contact Medtronic using the appropriate information provided on the blue and white contact information card packaged with each device; or contact your local distributor.

Medtronic Xomed, Inc.	US Help Line
6743 Southpoint Drive North	800 874 5797
Jacksonville, Florida 32216 USA	
www.medtronic.com	International Service
	International Customers should contact their local Medtronic Xomed office.

The following are trademarks or registered trademarks of Medtronic, Inc. in the United States and other countries: CD HORIZON®, Endo-Scrub®, Hydrodebrider®, Indigo™, Intelliflow®, IPC®, Legend®, Legend EHS®, Legend EHS Stylus®, Magnum®, Midas Rex®, Mednext®, NIM®, NIM-ECLIPSE®, POWEREASE®, Skeeter®, SOLERA®, StraightShot®, Stylus Touch®, Triton®, TSRH® 3Dx™, Visao® and XPS®. All other trademarks, service marks, registered trademarks, or registered service marks are the property of their respective owners in the United States and other countries.

Contents

- GLOSSARY..... 1-1
- INDICATIONS FOR USE 1-1
- DEVICE DESCRIPTION 1-1
- CONTRAINDICATIONS 1-1
- WARNINGS..... 1-1
 - System Warnings..... 1-1
 - Component Warnings..... 1-2
 - Disposable Warnings 1-2
- PRECAUTIONS 1-3
- SYSTEM REQUIREMENTS AND SPECIFICATIONS..... 1-4
 - Console Specifications 1-4
 - Power Cord Product Numbers 1-4
- SYSTEM SOUNDS AND FIGURES 1-5
- PRE-OPERATING INSTRUCTIONS 1-8
 - When the System Arrives 1-8
 - Set up the IPC 1-8
 - Install the Pump Cartridges or Irrigation Tubing 1-8
 - Prepare IPC for Use 1-8
 - Calibrate Touchscreen..... 1-8
 - Change System Settings 1-9
 - Set up and Prime Pumps 1-10
 - Confirm System Operation 1-10
- IPC COMPONENTS 1-11
 - Auxillary Power to Console..... 1-11
 - Multifunction Foot Pedal..... 1-11
 - Y-Splitter 1-11
 - IntelliFlow Irrigation Remote Control 1-11
- GUIDANCE AND MANUFACTURER’S DECLARATION – ELECTROMAGNETIC IMMUNITY 1-12
- LIMITED WARRANTY..... 1-14
- FOR ITEMS CONTAMINATED WITH TSE AGENTS 1-14
- TRITON ELECTRIC HIGH-TORQUE HANDPIECE..... 2-1
- SUCTION IRRIGATOR 3-1
- ENDO-SCRUB 2 4-1
- SPINE SHAVER (SC1) HANDPIECE 5-1
- STRAIGHTSHOT M5, STRAIGHTSHOT M4, STRAIGHTSHOT MAGNUM II, AND STRAIGHTSHOT III..... 6-1
- LEGEND EHS AND LEGEND EHS STYLUS 7-1
- STYLUS TOUCH..... 8-1
- LEGEND EHS, LEGEND EHS STYLUS AND STYLUS TOUCH ATTACHMENTS 9-1
- SKEETER ULTRA-LITE OTO-TOOL 10-1
- VISAO HIGH-SPEED DRILL 11-1
- INDIGO HIGH-SPEED OTOLOGIC DRILL..... 12-1
- MIDAS REX MICROSAWS 13-1
- POWEREASE DRIVER 14-1

TROUBLESHOOTING AND ERROR CODES..... A-1

CLEANING AND STERILIZATION..... B-1

 Post-Operative Instructions.....B-1

 Triton Electric High-Torque HandpieceM000030A322

 Endo-Scrub 2 Finger Switch.....68E4005

 StraightShot M5, StraightShot M4, StraightShot Magnum II, StraightShot III, or Midas Rex Spine Shaver 68E3282

 Legend EHS and Legend EHS StylusM000030A234

 Stylus Touch.....68E4132

 Legend Attachments/TubesM000030A235

 Skeeter Oto-Flex Burs68E3968

 Skeeter Handpiece68E3969

 Visao68E3281

 Indigo High-Speed Otologic Drill68E4187

 Indigo High-Speed Otologic Drill Attachments68E4188

 Microsaws.....M000030A231

 POWEREASE Driver68E4189

SYMBOLS C

GLOSSARY

The following words and acronyms may be used in this guide.

FCU	Foot Control Unit
FWD	Forward - Rotation is clockwise
I.V.	Intravenous
IPC	Integrated Power Console
NIM	Nerve Integrity Monitor - One or all of the Nerve Integrity Monitor units: NIM-Response 2.0, NIM-Neuro 2.0, NIM-Response 3.0 and NIM-Neuro 3.0
NIM-ECLIPSE	Nerve Integrity Monitor for spinal surgeries
OSC	Oscillate
REV	Reverse - Rotation is counter-clockwise
XPS	Xomed Power System

INDICATIONS FOR USE

The IPC is indicated for the incision/cutting, removal, drilling and sawing of soft and hard tissue and bone in Head & Neck/ENT (Otologic, Neurologic, Neurotologic, Sinus, Rhinologic, Nasopharyngeal/Laryngeal), Oral/Maxillofacial and Plastic/Reconstructive/Aesthetic surgical procedures.

DEVICE DESCRIPTION

The IPC System is a powered microdebrider, drill and saw system that will remove soft tissue, hard tissue and bone during surgical procedures. The system consists of a power control console, foot pedal, connection cables and assorted handpieces to drive various burs, blades, drills, rasps, cannulae and saws. It includes integrated irrigation pumps for irrigation of blades, burs and for motor coolant.

In addition to the handpieces and pumps there is a connection for continuous stimulation of the Visao straight burs that enables nerve integrity monitoring during surgical procedures. The Nerve Integrity Monitor (NIM) is a separate device that stimulates and monitors the nerve. This system has connections that allow the NIM to be connected with the Visao handpiece and Stimulating Bur Guard, enabling the NIM to stimulate and monitor the nerve at the surgical site.

This device is intended for use by physicians trained in the procedures described.

The system can be used to clear the end of a rigid rod endoscope in order to maintain good visualization of endoscopic procedures without having to remove the scope from the surgical site.

CONTRAINDICATIONS

None.

WARNINGS

System Warnings

- W1 It is important that the IPC system operator be familiar with the system User's Guide, its precautions, procedures and safety issues.
- W2 Do not use the IPC POWEREASE system in the presence of flammable anesthetics. Avoid potential ignition or explosion of gases.
- W3 When not operating handpiece, eliminate accidental foot control activation. Control energy to and through the handpiece to prevent unintended tissue, bone, or nerve resection.
- W4 Disconnect power to the IPC system before cleaning the unit to avoid electrical macro shock.
- W5 Do not attach unapproved components to the IPC system to avoid electrical macro shock.
- W6 To avoid the risk of electrical shock, achieve electrical grounding reliability with proper connections. Connect the IPC system to hospital grade receptacles only.
- W7 This medical device complies with EN60601-1-2 safety standard for electromagnetic compatibility, requirements and test. However, if this equipment is operated in the presence of high levels of electromagnetic interference (EMI) or highly sensitive equipment, interference may be encountered and the user should take whatever steps are necessary to eliminate or reduce the source of the interference. Diminished performance may lengthen operating time for the anesthetized patient.
- W8 Medical Electrical Equipment needs special Precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this Guide.
- W9 Portable and mobile RF communications equipment can affect Medical Electrical Equipment.
- W10 Do not operate the IPC POWEREASE system in the presence of Magnetic Resonance Imaging devices.
- W11 Use of accessories and cables other than those specified and sold by Medtronic may result in increased emissions and decreased immunity of this unit.
- W12 The IPC system should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the IPC system should be observed to verify normal operation in the configuration in which it will be used.
- W13 Do not attempt to run the IPC POWEREASE system handpiece immediately after autoclaving. Allow an adequate "cool down" period (Typically 1 hour).
- W14 Consult the Legend Bone Mill product insert before use with the IPC system.
- W15 For metal transection, observe the following safety precautions:
 - W15a Eye wear protection is essential.
 - W15b Irrigate well to cool the cutting surfaces.
 - W15c Protect the wound site from metal debris.
 - W15d Use a clamp or grasping device to control loose fragments during transection of any metal component.
- W16 Do not operate the IPC POWEREASE system without eye protection.
- W17 All service must be performed by Medtronic qualified personnel only.
- W18 Repair and/or modification to the IPC system by anyone other than qualified service personnel may significantly compromise the unit's ability to perform effectively and/or void the equipment warranty.

INTEGRATED POWER CONSOLE (IPC)

Component Warnings

- W19 Do not use any parts other than Medtronic system components as damage or substandard performance could result.
- W20 Always inspect the components before and after use for any damage. If damage is observed, do not use damaged part until it is repaired or replaced. Damaged parts may deposit metal shavings on surgical site.
- W21 When precise location of blade tip is required, engage the rotation lock on the handpiece, then calibrate and verify the blade tip on Image Guided Surgery (IGS) system. Always lock M4 handpiece when driving non-rotatable blades to maintain their IGS calibration.
- W22 Employ visualization, including use of imaging techniques (for example: fluoroscopy, image guided surgery) when using rotating powered accessories. Discontinue powered application in the event of lack of visualization of surgical site.
- W23 Midas Rex Variable Exposure attachments: Surgeons should familiarize themselves with the performance of dissecting tools before use, and should explore the effect of various levels of tool exposure on dissection stability. If the tool exhibits excessive chatter, vibration, or movement, decrease the tool exposure.
- W24 Motors and attachments may fail due to extended use and allow a component to detach and fall from the motor or attachment, causing patient injury.
- W25 Electrical contacts must be dry prior to use.
- W26 Heavy side loads and/or long operating periods may cause the device to overheat.
- W27 Do not use an overheated device, as it may cause thermal injury to the patient or operator.
- W28 Use adequate irrigation. The use of a tool without irrigation may cause an inordinate amount of heat buildup resulting in a thermal injury to tissue. Depending on the amount of irrigation used, the drill bits and saw blades can achieve temperatures in excess of 50 °C.
- W29 Do not attempt to change a dissecting tool, saw blade, or attachment while the motor is running, or when the motor or attachment is in an overheated state.
- W30 Do not immerse the system components, except as noted.
- W31 Do not place motor, attachment and tool on the patient or in an unsecured location during surgery.
- W32 A system that is not functioning properly should not be used until all necessary repairs have been made and the unit is tested to ensure that it is functioning in accordance with Medtronic specifications.
- W33 Match the nomenclature and color code on the tool packaging to the same nomenclature and color code on the attachment.
- W34 Make sure that the attachment is still in the locked position after each adjustment of the tool exposure, as attempting to increase the tool exposure too far, may result in the attachment accidentally being unlocked.
- W35 Midas Rex Legend EHS Motor and Midas Rex Legend EHS Stylus Motor should only be operated when the attachment is in the locked position. Attempting to operate the Midas Rex Legend EHS Stylus Motor when the attachment is in the unlocked position may result in the motor stalling.
- W36 Smoke and/or excessive heat may be generated if attachment is not in the fully locked position. This may result in thermal injury to the surgeon or staff.
- W37 The Indigo and Legend EHS motors will not run properly unless the attachment is in the locked position.
- W38 DO NOT change accessory with handpiece running to prevent laceration of user and cross-contamination through compromised glove.
- W39 Remove Legend Footed Attachments cautiously and slowly as per instructions to avoid injury to the operator.
- W40 DO NOT modify accessories used with the handpiece. Performance could be diminished with modified accessories.
- W41 The safe use of the Endo-Scrub 2 System in procedures where surgical lasers are also employed has not been clinically demonstrated.
- W42 In order to ensure compliance with requirements of IEC 60601-1, use a Medtronic approved power cable.
- W43 To avoid the risk of electric shock, this equipment must only be connected to a supply main with protective earth.
- W44 Keep NIM Muting Probe cable away from IPC system cables.
- W45 Verify reusable device was sterilized prior to use. If not sterilized, do not use.

Disposable Warnings

- W46 Tools are available for resection of soft tissue and bone for surgical procedures. Use of tools depends on the intended application and patient needs. Sharp-cutting powered tools induce bleeding and removal of significant tissue and bone.
- W47 Use methods at the operative site to control bleeding that do not compromise patient safety during at-risk surgery.
- W48 Always keep the cutting area of the tool/saw blade away from fingers and loose clothing. Prevent laceration of user and cross-contamination through compromised glove.
- W49 Operate the tool only after the appropriate anatomical landmarks and the intended surgical site have been confirmed.
- W50 Use care in application of the moving cutting end to only appropriate anatomical landmarks and the intended surgical site when using powered accessories.
- W51 Insertion of metal objects in accessory tip may cause the accessory to break leaving fragments in the wound. The fragments may be difficult to remove, causing irritation, inflammation and foreign-body response at surgical site.
- W52 Bending or prying may break the accessory, causing harm to patient or staff.
- W53 Do not use excessive force to pry or push bone with the attachment, tool or blade during dissection.
- W54 A tool's size and geometry may create excessive vibration at certain speeds. Increase or decrease speed on console. Change to a new tool to prevent unintended tissue removal from patient.
- W55 Test for wobble at desired speed prior to use. Discontinue use of accessory if tip begins to wobble and replace accessory to prevent unintended tissue removal from patient.
- W56 Eccentricity of the tool can cause tool vibration and may result in excess tissue and bone destruction and hearing damage.
- W57 Excessive noise from the tool when drilling close to the cochlea or ossicular chain may cause hearing damage.
- W58 CONSULT the cranial perforator device labeling for the recommended speed specifications.
- W59 Tools with "L" identification are longer tools intended for light bone dissection. The increased tool head/stem configuration may affect dissection stability.
- W60 Tool flutes and blade teeth are sharp and may perforate surgical gloves. Tools/blades may be grasped with a hemostat to aid in installation and removal.
- W61 DO NOT attempt to resharpen used tools. Worn tools should be replaced with new ones frequently to ensure effective cutting and control.
- W62 Carefully inspect tool both prior to and following each use for signs of excessive wear, fragmentation, eccentricities or other defects. Replace any suspicious tools with a new one prior to use.
- W63 Excessive pressure applied to bur may cause bur fracture. Should a tool fracture in use, extreme care must be exercised to ensure that all fragments of the tool are retrieved and removed from the patient. Unremoved tool fragments may cause tissue damage to the patient.
- W64 Do not use metal-cutting tools on bone.
- W65 Use only rotary tools specifically designed for use with this drill system.
- W66 When using non-rotatable tools, ensure rotation lock is engaged to prevent inadvertent rotation.
- W67 The use of powered reciprocating instruments may result in vibration / related injury.
- W68 Powered blades should be operated in the Oscillate mode only. Operating in the Forward mode may cause damage to the blade.
- W69 Do not attempt to sterilize disposable devices. The disposables are packed sterile and are not intended for repeat use. To prevent contamination, use only once.
- W70 Any tubing or other tip protectors used during shipping must be removed prior to cleaning and sterilization.
- W71 Do not use accessory if package is opened or damaged. Broken seal offers no protection against cross-contamination.

- W72 Properly dispose of single-use devices removed from sterile packages. Devices lose sterility upon removal from packaging.
- W73 Do not use dull, damaged or bent tools. Use of dull tools can reduce handpiece effectiveness and cause the handpiece temperature to increase.
- W74 T&A Blades: Gently remove the inner tube from the outer tube. The inner tube may elongate upon removal from the outer tube. If this occurs, the inner tube may not lock properly into the handpiece or the blade may not work properly.
- W75 T&A Blades: Rotate the inner tube when removing and inserting it in the outer tubes to prevent damage to the internal seal. If the seal is damaged, the blade will leak at the handpiece.
- W76 Always ensure that the drill is securely engaged into the handpiece prior to operating the system.
- W77 Always examine operation of each tool in a handpiece before use.
- W78 Powered burs and drills should be operated in the Forward mode only.
- W79 This system requires insulated connectors for the StraightShot M4 Microdebrider, StraightShot Magnum II Microdebrider, StraightShot III Microdebrider, Midas Rex SC1, Visao, or Skeeter handpieces and the Multi Function Foot Control Unit.
- W80 Sterilize and dry reusable device before storing the system. Decrease likelihood of cross-contamination with timely sterilization.
- W81 After each procedure, properly clean all reusable system components.
- W82 Auxiliary Power Outlet with protective cover is for use with the Hydrodebrider or Bone Mill consoles only.
- W83 Place Stylus Touch in safe mode while not in use.
- W84 Do not place Stylus Touch handpiece in the proximity of magnetic field, such as magnetic drape and MRI equipment, to avoid inadvertent handpiece activation.
- W85 Do not apply excessive side loading. Excessive side loading could cause angled attachments to unlock accidentally from motor.

PRECAUTIONS

- P1 PRIME/FLUSH Priming is a feature designed to purge air out of the tubing set(s) during setup. The first time a Prime or Flush button is pressed it will turn on pump 1 and/or 2 long enough to purge air out of the tubing set(s). Turning power Off and On resets the PRIME feature. Once pressed all Prime buttons will change to Flush buttons.
- P2 To prevent damage to curved tools, disconnect suction tube prior to changing tool during procedure.
- P3 When using an angled attachment, hold the handpiece assembly by the attachment so that the attachment does not inadvertently loosen from the handpiece.
- P4 For Legend tools only:
If a tool package is opened, but the tool is not used or contaminated, the tool can be re-sterilized. Remove tool from original packaging and place into an approved autoclave package. Steam sterilize as follows:
High-Vacuum Steam 132 °C for 5 minutes
Gravity Displacement 132 °C for 15 minutes
The re-sterilized tool must be used promptly following re-sterilization. If rust or corrosion is encountered after re-sterilization, do not use the re-sterilized tool.
- P5 DO NOT run the 16-MF attachment with operating speed above 62,000 rpm. This may cause over heating and damage to internal gears of attachment.
- P6 DO NOT use twist drill or Contra-Angle tool at an operating speed over 62,000 rpm.
- P7 Do not attempt to disconnect the cable from the Midas Rex Legend EHS Stylus Motor.
- P8 Do not kink cables. Inspect cables and pins for cracks, tears or corrosion.
- P9 Do not use anti-fog on scope or sheath, as weeping or leaking may result.
- P10 Disconnect cable from Midas Rex Legend EHS motor prior to sterilization.
- P11 The use of a washer-disinfector for cleaning may cause a pre-mature degradation in performance.
- P12 Remove devices from instrument case before placing into washer disinfector and allow devices to drain.
- P13 Orient devices in the washer-disinfector by following manufacturer recommendations.
- P14 DO NOT use low-temperature hydrogen peroxide gas plasma sterilization due to the lumen internal diameter and length restrictions.
- P15 DO NOT use low-temperature liquid peracetic acid sterilization due to immersion procedure.
- P16 DO NOT steam or EO sterilize the Legend Attachment Cleaning Nozzle.
- P17 Remove and discard accessories following local regulations for proper disposal of contaminated materials.
- P18 Disposable devices are for single-use only.
- P19 Clean the motor and cable while still connected together. This will help to reduce ingress of debris.
- P20 Use ONLY recommended cleaning agents.
- P21 Do not use excessive force to insert the endoscope into the Endo-Scrub 2 sheath. This will damage the endoscope as well as the Endo-Scrub 2 sheath.
- P22 If the endoscope tip can be seen extending beyond the tip of the Endo-Scrub 2 sheath, then the sheath has been damaged. Damaged product must be immediately discarded.
- P23 When using a Y-Splitter, only one Multifunction Foot Pedal shall be active at a time.
- P24 Only one Y-Splitter shall be used at a time.
- P25 Ensure the blade or bur is fully engaged in the microdebrider. For blades, verify the tip is fully engaged with the outer cannula prior to use.

SYSTEM REQUIREMENTS AND SPECIFICATIONS

Console Specifications

Functional Standards for Electrical Systems

ANSI/AAMI ES60601-1	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance	2005
IEC 60601-1	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance	2005
EN 60601-1	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance (IEC 60601-1:2005)	2006
IEC 60601-1-4	Medical electrical equipment - Part 1: General Requirements for Safety, Part 4: Programmable Electrical Medical Systems	2000
EN 60601-1-2	Medical electrical equipment - Part 1-2: General Requirements for Safety - Collateral Standard: Electromagnetic Compatibility - Requirements and Tests	2001/ A1: 2006
CAN/CSA C22.2 No. 60601-1	Medical Electrical Equipment - Part 1: General Requirements for basic safety and essential performance	2005

Physical Dimensions

Size	277 mm Width x 353 mm Height x 267 mm Depth
Weight	7.3 kg

Operational Environment

Temperature	+10 °C to +33 °C
Humidity	30 % to 75 % RH
Barometric Pressure	700 - 1060 hPa

Transport and Storage Environment

Temperature	-40 °C to +70 °C
Humidity	10 % to 95 % RH
Barometric Pressure	500 to 1060 hPa

Display / Touchscreen

Type	High contrast, digital, graphic color, visible in complete darkness
Resolution	Display 21 cm diagonal, resolution 480 X 640 pixels

Audio Output

Baseline Audio Sound Level 60 dBA minimum SPL (1 m)

Electrical

Input Voltage	100 V-240 V ± 10%
Frequency	50/60 Hz
Power Consumption	500 VA
Auxiliary AC output	200 VA Max.
Internal Fuse	5 x 20 mm T. L. 5 A, 250 V Medtronic Xomed P/N 11270066
Duty Cycle for Applied Part	Maximum On Time 120 Seconds Minimum Off Time 180 Seconds

Power Cord Product Numbers

North America: USA, Barbados, Belize, Bolivia, Canada, Colombia, Ecuador, Venezuela Standard P/N EA600 or 1895820 6 meter P/N EA650 or 189721	United Kingdom, Ireland, Hong Kong, Malaysia, Singapore P/N EA606 or 1895821	Continental Europe: Austria, Belgium, Finland, France, Germany, Greece, Korea, Luxembourg, Netherlands, Norway, Portugal, Spain, Sweden P/N EA602 or 1895822
China P/N EA604	India, South Africa P/N EA607	Switzerland P/N EA601
Argentina P/N EA608	Israel P/N EA609	Denmark P/N EA610
Australia, New Zealand P/N EA605	Japan P/N EA603 or 1895823	Italy, Chile P/N EA611

SYSTEM SOUNDS AND FIGURES

Audible Alarms and Tones

The following alarms and tones can sound while using the IPC Console.

Audible Alarm


When the system detects an error, a message appears on the touchscreen and the system emits a sequence of three tones.

Audible Tones

IPC Tone	Cause(s)
1 Tone	• Confirmation of change button pressed. • Change from Forward to Oscillate. • Change of active handpiece.
2 Tones	Change from Oscillate to Forward.
3 Tones	• Audible Alarm. Error detected. See screen for error message. • Active handpiece is in Reverse and foot pedal pressed. • First time accessory changes from Forward to Reverse.
Long Tone	Change from handpiece to drill.


System Figures

Figure 1-1. IPC Console Front

- 1 Pump 1: Coolant, lens cleaning or irrigation
- 2 Touchscreen
- 3 Power on/off
- 4 Pump 2: Irrigation or lens cleaning
- 5 Console connector panel for peripheral devices

Figure 1-2. IPC Console Back

- 1 Pole clamp
- 2 Compact flash card port (Medtronic Use)
- 3 Manual start/stop
- 4 Auxiliary power outlet
- 5 Endo-Scrub 2 connector
- 6 Fuse access
- 7 Equipotential Ground Connector. Apply potential equalization conductor.
- 8 Hospital grade power cord connector

Figure 1-3. IPC Console Connector Panel


Figure 1-4. Multifunction Foot Pedal & Y-Splitter


Figure 1-5. IPC Touchscreen


Figure 1-6. IPC Pumps Screen


Figure 1-7. Operating Room Setup




Figure 1-9. IntelliFlow Remote Control

- | | |
|---|--|
| 1 | Pause/On-Off |
| 2 | Increase/Decrease Fine Adjustment |
| 3 | Increase/Decrease Coarse Adjustment OR
Select stainless steel tubing size (French size) for
suction irrigator. |

Figure 1-8. IPC System Configuration

- | | | | |
|---|-----------------------------|----|-------------------------------|
| 1 | Irrigation and coolant bags | 7 | Irrigation pole base diameter |
| 2 | Irrigation pole | 8 | Irrigation pole basket |
| 3 | IPC console | 9 | Power cord |
| 4 | Console connector panel | 10 | Pump 2 |
| 5 | Accessory cables | 11 | Pump 1 |
| 6 | Console height | | |

PRE-OPERATING INSTRUCTIONS

The following are general IPC pre-operating instructions. Refer to other sections of this IFU for operating instructions specific to individual handpieces or accessories.

When the System Arrives

- Verify the contents of the box match the packing slip. If incomplete or damaged, notify Medtronic Customer Service.
- If container is damaged, or cushioning material shows stress, notify carrier and Medtronic Customer Service. Keep shipping materials for carrier inspection.
- Save the cartons and packing material. If the instrument is to be shipped the shipping package will provide proper protection.

Set up the IPC

Refer the related topics for detailed instruction.

1. Install pump cartridges or irrigation tubing.
2. Prepare IPC for use.
3. Calibrate touchscreen, if necessary.
4. Change system settings, if necessary.
5. Set up and prime pumps.
6. Confirm system operation.
7. Press the manual start/stop button on the back of the console (Figure 1-2) and verify you can start/stop the handpiece, irrigation and/or coolant flow.

Install the Pump Cartridges or Irrigation Tubing

1. Locate the correct pump and lift up the lock (Figure 1-10).
Pump 1: Coolant, lens cleaning or irrigation
Pump 2: Irrigation
Important: The number on the pump must match the number on the cartridge (either 1/1 or 2/2). If the cartridge does not have a pump designator number, use the Pump Setup Screen to install the pump cartridge.
2. Insert the pump cartridge.
3. Snap the pump lock shut.
Warning: Ensure the pump cartridge does not crimp the tubing.

Figure 1-10. Install Pump Cartridge

Prepare IPC for Use

1. Verify Operation Room set up (Figure 1-7). The surgeon may have preferences to the location and visibility.
2. Verify the wheels are locked on the IPC cart.
3. Inspect all components for damage and determine if the system is ready for use.
4. Mount the IPC and irrigation/coolant bags on the I.V. pole (Figure 1-8).
Important:
 - Mount irrigant and coolant bags above the IPC to ensure adequate flow.
 - It is recommended to use an irrigation pole with minimum base diameter of 53 cm and to mount all items as low as possible to increase stability during use.
 - For transport or uneven floor conditions greater than 10 degrees, maximum height to mount the console is 38 cm if irrigation and coolant bags are at fully extended pole height.
5. Plug the IPC into the power source. Position the IPC so that it does not obstruct the power source for the purpose of disconnecting the Main voltage by the power cord.
6. Locate the correct foot pedal or accessory connection port on the connector panel (Figure 1-3), align the mark on the connector to the mark on the console, and then insert the connector.
7. Connect suction, cooling and/or irrigation tubing.
8. Turn on the IPC and verify the system passes the self-test.
Note: If the IPC does not detect a handpiece or foot pedal the Connect Handpiece/Connect Foot Switch screen appears. Do the following:
 - Verify the cable is connected to the correct connection port.
 - Press [OK] in the Connect Handpiece/Connect Foot Switch message window to continue use of the IPC without the handpiece or foot pedal.

Calibrate Touchscreen

Note: This step is optional.

1. Turn on the IPC console.
2. When the system starts, on the Splash screen, press [Settings].
3. On the Settings screen, press **Touch Screen Calibration** and follow the screen prompts.

Change System Settings

Note: During surgery, system settings can be overwritten.

1. Turn on the IPC console.
2. While the system starts, on the Splash screen, press [Settings].
3. To change the language, press the appropriate language.
4. To change the default settings, press [Default].
 - On the Default screen, press the Forward or Backward arrow to change the accessory.
 - Make changes to the default settings.
 - To confirm system settings and return to the Splash screen, press [OK].
5. For accessories with audible tones, press the REV Audible Tones button to control the following:
 - The system delivers one set of reverse beeps when the Reverse mode is activated.
 - The system delivers one set of reverse beeps the first time the drill is used in Reverse mode after the Reverse mode has been activated.
6. To confirm system settings and continue to the IPC touchscreen, press [OK].
7. To restore settings to factory default, press [Restore].

Handpiece Default Settings

The system configuration is dependent on the handpiece(s) connected to the console. The following table defines the default configurations, default settings (X) and default options (O).

Note: During use, update the default settings for the active handpiece by pressing the set active handpiece settings as default settings button (Figure 1-5).

Table 1. IPC Touchscreen Default Configurations												
Handpiece	Speed			Mode or Mode Select Switch								
	rpm	cpm	%	Forward	Oscillate	Reverse	Acceleration	Deceleration	Size	Flow	Irrigation	Control
Visao	80000			X		O					30	
Indigo	52000			X		O					30	
Midas Rex SC1		3400		O	X						60	
StraightShot M5	5000			O	X						05	
StraightShot M4	5000			O	X						05	
StraightShot III, Magnum II	5000			O	X						30	
Legend EHS Stylus	60000			X		O	45 %	45 %			0	
Legend EHS	70000			X		O					0	
Stylus Touch	60000			X		O	80 %	100 %			0	Finger
Skeeter	16000			X		O	v				0	
Endo-Scrub 2										3		
Suction Irrigator									8	50 %		
Triton			100	X								
			100			X						
Midas Rex Microsaws			100								0	
POWEREASE	120			X								
	120				X							
	250					X						

Set up and Prime Pumps

- The IPC turns on pump 1 and/or 2 long enough to purge air out of the tubing set(s) the first time the prime button is pressed.
- The IPC resets the prime feature when you turn IPC power Off and On.
- After you prime the pump, the prime button and functionality become flush functionality.
 1. Connect tubing from an IPC cartridge to irrigation or coolant tubing on an accessory.
 2. On the irrigation tubing, turn the clamp to OPEN.
 3. If an accessory uses the clear drip chamber (Visao), fill the clear drip chamber with coolant. To fill, squeeze and release the chamber until full.
 4. On the IPC touchscreen (Figure 1-5), press the pumps button:
Note: The IPC pumps screen is also available from the Connect Handpiece/Connect Foot Switch screen which the system displays during IPC preparation for use if a handpiece or foot switch is not detected by the system.
 5. On the IPC pumps screen (Figure 1-6), select the accessory for each pump.
 6. For each pump, press the prime button and verify the following:
 - Pump(s) run until air is completely purged from tubing.
 - Small amount of lubricant flows at the tip of the irrigation device.
 - Pump(s) turns off.
 7. Press the close button.

Pump Default Configurations

The pump configuration is dependent on the handpiece(s) connected to the console. The following table defines the pump default settings (X) and default options (O).

Table 2. IPC Pumps Screen Default Configurations							
Handpiece	Pump 1		Pump 2	Endo-Scrub 2		Suction Irrigator	
	Cooling	Irrigation	Irrigation	Pump 1	Pump 2	Pump 1	Pump 2
Visao	X		X		O		O
Indigo		O	X	O	O	O	O
Midas Rex SC1		O	X	O	O		
StraightShot M5		O	X*	X	O		
StraightShot M4		O	X*	X	O		
StraightShot III, Magum II		O	X*	X	O		
Legend EHS Stylus		X	O*	O	O	O	O
Legend EHS		X	O	O	O	O	O
Stylus Touch		X	O	O	O	O	O
Skeeter				O	O	O	O
Endo-Scrub 2		X	O	X	O		
Suction Irrigator		O	O			O	O
Midas Rex Microsaws		O	X	O	O	O	O
POWEREASE				O	O	O	O
* When the IPC detects both the StraightShot M4 or StraightShot M5 and the Legend EHS Stylus Touch handpiece, by default, the system sets pump 2 as a "shared" irrigation pump. You must manually connect the irrigation tubing to the active handpiece.							

Confirm System Operation

1. Confirm the irrigation pedal starts handpiece and irrigation flow. Verify the speed changes from white to yellow in the Speed box on the touchscreen.
2. Confirm the foot pedal buttons operate. Refer to "Multifunction Foot Pedal" for details.
3. On the touchscreen, verify you can do all of the following:
 - Adjust Speed: In the Speed box, press the plus and minus buttons.
 - Change Modes: In a Mode box, press any mode button.
 - Adjust Flow Rate: In the Irrigation box, press the plus and minus buttons.

IPC COMPONENTS

Auxillary Power to Console

Warning: The auxillary power outlet is available for use with the Hydrodebrider and Bone Mill IPC consoles only (see W82). The auxillary power outlet is for use at grid voltage ≤ 120 VAC only.

Multifunction Foot Pedal

You can use the multifunction foot pedal (Figure 1-4) to start/stop the handpiece, control handpiece speed, handpiece selection and mode of operation. Refer to the Multifunction Foot Pedal Controls topic for each handpiece for specific use and control.

Y-Splitter

Y-Splitter (Figure 1-4) allows using a maximum of two multifunction foot pedals connected to a single IPC (see P23 & P24). In this configuration, the Y-Splitter shall be connected to the IPC, and the Multifunction Foot Pedal(s) shall be connected to the Y-Splitter. When connecting a single foot pedal to the Y-Splitter, you may connect to either Port 1 or 2.

IntelliFlow Irrigation Remote Control

Use the IntelliFlow irrigation remote control (Figure 1-9) to start/stop and change irrigation flow while in the sterile field.

If you are using handpiece irrigation:

- To pause irrigation flow, press the Pause/On-Off button.
- To adjust flow rate, press the Fine Adjustment or Coarse Adjustment Increase/Decrease button.

If you are using the Suction Irrigator:

- To pause or turn on/off the Suction Irrigator, press the Pause/On-Off button.
- To adjust flow rate, press the Fine Adjustment Increase/Decrease button.
- To select the stainless steel tubing size (French size), press the Stainless Steel Tubing Size button.

GUIDANCE AND MANUFACTURER’S DECLARATION • ELECTROMAGNETIC IMMUNITY

Part I

The IPC is intended for use in the electromagnetic environment specified below. The customer or the user of the IPC should assure that it is used in such an environment.			
Immunity test	IEC/EN60601-1-2 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0.5 cycle	<5 % UT (>95 % dip in UT) for 0.5 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the IPC requires continuous operation during power mains interruptions, it is recommended that the IPC be powered from an uninterruptible power supply or a battery.
	40 % UT (60 % dip in UT) for 5 cycles	40 % UT (60 % dip in UT) for 5 cycles	
	70 % UT (30 % dip in UT) for 25 cycles	70 % UT (30 % dip in UT) for 25 cycles	
	<5 % UT (>95 % dip in UT) for 5 sec	<5 % UT (>95 % dip in UT) for 5 sec	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: UT is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic emissions		
The IPC is intended for use in the electromagnetic environment specified below. The customer or the user of the IPC should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The IPC uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The IPC is suitable for use in all establishments, other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings for domestic purpose.
RF emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations IEC 61000-3-3	Complies	

Recommended separation distances between portable and mobile RF communications equipment and the IPC			
The IPC is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the IPC can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the IPC as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum power of transmitter W	Separation distance according to frequency of transmitter meters		
	150 kHz to 80 MHz d = 1.2√P	80 MHz to 800 MHz d = 1.2√P	800 MHz to 2.5 GHz d = 2.3√P
0.01	0.12	0.12	0.23
0.10	0.38	0.38	0.73
1.00	1.20	1.20	2.30
10.00	3.80	3.80	7.30
100.00	12.00	12.00	23.0
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Part II

The IPC is intended for use in the electromagnetic environment specified below. The customer or the user of the IPC should assure that it is used in such an environment.			
Immunity test	IEC/EN60601-1-2 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the IPC , including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance
Conducted RF	3 Vrms	3 Vrms	$d = 1.2 \sqrt{P}$
IEC 61000-4-6	150 kHz to 80 MHz		
Radiated RF	3 V / m	3 V / m	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz
IEC 61000-4-3	80 MHz to 2.5 GHz		$d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz
			Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. NOTE 3 When operating the IPC with Stylus Touch, the compliance level is 3 V/m except from 88 MHz to 91 MHz where it is 1 V/m. The formula for separation distance for the IPC with Stylus Touch will be $d = 3.5 \sqrt{P}$ in that frequency range.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the IPC is used exceeds the applicable RF compliance level above, the IPC should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the IPC. ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

LIMITED WARRANTY

- A. This Limited Warranty provides the following assurance for the customer who purchases a Medtronic IPC System. This Limited Warranty is extended only to the buyer purchasing the IPC System directly from Medtronic or from its affiliate or its authorized distributor or representative. The IPC System includes the console, motor or handpiece, foot control, motor cables, instrumentation cases and trays (hereafter referred to as System Components), straight and angled motor attachments (hereinafter referred to as "Attachments"), bur guards and telescoping tubes (hereinafter referred to as Semi-reusable Components) and dissecting tools, irrigation and coolant tubing, and Intelliflow remote control (hereinafter referred to as Single Use Components) and jointly referred to as the IPC System, unless specifically noted.
- i. Should a System Component fail to function to Medtronic's published specifications during the term of this Limited Warranty (one [1] year from the date of sale of a new System Component or ninety [90] days from the date of sale of a refurbished or used System Component), Medtronic will either repair or replace the Motor Component or any portion thereof.
 - ii. Should an Attachment fail to function to Medtronic's published specifications during the term of this Limited Warranty (ninety [90] days from the date of sale of a new Attachment), Medtronic will either repair or replace the Attachment or any portion thereof.
 - iii. Should a Semi-reusable Component fail to function to Medtronic's published specifications during the term of this Limited Warranty (thirty [30] days from the date of sale of a new Semi-reusable Component), Medtronic will replace the Semi-reusable Component or any portion thereof.
 - iv. Should a Single Use Component fail to function to Medtronic's published specifications prior to its "use by" date Medtronic will replace the Single Use Component.
- B. To qualify for this Limited Warranty, the following conditions must be met:
- i. The Product must be used on or before its "Use By" or "Use Before" date, if applicable.
 - ii. The Product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.
 - iii. Medtronic must be notified in writing within thirty (30) days following discovery of a defect.
 - iv. The Product must be returned to Medtronic within thirty (30) days of Medtronic receiving notice as provided for in (3) above.
 - v. Upon examination of the Product by Medtronic, Medtronic shall have determined that: (i) the Product was not repaired or altered by anyone other than Medtronic or its authorized representative, (ii) the Product was not operated under conditions other than normal use, and (iii) the prescribed periodic maintenance and services, if applicable, have been performed on the Product.
- C. This Limited Warranty is limited to its express terms. THIS LIMITED WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED WHETHER STATUTORY OR OTHERWISE, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. In no event shall Medtronic be liable for any consequential, incidental, prospective or other similar damage resulting from a defect, failure, or malfunction of the IPC System, whether a claim for such damage is based upon the warranty, contract, negligence or otherwise.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. Users may benefit from statutory warranty rights under legislation governing the sale of consumer goods. If any part or term of this Limited Warranty is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the Limited Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Limited Warranty did not contain the particular part or term held to be invalid.

FOR ITEMS CONTAMINATED WITH TSE AGENTS

Medtronic ENT/NT Transmissible Spongiform Encephalopathy (TSE) Return Policy

Medtronic will not authorize or accept the return of products that directly contact patients or is contaminated with a patient's body fluids suspected or confirmed with a Transmissible Spongiform Encephalopathy / Creutzfeldt-Jakob Disease (TSE/CJD) diagnosis.

The following are recommended guidelines and may vary according to specific policy and procedures among hospitals. Hospital personnel should contact their infection control personnel for current procedures and policy for reusable equipment processing when suspected of contamination with Creutzfeldt-Jakob Disease (CJD) or other Transmissible Spongiform Encephalopathy (TSE) agent.

Medtronic dissecting tools, burs, or blades used on a patient suspected of a TSE/CJD diagnosis should be incinerated. Reusable equipment that has been used on patients with suspected Creutzfeldt-Jakob Disease (CJD) or other Transmissible Spongiform Encephalopathy (TSE) should be quarantined and not reused until diagnosis is confirmed or excluded. Reusable equipment should be quarantined after having been cleaned, decontaminated, sterilized and packed in a rigid sealed container until final diagnosis. If TSE/CJD is excluded as a diagnosis, the quarantined reusable equipment may be returned for use after appropriate cleaning, decontamination and sterilization.

Medtronic recommends that all Medtronic products used directly on a patient confirmed with a TSE diagnosis be incinerated. Contact your Sales Representative to purchase replacement products or secure loaner equipment.

For additional information contact your Customer Service Representative.

TRITON ELECTRIC HIGH-TORQUE HANDPIECE

DEVICE DESCRIPTION

The Triton Electric High-Torque Handpiece is capable of removing hard and soft tissue, drilling pilot holes, and driving screws, wires, and pins during spinal, cranial, and small-bone surgical procedures performed in an operating-room environment by surgeons trained in its use.

The following instructions for the Triton Electric High-Torque Handpiece are in addition to “Set up the IPC” general assembly instructions. Complete IPC setup, then continue to the instructions below.

BEFORE USE

Clean and sterilize the device, prior to first use. Refer to the Reprocessing Instructions contained in this manual for additional information.

TRITON SAGITTAL SAW ASSEMBLY

Warning: Triton saw attachments should only be used with Medtronic Triton saw blades. Refer to the Triton Quick Reference Saw Blade Guide (LIT200017) for additional information on Triton saw blades.

1. Insert the Sagittal Saw Attachment into the handpiece in a position to allow easy insertion of the Sagittal Saw Key (Figure 2-1).
Note: You can install the attachment in 12 different positions to facilitate proper surgical access.
2. Insert the Sagittal Saw Key into the attachment and turn counterclockwise until there is slight resistance.
3. Insert the blade into the space between the two jaws, ensuring that the blade is fully seated.
4. Turn the Sagittal Saw Key clockwise to lock the blade. Run briefly, then retighten blade.
Caution: Do not over-tighten.

Figure 2-1. Sagittal Saw Assembly

TRITON SAGITTAL SAW DISASSEMBLY

To remove the Triton saw blade, insert the Sagittal Saw Key into the Attachment and turn counterclockwise.

TRITON RECIPROCATING SAW ASSEMBLY

Warning: Triton saw attachments should only be used with Medtronic Triton saw blades. Refer to the Triton Quick Reference Saw Blade Guide (LIT200017) for additional information on Triton saw blades.

1. Loosen the collet nut then insert the blade until it is fully seated (Figure 2-2).
Note: You can install the attachment in different positions to facilitate proper surgical access.
2. Finger-tighten the collet nut. Run briefly, then retighten collet nut.

Figure 2-2. Reciprocating Saw Assembly

TRITON RECIPROCATING SAW DISASSEMBLY

To remove the Triton Reciprocating Saw blade, unscrew the collet nut.

TRITON AO/SYNTHES CHUCK AND TRINKLE CHUCK ASSEMBLY

1. To install a drill bit, pull back on the attachment collar (Figure 2-3).
2. Insert the drill bit and release the attachment collar.

Figure 2-3. AO/Synthes Chuck and Trinkle Chuck Assembly

TRITON AO/SYNTHES CHUCK AND TRINKLE CHUCK DISASSEMBLY

1. To remove a drill bit, pull back on the attachment collar.
2. Remove the drill bit and release the attachment collar.

TRITON JACOBS CHUCK ASSEMBLY

- 1. To install a drill bit, turn the key to open the chuck or spin the collar if using a keyless chuck attachment (Figure 2-4).
- 2. Insert the drill bit.

Figure 2-4. Jacobs Chuck Assembly

TRITON JACOBS CHUCK DISASSEMBLY

- 1. To remove a drill bit, turn the key to open the chuck or spin the collar if using a keyless chuck attachment.
- 2. Remove the drill bit.

TRITON HUDSON AND ZIMMER CHUCK ASSEMBLY

To install an instrument, pull back on the attachment collar, then insert the male end of the instrument into the chuck (Figure 2-5).

Figure 2-5. Hudson and Zimmer Chuck Assembly

TRITON HUDSON AND ZIMMER CHUCK DISASSEMBLY

To remove an instrument, pull back on the attachment collar, then remove the instrument from the chuck.

TRITON WIRE AND PIN COLLET ASSEMBLY


The Wire Collet accepts wires up to 1.6 mm (.062 inches) in diameter. The Pin Collet accepts pins up to 3.2 mm (.125 inches) in diameter.

- 1. Insert the Wire or Pin Collet (Figure 2-6) while the handpiece is in the SAFE position (Figure 2-7).
- 2. Screw the Cannulated Extension in the back of the handpiece to protect the operator from the point of the wire or pin, as necessary.
- 3. Insert the wire or pin into the front or back of the handpiece.
- 4. Put the instrument in the RUN position by positioning the trigger control vertically.
- 5. Turn the Mode select switch at the base of the handle to the FORWARD position.
- 6. Squeeze the Wire/Pin Advance lever and hold it down.
- 7. Press the trigger control to drive the wire/pin. The pressure-sensitive trigger allows variable speed operation.
- 8. To obtain additional wire/pin length, release the wire/pin.
- 9. Advance the lever and trigger control.
- 10. Pull back on the instrument.
- 11. Squeeze the Wire/Pin Advance lever and trigger control to drive the wire.

TRITON WIRE AND PIN COLLET DISASSEMBLY

To remove threaded wire/pin, put the Mode select switch in REVERSE, squeeze the Wire/Pin Advance lever and press the trigger control.

Figure 2-6. Wire and Pin Collet Assembly

TRITON ELECTRIC HIGH-TORQUE HANDPIECE OPERATION

You can preload attachments before insertion into the handpiece.

Caution: Insert all attachments into the handpiece with the handpiece in the SAFE position (Figure 2-7).

Note: The handpiece has an Extension Handle that screws into the back of the handpiece. The handle extension provides balance and two-handed control for various drilling and cutting applications.

- 1. With the handpiece in the SAFE position (Figure 2-7), insert a preloaded attachment by pressing the quick-release button on top of the handpiece. Snap the attachment into the handpiece with a slight twisting motion until it is seated.
- 2. Place the handpiece in the RUN position, with the trigger control vertical. The Mode select switch at the base of the handle should be in the FORWARD position.
- 3. After use, return the trigger control to the SAFE position prior to removing the attachment.
- 4. Remove the attachment by pressing the quick-release button on top of the handpiece.

Safe Position

The handpiece will not operate in the SAFE position. To operate the handpiece, activate and press the trigger control.

Figure 2-7. Triton Handpiece Safe and Run positions


CONNECT TRITON ELECTRIC HIGH-TORQUE HANDPIECE TO IPC

Locate the Triton Electric High-Torque Handpiece connection port on the connector panel (Figure 2-8) and insert the connector.

Note: To insert multi-pin connectors (indicated by a silver or red mark on the connector), align the mark on the connector to the mark on the console, then insert the connector.


Figure 2-8. Triton IPC Connection Ports

TRITON ELECTRIC HIGH-TORQUE HANDPIECE TOUCHSCREEN CONTROLS

To adjust Triton Electric High-Torque Handpiece variable speed, on the IPC touchscreen, in the FWD Speed or REV Speed control box (Figure 2-9), press the plus button to increase variable speed or the minus button to decrease variable speed.

Figure 2-9. Triton Touchscreen

TRITON ELECTRIC HIGH-TORQUE HANDPIECE MODE SELECT SWITCH

Use the mode select switch to change the handpiece from Forward to Reverse when the handpiece is the active handpiece. When the handpiece is the inactive handpiece, use the mode select switch to activate the handpiece.

TRITON ELECTRIC HIGH-TORQUE HANDPIECE MULTIFUNCTION FOOT PEDAL CONTROLS

Important: By default, press each button on the foot pedal for at least 100 mS for the selection to become active. Use the IPC touchscreen Settings screen to change the default value.

To use the multifunction foot pedal (Figure 2-10) to control the handpiece to do the following:

- To toggle between the start/stop mode and variable speed mode, press the control button. 
- To change the handpiece, press the handpiece button. 

Triton Electric High-Torque Handpiece Reverse Foot Pedal Control

When you connect the optional foot pedal to the IPC console, the pedal can be used as an alternative method of activating Reverse mode (without manipulating the Mode select switch at the bottom of the handpiece). When the pedal is connected, the Reverse Pedal control box appears on the screen (Figure 2-11).

Note: By default, pedal functionality is turned OFF.

1. To use the pedal, press ON in the Reverse Pedal control box on the IPC touch screen.
2. With the pedal turned on, step on the foot pedal to put the handpiece into Reverse mode (regardless of what position the Mode select switch is in).
NOTE: The trigger control must be pressed to activate the handpiece, even when you are stepping on the foot pedal.
3. Remove your foot from the pedal to return the handpiece to the mode currently defined by the Mode select switch on the handpiece.

Figure 2-10. Multifunction Foot Pedal & Y-Splitter

- | | |
|---------------------------|--------------|
| 1 Mode button | 5 Foot pedal |
| 2 Handpiece button | 6 Y-Splitter |
| 3 Control button | 7 Port 1 |
| 4 Slip-resistant foot pad | 8 Port 2 |

Figure 2-11. Reverse Pedal Control Box

TRITON ELECTRIC HIGH-TORQUE HANDPIECE CLEANING AND STERILIZATION INSTRUCTIONS

Refer to document M000030A322 in the Cleaning and Sterilization section.

TRITON ELECTRIC HIGH-TORQUE HANDPIECE TECHNICAL SPECIFICATIONS

Triton Electric High-Torque Handpiece ED500

Size	3.5 inches Length x 5.4 inches Height x 1.1 inches Width
Weight	2.1 lbs
Speed	180-1800 rpm (actual speed depends on attachment used)
Duty Cycle for Applied Part	Cycle Time: 20 seconds on maximum / 20 seconds off minimum Maximum number of cycles before resting handpiece: 6 Maximum number of cycles before resting attachment: 3 Minimum rest period: 25 minutes

SUCTION IRRIGATOR

The following instructions for the Suction Irrigator are in addition to "Set up the IPC" general assembly instructions. Complete IPC setup, then continue to the instructions below.


SUCTION IRRIGATOR ASSEMBLY

1. Connect the suction tubing from a suction source to the suction fitting on the Suction Irrigator (Figure 3-1).
2. On the IPC touchscreen (Figure 1-5), press the pumps button. 
3. On the IPC pumps screen (Figure 1-6), select a pump for the Suction Irrigator.
Note: When using a handpiece that connects to a pump, the IPC automatically incorporates the Suction Irrigator at the pump not in use by the handpiece.
4. Connect the irrigation tubing from the IPC cartridge (Figure 1-8) to the irrigation fitting on the Suction Irrigator (Figure 3-1).
5. On the irrigation tubing, turn the clamp to OPEN.

SUCTION IRRIGATOR ADAPTER KIT

1. Connect an adapter to the high-speed irrigation tubing (blue adapter) or the IPC tubing (white adapter).
2. Connect an adapter to the irrigation connector tube (Figure 3-2).
3. Connect an irrigation connector tube to the irrigation fitting on the Suction Irrigator.

Figure 3-1. Suction Irrigator

- | | |
|-------------------|----------------------|
| 1 Suction tube | 4 Irrigation fitting |
| 2 Irrigation tube | 5 Suction fitting |
| 3 Tube size | |

Figure 3-2. Suction Irrigator Adapter Kit

SUCTION IRRIGATOR IN STANDALONE MODE

The Suction Irrigator tool can operate as a standalone device by changing the IPC system defaults.

- 1. On the Defaults menu (Figure 3-3), in the Stand Alone section, select **Suction Irrigator** and then press [OK].
- 2. On the Handpiece Connection screen (Figure 3-4), in the Suction Irrigator box, press [OK].

Figure 3-3. IPC General System Default Menu


Figure 3-4. Handpiece Connection Screen

SUCTION IRRIGATOR TOUCHSCREEN CONTROLS

To set or adjust Suction Irrigator controls, on the IPC touchscreen, in the Suction Irrigator control box (Figure 3-5), do the following:

- To set the tubing size, in the Size control box, press the plus and minus buttons.
- **Note:** The system defaults to size 8.
- To enable or disable the irrigation flow, in the Flow control box, select the On/Off box.
- To adjust the flow rate, in the Flow control box, press the plus and minus buttons.

Figure 3-5. Suction Irrigator Touchscreen

ENDO-SCRUB 2


The following instructions for the Endo-Scrub 2 are in addition to “Set up the IPC” general assembly instructions. Complete IPC setup, then continue to the instructions below. Refer to the Endo-Scrub 2 System Instructions for Use, Endo-Scrub Sheaths Instructions for Use and Endo-Scrub 2 Finger Switch Instructions for Use for additional information.

The IPC System incorporates Endo-Scrub 2 functionality by using irrigation pump number one (1) and controlling operation with the touch screen and an external foot pedal or finger switch.

DO NOT use the Endo-Scrub 2 for infusion, for disinfection or sterilization of an endoscope, or for suction removal of blood and debris.


Use the Endo-Scrub 2 sheath only with an endoscope listed on the sheath product label, as malfunction or poor performance could result.

Figure 4-1. Endo-Scrub 2 Finger Switch

- | | | | |
|---|-----------------------|---|---------------------|
| 1 | Finger switch | 4 | Light connection |
| 2 | Endo-Scrub 2 sheath | 5 | Finger Switch cable |
| 3 | Irrigation connection | | |

Figure 4-2. Endo-Scrub 2 Foot Pedal

ENDO-SCRUB 2 ASSEMBLY

1. Wet the endoscope.
2. Slowly, slide the approved endoscope into the Endo-Scrub 2 sheath (Figure 4-3).
3. Connect the irrigation tubing and a light source (Figure 4-4).

Figure 4-3. Endo-Scrub 2 Assembly

Figure 4-4. Endo-Scrub 2 Assembly

ENDO-SCRUB 2 FINGER SWITCH ASSEMBLY

If using the Endo-Scrub 2 finger switch, complete the following:

1. Slide the finger switch onto the Endo-Scrub 2 sheath (Figure 4-1). Align the cutout section of the ring with the luer connector of the tubing set. The finger switch is properly installed when the cutout section of the ring is firmly seated against the luer connector.
2. Activate the pump by pressing the actuator button located on the finger switch.

ENDO-SCRUB 2 ACTIVATION

Note: The procedure below also applies if using the multifunction foot pedal.

1. To activate the Endo-Scrub wash cycle, press and release the finger switch.
2. To initiate a continuous flow of irrigant, press and hold the finger switch.

CONNECT ENDO-SCRUB 2 TO IPC CONSOLE

1. Locate the Endo-Scrub 2 connector cover on the back of the IPC console (Figure 1-2).
2. Insert a small screwdriver in the notch on the cable connector cover and pull.
3. Connect the control switch cable to the cable connector.
4. Connect the Endo-Scrub 2 finger switch (Figure 4-1) or the foot pedal (Figure 4-2) to the console (Figure 4-5).

Figure 4-5. IPC Endo-Scrub 2 Connection Port

- 1 Finger switch or foot pedal connection port

ENDO-SCRUB 2 TOUCHSCREEN CONTROLS

To set or adjust Endo-Scrub 2 controls, on the IPC touchscreen, in the Flow section of the Endo-Scrub control box (Figure 4-6), do the following:

- To enable the Endo-Scrub 2 , press the On/Off check-box.
- To adjust the flow rate, press the plus button to increase flow rate or the minus button to decrease flow rate.  
- To prime the pump, press the prime button. 

Figure 4-6. Endo-Scrub 2 Touchscreen

ENDO-SCRUB 2 IN STANDALONE MODE

The Endo-Scrub 2 can operate as a standalone device by changing the IPC system defaults.

1. On the Defaults menu (Figure 4-7), in the Stand Alone section, select **Endo-Scrub** and then press **[OK]**.
2. On the Handpiece Connection screen (Figure 4-8), in the Endo-Scrub box, press **[OK]**.

Figure 4-7. IPC General System Default Menu


Figure 4-8. Handpiece Connection Screen

ENDO-SCRUB 2 CLEANING AND STERILIZATION INSTRUCTIONS

Refer to document 68E4005 in the Cleaning and Sterilization section.

Endo-Scrub 2 Foot Pedal Cleaning

Important: If debris is found under the foot pedal boot, return for warranty service.

DO NOT immerse or sterilize the foot pedal unit.

DO NOT use alcohol, other solvents or abrasive cleaners.

1. Wipe down the Endo-Scrub 2 foot pedal with a cloth dampened with a neutral enzymatic detergent, pH 6.0-8.0 or phenol based disinfectant.
2. Dry the unit with a clean, non-abrasive cloth.


SPINE SHAVER (SC1) HANDPIECE

DEVICE DESCRIPTION

The following instructions for the Spine Shaver (SC1) Handpiece are in addition to “Set up the IPC” general assembly instructions. Complete IPC setup, then continue to the instructions below. Refer to the Midas Rex User’s Guide for additional information.

The IPC incorporates the Midas Rex Spine Shaver (SC1) at pump 2. Control operation of the Midas Rex Spine Shaver (SC1) with the IPC touchscreen and the multifunction foot pedal.

Figure 5-1. Spine Shaver (SC1) Handpiece

SPINE SHAVER (SC1) BLADE OR BUR ASSEMBLY

1. Insert the tool aligning the tabs with the notches (Figure 5-2). Orientate the irrigation barb to the left or right side.
Note: The StraightShot M4 uses a four-tab alignment system.
 - For rotating straight blades, orient the irrigation barb at the 3 o’clock position for right-handed surgeons and 9 o’clock for left-handed surgeons.
 - For rotating curved blades, orient the irrigation barb at 3 o’clock.
2. Press the locking collar (Figure 5-3).
3. Release the locking collar.
Note: If collar does not return to full out position adjust the finger wheel with small back-and-forth motions until collar pops out.
4. Pull on the blade or bur to ensure engagement and visually check to make sure the distal tip of the inner blade is in contact with the distal tip of the outer cannula (Figure 5-4).

Figure 5-2. Spine Shaver (SC1) Assembly

Figure 5-3. Spine Shaver (SC1) Assembly



Figure 5-4. Spine Shaver (SC1) Assembly

SPINE SHAVER (SC1) HANDPIECE

SPINE SHAVER (SC1) SUCTION AND IRRIGATION TUBE ASSEMBLY

- 1. Attach a suction tube to the suction source (Figure 5-5) and an irrigation tube on the irrigation barb (Figure 5-5).
- 2. Secure suction and irrigation in the irrigation groove on the handpiece (Figure 5-6).

Figure 5-5. Spine Shaver (SC1) Suction Irrigation Assembly



Figure 5-6. Spine Shaver (SC1) Suction Irrigation Assembly


CONNECT SPINE SHAVER (SC1) TO IPC CONSOLE

Locate the Spine Shaver (SC1) connection port on the connector panel (Figure 5-7) and insert the connector.

Note: To insert multi-pin connectors (indicated by a silver or red mark on the connector), align the mark on the connector to the mark on the console, then insert the connector.

Figure 5-7. Spine Shaver (SC1) IPC Connection Ports

- 1 Spine Shaver (SC1) handpiece connection port
- 2 Multifunction foot pedal connection port

SPINE SHAVER (SC1) TOUCHSCREEN CONTROLS

To set or adjust Spine Shaver (SC1) controls, on the IPC touchscreen, in the control box (Figure 5-8), do the following:

- To change rotation mode, in the Mode control box select OSC (oscillating) or FWD (forward).
- To adjust speed, in the Speed control box, press the plus button to increase speed and the minus button to decrease speed. 

Forward Mode: Default, 12000 rpm; variable adjustment from 50 to 12000 rpm.
Oscillate Mode: Default, 3400 cpm; variable adjustment from 50 to 5000 rpm.

- To adjust the irrigation flow rate, in the Pump control box, press the plus button to increase flow rate or the minus button to decrease flow rate. If intermittent flow is available, pressing the plus or minus button progresses the system through intermittent and continuous flow. The system displays **Intermittent** when in intermittent flow mode.
Forward Mode: Default, 30 cc per minute.
Oscillate Mode: Default, 60 cc per minute.

Note: To adjust flow rate, you can use the touchscreen or the IntelliFlow Irrigation remote control.

- In oscillating mode only, you can use the Blade Position control box to do any of the following:


To rotate the inner blade 180° press the delta button. 

To rotate the inner blade in small increments, press the clockwise or counter-clockwise buttons. 

Note: The motion indicator indicates rotation direction of the blade. 

- To rotate the outer blade, use the finger wheel (Figure 5-1).

Figure 5-8. Spine Shaver (SC1) Touchscreen


SPINE SHAVER (SC1) BLADE CONTROL

Note: Periodically submerge blade tip in sterile water, with suction on, to keep blades clear during the procedure.

- To rotate the outer blade (Figure 5-9), use the finger wheel (Figure 5-1).
- To rotate the inner blade, use the Blade Position control box on the IPC touchscreen.

Figure 5-9. Blade Dissection

- 1 Suction flow in through inner blade
Irrigation flow between inner and outer blades
- 2 Outer blade
- 3 Inner blade
- 4 Outer sleeve


SPINE SHAVER (SC1) MULTIFUNCTION FOOT PEDAL CONTROLS

Important: By default, press each button on the foot pedal for at least 100 mS for the selection to become active. Use the IPC touch screen Settings screen to change the default value.

Use the multifunction foot pedal (Figure 5-10) to control the following:

- To select forward or oscillating mode, press the mode button.
- To start or adjust the speed of a handpiece in variable mode, press the foot pedal.
- To rotate the inner blade (180°) press the control button.
- To change the handpiece, press the handpiece button.

Figure 5-10. Multifunction Foot Pedal & Y-Splitter

- 1 Mode button
- 2 Handpiece button
- 3 Control button
- 4 Slip-resistant foot pad
- 6 Foot pedal
- 7 Y-Splitter
- 8 Port 1
- 9 Port 2

SPINE SHAVER (SC1) CLEANING AND STERILIZATION INSTRUCTIONS

Refer to document 68E3282 in the Cleaning and Sterilization section.

SPINE SHAVER (SC1) TECHNICAL SPECIFICATIONS

Spine Shaver (SC1) ED100

Size	14.3 cm Length x 1.8 cm Width (ED100)
Weight	228 g ED100
Speed	50-5000 rpm oscillate 50-12000 rpm forward
Duty Cycle for Applied Part	Maximum on time = 60 seconds. Minimum off time = 30 seconds.

STRAIGHTSHOT M5, STRAIGHTSHOT M4, STRAIGHTSHOT MAGNUM II, AND STRAIGHTSHOT III

DEVICE DESCRIPTION

Note: The StraightShot III is not for sale in the United States.



The following instructions for the StraightShot M5, StraightShot M4, StraightShot Magnum II, and the StraightShot III are in addition to “Set up the IPC” general assembly instructions. Complete IPC setup, then continue to the instructions below.

The IPC incorporates the microdebrider at pump 2. Control operation of the microdebrider with the IPC touchscreen and the multifunction foot pedal.

When the IPC detects both the StraightShot M4 or StraightShot M5 and the Legend EHS Stylus handpieces, the system defaults Pump 2 to the **Shared** configuration. You must manually move the irrigation tubing from the inactive to the active handpiece. Use the pumps screen to override the **Shared** default by selecting the StraightShot M4 or the Legend EHS Stylus for Pump 1. Refer to “Set up and Prime Pumps” for more information.

Figure 6-1. StraightShot M5 Handpiece

Figure 6-2. StraightShot Magnum II and III Handpiece

- | | |
|----------------------------|---------------------|
| 1 Finger wheel | 4 Finger wheel lock |
| 2 Irrigation tubing groove | 5 Locking collar |
| 3 Suction barb | |

- | | |
|----------------|------------------|
| 1 Suction barb | 2 Locking collar |
|----------------|------------------|




STRAIGHTSHOT M4 AND STRAIGHTSHOT M5 BLADE OR BUR ASSEMBLY

1. Insert the tool aligning the tabs with the notches (Figure 6-3). Orientate the irrigation barb to the left or right side.
Note: The StraightShot M4 and the StraightShot M5 use a four-tab alignment system.
 - For rotating straight blades, orient the irrigation barb at the 3 o'clock position for right-handed surgeons and 9 o'clock for left-handed surgeons.
 - For rotating curved blades, orient the irrigation barb at 3 o'clock.
 - For rotatable blades, adjust the finger wheel with small back-and-forth motions until the blade is engaged.
2. Press the locking collar (Figure 6-4).
3. Release the locking collar.
Note: If the collar does not return to full out position adjust the finger wheel with small back-and-forth motions until collar pops out.
4. Pull on the blade or the bur to ensure engagement and visually check to make sure the distal tip of the inner blade is in contact with the distal tip of the outer cannula (Figure 6-5).

Figure 6-3. Blade or Bur Assembly

Figure 6-4. Blade or Bur Assembly

Figure 6-5. Blade or Bur Assembly

STRAIGHTSHOT M5, STRAIGHTSHOT M4, STRAIGHTSHOT MAGNUM II, and STRAIGHTSHOT III

STRAIGHTSHOT M4 AND STRAIGHTSHOT M5 SUCTION AND IRRIGATION TUBE ASSEMBLY

- 1. Attach a suction tube to the suction source and an irrigation tube on the irrigation barb (Figure 6-6).
- 2. Secure suction and irrigation in the irrigation groove on the handpiece (Figure 6-7).

Figure 6-6. Suction and Irrigation Assembly




Figure 6-7. Suction and Irrigation Assembly

STRAIGHTSHOT MAGNUM II AND STRAIGHTSHOT III BLADE AND BUR ASSEMBLY

- 1. Press the collet and insert blade in collet (Figure 6-8).
- 2. Release the collet (Figure 6-8).
- 3. Pull on the tool to ensure engagement and check distal tip of inner blade is in contact with the distal tip of the outer cannula.

Figure 6-8. StraightShot II and III Blade and Bur Assembly

STRAIGHTSHOT MAGNUM II AND STRAIGHTSHOT III SUCTION AND IRRIGATION ASSEMBLY

- 1. Attach a suction tube to the suction source and an irrigation tube on the irrigation barb (Figure 6-9).
- 2. Secure suction and irrigation tubing with tubing clips.


Figure 6-9. StraightShot II and III Suction and Irrigation Assembly

CONNECT STRAIGHTSHOT TO IPC CONSOLE

Locate the StraightShot connection port on the connector panel (Figure 6-10) and insert the connector.

Note: To insert multi-pin connectors (indicated by a silver or red mark on the connector), align the mark on the connector to the mark on the console, then insert the connector.

- 1 Handpiece connection port
- 2 Multifunction foot pedal connection port

STRAIGHTSHOT M4, STRAIGHTSHOT MAGNUM II AND STRAIGHTSHOT III TOUCHSCREEN CONTROLS

Note: The StraightShot M4, StraightShot Magnum II and StraightShot III handpiece screens feature the same controls as those shown on the StraightShot M4 Control Box.


Note: When you stop the blade, one of the following occurs:

- If the IPC button is visible on the touchscreen, the inner blade returns to the same position it began.
- If the XPS button is visible on the touchscreen, the inner blade stops the current position.

To set or adjust StraightShot controls, on the IPC touchscreen, in the StraightShot Handpiece control box (Figure 6-11), do the following:

- To change rotation mode, select OSC (oscillating) or FWD (forward).
Note: The system displays the default Oscillating or Forward mode speed.
- To adjust speed, in the Speed control box, press the plus to increase speed or the minus button to decrease speed.
Forward Mode: Default, 12000 rpm; variable adjustment from 50 to 12000 rpm.
Oscillate Mode: Default, 5000 rpm; variable adjustment from 50 to 5000 rpm.
- To adjust the irrigation flow rate, in the Pump control box, press the plus button to increase flow rate or the minus button to decrease flow rate. If intermittent flow is available, pressing the plus or minus button progresses the system through intermittent and continuous flow. The system displays **Intermittent** when in intermittent flow mode.
Forward Mode: Default, 30 cc per minute.
Oscillate Mode: Default, 5 cc per minute.
Note: To adjust flow rate, you can use the touchscreen or the IntelliFlow Irrigation remote control.
- To rotate outer blade, use the finger wheel (Figure 6-1).
- **In oscillating mode only**, you can use the Blade Position control box to do any of the following:
Note: The motion indicator indicates rotation direction of the blade.
To enable the multifunction foot pedal to change rotation displacement, press the delta button.
To rotate inner blade in small increments, press the counter-clockwise buttons.

Figure 6-11. StraightShot M4 Touchscreen

STRAIGHTSHOT M5 TOUCHSCREEN CONTROLS

Note: When you stop the blade, one of the following occurs:

- If the IPC button is visible on the touchscreen, the inner blade returns to the same position it began.
- If the XPS button is visible on the touchscreen, the inner blade stops the current position.

To set or adjust StraightShot controls, on the IPC touchscreen, in the StraightShot Handpiece control box (Figure 6-11), do the following:

- To change rotation mode, select OSC (oscillating) or FWD (forward).
Note: The system displays the default oscillating or forward mode speed.
- To adjust speed, in the Speed control box, press the plus to increase speed or the minus button to decrease speed.
Forward Mode: Default, 12000 rpm; variable adjustment from 50 to 30000 rpm.
Oscillate Mode: Default, 5000 rpm; variable adjustment from 50 to 7500 rpm.
- To adjust the irrigation flow rate, in the Pump control box, press the plus button to increase flow rate or the minus button to decrease flow rate. If intermittent flow is available, pressing the plus or minus button progresses the system through intermittent and continuous flow. The system displays **Intermittent** when in intermittent flow mode.
Forward Mode: Default, 30 cc per minute.
Oscillate Mode: Default, 5 cc per minute.
Note: To adjust flow rate, you can use the touchscreen or the IntelliFlow Irrigation remote control.
- To rotate outer blade, use the finger wheel (Figure 6-1).

In oscillating mode only, you can use the Blade Position control box to do any of the following:

Note: The motion indicator indicates rotation direction of the blade.

To enable the multifunction foot pedal to change rotation displacement, press the delta button.

To rotate inner blade in small increments, press the counter-clockwise buttons.

Figure 6-12. StraightShot M5 Touchscreen



Figure 6-13. Blade Dissection

- 1 Suction flow in through inner blade irrigation flow between inner and outer blades
- 2 Outer blade
- 3 Inner blade
- 4 Outer sleeve

MICRODEBRIDER BLADE CONTROL

Important: If airway blade becomes clogged during use, 1-5 cc of irrigant could be aspirated by the patient before you detect the clog.

Note: Periodically submerge blade tip in sterile water, with suction on, to keep blades clear during the procedure.

- To rotate the outer blade (Figure 6-13), use the finger wheel (Figure 6-1).
- To rotate the inner blade, use the Blade Position control box on the IPC touchscreen. Refer to the related accessory Controls topic for further information.


MULTIFUNCTION FOOT PEDAL CONTROLS

Important: By default, press each button on the foot pedal for at least 100 mS for the selection to become active. Use the IPC touch screen Settings screen to change the default value.

To use the multifunction foot pedal (Figure 6-14) to control the handpiece do the following:

- To select forward or oscillate mode, press the mode button.
- To start or adjust the speed of a handpiece in variable mode, press the foot pedal.
- To rotate the inner blade (60° or 180°), press the control button.
- To change the handpiece, press the handpiece button.

Figure 6-14. Multifunction Foot Pedal & Y-Splitter

- | | | | |
|---|-------------------------|---|------------|
| 1 | Mode button | 5 | Foot pedal |
| 2 | Handpiece button | 6 | Y-Splitter |
| 3 | Control button | 7 | Port 1 |
| 4 | Slip-resistant foot pad | 8 | Port 2 |

CLEANING AND STERILIZATION INSTRUCTIONS

Refer to document 68E3282 in the Cleaning and Sterilization section.

STRAIGHTSHOT M4 TECHNICAL SPECIFICATIONS

StraightShot M4 1898200T

Size	14.3 cm Length x 1.8 cm Width (1898200T)
Weight	228 g 1898200T
Speed	50-5000 rpm oscillate 50-12000 rpm forward
Duty Cycle for Applied Part	The StraightShot M4 handpiece under full load is rated for intermittent operation per the following: Maximum On Time: 60 seconds Minimum Off Time: 30 seconds

STRAIGHTSHOT M5 TECHNICAL SPECIFICATIONS

StraightShot M5 1899200

Size	12.5 cm Length x 2.0 cm Width (1899200)
Weight	203 g 1899200
Speed	50-7500 rpm oscillate 50-30000 rpm forward
Duty Cycle for Applied Part	The StraightShot M5 handpiece under full load is rated for intermittent operation per the following: Maximum On Time: 60 seconds Minimum Off Time: 30 seconds

STRAIGHTSHOT MAGNUM II AND STRAIGHTSHOT III TECHNICAL SPECIFICATIONS

StraightShot Magnum II 1897200

StraightShot III 1897201

Size	17 cm Length x 1.6 cm Width (1897200/1897201)
Weight	240 g 1897200, 1897201 254 g 1897200T
Speed	50-5000 rpm oscillate 50-12000 rpm forward
Duty Cycle for Applied Part	Under full load is rated for intermittent operation per the following: Maximum On Time: 60 seconds Minimum Off Time: 30 seconds

LEGEND EHS AND LEGEND EHS STYLUS

The Legend EHS motor (Figure 7-1) is a high speed, high torque, reversible electric motor used to dissect bone and biomaterial at selectable speeds from 200 to 75000 rpm.

The Legend EHS Stylus motor (Figure 7-2) is a smaller, compact, high speed, high torque, reversible electric motor used to dissect bone and biomaterials at selectable speeds from 200 to 75000 rpm.

The following instructions for the Legend EHS and Legend EHS Stylus are in addition to "Set up the IPC" general assembly instructions. Complete IPC setup, then continue to the instructions below.

Figure 7-1. Legend EHS Motor



Figure 7-2. Legend EHS Stylus Motor

LEGEND EHS AND LEGEND EHS STYLUS ATTACHMENT ASSEMBLY

Refer to "Legend EHS, Legend EHS Stylus and Stylus Touch Attachments" for attachment assembly instructions.

CONNECT LEGEND EHS CABLE TO MOTOR

Note: The Legend EHS Stylus motor cable is integrated in the handpiece and cannot be removed from the motor.

To connect the Legend EHS cable to the Legend EHS motor, align the mark on the cable to the red mark on the motor (Figure 7-3) and connect the two pieces.

REMOVE LEGEND EHS CABLE FROM MOTOR

Note: The Legend EHS Stylus motor cable is integrated in the handpiece and cannot be removed from the motor.

To remove the Legend EHS cable from the Legend EHS motor, (1) push the cable towards the motor, then (2) pull out the cable by the locking ring ONLY (Figure 7-4).

Figure 7-3. Legend EHS Motor Cable Connection



Figure 7-4. Disconnect Legend EHS Cable from Motor

CONNECT LEGEND EHS OR LEGEND EHS STYLUS TO IPC

Locate the Legend EHS or Legend EHS Stylus connection port on the connector panel (Figure 7-5) and insert the connector.

Note: To insert multi-pin connectors (indicated by a silver or red mark on the connector), align the mark on the connector to the mark on the console, then insert the connector. The IPC incorporates the Legend EHS or Legend EHS Stylus irrigation at pump 1. When the system detects the Legend EHS or Legend EHS Stylus, pump 2 defaults to **None**.

When the IPC detects both the Legend EHS Stylus handpieces and the StraightShot M4, the system defaults Pump 2 to the **Shared** configuration. You must manually move the irrigation tubing from the inactive to the active handpiece. Use the pumps screen to override the **Shared** default by selecting the StraightShot M4 or the Legend EHS Stylus for Pump 1. Refer to “Set up and Prime Pumps” for more information.

Control operation of the Legend EHS or Legend EHS Stylus with the IPC touchscreen and the multifunction foot pedal.

LEGEND EHS AND LEGEND EHS STYLUS TOUCHSCREEN CONTROLS

To set or adjust Legend EHS or Legend EHS Stylus controls, on the IPC touchscreen, in the control box (Figure 7-6), do the following:

- To change rotation mode, in the Mode control box, select FWD (forward) or REV (reverse).
Important: System configuration may be different from the default. If the REV (reverse) button appears raised (Figure 7-7) and does not have a selectable radio button, you cannot select the Reverse mode. If the REV button appears concave (Figure 7-7) and has a selectable radio button, you can select the Reverse mode via the touchscreen or the multifunction foot pedal.
- To adjust speed in Forward (FWD) or Reverse (REV) mode, in the Speed control box, press the plus button to increase speed or the minus button to decrease speed. 
Default, 70000 rpm; variable adjustment from 200 to 75000 rpm.
- To adjust the irrigation flow rate, in the Pump control box, press the plus button to increase flow rate or the minus button to decrease flow rate. If intermittent flow is available, pressing the plus or minus button progresses the system through intermittent and continuous flow. The system displays **Intermittent** when in intermittent flow mode.
Default, 0cc per minute in Forward or Reverse mode.
Note: To adjust flow rate, you can use the touchscreen or the IntelliFlow Irrigation remote control.
- To adjust acceleration and deceleration, on the Defaults menu (Figure 7-8), press the corresponding plus or minus button. Acceleration is the rate at which the motor speeds up to reach the target speed. Deceleration is the rate at which the motor slows down to reach the target speed or stop.
Note: While in the default menu, the motor will not be active to demonstrate the selected acceleration and deceleration. In order to determine the desired values, adjust the acceleration and deceleration during motor operation and note the preferred values.
- To enable adjustment of acceleration or deceleration during handpiece operation and to display the acceleration and deceleration adjustment options on the motor operational screen; On the Defaults menu, select **Show On Screen** and then press **[OK]**. To hide the acceleration and deceleration during handpiece operation, on the Defaults menu, deselect **Show On Screen**.

Figure 7-5. Legend EHS and Legend EHS Stylus Connection Ports

- 1 Legend EHS motor connection port
- 2 Legend EHS Stylus motor connection port
- 3 Multifunction foot pedal connection port

Figure 7-6. Legend EHS Stylus Touchscreen

Figure 7-7. Legend EHS or Legend EHS Stylus Mode



Figure 7-8. Stylus Default Setting Menu

LEGEND EHS AND LEGEND EHS STYLUS MULTIFUNCTION FOOT PEDAL CONTROLS

Important: By default, press each button on the foot pedal for at least 100 mS for the selection to become active. Use the IPC touch screen Settings screen to change the default value.


To use the multifunction foot pedal (Figure 7-9) to control the handpiece do the following:

- To select Forward or Reverse mode, press the mode button.
- To start or adjust the speed of a handpiece in variable mode, press the foot pedal.
- To toggle between the start/stop mode and variable speed mode, press the control button.

Note: Functionality of the control button may be changed in the handpiece Defaults menu to pause irrigation. Refer to Change System Settings in the Pre-Operating Instructions section.

- To change the handpiece, press the handpiece button.

Figure 7-9. Multifunction Foot Pedal & Y-Splitter

- | | |
|---------------------------|--------------|
| 1 Mode button | 5 Foot pedal |
| 2 Handpiece button | 6 Y-Splitter |
| 3 Control button | 7 Port 1 |
| 4 Slip-resistant foot pad | 8 Port 2 |

LEGEND EHS AND LEGEND EHS STYLUS CLEANING AND STERILIZATION INSTRUCTIONS

Refer to documents M000030A234 and M000030A235 in the Cleaning and Sterilization section.

LEGEND EHS TECHNICAL SPECIFICATIONS

Legend EHS EM100-A

Size	9.02 cm Length x 2.03 cm Diameter
Weight	180 g
Speed	75000 rpm forward/reverse
Duty Cycle for Applied Part	For use in operating room temperatures up to 40°C (104°F), the Legend EHS Motor is rated for a cutting time of 3 minutes at 70000 rpm.
	For normal operating room temperatures (typically 20°C/68°F), the Legend EHS Motor is rated for cutting time of 10 minutes followed by 25 minutes or rest.
	The Legend EHS Motor is rated for intermittent use of 20 seconds ON / 20 seconds OFF, indefinitely at 70000 rpm.

LEGEND EHS STYLUS TECHNICAL SPECIFICATIONS

Legend EHS Stylus EM200


Size	7.77 cm Length x 1.65 cm Diameter
Weight	90 g
Speed	75000 rpm forward/reverse
Duty Cycle for Applied Part	For use in operating room temperatures up to 40°C (104°F), the Legend EHS Stylus Motor is rated for 3 minutes at 60000 rpm followed by 25 minutes or rest.
	For normal operating room temperatures (typically 20°C/68°F), the Legend EHS Stylus Motor is rated for cutting indefinitely at 60000 rpm.

STYLUS TOUCH

The Stylus Touch (Figure 8-1) is a small, compact, high-speed, high-torque, reversible electric motor used to dissect bone and biomaterials at variable speeds from 200 to 75000 rpm. The Stylus Touch motor includes a rotating finger lever that emulates the functions of the multifunction foot pedal.

The following instructions for the Stylus Touch are in addition to “Set up the IPC” general assembly instructions. Complete IPC setup, then continue to the instructions below.

Figure 8-1. Stylus Touch Motor


STYLUS TOUCH ATTACHMENT ASSEMBLY

Refer to “Legend EHS, Legend EHS Stylus and Stylus Touch Attachments” for attachment assembly instructions.

ROTATE STYLUS TOUCH FINGER LEVER

1. Press the control lever ring forward (Figure 8-2).
2. Rotate the lever clockwise or counter clockwise until the lever locks in a new position.

Figure 8-2. Stylus Touch Finger Lever Rotation

CONNECT STYLUS TOUCH TO IPC CONSOLE


Locate the Stylus Touch connection port on the connector panel (Figure 8-3) and insert the connector.

Note: To insert multi-pin connectors (indicated by a silver or red mark on the connector), align the mark on the connector to the mark on the console, then insert the connector.

The IPC incorporates the Stylus Touch at pump 1. If you do not use irrigation for the Stylus Touch, manually change the Pump 1 to **None**. Refer to “Set up and Prime Pumps” for more information.

Control operation of the Stylus Touch with the IPC touchscreen and the multifunction foot pedal.

Figure 8-3. Stylus Touch Connection Ports

- | | |
|---|--|
| 1 | Stylus Touch connection port |
| 2 | Multifunction foot pedal connection port |

STYLUS TOUCH TOUCHSCREEN CONTROLS

To set or adjust Stylus Touch controls, on the IPC touchscreen, in the control box (Figure 8-4), do the following:

- To change rotation mode, in the Mode control box, select FWD (forward) or REV (reverse).
Important: System configuration may be different from the default. If the REV (reverse) button appears raised (Figure 8-7) and does not have a selectable radio button, you cannot select the Reverse mode. If the REV button appears concave (Figure 8-7) and has a selectable radio button, you can select the Reverse mode via the touchscreen or the multifunction foot pedal.
- To adjust speed, in the Speed control box, press the plus button to increase speed or the minus button to decrease speed. 
Default, 60000 rpm; variable adjustment from 200 to 75000 rpm.
Note: When the handpiece is in safe mode, the Speed control box displays SAFE (Figure 8-6). Refer to the Stylus Touch Safe Mode topic for additional information.
- To adjust the irrigation flow rate, in the Pump control box, press the plus button to increase flow rate or the minus button to decrease flow rate. If intermittent flow is available, pressing the plus or minus button progresses the system through intermittent and continuous flow. The system displays Intermittent when in intermittent mode.
Default, 0cc per minute.
Note: To adjust flow rate, you can use the touchscreen or the IntelliFlow Irrigation remote control.
- To adjust acceleration and deceleration, on the Defaults menu (Figure 8-5), press the corresponding plus or minus button. Acceleration is the rate at which the motor speeds up to reach the target speed. Deceleration is the rate at which the motor slows down to reach the target speed or stop.
Note: While in the default menu, the motor will not be active to demonstrate the selected acceleration and deceleration. In order to determine the desired values, adjust the acceleration and deceleration during motor operation and note the preferred values.
- To enable the adjustment of acceleration or deceleration during handpiece operation and to display the acceleration and deceleration adjustment options on the motor operational screen; On the Defaults menu, select **Show On Screen** and then press [OK]. To hide the acceleration and deceleration during handpiece operation, on the Defaults menu, deselect **Show On Screen**.

STYLUS TOUCH SAFE MODE

When the handpiece is in safe mode, it is inoperable until the safety is turned off. The Speed control box on the Stylus Touch touchscreen displays SAFE when the handpiece is the active handpiece and in safe mode (Figure 8-6).

Switch the device to safe mode any time it is attached to the console, but not currently being used. To set the Stylus Touch to Safe Mode, on the Stylus Touch handpiece, switch the Safe Mode finger lever (Figure 8-1) to on.

When more than one handpiece is attached to the console, use the safety switch of an inactive handpiece to activate that handpiece and make it ready for use.

Figure 8-4. Stylus Touch Touchscreen



Figure 8-5. Stylus Touch Default Setting Menu

Figure 8-6. Stylus Touch Safe Mode

Figure 8-7. Stylus Touch Rotation Mode


STYLUS TOUCH MULTIFUNCTION FOOT PEDAL CONTROLS

Important: By default, press each button on the foot pedal for at least 100 mS for the selection to become active. Use the IPC touch screen Settings screen to change the default value.

To use the multifunction foot pedal (Figure 8-8) to control the handpiece do the following:

- To select Forward or Reverse mode, press the mode button.
 - To start or adjust the speed of a handpiece in variable mode, press the foot pedal.
 - To toggle between the start/stop mode and variable speed mode, press the control button.
- Note:** Functionality of the control button may be changed in the handpiece Defaults menu to pause irrigation. Refer to Change System Settings in the Pre-Operating Instructions section.
- To change the handpiece, press the handpiece button.

Figure 8-8. Multifunction Foot Pedal & Y-Splitter

- | | |
|---------------------------|--------------|
| 1 Mode button | 5 Foot pedal |
| 2 Handpiece button | 6 Y-Splitter |
| 3 Control button | 7 Port 1 |
| 4 Slip-resistant foot pad | 8 Port 2 |

STYLUS TOUCH CLEANING AND STERILIZATION INSTRUCTIONS

Refer to documents 68E4132 and M000030A235 in the Cleaning and Sterilization section.

STYLUS TOUCH TECHNICAL SPECIFICATIONS

Stylus Touch EM210	
Size	15.26 cm Length x 1.65 cm Diameter
Weight	130 g
Speed	75000 rpm forward/reverse
Duty Cycle for Applied Part	For use in operating room temperatures up to 40°C, the Legend EHS Stylus Touch Motor is rated for 3 minutes at 60000 rpm followed by 25 minutes of rest.
	For normal operating room temperatures (typically 20°C), the Legend EHS Stylus Touch Motor is rated for cutting indefinitely at 60000 rpm.


LEGEND EHS, LEGEND EHS STYLUS AND STYLUS TOUCH ATTACHMENTS

The following instructions for the Legend EHS, Legend EHS Stylus and Stylus Touch attachments are in addition to “Set up the IPC” and motor-specific general instructions. Complete IPC setup, motor specific instructions and then continue to the instructions below.

DISSECTING TOOLS NOMENCLATURE

Note: Match the nomenclature and color code on the dissecting tool packaging to the same nomenclature and color code on the attachment.

Part numbers for dissecting tools follow a standard naming convention (Figure 9-1). A basic part number consists of five characters, representing the associated attachment length, the tool-head shape, and the tool-head diameter. Part numbers can also include a variety of prefixes to identify specific attachment types, as well as a variety of suffixes to provide additional information about the dissecting tool. Tools that use a design taken from the Mednext line are designated by an additional “-MN” suffix.

- | | | | |
|---|-------------------|---|--------------------------------------|
| 1 | Optional prefix | 4 | Tool head diameter (x.x millimeters) |
| 2 | Attachment length | 5 | Optional suffix |
| 3 | Tool head shape | 6 | Optional “-MN” suffix |

Tool Number Prefixes

- F For use with footed attachments.
- MC For use with metal-cutting attachments.
- T For use with telescoping attachments.

Tool Head Shapes

- | | | | |
|----|------------|----|---------------|
| AC | Acorn | MH | Match Head |
| BA | Ball | OV | Oval |
| CY | Cylinder | RT | Reverse Taper |
| HM | Hole Maker | TA | Tapered |
| HS | Hole Saw | TD | Twist Drill |

Tool Number Suffixes

Note: More than one suffix can be combined in a single part number.

- | | | | |
|---|---------|----|----------------------|
| L | Long | S | Spiral |
| D | Diamond | SH | Short |
| X | Extra | | |
| F | Fine | DC | Diamond Coarse |
| C | Carbide | DX | Diamond Extra Coarse |

ALIGN MOTOR COLLET

Note: An attachment will not seat on the Legend EHS, Legend EHS Stylus or Stylus Touch motor if the arrows on the collet flats are not aligned.

1. Verify alignment of arrows on motor collet flats. Prior to installation of an attachment and dissecting tool on the Legend EHS, Legend EHS Stylus or Stylus Touch motor, ensure that arrows on the motor collet flats are in proper alignment. (Figure 9-2).
2. If the arrows are not aligned (Figure 9-3), use the motor wrench (Figure 9-4) to turn the rotational collet until its arrow is aligned with the arrow on the stationary collet.

Note: DO NOT use any other components except for the motor wrench to align arrows on motor collet.

Figure 9-2. Correct alignment

Figure 9-3. Incorrect alignment

Figure 9-4. Motor wrench

STRAIGHT ATTACHMENT ASSEMBLY


Warnings: Refer to warnings W36 and W54.

Caution: Match the nomenclature and color code on the tool packaging to the same nomenclature and color code on the attachment.

Notes:

- An attachment will not seat on the motor if the arrows on the collet flats are not in alignment.
 - The Legend EHS motors will not run properly unless the attachment is in the locked position.
 - Smoke and/or excessive heat may be generated if attachment is not in the fully locked position. This may result in thermal injury to the surgeon or staff.
1. Slide a straight attachment over the motor collet aligning triangular arrows on the attachment and the motor case (Figure 9-5). An audible click, heard and perceptible by touch, confirms that the tool is fully seated.

Figure 9-5. Straight Attachment Assembly

2. Insert the tool into the attachment with a slight rotational motion (Figure 9-6). An audible click, heard and perceptible by touch, confirms that the tool is fully seated.

Figure 9-6. Straight Attachment Assembly

3. Rotate the attachment in the direction indicated by arrow until the attachment alignment mark is directly in line with the locked symbol (Figure 9-7). You will hear two clicks as the attachment is rotated.

Figure 9-7. Straight Attachment Assembly

4. Gently pull on the tool to ensure that it is locked into the handpiece.
Note: Tool should rotate freely. If not, unlock the attachment, re-seat the tool, and re-lock the attachment.

STRAIGHT ATTACHMENT DISASSEMBLY

1. Hold the motor in palm of hand. Rotate the attachment to the unlocked position.
2. Remove the dissecting tool from the attachment and discard the tool.
3. Use thumb and index finger to lift the attachment off of the motor.

ANGLED ATTACHMENT ASSEMBLY

Caution: When using an angled attachment, hold the handpiece assembly by the attachment so that the attachment does not inadvertently loosen from the handpiece.

Notes:


- A dissecting tool may be installed and locked in the attachment before the angled attachment is installed onto the motor.
 - Angled and straight attachments with the same length, marking and color band share the same dissecting tools.
 - The Legend EHS Motors will not run properly unless the attachment is in the locked position.
 - Smoke and/or excessive heat may be generated if attachment is not in the fully locked position. This may result in thermal injury to the surgeon or staff.
1. With the tool lock in the unlocked position, insert a tool into the angled attachment with a slight rotational motion (Figure 9-8). An audible click, heard and perceptible by touch, confirms that the tool is fully seated.

Figure 9-8. Angled Attachment Assembly


2. Rotate the tool lock in the direction indicated by arrow until the tool lock alignment mark is directly in line with the locked symbol (Figure 9-9).

Figure 9-9. Angled Attachment Assembly

3. Gently pull on the tool to ensure that it is locked into the handpiece.
Note: Tool should rotate freely. If not, unlock the attachment, re-seat the tool, and re-lock the attachment.
4. Slide the angled attachment over the motor collet aligning triangular arrows on the attachment and the motor case. An audible click, heard and perceptible by touch, confirms that the tool is fully seated.
5. Rotate the attachment in the direction indicated by the arrow until attachment alignment mark is directly in line with the locked symbol. You will hear two clicks as the attachment is rotated.
6. Verify that both the attachment to motor alignment mark and the tool lock alignment mark are directly in line with the locked symbol (Figure 9-10).

Figure 9-10. Angled Attachment Assembly

ANGLED ATTACHMENT DISASSEMBLY

1. Rotate the Tool Lock to the unlocked position to remove the tool from the attachment.
2. Rotate the attachment to the unlocked position and lift attachment off of the motor.

FIXED FOOTED ATTACHMENT ASSEMBLY

Warning: Refer to Warning W37.

1. Insert a dissecting tool into the motor collet with a slight rotational motion (Figure 9-11). An audible click, heard and perceptible by touch, confirms that the tool is fully seated.

Figure 9-11. Footed Attachment Assembly


2. Slide the footed attachment over the dissecting tool onto the motor aligning triangular arrows on the attachment and the motor case (Figure 9-12).

Figure 9-12. Footed Attachment Assembly

3. Pull the footed attachment towards the motor and rotate the attachment to the locked position on the motor case (Figure 9-13).

Figure 9-13. Footed Attachment Assembly

FIXED FOOTED ATTACHMENT DISASSEMBLY

1. To remove the footed attachment, hold the motor in the palm of your hand. Push the sleeve on the footed attachment distally while rotating the attachment to the unlocked position on the motor case. Release the sleeve (Figure 9-14).
2. To avoid injury from the dissecting tool, use thumb and index finger to cautiously and slowly lift the attachment off of the motor and away from the dissecting tool (Figure 9-15).

Figure 9-14. Footed Attachment Disassembly



Figure 9-15. Footed Attachment Disassembly

3. Pull the dissecting tool out of the motor collet and discard the tool (Figure 9-16).

Figure 9-16. Footed Attachment Disassembly

ROTATING FOOTED ATTACHMENT ASSEMBLY

Warnings: Refer to Warnings, W20 and W52.

- Notes:**
- Rotating and fixed footed attachments with the same length, marking and color band share the same dissecting tools.
 - The footed end of the attachment has 360 ° of unrestricted rotation.

Attach the Rotating Footed attachment using the Fixed Footed Attachment Assembly Instructions.


ROTATING FOOTED ATTACHMENT DISASSEMBLY

Remove the Rotating Footed attachment using the Fixed Footed Attachment Disassembly Instructions.

CONTRA-ANGLE ATTACHMENT AC-16 ASSEMBLY

1. On the IPC Touchscreen, adjust the speed setting to 62,000 rpm using the speed control buttons.
2. Slide the Contra-Angle Attachment over the motor collet aligning triangular arrows on the attachment and the motor case. An audible click, heard and perceptible by touch, confirms that the tool is fully seated.
3. Rotate the attachment to the locked position on the motor case (Figure 9-17).

Figure 9-17. Contra-Angle Attachment Assembly

4. Rotate the attachment head's lever laterally to the open position and insert a dissecting tool (Figure 9-18).
5. Close lever (Figure 9-19).

Figure 9-18. Contra-Angle Attachment Assembly

Figure 9-19. Contra-Angle Attachment Assembly

6. Gently pull on the dissecting tool shaft to ensure proper installation.

CONTRA-ANGLE ATTACHMENT AC-16 DISASSEMBLY

1. Rotate the lever on the attachment head laterally to remove the dissecting tool.
2. Discard the dissecting tool.
3. Rotate the 16-MF attachment to the unlocked position and lift the attachment off of the motor.

METAL CUTTING ATTACHMENT ASSEMBLY

Warning: Refer to Warning W15.

Important: The Metal Cutting attachment uses the tungsten carbide or diamond wheel dissecting tools. All metal cutting dissecting tools have an "MC" attachment prefix in their nomenclature (for example, MC254 and MC30). Metal cutting dissecting tools cannot be installed into any other Attachment.

1. Slide the Metal Cutting attachment over the motor collet aligning triangular markers on the attachment and the motor case. An audible click, heard and perceptible by touch, confirms that the tool is fully seated.
2. Rotate the attachment to the locked position on the motor case.
3. With the tool lock unscrewed several turns and holding the attachment UPRIGHT, insert a dissecting tool into attachment (Figure 9-20).
4. Rotate the dissecting tool until it drops into position and is fully seated. You will feel a tactile click indicating that the tool is fully seated (Figure 9-20).

Figure 9-20. Metal Cutting Attachment Assembly

METAL CUTTING ATTACHMENT DISASSEMBLY

1. Unscrew the tool lock with several turns, then withdraw the dissecting tool.
2. Rotate the attachment to the locked position on the motor case and lift attachment off of the motor.


VARIABLE EXPOSURE ATTACHMENT ASSEMBLY

Warnings: Refer to Warnings W23, W34, W53 and W57.

The Variable Exposure attachments can be distinguished from standard attachments by the dual color bands on the attachment. Match the color band on the attachment to the color code on the dissecting tool packaging.

1. Assemble the attachment using the Straight Attachment Assembly Instructions.
2. After assembly, use the TUBE adjustment ring to adjust the exposure of the dissecting tool (Figure 9-21).
With the tool pointing away from you, turn the ring to the right to increase the length of the tube, thereby decreasing the exposure of the tool. Turn the ring to the left to decrease the length of the tube, thereby increasing the exposure of the tool.

Figure 9-21. Variable Exposure Attachment Assembly


VARIABLE EXPOSURE ATTACHMENT DISASSEMBLY

Remove the Variable Exposure Attachment using the Straight Attachment Disassembly Instructions.

TELESCOPING STRAIGHT ATTACHMENT AT10 ASSEMBLY

- 1. Assemble the attachment using the Straight Attachment Assembly Instructions.
 - 2. Insert the base end of the selected telescoping tube into the attachment.
 - 3. To lock tube in place, turn the TUBE Locking Ring clockwise until finger tight (Figure 9-22).
- Note:** DO NOT over tighten.

Figure 9-22. Telescoping Straight Attachment Assembly


- 4. Be sure that the TOOL Locking Ring is in the unlocked position. Insert the dissecting tool into the telescoping tube (Figure 9-23). Tactile feedback is felt when the dissecting tool is fully seated.

Figure 9-23. Telescoping Straight Attachment Assembly

- 5. Turn the TOOL Locking Ring to the locked position (Figure 9-24).

Figure 9-24. Telescoping Straight Attachment Assembly


- 6. Verify that the tool is in place by gently pulling on the tool.
- 7. If the tube position needs to be changed, rotate the TUBE Locking Ring toward the unlocked position, re-position the tube, then rotate the TUBE Locking Ring toward locked.
- 8. Gently pull on the dissecting tool, then the tube, to ensure proper installation.

TELESCOPING CURVED ATTACHMENT AT10 ASSEMBLY

- 1. Assemble the attachment using the Straight Attachment Assembly Instructions.
 - 2. Insert the base end of the curved bur into the attachment until the hub is fully seated. To lock in place, turn the TUBE Locking Ring until finger tight (Figure 9-25).
- Note:** DO NOT over tighten.

Figure 9-25. Telescoping Curved Attachment Assembly

- 3. Verify that the hub is in place by gently pulling on the tool.
- 4. Seat the tool in the tool Locking Ring by applying a slight amount of inward pressure on the bur (Figure 9-26). An audible click, heard and perceptible by touch, confirms that the tool is fully seated.

Figure 9-26. Telescoping Curved Attachment Assembly

- 5. Rotate the Tool Locking Ring until the tool lock alignment mark is directly in line with the locked symbol (Figure 9-27).
- 6. Verify that the bur is in place by gently pulling on the bur.
- 7. Prior to initial use, soak the cooling sleeve by dipping it into a cup of saline or DI water (Figure 9-28).
- 8. During use, maintain copious irrigation of the cooling sleeve and bur by dribbling saline or DI water along the entire length of the cooling sleeve.

Figure 9-27. Telescoping Curved Attachment Assembly


Figure 9-28. Curved Bur Cooling

TELESCOPING TUBE ATTACHMENT ASSEMBLY

- 1. Assemble the attachment to the motor using the Straight Attachment Assembly Instructions.
 - 2. Insert the base end of the tool into the attachment until the hub is fully seated. To lock in place, turn the TUBE Locking Ring until finger tight (Figure 9-27).
- Note:** DO NOT over tighten.
- 3. Verify that the hub is in place by gently pulling on the tool.
 - 4. Seat the tool in the tool Locking Ring by applying a slight amount of inward pressure on the bur (Figure 9-27). An audible click, heard and perceptible by touch, confirms that the tool is fully seated.

TELESCOPING ATTACHMENT AT10 DISASSEMBLY

- 1. Rotate the TUBE Locking Ring towards the unlocked position.
- 2. Rotate the TOOL Locking Ring to the unlocked position.
- 3. Pull the telescoping tube out of the attachment.
- 4. Rotate the attachment to the unlocked position on the motor case and lift the attachment off the motor.

PERFORATOR ATTACHMENT AD01 & AD03 ASSEMBLY

Warning: Refer to Warning W56.

Note: A cranial perforator device may be installed in the attachment before the perforator attachment is installed on the motor.

Console Setting	Maximum Speed	
	AD01 Output Speed (MAX)	AD03 Output Speed (MAX)
60000 rpm	645 rpm	830 rpm
70000 rpm	745 rpm	965 rpm
72000 rpm	770 rpm	995 rpm
74000 rpm	790 rpm	1020 rpm
75000 rpm	805 rpm	1035 rpm

- 1. Assemble the attachment using the Straight Attachment Assembly Instructions.
- 2. To install a cranial perforator device with a Hudson shank, pull back proximally on the collar of the Perforator Attachment. Insert a device and release the collar to its original position (Figure 9-29).

Figure 9-29. Perforator Attachment Hudson Shank Assembly

PERFORATOR ATTACHMENT AD01 & AD03 DISASSEMBLY

Pull back proximally on the collar of the Perforator Attachment to remove the cranial perforator device. Rotate the Perforator Attachment to the unlocked position and lift the attachment off of the motor.

JACOBS CHUCK ATTACHMENT ASSEMBLY

- 1. Assemble the attachment using the Straight Attachment Assembly Instructions.
- 2. To install a drill bit, turn Jacobs key to open ridged collar. Insert drill bit and tighten collar with key (Figure 9-30).

Figure 9-30. Jacobs Chuck Tool Assembly

JACOBS CHUCK ATTACHMENT DISASSEMBLY

Use the Jacobs key to open the collar. Remove and discard the drill bit in an appropriate container. Rotate the Jacobs Chuck Attachment to the unlocked position and lift the attachment off of the motor.

ANGLED DOUBLE LOCK ATTACHMENT ASSEMBLY

Notes:

- Angled attachments with the same length, marking, and color band share the same dissecting tools.
 - You can insert and lock a tool in the attachment before the angled attachment is installed on the motor.
1. Assemble the attachment using the Straight Attachment Assembly Instructions.
 2. Insert the tool into the attachment with a slight rotational motion. An audible click, perceptible by touch, confirms that the tool is fully seated (Figure 9-31).
 3. Rotate the tool lock in direction indicated by arrow until the tool lock alignment mark is directly in line with the locked symbol (Figure 9-32).
 4. Pull on the tool to ensure that it is locked into the handpiece.
 5. The tool should rotate freely. If not, unlock the attachment, re-seat the tool, and re-lock the attachment.
 6. Verify that both the attachment to motor and the tool-lock alignment mark is directly in line with the locked symbol.

Figure 9-31. Angled Double Lock Attachment Assembly



Figure 9-32. Angled Double Lock Attachment Assembly

ANGLED DOUBLE LOCK ATTACHMENT DISASSEMBLY

- 1. To remove the attachment, hold the motor in the palm of your hand, and push the sleeve on the attachment distally while turning the attachment to the unlocked position (Figure 9-33).
- 2. Release the sleeve and remove the attachment.

Figure 9-33. Angled Double Lock Attachment Disassembly

TRANS-NASAL SKULL BASE BUR ATTACHMENT ASSEMBLY

- 1. Slide the attachment over the motor collet aligning triangular arrows on the attachment and the motor case. An audible click, heard and perceptible by touch, confirms that the tool is fully seated.
- 2. Rotate the attachment to the locked position on the motor case (Figure 9-34).
- 3. Attach irrigation tubing (Figure 9-35).

Figure 9-34. Trans-Nasal Skull Base Bur Assembly


Figure 9-35. Trans-Nasal Skull Base Bur Assembly

1 Irrigation tubing connection

TRANS-NASAL SKULL BASE BUR ATTACHMENT DISASSEMBLY


- 1. Remove irrigation tubing.
- 2. Rotate the attachment to the unlocked position on the motor case.
- 3. Remove the attachment from the motor.

IRRIGATION ASSEMBLY

NOTE: Clip may not fasten to small bore attachment after having been used on large bore attachment.

- 1. Adjust the plastic clip on the stainless-steel irrigation tube (Figure 9-36).
- 2. Bend the irrigation tube to the desired angle.
- 3. Snap the clip onto the handpiece near the tool.

Figure 9-36. Irrigation Assembly

LEGEND EHS, LEGEND EHS STYLUS AND STYLUS TOUCH CLEANING AND STERILIZATION INSTRUCTIONS

Refer to documents M000030A234, 68E4132 and M000030A235 in the Cleaning and Sterilization section.

Variable Exposure Attachment Cleaning

When cleaning, clean the attachment completely. First without adjusting the tube length, then with the tube fully extended, and finally with the tube fully retracted.


SKEETER ULTRA-LITE OTO-TOOL

The following instructions for the Skeeter Ultra-Lite Oto-Tool are in addition to “Set up the IPC” general assembly instructions. Complete the IPC setup, then continue to the instructions below. Refer to the Skeeter Ultra-Lite Oto-Tool System Instructions for Use for additional information.

SKEETER ASSEMBLY

1. Press the bur release button (Figure 10-1).
2. Load the desired bur for the procedure into the handpiece by inserting the bur shaft through the distal end of the handpiece with a slight twisting motion while simultaneously pressing the bur release button.
3. The bur is locked into place when a “click” is noted. Locking of the bur should be checked prior to use by firmly pulling on the bur after the “click” is noted.
4. Tug the bur to ensure it fits securely in the handpiece.
5. To remove the bur from the handpiece, press the bur release button on the handpiece and pull the bur out.

Figure 10-1. Skeeter Handpiece Assembly

CONNECT SKEETER HANDPIECE TO IPC CONSOLE

On the IPC Console, locate the Skeeter accessory connection port on the connector panel (Figure 10-2), align the mark on the connector to the mark on the console, then insert the connector.

Figure 10-2. IPC Skeeter Connection Port

The Skeeter does not use irrigation. By default, the system sets both pumps to **None**.

Control the operation of the Skeeter with the IPC touchscreen and the multifunction foot pedal.

SKEETER TOUCHSCREEN CONTROLS

To set or adjust Skeeter controls, on the IPC touchscreen, in the Skeeter control box (Figure 10-3), do the following:

- To change rotation mode, in the Mode control box, select FWD (forward) or REV (reverse).
 - To adjust speed, in the Speed control box, press the plus button to increase speed or the minus button to decrease speed.
- Default, 16000 rpm; variable adjustment from 1000 to 16000 rpm.

Figure 10-3. Skeeter Touchscreen


SKEETER MULTIFUNCTION FOOT PEDAL CONTROLS

Important: By default, press each button on the foot pedal for at least 100 mS for the selection to become active. Use the IPC touch screen Settings screen to change the default value.

To use the multifunction foot pedal (Figure 10-4) to control the handpiece do the following:

- To select Forward or Reverse mode, press the mode button.
- To start or adjust the speed of a handpiece in variable mode, press the foot pedal.
- To toggle between the start/stop mode and variable speed mode, press the control button.
- To change the handpiece, press the handpiece button.

Figure 10-4. Multifunction Foot Pedal & Y- Splitter

SKEETER

SKEETER CLEANING AND STERILIZATION INSTRUCTIONS

Refer to documents 68E3968 and 68E3969 in the Cleaning and Sterilization section.

SKEETER TECHNICAL SPECIFICATIONS

Skeeter 3055601

Size	17 cm Length x 1.6 cm Diameter
Weight	57 g
Speed	1000-16000 rpm forward/reverse
Duty Cycle for Applied Part	Continuous run
Storage Temperature	-40 °C to +70 °C
Humidity	10 % to 100 % RH
Barometric Pressure	500 to 1060 hPa

VISAO HIGH-SPEED DRILL

DEVICE DESCRIPTION

The IPC incorporates the Visao coolant at pump 1 and irrigation at pump 2. The Visao coolant pump cartridge uses both a pump tube and a drip chamber return tube. If you will not use irrigation, select None for pump 2. Refer to “Set up and Prime Pumps”. Control operation of the Visao with the IPC touchscreen and the multifunction foot pedal.

Important: During use, maintain free flow irrigation to the cooling sleeve and bur by dribbling saline or DI water along the entire length of the cooling sleeve.

The following instructions are specific for the use of the Visao with the IPC. They are in addition to “Set up the IPC” general accessory instructions.

SET UP VISAO PUMP

Fill the clear-drip chamber with coolant before you prime the coolant system. Refer to “Set up and Prime Pumps” for instruction.

VISAO BUR GUARD ASSEMBLY

Important: On the Visao, a bur guard (Figure 11-2) is required for use with all burs.

Slide bur guard over the front end of the Visao until fully seated (Figure 11-1).

Figure 11-1: Bur Guard Assembly



Figure 11-2: Visao Bur Guards

Visao Reusable Bur Guards with and without irrigation.

Visao Single Use Bur Guard with irrigation.


Visao Single Use STIM Bur Guard with and without irrigation. The STIM Bur Guard also provides nerve stimulation to standard burs in static and dynamic modes when used with both the Medtronic Nerve Integrity Monitor (NIM) and the IPC System.

VISAO STRAIGHT BUR ASSEMBLY

Important: On the Visao, a bur guard (Figure 11-2) is required for use with all burs.

1. Align the alignment point on the Visao locking collar with the unlock symbol (Figure 11-3).
2. Insert the tool until fully seated.
3. Move the locking collar so that the alignment point aligns with the lock symbol.
4. To ensure a secure fit, gently pull the tool.

Figure 11-3: Visao Straight Tool Assembly


VISAO CURVED BUR ASSEMBLY

Important: On the Visao, a bur guard (Figure 11-2) is required for use with all burs.

- 1. Align the alignment point on the Visao locking collar with the unlock symbol (Figure 11-4).
- 2. Align the notch on the tool with the notch on the Visao collar.
- 3. Gently press on the tool until full seated.
- 4. Align the alignment point on the Visao locking collar with the lock symbol.
- 5. To ensure a secure fit, gently pull the tool.

Figure 11-4. Visao Curved Bur Assembly


VISAO IRRIGANT AND COOLING TUBE ASSEMBLY

Caution: Do not confuse the coolant tube with irrigation tube.

- 1. Connect one coolant tube to each coolant port (Figure 11-5). The coolant can flow in the left or right port if you connect the return tubing to the opposite port.
 - 2. Connect the free end of the irrigation tube to the irrigation barb.
 - 3. Confirm the coolant pedal starts the handpiece and coolant flow.
- Note:** The coolant pump runs for 1 minute after you release the pedal
- 4. Prior to initial use, soak the cooling sleeve in a cup of saline or DI water.
 - 5. During use, maintain copious irrigation of the cooling sleeve and bur by dribbling saline or DI water along the entire length of the cooling sleeve.

Figure 11-5. Visao Irrigant and Cooling Assembly

CONNECT VISAO TO IPC CONSOLE

Locate the Visao connection port on the connector panel (Figure 11-6) and insert the connector.

Note: To insert multi-pin connectors (indicated by a silver or red mark on the connector), align the mark on the connector to the mark on the console, then insert the connector.

VISAO TOUCHSCREEN CONTROLS

To set or adjust Visao controls, on the IPC touchscreen, in the Visao control box (Figure 11-7), do the following:

- To change rotation mode, in the Mode control box, select FWD (forward) or REV (reverse).
 - To adjust speed, in the Speed control box, press the plus button to increase speed or the minus button to decrease speed. 
Default, 80000 rpm; variable adjustment from 200 to 80000 rpm.
 - To adjust the irrigation flow rate, in the Pump control box, press the plus button to increase flow rate or the minus button to decrease flow rate. If intermittent flow is available, pressing the plus or minus button progresses the system through intermittent and continuous flow. The system displays Intermittent when in intermittent mode.
Forward Mode: Default, 30 cc per minute.
- Note:** To adjust flow rate, you can use the touchscreen or the IntelliFlow Irrigation remote control.

Figure 11-6. IPC Visao Connection Port



Figure 11-7. Visao Touchscreen


VISAO MULTIFUNCTION FOOT PEDAL CONTROLS

Important: By default, press each button on the foot pedal for at least 100 mS for the selection to become active. Use the IPC touch screen Settings screen to change the default value.
To use the multifunction foot pedal (Figure 11-8) to control the handpiece do the following:

- To select Forward or Reverse mode, press the mode button.
 - To start or adjust the speed of a handpiece in variable mode, press the foot pedal.
 - To toggle between the start/stop mode and variable speed mode, press the control button.
- Note:** Functionality of the control button may be changed in the handpiece Defaults menu to pause irrigation. Refer to Change System Settings in the Pre-Operating Instructions section.
- To change the handpiece, press the handpiece button.

Figure 11-8. Multifunction Foot Pedal & Y-Splitter

- | | |
|---------------------------|--------------|
| 1 Mode button | 5 Foot pedal |
| 2 Handpiece button | 6 Y-Splitter |
| 3 Control button | 7 Port 1 |
| 4 Slip-resistant foot pad | 8 Port 2 |

VISAO CLEANING AND STERILIZATION INSTRUCTIONS

Refer to document 68E3281 in the Cleaning and Sterilization section.

VISAO TECHNICAL SPECIFICATIONS

Visao 3334800

Size	16.0 cm Length x 2.0 cm Diameter
Weight	148 g
Speed	200-80000 rpm forward/reverse. Visao High-Speed Drill, water cooled
Duty Cycle for Applied Part	The Visao High-Speed Drill under full load is rated for intermittent operation per the following:
	Maximum On Time: 60 seconds
	Maximum Off Time: 30 seconds

VISAO USE WITH IPC AND NIM

To use the NIM with the IPC use the Stim Bur Guard (Figure 11-2). The Stim Bur Guard connects the IPC to the NIM via the Stim bur stimulus output port (Figure 11-6). The Stim Bur Guard carries stimulating current to the tool's tip and provides nerve stimulation to standard burs in static and dynamic modes. Refer to the NIM User's Manual and Stim Bur Guard Product Information and Instructions for further information.


INDIGO HIGH-SPEED OTOLOGIC DRILL

DEVICE DESCRIPTION

The Indigo drill is a small, compact high-speed, high-torque, reversible electric drill that can be used to dissect bone and biomaterial at variable speeds from 200 to 60000 rpm. The cable cannot be removed from the drill. The following instructions for the Indigo drill are in addition to "Set up the IPC" general assembly instructions. Complete IPC setup, then continue to the instructions below.

INDIGO DRILL STRAIGHT ATTACHMENT ASSEMBLY

Figure 12-1. Indigo Alignment


Important: The Indigo High-speed Otologic drill is designed to work only with Medtronic Xomed burs and Indigo attachments. The use of other burs and attachments may result in sub-standard performance and will void the manufacturer's warranty. Do not use with Midas Rex attachments or burs.

1. Verify the alignment marks on the motor collet are in alignment (Figure 12-1).
- Note:** If the marks are misaligned, turn the collet until the marks are aligned.
2. Slide the attachment over the motor collet (Figure 12-2) so that the alignment mark on the attachment aligns with the alignment mark (unlocked symbol) on the motor collet.
3. Insert the bur with a slight twisting motion until you feel it seat into position.
4. Turn the attachment so the alignment mark aligns with the locked symbol on the motor collet. You will hear two clicks while rotating the attachment.
5. To ensure a secure fit, gently pull the tool.

Important: The Indigo motor will not run correctly unless the attachment is in the locked position.

Warning: Smoke may be generated if attachment is not in the locked position (W36).

Figure 12-2. Indigo Straight Attachment Assembly


INDIGO DRILL ANGLED ATTACHMENT ASSEMBLY

Important: The Indigo High-speed Otologic drill is designed to work only with Medtronic Xomed burs and Indigo attachments. The use of other burs and attachments may result in sub-standard performance and will void the manufacturer’s warranty. Do not use with Midas Rex attachments or burs.

- 1. Verify the alignment marks on the motor collet are in alignment (Figure 12-1).
- Note:** If the marks are misaligned, turn the collet until the marks are aligned.
- 2. Slide the attachment over the motor collet (Figure 12-3) so that the alignment mark on the attachment aligns with the alignment mark (unlocked symbol) on the motor collet.
- 3. Turn the attachment so the alignment mark aligns with the locked symbol on the motor collet. You will hear two clicks while rotating the attachment.
- 4. Verify the alignment marks on the tool lock ring align with the unlocked symbol on the attachment.
- 5. Insert the bur with a slight twisting motion until you feel it seat into position.
- 6. Turn the lock ring so the alignment mark aligns with the locked symbol on the attachment. You will hear two clicks while rotating the attachment.
- 7. To ensure a secure fit, gently pull the tool.

Important: The Indigo motor will not run correctly unless the attachment and lock ring is in the locked position.


Warning: Smoke may be generated if attachment is not in the locked position (W36).

INDIGO DRILL IRRIGATION AND TUBING CLIP ASSEMBLY

- 1. Snap or slide the irrigation clip on to the attachment (Figure 12-4).
- Note:** The irrigation clip is designed with two nodules that fit in the grooves on the attachment.
- 2. Bend the irrigation tube to a desirable angle and adjust the clip closer to or farther from the bur, as necessary.
- 3. Snap the tubing clip on the handpiece or angled attachment.

Figure 12-4. Irrigation and Tubing Clip Assembly

CONNECT INDIGO DRILL TO IPC


Locate the Indigo drill connection port on the connector panel (Figure 12-5) and insert the connector.

Note: To insert multi-pin connectors (indicated by a silver or red mark on the connector), align the mark on the connector to the mark on the console, then insert the connector.

The IPC incorporates Indigo on-drill irrigation at pump 2. If you do not use drill irrigation, select **None** for pump 2. If you will use a suction irrigator, you must manually select it for pump 1 or pump 2.

Control operation of the Indigo drill with the IPC touchscreen and the multifunction foot pedal.

Figure 12-5. IPC Indigo Connection Ports

INDIGO DRILL TOUCHSCREEN CONTROLS

To set or adjust Indigo drill controls, on the IPC touchscreen, in the control box (Figure 12-6), do the following:

- To change rotation mode, in the Mode control box, select FWD (forward) or REV (reverse).
Important: System configuration may be different from the default. If the REV (reverse) button appears raised and does not have a selectable radio button (Figure 12-7), you cannot select the Reverse mode. If the REV button appears concave (Figure 12-7) and has a selectable radio button, you can select the Reverse mode via the touchscreen or the multifunction foot pedal.

- To adjust speed, in the Speed control box, press the plus button to increase speed or the minus button to decrease speed.

Default, 52000 rpm; variable adjustment from 200 to 60000 rpm in increments of 200, 500 and then 1000 until 60000.

Note: The speed you set remains constant when you switch between modes.

Warning: A tool's size and geometry may create excessive vibration at certain speeds. Increase or decrease speed on console. Change to a new tool to prevent unintended tissue removal from patient (W54).


- To adjust the irrigation flow rate, in the Pump control box, press the plus button to increase flow rate or the minus button to decrease flow rate. If intermittent flow is available, pressing the plus or minus button progresses the system through intermittent and continuous flow of irrigation at 5, 10, 15 and 20 cc per minute. The system displays **Intermittent** when in intermittent mode.
Default, 30cc per minute.

Note: To adjust flow rate, you can use the touchscreen or the IntelliFlow Irrigation remote control.

Figure 12-6. Indigo Touchscreen

Figure 12-7. Indigo Mode

INDIGO DRILL MULTIFUNCTION FOOT PEDAL CONTROLS

Important: By default, press each button on the foot pedal for at least 100 mS for the selection to become active. Use the IPC touchscreen Settings screen to change the default value.

To use the multifunction foot pedal (Figure 12-8) to control the handpiece do the following:

- To select Forward or Reverse mode, press the mode button.
- To start or adjust the speed of a handpiece in variable mode, press the foot pedal.
- To toggle between the start/stop mode and variable speed mode, press the control button.

Note: Functionality of the control button may be changed in the handpiece Defaults menu to pause irrigation. Refer to Change System Settings in the Pre-Operating Instructions section.

- To change the handpiece, press the handpiece button.

Figure 12-8. Multifunction Foot Pedal & Y-Splitter

- | | | | |
|---|-------------------------|---|------------|
| 1 | Mode button | 5 | Foot pedal |
| 2 | Handpiece button | 6 | Y-Splitter |
| 3 | Control button | 7 | Port 1 |
| 4 | Slip-resistant foot pad | 8 | Port 2 |

INDIGO DRILL CLEANING AND STERILIZATION INSTRUCTIONS

Refer to documents 68E4187 and 68E4188 in the Cleaning and Sterilization section.

INDIGO DRILL TECHNICAL SPECIFICATIONS

Indigo Drill 1845000

Size	11.9 cm Length x 1.53 cm Diameter
Weight	102 g
Speed	60000 rpm forward/reverse
Duty Cycle for Applied Part	For continuous use in operating room temperatures up to 40°C, the Indigo motor is rated for 3 minutes at 60000 rpm, followed by 25 minutes of rest.
	For normal operating room temperatures (typically 20°C), the Indigo motor is rated for continuous cutting indefinitely at 60000 rpm.

MIDAS REX MICROSAWS

DEVICE DESCRIPTION

Use Midas Rex Microsaws to remove hard tissue and bone during surgical procedures. Both the foot and finger-controlled saws are powered by the IPC. Operate the foot-controlled saws via the multifunction foot pedal. Operate the finger-controlled saws by pressing the finger lever.


The IPC incorporates microsaws irrigation at pump 2. If you do not use irrigation, select **None** for pump 2. Control operation of the microsaws with the IPC touchscreen and the multifunction foot pedal.

Figure 13-1. Foot-Controlled Microsaw Handpieces

- 1 Reciprocating saw collet
- 2 Sagittal saw blade release button
- 3 Oscillating saw head
- 4 Oscillating saw collet
- 5 Color band

Figure 13-2. Finger-Controlled Microsaw (Microsaw Touch) Handpieces


- 1 Reciprocating saw collet
- 2 Sagittal saw blade release button
- 3 Oscillating saw head
- 4 Oscillating saw collet
- 5 Color band
- 6 Finger lever safe mode switch
- 7 Finger lever

BLADE AND RASP NOMENCLATURE

Important: You can only use reciprocating blades and rasps with the reciprocating saw. You can use sagittal and oscillating blades interchangeably.

Saw blade part numbers follow a standard naming convention.

1	Cut Depth Limit	XX.X mm*			
2	Saw Type	S	Sagittal / Oscillating		
		R	Reciprocating		
3	Cut Width / Length	X.X mm			
4	Blade Thickness	X.X mm			
5	Other Identifiers	F	Fine Teeth	PT	Pointed Blade
		ARC	Rounded Blade	ST	Straight Rasp
		BAY	Bayoneted Blade	CC	Crosscut Rasp


* Blade depth markings are approximate and should be used for reference only.

RECIPROCATING MICROSAW ASSEMBLY

Important: Verify that the color code on the blade package matches the color band on the handpiece (Figures 13-1 and 13-2).

1. On the reciprocating handpiece (red band), verify that the collet is in the unlocked position (Figure 13-3).
2. Insert the reciprocating blade into the collet.
3. Rotate the collet knob 90 ° to the locked position.
4. To ensure a secure fit, gently pull the blade.

Figure 13-3. Reciprocating Saw Assembly

SAGITTAL MICROSAW ASSEMBLY

Important: Verify that the color code on the blade package matches the color band on the handpiece (Figures 13-1 and 13-2).

- 1. On the sagittal handpiece (yellow band), press the blade release button (Figure 13-4).
- 2. Insert the blade at the desired angle (Figure 13-5).
- 3. Release the blade release button.
- 4. To ensure a secure fit, gently move the blade from side to side.

Figure 13-4. Sagittal Saw Assembly


Figure 13-5. Blade Installation Angles

OSCILLATING MICROSAW ASSEMBLY

Important: Verify that the color code on the blade package matches the color band on the handpiece (Figures 13-1 and 13-2).

- 1. On the oscillating handpiece (blue band), verify that the collet is in the unlocked position (Figure 13-6).
- 2. Insert the blade into the collet. The blade "clicks" into place.
Note: You can insert the blade every 45 ° in a 360 ° circle in up to 8 locations (Figure 13-7).
- 3. Rotate the collet to the locked position (Figure 13-6).
- 4. To ensure a secure fit, gently move the blade from side to side.

Figure 13-6. Oscillating Saw Assembly



Figure 13-7. Oscillating Pin Alignment

MICROSAW IRRIGATION AND TUBING CLIPS ASSEMBLY


Note: Saw irrigation is optional.

1. Snap the irrigation clip on to the handpiece (Figures 13-8 and 13-9).
Note: On finger-controlled handpieces, attach the rear portion of the irrigation clip to the rotatable finger lever base (Figure 13-9).
2. Disengage the front of the irrigation clip to adjust the irrigation length, if necessary.
3. Snap the tubing clip on to the handpiece.

Figure 13-8. Foot-Controlled Microsaw Irrigation and Tubing Clip Assembly

Figure 13-9. Finger-Controlled Microsaw Irrigation and Tubing Clip Assembly

CONNECT MICROSAW HANDPIECE TO IPC

On the IPC Console, locate a 12-pin port on the connector panel (Figure 13-10), align the mark on the connector to the mark on the console, then insert the connector.

Figure 13-10. Microsaws IPC Connection Ports

- 1 Microsaws handpiece connection port
- 2 Multifunction foot pedal connection port

MICROSAWS TOUCHSCREEN CONTROLS

Important: The IPC recognizes the connected saw, then displays the appropriate name and color ring.

- To adjust speed, in the Speed control box (Figure 13-11), press the plus button to increase maximum speed and the minus button to decrease maximum speed. Variable adjustment ranges from 0 % to 100 %.
Note: When the Finger-Controlled Microsaw is in safe mode, the Speed control box displays SAFE (Figure 13-12). Refer to the Finger-Controlled Microsaw Safe Mode topic for additional information.
- To adjust the irrigation flow rate, in the Pump control box, press the plus button to increase flow rate or the minus button to decrease flow rate. If intermittent flow is available, pressing the plus or minus button progresses the system through intermittent and continuous flow. The system displays **Intermittent** when in intermittent mode.
Default flow rate is 0 cc per minute. Maximum flow rate is 50 cc per minute.
Note: To adjust flow rate, you can use the touchscreen or the Intelliflow Irrigation remote control.
- To switch between foot and finger control or to use both, in the Control control box, select the appropriate option. Finger-controlled saws can be operated via the lever on the handpiece, or the foot pedal.
Note: This option is only available when using a finger-controlled handpiece.

FINGER-CONTROLLED MICROSAW SAFE MODE

When the handpiece is in safe mode, it is inoperable until the safety is turned off. The Speed control box on the Microsaws touchscreen displays SAFE when the handpiece is the active handpiece and in safe mode (Figure 13-12).

Switch the device to safe mode any time it is attached to the console, but not currently being used.

When more than one handpiece is attached to the console, use the safety switch of an inactive handpiece to activate that handpiece and make it ready for use.

Figure 13-11. Microsaws Touchscreen

Figure 13-12. Finger-Controlled Microsaw Safe Mode


MICROSAW MULTIFUNCTION FOOT PEDAL CONTROLS

Important: By default, press each button on the foot pedal for at least 100 mS for the selection to become active. Use the IPC Settings screen to change the default value.

To use the multifunction foot pedal (Figure 13-13) to control the handpiece, do the following:

- To start or adjust the speed of a handpiece in variable mode, press the pedal.
- To toggle between the start/stop mode and variable-speed mode, press the control button.
Note: Functionality of the control button may be changed in the handpiece Defaults menu to pause irrigation. Refer to Change System Settings in the Pre-Operating Instructions section.
- To change the active handpiece, press the handpiece button.
Note: You can switch between finger-controlled handpieces by pressing the finger lever on the handpiece.

Figure 13-13. Multifunction Foot Pedal & Y-Splitter

- | | |
|---------------------------|--------------|
| 1 Mode button | 5 Foot pedal |
| 2 Handpiece button | 6 Y-Splitter |
| 3 Control button | 7 Port 1 |
| 4 Slip-resistant foot pad | 8 Port 2 |

MICROSAWS CLEANING AND STERILIZATION INSTRUCTIONS

Refer to document M000030A231 in the Cleaning and Sterilization section.

MICROSAWS TECHNICAL SPECIFICATIONS

Size	Reciprocating (ES200): 2.3 cm Width x 20.3 cm Length Sagittal (ES300): 2.3 cm Width x 17.3 cm Length Oscillating (ES100): 2.3 cm Width x 19.8 cm Length Finger-Controlled Reciprocating (ES210): 2.3 cm Width x 3.2 cm Height x 20.3 cm Length Finger-Controlled Sagittal (ES310): 2.3 cm Width x 3.2 cm Height x 17.8 cm Length Finger-Controlled Oscillating (ES110): 2.3 cm Width x 3.2 cm Height x 21.0 cm Length
Weight	Reciprocating (ES200): 300 g Sagittal (ES300): 250 g Oscillating (ES100): 300 g Finger-Controlled Reciprocating (ES210): 320 g Finger-Controlled Sagittal (ES310): 300 g Finger-Controlled Oscillating (ES110): 320 g
Speed	All Reciprocating: 1400-14000 cpm All Sagittal: 2000-20000 cpm All Oscillating: 1600-16000 cpm

Duty Cycle

≤ 20°C Ambient (one cut period = 20 seconds on / 20 seconds off)

Reciprocating: ≤ 15 cut periods¹

Sagittal: ≤ 6 cut periods²

Oscillating: ≤ 4 cut periods³

> 20°C Ambient (one cut period = 20 seconds on / 40 seconds off)

Reciprocating: ≤ 5 cut periods⁴

Sagittal: ≤ 6 cut periods⁵

Oscillating: ≤ 3 cut periods⁶

Time for cool down after maximum number of cut periods listed above: 25 minutes.

¹5 periods with Light Force (0.8 Lbf), 5 periods with Medium Force (2.1 Lbf), 5 periods with Heavy Force (3.9 Lbf)

²2 periods with Light Force (0.8 Lbf), 2 periods with Medium Force (1.5 Lbf), 2 periods with Heavy Force (3.0 Lbf)

³1 period with Light Force (0.7 Lbf), 2 periods with Medium Force (1.5 Lbf), 1 period with Heavy Force (2.6 Lbf)

⁴1 period with Light Force (0.8 Lbf), 2 periods with Medium Force (2.1 Lbf), 2 periods with Heavy Force (3.9 Lbf)

⁵2 periods with Light Force (0.8 Lbf), 2 periods with Medium Force (1.5 Lbf), 2 periods with Heavy Force (3.0 Lbf)

⁶1 period with Light Force (0.7 Lbf), 1 period with Medium Force (1.5 Lbf), 1 period with Heavy Force (2.6 Lbf)

POWEEASE DRIVER

DEVICE DESCRIPTION

The IPC POWEEASE System is specifically designed for drilling, tapping, and driving screws during spine surgery.

The IPC POWEEASE system consists of a POWEEASE Driver (Figure 14-1) that is powered by the IPC, along with the following POWEEASE driver compatible working end instruments:

- Taps
- Drill Bits
- Drivers
- Rod Cutter
- Post Cutter
- Set Screw Break-off Instrument (SSBO)

The POWEEASE System and accessories are compatible with the following implant systems:

- CD HORIZON SOLERA Spinal System
- TSRH 3Dx Spinal System

ASSEMBLE THE IPC

Important: The manual start/stop button on the IPC console (Figure 1-2) does not function with the POWEEASE driver.

Note: The POWEEASE driver does not use the IPC pump system. The pumps are available for use with other accessories requiring irrigation or cooling.

1. Verify the wheels are locked on the IPC cart.
2. Inspect all components for damage and determine if the system is ready for use.
3. Mount the IPC and irrigation/collant bags on the I.V. pole.
Important: Mount the irrigant and coolant bags above the IPC to ensure adequate flow.
4. Plug the IPC into the power source. Position the IPC so that it does not obstruct the power source for the purpose of disconnecting the Main voltage by the power cord.
5. Locate the POWEEASE connection port on the IPC connector panel (Figure 14-2) and insert the connector from the POWEEASE driver.
6. If using the NIM-ECLIPSE with the IPC POWEEASE, assemble the IPC POWEEASE and the NIM-ECLIPSE system.
7. Turn on the IPC and do the following:
 - Verify the system passes the self-test.
 - Verify the POWEEASE screen appears on the IPC monitor.
8. Confirm the IPC POWEEASE System Operation.

IPC POWEEASE MODE SELECT SWITCH

Use the mode select switch to change the handpiece from Forward, Reverse or Oscillate when the handpiece is the active handpiece. When the handpiece is the inactive handpiece, use the mode select switch to activate the handpiece.

IPC POWEEASE SYSTEM OPERATION

1. On the POWEEASE driver, rotate the mode select switch (Figure 14-3) to Forward (F).
2. On the IPC touchscreen, verify the following:
 - The FWD Speed box is active and the default speed appears in white.
 - The Trigger Status is Ready For Use.

Note: If the trigger status is Locked, release and press the trigger to disengage the lock. Contact Customer Service if the trigger status does not change to Ready For Use.
3. On the POWEEASE driver, touch the trigger (Figure 14-4) and verify the speed in the FWD Speed box on the IPC touchscreen turns green.
4. On the POWEEASE driver, press the trigger and verify the following:
 - The speed increases as you increase pressure on the trigger.
 - The speed in the FWD Speed box on the IPC touchscreen turns yellow.
5. Remove finger from the POWEEASE driver trigger and verify the following:
 - The driver rotation stops.
 - The speed in the FWD Speed box on the IPC touchscreen turns white.
6. Repeat steps 1 through 5 for the Reverse (R) and Oscillate  modes.

Figure 14-1. POWEEASE Driver


Figure 14-2. IPC POWEEASE Connection Ports




Figure 14-3. POWEEASE Mode Select Switch

POWEEASE TOUCHSCREEN CONTROLS

To adjust the maximum speed of the driver using the IPC POWEEASE touchscreen (Figure 14-4), in the relevant speed control box, press the plus button to increase maximum speed or minus button to decrease maximum speed.

Figure 14-4. IPC POWEEASE Touchscreen


- Variable speed control:
FWD Speed: variable adjustment from 10 to 250 rpm in increments of 10 rpm.
REV Speed: variable adjustment from 10 to 250 rpm in increments of 10 rpm.
OSC Speed: variable adjustment from 10 to 200 rpm in increments of 10 rpm.
- IPC POWEEASE touchscreen color codes:
White: active speed at idle, finger off trigger.
Green: active speed at idle, finger on trigger.
Yellow: active speed in use.
Gray: inactive
- Trigger Status:
Ready For Use: Device is ready for use.
Locked: In the case of a trigger control anomaly, the trigger is locked until the anomaly has cleared and use of the POWEEASE may resume.
To disengage the lock, release and press the trigger. Contact Customer Service if the device does not return to Ready For Use.

CLASSIC WORKING END ASSEMBLY AND OPERATION

NOTE: Please see the POWEEASE System Working Ends User's Manual and IFU for a listing of compatible working ends and further details regarding their proper use. Contact Customer Service or your sales representative for the most up-to-date version of the package insert or manual.


Figure 14-5. Classic Working End

1 Four sided drive shaft
2 Flat Side
3 Drive shaft shoulder
4 Shaft

1. On the POWEEASE driver, pull back on and hold the quick disconnect to unlock the collet (Figure 14-1, 1).
Note: The quick disconnect must be held in the unlocked position when inserting or removing tools.
2. Insert the classic working end drive shaft (Figure 14-5, 1) into the collet until fully seated. Fully seated is on or near the drive shaft shoulder (Figure 14-6).
Note: If using a two sided shaft, align the flat sides on the drive shaft with the marks on the collet (Figure 14-1, 8).
3. Release the quick disconnect to place the working end in the locked position (Figure 14-1, 2).
4. To set the mode on the POWEEASE driver, rotate the mode select switch (Figure 14-3) to Forward (F), Reverse (R), or Oscillate 
 - Place mode select switch in Forward to drill, tap or insert screw.
 - Place mode select switch in Reverse to back out drill, back out tap or back out screw.
 - Place mode select switch in Oscillate to lock the tool in its current position. Refer to POWEEASE Driver Electronic Ratcheting for additional information.
5. To adjust the speed of the driver, increase pressure to increase speed or decrease pressure to decrease speed.
NOTE: The IPC touchscreen displays the maximum speed for the selected mode.
6. Reverse procedure to remove the classic working end accessory.

Figure 14-6. Classic Working End Seating

POWEEASE Driver Electronic Ratcheting

Ratcheting is intended to permit manual motion in one direction only.

- FWD Mode: With the trigger released, the handpiece functions as a ratchet allowing the handpiece to be turned in a clockwise direction while ratcheting in a counter clockwise direction.
- REV Mode: With the trigger released, the handpiece functions as a ratchet allowing the handpiece to be turned in a counter clockwise direction while ratcheting in a clockwise direction.
- OSC Mode: With trigger released, the handpiece locks the tool in its current position, and the surgeon can manually rotate the handpiece in a clockwise or counter-clockwise direction.

POWEEASE Ratchet Test


1. Load a drill, tap or screwdriver in the collet. Refer to Classic Working End Assembly and Operation for details.
2. On the POWEEASE driver, set the mode select switch to Forward (F), (Figure 14-3).
3. Hold the POWEEASE driver shaft firmly and verify the following:
Important: DO NOT hold by the tool. **DO NOT** touch the trigger sensor.
 - The driver ratchets when rotated counter clockwise.
 - The driver does not ratchet when rotated clockwise.
4. Set the mode select switch to Reverse (R) and verify the following:
 - The driver ratchets when rotated clockwise.
 - The driver does not ratchet when rotated counter clockwise.
5. Set the mode select switch to Oscillate  and verify the driver does not ratchet when rotated clockwise or counter clockwise.

MECHANIZED WORKING END ASSEMBLY AND OPERATION

NOTE: Please see the POWEREASE System Working Ends User’s Manual and IFU for a complete listing of compatible working ends and further details regarding their proper use. Contact Customer Service or your sales representative for the most up-to-date version of the package insert or manual.

1. On the POWEREASE driver, pull back on and hold the quick disconnect (Figure 14-1, 1) to unlock the collet.
Note: The quick disconnect must be held in the unlocked position when inserting or removing tools.
2. Align the slots on the working end collar with the mechanized working collar locking pins (Figure 14-7).
3. To set the mode on the POWEREASE driver, rotate the mode select switch (Figure 14-3) to Forward (F), Reverse (R), or Oscillate 
 - Place the mode select switch in Forward to cut/shear implantable rods, cut vertical posts, break set screws.
 - Place the mode select switch in Reverse to open the tool.
 - Oscillate is intended as a rotational lock for classic working ends and has no function with mechanized working ends. However, if the trigger is pulled, the tool shaft will oscillate 60 ° clockwise and counter-clockwise.
4. To adjust the speed of the driver, increase pressure to increase speed or decrease pressure to decrease speed.
NOTE: The IPC touchscreen displays the maximum speed for the selected mode.
5. Reverse this procedure to remove the classic working end accessory.

Figure 14-7. Mechanical Working End Seating

POWEREASE CLEANING AND STERILIZATION INSTRUCTIONS

Refer to document 68E4189 in the Cleaning and Sterilization section.

POWEREASE DRIVER TECHNICAL SPECIFICATIONS

POWEREASE Driver 2300000

Size	175 mm Width x 160 mm Height x 42 mm Length
Weight	900 g (1200 g including cable)
Speed	0 to 250 rpm (variable)
Torque	7 Nm
Cable Length (Driver)	4.5 m
Cable Length (IPC to Patient Interface Box)	4.5 m
Cable Length (Driver to IPC)	4.5 m
Duty Cycle for Applied Part	For continuous operating room temperatures up to 33 °C, the POWEREASE driver rates for continuous operation. For normal operating room temperatures (below 25 °C), the POWEREASE driver rates for continuous operation.

IPC POWEREASE WITH THE NIM-ECLIPSE SYSTEM

Important: Do not use the IPC POWEREASE - NIM-ECLIPSE System with multiple portable socket-outlet or extension cord.

The POWEREASE System is equipped with two cables that connect the IPC to the NIM-ECLIPSE. Wires within the POWEREASE make contact with the uncoated tool and carry stimulating current to the tool’s tip for nerve monitoring. Contact to your local NIM-ECLIPSE representative for more information.

Important: Only taps, screws and drills are available for nerve stimulation. Please see the POWEREASE System Working Ends User’s Manual and Instructions for Use for further details regarding the working ends available for use with nerve stimulation. When the IPC POWEREASE driver is used in combination with the NIM-ECLIPSE, the combination creates a medical electrical system. The system components are as follows:

- IPC
- POWEREASE with accessories
- NIM-ECLIPSE
- NIM cable set


Assemble IPC POWEREASE and NIM-ECLIPSE System

Important: Refer to the NIM-ECLIPSE User's Guide for set up and use of the NIM-ECLIPSE system, including the laptop computer.

Warnings:

- Only the POWEREASE handpiece with its accessories are intended to be used within the patient environment. The NIM-ECLIPSE cannot be used in the patient vicinity (1.5 meter radius from patient).
 - Do not touch the patient and the enclosure parts of the NIM-ECLIPSE simultaneously.
 - Dispose of the single-use NIM cables after each procedure.
 - Any peripheral devices connected to the system must meet IEC60601-1 requirements. Consult with your biomedical engineering department to determine if the external devices meet this requirement.
1. Assemble the IPC and confirm IPC POWEREASE system operation.
 2. Assemble the NIM-ECLIPSE according to the instructions in the NIM-ECLIPSE User's Guide. Verify the laptop computer power cord is plugged into the NIM-ECLIPSE 945ECLC Controller.
 3. On the NIM-ECLIPSE stimulator extender, connect one end of the black NIM cable to the Stimulus Output connection port (Figure 14-8).
 4. On the IPC connector panel, connect the other end of the black NIM cable to the Stimulus Input connection port (Figure 14-8).
Important: Route the wire so that it does not interfere with operating room personnel.
 5. Connect one end of the blue NIM cable to the Stimulus Output IPC connector panel (Figure 14-8).
 6. Connect the other end of the blue NIM cable to the NIM Stim Input on the POWEREASE driver.
Important: Route the wire so that it does not interfere with operating room personnel.

Figure 14-8. POWEREASE, IPC, NIM-ECLIPSE Connection

NIM-ECLIPSE

- | | |
|--|------------------------------|
| 1 Main POWERASE Cable | 7 Mouse |
| 2 NIM Cable 2 (blue) | 8 Portable Notebook Computer |
| 3 NIM Cable (black) | 9 USB Interface Cable |
| 4 Stimulator Extender | 10 Stimulator Cable |
| 5 Digital Preamplifier Interface Cable | 11 Controller |
| 6 Preamplifier Module | 12 Notebook Power Supply |

TROUBLESHOOTING AND ERROR CODES

For any troubleshooting items not corrected by the actions below, contact Customer Service.

General System Troubleshooting

Issue	Possible Cause	Action
Pump(s) does not run.	Failed internal components.	Contact Customer Service
Little or no irrigation flow.	Tubing Set improperly seated in pump.	Reposition tubing in pump, verify pump lid is fully closed with the fluid flow from left to right.
	Tubing is pinched or kinked.	Check tubing at side of pump, see Irrigation/Coolant Pumps. Check remaining tubing for pinched or kinked areas, if necessary replace tubing.
	Tubing clamps are restricting flow.	Set tubing clamps in "open" position.
	Irrigation flow rate setting low.	Adjust irrigation flow rate.
	Irrigator obstructed.	Replace irrigator.
Pump stall error.	Tubing set improperly placed in pump.	Reposition tubing in pump, verify pump lid is fully closed with the fluid flow from left to right.
	Tubing is pinched or kinked.	Check tubing is not pinched or kinked on side of pump (see section on "Irrigation/Coolant Pumps").
Handpiece connected but console reads "Connect Handpiece."	Moisture in cable conflicts with handpiece recognition.	Run a dry cycle when sterilizing.
Handpiece connected but console displays incorrect handpiece.		
Console does not power up.	Power cord not properly connected.	Connect power cord.
	No power.	Check power available (i.e. power strip is on, circuit breaker is closed, etc.).
	Power Inlet Fuses blown.	Replace fuses with 5.00 A, 250 V, time delayed fuses (P/N 11270066 Fuse Kit 1898125).
	Failed internal components.	Contact Customer Service.
Power switch light is on but Touchscreen does not come on.	Failed internal components.	Contact Customer Service.
Console does not power down.	Power switch failure.	Unplug power cord, Contact Customer Service.
Touchscreen does not respond.	Screen gasket displaced or failed internal components.	Contact Customer Service.
Touchscreen does not work properly.	Touchscreen not calibrated.	Calibrate Touchscreen.
Console displays wrong handpiece / motor type.	Console misidentified the handpiece / motor.	Disconnect and reconnect the motor cable. Turn console off then on. Change motor, motor cable, or console to isolate the problem.
	Moisture in cable conflicts with handpiece recognition.	Run a dry cycle when sterilizing.
Foot control unit buttons or pedal does not respond.	Faulty Y-Splitter.	Disconnect Y-Splitter and connect FCU to IPC directly.
	Incorrect use.	Press and hold buttons for at least 1 second, wait for console confirmation beep.
	Top button does not respond.	One (1) handpiece connected (top button has no function with 1 handpiece connected).
	Connector not fully inserted.	Disconnect and reconnect the FCU cable connector. Try different FCU or console to isolate the problem.
	Internal component failure.	Contact Customer Service.
Handpiece fails to rotate.	Failed foot switch.	Disconnect foot switch, use manual start/stop rocker switch on rear of console.
	Failed handpiece motor or motor driver.	Contact Customer Service.
	Cables are not property connected.	Ensure handpiece cable is properly connected.
	Cables are damaged.	Check cables for cracks, splits, or bent connector pins.

Issue	Possible Cause	Action
Handpiece is too hot to touch/hold.	Blade, bur, and/or attachment transferring heat to handpiece.	Remove and reinstall or change blade, bur, and/or attachment. Ensure the part is properly seated before resuming use.
	Inadequate irrigation and/or coolant.	Ensure adequate irrigation and/or coolant is present.

XOMED Blades or Burs

Issue	Possible Cause	Action
Appears to be damaged or defective.	Damaged or defective.	Remove and replace.
Tool vibrates excessively, abnormal noise movement.	Tool is not firmly seated.	Microdebriders, pull back locking collet and re-seat the tool.
		Visao, unlock collar, check/re-seat notch, lock collar.
No suction.	Blade opening is obstructed.	Use stylet to clear blade. Remove blade from surgical site and submerge the blade tip in sterile water with suction connected to the handpiece to evacuate the obstruction.
	Tubing obstructed.	Remove and inspect suction tubing, and if obstructed, remove obstruction, reconnect tubing.
Tool is leaking irrigant.	Tool not seated correctly in collet.	Check for proper tool insertion by pulling back locking collet, and re-seating tool.
	Low or no suction.	See Possible Cause, No Suction.
Tool wobbles in handpiece.	Tool wobbles in handpiece.	Reduce handpiece operating speeds.
		Use tools that are rated for the console speed selected.
		If necessary, use bur guard with burs medium, long and X-long.
		Operate handpiece at 50 % of full speed for medium, long and X-long burs.
Tool vibrates excessively, abnormal noise movement.	Tool is not firmly seated.	Select a new tool.
		Microdebriders, pull back locking collet and reseat the tool. If applicable for model, thumbwheel may be rotated slightly to check. Visao, unlock collar, check/reseat notch, lock collar.

Legend EHS, Legend EHS Stylus and Indigo Motors

Issue	Possible Cause	Action
Motor is too hot to touch/hold.	Inadequate cool down period following sterilization.	Motor must be allowed to cool down following steam sterilization.
	Attachment transferring heat to the motor.	Switch attachments to determine whether the heat is being generated by the motor or the attachment.
	Heavy side loading during dissection.	Discontinue use and rest the motor by using it intermittently or wrap the motor with a moist sterile towel.
	Inadequate irrigation.	Ensure adequate irrigation to surgical site during bone dissection.
Tool is difficult to remove from attachment.	Aging of attachment.	Contact Customer Service.
	Use of reprocessed tools.	
	Use of an unauthorized refurbisher.	
	Improper cleaning.	Clean the attachment thoroughly according to the instructions in this manual. Change tool.
Attachment will not seat properly on the motor.	Motor collet flats are not aligned.	Use the Legend motor wrench to rotate the flat closest to the motor case until its marker is aligned with the marker on the flat farthest away from the motor case.
Motor does not run.	Cables not properly connected.	Ensure motor and foot control cables are properly connected.
	Speed setting is too low.	Ensure that a speed greater than 10000 rpm (EHS) or 3000 rpm (Stylus) is selected.

Issue	Possible Cause	Action
Motor with attachment rotates, but an abnormal noise is heard.	Attachment not properly installed and locked onto the motor.	Remove and reinstall the attachment and dissecting tool to ensure proper installation.
	Internal failure of motor and/or console.	Change motor or console to isolate the problem.
	Foot control not properly functioning.	Check for obstruction under the foot pedal.
	Cables damaged.	Check cables for cracks, splits, or bent connector pins.
	Bearings are worn.	Change the attachment to isolate the location of the problem.
	Poor electrical connection.	Check all connections from electrical source to console. Ensure motor and foot control cables are properly connected.
	Internal failure of motor, console, or cable.	Change motor, console, or cable to isolate the failing component.
	Attachment not properly installed.	Remove and reinstall the attachment and dissecting tool.

Stylus Touch Motors

Issue	Possible Cause	Action
Motor does not run.	Finger switch not reaching maximum speed.	Check that the control lever ring is properly seated in one of the four possible positions.
	Finger switch not responding. Safety switch in safe mode.	Place switch in run mode.
	Finger control damaged.	Contact Customer Service.

Legend EHS, Legend EHS Stylus and Stylus Touch Attachments and Telescoping Tubes

Issue	Possible Cause	Action
Attachment or Telescoping Tube has uncomfortable temperature to touch/hold.	Heat from worn attachment/tube bearings.	DO NOT use. Try another attachment/tube. Telescoping Tubes are multi-use disposable. If problem is resolved with a new Telescoping Tube, discard the over-heated tube.
	Attachment/tube unclean due to improper cleaning procedures.	Check that appropriate cleaning procedures are being followed.
	Heavy side loading during dissection.	Discontinue use and rest the attachment by using intermittently, try another identical attachment or wrap the attachment interface with a moist sterile towel.
Attachment/telescoping tube is bent, loose, damaged or missing a component.	Attachment mishandled, failed due to extended use or excessive force applied during use.	DO NOT use. Dispose of telescoping tube. Telescoping Tubes are multi-use disposable. Contact Customer Service.
Color band on Attachment/Telescoping Tube fades or discolors.	Incorrect cleaning or sterilization method.	Use nomenclature markings on the attachment to match with a corresponding dissecting tool or Contact Customer Service.
	Use of chlorine based or corrosive agents.	Telescoping Tubes are multi-use disposable.
	Aging.	
Attachment has excess lubrication.	Over lubrication during cleaning process.	Visually inspect and wipe excess lubrication.
Footed attachment has a component missing from leg/foot area or foot is bent.	Attachment damaged by dissecting tool drilling out part or all of leg/foot area.	DO NOT use. Contact Customer Service.
	Bend caused by incorrect use.	
16-MF contra-angle attachment is overheating.	The contra-angle attachment operates by a set of internal gears to engage the drive shaft. It is normal for some heat to be generated approximately 2 cm from the distal end of the attachment and at the right of the angle head.	If heat continues or is excessive, Contact Customer Service.

Issue	Possible Cause	Action
Smoke is generated by the attachment or motor.	Attachment is not in the locked position.	Make sure the attachment is in the locked position.

Legend EHS, Legend EHS Stylus and Stylus Touch Tools

Issue	Possible Cause	Action
Tool wobbles.	A non-Legend tool is being used.	Replace with a Legend tool.
	Worn attachment or tube bearings.	Try another attachment or tube to isolate the location of the problem. If the attachment is failing, Contact Customer Service. If the tube is failing, dispose of it and use a new tube.
	Attachment/tube and tool are not compatible.	Match color code on the tool packaging to the color code on the attachment/tube.
	Motor is damaged.	Contact Customer Service.
	Tool's size and geometry may contribute to wobbling at certain speeds.	Adjust the speed by changing the pressure setting or foot/finger control. Do not use if wobbling persists. Change tool.
Tool vibrates excessively.	Tool's size and geometry may create excessive vibration at certain speeds.	Adjust the speed.
		Change tools.
Tool dull.	Extended use.	Change to a new tool.
	Reprocessed tool was used.	
	Incorrect geometry.	Contact Customer Service.
Tool will not seat properly in the motor or attachment collet.	Debris in collet of attachment or motor.	Clean the attachment or motor thoroughly according to the instructions in this manual. If cleaning does not correct the problem, Contact Customer Service.
	A non-Legend tool is being used.	Replace with a Legend tool.

ERROR CODES

Code #	Title	Cause	Description
1	MCB does not report that it booted within 5 seconds of AI telling it to start and subsequent reattempts fail.	System Error.	Power off. Wait 10 seconds. Power on. If error persists, call Customer Service.
2	NOT USED.	NOT USED.	NOT USED
3	UI-MCB Com Failure - Max resends exceeded.	System Error.	Power off. Wait 10 seconds. Power on. If error persists, call Customer Service.
4	UI-MCB Com Failure - Get answer failed.		
5	UI-MCB Com Failure - No status message received.		
6	UI-MCB Com Failure - Serialization ID error.		
7	UI-MCB Com Failure - Timeout exception.		
8	UI-MCB Com Failure - Variable not recognized.		
9	Pump 1 stalled (no transitions on opto sensor).	Pump #1 stalled.	Check tubing connection.
10	Pump 2 stalled (no transitions on opto sensor).	Pump #2 stalled.	Check tubing connection.
11	Unrecognized/damaged handpiece plugged in on port 1 (first 12 pin).	Handpiece.	Unplug handpiece and plug back in. If error persists, replace handpiece.
12	Unrecognized/damaged handpiece plugged in on port 2 (second 12 pin).		
13	Unrecognized/damaged handpiece plugged in on port 3 (4 pin).		
14	Unrecognized/damaged handpiece plugged in on port 4 (Skeeter).		
15	Handpiece stalled.	Handpiece stalled.	Check accessory.
16	MCB motor overcurrent detected.	Motor overcurrent.	Unplug handpiece and plug back in. If error persists, replace handpiece.
17	Unrecognized/damaged multifunction foot pedal plugged in.	Foot Pedal Connection error.	When using foot pedal only, unplug foot pedal and plug back in. If error persists, replace foot pedal or switch to manual control. When using Y-Splitter, disconnect foot pedal and Y-Splitter, then connect foot pedal directly to IPC. If after reconnecting foot pedal, the error goes away, replace Y-Splitter with another unit. If error persists, replace foot pedal and reconnect both Y-Splitter and foot pedal or switch to manual control.
18	Damaged handpiece or finger lever base out of position.	Finger Control error.	Stylus Touch - A finger control error has been detected. Check that the control lever ring is properly seated in one of the four possible positions. If error persists contact Customer Service. Press OK to use alternate control method. Triton/POWERAISE - Please check mode select switch to ensure a mode has been selected by rotating mechanism until it rests in a detent. If error persists, contact Customer Service.
19	UI self test failure - culture (language) registry entry.	Self Test Failed.	Power off. Wait 10 seconds. Power on. If error persists, call Customer Service.
20	UI self test failure - sector configuration registry entry.		
21	UI self test failure - corrupt usage data file or unable to create usage data file.		
22	NOT USED.	NOT USED.	NOT USED.

IPC TROUBLESHOOTING AND ERROR CODES

23	MCB non-specific self test failure.	Self Test Failed.	Power off. Wait 10 seconds. Power on. If error persists, call Customer Service.
24	MCB self test failure - port 1.		
25	MCB self test failure - port 2.		
26	MCB self test failure - port 3.		
27	MCB self test failure - port 4.		
28	MCB self test failure - bridge transitor 1 shorted.		
29	MCB self test failure - bridge transitor 2 shorted.		
30	MCB self test failure - bridge transitor 3 shorted.		
31	MCB self test failure - bridge transitor 4 shorted.		
32	MCB self test failure - bridge transitor 5 shorted.		
33	MCB self test failure - bridge transitor 6 shorted.		
34	MCB self test failure - A/D converter.		
35	MCB self test failure - motor error.		
36	MCB self test failure - 3.3 volt supply.		
37	MCB self test failure - 12 volt supply.		
38	MCB self test failure - 48 volt supply.		
39	MCB self test failure - FCU port.		

CLEANING AND STERILIZATION

Reprocessing Instructions are subject to change. Refer to manuals.medtronic.com for current Reprocessing Instructions.

POST-OPERATIVE INSTRUCTIONS

Disconnect Accessory Cable from Console

To disconnect non-silicone multi-pin cables from the console, push the cable toward the console and then pull out by the lock ring.

Note: Silicone insulated multi-pin and single pin cable connectors do not have a lock ring (1). Remove these types of cable connectors straight from the connector panel.

Warning: After disconnecting insulated connectors (see W79) from the console, connectors that have debris under the insulator must be

cleaned according to Cleaning and Sterilization instructions. If debris is still present after cleaning and sterilization, return for warranty servicing.

Clean the Multifunction Foot Pedal

Important: If debris is present under the foot pedal's boot, return for warranty service.

DO NOT immerse or sterilize the foot pedal.

DO NOT use alcohol, other solvents or abrasive cleaners.

1. On the slip resistant foot pad ONLY, spray a neutral enzymatic detergent, pH 6.0-8.0, or a phenol based disinfectant, mixed according to manufacturer's instructions.
2. Leave the solution on the foot pad for approximately 10 minutes.
3. Dampen a cloth with a neutral enzymatic detergent, pH 6.0-8.0, or a phenol based disinfectant, mixed according to manufacturer's instructions.
4. Wipe the foot pedal with the damp cloth until visually clean.
5. Dry the unit with a clean, non-abrasive cloth.

Clean the Y-Splitter

1. DO NOT immerse or sterilize the Y-Splitter.
2. DO NOT use alcohol, other solvents or abrasive cleaners.
3. Dampen a cloth with a neutral enzymatic detergent, pH 6.0 – 8.0, or a phenol based disinfectant, mixed according to manufacturer's instructions.
4. Wipe the Y-Splitter with a damp cloth until visually clean.
5. Dry the unit with a clean, non-abrasive cloth.

 Medtronic		Reprocessing Instructions Triton Electric High-Torque Handpiece		Reprocessing Instructions Cleaning and Sterilization instructions are subject to change. Up to date instructions are available online at manuals.medtronic.com . M000030A322 D	
Warnings and Precautions	- Do not soak/submerge devices. - Do not use ultrasound to clean Triton High-Torque Handpiece. - Do not use chlorine based or corrosive cleaning agents such as bleach, lye, acetone, sodium hydroxide, formic acid, or solutions containing glutaraldehyde. - Allow an adequate cooling period after steam sterilization. - The use of a washer-disinfector for cleaning may cause a pre-mature degradation in performance.				
Limitations	Verify functionality prior to re-use, by connecting the handpiece to the IPC and pressing the trigger. Verify that the drive shaft rotates.				
Point of Use	Disconnect device from the IPC. Remove and discard disposables.				
Containment and Transportation	It is recommended that devices are reprocessed as soon as is practical following use.				
Preparation for Decontamination	Disconnect device from the IPC. Remove and discard disposables.				
Cleaning: Automated (Do NOT use ultrasonic washer)	Handpiece and Attachments - Prior to placing devices in the automated washer, manually rinse under tap water, until no visible soil is noticed. - Move any movable parts back and forth to allow water to reach hard-to-rinse areas. - Transfer the devices into the washer for processing. Place the devices in the washer to facilitate drainage. - Verify that devices are visually clean after automated cleaning.				
	Recommended Washer Cycle	Phase	Recirculation Time	Water Temperature	Detergent Type and Concentration
		Pre-Wash	2 minutes	Cold tap water	Not applicable
		Wash	5 minutes	66 °C	Neutral enzymatic detergent, pH 6.0-8.0
		Rinse	1 minute	Hot tap water	Not applicable
Cleaning: Manual	Handpiece and Attachments 1. Rinse the entire device, including cable, thoroughly under running water. 2. Wipe all external surfaces of the device and cable with a cloth dampened with a neutral enzymatic detergent pH 6.0-8.0 (henceforth referred to as "cleaning solution"). 3. Scrub the device thoroughly with a soft bristled brush that has been dampened with the cleaning solution, paying close attention to rough surfaces, crevices, and difficult to reach areas. Move any moveable parts, to allow the solution to reach all areas. 4. Clean central cannulation shaft of the handpiece with a cleaning brush (P/N 120-028-1 or equivalent) dampened with the cleaning solution. 5. Insert the cleaning brush into the cannulation from the front and make at least 6 passes with the brush. If residual soil is seen on the cleaning brush in between strokes, rinse it off in the cleaning solution. 6. Repeat the step above from the back side of the handpiece. 7. Flush approximately 30 mL of cleaning solution from a syringe into the central cannulation shaft. 8. Flush approximately 30 mL of cleaning solution from a syringe into the trigger and Forward / Reverse switch, while moving moveable parts. 9. Repeat the flushing in steps 7 and 8, using water instead of cleaning solution. 10. Rinse the entire device under running water. Move any moveable parts, to allow the water to reach all areas. Also, allow water to thoroughly clean the central cannulation shaft. Rinse until there is no visible evidence of soil or debris. 11. Dry the entire device with a towel. Sagittal and Reciprocating Saw Attachments 1. Rinse the entire device thoroughly under running water. 2. Wipe all external surfaces of the device with a cloth dampened with a neutral enzymatic detergent pH 6.0-8.0 (henceforth referred to as "cleaning solution"). 3. Scrub the device thoroughly with a soft bristled brush that has been dampened with the cleaning solution, paying close attention to rough surfaces, crevices, and difficult to reach areas. Ensure to brush the blade opening, the gaps/details between and around the blade jaws. 4. Dip the distal end of the saw body (up to approximately 1.25 inches from the distal end) into cleaning solution and agitate it for 10-15 seconds. 5. Flush approximately 30 mL of cleaning solution from a syringe into the blade opening and the gaps and details between and around two jaws. Repeat the flushing using water, instead of cleaning solution. 6. Repeat the previous step one more time. 7. Rinse the device under running water. Allow water to thoroughly rinse the gaps and details between and around blade jaws. Rinse until there is no visible evidence of soil or debris. 8. Dry the device with a towel.				
Disinfection	Follow hospital procedures.				
Packaging	For sterilization, place devices in instrument tray. Devices may be unwrapped, or wrapped with up to two layers of 1-ply polypropylene wrap				
Sterilization (Temperatures are minimum required, times are minimum required)	Steam Sterilization				
	Cycle	Gravity	Pre-Vac	Pre-Vac (FR/WHO)	Pre-Vac (UK)
	Temperature	132 °C	132 °C	134 °C	134 °C
	Time	21 minutes	14 minutes	19 minutes	19 minutes
	Drying	12 minutes	10 minutes	30 minutes	22 minutes
	STERAD	Do not use low temperature hydrogen peroxide gas plasma sterilization due to lumen internal diameter and length restrictions.			
	STERIS	Not applicable.			
	100% EtO Sterilization Parameters				
	Preconditioning	51-59 °C, 70 ±5 % relative humidity, 30 min.			
	Sterilization	Temperature	51-59 °C		
		Relative Humidity	70 +/- 5 %		
		Ethylene oxide concentration	725 +/- 25 mg/L		
		Gas exposure time (full-cycle)	4 hours		
		Aeration	51-59 °C, 18 hours		
Maintenance, Inspection and Testing	- Inspect devices for any damage before and after each use. If damage is observed, do not use the device until it is repaired - Verify functionality prior to re-use.				
Storage	Store in a clean, dry area.				
Additional Information	None				

Note: The Pre-Vac (FR/WHO) and Pre-Vac (UK) sterilization cycles are not considered by the United States Food and Drug Administration (US FDA) to be standard sterilization cycles. Users should only use sterilizers and accessories (such as sterilization wraps, sterilization pouches, chemical indicators, biological indicators, and sterilization containers) that have been cleared by the US FDA for the selected sterilization cycle specifications (time and temperature).

Note: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing is actually performed using equipment, materials and personnel in the reprocessing facility to achieve the desired result. This normally requires validation and routine monitoring of the process.

Note: All validations performed per current AAMI TIR12, Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers.

Medtronic recommends incineration of devices that have directly contacted patients suspected or confirmed with Transmissible Spongiform Encephalopathy (TSE)/ CJD diagnosis.

 Reprocessing Instructions Endo-Scrub 2 Finger Switch		Reprocessing Instructions Cleaning and Sterilization instructions are subject to change. Up to date instructions are available online at manuals.medtronic.com . 68E4005 D
Warnings and Precautions	Disconnect the finger switch from the Endo-Scrub 2 pump before cleaning.	
Limitations	After cleaning and sterilization, verify functionality prior to re-use.	
Point of Use	- This product is provided non-sterile and must be cleaned and sterilized before the first use and any reuse. - After use, remove the finger switch from the sheath and disconnect the plug from the pump. - Thoroughly rinse with water following use.	
Containment and Transportation	It is recommended that devices are reprocessed as soon as is practical following use.	
Preparation for Decontamination	Promptly and thoroughly rinse with deionized water after each use.	
Cleaning: Automated (Do NOT use ultrasonic washer)	Not validated	
Cleaning: Manual	- Dip the finger switch housing in a diluted mixture of mild (pH 7.0 - 8.5) enzymatic detergent. (Follow detergent manufacturer's instructions for proper dilution.) - Thoroughly clean the housing with a soft instrument brush to remove any blood and tissue. - Rinse the housing thoroughly with tap water and wipe dry. - Note: If wiping the cord dry, be sure to hold the cord and not the housing to avoid stressing or breaking the electrical connections located inside the housing.	
Disinfection	Do not cold soak in glutaraldehyde, chlorine, or ammonium solutions, or dry heat sterilize, as damage to the instrument finish may occur.	
Packaging	- A standard, sterilization wrap may be used. In the US, an FDA approved surgical wrap must be used. Ensure that the pack is large enough to contain the instrument without stressing the seals. - In sets: Instruments may be loaded into dedicated instruments trays or general purpose sterilization trays. Ensure that cutting edges are protected. Wrap trays using appropriate method.	
Sterilization (Temperatures are minimum required, times are minimum required)	- Check the cleanliness and operation of the instrument. Clean again if debris is present and remove from use any damaged instrument. Close instruments with catches and racks on the first notch. Arrange the instruments in sterilization containers with perforations on the top and bottom, and on supports such as those used in microsurgery. Follow the appropriate cycle listed in the table below. - The sterilization parameters given below should be used for devices that are fully disassembled when disassembly is possible. Use basic aseptic technique during post-sterilization assembly to maintain the sterility of the instruments. - All steam cycles have been validated in the wrapped configuration and can be sterilized wrapped or unwrapped. These devices have only been validated for steam sterilization methods.	
	Steam Sterilization	
	Cycle	Gravity Pre-Vac Pre-vac (UK)
	Temperature	132 °C 132 °C 134 °C
	Time	10 minutes 4 minutes 3 minutes
	Drying	15-30 minutes, or until visibly dry.
	STERRAD	Not validated.
	100% EtO Sterilization Parameters	Not validated.
Maintenance, Inspection and Testing	- Inspect finger switch for any damage before and after each use. If damage is observed do not use the finger switch until it is repaired or replaced. - After cleaning and sterilization, verify functionality prior to re-use.	
Storage	Store in a clean, dry area.	
Additional Information	None	

Note: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing is actually performed using equipment, materials and personnel in the reprocessing facility to achieve the desired result. This normally requires validation and routine monitoring of the process.

Note: All validations performed per current AAMI TIR12, Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers.

Medtronic recommends incineration of devices that have directly contacted patients suspected or confirmed with Transmissible Spongiform Encephalopathy (TSE)/ CJD diagnosis.

Medtronic

Reprocessing Instructions Microdebriders

StraightShot M5, StraightShot M4, StraightShot Magnum II, StraightShot III or
Midas Rex Spine Shaver

Reprocessing Instructions
Cleaning and Sterilization instructions are subject to change. Up to date instructions are available online at **manuals.medtronic.com**.
68E3282 F

Warnings and Precautions	• Disconnect the power before cleaning. • Do not fully immerse, or ultrasonically clean, this instrument. • Do not cold soak sterilize this instrument in glutaraldehyde. This will void the warranty. • Do not use organic solvents to clean the bur chuck. • For Microdebrider handpiece cleaning, cover handpiece cable connector end with one of the following handpiece cable caps: • Handpiece Cleaning Cap, Universal - catalog number 3318520. • Handpiece Cable Cap, Small - catalog number 3318510. • Use 3318520 for StraightShot M4, Spine Shaver, StraightShot Magnum II, and StraightShot III. • Using 3318520 is optional for StraightShot M5. • Use Catalog number 3318510 for other Microdebriders. • After completion of the cleaning steps, remove Handpiece Cable Cap or other protective components installed prior to cleaning.																																																
Limitations	After cleaning and sterilization, verify functionality prior to re-use.																																																
Point of Use	• This product is provided non-sterile and must be cleaned and sterilized before the first use and any reuse. • To remove occasional residual buildup on handpiece cable connector, use a soft brush and isopropyl alcohol.																																																
Containment and Transportation	It is recommended that instruments are reprocessed as soon as is practical following use.																																																
Preparation for Decontamination	Remove the bur from the handpiece, otherwise disassembly is not required.																																																
Cleaning: Automated (Do NOT use ultrasonic washer)	• Remove instruments and equipment from any instrument trays. • Use alkaline (pH <10.5) or neutral pH detergent. • Do not immerse the handpiece. Auto Wash Pre-Cleaning Instructions - Neutral Detergent • Manually rinse the handpiece, cable, and tubing management clips with running tap water to remove all visible soil from the device. Using a soft bristled brush and a neutral detergent solution, gently and thoroughly clean the entire open lumen of the device. Actuate the collet while rinsing and turn the thumbwheel. Hold the handpiece with collet end down when rinsing and run stream of water through the suction port. • Place devices in the washer with the collet facing downward to facilitate drainage. <table><tr><th colspan="4">Auto wash Cycle - Neutral Enzymatic Detergent</th></tr><tr><th>Phase</th><th>Minimum Time</th><th>Temperature</th><th>Detergent</th></tr><tr><td>Pre-wash 1</td><td>02:00</td><td>Cold tap water</td><td>Not applicable</td></tr><tr><td>Enzyme wash</td><td>05:00</td><td>< 50°C</td><td>neutral pH enzymatic detergent, manufacturer's recommended concentration</td></tr><tr><td>Rinse 1</td><td>01:00</td><td>Hot tap water</td><td>Not applicable</td></tr><tr><td>Thermal rinse</td><td>01:00</td><td>82°C minimum</td><td>Not applicable</td></tr></table> Visually examine for cleanliness. Auto Wash Pre-Cleaning Instructions - Alkaline Detergent • Manually rinse the handpiece, cable, and tubing management clips with running tap water to remove all visible soil from the device. Using a soft bristled brush and a neutral detergent solution, gently and thoroughly brush the entire open lumen of the device. Actuate the collet while rinsing and turn the thumbwheel. Hold the handpiece with collet end down when rinsing and run stream of water through the suction port. • Place devices in the washer with the collet facing downward to facilitate drainage. <table><tr><th colspan="4">Auto wash Cycle - Alkaline Detergent</th></tr><tr><th>Phase</th><th>Minimum Time</th><th>Temperature</th><th>Detergent</th></tr><tr><td>Pre-wash 1</td><td>02:00</td><td>Cold tap water</td><td>Not applicable</td></tr><tr><td>Wash 1</td><td>05:00</td><td>66°C (set point)</td><td>alkaline detergent pH ≤ 10.5</td></tr><tr><td>Rinse 1</td><td>01:00</td><td>Hot tap water</td><td>Not applicable</td></tr><tr><td>Thermal rinse</td><td>01:00</td><td>82°C minimum</td><td>Not applicable</td></tr></table> Visually examine for cleanliness.	Auto wash Cycle - Neutral Enzymatic Detergent				Phase	Minimum Time	Temperature	Detergent	Pre-wash 1	02:00	Cold tap water	Not applicable	Enzyme wash	05:00	< 50°C	neutral pH enzymatic detergent, manufacturer's recommended concentration	Rinse 1	01:00	Hot tap water	Not applicable	Thermal rinse	01:00	82°C minimum	Not applicable	Auto wash Cycle - Alkaline Detergent				Phase	Minimum Time	Temperature	Detergent	Pre-wash 1	02:00	Cold tap water	Not applicable	Wash 1	05:00	66°C (set point)	alkaline detergent pH ≤ 10.5	Rinse 1	01:00	Hot tap water	Not applicable	Thermal rinse	01:00	82°C minimum	Not applicable
Auto wash Cycle - Neutral Enzymatic Detergent																																																	
Phase	Minimum Time	Temperature	Detergent																																														
Pre-wash 1	02:00	Cold tap water	Not applicable																																														
Enzyme wash	05:00	< 50°C	neutral pH enzymatic detergent, manufacturer's recommended concentration																																														
Rinse 1	01:00	Hot tap water	Not applicable																																														
Thermal rinse	01:00	82°C minimum	Not applicable																																														
Auto wash Cycle - Alkaline Detergent																																																	
Phase	Minimum Time	Temperature	Detergent																																														
Pre-wash 1	02:00	Cold tap water	Not applicable																																														
Wash 1	05:00	66°C (set point)	alkaline detergent pH ≤ 10.5																																														
Rinse 1	01:00	Hot tap water	Not applicable																																														
Thermal rinse	01:00	82°C minimum	Not applicable																																														
Cleaning: Manual	• Remove instruments and equipment from instrument tray. • Do not immerse the handpiece. • Wipe the handpiece and cable with neutral detergent applied to a clean, non-abrasive cloth. • Manually rinse the handpiece, cable, and tubing management clips with running tap water to remove all visible soil from the device. Hold the handpiece with collet end down when rinsing and run stream of water through the suction port. • Use an enzymatic detergent solution to loosen and remove collected tissues from the unit. • Using an enzymatic detergent solution and a soft cloth dampened with the detergent, wipe the motor and cable. Gently clean the handpiece, cable, and tubing management clips with a soft bristle brush wetted with the enzymatic detergent. Using a soft bristled brush, gently and thoroughly clean the entire open lumen of the device. Actuate the collet while cleaning and turn the thumbwheel. Use the brush to thoroughly clean the tubing management clips. • Manually rinse the handpiece, cable, and tubing management clips with running tap water to remove all visible soil from the device. • Hold the handpiece with collet end down when rinsing the lumen thoroughly.* * Additional Cleaning Instructions for XPS Straightshot M5/M4/Spine Shaver Microdebrider: • During the normal cleaning cycle, run a gentle stream of warm water into the collet (front end), and into the lock lever of the Straightshot M4/Spine Shaver handpiece. • While warm water is running into the collet, rotate the mechanism for several revolutions (rotate the wheel); and while water is running into the lock lever, actuate the lock lever several times (locking and unlocking). • Shake excess water from the handpiece. • PRECAUTION: Ensure the use of a very gentle stream of warm clean water during this additional cleaning step.																																																
Drying and Lubrication	• Dry the handpiece and cable with a lint-free towel. Make sure to dry off the electrical connection on the cable ends. • Apply a small amount of silicone spray into the front-end collet and outside of the handpiece. • Sterilize the handpiece immediately after cleaning.																																																

 Medtronic		Reprocessing Instructions Microdebriders StraightShot M5, StraightShot M4, StraightShot Magnum II, StraightShot II  Midas Rex Spine Shaver			Reprocessing Instructions Cleaning and Sterilization instructions are subject to change. Up to date instructions are available online at manuals.medtronic.com. 68E3282 F	
Disinfection	Do not cold soak in glutaraldehyde.					
Packaging	- A standard, sterilization wrap may be used. In the US, an FDA approved surgical wrap must be used. Ensure that the pack is large enough to contain the instrument without stressing the seals. - In sets: Instruments may be loaded into dedicated instrument trays. Wrap trays using appropriate method.					
Sterilization (Temperatures are minimum required, times are minimum required)	The sterilization parameters given below should be used for devices that are fully disassembled when disassembly is possible. Use basic aseptic technique during post-sterilization assembly to maintain the sterility of the instrument(s). All steam sterilization cycles have been validated in the wrapped configuration and instruments can be sterilized wrapped or unwrapped.					
	Steam Sterilization					
	Cycle	Gravity*	Gravity	Pre-Vac	Pre-Vac (FR/WHO)	Pre-vac (UK)
	Temperature	121°C	132°C	132°C	134°C	134°C
	Time	40 minutes	15 minutes	4 minutes	18 minutes	3 minutes
	Drying	25 minutes, or until visibly dry				
	100% EtO Sterilization Parameters (set points)					
	Preconditioning	54°C, 60% relative humidity, 60 minutes				
	Sterilization	Temperature			54°C	
		Relative Humidity			60%	
		Ethylene oxide concentration			725 mg/L	
		Gas exposure time (full-cycle)			120 minutes	
		Aeration			54°C, 8 hours	
Maintenance, Inspection and Testing	- Inspect components for any damage before and after each use. If damage is observed do not use the instrument until it is repaired. - After cleaning and sterilization, verify functionality prior to re-use.					
Storage	It is extremely important that the handpiece be rapidly and completely vacuum dried before storage to prevent corrosion and residue deposits in the bearing and motor.					
Additional Information	Increase temperatures higher than those stated when necessary to satisfy governmental or health care facility requirements so long as the temperature does not exceed 149° C. Heating above 149° C may damage the handpiece and will void the warranty.					

* This 40 minutes and 121°C gravity sterilization cycle is not considered by the United States Food and Drug Administration (USFDA) to be a standard sterilization cycle. Users should only use sterilizers and accessories (such as sterilization wraps, sterilization pouches, chemical indicators, biological indicators, and sterilization containers) that have been cleared by the USFDA for the selected sterilization cycle specifications (40 minutes and 121°C).

Note: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing is actually performed using equipment, materials and personnel in the reprocessing facility to achieve the desired result. This normally requires validation and routine monitoring of the process.

Note: All validations performed per current AAMI TIR12, Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers.

Medtronic recommends incineration of devices that have directly contacted patients suspected or confirmed with Transmissible Spongiform Encephalopathy (TSE)/ CJD diagnosis.

 Medtronic		Reprocessing Instructions Legend EHS and Legend EHS Stylus			Reprocessing Instructions Cleaning and Sterilization instructions are subject to change. Up to date instructions are available online at manuals.medtronic.com . M000030A234 C	
Warnings and Precautions		- Do not soak/submerge devices. - Do not use ultrasound to clean devices. - Do not use chlorine based or corrosive cleaning agents such as bleach, lye, acetone, sodium hypochlorite/bleach, sodium hydroxide, formic acid, or solutions containing glutaraldehyde. - The use of a washer-disinfector for cleaning may cause a degradation in performance. - Allow an adequate cooling period after steam sterilization.				
Limitations		Verify functionality prior to re-use.				
Point of Use		Follow hospital procedures.				
Containment and Transportation		It is recommended that devices are reprocessed as soon as is practical following use.				
Preparation for Decontamination		Follow hospital procedures.				
Cleaning: Automated (Do NOT use ultrasonic washer)		- Review the washer-disinfector warning above, before using this cleaning method. - Remove devices from instrument trays before placing into washer baskets. - Orient devices following recommendations of the washer/disinfector manufacturers.				
		Recommended Washer Cycle	Phase	Recirculation Time	Water Temperature	Detergent Type and Concentration
			Pre-Wash	2 minutes	Cold tap water	Not applicable
			Wash	5 minutes	66 °C (set point)	Neutral enzymatic detergent, pH 6.0-8.0
			Final Rinse	1 minute	Hot tap water	Not applicable
Cleaning: Manual		- Wipe all external surfaces of the motor and cable with a cloth dampened with a neutral enzymatic detergent, pH 6.0-8.0. - Brush motor case and collet with a nylon brush dampened with a neutral enzymatic detergent. - Rinse motor thoroughly under running water, collet end pointed down. Dry collet and motor with towel. - Verify that devices are visually clean after manual cleaning. - Rinse thoroughly with tap water.				
Disinfection		Follow hospital procedures.				
Packaging		For sterilization, place devices in instrument tray. Devices may be unwrapped, or wrapped with up to two layers of 1-ply polypropylene wrap.				
Sterilization (Temperatures are minimum required, times are minimum required)		Steam Sterilization				
		Cycle	Gravity	Pre-Vac	Pre-Vac	
		Temperature	132 °C	132 °C	134 °C	
		Time	25 minutes	4 minutes	3 minutes	
		Drying	15 minutes	15 minutes	10 minutes	
		STERAD	Do not use low temperature hydrogen peroxide gas plasma sterilization due to lumen internal diameter and length restrictions.			
		STERIS	Do not use liquid peracetic acid sterilization due to immersion procedure.			
		100% EtO Sterilization Parameters				
		Preconditioning	51-59 °C, 70 ±5 % relative humidity, 30 min.			
		Sterilization	Temperature	51-59 °C		
			Relative Humidity	70 +/- 5 %		
			Ethylene oxide concentration	725 +/- 25 mg/L		
			Gas exposure time (full-cycle)	4 hours		
			Aeration	51-59 °C, 18 hours		
Maintenance, Inspection and Testing		- Inspect devices for any damage before and after each use. If damage is observed, do not use the device until it is repaired. - Verify functionality prior to re-use.				
Storage		Store with other sterile devices.				
Additional Information		None				

Note: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing is actually performed using equipment, materials and personnel in the reprocessing facility to achieve the desired result. This normally requires validation and routine monitoring of the process.

Note: All validations performed per current AAMI TIR12, Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers.

Medtronic recommends incineration of devices that have directly contacted patients suspected or confirmed with Transmissible Spongiform Encephalopathy (TSE)/ CJD diagnosis.

**Medtronic**

Reprocessing Instructions

Stylus Touch

Reprocessing Instructions

Cleaning and Sterilization instructions are subject to change. Up to date instructions are available online at manuals.medtronic.com.
68E4132 D

Warnings and Precautions	- Do not soak/submerge devices. - Do not use ultrasound to clean devices. - Do not use chlorine based or corrosive cleaning agents such as bleach, lye, acetone, sodium hypochlorite/bleach, sodium hydroxide, formic acid, or solutions containing glutaraldehyde. - The use of a washer-disinfector for cleaning may cause a degradation in performance. - Allow an adequate cooling period after steam sterilization.					
Limitations	Verify functionality prior to re-use.					
Point of Use	This product is provided non-sterile and must be cleaned and sterilized before the first use and any reuse.					
Containment and Transportation	It is recommended that devices are reprocessed as soon as is practical following use.					
Preparation for Decontamination	Follow hospital procedures.					
Cleaning: Automated (Do NOT use ultrasonic washer)	Review the washer-disinfector warning above, before using this cleaning method. Remove devices from instrument trays before placing into washer baskets. Orient devices following recommendations of the washer/disinfector manufacturers. Verify that devices are visually clean after automated cleaning.					
	Recommended Washer Cycle	Phase	Recirculation Time	Water Temperature	Detergent Type and Concentration	
		Pre-Wash	2 minutes	Cold tap water	Not applicable	
		Wash	5 minutes	66 °C (set point)	Neutral enzymatic detergent, pH 6.0-8.0	
		Rinse	1 minute	Hot tap water	Not applicable	
Cleaning: Manual	- Wipe all external surfaces with a cloth dampened with a neutral enzymatic detergent, pH 6.0-8.0. - Brush motor case and collet with a nylon brush dampened with a neutral enzymatic detergent. - Rinse motor thoroughly under running water, collet end pointed down. - Dry with towel. - Verify that devices are visually clean after manual cleaning.					
Disinfection	Follow hospital procedures.					
Packaging	For sterilization, place devices in instrument tray. Devices may be unwrapped, or wrapped with up to two layers of 1-ply polypropylene wrap.					
Sterilization (Temperatures are minimum required, times are minimum required)	Steam Sterilization					
	Cycle	Pre-Vac	Gravity	Pre-Vac	Pre-Vac	
	Temperature	132 °C	132 °C	134 °C	134 °C	
	Time	4 minutes	25 minutes	3 minutes	18 minutes	
	Drying	15 minutes	15 minutes	10 minutes	20 minutes	
	STERRAD	Not validated				
	100% EtO Sterilization Parameters					
	Preconditioning	55 °C, 70 % relative humidity, Vacuum Set Point: 1.3 psia, Time: 30 minutes				
	Sterilization	Temperature	55 °C +/- 4° C			
		Relative Humidity	70 +/- 5 %			
		Ethylene oxide concentration	725 +/- 25 mg/L			
		Gas exposure time (full-cycle)	4 hours			
		Aeration	55 °C +/-4 °C, 12 hours			
	Maintenance, Inspection and Testing	- Inspect devices for any damage before and after each use. If damage is observed, do not use the device until it is repaired. - After cleaning and sterilization, verify functionality prior to re-use.				
	Storage	Store with other sterile devices.				
Additional Information	None					

Note: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing is actually performed using equipment, materials and personnel in the reprocessing facility to achieve the desired result. This normally requires validation and routine monitoring of the process.

Note: All validations performed per current AAMI TIR12, Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers.

Medtronic recommends incineration of devices that have directly contacted patients suspected or confirmed with Transmissible Spongiform Encephalopathy (TSE)/ CJD diagnosis.

 Medtronic		Reprocessing Instructions Legend Attachments/Tubes			Reprocessing Instructions Cleaning and Sterilization instructions are subject to change. Up to date instructions are available online at manuals.medtronic.com . M000030A235 D
Warnings and Precautions	- Do not soak/submerge devices. - Do not use ultrasound to clean devices. - Do not use chlorine based or corrosive cleaning agents such as bleach, lye, acetone, sodium hypochlorite/bleach, sodium hydroxide, formic acid, or solutions containing glutaraldehyde. - The use of a washer-disinfector for cleaning may cause a pre-mature degradation in performance. - Allow an adequate cooling period after steam sterilization. - Use ONLY nylon cleaning brushes. Non-nylon cleaning brushes leave residue that may prevent the tool from being secured properly in the handpiece.				
Limitations	Verify functionality prior to re-use.				
Point of Use	Follow hospital procedures.				
Containment and Transportation	It is recommended that devices are reprocessed as soon as is practical following use.				
Preparation for Decontamination	Follow hospital procedures.				
Cleaning: Automated (Do NOT use ultrasonic washer)	- Review the washer-disinfector warning above, before using this cleaning method. - Manually rinse attachments/tubes under tap water, until no visible soil is noticed, before placing them into the automatic washer. - Remove devices from instrument trays before placing into washer baskets. - Orient devices following recommendations of the washer/disinfector manufacturers.				
	Recommended Washer Cycle	Phase	Recirculation Time	Water Temperature	Detergent Type and Concentration
		Pre-Wash	2 minutes	Cold tap water	Not applicable
		Wash	5 minutes	66 °C (set point)	Neutral enzymatic detergent, pH 6.0-8.0
		Rinse	1 minute	Hot tap water	Not applicable
Cleaning: Manual	- Wipe all attachments and telescoping tubes with a cloth, dampened with a surgical instrument cleaning solution. - Immerse the head of Contra-Angle attachments in surgical instrument cleaning solution and run the motor for 1 minute. - Other attachments and tubes may be mechanically agitated in cleaning solution, but not soaked or immersed. - A nylon brush dampened with a surgical instrument cleaning solution may be used to clean the external surfaces and internal connecting surfaces of the attachments and tubes. - Straight attachments, footed attachments and telescoping straight tubes have special cleaning brushes sized to the attachment's or telescoping tube's internal diameter. Push the brush wet with surgical instrument cleaning solution through the attachment or telescoping tube from rear to front to loosen and remove debris trapped inside. - Move any moveable parts back and forth to allow solution to thoroughly clean attachment, for example, sleeve on footed attachment, perforator attachment. - Rinse thoroughly with tap water.				
Drying and Lubrication	- Thoroughly dry attachments. An air gun may be used to blow moisture out from rear to front of attachment. Note: Medtronic no longer recommends using the Legend attachment cleaning nozzle (PA120), as this may cause some attachments to overheat. - Using an aerosol spray lubricant (such as Pana Spray), perform the following steps to lubricate attachments: - Holding the can approximately 10-15 cm (3-6 inches) away from the attachment, spray all components that move, rotate, or slide with three quick squirts. - Articulate movable components to ensure proper lubrication. - Remove excess lubricant with a clean cloth.				
Disinfection	Follow hospital procedures.				
Packaging	For sterilization, place devices in instrument tray. Devices may be unwrapped, or wrapped with up to two layers of 1-ply polypropylene wrap.				
Sterilization (Temperatures are minimum required, times are minimum required)	Steam Sterilization				
	Cycle	Gravity	Pre-Vac	Pre-Vac (FR/WHO)	Pre-vac (UK)
	Temperature	132 °C	121 °C	134 °C	134 °C
	Time	25 minutes	4 minutes	18 minutes	3 minutes
	Drying	15 minutes	15 minutes	20 minutes	10 minutes
	STERRAD				
	Do not use low temperature hydrogen peroxide gas plasma sterilization due to lumen internal diameter and length restrictions.				
	STERIS				
	Do not use liquid peracetic acid sterilization due to immersion procedure.				
	100% EtO Sterilization Parameters				
	Preconditioning	51-59 °C, 70 ±5% relative humidity, 30 min			
	Sterilization	Temperature	51-59 °C		
		Relative Humidity	70 ±5 %		
		Ethylene oxide concentration	725 +/- 25 mg/L		
		Gas exposure time (full-cycle)	4 hours		
		Aeration	51-59 °C, 18 hours		
Maintenance, Inspection and Testing	- Inspect devices for any damage before and after each use. If damage is observed, do not use the device until it is repaired. - Verify functionality prior to re-use.				
Storage	Store with other sterile devices.				
Additional Information	None				

Note: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing is actually performed using equipment, materials and personnel in the reprocessing facility to achieve the desired result. This normally requires validation and routine monitoring of the process.

Note: All validations performed per current AAMI TIR12, Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers.

Medtronic recommends incineration of devices that have directly contacted patients suspected or confirmed with Transmissible Spongiform Encephalopathy (TSE)/ CJD diagnosis.

Reprocessing Instructions for Reusable Instruments Skeeter Oto-Flex Burs

Reprocessing Instructions

Cleaning and Sterilization instructions are subject to change. Up to date instructions are available online at manuals.medtronic.com.
68E3968 F

Warnings and Precautions	- Remove bur from handpiece before reprocessing of bur and/or handpiece. - Before sterilization, carefully inspect the bur tips. - Burs exhibiting the following conditions should be replaced: - nicks on cutting surfaces - noticeable wear on PTFE bearings - severe bends or crimps on bur shaft - Cold soak in glutaraldehyde, chlorine, or ammonium solutions, or dry heat sterilization is not recommended as damage to the bur may occur. - Store burs in a clean, dry area.																																																										
Limitations	Discard any burs that shows signs of damage or wear.																																																										
Point of Use	- Remove burs from the handpiece before cleaning and sterilizing. - Promptly and thoroughly rinse instruments with deionized water after each use.																																																										
Containment and Transportation	It is recommended that instruments are reprocessed as soon as is practical following use.																																																										
Preparation for Decontamination	Promptly and thoroughly rinse instruments with deionized water after each use.																																																										
Cleaning: Automated	- Remove burs from any instrument trays before placing into washer baskets. - Orient burs following recommendations of washer/disinfector manufacturers. - Use alkaline or neutral pH detergent recommended by washer/disinfector or detergent manufacturers. - Following cleaning, apply a light coating of silicone spray or Pana Spray in the following manner: grasp the PTFE bearing and rotate the bur to assure application of the spray inside the bearing. <table><tr><td>Phase</td><td>Recirculation Time</td><td>Water Temperature</td><td>Detergent Type</td></tr><tr><td>Pre Wash</td><td>3 minutes</td><td>Cold tap water</td><td>Not applicable</td></tr><tr><td>Enzyme Wash</td><td>5 minutes</td><td>Hot tap water</td><td>Neutral pH enzymatic detergent</td></tr><tr><td>Wash 1</td><td>3 minutes</td><td>66 °C set point</td><td>Neutral pH detergent</td></tr><tr><td>Rinse 1</td><td>3 minutes</td><td>Hot tap water</td><td>Not applicable</td></tr></table>					Phase	Recirculation Time	Water Temperature	Detergent Type	Pre Wash	3 minutes	Cold tap water	Not applicable	Enzyme Wash	5 minutes	Hot tap water	Neutral pH enzymatic detergent	Wash 1	3 minutes	66 °C set point	Neutral pH detergent	Rinse 1	3 minutes	Hot tap water	Not applicable																																		
Phase	Recirculation Time	Water Temperature	Detergent Type																																																								
Pre Wash	3 minutes	Cold tap water	Not applicable																																																								
Enzyme Wash	5 minutes	Hot tap water	Neutral pH enzymatic detergent																																																								
Wash 1	3 minutes	66 °C set point	Neutral pH detergent																																																								
Rinse 1	3 minutes	Hot tap water	Not applicable																																																								
Cleaning: Manual	- Soak in lukewarm (less than 43 °C), mild (pH 7.0 - 8.5), enzymatic detergent and deionized water for a minimum of 2 minutes. - Then clean ultrasonically in lukewarm solution of mild detergent (less than 43 °C; pH 7.0 - 8.5) and deionized water for at least 10 minutes. - Rinse thoroughly with deionized water and wipe dry. - Following cleaning, apply a light coating of silicone spray or Pana Spray in the following manner: grasp the PTFE bearing and rotate the bur to assure application of the spray inside the bearing. Note: When using an ultrasonic cleaner or a spray washing machine, follow the manufacturer's recommendations, particularly with regard to articulated instruments and positioning of instruments.																																																										
Disinfection	Do not cold soak in glutaraldehyde.																																																										
Packaging	- A standard sterilization wrap may be used. In the US, an FDA approved surgical wrap is required. Ensure that the pack is large enough to contain the instrument without stressing the seals. - In sets: Instruments may be loaded into dedicated Medtronic Xomed instrument trays. Ensure that cutting edges are protected. Wrap trays using appropriate method.																																																										
Sterilization (Temperatures are minimum required; times are minimum required)	- The sterilization parameters given below should be used for devices that are fully disassembled when disassembly is possible. Use basic aseptic technique during post-sterilization assembly to maintain the sterility of the instrument(s). - All steam cycles have been validated in the wrapped configuration and can be sterilized wrapped or unwrapped. <table><tr><td>Cycle</td><td>Gravity</td><td>Gravity</td><td>Pre-Vac</td><td>Pre-Vac (FR/WHO)</td><td>Pre-vac (UK)</td></tr><tr><td>Temperature</td><td>121 °C</td><td>132 °C</td><td>121 °C</td><td>134 °C</td><td>134 °C</td></tr><tr><td>Time</td><td>30 minutes</td><td>15 minutes</td><td>15 minutes</td><td>15 minutes</td><td>15 minutes</td></tr><tr><td>Drying</td><td colspan="5">15 minutes, or until visibly dry</td></tr><tr><td>STERAD</td><td colspan="5">Not validated.</td></tr></table> 100 % EtO Sterilization Parameters <table><tr><td colspan="2">Conditioning in Chamber</td><td colspan="2">Sterilization</td></tr><tr><td>Temperature</td><td>54 +/- 2 °C</td><td>Temperature</td><td>54 +/- 2 °C</td></tr><tr><td>Relative Humidity</td><td>60 +/- 5 %</td><td>Relative Humidity</td><td>60 +/- 5 %</td></tr><tr><td>Vacuum Set Point</td><td>1.3 psia</td><td>Ethylene Oxide Concentration</td><td>600 +/- 25 mg/L</td></tr><tr><td>Preconditioning Time</td><td>30 minutes</td><td>Gas Exposure Time (full cycle)</td><td>120 minutes</td></tr><tr><td>Aeration</td><td colspan="3">8 hours heated aeration at 48-52 °C or 24 hours aeration at ambient temperature</td></tr></table>					Cycle	Gravity	Gravity	Pre-Vac	Pre-Vac (FR/WHO)	Pre-vac (UK)	Temperature	121 °C	132 °C	121 °C	134 °C	134 °C	Time	30 minutes	15 minutes	15 minutes	15 minutes	15 minutes	Drying	15 minutes, or until visibly dry					STERAD	Not validated.					Conditioning in Chamber		Sterilization		Temperature	54 +/- 2 °C	Temperature	54 +/- 2 °C	Relative Humidity	60 +/- 5 %	Relative Humidity	60 +/- 5 %	Vacuum Set Point	1.3 psia	Ethylene Oxide Concentration	600 +/- 25 mg/L	Preconditioning Time	30 minutes	Gas Exposure Time (full cycle)	120 minutes	Aeration	8 hours heated aeration at 48-52 °C or 24 hours aeration at ambient temperature		
Cycle	Gravity	Gravity	Pre-Vac	Pre-Vac (FR/WHO)	Pre-vac (UK)																																																						
Temperature	121 °C	132 °C	121 °C	134 °C	134 °C																																																						
Time	30 minutes	15 minutes	15 minutes	15 minutes	15 minutes																																																						
Drying	15 minutes, or until visibly dry																																																										
STERAD	Not validated.																																																										
Conditioning in Chamber		Sterilization																																																									
Temperature	54 +/- 2 °C	Temperature	54 +/- 2 °C																																																								
Relative Humidity	60 +/- 5 %	Relative Humidity	60 +/- 5 %																																																								
Vacuum Set Point	1.3 psia	Ethylene Oxide Concentration	600 +/- 25 mg/L																																																								
Preconditioning Time	30 minutes	Gas Exposure Time (full cycle)	120 minutes																																																								
Aeration	8 hours heated aeration at 48-52 °C or 24 hours aeration at ambient temperature																																																										
Maintenance, Inspection and Testing	Discard any burs that shows signs of damage or wear.																																																										
Storage	Store burs in a clean, dry area.																																																										
Additional Information	None																																																										

Note: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing is performed using validated equipment to achieve the desired result. This normally requires validation and routine monitoring of the process.

Some devices have specific assembly instructions. In this case, refer to the assembly insert provided with the device for additional instructions.

All validations performed per current AAMI TIR12, "Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers."

Medtronic recommends incineration of devices that have directly contacted patients suspected or confirmed with Transmissible Spongiform Encephalopathy (TSE)/ CJD diagnosis.

 Reprocessing Instructions Skeeter Handpiece		Reprocessing Instructions Cleaning and Sterilization instructions are subject to change. Up to date instructions are available online at manuals.medtronic.com . 68E3969 G
Warnings and Precautions	- Disconnect the power before cleaning. - Do not fully immerse, or ultrasonically clean, this instrument. - Do not use any cleaning instruments in the cannulated shaft of the handpiece. - Do not cold soak sterilize this instrument in glutaraldehyde. This will void the warranty. - Do not use organic solvents to clean the bur chuck.	
Limitations	After cleaning and sterilization, verify functionality prior to re-use.	
Point of Use	- This product is provided non-sterile and must be cleaned and sterilized before the first use and any reuse. - To remove occasional residual buildup on handpiece cable connector, use a soft brush and isopropyl alcohol.	
Containment and Transportation	It is recommended that instruments are reprocessed as soon as is practical following use.	
Preparation for Decontamination	Disassembly not required, other than removal of the bur.	
Cleaning: Automated (Do NOT use ultrasonic washer)	- Remove instruments and equipment from any sterilization trays before placing into washer baskets. - Orient devices following recommendations of washer/disinfector manufacturers. - Use alkaline or neutral pH detergent recommended by washer/disinfector or detergent manufacturers. - These products have been validated for effective cleaning using an automatic washer/disinfector cycle consisting of a minimum 44 minutes total time, including a pre-wash, main wash & rinse, and thermal rinse. The thermal rinse shall be at least 10 minutes long at a minimum temperature of 60 °C.	
Cleaning: Manual	- Carefully clean with an enzymatic detergent. Do not fully immerse. - The cannulated needle nose should be cleaned by immersing in the detergent solution up to the level of the Bur Release button. Do not use any cleaning instruments in the cannulated shaft of the handpiece. - Rinse by immersing the distal end of the handpiece (up to the Bur Release button) in distilled water, using a gentle swirling motion to flush away residual cleaning solution.	
Drying and Lubrication	- Avoid water accumulation in the motor housing by shaking excess water out with a downward motion. - Silicone spray or Pana Spray should be sprayed into the cannulated shaft of the handpiece prior to sterilization. Apply silicone spray or Pana Spray until surplus lubricant is noted on the outside of the Bur Release Button. Wipe away excess lubricant from the handpiece. Following this procedure will ensure that the bur release mechanism is well lubricated for proper functioning. - Sterilize the handpieces immediately after cleaning.	
Disinfection	Do not cold soak in glutaraldehyde.	
Packaging	- A standard, sterilization wrap may be used. In the US, an FDA approved surgical wrap must be used. Ensure that the pack is large enough to contain the instrument without stressing the seals. - In sets: Instruments may be loaded into dedicated instrument trays or general purpose sterilization trays. Wrap trays using appropriate method.	
Sterilization (Temperatures are minimum required, times are minimum required)	- The sterilization parameters given below should be used for devices that are fully disassembled when disassembly is possible. Use basic aseptic technique during post-sterilization assembly to maintain the sterility of the instrument(s). - All steam sterilization cycles have been validated in the wrapped configuration and instruments can be sterilized wrapped or unwrapped.	
	Steam Sterilization	
	Cycle	Gravity Gravity Pre-vac Pre-vac (FR/WHO) Pre-vac (UK)
	Temperature	121 °C 132 °C 132 °C 134 °C 134 °C
	Time	30 minutes 10 minutes 4 minutes 18 minutes 3 minutes
	Drying	15 minutes, or until visibly dry.
	STERRAD	100S Compatible
	100% EtO Sterilization Parameters	
	Conditioning in Chamber	Sterilization
	Temperature	54 +/- 2 °C Temperature 54 +/- 2 °C
	Relative Humidity	60 +/- 5 % Relative Humidity 60 +/- 5 %
	Vacuum Set Point	1.3 psia Ethylene Oxide Concentration 600 +/- 25 mg/L
	Preconditioning Time	30 minutes Gas Exposure Time (full cycle) 120 minutes
	Aeration	8 hours heated aeration at 48-52 °C or 24 hours aeration at ambient temperature
Maintenance, Inspection and Testing	- Inspect components for any damage before and after each use. If damage is observed do not use the instrument until it is repaired. - After cleaning and sterilization, verify functionality prior to re-use.	
Storage	It is extremely important that the handpiece be rapidly and completely dried before storage to prevent corrosion and residue deposits in the bearing and motor.	
Additional Information	None	

Note: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing is actually performed using equipment, materials and personnel in the reprocessing facility to achieve the desired result. This normally requires validation and routine monitoring of the process.

Note: All validations performed per current AAMI TIR12, Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities. A guide for medical device manufacturers.

Medtronic recommends incineration of devices that have directly contacted patients suspected or confirmed with Transmissible Spongiform Encephalopathy (TSE)/ CJD diagnosis.

Medtronic

Reprocessing Instructions

Visao

Reprocessing Instructions

Cleaning and Sterilization instructions are subject to change.
Up to date instructions are available online at manuals.medtronic.com.
68E3281 F

Warnings and Precautions	• Disconnect the power before cleaning. • Do not fully immerse, or ultrasonically clean, this instrument. • Do not cold soak sterilize this instrument in glutaraldehyde. This will void the warranty. • Do not use organic solvents to clean the bur chuck. • For Visao handpiece cleaning (and Xcalibur Hi-Speed with angled cable), cover handpiece cable connector end with Handpiece Cleaning Cap, Universal, catalog number 3318520. • After completion of the cleaning steps, remove Handpiece Cable Cap or other protective components installed prior to cleaning.				
Limitations	After cleaning and sterilization, verify functionality prior to re-use.				
Point of Use	• This product is provided non-sterile and must be cleaned and sterilized before the first use and any reuse. • To remove occasional residual buildup on handpiece cable connector, use a soft brush and isopropyl alcohol.				
Containment and Transportation	It is recommended that instruments are reprocessed as soon as is practical following use.				
Preparation for Decontamination	Remove the bur from the handpiece, otherwise disassembly is not required.				
Cleaning: Automated	• Remove instruments and equipment from any sterilization trays before placing into washer baskets. Orient devices following recommendations of washer/ disinfectant manufacturers. • Use alkaline or neutral pH detergent recommended by washer/disinfectant or detergent manufacturers. • These products have been validated for effective cleaning using an automatic washer/disinfectant cycle consisting of a minimum 44 minutes total time, including a pre-wash, main wash & rinse, and thermal rinse. The thermal rinse shall be at least 10 minutes long at a minimum temperature of 60°C.				
Cleaning: Manual	• After surgery, clean the irrigation sleeves and bur guards with an enzymatic detergent solution. Wipe the handpiece and cable with disinfectant applied to a clean, non-abrasive cloth. • A chuck brush cleaner (REF 3112500) or an appropriately sized small (plastic bristle) bore brush may be inserted into the distal end of the Visao handpiece, irrigation sleeves and bur guards to assist in removing fluids, tissue, or bone fragments, making sure to clean all passages. Use an enzymatic detergent solution to loosen and remove collected tissues from the unit. • Rinse out the distal end of the handpiece. Shake excess water from the handpiece. • Ensure all water is drained from the cooling housing. If saline was used for cooling during surgery, use distilled water to rinse the housing prior to draining. • Using distilled water, rinse saline from the irrigation nozzles. Drain the nozzle of all water. • Sterilize the handpiece immediately after cleaning.				
Disinfection	Do not cold soak in glutaraldehyde.				
Packaging	• A standard, sterilization wrap may be used. In the US, an FDA approved surgical wrap must be used. Ensure that the pack is large enough to contain the instrument without stressing the seals. • In sets: Instruments may be loaded into dedicated instrument trays or general purpose sterilization trays. Wrap trays using appropriate method.				
Sterilization (Temperatures are minimum required, times are minimum required)	The sterilization parameters given below should be used for devices that are fully disassembled when disassembly is possible. Use basic aseptic technique during post-sterilization assembly to maintain the sterility of the instrument(s). All steam sterilization cycles have been validated in the wrapped configuration and instruments can be sterilized wrapped or unwrapped.				
	Steam Sterilization				
	Cycle	Gravity	Pre-Vac	Pre-Vac (FR/WHO)	Pre-vac (UK)
	Temperature	121 °C	132 °C	134 °C	134 °C
	Time	40 minutes	4 minutes	18 minutes	3 minutes
	Drying	8 minutes, or until visibly dry			
	STERRAD	100S Compatible (Handpiece Only)			
	100% EtO Sterilization Parameters				
	Preconditioning				
	Sterilization	Temperature	54-55 °C		
		Relative Humidity	60 +/- 5 %		
		Ethylene oxide concentration	600 +/- 25 mg/L		
		Gas exposure time (full-cycle)	120 minutes		
		Aeration	48-52 °C, 8 hours		
Maintenance, Inspection and Testing	• Inspect components for any damage before and after each use. If damage is observed do not use the instrument until it is repaired. • After cleaning and sterilization, verify functionality prior to re-use.				
Storage	It is extremely important that the handpiece be rapidly and completely vacuum dried before storage to prevent corrosion and residue deposits in the bearing and motor.				
Additional Information	Increase temperatures higher than those stated when necessary to satisfy governmental or health care facility requirements so long as the temperature does not exceed 149 °C. Heating above 149 °C may damage the handpiece and will void the warranty.				

Note: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing is actually performed using equipment, materials and personnel in the reprocessing facility to achieve the desired result. This normally requires validation and routine monitoring of the process.

Note: All validations performed per current AAMI TIR12, Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers.

Medtronic recommends incineration of devices that have directly contacted patients suspected or confirmed with Transmissible Spongiform Encephalopathy (TSE)/ CJD diagnosis.

 Medtronic		Reprocessing Instructions Indigo High-Speed Otologic Drill			Reprocessing Instructions Cleaning and Sterilization instructions are subject to change. Up to date instructions are available online at manuals.medtronic.com . 68E4187 D	
Warnings and Precautions		- Do not soak/submerge devices. - Do not use ultrasound to clean devices. - Do not use chlorine based or corrosive cleaning agents such as bleach, lye, acetone, sodium hypochlorite/bleach, sodium hydroxide, formic acid, or solutions containing glutaraldehyde. - The use of a washer-disinfector for cleaning may cause a pre-mature degradation in performance. - Allow an adequate cooling period after steam sterilization.				
Limitations		Verify functionality prior to re-use.				
Point of Use		This product is provided non-sterile and must be cleaned and sterilized before the first use and any reuse.				
Containment and Transportation		It is recommended that devices are reprocessed as soon as is practical following use.				
Preparation for Decontamination		Follow hospital procedures.				
Cleaning: Detergent		Use a cleaning agent that is suitable for use on aluminum surfaces. A neutral enzymatic to mild alkaline agent (pH 6.0 to 10.5) is preferred. If a washer-disinfector is used, see the instructions supplied with the washer-disinfector machine to select the recommended cleaning agent.				
Cleaning: Automated (Do NOT use ultrasonic washer)		- Remove devices from instrument trays before placing into washer baskets. - Prior to cleaning, cover the drill cable connector end with Handpiece Cleaning Cap, catalog number 3318520. - Orient devices following recommendations of the washer/disinfector manufacturers. - After completion of the cleaning steps, remove Handpiece Cleaning Cap or other protective components prior to sterilization.				
		Recommended Washer Cycle	Phase	Recirculation Time	Water Temperature	Detergent Type and Concentration
			Pre-Wash	2 minutes	Cold tap water	Not applicable
			Wash	5 minutes	66 °C (set point)	Neutral enzymatic to mild alkaline detergent, pH 6.0-10.5
			Rinse	1 minute	Hot tap water	Not applicable
Cleaning: Manual		- Prior to cleaning, cover the drill cable connector end with Handpiece Cleaning Cap, catalog number 3318520. - Wipe all external surface of the motor and cable with a cloth dampened with the detergent prepared with lukewarm tap water. - Brush motor case and collet with a nylon brush dampened with a neutral enzymatic detergent. - Rinse motor thoroughly under running water, collet end pointed down. - Dry collet and motor with lint free towel - After completion of the cleaning steps, remove Handpiece Cleaning Cap or other protective components prior to sterilization. - Verify that devices are visually clean after manual cleaning				
Disinfection		Follow hospital procedures.				
Packaging		Place devices in instrument tray, and double wrap instrument case with 1-ply polypropylene wrap. In the US, an FDA approved surgical wrap must be used.				
Sterilization (Temperatures are minimum required, times are minimum required)		Steam Sterilization				
		Cycle	Gravity	Pre-Vac	Pre-Vac (FR/WHO)	Pre-vac (UK)
		Temperature	132 °C	132 °C	134 °C	134 °C
		Time	25 minutes	4 minutes	18 minutes	3 minutes
		Drying	20 minutes			
		STERRAD	Not validated			
		100% EtO Sterilization Parameters				
		Preconditioning	55 °C, 70 % relative humidity, 30 minutes			
		Sterilization	Temperature		55 °C	
			Relative Humidity		70 %	
			Ethylene oxide concentration		725 mg/L	
			Gas exposure time (full-cycle)		240 minutes	
			Aeration		53-57 °C, 18 hours	
Maintenance, Inspection and Testing		- Inspect components for any damage before and after each use. If damage is observed do not use the instrument until it is repaired. - After cleaning and sterilization, verify functionality prior to re-use.				
Storage		Store with other sterile devices.				
Additional Information		None				

Note: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing is actually performed using equipment, materials and personnel in the reprocessing facility to achieve the desired result. This normally requires validation and routine monitoring of the process.

Note: All validations performed per current AAMI TIR12, Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers.

Medtronic recommends incineration of devices that have directly contacted patients suspected or confirmed with Transmissible Spongiform Encephalopathy (TSE)/ CJD diagnosis.

Medtronic

Reprocessing Instructions

Indigo High-Speed Otologic Drill Attachments

Reprocessing Instructions
Cleaning and Sterilization instructions are subject to change. Up to date instructions are available online at manuals.medtronic.com.
68E4188 D

Warnings and Precautions	- Do not soak/submerge devices. - Do not use ultrasound to clean devices. - Do not use chlorine based or corrosive cleaning agents such as bleach, lye, acetone, sodium hypochlorite/bleach, sodium hydroxide, formic acid, or solutions containing glutaraldehyde. - The use of a washer-disinfector for cleaning may cause a pre-mature degradation in performance. - Allow an adequate cooling period after steam sterilization.				
Limitations	Verify functionality prior to re-use.				
Point of Use	This product is provided non-sterile and must be cleaned and sterilized before the first use and any reuse.				
Containment and Transportation	It is recommended that devices are reprocessed as soon as is practical following use.				
Preparation for Decontamination	Follow hospital procedures.				
Cleaning: Detergent	A neutral enzymatic to mild alkaline agent (pH 6.0 to 10.5) is preferred. If a washer-disinfector is used, see the instructions supplied with the washer-disinfector machine to select the recommended cleaning agent.				
Cleaning: Automated (Do NOT use ultrasonic washer)	- All attachments must be thoroughly rinsed manually with tap water, ensuring all hard to reach areas are rinsed, prior to transfer to automatic washer for processing. - Remove devices from instrument trays before placing into washer baskets. - Orient devices following recommendations of the washer/disinfector manufacturers.				
	Recommended Washer Cycle	Phase	Recirculation Time	Water Temperature	Detergent Type and Concentration
		Pre-Wash	2 minutes	Cold tap water	Not applicable
		Wash	5 minutes	66°C (set point)	Neutral enzymatic to mild alkaline detergent, pH 6.0-10.5
		Rinse	1 minute	Hot tap water	Not applicable
Cleaning: Manual	- Use the detergent, prepared with lukewarm tap water for the cleaning process. - Thoroughly wipe the straight and angled attachments with a cloth dampened with the detergent. - Wipe entire attachment until all gross soil has been removed. - Using a brush wetted with the detergent, brush the attachments to clean the external surfaces and internal connecting surfaces. - For the angled attachment: Wet an appropriately sized cleaning brush (Ø 2.4mm) with the detergent. Insert the brush into the bore at the front of the attachment. Brush the bore to loosen debris trapped inside. Rinse the bore to remove debris. - For the straight attachment: Wet an appropriately sized cleaning brush (Ø 2.4mm) with the detergent. Push the brush through the straight attachment from the rear to front to loosen and remove debris trapped inside. Rinse bore to remove debris. - Place one half of the straight or angled attachment into the detergent. Do not immerse the entire attachment. Gently agitate the attachments in the detergent and actuate any moveable parts. - Place the other half of the attachment into the detergent and repeat. Do not immerse the entire attachment. - Rinse the attachment thoroughly with tap water. Flush both ends to remove detergent.				
Drying and Lubrication	- Thoroughly dry the attachments. An air gun may be used on the straight attachment to blow moisture out from the rear to front. - Using an aerosol spray lubricant (such as Pana Spray), perform the following steps to lubricate attachments: - Holding the can approximately 10-15 cm (3-6 inches) away from the attachment, spray all components that move, rotate, or slide with three quick squirts. - Articulate movable components to ensure proper lubrication. - Remove excess lubricant with a clean cloth.				
Disinfection	Follow hospital procedures.				
Packaging	Place devices in instrument tray, and double wrap instrument case with 1-ply polypropylene wrap. In the US, an FDA approved surgical wrap must be used.				
Sterilization (Temperatures are minimum required, times are minimum required)	Steam Sterilization				
	Cycle	Gravity	Pre-Vac	Pre-Vac (FR/WHO)	Pre-vac (UK)
	Temperature	132 °C	132 °C	134 °C	134 °C
	Time	25 minutes	4 minutes	18 minutes	3 minutes
	Drying	20 minutes			
	STERRAD	Not validated			
	100% EtO Sterilization Parameters				
	Preconditioning	55 °C, 70 % relative humidity, 30 minutes			
	Sterilization	Temperature	55 °C		
		Relative Humidity	70 %		
		Ethylene oxide concentration	725 mg/L		
		Gas exposure time (full-cycle)	240 minutes		
		Aeration	53-57 °C, 18 hours		
Maintenance, Inspection and Testing	- Inspect components for any damage before and after each use. If damage is observed do not use the instrument until it is repaired. - After cleaning and sterilization, verify functionality prior to re-use.				
Storage	Store with other sterile devices.				
Additional Information	None				

Note: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing is actually performed using equipment, materials and personnel in the reprocessing facility to achieve the desired result. This normally requires validation and routine monitoring of the process.

Note: All validations performed per current AAMI TIR12, Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers.

Medtronic recommends incineration of devices that have directly contacted patients suspected or confirmed with Transmissible Spongiform Encephalopathy (TSE)/ CJD diagnosis.

 Medtronic Reprocessing Instructions Microsaws		Reprocessing Instructions Cleaning and Sterilization instructions are subject to change. Up to date instructions are available online at manuals.medtronic.com. M000030A231 D		
Warnings and Precautions	- Do not soak/submerge devices. - Do not use ultrasound to clean Midas Rex Microsaw devices. - Do not use chlorine-based or corrosive cleaning agents such as bleach, lye, acetone, sodium hypochlorite/bleach, sodium hydroxide, formic acid, or solutions containing glutaraldehyde. - Allow an adequate cooling period after steam sterilization. - The use of a washer-disinfector for cleaning may cause a premature degradation in performance.			
Limitations	Verify functionality prior to reuse.			
Point of Use	Follow hospital procedures.			
Containment and Transportation	It is recommended that devices are reprocessed as soon as is practical following use.			
Preparation for Decontamination	- Disconnect handpiece from console. - Remove irrigation tubing/clips, saw blade, and cable management clamps.			
Cleaning: Automated (Do NOT use ultrasonic washer)	- Prior to placing the saws in the automated washer, manually rinse the saw body with the distal end pointed down under tap water, until no visible soil is noticed. - Move any movable parts back and forth to allow water to reach hard-to-rinse areas - Transfer the devices into the washer for processing, making sure that the oscillating saw collet is open. Place the devices in the washer to facilitate drainage. - Verify that devices are visually clean after automated cleaning.			
	Recommended Washer Cycle	Phase	Recirculation Time	Water Temperature
		Pre-Wash	2 minutes	Cold tap water
		Wash	5 minutes	66 °C (set point)
		Rinse	1 minute	Hot tap water
Detergent Type and Concentration				
Cleaning: Manual	- Rinse the device thoroughly under running water, with the collet end pointed down. - Wipe the saw body and cable with a cloth dampened with a neutral enzymatic detergent, pH 6.0-8.0 (henceforth referred to as "cleaning solution"). - With the blade collet end pointed down, brush saw body, blade-release button, and extended lever tip (if applicable) with a nylon brush that has been dipped in cleaning solution. - Sagittal Only: - Dip the distal end of the saw in cleaning solution and agitate for 10-15 seconds. - Press and hold the blade-release button, then repeat step 1. - With the blade-release button pressed, flush approximately 15 mL of cleaning solution from a syringe into the blade opening. - Flush all ports with 15 mL of cleaning solution. - Repeat steps 1-4 using water instead of cleaning solution. - Reciprocating Only: - With the distal end facing up, squeeze 30 mL of cleaning solution from a syringe into the distal end of the device so that the solution thoroughly reaches the inside of the collet. - While the device is in the vertical position, rotate the collet knob back and forth three times. - Dump the solution from the device. - Repeat steps 1-3 a second time with cleaning solution - Repeat steps 1-3 a third time with water. - Rinse the device under running water, ensuring that water enters all ports/openings. Move any moveable parts to allow the solution to reach all areas. Rinse until there is no visible evidence of soil or debris.			
Drying and Lubrication	- Dry the entire device with a lint-free towel. - Reciprocating and Oscillating Only: If desired, complete the following steps using Pana Spray to lubricate the saw collet: - Holding the can approximately 10-15 cm (4-6 inches) away from the saw collet, spray into the hole of the shaft (reciprocating saw) or onto the proximal/bottom portion of the collet near the locked/unlocked indicators (oscillating saw) with three quick squirts - Articulate the collet to ensure proper lubrication. - Remove excess lubricant with a clean cloth			
Disinfection	Follow hospital procedures.			
Packaging	For sterilization, place devices in instrument tray. Devices may be unwrapped, or wrapped with up to two layers of 1-ply polypropylene wrap			
Sterilization (Temperatures are minimum required, times are minimum required)	Steam Sterilization			
	Cycle	Pre-Vac	Pre-Vac (FR/WHO)	Pre-Vac (UK)
	Temperature	132 °C	134 °C	136 °C
	Time	3 minutes	15 minutes	3 minutes
	Drying	30 minutes	30 minutes	30 minutes
	STERAD	Do not use low temperature hydrogen peroxide gas plasma sterilization due to lumen internal diameter and length restrictions.		
	STERRIS	Do not use liquid peracetic acid sterilization due to immersion procedure		
Maintenance, Inspection and Testing	- Inspect components for any damage before and after each use. If damage is observed, do not use the instrument until it is repaired. - After cleaning and sterilization, verify functionality prior to reuse.			
Storage	Store with other sterile devices			
Additional Information	None			

Note: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for reuse. It remains the responsibility of the processor to ensure that the reprocessing is actually performed using equipment, materials and personnel in the reprocessing facility to achieve the desired result. This normally requires validation and routine monitoring of the process.

Note: All validations performed per current AAMI TIR12, Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers.

Medtronic recommends incineration of devices that have directly contacted patients suspected or confirmed with Transmissible Spongiform Encephalopathy (TSE)/ CJD diagnosis.

Reprocessing Instructions

POWEEASE Driver

Reprocessing Instructions
Cleaning and Sterilization instructions are subject to change. Up to date instructions are available online at manuals.medtronic.com.
68E4189 C

Warnings and Precautions	- Do not soak/submerge devices. - Do not use ultrasound to clean devices. - Do not use chlorine based or corrosive cleaning agents such as bleach, lye, acetone, sodium hypochlorite/bleach, sodium hydroxide, formic acid, or solutions containing glutaraldehyde. - The use of a washer-disinfector for cleaning may cause a pre-mature degradation in performance. - Allow an adequate cooling period after steam sterilization.				
Limitations	Verify functionality prior to re-use.				
Point of Use	No particular requirements.				
Containment and Transportation	It is recommended that devices are reprocessed as soon as is practical following use.				
Preparation for Decontamination	Disconnect handpiece from console.				
Cleaning: Manual	- Rinse tap water through the open lumen of the handpiece, collet pointed down, until all evidence of soil and debris is removed. Actuate and rinse collet. - Wipe external surface of the motor and cable with a cloth dampened with neutral enzymatic detergent, pH 6.0-8.0. - Thoroughly clean the motor case, collet, and trigger area with a nylon brush dampened with neutral enzymatic detergent, pH 6.0-8.0. Actuate the collet and trigger while brushing. - Use a lumen brush to brush the entire length of the lumen. - Rinse entire device, with collet pointed down, under running tap water. - Immerse the collet portion of the device into neutral enzymatic detergent, pH 6.0-8.0 for a minimum of 2 minutes. Actuate and thoroughly brush the collet while immersed. A syringe should be used to flush hard to reach areas. - Rinse entire device, with collet pointed down, under running tap water. A syringe should be used to rinse hard to reach areas. - Dry with lint free cloth. - Verify that device is visually clean.				
Cleaning: Automated (Do NOT use ultrasonic washer)	The POWEREASE must undergo the manual cleaning process prior to processing in a washer/disinfector. 1. Rinse tap water through the open lumen of the handpiece, collet pointed down until all evidence of soil and debris is removed. Actuate and rinse collet. 2. Wipe external surface of the motor and cable with a cloth dampened with neutral enzymatic detergent. 3. Brush motor case, collet, and trigger area with a nylon brush dampened with neutral enzymatic detergent. Actuate the collet and trigger while brushing. 4. Use a lumen brush to brush the entire length of the lumen. 5. Rinse entire device, with the collet pointed down, under running tap water. 6. Immerse the collet portion of the device into neutral enzymatic detergent for a minimum of 2 minutes. Actuate and brush the collet while immersed. A syringe should be used to flush hard to reach areas. 7. Rinse entire device, with the collet pointed down, under running tap water. A syringe should be used to rinse hard to reach areas. 8. Place POWEREASE into a suitable washer/disinfector basket and process through a washer/disinfector cycle.				
	Recommended Washer Cycle	Phase	Recirculation Time	Water Temperature	Detergent Type and Concentration
		Pre-Wash	2 minutes	Cold tap water	Not applicable
		Wash	5 minutes	66 °C (set point)	Neutral enzymatic detergent, pH 6.0-8.0
		Rinse	1 minute	Hot tap water	Not applicable
Disinfection	Follow hospital procedures.				
Packaging	Steam sterilization may be performed wrapped or unwrapped. If wrapped, in the US, an FDA approved surgical wrap must be used.				
Sterilization (Temperatures are minimum required, times are minimum required)	Steam Sterilization				
	Cycle	Gravity	Pre-Vac	Pre-Vac (FR/WHO)	Pre-Vac (UK)
	Temperature	132 °C	132 °C	134 °C	134 °C
	Time	25 minutes	4 minutes	18 minutes	1 minute
	Drying	40 minutes			
	STERRAD	Not validated			
	STERIS	Not validated			
	100% EtO Sterilization Parameters	Not validated			
Maintenance, Inspection and Testing	- Inspect devices for any damage before and after each use. If damage is observed, do not use the device until it is repaired. - After cleaning and sterilization, verify functionality prior to re-use.				
Storage	Store with other sterile devices.				
Additional Information	None				

Note: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing is actually performed using equipment, materials and personnel in the reprocessing facility to achieve the desired result. This normally requires validation and routine monitoring of the process.

Note: All validations performed per current AAMI TIR12, Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers.

Medtronic recommends incineration of devices that have directly contacted patients suspected or confirmed with Transmissible Spongiform Encephalopathy (TSE)/ CJD diagnosis.

SYMBOLS

The following symbols can appear on this device and related packaging.

			Package Contents		Follow Instructions For Use		Conforms to ANSI/AAMI ES60601-1, Certified to CAN/CSA-C22.2 No.60601-1, IEC/EN 60601-1
	EMC Compliance Mark				Do Not Immerse		Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.
	Fuse	IPX1	Protected Against Vertical Water Drops				Authorized Representative in the European Community
	Use By Date	IPX7	Protected Against The Effects Of Temporary Immersion In Water				Precaution: Pinch Hazard. Keep Fingers Clear Of Rollers
	Accessory		Type BF Applied Part				Contains DEHP (di-2-ethyl hexyl phthalate)
			Start/Stop		Foot Pedal Connector		Protective Earth
	Output		RF Transmitter (Interference May Occur)		Fine Irrigant Adjustment		Equipotential Ground Connector
	Is Approximately Equal To		Consult Instructions for Use		Left Foot Control Unit Button / Mode Button		
			Stim Bur Connector		Right Foot Control Unit Button / Control Button		RoHS - Environmental friendly use period - China (SJ/T11364-2006)
	Quantity		NIM Console Connector		Top Foot Control Unit Button / Handpiece Button		Do Not Dispose Of This Product In The Unsorted Municipal Waste Stream. Dispose Of This Product According to Local Regulations. See Recycling.Medtronic.com for Instructions On Proper Disposal Of This Product.
	Not Greater Than 120vac		Electrical High Speed Handpiece Connector				If the single use symbol is on the device label then this device is designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.
			World Wide Standard for Medical Tubing Diameter				
	Applied Part Duty Cycle		Date of Manufacture		Handpiece		Accessory
	Catalog Number						Adapter
	Lot Number			1	Pump Head 1		Attachment
	Serial Number		Temperature Limitation Outer Limits = Transport Inner Limits = Use	2	Pump Head 2		Bone Mill
	USA Only		Ergonomically Engineered		Refurbished		Brush
							Control Unit
							Dissecting Tool
							Instrument Case
							Lubricant/Diffuser
							Motor
							Multi-Use Disposable Attachment
							Regulator

The information contained in this document is accurate at time of publication. Medtronic reserves the right to make changes to the product described in this manual. Refer to [manuals.medtronic.com](#) for the current version.

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© 2016 Medtronic, Inc.