The Swiss Instruments

EN Operation Manual

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
Congratulations on your purchase of a NOUVAG AG the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories. Thank you for the confidence shown in our products. Please consult the instruction manual for the use and maintenance of the device in order to ensure that it will function properly and efficiently for many years.


- Please read instructions carefully before operating!

Table of Contents

Introduction	
1. INTENDED USE	5
2. MEDICAL GROUNDS	5
3. CONTRAINDICATIONS	6
4. SIDE EFFECTS	6
5. OPERATION CONDITIONS	6
6. DEVICE OVERVIEW	8
7. TECHNICAL DATA	15
8. DESIGN AND OPERATION	20
1. Operating principle	20
2. Method of application	20
3. Startup of the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4 with accessories	21
4. Operation. of the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4 with accessories	28
5. Startup of the Dental surgery portable motor device models: MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories	41
6. Operation. of the Dental surgery portable motor device models: MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories	48
7. How to use handpieces and contra angles	63
8. How to use SANICAR Mobile Device Trolley	69
9. Suitcase	73
9. SCOPE OF DELIVERY	74
10. CLEANING OF DEVICE COMPONENTS	79
11. MAINTENANCE	86
1. Replacing the control unit fuse	86
2. Safety inspections	87
12. SAFETY REQUIREMENTS	89
1. Manufacturer's EC declaration	89
2. Integrated peristaltic pump	93
3. Modification and misuse	93
4. Essential requirements	93
5. During use	94
6. Safety inspections	94
13. MEDICAL DEVICE CLASSIFICATION	95
14. DEVICE PACKAGING AND MARKING	96
1. Marking	96
2. Packaging	99
15. TRANSPORTATION AND STORAGE CONDITIONS	112
16. SERVICE LIFE	112
17. DISPOSAL	113
18. WARRANTY COVERAGE	113
APPENDIX 1. Device components	115
APPENDIX 2. Tables of device components	121

INTRODUCTION

Dental implantology is a relatively new branch of dental surgery dealing with dentofacial system reconstruction using alloplastic materials. Dental arches, jaws, jaw joints and face shape are reconstructed using various implants. Surgeon performing implantation must be able to use filler materials, membranes, he must know the basics of guided tissue regeneration because very often when the implant is installed the surgeon has to face the problems of bone and mucosa defect closure. In this respect it is important to offer maximum operational comfort to the surgeon. The Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories was designed to make dental surgeon operation most comfortable. MD design is simple and understandable; it consists of control unit, electronic motor, pedal and tubing set for irrigation fluid. When designing the MD the main effort was focused on further development of the electronic motor and control unit. Additionally, special attention was paid to ergonomic component of the device. All protruding parts and sharp edges of the unit body were smoothed and flattened for easy cleaning; display is designed in such a way that one can easily read all the information at a glance. Setting up is done with single function keys and accompanied by audible feedback which guarantee accurate setup. Inserting irrigation tube set is effortless and simple. A broad range of handpieces allows for a great number of various applications. Taking all this into consideration one can say that any Implantology or microsurgery procedures can be done effectively with these devices.

Pic. 1. The Dental surgery portable motor device MD 11 model.

Pic. 2. The Dental surgery portable motor device MD 30 model.

1. INTENDED USE

The Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories are used primarily in dental implantology, microsurgical applications, oral and maxillofacial surgical procedures.

2. MEDICAL GROUNDS

The Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories is used primarily in dental implantology in following cases:

- loss of tooth due to various causes
- soreness and loss of tooth and bone tissue
- jaw bone condensing after loss of tooth

it is used specifically in implantology for

- milling and drilling the implant bed
- cutting the thread for the implant
- screwing in the implant
- removing the fixture mount
- placing the cover screw

Use in microsurgery

- medical grounds for drilling, milling and sawing bone as well as screw insertion into bone

3. CONTRAINDICATIONS

Relative or absolute contraindications may result from specific cases when use of cutting surgical electric powered instruments hazards the patient's health considerably. For example: acute or chronic inflammation or tumors in surgical area.

4. SIDE EFFECT

Implantology and microsurgery are conventional medical procedures. There are no known non-tolerable side effects. Acceptable side effects are sensation of pain and adverse anesthesia reaction.

X-ray film is required to undergo full course of treatment. This procedure must be conducted under supervision of trained medical personnel.

5. OPERATING CONDITIONS

For correct operation of the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories, the following requirements must be met:

Ambient conditions

	Transport and storage:	Operation:
Relative humidity:	Max 90%	30- 80 %
Temperature:	0°C to 60 °C	10 °C to 45 °C
Atmospheric pressure:	700 hPa to 1060 hPa	800 hPa to 1060 hPa

	Please carefully read these instructions and instructions for other devices used during surgery procedures.		The Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories may only be operated by qualified and trained personnel. All requirements and directions of this manual must be followed.
	Device preparation must be conducted according to operation manual to assure safety of patients, users and third parties.		Improper use or repair of the device and failure to observe these instructions relieve us from any obligation arising from warranty provisions or other claims.
	The use of third party products is the responsibility of the operator. Functionality and patient safety cannot be guaranteed with third party accessories.		Prior to using this device, before startup and before operating, the user must always ensure that the device and accessories are in good working order and are clean, sterile and operational.

	Keep the device out of direct sunlight, central heating radiators, excessive dust, humidity, vibration and accidental mechanical damage.		Repairs may only be performed by authorized NOUVAG AG service technicians.
	The integrated peristaltic pump is used to cool tissue in order to prevent damage to the tissue. It may only be operated with watery solutions, such as 0.9% Sodium Chloride irrigation solution or "Ringer" solution (not included in delivery). Supplying medication using the integrated pump is expressly prohibited.		Use Nou-Clean spray (not included in delivery) for maintenance and care of the motor, handpiece and contra angles. Using other care products can result in malfunction and/ or cause the warranty to be revoked.
	The device is not sterile on delivery. All sterilizable parts must be sterilized before use (see Chapter "Cleaning, disinfection and sterilization").		Do not use device in the vicinity of flammable mixtures!
	Ventilation slots of the devices must not be obstructed in order to prevent overheating.		Check the voltage settings of the device. Make sure they match the supply voltage.
	Never operate the clamping mechanism of the handpiece or contra angles while the system is running. This could result in instrument damage.		Operating of the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories must be conducted under constant surveillance of qualified medical personnel.
	Handpieces and contra angles may only be attached when the electric motor is not running.		To prevent injury, never touch drill bits or burs while they are still running.
	At operation at a patient, make sure to produce as less as possible frictional heat. By excessive thermal influence on the tissue, the tissue easily gets necrotized. The heat generation is directly proportional to the speed delivery of the instrument and the contact pressure delivered by the instrument.		

Modification and misuse

- Modification or misuse of the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories are prohibited. The manufacturer is not liable for any damages resulting from unauthorized modifications or manipulations. The warranty will be cancelled.
- Use of the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories outside the indications described in Section "Intended Use" is prohibited. The user or operator is solely responsible for any such use.

Operator skills:**Education:**

- medical education at university level

Operator's health:

- surgeon should be with normal physical and mental abilities

Language command:

- English or
- language of the Operation Manual


Medical use:

The Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD **11-4**; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories must be operated by qualified dental surgeons in appropriately equipped surgical rooms in dental departments of medical institutions. MD operation is stable at ambient temperature from 10°C to 45°C, relative humidity 30% - 80%. MD is resistant to exposure to biological liquids and body tissue contacting device during surgical operation. MD is resistant to pre-sterilization cleaning, disinfection and sterilization. **When operating medical devices, the user must always wear gloves, protective goggles and medical mask.**

6. DEVICE OVERVIEW


Overview The Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4 with accessories

Front view


- | | |
|---|--|
| 1. Pedal socket, on the device rear panel | 13. Tubing set (tube 2m. sterile
(Marketing Authorization №FSZ (ΦC3) 2010/07320 dated
05.07.2010) |
| 2. Vario pedal | 14. Stand for irrigation fluid bottle with irrigation liquid
(Marketing Authorization №FSZ (ΦC3) 2010/07320 dated
05.07.2010) |
| 3. Contra angle | 15. Drip chamber |
| 4. Clip for tubing set attachment at handpiece and
contra angles | 16. Bleed valve |
| 5. Handpiece cradle
(Marketing Authorization №FSZ (ΦC3) 2010/07320 dated
05.07.2010) | 17. Bottle with irrigation fluid |
| 6. Electronic motor 21 | 18. Type plate with type designation, reference number,
serial number, information on power supply and
device fuse |
| 7. Clip for tubing set attachment at motor 21 cable | 19. Power plug socket |
| 8. Motor socket | 20. Main switch |
| 9. Control unit | 21. Fuse compartment |
| 10. Display | 22. Nou-Clean spray nozzle for electronic motor
maintenance (Marketing Authorization №FSZ (ΦC3)
2010/07320 dated 05.07.2010) |
| 11. Release key for irrigation tubing set bracket
(1706 tube 2m., sterile
(Marketing Authorization №FSZ (ΦC3) 2010/07320 dated
05.07.2010) | 23. Nou-Clean spray nozzle for handpiece and contra
angles maintenance and care (Marketing Authorization №FSZ
(ΦC3) 2010/07320 dated 05.07.2010) |
| 12. Peristaltic pump (built-in cooling system) | |

Rear view

Overview The Dental surgery portable motor device models: MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories

Front view


1. Pedal socket, on the device rear panel
2. Vario pedal
3. Contra angle and handpiece
4. Clip for tubing sett attachment at handpieces and contra angles

14. Irrigation tubing set (1706 tube 2m, sterile (Marketing Authorization №FSZ (ФЦ3) 2010/07320 dated 05.07.2010)
15. Stand for irrigation fluid bottle with irrigation liquid (Marketing Authorization №FSZ (ФЦ3) 2010/07320 dated 05.07.2010)
16. Drip chamber

- | | | | |
|-----|--|-----|--|
| 5. | Handpiece cradle
(Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010) | 17. | Bleed valve |
| 6. | Electronic motor 21 | 18. | Bottle with irrigation fluid |
| 7. | Clip for tubing set attachment at motor cable | 19. | Type plate with type designation, reference number, serial number, information on power supply and device fuse |
| 8. | 1 indicator light for each motor | 20. | Power plug socket |
| 9. | 2 motor sockets | 21. | Main switch |
| 10. | Control unit | 22. | Fuse compartment |
| 11. | Control unit display | 23. | Nou-Clean spray nozzle for electronic motor maintenance
(Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010) |
| 12. | Release key for irrigation tubing set bracket
(1706 tube 2m, sterile
(Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010) | 24. | Nou-Clean spray nozzle for handpiece maintenance and care (Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010) |
| 13. | Peristaltic pump (built-in cooling pump) | | |


Rear view

On the front panel [1] of MD there are electronic motor sockets (2 sockets for MD 30 model and 1 socket for MD 11 model). The electronic motors [2] are connected to the handpiece with Intra Coupling (ISO3964) connector. Depending on the requirements, various handpieces [3] can be connected to the electronic motor. On the rear panel of the device there is a pedal socket. Vario pedal connected to the unit with Intra Coupling (ISO3964). Electronic motor rotation speed, rotation direction and switching on the pump (built-in pump) of cooling system are controlled with the pedal treadplate. Bottle with irrigation fluid (not included in delivery) is hung on a bottle stand. Irrigation fluid is supplied to the handpiece during oration through the tubing set (Ref. 1706 2m, sterile (ISO3964). The Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories are equipped with control panel with all necessary control functions on it. On the control panel of the devise there is an LCD display featuring parameters related to device operation. Operating parameters are set with the keys on the panel. On the rear panel there is a main switch ON/OFF "O/I". Pushing "OFF" the device can be switched off at any moment irrespective of routine shut down procedure.


To prepare MD for operation do the following:

1. Plug in the electronic motor cable into the motor socket.


Pic. 3. 21 electric motor socket on the control unit.

2. If required, plug in the electronic motor 2 cable into the motor socket 2 on the control unit panel (applicable to the Dental surgery portable motor device models: MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories).

Pic. 4. Two 21 electric motor sockets on the control unit.

3. Plug in the pedal cable into the pedal socket on the rear panel of the control unit.


Pic. 5. Pedal socket on the rear panel of the control unit.

Pic. 6. The MD with power and pedal cables plugged in.

4. Insert the stand for irrigation fluid bottle into the stand holder on the control unit.

Pic. 7. Insertion of the stand for irrigation fluid bottle.

Detailed MD preparation procedures are given in chapter "Startup".

The MD parts and accessories description is given in the Appendix 2 of this Operation Manual

Table. 1. Intended use of parts and accessories of the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories.

No	Part name	Intended use
1.	2097nou Electronic Motor 21	torque transmission to handpiece
2.	Contra Angles	used to treat tooth surface during implantology manipulations
3.	5052nou Drilling Contra Angle 20:1 LED, shank $\varnothing 2.35$ mm with internal/external cooling with push button system	contra angle providing shadowless instrument surface and surgical area

4.	1706 Tube set, sterile, 2m, single-use (Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010)	irrigation fluid supply to handpiece and surgical area
5.	1770 Stand for irrigation fluid (Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010)	irrigation liquid placement (not included in delivery)
6.	1170 Handpiece cradle (Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010)	electronic motor and handpiece placement during surgical manipulations
7.	Clip set for tube set attachment	tube set attachment to handpiece and motor
8.	Tube set for internal cooling	handpiece internal irrigation coolant supply
9.	1712 Internal cooling nozzle for drilling handpiece with latch system	securing internal cooling tube on the handpiece
10.	Spray nozzle clips for handpiece internal cooling	setting and attachment of coolant tube of handpiece
11.	1958 Spray (Nou-Clean) nozzle attachment for surgical instrument e- coupling (Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010) 1974 Spray (Nou-Clean) nozzle attachment for electronic motor 21 (Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010)	adapters for cleaning handpieces and electronic motor with Nou-Clean spray (not included in delivery)
12.	1029 Suitcase for MD11/MD30	The MD11/MD30 portable implantology motor system storage and transporting
13.	1898d SANICAR Device Trolley	MD transport and placement when operating

When operating this medical devices, the user must always wear gloves, goggles and protective mask (not included in delivery, purchased separately). Handpieces are in direct contact with tissues and biological liquids of patients (saliva, blood, tissue particles, purulent effluent). Tubing set for irrigation (1773, 39150) with attached (1777, 39152) Y-connectors, handpiece internal cooling tubing set (1712, 39116), spray nozzle clips for handpiece internal cooling (39149, 39117) and clips for attaching tubes (1873, 1881) to the handpiece contact the patients indirectly.

Materials used:

The tubing set for handpiece cooling are made of silicon rubber COH-Rlastic; Y-connectors (1777, 39116) and internal cooling nozzle for contra angle with push button system (39116) are made of medical grade stainless steel 1.4301 and internal cooling nozzle for contra angle with latch system is made of medical grade stainless steel 1.4301 1.4035; spray clips for handpiece internal cooling (39149, 39117) are made of LCP Vectra A130 polymer; clips for tube set attachment to handpieces and electronic motor (1873, 1881) are made of polyphenylsulfone (PPSU) Radel R5100; control units (3330 and 3335) are made of ABS UL 94-VO plastic and control unit buttons are made of polyacetal (POM); electronic motor 21 (2097nou) is made of medical grade stainless steel 1.4305; electronic motor cable is made of polyetherketone (PEEK).

7. TECHNICAL DATA

Essential technical data on The Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories is given bellow.

Table 2 MD 11 Technical data: MD 11 control unit.

Voltage AC	100 B ~ / 115 V ~ / 230 V ~, 50-60 Hz
Fuse	2 fuses T3.15 AL, 250 VAC
Power consumption	120 VA
Number of electronic motors	1
Number of electronic motor connectors	1
Number of pedals	1
Motor speed	300 – 40'000 rpm
Transmission ratio	1:1 speed 300-40'000 rpm torque 1-6 Ncm (not included in delivery) 16:1 speed 20-2'400 rpm torque 5-27 Ncm (not included in delivery) 20:1 speed 15-1'700 rpm torque 10-70 Ncm 32:1 speed 15-1'000 rpm torque 10-55 Ncm (drilling function for all numbers)
Drive coupling	Intra coupling ISO3964
Instrument connector	ISO 3964 connector (DIN 13940)
Max. Electric motor cable length	2 m
Max. Pedal cable length	3.15 m
Cooling system	Integrated peristaltic pump (built-in cooling pump) adjustable 1% - 100%
Volumetric displacement of pump	1 -100 ml/min
Pump continuous sound pressure level (low)	30 dB
Electrical protection class	I
Method of protection against electric shock	Type B
Operating time rating, Max.	1 min (ON), 3 min (OFF)

Device type	hospital
User's interface	LCD 95x52 mm Display
Electronic control	Microprocessor: input voltage– 230/115/110 V, programming voltage – 5 V
Software	Built-in: Version 2.06/Print V4.0, date 08.04.16, Class I (A)
Response time	2 sec
Software security level	Class B
Overall program size	17'584 lines of source code
CPU frequency	8MHz
Pedal	IPX8
Dimensions (W x L x H)	260 x 250 x 110 mm
Net weight MD11 control unit	3.3 kg

Table 3. Technical data: 21 electric motor for MD11

Rotation speed	300 – 40'000 rpm
Electronic motor type	BLDC-Motor with Hall sensor
Max. Motor torque	6 Ncm
Weight, Max.	300/280g with/without cable
Continuous sound pressure level	50 dB
Rotor power output	120 W
Connection	E-Type (ISO 3964)
Dimensions (D x L), Max.	21 x 106 mm
Cable length Max.	2 m
Sterilization	Steam Max. temperature 135°C

Table 4. Technical data: MD11 control pedal.

Type	Vario pedal with 2 keys
Protection against contact and continual submerging	IPX8
Weight, Max.	1 kg
Dimensions (W x L x H), Max.	140 x 240 x 90
Pedal press force, not exceeding	40 N
Cable length, Max.	3,15 m
Keys functions	Yellow: sets rotation direction right/ left Blue: ON/OFF peristaltic pump

Table 5. Technical data: MD30 control unit.

Voltage AC	100 V ~ / 115 V ~ / 230 V ~, 50-60 Hz
Fuse	2 fuses T3.15 AL, 250 VAC
Power consumption	120 VA
Number of electronic motors	1
Number of motor connectors	2
Number of pedals	1
Motor speed	300 – 50'000 rpm
Transmission ratio	1:5 speed 1'500-250'000 rpm torque 1-1 Ncm (not included in

	delivery) 1:3 speed 900-150'000 rpm torque 1-2 Nm (not included in delivery) 1:2 speed 600-100'000 rpm torque 1-2 Nm (not included in delivery) 1:1 speed 300-50'000 rpm torque 1-6 (not included in delivery) 4:1 speed 75-12'500 rpm torque 1-20 Nm (purchased separately) 16:1 speed 20-3'000 rpm torque 5-27 Nm (not included in delivery) 20:1 скорость 15-2'500 rpm torque 10-70 Nm 32:1 speed 10-1'500 rpm torque 10-55 Nm Handpiece for saws 15'000 rpm torque 6 Nm (not included in delivery) Mucotome speed 8'000 rpm torque 6 (not included in delivery)
Electronic motor coupling	internal coupling ISO3964
Instrument connector	ISO 3964 Connector (DIN 13940)
Electronic motor cable length, Max.	2 m
Pedal cable length, Max.	3,15 m
Cooling system	Integrated peristaltic pump (built-in cooling pump) adjustable 1% - 100%
Volumetric displacement of pump	1 -100 ml/min
Operating pump noise level (low)	30 dB
Electrical protection class	I
Method of protection against electric shock	Typa B
Operating time rating, Max.	1 min (ON), 3 min (OFF)
Device type	hospital
User interface	LCD 95x52 mm Display
Electronic control	Microprocessor: input voltage– 230/115/110 V, programming voltage – 5 V
Software	Built-in: Version 1.13/Print V4.0, date 13.04.16, Class I (A)
Response time	2 sec
Software protection level	Class B
Overall program size	18'682 lines of source code
Bit-timing frequency	8 MHz
Pedal	IPX8
Dimensions (W x L x H), Max.	260 x 250 x 110 mm
MD 30 control unit net weight, Max.	3.7 kg

Table 6. Technical data: 21 electric motor for MD30.

Motor speed	from 300 to 50'000 rpm
Motor type	BLDC-Motor with Hall sensor
Electronic motor torque, Max.	6 Nm
Weight, Max.	300/280g with/without cable
Continuous sound pressure level	50dB
Rotor power output	120 W
Connection	E-Type (ISO 3964)
Dimensions (D x L), Max.	21 x 106

Cable length, Max.	2 m
Sterilization	Steam Max. temperature 135°C

Table 7. Technical data: MD30 control pedal.

Type	Vario pedal with 3 keys
Number of speed modes	3 (low (100 – 400 rpm) intermediate (300 - 700) high (500 – 1000 rpm))
Protection against contact and continual submerging	IPX8
Weight, Max.	1.68 kg
Dimensions (W x L x H), mm, Max.	255 x 235 x 60
Pedal press force, not exceeding	40 N
Cable length, Max.	3,15 m
Keys functions	Yellow: sets rotation direction right/ left Green: program selection Blue: ON/OFF peristaltic pump

Table 8. Technical data: Handpieces.

Material	Stainless Steel 1.4305.
Surface roughness	Ra 0.8
Hardness of material	200MPa
Surface finish	Non-gloss
Burr pulling out force kgf	Min. 30
Instrument extracting/inserting force N	Min. 45
Thrust kgf	Max. 3
Force of friction btw a motor and a handpiece	25 N
Continuous sound pressure level	50 dB
Burr shaft diameter	2,35 mm
Rate of water supply	Min 50 ml/ min. at 200kPa (2 Bar); 1 – 100 ml/ min.
Surgical area irrigation	coolant one channel supply coolant two channel supply
Number of sterilization cycles without quality losses	250
Light source design:	light source in handpiece
Light source brightness	22 000 lx
Light source voltage	Max. 25V
Sterilization	Steam Max. temperature 135°C
Average time of failure (hour) Min.	For contra angles - 300

Table 9. Technical data: SANICAR device trolley

Setting in motion	manhandled
Intended use	MD transport and placement when operating
Classification	Internal use, multipurpose, repairable
Number of panels	2
Basic dimensions mm, Max.	Height: 1125 mm Top Panel Size: 450 mm x 360 mm Middle Panel Size: 344 mm x Base size: 470mm x 445mm Slope angle of panels - leveled
Weight, Max.	18.5kg
Weight capacity, Max	50+/-0,5 kg
Voltage	230 V/50 Hz
Input current	Max. 16A
Number of sockets/ type	3 / Schuko
Cable length, Max.	5 m
Electrical protection class	I
Clearance between one of the rolls and the floor on horizontal surface	Max. 5 mm
Propel effort	Max. 100 N (10 kgf)
Rotary force free of roll jam	Max. 0.35 N
Pressure force applied to the pedal to lock the rolls	Max. 150N (15 kgf)
Disinfection resistance	Resistant to disinfection (80% alcohol or microbiologically active disinfecting solutions)
Average service life	7 years
Shelf life	3 years
Limit state criteria	Restoration impossible or unpractical due to technical-and-economic reasons
Material	Steel
Coating	Powder spraying
Colour	White
Number of rolls	4, 2 rear rolls are lockable
Roll diameter, Max.	80 mm

Table 10. Suitcase for MD11/MD30

Dimensions (W x L x H) cm, Max.	52 x 44 x 24
Weight, Max.	5
Body material	Aluminum
Inside lining	Synthetic Rubber Foam
Handle material	Aluminum and Plastic
Handle load to failure	Max. 400 N
Acceptable clearance between led and body	Max. 0.4 mm
Angle between led and body when opened	90° -120°
Number of locks	2 locks 5 cm x 5 cm each

8. DESIGN AND OPERATION

1. Operating principle

These Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories consist of control unit, electronic motor, control pedal and tube set for irrigation liquid (Ref. 1706 tube 2m., sterile (Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010). MD 11 model is provided with one and MD 30 model with two motor sockets which makes operator's work easier and bring better results. The control unit runs all other components of the MD, it is started by pressing the pedal. Operation parameters are set using control unit display. Rotating instruments are latched in the handpiece, their rotation and required speed are rendered by electronic motor power. Two motor versions of MD has separate program for each motor, it provides continuous operation. Optimized electronic motor control system allows for precision manipulations in compliance with individually selected parameters to make every use of the device safe. The MD is used for Drilling, grinding, polishing and screwing in implants (implant installation). There are also special programs for saws and Mucotome (not included in delivery). Built-in irrigation system (pump) for sterile coolant supply reduces instrument heating and minimizes risk of tissue damage and inflectional complications.

2 Application method

Dental implantology:

Surgical procedure under local anesthesia.

Implantation procedure is basically the same regardless of number of teeth to be replaced.

Gingiva incision is made to expose the bone where the implant will be placed. If there is bone deficiency, bone grafting is recommended. When healthy bone tissue is properly treated using the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories, a special burr is used to prepare the bone for implantation. When the bone preparations are done, the implant is placed and screwed in the implant bed using the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories. Then the gingival cut is sutured.

Healing time varies from 3 to 6 months. That's the period of time, normally required for the implant to become part of lower jaw, this process is often called osseointegration. Surgical suture is usually removed within 7 to 14 days after operation. When tissue is ready for further manipulations the dental support for a new crown is attached to the implant. To achieve natural tooth transparency degree, the porcelain crowns are often used with zirconium abutment attached to the implant.

When the cast dental model is done a natural looking implant crown is manufactured. The crown is attached and adhered to the implant fixture. The finished crown looks like a natural tooth.

In many cases dental implants are included in general treatment regimen, which can be focused on replacement of several teeth. In some cases all compromised teeth can be replaced with implants depending on bone tissue condition. Other treatment regimens may recommend implants for framework and Implant-supported denture attachment.

Microsurgery:


Instruments driven by the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories are used for

bone and tissue preparation in compliance with medical grounds. During operation the cooling system irrigates surgical area to prevent damage (necrotization) to tissue caused by excessive heat.

Operation safety considerations:

Operating the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories means certain risk due to such parts of the devices as electronic motors with handpieces and rotating instruments with sharp edges. Their use includes possible risks such as damage to blood vessels and other organs. Some of these risks are operator's (user's) responsibility. Operations performed under constant surveillance of qualified personnel trained in dental surgery. Thus, close control of surgical site during operations allows to take urgent measures to avoid undesirable side effects and complications. Other hazards are included in Novag AG Risk Management System.

3. Startup of the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4 with accessories

Pic. 8. The Dental surgery portable motor device MD 11 model.

1 Device setup

- Place the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4 with accessories and instruments on an even, non slip surface and make sure you have good access to all controls.
 - Do not allow the operating range of the device (including the cable and contra angle) to be compromised by limiting factors.
 - The system display must be fully visible at all times.
 - The pedal must be placed within stepping distance between the patient and the surgeon.
 - It must be explicitly ensured that no objects can fall on the pedal.
 - The power plug at the rear of the device must be accessible at all times.

2 Connecting to the power supply

Before plugging the power cable into the power socket for the first time, you must check the supply voltage setting next to the power switch.

If the voltage shown does not correspond to the local mains voltage, the grey fuse holder must be set to the correct voltage:


1. Switch off the device.
2. Unplug the power cable.
3. Use a screwdriver to open the fuse slot.
4. Remove the fuse holder.
5. Remove the grey fuse holder and reinsert it so that the local mains voltage setting is shown in the small window.
6. Slide the fuse holder back in and close the fuse slot.
7. Check the mains voltage shown on the fuse slot.
8. Plug the power cable back into the device.

In order to prevent the risk of electric shock, the device may only be connected to a power network with a PE protective ground conductor.


3 Device preparation

1. Sterilize the motor 21 (the motor is not sterile on delivery). If the motor has already been sterilized: when removing the motor from the sterile packaging, ensure that the sterile packaging is not damaged and that the sterility indicator confirms sterility (if no sterility indicator is provided, the sterile packaging must at least show the date on which the shelf life of the sterile item is due to expire). 
2. Insert the stand for the irrigation fluid into the stand holder.


Pic. 9. Insertion of the stand for irrigation fluid bottle.

3. Plug the motor plug of the electronic motor into the motor socket.


Pic. 10. 21 electronic motor socket on the control unit.

4. Plug the Vario pedal plug into the pedal socket at the rear of the control unit.

Pic. 11. Pedal socket on the rear panel of the control unit.

Pic. 12. Device with plugged in the electric supply and pedal cables.

5. Attach the sterilized handpiece to the electronic motor. Press the handpiece firmly onto the electronic motor unit until it clicks into position and make sure it is secure by moving it slightly in the opposite direction.

Pic. 13. The handpiece attached to the electronic motor resting on the handpiece cradle.

6. Assemble the tubing set. Attach tubing set Ref. 1706 2m., sterile (Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010) for cooling of the contra angle, like described in the following article.

Use only Nouvag AG tube set Ref. 1706 2m. sterile (Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010) for handpiece cooling. Otherwise the function cannot be guaranteed.

Check the integrity of the tubing set Ref. 1706 2m. sterile (Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010) and ensure that the packaging is not damaged. Using non-sterile tubing set can result in serious infection and, in extreme cases, can be fatal.


When inserting the tubing set Ref. 1706 2m. sterile (Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010), watch the arrow on the cover of the pump compartment. It indicates the flow direction of the cooling liquid.


Do not regulate the amount of irrigation liquid using the roller clamp on the tubing set, with the MD 11, this is regulated instead using the integrated pump. For this reason, make sure to open the roller clamp as far as it will go.

- A) Press release key for set bracket (on top of the control unit) to open the pump.
- B) The compartment with the integrated tubing bracket opens.
- C) Place the tubing set into the tubing bracket provided in such a way that the part of the tubing set with the spike exits the pump towards the rear of the device. Make sure the tubing is secure.
- D) With the tubing set inserted, press the compartment downwards until it clicks into place.

Picture 14. To access the pump, press release key by the display on the front of the control unit. It can be easily done with one hand, without inadvertent moving of the device


- E) Insert the spike at the end of the tubing set into the irrigation fluid bottle and hang the bottle onto the stand.
- F) Open the roller clamp on the tubing set as far as it will go.
- G) Open the bleed valve beneath the drip chamber.
- H) Connect the control unit to the power socket.

Ensure that the operating voltage setting corresponds to the local mains voltage.


4 Assembly of external irrigation system


- A) Connect the open end of the Ref. 1706 2m. sterile (Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010) irrigation tube set with the cooling pipe of the Contra Angle.
- B) Attach Clip Ref. 1881 to the tube set Ref. 1706 2m.
- C) Attach Clip (Ref. 1881) with tube set Ref. 1706 2m. to the Contra Angle.
- D) Connect the electronic motor with the Contra Angle.
- E) Attach Clip Ref. 1873 to the tube set Ref. 1706 2m. sterile (Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010).
- F) Attach Clip (Ref. 1873) with tube set to the cable of the motor.
- G) Ready to use Drill unit with attached external cooling system.

If required, secure additional clips to the motor cable.


5 Assembly of internal and external irrigation system

- A) Attach open end of the tube set Ref. 1706 2m. sterile (PY № ΦC3 2010/07320 or 05.07.2010 г.) to the single out branch of the Y-connector Ref.1777.
- B) Connect a 6 cm piece of tube Ref, 1773 to each of the branching tubes off the Y-connector.
- C) Connect one of the 6 cm pieces of tube, that branch out of Y-connector, with the internal cooling tube of the contra angle (Ref. 39116).
- D) Connect the second piece of tube that branches out of the Y-connector with the external cooling pipe of the contra angle (welded on laterally to the contra angle head).

- E) Both 6 cm pieces of tube are now connected to the cooling system of the contra angle.
- F) Secure a tube clip Ref. 1881 to each of the short tube pieces.
- G) Secure tube clips (Ref. 1881) to the contra angle/
- H) Contra angle with attached tube clip (Ref. 1881).

- I) Fix a gray clip Ref.1873 to the tube set leading to the Y-connector.
- J) Fix tube set with the clip (Ref. 1873) to the motor cable.
- K) Motor cable with ready attached tube set to the motor cable. Secure additional clips to the motor cable if needed.
- L) Tube set routing with Y-connector for internal and external cooling of the instrument.

Essential accessories for internal cooling:

Ref. 1777
Y-connector


Ref.1873
Motor cable clips

16 cm connection
tubes

4. Operation. The Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4 with accessories.

1 Switching the device on and off


Picture 15. Main switch on the rear panel.

The power switch I/O (at the rear) is used to switch the control unit on and off.


The device can be switched off at any time irrespective of any procedure for switching off the device.

2 Overview: control elements on operating panel

- A) Display: Shows the operating values (See chapter "Overview: Standard display").
- B) "Handpiece" key: Selection of the required handpiece/ contra angle.
- The factory default settings contain the handpiece/ contra angles 1:1 (not included in delivery), 16:1 (not included in delivery), 20:1 and 32:1. If you want to use 70:1 contra angle (not included in delivery) it has to be activated in the configuration menu (refer to "Configuration menu/ handpiece existing").
- C) "Program" key: Selection of program 1 to 4.
- The number of usable programs can be altered from 2 to 6 available programs (refer to "Configuration menu/ handpiece existing").
 - By pressing both "HP + PR" keys at the same time the programs can be reset to the factory default setting.
- D) "Speed" keys:
- Restrict the maximum speed that can be selected using the pedal.
- « ▲ » increases the maximum speed « ▼ » reduces the maximum speed
- By pressing both "Speed" keys at the same time "Speed ▲ + ▼" the handpiece calibration will be started (refer to "Calibrating handpieces").
- E) "Torque" keys:
- Restrict the maximum torque.
 - "▲" increases the maximum torque "▼" reduces the maximum torque.
 - By pressing both "Torque" keys at the same time "Torque ▲ + ▼" the torque models AL or AS are selected (see chapter "Torque limit function, AL mode (Automatic Limiter)" and chapter "Torque limit function, AS mode (Automatic Stopper)").
 - If a speed of 20 rpm is set, irrespective of the selected handpiece/ contra angle, the pressing of both torque keys, "Torque ▲ + ▼" causes the activation of the thread cutting function.
- F) "Pump" keys:
- Adjusting the maximum pump flow rate, controlled by the pedal.
 - "▲" increases the maximum supply quantity "▼" reduces the maximum supply quantity
 - By pressing both "Pump" keys at the same time "Pump ▲ + ▼" the pump will be put on call, pressing again will switch it off.

3 Overview: Standard display

- A) Information line
Information and error messages are displayed here. Display is illuminated red at error messages. At standby the display shows "Nouvag AG - MD11 V2.0" including the indication of device history.
- B) Program
Shows the selected program number 1 to 4.
In the configuration menu the number of available programs can be set from 2 to 6 (refer to "Configuration menu/ number of programs").
- C) Pump
The numerical value shows the pump flow rate in percent and the drop symbol together with the on/off indication shows if the pump is in standby mode or switched off.
- D) Rotation direction of the motor
The arrow indicates the rotation direction set for the motor. Change the rotation direction by pressing the button "FOR/REV" on the pedal.
- E) Current speed
At stand still of the device the maximum speed is displayed. If pedal is activated and motor begins to run the current speed is displayed in real time.
- F) Transmission ratio of contra angle
Shows transmission ratio of used contra angle, e.g. 1:1 (not included in delivery), 16:1 (not included in delivery), 20:1, 32:1 or 70:1 (not included in delivery).
- G) Speed range
Shows speed range of the handpiece used.
- H) Maximum torque
Shows maximum torque setting.
- I) Bar graph (current torque)
The bar graph is providing a graphical representation of the current torque. All bars active means max. torque reached.

J) AS/AL

Indication of the Auto Stopper (AS) mode or the Automatic Limiter (AL) mode (see chapter "Torque limit function, AL mode" and chapter "Torque limit function, AS mode").

The pump does not begin to operate until the motor has been activated by pressing the pedal.


4 Adjusting the programs

The values for operation settings depend on the connected handpiece/ contra angle as well as the task to be performed.

4.1 1: Selecting the handpiece/ contra angle to be used

The selected handpiece/ contra angle must comply with the factual used handpiece/ contra angle.

Press the "HP" key to select the handpiece/ contra angle.
Press the "PR" key to choose between programs 1 to 4.

In the configuration menu the number of available programs can be set from 2 to 6 (refer to "Configuration menu/ number of programs").

Table 11. Appropriate handpiece/ contra angles to be used

Name of the handpiece/ contra angles with transmission ratio	Display	Speed min. rpm	Speed max. rpm	Torque min. Ncm	Torque max. Ncm	AS-range (Factory def.) rpm	Limit AS- range rpm
Drill handpiece, 1:1 (not included in delivery)	1: 1	300	40000	1	6	-	-
Drill contra angle, 16:1 (not included in delivery)	16: 1	20	2400	5	27	up to 20	20 - 45
Drill contra angle, 20:1	20: 1	15	1 700	10	70	up to 20	15 - 45
Drill contra angle, 32:1	32: 1	10	1000	10	55	up to 20	10 - 45
Drill contra angle, 70:1 (not included in delivery)	70: 1	5	600	10	55	up to 20	5 - 45

4.2 Step 2: Calibrating handpiece or contra angles

To make sure the displayed parameters correspond with the actual, measurable parameters of the handpiece or contra angle, it is recommended to calibrate each handpiece or contra angle on a regular basis.


It is a procedure as simple as it is important to guarantee safety and precision of each handpiece or contra angle being used.

After you take care of all prior preparations such as sterilization, maintenance and care of handpiece, device preparation and the selection of the handpiece of use, the calibration procedure can be performed.

The calibration of handpieces or contra angles guarantees accurate torque. Due to wear as well as varying lubrication of the handpiece or contra angles and lack of maintenance and care, the distribution of torque can vary widely.

1. Select the handpiece, attached to the motor, by pressing the handpiece “**HP**” key. Make sure that this handpiece attached to the motor is also on the display.
2. Hold motor with mounted handpiece or contra angle in your hand at a safe distance from your body.
3. Press both “Speed” keys at the same time (Speed ▲ + Speed ▼). Display indicates “Calibrating handpiece XX”.

4. Motor and handpiece start running and pass several speed cycles up to maximum speed.

5. After a tone emitted and a message displayed the calibration is finished. The display returns to normal mode.

If a handpiece is not working within the range of the calibration algorithm, even when it was cleaned and lubricated before, the display shows an error with a red lit display "Handpiece XX is faulty". This indicates soiling, wear or technical defect of the handpiece. These handpieces have to undergo a maintenance and care cycle or have to be repaired or replaced.

At the calibration procedure, the handpiece are tested on their torque delivery behavior. With the 1:1 handpiece (not included in delivery) the control unit is additionally adjusted to altered conditions of the handpiece to keep the 1:1 handpiece delivering the appropriate torque.

4.3 Step 3: Setting the speed

The possible speed range depends on contra angle to be used. The maximum speed within this speed range can be restricted to the required value. Using the pedal, the speed can be varied from the minimum speed up to the maximum speed as set.

Setting the speed:

Press the "Speed" keys "▲" to increase or "▼" to decrease the maximum speed. When key is pressed continuously the speed will be shown in fast forward mode.

4.4 Step 4: Setting the torque

Once the speed has been selected, the torque can be limited from the torque range of the corresponding handpiece or contra angle.

Depending on the speed, the torque mode AL or AS is applied.

For information on AL and AS mode, see chapter "Speed reduction AL (Automatic Limited)" and chapter "Torque limit function AS (Automatic Stopper)".

Press the “Torque” keys “▲” to increase or “▼” to decrease the maximum torque. By pressing the key continuously the torques are changed in fast forward mode.


The buildup of torque is displayed as a bar graph. When maximum torque is reached all elements of the bar graph are fully displayed.

4.5 Step 5: Setting the pump supply quantity (built-in cooling pump)


Press the “Pump” keys “▲” to increase or “▼” to decrease the pump supply quantity. By pressing the key continuously the values are fast forwarded.

Minimum and maximum of pump (built-in cooling pump) supply quantity as well as incremental steps can be adjusted in the Configuration Menu (refer to “Configuration Menu, Parameter Level 2, Pump”).

To activate or deactivate the pump, press both “Pump” keys at the same time, “Pump ▲ + ▼” or push the foot switch “Pump” shortly.

5. Torque limit function, AL mode (Automatic Limiter)

The AL function limits the torque applied to the instrument, for example to prevent crack initiation and bone fracture.

The speed of the instrument remains constant until the preset torque is reached. If the applied force gets over this limit the speed will be reduced, if necessary down to a stop, but the torque remains constant. If the applied force is reduced the speed picks up again.

On the display this procedure is shown as a bar graph. The bars in the graph fill up to their full range as the torque increases to its preset value. When the bar graph shows full capacity the speed will be reduced. As soon as the force on the instrument is reduced the torque decreases and the speed starts to pick up again as shown on the display.

AL mode is active at all speeds, except at the speeds where the AS mode takes effect.

6. Torque limit function, AS mode (Automatic Stopper)

The AS function limits the torque applied to the device. As soon as the preconfigured torque is reached, the electric motor stops immediately. The electronic motor no longer generates any force. In order to restart the electronic motor, the pedal must be released and pressed again.

On the display, the bar graph completes up to its full range, until the maximum torque is reached, then drops to zero.

The function "AS" is just active for certain contra angles and only within a certain speed range.

From the minimal speed up to the speeds shown at the right the AS mode is activated automatically (refer to "Configuration Menu, Parameter Level 2, AS-Zone").

Handpiece	16: 1	20: 1	32: 1	70: 1
Speed rpm	20	20	20	20

The upper limit of the AS range can be adjusted in the configuration menu.

In the following speed ranges the AS-/AL- mode can also be switched manually by pressing both torque keys "Torque ▲ + ▼" at the same time.

Handpiece	16: 1	20: 1	32: 1	70: 1
From rpm:	20	15	10	5
To rpm:	45	45	45	45

7. Storing various programs

With the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4 with accessories up to 6 different settings can be set as fixed program (program 1 to program 6). Which program is currently active is shown on the display.

At switch-off, the settings made by user are automatically saved. This includes the following parameters:

- | | |
|--|---|
| <ul style="list-style-type: none"> • handpiece/ transmission ratio • maximal speed • maximal torque | <ul style="list-style-type: none"> • pump On/Off • pump performance • AS/AL mode |
|--|---|

In the configuration menu the number of available programs can be set from 2 to 6 (refer to "Configuration Menu/ number of programs").

To change a program go to the specific parameter and change it. When the device is switched off all parameters are saved in that program.

When the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4 with accessories is switched on, the display shows the prior used program with the prior used handpiece or contra angle.

8. Configuration menu

In the configuration menu the user can customize the device for one's purposes. The parameters are organized in two levels. The following information and parameters are stored in there.

- | | |
|---|---|
| <ul style="list-style-type: none"> • Software version • Mainboard serial number • Language selection, DE/EN • Display illumination • Number of programs • Power on, at last program • Operating hours MD 11 • Operating hours Motor • Operating hours of irrigation pump | <ul style="list-style-type: none"> • Error log • Activation of available handpiece • Speed limitation for each handpiece • Range of action of AS mode • Pump behavior • Motor behavior • Behavior of torque at reverse mode • Reset to factory default settings |
|---|---|

Be cautious when changing parameters. Unusual behavior of instruments while operating may provoke false reaction and jeopardize the patient. Every setting and the new behavior of the instrument has to be verified and tested.

1. Access to configuration menu:

- Press **"HP"** (handpiece) and **"Speed ▼"** (speed) key for 3 seconds until you hear a tone. On the display the first position of the configuration menu shows up.

The arrows at the beginning of the lines indicate the operation in the Configuration Menu.

2. Selection of parameters:

- Press **"HP"** or **"PR"** to choose the desired parameter.
- Press **"Speed ▼"** to select the parameter.
- Change the value by pressing **"HP"** for up or **"PR"** for down.
- To confirm the changes made, press **"Speed ▼"** for about one second, until you hear a short tone.
- To abort the settings made, press **"Speed ▼"** briefly. The setting returns to its previous value.

3. Exit from Configuration Menu:

- To leave the Configuration Menu press **“HP”** and **“PR”** for **3 seconds** until you hear a tone.

Table 12. Parameter Level 1

Group/Parameter	Rights	Factory	Definition
Software/ Version	read	VX.XXXX	Shows current software version
Hardware/ Serial number MB	read	XXXXXXXXXXXX	Shows serial number of the main board
Language, 0=DE, 1=EN	read/write	1	Select language between German and English
Backlight/brightness (0..10)	read/write	9	Display brightness, changeable :0-10
Number of programs	read/write	4	Select between 2 and 6 available programs
Power on at last program	read/write	Yes	Display of the last program after switching device on
Operating hours/Control Unit	read	0	Shows operating hours of Control Unit
Operating hours Motor	read	0	Shows operating hours of motor 1
Operating hours Pump	read	0	Shows operating hours of pump
Error memory/ 1-8	read	0	8 Error messages in chronological order

Handpiece activation	Name of handpiece on display	Choice	Factory default	Definition
Handpiece existing/ HP 02 (not included in delivery)	16: 1	yes/no	no	Only one handpiece or contra angle can be selected
Handpiece existing/ HP 03	20: 1	yes/no	yes	
Handpiece existing/ HP 04	32: 1	yes/no	no	
Handpiece existing/ HP 05 (not included in delivery)	70: 1	yes/no	no	

Table 13. Parameter Level 2

To change the values in Level 2 the password **“9403”** must be entered. *The password can not be changed.*

- Enter the access code: Press **Speed** ▼ .
- Select the access code by pressing **“HP”** (handpiece) to increase, or **“PR”** (program) to decrease the value.
- To confirm the access code press **Speed** ▼ for one second, until sound is emitted.
(for fast forwarding or backwarding keep keys pressed)

Handpiece Max. Speed	Name of handpiece on the display	Speed range rpm	Factory default	Definition
Handpiece max speed/ HP 01 (not included in delivery)	1:1	300 – 40'000	40'000	Limit the maximum speed of you handpieces according to your own experience.
Handpiece max speed/ HP 02 (not included in delivery)	16:1	20 – 2'400	2'400	
Handpiece max speed/ HP 03	20:1	15 – 1'700	1'700	
Handpiece max speed/ HP 04	32:1	10 – 1'000	1'000	

Handpiece max speed/ HP 05 (not included in delivery)	70:1	5 - 600	600
---	------	---------	-----

AS-Zone for handpieces	Handpiece on display	Selection ..to ..rpm	Factory default	Definition
Handpiece AS-Mode/HP 02 (not included in delivery)	16:1	20 - 45	20	Selection range for max. AS speed
Handpiece AS-Mode/HP 03	20:1	15 - 45	20	Selection range for max. AS speed
Handpiece AS-Mode/HP 04	32:1	10 - 45	20	Selection range for max. AS speed
Handpiece AS-Mode/HP 05 (not included in delivery)	70:1	5 - 45	20	Selection range for max. AS speed

Pump parameters (built-in cooling pump)	Range	Factory defaults	Definition
Pump/Backwards turn mode variable	No/Yes	Yes	The pressures in the tube vary according to the pump speed. In "variable mode" this is considered, to prevent spilling when pump is switched off in backwards mode.
Pump/Way backwards	1 – 100 %	25 %	Specify how far the pump turns backwards.
Pump/Speed backwards	10 – 50 %	33 %	Specify how fast the pump has to turn backwards to prevent spilling after the switch off of the motor.
Pump/Range 1 increment	1 – 10 %	5 %	Adjustment steps in section 1.
Pump/Range 1 end	5 – 50 %	10 %	Set the range where section 1 is active.
Pump/Range 2 increment	1 – 10 %	5 %	Adjustment steps in section 2.
Pump/Range 2 end	10 – 90 %	50 %	Set the range where section 2 is active.
Pump/Range 3 increment	1 – 10 %	10 %	Adjustment steps in section 3.
Pump/Range 3 end	20 – 100 %	100 %	Set the range where section 3 is active.

Motor type 21, 40'000 rpm	Range	Factory default	Definition
Motor type 21, 40'000 rpm/Min. speed	300 – 5'000 rpm	300 rpm	Set the min. speed of the motor.
Motor type 21, 40'000 rpm/Max. speed	5'000 – 40'000 rpm	40'000 rpm	Set the max. speed of the motor.
Motor type 21, 40'000 rpm/ Start ramp	1 – 1'000 ms / 10'000 rpm	100 ms	Set acceleration time to 10'000 rpm.
Motor type 21, 40'000 rpm/ Stop ramp	1 – 1'000 ms / 10'000 rpm	50 ms	Set breaking time from 10'000-0 rpm.

The maximum torque of the motor is boosted while starting the reverse rotation to disengage stuck implants or screws.

Reverse torque	Range	Factory default	Definition
Reverse torque/Increase	5 – 30 %	25 %	Increase of the selected torque in reverse rotation.
Reverse torque/Increase time	100 – 2'000 ms	500 ms	Time during which the torque is raised.


Attention: Upon resetting of factory default all parameters will appear with factory default values (except date and time).


Resetting to factory default	Choice	Factory default	Definition
Default value/ Set default value	Yes/ No	-	Resetting all parameters of the Configuration Menu to factory default.

Exit from Configuration Menu

- To leave the Configuration Menu press "HP" (handpiece) and "PR" (program) for 3 seconds, until you hear a long tone.

9. Operation using the Vario pedal

1. Key "Pump":

Pressing the key briefly: switching the pump on or off (see information on display). Pressing the key longer: increases the pump speed (see information on display).

2. Key "FOR/REV":

Pressing the key briefly: switching the direction of revolution (see information on display).

3. Treadplate:

With the treadplate pedal the motor speed is variably adjusted and pump is activated.

Table 14.

Treadplate...	Motor...	Pump...
... not pressed	Motor off	Pump off
... pressed gently	Motor runs slowly	Pump on, if pump "On" displayed (speed set on the control unit)
... pressed all the way down	Motor runs at maximum speed (speed as set on the control unit)	Pump on, if pump "On" displayed (speed set on the control unit)

For safety reasons, the unit can only be operated by pedal.

10. Functional check

Prior to the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4 with accessories startup or use of accessory equipment, the user must always ensure that each individual component is in good working order, free from defects and in clean, sterile and operational. All inscriptions on the device and its accessories must be readable and there must be no loose parts on the device. Once the device is switched on, the most recent settings entered appear on the display.

10.1 Electronic motor 21

Use the "Speed ▲" key to set the speed of the electronic motor to its maximum. With the 1:1 Contra Angle attached, it is 40'000 rpm. With all the other Contra Angles the maximum speed is lower, according to the transmission ratio of the Contra Angle. Press the pedal treadplate; the electronic motor starts running and accelerates to 40'000 rpm. When the treadplate is released, the electronic motor slows down again.

- The electronic motor 21 is designed for an intermittent duty operation at maximum speed of "1 min ON/3 min OFF" at 4 cycles. 15 min brake. Otherwise the electronic motor can be damaged by excessive heat generation, and by touching the motor hull it could cause serious burnings.
- The motor ventilation slots must be kept clear in order to prevent the motor temperature from becoming excessive.

10.2 Pump

Press the key "Pump" on the pedal briefly; the peristaltic pump is switched on, which is shown on the display by the symbol of a drop. Press the pedal treadplate; the peristaltic pump (built-in cooling pump) and the electronic motor start up. Water sprays from the irrigation needle on the contra angle.

10.3 Rotation direction of the electronic motor 21

Press the key "FOR/REV" on the pedal briefly, the rotation direction of the electronic motor changes. Press the pedal treadplate; the electronic motor rotates to the left and a continuous tone is emitted. Release the treadplate; the electronic motor ceases to operate and the tone is no longer heard. By pressing the motor key again, the rotation direction is switched back to right rotation, which is shown on the display by the symbol of a changing arrow.

5. Startup of the Dental surgery portable motor device models: MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories

Pic. 16. The Dental surgery portable motor device MD 30 model.

1 Device setup

- Place the Dental surgery portable motor device models: MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories and all required accessories and instruments on an even, non slip surface and make sure you have good access to all controls.
- Do not allow the operating range of the device (including the cable and contra angle) to be compromised by limiting factors.
- The system display must be fully visible at all times.
- The pedal must be placed within stepping distance between the patient and the surgeon.
- It must be explicitly ensured that no objects can fall on the pedal
- The power plug at the rear of the device must be accessible at all times.

2 Connecting to the power supply

Before plugging the power cable into the power socket for the first time, you must check the supply voltage setting next to the power switch.

If the voltage shown does not correspond to the local mains voltage, the grey fuse holder must be set to the correct voltage:


- A) Unplug the power cable.
- B) a screwdriver to open the fuse slot.
- C) Remove the fuse holder.
- D) Remove the grey fuse holder and reinsert it so that the local mains voltage setting is shown in the small window.
- E) Slide the fuse holder back in and close the fuse slot.
- F) Check the mains voltage shown on the fuse slot.
- G) Plug the power cable back into the device.

In order to prevent the risk of electric shock, the device may only be connected to a power network with a PE protective ground conductor.


3 Device preparation

1. Sterilize the motor (the motor is not sterile on delivery). If the motor has already been sterilized: when removing the motor from the sterile packaging, ensure that the sterile packaging is not damaged and that the sterility indicator confirms sterility (if no sterility indicator is provided, the sterile packaging must at least show the date on which the shelf life of the sterile item is due to expire). 
2. Insert the stand for the irrigation fluid into the stand holder.


Pic. 17. Insertion of the stand for irrigation fluid bottle

3. Plug the motor plug of the electronic motor into the motor socket.

Picture 18. Two 21 electronic motor sockets on the control unit

4. If required, plug in the second electronic motor into extra socket on the control unit.
5. Plug the pedal plug into the pedal socket at the rear of the control unit.

Pic. 19. Pedal socket on the rear panel of the control unit

Pic. 20. Device with plugged in the electric supply and pedal cables

6. Attach the sterilized handpiece to the electronic motor. Press the handpiece firmly onto the electronic motor unit it clicks into position and make sure it is secure by moving it slightly in the opposite direction.

Picture 21. The handpiece attached to the electronic motor resting on the handpiece cradle.

7. Assemble the tubing set. Attach tubing set Ref. 1706 2m, sterile (Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010) for cooling of the contra angle. For cooling two handpieces at the same time, use Ref. 6025 tubing set (purchased separately).


Check the expiry date of the tubing set and ensure that the packaging is not damaged. Using non-sterile tubing set can result in serious infection and, in extreme cases, can be fatal.

When inserting the tubing set, watch the arrow on the cover of the pump compartment. It indicates the flow direction of the cooling liquid.

Do not regulate the amount of irrigation liquid using the roller clamp on the tubing set, with the MD 30 model, this is regulated instead using the integrated pump. For this reason, make sure to open the roller clamp as far as it will go (refer to “Setting the pump supply quantity”).

- A) Press release key for set bracket (on top of the control unit) to open the pump.
- B) The compartment with the integrated tubing bracket opens.
- C) Place the tubing set into the tubing bracket provided in such a way that the part of the tubing set with the spike exits the pump towards the rear of the device. Make sure the tubing is secure
- D) With the tubing set inserted, press the compartment downwards until it clicks into place.


- E) Insert the spike at the end of the tubing set into the irrigation fluid bottle and hang the bottle onto the stand.
- F) Open the roller clamp on the tubing set as far as it will go.
- G) Open the bleed valve beneath the drip chamber.

H) Connect the control unit to the power socket.


Ensure that the operating voltage setting corresponds to the local mains voltage.

4 Assembly of external irrigation system

- A). Secure motor cable clip (Ref. 1873) and handpiece clip (Ref. 1881) to the irrigation tube.
- B). Connect the 6cm section of tube (Ref. 1773) for the Y-connector with the external cooling tube of the contra angle.
- C). Secure clips to the contra angle (Ref. 1881) and to the motor (Ref. 1873) cable.
- D). If required, secure additional clips (Ref. 1873) to the motor cable.

5. Assembly of internal and external irrigation system

A). Graft a 6 cm tube piece (Ref. 1773) on each of the branching tubes of the Y-connector (Ref. 1777).

B). Connect one of the 6 cm tube pieces that branch out of the Y-connector with the handpiece internal cooling tube (Ref. 1712).

C). Secure clips (Ref. 1873) and clips (Ref. 1881) to the tube set.


D). Secure clips (Ref. 1881) to the contra angle and clips (Ref. 1873) to the motor cable.

E). Put the internal cooling tube with attached tube piece into the hole in the head of the contra angle.


F). Insert the other end of the internal cooling tube into the hole in the hole of the handpiece.

G). Connect the 6 cm tube piece for Y-connector (Ref. 1773) with the tube of the external cooling.

H). If needed, secure more clips to the cable of the motor or the contra angle.


Pic. 22. How to insert coolant tube into the hole on the handpiece head.

Pic. 23. Assembling cooling systems of LED contra angle (top) and contra angle with push button system (bottom).

6. Operation. The Dental surgery portable motor device models: MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories

1 Switching the device on and off


Pic. 24. Main switch on the rear panel.

The power switch I/O (at the rear) is used to switch the control unit on and off.

The device can be switched off at any time irrespective of any procedure for switching off the device.

2 Overview: control elements on operating panel

- A) Display: Shows the operating values (See chapter "Overview: Standard display").
- B) "Program" key: Selection of program 1 to 10. There are 10 selectable programs for each motor. (To facilitate the selection of programs, the number of selectable programs can be limited to the ones you need. Refer to "Configuration Menu", Level 1, Programs/Number of programs).
- C) "Motor 1/2" key: Switching between the connected motors (the green indicator lights beside the motor sockets indicate the active motor). Pressing the Motor key longer changes direction of rotation.
- D) "Handpiece" key: Choice of handpieces or contra angle (can be individually deactivated, refer to "Configuration Menu").
- E) "Enter" key: Enters the configuration menu by pressing long (refer to "Configuration Menu").
- F) "Select" keys:
- By pressing the left "Select" key "Select ▲" the software version is shown.
 - By pressing both "Select" keys at the same time "Select ▲ + ▼" the programs will be reset to factory default settings.
 - The "Select" keys serve in the configuration menu to change values and parameters
 "▲": value adjustment (upwards) "▼": value adjustment (downwards)
- G) "Speed" keys:
- Restrict the maximum speed that can be selected using the pedal.
 "▲": increases the maximum speed "▼": reduces the maximum speed
 - By pressing both "Speed" keys at the same time "Speed ▲ + ▼" the handpiece calibration will be started (refer to "Calibrating handpieces").
- H) "Torque" keys:
- Restricts the maximum torque.
 "▲": increases the maximum torque "▼": reduces the maximum torque
 - By pressing both "Torque" keys at the same time "Torque ▲ + ▼" the torque modes AL and AS are exchanged (see chapter "Torque limit function, AL mode (Automatic Limiter)" and chapter "Torque limit function, AS mode (Automatic Stopper)").
- I) "Pump" keys:
- Adjusting the maximum pump (built-in cooling pump) flow rate, controlled by the pedal.
 "▲" increases the maximum supply quantity "▼" reduces the maximum supply quantity

- By pressing both “Pump” keys at the same time “Pump ▲ + ▼” the pump (built-in cooling pump) will be put on call, pressing again will switch it off.

3 Overview: Standard display


- A) Information line
Information and error messages are displayed here.
- B) Program
Shows the selected program for the active motor.
- C) Pump
The numerical value shows the pump (built-in cooling pump) flow rate in percent and the drop symbol together with the on/off indication shows if the pump (built-in cooling pump) is in standby mode or switched off.
- D) Rotation direction of the motor
The arrow indicates the rotation direction set for the motor. The rotation direction can be changed by pressing the button  on the pedal or by long pressing of the “Motor 1/2” key at the operating panel.
- E) AS/AL
Indication of the Auto Stopper (AS) mode or the Automatic Limiter (AL) mode (see chapter “Torque limit function, AL mode” and chapter “Torque limit function, AS mode”).
- F) Current speed
Shows current revolutions per minute of the motor.
- G) Clock
- H) Date
- I) Motor Shows selected motor → see also green indicator light beside the motor socket.
- J) Name of the handpiece or corresponding transmission ratio.
Shows name of the handpiece used or the selected transmission ratio (see also “Selecting the contra angle or transmission ratio”).
- K) Speed range
Shows speed range of the handpiece used.
- L) Maximum torque

Shows maximum torque setting.

M) Current torque

Bar graph providing a graphical representation of the current torque. All bars active means max. torque reached.

The pump does not begin to operate until the motor has been activated by pressing the pedal


4 Adjusting the programs

The values for operation settings depend on the connected handpiece/ contra angle as well as the task to be performed.

4.1 Step 1: Selecting the handpiece or transmission ratio

Depending on the handpiece or contra angle attached to the motor, the corresponding transmission ratio has to be adapted accordingly by using the handpiece key (HP).

Press the "Handpiece" key repeatedly until the name of the required handpiece or contra angle with the corresponding transmission ratio is shown on the display. When the key is pressed constantly the handpiece or contra angles will be shown in fast forward mode.

Table 15. Appropriate handpiece/ contra angles to be used

Name of the handpiece/ contra angles with transmission ratio	Display	Speed min. rpm	Speed max. rpm	Torque min. Ncm	Torque max. Ncm	AS-range (Factory def.) rpm	Limit AS- range rpm
Drill contra angle, 16:1 (not included in delivery)	1:5	1500	240,000	1	1	No	-
Drill contra angle, 16:1 (not included in delivery)	1:3	900	150,000	1	2	No	-
Drill contra angle, 16:1 (not included in delivery)	1:2	600	100,000	1	2	No	-
Drill contra angle, 16:1 (not included in delivery)	1:1	300	50,000	1	6	No	-
Drill contra angle, 16:1	4:1	75	12,000	1	20	up to 75	75-80*
Drill contra angle, 16:1 (not included in delivery)	16:1	20	2800	5	27	up to 20	20-45*
Drill contra angle, 16:1	20: 1	15	2100	10	70	up to 20	15 - 45 *
Drill contra angle, 16:1 (not included in delivery)	32: 1	10	1200	10	55	up to 20	10 - 45 *
Drill contra angle, 16:1 (not included in delivery)	70: 1	5	700	10	55	up to 20	5 - 45 *
MSS 5000, Micro Compass Saw (not included in delivery)	Com. Saw	fixed 15,000		fixed 6		No	-
OMS 5000, Micro Oscillating Saw (not included in delivery)	Osc.. Saw	fixed 15,000		fixed 6		No	-
MOS 5000, Micro Sagittal Saw (not included in delivery)	Sag. Saw	fixed 15,000		fixed 6		No	-
Mucotome (not included in delivery)	Mucotome	fixed 8,000		fixed 6		No	-
Kirschner handpiece (not included in delivery)	Kirschner	500	2800	fixed, 6		No	-

* This parameter can be adjusted in the Configuration Menu.

4.2 Step 2: Calibrating handpiece

To make sure the displayed parameters correspond with the actual, measurable parameters of the handpiece or contra angle, it is recommended to calibrate each handpiece or contra angle on a regular basis.

It is a procedure as simple as it is important to guarantee safety and precision of each handpiece or contra angle being used.

After you take care of all prior preparations such as sterilization, maintenance and care of handpiece, device preparation and the selection of the handpiece of use, the calibration procedure can be performed.

The calibration of handpieces or contra angles guarantees accurate torque. Due to wear as well as varying lubrication of the handpiece or contra angles and lack of maintenance and care, the distribution of torque can vary widely.

1. Hold motor with mounted handpiece or contra angle in your hand at a safe distance from your body.
2. Press both "Speed" keys at the same time (Speed ▲ + Speed ▼).
3. Motor and handpiece start running and pass several speed cycles up to maximum speed.
4. After a tone emitted and a message displayed the calibration is finished.


4.3 Step 3: Setting the speed

The possible speed range depends on contra angle to be used. The maximum speed within this speed range can be restricted to the required value.

Using the pedal, the speed can be varied from the minimum speed up to the maximum speed as set.

Setting the speed:

Press the "Speed" keys "▲" to increase or "▼" to decrease the maximum speed. When key is pressed continuously the speed will be shown in fast forward mode.

The following handpieces run only with one specific speed, it cannot be changed. The specific rpm values of these handpieces are shown in Table 15.

- MSS 5000, Micro Compass Saw (not included in delivery)
- OMS 5000, Micro Oscillating Saw (not included in delivery)
- MOS 5000, Micro Sagittal Saw (not included in delivery)
- Mucotome (not included in delivery)

4.4 Step 4: Setting the torque

Once the speed has been selected, the torque can be limited from the torque range of the corresponding handpiece or contra angle. AS. Depending on the speed, the torque mode AL (Automatic Limited) or AS (Automatic Stopper) is applied.

For information on AL and AS made, see chapter “Speed reduction AL (Automatic Limited)” and chapter “Torque limit function AS (Automatic Stopper)”.

Press the “**Torque**” keys “▲” to increase or “▼” to decrease the maximum torque. By pressing the key continuously the torques are changed in fast forward mode.


The following handpieces run on fixed torque only (6Ncm) and it cannot be changed:

- MSS 5000, Micro Compass Saw (not included in delivery)
- OMS 5000, Micro Oscillating Saw (not included in delivery)
- MOS 5000, Micro Sagittal Saw (not included in delivery)
- Mucotome (not included in delivery)
- Kirschner handpiece (not included in delivery)

4.5 Step 5: Setting the pump supply quantity

Press the “**Pump**” keys “▲” to increase or “▼” to decrease the pump supply quantity. By pressing the key continuously the values are fast forwarded.

Minimum and maximum of pump (built-in cooling pump) supply quantity as well as incremental steps can be adjusted in the Configuration Menu (refer to “Configuration Menu, Parameter Level 2, Pump”).


To activate or deactivate the pump, press both “**Pump**” keys at the same time, “**Pump** ▲ + ▼” or use the footswitch .

5 Torque limit function, AL mode (Automatic Limiter)

The AL function limits the torque applied to the instrument, for example to prevent crack initiation and bone fracture.

The speed of the instrument remains constant until the preset torque is reached. If the applied force gets over this limit the speed will be reduced, if necessary down to a stop, but the torque remains constant. If the applied force is reduced the speed picks up again.

On the display this procedure is shown as a bar graph. The bars in the graph fill up to their full range as the torque increases to its preset value. When the bar graph shows full capacity the speed will be reduced. As soon as the force on the instrument is reduced the torque decreases and the speed starts to pick up again as shown on the display.

The AL mode is active at all speeds, except at the speeds where the AS mode takes effect.

6 Torque limit function, AS mode (Automatic Stopper)

The AS function limits the torque applied to the device. As soon as the preconfigured torque is reached, the electric motor stops immediately. The electronic motor no longer generates any force. In order to restart the electronic motor, the pedal must be released and pressed again.

On the display, the bar graph completes up to its full range, until the maximum torque is reached, then drops to zero.

The function “AS” is just active for certain contra angles and only within a certain speed range.

Handpiece	4:1	16: 1	20: 1	32: 1	70: 1
Speed rpm	75	20	20	20	20

From the minimal speed up to the speeds shown at the right the AS mode is activated automatically (refer to “Configuration Menu, Parameter Level 2, AS-Zone”).

The upper limit of the AS range can be adjusted in the configuration menu.

In the following speed ranges the AS-/AL- mode can also be switched manually by pressing both torque keys “Torque ▲ + ▼” at the same time.

Handpiece	4: 1	16: 1	20: 1	32: 1	70: 1
From rpm:	75	20	15	10	5
To rpm:	80	45	45	45	45

7 Storing various programs

With the Dental surgery portable motor device models: MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories up to 10 different settings can be set as fixed program (program 1 to program 10). Which program is currently active is shown on the display.

At switch-off, the settings made by user are automatically saved. This includes the following parameters:

- | | |
|--|---|
| <ul style="list-style-type: none"> • handpiece/ transmission ratio • maximal speed • maximal torque | <ul style="list-style-type: none"> • pump On/Off • pump performance • AS/AL mode |
|--|---|

In the configuration menu the number of available programs can be set from 3 to 10.

When the MD 30 model is switched on, the display shows the prior used program and motor.

8 Configuration menu

In the configuration menu the user can customize the device for one's purposes. The parameters are organized in two levels. The following information and parameters are stored in there.

- | | |
|---|---|
| <ul style="list-style-type: none"> • Software version • Mainboard serial number • Date and Time • Display illumination • Number of programs • Operating hours MD 30 • Operating hours Motor 1 • Operating hours Motor 2 • Operating hours of irrigation pump | <ul style="list-style-type: none"> • Error log (the last 8) • Activation of available handpiece • Speed limitation for each handpiece • Range of action of AS mode • Pump behavior • Motor behavior for 4 different motors • Behavior of torque at reverse mode • Reset to factory default settings |
|---|---|

Be cautious when changing parameters. Unusual behavior of instruments while operating may provoke false reaction and jeopardize the patient. Every setting and the new behavior of the instrument has to be verified and tested.

Access to configuration menu:

- Press “Enter” key for 3 seconds until hear a long tone. The right arrow “>” in the information line indicates that you are in the Configuration Menu.
- With “Select ▲” or “Select ▼” choose the desired parameter.
 - Push “Enter” to go to the sub menu.
 - To change the value press “Select ▲” or “Select ▼”.
 - To confirm the changes made, press “Enter” for one second, until you hear a short tone.
 - To abort the settings made press “Enter” briefly, the returns to its previous value.

Exit from Configuration Menu:

- To leave the Configuration Menu press “Enter” for 3 seconds, until you hear a long tone.

Table 16. Parameter Level 1

Group/Parameter	Rights	Factory	Definition
Software/ Version	read	VX.XXXX	Shows current software version
Hardware/ Serial number MB	read	XXXXXXXXXXXX	Shows serial number of the main board
Date-Time/ Date format US	read/write	No	Set to US Date format
Date-Time/ Date	read/write	-	Change current date
Date-Time/ Time	read/write	-	Change current time
Backlight/brightness (0 - 10)	read/write	9	Display brightness, changeable :0-10
Programs/ Number of programs	read/write	10	Number of activated programs 3-10
Operating hours/ MD 30	read	0	Shows operating hours of MD 30
Operating hours Motor 1	read	0	Shows operating hours of motor 1
Operating hours Motor 2	read	0	Shows operating hours of motor 2
Error memory/ 1-8	read	0	8 Error messages in chronological order

Table 17. Parameter Level 2

To adjust the values in Level 2 the password “9403” must be entered. *The password can not be changed.*
 Entering password: press “Select ▲ or ▼” (for fast forward or backward keep the **key** pressed).

Handpiece activation	Name of handpiece on display	Choice	Factory default	Definition
Handpiece existing/ HP 01 (not included in delivery)	1:5	yes/no	no	Deselect handpieces that do not belong to your product range by switching to “no”. This will later shorten the scrolling list for “Handpieces”. Otherwise you will later have to scroll down the complete list of available handpieces, each time time you select a handpiece. <i>*Handpiece 1:1 cannot be deactivated.</i>
Handpiece existing/ HP 02 (not included in delivery)	1:3	yes/no	no	
Handpiece existing/ HP 03 (not included in delivery)	1:2	yes/no	no	
Handpiece existing/ HP 04 (not included in delivery)	1:1	yes/no	yes	
Handpiece existing/ HP 05 (not included in delivery)	4:1	yes/no	yes	
Handpiece existing/ HP 02 (not included in delivery)	16:1	yes/no	yes	
Handpiece existing/ HP 07	20:1	yes/no	yes	
Handpiece existing/ HP 08	32:1	yes/no	yes	
Handpiece existing/ HP 09 (not included in delivery)	70:1	yes/no	no	
Handpiece existing/ HP 10 (not included in delivery)	Com. Saw	yes/no	yes	
Handpiece existing/ HP 11 (not included in delivery)	Osc. Saw	yes/no	yes	
Handpiece existing/ HP 12 (not included in delivery)	Sag. Saw	yes/no	yes	
Handpiece existing/ HP 13 (not included in delivery)	Mucotome	yes/no	yes	
Handpiece existing/ HP 14 (not included in delivery)	Kirschner	yes/no	yes	

Handpiece Max. Speed	Name of handpiece on the display	Speed range rpm	Factory default	Definition
Handpiece max speed/ HP 01 (not included in delivery)	1:5	1500-240'000	240'000	Limit the maximum speed of your handpieces according to your own experience.
Handpiece max speed/ HP 02 (not included in delivery)	1:3	900-150'000	150'000	
Handpiece max speed/ HP 03 (not included in delivery)	1:2	600-100'000	100'000	
Handpiece max speed/ HP 04 (not included in delivery)	1:1	300-50'000	50'000	
Handpiece max speed/ HP 05 (not included in delivery)	4:1	75-12'000	12'000	

Handpiece max speed/ HP 06 (not included in delivery)	16:1	20-2'800	2'800	
Handpiece max speed/ HP 07	20:1	15-2'100	2'100	
Handpiece max speed/ HP 08	32:1	10-1'200	1'200	
Handpiece max speed/ HP 09 (not included in delivery)	70:1	5-700	700	
Handpiece max speed/ HP 10 (purchased separately)	Com. Saw	fix 15'000	15'000	
Handpiece max speed/ HP 11 (not included in delivery)	Osc. Saw	fix 15'000	15'000	
Handpiece max speed/ HP 12 (not included in delivery)	Sag. Saw	fix 15'000	15'000	
Handpiece max speed/ HP 13 (purchased separately)	Mucotome	fix. 8'000	8'000	
Handpiece max speed/ HP 14 (not included in delivery)	Kirschner	500-2'800	2'800	

AS-Zone for handpieces	Handpiece on display	Selection ..to ..rpm	Factory default	Definition
Handpiece AS-Mode/HP 05 (not included in delivery)	4:1	75-80	75	Effective AS range 75-80 rpm
Handpiece AS-Mode/HP 06 (not included in delivery)	16:1	20-45	20	Effective AS range 20-45 rpm
Handpiece AS-Mode/HP 07	20:1	15 - 45	20	Effective AS range 15-45 rpm
Handpiece AS-Mode/HP 08	32:1	10 - 45	20	Effective AS range 10-45 rpm
Handpiece AS-Mode/HP 09 (not included in delivery)	70:1	5-45	20	Effective AS range 5-45 rpm

Pump parameters	Range	Factory defaults	Definition
Pump/Backwards turn mode variable	yes/no	yes	The pressures in the tube vary according to the pump speed. In "variable mode" this is considered, to prevent spilling when pump is switched off in backwards mode.
Pump/Way backwards	1 – 100 %	25 %	Specify how far the pump turns backwards.
Pump/Speed backwards	10 – 50 %	33 %	Specify how fast the pump has to turn backwards to prevent spilling after the switch off of the motor.
Pump/Range 1 increment	1 – 10 %	5 %	Adjustment steps in section 1.
Pump/Range 1 end	5 – 50 %	10 %	Set the range where section 1 is active.
Pump/Range 2 increment	1 – 10 %	5 %	Adjustment steps in section 2.

Pump/Range 2 end	10 – 90 %	50 %	Set the range where section 2 is active.
Pump/Range 3 increment	1 – 10 %	10 %	Adjustment steps in section 3.
Pump/Range 3 end	20 – 100 %	100 %	Set the range where section 3 is active.

The Dental surgery portable motor device models: MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories can recognize the type of plugged-in motor. This enables the adaption of future motors and their safe operation.

Motor type 1	Range	Factory default	Definition
Motor type 1/Min. speed	300 – 5'000 rpm	300 rpm	Set the min. speed of the motor 2.
Motor type 1/Max. speed	5 000 – 50'000 rpm	50'000 rpm	Set the max. speed of the motor 2.
Motor type 1/ Start ramp	1 – 1'000 ms / 10'000 rpm	100 ms	Set acceleration time to 10'000 rpm.
Motor type 1/ Stop ramp	1 – 1'000 ms / 10'000 rpm	50 ms	Set breaking time from 10'000-0 rpm.

The maximum torque of the motor is boosted at starts and in reverse rotation.

Reverse torque	Range	Factory default	Definition
Reverse torque/Increase	5 – 30 %	25 %	Increase of the selected torque in reverse rotation.
Reverse torque/Increase time	100 – 2'000 ms	500 ms	Time during which the torque is raised.

Resetting to factory default	Choice	Factory default	Definition
Default value/ Set default value	yes/ no	no-	Resetting all parameters of the Configuration Menu to factory default.

- Attention: Upon resetting of factory default all parameters will appear with factory default values (except date and time).
- Note: For resetting the “programs” to factory default press both “Select” keys at the same time (“Select ▲ + ▼”). For this procedure you have to be out of the Configuration Menu.

a. Operation using the Vario pedal

1. Carrying bracket:

The carrying bracket can be operated by foot.

2. Key :

Pressing the key briefly: switches the pump on or off (see information on display).

Pressing the key longer, increases the pump speed (see information on display).

3. Key :

Pressing the key briefly: switches the program (+1, see information on display).

Pressing the key longer: switches the program (-1, see information on display).

4. Key :

Pressing the key briefly: switches the rotation direction (see information on display).

Pressing the key longer: switches the motor (see also the indicator lights beside the motor socket on the control unit).

5. Treadplate:

With the treadplate pedal the motor speed is variably adjusted and pump is activated.

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Treadplate...	Motor...	Pump...
... not pressed	Motor off	Pump off
... pressed gently	Motor runs slowly	Pump on, if pump "On" displayed (speed set on the control unit)
... pressed all the way down	Motor runs at maximum speed (speed as set on the control unit)	Pump on, if pump "On" displayed (speed set on the control unit)

For safety reasons, the unit can only be operated by pedal.

The following handpieces run only at one specific speed, it cannot be changed.

- MSS 5000, Micro Compass Saw (not included in delivery)
- OMS 5000, Micro Oscillating Saw (not included in delivery)
- MOS 5000, Micro Sagittal Saw (not included in delivery)
- Muotome (not included in delivery)

10 Functional check

Prior to the Dental surgery portable motor device models: MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories startup or use of accessory equipment, the user must always ensure that each individual component is in good working order, free from defects and in clean, sterile and operational. All inscriptions on the device and its accessories must be readable and there must be no loose parts on the device. Once the device is switched on, the most recent settings entered appear on the display and the green LED for motor 1 lights up.

10.1 Electronic motor 21

Use the "Speed" selection key to set the speed of the electronic motor to 40'000 or 50'000 rpm (depends on the version) . Press the pedal treadplate; the electronic motor starts running and accelerates to 40'000 or 50'000 rpm. When the treadplate is released, the electronic motor slows down again

- The Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories is designed for intermittent duty operation at maximum speed of "1 min ON/ 3 min OFF" at 4 cycles. Break after 15 min.
- The motor ventilation slots must be kept clear in order to prevent the motor temperature from becoming excessive.

10.2 Pump

Press the key "💧" on the pedal briefly, the peristaltic pump (built-in cooling pump) is switched on, which is shown on the display by the symbol of a drop. Press the pedal treadplate; the peristaltic pump and the electronic motor start up. Water sprays from the irrigation needle on the contra angle.

10.3 Rotation direction of the electronic motor 21

Press the key "Ⓜ" on the pedal briefly, the rotation direction of the electronic motor changes. Press the pedal treadplate; the electronic motor rotates to the left and a continuous tone is emitted. Release the treadplate; the electronic motor ceases to operate and the tone is no longer heard. By pressing the motor key

again, the rotation direction is switched back to right rotation, which is shown on the display by the symbol of a changing arrow.

10.4 Program

The required program can be set by repeatedly pressing the key “ P ” on the pedal.

7. HOW TO USE HANDPIECES AND CONTRA ANGLES

HANDPIECES AND CONTRA ANGLES

The handpieces used with the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories may have straight or angular design, may be equipped with torque increasing/decreasing system and various heads suitable for various functions as well as the coolant supply system. They may have an external or an internal coolant supply system as well as fiber-lit lighting.

Non-working part of the handpiece or the contra angle is marked using an electrical discharge machine ensuring marking legibility. The marking stays lasting and legible during operation, storage and transportation. The marking of the handpieces and contra angles includes name of the manufacturer (NOUVAG AG), serial number (SN) and reference number (REF), model type (CA 20:1 LED) and the Steam Sterilization symbol.

Picture. 25. Handpiece and contra angle marking.

Precautions while operating the handpieces and contra angles

- The contra angles are not sterile on delivery. Before first use and immediately after every operation they must be cleaned, disinfected and sterilized!
- Do not clean the handpiece and the contra angle with compressed air!
- Speed limit for the handpiece and contra angle is 50.000 revolutions per minute!
- Handpieces and contra angles may only be attached when the electronic motor is not running!
- Never operate the clamping mechanism of the handpieces and contra angles while the electronic

motor is running! This could result in instrument damage.

- The handpieces and contra angles may only be operated by qualified and trained personnel!
- Improper use or repair of the device and failure to observe these instructions relieve NOUVAG AG from any obligations arising from warranty provisions or any other claims.

Explanation of symbols

	CE mark, confirming CE certificate being available, with reference number of the notified body		Autoclavable at max. 135 °C		Suitable for thermal disinfection		Date of expiry
	Warning		Manufacturer		Observe the instructions for use		Symbol indicating the serial number
	Symbol indicating the order number						

Contra angles

Technical data

Transmission ratio	20:1 (5053nou)	32:1 (5201nou)	20:1 (5052nou LED)
Conformity standard	ISO 3964	ISO 3964	ISO 3964
File - Ø, mm	2,35	2,35	2,35
Max. speed, revolutions/minute	2 100	1 200	2 100
Max. torque, Nm	70	55	70
File length, mm	75	21	75
Total length of the instrument, mm	116	112	100
Weight, g	90	90	100
Internal/external irrigation system	• / •	• / •	• / •


Internal cooling pipe

USE

Mounting / Replacement


Using handpieces with a button clamp press the button on the head of the handpiece and insert the instrument (file) into the head hole until it clicks into position. Release the button. The instrument is attached. Make sure it is secure.

Pic. 26.

Connection to the electronic motor

Attach the handpiece to the electronic motor. Press the handpiece firmly onto the electronic motor until it clicks into position and make sure it is secure by moving it slightly in the opposite direction. Reliable connection of the electronic motor and the handpiece means that you cannot remove the handpiece from the electronic motor by pulling the handpiece forward, while the handpiece can be rotated easily (360 degrees).

Picture 27.

Assembly of the tube set

For handpieces with external or internal irrigation system, you must connect cooling pipes and attach them to the handpiece and cable of the electronic motor using special clips before operation.

Assembly of external irrigation system:

A. Connect the open end (1706) of the 2 m long, sterile irrigation tube set (Marketing Authorization № FSZ (ФСЗ) 2010/07320 dated 05.07.2010) with the cooling pipe of the contra angle.

B. Attach a clip (1881) to the tube set.

C. Attach the clip (1881) with the tube set to the contra angle.

D. Connect the electronic motor with the contra angle.

E. Attach a clip (1873) to the tube set.

F. Attach the clip (1873) with the tube set to the cable of the electronic motor.

G. The Drill unit with the attached external cooling system is ready to use.

Assembly of internal and external irrigation system:

A. Attach the open end (1706) of the 2 m long, sterile irrigation tube set (Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010) to one of the branches of the Y-connector (1777).

B. Connect a piece of a 6 cm long connection tube for Y-connector (1773) to each of the branching tubes of the Y-connector.

C. Connect one piece of the 6 cm long connection tube for Y-connector, that branches out of the Y-connector, with the internal cooling pipe of the contra angle (39116).

D. Connect the second piece of the connection tube that branches out of the Y-connector, with the external cooling pipe of the contra angle (welded on laterally to the contra angle's head).

E. Both pieces of the 6 cm long connection tube for Y-connector are now connected to the cooling system of the contra angle.

F. Secure a tube clip (1881) to each of the short tube pieces.

G. Secure the tube clip (1881) to the contra angle.



H. The contra angle with the attached tube clip (1881).

I. Fix a clip (1873) to the tube set leading to the Y-connector.

J. Fix the tube set with the clip (1873) to the cable of the electronic motor.

K. The cable of the electronic motor with the attached tube set. Secure further clips to the motor cable if needed.

L. The tube set routing with the Y-connector and external cooling of the instrument.

Cooling system being attached to the contra angle and the motor

1777
Y-connector

1873
Motor cable
clips

1773
16 cm long
connection tubes


Preparation instructions

Reprocessing restrictions	Frequent reprocessing has only a limited impact on the life period of the handpiece. The end of the product service life is normally determined by wear and damage through use. The handpiece is designed for 250 sterilization cycles.
----------------------------------	---

INSTRUCTIONS

At the location of use	Remove surface soiling with a disposable cloth / tissue paper.
Storage and transportation	No special requirements. Due to the risk of blood drying and further corrosion reprocessing must be performed without undue delay (longer than 30 minutes).
Preparation for cleaning	Remove the instrument, the cooling pipe, the fixating clips and release the handpiece.
Automatic cleaning and disinfection	Equipment: Washer-disinfector with a special load carrier that ensures handpieces are connected to the washer-disinfector (for example, Miele G7835 CD Thermal Disinfector) and channels are rinsed. For small parts use a fine-mesh basket. Rinse them from every side. Use neutral or alkaline cleaning agents (for example, Schülke perform or 3% DESOMEDAN ID) in the manufacturer's recommended concentration. 1. Place the handpiece in the load carrier (ensure that the channels can be rinsed). 2. Place the small parts in the fine-mesh basket and close it. 3. Set a cleaning cycle that offers sufficient cleaning and rinsing. Perform the final rinse with demineralized water. 4. Perform a 10-minutes rinse cycle at 93°C to facilitate thermal disinfection. 4. When removing the handpiece, inspect it visually to verify whether soiling is still visible in the gaps and the grooves. If necessary, repeat the cycle or clean manually.
Manual cleaning	Equipment: Neutral or alkaline cleaning agent (for example, Schülke perform or 3% DESOMEDAN ID) in the manufacturer's recommended concentration, a soft brush, running water. <i>Procedure:</i> 1. Rinse off surface soiling on the handpiece and small parts under running water. 2. Use a brush to apply the cleaning agent to all surfaces and gaps. 3. Rinse the handpiece thoroughly under running water.
Manual disinfection	For manual disinfection, wipe the handpiece and the small parts with a cloth soaked in an ethanol-based disinfecting agent.
Drying	If a drying program is not provided by the washer-disinfector, the handpiece and the small parts must be dried manually or in a drying cabinet until there is no visible moisture left on them.
Inspection and maintenance  № 1958	Perform a visual inspection to check for damage, corrosion and wear. After disinfection thoroughly and briefly spray the handpiece with NouClean Spray (optional) and wipe it with a moist cloth (see the instructions on the spray can). 
Packaging	Individual: Pack the handpiece in individual packaging for sterile items. The bag must be large enough to ensure that the seal is not subject to strain. Sets: Place the handpieces in the load carrier or on all-purpose sterilization trays.
Sterilization	Autoclave in the vacuum autoclave at 135° C for at least 5 minutes*. When sterilizing several instruments during one sterilization cycle, do not exceed the maximum sterilizer load. A drying cycle must be added in the case of autoclaves without a post-vacuum function. Allow the handpieces to dry in the bags with the paper side facing upwards. * Temperature exposure times are based on country-specific guidelines and standard.
Storage	No special requirements. If the sterilized handpiece is not used immediately after

sterilization, the material packaging must be labeled with the sterilization date. Nouvag AG recommends to use a sterility indicator on the packaging.

In case of deterioration of LED contra angle lightning you shall clean the optical port with an alcohol-soaked cotton swab. Clean all the soiling and dirt. Do not use any sharp instruments to clean the optical port!

Malfunctions and troubleshooting

Malfunction	Cause	Solution
The electronic motor is running but the instrument doesn't switch on	The handpiece is not correctly connected with the electronic motor	Press the handpiece until it clicks into position and make sure it is secure.
Incorrect operation of the instrument	Incorrect attachment of the instrument to the contra angle	Fix the instrument and tighten the clamp
The handpiece causes excessive noise during operation	Dirty or badly oiled ball bearings	Spray the insides of the contra angle with Nou-Clean spray (not included into delivery)


8. HOW TO USE A MOBILE DEVICE TROLLEY

SANICAR is a mobile device trolley powered by caretaking hospital staff. It is designed for transporting or carrying this medical product. SANICAR trolley is equipped with three bottom trays to place 3 devices (picture 28) to perform surgery. The trolley comes with 4 rolls, from which two are lockable to prevent random movements of the trolley. A 5 m long cable is used to connect the trolley to the power network. There are three sockets on the back panel (picture 29) to connect the devices used.


Pic. 28. How to use SANICAR trolley.

..

Pic. 29. SANICAR Mobile Device Trolley.

Pic. 30. Labeling of the SANICAR Mobile Device Trolley.

The labeling includes the following data:

**Manufacturer:
NOUVAG AG
St. Gallerstrasse, 23-25
CH - 9403 GOLDACH**

Name and address of the manufacturer

Type: Sanicar fahrbar

Product description

Serie No. 01/15

Serial number

Power: ~230V 50Hz max. 16A

Power consumption

CE mark

Warning

Art. Nr. 1898d

Reference number №

9. SUITCASE

Picture 31. Suitcase

The suitcase is designed to house the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories. There is enough space for the complete set with handpieces and contra angles as well as for a spare electronic motor. The suitcase is good for both storage and transportation of the medical product as it ensures product protection against mechanical damages and environmental hazards. Your device will be protected by a rigid outer shell made from aluminum and foamed rubber at the inside. The ergonomic handle and reliable locks eliminate the possibility of suitcase self-opening. The suitcase is suitable for disinfection (for example, with an alcohol-based solution or disinfecting liquid).

9. SCOPE OF DELIVERY

I. Dental surgery portable motor device MD 11-1, including:

1.	3335	MD11 control unit with an integrated cooling pump system, single-motor,	1 unit;
2.	2097nou	Electronic motor 21, speed 40,000 revolutions/minute, including 2 m long motor cable, (if necessary)	2 units;
3.	1866nou	Vario IPX8 foot pedal for MD 11 with 3.15 m long cable,	1 unit;
4.	1706	Tube set , sterile, 2 m long (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
5.	1770	Stand for irrigation fluid bottle (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
6.	1873	Clip set (10 units) for tube set attachment at motor, (if necessary)	2 sets;
7.	1881	Clip set (3 units) for tube set attachment at handpieces and contra angles, (if necessary)	2 sets;
8.	1958	Spray nozzle for Nou-Clean spray, for care of electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
9.	1974	Spray nozzle for Nou-Clean spray, for care of instruments (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
10.	1170	Stand for electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
11.	28879	Fuse T 3.15 AL 250 V AC,	2 units;
12.	22262	Power cable, 3 m long,	1 unit;
13.	31686	User Manual on CD-ROM,	1 unit.

II. Dental surgery portable motor device MD 11-2, including:

1.	3335	MD11 control unit with an integrated cooling pump system, single-motor,	1 unit;
2.	2097nou	Electronic motor 21, speed 40,000 revolutions/minute, including 2 m long motor cable, (if necessary)	2 units;
3.	1866nou	Vario IPX8 foot pedal for MD 11 with 3.15 m long cable,	1 unit;
4.	1706	Tube set , sterile, 2 m long (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
5.	1770	Stand for irrigation fluid bottle (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
6.	1873	Clip set (10 units) for tube set attachment at motor, (if necessary)	2 sets;
7.	1881	Clip set (3 units) for tube set attachment at handpieces and contra angles, (if necessary)	2 sets;
8.	1958	Spray nozzle for Nou-Clean spray, for care of electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
9.	1974	Spray nozzle for Nou-Clean spray, for care of instruments (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
10.	1170	Stand for electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
11.	5201nou	Drill contra angle, 32:1, shank diameter of $\varnothing$ 2.35 mm, with internal/external cooling, and button clamp,	1 unit;

12.	28879	Fuse T 3.15 AL 250 V AC,	2 units;
13.	22262	Power cable, 3 m long,	1 unit;
14.	31686	User Manual on CD-ROM,	1 unit.

III. Dental surgery portable motor device MD 11-3, including:

1.	3335	MD11 control unit with an integrated cooling pump system, single-motor,	1 unit;
2.	2097nou	Electronic motor 21, speed 40,000 revolutions/minute, including 2 m long motor cable, (if necessary)	2 units;
3.	1866nou	Vario IPX8 foot pedal for MD 11 with 3.15 m long cable,	1 unit;
4.	1706	Tube set , sterile, 2 m long (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
5.	1770	Stand for irrigation fluid bottle (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
6.	1873	Clip set (10 units) for tube set attachment at motor, (if necessary)	2 sets;
7.	1881	Clip set (3 units) for tube set attachment at handpieces and contra angles, (if necessary)	2 sets;
8.	1958	Spray nozzle for Nou-Clean spray, for care of electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
9.	1974	Spray nozzle for Nou-Clean spray, for care of instruments (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
10.	1170	Stand for electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
11.	5053nou	Drill contra angle, 20:1, shank diameter of $\varnothing$ 2.35 mm, with internal/external cooling, and button clamp,	1 unit;
12.	28879	Fuse T 3.15 AL 250 V AC,	2 units;
13.	22262	Power cable, 3 m long,	1 unit;
14.	31686	User Manual on CD-ROM,	1 unit.

IV. Dental surgery portable motor device MD 11-4, including:

1.	3335	MD11 control unit with an integrated cooling pump system, single-motor,	1 unit;
2.	2097nou	Electronic motor 21, speed 40,000 revolutions/minute, including 2 m long motor cable, (if necessary)	2 units;
3.	1866nou	Vario IPX8 foot pedal for MD 11 with 3.15 m long cable,	1 unit;
4.	1706	Tube set , sterile, 2 m long (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
5.	1770	Stand for irrigation fluid bottle (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
6.	1873	Clip set (10 units) for tube set attachment at motor, (if necessary)	2 sets;
7.	1881	Clip set (3 units) for tube set attachment at handpieces and contra angles, (if necessary)	2 sets;
8.	1958	Spray nozzle for Nou-Clean spray, for care of electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
9.	1974	Spray nozzle for Nou-Clean spray, for care of instruments (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
10.	1170	Stand for electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;

11.	5052nou	LED Drill contra angle, 20:1, shank diameter of $\varnothing$ 2.35 mm, with internal/external cooling, and button clamp,	1 unit;
12.	28879	Fuse T 3.15 AL 250 V AC,	2 units;
13.	22262	Power cable, 3 m long,	1 unit;
14.	31686	User Manual on CD-ROM,	1 unit.

V. Dental surgery portable motor device MD 30-1, including:

1.	3330	MD30 control unit with an integrated cooling pump system, two-motor,	1 unit;
2.	2097nou	Electronic motor 21, speed 50,000 revolutions/minute, including 2 m long motor cable, (if necessary)	2 units;
3.	1510nou	Vario IPX8 foot pedal for MD 30 with 3.15 m long cable,	1 unit;
4.	1706	Tube set , sterile, 2 m long (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
5.	1770	Stand for irrigation fluid bottle (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
6.	1873	Clip set (10 units) for tube set attachment at motor, (if necessary)	2 sets;
7.	1881	Clip set (3 units) for tube set attachment at handpieces and contra angles, (if necessary)	2 sets;
8.	1170	Stand for electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
9.	1958	Spray nozzle for Nou-Clean spray, for care of electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
10.	1974	Spray nozzle for Nou-Clean spray, for care of instruments (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
11.	28879	Fuse T 3.15 AL 250 V AC,	2 units;
12.	22262	Power cable, 3 m long,	1 unit;
13.	31686	User Manual on CD-ROM,	1 unit.

VI. Dental surgery portable motor device MD 30-2, including:

1.	3330	MD30 control unit with an integrated cooling pump system, two-motor,	1 unit;
2.	2097nou	Electronic motor 21, speed 50,000 revolutions/minute, including 2 m long motor cable, (if necessary)	2 units;
3.	1510nou	Vario IPX8 foot pedal for MD 30 with 3.15 m long cable,	1 unit;
4.	1706	Tube set , sterile, 2 m long (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
5.	1770	Stand for irrigation fluid bottle (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
6.	1873	Clip set (10 units) for tube set attachment at motor, (if necessary)	2 sets;
7.	1881	Clip set (3 units) for tube set attachment at handpieces and contra angles, (if necessary)	2 sets;
8.	1170	Stand for electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
9.	5201nou	Drill contra angle, 32:1, shank diameter of $\varnothing$ 2.35 mm, with internal/external cooling, and button clamp,	1 unit;
10.	1958	Spray nozzle for Nou-Clean spray, for care of electronic motor (Marketing	1 unit;

	Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	
11. 1974	Spray nozzle for Nou-Clean spray, for care of instruments (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
12. 28879	Fuse T 3.15 AL 250 V AC,	2 units;
13. 22262	Power cable, 3 m long,	1 unit;
14. 31686	User Manual on CD-ROM,	1 unit.

VII. Dental surgery portable motor device MD 30-3, including:

1. 3330	MD30 control unit with an integrated cooling pump system, two-motor,	1 unit;
2. 2097nou	Electronic motor 21, speed 50,000 revolutions/minute, including 2 m long motor cable, (if necessary)	2 units;
3. 1510nou	Vario IPX8 foot pedal for MD 30 with 3.15 m long cable,	1 unit;
4. 1706	Tube set , sterile, 2 m long (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
5. 1770	Stand for irrigation fluid bottle (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
6. 1873	Clip set (10 units) for tube set attachment at motor, (if necessary)	2 sets;
7. 1881	Clip set (3 units) for tube set attachment at handpieces and contra angles, (if necessary)	2 sets;
8. 1170	Stand for electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
9. 5053nou	Drill contra angle, 20:1, shank diameter of $\varnothing$ 2.35 mm, with internal/external cooling, and button clamp ,	1 unit;
10. 1958	Spray nozzle for Nou-Clean spray, for care of electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
11. 1974	Spray nozzle for Nou-Clean spray, for care of instruments (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
12. 28879	Fuse T 3.15 AL 250 V AC,	2 units;
13. 22262	Power cable, 3 m long,	1 unit;
14. 31686	User Manual on CD-ROM,	1 unit.

VIII. Dental surgery portable motor device MD 30-4, including:

1. 3330	MD30 control unit with an integrated cooling pump system, two-motor,	1 unit;
2. 2097nou	Electronic motor 21, speed 50,000 revolutions/minute, including 2 m long motor cable, (if necessary)	2 units;
3. 1510nou	Vario IPX8 foot pedal for MD 30 with 3.15 m long cable,	1 unit;
4. 1706	Tube set , sterile, 2 m long (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
5. 1770	Stand for irrigation fluid bottle (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
6. 1873	Clip set (10 units) for tube set attachment at motor, (if necessary)	2 sets;
7. 1881	Clip set (3 units) for tube set attachment at handpieces and contra angles, (if necessary)	2 sets;
8. 1170	Stand for electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;


9.	5052nou	LED Drill contra angle, 20:1, shank diameter of $\varnothing$ 2.35 mm, with internal/external cooling, and button clamp ,	1 unit;
10.	1958	Spray nozzle for Nou-Clean spray, for care of electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
11.	1974	Spray nozzle for Nou-Clean spray, for care of instruments (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
12.	28879	Fuse T 3.15 AL 250 V AC,	2 units;
13.	22262	Power cable, 3 m long,	1 unit;
14.	31686	User Manual on CD-ROM,	1 unit.

In line with regulations pertaining to hazardous materials, 1984 Nou-Clean Maintenance and Care Spray is not delivered with the control unit and must be ordered separately. Delivered separately!

10. CLEANING OF DEVICE COMPONENTS.

Schematic flow chart for resterilization as per EN ISO 17664

- Reusable medical devices are subject to repeated sterilization and further storage under the conditions excluding secondary microbial contamination.
- All device components being in direct contact with patients' oral fluid, mucosa and oral tissues are subject to pre-sterilization treatment and sterilization. Other device components shall be disinfected without sterilization.
- The above medical devices shall be sterilized immediately before their first use as well as after every patient in order to avoid cross-contamination.
- The medical devices contacting patients' blood and serum shall be disassembled, cleaned and rinsed by medical staff wearing rubber gloves, protective goggles, a facemask, a gown and a waterproof apron covering the gown. All the staff shall carefully follow personal hygiene rules.
- All the containers containing working solutions of cleaning and disinfecting agents shall be equipped with lids and legibly labeled with the following data: description of the agent, its concentration, intended use, date of preparation and date of expiry. Containers with ready-to-use agents shall be labeled with the description of the agent and its intended use, if the agent may be used repeatedly, as well as with the date of the first use.
- For disinfection and pre-sterilization cleaning with chemical solutions immerse the medical devices in the working chemical solution so that all the gaps and grooves are filled. Before immersing disassemble all the parts that need assembling, open all the instruments with locking mechanisms and make some operational movements with them in the solution.
- Disinfection by wiping is allowed only for those medical devices that do not touch patients directly or for those devices which can not be immersed into the disinfecting solution because of their design features.
- After disinfection, wash away any residues of the disinfecting agent from reusable medical devices following instructions for the agent use.
- Cleaning efficiency is confirmed by several tests: amidopyrine-based test is used to check for any blood residues left and phenolphthalein-based test is used to check for any residue of cleaning agent left. Quality inspection shall be performed daily. It is to cover 1% simultaneously treated instruments of every item, but at least three pieces. Both internal and external device surfaces are subject to inspection! The pre-sterilization treatment log is the document confirming such quality inspection.
- This medical device is subject for physical sterilization (autoclaving).
- Sterilization shall be performed in the clean area.
- To transport medical device into the clean area use special closed containers for sterilization.
- The hospital staff involved in pre-cleaning and disinfecting medical devices must wear special protection clothes as well as a facemask, protection goggles and rubber gloves as well as a water-proof apron covering the gown!
- After sterilization, store the medical devices under the conditions to exclude secondary microbial contamination in a special cabinet protected against dust, moisture and temperature fluctuations.

- Storage period is marked on the packaging and depends on the type of packaging material as per the instruction for its use.
- Do not use sterilized medical devices if their storage period after sterilization is expired.
- Information on device sterilization is registered in a special log.

General recommendations on device preparation

- All the devices before being sterilized are subject to pre-sterilization cleaning in order to remove protein, fat and mechanical soiling as well as various medicines.
- For pre-sterilization cleaning disassemble all the parts that need assembling.
- The instruments are to be cleaned in a cleaning solution (for example, Schülke perform or 3% DESOMEDAN ID) while meeting all the personal safety requirements. While using solutions carefully follow all the manufacturer's instructions for use!
- To rinse the devices use demineralized water only (VE-water).
- Use only fresh daily prepared working solutions. If you have any gross soiling, replace the solution more often.
- For cleaning we recommend to use soft napless cloth, plastic brushes or cleaning pistols. Use of metal instruments is prohibited. For treatment use plastic or enameled (with intact enamel) containers covered with lids.
- For manual cleaning from the outside use a soft disposable napkin and disinfected brushes. Always clean the instruments with a brush under the surface of the cleaning solution in order to avoid dirty liquid splashing.
- The electronic motor must not be cleaned with ultrasound or compressed air.
- The medical devices shall be disinfected and cleaned simultaneously.
- For automatic cleaning and disinfection use a washer/disinfector (for example, Miele G7835 CD Thermal Disinfector) with thermal disinfection function.
- To avoid damages all the device parts must be safely secured in the washer/disinfector. The load carriers must ensure thorough rinsing of internal instrument cavities by proper nozzles.
- Instrument drying is an integral stage of automatic cleaning and disinfection process.
- After cleaning and disinfection is over instruments shall be removed by hands in antiseptic-treated disposable rubber gloves.
- All the cleaned and disinfected devices are subject to visual inspection. Critical areas need more close attention. All the devices and instruments with inner channels must be checked for any obstructions. Should any instruments fail to pass the inspection you must replace them without delay.
- Any instruments that are not clean enough are subject to the second cleaning cycle with thorough rinsing afterwards.
- Any bent instruments as well as instruments with micro cracks or other wear pattern are subject to immediate replacement as they are not fully able to perform their functions.
- If any instruments need some maintenance procedures, those maintenance procedures

shall come before the performance inspection.

- Before the performance inspection please lubricate all the hinged and screwed instruments. Reliable operation is ensured by the inspection. Thus, all the disassembled instruments must be assembled. If it is required, you may disassemble the instruments again after the inspection to sterilize them. While assembling/disassembling the devices follow applicable manufacturer's instructions.
- Any damaged or defective medical devices shall be sent for repairing or disposed.
- It is recommended to place the medical devices into the sterilization packaging. Follow the manufacturer's instructions regarding operation and the period while devices stay sterile!
- The packaging shall be appropriate and big enough for the below data to be labeled on it:
 - date of sterilization,
 - packing staff,
 - date of expiry,
 - description of the contents.
- For sterilization without using sterilization packaging the period for storing sterilized devices is 6 hours for sterilization in skeleton boxes and 20 days for sterilization in filter boxes, craft bags may be stored for three days maximum under sterile conditions. Any devices being sterilized without packaging must be used immediately after sterilization.
- Sterilize in vacuum autoclave (class B or S as per EN 13060) at 135°C. Sterilization specifications:

Sterilization mode					
Steam pressure inside the sterilization chamber, MPa (kgf/cm ²)		Operational temperature inside the sterilization chamber, °C		Duration of sterilization, minutes	
Rated value	Limit	Rated value	Limit	Rated value	Limit
0,20	±0,02	135	±2	5-25	±2

- To check for sterility the manufacturer recommends to use sterility indicator tapes (for example, Comply™ indicator tapes, manufactured by 3M).
- Returned instruments are those which are returned to the manufacturer after being used or without being used.
- Possible reasons for return include need for repairs or maintenance. Should any medical devices be returned to the manufacturer or its authorized representative, the latter will return the devices to the user only if they (the devices):
 - Have been subject to cleaning and disinfecting procedures as per manufacturer's instructions, are dry and declared to be hygienically clean, and/or
 - are expressly designed as being not soiled and have reliable packaging.
- The returned devices shall be decontaminated as in any case without delay in order to avoid device damage. Decontamination is prohibited if it is able to change or destroy the device, compromise analysis results or make the analysis impossible. If you have any doubts please consult the manufacturer.
- Every medical device being returned to the manufacturer shall be accompanied by the

relative documents containing the below information:

- Intended use of the medical device;
- Decontamination method;
- Date of decontamination;
- Last name of the person responsible for contamination;
- Reason for return.

All the instructions were validated by the device manufacturer as suitable for device reprocessing. Any reprocessing facilities are liable for reprocessing as well as for the equipment and the materials used and for the staff involved in reprocessing. The procedure must be validated and verified. Any deviations in the procedure stated in the instruction must be evaluated for their efficiency and possible adverse consequences.

Cleaning, disinfection and sterilization

The following points in particular are important with regard to caring for the material:

- Perform cleaning, disinfection and sterilization after every treatment.
- Always autoclave the material in sterilization packaging.
- Make sure that sterilization packaging is no more than 80% full.
- Always autoclave the material at 135°C for at least 5 minutes.
- If sterilized material is not used immediately, the material packaging must be labeled with the sterilization date.
- Nouvag AG recommends including a sterility indicator.


1. Control unit and pedal

Wipe the outside using micro-biologically tested surface disinfectants or a 70% isopropyl alcohol solution. The front plate of the control unit is sealed for this purpose and can be wiped clean.

2 Electronic motor 21

Reprocessing restrictions	Frequent reprocessing has only a limited impact on the electronic motor. The end of the product service life is normally determined by wear and damage through use. Electronic motor 21 is designed for 250 sterilization cycles.
---------------------------	---

INSTRUCTIONS	
At the location of use	Remove surface soiling with a disposable cloth/paper towel.
Storage and transport	No special requirements. Due to the risk of drying and corrosion, reprocessing must be performed without undue delay.
Preparation for cleaning	Remove soiling from the electronic motor with a disposable cloth/paper towel. Unscrew the motor cap, remove the cable and unscrew the handpiece carrier. 
Automatic cleaning and disinfection	Equipment: Washer-disinfector with a special load carrier that ensures motors are connected to the washer-disinfector (for example, Miele G7835 CD Thermal Disinfector) and channels are rinsed. Rinse the motor from the front. Use neutral or alkaline cleaning agents (for example, Schülke perform or 3% DESOMEDAN ID) in the manufacturer's recommended concentration. Place the electronic motor in the load carrier (ensure that the channels can be rinsed). Place the motor cap and cable, and handpiece holder, in the basket. Set a cleaning cycle that offers sufficient cleaning and rinsing. Perform the final rinse with fully demineralized water. Perform a 10-minute rinse cycle at 93°C to facilitate thermal disinfection. When removing the electronic motor, check the motor cap and cable, and handpiece holder, to verify whether soiling is still visible in the gaps and the grooves. If necessary, repeat the cycle or clean manually.
Manual cleaning	Equipment: neutral or alkaline cleaning agent (for example, Schülke perform or 3% DESOMEDAN ID), soft brush, running water (< 38°C). Procedure: Rinse off and brush away surface soiling on the electronic motor, motor cap and cable, and handpiece holder. Use a brush to apply cleaning agent to all surfaces and gaps. Rinse the electronic motor, motor cap and cable, and handpiece holder thoroughly under running water.
Manual disinfection	For manual disinfection, wipe the electronic motor, the motor cap and the plug and cable underneath it, and the handpiece holder, with a certified non-chlorine disinfectant.
Drying	If a drying program is not provided by the washer-disinfector, the electronic motor must be dried manually or in a drying cabinet until no visible moisture is left. The handpiece holder must then be screwed back onto the motor,

Inspection and maintenance  Art. No. 1974	Perform a visual inspection to check for damage, corrosion and wear. Connect the electronic motor with the Nou Clean spray nozzle (Art. No. 1974) by screwing it on and spray electronic motor for about 3 seconds and briefly wipe with a moist cloth (see instructions on the spray can). Once thorough spraying has been completed, screw the cable and motor cap back onto the electronic motor. 
Packaging 	Individual: Pack the electronic motor in individual packaging for sterile items. The bag must be large enough to ensure that the seal is not subject to strain. Nouvag AG recommends including a sterility indicator. Sets: Sort electronic motors on trays intended for this purpose or place them on allpurpose sterilization trays.
Sterilization	Autoclave in the vacuum autoclave (Class B or S as per EN 13060) at 135°C for at least 5 minutes*. When sterilizing several instruments during one sterilization cycle, do not exceed the maximum sterilizer load. A drying cycle must be added in the case of autoclaves without a post-vacuum function. Allow the electronic motor to dry in the bag for at least 1 hour at room temperature with the paper side facing upwards. * Temperature exposure times are based on country-specific guidelines and standards. Maximum exposure time: 25 minutes.
Storage	If the sterilized electronic motor is not used immediately after sterilization, the material packaging must be labeled with the sterilization date. Including a sterility indicator is recommended.

3. Tubing set, 2 m long (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.20102) Art. No. 1706 and tubing set Art. No. 6024/6025 (optional)

- Tubing set, 2 m long, sterile (Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.20102) Art. No. 1706 and tubing set Art. No. 6024/6025 (optional) are single-use only and may not be reused.
- Tubing sets must be disposed of properly after use!
- Do not use tubing sets when packaging is already open or damaged!
- Do not use tubing sets when expiry date has passed.
- Use only Nouvag sterile tubing sets with Art. No. 1706 and Art. No. 6024/6025 (optional)!
- The period for the tubing sets to stay sterile, apyrogenic and non-toxic is 5 years.

Sterility cannot be guaranteed by reusing and re-sterilization of tubing sets. The characteristics of the material change in a manner that can result in failure of the system. This may result in serious infections or even patient death worst-case.

4. Sterilizable connecting tubes (20 units) for Y-connector for contra angles with latch system Art. No. 1773

For sterilizable connecting tubes it is essential that all their pieces contacting with solutions have to be rinsed after every treatment. Pre-cleaning procedure includes cleaning with a cleaning agent and air flushing of the connecting tube channel. Hollow tubes to be re-sterilized shall be cleaned with a 5 cm³ (5 ml) syringe. After preliminary cleaning inspect the tubes for any leakages and make sure all of them are waterproof. Tubes must be also cleaned with compressed air. After that inspect the tubes once again for any damages and perform standard cleaning/disinfection cycle in the washer/disinfector. Sterilizable connecting irrigation tubes (20 units) for Y-connector for contra angles with latch system are resistant to autoclaving and designed for up to 20 cleaning, disinfection and sterilization cycles. After cleaning, disinfection and sterilization they keep being wear-resistant.

5. Handpieces and contra angles

Sterilization procedure for handpieces and contra angles is specified in Chapter "HOW TO USE HANDPIECES AND CONTRA ANGLES".

6. Stand for the electronic motor


Soiled stands for the electronic motor are cleaned using a neutral cleaning agent and then sterilized in accordance with the same instructions as for the electronic motor.

11. Maintenance

1. Replacing the control unit fuse

Users can replace faulty control unit fuses themselves. These are located at the rear of the device in the fuse slot beside the power switch (see picture 5):

- Unplug the power plug.
- Open the fuse slot using a screwdriver.
- Replace the faulty fuse T 3.15 AL 250 V AC.
- Slide the fuse holder back in and close the fuse slot.
- Check the mains voltage shown on the fuse slot.
- Plug in the power plug again.

1. Fuse slot locking mechanism
2. Display window for voltage setting
3. Fuse slot
4. Fuse 1
5. Fuse 2

2. Safety inspections

The safety inspection is a regular safety check that is compulsory for those operating medical devices. The objective of this measure is to ensure that device defects and risks to patients, users or third parties are identified in time.

An inspection interval of 2 years applies to the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories. NOUVAG AG offers a safety inspection service for its customers. Inspection results shall be registered in a special maintenance log.

If you have any questions regarding device maintenance and repair as well as after-sales servicing please contact our authorized representative in the Russian Federation: MedTechConsult, Ltd, 1 Nemanskiy proezd, building 1, office 129, Moscow, 123181, Russia, phone/fax +7 (495) 757-44-48, E-mail: medkonsult@mail.ru.

Table 19. Malfunctions and troubleshooting

Malfunction	Cause	Solution	Refer to operating instructions
Device is not operational	Control unit is not switched on	Set the power switch "I/O" to "I"	Chapter "Switching the device on and off"
	Power connection not established	Connect the control unit to the mains power supply	Chapter "Connection to the power supply"
	Incorrect operating voltage	Check the mains voltage	Chapter "Connection to the power supply"
	Faulty fuse	Replace the fuse	Chapter "Replacing the control unit fuse"
Motor does not run	Motor not switched on	Switch on the motor using the Vario pedal	Chapter "Operation using the Vario pedal"
	Motor not connected	Connect the motor cable to the control unit	Chapter "Device overview" Chapter "Connection to the power supply"
	Handpiece or contra angle not correctly assembled	Press the handpiece firmly onto the electronic motor until it clicks into position and check that it is secure by moving it slightly in the opposite direction.	Chapter "How to use handpieces and contra angles"
No irrigation fluid for instrument	Peristaltic pump (cooling pump) not switched on	Switch on the peristaltic pump.	Chapter "Operation using the Vario pedal"
	Tubing set incorrectly inserted	Insert tubing set correctly (note the direction)	Chapter "Device preparation"
	Tubing set clogged/ crusted matter visible	Replace the tubing set	Chapter "Device preparation"

	Bottle with sodium chloride solution not ventilated	Open the ventilation filter at the drip chamber	Chapter "Device preparation"
	Tubing set is dripping	Replace the tubing set	Chapter "Device preparation"
	Roller clamp of tubing set is closed	Open roller clamp all the way	Chapter "Device preparation"
Pedal is not operational	Pedal not connected	Connect the pedal to the control unit	Chapter "Device preparation"
	Incorrect operation	Check Operating Instructions	Chapter "Operation using the Vario pedal"

If a fault cannot be rectified, please contact your supplier or an authorized service centre.

Table 20. Error messages on display screen		
Error message/error code	Cause	Solution
Basic initialization/ W00	First initialization	
Reset to factory default values/ W01	Resetting to factory default values	
Memory failure/ E02	System failure	Send the control unit to the service centre
Operation failure/ E03	System failure	Send the control unit to the service centre
Program failure SW/ E04	System failure	Send the control unit to the service centre
Failure of user configuration SW/ E05	System failure	Send the control unit to the service centre
Display failure/ E06	System failure	Send the control unit to the service centre
Pump failure/ E07	System failure	Send the control unit to the service centre
Save factory default settings / User Configuration & Program	The message is displayed when standard program settings and parameter values are being saved using Nou-dongle.	
Factory default settings for saving / Program	The message is displayed when standard program settings are being saved.	
Pedal not connected/ E10	a) Pedal not connected. b) Faulty connector or cable	a) Connect the pedal b) Send the control unit and the pedal to the service centre.
Pedal test mode/ W11	Pedal test mode activated	Switch off the device for 5 seconds and after that switch it on again.
Keyboard test mode/ W12	Keyboard test mode activated	Switch off the device for 5 seconds and after that switch it on again.

Motor not connected/ E13	a) Motor not connected b) Faulty motor, motor cable, motor connector or control unit	a) Connect the motor b) Send the motor and the control unit to the service centre.
Unknown motor/ E16	a) Incorrect motor active. b) Correct motor active, but faulty motor, motor cable, motor connector or control unit.	a) Connect the correct motor b) Send the motor and the control unit to the service centre.
Pump opened/ E20	Motor is not operational when the pump compartment is opened in order to prevent injuries	Close the cooling pump compartment
Motor or pump test mode/ W21	Motor or pump test mode activated.	Switch off the device for 5 seconds and after that switch it on again.
AS-mode torque reached	The message is displayed when maximum AS-mode torque is reached	Release the pedal and re-start the motor by pressing the pedal once more.
Pedal locked W26, release the pedal	If the pedal is pressed while the device is being switched on, the pedal will not be operational.	Release the pedal for one second.
Faulty handpiece XX / E29	The handpiece has been over-exposed to torque during calibration or testing.	– Clean the handpiece and thoroughly spray it with Nou-Clean spray (not included in the delivery). – If the message keeps appearing after testing send the handpiece to the service centre.
Handpiece XX is Ok!	Tested handpiece is properly functioning	
Calibrating HPXX is Ok!	Calibrated handpiece is properly functioning	
Testing the handpiece XX	The handpiece is being tested	
Calibrating handpiece XX!	The handpiece is being calibrated	
Nou-dongle is plugged in	The message is displayed for 1 second after Nou-dongle is plugged in	

12. SAFETY REQUIREMENTS

Your own safety, safety of your team and of course of your patients is very important for us. Thus, it is essential to pay attention to the following information.

1. Manufacturer's EC Declaration

Please take into notice the below information on electromagnetic compatibility

Electromagnetic compatibility (EC)

Remark:

The **Product** subsequently referred to herein always denotes the **Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories.**

Changes or modifications to this product are prohibited as they may result in increased emissions or decreased immunity performance of the product and could cause EMC issues with this or other equipment. This product is designed and tested to comply with applicable regulations regarding electromagnetic compatibility and shall be installed and put into service according to the EMC information stated as follows.

WARNING

Use of portable phones or other radio frequency (RF) emitting equipment near the product may cause unexpected or adverse operation.

WARNING

The product shall not be used adjacent to, or stacked with, other equipment. If adjacent use or storage is necessary, the product shall be tested to verify normal operation in the configuration in which it is being used.

WARNING

Compliant cables and accessories.

The use of accessories, transducers and cables other than those specified may result in increased emissions or decreased immunity performance of the product.

The table below lists cables, transducers and other applicable accessories to which the manufacturer claims EMC compliance.

NOTE: Any supplied accessories that do not affect EMC compliance are not listed.

Table 21.

Description	Cable length max.
Electronic motor 21	2.0 m
Vario Footpedal IPX8 for MD 11	3.15 m
Vario Footpedal IPX8 for MD 30	3.15 m

Table 22.

Guidance and manufacturer's declaration – electromagnetic emissions		
The Product is intended for use in the electromagnetic environment specified below. The customer or the user of the Product should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Product uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Product is suitable for use in all establishments, including
Harmonic emissions IEC 61000-3-2	Class A	independent establishments and those directly connected to the

Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	public low-voltage power supply network that supplies buildings used for domestic purposes.
---	----------	---

Guidance and manufacturer's declaration – electromagnetic immunity			
The Product is intended for use in the electromagnetic environment specified below. The customer or the user of the Product should assure that it is used in such an environment.			
Immunity tests	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 6 kV contact +/- 8 kV air	+/- 6 kV contact +/- 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %
Electrical fast transient/burst IEC 61000-4-4	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+/- 1 kV differential mode +/- 2 kV common mode	+/- 1 kV differential mode +/- 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5 % UT (> 95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles < 5 % UT (> 95 % dip in UT) for 5 seconds	< 5 % UT (> 95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles < 5 % UT (> 95 % dip in UT) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Product requires continued operation during power mains interruptions, it is recommended that the Product be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment..
Note: UT is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity for not life support equipment			
The Product is intended for use in the electromagnetic environment specified below. The customer or the user of the Product should assure that it is used in such an environment.			
Immunity tests	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance

			Portable and mobile RF communications equipment should be used no closer to any part of the Product, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF IEC 61000-4-6	3 V rms 150 kHz to 80 MHz outside ISM bands	10 V rms 150 kHz to 80 MHz outside ISM bands	$d = 1,2 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m 80 MHz to 2.5 GHz	$d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800 MHz to 2.5 GHz
			where P is the maximum output power rating in the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a, should be less than the compliance level in each frequency range ^b. Interference may occur in the vicinity of equipment marked with the following symbol: 

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- Fixed strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Product is used exceeds the applicable RF compliance level above, the Product should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Product.

^a over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the not life support equipment			
The Product is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Product can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Product as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power m of transmitter, W	Separation distance according to frequency of transmitter, m		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the higher frequency range applies.			
Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

2. Integrated peristaltic pump

The integrated peristaltic pump is used to cool tissue in order to prevent damage to the tissue. It may only be operated with watery solutions, such as 0.9% Sodium Chloride irrigation solution or "Ringer" solution (not included in delivery). **Medication supply using the integrated pump is expressly prohibited.**

3. Modification and misuse

- If modifications are made by the user/operator or a third party to the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories provided, or if these items are used by these individuals for a purpose other than that for which the items were designed (according to chapter "Intended Use"), the manufacturer assumes no responsibility and the guarantee is void.

4. Essential requirements

The Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories may only be operated by qualified and trained

Improper use or repair of the device and failure to observe these instructions relieve us from any obligation arising from warranty provisions or other claims.

personnel.

The use of third-party products is the responsibility of the operator. Functionality and patient safety cannot be guaranteed with third-party accessories.

Prior to using the device, before startup, and before operation, the user must always ensure that the device and accessories are in good working order and are clean, sterile (if it is necessary) and operational.

Repairs may only be performed by authorized NOUVAG service technicians.

Use Nou-Clean spray (is not included into the scope of delivery; optional) for maintenance and care of the electronic motor, handpieces and contra angles. Using other care products can result in malfunctions and/or cause the warranty to be revoked.

5. During use

The device is not sterile on delivery. All sterilizable parts must be sterilized before use (see Chapter "Cleaning, disinfection and sterilization").

To prevent injury, never touch drill bits or burrs while they are still running.

Handpieces and contra angles may only be attached when the electronic motor is not running.

Never operate the clamping mechanism of the handpieces or contra angles while the system is running. This could result in instrument damage.

Make sure that the supply voltage setting corresponds to the local mains voltage.

The manufacturer runs the quality control and risk analysis system. The risk management procedures performed meet the requirements of EN ISO 14971:2012. Risk control is effected at device production and post-production stages. Risk control is also performed during post marketing observance of device use results. Risk analysis is an integral part of the technical documents of the device.

6. Safety inspections

The safety inspection is a regular safety check that is compulsory for those operating medical devices. The objective of this measure is to ensure that device defects and risks to patients, users or third parties are identified in time.

An inspection interval of 2 years applies to the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories. NOUVAG AG offers a safety inspection service for its customers. Inspection results shall be registered in a special maintenance log.

If you have any questions regarding device maintenance and repair as well as after-sales servicing please contact our authorized representative in the Russian Federation: MedTechConsult, Ltd, 1 Nemanskiy proesd, building 1, office 129, Moscow, 123181, Russia, phone/fax +7 (495) 757-44-48, E-mail: medkonsult@mail.ru.

The medical device is entirely safe for the environment during its operation.

13. Medical device classification

Electrical protection class:	I
Method of protection against electric shock:	B Type
Intended use:	Power instruments
Frequency of application:	Multiple application
Mode of operation:	Intermittent mode of operation
Type of perceived mechanical impact	Stationary
Type of human body contact	Short-term contact (less than 24 hours) with oral mucosa and mouth bone tissue

Compliance with international standards

The Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories comply with Council Directive 93/42/EC on Medical Devices and with the below standards:

EN60601.1:2006	General electrical safety standards
EN60601-1-2/07.2007	Electromagnetic compatibility. Requirements and testing. Medical equipment. Biological evaluation of medical devices. Part 1. Evaluation and studies.
EN 62366-2008	Medical devices. Application of usability engineering to medical devices.
EN61000.4.2	(EMC). Part 4: Testing and measuring techniques.– Section 2: Electrostatic discharge immunity test. Basic EMC Publication .
EN61000.4.3	(EMC). Part 4: Testing and measuring techniques.– Section 3: Radiated, radio-frequency, electromagnetic field immunity test. Basic EMC Publication.
EN61000.4.4	(EMC). Part 4: Testing and measuring techniques.– Section 4: Electrical fast transient/burst immunity test. Basic EMC Publication.
EN61000.4.5	(EMC). Electromagnetic compatibility (EMC). Part 4: Testing and measuring techniques. Section 5. Surge immunity test.
EN61000.4.6	Electromagnetic compatibility (EMC). Part 4: Testing and measuring techniques. Section 6. Immunity to conducted disturbances, induced by radio-frequency fields.
EN61000.4.8	(Electromagnetic compatibility (EMC). Part 4: Testing and measuring techniques. Section 8. Power frequency magnetic field immunity test.
EN61000.4.11	Electromagnetic compatibility (EMC). Part 4: Testing and measuring techniques. Section 11. Voltage dips, short interruption and voltage variations immunity test.
ISO 3964	Dental handpieces. Coupling dimensions.

EN 10993-1:2009

Biological evaluation of medical devices. Part 1: Evaluation and testing within a risk management.

EN ISO 17664-2004

Sterilization of medical devices – Information to be provided by the manufacturer for the processing of resterilizable medical devices

ISO 13485:2012

Design and development, production, distribution and servicing of medical devices.

CE certificate

Compliance with Council Directive 93/42/EEC on Medical Devices

14. DEVICE PACKAGING AND MARKING

1. MARKING

The medical device is marked as follows:

Table 23. Medical device marking.


	Important information		Autoclavable at max. 135° for at least 5 minutes
	Do not use if the packaging is damaged		Suitable for thermal disinfection
	Warning		Sterilized using ethylene oxide
	Manufacturer		Observe the operation manual
	Certified by the Canadian Standards Association (CSA)		Electrical and electronic devices that have reached the end of their service life comprise hazardous waste and may not be disposed of together with household waste. Valid local disposal regulations apply.
	Type B method of protection against electric shock		Symbol indicating the serial number with the date of manufacturing (year/month).
	Do not reuse		Symbol indicating the order number.
	Biohazard		Symbol indicating the lot number.
	Motor 1		Motor 2
	Pedal		Protective ground
IPX8	Protection against contact and continual submerging		The device is designed for intermittent duty operation at "1 min ON/3 min OFF" at 4 cycles, after 15 min. brake
	Date of manufacturing		CE mark with reference number of the notified body
	Date of expiry		Refer to the instructions for use

Pic. 32. The plate on the rear panel of the device.


Name and address of the manufacturer


Description of the medical device


Reference number of the medical device


Serial number of the medical device

Symbol indicating power consumption


Fuses used

Symbol indicating that the CE certificate has been obtained and that the medical device is certified by the Canadian Standards Association (CSA).

Symbol referring to the Operation Manual


Symbol indicating hazardous waste which may not be disposed of together with household waste

Symbol indicating the country of origin (Switzerland)

2D datamatrix (QR-code)


The device is designed for intermittent duty operation at "1 min ON/3 min OFF" at 4 cycles, after 15 min. brake

Picture 33. Symbol indicating Type B method of protection against electric shock is located on the front panel of the control unit.

2. PACKAGING

Individual packaging of the handpiece or contra angle


Pic. 34. Individual packaging of the handpiece or contra angle

Contra angle packaging made of white carton; its dimensions are 205x075x035 mm

Contra angle packaging made of blue plastic; its dimensions are 200x070x030 mm

Picture 35. Marking on the individual packaging

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Description of the medical device

REF 5052nou

Reference number of the medical device

SN 0995H1601R

Serial number of the medical device

CE mark

Refer to the Operation Manual

Autoclavable at max. 135° for at least 5 minutes.

Suitable for thermal disinfection

2D datamatrix (QR-code)

Sterile products

1706 Tube set, sterile, 2 m long (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010).


Pic. 36. Individual packaging of the sterile product.


The individual commercial sterile packaging (see picture 36) is made of PP50/PE30 PEEL plastic film, packaging dimensions are 27cm X 12.5 cm. The packaging guarantees the products stay sterile, nontoxic and apyrogenic. Design of the commercial packaging ensures product integrity, prevents microbial penetration and protects the product against mechanical damages. The individual commercial packaging does not affect toxicological properties of the packed items. The individual packaging ensures the products stay sterile for at least 5 years.

Pic. 37. Marking on the individual packaging of the sterile product

assembly instructions

- 1 – description of the medical device
- 2 – catalogue number
- 3 – do not use if the packaging is damaged
- 4 – do not reuse
- 5 – do not use if the packaging is damaged
- 6 – address of the manufacturer
- 7 – CE mark
- 8 – contact details of the manufacturer
- 9 – sterilized using ethylene oxide
- 10 – latex-free
- 11 – necessary to get acknowledged with the instructions
- 12 – date of sterilization

Pic. 38. Shipping packaging of the sterile product

The shipping packaging of the sterile product is shown in picture 38. It is a white carton, dimensions 310x215x80 mm with 10 individual packaging's packed into it.

Pic. 39. Marking on the shipping packaging of the sterile product.

Einweg-Schlauchset 2m
Tubing set 2m

description of the medical device

REF 1706

catalogue number

 **QTY 10**

10 items per a packaging


warning against product use if the packaging is damaged

contains phthalates

STERILE EO

sterilized using ethylene oxide

do not reuse

 **P0803-1510**

lot code number

 **2015-10**

date of manufacturing

2020-10

date of expiry

1275

CE mark

2D datamatrix (QR-code)

Manufacturer:
NOUVAG AG
 St. Gallerstr. 23-25
 CH-9403 Goldach
 Tel. : 0041 71 846 66 00
 Fax. : 0041 71 845 35 36
 E-mail : info@nouvag.com
 Homepage: www.nouvag.com

contact details of the manufacturer


EU-Mandatory:
NOUVAG GmbH
 Schulthaissstr. 15
 D-78462 Konstanz
 Tel. : 0049 7531 1290 0
 Fax. : 0049 7531 1290 12
 E-mail : info-de@nouvag.com

contact details of the distributor within the EU

Pic. 40. Shipping packaging of SANICAR Mobile Device Trolley

The shipping packaging of SANICAR Mobile Device Trolley is a 103x65x15 cm carton. It ensures vibration- and shock- resistance (proper operation after bouncing during transportation) of SANICAR mobile device trolley. The marking on the shipping packaging is as follows.

Pic. 41. Marking on the shipping packaging

The marking contains the below information:

Manufacturer:
NOUVAG AG
St. Gallerstrasse, 23-25
CH - 9403 GOLDACH

name and address of the manufacturer

Type: Sanicar fahrbar

description of the product

Serie No. 01/15

serial number

Power: ~230V 50Hz max. 16A

power consumption

CE mark

warning

Art. Nr. 1898d

Reference number

Shipping packaging

Inside the shipping package the medical device is packed as follows:

Pic. 42. Packing of the medical device.

The device is packed into a 400x400 mm polyethylene flat antistatic bag and after that in a 308x300x385 mm foamed polystyrene packaging (2 units). The shipping box contains a CD-ROM with the Operation Manual on it, a warranty card and that double packaging.

The shipping packaging (picture 43) protects the medical device against any mechanical shocks caused by transportation. Being packed in the shipping packaging the medical device is vibration- and shock-resistant. The packaging also protects the device during handling operations. The packaging is a 395x395x405 mm carton marked and labeled as follows (pictures 43, 44, 45, 46, 47).

Pic. 43. Shipping packaging

Pic. 44. Shipping packaging labeling

Pic. 45. Marking on the shipping packaging

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medical device

description of the

REF 2005

Reference number

SET 3N 2482E1602R


Serial number

2018-08

Date of expiry

NOUVAG AG
St. Gallerstrasse 13-25
CH-9 73 Goldsch
www.nouvag.ch


name and address of the manufacturer

CE mark

2D datamatrix (QR-code)

Pic. 46. Information labeled on the shipping packaging.

The information labeled on the shipping packaging includes:

1. Company logo with address and phone numbers of the manufacturer.
2. The following note: "The operation manual is not available in printed form with this product. The operation manual as PDF is on the CD-ROM that was included in the delivery. A CD-ROM drive and the Adobe Acrobat Reader software is needed to read or print this operation manual. To order the

operation manual in printed form please check on the type plate for the address or visit our website (www.nouvag.com). If you send us the article number (31686 EN) with the language code (RU for Russian) and the note "Printed copy", a printed manual will be sent to you free of charge.


The shipping packaging is marked with the following handling symbols.

Fig 47. Handling symbols on the shipping packaging.


Table 24. Handling symbols on the packaging.

	Keep dry The device must not be located near moisture sources or in a humid facility
	Handle with care. Fragile goods. Goods shall be handled with care. Such a symbol is used to mark goods which may be damaged in the result of tossing or shock.

Transportation conditions

The symbol is used to mark the transportation requirements (humidity and atmospheric pressure).

Transportation conditions

The symbol is used to mark the transportation requirements (temperature).

15. TRANSPORTATION AND STORAGE CONDITIONS.

There are no special storage requirements. The device may be stored at the below conditions. The device packed in the manufacturer's packaging should be stored in the warehouse and put on the racks. Avoid putting any other equipment or goods on the device.

There are no special transportation requirements as well. The device doesn't require disassembling for transportation. The device may be shipped by any covered vehicles or in multi-purpose containers as per the relevant shipping regulations at the below ambient conditions. Arrangement and attachment of shipping packagings must ensure their stable position preventing packagings from moving and bumping into each other.

Transportation and storage conditions:

Relative humidity:	Maximum 90 %
Temperature:	0 °C up to 60 °C
Atmospheric pressure:	700 hPa up to 1060 hPa

16. SERVICE LIFE

Shelf life of the medical device is 5 years. Average service life is **7 years from the date of manufacturing provided the device is operated according to the Operation Manual**. After the service life has been expired this medical device must be disposed of. After exploitation resource is over the medical device must be disposed of as per national regulations of the importing country as well as recommendations of the manufacturer.

For sterile disposable devices:

1706 Tube set, sterile, 2 m long (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010).

The service life after the packaging is opened is 3 years. The guaranteed storage period is 5 years.

17. DISPOSAL

The medical device shall be disposed of and destroyed as per the classification and rules of assembling, operation, decontamination, handling, storage, transportation, recording and disposal enforced by the competent authority of the importing country (Class B of medical waste hazard as per the Sanitary Regulations SanPin 2.1.7.2790-10 according to Russian Classification).

Do not dispose of the device with household waste!

To ensure environmental protection, old devices can be returned to the dealer or manufacturer.

Electrical and electronic devices that have reached the end of their service life comprise hazardous waste and may not be disposed of with household waste. Prevailing national and local disposal regulations apply.

Electronic motors that have reached the end of their service life must not be disposed of with household waste. Electronic motors must be sterilized before disposal. Please observe currently valid national disposal regulations for infectious waste (Class B of medical waste hazard as per the Sanitary Regulations SanPin 2.1.7.2790-10 according to Russian Classification).

When disposing of the device, device parts and accessories, the requirements specified in legislation must be followed.

To ensure environmental protection, old devices can be returned to the dealer or manufacturer..

Contaminated single-use tubing sets are subject to specific disposal requirements. Please observe currently valid national disposal regulations for infectious waste (Class B of medical waste hazard as per the Sanitary Regulations SanPin 2.1.7.2790-10 according to Russian Classification).

The user is liable for proper disposal of this equipment at specially designated collection points for electrical and electronic waste recycling. Separate collection and recycling of waste electrical and electronic equipment during disposal contributes to preservation of natural resources. Besides, use of breakthrough technologies of waste recycling encourages protection of staff health and environment. If you have any questions regarding disposal of medical equipment with expired service life please refer to its supplier.

18. WARRANTY COVERAGE

The manufacturer Nuvag AG provides warranty coverage of the medical device. Manufacturer's address is St.Gallerstrasse, 23-25, CH 9403 Goldach, Switzerland. Purchasing the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories entitles you to a 1-year warranty. If you return the warranty card for registration within four weeks of the date of purchase, warranty coverage will be extended for a further 6 months. Warranty coverage doesn't extend to the consumables. Improper use or repair, or the device breaking because of failure to observe these instructions, will relieve the manufacturer from any obligation arising from warranty provisions or other claims.

Modification and misuse

- If modifications are made by the user/operator or a third party to the Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories provided, or if these items are used by these individuals for a purpose other than that for which the items were designed (according to chapter "INTENDED USE"), the manufacturer assumes no responsibility and the guarantee is void.

Any claims or questions regarding handling of this medical device please send to:

MedTechConsult, Ltd, 1 Nemanskiy prosed, building 1, office 129, Moscow, 123181, Russia, phone/fax +7 (495) 757-44-48, E-mail: medkonsult@mail.ru.

Appendix 1. Medical device components

The Dental surgery portable motor device models: MD 11 in versions MD 11-1; MD 11-2; MD 11-3; MD 11-4; MD 30 in versions MD 30-1; MD 30-2; MD 30-3; MD 30-4 with accessories

I. Dental surgery portable motor device MD 11-1, including:

1.	3335	MD11 control unit with an integrated cooling pump system, single-motor,	1 unit;
2.	2097nou	Electronic motor 21, speed 40,000 revolutions/minute, including 2 m long motor cable, (if necessary)	2 units;
3.	1866nou	Vario IPX8 foot pedal for MD 11 with 3.15 m long cable,	1 unit;
4.	1706	Tube set , sterile, 2 m long (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
5.	1770	Stand for irrigation fluid bottle (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
6.	1873	Clip set (10 units) for tube set attachment at motor, (if necessary)	2 sets;
7.	1881	Clip set (3 units) for tube set attachment at handpieces and contra angles, (if necessary)	2 sets;
8.	1958	Spray nozzle for Nou-Clean spray, for care of electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
9.	1974	Spray nozzle for Nou-Clean spray, for care of instruments (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
10.	1170	Stand for electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
11.	28879	Fuse T 3.15 AL 250 V AC,	2 units;
12.	22262	Power cable, 3 m long,	1 unit;
13.	31686	User Manual on CD-ROM,	1 unit.

II. Dental surgery portable motor device MD 11-2, including:

1.	3335	MD11 control unit with an integrated cooling pump system, single-motor,	1 unit;
2.	2097nou	Electronic motor 21, speed 40,000 revolutions/minute, including 2 m long motor cable, (if necessary)	2 units;
3.	1866nou	Vario IPX8 foot pedal for MD 11 with 3.15 m long cable,	1 unit;
4.	1706	Tube set , sterile, 2 m long (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
5.	1770	Stand for irrigation fluid bottle (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
6.	1873	Clip set (10 units) for tube set attachment at motor, (if necessary)	2 sets;
7.	1881	Clip set (3 units) for tube set attachment at handpieces and contra angles, (if necessary)	2 sets;
8.	1958	Spray nozzle for Nou-Clean spray, for care of electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
9.	1974	Spray nozzle for Nou-Clean spray, for care of instruments (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
10.	1170	Stand for electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
11.	5201nou	Drill contra angle, 32:1, shank diameter of $\varnothing$ 2.35 mm, with internal/external	1 unit;

	cooling, and button clamp,	
12. 28879	Fuse T 3.15 AL 250 V AC,	2 units;
13. 22262	Power cable, 3 m long,	1 unit;
14. 31686	User Manual on CD-ROM,	1 unit.

III. Dental surgery portable motor device MD 11-3, including:

1. 3335	MD11 control unit with an integrated cooling pump system, single-motor,	1 unit;
2. 2097nou	Electronic motor 21, speed 40,000 revolutions/minute, including 2 m long motor cable, (if necessary)	2 units;
3. 1866nou	Vario IPX8 foot pedal for MD 11 with 3.15 m long cable,	1 unit;
4. 1706	Tube set , sterile, 2 m long (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
5. 1770	Stand for irrigation fluid bottle (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
6. 1873	Clip set (10 units) for tube set attachment at motor, (if necessary)	2 sets;
7. 1881	Clip set (3 units) for tube set attachment at handpieces and contra angles, (if necessary)	2 sets;
8. 1958	Spray nozzle for Nou-Clean spray, for care of electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
9. 1974	Spray nozzle for Nou-Clean spray, for care of instruments (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
10. 1170	Stand for electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
11. 5053nou	Drill contra angle, 20:1, shank diameter of $\varnothing$ 2.35 mm, with internal/external cooling, and button clamp,	1 unit;
12. 28879	Fuse T 3.15 AL 250 V AC,	2 units;
13. 22262	Power cable, 3 m long,	1 unit;
14. 31686	User Manual on CD-ROM,	1 unit.

IV. Dental surgery portable motor device MD 11-4, including:

1. 3335	MD11 control unit with an integrated cooling pump system, single-motor,	1 unit;
2. 2097nou	Electronic motor 21, speed 40,000 revolutions/minute, including 2 m long motor cable, (if necessary)	2 units;
3. 1866nou	Vario IPX8 foot pedal for MD 11 with 3.15 m long cable,	1 unit;
4. 1706	Tube set , sterile, 2 m long (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
5. 1770	Stand for irrigation fluid bottle (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
6. 1873	Clip set (10 units) for tube set attachment at motor, (if necessary)	2 sets;
7. 1881	Clip set (3 units) for tube set attachment at handpieces and contra angles, (if necessary)	2 sets;
8. 1958	Spray nozzle for Nou-Clean spray, for care of electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
9. 1974	Spray nozzle for Nou-Clean spray, for care of instruments (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
10. 1170	Stand for electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;

	2010/07320 dated 05.07.2010),	
11.	5052nou LED Drill contra angle, 20:1, shank diameter of $\varnothing$ 2.35 mm, with internal/external cooling, and button clamp,	1 unit;
12.	28879 Fuse T 3.15 AL 250 V AC,	2 units;
13.	22262 Power cable, 3 m long,	1 unit;
14.	31686 User Manual on CD-ROM,	1 unit.

Accessories for dental surgery portable motor device MD 11 all versions:

1.	1712	Internal cooling nozzle for contra angles with latch system	up to 5 units
2.	1773	Sterilizable connecting tubes (20 units) for Y-connector for contra angles with latch system	up to 5 units
3.	1777	Y-connector for internal/external cooling of contra angles	up to 5 units
4.	39117	Internal cooling nozzle clip for contra angles with push button system	up to 5 units
5.	39116	Internal cooling pipe for contra angles with push button system	up to 5 units
6.	39149	Internal cooling nozzle clip for contra angles 20:1 LED with push button system	up to 5 units
7.	39150	Internal cooling pipe for contra angles 20:1 LED with push button system	up to 5 units
8.	39152	Y-connector with three tubes for internal/external cooling of contra angles 20:1 LED with push button system	up to 5 units
9.	1029	Suitcase	1 unit;
10.	1898d	SANICAR Mobile Device Trolley with two shelves and 5 m long cable	1 unit

V. Dental surgery portable motor device MD 30-1, including:

1.	3330	MD30 control unit with an integrated cooling pump system, two-motor,	1 unit;
2.	2097nou	Electronic motor 21, speed 50,000 revolutions/minute, including 2 m long motor cable, (if necessary)	2 units;
3.	1510nou	Vario IPX8 foot pedal for MD 30 with 3.15 m long cable,	1 unit;
4.	1706	Tube set , sterile, 2 m long (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
5.	1770	Stand for irrigation fluid bottle (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
6.	1873	Clip set (10 units) for tube set attachment at motor, (if necessary)	2 sets;
7.	1881	Clip set (3 units) for tube set attachment at handpieces and contra angles, (if necessary)	2 sets;
8.	1170	Stand for electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
9.	1958	Spray nozzle for Nou-Clean spray, for care of electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
10.	1974	Spray nozzle for Nou-Clean spray, for care of instruments (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;

- | | | | |
|-----|-------|--------------------------|----------|
| 11. | 28879 | Fuse T 3.15 AL 250 V AC, | 2 units; |
| 12. | 22262 | Power cable, 3 m long, | 1 unit; |
| 13. | 31686 | User Manual on CD-ROM, | 1 unit. |

VI. Dental surgery portable motor device MD 30-2, including:

- | | | | |
|-----|---------|--|----------|
| 1. | 3330 | MD30 control unit with an integrated cooling pump system, two-motor, | 1 unit; |
| 2. | 2097nou | Electronic motor 21, speed 50,000 revolutions/minute, including 2 m long motor cable, (if necessary) | 2 units; |
| 3. | 1510nou | Vario IPX8 foot pedal for MD 30 with 3.15 m long cable, | 1 unit; |
| 4. | 1706 | Tube set , sterile, 2 m long (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010), | 1 unit; |
| 5. | 1770 | Stand for irrigation fluid bottle (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010), | 1 unit; |
| 6. | 1873 | Clip set (10 units) for tube set attachment at motor, (if necessary) | 2 sets; |
| 7. | 1881 | Clip set (3 units) for tube set attachment at handpieces and contra angles, (if necessary) | 2 sets; |
| 8. | 1170 | Stand for electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010), | 1 unit; |
| 9. | 5201nou | Drill contra angle, 32:1, shank diameter of $\varnothing$ 2.35 mm, with internal/external cooling, and button clamp, | 1 unit; |
| 10. | 1958 | Spray nozzle for Nou-Clean spray, for care of electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010), | 1 unit; |
| 11. | 1974 | Spray nozzle for Nou-Clean spray, for care of instruments (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010), | 1 unit; |
| 12. | 28879 | Fuse T 3.15 AL 250 V AC, | 2 units; |
| 13. | 22262 | Power cable, 3 m long, | 1 unit; |
| 14. | 31686 | User Manual on CD-ROM, | 1 unit. |

VII. Dental surgery portable motor device MD 30-3, including:

- | | | | |
|-----|---------|---|----------|
| 1. | 3330 | MD30 control unit with an integrated cooling pump system, two-motor, | 1 unit; |
| 2. | 2097nou | Electronic motor 21, speed 50,000 revolutions/minute, including 2 m long motor cable, (if necessary) | 2 units; |
| 3. | 1510nou | Vario IPX8 foot pedal for MD 30 with 3.15 m long cable, | 1 unit; |
| 4. | 1706 | Tube set , sterile, 2 m long (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010), | 1 unit; |
| 5. | 1770 | Stand for irrigation fluid bottle (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010), | 1 unit; |
| 6. | 1873 | Clip set (10 units) for tube set attachment at motor, (if necessary) | 2 sets; |
| 7. | 1881 | Clip set (3 units) for tube set attachment at handpieces and contra angles, (if necessary) | 2 sets; |
| 8. | 1170 | Stand for electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010), | 1 unit; |
| 9. | 5053nou | Drill contra angle, 20:1, shank diameter of $\varnothing$ 2.35 mm, with internal/external cooling, and button clamp , | 1 unit; |
| 10. | 1958 | Spray nozzle for Nou-Clean spray, for care of electronic motor (Marketing | 1 unit; |

	Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	
11. 1974	Spray nozzle for Nou-Clean spray, for care of instruments (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
12. 28879	Fuse T 3.15 AL 250 V AC,	2 units;
13. 22262	Power cable, 3 m long,	1 unit;
14. 31686	User Manual on CD-ROM,	1 unit.

VIII. Dental surgery portable motor device MD 30-4, including:

1. 3330	MD30 control unit with an integrated cooling pump system, two-motor,	1 unit;
2. 2097nou	Electronic motor 21, speed 50,000 revolutions/minute, including 2 m long motor cable, (if necessary)	2 units;
3. 1510nou	Vario IPX8 foot pedal for MD 30 with 3.15 m long cable,	1 unit;
4. 1706	Tube set , sterile, 2 m long (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
5. 1770	Stand for irrigation fluid bottle (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
6. 1873	Clip set (10 units) for tube set attachment at motor, (if necessary)	2 sets;
7. 1881	Clip set (3 units) for tube set attachment at handpieces and contra angles, (if necessary)	2 sets;
8. 1170	Stand for electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
9. 5052nou	LED Drill contra angle, 20:1, shank diameter of $\varnothing$ 2.35 mm, with internal/external cooling, and button clamp ,	1 unit;
10. 1958	Spray nozzle for Nou-Clean spray, for care of electronic motor (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
11. 1974	Spray nozzle for Nou-Clean spray, for care of instruments (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010),	1 unit;
12. 28879	Fuse T 3.15 AL 250 V AC,	2 units;
13. 22262	Power cable, 3 m long,	1 unit;
14. 31686	User Manual on CD-ROM,	1 unit.

Accessories for dental surgery portable motor device MD 30 all versions:

1. 1712	Internal cooling nozzle for contra angles with latch system	up to 5 units
2. 1773	Sterilizable connecting tubes (20 units) for Y-connector for contra angles with latch system	up to 5 units
3. 1777	Y-connector for internal/external cooling of contra angles	up to 5 units
4. 39117	Internal cooling nozzle clip for contra angles with push button system	up to 5 units
5. 39116	Internal cooling pipe for contra angles with push button system	up to 5 units
6. 39149	Internal cooling nozzle clip for contra angles 20:1 LED with push button system	up to 5 units
7. 39150	Internal cooling pipe for contra angles 20:1 LED with push button system	up to 5 units


8.	39152	Y-connector with three tubes for internal/external cooling of contra angles	units
		20:1 LED with push button system	up to 5
9.	1029	Suitcase	units
			1 unit;
10.	1898d	SANICAR Mobile Device Trolley with two shelves and 5 m long cable	1 unit

Appendix 2. Tables of device components.

Table 1. Device components and accessories.

No	Description of the component	Appearance	Weight (grams), not exceeding	Dimensions (mm), not exceeding
1	3335 MD11 control unit with an integrated cooling pump system, single-motor	 A white, rectangular control unit with a digital display and several buttons. The text "MD 11" is visible on the front panel.	3300	250x260x110
2	3330 MD30 control unit with an integrated cooling pump system, two-motor		3700	250x260x110

				
3	2097nou Electronic motor 21, speed 50,000 revolutions/minute, including 2 m long motor cable		300	Electronic motor: length – 80, Diameter – 20, Length of motor coupling – 30, diameter – 8 Cable length 2000
4	1866 nou Vario IPX8 foot pedal for MD 11 with 3.15 m long cable		1000	140 x 240 x 90 (Width x Depth x Height)

				
5	1510 nou Vario IPX8 foot pedal for MD 30 with 3.15 m long cable		1680	255 x 235 x 60 (Width x Depth x Height)
6	1770 Stand for irrigation fluid bottle (Marketing Authorization №FSZ (ФСЗ))		140	Height 46 cm Length under the hook 50, hook 40

	2010/07320 dated 05.07.2010)			
7.	1170 Stand for electronic motor holder (Marketing Authorization №FSZ (ФСЗ) 2010/07320 dated 05.07.2010)		50	Length 100 Width 60 Height 35
8.	28879 Fuse T 3.15 AL 250 V AC		7	Length – 20 Diameter -5

9.	22262 Power cable, 3 m long		320	Length 3000
10	31686 User Manual on CD-ROM		10	Diameter 120
11	1706 Tube set , sterile, 2 m long (Marketing Authorization №FSZ (ΦСЗ) 2010/07320 dated 05.07.2010)		40	Length 2000

				
12	1958 Spray nozzle for Nou-Clean spray, for care of electronic motor (Marketing Authorization №FSZ (ΦСЗ) 2010/07320 dated 05.07.2010), 1974 Spray nozzle for Nou-Clean spray, for care of instruments (Marketing Authorization №FSZ (ΦСЗ) 2010/07320 dated 05.07.2010)		20	Nozzle length 35, Nozzle diameter 10, Connection diameter 25, Connection height 20.

13	1773 Sterilizable connecting tubes (20 units) for Y-connector for contra angles with latch system		4	Length -159 Diameter – 3
14	1873 Clip set (10 units) for tube set attachment at motor		5	Diameter 8 Height 8 Rib 10

				
15	1881 Clip set (3 units) for tube set attachment at handpieces and contra angles		3	Diameter 18, Height 20,
16	1777 Y-connector for internal/external cooling of contra		20	Length of the fastener 8, Diameter 5

	angles			
17	1712 Internal cooling nozzle for contra angles with latch system		10	Length of the base 25, Length of the fastener 20, Length of the latch 5 x 3
18	39116 Internal cooling pipe for contra angles with push button system		10	Length of the base 23, Fastener 25
19	39117 Internal cooling nozzle clip for		10	Diameter 10,

	contra angles with push button system			Height 10, Rib 5
20	39149 Internal cooling nozzle clip for contra angles 20:1 LED with push button system		50	Diameter 10, Height 20, Rib 6
21	39150 Internal cooling pipe for contra angles 20:1 LED with push button system		10	Base 23 Fastener 22

22	39152 Y-connector with three tubes for internal/external cooling of contra angles 20:1 LED with push button system		3	40x5x1 Length of the tube 40 Length of the connector 5 Diameter of the connector 1
23	1029 Suitcase		5000	Length x Width x Height 520x440x240
24	1898d SANICAR Mobile Device Trolley with two shelves and 5 m long cable		18500	Height 1125 The upper panel is 450 mm X 360 The middle panel is 344 mm X 360 The foundation is: 470 X 445 Angle of panel slope: No



Table. 2. Contra angles

Contra angles with dental handle, push button collet and INTRA type connection	Reference number	Transmission ratio	Length (mm) up to	Shank diameter	Connection length (for drill/file)	Speed range (revolutions/min)	Max. torque (Ncm)	Weight (grams) up to	Material type
	5053nou Drill contra angle, 20:1, shank diameter of $\varnothing 2.35$ mm, with internal/external cooling, and pushbutton clamp	20:1	116	$\varnothing 2,35$ mm	13 mm	15-2100	70	90	Stainless steel 1.4305.
	5201nou Drill contra angle, 32:1, shank diameter of $\varnothing 2.35$ mm, with internal/external cooling, and pushbutton clamp	32:1	112	$\varnothing 2,35$ mm	13 mm	10-1200	55	90	Stainless steel 1.4305.
	5052nou Drill contra angle, 20:1 LED, shank diameter	20:1	100	$\varnothing 2,35$ mm	13 mm	15-2100	70	100	Stainless steel 1.4305.

	of $\varnothing$ 2.35 mm, with internal/ external cooling, and pushbutton clamp								
---	---	--	--	--	--	--	--	--	--

Date: 15.12.2017
Leszek Kazana, Area Sales Manager

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