Lesley Gjelsmup Kristensen
Regulatory Affairs Manager
Ferring Medical Devices AS
Lesley Gjelsmup Kristensen

ETHICON™

SURGIFLO™

HAEMOSTATIC MATRIX

<i>Instructions for Use</i>	<i>Instrukcja użytkowania</i>
<i>Návod k použití</i>	<i>Instruções de utilização</i>
<i>Brugsanvisning</i>	<i>Instrucțiuni de utilizare</i>
<i>Gebruiksaanwijzing</i>	<i>Návod na použitie</i>
<i>Kasutusjuhend</i>	<i>Navodila za uporabo</i>
<i>Käyttöohjeet</i>	<i>Instrucciones de uso</i>
<i>Notice d'utilisation</i>	<i>Bruksanvisning</i>
<i>Gebrauchsanweisung</i>	<i>Kullanım Talimatları</i>
<i>Οδηγίες χρήσης</i>	<i>Инструкция по применению</i>
<i>Használati utasítás</i>	<i>Upute za upotrebu</i>
<i>Istruzioni per l'uso</i>	<i>تعليمات الاستخدام</i>
<i>Lietošanas instrukcija</i>	<i>Hướng dẫn sử dụng</i>
<i>Naudojimo instrukcija</i>	<i>使用説明</i>
<i>Bruksanvisning</i>	<i>使用说明</i>

4547627

REF **MS0010**

4547627_MASTER_V1_2016.04.21.indd 1

4/21/16 12:49 PM

BNO
FEMD MATERIAL NO: 4547627
PRODUCT CODE / REORDER NUMBER: MS0010
VERSION: MASTER V1
DATE: 2016.04.21
PRODUCT NAME: SURGIFLO™ Haemostatic Matrix
PACKAGE TYPE: IFU
LANGUAGES: English, Czech, Danish, Dutch, Estonian, Finnish, French, German, Greek, Hungarian, Italian, Latvian, Lithuanian, Norwegian, Polish, Portuguese, Romanian, Slovak, Slovene, Spanish, Swedish, Turkish, Russian, Croatian, Arabic, Vietnamese, Traditional Chinese, Simplified Chinese
SIZE: Folded: 114.3 mm (W) x 241.3 mm (H) / 4.50 in (W) x 9.50 in (H)
COLORS: PMS 293 (Blue)
BARCODE(S): Type: 2D DATAMATRIX Content: 4547627

2D Datamatrix Scan:

Scan: **4547627**

Note: Scan produced in this document has been produced by Handheld Products 4600g scan gun and reflects the exact data encoded.

MASTER
2016.04.21

SURGIFLO™ Haemostatic Matrix

Do not inject into blood vessels.

PRODUCT DESCRIPTION

SURGIFLO™ Haemostatic Matrix (SURGIFLO™) is intended for haemostatic use by applying to a bleeding surface.

SURGIFLO™ comes in a sterile tray with *all* sterile components for preparing the Flowable Gelatin Matrix:

- A sterile pre-filled blue plunger syringe containing the porcine gelatin matrix that is off-white in appearance
- A sterile empty syringe
- A sterile liquid transfer cup
- A sterile blue flexible applicator tip that is bendable in all directions
- A sterile white applicator tip that can be trimmed to desired length

2 mL of sterile sodium chloride (sterile saline) solution must be added to the Gelatin Matrix prior to use.

Once the haemostatic matrix is mixed with sterile saline, the appropriate applicator tip should be attached to the syringe for product delivery onto the bleeding site.

ACTIONS

SURGIFLO™ has haemostatic properties. The Flowable Gelatin Matrix provides an environment for platelets to adhere and aggregate within, building on the patient's natural coagulation cascade. When used in appropriate amounts, SURGIFLO™ is absorbed completely within 4-6 weeks. In an animal implantation study, tissue reactions were classified as minimal.

INTENDED USE / INDICATIONS

SURGIFLO™ is indicated for surgical procedures (except ophthalmic) for haemostasis, when control of capillary, venous and arteriolar bleeding by pressure, ligature, and other conventional procedures is ineffective or impractical.

CONTRAINDICATIONS

- Do not inject or compress SURGIFLO™ into blood vessels. Do not use SURGIFLO™ in intravascular compartments because of the risk of thromboembolism, disseminated intravascular coagulation, and increased risk of anaphylactic reaction.
- Do not use SURGIFLO™ in patients with known allergies to porcine gelatin.
- Do not use SURGIFLO™ in closure of skin incisions because it may interfere with the healing of skin edges. This interference is due to mechanical interposition of gelatin and is not secondary to intrinsic interference with wound healing.

WARNINGS

- Do not inject or compress SURGIFLO™ into blood vessels. For episisional use only.
- Do not apply SURGIFLO™ in the absence of active blood flow, e.g. while the vessel is clamped or bypassed due to the risk of intravascular clotting from intravascular injection.
- SURGIFLO™ is not intended as a substitute for meticulous surgical technique and the proper application of ligatures or other conventional procedures for haemostasis.
- SURGIFLO™ should not be used in the presence of infection. SURGIFLO™ should be used with caution in contaminated areas of the body. If signs of infection or abscess develop where SURGIFLO™ has been positioned, reoperation may be necessary in order to remove or drain the infected material.
- SURGIFLO™ should not be used in instances of pumping arterial haemorrhage. It should not be used where blood or other fluids have pooled or in cases where the point of haemorrhage is submerged. SURGIFLO™ will not act as a tampon or plug in a bleeding site.
- SURGIFLO™ should be removed from the site of application when used in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm. Care should be exercised to avoid overpacking. SURGIFLO™ may swell creating the potential for nerve damage.
- Excess SURGIFLO™ should be removed once haemostasis has been achieved because of the possibility of dislodgement of the device or compression of other nearby anatomic structures.
- The safety and effectiveness of SURGIFLO™ for use in ophthalmic procedures has not been established.
- SURGIFLO™ should not be used for controlling post-partum intrauterine bleeding or menorrhagia.

- The safety and effectiveness of SURGIFLO™ has not been established in children and pregnant women.
- The blue flexible applicator tip should not be trimmed to avoid exposing internal guidewire.
- The white straight applicator tip should be trimmed away from the surgical area. Cut a square angle to avoid creating a sharp tip.

PRECAUTIONS

- SURGIFLO™ is for single use only. Do not resterilize. Unused opened SURGIFLO™ should be discarded.
- While packing a cavity for haemostasis is sometimes surgically indicated, SURGIFLO™ should not be used in this manner unless excess product that is not needed to maintain haemostasis is removed. When confined in a clot, SURGIFLO™ may swell approximately 20% upon contact with additional fluid.
- Only the minimum amount of SURGIFLO™ needed to achieve haemostasis should be used. Once haemostasis is achieved, any excess SURGIFLO™ should be carefully removed.
- SURGIFLO™ should not be used in conjunction with autologous blood salvage circuits. It has been demonstrated that fragments of collagen-based haemostatic agents may pass through 40 µm transfusion filters of blood scavenging systems.
- SURGIFLO™ should not be used in conjunction with methylmethacrylate adhesives. Microfibrillar collagen has been reported to reduce the strength of methylmethacrylate adhesives used to attach prosthetic devices to bone surfaces.
- SURGIFLO™ should not be used for the primary treatment of coagulation disorders.
- SURGIFLO™ should be used with caution in neurosurgery as its use has not been studied through randomised, controlled clinical studies.
- Although the safety and effectiveness of the combined use of SURGIFLO™ with other agents has not been evaluated in controlled clinical trials, if in the physician's judgement, concurrent use of other agents is medically advisable, the product literature for that agent should be consulted for complete Prescribing Information.
- The safety and effectiveness of the combined use of SURGIFLO™ with antibiotic solutions or powders has not been established.
- The safety and effectiveness for use in urological procedures have not been established through a randomized clinical study.
- In urological procedures, SURGIFLO™ should not be left in the renal pelvis, renal calyces, bladder, urethra, or ureters to eliminate the potential foci for calculus formation.

Gelatin-Based Haemostatic Agents: Reported Adverse Events

In general, the following adverse events have been reported with the use of absorbable porcine gelatin-based haemostatic agents:

- Gelatin-based haemostatic agents may serve as a nidus for infection and abscess formation and have been reported to potentiate bacterial growth.
- Giant cell granulomas have been observed at implant sites when used in the brain.
- Compression of the brain and spinal cord resulting from the accumulation of sterile fluid has been observed.
- Multiple neurologic events were reported when absorbable gelatin-based haemostatic agents were used in laminectomy operations, including *cauda equina* syndrome, spinal stenosis, meningitis, arachnoiditis, headaches, paresthesias, pain, bladder and bowel dysfunction, and impotence.
- The use of absorbable gelatin-based haemostatic agents during the repair of dural defects associated with laminectomy and craniotomy operations has been associated with fever, infection, leg paresthesias, neck and back pain, bladder and bowel incontinence, *cauda equina* syndrome, neurogenic bladder, impotence, and paresis.
- The use of absorbable gelatin-based haemostatic agents has been associated with paralysis, due to device migration into foramina in the bone around the spinal cord, and blindness due to device migration in the orbit of the eye, during lobectomy, laminectomy, and repair of a frontal skull fracture and lacerated lobe.
- Foreign body reactions, "encapsulation" of fluid, and haematoma have been observed at implant sites.
- Excessive fibrosis and prolonged fixation of a tendon have been reported when absorbable gelatin-based sponges were used in severed tendon repair.
- Toxic shock syndrome was reported in association with the use of absorbable gelatin-based haemostats in nasal surgery.
- Fever, failure of absorption, and hearing loss have been observed when absorbable haemostatic agents were used during tympanoplasty.

HOW SUPPLIED

SURGIFLO™ comes in a sterile tray with all sterile components:

- A sterile pre-filled syringe with blue plunger containing the sterile Gelatin Matrix. The syringe is labeled **SURGIFLO™ Hemostatic Matrix**
- A sterile empty syringe
- A sterile liquid transfer cup
- A sterile blue applicator tip that is bendable in all directions
- A sterile white applicator tip that can be trimmed to desired length

The sterile Flowable Gelatin Matrix and all the accessories in the sterile tray are sterilized by gamma irradiation.

SURGIFLO™ also contains Instructions for Use and tracking labels.

STORAGE AND HANDLING

- SURGIFLO™ should be stored dry at controlled room temperature (2°C-25°C).
- SURGIFLO™ is for single use only.
- The Flowable Gelatin Matrix may be used up to eight (8) hours after mixing with the sterile saline.

DIRECTIONS FOR USE

Before use:

Inspect the SURGIFLO™ packaging for signs of damage. If the packaging is damaged, opened, or wet, sterility cannot be assured and the contents should not be used.

Unused opened packages of SURGIFLO™ should be discarded since they are not intended for reuse and/or re-sterilization.


Opening the Flowable Gelatin Matrix tray:

Open the outer package and deliver the sterile inner tray to the sterile field using aseptic technique. Once placed in the sterile field, the sterile inner tray may be opened.

Preparing the Flowable Gelatin Matrix with the sterile saline in the sterile field:


1) Draw 2 mL of sterile saline from the sterile liquid transfer cup into the empty sterile syringe.

2) Connect syringes

Remove the blue cap from the end of the sterile pre-filled syringe with the blue plunger containing the Gelatin Matrix. Attach this syringe to the sterile syringe containing the sterile saline.

3) Mix contents of the 2 syringes

Begin mixing by transferring the sterile saline into the sterile pre-filled syringe containing the Gelatin Matrix. Push the combined material back and forth 6 times until the consistency is even.

Once mixed, the haemostatic matrix should reside completely in the syringe with the blue plunger that is labeled **SURGIFLO™ Hemostatic Matrix**. Remove the empty syringe and discard.

4) Attach applicator tip

The product is now ready for clinical use.

- The blue flexible applicator tip is bendable in all directions. (Do not trim the flexible applicator tip to avoid exposing the internal guidewire.)
- The white applicator tip can be trimmed to desired length. The tip should be trimmed away from the surgical area. Cut a square angle to avoid creating a sharp tip. The tray can be used to contain the excess piece(s) for discarding.

Do not inject SURGIFLO™ into blood vessels. See the Contraindications, Warnings, and Precautions.

For open procedures

- Identify the source of bleeding.
- Deliver SURGIFLO™ to the source of bleeding. SURGIFLO™ can be used with or without one of the applicator tips attached to the syringe labeled **SURGIFLO™ Hemostatic Matrix**. Apply sufficient SURGIFLO™ to cover the entire bleeding surface.
- For tissue defects (cavities, divots, or craters), apply SURGIFLO™ at the deepest part of the lesion and continue applying material as the syringe (or applicator tip) is withdrawn from the lesion.
- Apply a sterile saline moistened gauze over the SURGIFLO™ to ensure the material remains in contact with the bleeding tissue.
- After 1-2 minutes, lift and remove the gauze and inspect the wound site. Once bleeding has ceased, irrigate excess SURGIFLO™ away gently so as not to disturb the new clot.
- In cases of persistent bleeding indicated by saturation and bleeding through the material, repeat application of SURGIFLO™.

For endoscopic sinus surgery and epistaxis:

- Deliver SURGIFLO™ to the source of bleeding using the selected applicator tip attached to the **SURGIFLO™ Hemostatic Matrix** syringe.
- Apply sufficient SURGIFLO™ to cover the entire bleeding surface.
- Using forceps or an appropriate instrument, carefully layer a sterile saline moistened gauze over the SURGIFLO™ for 1-2 minutes to ensure the material remains in contact with the bleeding tissue.
- In cases of persistent bleeding, indicated by saturation and bleeding through the material, insert the applicator tip through the center of the mass of previously placed SURGIFLO™ to deliver fresh material as close as possible to the tissue surface. After re-application of SURGIFLO™, use a sterile saline moistened gauze to approximate the material to the tissue for another minute, and then inspect the site. Repeat re-application, if necessary.
- Once haemostasis has been achieved, remove the gauze. If possible, excess SURGIFLO™ should be removed with gentle irrigation or careful suction. Avoid disrupting the clot complex. The remaining SURGIFLO™ does not have to be removed, as it will be bioresorbed.
- Use of nasal packing is not necessary when satisfactory haemostasis is achieved.
- If necessary, gentle irrigation and/or careful suction can be used in the postoperative period to remove the remaining SURGIFLO™.

SYMBOLS USED ON LABELING

	Store at controlled room temperature 2°C- 25°C (36°F- 77°F)		Sterilized using irradiation.
	See instructions for use.		Distributed by.
	Do not reuse.		Manufacturer.
	Do not resterilize.		Re-order number.
	Do not use if package is damaged or opened.		Use by: year, month and day.
	Do not inject into blood vessels.		Date of manufacture: year, month and day.
	Designates the packaging material to which it is applied is recyclable. Recycling programs may not exist in your area.		Batch number.
			CE-mark and identification number of notified body.
	Device components are not made with natural rubber latex.		

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SURGIFLO™ Haemostatic Matrix

Area of application

It is used only in medical and preventive health care institutions.

Disposal

Disposal must be carried out in accordance with accepted medical practice standards, and relevant local, state, and Federal legislative acts and regulations according to SanPiN 2.1.7.2790-10, class B.

Sterilization

See paragraph «How Supplied».

Kit sterility term covers the shelf life of each constituent of the kit

Note: For single patient use only. Re-sterilization and re-use are prohibited.

Shelf life

Shelf life for kit is determined according to the shortest constituent expiration date, see kit box.

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FerroSan Medical Devices A/S
Lestley Gjesbump Kristensen
Regulatory Affairs Manager
[Signature]
2017.02.27

[Handwritten signature]

Manufacturer's warranty

Manufacturer is not liable for any damage or accidents resulting from improper use, human error or external emergencies. Regarding the product quality, please contact an authorized manufacturer's representative.

Authorized manufacturer's representative

LLC Johnson & Johnson
Krylatskaya 17, Building 2, 121614 Moscow, Russia
tel. (495)5807777, fax (495)5807878.

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INTERPRETATION OF ADDITIONAL SYMBOLS ARE NOT NOTICED IN ORIGINAL IFU SYMBOLS

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Product passed procedure of declaration of conformity

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Confederation of Danish Industry
I hereby certify that the above
company is a member of our
organization and known to us as
worthy of confidence.

CONFEDERATION OF DANISH INDUSTRY
(DI)
MEMBER OF INDUSTRY AND COMMERCE SECTION

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© Ethicon, Inc. 2011