Ferrosan Medical Devices A/S  
Lestey Gelstrup Kristensen  
Regulatory Affairs Manager

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4547974\_DRAFT\_V2\_2016.12.16.indd 1

12/19/16 11:10 AM

| BNO                            |                                                                                                                                                                                                                                     |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| FEMD MATERIAL NO:              | 4547974                                                                                                                                                                                                                             |
| PRODUCT CODE / REORDER NUMBER: | MS0012                                                                                                                                                                                                                              |
| VERSION:                       | V2                                                                                                                                                                                                                                  |
| DATE:                          | 2016.12.16                                                                                                                                                                                                                          |
| PRODUCT NAME:                  | SURGIFLO™ Haemostatic Matrix Kit with Thrombin                                                                                                                                                                                      |
| PACKAGE TYPE:                  | IFU                                                                                                                                                                                                                                 |
| LANGUAGES:                     | English, Croatian, Czech, Danish, Dutch, Estonian, Finnish, French, German, Greek, Hungarian, Italian, Latvian, Lithuanian, Norwegian, Polish, Portuguese, Romanian, Russian, Slovak, Slovene/Slovenian, Spanish, Swedish & Turkish |
| SIZE:                          | 250 mm (H) x 167 mm (W)                                                                                                                                                                                                             |
| COLOURS:                       | PMS 293 (Blue)                                                                                                                                                                                                                      |
| BARCODE(S):                    | Type: 2D DATAMATRIX<br>Content: 4547974<br>Size: 8.127 mm x 8.127 mm                                                                                                                                                                |

2D Datamatrix Scan:

Scan: **4547974**

Note: Scan produced in this document has been produced by Handheld Products 4600g scan gun and reflects the exact data encoded.

**DRAFT**  
2016.12.16

|             |            |
|-------------|------------|
| English     | Page 3     |
| Hrvatski    | Stranica 8 |
| Česky       | Strana 14  |
| Dansk       | Side 20    |
| Nederlands  | Pagina 26  |
| Eesti       | Lk 32      |
| Suomi       | Sivu 38    |
| Français    | Page 44    |
| Deutsch     | Seite 50   |
| Ελληνικά    | Σελίδα 56  |
| Magyar      | Oldal 62   |
| Italiano    | Pagina 68  |
| Latviešu    | 74.lpp     |
| Lietuvių    | Psl. 80    |
| Norsk       | Side 86    |
| Polski      | Strona 92  |
| Português   | Página 98  |
| Română      | Pagina 104 |
| Русский     | Стр. 110   |
| Slovenský   | Strana 116 |
| Slovenščina | Stran 122  |
| Español     | Página 128 |
| Svenska     | Sida 134   |
| Türkçe      | Sayfa 140  |

## SURGIFLO™ Haemostatic Matrix Kit with Thrombin

**Do not inject into blood vessels.**

### PRODUCT DESCRIPTION

SURGIFLO™ Haemostatic Matrix Kit with Thrombin (SURGIFLO™) is intended for haemostatic use by applying to a bleeding surface.

The kit contains:

1. A sterile tray with *all* sterile components to prepare the Flowable Gelatin Matrix
2. A sterile tray with *all* surface sterilized Thrombin kit components to prepare the Thrombin Solution
  1. The Flowable Gelatin Matrix comes in a tray with *all* sterile components:
    - A sterile pre-filled blue plunger syringe containing the porcine Gelatin Matrix that is off-white in appearance
    - A sterile empty syringe
    - A sterile liquid transfer cup
    - A sterile blue flexible applicator tip that is bendable in all directions
    - A sterile white applicator tip that can be trimmed to desired length
  2. The surface sterilized components to prepare the Thrombin Solution:
    - A Thrombin vial containing 2000 International Units (IU) of sterile lyophilized human thrombin
    - A needle-free syringe containing 2 mL of Sterile Water for Injection (Sterile WFI)
    - A sterile vial adapter

Thrombin should be reconstituted using the vial adapter and the needle-free syringe with Sterile WFI.

The Thrombin Solution must be added to the Flowable Gelatin Matrix prior to use.

Once the haemostatic matrix is mixed with the Thrombin Solution, the appropriate applicator tip should be attached to the syringe for product delivery onto the bleeding site.

### ACTIONS

SURGIFLO™ has haemostatic properties. The Flowable Gelatin Matrix provides an environment for platelets to adhere and aggregate within, building on the patient's natural coagulation cascade.

The patient's endogenous thrombin is activated, and the patient's thrombin converts fibrinogen into an insoluble fibrin clot. The thrombin component of SURGIFLO™ provides an ancillary effect to the innate haemostatic property of the Flowable Gelatin Matrix.

When used in appropriate amounts, SURGIFLO™ is absorbed completely within 4-6 weeks. In an animal implantation study, tissue reactions were classified as minimal.

### INTENDED USE / INDICATIONS

SURGIFLO™ is indicated in surgical procedures (other than ophthalmic) as an adjunct to haemostasis when control of bleeding ranging from oozing to spurting by ligature or other conventional methods is ineffective or impractical.

### CONTRAINDICATIONS

- Do not inject or compress SURGIFLO™ into blood vessels. Do not use SURGIFLO™ in intravascular compartments because of the risk of thromboembolism, disseminated intravascular coagulation, and increased risk of anaphylactic reaction.
- Do not use SURGIFLO™ in patients known to have anaphylactic or severe systemic reaction to human blood products.
- Do not use SURGIFLO™ in patients with known allergies to porcine gelatin.
- Do not use SURGIFLO™ in closure of skin incisions because it may interfere with the healing of skin edges. This interference is due to mechanical interposition of gelatin and is not secondary to intrinsic interference with wound healing.

### WARNINGS

- Do not inject or compress SURGIFLO™ into blood vessels. For epilesional use only.
- Do not apply SURGIFLO™ in the absence of active blood flow, e.g. while the vessel is clamped or bypassed due to the risk of intravascular clotting from intravascular injection.
- SURGIFLO™ contains thrombin made from human plasma. Products made from human plasma may carry a risk of transmitting infectious agents, such as viruses, and theoretically, the Creutzfeldt-Jakob disease (CJD) agent. The risk of transmitting an infectious agent has been reduced by screening plasma donors for prior exposure to certain viruses, by testing for the presence of certain current virus infections, and by inactivating and removing certain viruses. Despite these measures, such products can still potentially transmit disease. There is also the possibility that unknown infectious agents may be present in such products. The physician should discuss the risks and benefits of this product with the patient. The measures taken are considered effective for enveloped viruses such as HIV, HCV, HBV, and for the non-enveloped virus HAV. The measures taken may be of limited value against non-enveloped viruses such as parvovirus B19. Parvovirus B19 infection may be serious for pregnant women (foetal infection) and for individuals with immunodeficiency or increased erythropoiesis (e.g. haemolytic anaemia).
- SURGIFLO™ is not intended as a substitute for meticulous surgical technique and the proper application of ligatures or other conventional procedures for haemostasis.

- SURGIFLO™ should not be used in the presence of infection. SURGIFLO™ should be used with caution in contaminated areas of the body. If signs of infection or abscess develop where SURGIFLO™ has been positioned, reoperation may be necessary in order to remove or drain the infected material.
- SURGIFLO™ should not be used in instances of pumping arterial haemorrhage. It should not be used where blood or other fluids have pooled or in cases where the point of haemorrhage is submerged. SURGIFLO™ will not act as a tampon or plug in a bleeding site.
- SURGIFLO™ should be removed from the site of application when used in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm. Care should be exercised to avoid overpacking. SURGIFLO™ may swell creating the potential for nerve damage.
- Excess SURGIFLO™ should be removed once haemostasis has been achieved because of the possibility of dislodgement of the device or compression of other nearby anatomic structures.
- The safety and effectiveness of SURGIFLO™ for use in ophthalmic procedures has not been established.
- SURGIFLO™ should not be used for controlling post-partum intrauterine bleeding or menorrhagia.
- The safety and effectiveness of SURGIFLO™ has not been established in children.
- The safety and effectiveness of SURGIFLO™ has not been established in pregnant women and nursing mothers.
- The blue flexible applicator tip should not be trimmed to avoid exposing the internal guidewire.
- The white straight applicator tip should be trimmed away from the surgical area. Cut a square angle to avoid creating a sharp tip.

#### PRECAUTIONS

- SURGIFLO™ is for single use only. Do not resterilize. Unused opened SURGIFLO™ should be discarded.
- While packing a cavity for haemostasis is sometimes surgically indicated, SURGIFLO™ should not be used in this manner unless excess product that is not needed to maintain haemostasis is removed. When confined in a clot, SURGIFLO™ may swell approximately 20% upon contact with additional fluid.
- Only the minimum amount of SURGIFLO™ needed to achieve haemostasis should be used. Once haemostasis is achieved, any excess SURGIFLO™ should be carefully removed.
- SURGIFLO™ should not be used in conjunction with autologous blood salvage circuits. It has been demonstrated that fragments of collagen-based haemostatic agents may pass through 40 µm transfusion filters of blood scavenging systems.
- SURGIFLO™ should not be used in conjunction with methylmethacrylate adhesives. Microfibrillar collagen has been reported to reduce the strength of methylmethacrylate adhesives used to attach prosthetic devices to bone surfaces.
- Similar to comparable products containing thrombin, the Thrombin Solution may be denatured after exposure to solutions containing alcohol, iodine or heavy metals (e.g. antiseptic solutions). Such substances should be removed to the greatest possible extent before applying the product.
- SURGIFLO™ should not be used for the primary treatment of coagulation disorders.
- As with other gelatin haemostatic agents, in urological procedures SURGIFLO™ should not be left in the renal pelvis, renal calyces, bladder, urethra or ureters to eliminate the potential foci for calculus formation. The safety and effectiveness of SURGIFLO™ for use in urological procedures has not been established through a randomised clinical trial.
- As with other gelatin haemostatic agents that swell, SURGIFLO™ should be used with caution in neurosurgery. Safe and effective use of SURGIFLO™ in neurosurgery has not been proven through randomised, controlled studies.
- Although the safety and effectiveness of the combined use of SURGIFLO™ with other agents have not been evaluated in controlled clinical trials, if in the physician's judgement, concurrent use of other agents is medically advisable, the product literature for that agent should be consulted for complete Prescribing Information.
- The safety and effectiveness of the combined use of SURGIFLO™ with antibiotic solutions or powders has not been established.

#### Gelatin-Based Haemostatic Agents: Reported Adverse Events

In general, the following adverse events have been reported with the use of absorbable porcine gelatin-based haemostatic agents:

- Gelatin-based haemostatic agents may serve as a nidus for infection and abscess formation and have been reported to potentiate bacterial growth.
- Giant cell granulomas have been observed at implant sites when used in the brain.
- Compression of the brain and spinal cord resulting from the accumulation of sterile fluid has been observed.
- Multiple neurologic events were reported when absorbable gelatin-based haemostatic agents were used in laminectomy operations, including *cauda equina* syndrome, spinal stenosis, meningitis, arachnoiditis, headaches, paresthesias, pain, bladder and bowel dysfunction, and impotence.
- The use of absorbable gelatin-based haemostatic agents during the repair of dural defects associated with laminectomy and craniotomy operations has been associated with fever, infection, leg paresthesias, neck and back pain, bladder and bowel incontinence, *cauda equina* syndrome, neurogenic bladder, impotence, and paresis.
- The use of absorbable gelatin-based haemostatic agents has been associated with paralysis, due to device migration into foramina in the bone around the spinal cord, and blindness due to device migration in the orbit of the eye, during lobectomy, laminectomy, and repair of a frontal skull fracture and lacerated lobe.
- Foreign body reactions, "encapsulation" of fluid, and haematoma have been observed at implant sites.
- Excessive fibrosis and prolonged fixation of a tendon have been reported when absorbable gelatin-based sponges were used in severed tendon repair.
- Toxic shock syndrome was reported in association with the use of absorbable gelatin-based haemostats in nasal surgery.
- Fever, failure of absorption, and hearing loss have been observed when absorbable haemostatic agents were used during tympanoplasty.

#### Adverse Reactions to Human Thrombin

As with any other plasma derivatives, hypersensitivity or allergic reactions may occur in rare cases. In isolated cases, these reactions may progress to severe anaphylaxis. Other adverse events reported in a clinical trial were abnormal laboratory tests (prolonged activated partial thromboplastin time, prolonged prothrombin time, increased INR, decreased lymphocyte count, increased neutrophil count) and haematoma.

### Adverse Reactions to Gelatin-Based Haemostatic Agents with Thrombin

Adhesion formation and small bowel obstruction are well-known and common complications following abdominal and gynaecological surgeries. Adverse events such as these, as well as inflammation and foreign body reaction including giant cell granulomas, have been reported within abdominal and gynaecological surgeries in which gelatin-based and collagen-based haemostatic agents with thrombin have been used in excessive amount.

Therefore, as with other gelatin haemostatic agents, only the minimum amount of SURGIFLO™ needed to achieve haemostasis should be used. Once haemostasis is achieved, any excess SURGIFLO™ should be carefully removed.

### HOW SUPPLIED

SURGIFLO™ consists of:

1. A sterile tray with *all* sterile components to prepare the Flowable Gelatin Matrix
2. A sterile tray with *all* surface sterilized Thrombin kit components to prepare the Thrombin Solution

SURGIFLO™ is provided in the configuration shown in the table below.

| SURGIFLO™ Haemostatic Matrix Kit with Thrombin                                                                                                                                                                                                                                       |                                                                                                                                                                                                                                                                                |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Flowable Gelatin Matrix Components                                                                                                                                                                                                                                                   | Thrombin Components                                                                                                                                                                                                                                                            |
| <ul style="list-style-type: none"><li>• A sterile pre-filled blue plunger syringe with the porcine gelatin matrix</li><li>• A sterile empty syringe</li><li>• A sterile liquid transfer cup</li><li>• A sterile blue flexible tip</li><li>• A sterile white applicator tip</li></ul> | <ul style="list-style-type: none"><li>• A Thrombin vial containing 2000 International Units (IU) of sterile lyophilized human thrombin</li><li>• A needle-free syringe containing 2 mL of Sterile Water for Injection (Sterile WFI)</li><li>• A sterile vial adapter</li></ul> |

#### The sterile Flowable Gelatin Matrix and the accessories:

The tray is sterilized by gamma irradiation.

#### The Thrombin kit components:

The tray is surface sterilized with ethylene oxide.

- The lyophilized thrombin (human) is sterilized through an autoclaved sterile filter cartridge.
- The Sterile Water for Injection (Sterile WFI) in the needle-free syringe is steam sterilized.
- The vial adapter is sterilized by gamma irradiation.

SURGIFLO™ contains Instructions for Use. Tracking labels are provided in the kit to record name and batch number of the product to link use to the patient record. It is strongly recommended that every time that SURGIFLO™ is administered to a patient, the name and batch number of the product are recorded in order to maintain a link between the patient and the batch of the product.

### STORAGE AND HANDLING

- SURGIFLO™ should be stored dry at controlled room temperature (2°C–25°C).
- SURGIFLO™ is for single use only.
- The Thrombin vial should be kept away from light.
- Do not freeze or refrigerate Thrombin once it has been reconstituted.
- The Thrombin Solution should be used together with the Flowable Gelatin Matrix and used only as indicated.
- The Flowable Gelatin Matrix may be used up to 8 hours after mixing with the Thrombin Solution.

### DIRECTIONS FOR USE

#### Before use:

Inspect the SURGIFLO™ trays for signs of damage. If the packaging has been damaged, opened, or wet, sterility cannot be assured and the contents should not be used.

Unused opened packages of SURGIFLO™ should be discarded since they are not intended for reuse and/or re-sterilization.

#### Opening the tray with Flowable Gelatin Matrix and the tray with Thrombin kit components:


Open the outer packages and deliver the sterile inner trays to the sterile field using aseptic technique. Once placed in the sterile field, the sterile inner tray may be opened.

#### Preparing the Thrombin Solution inside the sterile field:


Flip off the cap from the Thrombin vial, leaving the aluminum ring and the rubber stopper in place. Peel off the lid from the vial adapter package.

1. Place the Thrombin vial on a flat surface, seat the vial adapter on the centre of the rubber stopper and push down until the spike penetrates the rubber stopper and the vial adapter snaps into place. Remove the blister package.


2. Snap off the tamper cap on the needle-free syringe containing the Sterile Water for Injection (Sterile WFI).

3. Connect and screw on the *needle-free syringe* to the *vial adapter*. Transfer the entire Sterile WFI into the *Thrombin vial*.

4. Gently swirl the *Thrombin vial* until the *Thrombin Solution* is clear.

5. Draw up the *Thrombin Solution* into the *needle-free syringe*. Label the *needle-free syringe*: "Thrombin 2000 IU".

6. Disconnect the *needle-free syringe* from the *vial adapter* and transfer the *Thrombin Solution* into the sterile liquid transfer cup as shown in the next section (Figure 1).

After reconstitution, discard the components used for the thrombin reconstitution.


Alternatively, the Thrombin may be reconstituted **outside** the sterile field. Be careful not to touch the rubber stopper of the vial. After reconstitution the Thrombin Solution should be transferred into the sterile liquid transfer cup using aseptic technique.

Place the sterile liquid transfer cup near the edge of the sterile field to receive the Thrombin Solution transfer without contaminating the sterile field.

**Preparing the Flowable Gelatin Matrix with the Thrombin Solution inside the sterile field:**


1) Draw *Thrombin Solution* from sterile liquid transfer cup into the empty sterile syringe.

2) **Connect syringes**


Remove the blue cap from the end of the sterile pre-filled syringe with the blue plunger containing the Gelatin Matrix. Attach this syringe to the sterile syringe containing the Thrombin Solution.

3) **Mix contents of the 2 syringes**

Begin mixing by transferring the sterile Thrombin Solution into the sterile pre-filled syringe containing the Gelatin Matrix. Push the combined material back and forth 6 times until the consistency is even.

Once mixed, the haemostatic matrix should reside completely in the syringe with the blue plunger that is labeled **SURGIFLO™ Hemostatic Matrix**. Remove the empty syringe and discard.

4) **Attach applicator tip**

The product is now ready for clinical use.

- The blue flexible applicator tip is bendable in all directions. (Do not trim the flexible applicator tip to avoid exposing the internal guidewire).
- The white applicator tip can be trimmed to desired length. The tip should be trimmed away from the surgical area. Cut a square angle to avoid creating a sharp tip. The tray can be used to contain the excess piece(s) for discarding.

Do not inject SURGIFLO™ into blood vessels. See the Contraindications, Warnings and Precautions.

**For open procedures:**

- Identify the source of bleeding.
- Deliver SURGIFLO™ to the source of bleeding. SURGIFLO™ can be used with or without one of the applicator tips attached to the syringe labeled **SURGIFLO™ Hemostatic Matrix**. Apply sufficient SURGIFLO™ to cover the entire bleeding surface.
- For tissue defects (cavities, divots or craters) apply SURGIFLO™ at the deepest part of the lesion, and continue applying material as the syringe (or applicator tip) is withdrawn from the lesion.
- Apply a sterile saline moistened gauze over the SURGIFLO™ to ensure the material remains in contact with the bleeding tissue.
- After 1-2 minutes, lift and remove the gauze and inspect the wound site. Once bleeding has ceased, irrigate excess SURGIFLO™ away gently so as not to disturb the new clot.
- In cases of persistent bleeding indicated by saturation and bleeding through the material, repeat application of SURGIFLO™.

**For endoscopic sinus surgery and epistaxis:**

- Deliver SURGIFLO™ to the source of bleeding using the selected applicator tip attached to the **SURGIFLO™ Hemostatic Matrix** syringe.
- Apply sufficient SURGIFLO™ to cover the entire bleeding surface.
- Using forceps or an appropriate instrument, carefully layer a sterile saline moistened gauze over the SURGIFLO™ for 1-2 minutes to ensure the material remains in contact with the bleeding tissue.
- In cases of persistent bleeding, indicated by saturation and bleeding through the material, insert the applicator tip through the center of the mass of previously placed SURGIFLO™ to deliver fresh material as close as possible to the tissue surface. After re-application of SURGIFLO™, use a sterile saline moistened gauze to approximate the material to the tissue for another minute, and then inspect the site. Repeat re-application, if necessary.
- Once haemostasis has been achieved, remove the gauze. If possible, excess SURGIFLO™ should be removed with gentle irrigation or careful suction. Avoid disrupting the clot complex. The remaining SURGIFLO™ does not have to be removed, as it will be bioresorbed.
- Use of nasal packing is not necessary when satisfactory haemostasis is achieved.
- If necessary, gentle irrigation and/or careful suction can be used in the postoperative period to remove the remaining SURGIFLO™.

**SYMBOLS USED ON LABELING**

|                                                                                     |                                                                                                                        |                                                                                     |                                                     |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-----------------------------------------------------|
|    | Store at controlled room temperature 2°C- 25°C (36°F- 77°F).                                                           |    | Sterilized using irradiation.                       |
|    | See instructions for use.                                                                                              |    | Sterilized using aseptic processing techniques.     |
|   | Do not reuse.                                                                                                          |   | Sterilized using steam or dry heat.                 |
|  | Distributed by                                                                                                         |  | Sterilized using ethylene oxide.                    |
|  | Do not resterilize.                                                                                                    |  | Date of manufacture: year, month and day.           |
|  | Do not use if package is damaged or opened.                                                                            |  | Manufacturer.                                       |
|  | Do not inject into blood vessels.                                                                                      |  | Re-order number.                                    |
|  | Device components are not made with natural rubber latex.                                                              |  | Use by: year, month and day.                        |
|  | Designates the packaging material to which it is applied is recyclable. Recycling programs may not exist in your area. |  | Batch number.                                       |
|                                                                                     |                                                                                                                        |  | CE-mark and identification number of notified body. |

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*[Handwritten signature]*

**Addendum to IFU for medical device**

**SURGIFLO™ Haemostatic Matrix Kit with  
Thrombin**

**Area of application**

It is used only in medical and preventive health care institutions.

**Disposal**

After use the medical device is a potential biohazard because after use it can be a potential biological hazard due to contamination with blood and other biological materials.

Disposal must be carried out in accordance with accepted medical practice standards, and relevant local, state, and Federal legislative acts and regulations according to SanPiN 2.1.7.2790-10, class B.

**Sterilization**

See paragraph «How Supplied».

Kit sterility term covers the shelf life of each constituent of the kit.

Note: For single patient use only. Re-sterilization and re-use are prohibited.

**Shelf life for kit**

Shelf life for kit is determined according to the shortest constituent expiration date, see kit box.

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Ferrosan Medical Devices A/S  
Lestry Gjelstrup Kristensen  
Regulatory Affairs Manager  
*Lestry Gjelstrup Kristensen*

2017.02.27

**Ferrosan Medical Devices A/S**  
Sydmarken, 5, DK-2860, Soeborg, Denmark

**Manufacturer's warranty**

Manufacturer is not liable for any damage or accidents resulting from improper use, human error or external emergencies. Regarding the product quality, please contact an authorized manufacturer's representative.

**Authorized manufacturer's representative**

LLC Johnson & Johnson  
Krylatskaya 17, Building 2, 121614 Moscow, Russia  
tel. (495)5807777, fax (495)5807878.

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**INTERPRETATION OF ADDITIONAL SYMBOLS ARE NOT NOTICED IN ORIGINAL IFU SYMBOLS**

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Product passed procedure of declaration of conformity

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