SECRETARY OF STATE
STATE OF INDIANA

Connie Lawson
Secretary of State


APOSTILLE

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1. Country: United States of America
2. This Public Document has been signed *Vicky Chambers*
3. Acting in the capacity of notary public in & for *Owen County*
4. and bears the seal/stamp of certification clerk in & for Sec. of State Office

CERTIFIED

5. at Indianapolis, Indiana
6. this *Seventh* of *June*, 2017
7. by the Secretary of State of Indiana
8. No. A2017-0607105242
9. Seal/Stamp:
10. Signature:

Connie Lawson

Connie Lawson
Indiana Secretary of State

*Effective May 1, 2003 all apostilles from the Indiana Secretary of State will have an electronically printed seal.
This document was revised March 17, 2012.*

COOK INCORPORATED
750 DANIELS WAY
BLOOMINGTON, IN 47404 USA
PHONE: 812.339.2235 TOLL FREE: 800.457.4500
WWW.COOKMEDICAL.COM

CERTIFICATION STATEMENT

As a responsible representative of Cook Incorporated, I confirm this is an accurate version of the Instructions for Use for Ureteral Dilation Sets.

Kris Weathers
Regulatory Affairs Manager
Cook Incorporated

State of Indiana)

ss

County of Owen)

Subscribed and sworn to before me the 6th day of June month 2017 year.

My Commission expires October 13, 2024

Vicky Chambers

INSTRUCTIONS FOR USE

Ureteral Dilation Sets

INTENDED USE


Used for the dilation of various parts of the urinary system (see Instructions for Use).


CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner.)

DEVICE DESCRIPTION

Parameter	Conversion
1 French	1 mm = 1Fr / 3
1 Inch	25,4 mm

Ascend AO Ureteral Dilation Balloon Catheter Set (AUBS)


The Ascend AQ Ureteral Dilation Balloon Catheter Sets are used for transluminal dilation of ureteral strictures or ureteral dilation prior to ureteroscopy or stone manipulation. The sets are comprised of a balloon catheter with batwing manifold, a stop-cock and an inflation device. The balloon catheter is manufactured in diameter 6.0 Fr/ 2.0 mm with a working length of 65 centimeters, the distal tip of the catheter has a diameter of 5 Fr/1.66 mm. The balloon is manufactured in lengths of 4.0, 6.0 and 10 centimeters and with diameters of 4.0, 5.0, 6.0, 7.0 and 10.0 millimeter diameters. The Ascend AQ Ureteral Dilation Balloon Catheter has a maximum inflation pressure of 20 atmospheres ($\pm 10\%$). Nominal balloon pressure – 14 atm. ($\pm 10\%$).

The distal tip of the balloon catheter is characterized by a 0.040 inch/ 1.02 mm (± 0.005 inch/ 0.127 mm) diameter endhole rounded to a smooth tip. The distal tip of the balloon catheter is characterized by a 0.040 inch/ 0.02 mm (± 0.005 inch/ 0.127 mm) diameter endhole rounded to a smooth finish. The distal lumen of the balloon catheter will accept a 0.038 inch/ 0.97 mm diameter wire guide.

Two radiopaque markers are placed on the catheter shaft beneath the balloon. Depending upon the specific balloon catheter, the radiopaque marker bands are spaced 4.0, 6.0 or 10 centimeters apart.

The proximal end of the balloon catheter is characterized by an inflation lumen and a batwing manifold. A high pressure one-way stopcock is attached to the inflation lumen of the batwing manifold.

The inflation device is a one-piece, plastic, disposable, inflation device with a lock lever design that controls the manometer, and a connecting tube with a male rotating adapter. The manometer measures pressures ranging from vacuum to gauge capacity. The gauge is marked in 1 atmosphere (atm) increments. The gauge also has an inner scale of comparable PSI measurements. The inflation device has a 20 cc/30 atm capacity.

Stop-cock parameters:

Length «A» - 1.594 inch/ 40.49mm (± 0.005 inch/0.127mm);

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
Inner diameter «B» - 0.098 inch/ 2.49 mm (± 0.005 inch/0.127mm);


Outer diameter: «D» – 0.430inch/10.92mm (± 0.005 inch/0.127mm);

Luer-lock «Female»:

Inner diameter «C» - 0.169 inch/ 4.3 mm (± 0.005 inch/0.127mm);

Outer diameter «E» - 0.300 inch/7.62mm (± 0.005 inch/0.127mm).


Parameters of batwing manifold:

Lengths: "A" – 2.275 inch/ 57.785mm (± 0.005 inch/0.127mm); "B" – 1.675 inch/ 42.545 mm (± 0.005 inch/0.127mm); "C" – 0.600 inch/15.24mm (± 0.005 inch/0.127mm); "D" – 1.575 inch/40.005 mm (± 0.005 inch/0.127mm).

Diameter: "E" – 0.136 inch/3.454 mm (± 0.005 inch/0.127mm);

Angle "F" -- 40°

Order Number	Catheter			Inflated Balloon	
	Catheter Shaft Diameter Fr/mm	Catheter Tip Diameter Fr/mm	Length (cm)	Diameter (mm)	Length (cm)
Tolerance	$\pm 10\%$	$\pm 10\%$	$\pm 10\%$	$\pm 10\%$	$\pm 10\%$
AUBS-10-4	6,0 / 2,0	5,0/1.66	65	10,0	4,0
AUBS-4-10	6,0 / 2,0	5,0/1.66	65	4,0	10,0
AUBS-4-4	6,0 / 2,0	5,0/1.66	65	4,0	4,0
AUBS-4-6	6,0 / 2,0	5,0/1.66	65	4,0	6,0
AUBS-5-10	6,0 / 2,0	5,0/1.66	65	5,0	10,0

Order Number	Catheter			Inflated Balloon	
	Catheter Shaft Diameter Fr/mm	Catheter Tip Diameter Fr/mm	Length (cm)	Diameter (mm)	Length (cm)
AUBS-5-4	6,0 / 2,0	5,0/1.66	65	5,0	4,0
AUBS-5-6	6,0 / 2,0	5,0/1.66	65	5,0	6,0
AUBS-6-10	6,0 / 2,0	5,0/1.66	65	6,0	10,0
AUBS-6-4	6,0 / 2,0	5,0/1.66	65	6,0	4,0
AUBS-6-6	6,0 / 2,0	5,0/1.66	65	6,0	6,0
AUBS-7-10	6,0 / 2,0	5,0/1.66	65	7,0	10,0
AUBS-7-4	6,0 / 2,0	5,0/1.66	65	7,0	4,0
AUBS-7-6	6,0 / 2,0	5,0/1.66	65	7,0	6,0

Cook Sphere Inflation Device (CID-20-30)

The Cook Sphere Inflation Device is recommended for use while performing balloon dilation procedures to inflate the balloon, monitor the pressure within the balloon and deflate the balloon. This device is supplied sterile, disposable and for single use only.

Parameters of Cook Sphere Inflation Device:

	Volume, cc	Length, mm	Mass, g	Manometer value, ATM	Manometer value of division, ATM	Value of division Syringe, ml	Adapter type
Tolerance	±10%	±10%	±10%	±10%	±10%	±10%	-
Syringe-manometer CID-20-30	20	210	113	30	1	2	Male rotating Luer lock

INTENDED USE/INDICATION FOR USE

Ascend AO Ureteral Dilation Balloon Catheter Set (AUBS-) -used for transluminal dilation of ureteral strictures or ureteral dilation prior to ureteroscopy or stone manipulation. Intended for one-time use.

CONTRAINDICATIONS

1. Acute urethritis
2. Acute prostatitis
3. Acute cystitis
4. Urethral abnormalities (diverticula, doubling, etc.) because of the risk of perforation
5. Uncorrected coagulopathy

Component	Material	Patient Contacting
Batwing manifold	Polyamide 12-Compound # 40960 Red CI Pigment Red 149 CI Pigment Yellow 138 Calcium Carbonate Carbon Black	Indirect contact
Stopcock	Polycarbonate № 80020 Stem: Celcon Grey	Indirect contact
Inflation device	Syringe: Polycarbonate # 111211-39-3 Black Nylon # 32131-17-2 Plunger: Polycarbonate # 111211-39-3 Black Nylon # 32131-17-2 Manometer: ABS Acrylonitrile-Butadiene-Styrene # 9003-56-9 Connecting tube: Polyurethane # 037383-28-1 Connecting tube hub: Polycarbonate # 111211-39-3 Lock Lever : Black Nylon # 26098-55-5 (Half Nut) # 32131-17-2 (Links) Glass Fibers # 65997-17-3 Polytetrafluoroethylene (PTFE) # 9002-84-0	Non-contact

Individual packing

Component	Material	Patient Contacting
Packing	Part of the paper - medical paper TYVEK: 1073B/109P Film part – Nylon BON	Non-contact

INTENDED USE/INDICATION FOR USE

Ascend AQ Ureteral Dilation Balloon Catheter Set (AUBS-) -used for transluminal dilation of ureteral strictures or ureteral dilation prior to ureteroscopy or stone manipulation. Intended for one-time use.

CONTRAINDICATIONS

1. Acute urethritis
2. Acute prostatitis
3. Acute cystitis
4. Urethral abnormalities (diverticula, doubling, etc.) because of the risk of perforation
5. Uncorrected coagulopathy
6. Untreated urinary infection

COMPLICATIONS /POTENTIAL ADVERSE EVENTS

1. Urinary tract infections
2. Pyelonephritis
3. Hemorrhage
4. Trauma to the surrounding tissue (due to excessive dilation of the balloon)
5. Bleeding
6. Dysuria
7. Fever
8. Postoperative pains

INSTRUCTIONS FOR USE

Ascend® AQ® Ureteral Dilation Balloon Catheter Set (AUBS)

BALLOON PREPARATION

1. Upon removal from the package, inspect the balloon catheter to ensure no damage has occurred during shipment.
2. Remove the protective sleeve from the balloon and discard. Lock a 10 ml syringe filled with 5 ml of dilute contrast medium onto the balloon lumen and apply negative pressure for 20 to 30 seconds. Release negative pressure on syringe allowing medium to draw into balloon. Detach the syringe leaving medium in the balloon lumen.

BALLOON CATHETER INTRODUCTION AND INFLATION

1. Under fluoroscopic control, cystoscopically pass a wire guide (up to .038 inch in diameter) the desired distance into the ureter beyond the area of planned dilation.
2. Activate the AQ hydrophilic coating by immersing the balloon and catheter in a sterile physiological solution for approximately 30 seconds.
3. Carefully advance the balloon catheter over the previously placed wire guide under fluoroscopic guidance utilizing the radiopaque markers to ensure proper positioning. Advance catheter so the entire balloon extends beyond the end of the cystoscope before inflating the balloon.
4. Prepare inflation device in standard manner and purge all air from syringe and tubing. To ensure proper regulation of balloon pressure, use of a balloon inflation device and pressure gauge is recommended.
5. Attach the prepared inflation device to the balloon lumen, open the stopcock on the inflation device and inflate the balloon.

6. Close the stopcock on the inflation device for at least 30 seconds to maintain pressure and allow adequate dilation. Time, rather than excessive pressure, is the key factor in transluminal dilation.

BALLOON DEFLATION AND WITHDRAWAL

1. Aspirate the balloon completely.
2. When the balloon is completely deflated, remove the device by withdrawing the catheter slowly. Removal of the catheter is facilitated by rotating the shaft counterclockwise during withdrawal. Using excessive force to withdraw the balloon can inflict trauma to tissue and/or damage the cystoscope.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred. If your device did not come with a Cook Sphere Inflation Device then the below instructions are not applicable, if your device did come with a Cook Sphere Inflation Device please refer to the instructions below for proper use.

INFLATION DEVICE

INSTRUCTIONS FOR USE

PREPARATION

Make all aspiration and injection maneuvers with the lock lever pushed left, i.e., unlocked. Unlock the piston by pushing the lock lever left. In this position, you can freely pull the piston back for aspiration, or push it forward for injection. To lock the piston in position, slide the lever right to the straight up position.

1. Prepare a solution of contrast medium and normal saline in a small sterile bowl. Check balloon catheter and contrast medium instructions for specific contrast mixture recommendations.
2. Orient the tubing downward into the contrast medium.
3. Push the release lever left and aspirate enough solution to fill the syringe. (Attach stopcock, if applicable.)
4. Hold the device upright to purge the air from the syringe and connecting tube. Tap the syringe lightly, if necessary, to remove all the air bubbles and fill the connecting tube completely.
5. Inspect the syringe and tubing (and stopcock, if applicable) to ensure that the device has been completely purged of air bubbles.
6. Adjust the syringe volume to the desired amount. If more contrast is needed, submerge the syringe tip into the basin of solution and aspirate. (Close stopcock, if applicable.)

ATTACHING THE INFLATION DEVICE TO THE BALLOON DILATION CATHETER

1. Prepare and test the balloon dilation catheter according to the manufacturer's directions for use.
2. If a separate syringe was used to prepare the balloon catheter, remove it. When a stopcock is installed on the end of the inflation device connecting tube, it should be opened and purged with contrast media from the inflation device to eliminate air. Create a fluid-fluid connection between the balloon and the stopcock or connecting tube (male rotating adapter) of the inflation device by placing a drop of contrast solution from the syringe into each hub.
3. Hand-tighten the hubs securely.

OPERATING INFLATION DEVICE

1. Release the lock lever and allow the piston to move forward into neutral position (0 atm).
2. To inflate the balloon, engage the lock lever, turn the palm grip on the piston clockwise slowly until the desired inflation pressure is reached. The lock lever maintains the increasing pressure.
3. To gradually deflate the balloon, turn the palm grip on the piston counterclockwise slowly until the desired deflation pressure is reached.
4. To rapidly deflate the balloon, push the lock lever left, releasing the piston, and pull back. Slide the lock lever back to lock, if desired.
5. After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable laws and regulations. For one time use only. Do not reuse. Do not re-sterilize. Read directions prior to use.

WARRANTY

An authorized representative of the manufacturer in the territory of the Russian Federation is the ZAO "SHAG".

Russia, 119002, Moscow, Karmanitskiy per., 9.office 501a.

phone: (495) 956-13-09;

Fax: (495) 956-1310,

e-mail: info@schag.ru

Manufacturer's warranty binding on quality products.

COMPLIANCE WITH RELEVANT REGULATIONS AND TECHNICAL STANDARDS

Cook Incorporated manufactures all medical products in accordance with Medical Devices Directive 93/42/EEC. Cook Incorporated manufacture all medical products in accordance with the United States Food and Drug Administration (FDA) Quality System Regulations, 21 CFR Part 820, *Quality System Regulations*.

All of Cook facilities are also certified to ISO 13485 (2003) *Quality Management System*.

SHELF LIFE

The "use by" date as denoted by symbols on the labeling. The shelf life is shown by indication the time of manufacture (factory symbol) to the stated expiration date (hourglass symbol). The difference between both dates is the shelf life.

The shelf-life of the products falling under the scope of this technical data summary - is 3 (three) years.

STERILIZATION DETAILS

These devices are supplied sterile in peel open packages utilizing an EtO (Ethylene Oxide) validated method to a sterility assurance level (SAL) of 10^{-6} . These devices are intended for one-time use, sterile if unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile.

STORAGE RULES

These devices should be stored in a dark, dry, cool place; avoid extended exposure to light.

RULES FOR USE

Urethral Dilation Catheters in Sets are intended for use in hospitals and health care facilities.

TRANSPORTATION RULES

Transportation shall account for the type of transportation employed, the destination, the vendor used for transportation and an estimation of the frequency and harshness of any handling.

PACKAGE

These devices are supplied with: balloon catheter, diameter 6.0 Fr, length 65 cm – 1 pc; one-way stopcock – 1 pc; inflation device – 1 pc; instructions for use – 1 pc.

TRANSPORTATION MARKINGS

Individual packages are packed in corrugated cardboard boxes. Boxes are sealed with a control tape around the side or around the perimeter. The transportation markings include information on the dimensions of the box, unique barcode, lot number, catalogue number and the global product number (GPN-Global Product Number), the number of single packages inside the box, expiration date, the Cook Medical logo, and the following symbols:

- symbol for "Fragile, handle with care"

- symbol for "Temperature Limitation"

- symbol for "Keep Dry"

PERFORMANCE CHARACTERISTICS

Ureteral Dilation Sets are stable under vibration loads in the range 10-55 Hz, amplitude displacement of 0.15 mm.

Ureteral Dilation Sets, in packaging are stable under vibration loads in the range 10-55 Hz, amplitude displacement of 0.35 mm and shock loads at peak acceleration of 100 m / s² (10g), and duration of shock acceleration of 16 ms.

Ureteral Dilation Sets are stable at temperatures in the range of +32 to + 42 ° C.

Ureteral Dilation Sets, in packaging are stable at a transportation and storage temperature in the in the range of -50 to + 50 ° C.

Bond tensile strength – no less than 10 N.

The ability to regulate the flow rate - when the one-way stop cock is "closed" position, the flow ceases completely, when opened, the flow is free, with no reduction in flow rate.

Smooth operation of the plunger (piston) of the inflation device – the movement of the plunger of the inflation device is smooth, doesn't require additional force along the entire length of the plunger.

Luer-Lock Male Adapter check - the adapter corresponds to the described size and is hermetically sealed when connected to the external tip and test article, the adapter corresponds to the Luer-Lock Male type.

The compatibility of the catheter with the wire guide – wire guide with a diameter of 0.97 mm (0.038 inch), passes through the catheter without getting stuck and does not require the application of additional force, catheter is compatible with the wire guide.

When the balloon catheter was inflated to a pressure equal to 20 atm. ($\pm 10\%$) no breaks or leakages occurred.

Ureteral Dilation Sets are resistant to corrosion.

Ureteral Dilation Sets are stable under sudden changes in temperature.

Ureteral Dilation Sets are supplied in sterile packages, do not attempt to clean, disinfect or re-sterilize.

PACKAGING PARAMETERS

RPNs	Height, cm	Width, cm	Weight, gr	Tolerance
AUBS-10-4	55.4	28.1	340.4	$\pm 10\%$
AUBS-4-10	55.5	28	341.8	$\pm 10\%$
AUBS-4-4	55.4	27.6	338.5	$\pm 10\%$
AUBS-4-6	55.4	28	340.6	$\pm 10\%$
AUBS-5-10	55.4	27.8	341.3	$\pm 10\%$
AUBS-5-4	55.1	27.9	339.2	$\pm 10\%$
AUBS-5-6	55	28	338.4	$\pm 10\%$
AUBS-6-10	55.1	28	341.4	$\pm 10\%$
AUBS-6-4	55.6	27.9	339.9	$\pm 10\%$
AUBS-6-6	55	28	339.6	$\pm 10\%$
AUBS-7-10	22.9	6.62	340	$\pm 10\%$
AUBS-7-4	55.2	28	339.8	$\pm 10\%$
AUBS-7-6	55.9	3.11	341	$\pm 10\%$

LABELING

LABELING

Urethral Dilation Catheters are supplied in the individual packaging, which should show no signs of damage or opening. Upon detecting the bends, breaks and other defects the device cannot be used. The device cannot be used if the integrity of the packaging has been compromised.

The label of the device shows the following device parameters:

- company-manufacturer and manufacturing address;
- device name (trade mark);
- marking about CE certification;
- marking about sterilization;
- address of the EU representative;
- trade mark of the manufacturer;
- manufacturing date;
- sizes;
- schematic picture;
- lot number;
- marking for "Single use only"
- marking for "Attention. See Instructions for Use"
- shelf life
- marking for "Keep Dry"
- marking for "Caution"
- marking for "Keep away from sunlight"
- catalogue number

Catheter Set

DISPOSAL RULES

Upon completion of procedure, the devices are to be disposed of as per institutional guidelines for biohazardous medical waste.

For EU countries -in accordance with the World Health Organization guidance document "Safe management of waste from Health Care Activities"

For US countries – in accordance with the US Department of Labor of Occupational Safety and Health Administration (OSHA)

For Russian Federation – in accordance with national requirements. Risk class of the medical waste: class B.

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