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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of WUXI JIAHAO GLOVES CO., LTD. on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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Authorized
Signature : XuePing

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(Date: JUN. 22, 2017)

OPERATIONAL DOCUMENTATION FOR RUSSIA

ARCHDALE Gloves PVC Examination (inspection) non-sterile powder-free plain polymer coated
ViniMAX
production


Wuxi Jiahao Gloves Co.,Ltd.

«Approve»

Wuxi Jiahao Gloves Co.,Ltd.

General director

Wang Xiaodi

Signature

Date

14.06.2017

Seal

	Table of contents	Page
	Product data	3
	Name and designation of the medical product	3
	Manufacturer of the medical product	3
	Description and operation	3
1	Brief information on the medical device including the intended use	3
2.2	Technical specifications and characteristics warranted by the manufacturer	3
2.3	Composition (delivery completeness) of the medical product	5
2.4	Warnings, precautions and limitations on the medical product use	5
2.5	Description of the packaging, labeling and warning signs on the product	5
3	Medical product application	6
3.1	Operating procedure (product intended use)	6
4	Maintenance and current repair	6
5	Disposal	6
6	Warranties	6
7	Indications and contraindications, side effects.	7

1. Product information

1.1. Name and designation MI

ARCHDALE vinyl diagnostic gloves (viewing) non-sterile powder-free powdered with a polymer coating ViniMAX: dimensions XS, S, M, L, XL

1.2. Manufacturer of medical device

Information about the manufacturer:

The name Wuxi Jiahao Gloves Co., Ltd.

Address: Xielidang Village, Houqiao, Xishan District, Wuxi City, Jiangsu Province, People's Republic of China

Phone: + 86-510-88721060

Data on the authorized representative of the manufacturer in the Russian Federation:

Limited Liability Company "Ardale-Impex" (LLC Ardale-Impex)

115201, Moscow, Kashirsky proezd, house 23, building 5, office.3

Tel.: 8 495 777 1170, 8 495 777 1171, e-mail: abx@ardl.ru, INN 7724781425

2. Description and work

2.1. Brief information about the medical device, including the purpose

The product is a non-sterile gloves made of polyvinyl resin as a main material and other ingredients intended to create a protective barrier for the medical worker's hands during examination of the patient to protect the patient and researchers from cross-contamination.

Smooth working surface of gloves

Glove classification:

- 1) depending on the material - made of polyvinyl resin as the main material and additional ingredients;
- 2) depending on the finish, smooth over the entire outer surface;
- 3) depending on the finish - non-powdered;
- 4) depending on the edge of the cuff - with a whisk on the edge;
- 5) depending on the color - natural (without dye);
- 6) depending on the dimensions - XS, S, M, L, XL.;
- 7) depending on the form - non-anatomical;

Gloves are universal (for both hands).

They have a polyurethane coating on the inside for easy donning.


2.2. Manufacturer's guaranteed technical parameters and characteristics

Materials used for manufacturing


Ingredient name	chemical name	Function	%
Polyvinyl resin	Polyvinyl brand SG-3	Raw material	58
DINEF	Diisononyl phthalate Brand Jayflex DINP Plasticizer	Plasticizer	41
Polyurethane coating	Polyurethane brand NPUA-1D	On the inside for easy donning	1

Main technical characteristics

The points for measuring the width and length of gloves are shown in the figure below.

The points for measuring the thickness of the glove wall are shown in the figure below.

The thickness of the glove wall at the cuff level is measured at a distance of 25 ± 5 from the edge of the cuff

Physical characteristics of gloves

Pipe size	Width	* Length, mm, not less	*** The thickness at the points shown in the figure, at least, mm	** The thickness in the center of the palm, not more, mm
XS	≤ 80 , but not less than 70	240	0.08	0.22
S	80 ± 10	240		
M	95 ± 10	240		
L	110 ± 10	240		
XL	≥ 110 but not more than 120	240		

- * The maximum permissible length (regardless of size) is not more than 270 mm
 - ** the minimum thickness in the center of the palm is not less than 0.08 mm
 - *** the maximum permissible thickness (in the points indicated in the figures) is not more than 0.20 mm
- Weight of one glove is not more than 5 g (for all sizes)

Physical properties of gloves

Before accelerated aging		After accelerated aging	
Effort at break, N, not less than	Elongation at break, %, not less than	Effort at break, N, not less than	Elongation at break, %, not less than
7.0	350	7.0	350

Gloves are hermetic, the level of sealing quality (AQI) is not more than 1.5

Gloves with a whisk of circular cross-section on the edge of the cuff. the diameter of the corolla is not more than 2.5 mm.

Gloves, not powdered, residual amount of dusting substance < 2.0 mg per glove

2.3. Composition (completeness of delivery) of a medical device

One group package (block) contains 100 pieces. Gloves. Transport packaging - 10 blocks.

2.4. Warnings, precautions and limitations when using MI

Precautionary measures

- 1) The hand skin should be well groomed, without wounds, microcracks and other damages.
- 2) The nails should be short and neatly trimmed. Do not use artificial nails.
- 3) Jewellery and wrist watches during operation of gloves should be removed.
- 4) Gloves are intended solely to provide a protective barrier for hands health worker during the examination of the patient in order to protect the patient and the investigator from cross-contamination.
- 5) Treatment with alcohol-containing solutions is prohibited, t. These solutions destroy the protective layer of gloves.
- 6) Gloves contain a plasticizer (diisononyl phthalate) which can be unsafe for the user

2.5. Description of packaging, labeling and warning labels on the product

Packaging

A consumer cardboard box, in one box there are 100 (50 pairs) pieces of gloves. The box allows you to maintain the quality of gloves during transportation and storage.

The weight of the package with gloves (100 pcs.) Is not more than 700 g (irrespective of the size)

Marking

Marking is reliable, verifiable and readable. Marking is applied to consumer packaging in any way (printing, embossing, stamp) ensuring its clarity, clarity and readability. When using the printed method of marking, paint can not be scraped off.

The marking on the consumer package of gloves contains:

- name and location of the manufacturer;
- name and location of the representative in the country (s) of the product circulation;
- a record of the material used or an equivalent and / or figure that clearly reflects this information;
- inscriptions describing the abbreviated and full name of the product: "VINYL GLASS GLOVES" and "ARCHDALE GLOVES VINYL DIAGNOSTIC (LIQUID) NONSTERILE UNSAFFECTED SMOOTHES WITH POLYMERIC COVERAGE VINMAX";
- an inscription: "UNSTERY" and / or a drawing that clearly reflects this indication;
- information: "ONE-TIME" and / or drawing, clearly reflecting this indication;
- information: "SMOOTH" and / or drawing, clearly reflecting this indication;
- the inscription: "UNSPECIFIED" and / or a drawing that clearly reflects this indication;
- information: "WITH POLYMERIC COATING" and / or a drawing that clearly reflects this indication;
- information: "The product contains a plasticizer (diisononyl phthalate), which can be unsafe for the user";
- lot number;
- data (month, year) of manufacture;
- shelf life;
- Size and number of pairs of gloves;
- information about storage and transportation conditions.

It is allowed to supplement the marking with other information, for example, to put symbols and drawings explaining the consumer properties of gloves and their application.

Do not apply information about the special properties of gloves without the appropriate documentary confirmation.

The marking shall be applied in the language of the country of circulation of the product, excluding the proper names and names.

3. Exploitation of MI

3.1. Procedure (use the product for its intended purpose)

Open the group pack of gloves, remove one pair of gloves and successively put the gloves on the right, then on the left hand.

Consistently on the left and right gloves turn the bent edge so that the glove completely covers the cuff of the dressing gown.

Removal of gloves

Take the glove in the palm base with the fingers of the other hand.

Pull off the glove and hold it with the other hand.

Place the fingers of the free hand in the cuff of the glove on the other hand and pull it together like the first.

Turning the glove inside out, use it as a cover for another glove.

Place the gloves in a suitable container.

After removing the gloves, wash your hands with soap and treat them with an alcohol-containing cutaneous antiseptic.

It is necessary to change gloves:

- every two hours of continuous operation;
- after detecting a visible defect;
- after a possible or actual puncture;
- if the glove begins to swell, stretch, stick, or loosely sits on the arm;
- in accordance with the rules adopted at the place of application of these gloves.

Do not use the same pair of gloves when changing from one patient to another, or from one manipulation to another.

4. Maintenance and repair

Gloves belong to single-use products and are not subject to maintenance and repair.

Gloves operated at a temperature (+10 to +35) °C and humidity (80-100)%.

Stored in original package in a cool, dry and well ventilated space and eliminating dust, exposure to sunlight, moisture, X-ray radiation at a temperature of (-50 to +40) °C and humidity (60-75)%.

The product is transported by all types of covered vehicles in accordance with the traffic regulations in force for this transport at a temperature of (-50 to +50) °C and humidity (60-75)%.

5. Disposal

Disposal and destruction of used products in Russia is carried out in accordance with the requirements of SanPiN 2.1.7.2790-10, the used gloves are classified as class B waste (epidemiologically hazardous waste). Collection of used gloves should be carried out in sealed packages of yellow color with the marking (inscription).

6. Warranties

Gloves belong to medical products of a single use.

The shelf life of gloves is 5 years from the date of manufacture.

Guaranteed shelf life of gloves is 5 years.

7. Indications and contraindications to use, side effects.

Indications for use: to create a protective barrier for hands health worker during the examination of the patient in order to protect the patient and the investigator from cross-contamination.

Contraindications: Gloves contain a plasticizer (diisononyl phthalate), which can be unsafe for the user.

Side effects: contact dermatitis (in individuals due to the presence of a plasticizer (diisononyl phthalate)).

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