MERIT MEDICAL SYSTEMS, INC
1600 WEST MERIT PARKWAY
SOUTH JORDAN UTAH 84095
PHONE 801-253-1600
FAX 801-253-1688
www.merit.com

ONE Snare & ONE Snare Microsnare Instructions for Use w/ Annex
Dated July 18, 2017

ACKNOWLEDGEMENT


State of Utah)

§

County of Salt Lake)

On this 18th day of July, 2017, before me, Jennifer Fordham, a notary public, personally appeared Dan Clark, Director, Regulatory Affairs, proved on the basis of satisfactory evidence to be the person whose name is subscribed to the attached instrument, and acknowledge he executed the same.

Witness my hand and official seal.


Notary Public

My Commission Expires: June 27, 2021

Handwritten signature and date: 5013, 18, 2017

ONE Snare

ONE Snare

STANDARD SNARE KIT

INSTRUCTIONS FOR USE:

DESCRIPTION:

English

The ONE Snare endovascular snare system consists of the snare, snare catheter, insertion tool and torque. The snare is constructed of nitinol cable and gold plated tungsten loop. The pre-formed snare loop can be introduced through catheters without risk of snare deformation because of the snare's super-elastic construction. The snare catheter is constructed of polyether block amide and contains a platinum/iridium radiopaque marker band.

INDICATION FOR USE:

The ONE Snare endovascular snare system is intended for use in the coronary and peripheral vascular system or hollow viscous to retrieve and manipulate foreign objects. Retrieval and manipulation procedures include indwelling venous catheter repositioning, indwelling venous catheter fibrin sheath stripping, and central venous access venipuncture procedure assistance.

CONTRAINDICATIONS:

1. This device is not intended for the removal of foreign objects entrapped by tissue growth.
2. This device should not be used for fibrin sheath stripping in the presence of septal defects of Persistent Foramen Ovale.
3. This device is not intended for removal of implanted pacing leads.

WARNINGS:

1. Pull forces applied to catheters during fibrin sheath stripping may damage, stretch, or break indwelling catheters 6 French or smaller in diameter. Do not use excessive pull force when attempting fibrin sheath stripping of catheters 6 French or smaller in diameter.
2. Do not use excessive force when manipulating the catheter through an introducer. Excessive force may damage the snare catheter.
3. This device has been sterilized utilizing ethylene oxide and is considered sterile if the package is not opened or damaged. Do not use a device that has been damaged or if the package is open or damaged.
4. For single patient use only. Do not reuse, reprocess or sterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
5. After use this device may be a potential biohazard. Handle in a manner that will prevent accidental contamination.
6. Nitinol is a nickel titanium alloy. Possible reaction may occur for those patients who exhibit sensitivity to nickel.
7. Care should be observed when using this device for removal of a large fibrin sheath in order to minimize risk of pulmonary embolism.

Caution: Only: CAUTION - Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

POTENTIAL COMPLICATIONS:

1. Potential complications associated with foreign body retrieval devices in arterial vasculature include, but are not limited to:
 - Embolization
 - Stroke
 - Myocardial Infarction (depending upon placement)
2. Potential complications associated with snare retrieval devices in venous vasculature include, but are not limited to:
 - Pulmonary embolism
3. Other potential complications associated with foreign body retrieval devices include, but are not limited to:
 - Vessel perforation
 - Device entrapment

Catheter damage can occur when attempting fibrin sheath stripping on small French size diameter catheters. (See WARNINGS)

Incidence of pulmonary embolism after fibrin sheath stripping may occur. (See WARNINGS)

DIRECTIONS FOR USE:

Prepare the ONE Snare System:

Select the appropriate snare diameter range for the site in which the foreign body is located. The snare diameter range should approximate the size of the vessel in which it will be used.

1. Remove the snare and snare catheter from their hoop holders and inspect for any damage.
2. Remove the insertion tool and torque device from the proximal end of the snare shaft.
3. Load the snare into the snare catheter by inserting the proximal end of the snare into the distal (non-hubbed) end of the snare catheter, until the proximal end of the snare shaft exits the hub and the loop can be retracted into the distal end of the snare catheter.
4. Inspect the device by extending and retracting the snare loop through the distal end of the snare catheter 2-3 times, while carefully examining the snare catheter, radiopaque band and the device for any damage or defects.
5. When appropriate, the system (snare and snare catheter) can be advanced to the desired site as a single unit assembled as described above.

Alternative Preparation of the ONE Snare System:

If the snare catheter is already positioned within the vasculature, the provided insertion tool (located on the proximal end of the snare and just distal to the torque device) may be used to position the snare in the indwelling snare catheter.

1. Remove the snare from the protective holder and inspect for any damage.
2. Move the provided insertion tool (located on the proximal end of the snare, just distal to the torque device) distally until the loop of the snare is enclosed within the tubing portion of the insertion tool.
3. Insert the distal end of the insertion tool into the hub of the indwelling snare catheter until resistance is felt. This will indicate the tip of the insertion tool is properly aligned with the inner lumen.
4. Hold the insertion tool as straight as possible, grasp the shaft of the snare just proximal to the hub of the insertion tool and advance the snare until it is secure within the lumen of the snare catheter. The insertion tool can be removed by grasping the blue tab and firmly peeling it away from the snare shaft.

Snare Assisted Retrieval and Manipulation Suggestions:

1. If present, remove the indwelling delivery catheter.
2. If a guidewire is in a patient at the location of a foreign body, advance a snare catheter over the guidewire to the desired location. Then remove the guidewire and advance the snare through the snare catheter. Alternatively, cinch one loop of the snare over the proximal end of the guidewire and advance the entire system (snare and snare catheter assembly) into a guide catheter or introducer sheath until the distal end of the snare catheter is positioned proximal to the foreign body.
3. If a guidewire is not present, pull the snare into the distal end of the snare catheter and advance through a guide catheter or introducer sheath until it is positioned proximal to the foreign body. Alternatively, collapse the snare loop by pulling the device into the distal end of the snare insertion tool. Place the tapered end of the snare insertion tool into the proximal (hub) end of the snare catheter, guide catheter or sheath and advance the snare forward maintaining constant contact between the insertion tool and snare catheter hub. NOTE: When attempting to utilize guide catheters or sheaths not specifically manufactured for use with the ONE Snare system, it is important to test product compatibility prior to use.
4. Gently push the snare shaft forward to completely open the loop. The loop is then slowly advanced forward, and around the proximal end of the foreign body. Alternatively, the snare may be advanced beyond the target location and the loops brought back around the distal end of the foreign body.
5. By advancing the snare catheter, the loop of the snare is closed to capture the foreign body. Note that attempting to close the loop by pulling the snare into the snare catheter will move the loop from its position around the foreign body.
6. To manipulate a foreign body, maintain tension on the snare catheter to retain the hold on the foreign body, and move the snare and snare catheter together to manipulate a foreign body to the desired position.
7. To retrieve a foreign body, maintain tension on the snare catheter and move the snare and snare catheter assembly together proximally to or into a guide catheter or sheath. The foreign body is then withdrawn through or together with the guiding catheter or introducer sheath. Withdrawal of large foreign bodies may require the insertion of larger sheaths, guiding catheters, or a cut-down at the peripheral site.

Snare Assisted Removal of Fibrin Sheaths from Indwelling Catheters:

1. Using standard technique, prepare a femoral vein approach, advance the selected snare to the inferior vena cava or right atrium.
2. Advance a .035" (0.89 mm) guidewire through the end port (distal or venous port if more than one lumen) of the indwelling catheter and into the inferior vena cava or right atrium.
3. Position the snare loop around the guidewire.
4. Advance the snare over the distal end of the catheter to a position proximal to the fibrin sheath.
5. Close the snare around the catheter and continue applying light traction while gently pulling the snare down toward the distal end of the catheter over the end ports.
6. Repeat steps 4 & 5 until the catheter is free of fibrin sheath.

Snare Assisted Venous Canalization:

1. Introduce the snare at a patent venous access site and position in the vasculature at the desired site.
2. Open the snare loop to provide a target to guide an entry needle into the desired venous access site.
3. Introduce a guidewire through the needle and through the snare loop. Remove the needle.
4. Close the snare over the guidewire by advancing the snare catheter.
5. Pull the guidewire into the inferior vena cava.

ONE Snare System contains:

(1) Snare, (1) Snare Catheter, (1) Insertion Tool and (1) Torque Device.

MERITMEDICAL

CE 0086

Merit Medical Systems, Inc.
1600 West Merit Parkway, South Jordan, Utah 84095 U.S.A. 1-801-253-1600
U.S.A. Customer Service 1-800-356-3748

EC REP

Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland

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EC REP

Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ирландия

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INSTRUCTIONS FOR USE:

DESCRIPTION:

The ONE Snare[®] Endovascular Microsnare System consists of the snare, snare catheter, insertion tool and torque device. The microsnare is constructed of nitinol cable and a gold plated tungsten loop. The pre-formed microsnare loop can be introduced through catheters without risk of snare deformation because of the snare's super-elastic construction. The microsnare catheter is constructed of polyether block amide (Pebax[®]) and contains a platinum/iridium radiopaque marker band.

INDICATIONS FOR USE:

The ONE Snare[®] endovascular microsnare system is intended for use in adult and pediatric populations for the retrieval and manipulation of traumatic foreign bodies located in the coronary and peripheral cardiovascular system and the extra-cranial anatomy.

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WARNINGS:

1. Do not use excessive force when manipulating the catheter through an introducer. Excessive force may damage the catheter.
2. This device has been sterilized utilizing ethylene oxide and is considered sterile if the package is not opened or damaged. Do not use a device that has been damaged or if the package is open or damaged.
3. For single patient use only. Do not reuse, reprocess or sterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
4. After use this device may be a potential biohazard. Handle in a manner that will prevent accidental contamination.
5. Nickel is a nickel-titanium alloy. Possible reaction may occur for those patients who exhibit sensitivity to nickel.

Only: CAUTION - Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

POTENTIAL COMPLICATIONS:

1. Potential complications associated with snare retrieval devices in arterial vasculature include, but are not limited to:
 - Embolization
 - Stroke
 - Myocardial infarction (depending upon placement)
2. Potential complications associated with snare retrieval devices in venous vasculature include, but are not limited to:
 - Pulmonary embolism
3. Other potential complications associated with foreign body retrieval devices include, but are not limited to:
 - Vessel perforation
 - Device entrapment

DIRECTIONS FOR USE:

Prepare the ONE Snare[®] Endovascular Microsnare System:

Select the appropriate microsnare diameter range for the site in which the foreign body is located. The microsnare diameter range should approximate the size of the vessel in which it will be used.

1. Remove the microsnare and microcatheter from their hoop holders and inspect for any damage.
2. Remove the insertion tool and torque device from the proximal end of the microsnare shaft.
3. Load the microsnare into the microcatheter by inserting the proximal end of the microsnare into the distal (non-hubbed) end of the microsnare catheter, until the proximal end of the microsnare shaft exits the hub and the loop can be retracted into the distal end of the snare microcatheter.

Alternative Preparation of the ONE Snare[®] Endovascular Microsnare System:

If the snare microcatheter is already positioned within the vasculature, the provided insertion tool (located on the proximal end of the microsnare and just distal to the torque device) may be used to position the microsnare in the indwelling snare microsnare catheter.

1. Remove the microsnare from the protective holder and inspect for any damage.
2. Move the provided insertion tool distally until the loop of the microsnare is enclosed within the tubing portion of the insertion tool.
3. Insert the distal end of the insertion tool into the hub of the indwelling snare microcatheter until resistance is felt. This will indicate the tip of the insertion tool is properly aligned with the inner lumen.
4. Hold the insertion tool as straight as possible, grasp the shaft of the microsnare just proximal to the hub of the insertion tool and advance the microsnare until it is secure within the lumen of the microcatheter. The insertion tool can be removed by grasping the blue tab and firmly peeling it away from the microsnare shaft.

Position Microsnare Proximal to Foreign Body (Microsnare Assembled Outside the Body/ Guidewire Present):

1. Remove the indwelling balloon or delivery catheter, if present. It may be necessary to change or extend the indwelling guidewire to facilitate balloon removal, and to upsize the indwelling guide catheter/sheath to accommodate the 2.3-3.0F (0.81-1.0 mm) microsnare catheter.
2. Place the microsnare loop over the proximal end of the wire, and cinch until just snug by pulling the proximal end of the microsnare shaft.
3. Advance the microsnare, distal end first, into the guide catheter or sheath until the distal end of the catheter is positioned just proximal to the foreign body.

Microsnare Assembled Outside The Body/Guidewire Not Present:

1. Pull the microsnare loop just into the distal end of the microcatheter.
2. Advance the microsnare in the guide catheter until it is positioned just proximal to the foreign body.

Microsnare Assembled With An Indwelling Catheter:

1. Advance the microsnare shaft in the microcatheter until the loop is just proximal to the microcatheter end hole.
2. Advance the microsnare in the guide catheter until it is positioned just proximal to the foreign body.

Microsnare Assisted Manipulation/Retrieval:

1. Push the microsnare shaft forward gently to completely open the loop. Slowly advance the loop forward, and around the proximal end of the foreign body.
2. Once the loop is around the foreign body, advance the catheter to cinch the loop closed and grab the foreign body. (Note that attempting to close the loop by pulling the microsnare within the catheter will move the loop from its position around the foreign body.)
3. To manipulate a foreign body, maintain tension on the catheter to retain the hold on the foreign body, and move the microsnare and microcatheter together to manipulate foreign body to the desired position.
4. To retrieve a foreign body, maintain tension on the microcatheter and move the microsnare and catheter together (proximally) to, or into the guide catheter or sheath. The foreign body is then withdrawn through or together with the guiding catheter or vascular sheath. Withdrawal of large foreign bodies may require the insertion of larger sheaths, or guiding catheters, or a cut down at the peripheral site.

ONE Snare[®] Endovascular Microsnare System contains:

- (1) Microsnare, (1) Microsnare Catheter, (1) Insertion Tool and (1) Torque Device.

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ANNEX TO INSTRUCTIONS FOR USE

MEDICAL DEVICE NAME:

- ONE Snare Endovascular Snare System, configuration options: ONE500, ONE1000, ONE1500, ONE2000, ONE2500, ONE3000, ONE3500, ONE1001, ONE2501
- ONE Snare Endovascular Microsnare System, configuration options: ONE200, ONE201, ONE400, ONE401, ONE700, ONE701.

ATTENTION – the device should be used only in specialized medical institutions.

MEDICAL DEVICE COMPOSITION:

1. Snare device with an insertion tool (introducer) and a torque device
2. Snare catheter
3. Instructions for use

INTENDED USE:

The device is intended for manipulations in the lumen of coronary and peripheral vessels or in the lumen of hollow viscus, as well as for retrieval of foreign bodies.

Endovascular snare system is an auxiliary device that is temporarily introduced into the lumen of vessels or hollow viscus at certain stages of intravascular intervention.

PRECAUTIONS:

Inspect the device for any damage prior to use.

Do not use if the package is damaged.

Do not reuse, reprocess or resterilize.

After use the device may pose threat of biological contamination. Handle the device in such a way as to prevent incidental contamination.

An allergic reaction may develop in patients sensitive to nickel. Do not use the device in patients allergic to device materials.


Prior to operation carefully read the instructions for use.

PRINCIPLE OF OPERATION:

ONE Snare Endovascular Snare, ONE Snare Endovascular Microsnare and EN Snare Endovascular Snare Systems are used by clinicians to manually capture, manipulate or retrieve objects in patient's cardiovascular system or hollow viscus, which is typically performed under fluoroscopy. The snare is loaded by inserting the proximal end of the snare into the distal end of the snare catheter, until the proximal end of the snare shaft exits the hub. After that the assembled system is advanced to target site. Slowly advance the loop forward and around the proximal end of the foreign body. Once the loop is around the foreign body, advance the catheter to cinch the loop closed and grab the foreign body. To manipulate a foreign body, maintain tension on the catheter to retain the hold on the foreign body. To retrieve a foreign body move the snare and catheter together proximally or inside the guiding catheter or sheath.

Prior to use the appropriate snare diameter range for the foreign body is located is selected. The snare diameter range should approximate the diameter of the vessel in which it will be used.

Inspect the device for any damage prior to use.

[Signature]
DAN CLARK
JULY 18, 2017

MEDICAL DEVICE SPECIFICATIONS:

Catalog numbers	Loop diameter (mm)	Snare length (mm)	Catheter OD (in)/(mm)	Catheter ID (in)/(mm)	Catheter length (mm)	Insertion tool OD (in)/(mm)	Insertion tool ID (in)/(mm)	Insertion tool length (mm)	Torque device diameter (in)/(mm)	Torque device length (mm)
ONE Snare Endovascular Snare System										
ONE500	5	1200	0.057/1.45	0.040/1.02	1000	0.081/2.06	0.040/1.02	46	(0,018 – 0,038) / (0.46 – 0.97)	40
ONE1000	10	1200	0.057/1.45	0.040/1.02	1000	0.081/2.06	0.040/1.02	46	(0,018 – 0,038) / (0.46 – 0.97)	40
ONE1500	15	1200	0.082/2.11	0.056/1.42	1000	0.081/2.06	0.054/1.37	71	(0,018 – 0,038) / (0.46 – 0.97)	40
ONE2000	20	1200	0.082/2.11	0.056/1.42	1000	0.081/2.06	0.054/1.37	71	(0,018 – 0,038) / (0.46 – 0.97)	40
ONE2500	25	1200	0.082/2.11	0.056/1.42	1000	0.081/2.06	0.054/1.37	71	(0,018 – 0,038) / (0.46 – 0.97)	40
ONE3000	30	1200	0.082/2.11	0.056/1.42	1000	0.081/2.06	0.054/1.37	71	(0,018 – 0,038) / (0.46 – 0.97)	40
ONE3500	35	1200	0.082/2.11	0.056/1.42	1000	0.081/2.06	0.054/1.37	71	(0,018 – 0,038) / (0.46 – 0.97)	40
ONE1001	10	650	0.057/1.45	0.040/1.02	480	0.081/2.06	0.041/1.04	46	(0,018 – 0,038) / (0.46 – 0.97)	40
ONE2501	25	650	0.082/2.11	0.056/1.42	480	0.081/2.06	0.054/1.37	71	(0,018 – 0,038) / (0.46 – 0.97)	40
ONE Snare Endovascular Microsnare System										
ONE200	2	1750	(0.032-0.041) / (0.81-1.00)	0.019/0.48	1500	0.081/2.06	0.020/0.51	104	(0,012 – 0,025) / (0.31 – 0.64)	40
ONE201	2	2000	(0.032-0.041) / (0.81-1.00)	0.019/0.48	1750	0.081/2.06	0.020/0.51	104	(0,012 – 0,025) / (0.31 – 0.64)	40
ONE400	4	1750	(0.032-0.041) / (0.81-1.00)	0.019/0.48	1500	0.081/2.06	0.020/0.51	104	(0,012 – 0,025) / (0.31 – 0.64)	40
ONE401	4	2000	(0.032-0.041) / (0.81-1.00)	0.019/0.48	1750	0.081/2.06	0.020/0.51	104	(0,012 – 0,025) / (0.31 – 0.64)	40
ONE700	7	1750	(0.032-0.041) / (0.81-1.00)	0.019/0.48	1500	0.081/2.06	0.020/0.51	104	(0,012 – 0,025) / (0.31 – 0.64)	40
ONE701	7	2000	(0.032-0.041) / (0.81-1.00)	0.019/0.48	1750	0.081/2.06	0.020/0.51	104	(0,012 – 0,025) / (0.31 – 0.64)	40

DEVICE MATERIALS:

The materials used in the manufacture of the components of ONE Snare Endovascular Snare System and ONE Snare Endovascular Microsnare System are identified below:

Component Material Specifications		
Snare Materials		Type and duration of the contact
Component	Material	
Guidewire (core)	Black oxide, nitinol SE508	Externally communicating device, contact with circulating blood, limited contact duration (≤ 24 hours)
Standard snare cable	Nitinol	
Crimp stem (PTFE shrink tubing)	PTFE, TFX	
Snare loop	Gold-plated tungsten wire / nitinol	
Bonding material (adhesive)	Loctite 4304	
Insertion tool (introducer) materials		No patient contact
Tubing	Polypropylene natural 13R9A	
Sleeve	Polypropylene natural P4G2K-152A	
Torque device materials		No patient contact
Cap	Polycarbonate, Dow PDX08115fc-03000460	
Body	Polypropylene natural 13R9A	
Collet sleeve	Brass (high cutting capacity D36000)	
Catheter materials		Externally communicating device, contact with circulating blood, limited contact duration (≤ 24 hours)
Tubing	Pebax 7233 SA01	
Marker band	Platinum – 90% Iridium – 10%	
Luer hub	Polycarbonate – injection molding resin, clear, Dow PDX08115fc-03000460	
Catheter sleeve (cannula)	Pebax 3533 SA01MED	

SAFETY REQUIREMENTS:

ONE Snare Endovascular Snare System and ONE Snare Endovascular Microsnare System are safe and effective for their intended use as labeled.

OPERATION, STORAGE AND TRANSPORTATION CONDITIONS:

Operation conditions: in medical institutions at temperature from +10°C to +35°C and relative humidity up to 85% at 25°C.

Transportation conditions: at temperature from -50°C to +50°C and relative humidity up to 100%. After long-term exposure to subzero temperatures the device should be conditioned in its packaging for several days at storage temperature.

Storage conditions: in a cool dry place at temperature from +5°C to +40°C and relative humidity up to 85% at +25°C.


PACKAGING:

Primary packaging (inner and outer):

The snare device, insertion tool and torque device (insertion tool is threaded onto the proximal end of the snare shaft to minimize damage during shipping/storage) are placed in polypropylene spiral


hoop (dispenser / protective retaining ring) and then inside a separate Tyvek/Nylon pouch that is sealed using a validated heat seal process. A label containing information on the device **is applied**. The catheter is placed in polypropylene spiral hoop (dispenser / protective retaining ring) and then inside a separate Tyvek/Nylon pouch that is sealed using a validated heat seal process. A label containing information on the device is applied.

Snare device / torque device / insertion tool and catheter placed in a spiral hoop (dispenser / protective retaining ring):

Retail packaging:

Both pouches and instructions for use are placed in a carton unit box.

Shipping package: Then twenty (20) cartons are placed into a shipper box.

LABELING:

The endovascular snare system is labeled as required by Annex I, Council **Medical Devices Directive 93/42/EEC**, and graphical symbols used in the labeling comply with **EN 980:2003 requirements**.


Inner packaging (pouch) labeling layout for Russia contains the following information:

- manufacturer's name
- device trade name
- device dimensions (loop diameter and snare length)
- graphical scheme of pouch contents
- bar code
- CE marking
- catalog number (REF)
- lot number
- symbol "Best before", shelf life

- symbol "Sterile, sterilized with ethylene oxide"
- symbol "Do not use if package is damaged"
- symbol "Single use only"
- symbol "Non-pyrogenic"
- symbol "Consult accompanying documents"
- manufactured by... name and full address of the manufacturer
- EC representative address
- manufacturing date

Carton unit box labeling contains the following information:

- manufacturer's name
- device trade name
- device dimensions (loop diameter and snare catheter length)
- graphical scheme of pouch contents
- bar code
- CE marking
- catalog number (REF)
- lot number
- symbol "Best before", shelf life
- symbol "Sterile, sterilized with ethylene oxide"
- symbol "Do not use if package is damaged"

- symbol "Single use only"
- symbol "Non-pyrogenic"
- symbol "Consult accompanying documents"
- manufactured by... name and full address of the manufacturer
- EC representative address
- manufacturing date

Table of symbols used in packaging:

Symbol	Meaning
	This symbol means that the device meets essential requirements of EC Directive 93/42/EEC concerning medical devices
REF	This symbol means device number according to manufacturer's catalog. The number itself is stated to the right of the symbol
	This symbol means device lot number. A letter and number code is stated to the right of the symbol
	This symbol means device shelf life. Shelf life (as a date) is stated to the right of the symbol
	This symbol means that the device is sterile
	This symbol means that the device has been sterilized using ethylene oxide
	This symbol prohibits to use the device if the package is not intact
	This symbol warns to operate the device without reading instructions for use
	This symbol indicates that the medical device is non-pyrogenic
	This symbol prohibits to reuse the device
	Manufacturer's address is stated next to this symbol
	EC representative's address is stated next to this symbol

METHODS AND MEANS OF DISINFECTION AND PRE-STERILIZATION CLEANING:
Not applicable. ONE Snare Endovascular Snare System and ONE Snare Endovascular Microsnare System are sterile devices for single use only.
Used devices are subject to disinfection prior to disposal.

STERILITY:

ONE Snare Endovascular Snare System and ONE Snare Endovascular Microsnare System are sterilized using a traditional ethylene oxide sterilization process. Prior to sterilization the device is assembled and packaged into dispensers and pouches (inner packaging). After that the system is sterilized. The device is not intended to be sterilized by the user.

DISPOSAL AND DESTRUCTION PROCEDURE

Used devices are gathered in a waste container and disposed in accordance with applicable regulations and local rules for disposal.

Devices are subject to disinfection prior to disposal. For single-use devices decontamination apply: 6% peroxide, 3% chloramine solution, 1.5% calcium hypochlorite solution, 2% Lysoformin 3000 solution. Exposure to disinfectants – 60 minutes or in 0.05% anolyte – 2 hours.

MEDICAL DEVICE CLASSIFICATION:

Potential risk of medical device use classified in accordance with nomenclature of medical devices – 3.

Complies with essential requirements of EC Medical Devices Directive 93/42/EEC.
Products of Merit Company meet ISO 13485:2003 requirements.

WARRANTY:

The manufacturer guarantees medical device sterility, as well as its parameters and performance within the stated shelf life, i.e. three years after manufacturing date, provided individual packaging remains intact and storage, transportation and operation conditions are observed according to the instructions for use.

MEDICAL DEVICE MAINTENANCE AND REPAIR. AFTER-SALES SERVICE:

ONE Snare Endovascular Snare System and ONE Snare Endovascular Microsnare System contain no parts requiring maintenance during operation. Thus, field repair is not performed.

SHELF LIFE:

3 years.

RECLAMATIONS:

Reclamations should be sent to the manufacturer's authorized representative in Russia:
Merit Technologies, LLC
115054, Moscow, Bakhrushina St. 32, bldg. 1
Tel: + 7 495 221 89 02

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EN Snare Instructions for Use w/ Annex
Dated July 18, 2017

MERIT MEDICAL SYSTEMS, INC
1600 WEST MERIT PARKWAY
SOUTH JORDAN UTAH 84095
PHONE 801-253-1600
FAX 801-253-1688
www.merit.com

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
INSTRUCTIONS FOR USE:

INDICATION FOR USE:

The EN Snare Endovascular Snare and Catheter is intended for use in the cardiovascular system or hollow viscous to retrieve and manipulate foreign objects. Manipulation procedures include indwelling venous catheter repositioning, indwelling venous catheter fibrin sheath stripping, and central venous access venipuncture procedure assistance.

DESCRIPTION:

The EN Snare system consists of three interlaced, cabled, super-elastic Nitinol, preformed loops. The super-elastic Nitinol construction enables the loops to be introduced through catheters without the risk of device deformation.

WARNINGS:

1. This device is not intended for the removal of foreign objects entrapped by tissue growth.
2. This device should not be used for fibrin sheath stripping in the presence of septal defects of Persistent Foramen Ovale.
3. This device is not intended for removal of implanted pacing leads.
4. Pull forces applied to catheters during fibrin sheath stripping may damage, stretch, or break indwelling catheters 6 French or smaller in diameter. Do not use excessive pull force when attempting fibrin sheath stripping of catheters 6 French or smaller in diameter.
5. Do not use excessive force when manipulating the catheter through an introducer. Excessive force may damage the catheter.
6. This device has been sterilized utilizing Ethylene Oxide and is considered sterile if the package is not opened or damaged. It is intended for Single Patient Use Only. Do not attempt to clean or re-sterilize the device. After use this device may be a potential biohazard. Handle in a manner that will prevent accidental contamination. Do not use a device that has been damaged or if the package is open or damaged.
7. Nitinol is a nickel titanium alloy. Possible reaction may occur for those patients who exhibit sensitivity to nickel.

Rx Only: CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

PRECAUTION: Care should be observed when using this device for removal of a large fibrin sheath in order to minimize risk of pulmonary embolism.

POTENTIAL COMPLICATIONS:

1. Potential complications associated with foreign body retrieval devices in arterial vasculature include, but are not limited to:
 - Embolization
 - Stroke
 - Myocardial infarction (depending upon placement)
2. Potential complications associated with snare retrieval devices in venous vasculature include, but are not limited to:
 - Pulmonary embolism
3. Other potential complications associated with foreign body retrieval devices include, but are not limited to:
 - Vessel perforation
 - Device entrapment

Catheter damage can occur when attempting fibrin sheath stripping on small French size diameter catheters. (See WARNINGS) Incidence of pulmonary embolism after fibrin sheath stripping may occur. (See PRECAUTION).

Prepare the EN Snare System:

- Select the appropriate Snare diameter range for the site in which the foreign body is located. The Snare diameter range should approximate the size of the vessel in which it will be used.
1. Remove the Snare and Snare Catheter from their hoop holders and inspect for any damage.
 2. Remove the Introducer and Torque Device from the proximal end of the snare shaft.
 3. Load the Snare into the Snare Catheter by inserting the proximal end of the snare into the distal (non-hubbed) end of the Snare Catheter, until the proximal end of the Snare shaft exits the hub and the loops can be retracted into the distal end of the Snare Catheter.
 4. Test and inspect the device by extending and retracting the snare loops through the distal end of the snare catheter 2-3 times, while carefully examining the Snare Catheter, radiopaque band and the device for any damage or defects.
 5. When appropriate, the system (Snare and Snare Catheter) can be advanced to the desired site as a single unit assembled as described above.

Alternative Preparation of the EN Snare :

If the Snare Catheter is already positioned within the vasculature, the provided Introducer (located on the proximal end of the Snare and just distal to the Torque Device) may be used to position the Snare in the indwelling Snare Catheter.

1. Remove the Snare from the protective holder and inspect for any damage.
2. Move the provided Introducer (located on the proximal end of the Snare, just distal to the Torque Device) distally until the loops of the Snare are enclosed within the tubing portion of the Introducer.
3. Insert the distal end of the Introducer into the hub of the indwelling Snare Catheter until resistance is felt. This will indicate the tip of the Introducer is properly aligned with the inner lumen.
4. Hold the Introducer as straight as possible, grasp the shaft of the Snare just proximal to the hub of the Introducer and advance the Snare until it is secure within the lumen of the Snare Catheter. The Introducer can be removed by first removing the Torque Device and pulling the Introducer off the proximal end of the Snare's shaft.

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1. If present, remove the indwelling delivery catheter.
2. If a guidewire is in a patient at the location of a foreign body, advance a Snare Catheter over the guidewire to the desired location. Then remove the guidewire and advance the Snare through the Snare Catheter. Alternatively, cinch one loop of the Snare over the proximal end of the guidewire and advance the entire system (Snare and Snare Catheter assembly) into a guide catheter or introducer sheath until the distal end of the Snare Catheter is positioned proximal to the foreign body.
3. If a guidewire is not present, pull the Snare into the distal end of the Snare Catheter and advance through a guide catheter or introducer sheath until it is positioned proximal to the foreign body. Alternatively, collapse the Snare loops by pulling the device into the distal end of the Snare Introducer. Place the tapered end of the Snare Introducer into the proximal (hub) end of the Snare Catheter, guide catheter or sheath and advance the Snare forward maintaining constant contact between the Introducer and Snare Catheter hub.
NOTE: When attempting to utilize guide catheters or sheaths not specifically manufactured for use with the EN Snare System, it is important to test product compatibility prior to use.
4. Gently push the Snare shaft forward to completely open the loops. The loops are then slowly advanced forward, and may be rotated if desired, around the proximal end of the foreign body. Alternatively, the Snare may be advanced beyond the target location and the loops brought back around the distal end of the foreign body.
5. By advancing the Snare Catheter, the loops of the device are closed to capture the foreign body. (Note that attempting to close the loops by pulling the Snare into the Snare Catheter will move the loops from their position around the foreign body.)
6. To manipulate a foreign body, maintain tension on the Snare Catheter to retain the hold on the foreign body, and move the Snare and Snare Catheter together to manipulate a foreign body to the desired position.
7. To retrieve a foreign body, maintain tension on the Snare Catheter and move the Snare and Snare Catheter assembly together proximally to, or into a guide catheter or sheath. The foreign body is then withdrawn through or together with the guiding catheter or introducer sheath. Withdrawal of large foreign bodies may require the insertion of larger sheaths, guiding catheters, or a cut-down at the peripheral site.

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1. Using standard technique, prepare a femoral vein approach, advance the selected Snare to the inferior vena cava or right atrium.
2. Advance a .035" guidewire through the end port (distal or venous port if more than one lumen) of the indwelling catheter and into the inferior vena cava or right atrium.
3. Position one of the Snare loops around the guidewire.
4. Advance the Snare over the distal end of the catheter to a position proximal to the fibrin sheath.
5. Close the Snare around the catheter and continue applying light traction while gently pulling the Snare down toward the distal end of the catheter over the end ports.
6. Repeat steps 4 & 5 until the catheter is free of fibrin sheath.

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1. Introduce the Snare at a patent venous access site and position in the vasculature at the desired site.
2. Open the Snare loops to provide a target to guide an entry needle into the desired venous access site.
3. Introduce a guidewire through the needle and through the Snare loops. Remove the needle.
4. Close the Snare over the guidewire by advancing the Snare Catheter.
5. Pull the guidewire into the desired location.

EN Snare FOREIGN BODY RETRIEVAL DEVICE SYSTEM*

Description	Snare Diameter Range	Snare Length	Catheter Size	Catheter Length
Mini EN Snare System	2-4 MM	175 CM	3.2 F	150 CM
Mini EN Snare System	4-8 MM	175 CM	3.2 F	150 CM
Description	Snare Diameter Range	Snare Length	Catheter Size	Catheter Length
EN Snare System	6-10 MM	120 CM	6 F	100 CM
EN Snare System	9-15 MM	120 CM	6 F	100 CM
EN Snare System	12-20 MM	120 CM	6 F	100 CM
EN Snare System	18-30 MM	120 CM	7 F	100 CM
EN Snare System	27-45 MM	120 CM	7 F	100 CM

* Every EN Snare System contains: (1) Snare, (1) Snare Catheter, (1) Introducer and (1) Torque Device.

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INTENDED USE:

EN Snare Endovascular Snare System is designed for use in the cardiovascular system or hollow viscus to retrieve and manipulate foreign bodies. Manipulation procedures include indwelling venous catheter repositioning, indwelling venous catheter fibrin sheath stripping and central venous access venipuncture procedure assistance.

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PRINCIPLE OF OPERATION:

EN Snare Endovascular Snare Systems is used by clinicians to manually capture, manipulate or retrieve objects in patient's cardiovascular system or hollow viscus, which is typically performed under fluoroscopy. The snare is loaded by inserting the proximal end of the snare into the distal end of the snare catheter, until the proximal end of the snare shaft exits the hub. After that the assembled system is advanced to target site. Slowly advance the loop forward, and around the proximal end of the foreign body. Once the loop is around the foreign body, advance the catheter to cinch the loop closed and grab the foreign body. To manipulate a foreign body, maintain tension on the catheter to retain the hold on the foreign body. To retrieve a foreign body move the snare and catheter together proximally or inside the guiding catheter or sheath.

Prior to use the appropriate snare diameter range for the site in which the foreign body is located is selected. The snare diameter range should approximate the size of the vessel in which it will be used.

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Handwritten signature and date:
DAN CLARK
JUL 18, 2017

MEDICAL DEVICE SPECIFICATIONS:

Catalog numbers	Loop diameter (mm)	Snare length (mm)	Catheter OD (in)/(mm)	Catheter ID (in)/(mm)	Catheter length (mm)	Insertion tool OD (in)/(mm)	Insertion tool ID (in)/(mm)	Insertion tool length (mm)	Torque device diameter (in)/(mm)	Torque device length (mm)
EN Snare Endovascular Snare System										
EN1003004	2-4	1750	0.045/1.14	0.030/0.76	1500	0.081/2.06	0.050/1.27	46	(0,018 – 0,038) / (0.46 – 0.97)	40
EN1003008	4-8	1750	0.045/1.14	0.030/0.76	1500	0.081/2.06	0.050/1.27	46	(0,018 – 0,038) / (0.46 – 0.97)	40
EN2006010	6-10	1200	0.083/2.11	0.062/1.57	1000	0.090/2.29	0.050/1.27	84	(0,018 – 0,038) / (0.46 – 0.97)	40
EN2006015	9-15	1200	0.083/2.11	0.062/1.57	1000	0.090/2.29	0.050/1.27	84	(0,018 – 0,038) / (0.46 – 0.97)	40
EN2006020	12-20	1200	0.083/2.11	0.062/1.57	1000	0.090/2.29	0.050/1.27	84	(0,018 – 0,038) / (0.46 – 0.97)	40
EN2007030	18-30	1200	0.095/2.41	0.074/1.88	1000	0.090/2.29	0.050/1.27	84	(0,018 – 0,038) / (0.46 – 0.97)	40
EN2007045	27-45	1200	0.095/2.41	0.074/1.88	1000	0.090/2.29	0.050/1.27	84	(0,018 – 0,038) / (0.46 – 0.97)	40

DEVICE MATERIALS:

The materials used in the manufacture of the components of EN Snare Endovascular Snare System are identified below:

Component Material Specifications		
Snare Materials		Type and duration of the contact
Component	Material	
Wire core (guidewire)	Nitinol (titanium mononickelide NiTi) SE508	Externally communicating device, contact with circulating blood, limited contact duration (≤ 24 hours)
Loop cables of (mini) snare (loop diameter 2-4 mm, 4-8 mm)	Nitinol SE508 and platinum / iridium (90% / 10%)	
Loop cables of (standard) snare (diameter 6 up to 45 mm)	Nitinol SE508 and platinum	
Crimp stem	Stainless steel 304 (UNS S30400)	
Wire solder	Tin / silver (95% / 5%)	
Bonding material (adhesive)	Loctite 4014	
Insertion tool (introducer) materials		
Tubing	Polypropylene, Flint Hills III P4G3Z-039	No patient contact
Sleeve with Luer hub	Polypropylene, PolyOne #CC10200699 polymer	
Torque device materials		
Cap	Polycarbonate, Polystyrol 2081-15-FC030004	No patient contact
Body	Polypropylene III, FHR 13R9A	
Collet sleeve	Brass	
Catheter materials		
Catheter tubing		Externally communicating device, contact with circulating blood, limited contact duration (≤ 24 hours)
Outer layer	Fluorinated ethylene propylene, Pebax 7233-SA01-MED	
Inner layer	Polytetrafluoroethylene (Zeus)	
Marker band	Platinum / iridium (90%/10%)	
Luer hub	Polycarbonate, Polystyrol 2081-15-FC 030004	Externally communicating device, contact with circulating blood, limited contact duration (≤ 24 hours)
Catheter sleeve (cannula)	Pebax 3533-SA01-MED	

SAFETY REQUIREMENTS:

EN Snare Endovascular Snare System is safe and effective for its intended use as labeled.

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Transportation conditions: at temperature from -50°C to +50°C and relative humidity up to 100%.

After long-term exposure to subzero temperatures the device should be conditioned in its packaging for several days at storage temperature.

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
PACKAGING:

Primary packaging (inner and outer):


The snare device, insertion tool and torque device (insertion tool is threaded onto the proximal end of the snare shaft to minimize damage during shipping/storage) are placed in polypropylene spiral hoop (dispenser / protective retaining ring) and then inside a separate Tyvek/Nylon pouch that is sealed using a validated heat seal process. A label containing information on the device is applied. The catheter is placed in polypropylene spiral hoop (dispenser / protective retaining ring) and then inside a separate Tyvek/Nylon pouch that is sealed using a validated heat seal process. A label containing information on the device is applied.

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Shipping package:

Then several (up to five) cartons are placed into a large corrugated shipper box delivered to customers.

LABELING:

The endovascular snare system is labeled as required by Annex I, Council Medical Devices Directive 93/42/EEC, and graphical symbols used in the labeling comply with EN 980:2003 requirements.

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
- manufacturer's name
- device trade name
- device dimensions (loop diameter and snare length)
- graphical scheme of pouch contents
- bar code
- CE marking
- catalog number (REF)
- lot number
- symbol "Best before", shelf life
- symbol "Sterile, sterilized with ethylene oxide"
- symbol "Do not use if package is damaged"
- symbol "Single use only"
- symbol "Non-pyrogenic"
- symbol "Consult accompanying documents"
- manufactured by... name and full address of the manufacturer

- EC representative address
- manufacturing date

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- manufacturer's name

- device trade name
- device dimensions (loop diameter and snare catheter length)
- graphical scheme of pouch contents
- bar code
- CE marking
- catalog number (REF)
- lot number
- symbol "Best before", shelf life
- symbol "Sterile, sterilized with ethylene oxide"
- symbol "Do not use if package is damaged"
- symbol "Single use only"

- symbol "Non-pyrogenic"
- symbol "Consult accompanying documents"
- manufactured by... name and full address of the manufacturer
- EC representative address

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METHODS AND MEANS OF DISINFECTION AND PRE-STERILIZATION CLEANING:
 Not applicable. EN Snare Endovascular Snare System is a sterile device for single use only.
 Used devices are subject to disinfection prior to disposal.

STERILITY:

EN Snare Endovascular Snare System is sterilized using a traditional ethylene oxide sterilization process. Prior to sterilization the device is assembled and packaged into dispensers and pouches (inner packaging). After that the system is sterilized. The device is not intended to be sterilized by the user.

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Complies with essential requirements of EC Medical Devices Directive 93/42/EEC.

Products of Merit Company meet ISO 13485:2003 requirements.

WARRANTY:

The manufacturer guarantees medical device sterility, as well as its parameters and performance **within the stated shelf life**, i.e. three years after manufacturing date, provided individual packaging remains intact and storage, transportation and operation conditions are observed according to the instructions for use.

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EN Snare Endovascular Snare System contains no parts requiring maintenance during operation. Thus, field repair is not performed.

SHELF LIFE:

3 years.

RECLAMATIONS:

Reclamations should be sent to the manufacturer's authorized representative in Russia:

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Instructions for Use Supplement
ONE Snare
Consisting of 3 pages
Dated June 15, 2017


MERIT MEDICAL SYSTEMS, INC
1600 WEST MERIT PARKWAY
SOUTH JORDAN UTAH 84095
PHONE 801-253-1600
FAX 801-253-1688
www.merit.com

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Notary Public

My Commission Expires: June 27, 2021

Supplement (specification) to the instructions for use No. 1 dated June 15, 2017

Endovascular snare system in the following implementations:

ONE Snare endovascular snare system

ONE Snare Microsnare endovascular microsnare system

John Clark
JUL 13, 2017


1) LABELING section: interpret the symbol as:

	<i>Indicates that the medical device is non-pyrogenic</i>
---	--


Add the following information:

The shipping package labeling contains the following information:

- name or trademark of the manufacturer;
- commercial name of the device;
- dimensions of the device (loop diameter and snare length);
- pouch contents as a graphical scheme;
- barcode;
- catalog number (REF);
- lot number (LOT);
- "Use before" symbol, expiry date;
- "For single use only" symbol;
- "Do not use if package is damaged" symbol;
- "Sterile, sterilized with ethylene oxide" symbol;
- "Consult supporting documents" symbol;
- "Non-pyrogenic" symbol;
- "RX only" symbol;
- international CE mark;
- number of pieces in the package;
- manufactured in... name and full address of the manufacturer;
- address of the EU representative of the manufacturer.

Label layout for Russia

2) TECHNICAL SPECIFICATIONS OF THE MEDICAL DEVICE section: add the following information:

	<i>Snare loop mass Not more than (g)</i>	<i>Catheter mass Not more than (g)</i>	<i>Insertion tool mass Not more than (g)</i>	<i>Torque device mass Not more than (g)</i>
ONE500	2.25	3.20	0.45	1.90
ONE1000	2.80	3.20	0.45	1.90
ONE1500	3.00	3.20	0.45	1.90
ONE2000	3.00	3.20	0.45	1.90
ONE2500	3.05	3.20	0.45	1.90
ONE3000	3.10	3.20	0.45	1.90
ONE3500	3.20	3.20	0.45	1.90
ONE1001	1,75	2,20	0,45	1,90
ONE2501	1,75	2,20	0,45	1,90
ONE200	1,00	2,20	0,45	1,90
ONE201	1,20	2,20	0,45	1,90
ONE400	1,20	2,40	0,45	1,90
ONE401	1,20	2,40	0,45	1,90
ONE700	1,25	2,40	0,45	1,90
ONE701	1,25	2,40	0,45	1,90

Strength characteristics of ONE Snare and ONE Snare Microsnare implementations:

Microsnare 4 mm and 7 mm:

- *Minimum pull-off force of loop wires from the crimping shaft – 2.25 pounds-force (10.01 N).*
- *The crimping shaft must remain attached to the nitinol wire core.*

Snare (standard) 10 mm:

- *Minimum pull-off force of loop wires from the crimping shaft – 5.0 pounds-force (22.24 N).*
- *The crimping shaft must remain attached to the nitinol wire core.*

Snare (standard) 15 mm, 20 mm, 25 mm, 30 mm and 35 mm:

- *Minimum pull-off force of loop wires from the crimping shaft – 8.0 pounds-force (35.59 N).*
- *The crimping shaft must remain attached to the nitinol wire core.*

Introducer/snare pull-off force:

ONE Snare endovascular snare system 4 F, 6 F – < 1.0 pound-force (4.45 N)

ONE Snare Microsnare endovascular microsnare system 2.3-3.0 F – < 0.5 pounds-force (2.22 N)

3) Remove the following:

The introducer is made of the same materials as the catheter, but does not feature a marker band...

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