
Carestream

CS 3500 INTRAORAL SCANNER with Accessories

User and Installation Guide



Robert C. Meagher
Senior Director
World Wide Regulatory Affairs
Carestream Health, Inc.
150 Verona Street
Rochester, New York 14608
Telephone 585-627-6528
STATE OF NEW YORK

COUNTY OF MONROE

Signed and attested before me on February 20, 2017 by

CHRISTINE A. BENOIT-AVERILL
NOTARY PUBLIC-STATE OF NEW YORK
No. 01BE6320594
Qualified in Monroe County
My Commission Expires March 09, 2019

Notary Public

User and Installation Guide

CS 3500 INTRAORAL SCANNER with Accessories

Notice

No part of this publication may be reproduced, stored in a retrieval system, translated to another language, or transmitted in any form by any means, electronic, mechanical, photocopied, recorded, or otherwise, without prior written permission.

The information in this document is subject to change. Neither Carestream Health, Inc. nor any of its subsidiaries shall be liable for errors contained herein or for incidental damages in conjunction with the furnishing, performance, or use of this material.

We recommend that you thoroughly familiarize yourself with this Guide in order to make the most effective use of your system.

WARNING: We recommend that you consult the "Technical information" before using the CS 3500 Systems.

Carestream Health is a trademark of Carestream Health, Inc.

All other trademarks and registered trademarks are the property of their respective holders.

US Federal law restricts this device to sale by or on the order of a dentist.

I. CS 3500 INTRAORAL SCANNER with Accessories:

1 CS 3500 INTRAORAL SCANNER - 1 pc.

2 Operational documentation - 1 pc.

II. Accessories:

1 REGULAR TIPS CS 3500, standard size – 1 package (5 pieces)

2 REGULAR TIPS CS 3500, small size - 1 package (5 pieces)

3 TOOTH PLASTER MODELCS 3500 - 1pc.

4 COLLAR PACKING CS 3500 - 1 package (10 pieces)

5 Adapter CS 3500 - 1pc.

6 WALL HOLDER CS 3500 - 1pc.

7 DESK HOLDER CS 3500 - 1pc.

Manual Name: *CS 3500 User and Installation Guide*

Part Number: 9H4853

Revision Number: 04

Print Date: 2014 - 07

Contents

Chapter 1	Conventions in This Guide.....	1
Conventions in This Guide		
Chapter 2	CS 3500 Overview	3
CS 3500 Overview	CS 3500 Holders Overview	5
Chapter 3	Minimum Computer System Requirements.....	7
CS 3500 Software Overview	General Software Overview	7
	Imaging Software	7
	CS 3500 Acquisition Interface Overview	8
	Toolbar Overview	9
	Scanning Area Selection	12
	Preparation Check	12
	Arch Display	13
	Preferences Dialog Boxes Overview.....	13
	General Preferences	13
	Scanner Preferences	14
	Tools Preferences.....	15
Chapter 4	Setting Up the CS 3500	17
Setting Up the CS 3500	Testing Image Quality and Accuracy.....	23
	Installing the CS 3500 Desktop Holder	24
	Installing the CS 3500 Wall Mount Holder	25
Chapter 5	Using the Imaging Software.....	27
Getting Started	Accessing the Acquisition Interface.....	27
	Accessing the Acquisition Interface Using the Imaging Software	27
	Accessing the Acquisition Interface from Your Practice Management Software	27
	Setting Up Preferences.....	28
	Running the CS 3500 Demo	29

	CS 3500 Tools and Tips	30
	Tooth Preparation	30
	CS 3500 Position	30
	Acquisition Modes	30
	Important Scanning Procedures—Restoration	
	Acquisition.....	31
	Scanning One Surface at a Time.....	31
	Scanning Each Tooth Individually	31
	Important Scanning Procedures—Orthodontic	
	Acquisition.....	31
	Scanning One Surface at a Time.....	31
	Scanning Each Tooth Individually	31
	Scanning Metal, Crowns, or Holes	32
	Preparing the CS 3500	33
	Using the Collar to Scan Anterior Teeth	34
Chapter 6	Scanning Teeth on the Upper and Lower Jaw	35
Acquiring a	Scanning the Buccal Bite Registration	39
3D Model for	Refinement.....	41
Restoration	Preparation Check.....	44
	Drawing Margin Lines	46
	Acquiring 2D Images.....	48
Chapter 7	Scanning Teeth on the Upper and Lower Jaw	49
Acquiring a	Scanning the Buccal Bite Registration	53
3D Model for	Refinement.....	55
Orthodontics	Preparation Check.....	57
	Acquiring 2D Images.....	58
Chapter 8	Cleaning, Disinfecting, and Sterilizing	59
Maintenance		
Chapter 9	Manufacturer's Address.....	61
Contact Information	Authorized Representatives	61

1 Conventions in This Guide

Conventions in This Guide

The following special messages emphasize information or indicate potential risk to personnel or equipment:

WARNING: Warns you to avoid injury to yourself or others by following the safety instructions precisely.

CAUTION: Alerts you to a condition that might cause serious damage.

Important: Alerts you to a condition that might cause problems.

Note: Emphasizes important information.

Tip: Provides extra information and hints.

2 CS 3500 Overview


The CS 3500 INTRAORAL SCANNER with Accessories is designed for non-contact measurement of 3D surface of the tooth directly into the mouth of the patient, and allows dentists to quickly scan the mouth and send 3D model in the laboratory.

The CS 3500 is designed to acquire 2D and 3D still images in the following modes:

- Lower jaw
- Upper jaw
- Buccal bite registration

CS 3500 Overview

Figure 1 CS 3500 Components Overview

- 1 Reusable tip
- 2 Acquisition buttons

The tip can be installed facing upward or downward. Tips come in two sizes: regular and small. Tips are autoclavable for up to 20 times.

Press to acquire a 3D image.

- The CS 3500 is in auto acquisition mode by default. When the CS 3500 is stable and the image is clear, the acquisition begins.
- To make a manual acquisition, press the acquisition button once to acquire a single image. To keep the CS 3500 in manual acquisition mode, deselect Auto Capture in the Preferences window. See "Scanner Preferences" for more information.
- When the guidance system is enabled in the Preferences window, the LED ring surrounding the buttons changes color.

- -- Indicates an image was acquired correctly (Feedback option).
- -- Indicates there is a proper overlap, and you can make an acquisition (Guide option).
- -- Indicates an image was not acquired correctly (Feedback option).
- -- Indicates there is not enough overlap or the angle of the CS 3500 is incorrect (Guide option).

If you do not enable the guidance system, the LED ring always displays in green.

3 Mode indicators

•  Lower jaw scan mode

•  Upper jaw scan mode

•  Buccal bite registration mode

4 Mode button

Press the button to switch between different modes.

5 Power button

- Press one second to power ON.
- Press three seconds to power OFF.

Power indicator LED:

-  Indicates the CS 3500 is active.
-  Indicates the CS 3500 is in sleep mode, holder mode, standby mode, or other inactive mode.
-  Indicates the power is OFF.

6 USB connection indicator

-  Indicates the connection is established.
-  Indicates the connection is disconnected.

7 USB cable

One end of the cable is hard-wired to the bottom of the CS 3500, and the other end is inserted in the computer's USB port.

8 Power adapter

One end of the adapter is inserted in the jack on the USB plug, and the other end is inserted in the power outlet.

9 Service Calibration Diagnostic Tip

The service calibration diagnostic tip enables you to verify the accuracy of the CS 3500.

10 Reusable Collar

The collar is placed over the CS 3500 tip and can help you acquire images of the anterior teeth. The collar is autoclavable for up to 20 times.

CS 3500 Holders Overview

The CS 3500 has both a desktop and a wall mount holder. Place the CS 3500 in the holder when you are not using it.

Figure 2 CS 3500 and Desktop Holder



Figure 3 CS 3500 and Wall Mount Holder

Note: The CS 3500 goes into holder mode when inserted into the holder. To use it again, take it out of the holder.

Note: The CS 3500 goes into sleep mode when left idle for one minute (if set on a desk, for example). To use it again, pick it up or press the power button.

Note: The CS 3500 goes into standby mode when the tip is removed. To use it again, put on the tip.

Note: When the CS 3500 is idle for more than an hour, it powers off.

3 CS 3500 Software Overview

Minimum Computer System Requirements

For the minimum computer system requirements, see the CS 3500 «**Technical information**»

Important: It is MANDATORY to check that your computer system configuration is compatible with the computer system requirements for the CS 3500 software.

General Software Overview

The CS 3500 system operates with the following software:

- Imaging software
- CS 3500 Acquisition interface

Imaging Software

Carestream's imaging software is a user-friendly working interface that was designed and developed as a common imaging platform for all our dental digital systems.

For more details, see:

- Dental Imaging Software online help and the Dental Imaging Software Quick User Guide
- OR
- CS Orthodontic Imaging software online help and the CS Orthodontic Imaging Quick User Guide
- OR
- CS OMS Imaging software online help and the CS OMS Imaging Quick User Guide

CS 3500 Acquisition Interface Overview


The CS 3500 Acquisition interface enables you to acquire images in two ways:

- Partial arch scan: Several teeth in the preparation area on both the upper and lower jaws, and a buccal bite registration.
- Full arch scan: Upper jaw, lower jaw, and buccal bite registration.

Note: In some cases, you can acquire images of a single arch (partial or full) and not obtain a buccal bite registration (for example, if there are no teeth in the opposing arch), but it is recommended that you acquire both arches and a buccal bite registration when possible.

Figure 4 CS 3500 Acquisition Interface Overview

1 Title Bar: Displays:

- Patient's name and type of acquisition
-  Minimizes the Acquisition interface.
-  Maximizes or restores down the Acquisition interface.
-  Closes the Acquisition interface.

- 2 CS 3500 Status Bar: Indicates the current status of the CS 3500 and enables you to change the current acquisition type.
- 3 Video Preview Screen: Displays live video of the scan of the patient's teeth, as well as the alignment guide, when enabled.
- 4 Tutorial Animation: Demonstrates the image acquisition process. Float your mouse over the image of the jaw, and click Show Scanning Sample to launch the demo. Click OK to stop the demo and close the window.
- 5 Progress Bar: Indicates the status of the image acquisition process.
- 6 3D Model Display Screen: Displays the 3D model created from the scan.
- 7 Image Toolbar: Enables you to select and manage 3D images.
- 8 Settings Toolbar: Enables you to import and export files, set preferences, and access system information and online help.
- 9 Acquisition Feedback Indicator: Indicates if the acquisition is successful or if it failed; or if the CS 3500 was not steady.
- 10 Acquisition Mode Toolbar: Enables you to select the acquisition mode.

Toolbar Overview

Acquisition Mode Toolbar

The acquisition mode toolbar enables you to select the acquisition mode.

Lower Jaw button: Acquires a 3D image of the lower jaw.

Upper Jaw button: Acquires a 3D image of the upper jaw.

Buccal Bite Registration button: Acquires a 3D image of the bite.

Image Toolbar

The image toolbar enables you to select and manage images.

Intraoral Mode button. Acquires two-dimensional intraoral images. You can acquire a maximum of 15 2D images.

Switch button: Changes the acquisition mode from upper to lower or vice versa, if you accidentally scan teeth on the wrong jaw.

True Color button: Displays the 3D model in the actual color of the patient's hard and soft tissues.

Freeze button: Prevents the 3D model from being resized and rotated.

Zoom Fit button: Scales the 3D model to the best size to fit the display region.

Draw Margin Line button: Draw and save the margin line on a snapshot of the 3D model. This button is available after the model has been refined.

Margin Line Viewer button: View and delete previously saved margin line images. This button is available after the model has been refined.

Select by window button: Selects the view or views by dragging the mouse to draw a selection box.

Free Cut button: Selects and deletes unwanted data, such as soft tissue, in the image using a free-form shape. This button is available after the model has been refined.

Cut button: Selects and deletes unwanted data, such as soft tissue, in the image. This button is available after the model has been refined.

Delete button: Deletes the selected view/views.

Delete All button: Deletes all views of the current mode.

Restore button: Restores the last view or views that were deleted.

Reset button: Reverses all deletions and returns to initial acquisition. This button is available after the model has been refined.

Next button: Moves to the next step in the acquisition workflow.

Previous button: Moves to the previous step in the acquisition workflow.

Done button: Sends the completed image to the Dental Imaging Software, CS Orthodontic Imaging software, or CS OMS Imaging software.

Settings Toolbar

The settings toolbar enables you to import and export files, set preferences, and access system information and online help.

Export button: Accesses the Documents library, where you can select the folder to which you want to export the image file. Use the export function if you want to save the acquired images for later use.

Import button: Accesses the Documents library, where you can select the folder from which to import an image file. Use the import feature to access an acquisition you have previously saved.

Preference Setting button: Accesses the Preferences dialog box.

About button: Accesses system information, including the firmware and software version, and CS 3500 serial number.

Help button: Accesses the online help.

CS 3500 Status Icons

The CS 3500 status icons display the current status of the CS 3500.

Select Acquisition Type button: Accesses the orthodontic acquisition mode. You can toggle between Acquisition for Orthodontics and Acquisition for Restoration by clicking the button.

Select Acquisition Type button: Accesses the restoration acquisition mode. You can toggle between Acquisition for Restoration and Acquisition for Orthodontics by clicking the button.

In holder icon: Indicates the CS 3500 is in the holder.

Sleep mode icon: Indicates the CS 3500 is in sleep mode.

Standby icon: Indicates the CS 3500 is in standby mode.

Overheating icon: Indicates the CS 3500 is overheating. If this icon is displayed, place the CS 3500 in the holder for 5 to 10 minutes. The CS 3500 will become inactive and cool down.

Acquisition Feedback Status Bar

The acquisition feedback status bar enables you to see if the acquisition is successful.

Acquisition successful icon: Indicates the image was acquired successfully.

Acquisition failed due to movement icon: Indicates the CS 3500 was not stable. Reacquire the image.

Acquisition failed icon: Indicates that there is not enough overlap between the last image and the newest image. Reacquire the image, being sure to overlap with the previous tooth by approximately 30%.

Scanning Area Selection

The Scanning Area Selection area is displayed after the model has been refined and enables you to indicate the teeth that are part of the model.

Note: The Scanning Area Selection area is displayed only during a restoration acquisition.


Scanning Area Selection: Click  to display the arch diagram.

Click the Adult or Child icon to display the appropriate arch, and select the teeth that are in the model, including the upper and lower jaw. Click the Reset button to clear the selection.

Preparation Check

The Preparation Check section is displayed after the model has been refined and enables you to verify important restoration data before continuing with the process.


Restoration Preparation Check: A list of items you must check on the restorative 3D model to verify the bite registration and preparation design.

You must examine the appropriate area of the 3D model and select the corresponding option from the list before you can move to the next step.

Verify the following:

- Occlusal space
- Undercut
- Margin line

Orthodontic Preparation Check: A list of items you must check on the orthodontic 3D model to verify the accuracy of the model before sending it to the imaging software.

You must examine the appropriate area of the 3D model and select the corresponding option from the list before you can move to the next step.

Verify the following:

- No holes in the image
- Jagged boundaries removed
- Irregular soft tissue removed

Arch Display

The Arch Display section is displayed during refinement selection and after the model has been refined and enables you to view or hide the upper and lower jaws.

Upper: Click the upper jaw button to view or hide the upper jaw.


Lower: Click the lower jaw button to view or hide the lower jaw.

Preferences Dialog Boxes Overview

The Preferences dialog boxes enable you to select the CS 3500 Acquisition interface settings. It is recommended that you set up your preferences before using the CS 3500.

General Preferences

Figure 5 General Preferences


The General Preferences window enables you to select these options:

- Language: Select the menu language.
- Help language: Select the online help system language.
- Enable Demo: Select to enable the software demo mode the next time the CS 3500 Acquisition interface is opened.

Scanner Preferences

Figure 6 Scanner Preferences

The Scanner Preferences window enables you to select these options:

- Auto capture: Select to enable auto acquisition. Click the slider to adjust the acquisition sensitivity.
- Fast capture: Select to enable the fastest capture speed. Fast capture is recommended for experienced users who are familiar with using the CS 3500 and maintaining a 30% overlap with each image.

Note: When Fast capture is selected, the guidance systems are disabled.

- Guidance: Select to enable the guidance systems.
 - Feedback: Select to enable the projection light and the LED ring surrounding the acquisition button. These indicators tell you whether the *last image* was acquired correctly (green) or if you need to reacquire the image (orange).
 - Guide: Select to enable the projection light, LED ring, and the alignment guide. These indicators guide you in acquiring your *next image* by helping to ensure the correct amount of overlap with each image. A green indicator means you can proceed with the acquisition; an orange indicator means there is not enough overlap or the angle of the CS 3500 is incorrect.

Note: The guidance systems are disabled when Fast capture is selected.

- Capture beep: Select to enable a beeping sound each time a scan is attempted.
- Enable image review: Select to display color images immediately after an acquisition.


Note: When Enable image review is selected, acquisition speed is reduced.

- Default acquisition mode: Select the jaw you will always scan first.

Tools Preferences

Figure 7 Tools Preferences

The Tools Preferences window enables you to select these options:

- Create Support File: Click to select the folder where the log files for service information are saved.
- Diagnose Scanner: After attaching the service calibration diagnostic tip, click to run a diagnostic test to verify the accuracy of the CS 3500. See "Testing Image Quality and Accuracy" for more information.

4 Setting Up the CS 3500

Setting Up the CS 3500

To set up the CS 3500, follow these steps:


- 1 Install the imaging software on the computer. To install the Dental Imaging Software, CS Orthodontic Imaging software, or CS OMS Imaging software, see the relevant installation and quick user guides.

- 2 Insert the CS 3500 Installation Wizard flash drive  into the USB drive. The AutoPlay window is displayed.

Note: If the AutoPlay window does not automatically display, navigate to the CS 3500 flash drive in Windows Explorer and double-click wizard.exe.


- 3 Select Run wizard.exe. The Autorun window is displayed.

Important: To view CS 3500 documentation, click Documentation. To view training videos about the CS 3500, click Training Videos.


- 4 Select a language from the drop-down list, and click Run Installation Wizard. The Installation Wizard window is displayed.

The wizard verifies that your computer meets the minimum requirements to run the CS 3500. Keep the following in mind:

-  — Indicates that the recommended requirements are met for that area.
-  — Indicates that the minimum requirements are met, but not the recommended requirements. You can continue with the installation, but the CS 3500 might not run as well as it could if the recommended requirements were met.
-  — Indicates that the CS 3500 might not run at all. Stop the installation and upgrade your computer before continuing. Click the Help link [Help?](#) for information about the requirement that is not met.


- 5 If there are no **X**, click Next. The wizard begins installing the software.


- 6 When the software is installed, you are prompted to copy the documentation and training videos.

- 7 Click Next. The Browse for Folder window is displayed.


- 8 Select a location for the documentation and training video folders. Click Make New Folder, type a folder name, and click OK. Shortcuts to this folder are added to your desktop.
- 9 When the documentation and training videos have been copied, the Installation Wizard Completed message is displayed.


- 10 Click Finish.


- 11 Firmly slide the tip onto the end of the CS 3500.


- 12 Insert the DC power supply (A) into the jack on the USB plug (B), and insert the power adapter (C) into an outlet.

- 13 Connect the USB cable to the computer.
- 14 Press the power button one second to power on the CS 3500. Ensure the power indicator turns orange.

15 Wait until the USB connection indicator turns blue.

Note: The power indicator turns green when the acquisition interface is open and the CS 3500 is ready to acquire images.

Testing Image Quality and Accuracy


After setting up the CS 3500, you should verify the accuracy of the CS 3500, using the service calibration diagnostic tip. We recommend that you run this test once per month.

In addition, you should verify the accuracy of the CS 3500:

- If the CS 3500 is dropped.
- If a restoration does not accurately fit in a patient's mouth.

To test the accuracy of the CS 3500, follow these steps:

- 1 Firmly slide the service calibration diagnostic tip onto the end of the CS 3500. Ensure that the arrow on the tip is pointing toward the acquisition button on the top of the CS 3500.

- 2 On your desktop, double-click  (Carestream's Dental Imaging Software) or  (CS Orthodontic Imaging software or CS OMS Imaging software), or open the imaging software from within your practice management software.
- 3 Find or create the patient record.
- 4 Double-click the patient record to access the imaging window or tab.
- 5 In the imaging window or tab, click  to open the CS 3500 Acquisition interface.
- 6 On the Acquisition interface, click . The Preferences dialog box is displayed.
- 7 Click . The Tools Preferences dialog box is displayed.
- 8 Press the power button for one second to power on the CS 3500.


Tip: If the CS 3500 is already powered on, verify that the Sleep or In Holder icons are not displayed in the upper-left corner of the Acquisition interface.

- 9 Click Diagnose Scanner. The test begins.
- 10 When the test is finished, a message is displayed, stating whether the CS 3500 has passed or failed the diagnostic test. If the CS 3500 has failed the test, click Create Support File, and save the .zip file to a folder you select. Contact your local service provider for assistance before using the CS 3500 again.

Installing the CS 3500 Desktop Holder

To install the CS 3500 desktop holder, follow these steps:

- 1 Select a clean surface area that you can access easily.
- 2 Clean the surface on which you are going to fasten the CS 3500 holder.
- 3 Attach the adhesive tape to the indentations in the holder (A).
- 4 Remove the adhesive tape protection (B).
- 5 Place the taped side of the CS 3500 holder on the clean surface (C) and press firmly several times for correct adherence. The maximum adherence is obtained after two hours.

Important: You must avoid all stress on the CS 3500 holder during the first two hours.

Important: Adhesive tape can lose its adherence capability due to environmental and other factors. Check the adherence regularly to ensure that it is secure.


Important: Carestream Health is not responsible for accidents or damage involving the CS 3500 falling due to loosening of the adhesive tape.

Installing the CS 3500 Wall Mount Holder

To install the CS 3500 wall mount holder, follow these steps:

- 1 Select an area that you can access easily.
- 2 Insert screws through the appropriate holes in the holder to affix it to a solid surface.

Important: If the holder is not properly installed, there is a risk that the holder can fall off the wall, resulting in damage to the CS 3500.

5 Getting Started

Using the Imaging Software

To use the Dental Imaging Software, see the Dental Imaging Software Quick User Guide and the online help delivered with the software.

To use the CS Orthodontic Imaging software, see the CS Orthodontic Imaging Quick User Guide and the online help delivered with the software.

To use the CS OMS Imaging software, see the CS OMS Imaging Quick User Guide and the online help delivered with the software.

Accessing the Acquisition Interface

You can access the CS 3500 Acquisition interface from inside the Dental Imaging Software, the CS Orthodontic Imaging software, the CS OMS Imaging software, or from inside your practice management software.

Accessing the Acquisition Interface Using the Imaging Software

To access the CS 3500 Acquisition interface from the Dental Imaging Software, the CS Orthodontic Imaging software, or the CS OMS Imaging software, follow these steps:

- 1 On your desktop, double-click  (Carestream's Dental Imaging Software) or  (CS Orthodontic Imaging software or CS OMS Imaging software).
- 2 Find or create the patient record.
- 3 Double-click the patient record to access the imaging window or tab.
- 4 In the imaging window or tab, click  to open the CS 3500 Acquisition interface.
- 5 Select your preference settings.

Accessing the Acquisition Interface from Your Practice Management Software

To access the CS 3500 Acquisition interface from your practice management software, follow these steps:

- 1 Open your practice management software.
- 2 Find the patient record and open it.
- 3 Open the imaging software.
- 4 Click  to open the CS 3500 Acquisition interface.
- 5 Select your preference settings.

Setting Up Preferences

Set up your preferences before using the CS 3500. To set up the preferences, follow these steps:

1 On the Acquisition interface, click . The Preferences dialog box is displayed.

2 Click  to set the general preferences.

3 Click  to set the CS 3500 preferences.

4 Click  to set the tools preferences.

You can now start acquiring 3D and 2D still images.

Running the CS 3500 Demo

To help you become familiar with the CS 3500 Acquisition interface, you can run a demo version that does not require the CS 3500 to be connected.

To run the demo version of the CS 3500, follow these steps:

- 1 In the Acquisition interface, click . The Preferences dialog box is displayed.
 - 2 Click Enable Demo and click OK.
 - 3 Close the Acquisition interface.
 - 4 Re-open the Acquisition interface. The Select a Demo Data File window is displayed.
 - 5 Select the demo type and click Open. Wait for the demo data to finish loading.
 - 6 Click  once to simulate a single image acquisition on the upper jaw. Press Shift and click  to acquire multiple images for the demo, and continue with the acquisition process as described in "Scanning Teeth on the Upper and Lower Jaw" for restoration acquisitions or "Scanning Teeth on the Upper and Lower Jaw" for orthodontic acquisitions.
 - 7 When you are finished with the demo, do one of the following:
 - Click , deselect Enable Demo, and click OK.
 - Close the Acquisition interface.
- OR
- Close the Acquisition interface.
 - When you open the Acquisition interface next time, the Select a Demo Data File window is displayed. Close the window, and a prompt asking if you want to close the Acquisition interface is displayed.
 - Select Run live mode next time, and click Yes.

Note: Any time you open the Acquisition interface without the CS 3500 connected, a prompt asking if you want to run in demo mode is displayed.

- Re-open the Acquisition interface.

CS 3500 Tools and Tips

Use the following tools and tips to help you acquire quality images:

- **Projection Light:** When the guidance system is enabled with Feedback selected, a light is emitted from the tip of the CS 3500 to indicate if the acquisition is correct. The light is green if the image is scanned correctly and orange if the image is not scanned correctly.

When the guidance system is enabled with Guide selected, a green light means you can proceed with the acquisition; an orange light means there is not enough overlap or the angle of the CS 3500 is incorrect.

- **LED Ring Surrounding Acquisition Button:** When the guidance system is enabled with Feedback selected, the LED ring surrounding the acquisition button changes color to indicate if the acquisition is correct. The LED ring is green if the image is scanned correctly and orange if the image is not scanned correctly.

When the guidance system is enabled with Guide selected, a green LED ring means you can proceed with the acquisition; an orange LED ring means there is not enough overlap or the angle of the CS 3500 is incorrect.

- **Alignment Guide:** When the guidance system is enabled with Guide selected, alignment indicators guide you in acquiring your *next image* by helping to ensure the correct amount of overlap with each image. A green indicator means you can proceed with the acquisition; an orange indicator means there is not enough overlap or the angle of the CS 3500 is incorrect.
- **Capture Beep:** When the capture beep is enabled, the CS 3500 beeps each time a scan is attempted.
- **Reusable Collar:** Use the collar when scanning anterior teeth to improve the image quality.

Tooth Preparation

- Dry the teeth thoroughly before starting an acquisition.
- Re-dry the teeth as appropriate throughout the acquisition process.

CS 3500 Position

- To enable the antifog feature, place the tip on the CS 3500 and allow the CS 3500 to warm up for approximately three minutes before starting the acquisition.
- Hold the CS 3500 at a 90-degree angle to the occlusal surface of the teeth; hold the CS 3500 at a 45-degree angle to the buccal and lingual surfaces of the teeth.
- Rest the tip of the CS 3500 on the tooth surface to steady the CS 3500.
- When scanning a tooth, keep approximately a 30% overlap with the last image.

Acquisition Modes

- **Auto acquisition mode (default):** Hold the CS 3500 at a 90-degree angle to the occlusal surface of the teeth. When the CS 3500 is stable and the image is clear, the acquisition begins. You can rest the tip on the tooth surface during acquisition to steady the CS 3500. Live video is displayed on the video preview screen, and images are automatically acquired and transferred for reconstruction.
- **Manual acquisition mode:** Hold the CS 3500 at a 90-degree angle to the occlusal surface of the teeth, and press the acquisition button on the CS 3500 once to acquire a 3D image. You can rest the tip on the tooth surface during acquisition to steady the CS 3500.

Important Scanning Procedures—Restoration Acquisition

Scanning One Surface at a Time

- Start with the occlusal surface, beginning with the preparation area and then scanning the surrounding teeth on the jaw.
- Next, scan the buccal surface, starting with the preparation area and then scanning the surrounding teeth on the jaw.
- Scan the lingual surface, starting with the preparation area and then scanning the surrounding teeth on the jaw.
- Finally, acquire the buccal bite registration.

Scanning Each Tooth Individually

- Scan one tooth at a time, following the occlusal/buccal/lingual, buccal/occlusal/lingual, or lingual/occlusal/buccal scan path.
- Once the entire tooth is scanned, move to the next tooth and repeat the scanning process until you have scanned all the teeth in the preparation area, as well as the surrounding teeth on the jaw.
- Finally, acquire the buccal bite registration.

Important Scanning Procedures—Orthodontic Acquisition

Scanning One Surface at a Time

- Start with the occlusal surface, and scan the teeth on the jaw.
- Next, scan the buccal surface of the teeth and soft-tissue areas. Be sure to include several millimeters of gingival tissue in the buccal scan.
- Scan the lingual surface of the teeth, palate, and soft-tissue areas. Be sure to include several millimeters of gingival tissue in the lingual scan.
- Finally, acquire the buccal bite registration.

Scanning Each Tooth Individually

- Scan one tooth at a time—following the occlusal/buccal/lingual, buccal/occlusal/lingual, or lingual/occlusal/buccal scan path—and include the palate and soft-tissue areas. Be sure to include several millimeters of gingival tissue in the buccal and lingual scans.
- Once the entire tooth and surrounding areas are scanned, move to the next tooth and repeat the scanning process until you have scanned all the teeth, the palate, and the soft-tissue areas. Be sure to include several millimeters of gingival tissue in the buccal and lingual scans.
- Finally, acquire the buccal bite registration.

Scanning Metal, Crowns, or Holes

- When scanning metal amalgams and crowns, thoroughly dry the teeth before scanning.
- Adjust the operatory light away from the patient to decrease light scatter.
- If holes are displayed in the scanned image near the restoration area, re-scan the area until the holes are filled. Holes away from the restoration area are not of great concern. To better view

any holes, turn off True Color by clicking . Then, place the cursor over the 3D image, left-click and hold to rotate the model. Use the mouse wheel to zoom in on the 3D model.

- If you detect soft-tissue artifacts in an image during acquisition, click on the soft tissue to

highlight that scan, and click . Re-scan the area to fill any holes.

Preparing the CS 3500


The reusable tip attaches to the body of the CS 3500 and provides a sanitary shield for the patient. Always disinfect the body of the CS 3500, and clean and sterilize the tip after each use.

Important: The removable scanner tip is autoclavable up to 20 cycles. After 20 cycles, discard the tip. See the «Technical information» for more information.

To prepare the CS 3500, follow these steps:

- 1 Make sure the lens window at the base of the CS 3500 is clean by wiping it with a moist, lint-free cloth or lens tissue.
- 2 Slide the tip onto the CS 3500 as illustrated, with the lens facing downward (A) for a lower jaw view or upward (B) for an upper jaw view.


WARNING: The base of the CS 3500 becomes heated when the CS 3500 is turned on. Do not touch the heating element at the base of the CS 3500.

Note: Always begin scanning with the regular tip, and scan as much of the mouth as possible. If necessary, change to the small tip to reach areas that are not accessible with the regular tip.


- 3 Press the power button for one second to power on the CS 3500.

- 4 Let the CS 3500 warm up for approximately three minutes to enable the antifog feature on the tip.

Using the Collar to Scan Anterior Teeth

It is recommended that you use the collar when scanning the anterior teeth. Slide the collar onto the CS 3500 tip as illustrated.

Scan the anterior teeth, and remove the collar to scan the remaining teeth. Always disinfect the body of the CS 3500, and clean and sterilize the tip and collar after each use.

Important: The removable scanner tip and collar are autoclavable up to 20 cycles. After 20 cycles, discard the tip and collar. See the «Technical information» for more information.

6 Acquiring a 3D Model for Restoration

You can use the CS 3500 to scan either a partial or full arch. For either type of scan, you should acquire images of the upper jaw, lower jaw, and the buccal bite registration. The software combines these images to create a 3D model.

To acquire a 3D model, follow these procedures:

- Scan the upper and lower jaw.
- Scan the buccal bite registration.
- Refine the image.
- Complete the preparation check and export the image.

Note: In some cases, you can acquire images of a single arch (partial or full) and not obtain a buccal bite registration (for example, if there are no teeth in the opposing arch), but it is recommended that you acquire both arches and a buccal bite registration when possible.

Scanning Teeth on the Upper and Lower Jaw

To acquire images using auto acquisition, follow these steps:

- 1 Dry the teeth thoroughly before starting an acquisition.
- 2 Access the CS 3500 Acquisition interface by clicking  in the imaging software.
- 3 Select  in the Select Acquisition Type window.


Note: To make Acquisition for Restoration the default acquisition type, deselect Always prompt to select acquisition type when the interface is started. To perform an orthodontic acquisition, click the Select Acquisition Type icon in the upper-left corner of the Acquisition window and select Acquisition for Orthodontics.

- 4 On the Acquisition interface, select the Upper Jaw acquisition mode.

OR

On the CS 3500, press the mode button for one second to select the acquisition mode. The mode indicator LED turns green.


i Note: Upper jaw is the default setting until you change it. For information on how to change the default setting, see “Scanner Preferences”.

- 5 Hold the CS 3500 at a 90-degree angle to the occlusal surface of the preparation area. Rest the tip on the tooth surface to steady the CS 3500. Live video is displayed on the video preview screen.

When the CS 3500 is stable and the image is clear, the acquisition begins. Images are automatically acquired and transferred to the 3D model display screen.

The image below shows several acquisitions of the upper jaw.

💡 Tip: Watch the demo at the lower-left of the Acquisition interface for information on how to scan the teeth. Float your mouse over the image of the jaw, and click Show Scanning Sample to launch the demo.

After an acquisition, the feedback indicator displays the following:

-  — Indicates the image was acquired successfully.
-  — Indicates that there is not enough overlap between the last image and the newest image. Reacquire the image, being sure to overlap with the previous tooth by approximately 30%.
-  — Indicates the CS 3500 was not stable. Reacquire the image.

Important: When scanning a tooth, keep approximately a 30% overlap with the last image.

- 6 Slowly move the CS 3500 tip along the occlusal surface to scan the remaining teeth in the preparation area, keeping approximately a 30% overlap with the last image.


Tip: It is recommended that you use the collar when scanning anterior teeth.

- 7 When the occlusal surface scan is complete, scan the lingual surface of the teeth in the preparation area. Hold the CS 3500 at a 45-degree angle to the lingual surface of the teeth.
- 8 When the lingual surface scan is complete, scan the buccal surface of the preparation area. Hold the CS 3500 at a 45-degree angle to the buccal surface of the teeth.

Important: Re-dry the teeth as appropriate throughout the acquisition process.

The example below shows a 3D model of the upper jaw when the occlusal, lingual, and buccal surfaces have been completely scanned.

Important: If holes are displayed in the scanned image on the preparation area, re-scan the area until the holes are filled. Use the mouse wheel to zoom in on the preparation area for a closer look.

- 9 Once the upper jaw has been scanned, you can begin scanning the lower jaw. Click  and repeat steps 5 through 8 until the teeth in the area of restoration on the lower jaw are scanned.

The example below shows a 3D model of the lower jaw when the occlusal, lingual, and buccal surfaces have been completely scanned.

- 10 Visually inspect the 3D model for any holes. If holes are present near the restoration area, re-scan the area until the holes are filled.
- 11 When you are satisfied with the 3D model, continue with the buccal bite registration.

Important: If you acquire images of teeth from both the upper and lower jaw, you must also acquire the buccal bite registration.


Scanning the Buccal Bite Registration

To acquire a buccal bite registration using auto acquisition, follow these steps:

- 1 Click .
- 2 Have the patient bite down.
- 3 Position the CS 3500 at a 90-degree angle to the buccal surface of the preparation area, and align the point where the upper and lower teeth meet in the middle of the video preview screen. Rest the tip on the tooth surface to help steady the CS 3500.


When the CS 3500 is stable and the image is clear, the acquisition begins.

The example below shows a buccal bite registration.

A dot  is displayed at the bottom of the window to indicate the capture was successful. A successful bite image includes both the upper and lower arch.

Half of a dot (only one arch, upper or lower) indicates a partial buccal bite. The example below shows that one complete buccal bite image and one-half of a buccal bite image have been acquired. Note the green box on the teeth on the lower jaw. To complete the bite, acquire the teeth on the opposing jaw, directly opposite of the green box.

To complete a partial buccal bite image, move the CS 3500 slightly more gingival, *toward the missing arch*, and ensure that there is at least a 30% overlap between the second buccal bite acquisition and the previous one.

Tip: If you need to acquire a second buccal bite image because an opposing arch is missing, always move toward the opposing arch with the CS 3500.

Once both arches are displayed in the Acquisition interface, the bite image is successful, and you can acquire additional bite images, which will increase the accuracy of the bite.

Important: For a partial arch, you should acquire at least three buccal bite images—one on the preparation, one mesial to the preparation, and one distal to the preparation. For a full arch, you should acquire at least four buccal bite images—one on each side of the mouth at the molars, and one on each side of the mouth at the canines.

- 4 You should acquire at least two more buccal bite images, one on either side of the preparation area. You can take up to six bite images for a full arch.
- 5 Once the bite has been registered, rotate the model and zoom the view to ensure that the bite is accurate and that there are no areas where the bite is mismatched. Click a dot to view the bite for that acquisition.

The example below shows several buccal bite registrations.

- 6 When you have finished acquiring the buccal bite registration, proceed to the refinement selection process.

Refinement

The refinement step enables you to select specific teeth in the 3D model display area and automatically refine them so that these teeth have the highest accuracy.


Important: If you have acquired fewer than five teeth, the software skips refinement selection, and the entire image is refined. You would then continue with Step 5.

To select specific teeth on the 3D model for refinement, follow these steps:

- 1 Click  to perform refinement selection.
- 2 Click .

- 3 Place the cursor over the area you want to refine, and left-click and drag the mouse to select the teeth for refinement. You should include the preparation, the adjacent surfaces of neighboring teeth, the teeth on the opposite arch, and at least one buccal bite registration acquisition in this selection.


Important: Use the mouse to angle the image, and select teeth from above -- not from the side. This prevents teeth behind the target teeth from being refined.


Important: You cannot select more than half an arch, or you will not be able to continue with the refinement process.

Release the mouse button. The selected teeth are highlighted in blue.

Tip: If necessary, you can select additional areas for refinement. Press Ctrl, place the cursor over the additional area to refine, and left-click and drag the mouse to select the additional teeth.

- 4 Click  to refine the image.

- 5 Manipulate the refined 3D model using the following methods:
- Right-click and hold on the 3D model to move it in the window.
 - Left-click and hold on the 3D model to rotate it.
 - If your mouse has a scroll wheel, use the wheel to zoom in or zoom out on the 3D model.
 - Click  to scale the model to its best view.

- Click  or  to hide the upper or lower jaw. Click the button again to restore the view of the jaw.
- Click  or  to select and delete excess soft tissue in the image.

- If you find holes in the preparation area, click  and re-scan the teeth that are missing data.
- Repeat steps 1 through 6 until you are satisfied with the 3D model.
- When you have successfully finished the refinement process, proceed to the preparation check process.

Preparation Check

The preparation check step enables you to verify that the 3D model is satisfactory before the 3D model file is sent for the final restoration process.

When performing the preparation check, you must examine the model for each of the following:

- Occlusal space
- Undercut
- Margin line
- Identify the appropriate arch
- Identify the appropriate teeth in the preparation area

To complete the preparation check, follow these steps:

- Examine the occlusal space on the 3D model and do one of the following:
 - If the occlusal space is acceptable, click the corresponding option in the Preparation Check section.
 - If the occlusal space is not acceptable, click  and re-scan the affected teeth.
- Examine the undercut on the 3D model and do one of the following:
 - If the undercut is acceptable, click the corresponding option in the Preparation Check section.
 - If the undercut is not acceptable, make the necessary adjustments to the preparation and then re-scan the preparation area and surrounding teeth.
- Examine the margin line on the 3D model and do one of the following:
 - If the margin line is acceptable, click the corresponding option in the Preparation Check section.
 - If the margin line is not acceptable, click  and re-scan the affected teeth.

- 4 Click . The arch diagram is displayed.

- 5 Click the Adult or Child icon to display the appropriate arch.
- 6 Select the teeth on the arch that were acquired for the 3D model, including those on the upper and lower jaw.
- 7 If you are exporting the 3D model to a dental lab, you can draw a margin line reference on a snapshot of the 3D model. See "Drawing Margin Lines".
- 8 When satisfied with the 3D model, click  to send the 3D model to the imaging software.


Important: If your office is sending files to a dental lab, send the dental lab a copy of the "Installing and Using the Orientation Converter" document and the setupCSOrientationConverter.zip file. Both files are found in the Orientation Converter folder, located in the Documentation folder you downloaded during the installation of the CS 3500.

Drawing Margin Lines

If you are exporting the 3D model to a dental lab, you can draw a margin line reference on a snapshot of the 3D model that can be exported with the model.


To draw a margin line, follow these steps:



- 1 Click . The Margin Line window containing the 3D model is displayed.
- 2 Left-click and drag to draw the margin line around the preparation.

- 3 When you are finished, click Save. You can draw up to five margin line references, which are saved in JPEG format.

- 4 To review the margin line references, click  to open the Margin Line Viewer window. The images below show several margin line references for the preparation area.

Acquiring 2D Images

You can use the CS 3500 to acquire 2D intraoral images that you can send to the dental lab along with the 3D model. You can acquire a maximum of 15 2D images.

Note: The Intraoral Mode button is available only during acquisition. If you need to acquire 2D images after refinement, you must repeat the steps for acquiring a 3D model.

To acquire 2D images, follow these steps:

- 1 Prepare the CS 3500 for use as described under “Preparing the CS 3500”.
- 2 Access the CS 3500 Acquisition interface by clicking  in the imaging software.
- 3 In the Acquisition interface, click .
- 4 Hold the CS 3500 over the area you want to view. Live video is displayed in the window.
- 5 Press the acquisition button to acquire an image.
- 6 When you are finished acquiring 2D images, click  to return to the Acquisition interface.

7

Acquiring a 3D Model for Orthodontics

You can use the CS 3500 to scan a full arch to create a 3D model. You then export this model to the CS Model software, where you can use the electronic model to take basic distance measurements, eliminating the need for taking impressions and working with or storing stone models. You should acquire images of the upper jaw, lower jaw, and the buccal bite registration. The software combines these images to create the 3D model.

To acquire a 3D model for orthodontic use, follow these procedures:

- Scan the upper and lower jaw.
- Scan the buccal bite registration.
- Refine the image.
- Complete the preparation check and export the image.

Note: In some cases, you can acquire images of a single arch (partial or full) and not obtain a buccal bite registration (for example, if there are no teeth in the opposing arch), but it is recommended that you acquire both arches and a buccal bite registration when possible.

Scanning Teeth on the Upper and Lower Jaw

To acquire images using auto acquisition, follow these steps:

- 1 Dry the teeth thoroughly before starting an acquisition.
- 2 Access the CS 3500 Acquisition interface by clicking  in the imaging software.
- 3 Select  in the Select Acquisition Type window.


Note: To make Acquisition for Orthodontics the default acquisition type, deselect Always prompt to select acquisition type when the interface is started. To perform a restoration acquisition, click the Select Acquisition Type icon in the upper-left corner of the Acquisition window and select Acquisition for Restoration.

- 4 On the Acquisition interface, select the Upper Jaw acquisition mode.

OR

On the CS 3500, press the mode button for one second to select the acquisition mode. The mode indicator LED turns green.


Note: Upper jaw is the default setting until you change it. For information on how to change the default setting, see "Scanner Preferences".

- 5 Hold the CS 3500 at a 90-degree angle to the occlusal surface of the teeth. Rest the tip on the tooth surface to steady the CS 3500. Live video is displayed on the video preview screen.

When the CS 3500 is stable and the image is clear, the acquisition begins. Images are automatically acquired and transferred to the 3D model display screen.

The image below shows several acquisitions of the upper jaw.

Tip: Watch the demo at the lower-left of the Acquisition interface for information on how to scan the teeth. Float your mouse over the image of the jaw, and click Show Scanning Sample to launch the demo.

After an acquisition, the feedback indicator displays the following:

-  — Indicates the image was acquired successfully.
-  — Indicates that there is not enough overlap between the last image and the newest image. Reacquire the image, being sure to overlap with the previous tooth by approximately 30%.
-  — Indicates the CS 3500 was not stable. Reacquire the image.

Important: When scanning a tooth, keep approximately a 30% overlap with the last image.

- 6 Slowly move the CS 3500 tip along the occlusal surface to scan the remaining teeth, keeping approximately a 30% overlap with the last image.


Tip: It is recommended that you use the collar when scanning anterior teeth.

- 7 When the occlusal surface scan is complete, scan the lingual surface of the teeth. Hold the CS 3500 at a 45-degree angle to the lingual surface of the teeth. Be sure to include several millimeters of gingival tissue in the lingual scan.
- 8 When the lingual surface scan is complete, scan the buccal surface of the teeth. Hold the CS 3500 at a 45-degree angle to the buccal surface of the teeth. Be sure to include several millimeters of gingival tissue in the buccal scan.

Important: Re-dry the teeth as appropriate throughout the acquisition process.


The example below shows a 3D model of the upper jaw when the occlusal, lingual, and buccal surfaces have been completely scanned, as well as much of the palate.

Important: If holes are displayed in the scanned image, re-scan the area until the holes are filled. Use the mouse wheel to zoom in on the image for a closer look.

- 9 Once the upper jaw has been scanned, you can begin scanning the lower jaw. Click  and repeat steps 5 through 8 until the teeth on the lower jaw are scanned.

The example below shows a 3D model of the lower jaw when the occlusal, lingual, and buccal surfaces have been completely scanned, including soft tissue.

- 10 Visually inspect the 3D model for any holes. If holes are present, re-scan the area until the holes are filled.
- 11 When you are satisfied with the 3D model, continue with the buccal bite registration acquisition.

Scanning the Buccal Bite Registration

To acquire a buccal bite registration using auto acquisition, follow these steps:

- 1 Click .
- 2 Have the patient bite down.
- 3 Position the CS 3500 at a 90-degree angle to the buccal surface of the teeth, and align the point where the upper and lower teeth meet in the middle of the video preview screen. Rest the tip on the tooth surface to help steady the CS 3500.


When the CS 3500 is stable and the image is clear, the acquisition begins.

The example below shows a buccal bite registration.

A dot  is displayed at the bottom of the window to indicate the capture was successful. A successful bite image includes both the upper and lower arch.

Half of a dot (only one arch, upper or lower) indicates a partial buccal bite. The example below shows that one complete buccal bite image and one-half of a buccal bite image have been acquired. Note the green box on the teeth on the lower jaw. To complete the bite, acquire the teeth on the opposing jaw, directly opposite of the green box.

To complete a partial buccal bite image, move the CS 3500 slightly more gingival, *toward the missing arch*, and ensure that there is at least a 30% overlap between the second buccal bite acquisition and the previous one.

Tip: If you need to acquire a second buccal bite image because an opposing arch is missing, always move toward the opposing arch with the CS 3500.


Once both arches are displayed in the Acquisition interface, the bite image is successful, and you can acquire additional bite images, which will increase the accuracy of the bite.

Important: For a full arch, you should acquire at least four buccal bite images—one on each side of the mouth at the molars, and one on each side of the mouth at the canines.

- 4 Acquire several additional buccal bite images. You can take up to six bite images for a full arch.
- 5 Once the bite has been registered, rotate the model and zoom the view to ensure that the bite is accurate and that there are no areas where the bite is mismatched. Click a dot to view the bite for that acquisition.

The example below shows several buccal bite registrations.


- 6 When you have finished acquiring the buccal bite registration, proceed to the refinement step.

Refinement

The refinement step enables you to further process the 3D image to obtain the highest accuracy.

To refine the 3D model and further examine it, follow these steps:

- 1 Click  to refine the image.

2 Manipulate the refined 3D model using the following methods:

- Right-click and hold on the 3D model to move it in the window.
- Left-click and hold on the 3D model to rotate it.
- If your mouse has a scroll wheel, use the wheel to zoom in or zoom out on the 3D model.
- Click  to scale the model to its best view.
- Click  or  to hide the upper or lower jaw. Click the button again to restore the view of the jaw.

- Click  or  to select and delete excess soft tissue in the image.
- 3 If you find holes in the image, click  and re-scan the teeth that are missing data.
 - 4 Repeat steps 1 through 3 until you are satisfied with the 3D model.
 - 5 When you have successfully finished the refinement process, proceed to the preparation check process.

Preparation Check

The preparation check step enables you to verify that the 3D model is satisfactory before you export the 3D model file to the imaging software.

When performing the preparation check, you must examine the model for each of the following:

- Holes in the image
- Jagged boundaries
- Irregular soft tissue

To complete the preparation check, follow these steps:

- 1 Examine the 3D model for holes and do one of the following:
 - If there are no holes in the image, click the corresponding option in the Preparation Check section.
 - If there are holes that need to be filled, click  and re-scan the affected teeth.
- 2 Examine the 3D model for jagged boundaries and do one of the following:
 - If there are no jagged boundaries in the image, click the corresponding option in the Preparation Check section.
 - If you find jagged boundaries in the image, remove the boundaries using the  or  tool.
- 3 Examine the 3D model for irregular soft tissue and do one of the following:
 - If there is no irregular soft tissue in the image, click the corresponding option in the Preparation Check section.
 - If you find irregular soft tissue, remove the irregular soft tissue using the  or  tool.
- 4 When satisfied with the 3D model, click  to send the 3D model to the imaging software.

Important: If your office is sending files to a dental lab, send the dental lab a copy of the "Installing and Using the Orientation Converter" document and the setupCSOrientationConverter.zip file. Both files are found in the Orientation Converter folder, located in the Documentation folder you downloaded during the installation of the CS 3500.

Acquiring 2D Images

You can use the CS 3500 to acquire 2D intraoral images that you can send to the dental lab along with the 3D model. You can acquire a maximum of 15 2D images.

Note: The Intraoral Mode button is available only during acquisition. If you need to acquire 2D images after refinement, you must repeat the steps for acquiring a 3D model.

To acquire 2D images, follow these steps:

- 1 Prepare the CS 3500 for use as described under "Preparing the CS 3500".
- 2 Access the CS 3500 Acquisition interface by clicking  in the imaging software.
- 3 In the Acquisition interface, click .
- 4 Hold the CS 3500 over the area you want to view. Live video is displayed in the window.
- 5 Press the acquisition button to acquire an image.
- 6 When you are finished acquiring 2D images, click  to return to the Acquisition interface.

8 Maintenance

The user is responsible for the operation and maintenance of the device. Only qualified professionals can manage data device. Before working with the device, must necessarily be trained. At need a qualified service technician for inspection and maintenance.

Warranty: 1 year

During the warranty period the manufacturer provides warranty repair or replacement of defective product or of its individual element (component) free of charge.

Cleaning, Disinfecting, and Sterilizing

You must clean, disinfect, and sterilize the CS 3500 and accessories regularly.

The removable CS 3500 scanner tips and collars are autoclavable up to 20 cycles. After 20 cycles, discard the tip and collar.

Important: For information on cleaning, disinfecting, and sterilizing, see the Technical information.

9 Contact Information

Manufacturer's Address

Carestream Health, Inc.
150 Verona Street
Rochester, NY 14608, USA

Authorized Representatives

Authorized Representative in the European Community

Carestream Health France
1, rue Galilée
93192 Noisy-Le-Grand Cedex, France

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Phone: 8(495) 660-56-90
Fax: 8 (495) 660-56-91

A Division of Carestream Health, Inc.
150 Verona St.
Rochester, NY 14608
USA

For more information, visit: www.carestreamdental.com

To give documentation feedback, visit: www.carestreamdental.com/documentationfeedback

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