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NIALL SHEEHY

NOTARY PUBLIC


6/7 Glentworth Street,
Limerick.
Ireland.

Telephone + 353 61 414355
Facsimile + 353 61 414738
Email nsheehy@csg.ie

TO ALL TO WHOM these Presents shall come

I, **NIALL SHEEHY**, Notary Public of 6/7 Glentworth Street, Limerick duly appointed by the Chief Justice of Ireland to practise for life as a Notary Public in the County Borough of Limerick, in the Counties of Limerick and Clare and in the adjoining counties, do hereby certify that on the 1st of August 2017 **EIMEAR CARR**, known to me as Regulatory Affairs Specialist at Beckman Coulter Ireland Inc. of Lismeehan, O'Callaghan's Mills, Co. Clare, Ireland, acknowledged to me that on the 10th of August 2017 she did validly duly sign on behalf of Beckman Coulter Ireland Inc. the attached **Declaration of Product Instructions for Use** for Registration in the Russian Federation by signing her name thereto as the same now appears, and I certify that the signature subscribed thereto is the true and proper signature of the said **EIMEAR CARR**.

IN FAITH AND TESTIMONY whereby I have set my hand and affixed my seal of Office this 11th day of August 2017.

Niall Sheehy
Notary Public

**BECKMAN
COULTER**

Beckman Coulter Ireland Inc.,
Lismeehan,
O' Callaghan's Mills,
Co. Clare,
Ireland.
Tel: +353-65-6831100.
Fax: +353-65-6831122.

Declaration of Product Instructions for Use

To Whom it My Concern

I, Eimear Carr, as duly authorized Regulatory Affairs Specialist at, on behalf of Beckman Coulter Ireland Inc., to the best of my knowledge hereby declare that the attached document is a true copy of the Russian Instructions for Use (IFU) which are included with this registration dossier to provide supporting information for product registration purposes for the below listed products.

Product Name

Urine/CSF Albumin calibrator

Catalogue No.

B38859

Dated this 10th day of August 2017

Signed on behalf of Beckman Coulter Ireland Inc.

Eimear Carr

Eimear Carr

Regulatory Affairs Specialist

Beckman Coulter Ireland Inc.

Registered office: Lismeehan,
O' Callaghan's Mills,
Co. Clare.

Directors: Piers Devereux
Paul Henry

Incorporated In Republic of Ireland
with limited liability
Co. Reg. No. 501471.

AU

Instructions For Use

© 2017 Beckman Coulter, Inc. All rights reserved.

UALB CAL
Urine/CSF Albumin Calibrator

REF B38859 Low 1 x 2 mL Calibrator # 1 White Cap
1 x 2 mL Calibrator # 2 Yellow Cap
1 x 2 mL Calibrator # 3 Orange Cap
1 x 2 mL Calibrator # 4 Red Cap
High 1 x 2 mL Calibrator # 5 Black Cap

For *in vitro* diagnostic use only.

ANNUAL REVIEW

Reviewed by	Date	Reviewed by	Date

PRINCIPLE

INTENDED USE

The Urine/CSF Albumin calibrator is intended to be used to calibrate the Urine/CSF Albumin reagent on the Beckman Coulter AU clinical chemistry systems.

REAGENTS

CONTENTS

The Urine/CSF Albumin Calibrator consists of known quantities of human albumin in a buffered solution.

WARNING AND PRECAUTIONS

WARNING - POTENTIAL BIOHAZARDOUS MATERIAL

The calibrator is manufactured from human material; therefore it should be handled as though capable of transmitting infectious disease. Each serum or plasma donor unit used in the preparation of this material was tested by United States Food and Drug Administration (FDA) approved methods and found to be negative for the antibodies to HIV and HCV and nonreactive for HBsAg. Because no test method can offer complete assurance that HIV, hepatitis B virus, and hepatitis C virus or other infectious agents are absent, this material and all patient samples should be handled as though capable of transmitting infectious disease. This product may also contain other human source material for which there is no approved test. The FDA recommends such samples be handled as specified in the Centers for Disease Control's Biosafety Level 2 guidelines¹.

REACTIVE INGREDIENTS

Final concentration of reactive ingredients.

Urine/CSF Albumin Calibrator Levels 1 - 5

Human Albumin	Various
Sodium Azide (used as a preservative)	< 0.1% (w/w)

CAUTION

Sodium azide preservative may form explosive compounds in metal drain lines. See NIOSH Bulletin: Explosive Azide Hazard (8/16/76). To avoid the possible build-up of azide compounds, flush wastepipes with water after the disposal of undiluted reagent. Sodium azide disposal must be in accordance with appropriate local regulations.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

Safety Data Sheet is available at techdocs.beckmancoulter.com

INDICATIONS OF DETERIORATION

Visible signs of microbial growth, turbidity, precipitate or any change in color in the Urine/CSF Albumin calibrator may indicate degradation and warrant discontinuation of use.

CALIBRATION

CALIBRATOR PREPARATION

The calibrators are ready for use.

No preparation is required.

Recommended Procedure

1. Allow contents of the bottles to come to room temperature.
2. Mix the contents by gently inverting the vials several times to achieve a homogenous mixture.
3. Transfer sufficient volume of the calibrator to sample cups. Handle this calibrator with the same care used for patient samples.
4. Replace the cap immediately and store unused calibrator at 2... 8°C. Do not mix caps among the various bottles in the calibrator set.
5. For detailed calibration instructions, refer to the AU Instrument User's Guide.
6. Good Quality Control Practices should be observed to assure proper system performance.

CALIBRATOR STORAGE AND STABILITY

1. The calibrator temperature requirements should be met during transport for all transport methods.
2. Calibrators should be stored in the original package at +2°C to +8°C in a temperature controlled area e.g. Refrigerator.
3. Do Not Freeze

- Opened bottles of calibrators are stable for 30 days, provided they are free from contamination, tightly capped immediately after each use and stored at 2...8°C.
- Calibrators shelf life is 18 months from date of manufacture. The expiration date and date of manufacture is indicated on the product packaging.
- The calibrator is stable, unopened, up to the stated expiry date when stored at 2...8°C.
- Do not use the calibrators after the expiry date on the label.

CALIBRATOR LIMITATIONS

This calibrator can only be used with the application of Urine/CSF Albumin reagent B38858/B46435.

The results obtained with this calibrator are dependent on several factors, including proper storage of the calibrator, and proper technique and use of Beckman Coulter AU Analyzers and their respective reagents.

The albumin result of a urine sample may be elevated when it immediately follows a serum sample. In order to eliminate this effect, it is recommended to calibrate the Urine/CSF albumin assay separately to serum assays.

TRACEABILITY

The Urine/CSF Albumin calibrator values are traceable to the International Federation of Clinical Chemistry Certified Reference Material CRM470.

Value Assignment

Assigned values are listed in the Value Assignment Sheet provided in Calibrator pack.

Calibration values have been determined using standardised procedures and Beckman Coulter System Reagents on Beckman Coulter analysers. Please ensure that the lot number on the calibrator vial is the same as the one listed in the table of assigned values in the package insert.

Note: Following calibration, the resulting curve should be visually reviewed, on the Beckman Coulter analyser, for acceptability using the software options - Routine, Calibration Monitor, Calibration Curve. Quality control procedures should be undertaken immediately following calibration in accordance with good laboratory practice.

ADDITIONAL INFORMATION

The manufacturer guarantees that the product meets the requirements of regulatory and technical documents. Safety and quality of the product are guaranteed until the expiry date. The Manufacturer is responsible for product defects, excluding defects that appeared because of incorrect usage, shipment and storage, or third persons or compelling force. The manufacturer assures that they will substitute the product at their own expense, if technical and functional features (consumer properties) don't meet the requirements of regulatory and technical documents, if stated defects were the subsequence of a hidden defect of materials or low-quality production by the manufacturer. If you have questions concerning quality of the product please address the authorized representative of the manufacturer in the Russian Federation.

The lot number on the vial is the same as the one listed in the table on the value assign sheet.

The selected value is appropriate for the units on the analyzer parameter settings.

REVISION HISTORY

Russian market specific IFU

REFERENCES

1. CDC-NIH manual, Biosafety in Microbiological and Biomedical Laboratories, 5th Edition (CDC 21-1112), U.S. Government Printing Office, Washington, D.C. (2009).

Beckman Coulter Ireland Inc., Lismeehan, O'Callaghan's Mill's Co. Clare, Ireland (001) 703-527-3887
www.beckmancoulter.com

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