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NIALL SHEEHY

NOTARY PUBLIC

6/7 Glentworth Street,
Limerick,
Ireland.

Telephone + 353 61 414355
Facsimile + 353 61 414738
Email nsheehy@csg.ie

TO ALL TO WHOM these Presents shall come

I, **NIALL SHEEHY**, Notary Public of 6/7 Glentworth Street, Limerick duly appointed by the Chief Justice of Ireland to practise for life as a Notary Public in the County Borough of Limerick, in the Counties of Limerick and Clare and in the adjoining counties, do hereby certify that on the 11th of August 2017 **EIMEAR CARR**, known to me as Regulatory Affairs Specialist at Beckman Coulter Ireland Inc. of Lismeehan, O'Callaghan's Mills, Co. Clare, Ireland, acknowledged to me that on the 11th of August 2017 she did validly duly sign on behalf of Beckman Coulter Ireland Inc. the attached **Declaration of Product Instructions for Use** for Registration in the Russian Federation by signing her name thereto as the same now appears, and I certify that the signature subscribed thereto is the true and proper signature of the said **EIMEAR CARR**.

IN FAITH AND TESTIMONY whereby I have set my hand and affixed my seal of Office this 11th day of August 2017.

Niall Sheehy
Notary Public
NIALL SHEEHY
6/7 Glentworth Street

Beckman Coulter Ireland Inc.
Lismeehan,
O' Callaghan's Mills,
Co. Clare,
Ireland
Tel: +353-65-6831100.
Fax: +353-65-6831122.

Declaration of Product Instructions for Use

To Whom it My Concern,

I, Eimear Carr, as duly authorized Regulatory Affairs Specialist at, on behalf of Beckman Coulter Ireland Inc., to the best of my knowledge hereby declare that the attached document is a true copy of the Russian Instructions for Use (IFU) which are included with this registration dossier to provide supporting information for product registration purposes for the below listed products.

Product Name

Urine/CSF Albumin

Urine/CSF Albumin

Catalogue No.

B38858

B46435

Dated this 11th day of August 2017

Signed on behalf of Beckman Coulter Ireland Inc.

Eimear Carr
Eimear Carr

Eimear Carr

Regulatory Affairs Specialist

Beckman Coulter Ireland Inc.

Niall Sheehy
NIALL SHEEHY
11-08-2017

11-08-2017

Registered office: Lismeehan,
O' Callaghan's Mills,
Co. Clare.

Directors: Piers Devereux
Paul Henry

Incorporated in Republic of Ireland
with limited liability
Co. Reg. No. 501471.

AU

Instructions For Use

© 2017 Beckman Coulter, Inc. All rights reserved.

UALB

URINE CSF ALBUMIN

REF

B38858 4 x 32.6 mL R1, 4 x 4.4 mL R2
B46435 4 x 54.8 mL R1, 4 x 7.3 mL R2

For *in vitro* diagnostic use only.

ANNUAL REVIEW

Reviewed by	Date	Reviewed by	Date

PRINCIPLE

INTENDED USE

The Urine/CSF Albumin reagent is intended to be used for the quantitation of albumin concentration in human urine and cerebrospinal fluid (CSF) on the Beckman Coulter AU clinical chemistry systems as an aid in the diagnosis of kidney diseases.

SUMMARY AND EXPLANATION

Reference^{1,2,3,4,5}

The earliest clinical evidence of nephropathy is the appearance of low, but abnormal levels (> 30 mg/day or $20 \mu\text{g}/\text{min}$) of albumin in the urine, and patients with moderately increased urine albumin are referred to as having incipient nephropathy. Conventional qualitative tests (chemical strips or dipsticks) for albuminuria do not detect the small increases in urinary albumin excretion seen in early stages of nephropathy. For this purpose, tests for moderately increased urine albumin are used. Moderately increased urine albumin is defined as an albumin excretion rate between $30\text{--}300$ mg/24h on 2 of 3 urine collections.

Increased urinary albumin loss is considered a clinically important indicator of deteriorating renal function in diabetic subjects and regular screening of urinary albumin loss is valuable in monitoring type 1 and type 2 diabetes. Prospective studies have demonstrated that increased urinary albumin excretion precedes and is highly predictive of diabetic nephropathy, end stage renal disease, cardiovascular mortality, and total mortality in patients with diabetes mellitus. In addition, increased urinary albumin excretion identifies a group of non-diabetic subjects at increased risk for coronary artery disease.

The degree of permeability of the blood-CSF barrier may be evaluated by the simultaneous measurement of serum and CSF albumin. CSF albumin measurements may also be used to determine the IgG CSF/albumin CSF ratio, which is an important factor in differentiating between intrathecal and localized synthesis of IgG.

Instructions for Use BLOS836058AU 08
AUGUST 2017

English

UALB Specific Protein Special Chemistry
Page 1 of 10

METHODOLOGY

Urine/CSF Albumin reagent is used to measure albumin concentration by a turbidimetric method. In the reaction, anti-human serum albumin antibodies combine with albumin from the sample to form immune complexes that scatter light in proportion to their size, shape and concentration. The absorbance of these aggregates is proportional to the albumin concentration in the sample. Change in absorbance is measured at 380nm with subtraction of a reference wavelength at 800nm.

SPECIMEN

TYPE OF SPECIMEN

Urine and CSF may be used. The urine specimen should be a fresh or a 24 hour urine.

Urine and CSF may be diluted 1 in 10 on-board the analyzer with either 0.9% saline or deionized water using automatic dilution function. Refer to the Beckman Coulter AU analyzer User Guide/Instructions For Use (IFU) for detailed instructions.

SPECIMEN STORAGE AND STABILITY

Urine^{6,7}: Stable for 1 month when stored refrigerated (2...8°C). Stable for 6 months when stored frozen at -20°C.

Cerebrospinal fluid^{6,7}: Stable for 72 hours when stored refrigerated (2...8°C). Stable for 6 months when stored frozen at -20°C.

REAGENTS

CONTENTS

Urine/CSF Albumin reagent

WARNING AND PRECAUTIONS

Exercise the normal precautions required for handling all laboratory reagents.

Dispose of all waste material in accordance with local guidelines.

Antisera was produced in healthy animals in facilities free from rinderpest, foot and mouth disease, peste des petits ruminants, Rift Valley fever, bovine spongiform encephalopathy and blue tongue disease.

Do not ingest. Harmful if swallowed.

Sodium azide preservative may form explosive compounds in metal drain lines. See NIOSH Bulletin 13: Explosive Azide Hazard (8/16/76)⁸. To avoid the possible build-up of azide compounds, flush wastepipes with water after the disposal of undiluted reagent.

This product contains material of animal origin. The product should be considered as potentially capable of transmitting infectious diseases.

REACTIVE INGREDIENTS

Final concentration of reactive ingredients:

Phosphate buffer	18 mmol/L
Goat Anti-human Albumin Antibody	Variable
Polyethylene glycol 8000	3.6%
Sodium Azide (used as a preservative)	< 0.1% (w/w) ⁹

The concentrations of the reactive components of the reagents shown on the kit label are the actual concentrations in the individual R1/R2 vials. The reagent composition which is shown in the Instructions For Use is the final concentration of these components in the reaction cuvette after addition of R1, Sample, and R2.

 CAUTION

Sodium azide preservative may form explosive compounds in metal drain lines. See NIOSH Bulletin: Explosive Azide Hazard (8/16/76). To avoid the possible build-up of azide compounds, flush wastepipes with water after the disposal of undiluted reagent. Sodium azide disposal must be in accordance with appropriate local regulations.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

	Safety Data Sheet is available at techdocs.beckmancoulter.com
---	--

MATERIALS NEEDED BUT NOT SUPPLIED WITH REAGENT KIT

Urine/CSF Albumin Calibrator (B38859)

Deionised water for performing reagent blank.

EQUIPMENT AND MATERIALS

For AU400/400[®]/480, AU640/640[®]/680, AU2700/5400/AU5800 and DxC 700 AU Beckman Coulter Analyzers.

REAGENT PREPARATION

The reagents are ready for use and can be placed directly on board the instrument.

REAGENT STORAGE AND STABILITY

Shipping and Storage should be carried out in the original package at 2...8°C

Do Not Freeze

The expiry date is printed on the package label.

The reagents are stable, unopened, up to the stated expiry date when stored at 2...8°C.

Opened bottles of reagent are stable for 60 days when stored in the refrigerated compartment of the analyzer.

Do not use the reagents after the expiry date on the label.

The reagent temperature requirements should be met during transport for all transport methods.

Reagents should be stored in the original package at +2°C to +8°C in a temperature controlled area e.g. Refrigerator.

Reagents shelf life is 18 months from date of manufacture. The expiration date and date of manufacture is indicated on the product packaging.

INDICATIONS OF DETERIORATION

Visible signs of microbial growth, turbidity or precipitate, or any change in reagent color may indicate degradation and warrant discontinuance of use.

CALIBRATION

CALIBRATION INFORMATION

Urine/CSF Albumin Calibrator Reference No. B38859.

The frequency of calibration for this procedure is every 60 days.

Recalibration of this test is required when any of these conditions exist:

Change in reagent lot or significant shift in control values;

Major preventative maintenance was performed on the analyser or a critical part was replaced.

Following calibration, the resulting curve should be visually reviewed, on the Beckman Coulter analyser, for acceptability using the software options - Routine, Calibration Monitor, Calibration Curve. Quality control procedures should be undertaken immediately following calibration in accordance with good laboratory practice.

TRACEABILITY

The Urine/CSF Albumin calibrator values are traceable to the International Federation of Clinical Chemistry Certified Reference Material CRM470.

QUALITY CONTROL

During operation of the Beckman Coulter AU analyzer at least two levels of an appropriate quality control material should be tested a minimum of once a day.

Only control materials, of human origin, with values determined by this Beckman Coulter system should be used.

Each laboratory should establish its own control frequency however good laboratory practice suggests that controls be tested each day patient samples are tested and each time calibration/blanking is performed.

Quality Controls should be used in accordance with the applicable local guidelines.

The results obtained by any individual laboratory may vary from the given mean value. It is therefore recommended that each laboratory generates analyte specific control target values and intervals based on multiple runs according to their requirements. These target values should fall within the corresponding acceptable ranges given in the relevant product literature.

If any trends or sudden shifts in values are detected, review all operating parameters.

Each laboratory should establish guidelines for corrective action to be taken if controls do not recover within the specified limits.

Please note that the recovery of some controls may vary with reagent lots due to the use of non-human materials in the controls.

TESTING PROCEDURE(S)

Refer to the appropriate Beckman Coulter AU analyser User Guide/Instructions For Use (IFU) for analyser-specific assay instructions for the sample type as listed in the Intended Use statement.

CALCULATIONS

The Beckman Coulter AU analyzers automatically compute every determination at the same time interval.

RESULTS INTERPRETATION

Automatically generated for each sample in mg/L (mg/dL).

REPORTING RESULTS

REFERENCE INTERVALS

American Diabetes Association Definition of Moderately increased urine albumin ¹			
Urine:	24-h collection (mg/24h)	Timed collection (µg/min)	Spot collection (µg/mg creatinine)
Normal	< 30	< 20	< 30
Moderately increased	30 - 299	20 - 199	30 - 299
Clinical albuminuria	≥ 300	≥ 200	≥ 300
CSF ⁹ :			
	3 mo – 4y	0 - 450 mg/L (0 - 45 mg/dL)	
	> 4 y	100 - 300 mg/L (10 - 30 mg/dL)	

Reference Intervals shown above were taken from the literature.

Expected values may vary with age, sex, sample type, diet and geographical location. Each laboratory should verify the transferability of the expected values to its own population, and if necessary determine its own reference interval according to good laboratory practice. For diagnostic purposes, results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

PROCEDURAL NOTES

LIMITATIONS

1. Samples with Albumin concentrations > 20,000 mg/L (2,000 mg/dL) can generate false low results without appropriate "Z" flags due to excess antigen in the sample.
2. The albumin result of a urine sample may be elevated when it immediately follows a serum sample. In order to eliminate this effect it is recommended to: Calibrate the Urine/CSF albumin assay separately to serum assays. Avoid requisition of a urine sample following a serum sample. For the AU2700, AU5400, AU480, AU680, AU5800 and DxC 700 AU Analyzers sample contamination parameters are available on the Beckman Coulter website. If sample contamination parameters are not programmed, place a sample cup containing 2% AU wash solution prior to urine samples and requisition a test for this sample.
3. Samples with extremely abnormal optical characteristics, including turbidity, interfere with test results. Extremely turbid samples should not be run.

INTERFERENCES

Reference¹⁰

Results of studies conducted show that the following substances may interfere with this procedure:

The criteria for no significant interference (NSI) is recovery within 10% or 2 mg/L (0.2 mg/dL) of the initial value.

Urine:

Creatinine:	NSI up to 300 mg/dL Creatinine
Glucose:	NSI up to 3,000 mg/dL Glucose
Urea:	NSI up to 5,000 mg/dL Urea

Ascorbate:	NSI up to 500 mg/dL Ascorbate
Citrate:	NSI up to 50 mg/dL Citrate
Magnesium:	NSI up to 400 mg/dL Magnesium
Oxalate:	NSI up to 30 mg/dL Oxalate
Conjugated Bilirubin:	NSI up to 40 mg/dL Conjugated Bilirubin
Hemoglobin:	NSI up to 500 mg/dL Hemoglobin
Acetone:	NSI up to 350 mg/dL Acetone
Uric Acid:	NSI up to 10 mg/dL Uric Acid
Urobilinogen:	NSI up to 2.25 mg/dL Urobilinogen
Acetaminophen:	NSI up to 300 mg/dL Acetaminophen
Ibuprofen:	NSI up to 400 mg/dL Ibuprofen
Metronidazole:	NSI up to 600 mg/dL Metronidazole
5-Aminosalicylate:	NSI up to 150 mg/dL 5-Aminosalicylate
Calcium:	NSI up to 78 mg/dL Calcium. Concentrations above this interfere with the test result

CSF:

Conjugated Bilirubin:	NSI up to 40 mg/dL Conjugated Bilirubin
Hemoglobin:	NSI up to 500 mg/dL Hemoglobin

The information presented is based on results from Beckman Coulter studies and is current at the date of publication. Beckman Coulter Ireland Inc. makes no representation about the completeness or accuracy of results generated by future studies. For further information on interfering substances, refer to Young¹¹ for a compilation of reported interferences with this test.

PERFORMANCE CHARACTERISTICS

PERFORMANCE CHARACTERISTICS

Data contained within this section is representative of performance on Beckman Coulter systems. Data obtained in your laboratory may differ from these values.

ANALYTIC RANGE

Urine: The analytical range for the Urine application is 7 – 450 mg/L. (0.7 – 45 mg/dL).

CSF: The analytical range for the CSF application is 10 – 450 mg/L. (1.0 – 45 mg/dL).

SENSITIVITY

The limit of detection (LoD) and Limit of Quantitation (LoQ) were determined in accordance with the CLSI EP17-A2 guideline¹². Correctly operating AU Systems should exhibit sensitivity less than or equal to 7 mg/L (0.7 mg/dL) for Urine and less than or equal to 10 mg/L (1.0 mg/dL) for CSF.

The following data was obtained on an AU5800 analyzer using Urine & CSF patient samples:

Niall Sheehy
11-08-2017

Urine Application:

LoB

0.0mg/L (0.00mg/dL)

LoD

0.7 mg/L (0.07 mg/dL)

LoQ

7 mg/L (0.7 mg/dL)

CSF Application:

LoB

0.0mg/L (0.00mg/dL)

LoD

1.3 mg/L (0.13 mg/dL)

LoQ

7 mg/L (0.7 mg/dL)

Will Shelly

11-08-2017

LoB was calculated as the upper 95th percent confidence interval. LoD is defined as the lowest albumin concentration that can be detected with a probability of 95%. LoQ is defined as the lowest concentration with an inter-assay precision of < 20% CV.

METHODS COMPARISON

Reference¹³

Patient Urine samples were run to compare this Urine Albumin assay on the AU5800 analyzer against another commercially available assay. Results of Deming regression analysis were as follows:

Slope	Intercept	r	n	Sample Range
1.09	0.3 mg/L (0.03 mg/dL)	1.00	131	8.1 - 407 mg/L (0.81 - 40.7 mg/dL)

Patient CSF samples were run to compare this CSF Albumin assay on the AU5800 analyzer against another commercially available assay. Results of Deming regression analysis were as follows.

Slope	Intercept	r	n	Sample Range
1.05	-7.7 mg/L (-0.77 mg/dL)	0.99	131	17.2 - 412 mg/L (1.72 - 41.2 mg/dL)

PRECISION

Correctly operating AU Systems should exhibit precision values less than or equal to the following for both Urine & CSF applications:

Repeatability	5% CV or 1 mg/L (0.1 mg/dL)
Within Laboratory Precision	10% CV or 2 mg/L (0.2 mg/dL)

Estimates of imprecision based on CLSI recommendations¹⁴ are consistent with typical performance. Comparative performance data for the AU5800 system evaluated using the CLSI approved guideline EP05-A2 appears in the table below. Each laboratory should characterize their own instrument performance for comparison purposes.

The following data was obtained on an AU5800 analyzer using 3 Urine & 3 CSF pools analyzed over 20 days:

11-08-12

--	--

TYPE OF IMPRECISION	SAMPLE TYPE	No. Data Points	MEAN		SD		%CV
			mg/L	mg/dL	mg/L	mg/dL	
Repeatability	Urine Pool 1	80	16.8	1.68	0.2	0.02	0.9
	Urine Pool 2	80	32.8	3.28	0.2	0.02	0.6
	Urine Pool 3	80	200	20.0	3	0.3	1.4
Within Laboratory	Urine Pool 1	80	16.8	1.68	0.7	0.07	4.3
	Urine Pool 2	80	32.8	3.28	0.8	0.08	2.5
	Urine Pool 3	80	200	20.0	4	0.4	2.1
Repeatability	CSF Pool 1	80	63.0	6.30	0.7	0.07	1.0
	CSF Pool 2	80	265	26.5	4	0.4	1.6
	CSF Pool 3	80	350	35.0	7	0.7	1.9
Within Laboratory	CSF Pool 1	80	63.0	6.30	1.1	0.11	1.8
	CSF Pool 2	80	265	26.5	6	0.6	2.4
	CSF Pool 3	80	350	35.0	9	0.9	2.5

ADDITIONAL INFORMATION

DxC 700 AU requires that each reagent application has a standard format of abbreviated Closed Test Name. This Closed Test Name is required to allow automated loading of the calibrator information for each application as part of the DxC 700 AU Closed System. Refer to the table below for the Closed Test Name assigned to each application for this assay.

Test Name	Description
ULB1G	Urine/CSF Albumin (Urine)
ULB1G	Urine/CSF Albumin (CSF)

The manufacturer guarantees that the product meets the requirements of regulatory and technical documents. Safety and quality of the product are guaranteed until the expiry date. The Manufacturer is responsible for product defects, excluding defects that appeared because of incorrect usage, shipment and storage, or third persons or compelling force. The manufacturer assures that they will substitute the product at their own expense, if technical and functional features (consumer properties) don't meet the requirements of regulatory and technical documents, if stated defects were the subsequence of a hidden defect of materials or low-quality production by the manufacturer. If you have questions concerning quality of the product please address the authorized representative of the manufacturer in the Russian Federation:

Setting Sheet Footnotes

User defined

† Beckman Coulter Urine/CSF Albumin Calibrator Cat No.: B38859.

* Values set for working in mg/L. To work in mg/dL, divide by 10.

Niall Sherry

11-08-2017

REVISION HISTORY

Russian market specific IFU

The IFU code has been revised to standardise the document naming convention for each market.

Nall Shelly

11-08-2017

Amir Khan
11-08-17

REFERENCES

1. American Diabetes Association. Diabetic Nephropathy. Diabetes Care, 2002; 25:(Suppl. 1):S85-S89.
2. Sacks DB, Bruns DE, Goldstein DE, MacLaren NK, McDonald JM, Parrot M. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clin Chem 2002;48:436-472.
3. Hortin GL. Aminoacids, peptides and Proteins. In Burtis CA, Ashwood ER, Bruns DE, eds. Tietz textbook of clinical chemistry and molecular diagnostics, 2012:507pp.
4. Lamb EJ, Price CP. Kidney Function Tests. In Burtis CA, Ashwood ER, Bruns DE, eds. Tietz textbook of clinical chemistry and molecular diagnostics, 2012:669pp.
5. Sacks DB. Diabetes Mellitus In Burtis CA, Ashwood ER, Bruns DE, eds. Tietz textbook of clinical chemistry and molecular diagnostics, 2012:1415pp.
6. Ehret W, Heil W, Schmitt Y, Töpfer G, Wissner H, Zawta B, et al. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 Rev.2: p46 & p50.
7. Tietz NW, ed Clinical guide to laboratory tests, 3rd ed. Philadelphia WB Saunders Company, 1995:p22 & p24
8. Current Intelligence Bulletin 13: Explosive Azide Hazard, August 16, 1976. DHHS (NIOSH) Publication Number 78-127.
9. Painter PC, Cope JY, Smith JL. Reference information for the clinical laboratory. In: Burtis CA, Ashwood ER, eds. Tietz textbook of clinical chemistry. Philadelphia:WB Saunders Company, 1999; 1800pp.
10. Clinical and Laboratory Standards Institute. Interference Testing in Clinical Chemistry; Approved Guideline - Second Edition. CLSI document EP07-A2, Wayne PA 2005.
11. Young DS. Effects of Drugs on Clinical Laboratory Tests, AACC, 5th ed. AACC Press, 2000.
12. Clinical and Laboratory Standards Institute. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline - Second Edition. CLSI document EP17-A2, Wayne PA 2012.
13. Clinical and Laboratory Standards Institute. Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline - Third Edition. CLSI document EP09-A3, Wayne PA 2013.
14. Clinical and Laboratory Standards Institute. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition. CLSI document EP05-A2, Wayne PA 2004.

Beckman Coulter Ireland Inc., Lismeehan, O'Callaghan's Mill's Co. Clare, Ireland (001) 703-527-3887
www.beckmancoulter.com

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1. American Diabetes Association. Diabetic Nephropathy. Diabetes Care, 2002; 25:(Suppl. 1):S85-S89.
2. Sacks DB, Bruns DE, Goldstein DE, MacLaren NK, McDonald JM, Parrot M. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clin Chem 2002;48:436-472.
3. Hortin GL. Aminoacids, peptides and Proteins. In Burtis CA, Ashwood ER, Bruns DE, eds. Tietz textbook of clinical chemistry and molecular diagnostics, 2012:507pp.
4. Lamb EJ, Price CP. Kidney Function Tests. In Burtis CA, Ashwood ER, Bruns DE, eds. Tietz textbook of clinical chemistry and molecular diagnostics, 2012:669pp.
5. Sacks DB. Diabetes Mellitus In Burtis CA, Ashwood ER, Bruns DE, eds. Tietz textbook of clinical chemistry and molecular diagnostics, 2012:1415pp.
6. Ehret W, Heil W, Schmitt Y, Topfer G, Wisser H, Zawta B, et al. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 Rev.2: p46 & p50.
7. Tietz NW, ed Clinical guide to laboratory tests, 3rd ed. Philadelphia WB Saunders Company, 1995:p22 & p24
8. Current Intelligence Bulletin 13: Explosive Azide Hazard, August 16, 1976. DHHS (NIOSH) Publication Number 78-127.
9. Painter PC, Cope JY, Smith JL. Reference information for the clinical laboratory. In: Burtis CA, Ashwood ER, eds. Tietz textbook of clinical chemistry. Philadelphia:WB Saunders Company, 1999; 1800pp.
10. Clinical and Laboratory Standards Institute. Interference Testing in Clinical Chemistry; Approved Guideline - Second Edition. CLSI document EP07-A2, Wayne PA 2005.
11. Young DS. Effects of Drugs on Clinical Laboratory Tests, AACC, 5th ed. AACC Press, 2000.
12. Clinical and Laboratory Standards Institute. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline - Second Edition. CLSI document EP17-A2, Wayne PA 2012.
13. Clinical and Laboratory Standards Institute. Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline - Third Edition. CLSI document EP09-A3, Wayne PA 2013.
14. Clinical and Laboratory Standards Institute. Evaluation of Precision Performance of Quantitative Measurement Methods, Approved Guideline - Second Edition. CLSI document EP05-A2, Wayne PA 2004.

 Beckman Coulter Ireland Inc., Lismeehan, O'Callaghan's Mill's Co. Clare, Ireland (001) 703-527-3887
www.beckmancoulter.com

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