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State of Minnesota

SECRETARY OF STATE

APOSTILLE

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1. Country: United States of America

This Public Document: Instruction for Use for 3M ESPE RelyX Luting Cement Kit

2. has been signed by Lori J Sheetz

3. acting in the capacity of Notary Public, State of Minnesota

4. bears the seal/stamp of Lori J Sheetz, Notary Public, State of Minnesota

CERTIFIED

5. at: St. Paul, Minnesota

6. Date: May 03, 2017

7. by: Secretary of State, State of Minnesota

8. File No: 9488836-3

9. Seal/Stamp:

10. Signature:

Steve Simon

Steve Simon
Secretary of State

3M ESPE
Dental Products

2510 Conway Avenue
St. Paul, MN 55144-1000

To Whom It May Concern:

In my capacity as Regulatory Affairs Specialist, 3M ESPE Dental Products, St. Paul, MN, USA, I certify that the attached English version of the Instructions for Use (IFU) for 3M™ ESPE™ RelyX™ Luting Cement Kit intended for Russia, is truthful and accurate.

Note: 3M ESPE is responsible for the information contained in the English version of the IFU. 3M Russia is responsible for the information contained in the Russian version of the IFU.

Lin Ma, Ph.D., CQE
Regulatory Affairs Specialist
3M ESPE Dental Products

02 May 2017
Date

STATE OF MINNESOTA)
COUNTY OF RAMSEY) ss

On May 2, 2017, before me, a Notary Public within and for said county, personally appeared Lin Ma, known by me to be the person described in and who executed the foregoing instrument and acknowledged that she executed the same as her own free act and deed.

Lori J. Sheetz, Notary Public
My Commission Expires, January 31, 2020

ENGLISH

RelyX Luting Cement

Instruction for Use

General Information

3M ESPE RelyX Luting cement is a glass ionomer system comprised of two parts: a powder and a liquid. The powder is a radiopaque, fluoroaluminosilicate glass. The liquid is an aqueous solution of a modified polyalkenoic acid. RelyX Luting Cement provides the major benefits of glass ionomer cements – adhesion to tooth structure, fluoride release along with additional benefits of low solubility, improved fracture toughness and a low viscosity, non-stringy, slump-resistant mix.

RelyX Luting Cement is a self-cure system using two dark-cure reactions to provide the improved properties. RelyX Luting cement is radiopaque and is available in a white shade. The powder jar contains 16gm of material and the liquid bottle contains 9ml (11gm) of material.

Indications

RelyX Luting Cement indicates for use when luting:

- Porcelain fused to metal crowns and bridges to tooth structure, amalgam composite or glass ionomer core buildups
- Metal inlays, onlays or crowns
- Crowns made with all-alumina or all-zirconia cores such as 3M ESPE Lava or Procera AllCeram
- Prefabricated and cast post cementation
- Orthodontic bands and appliances
- Porcelain fused to metal (PFM), metal, all-alumina or all-zirconia core restorations on implant abutments

Contraindications: There is no known contraindication

Scope of use: Dentistry. Only for professional use

Precautionary Information for Patients:

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of this product in patients with known acrylate allergies. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If allergic reaction occurs, seek medical attention as needed, remove the product if necessary and discontinue future use of the product.

Contains potassium persulfate that may trigger an allergic respiratory reaction in certain individuals. This product may not be appropriate for use in those individuals with known sensitivity to sulfites, since a cross-reaction may occur with potassium sulfate.

Precautionary Information for Dental Personnel:

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If product contact glove, remove and discard glove, wash hands immediately with soap and water and then re-glove. If allergic reaction occurs, seek medical attention as needed.

3M ESPE MSDS information can be obtained from www.3MESPE.com or contact your local subsidiary.

Instructions for Use

Preparation

1. Remove the provisional restoration and all residual temporary cement.
2. Thoroughly clean the preparation with an oil-free pumice paste. Rinse with water and let dry.

Note: Do not desiccate the tooth. Desiccation of tooth structure is believed to cause post-operative sensitivity in some individuals.

3. Pulp protection: Use a hand setting calcium hydroxide material with near exposures. Note: Use of some tooth preparation treatments is contraindicated as their residue may adversely affect bond strength or the setting reactions (e.g., Tubulicid Red). Contact your local 3M ESPE representative or the 3M ESPE Customer Hotline for specific details.
4. Casting preparation: Thoroughly clean interior surfaces of cast crown, inlay or onlay.

Directions

1. Dispense powder and liquid: The standard powder/liquid ratio of 1.6:1 by weight can be obtained by using an equal number of level scoops of powder and liquid drops. Three scoops of powder and three drops of liquid will provide an adequate amount of material to seat one typical crown.
 - 1.1. Shake the jar to fluff the powder before dispensing.
 - 1.2. Insert the scoop into the jar, overfill it with loosely packed powder and withdraw it against the plastic leveler to remove excess powder and obtain a level scoop.
 - 1.3. Dispense the desired number of powder scoops onto the mixing pad.
Note: The glass ionomer powders are sensitive to high humidity. Store with jar caps securely tightened and away from high humidity. Replace cap immediately after use.
 - 1.4. To obtain a proper liquid drop size, hold the RelyX Luting Cement liquid vial vertically with the dropper tip down and without the tip contacting the mixing pad.
 - 1.5. Squeeze the vial to dispense the desired number of liquid drops onto the mixing pad. Replace cap immediately after use.
2. Mixing: Using a cement spatula, mix the powder into the liquid. To minimize water evaporation and maximize working time, confine spatulation of the powder and liquid to a small area of the mixing pad. All of the powder should be incorporated into the liquid within 30 seconds.
3. Working time of the standard powder/liquid ratio is at least 2.5 minutes from the start of mix at a room temperature. Higher temperatures and vigorous spatulation will shorten working time. Lower temperatures will lengthen working time. **Note:** Working time can be lengthened by using refrigerated liquid or by mixing on a cold slab.
4. Crown placement: Load the crown by spreading a layer of the cement over all the interior surfaces of the crown. Seat crown. Maintain pressure on crown to maintain position during setting process.
5. Clean up express. Excess material can be removed when cement reaches a waxy stage after a minimum of 3 minutes from placement in the mouth (37°C or 98°F). Use a suitable instrument for this process. Cement will set hard within 10 minutes after placement in the mouth.

Storage and Use

This product is designed to be stored and used at room temperature.
Values of humidity and pressure are not regulated.

Shelf life at room temperature is 24 months from the date of manufacture.

Ambient temperatures routinely higher than +27°C/80°F, or lower than +10°C/50°F may reduce shelf life. See outer package for expiration data. Do not use after expiration date.

Glass ionomer powders are sensitive to high humidity. Store with jar caps securely tightened and away from high humidity.

- Close Jar Lid tightly after use – Moisture Sensitive Materials.
- Use an appropriate technique that minimizes jar and scoop contact with liquid disinfectants.
- Do not refrigerate.

Transportation conditions: No special transportation condition required

Do not substitute 3M ESPE Vitremer Core Buildup/Restorative powder or liquid for RelyX Luting Cement powder or liquid.

Cleaning and Disinfection

Disinfect this product using an intermediate level disinfection process (liquid contact) as recommended by the Centers for Disease Control and endorsed by the American Dental Association. Guidelines for Infection Control in Dental Health-Care Settings – *MMWR*, December 19, 2003:52(RR-17), Centers for Disease Control and Prevention.

Disposal

Packaging/contents should be disposed in accordance with the local/regional/national/ international regulations.

Medical products should be disposed by the local public law organizations in accordance with SanPiN 2.1.7.2790-10 "Sanitary and Epidemiological Requirements for Medical Waste Treatment" waste Class B).

See the Material Safety Data Sheet (available at www.3MESPE.com or through your local subsidiary) for additional disposal information.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Caution: U.S. Federal Law restricts this device to sale or use on the order of a dental professional.

Contents

The outer pack (carton box) contents: 1 Powder Jar, 1 Liquid Vial, 1 Spoon, 2 Mixing pads, 1 Instruction for Use.

Technical Characteristics

Mixing Time: 30 sec. (by hand); Working Time: ≥ 1 min. 30 sec.; Setting Time ≤ 8 min.; Film Thickness ≤ 25 μ m; Flexural Strength ≥ 10 MPa; Radiopacity ≥ 1 mm Al.

Powder: 16 g (in one Jar); Liquid: 9 ml (in one Vial); Spoon: Length 8.38 ± 0.08 cm, Width 0.86 ± 0.08 cm, Volume 86.85 mm³, Weight 1.34 ± 0.10 g; Mixing Pad: Length 7.47 ± 0.15 cm, Width 6.20 ± 0.15 cm, Number of sheets: 60.

Composition

Powder: Silanized glass, encapsulated potassium persulfate and ascorbic acid, and white pigment.

Liquid: Methacrylates, water, tartaric acid, and BHT.

Description of packaging

Inner package: Liquid Vial (Oval Vial: Low Density Polyethylene; parameters: $1.70 \pm 0.04 \times 2.72 \pm 0.04$ cm, Height: 5.05 ± 0.08 cm, Volume: 11.2 ± 1.0 ml, Weight: 3.5 ± 0.6 g, Neck Finish: 13-425. Cap for Vial: Low Density Polyethylene; parameters: Height: 1.98 ± 0.25 cm, Thread pitch: 2.0 mm, Thread length: 13.0 mm, Outer part diameter: 1.62 ± 0.25 cm, Outer diameter of the cap top: 4.0 ± 0.1 mm), Glass Jar for Powder (Glass Jar: Amber glass; parameters: Height: $49.2 \div 50.8$ mm, Diameter: $30.5 \div 32.3$ mm, Volume: 17.0 ml, Weight: 27.0 ± 1.5 g, Neck Finish: R3-24. Cap for Glass Jar: Low Density Polyethylene; parameters: Height: 1.32 ± 0.05 cm, Outer diameter: 2.74 ± 0.05 cm, Inner diameter: 1.63 ± 0.05 cm, Thread pitch: 8 threads per 25.4 mm

Outer package: Carton box (Fiberboard; box size (LxWxH) 13.97 cm x 9.84 cm x 4.29 cm)

Labelling

Products are supplied in packaging labeled with the following information: the name and REF number; the name and address of the manufacturer; content; LOT number; storage conditions; expiration date; information about trademarks the name and address of EU Authorized Representative, "Caution: U.S. Federal Law restricts this device to sale by or on the order of a dental professional"; ISO 9917-2 standard compliance; and symbols - Caution, refer to IFU, CE symbol, environmental packaging, recycling sign.

Maintenance and repair

This medical device isn't subject to maintenance and repair.

Warranty

3M ESPE warrants this product to be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Name and address of the manufacturer:

3M ESPE Dental Products, 2510 Conway Avenue, St. Paul MN 55144-1000, USA

Authorized representative on the territory of the Russian Federation:

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3M, ESPE, RelyX, Lava, and Vitremer are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. © 3M 2016. All rights reserved. 34-8718-3422-1-A