**State of Colorado
Secretary of State****APOSTILLE**

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1. Country: Pays / País:	United States of America		
This public document Le présent acte public / El presente documento público			
2. has been signed by a été signé par ha sido firmado por	Deanna Jacobs		
3. acting in the capacity of agissant en qualité de quien actúa en calidad de	Notary Public, State of Colorado		
4. bears the seal / stamp of est revêtu du sceau / timbre de y está revestido del sello / timbre de	Deanna Jacobs, Notary Public State of Colorado		
Certified Attesté / Certificado			
5. at à / en	Denver, Colorado	6. the le / el día	17th day of August, 2017
7. by par / por	Wayne W. Williams, Secretary of State of the State of Colorado		
8. Number sous n° bajo el número	5868385429		
9. Seal / stamp: Sceau / timbre : Sello / timbre:		10. Signature: Signature : Firma:	

This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears.
This Apostille does not certify the content of the document for which it was issued.
This Apostille is not valid for use anywhere within the United States of America, its territories or possessions.
To verify the issuance of this Apostille, see: www.sos.state.co.us/auth

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Date: 9/1/17

Signature: *[Handwritten Signature]*

Position: Associate RegOps Specialist

Name: Deanna Jacobs

Addition to Instruction for use

1. MEDICAL DEVICE NAME

Spectra Optia platelet, plasma set

2. INFORMATION ON MANUFACTURERS OF THE MEDICAL DEVICE

Developer

Terumo BCT, Inc.

10811 W. Collins Avenue

Lakewood, Colorado 80215

USA

Tel: +1 303 231 4357

Fax: +1 303 542 5215

E-mail: Pascale.VanHoydonck@terumobct.com

www.terumobct.com

Manufacturer

Terumo BCT, Inc.

10811 W. Collins Avenue

Lakewood, Colorado 80215

USA

Tel: +1 303 231 4357

Fax: +1 303 542 5215

E-mail: Pascale.VanHoydonck@terumobct.com

www.terumobct.com

Manufacturing Facility

Terumo BCT, Inc.

10811 W. Collins Avenue

Lakewood, Colorado 80215

USA

Tel: +1 303 231 4357

Fax: +1 303 542 5215

E-mail: Pascale.VanHoydonck@terumobct.com

www.terumobct.com

Authorize Representative of Manufacturer in Russia

Terumo Russia LLC

BC Northern Tower, 13th Floor,

10 Testovskaya Street, Moscow 123317, Russia

Tel: +7-495-988-4740

Fax: +7-495-988-4739

E-mail: andrey.shevelev@terumobct.com

E-mail: Quality.EMEA@terumobct.com

www.terumobct.com

3. THE USE OF THE MEDICAL DEVICE FOR WHICH INTENDED BY THE MANUFACTURER

Indications

The Spectra Optia platelet and plasma set is intended for use in platelet and plasma collections. It may collect the following products alone or in combination:

- Platelets Pheresis (single, double, or triple units)
- Platelets Pheresis, Leukocytes Reduced (single, double, or triple units)
- Plasma
- Plasma, Leukocytes Reduced

Counter-indications

There are no known contraindications for the use of the Spectra Optia platelet and plasma set except those associated with the infusion of solutions and replacement fluids as required by the apheresis procedure, and those associated with all types of automated apheresis systems.

The platelet, plasma set is used on healthy donor for the collection of blood component products.

The Platelet, plasma disposable by itself does not interact with implantable devices.

The Disposable Kits described here contain phthalates (DEHP). Patient groups including pregnant or nursing women and children are considered to be the most at risk of harmful effects due to DEHP exposure. Based on five decades of successful use in medical devices without a reported adverse event related to DEHP, the risk of not doing a needed procedure is far greater than the risk associated with exposure to DEHP.

There is no influence on feature of driving, mechanisms.

Possible side effects are those associated with therapeutic aphaeresis procedure

Possible side effects are those associated with therapeutic aphaeresis procedure.

Platelets are collected with target cells during continuous mononuclear cell collection (CMNC) procedures, granulocyte (PMN) collection procedures, and white blood cell depletion (WBCD) procedures. If the initial platelet count is low and the patient is not properly managed, excessive platelet depletion can occur.

Red blood cells (RBC) are collected with target cells during granulocyte (PMN) collection procedures, white blood cell depletion (WBCD) procedures, and platelet depletion (PLTD) procedures. If the procedure is not properly monitored, excessive RBC depletion may occur. To monitor excessive RBC depletion, consider testing the patient's hematocrit throughout the run.


Do not use a blood warmer on the replace line during a white blood cell depletion (WBCD) procedure or a platelet depletion (PLTD) procedure. Using a blood warmer on the replace line slows the flow of replacement fluid into the reservoir, affecting the system's management of the fluid balance.

During collection and depletion procedures, the operator must monitor the plasma line for hemolysis. The RBC detector is not used to monitor the plasma line during these procedures. If hemolysis occurs that is not related to the patient's condition, the operator should consider consulting the physician to determine whether or not to continue the procedure and perform rinseback.

Before using a granulocyte-colony stimulating factor (G-CSF) or other agents to prepare the patient for a procedure performed on the Spectra Optia system, consult the manufacturer's package insert for potential contraindications and adverse events.

When transferring plasma from the plasma bag to the collection bag during an MNC collection procedure, you must put the collection bag lower than the plasma bag as indicated in the instructions on the screen. Otherwise, cells in the collection bag can flow into the plasma bag.

When performing a collection procedure using the IDL Set, seal the collect line above the T-shaped connector to the plasma line as shown in the illustration below. If you seal the collect line below the connector, the cells in the collection bag could flow back into the tubing set when the system raises the cassette.

If you use a low inlet pump flow rate during a CMNC procedure, it will take the system longer to process the required volume to establish the interface and to collect the target cells. Be aware of how this is relevant to a patient with a small total blood volume. Also, try to avoid conditions that will cause a delay in establishing the interface, such as pump pauses.

Indication of possibility and peculiarities of using the medical device by people with implants medical devices, pregnant women, breastfeeding women, children, adults with chronic diseases

The plasma, platelet disposable by itself does not interact with implantable devices.

The Disposable Kits described here contain phthalates (DEHP). Patient groups including pregnant or nursing women and children are considered to be the most at risk of harmful effects due to DEHP exposure. Based on five decades of successful use in medical devices without a reported adverse event related to DEHP, the risk of not doing a needed procedure is far greater than the risk associated with exposure to DEHP.

Before undergoing apheresis, patients with active implantable medical devices should review the implant manufacturer's instructions for any cautions or contraindications concerning the use of and/or proximity to devices that emit electromagnetic energy. Examples of active implantable devices include cardiac pacemakers and cardioverter defibrillators, cochlear implants, vagus nerve stimulators, and devices that stimulate the peripheral or central nervous systems. This caution also extends to other devices that may be connected to the patient such as ECG equipment and infusion systems. Care should be taken to prevent the Seal Safe system from touching the patient or any cables of the devices that may be connected to the patient.

Conditions of use

Use of the Spectra Optia platelet and plasma set is restricted to prescription only.

Each operator should be thoroughly familiar with the Spectra Optia set's operating instructions before using the system. All procedures should be performed by qualified medical personnel under the supervision of a physician. Spectra Optia platelet, plasma set should be used indoor only and is typically used in a medical institution.

Characteristics	Performance	Conditions
Temperature for use	10°C to 40°C (50°F to 104°F) ○	
Humidity for use	10% to 90%	Relative humidity (RH), non-condensing
Storage and transportation pressure	700 - 1060 mbar (from -400 m to 3000 m)	

4. METHOD OF APPLICATION

- 1) Connect the patient or donor to the sterile, non-pyrogenic platelet, plasma disposable tubing set using needle. The needle connected by donor line
- 2) Connect the tubing set to an anticoagulant (ACD-A), saline, and any needed replacement fluid. The set is loaded onto the equipment so that the flow through various portions of the disposable set can be controlled by pumps and valves in the equipment.

The requirements for assembly and installation

The blood and fluid pathways of the tubing set are sterilized with ethylene oxide and are nonpyrogenic. Do not use the set if any of the following conditions are true:

- Severe kinks in the tubing are apparent, such that the interior diameter of the tubing is smaller than the interior diameter of the patient's inlet or return access.
- Tubing set appears to be incorrectly assembled.
- Tubing set is damaged.
- Any clamps are closed on the lines.
- The end caps on the lines of the tubing set are missing when you open the package.

You may load the tubing set up to 24 hours before the procedure as long as you do not lower the cassette or prime the set. Once you lower the cassette or prime the set, you must use the set during the same work shift.

The instructions below to unpack the tubing set needs to be followed.


Load the set on the system, and test the set. After you load the set, the system performs a test to ensure that you have loaded the correct set for the procedure.

- Verify that the tubing set has not expired by checking the expiration date on the cover of the package.
- Document the lot number of the tubing set to use for future reference.
- Take the product bags and the vent bag out of the package and hang the bags from left to right on the IV pole in the following order: vent bag, plasma bag, collection bag.
- Take out the replace line and put the line over the right side of the front panel between the two sides of the IV pole.
- Take out the coiled inlet line (red clamps) and remove the paper tape from the coil. Hang the inlet connection on the left end of the IV pole. Repeat this step with the return line (blue clamps).
- Take out the AC line (orange) and the saline line (green) and hang the lines over the left side of the front panel.

Loading the Cassette

Remove the cassette from the package.

Using both hands, position the bottom of the cassette inside the bottom edge of the cassette tray, aligning the pump headers with the pumps. Ensure that the bottom of the cassette is captured on all three points along the bottom rail (Figure below - Spectra Optia Cassette Loading Points).

Using both hands, press the upper corners of the cassette to clamp it into the cassette tray until it "clicks" into fully loaded position.

Loading the Channel into the Centrifuge

Take the channel out of the package and put it on top of the centrifuge cover. Set the empty package aside.

Open the centrifuge door.


Locate the pin on the filler latch. Raise the latch by pushing the pin toward the center of the filler while pulling up the latch.

Turn the centrifuge so that the loading port faces you.

Extend the lines between the cassette and the channel, and ensure that they are not twisted.

Pull the channel up through the loading port and the through opening in the center of the filler. To avoid damaging the chamber and the collect line when using a Collection Set, guide the chamber through the loading port with one hand as you pull the channel through the loading port with the other hand (Figure below). Lower the filler latch and lock it in place.



Figure - Pushing the channel through the loading port

Position the lower collar in the collar holder on the filler latch so that the pink line is not obstructed by the other lines. Ensure that the base of the pink line aligns with the space between the two screws on the filler latch (see the figure - Lower collar correctly positioned and locked with the notch visible).



Grasp the centrifuge loop below the lower collar and gently pull down the collar until you hear the "click" of the locking pin as it pops out and locks the collar in the collar holder. Ensure that the notch at the base of the locking pin is visible. If the collar is locked, you can see the notch (see the figure - Lower collar correctly positioned and locked with the notch visible).

Starting with the connector, insert the channel into the groove in the filler, finishing with the inlet port (Figure below - Pressing the channel into the groove). Run your finger over the groove and push down any section of the channel that is not completely inserted in the groove. The channel must sit flush with the groove.


bearing
upper
the lower
upper
the loop

Insert the narrow part of the lower bearing into the lower holder, and the narrow part of the upper bearing into the bearing holder (Figures below - Lower bearing inserted in bearing holder (left) and the upper bearing inserted in the bearing holder (right)). Ensure that the braided section of is not twisted.

Position the upper collar below the upper collar holder and insert the line into the holder. Pull the line up to secure the upper collar in the holder (Figure below - Inserting the upper collar into the holder).

5. CLASSIFICATION OF MEDICAL DEVICE

1. The Spectra Optia platelet, plasma set is classified as Class IIb under Rule 6 of Annex IX of the Council Directive 93/42/EEC concerning Medical Devices (MDD);
2. It is invasive medical device;
3. For single use only;
4. The Optia platelet, plasma set is sterile device. Sterilization by ethylene oxide;
5. The contact with human body is short (less than 24 hours)

6. PRECAUTION MEASURES APPLIED WHEN USING THE MEDICAL DEVICE

Warnings

1. If the integrity of the tubing set is compromised for any reason, the tubing set is no longer functionally closed, and the product should not be stored beyond 24 hours.
2. Tubing sets manufactured by Terumo BCT contain DEHP. Results from selected studies of rodents exposed to DEHP (typically via an oral route) have shown potentially harmful effects to male reproductive organs during fetal development. Because of this, patient groups that include pregnant or nursing women and children are considered to be most at risk to potential harmful effects of exposure to DEHP. However, no studies performed to date involving other mammals and humans exposed to intravenous DEHP have shown similar results. Regulatory bodies have noted that the risk of not doing a needed procedure is far greater than the risk associated with exposure to DEHP. It is the responsibility of treating physicians to balance this risk for their patients.

Cautions

1. **Not for Reuse:** Tubing sets manufactured by Terumo BCT bearing the not for reuse symbol are intended for single use only and are not intended to be reprocessed or reused. Terumo BCT cannot ensure the functionality or sterility of the set if it is reprocessed or reused.
2. Each operator should be thoroughly familiar with the Spectra Optia system's operating instructions before using the system. All procedures should be performed by qualified medical personnel under the supervision of a physician.
3. Use aseptic technique throughout all procedures to ensure donor safety and product quality.
4. The blood and fluid pathways of the tubing set are sterilized with ethylene oxide and are nonpyrogenic. Under any of the following conditions, DO NOT use the set:
 - The end caps are not in place.
 - Severe kinks in the tubing are apparent.
 - The tubing set is incorrectly assembled.
 - The tubing set is damaged.
5. When loading the channel and the LRS chamber into the centrifuge filler, be careful not to stretch the tubes, as tubing leaks could develop.
6. Use only your fingers to load the channel into the filler. To avoid puncturing the channel, never use a sharp object.

7. PRINCIPLE OF MEDICAL DEVICE OPERATION

The Spectra Optia platelet, plasma set is used with the Spectra Optia System. It uses well established technologies and methods to complete the apheresis procedures.

The patient or donor is connected to the sterile, non-pyrogenic platelet, plasma disposable tubing set. The tubing set is connected to an anticoagulant (ACD-A), saline, and any needed replacement fluid. The set is loaded onto the equipment so that the flow through various portions of the disposable set can be controlled by pumps and valves in the equipment.

The system automatically collects the products selected. The system first completes a blood prime, during which the Status Line message "Blood prime in progress" displays.

Anticoagulated blood is pumped up the draw line until the blood reaches the centrifuge.

The return pump then reverses and draws anticoagulated blood up the return line until the lower level sensor in the return reservoir detects fluid. At this point, blood prime is complete and the draw cycle begins.

The system generates an audible tone to alert you to the start of the first return cycle.

The system first begins collecting plasma. It then begins collecting platelets concurrently with plasma, pausing plasma collection occasionally and only collecting platelets. If necessary, the system completes plasma collection after platelet collection is done.

When all products are collected, the system automatically continues into Rinseback.

This procedure requires the use of the Spectra Optia plasma, platelet set.


The patient's blood is continuously pumped through the set.

The inlet pump pumps the anticoagulated whole blood through the inlet line into the channel.

The buffy coat is separated from the RBC and plasma in the connector.

The collect pump pumps the buffy coat from the channel into the collection bag. At the beginning of the procedure the collect valve is in the return position. As soon as cells are detected in the collect port by the AIM system, the collect valve moves to the collect position and the cells are pumped into the collection bag.


Cell Separation: Channel

1. Whole blood enters the channel.
2. Red blood cells (RBC) flow passively to the reservoir.
3. Plasma is pumped into the reservoir or it is collected.
4. Cells are continuously pumped into the collection bag.


8. TECHNICAL SPECIFICATIONS OF THE MEDICAL DEVICE

Structure (basic) and list of accessories (MD complete set options).

- | | | |
|-------------------------------|----------------------------------|-----------------------------|
| 1 Donor line | 13 Return reservoir | 25 Inlet port |
| 2 Sample bag | 14 Return pump header | 26 Collection chamber |
| 3 AC/draw/return manifold | 15 Centrifuge pressure sensor | 27 LRS chamber |
| 4 Anticoagulant (AC) line | 16 Inlet line to centrifuge | 28 Plasma collect line |
| 5 Sterile barrier filter | 17 RBC line from centrifuge | 29 Plasma bag |
| 6 Cassette | 18 Platelet line from centrifuge | 30 Platelet collect line |
| 7 AC pump header | 19 Plasma line from centrifuge | 31 Platelet bag |
| 8 Draw/return pressure sensor | 20 Plasma pump header | 32 Platelet product sampler |
| 9 AC Line | 21 Cassette label | 33 Vent bag line |
| 10 Return line | 22 Collect pump header | 34 Vent bag |
| 11 Draw line | 23 Centrifuge loop | |
| 12 Inlet pump header | 24 Channel | |

Photo or image

This set consists of a preconnected channel, cassette, product bags, and blood tubing. The Spectra Optia Platelet, Plasma Set features a cassette that incorporates self-loading pump headers, filters, pressure diaphragms, sensors, and a return reservoir. The cassette provides the hardware interface for sensors, valves, and an RBC detector.

The tubing set is designed for ease of loading and disposal. It routes blood components from the donor needle site through the Spectra Optia system and back to the donor. The donor lines, the channel, and the product bags are preconnected. The blood separator channel is a sterile, nonpyrogenic device that spins in a centrifuge to separate blood components.

1. Donor line: provides access to and return to the donor.
2. Sample bag: used to collect blood from the donor for sampling purposes
3. AC/draw/return manifold: consists of the access injection site and the connections for the anticoagulant line (4), the draw line (11), and the return line (10).
4. Anticoagulant (AC) line (with orange spike): carries anticoagulant from the AC bag to the cassette (6).
5. Sterile barrier filter: prevents bacteria from entering the system to maintain a functionally closed system.
6. Cassette: guides the flow of blood and products through the tubing set.
7. AC pump header: the tubing segment that fits into the AC pump.
8. Draw/return pressure sensor: monitors the draw and return pressures at the donor access site.
9. AC line: carries AC from the cassette (6) to the AC/draw/return manifold (3).
10. Return line: carries blood components back to the donor.
11. Draw line: carries anticoagulated whole blood into the tubing set.
12. Inlet pump header: the tubing segment that fits into the inlet pump.
13. Return reservoir: receives uncollected components for return to the donor. Also provides a return filter (200 micron) and air detection.
14. Return pump header: the tubing segment that fits into the return pump.
15. Centrifuge pressure sensor: detects high pressure in the centrifuge.
16. Inlet line to centrifuge: carries blood to the centrifuge.
17. RBC line from centrifuge: carries RBC (red blood cells) from the centrifuge for collection or return to the donor.
18. Platelet line from centrifuge: carries platelets from the centrifuge for collection or return to the donor.
19. Plasma line from centrifuge: carries plasma from the centrifuge for collection or return to the donor.
20. Plasma pump header: the tubing segment that fits into the plasma pump.
21. Cassette label: used by the Spectra Optia system RBC detector to identify a tubing set as capable of collecting platelet products.
22. Collect pump header: the tubing segment that fits into the collect pump.
23. Centrifuge loop: consists of lines that are loaded in the centrifuge. Includes: Four-lumen tubing: used to carry fluid into and out of the channel.
Sleeves: used to reinforce the tubing.

Collars: used to secure the ends of the loop in the centrifuge.

Bearings: used to connect the loop to the centrifuge.

24. Channel: used during centrifugation to separate the donor's blood into cellular components.

25. Inlet port: routes incoming blood into the channel.

26. Collection chamber: routes separated blood products to the appropriate collection lines.

27. LRS chamber: provides leukoreduction of collected platelets.

28. Plasma collect line: carries the collected plasma to the plasma bag (29).

29. Plasma bag: 1 L bag that holds concurrently collected plasma.

30. Platelet collect line: carries the collected platelets to the platelet bags (31).

31. Platelet bag: the bag where the collected platelets are stored. The bag is made from PVC with a citrate plasticizer.




32. Platelet product sampler: used for product quality testing





33. Vent bag line: carries displaced air to and from the vent bag (34)









34. Vent bag: holds displaced air from the system.




There is no accessory related to the Spectra Optia platelet, plasma set.

Technical characteristics of the components

	Donor line	Dimensions: Needle - 17 gauge (Nominal outer diameter – 1.473 mm $\pm$ 0.013 mm; Nominal inner diameter - 1.067mm $\pm$ 0.038mm; Nominal wall thickness - 0.203 $\pm$ 0.013 mm Length of cannula: 24.3 $\pm$ 0.0127 mm Tubing: Wall Thickness: 1.00 mm Outer Diameter: 5.00 +0.05/-0.1 mm Inner Diameter: 3 $\pm$ 0,1 mm
	Sample bag	Dimensions – 7,95 x 11,63 cm
	AC/draw/return manifold	Wall Thickness: 0.94 $\pm$ 0,05 mm Outer Diameter: 4.75 $\pm$ 0,1 mm Inner Diameter: 2.87 mm
	Anticoagulant (AC) line (with orange spike)	Dimensions – Length: 43,18 cm + 7,62 cm; Wall Thickness: 0.94 $\pm$ 0,05 mm Outer Diameter: 4.75 mm Inner Diameter: 2.87 $\pm$ 0,05 mm
	Sterile barrier filter	Filter Porous Dimensions: 0.2 microns Filtration rate: 0,45 liters/hour at 914.4 mm (36 inches) head height
	Cassette	Length 253.1 mm to 242.2 mm; Height 105.2 mm plus 0.46mm minus 0.38mm Depth at the corners 17.3 $\pm$ 0.2 mm.

	AC pump header	Dimensions: Wall Thickness: $0.94 \pm 0,05$ mm Outer Diameter: 4.75 mm Inner Diameter: $2.87 \pm 0,05$ mm
	Draw/return pressure sensor	Diameter – $16 \pm 0,05$ mm (0.590 Inch)
	AC line	Dimensions: Length: $182,88 \pm 5.08$ cm (72 in.); Outer diameter: $4,75 \pm 0,05$ mm (0,187 in.) Inner diameter: 2,89 mm (0,113 in.), Thickness of wall: 1,8 mm (0,071 inch)
	Return line	Dimensions: Length: $182,88 \pm 5.08$ cm (72 in.); Outer diameter: $4,75 \pm 0,05$ mm (0,187 in.) Inner diameter: 2,89 mm (0,113 in.), Thickness of wall: 1,8 mm (0,071 inch)
	Draw line	Dimensions: Length: $182,88 \pm 5.08$ cm (72 in.); Outer diameter: $4,75 \pm 0,05$ mm (0,187 in.) Inner diameter: 2,89 mm (0,113 in.), Thickness of wall: 1,8 mm (0,071 inch)
	Inlet pump header	Dimensions - 43,18 cm (length 17 in.) + 7,62 cm (3 in.); Wall Thickness: $0.94 \pm 0,05$ mm Outer Diameter: 4.75 mm Inner Diameter: $2.87 \pm 0,05$ mm
	Return reservoir	Dimensions: Top Width: 48 ± 5 mm Bottom Width: 30.5 ± 5 mm Height: 63.5 ± 2.54 mm Depth: 41.91 ± 2.54 mm
	Return pump header	Dimensions: Wall Thickness: 0.94 ± 0.05 mm Outer Diameter: 6.63 mm Inner Diameter: 4.8 ± 0.076 mm

	Centrifuge pressure sensor	Diameter: 14.48 ± 0.05 mm
	Inlet line to centrifuge	Dimensions: Outer Diameter: 3.98 ± 0.076 mm Inner Diameter: 2.38 ± 0.1 mm
	RBC line from centrifuge	Dimensions: Outer Diameter: 3.98 ± 0.076 mm Inner Diameter: 2.38 ± 0.1 mm
	Platelet line from centrifuge	Dimensions: Outer Diameter: 2.38 ± 0.076 mm Inner Diameter: 1.09 ± 0.1 mm
	Plasma line from centrifuge	Dimensions: Outer Diameter: 3.07 ± 0.076 mm Inner Diameter: 1.49 ± 0.05 mm
	Plasma pump header	Dimensions: Wall thickness: 0.93 ± 0.05 mm Outer Diameter: 4.75 mm Inner Diameter: 1.5 ± 0.05 mm
	Cassette label	Dimensions: Width: 25.4 ± 2.54 mm Height: 12.7 ± 2.54 mm
	Collect pump header	Dimensions: Wall thickness: 0.93 ± 0.05 mm Outer Diameter: 4.75 mm Inner Diameter: 1.5 ± 0.05 mm
	Centrifuge loop	Length – 49.83 ± 1.27 cm (19.1 inch);
	Channel	Length – 65.24 ± 0.254 cm Width – 3.81 ± 0.005 cm
	Inlet port	Dimensions: Outer Diameter: 7.45 ± 0.762 mm Inner Diameter: 3.7 ± 0.127 mm

	Collection chamber	Length: 52.76 ± 0.127 mm Height: 35.48 ± 2.54 mm
	LRS chamber	Lower diameter: 2.54 ± 0.05 mm Upper diameter: 2.34 ± 0.05 mm
	Plasma collect line	Dimensions: Length: 76.2 ± 5.08 cm Outer diameter: 4.06 ± 0.076 mm; Inner diameter: 2.95 mm; Wall thickness 0.47 ± 0.05 mm
	Plasma bag	Dimensions Height: 278 ± 0.1 mm, Width: 128 ± 0.05 mm, thickness – 0,47 mm thickness of wall – 0,112 cm Weight 0,045 kg
	Platelet collect line	Dimensions Length: 30.48 ± 0.64 cm Outer diameter: 4.06 ± 0.0762 mm Inner diameter: 2.95 mm Thickness of wall – 0.7 ± 0.05 mm
	Platelet bag	Dimensions: Height: 250 ± 4.8 mm Width: 175.5 ± 4.06 mm bond socket: inner diameter 4.01 mm, tubing: outer diameter: 4.06 mm thickness wall tubing: 0,56 mm thickness – 0,38 mm Weight 0,43kg
	Platelet diameteruct sampler	Dimensions: Length: 77.3 ± 0.762 mm Diameter: 1.27 ± 0.02 mm
	Vent bag line	Dimensions Length: 76.2 ± 5.08 cm; Outer diameter: 4.75 mm Inner diameter 2.89 ± 0.05 mm thickness of wall – 1.8 ± 0.05 mm (0,071 inch)

	Vent bag	Dimensions Height: 204.2 ± 4.57 mm Width: 117.07 ± 5.1 mm bond socket: Inner diameter: 4.75 mm Thickness – 0,305 mm Weight 0,022 kg
---	----------	---

General Specifications

Extracorporeal Volume	- Approximately 196 ml
Possibility to collect products	- Platelets Pheresis (single, double, or triple units) - Platelets Pheresis, Leukocytes Reduced (single, double, or triple units) - Plasma - Plasma, Leukocytes Reduced
Weight	0.983kg (2.16lb)

Shelf life (expiry date)

The product expires 2 years from date of manufacture. Refer to product labelling for expiry date.

Furnishing options of the medical device

There's only 1 furnishing option of the Spectra Optia platelet, plasma set

Opportunities and ways of integration with other MD

The Spectra Optia Platelet, plasma set has to be used with the Spectra Optia System.

9. INFORMATION ON TECHNICAL MAINTENANCE

Spectra Optia platelet, plasma set is a single-use medical device and it does not need technical maintenance

10. PROCEDURE OF DISPOSAL (UTILIZATION) AND DESTRUCTION OF THE USED MEDICAL DEVICE (indication (if necessary) special precautionary measures in the destruction unused MD)

The disposable Set should be disposed of in accordance with local regulations.

11. CONDITIONS OF TRANSPORTATION AND STORAGE

Environmental: Tubing sets

The product needs to be kept away from the sun and be used in a dry environment.

Characteristics	Performance	Conditions
Storage and transportation temperature	- Long term storage range: 0°C to 35°C (32°F to 95°F) - Permitted excursions: -29°C to 0°C (-20°F to 32°F) for up to 72 hours 35°C to 50°C (95°F to 122°F) for up to 6	

	weeks	
Storage and transportation humidity	0% to 75%	Relative humidity (RH), non-condensing
Storage and transportation pressure	500 - 1060 mbar (from -400 m to 3000 m)	

The device can be transported by any roofed transport when following conditions of transportation

Packaging specification

Quantity of tubing sets per case	6
External Case dimension	750 mm x 350 mm x 240 mm (29.53 in x 13.78 in x 9.45 in) (LxWxH)
Case weight	- Containing disposables: 6.27 kg (13.82 lb) - Disposable Tray Lid: 0.009 kg (0.02 lb) - Disposable Tray: 0.15 kg (0.33 lb) - Box: Cardboard 0.83 kg (1.83 lb) - Box Divider: Cardboard 0.1 kg (0.22 lb)
Tubing Set Primary Dimension	370 mm x 350 mm x 76 mm (17.57 in x 13.78 in x 2.99 in) (LxWxH)

12. WARRANTY LIABILITIES

When purchasing a Terumo BCT product, the general terms and conditions are applicable. This includes the following warranty statement:

LIMITED WARRANTY. Warranty information applicable to the equipment or disposables specified on this invoice is set forth in the price list and any contract for such equipment or disposables and shall be provided by Terumo BCT upon request. For a period of ninety (90) days from the date of the invoice, Terumo BCT will repair or replace Terumo BCT spare parts or optional feature components that Terumo BCT finds defective in material or workmanship. Terumo BCT is not responsible for, and shall have no obligation with respect to any failure caused by normal wear and tear, misuse, unauthorized alterations, accident, neglect, use of non-standard accessory attachments, and/or improper customer maintenance. Terumo BCT does not verify the safety or efficacy of, and makes no warranties, express or implied, with respect to any non-Terumo BCT components or non-Terumo BCT spare parts included at Purchaser's request or any components or products used other than as expressly intended by their manufacturer. Failure to perform preventive maintenance procedures will invalidate this warranty. Terumo BCT is not responsible for any events related to non-Terumo BCT parts installed in the equipment. **THIS LIMITED WARRANTY CONTAINS THE PURCHASER'S EXCLUSIVE WARRANTIES AND REMEDIES.**

Return of used product

If for any reason this product must be returned to Terumo BCT, Inc., a returned goods authorization (an RGA number) is required from Terumo BCT prior to shipping. Instructions for cleaning and materials, including appropriate shipping containers, proper labeling, and an RGA number, may be obtained from the Terumo BCT Quality Assurance Department. **IT IS THE RESPONSIBILITY OF THE HEALTH CARE INSTITUTION TO ADEQUATELY PREPARE AND IDENTIFY THE PRODUCT FOR RETURN SHIPMENT.** Please contact your local representative for information regarding returned goods and product complaints.

CLAIM

Representative office in Russia:

Terumo Russia LLC
 BC Northern Tower, 13th Floor,
 10 Testovskaya Street, Moscow 123317, Russia
 Tel: +7-495-988-4740
 Fax: +7-495-988-4739
 E-mail: Andrey.Shevelev@terumobct.com
 E-mail: Quality.EMEA@terumobct.com

www.terumobct.com

13 METHODS AND MEANS OF DISINFECTION AND PRESTERILIZING CLEANING

The Spectra Optia plasma, platelet set is for single use. This medical device is provided sterile and does not need in additional sterilization and cleaning.

14 METHODS AND CONDITIONS OF STERILIZATION, THE TERM OF MAINTAINING STERILITY

Not for Reuse: Tubing sets manufactured by Terumo BCT bearing the not for reuse symbol are intended for single use only and are not intended to be reprocessed or reused. Terumo BCT cannot ensure the functionality or sterility of the set if it is reprocessed or reused. Spectra Optia plasma, platelet set (the blood and fluid pathways of the tubing set) is sterilized with ethylene oxide and are nonpyrogenic. Spectra Optia plasma, platelet set is sterilized under the Terumo BCT Sterilization Processing Group 1 which is sterilized using the Terumo

BCT cycle P/N 377200-091. The process used to sterilize the Spectra Optia plasma, platelet set will result in a sterility assurance level (SAL) of $< 10^{-6}$ in accordance with ANSI/AAMI/ISO 11135-1:2007.

The Spectra Optia product family is outgassed under the Trima Aeration Processing group. These products are outgassed according to MOP P/N 377200-102. The process used to outgas the Spectra Optia plasma, platelet set will result in residuals ~ 6 mg EtO and ≤ 6 mg ethylene chlorohydrin (ECH) which meet the requirements listed in ANSI/AAMI/ISO 10993-7:2009.

Spectra Optia plasma, platelet set is outgassed using the heated aeration room at Terumo BCT Sterilization Services. Presently there are 3 different heated aeration processes used for outgassing. All three processes ensure products are released from heated aeration meeting the Terumo BCT requirements of residual levels for EtO and ethylene chlorohydrin (ECH) which are 6 mg or less. These limits meet the requirements listed in ANSI/AAMI/ISO 10993-7:2009

- *Biological Evaluation of Medical Devices - Part 7 Ethylene Oxide Sterilization Residuals*. Each process was validated using a master product to represent all other products in the processing group. The master product pallet configuration described above was also used for performing the outgassing validations. The following Table describes the 3 heated aeration processes.

Heated Aeration Processes

Heated Aeration Processing Group / MOP P/N	Minimum and Maximum Heated Aeration Times	Heated Aeration Processing Group Master Product	Quality Labs Validation Report
Trima Aeration Processing Group – MOP P/N 377200-102	1 day and 16 hours – 6 days and 8 hours	777800-550	R/N 09-0012
Trima Aeration Processing Group – MOP P/N 377200-103	3 days and 1 hour – 6 days and 8 hours	70500	R/N 09-0013
Trima Aeration Processing Group – MOP P/N 377200-106	2 days and 3 hours – 6 days and 8 hours	N/A (only 1 product configuration)	R/N 12-0053

Spectra Optia plasma, platelet set maintains the sterility throughout the shelf-life (2 years).

The product is sterilized by ethylene oxide and has no other methods of sterilization

15 INFORMATION OF MARKING AND PACKING

Products supplied in packing marked with the following information:

- device name and reference number;
- LOT of the product;
- trademark and address of the manufacturer;
- Information about sterile;
- Number of products;
- expiry date;
- and other information and symbols.

Users can find the lot number adjacent to this symbol in the following locations on the packaging: box label, disposable package label, collection bag label, and plasma bag label. (Note: Lot numbers are structured for Terumo BCT internal use only. Refer to the Terumo BCT Support Center (TERUMOBCT.COM) for certificates of compliance.

Users can find the expiration date adjacent to this symbol in the following locations on the packaging: box label, disposable package label, collection label, and plasma bag label

The carried out tests of packaging of medical devices subject to final sterilization:

1 Objective: Test 1 / Fit on Spectra Optia System

- 1.1 The objective of this verification was to verify that the disposable set packaged in the deep tray meets its requirement related to fitting on the Optia system during loading.

7 Conclusion

- 7.1 All acceptance criteria were met. Test 1 resulted in a 'Pass.'

8 Objective: Test 2 / Open without tools

- 8.1 The objective of this verification was to verify that the Spectra Optia tray sealed with Tyvek can be opened without tools.

14 Conclusion

- 14.1 All acceptance criteria were met. Test 2 resulted in a 'Pass.'

15 Objective: Test 3 / Case Weight

- 15.1 The objective of this verification was to verify that the weight of a case of product does not exceed 35 pounds.

21 Conclusion

- 21.1 All acceptance criteria were met. Test 3 resulted in a 'Pass.'

22 Objective: Test 4 / Recyclable materials

- 22.1 The objective of this verification by inspection was to verify that the Spectra Optia tray and carton are made of recyclable materials.

28 Conclusion

- 28.1 All acceptance criteria were met. Test 4 resulted in a 'Pass.'

29 Objective: Test 5 / Breaking down and nesting of materials

- 29.1 The objective of this verification by inspection was to verify that the Spectra Optia packaging materials can be broken down and nested.

35 Conclusion

- 35.1 All acceptance criteria were met. Test 5 resulted in a 'Pass.'

36 Objective: Test 6 / Sharp edges and corners

- 36.1 The objective of this verification by inspection was to verify that the Spectra Optia tray does not have any sharp edges or corners.

42 Conclusion

- 42.1 All acceptance criteria were met. Test 6 resulted in a 'Pass.'

43 Objective: Test 7 / Tyvek color

- 43.1 The objective of this verification by inspection was to verify that the Spectra Optia primary packaging (tray and Tyvek) includes a visual indicator to identify set type (in this case, differing printing colors on the Tyvek lidstock).

49 Conclusion

- 49.1 All acceptance criteria were met. Test 7 resulted in a 'Pass.'

50 Objective: Test 8 / Tamper evidence


- 50.1 The objective of this verification by inspection was to verify that the Spectra Optia primary packaging provides tamper evidence.

56 Conclusion

- 56.1 All acceptance criteria were met. Test 8 resulted in a 'Pass.'

Results:
Test 1:
☒ **PASS** ☐ **PASS WITH DEVIATION** ☐ **FAIL**
Test 2:
☒ **PASS** ☐ **PASS WITH DEVIATION** ☐ **FAIL**
Test 3:
☒ **PASS** ☐ **PASS WITH DEVIATION** ☐ **FAIL**
Test 4:
☒ **PASS** ☐ **PASS WITH DEVIATION** ☐ **FAIL**
Test 5:
☒ **PASS** ☐ **PASS WITH DEVIATION** ☐ **FAIL**
Test 6:
☒ **PASS** ☐ **PASS WITH DEVIATION** ☐ **FAIL**
Test 7:
☒ **PASS** ☐ **PASS WITH DEVIATION** ☐ **FAIL**
Test 8:
☒ **PASS** ☐ **PASS WITH DEVIATION** ☐ **FAIL**
Symbols:

	Indicates the product quantity when the quantity is placed in the square
	Indicates the expiration date of the product when accompanied by a specific date
	Indicates the date of manufacture (or sterilization date, if the product is sterile) when accompanied by a specific date
	Indicates the product catalog number when accompanied by a number
	Indicates the product lot number when accompanied by a number
	Indicates the manufacturer of the product when accompanied by the name of the manufacturer.
	Indicates the product was manufactured in accordance with Annex II of the European Council Directive 93/42/EEC, as amended.
	Indicates a sterile fluid pathway using ethylene oxide (EtO)
	Indicates that the user must read the instructions for use before application
	Not for Reuse: Tubing sets manufactured by Terumo BCT bearing the not for reuse symbol are intended for single use only and are not intended to be reprocessed or reused. Terumo BCT cannot ensure the functionality or sterility of the set if it is reprocessed or reused.

	Indicates that the product is non-pyrogenic.
Rx Only	Indicates that the use of product is by prescription only
	Indicates that the product must be kept away from sunlight
	Indicates that the product must be kept dry
	Indicates that the product contains phthalates, specifically Di(2-ethylhexyl)phthalate (DEHP)
	Indicates that the product is made from High Density Polyethylene.
	Indicates that the product packaging complies with European Directive 94/62/EC for packaging and packaging waste.

Labels on the Platelet, plasma set

Case label

Spectra Optia® Platelet, Plasma Sets
 Kits collecte Plaquettes + Plasma Spectra Optia®
 Circuit per plaquettes e plasma Spectra Optia®
 Spectra Optia® Trombozyten-, Plasma-Sets
 Equipos Spectra Optia® para plaquetas y plasma
 Conjuntos de plasma e plaquetas Spectra Optia®
 Sedy Spectra Optia® pro trombocyty a plazmu
 Spectra Optia® trombocyt, plasmaset
 Spectra Optia® -trombocytyl, plasmaset
 Zestawy Spectra Optia® do pobierania krwi i osocza do systemu Spectra Optia®
 Колекторы Spectra Optia® для тромбоцитов и плазмы
 Spectra Optia® Trombocyt- och plasmaset
 Spectra Optia® Trombosit, Plasma Setten


Sterile and non-pyrogenic fluid pathway. Do not use the set if the end caps are not in place. Read instructions for use before application. Sterilized with ethylene oxide. Intended for single use only. Do not re-sterilize.

EXPIRY DATE  2016-05-01

MANUFACTURE DATE  2014-05-01

CATALOG NO. **REF** 10400

LOT NO. 05W1234

CE 0086

TERUMOBCT
10011 W. Collins Ave
 Laboratory, CO 80215
 MADE IN USA

(01)3502058 3104007(17)100501(30)5(10)05W1234
 TPN - 777815-168

Set label – plasma label and platelet labels.

Spectra Optia® Platelet, Plasma Set
 Kit de collecte Plaquettes + Plasma Spectra Optia®
 Circuito per piastrine e plasma Spectra Optia®
 Spectra Optia® Thrombozyten + Plasma-Set
 Equipo Spectra Optia® para plaquetas y plasma
 Conjunto de plasma e plaquetas Spectra Optia®
 Set Spectra Optia® pro trombocyty a plazmu
 Spectra Optia® trombocyt, plasmaset
 Spectra Optia® -trombozytti, plasmasetti
 Zestaw do pobierania płytek krwi i osocza do systemu Spectra Optia®
 Комплект Spectra Optia для тромбоцитов и плазмы
 Spectra Optia® Trombocyt- och plasmaset
 Spectra Optia® Trombositi, Plasma Seti

1

Sterile and non-pyrogenic fluid pathway. Do not use the set if the end caps are not in place. Read instructions for use before application. Sterilized with ethylene oxide. Intended for single use only. Do not re-sterilize.

EXPIRY DATE 2016-05-01

MANUFACTURE DATE 2014-05-01

CATALOG NO. REF 10400

LOT NO. 05W1234

TERUMOBCT

10011 W. Collins Ave
 Lakewood, CO 80215
 MADE IN USA

(01)05020583104006

(1 7)1 60501 (10)05W1234
 TPN - 777815-187

2016-05-01

+

TPN - 777815-127

REF 10400


LOT 05W1234

BC00010400

00005W1234

16 INFORMATION OF PACKING AND COMPLETENESS OF DELIVERY

The platelet and plasma disposable is packed by 6 in a box of corrugated fiberboard. Each set is packaged in a disposable tray covered with a disposable tray lid and contain Instruction for use.

Attention! Do not use the Set if its packaging is damaged. The sterility of the Set may be impaired.

17 MATERIALS OF WHICH THE MEDICAL DEVICE OR ITS MAIN COMPONENTS ARE MADE

Part of device	Contents of parts of device	Type of material	Grade (Trade name) of material	Manufacturer
1. Donor line	Luer connector (male)	Polycarbonate	2558-550115 clear	Bayer
	Tubing	PVC-DEHP EP 3.1.1.2	2222-NK-K68	AlphaGARY

	Needle	Stainless steel	SUS304 (C - 0,08%, Cr - 19%, Fe - 68.85%, Mn 1,5%, Ni - 9,5%, P - 0,045%, S - 0,025, Cu - 1%)	Exel International Inc
2. Sample bag	Sample system	PVC- DEHP 3.1.1.1.	30-2116FT	RENOLIT (SOLVA Y- DRAKA)
	Luer adapter vacutainer	PVC-DEHP EP 3.1.1.2	2222-NK-K75	AlphaGARY
3 AC/draw/return manifold	Manifold 3-1	COPOLYESTER	PCTG DN003	EASTMAN
	Injection Plug	SEBS	C-FLEX 35A R70-028	CPT
	3/16 Y connector	POLYESTER	PETG MN211	EASTMAN
	clear tube	PVC - DEHP EP 3.1.1.2	unichem 6711-00kel	COLORITE
4 Anticoagulant (AC) line	Luer connector (male)	Polycarbonate	2558-550115 clr	bayer
	Spike (orange)	PVC-DEHP EP 3.1.1.2	VJ 373/1 CLEAR 124	IENOS (HYDRO POLYMERS)
		ABS	GUARDIAN (LUSTRA N) ABS 348-201348	IENOS (BAYER)
	Tubing	PVC-DEHP EP 3.1.1.2	ALPHA 2222-NK-K68	AlphaGARY

5 Sterile barrier filter	Sterile barrier filter (0.2 microns)	ACRYLIC	CRYOLITE G-20	GELMAN/PALL
		POLYETHERSULFONE	SUPOR-200	
6. Cassette		Polyethylene Terephthalate (PET- G)	PETG MN21	EASTMAN
7 AC pump header	Pump Header tubing	PVC-DEHP EP 3.1.1.2	6711-00KEL	COLORITE
8 Draw/return pressure sensor		POLYESTER	PETG MN211	EASTMAN
		SILICONE	SILASTIC Q7-4840	DOW CORNING
		SILICONE	1950-40 CLEAR	SHIN-ETSU
9. AC line	Clear tubing	PVC-DEHP EP 3.1.1.2	6711-00KEL	COLORITE
10 Return line	Blue stripe tubing	PVC-DEHP EP 3.1.1.2	TEKNOR APEX 89-P057A-NT BLUE #1236	COLORITE
11 Draw line	Red stripe tubing	PVC-DEHP EP 3.1.1.2	TEKNOR APEX 89-P057A-NT BLUE	COLORITE
12 Inlet pump header	Pump Header tubing	PVC-DEHP EP 3.1.1.2	6711-00KEL	COLORITE
13 Return reservoir	Return reservoir	Polyethylene Terephthalate (PET- G)	PETG MN21	EASTMAN
	Return filter	POLYESTER	PE 200 HC	ZBF
		PVC-DEHP EP 3.1.1.2	7811BG015	COLORITE
14 Return pump header		PVC-DEHP EP 3.1.1.2	6711-00KEL	COLORITE
15 Centrifuge pressure sensor		POLYESTER	PETG MN211	EASTMAN

		SILICONE	SILAST IC Q7- 4840	DOW CORNING
		SILICONE	1950-40 CLEAR	SHIN-ETSU
16 Inlet line to centrifuge	red tubing	PVC-DEHP EP 3.1.1.2	2222-nk- k75 w/apx89- pO57a- ntred 889 4%	alpha/teknor
17 RBC line from centrifuge	Clear tubing	PVC-DEHP EP 3.1.1.2	2222-nk- k75	AlphaGARY
18 Platelet line from centrifuge	Clear tubing	PVC-DEHP EP 3.1.1.2	2222-nk- k75	AlphaGARY
19 Plasma line from centrifuge	yellow tubing	PVC-DEHP EP 3.1.1.2	2222-nk- k75 w/apx89- pO57a- nt	alpha/teknor
20 Plasma pump header	Pump Header tubing	PVC- DEHP EP 3.1.1.2	6711- 00KEL	COLORITE
21 Cassette label		Paper white	TransTherm 0TTFP	Avery Dennison
		Adhesive material – Glue	758M	Green Bay
		Film – Wax/Resin	ITW Thermal B128165EOS	ITW Thermal Films
22 Collect pump header	Pump Header tubing	PVC- DEHP EP 3.1.1.2	UNICHEM 6711 B or A w/ 0.5% TEGOME R H-Si6440	COLORITE
23 Centrifuge loop	Four-lumen tubing	PVC-DEGF EP 3.1.1.2 -	ALPHA 2222-NK- K75	AlphaGARY
	Sleeves	Polyesterpolycaprolactone copolymer white	Merquinsa D15N85Uv	Merquinsa

	Collars	Acetate copolymer white -	Acetal Celanese Celcon LW90	Celanese
	Bearings	Acetate copolymer white	Acetal Celanese Celcon LW90	Celanese
24 Channel		PVC-DEHP EP 3.1.1.2.	0-4139 (EC 93- 710)	RENOLIT
25 Inlet port		PVC-Rigid	2212C2B /H-118 CLEAR 003	ALPHAGAR Y
26 Collection chamber	Collection chamber	PVC-Rigid	2212C2B /H-118 CLEAR 003	ALPHAGAR Y
	Control "Y"	PVC-Rigid	2212C2B /H-118 CLEAR 004	ALPHAGAR Y
	Tubing, control, single stage blue stripe	PVC-DEHP EP 3.1.1.2.	2222-NK- K75 Blue Stripe: Tecnor APEX 89-	ALPHAGAR Y
	Clear tubing	PVC-DEHP EP 3.1.1.2	2222-NK- K75	AlphaGARY
27 LRS chamber		POLYESTER	PETG MN211	EASTMAN
28 Plasma collect line	Tubing 0.158 OD x 0.022 wall	PVC-DEHP EP 3.1.1.2	2222-NK- K75	AlphaGARY
	Tubing 0.113 x .037 wall	PVC-DEHP EP 3.1.1.2	2222-NK- K68	AlphaGARY
29 Plasma bag	Plasma bag	PVC plasticized with DEHP	ES- 2967 ZPH (ES300 0)	RENOLIT

	Tubing 0.158 OD x 0.022 wall	PVC-DEHP EP 3.1.1.2	2222-NK-K75	AlphaGARY
	0.158" ID bond socket	PVC-DEHP EP 3.1.1.2	2222-NK-K75	ALPHAGARY
	Spike receptor	PVC-DEHP EP 3.1.1.1	3250	RENOLIT
30 Platelet collect line	Tubing 0.113 x .037wall	PVC-DEHP EP 3.1.1.2	2222-NK-K68	AlphaGARY
	Tubing 0.158 OD x 0.022 wall	PVC-DEHP EP 3.1.1.2	2222-NK-K75	AlphaGARY
	Reducer 0.187 to 0.158 00	POLYESTER	PETG MN211	EASTMAN
	"Y" connector, offset 0.158	POLYESTER	MN211	EASTMAN
31 Platelet bag	Platelet bag film	PVC-CITRATE NO DEHP	ELLA Y ES 3049	RENOLIT
	Platelet bag tubing	PVC-DEHP EP 3.1.1.2	2222-NK-K75	KELCOURT
	Spike receptor	PVC-DEHP EP 3.1.1.1	3250	RENOLIT
32 Platelet product sampler	Sample bulb	PVC-DEHP EP 3.1.1.2	WELVIC WE 9075	EVC
33 Vent bag line	Tubing 0.158 OD x 0.022 wall	PVC-DEHP EP 3.1.1.2	2222-NK-K75	AlphaGARY
34 Vent bag	Vent bag, film	PVC-DEHP EP 3.1.1.1 & EP 3.1.1.2	30-2116 FT	RENOLIT

18 LIST OF NATIONAL AND INTERNATIONAL REGULATORY DOCUMENTS/STANDARDS, THE MEDICAL DEVICE COMPLIES WITH

- ISO 13485:2011, Medical Devices – Quality Management Systems – Requirements for Regulatory Purposes

- ISO 1135-4, 2012, (section 5.4) Transfusion Equipment for Medical Use -- Part 4: Transfusion Sets for Single Use
 - ISO 3826-1:2013, Plastics Collapsible Containers for Human Blood and Blood Components — Part 1: Conventional Containers
 - ISO 3826-2: 2008, Plastics Collapsible Containers for Human Blood and Blood Components — Part 2: Graphical Symbols for Use on Labels and Instruction Leaflets
 - EN ISO 14155: 2011, Clinical Investigation of Medical Devices for Human Subjects Good Clinical Practice
 - ISO 594-1:1986, Conical Fittings with a 6% (Luer) Taper for Syringes, Needles and Certain Other Medical Equipment – Part 1: General Requirements
 - ISO 8536-4: 2010, (section 6.4) Infusion Equipment for Medical Use -- Part 4: Infusion Sets for Single use, Gravity Feed
 - ISO 14971:2012, Medical Devices – Application of Risk Management to Medical Devices
 - ISO/IEC 62366:2015 Medical Devices – Application of Usability Engineering to Medical Devices
 - ISO 23908: 2013, Sharps injury protection — Requirements and test methods — Sharps protection features for single-use hypodermic needles, introducers for catheters and needles used for blood sampling Labeling
 - EN 1041:2008, Information Supplied by the Manufacturer with Medical Devices
 - EN ISO 15223-1: 2014, Medical devices —Symbols to be used with medical device labels, labeling, and information to be supplied —Part 1: General requirements
 - EN ISO 15986: 2011, Symbol for use in the labeling of medical devices - requirements for labeling of medical devices containing phthalate
 - EN ISO 10993-1:2011, Overall Biological Evaluation as Part of Risk Management
 - EN ISO 10993-2:2009, Animal Welfare
 - EN ISO 10993-3:2014, Genotoxicity (Reaffirmed by AMMI in 2009)
 - EN ISO 10993-4:2011, Biological evaluation of Medical Devices – Part 4: Selection of tests for interactions with Blood.
 - EN ISO 10993-5:2011, Cytotoxicity
 - EN ISO 10993-7: 2011, EtO Residuals
 - EN ISO 10993-10: 2011, Irritation and Skin sensitization
 - EN ISO 10993-11:2011, Systemic Toxicity
 - EN ISO 10993-12: 2015, Sample Preparation
 - EN ISO 10993-17: 2011, Establishment of Allowable Limits for Leachable Substances
 - EN ISO 10993-18: 2011, Chemical Characterization
- ISO 11135-1:2007, Sterilization of Health Care Products – Ethylene Oxide – Part 1: Requirements for Development, Validation and Routine Control of a Sterilization Process for Medical Devices

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Sterile and non-pyrogenic fluid pathway. Do not use the set if the end caps are not in place. Read instructions for use before application. Sterilized with ethylene oxide. Intended for single use only. Do not resterilize.

EXPIRY DATE  2016-05-01

MANUFACTURE DATE  2014-05-01

CATALOG NO. **REF** 10400

LOT NO. 05W1234

CE 0086

TERUMOBCT
 10811 W. Collins Ave
 Lakewood, CO 80215
 MADE IN USA

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Sterile and non-pyrogenic fluid pathway. Do not use the set if the end caps are not in place. Read instructions for use before application. Sterilized with ethylene oxide. Intended for single use only. Do not resterilize.

EXPIRY DATE  2016-05-01

MANUFACTURE DATE  2014-05-01

CATALOG NO. **REF** 10400

LOT NO. 05W1234

 0001


 (01)05020583104006

 (17)160501(10)05W1234
 TPN - 777815-167

   0086

 **TERUMOBCT**
 10811 W. Collins Ave.
 Lakewood, CO 80215
 MADE IN USA

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