State of Minnesota

SECRETARY OF STATE

APOSTILLE

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1. Country: United States of America

This Public Document: Documents for 3M Unitek

2. has been signed by Lori J Sheetz

3. acting in the capacity of Notary Public, State of Minnesota

4. bears the seal/stamp of Lori J Sheetz, Notary Public, State of Minnesota

CERTIFIED

5. at: St. Paul, Minnesota

6. Date: April 24, 2017

7. by: Secretary of State, State of Minnesota

8. File No: 9473893-3

9. Seal/Stamp:

10. Signature:

Steve Simon

Steve Simon
Secretary of State

To Whom It May Concern:

In my capacity as Regulatory Affairs Manager, 3M Unitek, I certify that the attached English version of the Instruction for Use (IFU) for Orthodontics Hand Instruments intended for Russia is truthful and accurate.

Sincerely,


Mark Burville
Regulatory Affairs Manager
3M Unitek

21 APR 2017
Date

STATE OF MINNESOTA)
COUNTY OF RAMSEY) ss

On April 21, 2017, before me, a Notary Public within and for said county, personally appeared Mark Burville, known by me to be the person described in and who executed the foregoing instrument and acknowledged that he executed the same as his own free act and deed.

Lori J. Sheetz, Notary Public
My Commission Expires, January 31, 2020

3M Unitek Orthodontics Hand Instruments

Instructions for Use

Intended use

The instruments are designed for the positioning, disengagement, and debonding of devices for orthodontic applications.

Indication(s):

- Cutters are used to cut various sizes of orthodontic wires, pins, and auxiliaries.
- Pliers are used to adjust, place, bend, loop and twist wires, bands, auxiliaries, ligatures, and appliances.
- SmartClip Appliance Wire Disengagement Instrument is designed for removing archwires.
- Unitek Self-Ligating Bracket Debonding Instrument is designed to debond 3M Unitek Self-Ligating brackets.
- Unitek Bracket Debonding Instrument debonds Clarity brackets with a squeeze.
- Unitek Bracket Positioning Gauge is used to measure from bracket slot to incisal/occlusal surface of tooth.
- Unitek Bracket Placement Instrument is used to pick up Self-Ligating Brackets and Standard Twin Brackets and place the bracket onto the tooth.

Scope of use: Orthodontia. For professional use only

Contraindications

No known contraindications

Warning

1. Use instruments for indicated use only. Failure to do so may result in breakage of the instrument and injury to patient or orthodontic professional. Safety glasses recommended for patient and orthodontic staff when using these instruments.
2. Before cutting wire in the mouth with any non-safety hold cutter, place a gauze pad or cotton roll next to the wire end to catch wire fragments and to prevent injury to patient.
3. Stainless steel contains nickel and chromium. A small percentage of the population is known to be allergic to these metals.

If an allergic reaction occurs, direct patient to consult a physician.

Safety glasses recommended when using these instruments.

Authorized Service Center: Precision Plier Service. All instruments should be sent to the 3M Unitek Authorized Service Center for evaluation to determine credit, replacement, or repair. The Authorized Service Center will determine, upon receipt of instrument, if the instrument is to be covered by warranty. Materials, workmanship, and corrosion resistance are the only items covered by the warranty. Misuse or abuse will void warranty. Customer will bear any charges regarding repair of out of warranty instruments.

For technical assistance, call the 3M Unitek Technical Helpline at 1-800-265-1943, 1-626-574-4577.

For regular maintenance (cleaning, re-sharpening, re-conditioning), please send instruments to the Precision Plier Service.

Disinfection Information

These instruments are not supplied in a sterile state. Cleaning, disinfection and sterilization is required for the device's first use after removal of the protective packaging and prior to each use.

Effective cleaning and disinfection is an indispensable requirement for proper device sterilization. Three cleaning and disinfection procedures are provided; an automatic procedure using a thermal disinfectant unit (washer-disinfectant), as well as two manual procedures (ultrasonic cleaning/disinfection and manual cleaning/disinfection).

The user is responsible for the sterility of the devices. Ensure that only validated procedures are used for cleaning, disinfection and sterilization. Sterilization equipment must also be maintained and checked regularly, as well as the validated parameters for the cleaning, disinfection and sterilization cycles.

Cautions

- Do not use a metal brush or steel wool to clean instruments, as this will damage the finish of the instrument.
- Do not apply mineral oil or silicone lubricants to instrument hinges because they introduce microorganisms, prevent penetration of the steam to the instrument's surface, and are difficult to remove.
- Do not expose devices, cassettes, trays or sterilization containers to temperatures higher than 141°C / 286°F. Exposure to higher temperatures is the responsibility of the user.
- Do not use detergents or disinfectants containing the following substances:

Do not use detergents or disinfectants containing the following substances:

- (a) Strong alkalines (>pH9)
- (b) Strong acids (<pH4)
- (c) Phenols or iodophors
- (d) Hydrogen peroxide (H₂O₂)
- (e) Interhalogenic agents / halogenic hydrocarbons / iodophors
- (f) Strong oxidizing agents / peroxides
- (g) Organic solvents

- Water quality may influence the result of the cleaning and disinfection of the devices. Therefore, use only sterile or low contaminated water (maximum 10 germs/ml, maximum 0.25 endotoxin units/ml), for example, purified water/highly-purified water (PW/HPW according to the pharmacopeias) for post-rinsing, as well as for other processing steps that require water. Corrosion of instruments may be caused by high contents of chloride or other minerals in tap water. If stains and corrosion occur and other reasons can be excluded, it may be necessary to test the quality of the tap water in your area. Most water quality problems can be avoided with the use of deionized or distilled water.

- Please refer to subsection 12 for additional instructions for specific instruments / devices.

Limitations on Reprocessing

3M Unitek warrants that these instruments will be free from defects in materials, workmanship or separation of brazed cutting edges for a period of 1 year from the date of purchase. 3M Unitek further warrants that this instrument will be resistant to corrosion for 5 years from the date of purchase. Failure to follow recommended use, maintenance and sterilization procedures will void all warranty claims. Recommended maintenance includes cutter tip blade re-sharpening and reconditioning (cleaning, polishing, adjusting screw as necessary) every 6 months.

The user is responsible for inspecting instruments prior to each use.

Preparation at Point of Use

- Handle used and contaminated devices with protective gloves that fulfill the requirements of Directive 89/686/EEC.
- Remove excess soil with disposable cloth/paper wipes.
- Avoid drying and fixation of tissue debris onto the surface of devices.

Containment & Transportation

- Contain contaminated devices in an approved sealed container during transport from the point of use to the decontamination area.
- Reprocess devices as soon as is reasonably practical following use.

Preparation

Selection of cleaning detergents – consider the following in the selection of detergents:

- Suitability of the cleaning detergent for cleaning devices.
- Compatibility of cleaning detergent with the devices (see section 1, Cautions).
- Compatibility of cleaning detergent with thermal disinfectant unit (washer-disinfectant) or with ultrasonic treatment.
- Compatibility of the cleaning detergent with the disinfectant.
- Instructions of the cleaning detergent manufacturer regarding concentration and soaking time.

Selection of disinfectant – Consider the following in the selection of disinfectants:

- Suitability of the disinfectant for disinfection of devices.
- Compatibility of the disinfectant with the devices (see section 1, Cautions).
- Compatibility of the disinfectant with the cleaning detergent.
- Approved efficiency of disinfectant, i.e., DGHM/VAH, RKI approval or CE marking or at least intermediate level disinfectant.
- Instructions of the disinfectant manufacturer regarding concentration and soaking time.

Preparation (required for all three cleaning and disinfection procedures listed in section 6)

- Remove coarse impurities on the devices immediately after application (within a maximum of 2 hours).
- Use a soft brush for manual removal of coarse impurities—do not use metal brushes or steel wool.
- For hinged instruments, open and close the hinge at least 5 times with the instrument submerged in an enzymatic cleaner or disinfectant solution.
- For instruments with blind holes (i.e., a hole that does not go through the material), rinse the blind hole intensively with a syringe (use a syringe with a needle if the size of the hole is too small to adequately rinse with just the syringe).
- The disinfectant should be aldehyde-free (otherwise fixation of blood impurities could occur), possess an approved efficiency (i.e., DGHM, RKI approval or CE marking or at least intermediate level disinfectant), be suitable for disinfection of instruments and be compatible with the instruments (see section 1, Cautions).
- Follow the instructions of the manufacturers for the cleaning detergents and the disinfectants.

Cleaning and Disinfection

Method 1 – Automated Cleaning & Disinfection

Equipment: Thermal Disinfectant Unit (Washer-Disinfectant)

Consider the following when using a thermal disinfectant unit (washer-disinfectant):

- Efficiency of the thermal disinfectant unit (washer-disinfectant) (e.g. EN ISO 15883-1, HTM 2030 if required by local regulations and guidelines).
- Thermal disinfection program that is approved [e.g. according to EN/ISO 15883-1 with A_0 value > 3000 (for example to be established by application of 90°C / 194°F for 5 minutes) or by application of formerly common programs with “overkill” (for example: 93°C / 200°F for 10 minutes)] and suitability for devices as well as sufficient rinsing steps.[The application of A_0 values > 600 (as required in some European countries) is not recommended. If an A_0 value > 600 is applied, the user is responsible for demonstrating sufficient efficiency of the cleaning and disinfection process by performance of a product, process, and device specific validation.]

- Thoroughly post-rinse, preferably with high quality water, e.g., PW/HPW (see section 1, Cautions).

Washer-Disinfector Procedure

1. Disassemble instruments if applicable (not required for hinged instruments) or open hinged instruments completely.
2. Place the devices in the thermal disinfector unit (washer-disinfector), assuring no contact between the instruments.
3. Initiate the cycle.
4. Remove the devices from the thermal disinfector unit (washer-disinfector) after the end of cycle.

The fundamental suitability of the devices for an effective automated cleaning and disinfection was demonstrated by an independent accredited test laboratory by application of the disinfector G 7836 GD, Miele & Cie, GmbH & Co., Gütersloh, Germany, (thermal disinfection) and the cleaning detergent Neodisher medizym (Dr. Weigert GmbH & Co. KG, Hamburg, Germany) considering the specified procedure.

Method 2 – Ultrasonic Cleaning & Disinfection

Equipment: Ultrasonic Cleaning Unit.

Cleaning Procedure:

1. Disassemble instruments if applicable (not required for hinged instruments) or open hinged instruments completely.
2. Completely submerge devices in the cleaning solution for the time that the manufacturer of the cleaning detergent recommends. Note: There should be no contact between the instruments; do not overload the ultrasonic cleaning unit.
3. Use the processing time recommended by the manufacturer of the ultrasonic unit. If using a cassette system, the ultrasonic cleaning time has to be at least 16 minutes, unless the manufacturer of the cleaning detergent requires a longer exposure time.
4. Remove devices from the cleaning solution and post-rinse thoroughly at least 3 times for at least one-half minute each time with high quality water, e.g., PW/HPW (see section 1, Cautions).
5. Inspect devices for good cleaning result and repeat procedure if necessary.

Disinfection Procedure:

1. Completely submerge devices in the disinfectant solution and soak for the recommended soaking time of the disinfectant manufacturer. There should be no contact between the devices.
2. Remove devices from the disinfectant solution and post-rinse them thoroughly at least 5 times for at least one-half minute each time with high quality water, e.g., PW/HPW (see section 1, Cautions).

The fundamental suitability of the devices for an effective cleaning and disinfection was demonstrated by an independent accredited test laboratory by application of the ultrasonic bath SONOREX SUPER RK 514 H (BANDELIN electronic GmbH & Co. KG, Berlin, Germany) as well as of the cleaning detergent Cidezyme/Enzol and the disinfectant Cidex OPA (ASP, Irvine, USA / Johnson & Johnson GmbH, Norderstedt, Germany) considering the specified procedure.

Method 3 – Manual Cleaning and Disinfection

Equipment: Not applicable.

Cleaning procedure:

1. Disassemble instruments if applicable (not required for hinged instruments) or open hinged instruments completely.
2. Completely submerge devices in the cleaning solution and soak for the recommended soaking time of the cleaning detergent manufacturer. If necessary, use a soft brush to carefully brush instrument.
3. Remove devices from cleaning solution and rinse thoroughly at least 3 times for at least one-half minute each time with high quality water, e.g., PW/HPW (see section 1, Cautions).
4. Inspect the devices for good cleaning result and repeat procedure if necessary.

Disinfection procedure:

1. Completely submerge the devices in the disinfectant solution and soak for the recommended soaking time of the disinfectant manufacturer. There should be no contact between the devices.
 2. Remove devices from the disinfectant solution and post-rinse thoroughly at least 5 times for at least one-half minute each time with high quality water, e.g., PW/HPW (see section 1, Cautions).
- The fundamental suitability of the instruments for an effective cleaning and disinfection was demonstrated by an independent accredit test laboratory by application of the cleaning detergent Cidezyme/Enzol and the disinfectant Cidex OPA (ASP, Irvine, USA / Johnson & Johnson GmbH, Norderstedt, Germany) considering the specified procedure.*

Drying

- Allow post-drying step in a clean place, or use filtered air for drying to prevent recontamination.
- If drying is achieved as part of a washer-disinfector cycle, do not exceed 120°C / 248°F.
- Completely dry devices before packaging.

Inspection, Maintenance, Testing

- Inspect devices for impurities, damage and wear. Instrument cutting edges should be free of nicks. If devices are still dirty, repeat cleaning and disinfection procedures.
- Reassemble disassembled instruments. Check for smooth movement of hinge without excessive "play". Locking (ratchet) mechanisms should be checked for action.
- Lubricate hinges slightly with sterilization-capable, vapor-permeable lubricant (e.g., a paraffinic white oil) suitable for steam sterilization with manufacturer's certification that the lubricant will not affect the sterilization efficiency of the oiled instruments, and that the biocompatibility of the lubricant will be maintained during the sterilization cycle.

Packaging

Consider the following in the selection of suitable sterilization containers:

- Conformity with ANSI AAMI (EN) ISO 11607 (EN 868-8)
 - Suitable for steam sterilization (temperature resistance up to 141°C (286°F))
 - Sufficient protection of the instruments and the sterilization packaging against mechanical damage
 - Sterilization container - regular maintenance according to the manufacturer's instructions
- Alternatively, standard sterilization packaging (e.g., paper/foil, single or double packaging) can be used. Consider the following during the selection of suitable sterilization packaging:
- Conformity with ANSI AAMI (EN) ISO 11607 (EN 868-2...)
 - Suitable for steam sterilization (temperature resistance up to 141°C (286°F))
 - Sufficient protection of the instruments and the sterilization packaging against mechanical damage

Sterilization

Use only the procedure recommended below. Other sterilization procedures are the responsibility of the user.

- Use a fractionated vacuum autoclave or a gravity displacement autoclave (with sufficient product drying). If possible, a fractionated vacuum procedure should be used for sterilization. A gravity procedure should be used only if a fractionated vacuum autoclave is not available; in this case, the significantly lower efficiency of a gravity displacement autoclave has to be considered.
- The sterilizer should be validated according to ANSI AAMI (EN) ISO 17665 (and HTM 2010 if required by local regulations and guidelines).
- Do not use the flash sterilization procedure.
- Do not exceed a sterilization temperature of 138°C / 280°F.

- Do not exceed sterilizer's maximum load when sterilizing multiple devices in one autoclave cycle
- Steam sterilize according to EN 13060, and EN 285 if required by local regulations and guidelines.
- Use one of the following sterilization exposure time at the sterilization temperature: See section 12, Special Procedures, for additional sterilization exposure time and temperature for selected instruments.
- Fractionated vacuum autoclave
 - three minutes at 132°C / 134°C (270°F / 273°F), or
 - twenty minutes at 121°C (250°F)
- Gravity displacement autoclave
 - five minutes at 132°C / 134°C (270°F / 273°F), or
 - twenty minutes at 121°C (250°F)

The fundamental suitability of the devices for an effective steam sterilization was demonstrated by an independent accredited test laboratory by application of the steam sterilizer Systec V-150 (Systec GmbH Labor-Sytemtechnik, Wettenberg) and of the fractionated vacuum procedure as well as of the gravity procedure. For this, typical conditions in the clinic and doctor's practice as well as the specified procedure were considered.

Storage

Store devices after sterilization in a dry and dust-free place. Sterilization can only be maintained if the devices remain packaged or wrapped, impermeable to microorganisms, following a validated sterilization. The status of the sterilization has to be clearly indicated on the wrapped packages or the containers. For safety reasons, keep sterile and non-sterile devices separated.

Storage Conditions

Unitek Orthodontics Hand Instruments should be stored at room temperature in a clean, ventilated room without corrosive air.

Transportation conditions

Values of temperature, humidity and pressure are not regulated.

Conditions of use

This product is designed to be used at room temperature. The humidity and pressure are not regulated.

Special Procedures

Hinged Instruments	Open and close instrument several times during cleaning process and disinfection processes. Lubricate hinge before sterilization and process in an open position during sterilization.
Hinged Instruments REF 804-171, Unitek Bracket Placement Instrument	Use <u>only the following to clean/disinfect these devices</u> . Method 1, Automated Cleaning & Disinfection, or Method 3, Manual Cleaning & disinfection
REF 804-160, SmartClip Appliance Wire Disengagement Instrument	Remove plastic slide and relax the lateral screw to enable full opening of the instrument prior to pretreatment, cleaning and disinfection (only packaging and sterilization in the mounted state). After sterilization, tighten the lateral screw with the allen wrench that is supplied with the instrument. Minimum number of sterilization cycles for plastic slide is 20 cycles. Replacement plastic slides are available (REF 804-161).
Prestige Series Instruments	Prior to use, use a lubricant such as DSX Medic 365 and work

	the joints gently open and close for a few minutes. Proceed with sterilization instructions according to the manufacturer's instructions. Lubricate again after instruments have been cooled and dried.
--	---

Warranty Information

3M Unitek warrants that this instrument will be free from defects in materials, workmanship or separation of brazed cutting edges for a period of 1 year from the date of purchase. 3M Unitek further warrants that this instrument will be resistant to corrosion for 5 years from the date of purchase. Failure to follow recommended use, maintenance and sterilization procedures will void all warranty claims.

Recommended maintenance includes cutter tip/blade re-sharpening and reconditioning (cleaning, polishing, adjusting screw as necessary) every six months.

3M Unitek's sole obligation and customer's sole remedy in the event of a claimed defect shall be limited to, at 3M Unitek's option, repair, replacement of merchandise, or refund of the purchase price. ALL IMPLIED WARRANTIES, INCLUDING IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE OR USE ARE DISCLAIMED.

Limitation of Liability

Except where prohibited by law, 3M Unitek will not be liable for any loss or damages arising from this Product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability. All instruments should be sent to the 3M Unitek Authorized Service Center at the address below for evaluation to determine credit, replacement, or repair. The Authorized Service Center will determine, upon receipt of instrument, if the instrument is to be covered by warranty. Materials, workmanship, and corrosion resistance are the only items covered by this warranty. Misuse or abuse will void warranty. Customer will bear any charges regarding repair of out of warranty instruments.

Materials the medical devices are made of

Hard Wire Cutter - Stainless Steel BG42, 410, 430

Universal Cutter Long Handle - Stainless Steel BG42, 410, 430

Slimnose Weingart Pliers - Stainless Steel 425M or DIN 1.4116, 430

SmartClip Appliance Wire Disengagement Hand Instrument - Plastic and Stainless Steel 430F or 303, 425M or 420, 18-8

Unitek Self-Ligating Bracket Debonding Instrument - Nickel-Titanium (Nitinol) Alloy and Stainless Steel 32Cr13Mo GB/T 20878-2007, 430F or 303

Unitek Debonding Instrument - Stainless Steel 425M or DIN 1.4116, 430

Unitek Bracket Positioning Gauges - Stainless Steel 17-7PH, 18-8, and Brass

Unitek Bracket Placement Instrument - Stainless Steel 420 or DIN 1.4024

Labelling of the medical devices contains the following information:

Products are supplied in a package with labelling containing the following information: product name and reference number; name and address of the manufacturer; contents; batch number; information about trade mark; and symbols – CE conformity mark; Attention! Refer to IFU.

Description of Packaging

Primary package: Moulding Foam (Plastic Case or Poly Bag). Secondary package: Paperboard Box and/or Poly Bag.

Shelf Life

There is no expiration date for the Unitek Orthodontics Hand Instruments.

Disposal and recycling

Packaging/contents should be disposed in accordance with the local/regional/national/international regulations.

Medical devices should be disposed by the local public law organizations in accordance with SanPiN 2.1.7.2790-10 "Sanitary and Epidemiological Requirements for Medical Waste Treatment" waste Class B).

Legal manufacturer

3M Unitek Corporation, USA

2724 South Peck Road, Monrovia, California, 91016, USA

Authorized representative on the territory of the Russian Federation:

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