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Möller Medical GmbH

Wasserkuppenstr. 29-31, D-36043 Fulda · Postfach 1737, D-36007 Fulda
Tel.: +49 (0) 661-94195-0 · Fax: +49 (0) 661-94195-850
E-Mail: info@moeller-medical.com · Web: www.moeller-medical.com

MöllerMedical

Peter Schrempp, General Manager

MöllerMedical
GmbH

D-36043 Fulda
Wasserkuppenstr. 29-31
Tel.: +49 (0) 661-94195-0
Fax: +49 (0) 661-94195-850
www.moeller-medical.com

04.09.2017

Date

INSTRUCTION

for use of the medical device

“NEEDLE KITS FOR BONE CEMENT”

manufactured by Möller Medical GmbH, Germany

TABLE OF CONTENTS

GENERAL INFORMATION	3
PRESCRIPTION, INDICATIONS FOR USE	4
CONTRAINDICATIONS, SIDE EFFECTS, COMPLICATIONS	4
PRECAUTIONS	5
PREPARATION FOR OPERATION AND OPERATION PROCEDURE	6
RULES OF USE AND UTILIZATION	12
TRANSPORTATION AND STORAGE	12
PACKAGING AND LABELING	12
DELIVERY PACKAGE	15
10. GENERAL TECHNICAL CHARACTERISTICS	18
11. REQUIREMENTS FOR MAINTENANCE AND REPAIR	36
12. REQUIREMENTS FOR ENVIRONMENTAL PROTECTION	36
13. WARRANTY, POSTMARKETING SURVEILLANCE AND COMPLAINTS	36

I. GENERAL INFORMATION

The medical device "NEEDLE KITS FOR BONE CEMENT" includes kits and separately packaged needles of different types for necessary surgical techniques when performing vertebroplasty and biopsy.

The kit also includes wires with depth marks, trocars and stylets for easier biopsy and more precise injection of the bone cement into the vertebral body. There are square and bevelled end needles for vertebroplasty which allow to direct the bone cement in the desired way.

Vertebroplasty needle, gauge 8G, with a side-hole, a square end, a Luer taper

Vertebroplasty needle, gauge 10G, with a side-hole, a square end, a Luer taper

Biopsy needle for the Vertebroplasty Needle Kit, gauge 8G, with a Luer taper

Wire with depth marks

Cannulated trocar

Trocar

Injection needle with a Luer taper

2. PRESCRIPTION, INDICATIONS FOR USE

The medical device "Needle kits and separately packaged needles for bone cement" is intended for the percutaneous treatment of the vertebral body, injections of the PMMA cement in the vertebral bodies while vertebroplasty.

INDICATIONS.

Vertebroplasty Needle Kits and Biopsy Kits are used only for injections of the bone cement VertecemV+ manufactured by Synthes GmbH, Switzerland (not included in this medical device) into the vertebral body while vertebroplasty and bone biopsy.

The Bone Access Needle is used to provide access to the bone.

The Biopsy Needle is used for bone biopsy.

The medical device "Needle Kits for Bone Cement" is introduced into the vertebral body or the intervertebral space during surgical interventions along with simultaneous application of the medical imaging techniques (X-rays, computer tomography, bone scanning).

Biopsy Kits are intended for use only with Vertebroplasty Needle Kits.

Vertebroplasty Needle Kits and Biopsy Kits are used only for injections of the bone cement Vertecem V+ manufactured by Synthes GmbH, Switzerland (not included in this medical device).

The bone cement Vertecem V+ manufactured by Synthes GmbH, Switzerland is a radiopaque injectable cement based on polymethylmethacrylate (PMMA) used in case of spine pathologies in order to augment the spongy bone tissue. The bone cement belongs to self-polymerizing two-component systems (powder and liquid components), starting to polymerize at room temperature immediately after mixing.

Bone Access Needles are used along with Vertebroplasty Needle Kits. Biopsy Needles are used with Biopsy Kits manufactured by Möller Medical.

3. CONTRAINDICATIONS, SIDE EFFECTS, COMPLICATIONS

The medical device can not be used in patients with coagulation disorders or infections, or in patients for whom who vertebroplasty and biopsy are contraindicated.

There is no information about side effects.

Information about complications: possible complications are related to incorrect techniques of the application. If some PMMA is in the vertebral canal, it can lead to a compression of nerve structures resulting in a pain syndrome or weakness in the extremities. In this case, a re-surgery is required in order to remove the excessive cement from the vertebral canal.

Infections of the surgical wound caused by such procedures are rare cases, since the surgical procedures are performed through a small incision, and during polymerization the bone cement is heated up to 65-70°C

Information about complications: possible complications are related to incorrect techniques of the MD application. If some PMMA is in the vertebral canal, it can lead to a compression of nerve structures resulting in a pain syndrome or weakness in the extremities. In this case, a re-surgery is required in order to remove the excessive cement from the vertebral canal.

Infections of the surgical wound caused by such procedures are rare cases, since the surgical procedures are performed through a small incision, and during polymerization the bone cement is heated up to 65-70°C

4. PRECAUTIONS

Vertebroplasty Needle Kit is used in technically complex operations, so the medical specialist should know how to use the kit and other tools.

Vertebroplasty needle (for instance, of gauge 8G, with a side hole, a square end, a Luer taper) may be introduced into the bone by careful hammer tapping. It is necessary to avoid hard hammer blows even in case of a sclerotic or thickened cortical bone. Under certain circumstances, the risk of bone breakdown can develop.

Since the Bone Access Needle and the Biopsy Needle are used in technically complex operations, the medical specialist should know how to handle the kit, to use access techniques and therapeutic materials.

Safe use of the Bone Access Needle and the Biopsy Needle can not be guaranteed when used along with MRI that involves certain risks, including:

- heating or dislocation;
- artefacts in the MRI image.

5. PREPARATION FOR OPERATION AND OPERATION PROCEDURE

Place the patient in an inclined position on a lumbar support. The table should be radiolucent in all planes. Through an injection needle with a Luer taper the bone cement is introduced in the patient's body. Transpedicular (Fig.1) or parapedicular approach (Fig.2) is used to access the vertebral body:


Fig.1


Fig.2

A wide range of vertebroplasty needles of necessary gauge with a bevelled or square end and a Luer taper provides a choice of procedures. For optimal control of filling of the vertebral cavity with the cement, the needles with side and direct holes are used.

Direction of the cement introduction through a side hole (Fig.3):


Fig.3

Types of approaches:

1. Cement introduction with a wire with depth marks
2. Direct approach for the cement introduction

Biopsy method

method of approach with a wire with depth marks and the direct approach method described in the following sections are based on a minimally invasive method. In both cases, the access is provided by a combination of the following types: trans- or parapedicular approaches; uni- or bipedicular approaches.

The cement introduction with a wire with depth marks (Vertebroplasty Needle Kits) 1,2), a minimally invasive method with the transpedicular approach:

The C-type X-ray machine is set for a front projection, a precise place for incision is determined under fluoroscopy. Under control of the C-type X-ray machine, the skin is punctured. A wire with depth marks (hereinafter – the wire) is directed through soft tissues until the vertebral body surface (under X-ray control in the front plane). After reaching the bone surface, the wire tip is directed aside from the pedicle of arch of the vertebra in the upper third. In case of a thoracic surgery, the wire is placed on the rib transverse process, and in case of a lumbar spine surgery is placed on the edge of the lateral joint surface and transverse process. Then the wire is moved along the pedicle projection. Some careful hammer taps are needed to pass the bone surface. If necessary, the wire direction can be adjusted, and, under X-ray control, it is moved further in the front projection. After the tip reaches the pedical medial border, the wire depth is monitored in the lateral projection.

The wire is previously introduced to all levels where the cement introduction is planned. It allows to save time and avoid contamination after changing the X-ray projections. Using the image intensifier, save images of the previously introduced wire for each vertebra.

When changing the X-ray projection, the wire tip should be positioned at least at the level of the vertebral posterior wall; otherwise, it is necessary to change the wire position, returning to the front projection of the X-ray image. Then the needle is carefully moved further with hammer taps, changing the direction, if necessary, in order to reach the vertebral body centre. The wire introduction is planned in accordance with the front image, then, a local anesthetic is administered, and the wire is moved further under X-ray control. The wire has 1 cm marks. It allows to assess the wire moving progress.

roduction of the needle.

needle of the required gauge with a side hole, and together with a cannulated trocar, is moved along the wire by rotational movements. This procedure can be painful, it is necessary to warn the anesthesiologist about the need for appropriate anesthesia. The final X-ray control takes place after the final installation of the needle. The tip of the needle is moved until it reaches the anterior half of the vertebral body. After reaching the final destination, the wire and the cannulated trocar are removed (Fig.4).

Fig.4


Fig.5

Fig.6

Before the introduction of an injection needle with a Luer taper, make sure that the whole bone tissue has been removed from the outer cannula (the vertebroplasty needle of the required gauge) using the trocar. The injection needle with a Luer taper is introduced in order to close the front hole of the outer cannula. The side hole of the injection needle with a Luer taper allows to direct the cement in the desired way. The injection needle with a Luer taper is placed correctly that is monitored by an opening mechanism on the handle (Fig.5). An arrow on the handle indicates the needle side hole (Fig.6). Before operation, make sure that the mark on the injection needle with a Luer taper and the arrow on the handle coincide. It opens the side hole for the cement introduction. The needle has 1 cm marks, that fact allows to control the needle introduction. If there are any doubts about the needle position, return to the X-ray image in the front projection.

1 The direct approach for the cement introduction (Vertebroplasty Needle Kits 3, 4, 5), a minimally invasive method with the transpedicular approach:

The C-type X-ray machine is set for a front projection, a precise place for incision is determined under fluoroscopy.

Under control of the C-type X-ray machine, the skin is punctured. A needle of the required gauge with a side hole, a Luer taper and a selected tip is directed through soft tissues until the vertebral body surface. Under X-ray control, monitor whether the injection needle with a Luer taper is positioned correctly. To avoid fingers radiation, use a radiolucent needle holder (separate MD).

Under X-ray control, the needle of the required gauge with a side hole, a Luer taper and a selected tip (hereinafter – the needle) is introduced in the front plane. After reaching the bone surface, the needle tip is directed aside from the pedicle of arch of the vertebra in the upper third. Introduce a wire with depth marks previously to all levels where the cement introduction is planned. In case of a thoracic surgery, the wire is placed on the rib transverse process, and in case of a lumbar surgery it is placed on the edge of the lateral joint surface and transverse process. Then the wire is moved along the pedicle projection.

After passing through the bone surface, careful pressing combined with rotational movements is necessary (careful hammer tapping is acceptable as well). To correct the needle direction, the tip is turned to the desired direction, and it is moved further under X-ray control in the front projection. As the procedure can be painful, it is necessary to warn the anesthesiologist about the need for appropriate anesthesia. The final X-ray control takes place after the final installation of the needle.

After the tip reaches the pedicle medial border, the needle depth is monitored in the lateral projection. Under X-ray control, the needle is previously introduced to all levels where the introduction of the bone cement VertecemV+ is planned. It allows to save time and avoid examination after changing the X-ray projections. Using the image intensifier, save images of the previously introduced wire for each vertebra. Change the front projection to the lateral one. In the lateral projection the needle tip should be positioned at least at the level of the vertebral posterior wall; otherwise, it is necessary to change the needle position, returning to the front projection of the X-ray image.

To provide room for the injection needle with a Luer taper, it is necessary to compress the bone tissue at the needle tip. For this purpose, rotate the needle and the trocar by 360 ° (Fig.7). Then, remove the trocar. It is necessary if the injection needle with a Luer taper is used to introduce the cement.


Fig. 7

Introduction of the injection needle with a Luer taper

Before the introduction of the injection needle with a Luer taper, make sure that the whole bone tissue has been removed from the outer needle using the trocar. The vertebroplasty needle with a side hole is introduced so as to close the front hole of the outer cannula. The opening side hole allows to direct the cement in the desired way. Position the injection needle with a Luer taper correctly and monitor by the opening handle mechanism.

3 The biopsy method (Vertebroplasty Needle Kits 6, 7):

When the Biopsy Needle (to the Vertebroplasty Needle Kit) of the required gauge has passed through the posterior wall, use a trocar with a bevelled tip. This technique is also used for the vertebroplasty needles of gauge 8 or 10 from Vertebroplasty Needle Kits 1, 2 or 3, but the trocar with a bevelled tip or a cannulated one are removed as described above: together with the wire with depth marks. Introduce the biopsy needle and block it on the handle of the vertebroplasty needle of the desired gauge, with a side hole, a selected tip and a Luer taper (hereinafter - the needle, it may be used as the outer cannula) (see Fig. 8). Move the needle until the final position in the vertebral body (Fig. 8). At the same time, take the bone tissue with the biopsy needle. Do not use the hammer and avoid punctures of the vertebral front wall. Rotate the whole device by 360° at least once. It will help to expose the biopsy material. Next, connect a syringe (not included in the MD), with a Luer tip (1-3 ml), with the biopsy needle, and pull back the plunger for biopsy sampling (Fig. 10). Disconnect the biopsy needle with the syringe from the outer needle. Remove the syringe and use a stylet to remove the sampled bone tissue from the biopsy needle (Fig. 11).

Fig. 8


Fig. 9

Fig. 10


Fig. 11

If the cement is introduced through an injection needle with a Luer taper, the force necessary for it increases with time. Smaller syringes require less force. It is recommended to use 2 ml

ranges first. 1 ml syringes are used when the force applied to a 2 ml syringe becomes sufficient for the cement introduction. At the end of the introduction it is possible to use a trocar further the cement carefully from the needle. Applying force to the trocar, avoid moving the needle further. This method allows to use cement of high viscosity. Viscosity is a key indicator of safety, and the risk of leakage is reduced with increasing viscosity. After starting the introduction, it is recommended to monitor the filling character with X-ray control in the front projection.

Removal of the needle.

Remove the needle only after the cement introduction. Close the needle side hole turning the inner sleeve (Fig.12). It ensures no contact of the cement with the vertebral body. This unique feature ensures no "spikes" when removing the needle. A similar rotary movement releases the needle, and it can be pulled out. The wound is strongly sutured.


Fig. 12

Postsurgical procedures.

After the procedure, the patient is placed on his back for 1 hour to compress the wound. A hematoma may occur at the puncture point. After this, mobilization is possible if permitted by the surgeon. The procedure effect can be assessed immediately as the pain is reduced below the usual level, and there is only discomfort associated with the

puncture.

Method of application of Bone Access Needles

Bone Access Needle is intended for procedures that require access to the bone, as well as to provide access to wires with depth marks, to obtain bone marrow through aspiration techniques, and to administer therapeutic materials, including bone cements. Using Bone Access Needle along with other tools and/or therapeutic, study the surgical technique and/or instructions for use provided by the manufacturers.

6. RULES OF USE AND UTILIZATION

The medical device "Needle Kits for Bone Cement" is intended for use in health care institutions. The device should be used at temperature of 15-25°C, at relative humidity of 30-70%. The personnel, using the product, should be well trained in the biopsy techniques, the kit handling techniques providing approach to the vertebrae and methods of the therapeutic materials application.

The medical device "Needle Kits for Bone Cement" is introduced into the bone during surgical procedures along with simultaneous application of the medical imaging facilities (X-rays, computer tomography, bone scanning).

The medical device "Needle Kits for Bone Cement" belongs to class B of the medical waste classification. The product is intended for single use. Re-use of the device leads to a potential risk of infection. The product is intended for single use only. The product consists of a single-use pieces that should not and can not be re-used, recycled or re-sterilized. They are used for a single surgical procedure. One unit of the product can not be used to treat more than one patient. At the end of the shelf life the device is utilized in accordance with the applicable national legislation.

7. TRANSPORTATION AND STORAGE


The "Needle Kits for Bone Cement" is transported by all types of covered transport in accordance with the current transportation rules at temperature from 0 to + 40°C and at humidity from 30 to 95%.

Do not use after the expiration date on the label. The minimum shelf life is 5 years.

8. PACKAGING AND LABELING

8.1. Needle Kits

The simulated trocar is placed in the vertebroplasty needle. The kits are packed in an individual sealed blister with cover of Tyvek, 364x210x20 mm, with an envelope of Tyvek (Steriking), 364x210 mm. The blister with instructions for use is placed in a cardboard box (consumer package Duplex 230 g/m², Mikrovell, 510x255x30 mm). The cardboard box is closed with two seals with protection from opening.

Biopsy Kits

The stylet should be placed in the biopsy needle. The stylet and the biopsy needle are placed in a protective sheath (polyethylene RAU - PE 140 natur). The biopsy kit is packed in a double pack tubular film (individual package) of Tyvek (Steriking), 75×320 mm internal, 100×385 mm external. This package is placed into a consumer package – a cardboard box (Duplex 230 g/m², Mikrovell, 409×90.5×30 mm). The cardboard box is closed with two labels with protection from opening.

Vertebroplasty Access Needle Kits

The stylet should be placed into the vertebroplasty bone access needle. The stylet and the bone access needle are placed in a protective sheath (polyethylene RAU - PE 140 natur). The vertebroplasty bone access needle is packed in a with cover of Tyvek, 260×85×22 mm. The blister with instructions for use is placed in a cardboard box (Duplex 230 g/m², Mikrovell, 258×82×20 mm) with a central hole for easy viewing.

Biopsy Needles 03.702.250S -252S

The stylet should be placed into the biopsy needle. The stylet and the biopsy needle are placed in a protective sheath (polyethylene RAU - PE 140 natur). The biopsy needles are packed in a blister with cover of Tyvek, 290×90 mm. The blister is placed in an envelope of Tyvek (Steriking), 335 mm (blister and envelope – an individual package). The envelope with instructions for use is placed in a cardboard box (consumer package of Duplex 230 g/m², Mikrovell, 409×90.5×30 mm). The cardboard box is closed with two labels with protection from opening.

The gross weight of the devices in the consumer package is given in the table of technical characteristics.

Transportation packaging: the medical devices are packed in a corrugated cardboard for transportation (585×295×310 mm) which protects from external damage and effects, and are sent to the destination.


Individual package. On the cover of Tyvek for blister and the case of Tyvek there are the following labels: "Item No."; "MD name"; "Lot No."; symbols "Use before"; "Sterilized by ethylene oxide"; manufacturer's name and address; "Product certification code in the UE"; "Used material"; company-distributor's logo; "Do not re-use"; "Do not re-sterilize"; "Refer to the instruction"; "Keep out of sunlight"; "Keep dry"; text "Contents are sterile until opening"; "Amount of 1 pc". Name and address of the manufacturer; Authorized representative in Russia, Russian registration certificate, label version number.

In the box there is also a label which protects from opening with the text "Warranty is annulled if seal is damaged".

Transportation packaging labeling: "Product certification code in the UE"; "Manufacturer's name and address"; "Lot No."; "Production date"; "Item No."; symbols "Keep dry"; "Please, refer to accompanying documents"; "Keep out of sunlight".

Explanation of symbols used in the medical device labeling:

	Refer to the instruction for use		Lot number
	Do not re-use		Item number
	Sterilized by ethylene oxide		Keep out of sunlight
	Manufacturer		Keep dry
	Use before year-month		Do not re-sterilize
	Do not use if the package is damaged		Product certification code in the EU
Content sterile leakless inner package open or damaged			

DELIVERY PACKAGE

Medical device "Needle Kits and Separately Packaged Needles for Bone Cement" includes:

Instruction for use - 1 pc.

Vertebroplasty Needle Kits, variants:

Vertebroplasty Needle Kit includes:

Needle of gauge 8 G, with a side hole, a square end, a Luer taper – 2 pcs, cannulated

trocars – 2 pcs.

Injection needle with a Luer taper – 2 pcs.

Wires with depth marks – 2 pcs.

Instruction for use – 1 pc.

Vertebroplasty Needle Kit includes:

Needle of gauge 8 G, with a side hole, a square end, a Luer taper – 2 pcs.

Cannulated trocars – 2 pcs.

Trocar – 2 pcs.

Injection needle with a Luer taper – 2 pcs.

Wires with depth marks – 2 pcs.

Instruction for use – 1 pc.

Vertebroplasty Needle Kit includes:

Needle of gauge 8 G, with a side hole, a bevelled end, a Luer taper – 2 pcs.

Trocar with a bevelled end – 2 pcs.

Injection needle with a Luer taper – 2 pcs.

Wires – 2 pcs.

Instruction for use – 1 pc.

Vertebroplasty Needle Kit includes:

Needle of gauge 12 G, with a side hole, a square end, a Luer taper – 2 pcs.

er Medical GmbH

Postfach 1737, D-36007 Fulda
Tel: +49 (0) 661-94195-0 Fax: +49 (0) 661-94195-850
info@moeller-medical.com Web: www.moeller-medical.com

car – 2 pcs.

ection needle with a Luer taper – 2 pcs.

ylet – 2 pcs.

struction for use – 1 pc.

Vertebroplasty Needle Kit includes:

eedle of gauge 12 G, with a side hole, a bevelled end, a Luer taper – 2 pcs.

rocar with a bevelled end – 2 pcs.

jection needle with a Luer taper – 2 pcs.

stylet – 2 pcs.

struction for use – 1 pc.

Biopsy Kits, variants:

1. Biopsy Kit includes:

biopsy needle for Vertebroplasty Needle Kit of gauge 8 G, with a Luer taper – 1 pc.

stylet – 1 pc.

struction for use – 1 pc.

2. Biopsy Kit includes:

biopsy needle for Vertebroplasty Needle Kit of gauge 10 G, with a Luer taper – 1 pc.

stylet – 1 pc.

struction for use – 1 pc.

Core Access Needle Kits, variants:

1. Core Access Needle Kit includes:

core access needle with a square end, gauge 8 G, length 140 mm – 1 pc.

stylet – 1 pc.

struction for use – 1 pc.

2. Core Access Needle Kit includes:

core access needle with a square end, gauge 11 G, length 100 mm – 1 pc.

Set – 1 pc.

Instruction for use – 1 pc.

Bone Access Needle Kit includes:

Bone access needle with a bevelled end, gauge 11 G, length 100 mm – 1 pc.

Set – 1 pc.

Instruction for use – 1 pc.

Bone Access Needle Kit includes:

Bone access needle with a square end, gauge 11 G, length 150 mm – 1 pc.

Set – 1 pc.

Instruction for use – 1 pc.

Bone Access Needle Kit includes:

Bone access needle with a bevelled end, gauge 11 G, length 150 mm – 1 pc.

Set – 1 pc.

Instruction for use – 1 pc.

Bone Access Needle Kit includes:

Bone access needle with a square end, gauge 13 G, length 100 mm – 1 pc.

Set – 1 pc.

Instruction for use – 1 pc.

Bone Access Needle Kit includes:

Bone access needle with a bevelled end, gauge 13 G, length 100 mm – 1 pc.

Set – 1 pc.

Instruction for use – 1 pc.

Spinal Needles, variants:

Spinal Needles include:

Spinal needle for a needle of gauge 8 G, length 140 mm – 1 pc.

er Medical GmbH

29-31, D-36043 Fulda · Postfach 1737, D-36007 Fulda
Tel: +49 (0) 661-94195-0 · Fax: +49 (0) 661-94195-850
Email: moeller-medical.com · Web: www.moeller-medical.com

set – 1 pc.

Biopsy Needles include:

Biopsy needle for a needle of gauge 11 G, length 100 mm – 1 pc.

set – 1 pc.

Instruction for use – 1 pc.

Biopsy Needles include:

Biopsy needle for a needle of gauge 11 G, length 150 mm – 1 pc.

set – 1 pc.

Instruction for use – 1 pc.

GENERAL TECHNICAL CHARACTERISTICS

Needles, trocars, stylet, wires with depth marks, injection needles with a Luer taper are made of chrome nickel austenitic stainless steel of 1.4301 mark in accordance with the standard AISI 304, with approved corrosion resistance, wear resistance, stability in heat treatment and durability of sharpening. The devices are resistant to climatic effect while operation, transportation and storage.

Surface finish: electropolishing.

Plastic parts (needle handle, stylet and all trocars cover) are made of polybutylene terephthalate (PBT) of high density (Celanex 2401 MT); adhesive (connection of plastic parts with steel Monopox the MG 017).

Quality criteria for the glued parts resistance to tensile forces: ≥ 450 N

Lock screw connection: low pressure polyethylene;

Protective sheathe for needles: polyethylene RAU - PE 140 natur;

There are needles with blue (gauge 8), purple (gauge 11), yellow (gauge 10), orange (gauge 13) and green (gauge 12) handles.

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2.0-3.0% RAL Design 0606080;

Technical characteristics:

Material	Stainless steel	Chrome nickel austenitic stainless steel in accordance with the standard AISI 304
Mark		1.4301
Density		7900 kg/m ³
Rockwell hardness		195 HB
Surface roughness, µm		0.8 µm Ra
Beveling part sharpness		Not more than 0.025 mm
Corrosion resistance		High grade of resistance
Strength		360 N/mm ²

The medical device "Needle Kit for Bone Cement" includes the following types of needles:


Fig. 1


Fig. 2

Vertebroplasty needle from Vertebroplasty Needle Kit (Fig. 1). Needle of gauge 10 G, with a square end) a side hole (Fig. 2), a square (Fig. 4)/bevelled (Fig. 5) end, a Luer taper (Fig. 3).

The needle has 1 cm marks. Used as an outer cannula; an injection needle with a Luer taper, or trocar, or a biopsy needle is inserted into it. Provides access to the vertebral body.



Fig.3

The arrow on the handle points to a side hole for the injected cement direction.

The Luer taper of the needle is connected with a standard syringe with a Luer taper, which is filled with the bone cement, to be introduced into the vertebral body.


Fig. 4


Fig.5

Trocar, cannulated trocar, trocar with a bevelled end help to compress the bone tissue at the needle tip and remove it from the outer needle.



Fig. 1


Fig. 2

An injection needle with a Luer taper (Fig.1) is used to introduce the bone cement Vertecem V+. A side hole (Fig.2) allows to direct the cement to the desired way. The injection needle with a Luer taper is inserted in a needle of the required gauge from Vertebroplasty Needle Kit (used as an outer cannula). The arrow on the needle handle points the side hole direction (Fig.4). The Luer taper (Fig.3) is connected with a syringe with a Luer taper for introduction of the bone cement. Fig.3 shows an inner sleeve. If turn the sleeve counterclockwise, the side hole of the outer cannula used for the cement flow closes.



Fig. 3





Fig. 4

A stylet is used to remove the
sampled bone tissue from the
biopsy needle.


A wire with depth marks. Limits the depth of the vertebroplasty or biopsy needle injection. Has 1 cm marks.

A biopsy needle for Vertebroplasty Needle Kit, with a Luer taper. Used for the bone biopsy. Inserted in the vertebroplasty needle. The Luer taper is connected with a Luer taper of a standard syringe.


Bone access needle with a square/ bevelled end.


Used to provide access to the bone, to obtain the bone marrow through aspiration techniques and to introduce therapeutic materials including bone cements. Has 1 cm marks. An injection needle with a Luer taper inserts into the bone access needle with a square/ bevelled end.

er Medical GmbH

ppenstr. 29-31, D-36043 Fulda · Postfach 1737, D-36007 Fulda
(0) 661-94195-0 · Fax: +49 (0) 661-94195-850
@moeller-medical.com · Web: www.moeller-medical.com

MöllerMedical

Biopsy needle for needles of gauge 8 G and 11 G. Used to remove the bone samples together with the bone access needle.

Table 2. Technical characteristics of needle kits and separately packaged needles.(characteristics can be provided in shortened form)

No.	Item/Device name	Number of pcs. in the kit	Used material	Length of needle with handle/length of steel cannula**, mm	Diameter**, mm	Weight**, g
Acceptable deviations from dimensions and weight $\pm 5\%$.						
Vertebroplasty Needle Kits include:						
(there are needles with blue (gauge 8), yellow (gauge 10) or green (gauge 12) handles).						
1.	03.702.216S Vertebroplasty Needle Kit includes:					Packed*: 213.90
	Needle of gauge 8 G, with a side hole, a square end, a Luer taper	2	Stainless steel Cutting angle 15°	180/140	4.20	19.42
	Cannulated trocar, blue	2	Stainless steel	187/179	3.40	8.82
	Trocar	2	Stainless steel	189/180	3.00	12.07
	Injection needle with a Luer taper	2	Stainless steel	195/176	3.40	3.51
	Wire with depth marks	2	Stainless steel	230	2.00	5.65


						Packed*: 310,38
2.	03.702.218S	Vertebroplasty Needle Kit includes:				
	Needle of gauge 10 G, with a side hole, a square end, a Luer taper	2	Stainless steel Cutting angle 20°	160/120	3.40	17.88
	Cannulated trocar, yellow	2	Stainless steel	167/159	2.60	3.50
	Trocar	2	Stainless steel	168/159	2.30	7.23
	Injection needle with a Luer taper	2	Stainless steel	175/156	2.60	3.03
	Wire with depth marks	2	Stainless steel	210	2.00	5.15
						Packed*: 308.72
3.	03.702.219S	Vertebroplasty Needle Kit includes:				
	Needle of gauge 10 G, with a side hole, a square end, a Luer taper	2	Stainless steel Cutting angle 20°	160/120	3.40	17.77
	Trocar with a bevelled end, yellow	2	Stainless steel	163/154	2.60	8.51
	Injection needle with a Luer taper	2	Stainless steel	174/155	2.60	3.05

	Stylet	2	Stainless steel	160/151	2.30	6.96
4.	03.702.220S Vertebroplasty Needle Kit includes:					Packed*: 297.06
	Needle of gauge 12 G, with a side hole, a square end, a Luer taper	2	Stainless steel Cutting angle 20°	160/120	2.70	17.03
	Trocar	2	Stainless steel	168/159	1.90	5.56
	Injection needle with a Luer taper	2	Stainless steel	174/155	1.90	2.75
	Stylet	2	Stainless steel	153/144	1.60	4.21
5.	03.702.221S Vertebroplasty Needle Kit includes:					Packed*: 297.20
	Needle of gauge 12 G, with a side hole, a bevelled end, a Luer taper	2	Stainless steel Cutting angle 20°	160/120	2.70	17.01
	Trocar with a bevelled end	2	Stainless steel	163/154	1.90	5.40
	Injection needle with a Luer taper	2	Stainless steel	174/155	1.90	2.76
	Stylet	2	Stainless steel	160/151	1.60	4.27

Biopsy Kits, variants: (there are needles with blue (gauge 8) and yellow (gauge 10) handles)						
1.	03.702.222S Biopsy Kit includes:					Packed*: 114.0
	Biopsy needle for Vertebroplasty Needle Kit, gauge 8 G, with a Luer taper	1	Stainless steel Cutting angle 15°	194/175	3.40	In a protective sheath 7.06 /without a protective sheath 4.92
	Stylet	1	Stainless steel	204.5/196	2.90	11.36
2.	03.702.223S Biopsy Kit includes:					Packed*: 107.40
	Biopsy needle for Vertebroplasty Needle Kit, gauge 10 G, with a Luer taper	1	Stainless steel Cutting angle 20°	174/155	2.60	In a protective sheath 5.91 / without a protective sheath 4.03

Möller Medical GmbH

Wasserkuppenstr. 29-31, D-36043 Fulda · Postfach 1737, D-36007 Fulda
 Tel.: +49 (0) 661-94195-0 · Fax: +49 (0) 661-94195-850
 E-Mail: info@moeller-medical.com · Web: www.moeller-medical.com


**MöllerMedical**

	Stylet	1	Stainless steel	184.5/176	2.10	5.84
Bone Access Needle Kit, variants: (there are needles with blue (8G), purple (11G) or orange (13G) handles)						
1.	03.702.240S Bone Access Needle Kit includes:					
	Bone access needle with a square end, gauge 8 G, length 140 mm	1	Stainless steel. Needle with wave cutting and a Luer taper. Cutting angle 20°.	180/140	4.20	In a protective sheath 40.1/ without a protective sheath 34.6 Packed*: 92.1 g
	Stylet	1	Stainless steel			
2.	03.702.241S Bone access needle with a square end, gauge 11 G, length 100 mm					
		1	Stainless steel. Needle with wave cutting and a Luer	140/100	3.15	In a protective sheath 25.60/ without a

			taper. Cutting angle - 20°			protective sheath 24.24 Packed*: 77.60 g
	Stylet	1	Stainless steel	189,5/180	3,4	13,3g
3.	Bone access needle with a bevelled end, gauge 11 G, length 100 mm 03.702.243S	1	Stainless steel. Needle with wave cutting and a Luer taper. Cutting angle - 20°	140/100	3.15	In a protective sheath 25.5/ without a protective sheath 24.2 Packed*: 77.5 g
	Stylet	1		148,5/139	2,63	3,2g
4.	03.702.244S Bone access needle with a square end, gauge 11 G, length 150 mm	1	Stainless steel. Needle with wave cutting and a Luer taper. Cutting angle- 20°	190/150	3.15	In a protective sheath 28.98/ without a protective sheath 27.07 Packed*: 81.1 g

Möller Medical GmbH

Wasserkuppenstr. 29-31, D-36043 Fulda · Postfach 1737, D-36007 Fulda
 Tel.: +49 (0) 661-94195-0 · Fax: +49 (0) 661-94195-850
 E-Mail: info@moeller-medical.com · Web: www.moeller-medical.com

**MöllerMedical**

	Stylet	1	Stainless steel.	198,5/189	2,63	3,5g
5.	03.702.245S Bone access needle with a bevelled end, gauge 11 G, length 150 mm Packed*: 80.95 g	1	Stainless steel. Needle with wave cutting and a Luer taper.	190/150	3.15	In a protective sheath 20.94/ without a protective sheath 27.04
	Stylet	1	Stainless steel.	193,83/184,33	2,63	8,3g
6.	03.702.246S Bone access needle with a square end, gauge 13 G, length 100 mm	1	Stainless steel. Needle with wave cutting and a Luer taper. Cutting angle - 20°	140/100	3.15	In a protective sheath 22.51/ without a protective sheath 21.2 Packed*: 74.50 g
	Stylet	1	Stainless steel.	147,5/ 38	1,98	3,5g
7.	03.702.247S Bone access needle with a bevelled end, gauge 13 G, length 100 mm	1	Stainless steel. Needle with wave cutting and a Luer	140/100	2.50	In a protective sheath 22.5/ without a


			taper. Cutting angle - 20°			protective sheath 21.2 Packed*: 74.50 g
	Stylet	1	Stainless steel.	143,75/134,25	1,98	1,5g
Biopsy Needles, variants: (handles differ by colors – blue handles are intended for bone access needles of gauge 8, purple handles are intended for needles of gauge 11)						
1.	03.702.250S Biopsy needle for the needle of gauge 8 G, length 140 mm	1	Stainless steel. Needle tube cutting - 15°.	195/176, for length of 140 mm	3.40	In a protective sheath 21.4/ without a protective sheath 16.8 Packed*: 77.8
	Stylet	1	Stainless steel	205,4/197	2,9	11,0g

2.	03.702.251S Biopsy needle for the needle of gauge 11 G, length 100 mm	1	Stainless steel. Needle tube cutting 15°.	155/136, for length of 100 mm	2.60	In a protective sheath 13.11/ without a protective sheath 9.53 Packed*: 69.55
	Stylet	1	Stainless steel	165,4/157	2,1	4,7g
3.	03.702.252S Biopsy needle for the needle of gauge 11 G, length 150 mm	1	Stainless steel. Needle tube cutting - 15°.	205/186, for length of 150 mm	2,60	In a protective sheath 16.38/ without a protective sheath 11.49 Packed*: 72.8
	Stylet	1	Stainless steel	215,4/207	2,1	5,0g
Technical characteristics:						
The radius of the working parts of blunting needles		0,02 mm				
bending strength		0,30 mm				

fracture strength	They do not have signs of fracture after testing
corrosion resistance	Corrosion
The severity of the stitching	0,025 mm

Table 3. Maximum puncture force parameter information

Teilenummer	Synthes	Ø + tip type	Test in PU-film		Test in bone substitute material		Basis	Adoption from
	Code		Needle + Trokar	K-wire	Needle + Trokar	K-wire		
10001541	03.702.216S	4,2 Wavy	17,3 N	3,4 N	202,0 N	n/a	Test result	n/a
10001542	03.702.218S	3,4 Wavy	27,0 N	7,3 N	172,0 N	n/a	Test result	n/a
10001543	03.702.220S	2,7 Wavy	14,2 N	n/a	132,0 N	n/a	Test result	n/a
10001544	03.702.219S	3,4 Diagonal	13,1 N	n/a	164,0 N	n/a	Test result	n/a
10001545	03.702.221S	2,7 Diagonal	17,3 N	n/a	123,0 N	n/a	Test result	n/a
10001853	03.702.241S	3,15 Wavy	12,8 N	n/a	140,0 N	n/a	Test result	n/a
10001855	03.702.243S	2,5 Wavy	9,4 N	n/a	113,0 N	n/a	Test result	n/a

10001857	03.702.245S	3,15 Diagonal	11,4 N	n/a	143,0 N	n/a	Test result	n/a
10001859	03.702.250S	3,4 Biopsy	13,1 N	n/a	164,0 N	n/a	Adoption	03.702.219S
10001861	03.702.251S	2,6 Biopsy	17,3 N	n/a	123,0 N	n/a	Adoption	03.702.221S
10001862	03.702.252S	2,6 Biopsy	17,3 N	n/a	123,0 N	n/a	Adoption	03.702.221S
10001574	03.702.222S	3,4 Biopsy	13,1 N	n/a	164,0 N	n/a	Adoption	03.702.219S
10001573	03.702.223S	2,6 Biopsy	17,3 N	n/a	123,0 N	n/a	Adoption	03.702.221S
10001852	03.702.240S	4,2 Wavy	17,3 N	n/a	202,0 N	n/a	Adoption	03.702.216S
10001854	03.702.244S	3,15 Wavy	12,8 N	n/a	140,0 N	n/a	Adoption	03.702.241S
10001855	03.702.246S	2,5 Wavy	9,4 N	n/a	113,0 N	n/a	Adoption	03.702.243S
10001858	03.702.247S	2,5 Diagonal	17,3 N	n/a	123,0 N	n/a	Adoption	03.702.221S

Möller Medical GmbH

Wasserkuppenstr. 29-31, D-36043 Fulda · Postfach 1737, D-36007 Fulda
Tel.: +49 (0) 661-94195-0 · Fax: +49 (0) 661-94195-850
E-Mail: info@moeller-medical.com · Web: www.moeller-medical.com

11. REQUIREMENTS FOR MAINTENANCE AND REPAIR

The device is intended for single use, it is not subject to maintenance or repair.

12. REQUIREMENTS FOR ENVIRONMENTAL PROTECTION

The device "Needle Kits for Bone Cement" manufactured by Möller Medical GmbH, Germany, is safe for the environment.

13. WARRANTY, POSTMARKETING SURVEILLANCE AND COMPLAINTS

The manufacturer guarantees that the product complies with the technical documentation requirements if the the rules of transportation, storage and use are observed and if the packaging is not opened or damaged. The warranty covers the entire shelf life of the medical device. Möller Medical GmbH provides postmarketing surveillance in accordance with Council Directive 93/42/EEC in order to monitor any signal related to safety, new side effects, emerging issues, new category of complaint/complaints regarding the original parameters, or additional risks.

Manufacturer:

Möller Medical GmbH, Wasserkuppenstrasse 29-31, 36043 Fulda, Germany

Telephone/fax: +49 (0) 661 941 95 0

E-mail: info@moeller-medical.com

Website: www.moeller-medical.com

Manufacturing site: Möller Medical GmbH Wasserkuppenstrasse 29-31, 36043 Fulda, Germany

Telephone/fax: +49 (0) 661 941 95 0

E-mail: info@moeller-medical.com


For any questions on device operation, quality, etc., within the territory o the Russian Federation, please refer to Johnson&Johnson LLC:


Legal address: Krylatskaya str., 17, build.2, Moscow 121614, Russia

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