INSTRUCTION

For ophtalmic spectacle frames

1. **PURPOSE:** Ophthalmic spectacle frames made of metal, plastic, combined produced by "OWP Brillen GmbH", meant for holding the lenses in a correct fixed position over the eyes.

2. MAIN TECHNICAL CHARACTERISTICS AND CONSTRUCTION

Ophthalmic spectacle frames, depending on design features and manufacturing technology, are divided into types indicated in Table 1.

Table 1.

Designation of frame type	Name of frame type	Type and designation of temples
FP	Plastic optical frames	Plastic rigid PR; Plastic elastic PE.
FM	Metal optical frames	Metal rigid PR; Metal elastic PE.
FC	Combined optical frames	Metal rigid MR; Metal elastic PE; Plastic rigid PR; Plastic elastic PE

Figure. 1. Metal optical frames (FM)

Figure. 2. Plastic optical frames (FP)

Figure. 3. Combined optical frames (FC)

Ophthalmic spectacle frames consist of a front, temples, nose pads, hinges. Hinges may be simple and spring-based.


Figure. 4. Structure of ophthalmic spectacle frames

Ophthalmic spectacle frames are characterized by 3 main geometric parameters: width (a), bridge width (b) and temple length to the curve (L).

The sizes used for spectacle frames by “OWP Brillen GmbH” are given in the table 2:

Ophthalmic spectacle frames		Distance between the vertical tangents to the facet groove of the rim (a). (mm)	Minimum distance between the vertical tangents to the bottom of the facet groove of the rim (b). (mm)	Temple length from the hinge axis to the curve (L) (mm)	Temple mass without lenses (g)
Ophthalmic spectacle frames (metal)	men	55 ± 0,2	16 ± 0,2	100 ± 1,0	22 ± 1,0
	men	56 ± 0,2	20 ± 0,2	95 ± 1,0	22 ± 1,0
	women	51 ± 0,2	17 ± 0,2	95 ± 1,0	12 ± 1,0
	women	54 ± 0,2	18 ± 0,2	100 ± 1,0	20 ± 1,0
	unisex	50 ± 0,2	18 ± 0,2	95 ± 1,0	14 ± 1,0
	unisex	56 ± 0,2	20 ± 0,2	90 ± 1,0	18 ± 1,0
Ophthalmic spectacle frames (plastic)	men	47 ± 0,2	18 ± 0,2	95 ± 1,0	26 ± 1,0
	men	58 ± 0,2	20 ± 0,2	95 ± 1,0	24 ± 1,0
	women	54 ± 0,2	15 ± 0,2	95 ± 1,0	24 ± 1,0
	women	54 ± 0,2	13 ± 0,2	95 ± 1,0	22 ± 1,0
	unisex	47 ± 0,2	18 ± 0,2	90 ± 1,0	20 ± 1,0
	unisex	55 ± 0,2	18 ± 0,2	95 ± 1,0	24 ± 1,0
Ophthalmic spectacle frames (combined)	men	58 ± 0,2	20 ± 0,2	90 ± 1,0	24 ± 1,0
	men	53 ± 0,2	16 ± 0,2	95 ± 1,0	20 ± 1,0
	women	54 ± 0,2	18 ± 0,2	95 ± 1,0	22 ± 1,0
	women	51 ± 0,2	17 ± 0,2	100 ± 1,0	20 ± 1,0
	unisex	54 ± 0,2	18 ± 0,2	95 ± 1,0	24 ± 1,0
	unisex	50 ± 0,2	16 ± 0,2	95 ± 1,0	22 ± 1,0

On the surface of the frames there may be no sharp edges or burrs, shells, scratches, bubbles, sagging solder visible to the naked eye.
Temples should rotate on the axis smoothly, without jamming.

Decorative trim (if any) is from nickel silver (N18 CUNI18ZN20).

Ophthalmic spectacle frames correspond to the standards: EN ISO 12870:2014, EN 12472:2009, EN 16128:2015, to the regulations EU-Directive 93/42/EEC concerning medical devices; Regulation (EC) No 1907/2006 on the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH).

Ophthalmic spectacle frames are chosen according to the anthropometric parameters of the patient.

3. INDICATIONS

Application on doctor's orders.

4. The following may be considered CONTRAINDICATIONS:
- individual intolerance to sight correction by means of spectacles;

- impossibility of sight correction by means of spectacles.

5. MATERIALS:

Frames are made of materials that do not cause irritation, allergic or toxic reactions during their exploitation upon contact with human skin.

Plastic frames are produced of cellulose acetate (cellulose acetate polymer (69-75 %) and diethyl phthalate (25-31%)), hinges and screws nickel silver (N18 CUNI18ZN20)

Metal and combined frames are produced of nickel silver (N18 CUNI18ZN20) and cellulose acetate (cellulose acetate polymer (69-75 %) and diethyl phthalate (25-31%)), or clean titanium (number CAS 7740-32-6, number EINECS 241-036-9), hinges and screws are produced of nickel silver (N18 CUNI18ZN20). Tips of temples are produced of cellulose acetate (cellulose acetate polymer (69-75 %) and diethyl phthalate (25-31%)). Nose pads are produced of silicone LSR50.

The transport insert made of acrylic Polycryl Clear.

Materials comply with the requirements of biological safety.

6. MARKING:

The right temple indicates the manufacture's name, conformity mark and country of manufacture.


Figure 5. Frame marking. Right temple.

Left temple shall have factory marking with indication of:

model number (1372), technical size (54-14-135), colour or coating code (300).


Figure 6. Frame marking. Left temple.

Each frame is packed into an individual consumer plastic packaging. Each consumer packaging (plastic bag) is marked with a conformity mark.

Figure 7. Individual plastic packaging.

Plastic packaging of each frame or the label placed into plastic packaging should contain:

- nominal value of basic sizes;
- model number;
- model colour;
- name of product;
- manufacturer;
- date of manufacturing;
- storage life;
- conformity mark.

Frames in plastic packaging should be compactly placed into group packaging box. Group packaging box should contain frames of the same type with the same basic sizes, in quantity of 8 - 10 pcs.

Each group packaging box should have a sticker of the manufacturer on it, which should contain following inscriptions:

- model number;
- model colour;
- nominal value of basic sizes;
- quantity of frames

Dimensions of the packaging box, in cm (length x width x height): 23 x 15 x 5,8.

Maximum load on the group packaging box is 10 kg.

Transport packaging is made only by agreement with the consumer.

Frames packed in transport packaging should be resistant to transport shaking with a maximum acceleration of 30 m/s² and frequency of 2 to 3 Hz.

Marking of transport packaging with application of the handling mark "Caution, fragile".

7. PREPARING THE ITEM FOR USE: The spectacle frame should not lose its original shape after fitting the lenses. The rims should be symmetrical. No dents, scratches or gaps between spectacle lenses and rims, nor any other assembling defects are allowed.

Spectacle lenses fixed in the frame should not move or fall out of the rims during normal spectacle use.

8. TERMS OF STORAGE AND USE:

- It is not recommended to leave the rim in direct sunlight at high temperatures (eg, car control panel);
- Need to avoid the fall of the rim on hard surfaces, and mechanical loads that can damage the rim;
- Storage frames need to use a case or pouch, is not permitted the wearing of eyeglasses in the pockets of clothing or bags in order to avoid scratches on the cover and mechanical damage;
- To put on and take off glasses with both hands is necessary - this prevents their deformation;
- Need to avoid wearing eyeglasses during sleep;
- Are not allowed to apply to the rim significant mechanical forces;
- Caution should be exercised when doing sports and the production of works involving sudden movements of the head, the slopes, jumps and bumps.
 - Do not apply significant mechanical efforts to the frame;
 - Keep the frame in a case or any other package to protect it from mechanical damage;
 - Full average term of use is no less than 3 years..
 - Terms of storage: in heated storages under the environmental temperature from + 5°C to + 50°C, relative air moisture 80%, atmospheric pressure 700-1060 hPa
 - No impact of aggressive substances (acids, alkali) and organic solvents and their fumes is allowed during storage.

9. TERMS OF USE

Ophtalmic spectacle frames are used in all climate zones at any time of the day.

10. PRECAUTIONS FOR USE

While using ophthalmic spectacle frames it is necessary to exercise caution practising sports and activities involving sudden head movements, bends, jumps and bumps.

11. POSSIBLE SIDE EFFECTS

Allergic reaction in case of individual intolerance of frames material is possible.

12. INFORMATION ON THE POSSIBLE IMPACT OF THE USE ON ABILITY TO DRIVE VEHICLES, MACHINERY

Drivers with vision problems may drive vehicles and machinery using spectacle frames (in the manufacture of frames for corrective glasses).

13. STERILITY

Unsterile.

14. MAINTENANCE AND REPAIRS

Ophthalmic spectacle frames should be repaired in special laboratories licensed for production, maintenance and repairs.

15. WARRANTY OBLIGATIONS

Storage life: plastic and combined frames – 5 years, metal frames – 10 years.

Full storage life – not less than 3 years with temple bends not more than 30 000.

16. TRANSPORTATION

Ophthalmic spectacle frames are transported by all kinds of transport in covered vehicles in accordance with the rules of transportation of goods applicable to this kind of transport.

Temperature from - 45°C to + 50°C

Relative air moisture up to 98%

Atmospheric pressure 700-1060 hPa

17. DISPOSAL

Spectacle frames do not contain harmful substances and components that pose a threat to human health and the environment during and after the end of life and disposal. Ophthalmic spectacle frames should be disposed of as household waste.

18. CERTIFICATES

The manufacturer uses its own quality system that is subject to verification and updating on an annual basis.

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
Legal address: 111024, Moscow, Andronovskoe highway, 26, build. 6, off. 11B
Phone number 8 (495) 660-01-25

Klaus Lippert
Chief financial officer
OWP Brillen GmbH

Signature

Date 16th of March 2017

B 0451/17

Protocol No. B 0451/17

I, Gottfried Braun, notary public in Passau, hereby certify,
that the above is the true signature, subscribed in my pres-
ence, of

Mr. Klaus Lippert,
born on the 10th of March 1974,
address: Mühlbergstr. 16, 94535 Eging am See
who is personally known to me,
acting


on behalf of


OWP Brillen GmbH

a company with registered office at Passau, Germany.

Upon my inspection of the electronic commercial register on
the 15th day of March 2017, I certify that OWP Brillen GmbH
is registered with the district court of Passau under No. HRB
4477 and that Mr. Klaus Lippert is entitled to act individu-
ally as this company's legal representative and authorized
procurator.

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