VEINO PLUS®

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Ad Rem Technology
ISO 13485 : 2012 certified
VEINOPLUS® is a registered
trademark of Ad Rem Technology

Ad Rem Technology
162, rue du Faubourg Saint-Honoré
75008 Paris - France
RCS 453 317 265
Tel : +33 (0)1 42 60 00 22
Fax : +33 (0)1 42 60 00 63

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


Distribution, whether in whole or in part, in any manner whatsoever and by any means whatsoever (printed or electronic), unless approved in advance and in writing by AD REM TECHNOLOGY, is prohibited and deemed to be infringement of intellectual property rights which can be prosecuted under articles 335-2 et L. 335-3 of the French Intellectual Property Code (Code de la Propriété Intellectuelle).

DEFINITION OF SYMBOLS ON THE DEVICE

Follow instructions for use

The device has one applied part type BF (electrodes)

Adjustments of intensity buttons

ON/OFF switch button

This symbol means that VEINOPLUS® is in conformity with the requirements of the Medical Device Directive 93/42/EEC, 0535 is the number of notified body.

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INTERNATIONAL WARRANTY

VEINOPLUS® BACK is manufactured according to high-quality standards.

VEINOPLUS® BACK is guaranteed for 12 months from the date of the purchase against any manufacturing defect, except for the battery, glass breakage, abnormal wear and tear of electrodes and cable, as well as any damage due to abusive handling or removing parts of VEINOPLUS® BACK by unauthorised third parties.

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Solution: Check that **both** electrodes are correctly applied to the skin and that both are connected to the device. Increase the intensity of stimulation by single short pushes on the key ⊕.

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Solution: Check that your skin is clean and free from grease and lotion. Improve stimulation effectiveness and comfort by gently rubbing one or two drops of water onto the gel surface. If the problem persists, you need to change the electrodes.

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KEY TO DIAGRAM ON INTERNAL COVER PAGE

- | | |
|----------------------------|-------------------------------|
| 1 ON / OFF | 5 Battery replacement display |
| 2 Intensity control | 6 Usage counter |
| 3 Intensity of stimulation | 7 Cable jack |
| 4 Remaining time | 8 Electrodes |

ENVIRONMENTAL CONDITIONS FOR USE, STORAGE AND TRANSPORT

■ VEINOPLUS® BACK need special precautions regarding EMC and need to be installed and put into service according to the EMC information provided in the Instruction Manual (see EMC tables).

■ Portable and mobile radio-frequency (RF) communications equipment can affect medical electrical equipment.

■ The use of accessories, electrodes and cables other than recommended by Ad Rem Technology may result in increased Emissions or decreased Immunity of VEINOPLUS® BACK.

■ VEINOPLUS® BACK does not produce excessive electromagnetic interference (EMI) in test according to IEC 60601-1-2 standard. **It may be used during air travel but not during take-off and landing.** (see EMC tables).

■ VEINOPLUS® BACK should not be used adjacent to or stacked with other equipment.

■ Operation in close proximity (e.g. 1 m) to a shortwave or microwave therapy equipment may produce instability in the stimulator output (see EMC tables).

■ VEINOPLUS® BACK is not intended for use in sterile environments.

■ VEINOPLUS® BACK and accessories must be used/ stored/ transported in temperatures from +5 to +27 °C and with Relative Humidity (RH) 0~50%, pressure 850 to 1060 hPa and altitude <3000m.

■ At the end of their life the VEINOPLUS® BACK devices should be disposed of properly with respect to environmental regulations.

SERVICE

■ VEINOPLUS® BACK requires neither routine maintenance nor calibration. If necessary, the device can be cleaned by alcohol impregnated wipes.

■ The VEINOPLUS® Pack electrodes and the cable supplied with VEINOPLUS® BACK can be replaced through your local VEINOPLUS® distributor.

■ All services are handled either by Ad Rem Technology or by authorised representatives. To obtain this service please contact your vendor.

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SAFETY

To ensure safe and effective usage of VEINOPLUS® BACK, please follow the instructions of this manual.

VEINOPLUS® BACK complies with the IEC 60601-1 standard on general safety requirements for electro-medical devices, the IEC 60601-1-2 standard on Electromagnetic Compatibility (EMC) and the IEC 60601-2-10 standard on particular safety requirements for nerve and muscle stimulators. The impulses from VEINOPLUS® BACK carry a very low quantity of electricity. This quantity is below the limits established by any international safety standards including United States AAMI/ANSI NS-4 Standard (ANSI = American National Standard Institute).

TECHNICAL DESCRIPTION

VEINOPLUS® BACK is an internally powered (9 V battery) electrical muscle stimulator that delivers specific low frequency and low voltage stimulation impulses. The device is intended for continuous use and the degree of protection against ingress of water is IPX0. Input power is below 0.5 W and output power of VEINOPLUS® BACK is below 0.05 W.

Output parameters of VEINOPLUS® BACK:

- The voltage waveform of the VEINOPLUS® BACK is rectangular monopolar and balanced with ZERO net d.c. component,
- The time duration of output impulses ranges from 25 to 240 microseconds,
- During impulses, the instantaneous output voltage reaches a peak value of 50Vp $\pm$ 5%, when loaded with impedance in a range from 250 Ω to 10,000 Ω .

Maximum ratings of the VEINOPLUS® BACK output when loaded with 500 Ω resistance:

The current intensity may reach a peak value of 100 mA $\pm$ 5% during impulses. The waveform is modulated in a special pattern resulting in muscle contractions ranging from 2 to 120 beats-per-minute.

VEINOPLUS® BACK is supplied with all the accessories required for an easy, safe and effective treatment.

These include self-adhesive skin electrodes (VEINOPLUS® Pack), one cable, one 9 V alkaline battery and one instruction manual.

Electrodes:

The VEINOPLUS® BACK produces electrical impulses and delivers energy to the surface of the skin by means of a pair of bio-compatible electrodes (VEINOPLUS® Pack).

For effective and safe treatment, use only VEINOPLUS® Pack electrodes in their sealed original package. It is possible to use the electrodes several times. It is recommended that only one person uses the same set of electrodes. Before applying the electrodes, the skin must be clean and free from grease and lotion.

Used electrodes should be replaced as soon as they become dry or are no longer adhesive. When electrodes begin losing adhesion, you can gently rub one or two drops of water onto the gel surface to extend usage. Specific instructions for use and maintenance of the electrodes are included in their packages.

Batteries:

The battery should be changed when the icon on the screen blinks intermittently. To change the battery, follow the drawings from the cover page. Use only **alkaline or rechargeable** (contact your dealer for the availability of rechargeable option) **9 V batteries** type 6LR61 (or equivalent e.g. 6LF22, AM6, MN1604 or GP1604). Don't use rechargeable lithium batteries. If other batteries are used, the VEINOPLUS® BACK may temporarily stop functioning.

Precautions concerning the handling of batteries:

- Do not recharge alkaline batteries,
- Do not short-circuit the batteries by contact with metallic objects,
- Do not dispose of batteries in fire,
- Store batteries in temperatures ranging from +5 to +27 °C,
- When your device is not used for a long period of time, remove the battery,
- Waste batteries must be disposed properly, to preserve the environment.

Cable:

Use only cable supplied with VEINOPLUS® BACK device.
Cable is 135 cm long.

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HOW TO USE VEINOPLUS® BACK?

Please read the instructions and refer to the drawings from the cover page.

- 1 VEINOPLUS® BACK can be applied on the upper, middle or lower back, depending on the painful area.
- 2 The battery compartment is located on the back of the device. Do not use any tools to open the battery compartment.
- 3 Insert the battery into the device according to the +/- signs and close the battery compartment firmly.
- 4 Remove the electrodes from the plastic shields without pulling on the cable.
- 5 Place both electrodes on the painful area: neck/shoulder(5a), middle of the back(5b) or lumbar area(5c).
- 6 DO NOT switch on the device before connecting the electrodes.
- 7 Plug the electrodes in the cable and then plug the cable in the device.
- 8 Turn on VEINOPLUS® BACK by pressing the central key marked ① once.
SLOWLY increase the intensity of stimulation by pushing the ⊕ key.
- 9 If sensations become too strong and/or uncomfortable you may reduce the intensity by short pushes on the ⊖ key.
- 10 VEINOPLUS® BACK will switch OFF automatically after 45 minutes but you can stop the stimulation at any time by pushing "ON/OFF" button.
- 11 Always disconnect the electrodes after the device has been turned "OFF".
- 12 Remove the electrodes from the skin without pulling on the cable.
- 13 Apply the electrodes onto the side "ON" of the plastic shields and store them in the sealed bag provided.

ELECTRODE PLACEMENT

Please refer to the drawings from the cover page (5 a. b. c.) regarding the possible positions of electrodes.

CONTRA-INDICATIONS AND WARNINGS

Contra-indications:

VEINOPLUS® BACK must not be used on patients with implanted cardiac pacemakers or defibrillators.

Warnings:

- Do not apply electrodes on skin with a lesion of any kind (cancerous lesions, etc), inflammation, infection, or severe varicose veins.
- Should a skin rash or a skin burn occur, immediately discontinue use and if this persists, contact your health practitioner.
- Pregnant women should not use VEINOPLUS® BACK.

Precautions and limitations:

- Simultaneous connection from a patient to a high-frequency electro-surgical equipment may result in skin burns at the site of electrodes and possible damage to the stimulator.
- Application of the electrodes across the thorax may increase the risk of cardiac fibrillation.
- Do not use VEINOPLUS® BACK while operating powered machinery or operating a vehicle.
- Keep out of reach of children.
- Use only VEINOPLUS® accessories (electrodes, cable, etc...) with your VEINOPLUS® BACK device.

WHEN TO USE VEINOPLUS® BACK?

Use VEINOPLUS® BACK on painful areas or zones needing relaxation, when needed.

The more you use VEINOPLUS® BACK, the more effective it becomes.

There is no risk of overuse.

Each session lasts 45 minutes and the device stops automatically.

You can have several sessions in a row. There is no minimum rest time between sessions.

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WELCOME

Dear Customer,

We thank you for your confidence in VEINOPLUS® stimulators.

VEINOPLUS® is a patented technology as a result of several years of Research & Development to offer you a very effective yet comfortable and safe stimulation.

VEINOPLUS® BACK is a Medical Device, according to the Medical Device Directive 93/42/EEC, that can relieve your pain and help you in your daily life.

Please read carefully this Instruction Manual before using the device.

Ad Rem Technology and its distributors decline any responsibility for any use of the VEINOPLUS® BACK that does not comply with this Instruction Manual.

HOW IT WORKS

VEINOPLUS® BACK combines a complex and unique stimulation pattern including a TENS component and a NMES component.

TENS (Transcutaneous Electrical Nerve Stimulation) is widely recognized to alleviate pain, whereas NMES (Neuro-muscular Electrical Stimulation) stimulates the muscles to cause their contractions and can relieve (relax) muscle spasms.

Specifically developed to bring comfort and ease of use, VEINOPLUS® BACK combines 4 effects aiming back pain relief into one single stimulation program.

INDICATIONS

VEINOPLUS® BACK is designed to:

- Relieve pain by local analgesia in the area of stimulation,
- Deliver a relaxing massage to the area of stimulation,
- Relieve muscle spasms in the area of stimulation,
- Allow the release of natural endorphins to relieve general pain.

Results are not guaranteed, and can vary from person to person. If symptoms persist, consult your doctor.

Table 1 : Guidance and manufacturer's declaration – electromagnetic emissions

The VEINOPLUS® 2.1 is suitable for use in the specified electromagnetic environment. The customer and/or the user of the equipment should assure that it is used in an electromagnetic environment as described below.		
Emission Test	Compliance	electromagnetic environment Guidance
RF emission CISPR 11	Group 1	The VEINOPLUS® 2.1 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The VEINOPLUS® 2.1 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emission CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	

Table 2 : Guidance and manufacturer's declaration – electromagnetic immunity

The VEINOPLUS® 2.1 is intended for use in the electromagnetic environment specified below. The customer or the user of the VEINOPLUS® 2.1 should assure that it is used in such an environment.			
Immunity Test	Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic discharge (ESD)	± (2, 4, 6) kV contact ± (2, 4, 8) kV air	± (2, 4, 6) kV contact ± (2, 4, 8) kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic materials, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec.	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the VEINOPLUS® 2.1 requires continued operation during power mains interruptions, it is recommended that the VEINOPLUS® 2.1 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE : UT is the a.c. mains voltage prior to application of the test level.			

Table 4 : Guidance and manufacturer's declaration – electromagnetic immunity

The VEINOPLUS® 2.1 is intended for use in the electromagnetic environment specified below. The customer or the user of the VEINOPLUS® 2.1 should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the VEINOPLUS® 2.1, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1,2\sqrt{P}$ $d = 1,2\sqrt{P}$ 80 MHz to 800 MHz $d = 1,2\sqrt{P}$ 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a, should be less than the compliance level in each frequency range ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	
NOTE 1 : At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 : These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the VEINOPLUS® 2.1 is used exceeds the applicable RF compliance level above, the VEINOPLUS® 2.1 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the VEINOPLUS® 2.1. b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

**Table 6 :
Recommended separation distances between
portable and mobile RF communications equipment and the VEINOPLUS® 2.1**

The VEINOPLUS® 2.1 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the VEINOPLUS® 2.1 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the VEINOPLUS® 2.1 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,5 GHz $d = 1,2\sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 : At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 : These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

VEINOPUS[®] BACK is manufactured according to high quality standards

VEINOPUS[®] BACK is guaranteed for 12 months from the date of the purchase against any manufacturing defect, except for the battery, glass breakage, abnormal wear and tear of electrodes and cable, as well as any damage due to abusive handling or removing parts of VEINOPUS[®] BACK by unauthorised third parties.

This International Warranty does not affect your consumer statutory rights under applicable local Laws, nor your rights against the seller arising from your purchase transaction. Please contact Ad Rem Technology or your local distributor.

Problem: Device non-functioning, blank screen.

Solution: Insert a fresh alkaline battery.

Problem: Screen functioning but there is no feeling of stimulation.

Solution: Check that **both** electrodes are correctly applied to the skin and that both are connected to the device. Increase the intensity of stimulation by single short pushes on the key 

Problem: Stimulation is less effective or less pleasant.

Solution: Check that your skin is clean and free from grease and lotion.

Improve stimulation effectiveness and comfort by gently rubbing one or two drops of water onto the gel surface. If the problem persists, you need to change the electrodes.

If your problem is not solved by following the above instructions, please contact Ad Rem Technology or an authorised representative.

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© Copyright 2013 Paris - Ad Rem Technology