InQwire Diagnostic Guide Wires with accessories
Instructions For Use
Consisting of 8 pages
Dated September 13, 2017

ACKNOWLEDGEMENT


State of Utah)

208

County of Salt Lake)

On this 13th day of September, 2017, before me, Jennifer Wojciechowski, a notary public, personally appeared Dan Clark, Director, Regulatory Affairs, proved on the basis of satisfactory evidence to be the person whose name is subscribed to the attached instrument, and acknowledge he executed the same.

Witness my hand an official seal.

Notary Public

My Commission Expires: March 6, 2019

InQwire Diagnostic Guide Wires with accessories

INSTRUCTIONS FOR USE

Dan Clark
Director, Regulatory Affairs
Merit Medical Systems, Inc.
September 13, 2017

Indications: Merit Medical guide wires are fabricated from high quality stainless steel utilizing sophisticated construction process and are available with or without a PTFE coating. Guide wires are supplied sterile, non-pyrogenic, and are intended for single use only.

Medical guide wires are packaged in a plastic hoop, which is fitted with a luer hub. Packaging is provided to facilitate compliance with the manufacturer recommended instructions that the wire be flushed with saline or heparinized saline prior to use. (See directions for use - Note).

Intended Use: The Merit Diagnostic Guide Wires are intended in the peripheral vascular system to facilitate the placement of devices during diagnostic and interventional procedures.

Contraindications: InQwire diagnostic guide wires are contraindicated for use in the coronary and cerebral vasculature.

Precautions: Angiography should be undertaken only by an experienced angiographer. **One Time Use Only.** Guide wires will collect blood and other foreign material in their lumens; neither autoclaving nor ultrasonic cleaning will completely remove foreign material, therefore guide wires are recommended for one time use.


PRECAUTION STATEMENT

Single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or reuse may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or reuse may also create a risk of contamination of the device and/or cause patient injury or cross-infection, including, but not limited to, the transmission of infectious agents from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Inspect all guide wires prior to use. Do not use any unit if the package is open or damaged.

Use aseptic technique during removal from the package and during use.

Guide wires are secured in the hoop dispenser by the J-tip straightener. To avoid damaging the guide wire during removal from the flush hoop, grasp the J-tip straightener near the base and slide it forward approximately 10 cm or until the J-tip straightener is no longer attached to the flush hoop adaptor. Holding both the guide wire and J-tip straightener, continue to dispense guide wire from the hoop.

Do not use excessive force to advance the moveable core while the guide wire is in a vessel. Placement with excessive force may cause coil penetration and vessel damage.

When withdrawing PTFE coated guide wires back through a metal needle. The sharp edge of the needle may scrape the coating. It is suggested that a catheter or PTFE vessel dilator be used to remove the needle as soon as the guide wire has reached the appropriate position.



During advancement of the catheter and guide wire within the aorta, it is recommended that the guide wire be removed at the appropriate level of the aorta.

Caution should be taken when manipulating a catheter during placement and withdrawal to prevent possible intravascular tissue damage. If resistance is felt when removing a guide wire from a catheter, the guide wire and catheter should be removed as a unit to prevent potential damage to the vessel wall.

The guide wire is a delicate instrument and remains the most fallible instrument used in a percutaneous procedure. Any time that a guide wire is used there is a possibility of thrombosis, embolism, vessel wall damage, and plaque dislodgement, which could result in myocardial infarction, cardiac arrhythmia, stroke or death. The physician should be familiar with the use of angiography products and the literature concerning the complications of angiography.

Flush Hoop Instruction Guide:

1. Attach flush filled syringe to flush port luer
2. Rotate syringe clockwise (as pictured)
3. Inject saline into hoop
4. Detach syringe from hoop
5. Dispense guide wire

Note: In order to reduce the potential of clot formation, it is recommended that the guide wire be flushed with saline or heparinized saline prior to use. Attach a filled syringe to the luer hub located at the end of the plastic hoop, flush several times. After flushing, remove guide wire from hoop and use as described above.


Directions for use: The following schematic shows a typical procedure for percutaneous entry utilizing the Seldinger technique. Variations in individual patient anatomy may preclude the utilization of this technique.

1. Vessel puncture with a two part needle
2. When using a two part needle remove stylet leaving cannula in place, insert flexible (distal) end of guide wire through cannula and into vessel.
3. "J" guide wires are shipped with "J" straightener to aid in the insertion of the wire into the puncture needle. Advance the straightener until 2-3 mm of the tip extends from the tip. Insert wire into hub and through needle. Remove "J" straightener proximally and discard.
4. Remove needle cannula leaving the guide wire within the lumen of the vessel.
5. Pass dilator or catheter over the guide wire directly into the vessel.
6. Carefully remove the guide wire leaving the catheter in place.

Note: Device is sterilized as stated on the package label.

How Supplied: Individually packaged 5-10 per box, refer to catalog for ordering information.

DIAGNOSTIC GUIDE WIRES

DEVICE SPECIFICATION

The following product configurations are available in the InQwire product family:

Table 1: InQwire Configurations

Characteristic	Description
Outer Diameter	0.018" (.46mm), 0.021" (.53mm), 0.025" (.64mm), 0.032" (.81mm) 0.035" (.89mm) or 0.038" (.97mm)
Length	50cm to 260cm
Mass/Weight	≤9.5g
Core Wire Distal Tip Configuration	Straight Configurations: Straight(S) 9cm distal flex segment Bentson(B) 23cm distal flex segment Bentson C(BC) 15cm distal flex segment Bentson Short Taper(BST) 10cm distal flex segment Bentson Exchange Amplatz Stiff Super Stiff Super Stiff Exchange Standard Exchange Rosen Exchange Rosen (RS) Newton Long Taper(NLT) 13.5cm distal flex segment Newton Long Long Taper(NLLT) 18.5cm distal flex segment Newton Long Long Long Taper(NLLLT) 23.5cm distal flex segment Curved Configurations: J-tip(J) Curved and Straight Configurations: J/S
Distal Tip radius Size	1.5mm, 3mm, 6mm, 15mm
Wire Body	Standard

Characteristic	
	Firm(F) Heavy Duty(HD)

The materials used in the manufacture of the InQwire guide wire are identified below.

Table 2: Component Material Specifications	
Component	Material
Core Wire	304V or 304 Stainless Steel
Wire Coil	304 Stainless Steel
Wire Coil Coating	Polytetrafluoroethylene (PTFE)
Safety Wire	304V or 304 Stainless Steel

POTENTIAL ADVERSE REACTIONS

Potential adverse reactions which may result from the improper use of a guide wire include, but are not limited to:

- Thrombus
- Emboli
- Arterial or venous vessel wall damage
- Plaque dislodgment
- Hematoma at the puncture site
- Infection
- Vessel perforation
- Vessel spasm
- Hemorrhage
- Vascular thrombosis

DEVICE COMPATIBILITY

The InQwire guide wire may be used with any device with an inner diameter appropriately sized to accommodate a 0.018" (.46mm), 0.021" (.53mm), 0.025" (.64mm), 0.032 inch (0,81mm), 0.035 inch (.89mm) or 0.038 (.97mm) guide wire.

HOW SUPPLIED

Individually packaged 5-10 per box with instruction for use.

Following separately packaged accessories are available to InQwire guide wire:

- RingMASTER guide wire basin.
- SeaDragon guide wire torque device.

ACCESSORIES:

1. SeaDragon guide wire torque device

DEVICE DESCRIPTION

Torque Devices are a guide wire accessory developed to assist in the manipulation and placement of the guide wire during interventional procedures. Interventional as used in this discussion includes all diagnostic and interventional procedures. They are designed for hydrophilic and other typical guide wires. The SeaDragon guide wire torque device consists of an ABS body, a polypropylene core, and a stainless steel spring. The core of the SeaDragon is a spring-loaded compression device. A tear dropshaped tunnel going through the core and the body of the SeaDragon align

when the core button is depressed and held. In this state, the guide wire is inserted through the tunnel and when the desired length of guide wire has been passed through the device, the core button is released. The alignment of the tunnel is broken and the core is displaced relative to the body. The spring provides sufficient pressure on the guide wire to prevent slipping and allows the clinician to use the body of the SeaDragon as a handle to manipulate the guide wire. The SeaDragon accepts guide wires with diameters ranging from 0.018" to 0.038".

INTENDED USE / INDICATION FOR USE



The SeaDragon torque device is intended for facilitating guide wire manipulation during interventional procedures.

Figure 1. SeaDragon Torque Device Outside perspective (A) and inner mechanism (B).

PERFORMANCE SPECIFICATIONS

Device	Characteristic	Specification
SeaDragon Torque Device	Length	1.378" < length < 2.363"
	Diameter of grip area	0.295" < diameter of grip area < 0.394"
	Device accommodation	Guide wires 0.018"-0.038" in diameter

rs, and other associated
es the customer with a cost
ng a procedure. The

Outer Diameter 7.83 inches (20cm).

INTENDED USE / INDICATION FOR USE

The RingMASTER guide wire basin is a containment device that holds guide wires, catheters, and other associated medical devices such as balloons and stents.

MATERIALS OF CONSTRUCTION: basin - Polypropylene Base Resin and Polypropylene Colorant Orange.

TRANSPORTATION, STORAGE AND OPERATION CONDITIONS

Transportation temperature: -35 °C to +55 °C

Storage temperature: -35 °C to +55 °C

Usage: 15 °C to 35 °C, relative humidity ≤80%, atm. Pressure 84.0 – 106.7 kPa (630-800 mm Hg).

Contact with the internal environment of the body at a temperature of 32 ° C - 42 ° C

SHELF LIFE

Shelf life - 3 year.

STERILIZATION


Ethylene oxide sterilization

DISPOSAL

Dispose the devices according to local environmental requirements.

Class of hazard for medical wastes by reference to specification of morphological composition (in accordance with Russian Sanitary Regulations and Standards 2.1.7.2790-10) Class B

Applicable symbols for InQwire guidewires and accessories

	This symbol means that the device fully complies with European Directive 93/42/EEC.
	Manufacturer
	Authorized Representative
	Single use
	Catalog Number
	Lot Number
	Sterilization: Ethylene-Oxide Gas
	Caution: U. S. federal law restricts this device to sale by or on the order of a physician.
	Use by

	Caution: Consult Accompanying Documents
	Do not use if package is damaged.
	Non-pyrogenic
Нетоксично	Нетоксично

WARRANTY

Limited Warranty and Disclaimer

Merit represents and warrants that the products shall materially conform to the specifications at the time of shipment to Customer and will remain materially free of defects in materials and workmanship for the period specified on the written limited warranty separately provided with each product. This warranty shall not apply to products that have been (i) modified, changed, repaired, refurbished, reprocessed, or altered by anyone other than Merit; (ii) subjected to misuse, mishandling, re-use, re-sterilization, accident, abuse, neglect or tampering; (iii) damaged by excessive physical, environmental or electrical stress; or that have had a serial number altered, defaced or removed; (iv) combined with any other product; (v) used outside the approved "Indications for Use" as cleared by the relevant competent authority, used contrary to the use outlined in the product specifications or in an application or environment for which such product was not designed or contemplated. Any technical product advice furnished by Merit is provided without duty or compensation and Merit assumes no obligation or liability, all such advice being given and accepted at Customer's risk. EXCEPT AS EXPRESSLY SET FORTH HEREIN MERIT MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, TO CUSTOMER, ITS END-USERS OR ANY THIRD PARTIES WITH RESPECT TO THE PRODUCTS. MERIT DISCLAIMS AND EXCLUDES ANY AND ALL OTHER EXPRESS, IMPLIED, AND STATUTORY WARRANTIES, INCLUDING WITHOUT LIMITATION, WARRANTIES AGAINST INFRINGEMENT, AND THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

Limitation of Liability

UNDER NO CIRCUMSTANCES SHALL MERIT BE LIABLE TO CUSTOMER OR ANY OTHER PERSON OR ENTITY FOR ANY PUNITIVE, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, WHETHER BASED UPON LOST GOODWILL, LOST RESALE PROFITS, WORK STOPPAGE, PRODUCT FAILURE, IMPAIRMENT OF OTHER PRODUCTS OR OTHERWISE AND WHETHER IN A CONTRACT, TORT OR OTHER ACTION FOR OR ARISING OUT OF BREACH OF WARRANTY, BREACH OF CONTRACT, DELAY, NEGLIGENCE, STRICT LIABILITY, OR OTHERWISE.

Customer Service

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095 U.S.A.
+ 1 801 253 1600

Merit Technologies LLC
115054, Moscow
Bakhrushina Str., 32, bldg. 1, Russia
+ 7 495 221 89 02

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