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
This Apostille is not valid for use anywhere within the United States of America, its territories or possessions.

This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears.

This Apostille does not certify the content of the document for which it was issued.


Apostille

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1. District of Columbia, United States of America
AHMAD MISKI
2. This public document has been signed by _____
NOTARY PUBLIC IN AND FOR THE DISTRICT OF COLUMBIA
3. acting in the capacity of _____
AHMAD MISKI, NOTARY PUBLIC IN AND FOR THE
DISTRICT OF COLUMBIA
4. bears the seal/stamp of _____

CERTIFIED

5. at Washington, D.C.
6. the 1 day of JUNE 2018
7. by Secretary of the District of Columbia
8. No. 482529
9. Seal/Stamp
10. Signature:

LAUREN C. VAUGHAN

Lauren C. Vaughan
Secretary of the District of Columbia

IS
COLOR MATCH BLUE FRONT
UNAUTHORIZED COPY PANTOGRAPH
COLOR BACKGROUND IS
ERASURE RESISTANT

IS
COLOR MATCH BLUE FRONT
UNAUTHORIZED COPY PANTOGRAPH
COLOR BACKGROUND IS
ERASURE RESISTANT

SECURITY PAPER
COLOR MATCH BLUE FRONT/WHITE BACK
UNAUTHORIZED COPY PANTOGRAPH
COLOR BACKGROUND IS
ERASURE RESISTANT

SECURITY PAPER
COLOR MATCH BLUE FRONT/WHITE BACK
UNAUTHORIZED COPY PANTOGRAPH
COLOR BACKGROUND IS
ERASURE RESISTANT


May 29, 2018

To Whom it May Concern:


This is to certify that the attached document is the Dimension assay Instructions for Use with Addendum for:

Products: Dimension CHK Flex reagent cartridge

Sincerely,

Karlyn Kellogg, Regulatory Technical Specialist,
Siemens Healthcare Diagnostics Inc.

Vickie M. Coulter-Phillips,
Notary Public

Vickie M. Coulter-Phillips
Notary Public
State of Delaware
My Commission Expires
February 10, 2020

Only for registration in Russia

Disclosure and Confidential Statement as of 5/29/2018

This document contains information which is strictly confidential and proprietary. It is to be used exclusively for registration by Siemens Regulatory Representatives, Local Agency Reviewers and Health Authorities. The distribution, copying or publication of any information in this document is prohibited.

Siemens Healthcare Diagnostics Inc.

P.O. Box 6101
Newark, DE 19714-6101

1-800-242-3233
www.siemens.com/diagnostics

REF DF179

SIEMENS

Biochemical analyzer Dimension

Cartridge with reagent CHK (for check-up) for biochemical analyzers series "Dimension" (Dimension CHK Flex reagent Cartridge (CHK))

Scope of delivery:

1. Cartridge with reagent (for check-up) (Dimension CHK Flex reagent Cartridge (CHK)) – 8 pcs
2. Package insert;
3. Instruction for use.

Application

The cartridge is designed to verify the proper functioning of sampling systems and liquid reagents as well as photometric measurement system, which are a part of the analyzers of the "Dimension" series. It is used for daily checks of the analyzers of the "Dimension" series.

Scope

Clinical laboratory diagnostics.

Contraindications

Not applicable for this medical device.

Possible side effects when using the medical device

Not applicable for this medical device.

Operator (user) requirements

Medical device should be used only by a qualified, trained (in the field of clinical laboratory diagnostics) personnel.

Description

In the test for optical density measurement, red organic dye is used as the sample and reagent to check for reliable functioning of the reagents and samples dosing system as well as photometric system of measurement on the analyzer for clinical biochemistry Dimension.

Medical device content

Reagent wells ^a	Ingredient	Concentration ^b
1–6	Ponceau S	0.14 mg/ml
7–8	empty	-

^a Cells are numbered sequentially starting from the broad side of the cartridge.

^b Nominal value of the test during manufacturing. Only for in vitro diagnostics.

Cartridge CHK Flex contains stable liquid reagent consisting of red organic dye (Ponceau S), buffer solution, and preservative.

Product/ingredient name	Identifiers	Siemens Healthineers Diagnostics I
2-methyl-4-isothiazolin-3-one	EC: 220-239-6 CAS: 2682-20-4	500 GBC Drive PO Box 6101 Newark, DE 19714 <0.1

Safety measures

Only for in vitro diagnostics. Only for use by professionals.

Cartridge should only be used by trained personnel. Special protective clothing, gloves, and eye/face protection is required, operation must be carried out in accordance with the requirements of good laboratory practice.

Precautions associated with the procedure

Prior to using, read the entire manual. Follow the instructions carefully, failure to follow the instructions may lead to deviation of test results.

Do not use the cartridge if packaging is damaged.

The use of the cartridge for analysis of samples different from those which are permitted for use in this test system may lead to deviation of test results.

Do not use the cartridge after the expiration date.

Always check the expiration date prior to testing.

Used cuvettes contain body fluids, so treat them with due caution to avoid contact with skin and ingestion.

Environmental protection

When handled, transported, or stored properly, this device causes no harm to the environment.

Waste disposal procedure

If possible, avoid waste generation or minimize its quantity. When disposing of this product, solutions and derivative products requirements of legislation on environmental protection and waste disposal as well as any requirements of local authorities should be followed. Disposal of excess product and products not subject to recycling should be carried out by a contractor organization licensed to recycle waste. Dumping of untreated waste into the sewer is prohibited, except in full conformity with the requirements of all authorities with appropriate authorization.

Product classification may meet hazardous waste criteria.

If possible, avoid waste generation or minimize its quantity. Leftover packaging is subject to recycling. Incineration or waste burial should only be considered when recycling is not feasible.

This material and its container must be disposed of in a safe way. Caution should be taken when handling empty containers that have not been cleaned or rinsed. Empty containers and packing materials may contain product residues. Avoid spreading of spilled material as well as its entering into soil, watercourses, sinks, or

drains.

Preparation of the reagent working solution

The solution is liquid and ready for use.

Expiry date:

12 months

Expiry date is indicated on the packaging of individual sealed reagent cartridges. After opening, the cartridges are stable for 30 days.

Storage conditions

Storage temperature: (+2...+8) °C

Relative humidity: 10% to 80% noncondensing.

Pressure: atmospheric.

To be stored in accordance with local regulations. Store in original container protected from direct sunlight, in a cool dry place with good ventilation, away from incompatible materials, food, and beverages. Prior to using, the container should be tightly sealed. Opened containers should be resealed and kept upright to prevent leakage. Do not store in containers without labeling. Use appropriate measures to prevent environment pollution.

Do not ingest. Avoid contact with eyes, skin, or clothing. Avoid inhalation of vapor or mist. Avoid release to the environment. Store in the original container or an approved alternative container made of compatible material; seal tightly when not in use. Empty containers contain product residues and may be hazardous. Do not reuse container.

It is prohibited to eat, drink, or smoke in areas where this material is being handled or in places of its storage and processing. Before eating, drinking, or smoking, workers should wash their face and hands. Before entering the eating area, remove contaminated clothing and protective equipment.

Sealed cartridges are stable until the expiry date stated on the product label at a storage temperature of (+2...+8) °C. After opening, store in a dark place.

Transportation conditions

Transportation should be performed all types of the covered transport, observing temperature condition of storage of the reagents.

Transportation within customer premises: transportation should always be carried out in sealed and protected containers in upright position. Ensure that persons transporting the product know what to do in case of damage or leakage of the product.

Stability during transportation: during transportation, the cartridge is stable under storage conditions until the expiration date stated on the packaging.

Relative humidity: 10% to 80% noncondensing.

Pressure: atmospheric.

The stability of an open cell

3 days for cells 1–6

Collection of samples (in the system without NM module)

Before removing the contents of one cell with the pipette and placing it in the tube for samples, gently invert CHK Flex reagent cartridge.

CHK cartridge is intended only to provide a sample—it should not be installed on the analyzer for testing of the reagents.

Procedure

Provided materials

CHK cartridge with reagent (for check-up) for biochemical analyzers series "Dimension" (Dimension CHK Flex reagent Cartridge (CHK)) Cat. No. DF179

Required (but not provided) materials

None.

Testing stages

Sample collection, addition of the reagent, mixing, separation, processing, and printing of results are carried out automatically by the Dimension system. Further details are provided in the Dimension analyzer's operator manual. Container for the collection of samples should be of appropriate size so that it could contain a control volume and have some free space; there is no need to completely fill the container.

Calibration

Test for optical density measurement is used to verify performance of the analyzer. Calibration is not needed.

Quality control

Use of CHK test should be conducted as part of system check-up during daily maintenance. Check-up consists of 5 repetitions of various tests that allow ensuring that the sampler and the subsystem for work with liquid reagents are functioning normally.

For the system without the NM module, to start the system check-up, it is necessary for the operator to manually prepare a CHK reagent sample cup. Instruments with the NM module are configured to automatically check the system and do not require manual preparation of the CHK fluid.

Results

Many different tests are used as part of the system check-up, and each of them uses a different system for adding fluid. In each test, the difference of optical density between the two wavelengths of Ponceau S, diluted as a sample or a reagent, is measured by means of the endpoint method with two filters (510 nm

and 600 nm). The results of these measurements must meet the parameters set out in the operator's manual. If the parameters are not met, this can serve as an indicator of system malfunction.

Test	Condition
Temperature	37 °C
Wavelength	510 nm, 600 nm
Measurement type	Bichromatic, endpoint

Specification of system check-up for CHK

System check-up		Specifications
Photometer		-2.5 mAU to +2.5 mAU only for 293 nm filter -1.5 mAU to +1.5 mAU for all other filters
LOCI		All values are in the range from 0.001 kE to 0.144 kE
Mean		Std. dev.
Reagent holder No. 1	Certified value provided on the CHK package, ±15 mAU	≤ 3.8
Reagent holder No. 2	Certified value provided on the CHK package, ±15 mAU	≤ 3.8
Sampler	10% from the certified value provided on the CHK package, ±2 mAU	≤ 1.6
GIA wash system	10% from the certified value provided on the CHK package, ±4 mAU	≤ 1.6
RMS	Certified value provided on the CHK package, ±15 mAU	≤ 3.8
Sampler probe cleaner	True	
Reagent intake probe cleaner	True	

Limitations of the procedure

Information system of the instrument contains error messages to warn the operator about certain malfunctions. Prior to proceeding, read the Dimension analyzer's operator manual to determine the malfunction.

Warranty

This device has been designed, manufactured, tested, and packaged in compliance with all appropriate precautions. Siemens Healthcare Diagnostics Inc. provides warranty for the device until the expiry date, provided that its packing is not damaged. Manufacturer guarantees safety of the medical device, absence of unacceptable risk of harm to the human life,

health, and the environment when the cartridge is used as intended in conditions stipulated by the manufacturer.

In case of any issues related to manufacturer's warranty obligations, contact the official representative of the manufacturer in Russia (Siemens Healthcare LLC).

Legend

Symbol	Description
	Use before
	Batch code
	Catalog number
	Manufacturer
	Authorized representative in the European Community
	Contains enough for conducting <n> tests
	Medical device for in vitro diagnostics
	Temperature range
	Refer to instructions for use

Dimension and Flex are brand names of the Siemens Healthcare Diagnostics Inc.

©2008 Siemens Healthcare Diagnostics Inc.
All rights reserved.


Siemens Healthcare Diagnostics Inc. 500 GBC Drive, P.O. Box 6101, Newark, Delaware 19714-6101, USA

Authorized representative in Russia:
Siemens Healthcare LLC
96 Dubinskaya St.,
Moscow 115093, Russia
Tel.: (495) 737-12-52
Fax: (495) 737-13-20
<https://www.healthcare.siemens.ru/>


Based upon available information

US Legalization

I, Hamadi Hamu, Vice President, hereby
certify that is an original document.

Authorized signature

District of Columbia: SS
Subscribed and sworn to
before me,
this 31 day of May, 2018

Ahmad Miski, Notary Public,
District of Columbia

Notary Public - District of Columbia

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