certified copy

DR. CHRISTINE VON HERTZBERG
NOTARIN IN NEUSS

I hereby certify that the below is a true copy of the copy.

Neuss, 18 May 2017

Dr. Christine von Hertzberg, Notary Public,
with official seat in Neuss, Germany

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3M™ Durapore™ Surgical Tape

Description:

3M™ Micropore™ Surgical Tape is general-purpose gentle tape used to secure dressings, tubing and devices to patients.

3M™ Transpore™ Surgical Tape is a transparent, perforated, latex-free medical adhesive tape. 3M™ Transpore™ Surgical Tape can be used in combination with other devices such as tubes, primarily for secural purpose. These uses are typically in areas with oversight by medical professionals.

3M™ Durapore™ Surgical Tape is a general purpose tape with high strength and adhesion. The backing is coated with a skin friendly acrylic based adhesive.

Indications for Use

The surgical tapes are indicated for use as non-invasive devices which come into contact with intact skin and are intended to be used to secure most dressings, tubing and appliances to patient.

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Directions for Use

1. Always apply tape to skin that is dry, clean and free from soaps, chemicals and oils. If needed, apply a skin protector and allow it to dry before using tape. Avoid using tincture of benzoin, as it can be a skin irritant.
2. Always apply the tape without tension and smooth from the center out. Allow at least one inch of tape on the skin around dressing edges.
3. Do not secure one end of the tape and apply tension while securing the other end of the tape. This may induce skin trauma in the form of blistering or skin tearing.
4. Maximize adhesion by firmly pressing or rubbing the tape into place on skin, tubing or appliances.

Tape Removal Technique:

1. Loosen all ends of the tape strips, or edges of dressings. Grasp the full width of the tape and slowly and gently pull the tape back over itself toward the wound to minimize disruption of healing tissue.
2. Tape should be removed in the direction of hair growth whenever possible, if this will not disturb the wound area. Keep the tape close to the skin surface as you pull it back, and support the skin immediately adjacent to the strip being removed.

Storage/Shelf Life

3M™ Micropore™ Surgical Tape (US-made), 3M™ Transpore™ Surgical Tape,

These products are expected to perform for a minimum of 5 years, and after that the product should be discarded.

3M™ Micropore™ Surgical Tape (German-made), 3M™ Durapore™ Surgical Tape

German-made Micropore has a shelf-life of 3 (three) years.

These medical devices aren't subject to maintenance and repair.

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Cat Nos. 1037, 1040, 1050, 1051,

2037, 2040, 2045, 2050, 2051

Authorised representative of the
manufacturer in Russia:

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3M™ Coban™ Self-Adherent Wrap

Description: 3M™ Coban™ Self-Adherent Wrap is a cohesive elastic wrap constructed from a nonwoven material and elastic fibers. The cohesive properties allow the wrap to stick to itself but not to other materials or skin. The product is supplied non-sterile.

Indications for Use

3M™ Coban™ Self-Adherent Wrap is intended for use as an elastic wrap to provide compression, support, or to secure dressings or devices.

Contraindications

No contraindications have been identified for these products.

Potential complications

No potential complications have been identified for these products.

Warnings

1. Coban self-adherent wrap will not slip or loosen after application. DO NOT wrap too tightly to avoid impairing circulation and causing serious injury. The risk of impaired circulation is higher for children.
2. Monitor wrapped area frequently (no less than every 24 hours) for signs of swelling, discoloration, pain, numbness, tingling, or other changes in sensation. If these symptoms occur, remove wrap immediately and contact your health care provider.
3. DO NOT place Coban self-adherent wrap in direct contact with an open wound because it may lead to infection. Breached skin shall be covered with an wound appropriate dressing before applying the Coban self-adherent wrap.

Caution


This product contains natural rubber latex which may cause allergic reactions.

Precautions

1. When used as a component in a compression wrap system, it is important to ensure adequate arterial blood flow to avoid injury secondary to decreased circulation. Coban self-adherent wrap should be used as a component of a compression wrap system only under the supervision of a wound care specialist.
2. Do Not Reuse. Reuse may result in compromising product integrity and / or lead to device failure.

Directions for Use

1. Unwind approximately 12 inches of Coban self-adherent wrap from the roll and allow the wrap to relax.

2. Without stretching, apply one full wrap and overlap. Press overlapped area lightly to keep the end in place. See Figure 1.

3. Proceed with wrapping, applying tension to the wrap only to obtain desired compression. Always unwind more Coban self-adherent wrap from the roll before making a wrap so that the product is not applied with too much tension. See Figure 2.

4. When application is complete, cut off excess Coban self-adherent wrap. Secure the end of the application by gently pressing.

Storage/Shelf Life

3M™ Coban™ Self-Adherent Wrap has a shelf-life of 3 (three) years.

This product should be stored at room temperature. Avoid excessive heat and humidity. For best performance, store Coban self-adherent wrap in its individual package until use. For shelf life, refer to the date printed on each primary carton.

This medical device isn't subject to maintenance and repair.

Utilization

Dispose of the device in accordance with applicable local/regional/national/international regulations. Medical devices in Russian Federation must be disposed of in the local public entities in accordance with SanPiN 2.1.7.2790-10 "Sanitary-epidemiological requirements for the handling of medical waste" to the waste class B).

Transport conditions

Shelf life of the products is guaranteed as labeled, if stored according to the recommended storage conditions. During transportation the products may experience short term temperature variations without affecting product safety and performance.

Transportation of medical device is carried out by all type of transport.

Explanation of Symbols

Caution: This product contains natural rubber latex which may cause allergic reactions.

3M Deutschland GmbH
Health Care Business
Carl-Schurz-Str. 1
41453 Neuss
Germany

3M Health Care
2510 Conway Ave.
St. Paul, MN 55144 USA
1-800-228-3957
www.3M.com/C3CD
Cat Nos. 1037, 1040, 1050, 1051,
2037, 2040, 2045, 2050, 2051

Authorised representative of the
manufacturer in Russia:
ZAO "3M Russia", Russia 121614,
Moscow, Krylatskaya str. bldg.17/3
T.+74957847474

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