The Swiss Instruments

- DE** Bedienungsanleitung
- EN** Operation Manual
- FR** Mode d'emploi
- IT** Istruzione d'uso
- ES** Instrucciones de uso

HighSurg 30

Motor System for Medical Surgery

Date: 1.6.2017

Leszek Kazana, Area Sales Manager

L. Kazana

Nouvag AG St.Gallerstr. 23-25 CH-9403 Goldach
Phone +41 (0)71 846 66 00 Fax +41 (0)71 846 66 70
info@nouvag.com www.nouvag.com

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Congratulations on your purchase of a NOUVAG AG the "HighSurg 30 Motor System for Medical Surgery". Thank you for the confidence shown in our products. Please consult the instruction manual for the use and maintenance of the device in order to ensure that it will function properly and efficiently for many years.

• **Please read instructions carefully before operating!**

Date: 1.6.2017

Leszek Kazana, Area Sales Manager

Table of Contents

Introduction	3
1. INTENDED USE AND OPERATION	4
2. INDICATIONS	4
3. CONTRAINDICATIONS	4
4. SIDE EFFECTS	4
5. OPERATING CONDITIONS	4
6. DEVISE OVERVIEW	7
7. TECHNICAL DATA	12
8. DESIGN AND OPERATION	15
1. Operating principle	15
2. Method of application	15
3. Startup of HighSurg 30 Motor System for Medical Surgery	16
4. Operation. HighSurg 30 Motor System for Medical Surgery	19
5. How to use handpieces	34
9. SCOPE OF DELIVERY	39
10. CLEANING OF DEVICE COMPONENTS	40
11. MAINTENANCE	47
1. Replacing the control unit fuse	47
2. Safety inspections	51
12. SAFETY REQUIREMENTS	51
1. Manufacturer's EC declaration	51
2. Integrated peristaltic tube pump	52
3. Modification and misuse	52
4. Essential requirements	52
5. During use	52
6. Safety inspections	53
7. Environment protection requirements during MD operation	53
13. MEDICAL DEVICE CLASSIFICATION	53
14. DEVICE PACKAGING AND MARKING	54
1. Marking	54
2. Packaging	57
15. TRANSPORTATION AND STORAGE CONDITIONS	71
16. SERVICE LIFE	71
17. DISPOSAL	71
18. WARRANTY COVERAGE	72
APPENDIX 1. Device components	73
APPENDIX 2. Tables of device components	74
APPENDIX 3. EMC Manufacturer's Declaration of Conformity	82

INTRODUCTION

Surgery is one of the oldest branches of medicine which is constantly developing and improving. Expertise and specific skills of the surgeon, however important in invasive methods of treatment are sometimes useless without sophisticated surgical devices. As part of new generation of Nuvag AG surgery motor systems, the "HighSurg 30 Motor System for Medical Surgery" is designed for a wide range of applications such as ENT surgery (e.g. Otorrhinology, Head Surgery) and Maxillofacial Surgery. The main idea of developing and designing of this device was patients treatment using minimally invasive surgery techniques. The result of this development is an MD with electronic motor with control unit set operation mode and smooth instrument operation with a speed up to 80'000 rpm.

The 2 motor device allows the surgeon to concentrate on surgical procedures saving him change of instruments. High precision and accuracy of incisions and possibility to use continuous or oscillating shaver are specific features of the device. There is a broad selection of shaver instruments for various applications. The device has an integrated peristaltic pump to cool the treated tissue and used instruments. Wide LCD display of the control unit is designed in such a way that one can easily read all the information at a glance. The display provides real time rotation speed and torque (graphically) values. Pedal, included in scope of delivery, controls the device switch on/off, peristaltic pump and electronic motor operation.

Pic. 1. The "HighSurg 30 Motor System for Medical Surgery".

1. INTENDED USE

The "HighSurg 30 Motor System for Medical Surgery" is designed for drilling, milling, screwing and sawing of the bone as well as for shaver procedures in ENT surgery.

The MD is used in ENT surgery (e.g. Otosurgery, Head Surgery) and in oral and Maxillofacial Surgery.

2. INDICATIONS

HighSurg 30 Motor System for Medical Surgery is specially developed for modeling bone tissue and soft and hard tissue resection during surgery procedures in different fields such as:

- Surgical Restoration and Plastic Surgery
- Spinal Surgery
- ENT Surgery
- Neurological Surgery
- Orthopedic/Arthroscopy

3. CONTRAINDICATIONS

Relative or absolute contraindications may result from specific cases when use of cutting surgical electric powered instruments hazards the patient's health considerably.

4. SIDE EFFECT

HighSurg 30 Motor System for Medical Surgery operating must be conducted under supervision of trained medical personnel. Constant treatment and procedure surveillance makes it possible to duly answer side effects and complications should such occur.

Electric powered medical devices with integrated cooling systems are widely used in surgical procedures. Acceptable side effects are sensation of pain, swelling and irritation.

Adverse anesthesia reaction is possible. This procedure must be conducted under supervision of trained medical personnel.

5. OPERATING CONDITIONS

For correct "HighSurg 30 Motor System for Medical Surgery" operation, the following requirements must be met:

Ambient conditions

	Transport and storage:	Operation:
Relative humidity:	Max 90%	30- 80 %
Temperature:	0 °C to 60 °C	10 °C to 45 °C
Atmospheric pressure:	700 hPa to 1060 hPa	800 hPa to 1060 hPa

	Please carefully read these instructions and instructions for other devices used during surgery procedures.		The HighSurg 30 Motor System for Medical Surgery may only be operated by qualified and trained personnel. All requirements and directions of this manual must be followed.
	Device preparation must be conducted according to operation manual to assure safety of patients, users and third parties.		Improper use or repair of the device and failure to observe these instructions relieve us from any obligation arising from warranty provisions or other claims.
	The use of third party products is the responsibility of the operator. Functionality and patient safety cannot be guaranteed with third party accessories.		Prior to using this device, before startup and before operating, the user must always ensure that the device and accessories are in good working order and are clean, sterile and operational.
	Keep the device out of direct sunlight, central heating radiators, excessive dust, humidity, vibration and accidental mechanical damage.		Repairs may only be performed by authorized NOUVAG AG service technicians.
	The integrated peristaltic pump is used to cool tissue in order to prevent damage to the tissue. It may only be operated with watery solutions, such as 0.9% Sodium Chloride irrigation solution or "Ringer" solution (not included in delivery). Supplying medication using the integrated pump is expressly prohibited.		Use Nou-Clean spray (not included in delivery) for maintenance and care of the motor, handpiece and contra angles. Using other care products can result in malfunction and/ or cause the warranty to be revoked.
	The device is not sterile on delivery. All sterilizable parts must be sterilized before use (see Chapter "Cleaning, disinfection and sterilization").		Do not use device in the vicinity of flammable mixtures!
	Ventilation slots of the devices must not be obstructed in order to prevent overheating.		Check the voltage settings of the device. Make sure they match the supply voltage.

	Never operate the damping mechanism of the handpiece or contra angles while the system is running. This could result in instrument damage.		Operating of the HighSurg 30 Motor System for Medical Surgery must be conducted under constant surveillance of qualified medical personnel.
	Handpieces and contra angles may only be attached when the electric motor is not running.		To prevent injury, never touch drill bits or burs while they are still running.
	At operation at a patient, make sure to produce as less as possible frictional heat. By excessive thermal influence on the tissue, the tissue easily gets necrotized. The heat generation is directly proportional to the speed delivery of the instrument and the contact pressure delivered by the instrument		

Modification and misuse

- Modification or misuse of the HighSurg 30 Motor System for Medical Surgery are prohibited. The manufacturer is not liable for any damages resulting from unauthorized modifications or manipulations. The warranty will be cancelled.
- Use of the HighSurg 30 Motor System for Medical Surgery outside the indications described in Section "Intended Use" is prohibited. The user or operator is solely responsible for any such use.

Operator skills:

Education:

- medical education at university level

Operator's health:

- surgeon should be with normal physical and mental abilities

Language command:


- English or
- language of the Operation Manual

Medical use:

The HighSurg 30 Motor System for Medical Surgery must be operated by qualified dental surgeons in appropriately equipped surgical rooms in dental departments of medical institutions. MD operation is stable at ambient temperature from 10°C to 45°C, relative humidity 30% - 80%. MD is resistant to exposure to biological liquids and body tissue contacting device during surgical operation. MD is resistant to pre-sterilization cleaning, disinfection and sterilization. **When operating medical devices, the user must always wear gloves, protective goggles and medical mask (not included in delivery).**

6. DEVICE OVERVIEW

Front view


- | | |
|--|---|
| 1. Pedal socket, on the device rear panel | 14. Irrigation tubing set (1706 tube 2m., sterile (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010)) |
| 2. Vario pedal | |
| 3. Handpiece (not included in delivery) | 15. Stand for irrigation fluid bottle with irrigation liquid (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010) |
| 4. Clip for tubing sett attachment at handpieces and contra angles | 16. Drip chamber |
| | 17. Bleed valve |
| 5. Handpiece cradle (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010) | 18. Bottle with irrigation fluid |
| 6. Electronic motor 21 | |
| 7. Clip for tubing sett attachment at motor 21 cable | 19. Type plate with type designation, reference number, serial number, information on power supply and device fuse |
| 8.1 1 indicator light for each motor | |

HighSurg 30 Operation manual, Ref. 31666, V13/15

- | | | | |
|-----|--|-----|--|
| 9.2 | 2 motor sockets | 20. | Power plug socket |
| 10. | Operating panel | 21. | Main switch |
| 11. | Control unit display | 22. | Fuse compartment |
| 12. | Release key for irrigation tubing set bracket(
(Marketing Authorization № FSZ (ФЦ3) 2010/07320 dated
05.07.2010) | 23. | Sspray nozzle for electronic motor
maintenance |
| 13. | Peristaltic pump | 24. | Nou-Clean spray nozzle for handpiece
maintenance and care |

Rear view

The operation procedure is the following

On the front panel of MD there are 2 electronic motor sockets. The electronic motors are connected to the handpiece with Intra Coupling (ISO3964) connector. Depending on the requirements, various handpieces can be connected to the electronic motor. On the rear panel of the device there is a pedal socket. Vario pedal connected to the unit with Intra Coupling (ISO3964). Electronic motor rotation speed, rotation direction and switching on the pump (built-in pump) of cooling system are controlled with the pedal treadplate. Bottle with irrigation fluid (not included in delivery) is hung on a bottle stand. Irrigation fluid is supplied to the handpiece during operation through the tubing set (2m., sterile). The HighSurg 30 Motor System for Medical Surgery is equipped with control panel with all necessary control functions on it. On the control panel of the device there is an LCD display featuring parameters related to device operation. Operating parameters are set with the keys on the panel. On the rear panel there is a main switch ON/OFF "O/I". Pushing "OFF" the device can be switched off at any moment irrespective of routine shut down procedure.

To prepare MD for operation do the following:

1. Plug in the electronic motor cable into the motor socket.
2. If required, plug in the electronic motor 2 cable into the motor socket 2 on the control unit panel

Pic. 2. Two 21 electric motor sockets on the control unit.

3. Plug in the pedal cable into the pedal socket on the rear panel of the control unit.

Pic. 3. Pedal socket on the rear panel of the control unit.

Pic. 4. The MD with power and pedal cables plugged in.

4. Insert the stand for irrigation fluid bottle into the stand holder on the control unit.

Pic. 5. Insertion of the stand for irrigation fluid bottle.

Detailed MD preparation procedures are given in chapter "Startup".

The MD parts and accessories description is given in the Appendix 2 of this Operation Manual

Table. 1. Intended use of parts and accessories of the HighSurg 30 Motor System for Medical Surgery.

No	Part name	Intended use
1.	Electronic Motor 21	torque transmission to handpiece
2.	Surgical Handpiece	used for drilling, milling, screwing and sawing of the bone for nasal and lower jaw surgery
3.	Surgical Contra Angle	used for drilling, milling, screwing and sawing of the bone

HighSurg 30 Operation manual, Ref. 31666, V13/15

		in frontal sinus, upper jaw and ear
4.	1510nou Vario IPX8 pedal	MD control (ON/OFF, built-in pump, electronic motor 21)
5.	6024 Tube set, sterile, 3m, (Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010)	irrigation fluid supply to handpiece and surgical area
6.	1770 Stand for irrigation fluid	irrigation liquid placement (not included in delivery)
7.	1170 Handpiece cradle	electronic motor and handpiece placement during surgical manipulations
8.	1873 Clip set for tube set attachment to the motor 10 pcs	tube set (6024 tube set, sterile, 3m)attachment to the motor
9.	1881 Clip set for tube set attachment to the handpiece 3 pcs	tube set (6024 tube set, sterile, 3m)attachment to the handpiece
10.	1958 Spray nozzle attachment for handpiece cleaning	adapters for cleaning handpieces with Nou-Clean spray (not included in delivery)
11.	1974 Spray nozzle attachment for electronic motor 21 cleaning	adapters for cleaning electronic motor with Nou-Clean spray (not included in delivery)

When operating this medical devices, the user must always wear gloves, goggles and protective mask. (not included in delivery, purchased separately). Handpieces are in direct contact with tissues and biological liquids of patients (saliva, blood, tissue particles, purulent effluent). Tubing set 6024 (sterile, 3m, ((Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010)) 10pcs/pack contact the patients indirectly. Also possible inadvertent contact with electronic motors 21 2098nou and 2099nou, 1881 clip set for tube set attachment to the handpiece 3 pcs/pack and 1873 Clip set for tube set attachment to the motor 10 pcs/pack.

Materials used:

1881 clip set for tube set attachment to the handpiece 3 pcs/pack is made of polyphenylsulfone (PPSU) Radel R5100; 1873 Clip set for tube set attachment to the motor 10 pcs/pack is made of polyacetal POM Delrin 500P; 1958 Spray nozzle attachment for handpiece cleaning and 1974 Spray nozzle attachment for electronic motor 21 cleaning are made of EN AW-6082 aluminum; 1170 Handpiece cradle is PP30%GF (polypropylene + 30% glass fiber); electronic motors 21 2098nou and 2099nou are made of medical grade stainless steel 1.4305; electronic motor cable is made of perfluoroalkoxy (PFA-silicon); tow- motor 3360 the

HighSurg 30 control unit with integrated peristaltic pump is made of ABS UL 94-VO plastic and the control unit buttons are polyacetal (POM).

7. TECHNICAL DATA

Essential technical data on The HighSurg 30 Motor System for Medical Surgery are given bellow.

Table 2. 3360 HighSurg 30 Control Unit Technical Data:

Voltage AC	100 V ~ / 115 V ~ / 230 V ~, 50-60 Hz
Fuse	2 fuses T3.15 AL, 250 VAC AC
Power consumption	120 VA
Number of electronic motors	2
Number of motor connectors	2
Number of pedals	1
Motor speed	300 – 50'000 rpm 300 – 80'000 rpm for "HighSpeed"
Handpieces for 2098nou electronic motor 21 "HighSpeed" 80'000 rpm	1:5 speed 1'500-240'000 rpm, torque 1 Ncm (not included in delivery) 1:3 speed 900-150'000 rpm, torque 2 Ncm (not included in delivery) 1:2 speed 600-100'000 rpm, torque 2 Ncm (not included in delivery) 1:1 speed 300-80'000 rpm, torque 6 Ncm Continuous Shaver speed 300-6'000 rpm, torque 3 Ncm (not included in delivery) Oscillating Shaver speed 50-3'000 rpm, oscillation period 0.3-3 sec (not included in delivery) Jacobs Chuck speed 80-2'100 rpm, torque 60 Ncm (not included in delivery) Craniotome speed 1000-50'000 rpm, torque 6 Ncm (not included in delivery) Perforator speed 80-1'000 rpm, torque 55 Ncm (not included in delivery)
Handpieces for 2099nou electronic motor 21 50'000 rpm	1:5 speed 1'500-240'000 rpm, torque 1-1 Ncm (not included in delivery) 1:3 speed 900-150'000 rpm, torque 1-2 Ncm (not included in delivery) 1:2 speed 600-100'000 rpm, torque 1-2 Ncm (not included in delivery) 1:1 speed 300-50'000 rpm, torque 1-6 Ncm CompassSaw speed 15'000 rpm, torque 6 Ncm (not included in delivery) Oscillating Saw speed 15'000 rpm, torque 6 Ncm (not included in delivery) Sagittal Saw speed 15'000 rpm, torque 6 Ncm (not included in delivery) Dermatome speed 8'000 rpm, torque 6 Ncm (not included in delivery)

HighSurg 30 Operation manual, Ref. 31666, V13/15

	Kirschner HP speed 500-2'800 rpm, torque 48 Ncm (not included in delivery) Tattoo HP fixed speed 9'000/10'000/11'000/12'000 rpm, torque 6 Ncm (not included in delivery) Continuous Shaver speed 50-3'000 rpm, torque 6 Ncm (not included in delivery) Oscillating Shaver speed 50-3'000 rpm, oscillation period 0.3-3 sec (not included in delivery) Jacobs Chuck speed 80-2'100 rpm, torque 60 Ncm (not included in delivery) Craniotome speed 1000-50'000 rpm, torque 6 Ncm (not included in delivery) Perforator speed 80-1'000 rpm, torque 55 Ncm (not included in delivery)
Electronic motor coupling	internal coupling ISO3964
Instrument connector	ISO 3964 Connector (DIN 13940)
Electronic motor cable length, Max.	3 m
Pedal cable length, Max.	3,15 m
Cooling system	Integrated peristaltic pump (built-in cooling pump) adjustable 1% - 100%
Volumetric displacement of pump	1 -100 ml/min
Operating pump noise level (low)	30 dB
Electrical protection class	I
Method of protection against electric shock	Typa B
Operating time rating, Max.	1 min (ON), 3 min (OFF)
Device type	hospital
User interface	LCD 95x52 mm Display
Electronic control	Microprocessor: input voltage– 230/115/110 V, programming voltage – 5 V
Software	Built-in: Version 1.13/Print V4.0, date 15.04.16, Class I (A)
Response time	2 sec
Software protection level	Classc B
Overall program size	18'425 lines of source code
Bit-timing frequency	8 MHz
Battery:	Back-up battery CR 2032
Type of battery	Lithium battery
Battery capacity	225mAh
Operation temperature	-30°C to +70°C
Dimensions Max.	Diameter: 20 mm Hight: 3.2 mm
Pedal	IPX8
Dimensions (WxLxH), Max.	260 x 250 x 110 mm
HighSurg 30 Motor System for Medical Surgery control unit net weight, Max.	3.7 kg

Table 3. Technical data: 2099nou 21 Electric Motor:

Motor speed	from 300 to 50'000 rpm
Motor type	BLDC-Motor with Hall sensor
Electronic motor torque, Max.	6 Ncm
Weight, Max.	300/280g with/without cable
Continuous sound pressure level	50dB
Rotor power output	80 W
Connection	E-Type (ISO 3964)
Dimensions (D x L), Max.	20 x 110
Cable length, Max.	3 m
Sterilization	Steam Max. temperature 135°C

Table 4. Technical data: 2098nou 21 HighSpeed Electric Motor:

Motor speed	from 300 to 80'000 rpm
Motor type	BLDC-Motor with Hall sensor
Electronic motor torque, Max.	3 Ncm
Weight, Max.	300/280g with/without cable
Continuous sound pressure level	50dB
Rotor power output	120 W
Connection	E-Type (ISO 3964)
Dimensions (D x L), Max.	20 x 110
Cable length, Max.	3 m
Sterilization	Steam Max. temperature 135°C

Table 5. Technical data: 1510nou Vario Control Pedal.

Type	Vario pedal with 3 keys
Number of speed modes	3 (low (100 – 400 rpm) intermediate (300 - 700) high (500 – 1000 rpm))
Protection against contact and continual submerging	IPX8
Weight, Max.	1.68 kg
Dimensions (W x L x H), mm, Max.	255 x 235 x 60
Pedal press force	40 N
Cable length, Max.	3,15 m
Keys functions	Yellow: sets rotation direction right/ left Green: program selection Blue: ON/OFF peristaltic pump

Table 6. Technical data: Handpieces.

Material	Stainless Steel 1.4305.
Surface roughness	Ra 0.8
Hardness of material	200MPa
Surface finish	Non-gloss
Burr pulling out force kgf	Min. 30
Instrument extracting/inserting force N	Min. 45
Thrust kg	Max. 3
Frictional force between handpiece and	25N

HighSurg 30 Operation manual, Ref. 31666, V13/15

electric motor	
Continuous sound pressure level	50 dB
Burr shank diameter	2,35 mm
Rate of water supply	Min 50 ml/ min. at 200kPa (2 Bar); 1 – 100 ml/ min.
Surgical area irrigation	coolant one channel supply
Number of sterilization cycles without quality losses	250
Sterilization	Steam Max. temperature 135°C
Average time of failure (hour) Min.	For straight handpieces - 500 For contra angles - 300

8. DESIGN AND OPERATION

1. Operating principle

This HighSurg 30 Motor System for Medical Surgery consists of control unit, electronic motor, control pedal and tube set for irrigation liquid. The MD is provided with two motor sockets which makes operator's work easier and brings better results. The control unit runs all other components of the MD, it is started by pressing the pedal. Operation parameters (speed, torque) are set using control unit display. Rotating instruments are latched in the handpiece, their rotation and required speed are rendered by electronic motor power. The MD has separate program for each motor, it provides continuous operation. Optimized electronic motor control system allows for precision manipulations in compliance with individually selected parameters to make every use of the device safe. Cooling of rotating tools and tissues is provided with adjustable built-in peristaltic pump equipped with tubing sets for coolant. Vario pedal is used to control peristaltic cooling pump, switches programs on and off and controls electric motor function. The MD is used for drilling, grinding and sawing bones. There are also special programs for shavers, saws, Dermatome, Kirschner HP, perforator, Craniotome, Tattoo HP and Jacobs handpieces (not included in delivery, ordered and supplied separately). Built-in irrigation system (pump) for sterile coolant supply reduces instrument heating and minimizes risk of tissue damage and inflectional complications.

2 Method of application

The procedure is conducted under local or general anaesthesia.

The electronic motor is attached to the handpiece equipped with a cutting blade or modeling or drilling bit. The pump is activated for cooling. Tissue is prepared with instruments. There is no direct contact between the electric motor and the patient.

Operation safety considerations:

Operating the HighSurg 30 Motor System for Medical Surgery means certain risk such as damage to blood vessels and other organs. This information should be available to the operator, because some of these risks are operator's (user's) responsibility. Other hazards are included in Nouvag AG Risk Management System.

3. Startup of the HighSurg 30 Motor System for Medical Surgery

1. Device setup

- Place the HighSurg 30 Motor System for Medical Surgery with accessories and instruments on an even, non slip surface and make sure you have good access to all controls.
- Do not allow the operating range of the device (including the cable and contra angle) to be compromised by limiting factors.
- The system display must be fully visible at all times.
- The pedal must be placed within stepping distance between the patient and the surgeon.
- It must be explicitly ensured that no objects can fall on the pedal.
- The power plug at the rear of the device must be accessible at all times.
- The motor ventilation slots must be kept clear in order to prevent the motor temperature from becoming excessive.

2. Connection to the power supply

Before plugging the power cable into the power socket for the first time, you must check the supply voltage setting next to the power switch.

If the voltage shown does not correspond to the local mains voltage, the grey fuse holder must be set to the correct voltage:

1. Switch off the device.
2. Unplug the power cable.
3. Use a screwdriver to open the fuse slot.
4. Remove the fuse holder.
5. Remove the grey fuse holder and reinsert it so that the local mains voltage setting is shown in the small window.
6. Slide the fuse holder back in and close the fuse slot.
7. Check the mains voltage shown on the fuse slot.
8. Plug the power cable back into the device.

In order to prevent the risk of electric shock, the device may only be connected to a power network with a PE protective ground conductor.

3. Device preparation

1. Sterilize the motor 21 (the motor is not sterile on delivery). If the motor has already been sterilized: when removing the motor from the sterile packaging, ensure that the sterile packaging is not damaged and that the sterility indicator confirms sterility (if no sterility indicator is provided, the sterile packaging must at least show the date on which the shelf life of the sterile item is due to expire). 
2. Insert the stand for the irrigation fluid into the stand holder.
3. Plug the motor plug of the electronic motor into one of the motor sockets.
4. Where appropriate, plug the motor plug of the second electronic motor into a motor socket.
5. Plug the Vario pedal plug into the pedal socket at the rear of the control unit.
6. Attach the sterilized handpiece to the electronic motor. Press the handpiece firmly onto the electronic motor unit until it clicks into position and make sure it is secure by moving it slightly in the opposite direction.
7. Assemble the tubing set: Decide whether you will use tubing set Ref. 6024 (sterile, 3m, (Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010)), for cooling of a single contra angle or if you need to use tubing set Ref.6025 (not included in delivery, supplied separately) with an integrated 3 way cock for cooling two contra angles at the same time when two motors are used.

Use only Nouvag AG tube set Ref. 6024 3m. sterile (Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010) and 6025, otherwise the function cannot be guaranteed.


Check the expiry date of the tubing set and ensure that the packaging is not damaged. Using non-sterile tubing set can result in serious infection and, in extreme cases, can be fatal.

When inserting the tubing set, watch the arrow on the cover of the pump compartment. It indicates the flow direction of the cooling liquid.

Do not regulate the amount of irrigation liquid using the roller clamp on the tubing set, with the HighSurg 30, this is regulated instead using the integrated pump. For this reason, make sure to open the roller clamp as far as it will go.

- A) Press release key for set bracket (on top of the control unit) to open the pump.
- B) The compartment with the integrated tubing bracket opens.
- C) Place the tubing set into the tubing bracket provided in such a way that the part of the tubing set with the spike exits the pump towards the rear of the device. Make sure the tubing is secure.
- D) With the tubing set inserted, press the compartment downwards until it clicks into place.

Pic 6. To access the pump, press release key by the display on the front of the control unit. It can be easily done with one hand, without inadvertent moving of the device


- 8. Insert the spike at the end of the tubing set into the irrigation fluid bottle and hang the bottle onto the stand.
- 9. Open the roller clamp on the tubing set as far as it will go.
- 10. Open the bleed valve beneath the drip chamber.
- 11. Connect the control unit to the power socket.

Ensure that the operating voltage setting corresponds to the local mains voltage.

4. External cooling system assembly

- A) Secure motor cable clips (1873) and handpiece clips (1881) to the irrigation tube.
- B) Connect the tube set with the external cooling tube of the instrument (Example: Shaver blade (not included in delivery)).
- C) Secure clips to the handpiece and to the motor.
If required, secure additional clips to the motor.

4. Operation. The HighSurg 30 Motor System for Medical Surgery.


Switching the device on and off

Pic. 7. Rear panel of the control unit.


The power switch I/O (at the rear) is used to switch the control unit on and off.
The device can be switched off at any time irrespective of any procedure for switching off the device.

Overview: control elements on operating panel

- A) **Display:** Shows the operating values (See chapter "Overview: Standard display").
- B) **"Program" key:** Selection of program 1 to 10. There are 10 selectable programs for each motor. (To facilitate the selection of programs, the number of selectable programs can be limited to the ones you need. Refer to "Configuration Menu", Level 1, Programs/Number of programs).
- C) **"Motor 1/2" key:** Switching between the connected motors (the green indicator lights beside the motor sockets indicate the active motor). Pressing the Motor key longer changes direction of rotation.
- D) **"Handpiece" key:** Choice of handpieces or contra angle (can be individually deactivated, refer to "Configuration Menu").
- E) **"Enter" key:** Enters the configuration menu by pressing long (refer to "Configuration Menu").
- F) **"Select" keys:**
 - By pressing the left "Select" key "Select ▲" the software version is shown.
 - By pressing the right "Select" key "Select ▼" the software version is shown.
 - By pressing both "Select" keys at the same time "Select ▲ + ▼" the programs will be reset to factory default settings.
 - The "Select" keys serve in the configuration menu to change values and parameters
"▲": value adjustment (upwards) "▼": value adjustment (downwards)
- G) **"Speed" keys:**
 - Restrict the maximum speed that can be selected using the pedal.
"▲": increases the maximum speed "▼": reduces the maximum speed
 - By pressing both "Speed" keys at the same time "Speed ▲ + ▼" the handpiece calibration will be started (refer to "Calibrating handpieces").
- H) **"Torque/Frequency" keys:**
 - Restricts the maximum torque.
"▲": increases the maximum torque "▼": reduces the maximum torque
 - With the shaver handpiece (not included in delivery) the key pair is used to set oscillation period, if oscillation mode is selected.
- I) **"Pump" keys:**
 - Adjusting the maximum pump (built-in cooling pump) flow rate, controlled by the pedal.
"▲" increases the maximum supply quantity "▼" reduces the maximum supply quantity
 - By pressing both "Pump" keys at the same time "Pump ▲ + ▼" the pump (built-in cooling pump) will be put on call, pressing again will switch it off.

Overview: Standard display

A) Information line

Information and error messages are displayed here. Display is illuminated red when error messages are displayed.

B) Program

Shows the selected program for the active motor.

C) Pump

The numerical value shows the pump flow rate in percent and the drop symbol together with the on/off indication shows if the pump is in standby mode or switched off.

D) Rotation direction of the motor

The arrow indicates the rotation direction set for the motor. The rotation direction can be changed by pressing the button  on the pedal or by long pressing of the "Motor 1/2" key at the operating panel. In Shaver mode both arrows are displayed at the same time, which indicates the oscillation mode of the shaver.

E) Current speed

At stand still of the motor the set maximum speed is displayed. As the motor starts to run when the pedal is pressed, the evolving speed is shown in real time.

F) Date

G) Clock

H) Motor

Motor Shows selected motor → see also green indicator light beside the motor socket.

I) Name of the handpiece or corresponding transmission ratio

Shows name of the handpiece used or the selected transmission ratio (see also "Selecting the handpiece or transmission ratio").

- J) **Speed range**
Shows speed range of the handpiece used.
- K) **Maximum torque**
Shows maximum torque setting.
- L) **Current torque**
Bar graph providing a graphical representation of the current torque. All bars active means max. torque reached.

The pump does not begin to operate until the motor has been activated by pressing the pedal

Adjusting the programs

The values for operation settings depend on the connected handpiece/ contra angle as well as the task to be performed.

Step 1: Selecting the handpiece or transmission ratio

Depending on the handpiece or contra angle attached to the motor, the corresponding transmission ratio has to be adapted accordingly by using the handpiece key.

Handpiece

Press the “Handpiece” key repeatedly until the name of the required handpiece with the corresponding transmission ratio is shown on the display. When the key is pressed constantly the handpiece or contra angles will be shown in fast forward mode.

Table 7. Table of possible Handpieces, available for the 50'000 rpm motor

Name of the handpiece with transmission ratio	Handpiece name on	Speed min. rpm	Speed max. rpm	Torque min. Ncm	Torque max. Ncm
--	----------------------	-------------------	-------------------	--------------------	--------------------

HighSurg 30 Operation manual, Ref. 31666, V13/15

	display				
Handpiece, 1:5(not included in delivery)	1:5	1500	240 000	1	1
Handpiece, 1:3(not included in delivery)	1:3	900	150 000	1	2
Handpiece, 1:2(not included in delivery)	1:2	600	100 000	1	2
Handpiece, 1:1	1:1	300	50 000	1	6
Compass Saw, MSS 5000 (not included in delivery)	Com. Saw	fix.. 15 000		fix. 6	
Oscillating Saw, OMS 5000 (not included in delivery)	Osc. Saw	fix.. 15 000		fix. 6	
Sagittal Saw, MOS 5000 (not included in delivery)	Sag. Saw	fix.. 15 000		fix. 6	
Dermatome (not included in delivery)	Dermatome	fix.. 8 000		fix. 6	
Kirschner Handpiece (not included in delivery)	Kirschner	500	2800	fix. 48	
Tattooing Handpiece (not included in delivery)	Tattoo Hp	fix. 9000/10000/11000/12000		fix. 6	
Jacobs Chuck (not included in delivery)	Jacobs Chuck	80	2100	fix. 60	
Perforator (not included in delivery)	Perforator	80	1000	fix. 55	
Kraniotome (not included in delivery)	Kraniotome	1000	50000	fix. 6	

Table 8. Table of possible Handpieces, available for the 80'000 rpm motor

Name of the handpiece with transmission ratio	Handpiece name on display	Speed min. rpm	Speed max. rpm	Torque min. Ncm	Torque max. Ncm
Handpiece, 1:5(not included in delivery)	1:5	1500	240 000	fix. 1	
Handpiece, 1:3(not included in delivery)	1:3	900	150 000	fix. 2	
Handpiece, 1:2(not included in delivery)	1:2	600	100 000	fix. 2	
Handpiece, 1:1	1:1	300	80 000	fix. 6	
Jacobs Chuck (not included in delivery)	Jacobs Chuck	80	2100	fix. 60	
Perforator (not included in delivery)	Perforator	80	1200	fix. 55	
Kraniotome (not included in delivery)	Kraniotome	1000	50000	fix. 6	

Table 9. Shaver functions table, available for 50 000 and 80 000 rpm motor

Name of the handpiece with transmission ratio	Display	Speed min. rpm	Speed max. rpm	Oscillation period from seconds	Oscillation period to seconds
Continuous Shaver (not included in delivery)	Cont. Shaver	300	6000		

HighSurg 30 Operation manual, Ref. 31666, V13/15

Oscillating Shaver (not included in delivery)	Osc.Shaver	300	5000	0,30	3
---	------------	-----	------	------	---

The oscillation period for oscillating shavers (not included in delivery) is determined by the speed. If the oscillation period is set too low for the set speed, it switches automatically to the next possible level. This ensures that the shaver knife makes a full rotation for the speed setting.

Handpieces that don't belong to one's own assortment can be deactivated in the Configuration Menu. (refer to "Configuration Menu") That makes finding the right handpiece in the handpiece list easier, because the list of available handpieces becomes shorter.

Step 2: Calibrating handpieces

To make sure the displayed parameters correspond with the actual, measurable parameters of the handpiece, it is recommended to calibrate each handpiece or contra angle on a regular basis.


It's a procedure as simple as it is important to guarantee safety and precision of each handpiece being used.

After you take care of all prior preparations such as sterilization, maintenance and care of handpieces, device preparation and the selection of the handpiece of use, the calibration procedure can be performed.

The calibration of handpieces or contra angles guarantees accurate torque. Due to wear as well as varying lubrication of the handpieces or contra angles, and lack of maintenance and care, the distribution of torque can vary widely.

1. Select the correct handpiece, corresponding to the one mounted on the motor and double check on the display for correct selection.
2. Hold motor with mounted handpiece or contra angle in your hand at safe distance to your body.
3. Press both "Speed" keys at the same time (**Speed ▲ + ▼**). "Testing the handpiece **XX**" is displayed.

4. Motor and handpiece start running and pass several speed cycles up to maximum speed.
5. After a tone is emitted the calibration is finished. The display shows the message "Handpiece is OK".

If a handpiece does not work within the parameters established by calibration even after cleaning and lubrication, the device displays a red-illuminated error message, "Handpiece **XX** is faulty".

This indicates soiling, wear or a technical defect. Such handpieces must be repaired or replaced.

The torque response of all handpieces is checked during calibration.

With the 1:1 handpiece, control of the device is also adjusted to altered handpiece conditions to remain with the torque response tolerances.


Step 3: Setting the speed

The possible speed range depends on the attached handpiece. The maximum speed within this speed range can be restricted to the required value.

Using the pedal, the speed can be varied from the minimum speed up to the maximum speed as set.

Setting the speed:

Press the «Speed» keys “▲” to increase or “▼” to decrease the maximum speed. When key is pressed continuously the speeds will be shown in fast forward mode.


The following handpieces run only with one specific speed, it cannot be changed. The specific rpm values of all handpieces are shown in the tables 7, 8.

- Micro Compass Saw MSS 5000 (not included in delivery)
- Micro Oscillating Saw OMS 5000 (not included in delivery)
- Micro Sagittal Saw MOS 5000 (not included in delivery)
- Dermatome (not included in delivery)
- Tattooing handpiece (4 steps) (not included in delivery)

Step 4: Setting the torque

Once the speed has been selected, the torque can be determined from the corresponding torque range.

Press the “Torque” keys “▲” to increase or “▼” to decrease the maximum torque. By pressing the key continuously the torques are shown in fast forward mode.

When the shaver handpiece is used in oscillation mode, the number of oscillations per second of the shaver knife (number of changes of direction of the rotating knife per second) is set with these two keys. In this function, the arrow displays both directions.


The following handpieces run on a fixed torque.

- Handpiece, 1:5 (1 Ncm) (not included in delivery)
- Micro Compass Saw, MSS 5000 (6 Ncm) (not included in delivery)
- Oscillating Micro Saw, OMS 5000 (6 Ncm) (not included in delivery)
- Micro Sagittal Saw, MOS 5000 (6 Ncm) (not included in delivery)
- Dermatome (6 Ncm) (not included in delivery)
- Kirschner Handpiece (48 Ncm) (not included in delivery)
- Tattooing Handpiece (6 Ncm) (not included in delivery)
- Jacobs Chuck (60/60 Ncm)* (not included in delivery)
- Perforator (55/55 Ncm)* (not included in delivery)
- Craniotome (6/6 Ncm)* (not included in delivery)

* The value before the slash refers to the use with the 50,000 rpm motor, the value after the slash to the 80,000 rpm motor.

Step 5: Setting the pump supply quantity

Press the "Pump" keys "▲" to increase or "▼" to decrease the pump supply quantity. By pressing the key continuously the values are fast forwarded.

Minimum and maximum of pump supply quantity as well as operating steps can be adjusted in the Configuration Menu (refer to « Configuration Menu, Parameter level 2, Pump»).

To activate or deactivate the pump, press both "Pump" keys at the same time, "Pump ▲ + ▼", or use the foot switch .

Torque limiter

The HighSurg 30 Motor System for Medical Surgery has an automatic torque limiter, which functions like a torque wrench.

If the attached instrument meets resistance, the torque increases as far as the maximum level, and the speed then reduces, to standstill if necessary. The torque on the instrument is maintained. If the load on the instrument reduces, the speed increases again up to the set maximum level.

This process can easily be seen on the display by means of the bar chart. The segments of the bar chart fill up as the resistance on the instrument increases. When the torque has reached the set maximum, that is, when all segments are visible, the speed reduces. As soon as the pressure on the instrument is reduced, the torque falls again. The bar chart on the display regresses and the instrument speed increases again.

Storing various programs

With the HighSurg 30, up to 10 different settings per motor can be set as fixed programs (program 1 to program 10). Which program is currently active can be seen on the display.

At switch-off, the settings made by user are automatically saved. This includes the following parameters:

- Handpiece/Transmission ratio
- Maximal speed
- Maximal torque
- Pump On/Off
- Pump performance
- Shaver oscillation period

To change a program go to the specific parameter and change it. When the device is switched off all parameters are saved in that program.

The number of storable programs can be limited in the configuration menu from 3 to 10 programs.

In the configuration menu it can be chosen what the display shows at start-up of the device. Choose between Program 1 with Motor 1 or the last used setting before the device was switched off.

Configuration Menu

In the Configuration Menu the user can customize the device after his favor. The parameters are organized in two levels. "Parameter Level 1" comprises all basic information and possible settings that affect the device itself. Parameter Level 2 includes all information and settings that involve the surgeon's work. In the configuration menu the following information can be obtained or parameters can be set according to one's own needs:

- Software version
- Mainboard serial number
- Date format US?
- Date/Time
- Battery voltage
- Display brightness
- Number of programs
- Behavior after switching on
- Operating hours HighSurg 30
- Operating hours Motor 1
- Operating hours Motor 2
- Operating hours of irrigation pump
- Error-messages (the last 8)
- Activation of available handpieces (15 possibilities)
- Speed limitation for each handpiece
- Pump behavior
- Motor behavior for 2 different preconfigured motor types
- Reset to factory default settings

Be cautious when changing parameters. Unusual behavior of instruments while operating may provoke false reactions and jeopardize the patient.
Every setting and the new behavior of the instrument has to be verified and tested.

Access to Configuration Menu:

- Press "**Enter**" key for 3 seconds until a long tone is emitted.
The right arrow ">" in the information line on the display indicates that you're in the Configuration Menu.

Enter

>Software
>Version V1.06bP

Changing certain parameters:

- Choose the desired parameter by pressing "**Select ▲**" or "**Select ▼**".
- Press "**Enter**" to activate setting mode.
The value of the parameter is now displayed in parentheses. To change the value press "**Select ▲**" or "**Select ▼**".

>Backlight
>brightness (0..10) [02]

>Backlight
>brightness (0..10) 10

- To confirm the changes made, press "Enter" for one sec., until a short tone is emitted. - To abort the settings made, press "Enter" briefly, the setting returns to its previous value, without emitting a tone.	
---	--

Exit from Configuration Menu:

- To leave the Configuration Menu press **"Enter"** for 3 seconds, until a long tone is emitted.

Table 10. Parameter Level 1

Group/Parameter	Rights	Factory	Definition
Software/Version	read	VX_XXXX	Shows current software version
Hardware/Serial number MB	read	XXXXXXXXXX	Shows serial number of the main board
Date-Time/Date format US	read/write	no	Set to US Date format
Date-Time/Date	read/write	-	Change current date
Date-Time/Date	read/write	-	Change current time
Battery/Voltage	read	e.g. 3120 mB	State of charge of battery (Capacity about 3100 mV, Low 1800 mV)
Backlight/brightness (0..10)	read/write	9	Display brightness, variable: 0, ..., 10
Program/Number of programs	read/write	10	Number of activated Programs: 3, ..., 10
Program/Power on at last program	read/write	yes	No: Display of Motor 1 with Program 1 Yes: Display of the last used program/motor
Рабочие часы/Панель управления Operating hours/Control Unit	read	0	Shows operating hours of HighSurg 30
Operating hours/Motor 1	read	0	Shows operating hours of motor 1
Operating hours/Motor 2	read	0	Shows operating hours of motor 2
Operating hours/Pump	read	0	Shows operating hours of pump
Error memory/ 1 – 8	read	0	8 Error messages in chronological order. <i>For the list of error codes please refer to chapter Maintenance.</i>

Table 11. Parameter Level 2

Values in Parameter Level 2 can only be adjusted after the password **"9403"** was entered. *The password cannot be changed.*
 Entering password: Press **"Enter"**, press **"Select ▲ or ▼"**, (for fast forward or backward keep key pressed).
 To confirm press **"Enter"** for 1 second, until a tone is emitted.

Handpiece activation	Name of handpiece on display	Choice	Factory default	Definition
Handpiece existing/HP 01	1:5	yes/no	yes	

Handpiece existing/HP 02	1:3	yes/no	yes	Deselect handpieces that do not belong to your product range by switching to "no". This will later shorten the scrolling list for "Handpieces". Otherwise you will later have to scroll down the complete list of available handpieces, each time you select a handpiece.
Handpiece existing/HP 03	1:2	yes/no	yes	
Handpiece existing/HP 05	Corr.Saw	yes/no	yes	
Handpiece existing/HP 06	Osc.Saw	yes/no	yes	
Handpiece existing/HP 07	Sag.Saw	yes/no	yes	
Handpiece existing/HP 08	Dermatome	yes/no	yes	
Handpiece existing/HP 09	Kirschner Hp	yes/no	yes	
Handpiece existing/HP 10	Tattoo Hp	yes/no	yes	
Handpiece existing/HP 11	Cont.Shaver	yes/no	yes	
Handpiece existing/HP 12	Osc.Shaver	yes/no	yes	
Handpiece existing/HP 13	Jacobs Chuck	yes/no	no	Because HP 04, 1:1 is permanently active, it's not listed in this table.
Handpiece existing/HP 14	Perforator	yes/no	no	
Handpiece existing/HP 15	Craniotome	yes/no	no	

Handpiece Max. Speed	Name of handpiece on the display	Speed range rpm	Factory default	Definition
Handpiece max speed/HP 01	1:5	1500 – 240 000	240 000	Limit the maximum speed of your handpieces according to your own experience
Handpiece max speed/HP 02	1:3	900 – 150 000	150 000	
Handpiece max speed/HP 03	1:2	600 – 100 000	100 000	
Handpiece max speed/HP 04	1:1	300 – 80 000	80 000	
Handpiece max speed/HP 09	Kirschner Hp	500 – 2800	2800	
Handpiece max speed/HP 10	Tattoo Hp	fix. 12 000	12 000	
Handpiece max speed/HP 11	Cont. Shaver	fix. 6000	6000	
Handpiece max speed/HP 12	Osc.Shaver	fix. 5000	5000	
Handpiece max speed/HP 13	Jacobs Chuck	fix. 2100	2100	
Handpiece max speed/HP 14	Perforator	1200	1200	
Handpiece max speed/HP 15	Craniotome	50 000	50 000	

*Scope of delivery contains "Activated handpiece HP04 1:1". Other handpieces are ordered and supplied separately!

Pump parameters	Range	Factory default	Definition
Pump/Backwards turn mode variable	No/Yes	Yes	The pressures in the tube set vary according to the pump speed. In "variable mode" this is considered, to prevent of spilling when pump is switched off in backwards mode.
Pump/Way backwards	1 - 100 %	25 %	Specify how far the pump turns backwards.
Pump/Speed backwards	10 - 50 %	33 %	Specify how fast the pump has to turn backwards to prevent of spilling after the switch off of the motor.
Pump/Range 1 increment	1 – 10 %	1 %	Adjustment steps in section 1.
Pump/Range 1 end	5 – 50 %	20 %	Set the range where section 1 is active.
Pump/Range 2 increment	1 – 10 %	5 %	Adjustment steps in section 2.
Pump/Range 2 end	10 – 90 %	50 %	Set the range where section 2 is active.
Pump/Range 3 increment	1 – 10 %	10 %	Adjustment steps in section 3.
Pump/Range 3 end	20 – 100 %	100 %	Set the range where section 3 is active

The MD HighSurg 30 Motor System for Medical Surgery can recognize the type of a plugged-in motor. This enables the adaption of the motors and their safe operation.

Table 12. Activated motors 21


Motor 21, 50,000 rpm	Range	Factory Default	Definition
Motor 21, 50,000 rpm/Min. Speed	300 – 5000 rpm	300 rpm	Set the min. speed of motor 1
Motor 21, 50,000 rpm/Max. Speed	5000 – 50 000 rpm	50 000 rpm	Set max. speed of motor 1
Motor 21, 50,000 rpm/Ramp start	1 – 1000 mc/10 000 rpm	250 ms	Set acceleration time to 10,000 rpm
Motor 21, 50,000 rpm/Ramp stop	1 – 1000 mc/10 000 rpm	50 ms	Set breaking time from 10,000 – 0 rpm
Motor 21, 80,000 rpm	Range	Factory default	Definition
Motor 21, 80,000 rpm/Min. Speed	300 – 5000 rpm	300 rpm	Set the min. speed of motor 2
Motor 21, 80,000 rpm/Max. Speed	5000 – 50 000 rpm	50 000 rpm	Set max. speed of motor 2
Motor 21, 80,000 rpm/Ramp start	1 – 1000 mc/10 000 rpm	250 ms	Set acceleration time to 10,000 rpm
Motor 21, 80,000 rpm/Ramp stop	1 – 1000 mc/10 000 rpm	50 ms	Set breaking time from 10,000 – 0 rpm
Max. speed warning delay	0 -1000 ms	1000 ms	0: Warning is switched off 1 – 1000: Delay time of warning when exceeding 50,000 rpm.

Resetting to factory default	Choice	Factory default	Definition
Default value/Set default value	Да/Нет	-	Resetting all parameters of the Configuration Menu to factory default.

- **Attention:** Upon resetting to factory default all parameters will appear with factory default values (except date and time and the operating hours counter).
- **Note:** For resetting the “programs” to factory default press both “**Select**” keys at the same time (“**Select** ▲ + ▼”). For this procedure you have to be out of the Configuration Menu.

Operation using the Vario pedal

1. **Carrying bracket**
The carrying bracket can be operated by foot (collapsible).
2. **Key  :**
Pressing the key briefly: switches the pump on or off (check display).
Pressing the key longer: increases the pump speed (check display).
3. **Key  :**
Pressing the key briefly: switches the program (+1) (check display).
Pressing the key longer: switches the program (-1) (check display).
4. **Key  :**
Pressing the key briefly: switches the rotational direction (check display).
Pressing the key longer: switches the motor (check display and indicator lights beside the motor sockets on the control unit).
5. **Treadplate**
With the treadplate pedal the motor speed is variably adjusted and the pump is activated.

Treadplate...	Motor:	Pump:
... not pressed	Motor off	Pump off
... pressed gently	Motor runs slowly	Pump on, if pump "On" displayed (speed as set on the control unit)
... pressed all the way down	Motor runs at maximum speed (speed as set on the control unit)	Pump on, if pump "On" displayed (speed as set on the control unit)

For safety reasons, the unit can only be operated by pedal.

The speed of the following handpieces is fixed and cannot be changed.

- Micro Compass Saw MSS 5000 (not included in delivery)
- Micro Oscillating Saw OMS 5000 (not included in delivery)
- Tattooing handpiece (not included in delivery)
- Micro Sagittal Saw MOS 5000 (not included in delivery)
- Dermatome (not included in delivery)

Functional check

Prior to the HighSurg 30 Motor System for Medical Surgery 30 startup or use of accessory equipment, the user must always ensure that each individual component is in good working order, free from defects, and is clean, sterile (if required) and operational. All inscriptions on the device and its accessories must be readable and there must be no loose parts in the device. Once the device is switched on, depending on configuration the most recent settings entered appear on the display and the green LED for motor 1 lights up.

1. Electronic motor 21

The performance of the electronic motor is checked without a handpiece or contra angle. Nevertheless, the 1:1 handpiece must be activated via the «Handpiece» key to check the maximum speed. Use the “Speed” selection keys to set the speed of the electronic motor to 50,000 rpm. Press the pedal treadplate; the electronic motor starts and accelerates to up to 50,000 rpm. When the treadplate is released, the electronic motor slows down again. For motor types with a permitted speed of 80,000 rpm, the maximum speed must be set accordingly to 80,000 rpm.

- The HighSurg 30 Motor System for Medical Surgery is designed for intermittent duty operation at maximum speed of “1 min ON/3 min OFF” at 4 cycles. After 15 min brake. Disregarding the directions of intermittent duty operation can cause burns when touching the handpiece or the motor.
- The motor ventilation slots must be kept clear in order to prevent the motor temperature from becoming excessive.

2. Pump

Press the key  on the pedal briefly; the peristaltic pump is switched on, which is shown on the display by the symbol of a drop. Press the pedal treadplate; the peristaltic pump and the electronic motor start up. Water sprays from the irrigation needle on the handpiece.

3. Rotational direction of the electronic motor

Press the key  on the pedal briefly; the rotational direction of the electronic motor changes. Press the pedal treadplate; the electronic motor rotates to the left and a continuous tone is emitted. Release the treadplate; the electronic motor ceases to operate and the tone is no longer heard. By pressing the motor

key again, the rotational direction is switched back to right rotation, which is shown on the display by the symbol of a changing arrow.

4. Program

The required program can be set by repeatedly pressing the key  on the pedal.

5. HOW TO USE HANDPIECES

HANDPIECES

The handpieces used with the The HighSurg 30 Motor System for Medical Surgery may have straight or angular design and also special handpieces such as continuous and oscillating shavers (not included in delivery), Kirschner handpieces, Perforator HP, Craniotome, Dermatorne, Tattooing HP and Jacobs Chuck (not included in delivery, ordered and supplied separately) can be used. They may be equipped with speed increasing/decreasing system and have a coolant supply system.

Non-working part of the handpiece or the contra angle is marked using an electrical discharge machine ensuring marking legibility. The marking stays lasting and legible during operation, storage and transportation. The marking of the handpieces and contra angles includes name of the manufacturer (NOUVAG AG), serial number (SN), model type (HP:1) and the Steam Sterilization symbol.

Pic. 8. Handpiece marking.

- The handpieces are not sterile on delivery. Before first use and immediately after every operation they must be cleaned, disinfected and sterilized!
- Do not clean the handpiece with compressed air!
- Speed limit for the handpieces are 50'000 and 80'000 revolutions per minute!
- Handpieces may only be attached when the electronic motor is not running!
- Never operate the clamping mechanism of the handpieces while the electronic motor is running! This could result in instrument damage.

HighSurg 30 Operation manual, Ref. 31666, V13/15

- The handpieces may only be operated by qualified and trained personnel!
- Improper use or repair of the device and failure to observe these instructions relieve NOUVAG AG from any obligations arising from warranty provisions or any other claims.

Explanation of symbols

CE mark with
Notified Body

Autoclavable at
max. 135 °C

Can be thermally
disinfected

Specification
of Serial
number

Warning

Manufacturer

Observe the
instructions for use


Specification
of order
number

Technical data

				
Ref. number	1904nou	1907nou	1950nou	1960nou
Transmission ratio	1:1	1:1	1:1	1:1
Conformity standard	ISO 3964	ISO 3964	ISO 3964	ISO 3964
Shank - Ø, mm	2,35 mm	2,35 mm	2,35 mm	2,35 mm
Max. speed, rpm	80.000	80.000	50.000	50.000
Max. torque, Ncm	6 Ncm	6 Ncm	6 Ncm	6 Ncm
File length, mm	31 mm	31 mm	21 mm	21 mm
Total length of the instrument, mm	95 mm	95 mm	70 mm	70 mm
Weight (with cooling clamp)	82g	110 g	85 g	113 g



Handpiece overview

Pic. 9. Handpiece overview

Installing and exchanging instrument head

1. Hold the handpiece from underneath with one hand and turn quick release locking mechanism counter clockwise with the other till it clicks open.
2. Remove placeholder pin or instrument.
3. Insert new instrument head or placeholder pin.
4. Close quick release locking mechanism, turning it in the opposite direction.
5. Check for perfect fit by slightly pulling at the instrument head.

Pic. 10. Instrument head inserted in the handpiece

When the instrument is not used, the placeholder must be inserted in the handpiece.

Use only appropriate instruments. Follow manufacturer's operation instructions.


Hightsurg 30 Operation manual, Ref. 31666, V13/15

Prior to use check the MD settings on the control unit!

Connection to the electronic motor

Attach the handpiece to the electronic motor. Press the handpiece firmly onto the electronic motor until it clicks into position and make sure it is secure by moving it slightly in the opposite direction. Reliable connection of the electronic motor and the handpiece means that you cannot remove the handpiece from the electronic motor by pulling the handpiece forward, while the handpiece can be rotated easily (360 degrees).

Pic. 11.


Assembly of the tube set

Before operation you must connect cooling pipes and attach them to the handpiece and cable of the electronic motor using special clips.

Assembly of irrigation system:

2.

1. Secure the cooling nozzle on the handpiece
2. Connect the tubing set (6024, 3m, sterile Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.2010) or 6025) with the cooling pipe of the handpiece.
3. Attach the handpiece clips (1881) to the cooling tubing set. Attach the clips to the handpiece.

Preparation instructions

Reprocessing restrictions

Frequent reprocessing has only a limited impact on the life period of the handpiece. The end of the product service life is normally determined by wear and damage through use. The handpiece is designed for 250 sterilization cycles.

INSTRUCTIONS

At the location of use	Remove surface soiling with a disposable cloth / tissue paper.
Storage and transportation	No special requirements. Due to the risk of blood drying and further corrosion reprocessing must be performed without undue delay (longer than 30 minutes).
Preparation for cleaning	Remove the instrument, the cooling pipe, the fixating clips and release the handpiece.
Automatic cleaning and disinfection	Equipment: Washer-disinfector with a special load carrier that ensures handpieces are connected to the washer-disinfector (for example, Miele G7835 CD Thermal Disinfector) and channels are rinsed. For small parts use a fine-mesh basket. Rinse them from every side. Use neutral or alkaline cleaning agents (for example, Schülke perform or 3% DESOMEDAN ID) in the manufacturer's recommended concentration. 1. Place the handpiece in the load carrier (ensure that the channels can be rinsed). 2. Place the small parts in the fine-mesh basket and close it. 3. Set a cleaning cycle that offers sufficient cleaning and rinsing. Perform the final rinse with demineralized water. 4. Perform a 10-minutes rinse cycle at 93°C to facilitate thermal disinfection. When removing the handpiece, inspect it visually to verify whether soiling is still visible in the gaps and the grooves. If necessary, repeat the cycle or clean manually.
Manual cleaning	Equipment: Neutral or alkaline cleaning agent (for example, Schülke perform or 3% DESOMEDAN ID) in the manufacturer's recommended concentration, a soft brush, running water. <i>Procedure:</i> 1. Rinse off surface soiling on the handpiece and small parts under running water. 2. Use a brush to apply the cleaning agent to all surfaces and gaps. 3. Rinse the handpiece thoroughly under running water.
Manual disinfection	For manual disinfection, wipe the handpiece and the small parts with a cloth soaked in an ethanol-based disinfecting agent.
Drying	If a drying program is not provided by the washer-disinfector, the handpiece and the small parts must be dried manually or in a drying cabinet until there is no visible moisture left on them.
Inspection and maintenance  No 1958	Perform a visual inspection to check for damage, corrosion and wear. After disinfection thoroughly and briefly spray the handpiece with NouClean Spray (optional) and wipe it with a moist cloth (see the instructions on the spray can). 
Packaging	Individual: Pack the handpiece in individual packaging for sterile items. The bag must be large enough to ensure that the seal is not subject to strain. Sets: Place the handpieces in the load carrier or on all-purpose sterilization trays.
Sterilization	Autoclave in the vacuum autoclave at 135°C for at least 5 minutes*. When sterilizing several instruments during one sterilization cycle, do not exceed the maximum sterilizer load. A drying cycle must be added in the case of autoclaves without a post-vacuum function. Allow the handpieces to dry in the bags with the paper side facing upwards. * Temperature exposure times are based on country-specific guidelines and standard.
Storage	No special requirements. If the sterilized handpiece is not used immediately after sterilization, the material packaging must be labeled with the sterilization date. Nouvag AG recommends to use a sterility indicator on the packaging.

In line with regulations pertaining to hazardous materials, 1984 Nou-Clean Maintenance and Care Spray is not delivered with the control unit and must be ordered separately. Delivered separately!

All the instructions were validated by the device manufacturer Nouvag AG as compliant to existing regulations. Those in charge of treatment are fully responsible for the results of treatment conducted in equipped hospital conditions as well as for the equipment and the materials used and for the staff involved in operating the device. The procedure must be validated and verified. Any deviations in the procedure stated in the instruction must be evaluated for their efficiency and possible adverse consequences.

Table 13. Malfunctions and troubleshooting

Malfunction	Cause	Solution
The electronic motor is running but the instrument doesn't switch on	The handpiece is not correctly connected with the electronic motor	Press the handpiece until it clicks into position and make sure it is secure.
Incorrect operation of the instrument	Incorrect attachment of the instrument to the contra angle	Fix the instrument and tighten the clamp
The handpiece causes excessive noise during operation	Dirty or badly oiled ball bearings	Spray the insides of the contra angle with Nou-Clean spray (not included into delivery, ordered and supplied separately)

9. SCOPE OF DELIVERY

“The The HighSurg 30 Motor System for Medical Surgery”

1.	3360	2-motor HighSurg 30 control unit with built-in peristaltic pump.	1 unit
2.	2098nou	Electronic motor 21 incl. 3 m motor cable, up to 80,000 rpm	1 unit
3.	2099nou	Electronic motor 21 incl. 3 m motor cable, up to 50,000 rpm	1 unit
4.	1510nou	Vario pedal; IPX8; electronic incl. 3.15 m cable	1 unit
5.	1170	Handpiece cradle	2 units
6.	1770	Stand for irrigation fluid	1 unit
7.	22262	Power cable 3 m	1 unit
8.	1904nou	Surgery high speed handpiece 1:1 FG shank $\varnothing$ 2,35 mm, length 95 mm	1 unit
9.	1907nou	Surgery high speed angled handpiece 1:1 FG shank $\varnothing$ 2,35 mm, length 95 mm	1 unit
10.	1950nou	Surgery high speed handpiece 1:1 FG shank $\varnothing$ 2,35 mm, length 70 mm	1 unit
11.	1960nou	Surgery high speed angled handpiece 1:1 FG shank $\varnothing$ 2,35 mm, length 95 mm	1 unit
12.	6024	Tube set, sterile, 3 m, single-use (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010) 10 pcs/unit	1 unit
13.	1873	Clip set (10 pcs/unit) for tube set attachment at electronic motor	1 unit
14.	1881	Clip set (3 pcs/unit) for tube set attachment at handpiece	1 unit

HighSurg 30 Operation manual, Ref. 31666, V13/15


15.	1958	Spray nozzle for care and maintenance of handpieces	1 unit
16.	1974	Spray nozzle for care and maintenance of electronic motor 21	1 unit
17.	28879	Fuse T 3.15 AL 250 W, AC	2 pcs
18.	31666	User Manual on CD	1 pcs

In line with regulations pertaining to hazardous materials, 1984 Nou-Clean Maintenance and Care Spray is not delivered with the control unit and must be ordered separately. Delivered separately!

10. CLEANING OF DEVICE COMPONENTS.

Schematic flow chart for resterialization as per EN ISO 17664

General information on device cleaning, disinfection and sterilization

- Reusable medical devices are subject to repeated sterilization and further storage under the conditions excluding secondary microbial contamination.
- All device components being in direct contact with patients' oral fluid, mucosa and oral tissues are subject to pre-sterilization treatment and sterilization. Other device components shall be disinfected without sterilization.
- The above medical devices shall be sterilized immediately before their first use as well as after every patient in order to avoid cross-contamination.
- The medical devices contacting patients' blood and serum shall be disassembled, cleaned and rinsed by medical staff wearing rubber gloves, protective goggles, a facemask, a gown and a waterproof apron covering the gown. All the staff shall carefully follow personal hygiene rules.
- All the containers containing working solutions of cleaning and disinfecting agents shall be equipped with lids and legibly labeled with the following data: description of the agent, its concentration, intended use, date of preparation and date of expiry. Containers with ready-to-use agents shall be labeled with the description of the agent and its intended use, if the agent may be used repeatedly, as well as with the date of the first use.
- For disinfection and pre-sterilization cleaning with chemical solutions immerse the medical devices in the working chemical solution so that all the gaps and grooves are filled. Before immersing disassemble all the parts that need assembling, open all the instruments with locking mechanisms and make some operational movements with them in the solution.
- Disinfection by wiping is allowed only for those medical devices that do not touch patients directly or for those devices which cannot be immersed into the disinfecting solution because of their design features.
- After disinfection, wash away any residues of the disinfecting agent from reusable medical devices following instructions for the agent use.
- Cleaning efficiency is confirmed by several tests: amidopyrine-based test is used to check for any blood residues left and phenolphthalein-based test is used to check for any residue of cleaning agent left. Quality inspection shall be performed daily. It is to cover 1% simultaneously treated instruments of every item, but at least three pieces. Both internal and external device surfaces are subject to inspection! The pre-sterilization treatment log is the document confirming such quality inspection.
- This medical device is subject for physical sterilization (autoclaving).
- Sterilization shall be performed in the clean area.
- To transport medical device into the clean area use special closed containers for sterilization.
- The hospital staff involved in pre-cleaning and disinfecting medical devices must wear special protection clothes as well as a facemask, protection goggles and rubber gloves as well as a water-proof apron covering the gown!
- After sterilization, store the medical devices under the conditions to exclude secondary microbial contamination in a special cabinet protected against dust, moisture and temperature fluctuations.
- Storage period is marked on the packaging and depends on the type of packaging

material as per the instruction for its use.

- Do not use sterilized medical devices if their storage period after sterilization is expired.
- Information on device sterilization is registered in a special log.

General recommendations on device preparation

- All the devices before being sterilized are subject to pre-sterilization cleaning in order to remove protein, fat and mechanical soiling as well as various medicines.
- For pre-sterilization cleaning disassemble all the parts that need assembling.
- The instruments are to be cleaned in a cleaning solution (for example, Schülke perform or 3% DESOMEDAN ID) while meeting all the personal safety requirements. While using solutions carefully follow all the manufacturer's instructions for use!
- To rinse the devices use demineralized water only (VE-water).
- Use only fresh daily prepared working solutions. If you have any gross soiling, replace the solution more often.
- For cleaning we recommend to use soft napless cloth, plastic brushes or cleaning pistols. Use of metal instruments is prohibited. For treatment use plastic or enameled (with intact enamel) containers covered with lids.
- For manual cleaning from the outside use a soft disposable napkin and disinfected brushes. Always clean the instruments with a brush under the surface of the cleaning solution in order to avoid dirty liquid splashing.
- The electronic motor must not be cleaned with ultrasound or compressed air.
- The medical devices shall be disinfected and cleaned simultaneously.
- For automatic cleaning and disinfection use a washer/disinfector (for example, Miele G7835 CD Thermal Disinfector) with thermal disinfection function.
- To avoid damages all the device parts must be safely secured in the washer/disinfector. The load carriers must ensure thorough rinsing of internal instrument cavities by proper nozzles.
- Instrument drying is an integral stage of automatic cleaning and disinfection process.
- After cleaning and disinfection is over instruments shall be removed by hands in antiseptic-treated disposable rubber gloves.
- All the cleaned and disinfected devices are subject to visual inspection. Critical areas need more close attention. All the devices and instruments with inner channels must be checked for any obstructions. Should any instruments fail to pass the inspection you must replace them without delay.
- Any instruments that are not clean enough are subject to the second cleaning cycle with thorough rinsing afterwards.
- Any bent instruments as well as instruments with micro cracks or other wear pattern are subject to immediate replacement as they are not fully able to perform their functions.
- If any instruments need some maintenance procedures, those maintenance procedures shall come before the performance inspection.
- Before the performance inspection please lubricate all the hinged and screwed instruments. Reliable operation is ensured by the inspection. Thus, all the disassembled

instruments must be assembled. If it is required, you may disassemble the instruments again after the inspection to sterilize them. While assembling/disassembling the devices follow applicable manufacturer's instructions.

- Any damaged or defective medical devices shall be sent for repairing or disposed.
- It is recommended to place the medical devices into the sterilization packaging. Follow the manufacturer's instructions regarding operation and the period while devices stay sterile!
- The packaging shall be appropriate and big enough for the below data to be labeled on it:
 - date of sterilization,
 - packing staff,
 - date of expiry,
 - description of the contents.
- For sterilization without using sterilization packaging the period for storing sterilized devices is 6 hours for sterilization in skeleton boxes and 20 days for sterilization in filter boxes, craft bags may be stored for three days maximum under sterile conditions. Any devices being sterilized without packaging must be used immediately after sterilization.
- Sterilize in vacuum autoclave (class B or S as per EN 13060) at 135°C. Sterilization specifications:

Sterilization mode					
Steam pressure inside the sterilization chamber, MPa (kgf/cm ²)		Operational temperature inside the sterilization chamber, °C		Duration of sterilization, minutes	
Rated value	Limit	Rated value	Limit	Rated value	Limit
0,20	±0,02	135	±2	5-25	±2

- To check for sterility the manufacturer recommends to use sterility indicator tapes (for example, Comply™ indicator tapes, manufactured by 3M).
- Returned instruments are those which are returned to the manufacturer after being used or without being used.
- Possible reasons for return include need for repairs or maintenance. Should any medical devices be returned to the manufacturer or its authorized representative, the latter will return the devices to the user only if they (the devices):
 - Have been subject to cleaning and disinfecting procedures as per manufacturer's instructions, are dry and declared to be hygienically clean, and/or
 - are expressly designed as being not soiled and have reliable packaging.
- The returned devices shall be decontaminated as in any case without delay in order to avoid device damage. Decontamination is prohibited if it is able to change or destroy the device, compromise analysis results or make the analysis impossible. If you have any doubts please consult the manufacturer.
- Every medical device being returned to the manufacturer shall be accompanied by the relative documents containing the below information:
 - Intended use of the medical device;
 - Decontamination method;

- Date of decontamination;
- Last name of the person responsible for decontamination;
- Reason for return.

All the instructions were validated by the device manufacturer as suitable for device reprocessing. Any reprocessing facilities are liable for reprocessing as well as for the equipment and the materials used and for the staff involved in reprocessing. The procedure must be validated and verified. Any deviations in the procedure stated in the instruction must be evaluated for their efficiency and possible adverse consequences.

Cleaning, disinfection and sterilization

The following points in particular are important with regard to caring for the material:

- Perform cleaning, disinfection and sterilization after every treatment.
- Always autoclave the material in sterilization packaging.
- Make sure that sterilization packaging is no more than 80% full.
- Always autoclave the material at 135°C for at least 5 minutes and 25 minutes max.
- If sterilized material is not used immediately, the material packaging must be labeled with the sterilization date.
- Nouvag AG recommends including a sterility indicator.

Control unit and pedal

The control unit and the pedal do not contact the patient.
Wipe the outside using micro-biologically tested surface disinfectants or a 70% ethanol solution. The front plate of the control unit is sealed for this purpose and can be wiped clean.

Electronic motor 21

Reprocessing restrictions	Frequent reprocessing has only a limited impact on the electronic motor. The end of the product service life is normally determined by wear and damage through use. Electronic motor 21 is designed for 250 sterilization cycles.
---------------------------	---

INSTRUCTIONS	
At the location of use	Remove surface soiling with a disposable cloth/paper towel.
Storage and transport	No special requirements. Due to the risk of corrosion when drying, reprocessing must be performed without undue delay.
Preparation for cleaning	Remove soiling from the electronic motor with a disposable cloth/paper towel. Unscrew the motor cap, remove the cable and unscrew the handpiece carrier. 

Automatic cleaning and disinfection	Equipment: Washer-disinfector with a special load carrier that ensures motors are connected to the washer-disinfector (for example, Miele G7835 CD Thermal Disinfector) and channels are rinsed. Rinse the motor from the front. Use neutral or alkaline cleaning agents (for example, Schülke perform or 3% DESOMEDAN ID) in the manufacturer's recommended concentration. Place the electronic motor in the load carrier (ensure that the channels can be rinsed). Place the motor cap and cable, and handpiece holder, in the basket. Set a cleaning cycle that offers sufficient cleaning and rinsing. Perform the final rinse with fully demineralized water. Perform a 10-minute rinse cycle at 93°C to facilitate thermal disinfection. When removing the electronic motor, check the motor cap and cable, and handpiece holder, to verify whether soiling is still visible in the gaps and the grooves. If necessary, repeat the cycle or clean manually.
Manual cleaning	Equipment: neutral or alkaline cleaning agent (for example, Schülke perform or 3% DESOMEDAN ID), soft brush, demineralized water (< 38°C). Procedure: Rinse off and brush away surface soiling on the electronic motor, motor cap and cable, and handpiece holder. Use a brush to apply cleaning agent to all surfaces and gaps. Rinse the electronic motor, motor cap and cable, and handpiece holder thoroughly under running water.
Manual disinfection	For manual disinfection, wipe the electronic motor, the motor cap and the plug and cable underneath it, and the handpiece holder, with a certified non-chlorine disinfectant.
Drying	If a drying program is not provided by the washer-disinfector, the electronic motor must be dried manually or in a drying cabinet until no visible moisture is left. The handpiece holder must then be screwed back onto the motor.
Inspection and maintenance  Art. No. 1974	Perform a visual inspection to check for damage, corrosion and wear. Connect the electronic motor with the Nou Clean spray nozzle (Art. No. 1974) by screwing it on and spray electronic motor for about 3 seconds and briefly wipe with a moist cloth (see instructions on the spray can). Once thorough spraying has been completed, screw the cable and motor cap back onto the electronic motor. 
Packaging 	Individual: Pack the electronic motor in individual packaging for sterile items. The bag must be large enough to ensure that the seal is not subject to strain. Nouvag AG recommends including a sterility indicator. Sets: Sort electronic motors on trays intended for this purpose or place them on allpurpose sterilization trays.
Sterilization	Autoclave in the vacuum autoclave (Class B or S as per EN 13060) at 135°C for at least 5 minutes*. When sterilizing several instruments during one sterilization cycle, do not exceed the maximum sterilizer load. A drying cycle must be added in the case of autoclaves without a post-vacuum function. Allow the electronic motor to dry in the bag for at least 1 hour at room temperature with the paper side facing

HighSurg 30 Operation manual, Ref. 31666, V13/15

	upwards. * Temperature exposure times are based on country-specific guidelines and standards. Maximum exposure time: 25 minutes.
Storage	If the sterilized electronic motor is not used immediately after sterilization, the material packaging must be labeled with the sterilization date. Including a sterility indicator is recommended.

Tubing set, Art. No. 6024, 3 m long sterile (Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.20102) 10pcs/unit and tubing set Art. No. 6025 (optional)

- Tubing set, 3 m long, sterile (Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.20102) Art. No. 6024 10pcs/unit and tubing set Art. No. 6025 (optional) are single-use only and may not be reused.
- Tubing sets must be disposed of properly after use!
- Do not use tubing sets when packaging is already open or damaged!
- Do not use tubing sets when expiry date has passed.
- Use only Novvag sterile tubing sets with Art. No. 6024 3 m long, sterile (Marketing Authorization № FSZ (ΦC3) 2010/07320 dated 05.07.20102) 10pcs/unit and Art. No. 6025 (optional)!
- The period for the tubing sets to stay sterile, apyrogenic and non-toxic is 5 years.

Sterility cannot be guaranteed by reusing and re-sterilization of tubing sets. The characteristics of the material change in a manner that can result in failure of the system. This may result in serious infections or even patient death worst-case.

Handpieces

Sterilization procedure for handpieces is specified in Chapter "HOW TO USE HANDPIECES".

Handpiece cradle


Soiled handpiece cradles are cleaned using a neutral cleaning agent and then sterilized in accordance with the same instructions as for the electronic motor.

11. Maintenance

1. Replacing the control unit fuse

Users can replace faulty control unit fuses themselves. These are located at the rear of the device in the fuse slot beside the power switch (see picture 5):

- Unplug the power plug.
- Open the fuse slot using a screwdriver.
- Replace the faulty fuse T 3.15 AL 250 V AC.
- Slide the fuse holder back in and close the fuse slot.
- Check the mains voltage shown on the fuse slot.
- Plug in the power plug again.

- 1. Fuse slot locking mechanism
- 2. Display window for voltage setting
- 3. Fuse slot
- 4. Fuse 1
- 5. Fuse 2

Table 14. Malfunctions and troubleshooting

Malfunction	Cause	Solution	Refer to operating instructions
Device is not operational	Control unit is not switched on	Set the power switch "I/O" to "I"	Chapter "Switching the device on and off"
	Power connection not established	Connect the control unit to the mains power supply	Chapter "Connection to the power supply"
	Incorrect operating voltage	Check the mains voltage	Chapter "Connection to the power supply"
	Faulty fuse	Replace the fuse	Chapter "Replacing the control unit fuse"
Motor does not run	Motor not switched on	Switch on the motor using the Vario pedal	Chapter "Operation using the Vario pedal"
	Inappropriate active motor	Activate appropriate motor, using Vario pedal	Chapter "Operation using the Vario pedal"
	Motor not connected	Connect the motor cable to the control unit	Chapter "Device overview" Chapter "Connection to the power supply"
	Handpiece is not correctly assembled	Press the handpiece firmly onto the electronic motor until it clicks into position and check that it is secure by moving it slightly in the opposite direction.	Chapter "How to use handpieces"
No irrigation fluid for	Peristaltic pump (cooling	Switch on the peristaltic	Chapter "Operation using the

instrument	pump) not switched on	pump.	Vario pedal"
	Tubing set incorrectly inserted	Insert tubing set correctly (note the direction)	Chapter "Device preparation"
	Tubing set clogged/ crusted matter visible	Replace the tubing set	Chapter "Device preparation"
	Bottle with sodium 0.9% chloride solution not ventilated	Open the ventilation filter at the drip chamber	Chapter "Device preparation"
	Tubing set is dripping	Replace the tubing set	Chapter "Device preparation"
	Roller clamp of tubing set is closed	Open roller clamp all the way	Chapter "Device preparation"
Pedal is not operational	Pedal not connected	Connect the pedal to the control unit	Chapter "Device preparation"
	Incorrect operation	Check Operating Instructions	Chapter "Operation using the Vario pedal"

If a fault cannot be rectified, please contact your supplier or an authorized service centre. Contact details are given on the last page of these Instructions for Use

Table 15. Error messages on display screen

HighSurg 30, Error messages on display screen		
Сообщения об ошибках/ Код ошибки	Причина	Устранение
Базовая инициализация/ W00	Первая инициализация	
Reset to factory default values/ W01	Оброс на параметры по умолчанию	
Memory failure/ E02	System failure	Send the control unit to the service centre
Operation failure/ E03	System failure	Send the control unit to the service centre
Program failure SW/ E04	System failure	Send the control unit to the service centre
Failure of user configuration SW/ E05	System failure	Send the control unit to the service centre
Display failure/ E06	System failure	Send the control unit to the service centre
Pump failure/ E07	System failure	Send the control unit to the service centre
Storing factory settings/ User Config & Program	Message while default values of parameters and programs are saved to the NOU-dongle.	
Storing factory settings/ Program	Message while default values of programs are stored.	
Pedal not connected/ E10	a) Pedal is not plugged in. b) Plug or cable is defective.	a) Plug in Pedal b) Send Control Unit and pedal to Service Center.

HighSurg 30 Operation manual, Ref. 31666, V13/15

Pedal test mode/ W11	Pedal test mode switched on	Switch off device for 5 seconds, than switch on again.
Keyboard test mode/ W12	Keyboard test mode switched on	Switch off device for 5 seconds, than switch on again.
No motor connected/ E13	a) No motor connected b) Motor, motor cable, motor plug or Control Unit is defective	a) Plug in motor b) Send motor and Control Unit to Service Center.
Motor 2 not connected/ E14	a) Motor 2 is selected but no motor plugged in b) Motor connected to motor socket 2, but motor, motor cable, motor plug or Control Unit is defective.	a) Plug in motor b) Send motor and Control Unit to Service Center.
Motor 1 not connected/ E15	a) Motor 1 is selected but no motor plugged in b) Motor connected to motor socket 1, but motor, motor cable, motor plug or Control Unit is defective.	a) Plug in motor b) Send motor and Control Unit to Service Center.
Unknown motor 2/ E16	a) Motor 2 is selected but wrong motor is plugged in. b) Right motor is connected to motor socket 2, but motor, motor cable, motor plug or Control Unit is defective.	a) Plug in correct motor b) Send motor and Control Unit to Service Center.
Unknown motor 1/ E17	a) Motor 2 is selected but wrong motor is plugged in. b) Right motor is connected to motor socket 1, but motor, motor cable, motor plug or Control Unit is defective.	a) Plug in correct motor b) Send motor and Control Unit to Service Center.
Hp not allowed for motor 2/ E18	The Handpiece may not be operated with the motor plugged in at the motor socket 2.	Choose appropriate Handpiece or change to correct motor type.
Hp not allowed for motor 1/ E19	The Handpiece may not be operated with the motor plugged in at the motor socket 1.	Choose appropriate Handpiece or change to correct motor type.
Pump is open/ E20	Motor is not working when pump compartment is open to prevent of injuries.	Close pump compartment.
Motor or pump test mode/ W21	Motor or pump test mode is switched on.	Switch off device for 5 seconds, than switch on again.
Pedal locked/ W26, pedal let go	If pedal is pressed at switch on procedure, pedal will not work.	Release pedal for one second.
Battery is almost empty/ W27, continue with „Enter“	Battery is almost empty.	- After pressing „Enter“, work can be continued. - Send Control Unit as soon as possible to Service Center.

HighSurg 30 Operation manual, Ref. 31666, V13/15

Watch Error XX/ E28, continue with „Enter“	a) Clock at Control Unit is defective. b) Device was switched on after battery change but clock is not set yet.	a) After pressing „Enter“, work can be continued but Control Unit has to be sent to Service Center as soon as possible. b) Set date and time.
Handpiece XX is faulty/ E29	At calibrating or testing handpiece was overexposed to high torque.	- Clean handpiece and spray it thoroughly with Nou-Clean spray. - If message is still displayed after test procedure, handpiece has to be sent to Service Center.
Handpiece XX is Ok!	Tested handpiece is OK.	
Testing the handpiece XX	Handpiece is testing.	
Check max. Speed handpiece/ V30, wait 1 second	Message is displayed when motor with 80,000 rpm is selected and max. speed is set to 50,000 rpm. Check Handpiece for max speed clearance before you proceed.	Check handpiece and wait for the message to disappear. Motor starts if the pedal still is pressed.
NOU-Dongle is plugged in	Message is displayed for 1 second after NOU-dongle was plugged in.	

The red highlighted error messages are also illuminated red on the control device display. The other messages are for information and do not require any action by the user.

2. Safety inspections

The safety inspection is a regular safety check that is compulsory for those operating medical devices. The objective of this measure is to ensure that device defects and risks to patients, users or third parties are identified in time.

An inspection interval of **2 years** applies to the “HighSurg 30 Motor System for Medical Surgery” It is conducted by trained personnel. Inspection results shall be registered in a special maintenance log.

If you have any questions regarding device maintenance and repair as well as after-sales servicing please contact our authorized representative in the Russian Federation: **ООО MedTechKonsult, 1, Nemanskiy Proezd Str., bldg “1”, office 129, Moscow, 123181, Russia, phone/fax +7 (495) 757-44-48, +7 (495) 660-92-80, E-mail: medkonsult@mail.ru.**

12. SAFETY REQUIREMENTS

Your own safety, safety of your team and of course of your patients is very important for us. Thus, it is essential to pay attention to the following information.

12.1. Manufacturer’s EC Declaration

Please take into notice the below information on electromagnetic compatibility, provided with these Instructions for Users (see Appendix 3).

HighSurg 30 Operation manual, Ref. 31666, V13/15

12.2 Integrated peristaltic pump

The integrated peristaltic pump is used to cool tissue in order to prevent damage to the tissue (necrotization). It may only be operated with watery solutions, such as 0.9% Sodium Chloride irrigation solution (not included in delivery) or “Ringer” (not included in delivery) solution. **Medication supply using the integrated pump is expressly prohibited.**

12.3 Modification and misuse

If modifications are made by the user/operator or a third party to the “HighSurg 30 Motor System for Medical Surgery” provided, or if these items are used by these individuals for a purpose other than that for which the items were designed (according to chapter “Intended Use”), the manufacturer assumes no responsibility and the guarantee is void.

12.4 Essential requirements

The “HighSurg 30 Motor System for Medical Surgery” may only be operated by qualified and trained personnel.

Improper use or repair of the device and failure to observe these instructions relieve us from any obligation arising from warranty provisions or other claims.

The use of third-party products is the responsibility of the operator. Functionality and patient safety cannot be guaranteed with third-party accessories.

Prior to using the device, before startup, and before operation, the user must always ensure that the device and accessories are in good working order and are clean, sterile (if it is necessary) and operational.

Repairs may only be performed by authorized NOUVAG service technicians.

Use Nou-Clean spray (is not included into the scope of delivery; optional) for maintenance and care of the electronic motor, handpieces and contra angles. Using other care products can result in malfunctions and/or cause the warranty to be revoked.

12.5 During use

The device is not sterile on delivery. All sterilizable parts must be sterilized before use (see Chapter “Cleaning, disinfection and sterilization”).

To prevent injury, never touch drill bits or burrs while they are still running.

Handpieces and contra angles may only be attached when the electronic motor is not running.

Never operate the clamping mechanism of the handpieces or contra angles while the system is running. This could result in instrument damage.

Make sure that the supply voltage setting corresponds to the local mains voltage.

Do not use device in the vicinity of flammable mixtures!

The manufacturer runs the quality control and risk analysis system. The risk management procedures performed meet the requirements of EN ISO 14971:2012. Risk control is effected at device production and post-production stages. Risk control is also performed during post marketing observance of device use results. Risk analysis is an integral part of the technical documents of the device.

12.6. Safety inspections

The safety inspection is a regular safety check that is compulsory for those operating medical devices. The objective of this measure is to ensure that device defects and risks to patients, users or third parties are identified in time.

An inspection interval of 2 years applies to the “HighSurg 30 Motor System for Medical Surgery”. NOUVAG AG offers a safety inspection service for its customers. Inspection results shall be registered in a special maintenance log.

If you have any questions regarding device maintenance and repair as well as after-sales servicing please contact our authorized representative in the Russian Federation: OOO MedTechKonsult, 1, Nemanskiy Proezd Str., bldg “1”, office 129, Moscow, 123181, Russia, phone/fax +7 (495) 757-44-48, +7 (495) 660-92-80, E-mail: medkonsult@mail.ru.

12.7. Environment protection requirements during MD operation

The medical device is entirely safe for the environment during operation.

13. MEDICAL DEVICE CLASSIFICATION

Electrical protection class:	I
Method of protection against electric shock:	B Type
Intended use:	Power instruments
Frequency of application:	Multiple application
Application	Professional medical device
Mode of operation:	Intermittent mode of operation
Type of perceived mechanical impact	Stationary
Type of human body contact	Short-term contact (less than 24 hours) with oral mucosa and mouth bone tissue

The “HighSurg 30 Motor System for Medical Surgery” complies with Council Directive 93/42/EC on Medical Devices and with the below standards:

EN60601.1:2006	General electrical safety standards
EN60601-1-2/07.2007	Electromagnetic compatibility. Requirements and testing. Medical equipment. Biological evaluation of medical devices. Part 1. Evaluation and studies..
EN 62366-2008	Medical devices. Application of usability engineering to medical devices.

HighSurg 30 Operation manual, Ref. 31666, V13/15


EN61000-3-2	"Electromagnetic compatibility (EMC) - Part 3, Section 2: Limits - Limits for harmonic current emissions (equipment input current ≤ 16 A per phase)"
EN61000-3-3	Electromagnetic compatibility (EMC) - Part 3 Section 3: Limits - Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems, for equipment with rated current ≤ 16 A per phase and not subject to conditional connection
EN61000-4-2	(EMC). Part 4: Testing and measuring techniques.– Section 2: Electrostatic discharge immunity test. Basic EMC Publication .
EN61000-4-3	(EMC). Part 4: Testing and measuring techniques.– Section 3: Radiated, radio-frequency, electromagnetic field immunity test. Basic EMC Publication.
EN61000-4-4	(EMC). Part 4: Testing and measuring techniques.– Section 4: Electrical fast transient/burst immunity test. Basic EMC Publication.
EN61000-4-5	(EMC). Electromagnetic compatibility (EMC). Part 4: Testing and measuring techniques. Section 5. Surge immunity test.
EN61000-4-6	Electromagnetic compatibility (EMC). Part 4: Testing and measuring techniques. Section 6. Immunity to conducted disturbances, induced by radio-frequency fields.
EN61000-4-8	(Electromagnetic compatibility (EMC). Part 4: Testing and measuring techniques. Section 8. Power frequency magnetic field immunity test.
EN61000-4-11	Electromagnetic compatibility (EMC). Part 4: Testing and measuring techniques. Section 11. Voltage dips, short interruption and voltage variations immunity test.
IEC 60601-1-1-2011	Medical electrical equipment - Part 1-1. General requirements for basic safety. Safety requirements for medical electrical systems.
ISO 3964	Dental handpieces. Coupling dimensions.
EN 10993-1:2009	Biological evaluation of medical devices. Part 1: Evaluation and testing within a risk management.
EN ISO 17664-2004	Sterilization of medical devices – Information to be provided by the manufacturer for the processing of resterilizable medical devices
ISO 13485:2012	Design and development, production, distribution and servicing of medical devices.
CE certificate	Compliance with Council Directive 93/42/EEC on Medical Devices

14. DEVICE PACKAGING AND MARKING

1. MARKING


The medical device is marked as follows:

Table 16. Medical device marking.

	Important information		Autoclavable at max. 135° for at least 5 minutes
	Do not use if the packaging is damaged		Suitable for thermal disinfection
	Warning		Sterilized using ethylene oxide
	Manufacturer		Observe the operation manual
	The device is designed for intermittent duty operation at "1 min ON/3 min OFF"		Electrical and electronic devices that have reached the end of their service life comprise hazardous

HighSurg 30 Operation manual, Ref. 31666, V13/15

	at 4 cycles, after 15 min. brake		waste and may not be disposed of together with household waste. Valid local disposal regulations apply.
	Type B method of protection against electric shock	SN	Symbol indicating the serial number with the date of manufacturing (year/month).
	Do not reuse	REF	Symbol indicating the order number.
	Biohazard	LOT	Symbol indicating the lot number.
	Motor 1		Motor 2
	Pedal		contains phthalates
IPX8	Protection against contact and continual submerging		Certified by the Canadian Standards Association (CSA)
	Date of manufacturing		Date of expiry
	CE mark with reference number of the notified body		Refer to the instructions for use

Pic. 12. The plate on the rear panel of the device.

NOUVAG AG
St. Gallerstrasse 23-25
CH-9403 Goldach

Name and address of the manufacturer

Type : HighSurg 30

Description of the medical device

REF

3360

Reference number of the medical device

SN

7228U1511R

Serial number of the medical device


Power : 100V~/115V~/230V~ 50-60Hz 120VA

consumption

Symbol indicating power

Fuse : T 3.15 AL 250 VAC

Fuses used

Symbol indicating that the CE certificate has been obtained and that the medical device is certified by the Canadian Standards Association (CSA).

Symbol referring to the Operation Manual

Symbol indicating hazardous waste which may not be disposed of together with household waste


Symbol indicating the country of origin (Switzerland)

HighSurg 30 Operation manual, Ref. 31666, V13/15

2D datamatrix (QR-code)

The device is designed for intermittent duty operation at "1 min ON/3 min OFF" at 4 cycles, after 15 min. brake

Pic. 13. Symbol indicating Type B method of protection against electric shock.

2. PACKAGING

Individual packaging

Pic. 14. Individual packaging of the handpiece

HighSurg 30 Operation manual, Ref. 31666, V13/15

Contra angle packaging made of white carton; its dimensions are 205x075x035 mm
 Contra angle packaging made of blue plastic; its dimensions are 200x070x030 mm

Рис. 15. Marking on the individual packaging

Information on the manufacturer

Typ : High-Speed Handstück gerade 95mm

Type : High-Speed Handpiece straight 95mm

Description of the medical device

1904nou

Reference number of the medical device

3258H1508R

Serial number of the medical device

CE mark

Refer to the Operation Manual

Autoclavable at max. 135° for at least 5 minutes.

Suitable for thermal disinfection

2D datamatrix (QR-code)

Symbol indicating the country of origin (Switzerland)

Sterile products

6024 Tube set, sterile, 3 m long (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010)
10pcs/unit.



Рис. 16. Individual packaging of the sterile product.


The individual commercial sterile packaging (see picture 16) is made of PP50/PE30 PEEL plastic film, packaging dimensions are 210 mm X 100 mm. The packaging guarantees the products stay sterile, nontoxic and apyrogenic. Design of the commercial packaging ensures product integrity, prevents microbial penetration and protects the product against mechanical damages. The individual commercial packaging does not affect toxicological properties of the packed items. The individual packaging ensures the products stay sterile for at least 5 years.

Pic. 17. Marking on the individual packaging of the sterile product

The "HighSurg 30 Motor System for Medical Surgery" and the "Dental surgery portable motor device MD 20" version (Marketing Authorization №FSZ (ФСЗ) 2010/07320 dated 05.07.2010) assembly instructions


- 1 – description of the medical device
- 2 – catalogue number
- 3 – do not use if the packaging is damaged
- 4 – do not reuse
- 5 – do not use if the packaging is damaged
- 6 – address of the manufacturer
- 7 – CE mark
- 8 – contact details of the manufacturer
- 9 – sterilized using ethylene oxide
- 10 – latex-free
- 11 – necessary to get acknowledged with the instructions
- 12 – date of sterilization

Pic. 18. Shipping packaging of the sterile product

The shipping packaging of the sterile product is shown in picture 18. It is a white carton, dimensions 310x215x80 mm with 10 individual packagings packed into it.

Pic. 19. Marking on the shipping packaging of the sterile product.

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DEHPfree

phthalates free

warning against product use if the packaging is damaged

2015-10

date of manufacturing

2021-06-28

date of expiry

CE mark

QTY: 10

10 items per a packaging

2D datamatrix (QR-code)

Manufacturer:

NOUVAG AG

St. Gallerstr. 23-25

CH-9403 Goldach

Tel. : 0041 71 846 66 00

Fax. : 0041 71 845 35 36

E-mail : info@nouvag.com

Homepage: www.nouvag.com

contact details of the manufacturer


HighSurg 30 Operation manual, Ref. 31666, V13/15

NOUVAG GmbH
 Schulthaißstrasse 15
 D-78462 Konstanz
 Tel.: +49 7531 1290-0
 Fax: +49 7531 1290-12
 E-Mail: info-de@nouvag.com
 www.nouvag.com

contact details of the distributor within the EU

MD Shipping packaging

Inside the shipping package the medical device is packed as follows:

Pic. 20. Packing of the medical device.

HighSurg 30 Operation manual, Ref. 31666, V13/15

The device is packed into a 400x400 mm polyethylene flat antistatic bag and after that in a 308x300x385 mm foamed polystyrene packaging (2 units). The shipping box contains a CD with the Operation Manual on it, a warranty card and that double packaging.

The shipping packaging (picture 21) protects the medical device against any mechanical shocks caused by transportation. Being packed in the shipping packaging the medical device is vibration- and shock-resistant. The packaging also protects the device during handling operations. The packaging is a 395x395x405 mm carton marked and labeled as follows (pictures 21, 22, 23, 24, 25).


Pic. 21. Shipping packaging

Рис. 22. Shipping packaging labeling

Рис. 23. Marking on the shipping packaging (example)

the medical device

description of


HighSurg 30 Operation manual, Ref. 31666, V13/15


Reference number

Serial number


name and address of the manufacturer

CE mark

2D datamatrix (QR-code)

Pic. 24. Information labeled on the shipping packaging.

The information labeled on the shipping packaging includes:


1. Company logo with address and phone numbers of the manufacturer.
2. The following note: "The operation manual is not available in printed form with this product. The operation manual as PDF is on the CD-ROM that was included in the delivery. A CD-ROM drive and the Adobe Acrobat Reader software is needed to read or print this operation manual. To order the operation manual in printed form please check on the type plate for the address or visit our website (www.nouvag.com). If you send us the article number (31686 EN) with the language code (RU for Russian) and the note "Printed copy", a printed manual will be sent to you free of charge.


The shipping packaging is marked with the following handling symbols.

Pic 25. Handling symbols on the shipping packaging.

Table 17. Handling symbols on the packaging.

	Keep dry The device must not be located near moisture sources or in a humid facility
	Handle with care. Fragile goods. Goods shall be handled with care. Such a symbol is used to mark goods which may be damaged in the result of tossing or shock.
	Transportation conditions The symbol is used to mark the transportation requirements (humidity and atmospheric pressure).

Transportation conditions

The symbol is used to mark the transportation requirements (temperature).

15. TRANSPORTATION AND STORAGE CONDITIONS.

There are no special storage requirements. The device may be stored at the below conditions. The device packed in the manufacturer's packaging should be stored in the warehouse and put on the racks. Avoid putting any other equipment or goods on the device.

There are no special transportation requirements as well. The device doesn't require disassembling for transportation. The device may be shipped by any covered vehicles or in multi-purpose containers as per the relevant shipping regulations at the below ambient conditions. Arrangement and attachment of shipping packagings must ensure their stable position preventing packagings from moving and bumping into each other.

Transportation and storage conditions:

Relative humidity:	Maximum 90 %
Temperature:	0 °C up to 60 °C
Atmospheric pressure:	700 hPa up to 1060 hPa

16. SERVICE LIFE

Shelf life of the medical device is 5 years. Average service life is **5 years from the date of manufacturing provided the device is operated according to the Operation Manual**. The device is regularly submitted to safety checks, after which the service life of 5 years may be extended. After exploitation resource is over the medical device must be disposed of as per national regulations of the importing country as well as recommendations of the manufacturer.

For sterile disposable devices:

6024 Tube set, sterile, 3 m long, single use, 10 items/unit.

The service life after the packaging is opened is 3 hours. The guaranteed storage period is 5 years.

17. DISPOSAL

The medical device shall be disposed of and destroyed as per the classification and rules of assembling, operation, decontamination, handling, storage, transportation, recording and disposal enforced by the

competent authority of the importing country (Class B of medical waste hazard as per the Sanitary Regulations SanPin 2.1.7.2790-10 according to Russian Classification).

Do not dispose of the device with household waste!

To ensure environmental protection, old devices can be returned to the dealer or manufacturer.

Electronic motors that have reached the end of their service life comprise hazardous waste and may not be disposed of with household waste. Prevailing national and local disposal regulations apply.

Contaminated single-use tubing sets are subject to specific disposal requirements. Please observe currently valid national disposal regulations for infectious waste (Class B of medical waste hazard as per the Sanitary Regulations SanPin 2.1.7.2790-10 according to Russian Classification).

The user is liable for proper disposal of this equipment at specially designated collection points for electrical and electronic waste recycling. Separate collection and recycling of waste electrical and electronic equipment during disposal contributes to preservation of natural resources. Besides, use of breakthrough technologies of waste recycling encourages protection of staff health and environment. If you have any questions regarding disposal of medical equipment with expired service life please refer to its supplier.

18. WARRANTY COVERAGE

The manufacturer Nouvag AG provides warranty coverage of the medical device. Manufacturer's address is St.Gallerstrasse, 23-25, CH 9403 Goldach, Switzerland. Purchasing the "HighSurg 30 Motor System for Medical Surgery" entitles you to a 1-year warranty. If you return the warranty card for registration within four weeks of the date of purchase, warranty coverage will be extended for a further 6 months. Warranty coverage doesn't extend to the consumables. Improper use or repair, or the device breaking because of failure to observe these instructions, will relieve the manufacturer from any obligation arising from warranty provisions or other claims.

Modification and misuse

- If modifications are made by the user/operator or a third party to the "HighSurg 30 Motor System for Medical Surgery" provided, or if these items are used by these individuals for a purpose other than that for which the items were designed (according to chapter "INTENDED USE"), the manufacturer assumes no responsibility and the guarantee is void.

Any claims or questions regarding handling of this medical device please send to:

the Russian Federation: OOO MedTechKonsult, 1, Nemanskiy Proezd Str., bldg "1", office 129, Moscow, 123181, Russia, phone/fax +7 (495) 757-44-48, +7 (495) 660-92-80, E-mail: medkonsult@mail.ru.

HighSurg 30 Operation manual, Ref. 31666, V13/15

Appendix 1. Medical device components

I. The “HighSurg 30 Motor System for Medical Surgery”:


1.	3360	2-motor HighSurg 30 control unit with built-in peristaltic pump.	1 unit
2.	2098nou	Electronic motor 21 incl. 3 m motor cable, up to 80,000 rpm	1 unit
3.	2099nou	Electronic motor 21 incl. 3 m motor cable, up to 50,000 rpm	1 unit
4.	1510nou	Vario pedal; IPX8; electronic incl. 3.15 m cable	1 unit
5.	1170	Handpiece cradle	2 units
6.	1770	Stand for irrigation fluid	1 unit
7.	22262	Power cable 3 m	1 unit
8.	1904nou	Surgery high speed handpiece 1:1 FG shank $\varnothing$ 2,35 mm, length 95 mm	1 unit
9.	1907nou	Surgery high speed angled handpiece 1:1 FG shank $\varnothing$ 2,35 mm, length 95 mm	1 unit
10.	1950nou	Surgery high speed handpiece 1:1 FG shank $\varnothing$ 2,35 mm, length 70 mm	1 unit
11.	1960nou	Surgery high speed angled handpiece 1:1 FG shank $\varnothing$ 2,35 mm, length 95 mm	1 unit
12.	6024	Tube set, sterile, 3 m, single-use (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010) 10 pcs/unit	1 unit
13.	1873	Clip set (10 pcs/unit) for tube set attachment at electronic motor	1 unit
14.	1881	Clip set (3 pcs/unit) for tube set attachment at handpiece	1 unit
15.	1958	Spray nozzle for care and maintenance of handpieces	1 unit
16.	1974	Spray nozzle for care and maintenance of electronic motor 21	1 unit
17.	28879	Fuse T 3.15 AL 250 W, AC	2 pcs
18.	31666	User Manual on CD	1 pcs

Appendix 2. Tables of device components.

Table 1. Device components and accessories.

No	Description of the component	Appearance	Weight (grams), not exceeding	Dimensions (mm), not exceeding
1	3360 HighSurg 30 control unit		3700	250x260x110

2	2098nou Electronic motor 21 "HighSpeed", speed 80,000 revolutions/minute, including 3 m long motor cable		300	Electronic motor: length – 80, Diameter – 20, Length of motor coupling – 30, diameter – 8 Cable length 3000
3	2099nou Electronic motor 21, speed 50,000 revolutions/minute, including 3 m long motor cable		300	Electronic motor: length – 80, Diameter – 20, Length of motor coupling – 30, diameter – 8 Cable length 3000
4	1510 nou Vario IPX8 foot pedal for MD 30 with 3.15 m long cable		1680	255 x 235 x 60 (Width x Depth x Height)

5	1770 Stand for irrigation fluid bottle		140	Height 460 mm Length under the hook 50 mm, hook 40 mm
6	1170 Handpiece cradle		50	Length 100 mm Width 60 mm Height 35 mm
7	22262 Power cable, 3 m long		320	Length 3000 Diameter 75

8	31666 User Manual on CD		10	Diameter 120
9	28879 Fuse T 3.15 AL 250 V AC		7	Length – 20 Diameter -5

10	6024 Tube set , sterile, 3 m long (Marketing Authorization №FSZ (ΦC3) 2010/07320 dated 05.07.2010) 10items/unit		40	Length 3000
11	1958 Spray nozzle for handpiece cleaning		8	Nozzle diameter 8, Nozzle length 55, Nozzle width 20

12	1974 Spray nozzle for care and cleaning of electronic motor 21		11	Nozzle diameter 10, Nozzle length 64, Nozzle width 21 Thread M20
13	1873 Clip set (10 items) for tube set attachment at motor		5	Diameter 8 Height 8 Rib 10 10x8x8

14	1881 Clip set (3 items) for tube set attachment at handpieces and contra angles		3	Diameter 18, Height 20,
----	---	--	---	----------------------------

Table. 2. Handpieces

Surgery handpieces with spray nozzle and quick-lock heads	Reference number	Transmission ratio	Length (mm) up to	Shank diameter	Connection length (for drill/file)	Speed range (revolutions/min)	Max. torque (Nm)	Weight (grams) up to	Material type
	1907nou Surgery HighSpeed contra angle 1:1, FG shank $\varnothing$ 2,35 mm, length 95 mm	1:1	95	$\varnothing$ 2,35 mm	31 mm	300- 80000	1-3	110	Stainless steel 1.4305., Cooling nozzle 1.4401
	1904nou Surgery HighSpeed handpiece 1:1, FG shank $\varnothing$ 2,35 mm, length 95 mm	1:1	95	$\varnothing$ 2,35 mm	31 mm	300- 80000	1-3	82	Stainless steel 1.4305., Cooling nozzle 1.4401

	1950nou Surgery handpiece 1:1, FG shank $\varnothing$ 2,35 mm , length 70 mm	1:1	70	$\varnothing$ 2,35 mm	21 mm	300- 50000	1-6	85	Stainless steel 1.4305., Cooling nozzle 1.4401
	1960nou Surgery contra angle 1:1, FG shank $\varnothing$ 2,35 mm, length 70 mm	1:1	70	$\varnothing$ 2,35 mm	21 mm	300- 50000	1-6	113	Stainless steel 1.4305., Cooling nozzle 1.4401

APPENDIX 3. EMC Manufacturer's Declaration of Conformity

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
Emission and Immunity Tests

Test Item :

Control unit for electronic motor and pump
HighSurg 30

Manufacturer:

NOUVAG AG

Date: 1.6.2017

Leszek Kazana, Area Sales Manager

Without permission of the test centre this report is not permitted to be duplicated in extract.

Report No.: PLE150611

Test Item: Control unit for electronic motor and pump

Type or Model: HighSurg 30

Applicant: NOUVAG AG
St. Gallerstraße 23-25
CH-9403 Goldach

Arrival Date: 09. June 2015

Place of Testing: PRO EMV Labor Strausberg GmbH
Garzauer Chaussee Phone: +493341 - 335255
15344 Strausberg Fax: +493341 - 335368
e-mail: proemv@proemv.de www.proemv.de

Standards: EN 60601-1-2 / 07.2007

Date of Testing: 12. June to 16. June 2015

Procedure: Section 6 (Electromagnetic compatibility)

Test Result: The device meets the requirements for class B.

Tested by: A. Koch
Lab. Manager

Inspected : R. Erxleben
General Manager

2015-06-16

Date, Signature

2015-06-16

Date, Signature

Contents

1. Instruments And Equipment Used	4
2. Equipment under Test (EUT) Description	6
3. Test Set-up and Mode of Operation During the Tests	7
4. Test Program / Summary and Test Results	8
5. Performance Criteria	9
5.1. Description of EUT quality of performance	9
6. Measurements	10
6.1. Conducted Continuous Disturbance at Mains Port	10
6.1.1. Standards	10
6.1.2. Test Description	10
6.1.3. Limits	10
6.1.4. Measured Values	10
6.1.5. Test Result	10
6.2. Radiated Disturbance 30 MHz to 1 GHz	11
6.2.1. Standards	11
6.2.2. Test Description	11
6.2.3. Limits	11
6.2.4. Set-up and Operating Conditions During the Measurements	11
6.2.5. Measured Values	11
6.2.6. Test Result	12
6.3. Harmonic Currents in Low-voltage Supply Systems	13
6.3.1. Standards	13
6.3.2. Test Description	13
6.3.3. Measured Values	13
6.3.4. Test Result	13
6.4. Voltage Fluctuations and Flicker in Low-voltage Supply Systems	16
6.4.1. Standards	16
6.4.2. Description of the Test Level	16
6.4.3. Measured Values	16
6.4.4. Test Result	16
6.5. Immunity to Radiated Electromagnetic Fields	17
6.5.1. Standards	17
6.5.2. Set-up and Operating Conditions During the Measurements	17
6.5.3. Test Description	17
6.5.4. Test Result	17
6.6. Immunity to Conducted Disturbance by RF fields, injected Current	18
6.6.1. Standards	18
6.6.2. Test Description	18
6.6.3. Set-up and Operating Conditions During the Measurements	18
6.6.4. Test Result	18
6.7. Immunity to Fast Transient Disturbances (Bursts)	19
6.7.1. Standards	19
6.7.2. Test Description	19
6.7.3. Set-up and Operating Conditions During the Measurements	19
6.7.4. Test Result	19
6.8. Immunity to Electrostatic Discharge (ESD)	20
6.8.1. Standard	20
6.8.2. Test Description	20
6.8.3. Test Result	21
6.9. Immunity to Surge Pulses	22
6.9.1. Standards	22
6.9.2. Test Description	22
6.9.3. Measurements	22
6.9.4. Test Result	22

6.10. Voltage dips, interruptions and fluctuations immunity tests	23
6.10.1. Standards	23
6.10.2. Test levels	23
6.10.3. Test Result	23
6.11. Power frequency magnetic field immunity test	24
6.11.1. Standards	24
6.11.2. Test Description	24
6.11.3. Set-up and Operating Conditions During the Measurements	24
6.11.4. Test Result	24
6.12. Immunity to variations of power frequency	25
6.12.1. Standard	25
6.12.2. Test Description	25
6.12.3. Set-up and Operating Conditions During the Measurements	25
6.12.4. Test Result	25
7. Photo documentation	34

Other Aspects: Diagrams

Page 26 to 33

1. Instruments And Equipment Used

Conducted disturbance:

	Type	Make	Ser. / Inv. No.
EMI test receiver	ESIB 26	Rohde&Schwarz	SN 100172
EMI test software	ES-K1 V1.60	Rohde&Schwarz	
Network simulation	ESH2-Z5	Rohde&Schwarz	SN 832769/008

Radiated disturbance:

	Type	Make	Ser. / Inv. No.
EMI test receiver	ESIB 26	Rohde&Schwarz	SN 100172
EMI test software	ES-K1 V1.60	Rohde&Schwarz	
Bilog antenna	CBL 6111D	Chase	SN 35421

Harmonic current emissions and voltage fluctuations, flicker:

	Type	Make	Ser. / Inv. No.
Harmonics and flicker equipment	DPA 500	EMTEST	SN 0604-03
AC-voltage source	ACS 500	EMTEST	SN 0804-01
Harmonics and flicker software	ISM DPA V.3.11	EMTEST	

Immunity test (RF field):

	Type	Make	Ser. / Inv. No.
Power amplifier	30S1G3	amplifier research	SN 30778
Power amplifier	50W1000A	amplifier research	SN 14506
Signal generator	SMB100A	Rohde&Schwarz	SN 108799
Test software (RF field)	Cosori V02.03.06	CONFORMITAS	
Power meter	NRVD	Rohde&Schwarz	SN 836519/011
Directional coupler	DC-6180	amplifier research	SN 14375
Thermal power head	NRV-Z 51	Rohde&Schwarz	SN 834519/023
Field strength meter	PMM 8051	PMM	SN 0106
Probe	BA 01	PMM	SN 059
Bilog antenna	CBL 6144	Chase	SN 28870

Immunity test (Conducted disturbance by RF field, current injection):

	Type	Make	Ser. / Inv. No.
Power amplifier	50W1000A	amplifier research	SN 14506
Signal generator	SML03	Rohde&Schwarz	SN 100097
Test software (current injection)	OATS-sX V6.0.9	CONFORMITAS	
Power meter	NRVD	Rohde&Schwarz	SN 836519/011
HF-Voltage head	URV5-Z7	Rohde&Schwarz	SN 860620
Coupling networks	CDN801	MEB	
Coupling clamp	KEMZ 801	MEB	

Immunity test (Burst pulses):

	Type	Make	Ser. / Inv. No.
Burst generator	NSG 3040	Teseq	SN 1933
Software Burst Test	FW 0002.31	Teseq	

Immunity test (ESD):

	Type	Make	Ser. / Inv. No.
ESD simulator	NSG 435	Schaffner	SN 00000599

Immunity test (Surge pulses):

	Type	Make	Ser. / Inv. No.
Surge simulator	NSG 3040	Teseq	SN 1933
Test software (Surge)	FW 0002.31	Teseq	

Immunity test (Voltage dips, interruptions and fluctuations):

	Type	Make	Ser. / Inv. No.
Simulator	NSG 3040	Teseq	SN 1933
Software (Voltage fluctuations)	FW 0002.31	Teseq	

Immunity magnetic LF fields:

	Type	Make	Ser. / Inv. No.
AC-Source	ACS 500	EMTEST	SN 0804-01
Transformer	T 500	peco	SN 6018
Helmholtz coils	HHS1	PRO EMV GmbH	IN 007/96

Immunity to variations of power frequency:

	Type	Make	Ser. / Inv. No.
AC-Source	ACS 500	EMTEST	SN 0804-01

Measuring equipment:

	Type	Make	Ser. / Inv. No.
Absorbing chamber	8,4 x 7,2 x 5,3 [m]	Frankonia	
Antenna mast	MA 240	Deisel	
Rotary disk	DS 415	Deisel	
Controller turn table	HD 100	Deisel	
Controller antenna mast	CO3000	innco systems	

2. Equipment under Test (EUT) Description

The microprocessor-controlled drilling and shaver system HighSurg 30 is used in ENT surgery (e.g. Otosurgery, Head Surgery) and in oral and Maxillofacial Surgery. The device is designed for drilling, milling, screwing and sawing of the bone as well as for shaver procedures in ENT surgery. The HighSurg 30 has an integrated peristaltic pump (roller pump or tube pump) to cool the treated tissue and the used instruments in order to prevent tissue damage.

The test item is designed for use in medical areas (hospital, clinic, doctor's office).

Classification of the EUT into group and class based on CISPR 11:

The HighSurg 30 is an equipment of group1, class B.

Class A equipment is equipment suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Class B equipment is equipment suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes

Technical Data:

Mains voltage (V)	: 100 / 115 / 230 V AC (switchable)
Mains frequency (Hz)	: 50 - 60
Power consumption	: 120 VA
Protecting class	: I
Dimensions (cm)	: 26 x 25 x 11
Case design	: Plastic
Interfaces	: • Power supply "In" • Foot switch • Electronic Motor M1 • Electronic Motor M2 • Protective Earth (PE)
Firmware version	: V1.10dP
Serial number	: 1068U1412R

For any further information see the applicant's system documentation.

3. Test Set-up and Mode of Operation During the Tests

The measurements were performed in a typical test configuration providing a maximum interference capability.
The test set up was done according to the named EMC base standards.

Mode of operation during the tests	:	Normal operating conditions
Mains voltage at tests	:	• 230 V AC / 50 Hz • 100 V AC / 60 (50) Hz
Adjustments at the EUT	:	• Torque 3 Ncm • Pump 100 % • Speed 80000 rpm
Connected lines, Length (meters)	:	• Power supply "In", unshielded (3,0) • Foot switch, shielded (3,0) (Footcontrol vario , NOUVAG AG, SN: 1100F1412R) • Electronic motor M1, shielded (3,0) (NOUVAG AG, SN: 5145M1504S) • Protective Earth (PE), unshielded (1,2)
Indicators for the evaluation of the EUT behaviour	:	• Display • Motor driving • Pump driving

The test item was tested as a table top equipment.

Climatic conditions:

	Required	Actual
Ambient temperature	15 to 35°C	✓
Relative humidity	30 to 60 %	✓
Atmospheric pressure	86 to 106 kPa	✓

Unless otherwise noted these conditions are valid for all following measurements.

Details of device settings and test arrangements may also be seen in the photo documentation.

4. Test Program / Summary and Test Results

The device presented for testing was tested for compliance with the mentioned standards. The following tests (measurements) have been performed on the test item:

Test	Test level	Result
Conducted continuous disturbance at mains port	class B	passed
Radiated disturbance 30 MHz to 1 GHz	class B	passed
Harmonic currents in low-voltage supply systems	class A	passed
Voltage fluctuations and flicker in low-voltage supply systems		passed
Immunity to radiated electromagnetic fields	3 V/m	passed
Conducted immunity, injected currents	3 V	passed
Immunity to fast transient disturbances	1 / 2 kV	passed
Immunity to electrostatic discharge	6 / 8 kV	passed
Immunity to surge pulses	1 / 2 kV	passed
Voltage dips, interruptions and fluctuations		passed
Immunity magnetic LF fields	3 A/m	passed
Immunity to variations of power frequency	49 / 51 59 / 61	passed

5. Performance Criteria

Performance criteria corresponding with EN 60601-1-2:

Under the test the equipment (EUT) shall be able to provide the intended clinical benefit and remain safe. The EUT may exhibit degradation of performance (e.g. deviation from manufacturer's specification). However, the following degradations shall not be allowed if they are related to basic safety or essential performance.

- component failures
- changes in programmable parameters
- reset to factory defaults (manufacturer's presets)
- change of operating mode
- false alarms
- cessation of any intended operation, even if accompanied by an alarm
- initiation of any unintended operation, including unintended or uncontrolled motion, even if accompanied by an alarm
- error of a displayed numerical value sufficiently large to affect diagnosis or treatment
- noise on a waveform in which the noise is indistinguishable from physiologically - produced signals or the noise interferes with interpretation of physiologically - produced signals
- artefact or distortion in an image in which the artefact is indistinguishable from physiologically - produced signals or the distortion interferes with interpretation of physiologically - produced signals
- failure of automatic diagnosis or treatment EUT to diagnose or treat, even if accompanied by an alarm

5.1. Description of EUT quality of performance

The manufacturer has declared that the device has no essential performance based on the definition given in EN 60601-1:2006.

6. Measurements

Note:

The test results shall apply exclusively to the device under test. They shall not represent a generally valid opinion on the properties of the respective products from the running production process.

6.1. Conducted Continuous Disturbance at Mains Port

6.1.1. Standards

EN 60601-1-2 /07.2007

CISPR 11 /2003 + A1 /2004, mod. + A2 /2006

group 1

class B

6.1.2. Test Description

The conducted disturbance at the mains port was tested. It was measured in the frequency range from 0.15 MHz to 30 MHz.

6.1.3. Limits

The limiting values are:

Frequency range	Limit class B (Detector)
0,15 - 0,5 MHz	66 - 56 dB μ V (QP) 56 - 46 dB μ V (AV)
0,5 - 5 MHz	56 dB μ V (QP) 46 dB μ V (AV)
5 - 30 MHz	60 dB μ V (QP) 50 dB μ V (AV)

Remark: for example 66 to 56 dB μ V means, the limiting value falls linearly with the logarithm of the frequency.

6.1.4. Measured Values

Input voltage (U_{in}): 230 V AC / 50 Hz

Measurement	Diagram	Phase	Detector	Remark
1	page 26, 28	L1/N	PK/QP	Test ok
2	page 27, 28	L1/N	AV/AV	Test ok

6.1.5. Test Result

The test item meets the requirements.

Note: No exceeding of the limits was observed.

6.2. Radiated Disturbance 30 MHz to 1 GHz

6.2.1. Standards

EN 60601-1-2 /07.2007

CISPR 11 /2003 + A1 /2004, mod. + A2 /2006

group 1

class B

6.2.2. Test Description

The radiated disturbance was tested. It was measured in the frequency range from 30 MHz to 1000 MHz at a measuring distance of 10 m.

6.2.3. Limits

The QP limiting values are:

Frequency range	Limit class B
30 - 230 MHz	30 dB μ V/m
230 - 1000 MHz	37 dB μ V/m

6.2.4. Set-up and Operating Conditions During the Measurements

The disturbance levels indicated in the measurement diagrams (measured in an absorbing chamber) shall be reduced by approximately 10 dB before comparing them with the limit (conversion of the measuring distance from 3 m to 10 m). The correction values: test assembly attenuation, antenna factor, and cable attenuation, are taken into consideration in the measurement result.

The measurement diagrams with QP verification (decisive for the rating) represent the maximum which could be reached by turning the test item, by variations in the height of the antenna (1 - 4 m), and by changing the antenna polarisation.

Input voltage (U_{in}): 230 V AC / 50 Hz

6.2.5. Measured Values

6.2.5.1. Measured values of the pre-measurements in the absorbing chamber

Antenna height 2.0 m; horizontal and vertical polarization; antenna-to-EUT azimuth 0, 90, 180, 270 [dg]; MaxHold scan

Measurement	Diagram	Detector	Remark
1	page 29	PK	horizontal antenna
2	page 30	PK	vertical antenna

6.2.5.2. Measured values of the final measurement

Measurement	Diagram	Detector	Remark
3	page 31, 33	QP	horizontal antenna, Test ok
4	page 32, 33	QP	vertical antenna, Test ok

Remark:

Due to the radio interference of the environment a final measurement was not possible on the open area test site in 10 m and 3 m measuring distance. Therefore it was measured according to the current edition CISPR 11 /2009, mod. + A1 /2010 on a test setup according to the CISPR 16-1-4 with a measuring distance of 3 m and up to 10 dB increased limiting value.

6.2.6. Test Result

The test item meets the requirements.

Note: No exceeding of the limits was observed.

6.3. Harmonic Currents in Low-voltage Supply Systems

6.3.1. Standards

EN 60601-1-2 /07.2007

IEC 61000-3-2

6.3.2. Test Description

Harmonic Limits for classes A, B, C, D reproduced from IEC 61000-3-2, table 1 to 3.

Special regulation:

There are no limits for following kinds of devices:

- equipments with a power ≤ 75 W, which are not lighting equipments;
- professional equipments > 1 kW;
- Symmetrical controlled electric heating elements ≤ 200 W;
- Dimmer for electric light bulbs ≤ 1 kW.

6.3.3. Measured Values

Measurement	Diagram	Remarks
1	---	Apparent power (VA): 56 Real power (W): 36
2	page 14, 15	Test for class A.

Class D: Personal computers and personal computer monitors; Television receivers
 $75 \text{ W} < P < 600 \text{ W}$.

Class C: Lighting equipment.

Class B: Portable tools.

Class A: Class A apparatus are defined as all apparatus > 75 W that are not classified as Class B, C or D according to IEC 61000-3-2.

6.3.4. Test Result

The test item meets the requirements of class A.

Harmonic Currents (Average)

Hn	I _{eff} [A]	I _{eff} [%]	Limit [A]	Result
1	82.621E-3	100.000		
2	3.039E-3	3.678	1.08	passed
3	66.704E-3	80.735	2.30	passed
4	2.784E-3	3.369	430.00E-3	passed
5	56.725E-3	68.656	1.14	passed
6	1.773E-3	2.146	300.00E-3	passed
7	37.957E-3	45.941	770.00E-3	passed
8	1.165E-3	1.410	230.00E-3	passed
9	23.333E-3	28.241	400.00E-3	passed
10	647.519E-6	0.784	184.00E-3	passed
11	9.873E-3	11.949	330.00E-3	passed
12	274.599E-6	0.332	153.33E-3	passed
13	7.063E-3	8.549	210.00E-3	passed
14	228.707E-6	0.277	131.43E-3	passed
15	7.775E-3	9.410	150.00E-3	passed
16	160.214E-6	0.194	115.00E-3	passed
17	6.083E-3	7.363	132.35E-3	passed
18	212.387E-6	0.257	102.22E-3	passed
19	3.066E-3	3.711	118.42E-3	passed
20	139.162E-6	0.168	92.00E-3	passed
21	1.581E-3	1.914	160.71E-3	passed
22	156.333E-6	0.189	83.64E-3	passed
23	1.922E-3	2.326	146.74E-3	passed
24	136.776E-6	0.166	76.66E-3	passed
25	1.966E-3	2.380	135.00E-3	passed
26	126.932E-6	0.154	70.77E-3	passed
27	1.410E-3	1.706	124.99E-3	passed
28	122.251E-6	0.148	65.71E-3	passed
29	1.113E-3	1.347	116.39E-3	passed
30	119.087E-6	0.144	61.33E-3	passed
31	1.307E-3	1.582	108.87E-3	passed
32	127.754E-6	0.155	57.50E-3	passed
33	1.158E-3	1.402	102.27E-3	passed
34	114.842E-6	0.139	54.12E-3	passed
35	763.429E-6	0.924	96.44E-3	passed
36	129.647E-6	0.157	51.11E-3	passed
37	550.806E-6	0.667	91.21E-3	passed
38	108.404E-6	0.131	48.42E-3	passed
39	655.165E-6	0.793	86.53E-3	passed
40	142.551E-6	0.173	46.00E-3	passed

Harmonic Currents (Maximum)

Hn	I _{eff} [A]	I _{eff} [%]	Limit [A]	Result
1	153.512E-3	100.000		
2	3.318E-3	2.162	1.62	passed
3	129.376E-3	84.277	3.45	passed
4	3.225E-3	2.101	645.00E-3	passed
5	103.842E-3	67.644	1.71	passed
6	2.428E-3	1.581	450.00E-3	passed
7	66.464E-3	43.296	1.15	passed
8	1.769E-3	1.152	345.00E-3	passed
9	34.944E-3	22.763	600.00E-3	passed
10	1.342E-3	0.874	276.00E-3	passed
11	16.655E-3	10.849	495.00E-3	passed
12	763.467E-6	0.497	229.99E-3	passed
13	8.600E-3	5.602	315.00E-3	passed
14	474.977E-6	0.309	197.15E-3	passed
15	11.377E-3	7.411	225.00E-3	passed
16	397.739E-6	0.259	172.50E-3	passed
17	9.096E-3	5.925	198.52E-3	passed
18	425.707E-6	0.277	153.33E-3	passed
19	6.208E-3	4.044	177.63E-3	passed
20	343.746E-6	0.224	138.00E-3	passed
21	3.538E-3	2.305	160.71E-3	passed
22	377.001E-6	0.246	125.46E-3	passed
23	4.121E-3	2.684	146.74E-3	passed
24	338.330E-6	0.220	114.99E-3	passed
25	3.720E-3	2.423	135.00E-3	passed
26	317.354E-6	0.207	106.16E-3	passed
27	2.655E-3	1.730	124.99E-3	passed
28	268.844E-6	0.175	98.57E-3	passed
29	2.063E-3	1.344	116.39E-3	passed
30	250.406E-6	0.163	92.00E-3	passed
31	1.962E-3	1.278	108.87E-3	passed
32	248.625E-6	0.162	86.25E-3	passed
33	1.993E-3	1.298	102.27E-3	passed
34	249.864E-6	0.163	81.18E-3	passed
35	1.473E-3	0.960	96.44E-3	passed
36	238.411E-6	0.155	76.66E-3	passed
37	1.323E-3	0.862	91.21E-3	passed
38	227.880E-6	0.148	72.63E-3	passed
39	1.228E-3	0.800	86.53E-3	passed
40	323.820E-6	0.211	69.00E-3	passed

6.4. Voltage Fluctuations and Flicker in Low-voltage Supply Systems

6.4.1. Standards

EN 60601-1-2 /07.2007

IEC 61000-3-3

6.4.2. Description of the Test Level

Flicker and voltage limits as stated in the mentioned standard.

- P_{st} = Short term flicker, Limit: 1,0
- P_{lt} = Long term flicker, Limit: 0,65
- d_C = Relative steady state voltage change, Limit: 3.3 %
- $d(t)$ = The value $d(t)$ during a voltage deviation must not exceed 3 % for more than 500 ms.
- d_{max} = Maximum relative voltage change, Limit: 4 %

6.4.3. Measured Values

FREQ = 50.0 Hz	
$R_A = 0.24 \Omega$	$X_A = 0.15 \Omega$
$R_N = 0.16 \Omega$	$X_N = 0.10 \Omega$

Flicker Meter Results

	EUT Werte	Limit	Result
Pst	0.028	1.00	passed
Plt	0.028	0.65	passed
dc [%]	0.026	3.30	passed
dmax [%]	0.107	4.00	passed
dt [s]	0.000	0.50	passed

6.4.4. Test Result

The test item meets the requirements

6.5. Immunity to Radiated Electromagnetic Fields

6.5.1. Standards

EN 60601-1-2 /07.2007

IEC 61000-4-3

6.5.2. Set-up and Operating Conditions During the Measurements

The EUT was placed on a wooden rotation desk so that it was positioned in a 0.8 m distance from the chamber's floor and in a 3 m distance in front of the antenna.

The EUT was exposed to radiation from the front, from the right side, from the left side, and from behind.

The frequency range was tested once in each polarization.

Input voltage (U_{in}): 230 V AC / 50 Hz

6.5.3. Test Description

Frequency range	80 - 2500 MHz
Step size	1 % of previous step
Modulation	AM 80%, 1 kHz sinus
Field strength	3 V/m
Polarisation	horizontal and vertical
Performance criteria	see section 5 / 5.1
Dwell time	2 s for each step

6.5.4. Test Result

Test item reactions:

Frequency range [MHz]	Field strength [V/m] Polarisation	Reactions
80 to 90 303 to 370	3 horizontal	The motor speed varies up to two steps (1000 rpm).
80 to 86 170 to 215 250 to 370	3 vertical	The motor speed varies one steps (500 rpm).

No unacceptable loss of performance or loss of data was observed.

After the influence of the disturbance the test item fulfilled its normal function properly.

6.6. Immunity to Conducted Disturbance by RF fields, injected Current

6.6.1. Standards

EN 60601-1-2 /07.2007

IEC 61000-4-6

6.6.2. Test Description

Frequency range	0,15 - 80 MHz
Frequency steps	1%
Modulation	AM 80%, 1 kHz sinus
RF-voltage	3 V (unmodulated wave) *
Injected Signal lines	3 V (modulated wave)
Injected DC lines	3 V (modulated wave)
Injected AC-lines	3 V (modulated wave)
Rating criteria	see section 5 / 5.1
Dwell time	2 s for each frequency step

*) Deviating to IEC 61000-4-6 the test level was calibrated according to EN 60601-1-2 within: -0% to +25%

6.6.3. Set-up and Operating Conditions During the Measurements

For the tests the following connected lines were led through coupling networks (CDN).

Connected Line	Coupling network (0,15 - 80 MHz)
230 V AC-power supply	CDN 801-M3
Protective Earth (PE)	CDN 801-M1
Foot switch	KEMZ 801
Electronic motor M1	KEMZ 801

The coupling network were separately connected to a RF generator for each test run. One RF port of the inactive coupling networks was equipped with 50 R.

6.6.4. Test Result

Test item reactions: No unacceptable loss of performance or loss of data was observed.

During and after the influence of the disturbance the test item fulfilled its normal function properly.

6.7. Immunity to Fast Transient Disturbances (Bursts)

6.7.1. Standards

EN 60601-1-2 /07.2007

IEC 61000-4-4

6.7.2. Test Description

Pulse shape	5 / 50 ns
Repetition frequency	5 kHz
Test duration	> 1 min
Signal lines * (cap. coupling pliers)	1,0 kV (0,5 kV start, 0,5 kV steps to 1,0 kV)
DC lines (direct coupling)	2,0 kV (0,5 kV start, 0,5 kV steps to 2,0 kV)
AC-lines (direct coupling)	2,0 kV (0,5 kV start, 0,5 kV steps to 2,0 kV)
Rating criteria	see section 5 / 5.1

*) Applicable only to ports interfacing with cables whose total length according to the manufacturers functional specification may exceed 3 m.

Input voltage (U_{in}): 100 V AC / 50 Hz and 230 V AC / 50 Hz

6.7.3. Set-up and Operating Conditions During the Measurements

The burst pulses were coupled into the following lines:

Test voltage [kV]	Input coupling
+/- 2,0	AC-power supply (L, N, PE)
+/- 1,0	Foot switch
+/- 1,0	Electronic motor M1

The test voltage was increased in steps from the lowest up to the selected test level.

6.7.4. Test Result

Test item reactions: No unacceptable loss of performance or loss of data was observed.

During and after the influence of the disturbance the test item fulfilled its normal function properly.

6.8. Immunity to Electrostatic Discharge (ESD)

6.8.1. Standard

EN 60601-1-2 /07.2007

IEC 61000-4-2

6.8.2. Test Description

Air discharge	8 kV (2 kV start, 2 kV steps to 8 kV)
Contact discharge	6 kV (2 kV start, 2 kV steps to 6 kV)
Discharges per pole	> 10
Discharge-R	330 R
Discharge-C	150 pF
Performance criteria	see section 5 / 5.1

Test points:	C	A	Test level
horizontal / vertical coupling area	X		a.l.
Housing (Plastic parts)		X	a.l.
Primary power receptacle		X	a.l.
Main switch		X	a.l.
Fuse compartment		X	a.l.
AC power cordset		X	a.l.
Connector PE	X		a.l.
Foot pedal cable connection		X	a.l.
Foot pedal connector housing		X	a.l.
Foot pedal cable		X	a.l.
Motor cable connection (socket) M1	X		a.l.
Motor cable connector housing M1	X		a.l.
Motor cable M1		X	a.l.
Readiness LED M1		X	a.l.
LED socket M1	X		a.l.
Motor cable connection (socket) M2	X		a.l.
Readiness LED M2		X	a.l.
LED socket M2	X		a.l.
Display		X	a.l.
Push buttons		X	a.l.
Stand for irrigation fluid bottle	X		a.l.
Release key for tubing set bracket		X	a.l.
Peristaltic pump (compartment inside)		X	a.l.
Foot switch			
Housing (Plastic parts)		X	a.l.
Foot pedal cable		X	a.l.
Foot pedal		X	a.l.
Push buttons		X	a.l.
Housing bottom	X		a.l.
Electronic motor			
Housing	X		a.l.
Motor cable		X	a.l.

C: Contact discharge
a.l.: All test levels

A: Air discharge

Input voltage (U_{in}): 230 V AC / 50 Hz

6.8.3. Test Result

Test item reactions:

Test level	Test point	Following effects detected
- 8 kV (A)	Readiness LED M1 and M2	Motor and pump function is stopped. Operation locked. After the disturbance the normal operation can be recovered by off- and on-switching of the device.

No unacceptable loss of performance or loss of data was observed.

After the influence of the disturbance the test item fulfilled his normal function properly.

6.9. Immunity to Surge Pulses

6.9.1. Standards

EN 60601-1-2 /07.2007

IEC 61000-4-5

6.9.2. Test Description

Open circuit pulse shape	1,2μs / 50μs
Short circuit pulse shape	8μs / 20μs
Number of pulses	5 pos. and 5 neg.
Repetition pulse	all 30 seconds pulse
Angles	0°, 90°, 270°
AC lines	1 kV sym. 2 kV asym.
Performance criteria	see section 5 / 5.1

6.9.3. Measurements

Input voltage (U_{in}): 100 V AC / 50 Hz and 230 V AC / 50 Hz

AC power port:

sym. (S) = differential mode :L against N

$R_i = 2 \text{ Ohm}$ and $C = 18\mu\text{F}$ (LZ)

asym. (A) = common mode :L, N against PE

$R_i = 12 \text{ Ohm}$ and $C = 9\mu\text{F}$ (HZ)

Number of pulses	Form	Pol.	Up	Angle	Surge Tot.	Result
5 x (L against N)	Surge LZ	+	1000	S 0	5	passed
5 x (L against N)	Surge LZ	-	1000	S 0	10	passed
5 x (L against N)	Surge LZ	+	1000	S 90	15	passed
5 x (L against N)	Surge LZ	-	1000	S 90	20	passed
5 x (L against N)	Surge LZ	+	1000	S 270	25	passed
5 x (L against N)	Surge LZ	-	1000	S 270	30	passed
5 x (L against PE)	Surge HZ	+	2000	A ---	35	passed
5 x (L against PE)	Surge HZ	-	2000	A ---	40	passed
5 x (N against PE)	Surge HZ	+	2000	A ---	45	passed
5 x (N against PE)	Surge HZ	-	2000	A ---	50	passed

The test voltage was increased in steps from the lowest up to the selected test level.

6.9.4. Test Result

Test item reactions: No unacceptable loss of performance or loss of data was observed.

During and after the influence of the disturbance the test item fulfilled its normal function properly.

6.10. Voltage dips, interruptions and fluctuations immunity tests

6.10.1. Standards

EN 60601-1-2 /07.2007

IEC 61000-4-11

6.10.2. Test levels

Input voltage: $U_{in} = 230 \text{ V AC} / 50 \text{ Hz}$
 $U_{in} = 100 \text{ V AC} / 50 \text{ Hz}$

Basic standard	IEC 61000-4- 11			
	dips	dips	interruptions	interruptions
Input voltage U_{in}	U_{in}	U_{in}	U_{in}	U_{in}
Variable voltage U_{var}	$U_{in} -30 \%$	$U_{in} -60 \%$	0 V (-100 %)	0 V (-100%)
Duration of the events t_{event}	500 ms	100 ms	10 ms	5000 ms
Repetition time t_{rep}	20 s	20 s	20 s	20 s
Phase angle t_{phase} [grd]	0, 180	0, 180	0, 180	0, 180
Test duration	3 x per phase position	3 x per phase position	3 x per phase position	3 x per phase position
Performance criteria	see section 5 / 5.1	see section 5 / 5.1	see section 5 / 5.1	see section 5 / 5.1

6.10.3. Test Result

The test item meets the stipulated requirements.

Test item reactions 230 V AC:

Test	Parameter U_{var} / t_{event}	Remarks
interruptions	0 V / 10 ms	During and after the test the EUT functioned as intended.
dips	160 V / 500 ms	Motor and pump speed varies in the rate of the disturbance signal.
dips	90 V / 100 ms	
interruptions	0 V / 5000 ms	EUT shuts down. After the disturbance the normal operation can be recovered by off- and on-switching of the foot switch.

Test item reactions 100 V AC:

Test	Parameter U_{var} / t_{event}	Remarks
interruptions	0 V / 10 ms	During and after the test the EUT functioned as intended.
dips	70 V / 500 ms	Motor and pump speed varies in the rate of the disturbance signal.
dips	40 V / 100 ms	
interruptions	0 V / 5000 ms	EUT shuts down. After the disturbance the normal operation can be recovered by off- and on-switching of the foot switch.

6.11. Power frequency magnetic field immunity test (immunity magnetic LF fields)

6.11.1. Standards

EN 60601-1-2 /07.2007

IEC 61000-4-8

6.11.2. Test Description

Basic standard	IEC 61000-4-8
Mag. field strength	3 A/m
Frequency	50 / 60 Hz
Test duration	> 1 min (every direction)
Rating criteria	see section 5 / 5.1

6.11.3. Set-up and Operating Conditions During the Measurements

The EUT was arranged in the focal point of a Helmholtz coils and exposed to the magnetic field.

The EUT was rotated by 90° in order to expose it to the test field with different orientations

Input voltage (U_{in}): 230 V AC / 50 Hz

6.11.4. Test Result

Test item reactions: No unacceptable loss of performance or loss of data was observed.

During and after the influence of the disturbance the test item fulfilled its normal function properly.

6.12. Immunity to variations of power frequency

6.12.1. Standard

EN 60601-1-2 /07.2007

IEC 60601-1

6.12.2. Test Description

Basic standard	IEC 60601-1	IEC 60601-1
Nominal voltage	230 V	100 V
Nominal frequency	50 Hz	60 Hz
Test parameter	± 1 Hz of nominal frequency	± 1 Hz of nominal frequency
Test frequency	49 Hz (50 Hz - 1 Hz)	59 Hz (60 Hz - 1 Hz)
Test frequency	51 Hz (50 Hz + 1 Hz)	61 Hz (60 Hz + 1 Hz)
Duration of test	> 1 min	> 1 min
Rating criteria	see section 5 / 5.1	see section 5 / 5.1

6.12.3. Set-up and Operating Conditions During the Measurements

Nominal voltage: 100 V AC and 230 V AC

The test was carried out at the nominal voltage and $\pm 10\%$ of the nominal voltage.

6.12.4. Test Result

Test item reactions: During the test -10 % nominal voltage reduced the motor speed of 80000 on 76000 rpm.

No unacceptable loss of performance or loss of data was observed.

After the influence of the disturbance the test item fulfilled its normal function properly.

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





7. Photo documentation


Figure 1 : Test Set-up "Conducted continuous disturbance at mains port"

Figure 2 : Test Set-up "Radiated electromagnetic fields and immunity to radiated electromagnetic fields"


Figure 3 : Test Set-up "Radiated electromagnetic fields and immunity to radiated electromagnetic fields"


Figure 4 : Test Set-up "Radiated electromagnetic fields and immunity to radiated electromagnetic fields"


Figure 5 : Test Set-up "Harmonics and voltage fluctuations, flicker", "Immunity to variations of power frequency"


Figure 6 : Test Set-up "Immunity to conducted disturbance by RF field, current injection"


Figure 7 : Test Set-up "Immunity to fast transient disturbances (Bursts)"


Figure 8 : Test Set-up "Immunity to electrostatic discharge (ESD)"


Figure 9 : Test Set-up "Immunity to surge pulses (surge)"


Figure 10 : Test Set-up "Immunity to voltage dips, interruptions and fluctuations"


Figure 11 : Test Set-up "Immunity magnetic LF fields"

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
[Handwritten signature]

Date: 1.6.2017
Leszek Kazana, Area Sales Manager

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