Robert C. Meagher
Senior Director
World Wide Regulatory Affairs
Carestream Health, Inc.
150 Verona Street
Rochester, New York 14608
Telephone 585-627-6528

" 30 " July 2015 r

Carestream

INSTRUCTION FOR USE
CARESTREAM Vita Flex CR System

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Photo, diagram, drawing medical device

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Notice

All rights reserved. No part of this manual may be reproduced or copied in any form by any mean graphic, electronic, or mechanical, including photocopying, typing, or information retrieval systems without written permission.

Use of Manual

CARESTREAM CR Systems are designed to meet international safety and performance standards. Personnel operating the unit must have a thorough understanding of the proper operation of the system. This manual has been prepared to aid medical and technical personnel to understand and operate the system. Do not operate the system before reading this manual and gaining a clear understanding of the operation of the system. If any part of this manual is not clear, please contact your Carestream representative for clarification

Authorized European Representa-tive

Carestream Health France
1, rue Galilée
93192 NOISY-LE-GRAND CEDEX
France

© Carestream Health, 2014 г.
CARESTREAM is a trademark of 150 Verona St. Carestream Health, Inc.

Importer for European Union
Carestream Health Netherlands

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PRODUCT DESCRIPTION:

CARESTREAM Vita Flex CR System.

The CR System is a digital imaging system used for image acquisition and processing of static projection radiography images. The System utilizes a phosphor screen with energy storage capability as an X-ray image receptor.

After exposure, a laser beam stimulates the phosphor screen, creating luminescence (blue light) proportional to the local X-ray exposure. The luminescence signal is captured by a Photo Multiplier Tube (PMT), digitized, and processed for review.

The System enables you to read a storage phosphor screen quickly and erase it for the next exposure.

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INDICATION

The CR System is designed for medical professionals to read phosphor CR screens. The system includes:

- CARESTREAM IMAGE SUITE Software (or approved third party compatible software)
- Plug and play USB2 interface
- Optional screens and cassettes in the following standard sizes: 20 x 25 cm (8 x 10 in.); 25 x 30 cm (10 x 12 in.); 28 x 35 cm (11 x 14 in.); 35 x 35 cm (14 x 14 in.); 35 x 43 cm (14 x 17 in.); 24 x 30 cm.

Note

From this point on, the Carestream CR System will be referred to as "CR System."

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1. Safety and Regulatory Guide

Caution

To avoid physical injuries to yourself or damage to the scanner, do not open the scanner covers and do not try fixing or adjusting any of the internal parts. Perform only user service and maintenance procedures specified in this document.

Introduction

The information contained herein is based on the experience and knowledge relating to the subject matter gained by Carestream Health, Inc. prior to publication.

No patent license is granted by this information.

Carestream Health, Inc. reserves the right to change this information without notice, and makes no warranty, express or implied, with respect to this information. Carestream Health shall not be liable for any loss or damage, including consequential or special damages, resulting from any use of this information, even if loss or damage is caused by Carestream Health's negligence or other fault.

General Safety Guidelines

- This product is designed and manufactured to ensure maximum safety of operation. Operate and maintain it in strict compliance with the safety precautions and operating instructions contained in this manual.
- This product meets all the safety requirements applicable to medical equipment. However, anyone attempting to operate the system must be fully aware of potential safety hazards.
- Do not modify this product in whole or in part without prior written approval from Carestream Health.
- End users are allowed to perform only the operation and maintenance procedures specified in the User Guide (6K7304).
- To ensure safety, read all user manuals carefully before using the system and observe all Cautions, Important, and Notes located throughout the manual.

Keep this manual with the equipment.

- Reading this manual does not qualify you to operate, test, or calibrate this system.
- Unauthorized personnel are not allowed access to the system.
- If the product does not operate properly or fails to respond to the controls as

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described in this manual:

- Follow the safety precautions as specified in this manual.
- Stop using the unit and prevent any changes to it.
- Use only legally marketed cassettes. Check periodically the quality of the cassettes, and replace them if any defects are apparent.
- The images provided by this system are intended as tools for the trained user. They are explicitly not to be regarded as a sole incontrovertible basis for clinical diagnosis.
- Be aware of the product specifications and of system accuracy and stability limitations. Consider these limitations before making any decision based on quantitative values. If you have any doubts, consult your Sales Representative.
- This system is Class I continuous operated stationary equipment without applied parts and has one signal input/output part.
- The appliance coupler of the flexible cord is used as a disconnecting device.
- In order to guarantee medical-grade leakage current, the computer that is connected to the system should be a medical-grade computer or connected to the mains through a medical-grade isolation transformer.
- To achieve grounding reliability in USA installations, the equipment must be connected to an equivalent receptacle marked "Hospital Only" or "Hospital Grade".
- Product is provided with ordinary protection against the harmful ingress of water.
- Not suitable for use in the presence of a flammable anesthetics mixture with air or with oxygen or with nitrous oxide.
- The use of accessory equipment not complying with the equivalent safety requirements of this equipment may lead to a reduced level of safety of the resulting system. Consideration relating to the choice shall include:
 - Use of the accessory in the patient vicinity.
 - Evidence that the safety certification of the accessory has been performed in accordance with IEC 60601-1 or the system to IEC 60601-1-1 or local equivalent.

Important

When the System is not in use, disconnect the power cable from the electrical outlet.

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Electrical Hazards

- Do not touch the power supply or the power supply cable if it appears to be damaged or if any wire is exposed.
- Do not remove or open system covers or plugs. Internal circuits use high voltage capable of causing serious injury.
- Fuses blown within 36 hours of being replaced by a qualified technician may indicate malfunctioning electrical circuits within the system. Have the system checked by qualified service personnel. Do not attempt to replace any fuse.
- Fluids that seep into the active circuit components of the system may cause short circuits that can result in electrical fires. Therefore, do not place any liquid or food on any part of the system.
- To avoid risk of electric shock, this equipment must only be connected to mains supply with protective earth.

Electromagnetic Emissions/Immunity

- Electromagnetic Compatibility Precautions
 - Medical electrical equipment requires special precautions regarding electromagnetic compatibility (EMC). Medical equipment must be installed and put into service according to the EMC information provided in the following documentation.
- Communications Equipment
 - Portable and mobile radio frequency (RF) communications equipment can affect medical electrical equipment EMC performance.
 - The use of cables, accessories or transducers other than those specified below with the exception of transducers or cables sold by the manufacturer of the equipment as replacement parts for internal components, may result in increased emissions or decreased immunity of the medical equipment.Replacement of Cables, Accessories or Transducers
- Other Equipment
 - The Vita Flex CR should not be used adjacent to or stacked on other equipment. If adjacent or stacked use is necessary, verify normal operation in the configuration in which the Scanner will be used.
- Cable, Accessory and Transducer Information for the Vita Flex CR.

Port Type	Port Description	From	To	Cable Type	Cable Length
Power	AC Power	Vita Flex CR	AC Mains	Unshielded	3 m (9.84 ft.)

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Signal	USB	Vita Flex CR	PC	Shielded	3 m (9.84 ft.)
--------	-----	--------------	----	----------	----------------

Guidance and manufacturer's declaration - electromagnetic emissions

Vita Flex CR is intended for use in the electromagnetic environment specified below. The customer or the user of the Vita Flex CR should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Vita Flex CR uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The Vita Flex CR is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings for domestic purposes.
Harmonics emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

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Caution

- Do not operate the equipment in the presence of explosive liquids, vapors, or gases.
- Do not plug in or turn on the system if hazardous substances are detected in the environment.

If these substances are detected after the system has been turned on, do not attempt to turn off the unit or unplug it. Evacuate and ventilate the area before turning off the system.

Overheating

Caution

Do not touch the erase LED unit if hot. Allow the system ten minutes to cool down before touching any of the internal parts.

Vita Flex CR Laser Safety Instructions

The Vita Flex CR is classified as a Class 1 Laser Product (with the outer cover in place).

Laser Warning

- The Vita Flex CR contains a Class 1 laser. During normal operation, always keep the unit enclosed in its protective cover.
 - Do not attempt to remove the cover. Only a qualified technician may remove the cover to service this product.
- Complies with 21 CFR 1040.10 and 1040.11 except for deviations pursuant to Laser Notice No. 50, dated June 24, 2007.
Class 1 Laser Product, and IEC/EN 60825-1.

Device-Specific Safety Information

LIFTING HAZARD

The Vita Flex CR weighs 25 kg (55 lb). Do not try to lift the unit by yourself. Always seek assistance from another person. Lifting equipment that is too heavy may result in serious injury and/or damage to equipment.

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Health and Safety Compliance

Regulatory Information

The Vita Flex CR was examined for compliance and has classifications and licenses as follows:

Safety

USA	ANSI/AAMI ES60601-1:2005 - Part 1: Medical Electrical Equipment General Requirements for safety and essential performance
Canada	CAN/CSA C22.2 № 601.1-M90 - Medical Electrical Equipment — Part 1: General Requirements for safety and essential performance
European Union	EN 60601-1: 2006 — Medical Electrical Equipment — Part 1: General requirements for safety and essential performance EN 60825-1: 2007 — Safety of Laser Products
International	IEC 60601-1: 2005 — Medical Electrical Equipment — Part 1: General requirements for safety and essential performance IEC 60825-1: 2007 — Safety of Laser Products JIS 60601 - 1: 2012 — Medical Electrical Equipment UL 60601 -1: 2003 — Medical Electrical Equipment 1st edition.

EMC

IEC 60601-1-2: 2007 Medical Electrical Equipment - Electro-Magnetic Compatibility Requirements & Tests.
EN 60601-1-2: 2007 Medical Electrical Equipment - Electro-Magnetic Compatibility Requirements & Tests

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CE Conformity

This product conforms to the requirements of Council Directive 93/42/EEC. The Vita Flex CR bears the following mark of conformity.

0086

IEC Symbols Used

The system may have labels with one or more of the following symbols. These symbols indicate the IEC standards to which the system conforms.

Caution – consult accompanying documents

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
APPLICATION

2. Site Preparation

Before installing the CR system, make sure your site meets the following requirements:

Work Area

- Prepare a table with a flat stable surface and without wheels.
- Make sure you have enough free space around the work area.

Power Requirements

Make sure you have at least three wall power outlets to connect the following equipment:

- CR system
- Monitor
- Computer

Using extension cords or power strips is not recommended.

Computer

For the computer specifications, refer to the Specifications for the CARESTREAM Vita Flex CR System 6K9765 located on the documentation DVD that is included with the system.

Internet Access

It is recommended to have internet access to the computer of the system during installation for license activation. However, license activation can be completed from another computer, using email. Email instructions are provided in the Quick Reference Guide. It is required to have internet access to the computer of the system for remote access to and from the system, after the installation.

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
Electromagnetic Immunity

Guidance and manufacturer's declaration - electromagnetic immunity

The VITA FLEX CR is intended for use in the electromagnetic environment specified below. The customer or the user of the Scanner should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.

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Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: U_T is the AC mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration - electromagnetic immunity			
The Vita Flex system is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance

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			$d = 2.33 \cdot \sqrt{P}$ 800 MHz to 2.5 GHz
			where P is the maximum output rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).
			Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey*, should be less than the compliance level in each frequency range†.
			Interference may occur in the vicinity of equipment marked with the following symbol: 
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

* Field from fixed transmitters, such as base station for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Scanner is used exceeds the applicable RF compliance level above, the Scanner should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Scanner.

† Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Recommended separation distance between portable and mobile RF communications equipment and the VITA FLEX CR

The VITA FLEX CR is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the CR System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the CR System as recommended below, according to the maximum output of the communications equipment.

W Rated maximum output power of transmitter	m Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz d = 1.17-P	80 MHz to 800 MHz d = 1.17-P	800 MHz to 2.5 GHz d = 2.33-P
0.01	0.117	0.117	0.233
0.1	0.37	0.37	0.737
1	1.17	1.17	2.33
10	3.7	3.7	7.36
100	11.7	11.7	23.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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
- The use of mobile radio frequency communications equipment can affect the unit.
- Requires the use of special measures to ensure the electromagnetic compatibility (EMC) and must be installed and commissioned in accordance with the information relating to EMC, given in the documentation

3. CR System Components Overview


System LED Indicator

The below table shows status indicator colors and their meaning:

System Status	LED
Powering up and Homing	Orange
Ready for Scan	Steady Green
Scanning	Blinking Green
Error	Blinking Orange
No communication	Blinking Green and Orange

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
Installation Work Flow


Unpack the System

Unpack the box and verify that the following components are in the box:

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

Training and Documentation DVD package	8	USB cable
Image Suite SW DVD and Documentation DVD	9	Cassettes and screens
Quick Reference Guide for the CR system installation	10	Software License Card(s)
Safety and Regulatory Guide for the Vita Flex CR system	11	Screen extraction tool

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CR system	12	X-OMAT screen and rollers cleaning fluid
Power supply	13	ATP phantom
Power cord	14	Cover removal tool (for cable compartment cover)

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Remove the System from the Box

Two people are required to remove the system from the box.

Position the CR System

Place the CR system on a flat surface according to the following requirements:

- The CR System should be placed on a flat stable surface.
- The CR System can be placed in either horizontal or vertical orientation - it will scan in either orientation.

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

Hardware Installation

Prerequisites:

Make sure that the power source for the system is in the range of 100V AC- 240V AC 1.5 A

Note

For optimal results it is recommended that the light level in the room is less than 500 lux to minimize the imaging plates exposure to light. Avoid placing the scanner beneath windows or sources of direct, strong light. For more information on ambient operating condition requirements, see the **VITA FLEX Safety and Regulatory Guide**

Important

The CR system should be positioned so that there is always easy access to the mains power supply 1. Place the CR System on a stable flat surface and turn it onto its back to gain access to the service door at the rear of the scanner.

2. Open and remove the cable compartment cover located in the rear of the scanner. Use the cover removal tool.

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
3. Connect the main power cable and the USB cable to the scanner and close the cable compartment cover.

Note

This unit shall be powered by external dedicated AC/DC adapter, model name MENB1121A2451F02 by SL Power/AULT. This adapter is part of unit.


Note


You Must close the door to the power supply, if you do not close the door to the power supply properly, there is a risk of light penetration resulting in loss of image quality.

4. Plug the USB cable into your computer.

5. Plug the power supply into the Mains wall socket.

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
6. Press and hold the power button until the LED indicator is orange.

Note

The LED indicator will blink green until the system and software installation is complete. When the installation is complete and communication is established, the LED will be solid green.

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Image Suite Software Installation

To operate the system, it is required to install image acquisition software. This chapter provides explanations on how to install Image Suite version 4.0 software.

Note

For installation instructions of image acquisition software other than Image Suite, refer to relevant software company documentation.

1. Insert Image suite installation CD into the PC and follow the setup wizard:
2. In the Login window, type User Name **admin** and Password **1234**.
3. When logging into Image Suite for the first time, you are required to change the password. Change the password to **123456**.

Note

For further information refer to the Image Suite Installation Instructions for the Image Suite V4 Software, 6M1562 provided electronically on the Image Suite Documentation CD. You can also refer to the Quick Reference Guide for the CARESTREAM Vita Flex CR System Installation, 6M1492, that is included with the system.

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Configure Cassette Size

By default, the scanner is set up to work with the following cassette sizes:

- 20 x 25 cm (8 x 10 in.)
- 25 x 30 cm (10 x 12 in.)
- 28 x 35 cm (11 x 14 in.)
- 35 x 35 cm (14 x 14 in.)
- 35 x 43 cm (14 x 17 in.)


Therefore, no changes in scanner settings are required for these cassette sizes.

To use **24 x 30 cm** cassette size instead of 25 x 30 cm (10 x 12 in.) cassette size, scanner settings need to be updated as specified in the following steps:

1. Exit the **Image Suite** software.
2. On the desktop, click the **Image Suite Utility Center** icon.

3. In the Login window that appears, type **User name (admin)** and **Password (123456)**.

4. From the Utility Center side menu, select **CR Tools >>POC Vita System Settings >>POC Vita System Settings**.

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


The Vita CR Settings window with the default scanner settings appears.

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5. Clear the 10" x 12" (10" x 12" / 24cm x 30cm) check box to set up the scanner to work with 24 cm x 30 cm cassettes.

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
Note

To set up the scanner to work with 10" x 12" (25 cm x 30 cm) cassettes again, select the 10" x 12" (10" x 12" / 24cm x 30cm) check box.

6. Click **Save & Exit** to save the settings and close the **Vita CR Settings** window.


The scanner is now set up to work with 24 cm x 30 cm cassettes.

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Prepare the Cassettes

Screens and cassettes are supplied separately. It is required to insert a screen into the cassette before using it in the system. Use the following procedures to insert a screen into a cassette.

Prerequisites:

Remove the screens and cassettes from their packaging. Clean the screen, see Clean the Cassettes and Flexible Phosphor Screens


1. Slide the edge of the extraction tool (PN SK250056) over the back edge of the cassette.

2. Pull the tool down until it enters the slots in the back of the cassette.

3. With the extraction tool in place, fully insert the screen into the cassette.

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

Note

Make sure that the screen is inserted correctly into the cassette:

- The white side of the screen (phosphor side) should face the front side (gray) of the cassette.
- The black side of the screen should face back side (black) of the cassette (with label as shown above).

4. Pull up the release lever and remove the extraction tool.

5. Apply the screen type label, provided with the screen, to the back side of the cassette. Write down the installation date on the label.

Note


Applying the human skeleton label, provided with the screen, is optional

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

Acceptance Test Procedure

The Acceptance Test Procedure (from this point on will be referred as ATP) confirms that the CR system was not damaged during shipment and is ready to begin patient testing.

Note

This test is mandatory. If not performed, a warning message will appear requesting to perform the test.

Access ATP

1. From the desktop of the Image Suite computer, click the **Image Suite 4.0** icon to open the Image Suite application.

2. In the Login window that appears, type **User name (admin)** and **Password (123456)**.

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
3. On the header bar of the Image Suite application, click the **Help** menu button.

4. From the displayed Help menu, select **Launch Acceptance Test**.


Run ATP

1. The **Acceptance Test** window appears with the ATP instructions. Follow the displayed instructions.

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
2. Place the ATP phantom on one of the following cassettes: 14 x 17 in./10 x 12 in./24 x 30 cm.
3. Properly align the edges of the phantom with the edges of the cassette.


4. Expose the cassette with the phantom to a dose of 5 mR using x-ray settings: 60KVp, 1.5mAs, SID = 1.0 m. (40 in.).


Note

If a dosimeter is available, use it for fine-tuning the X-ray to achieve 5 +/- 1 mR, before exposing the cassette.

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5. Insert the exposed cassette into the scanner. Phosphor screen is automatically loaded and scanned.

6. After the scan is complete, the Acceptance Test results are displayed:

- If the Acceptance Test succeeds, the system is ready for use.


- If the Acceptance Test fails, follow the instructions displayed on screen and readjust the exposure level. Then, perform the Acceptance Test again.

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
3. Type your **User name (admin)** and **Password (123456)** and click **Login**. Carestream Image Suite opens


Note

The next chapters provide information only on basic operations using the Image Suite software, such as scan, abort and erase. For more information, refer to the software Help that can be accessed via the Help button on the Notification Bar.

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
Prerequisites:

Important

During the scan, do not disconnect the scanner, the USB or power supply cables and do not press the **On/Off** button. Doing so may result in loss of data.

1. Log in to Image Suite. The **Work list** screen opens.
2. Click **New Patient**.

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
Corestream Image Suite

Study Date: 12/10/2014 Study Time: 15:58 Patient ID: P01080002 Family Name: Given Name: Middle Name:

Date of Birth: Age: Gender:

Body Part List: New Previous

Current View: [Empty Box]

Body Part List: [Empty Box]

Current View: [Empty Box]

3. Fill out the following fields:

- Study Date
- Patient ID
- Family Name
- Given Name
- Middle Name

4. Select a body part.

Corestream Image Suite

Study Date: 12/10/2014 Study Time: 15:58 Patient ID: P01080002 Family Name: Given Name: Middle Name:

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5. In the **Views** drop down list, select a view.

6. Click  **Start**.

7. Expose a patient body part with the desired dose on any supported cassette size.


Note

If you are using Image Suite for the first time, it is recommended that you expose a phantom or any metal object in order to get familiar with the system before clinical use.

8. Insert the cassette into the scanner until the cassette locks into place with:

- The gray side up
- The open edge (with the red arrow) toward the scanner
- The cassette centered exactly in the scanner slot

The screen is then pulled automatically into the scanner and scanning begins. After the scan is complete, the screen is automatically erased and the scanned image appears on screen.

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Note

9. Click the green check mark to approve the image.
10. Click the **End Study** icon to end the study.

Abort Scan

Abort operation will stop the scan and eject the cassette. Perform the abort only if you are sure this operation is necessary (for example, an unexposed cassette was inserted by mistake).

Caution

Performing abort during the scan cycle will result in losing data. Be sure you are aware of the consequences.

1. While the scan is in progress, click the Status button

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

The **Scanner Status** window appears.

2. Click the **Abort** button.

3. In the **Image Suite** dialog box that appears, click **OK** to confirm the action.

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
4. Click OK to complete the abort operation.


Be aware that you may need to perform a manual erase of the cassette. Manually Erase the Cassette.


Manually Erase the Cassette

Normally, after the scan is complete, the screen is automatically erased and manual erase is not required. If, for any reason, you wish to perform manual erase (for example, the screen was not in use for a long period of time), perform the following procedure.

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
Note

The alternative for manual erase, is performing a normal scan cycle at the end of which the screen will be automatically erased.

1. On the Scanner Status dialog panel, make sure the displayed status is **Ready**


2. Click the status button (Ready).
The **Scanner Status** window appears.
3. Click the **Erase Mode** button.

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




4. Insert the cassette into the scanner.

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
6. Scanner performs the erasing of the cassette.

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
After the erasing is complete, scanner returns to the initial erasing mode state.

7. Two options are available to continue:

- If you want to continue with erasing process, insert another cassette and perform steps 4 - 6.
- To exit the Erase mode, click the exit button to return to the initial Scanner Status window, and then click exit button again to exit the Scanner Status window.

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System Specifications

Physical Specifications

Weight	25 kg (55 lb)
Dimensions (W x D x H)	68.0 x 34.0 x 45.0 cm (26.7 x 13.4 x 17.7 in.)
Screen Type	Flexible GP2 Phosphor Screen
Product Positioning	Horizontal or Vertical
Dynamic Range	0 < 40 mR
Electric Rating	100 - 240 V AC, 50 - 60 Hz
Maximum Power Consumption	120 W
Power Consumption While in "Stand-by"	42 W
Heat Dissipation While in Operation	330 BTU/H (97 W)

Erase Specifications

Erase Range 0 < 40 mR
Erase Time 7 < 80 sec (depending on screen exposure)

Icons Used in this Guide

Caution

Cautions point out procedures that you must follow precisely to avoid damage to the system or any of its components, yourself or others, loss of data or corruption of files in software applications.

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Important

Important highlights critical policy information that affects how you use this manual and this product.

Note

Notes provide additional information, such as expanded explanations, hints, or reminders.

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The materials from which it is made or its medical device basic parts

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Storage and Operation Environments

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6. Maintain and Clean the System (Cassette, Screen & Roller Cleaning)

Clean the Rollers

The rollers need cleaning every 1000 cycles or when there is visible dirt that effects the image quality.

Prerequisites:

Note

When cleaning rollers, it is recommended to clean also the screens.
Read and follow the instructions in the manufacturer's Material Safety Data Sheet (MSDS) prior to use.

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Cleaning Materials


- Several soft, lint-free cloths
- Carestream X-OMAT Screen Cleaner CAT number 845-4977 (473 mL/ 16 Fluid Ounces) or CAT number 103-0428 (177 mL/ 6 Fluid Ounces).


Caution

Cleaning materials other than those recommended can contain chemicals that cause visible or hidden damage to the rollers. Never use isopropyl alcohol (isopropanol, rubbing alcohol) to clean the rollers.

1. On the header bar of the Image Suite application, click the **Help** menu button
2. From the displayed Help menu, select **Access Rollers Cleaning**.

3. Click **OK** to shut down the scanner and continue with cleaning procedure.

4. Follow the instructions displayed on screen.

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


Roller Cleaning

Clean rollers as following:

1. Remove the Screen Transportation Unit.
2. Clean the rollers.
3. Insert the Screen Transportation Unit into the scanner.
4. Turn on the scanner.


- a. Disconnect the power supply cable from the Mains wall socket.
- b. Remove the Screen Transportation Unit, See Remove the Screen Transportation Unit

Caution

When you remove the Screen Transportation Unit, be careful not to touch the erase LED heat sink which may be very hot and have sharp edges. Do Not insert your fingers between the cassette opening and the door.

- c. Thoroughly wipe the rollers' surface with Carestream X-OMAT screen cleaner

- d. Rotate the rollers and wipe the visible surface, until the entire rollers surface is clean.

Carestream

Roller Cleaning

Clean rollers as following:

1. Remove the Screen Transportation Unit.
2. Clean the rollers.
3. Insert the Screen Transportation Unit into the scanner.
4. Turn on the scanner.

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
Note

Before inserting the Screen Transportation unit into the system, ensure that the rollers are dry.

- e. Install the Screen Transportation Unit into the system, See Install the Screen Transportation Unit
- f. Connect the power supply.
- g. Press and hold the power button for 1 second.

5. Click **OK** to approve the performed procedure and to reset the counters and then exit the **Utilities Center** software.

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
Note

If you click **Cancel**, the counters will not be reset and cleaning messages will appear for the next 100 cycles.


Confirm the Rollers Function

1. Switch on the system and verify that the LED on the Scanner is green and that your imaging software shows that the Scanner status is **OK**.
2. Insert an unexposed 35.6 x 43.2 cm (14 x 17 in.) cassette, or if unavailable, any other size cassette and perform a scan.
3. Verify that the cycle completes without error.

Clean the Cassettes and Flexible Phosphor Screens

Remove a Phosphor Screen from the Cassette

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Read and follow the instructions in the manufacturer's Material Safety Data Sheet (MSDS) prior to use.

Cleaning Materials

- Several soft, lint-free cloths
- Carestream X-OMAT Screen Cleaner CAT number 845-4977 (473 mL / 16 Fluid Ounces) or CAT number 103-0428 (177 mL / 6 Fluid Ounces).

Caution


Cleaning materials other than those recommended can contain chemicals that cause visible or hidden damage to the rollers.

Never use isopropyl alcohol (isopropanol, rubbing alcohol) to clean the rollers.

1. Slide the edge of the extraction tool (PN SK250056) over the back edge of the cassette.

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3. Push the slider on the cassette upwards.
A centimeter of the screen sticks out of the cassette.

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
4. Gently pull out the rest of the screen.

Clean the GP-2 Phosphor Screens

Under normal use, the GP-2 Phosphor Screens will eventually show wear. Screen wear can result in image artifacts. This wear may occur from abrasion of the protective overcoat or damage to the surface of the screen.

Prerequisites:

Read and follow the instructions in the manufacturer's Material Safety Data Sheet (MSDS) prior to use.

Cleaning Materials

- Several soft, lint-free cloths
- Carestream X-OMAT Screen Cleaner CAT number 845-4977 (473 mL/ 16 Fluid Ounces) or CAT number 103-0428 (177 mL/ 6 Fluid Ounces).

Caution

Cleaning materials other than those recommended can contain chemicals that cause visible or hidden damage to the screen or cassette and could result in immediate or future image artifacts. Never use isopropyl alcohol (isopropanol, rubbing alcohol) to clean screens or cassettes.

Important

Failure to routinely clean the screens with a recommended solution may create the following result:

- The screen may not move smoothly from the cassette into the reader, causing an error message.
- The screen may not fully return to the cassette at the end of the cycle, preventing the cassette from latching;

no error message appears.

• Image artifacts may appear.

• Rub the cloth with a small amount of screen cleaner.

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

Caution

Do not pour screen cleaner directly onto the screen.

2. Wipe the screen thoroughly, one side at a time, applying pressure to remove persistent contaminations if necessary.
3. Wipe the screen thoroughly dry on both sides with a dry cloth.

Caution

Do not leave the screen to air dry.

Never insert a phosphor screen into a cassette unless the screen and cassette are thoroughly dry.

Load a Phosphor Screen into the Cassette

1. With the extraction tool in place, hold the screen with the white (phosphor) side down and insert the screen so that the dark side is facing the cassette label.
2. Hold the cassette vertically and make sure the screen is fully inserted into the cassette.
3. Pull up the release lever and remove the extraction tool.

Note

If you insert the screen incorrectly:

- The screen may not feed properly from the cassette, causing an error message.
- An exposed image might not be processed.
- The screen may fall out of the cassette.

Disinfect the Equipment, Cassettes, and Screens

Prerequisites:

The equipment or cassette/screen is visibly contaminated with blood or bodily fluids.

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remove the blood/bodily fluids followed by intermediate-level disinfection as indicated prior to use.

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Several soft, lint-free cloths

- For scanner and cassettes—a 1:10 bleach and water solution (one part of 5.25% sodium hypochlorite to 10 parts water)
- For screens—a commercially-prepared equivalent solution of diluted bleach that contains no materials that can cause screen damage.

Cleaning and Disinfection Method: Scanner and Cassettes

1. Dampen the cloth thoroughly with a small amount of the bleach and water solution.

Caution

Do not pour solution directly onto the scanner or cassette.

2. Wipe the scanner or the cassette thoroughly, one side at a time, applying pressure to remove persistent contaminations if necessary.
3. Wipe the scanner, or the cassette thoroughly dry.

Caution

Do not leave to air dry.

Postrequisites:

Note

Be sure all external and internal surfaces are thoroughly dry before reassembling and returning the screens and cassettes to use.

Disinfection Method: Screens

1. Dampen the cloth with a small amount of diluted bleach solution.

Caution

Do not pour solution directly onto the screen.

2. Wipe the screen thoroughly, one side at a time, applying pressure to remove persistent contaminations if necessary.
3. Wipe the screen thoroughly dry.

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7. CR System Service Procedures

CR System has been designed to be very reliable and does not need servicing in normal use. The service procedures in this manual are for the user to clean and maintain the system. Perform only maintenance procedures intended for the user and specified in this chapter.

Important

When doing the procedures in this document, you must use safe work practices and wear the correct personal protective equipment (for example, safety eyewear) according to your company's standard operating procedures

Parts Replacements Table

The below table provides a list of replacement parts with their part numbers and pictures. In case you need to replace one of the listed parts, use the provided part number to order the required part.

Part Number	Picture	Description
SK000455		USB cable—See, Replace the Power Supply and USB cables The maximum cable length is not more than 3 m

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
SK000484		USB & PS Door (cable compartment cover)—See, Replace the Cable Compartment Cover
SK000482		Vita Flex CR Replacement System —Replacement instructions are provided with the replacement system.

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Caution

Use of accessories and cable converters, are not contained in the list, except converters and cables supplied by the device manufacturer as replacement parts for interior parts may lead to an increase in electromagnetic emissions or decreased immunity of the device

Replace the Screen Transportation Unit

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

Remove the Screen Transportation Unit

1. Press the On/Off button to switch off the scanner and ensure that the LED power indicators are off.
2. Disconnect the power supply cable from the Mains wall socket.

When removing the Screen Transportation Unit, after disconnecting the power supply cable, wait for 10 minutes to allow the scanner to cool down and be careful not to touch the erase LED heat sink which may be very hot and have sharp edges. Do Not insert your fingers between the cassette opening and the door. Do Not get your fingers caught between the loader and the screen transportation unit

3. Release the two tightening knobs.

4. Holding the plastic handles, remove the Screen Transportation Unit.

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

Caution

Do Not insert your fingers between the cassette opening and the door. Do Not get your fingers caught between the loader and the screen transportation unit


Install the Screen Transportation Unit

1. Insert the Screen Transportation Unit, and secure the two tightening knobs.
2. Connect the power supply cable to the Mains wall socket.
3. Switch on the system and verify that the LED on the Scanner is green and that your imaging software shows that the Scanner status is **OK**.

Confirm the Scanner Function after Replacing Screen Transportation Unit

1. Perform the System Self Test to make sure system functions properly. See Run System Self Test
2. Perform the Acceptance Test Procedure to make sure the image acquisition functions properly. See Acceptance Test Procedure
3. Open your imaging software.
4. Scan an unexposed 35.6 x 43.2 cm (14 x 17 in.) cassette, or if unavailable, any other size cassette.
5. Verify that the cycle completes without error.

Replace the Power Supply and USB cables

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Remove the Power Supply Cable

1. Press the On/Off button to switch off the scanner and ensure that the LED power indicators are off.

Note

You must make sure that the LED power indicators are off before any service or maintenance procedures.

2. Unplug the power supply from the Mains wall socket.

3. Turn the scanner onto its back.

4. Remove the cable compartment cover at the rear of the CR System. See Remove the Cable Compartment Cover

Caution

Check that the power supply is not hot before touching it. If the power supply malfunctions, it may get very hot.

5. Disconnect the power supply cable and the USB cable.

Note

Make sure that you release the clip before trying to remove the main power connector.

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
Install the Power Supply Cable

1. Connect the power supply cable and the USB cable to the CR System.

Note

- You Must close the door to the power supply, if you do not close the door to the power supply properly, there is a risk of light penetration resulting in loss of image quality
- 2. Install the cable compartment cover at the rear of the CR System. See Install the Cable Compartment Cover
- 3. Plug the power supply into the Mains wall socket
- 4. Switch on the CR System by pressing the On/Off button.

Important


Make sure that you install the door at the rear of the CR System. If you do not install the door at the rear of the CR System, light can penetrate the system and affect image quality.

Replace the Cable Compartment Cover

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
SK000484 Cable Compartment Cover

Tools Required

- Flat head screwdriver
- Cover removal tool

Remove the Cable Compartment Cover

1. Press the On/Off button to switch off the scanner and ensure that the LED power indicators are off.
2. Disconnect the power supply cable from the mains wall socket.
3. Remove the cable compartment cover using cover removal tool.

Install the Cable Compartment Cover

1. Install new cable compartment cover and make sure the Power Supply cable and USB cables are properly inserted into the cover groove

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

2. Connect the power supply cable to the Mains wall socket.
3. Press the On/Off button to switch on the scanner.

Replace the SD Card

SD card contains system licenses. It is very rare that the SD card malfunctions. In case it is lost or you are instructed by Service to replace the SD card, perform the below procedure.

Parts Required


SK000483 SD Card

Tools Required

- Flat head screwdriver

Remove the SD Card

1. Press the On/Off button to switch off the scanner and ensure that the LED power indicators are off.
2. Disconnect the power supply cable from the mains wall socket.
3. Remove the cable compartment cover. See Remove the Cable Compartment Cover
4. Using flat head screwdriver, slightly press on the SD card to remove it from slot, and then manually remove the SD card.

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

SD Card Location

Using flat head screwdriver
press on SD card

Manually remove SD card

Install the SD Card

1. Insert the new SD card into its slot in cable compartment. Pay attention to the SD card direction:


2. Connect the power supply cable and USB cable to the scanner
3. Install the cable compartment cover. See Install the Cable Compartment Cover
4. Connect the USB cable to the work station computer.
5. Connect the power supply cable to the Mains wall socket.
6. Press the On/Off button to switch on the scanner.

Confirm the Scanner Function after Replacing the SD Card

1. Perform the System Self Test to make sure system functions properly. See Run System Self Test
2. Perform the Acceptance Test Procedure to make sure the image acquisition functions

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properly. See Acceptance Test Procedure
3. Open your imaging software.
4. Scan an unexposed 35.6 x 43.2 cm (14 x 17 in.) cassette, or if unavailable, any other size cassette.
5. Verify that the cycle completes without error.

8. Troubleshooting the CR System
General Failures

Symptom	Possible Cause	Solution
Scanner does not turn on, front panel LEDs are off.	- Power supply is not connected to the scanner. - Power supply is not connected to the wall socket. - Power supply failure.	- Connect the power supply to the scanner and to the wall socket. - Replace the power supply.
Front panel LEDs are blinking in orange and the status in the Image Suite notification bar is Not Ready .	- USB cable is not connected. - Failure of the USB port on the PC. - Scanner driver is not recognized by Windows.	- Verify that the USB cable is connected to the PC. - Connect the USB cable to a different USB port on the PC, preferably on the rear of the PC. - Restart Windows.

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System Self Test Failures

Symptom	Possible Cause	Solution
System self test could not proceed since scanner is not connected.	USB cable is not connected to the scanner and/or PC, or the PC is turned off. • Power supply is not connected to the scanner and/or is not connected to the power supply. • Power supply failure.	• Verify that the USB cable and power supply cable are properly connected. • Verify that the scanner is turned on. • Verify that the power supply cable is connected to the wall socket.

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Acceptance Test Procedure Failures

Symptom	Possible Cause	Solution
Need to perform System Acceptance Test	- Software was installed/reinstalled, but the Acceptance Test Procedure was not performed. - Service procedure, such as part replacement and/or recalibration, was performed by qualified technician.	Perform the Acceptance Test Procedure.
Test failed due to "Wrong angle detected".	ATP phantom was not properly aligned with the edges of the cassette.	Properly align the edges of the ATP phantom with the cassette edges.
Test failed due to "Wrong exposure level detected".	Measured Average Pixel Value (APV) is not within the accepted range of 5 +/- 1 mR.	Readjust the X-ray settings according to the displayed instructions:

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Error Messages

Note

The first two digits of the error message indicate the error type. If the first two digits are **70XXX**, the occurred error is System error. If the first two digits are **80XXX**, the occurred error is related to Screen Transportation Unit.

Note

The below table provides possible causes and common solutions to the listed errors. After trying the provided solutions, perform the System Self Test to make sure the system functions properly. Run System Self Test

Symptom	Possible Cause	Solution

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Image Failures

Symptom	Possible Cause	Solution
Exposure index does not match the exposure dose.	- X ray generator does not function properly - System is not properly calibrated.	- X ray generator should be tested by the qualified service engineer. - Recalibrate the scanner (contact Service for further information)
Loss of data (full or partial image)	- Scanner is disconnected from the USB during the scan. - Scanner is disconnected from the power supply during the scan. - External devices other than the scanner or printer, connected to the PC (such as cell phone charger, mp3 player and etc.).	- Short-press the On/Off button to reset the scanner. - Properly secure the USB and power supply connectors. - Verify that in the notification bar, the status is Ready and the front panel LED is green.

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9. System Self Test

Test System Hardware using System Self Test

System Self-Test is an automatic diagnostics tool that can assist in troubleshooting the system malfunctions.

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Access System Self Test

Option 1: From the Image Suite Software

1. On the header bar of the Image Suite application, click the **Help** menu button.
2. From the displayed Help menu, select **Run Scanner Self-test**.

The System Self Test window appears.

Option 2: From Software Error Message

In the ERROR window, click **Run Scanner Self-Test**.
The System Self Test window appears.


Run System Self Test

1. In the System Self Test window that appears, click **Start**.

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System Self Test

Welcome to Vita Flex Self Test, during this procedure system will check functionality of Software and Hardware components of the system.

Before you proceed, please note the following:

- *) A standard size cassette is required to be exposed and scanned
- *) The scanned cassette will be automatically erased - all the data on the screen will be lost
- *) procedure takes a few minutes

Press Start to begin test, Cancel to abort

After the check is completed, the Screen Transportation Test window appears.

System Self Test V1.00

Screen Transportation Test

Please Insert any cassette.

The following procedure will test the scanner's transportation basic functions.

Now testing:

- Screen load
- Screen Erase and Unload
- Eject Cassette

System silently checks the internal components of the scanner. The check may take a few seconds.

2. Insert an unexposed any size cassette into the scanner.
System performs an automatic test of the screen transportation unit.
A green check mark appears next to each tested function.

3. System performs Motors test.
A green check mark appears next to each tested function.

System Self Test V1.00

System Self Test

Welcome to Vita Flex Self Test, during this procedure system will check functionality of Software and Hardware components of the system.

Before you proceed, please note the following:

- *) A standard size cassette is required to be exposed and scanned
- *) The scanned cassette will be automatically erased - all the data on the screen will be lost
- *) procedure takes a few minutes

Press Start to begin test, Cancel to abort

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


4. Expose one of the following cassettes: 14 x 17 in./10 x 12 in./24 x 30 cm. to a dose of 5 mR using x-ray settings: 60KVp, 1.5mAs, SID = 1.0 m. (40 in.). Then, insert the cassette into the system. The Image Readout Test starts.

5. When the progress bar reaches 100%, the test status is indicated by a green check mark (meaning the test succeeded) or by a red cross (meaning the test failed). Do one of the following:
- If the test succeeded, click **Exit** to view the System Self Test Report.

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
- If the test failed, click **Redo** to perform the Image Readout Test again, or exit to view the System Self Test Report.

6. After the Image Readout test is complete, the System Self Test Report window appears showing the test results.


If a success report is displayed (no detected errors), the system functions properly. Click **View Report File** to view/print/save the generated basic report.

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


10. Storage and Operation Environments

	Operation	Transport/Storage
Ambient temperature	+5 to + 45°C (+41 to + 113°F)	-23 to + 66°C (-10 to + 150 °F)
Relative humidity	25 - 81%	10 - 86%
Atmospheric pressure	700 - 1060 hPa	

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
Robert C. Meagher
Senior Director
International Regulatory Affairs
Carestream Health, Inc.
150 Verona Street
Rochester, New York 14608
Telephone 585-627-6528
STATE OF NEW YORK
COUNTY OF MONROE

Signed and attested before me on July 30, 2015 by

CHRISTINE A. BENOIT-AVERILL
NOTARY PUBLIC-STATE OF NEW YORK
No. 01BE6320594
Qualified in Monroe County
My Commission Expires March 09, 2019

Notary Public

Appendix A. Technical specification

Functional parameters:

The spatial resolution	5.8 lp / mm										
Number of screen withstand cycles	Not less than 8000										
The wavelength of the reading laser beam	665 nm										
The diameter of the laser beam	10 microns										
Bit of digitizing	16 bits										
The number of pixels horizontally and vertically for all screens formats	<table> <tr> <td>35x43</td><td>4048x4932</td></tr> <tr> <td>35x35</td><td>4048x4048</td></tr> <tr> <td>25x30</td><td>2868x3460</td></tr> <tr> <td>20x25</td><td>2280x2868</td></tr> <tr> <td>24x30</td><td>2708x3415</td></tr> </table>	35x43	4048x4932	35x35	4048x4048	25x30	2868x3460	20x25	2280x2868	24x30	2708x3415
35x43	4048x4932										
35x35	4048x4048										
25x30	2868x3460										
20x25	2280x2868										
24x30	2708x3415										
Contrast sensitivity at a dose of 2.5 mcg.	2%										
The quantum detection efficiency at near-zero frequency	62%										
Image processing time	from 46 to 85 seconds before the image										
System productivity	At least 30 shots per hour										

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Appendix B. Electromagnetic compatibility

ATTENTION!

The **The Vita Flex CR** requires special measures to ensure EMC and should be installed and put into service in accordance with the EMC-related information specified in operating instructions.

ATTENTION!

The **Vita Flex CR** could be affected by RF communications equipment.

ATTENTION!

Use of accessories, transducers and cables, not specified in the list, other than supplied with **The Vita Flex CR** as replacement parts for internal components can lead to increase of electromagnetic emissions or affect electromagnetic immunity of **The Vita Flex CR**

Table B.1

Guidance and manufacturer's declaration – Electromagnetic emissions		
The Vita Flex CR is intended for use in the electromagnetic environment specified below. The customer or the user of The Vita Flex CR should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions, CISPR 11	Group 1	The Vita Flex CR uses RF energy for internal functions only. RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

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RF emissions, CISPR 11	Class A	The Vita Flex CR is suitable for use in all establishments including domestic buildings and those directly connected to the power distribution network that supplies domestic buildings.
Harmonic emissions, IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions, IEC 61000-3-3	Complies	

Table B.2

Guidance and manufacturer's declaration – Electromagnetic immunity			
The VITA FLEX CR is intended for use in the electromagnetic environment specified below. The customer or the user of The VITA FLEX CR should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge, (ESD) IEC 61000-4-2	±6 kV contact	±6 kV contact	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
	±8 kV air	±8 kV air	

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Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
---	-------	-------	---

NOTE: U_T is the AC mains voltage prior to application of the test level.

Table B.3

Guidance and manufacturer's declaration – Electromagnetic immunity			
The VITA FLEX CR is intended for use in the electromagnetic environment specified below. The customer or the user of the Product should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance

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			Portable and mobile RF communications equipment should be used no closer to any part of The VITA FLEX CR, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF, IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz over the frequency range for RF equipment ^{a)}	3 V	$d=1,2 \sqrt{P}$

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a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the **Scanner** is used exceeds the applicable RF compliance level above, the **Scanner** should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the **Scanner**.

b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Notes:

1. At 80 MHz and 800 MHz, the higher frequency range applies.
2. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Table B.4

Recommended separation distances between portable and mobile RF communications equipment and the **VITA FLEX CR**

The **VITA FLEX CR** is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The buyer or the user of the **CR System** can prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the **CR System** as recommended below, according to the maximum output power of the communications equipment.

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Notes

1. At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
2. For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
3. These rates may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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