State of Minnesota

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SECRETARY OF STATE

APOSTILLE

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1. Country: United States of America

This Public Document: Ranger Blood/Fluid Warming Unit Model 245 Operator's Manual

2. has been signed by Alison Leigh Chulla

3. acting in the capacity of Notary Public, State of Minnesota

4. bears the seal/stamp of Alison Leigh Chulla, Notary Public, State of Minnesota

CERTIFIED

5. at: St. Paul, Minnesota

6. Date: August 24, 2017

7. by: Secretary of State, State of Minnesota

8. File No: 9623765-5

9. Seal/Stamp:

10. Signature:

Steve Simon

Steve Simon
Secretary of State

3M Health Care Business


3M Center
2510 Conway Ave., Bldg. 275-5W-06
St. Paul, MN 55144 U.S.A.
651 733 1110

3M

Russia
Apostille

To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M Health Care Business, St. Paul, MN, USA, I certify that the attached Instruction for Use (IFU) for 3M™ Ranger Irrigation Blood/Fluid Warming Unit model 245 and disposable warming sets in English is truthful and accurate.


Regina Feferman
International Regulatory Affairs Associate
3M Health Care

8/22/17


Date

State of Minnesota)
County of Dakota) ss


Subscribed and sworn to before me this 22nd day of August, 2017.

Alison Leigh Chulla, Notary Signature
My Commission Expires: January 31, 2021

3M

Ranger

Blood/Fluid Warming Unit

Model 245

Operator's Manual

Technical Service and Order Placement (USA)

TEL: 1-800-228-3957

OUTSIDE OF THE USA

Contact Technical Support Service of local authorized representative of the manufacturer.

When contacting Technical Support Service, please, provide serial number of your device. Marking with serial number is located in the lower side of the warming device.

In-warranty repair and exchange (USA)

If your RANGER Infusion Fluid Warming System 245 requires factory servicing, call 3M Customer Support Service. To return a device to 3M for service, first obtain a Return Authorization (RA) number from a customer service representative. Please use the (RA) number on all correspondence when returning a device for service. A shipping carton will be delivered to you at no charge, if needed. Call your local supplier or sales representative to inquire about loaner devices while your device is being serviced.

Disposable warming sets are single use and are not object of maintenance and repair.

OUTSIDE USA

For in-warranty repair and replacement, please contact your local authorized representative.

INTRODUCTION

The Ranger fluid warming system includes a warming unit, Model 245, and a disposable warming sets.

The Ranger infusion fluid warming system, Model 245, is designed to warm blood, blood products, and liquids and deliver these at flow rates from KVO (regime for "open vein maintenance" - from 8 to 15 drops per minute) to 500 mL/min (30,000 ml/hour). At these flow rates, the device maintains fluid output temperatures ranging from 33°C to 41°C (Note: Output temperatures are dependant on input fluid temperatures and flow rate). It takes less than 2 minutes to warm up to the 41°C set point temperature.

Disposable warming sets are sterile, latex-free, single-use-only items and are designed to be used with the reusable warming unit, Model 245. For information about using Ranger fluid warming sets with the Ranger infusion fluid warming unit, Model 245, refer to the "Instructions for Use" attached to each disposable set.

The Ranger fluid warming system should only be used in healthcare facilities by trained medical professionals.

The Ranger fluid warming unit is designed to be mounted to an I.V. pole. A handle located on the top of the unit makes transport easy. When mounted to the I.V. pole, the unit fits easily above the 3M Bair Hugger warming unit.

This manual includes operating instructions and unit specifications for the Ranger infusion fluid warming system.

Indications for use

Indication:

The RANGER fluid warming system is intended to warm blood, blood products, and irrigation and infusion solutions.






Indication for use:

- the RANGER infusion fluid warming system is intended to warm blood, blood products, and infusion solutions.


Contraindications

Not known

Definition of Symbols

Recycle to avoid environmental contamination. This product contains recyclable parts. For information on recycling - please contact your nearest 3M Service Center for advice.

Manufacturer

STERILE|EO Sterile, ethylene oxide

Store at room temperature

DEHP-free

Container quantity

Single use

Latex free

Keep dry

Temperature limits

Explanation of Signal Word Consequences

WARNING.

Indicates a hazardous situation, which, if not avoided, could result in death or serious injury.

CAUTION.

Indicates a hazardous situation, which, if not avoided, could result in minor or moderate injury.

NOTICE.

Indicates a situation, which, if not avoided, could result in property damage only.

WARNING.

1. To reduce the risks associated with hazardous voltage and fire and thermal energy hazards follow requirements below.

- Do not substitute other devices (i.e. Model 247) for the Ranger infusion fluid warming unit or Ranger disposable warming sets.
- Do not continue use of the unit if the over-temperature alert continues to sound and the temperature does not return to the set point temperature. Immediately stop fluid flow and discard the warming set. Have the fluid warming unit tested by a biomedical technician or call 3M Patient Warming technical service.

2. To reduce the risks associated with hazardous voltage and fire follow requirements below.

- Do not modify or service this device, and do not open the warming unit case as there are no user-serviceable parts in the unit.
- Connect power cord to receptacles marked "Hospital Only," "Hospital Grade," or a reliably grounded outlet.

- Use only the power cord specified for this product and certified for the country of use.
 - Do not allow the power cord to get wet.
 - Do not use the Ranger blood/fluid warming system when it appears the unit power cord, or warming set is damaged. Use only 3M specified replacement parts.
 - Keep power cord visible and accessible at all times. The plug on the power cord serves as the disconnect device. The wall socket outlet shall be as close as practical and shall be easily accessible.
3. To reduce the risks associated with air embolism follow requirements below.
- Never infuse fluids if air bubbles are present in the fluid line.

CAUTION.

1. To reduce the risks associated with cross-contamination perform the follow.
 - The cleaning tool provides only superficial cleaning, it does not disinfect or sterilize the interior of the unit.
2. To reduce the risks associated with impact and facility medical device damage follow recommendations below.
 - Clamp the Ranger infusion fluid warming unit to an I.V. pole with a minimum 14" (35.6 cm) radius wheelbase and at a height no higher than 44" (112 cm).
3. To reduce the risks associated with environmental contamination perform the follow.
 - Follow applicable regulations when disposing of this device or any of its electronic components.
4. Do not use for direct cardiac application. When using the Ranger infusion fluid warming system with a central venous catheter (CVC), ensure that the catheter tip does not have direct contact with the heart and ensure that all electrical devices connected to or near the patient have an appropriate leakage current rating for the application. If a CVC tip is found to be in direct contact with a patient's heart, the infusion fluid warming unit, Model 245 should be immediately disconnected from the CVC until the CVC is safely repositioned. Failure to follow these precautions may cause cardiac disturbance and/or patient injury.
5. Follow the AABB Guidelines for the use of blood warming device which caution against warming when administering platelets, cryoprecipitates, or granulocyte suspensions.

NOTICE.

1. Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare professional.
2. To avoid device damage follow the indications below:
 - Do not clean the infusion fluid warming unit with solvents. Damage to the case, label, and internal components may result.
 - Do not immerse the infusion fluid warming unit in cleaning or sterilizing solutions. The unit is not liquid proof.
 - Do not insert metallic instruments in the infusion fluid warming unit.
 - Do not use abrasive materials or solutions to clean the heater plates.
 - Do not allow spills to dry inside the unit, as this may make it more difficult to clean the unit.
3. The Ranger infusion fluid warming unit meets medical electronic interference requirements. If radio frequency interference with other equipment should occur, connect the unit to a different power source.

Product Description

The Ranger infusion fluid warming system consists of a Model 245 warming unit and a sterile disposable warming set.

The Ranger infusion fluid warming system


The infusion fluid warming unit is a compact, lightweight, liquid-resistant device with a clamp located on the side for attachment to an I.V. pole. A carrying handle on the top of the unit makes it easy to transport. Located on the front panel you will find:

- Alphanumeric display that indicates the heater temperature during normal operation.

In an over-temperature condition, the display alternately flashes a temperature of 43°C or higher and the word "HI."

An audible alert also sounds. In an undertemperature condition, the display alternately flashes a temperature of 33°C or lower and the word "LO."

- Alert indicator light that comes on when either an over- or undertemperature condition occurs.

Ranger infusion fluid warming disposable sets

Infusion fluid warming sets include a cassette, tubing and connectors compatible with hospital standard I.V. infusion fluid administration sets, bubble trap, and injection port. The high flow set also includes 2 infusion fluid bag spikes and a 150 micron filter located in a central drip chamber.

Ranger infusion fluid warming set specification

	Standard Flow Fluid Warming Set		High Flow Fluid Warming Set	Pediatric Fluid Warming Set
Model	24200	24250	24355	24450
Flow rate	KVO-9,000 ml/h	KVO-9,000 ml/h	KVO-3,000 ml/h	KVO-6,000 ml/h
Priming volume	39 ml	44 ml	150 ml	20 ml
Spike	None	None	2 spikes	None
Heat exchanger, dimensions	Length 12.662 inches Width 5.825 mm	Length 12.662 inches Width 5.825 mm	Length 12.662 inches Width 5.825 mm	Length 12.662 inches Width 5.825 mm
Length of patient line	76 cm	145 cm	152 cm	60 cm
Inner diameter of patient line	0.118±0.003 inches	0.118±0.003 inches	0.189±0.004 inches	0.098±0.003 inches
Outer diameter of patient line	0.161±0.003 inches	0.161±0.003 inches	0.268±0.005 inches	0.161±0.003 inches
Drip chamber	None	None	Yes	None
Bubble trap	None	None	Yes	None
Filter (Model 90029)	None	None	Yes	None
Pressure infusion	Up to 300 mm Hg	Up to 300 mm Hg	Up to 300 mm Hg	N/A
Sterile	EtO	EtO	EtO	EtO

Sterilization validity	3 years	3 years	3 years	3 years
Use	Single use	Single use	Single use	Single use

Standard Flow Fluid Warming Set (Model 24200)

Standard Flow Fluid Warming Set are designed for cases where a moderate flow of fluid is required

- For flow rates up to 150 ml/min (9000 ml/hour)
- Priming volume: 39 ml
- Patient line: 76 cm
- STERILE and nonpyrogenic fluid pathway. Do not use if sterile barrier system has been breached
- For use with Range infusion fluids warming unit
- Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare professional
- Not made with natural rubber latex

Instructions For Use

Do not use if package has been previously opened or is damaged. Refer to the Range infusion fluid warming unit Operator's Manual for set up and use of the warming system.


WARNING

- To reduce the risks associated with air embolism and incorrect routing of fluids follow recommendations below:
- Never infuse fluids if air bubbles are present in the fluid line.
- Ensure all luer connections are tightened.
- Do not mount the Ranger blood fluid warming system higher than patient level, as air embolism may result.
- This product is designed for single patient use only. To reduce the risk of cross contamination or infection do not reuse.
- To reduce the risk associated with patient exposure to leakage current that exceeds Type BF equipment limits, do not use in combination with other medical equipment without verifying that the total leakage current from the combined equipment does not exceed the safety limits for Type BF equipment.
- To reduce the risks associated with potential blood loss, do not use in combination with an extracorporeal circuit.

CAUTION

- Do not exceed 300 mmHg pressure.
- Follow the AABB Guidelines for the use of blood warming devices which caution against warming when administering platelets, cryoprecipitate, or granulocyte suspensions.
- To reduce the risks associated with potential material toxicity, do not use in combination with an extracorporeal circuit as the Ranger infusion fluid warming system is not designed for use with a circulating blood system; therefore, biocompatibility testing has not been completed to assess genotoxicity.

Illustration of Disposable Set

Instructions

Remove excess air from I.V. solution containers before spiking bags.

1. Before priming, slide the fluid warming cassette into the slot in the warming unit.
2. Connect the short inlet line (A) on the Ranger warming set to the fluid source.
3. Prime all tubing to purge air from the line. Invert the bubble trap (D) while priming and fill.
4. Turn the bubble trap right side up and prime the patient line (F).
5. Place the bubble trap into the holder on the warming unit.
6. Close all clamps.
7. Turn the warming unit on.

The fluid warming set is now ready for use.

Removing Air from the Disposable Set

1. Close the clamp (H) between the injection port (G) and the patient connection (I).
2. Invert the bubble trap.
3. Insert a syringe in the injection port and aspirate air until the bubble trap and I.V. line are free of air.
4. Place the bubble trap right side up in the holder.
5. Open the clamp and continue the infusion.

Removing the Disposable Set

1. Close the inlet clamp (B) and open the patient line clamps (E, H).
2. Disconnect the warming set from the fluid source (if applicable).
3. Allow fluid to flow into the patient (this may take 2-3 seconds).
4. Remove the fluid warming cassette from the warming unit and discard according to institutional protocol.
5. Reconnect the patient I.V. line to the fluid source to continue the infusion without warming.

Standard Flow Fluid Warming Set with Extension (Model 24250)

Standard Flow Fluid Warming Set with Extension is designed for RANGER infusion fluid warming device application in those cases where a moderate flow of fluid is required.

- For flow rates up to 150 ml/min (9000 ml/hour).
- Priming volume: 44 ml.
- Patient line: 76 cm + plus 76 cm extension.
- STERILE and nonpyrogenic fluid pathway. Do not use if sterile barrier system has been breached.
- For use with Range infusion fluids warming unit.
- Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare professional.
- Not made with natural rubber latex.

Instructions For Use

Do not use if package has been previously opened or is damaged. Refer to the Range infusion fluid warming unit Operator's Manual for set up and use of the warming system.

WARNING

- To reduce the risks associated with air embolism and incorrect routing of fluids follow recommendations listed below:
 - Never infuse fluids if air bubbles are present in the fluid line.
 - Ensure all luer connections are tightened.
- Do not mount the Ranger infusion fluid warming system higher than patient level, as air embolism may result.
- This product is designed for single patient use only. To reduce the risk of cross contamination or infection do not reuse.


- To reduce the risk associated with patient exposure to leakage current that exceeds Type BF equipment limits, do not use in combination with other medical equipment without verifying that the total leakage current from the combined equipment does not exceed the safety limits for Type BF equipment.
- To reduce the risks associated with potential blood loss, do not use in combination with an extracorporeal circuit.

CAUTION

- Do not exceed 300 mm Hg pressure.
- Follow the AABB Guidelines for the use of blood warming devices which caution against warming when administering platelets, cryoprecipitate, or granulocyte suspensions.
- To reduce the risks associated with potential material toxicity, do not use in combination with an extracorporeal circuit as the Ranger blood/fluid warming system is not designed for use with a circulating blood system; therefore, biocompatibility testing has not been completed to assess genotoxicity.

Illustration of Disposable Set

- A) Inlet line.
 B) B, E, H) Clamps.
 C) Fluid warming cassette.
 D) Bubble trap.
 F) Patient line.
 G) Needleless injection port.
 I) Patient connection.

Instructions

Remove excess air from I.V. solution containers before spiking bags.

1. Before priming, slide the fluid warming cassette into the slot in the warming unit.
2. Connect the short inlet line (A) to the fluid source.
3. Prime all tubing to purge air from the line. Invert the bubble trap (D) while priming and fill.
4. Turn the bubble trap right side up and prime the patient line (F).
5. Place the bubble trap into the holder on the warming unit.
6. Close all clamps.
7. Turn the warming unit on.

The fluid warming set is now ready for use.

Removing Air from the Disposable Set

1. Close the clamp (H) between the injection port (G) and the patient connection (I).
2. Invert the bubble trap.
3. Insert a syringe in the injection port and aspirate air until the bubble trap and I.V. line are free of air.
4. Place the bubble trap right side up in the holder.
5. Open the clamp and continue the infusion.

Removing the Disposable Set

1. Close the inlet clamp (B) and open the patient line clamps (E, H).
2. Disconnect the warming set from the fluid source (if applicable).
3. Allow fluid to flow into the patient (this may take 2-3 seconds).
4. Remove the fluid warming cassette from the warming unit and discard according to institutional protocol.
5. Reconnect the patient I.V. line to the fluid source to continue the infusion without warming.

Pediatric Fluid Warming Set with Fluid Aspiration Port (Model 24450)

Pediatric Fluid Warming Set with Fluid Aspiration Port is intended for use with RANGER Infusion Fluids Warming Unit in the patients who require less infusion volumes, in particular, a group of pediatric / neonatal patients.

- For flow rates up to 100 mL/min (6000 ml/hour)
- Priming volume: 20 mL
- Patient Line: 60 cm
- STERILE and nonpyrogenic fluid pathway. Do not use if sterile barrier system has been breached.
- For use with the Ranger Infusion Fluid Warming System
- Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare professional
- Not made with natural rubber latex

Instructions for Use

Do not use if package has been previously opened or is damaged.

Refer to Ranger Infusion Fluid Warming Unit Operator's Manual for set up and use of the warming system.


WARNING:

- To reduce the risks associated with contamination, biohazard exposure and delay of therapy, do not use if caps are not in place and ensure that all luer connections are securely tightened prior to priming set.
- To reduce the risks associated with air embolism and incorrect routing of fluids:
 - Never infuse fluids if air bubbles are present in the fluid line.
 - Ensure all luer connections are tightened.
- Do not mount the Ranger infusion fluid warming system higher than patient level, as air embolism may result.
- This product is designed for single patient use only.
- To reduce the risk of cross contamination or infection do not reuse.
- To reduce the risk associated with patient exposure to leakage current that exceeds Type BF equipment limits, do not use in combination with other medical equipment without verifying that the total leakage current from the combined equipment does not exceed the safety limits for Type BF equipment.
- To reduce the risks associated with potential blood loss, do not use in combination with an extracorporeal circuit.

CAUTION:


- Do not use with pressure, use only gravity.
- Do not inject fluids or medication through aspiration ports.
- Follow the AABB Guidelines for the use of blood warming devices which caution against warming when administering platelets, cryoprecipitate, or granulocyte suspensions.
- To reduce the risks associated with potential material toxicity, do not use in combination with an extracorporeal circuit as the Ranger blood/fluid warming system is not designed for use with a circulating blood system; therefore, biocompatibility testing has not been completed to assess genotoxicity.

Illustration of Disposable Set

- A. Fluid warming cassette
- B. Inlet clamp
- C. Inlet line
- D. Patient line
- E. Fluid aspiration port
- F. Air aspiration port
- G. Bubble trap
- H. Outlet clamp

Instructions:

See corresponding illustrations at left

Remove excess air from I.V. solution containers before spiking bags.

1. Before priming, slide the fluid warming cassette (A) into the slot in the warming unit.

2. Connect the inlet line (B) on the Ranger warming set to the fluid source.

3. Prime all tubing to purge air from the line. Invert the bubble trap (G) while priming and fill. Turn the bubble trap right side up and prime the patient line (D).

4. Place the bubble trap into the holder on the warming unit.

5. Close the outlet line clamp and turn the warming unit on. The pediatric infusion disposable cassettes with Aspiration Port is now ready for use.

Removing the Disposable Set

1. Close the inlet clamp proximal to the fluid warming cassette and open the outlet clamp distal to the cassette.

2. Allow fluid to flow into the patient (this may take 2-3 seconds).

3. Remove the fluid warming cassette from the warming unit and discard according to institutional protocol.

High Flow Fluid Warming Set (Model 24355)

High Flow Fluid Warming Set is designed for use with RANGER infusion fluids warming device for those cases when intensive fluid flow is required.

- For flow rates up to 500 mL/min (30000 ml/h)
- Priming volume: 150 mL
- Patient line: 76 cm and 76 cm extension
- Drip rate: 20 drips per mL
- STERILE and nonpyrogenic fluid pathway. Do not use if sterile barrier system has been breached.
- For use with a Ranger Infusion Fluid Warming System.
- Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare professional.
- Not made with natural rubber latex

Instructions for Use

Do not use if package has been previously opened or damaged. Refer to the Ranger Infusion Fluid Warming System Operator's Manual for set up and use of the warming system.

WARNING


- To reduce the risks associated with contamination, biohazard exposure and delay of therapy, do not use if caps are not in place and ensure that all luer connections are securely tightened prior to priming set.
- To reduce the risks associated with air embolism and incorrect routing of fluids:
 - Never infuse fluids if air bubbles are present in the fluid line.

- Ensure all luer connections are tightened.
- Do not mount the Ranger infusion fluid warming system higher than patient level, as air embolism may result.
- This product is designed for single patient use only. To reduce the risk of cross contamination or infection do not reuse.
- To reduce the risk associated with patient exposure to leakage current that exceeds Type BF equipment limits, do not use in combination with other medical equipment without verifying that the total leakage current from the combined equipment does not exceed the safety limits for Type BF equipment.
- To reduce the risks associated with potential blood loss, do not use in combination with an extracorporeal circuit.

CAUTION

- Do not exceed 300 mmHg pressure.
- Follow the AABB Guidelines for the use of blood warming devices which caution against warming when administering platelets, cryoprecipitate, or granulocyte suspensions.
- To reduce the risks associated with potential material toxicity, do not use in combination with an extracorporeal circuit as the Ranger blood/fluid warming system is not designed for use with a circulating blood system; therefore, biocompatibility testing has not been completed to assess genotoxicity.

Illustration of Disposable Set

- A) Bag spikes
- C) Drip chamber with filter
- D) Inlet clamp
- F) Bubble trap
- H) Outlet clamp

NOTE: Change the filter when it becomes clogged or according to institutional protocol. Additional filters may be obtained from 3M Patient Warming.

Instructions

Remove excess air from I.V. solution containers before spiking bags.

1. Tighten luer connections (B).
2. Before priming, slide the fluid warming cassette (E) into the slot in the warming unit.
3. Adjust warming unit height to patient level or lower.
4. Insert the Ranger set spikes (A) into the I.V. bag ports.
5. Invert the drip chamber (C) and fill 2/3 full. Turn upright.
6. Close the clamp below the bubble trap (F) and allow the bubble trap to fill. Flow will stop when it is properly filled. Gently tap the bubble trap against the warming unit to expel any trapped air, if necessary.
7. Place the bubble trap into the holder (I) on the warming unit.
8. Open the clamp and continue to prime the patient line.
9. Close any clamp to stop flow.
10. Turn the warming unit on.

The fluid warming set is now ready for use.

Removing Air from the Bubble Trap

1. The bubble trap is designed to continuously vent air from the patient line. Monitor fluid lines to ensure that they are air-free.
2. Gently tap the bubble trap against the warming unit to expel any trapped air.

Removing the Disposable Set

1. Close the inlet clamp (D) proximal to the fluid warming cassette and open all clamps distal to the cassette.
2. Disconnect the Ranger set from the fluid source (if applicable).
3. Allow fluid to flow into the patient (this may take 2 to 3 seconds) and close outlet clamp (H).
4. Remove the fluid warming cassette from the warming unit and discard according to institutional protocol.
5. Reconnect the patient I.V. line to the fluid source to continue the infusion without warming.

Product safety features, Model 245

The following chart describes the safety alert features of the Ranger infusion fluid warming unit.

ALERT TYPE	WHAT TO LOOK FOR	DESCRIPTION/CAUSE	ACTION
Over-temperature alert - 43°C	Alert indicator light illuminates and alert sounds, alphanumeric display alternately flashes a temperature of 43°C or higher and the word "HI."	Heater temperature rose to 43°C because of transient conditions.*	Observe alphanumeric display. If temperature does not drop to 41°C (may take a few minutes), discontinue use of unit. Call 3M Patient Warming Customer Service.
Under-temperature alert - 33°C	Alert indicator light illuminates and alert sounds, alphanumeric display alternately flashes a temperature of 33°C or lower and the word "LO."	Heater temperature has fallen to 33°C.	Alerts stop when temperature rises above 33°C. Continue use of unit. If the temperature does not rise above 33°C, unplug the unit and call 3M Patient Warming Customer Service.

* Transient conditions may cause an over-temperature alert condition. These conditions include:

- There was an extreme change in flow rates (e.g., from 500 mL/min to stop flow).
- The unit was turned on and reached the set point temperature of 41°C before the warming cassette was inserted in the device.
- Fluids were prewarmed to above 42°C before being infused.

Warming units Rev A to M

ALERT TYPE	WHAT TO LOOK FOR	DESCRIPTION/CAUSE	ACTION
Independent back-up safety alert - 46°C	• Alert indicator light illuminates and alert sounds, alphanumeric display alternately flashes a temperature of 43°C or higher and the word "HI." • Alphanumeric display is dark, alert sounds (backup safety alert still works even if display is dark).	Heater temperature rose to 46°C. Safety backup system is activated at 46°C and unit shuts off power to the heating plates.	TURN OFF THE UNIT AND UNPLUG IT. Do not use the unit. Discard the disposable set. Call 3M Patient Warming Customer Service.

Warming units Rev N and newer

ALERT TYPE	WHAT TO LOOK FOR	DESCRIPTION/CAUSE	ACTION
Independent back-up safety alert - 44°C	- Alert indicator light illuminates and alert sounds, alphanumeric display alternately flashes a temperature of 43°C or higher and the word "HI." - Alphanumeric display is dark, alert sounds (backup safety alert still works even if display is dark).	Heater temperature rose to 44°C. Safety backup system is activated at 44°C and unit shuts off power to the heating plates.	TURN OFF THE UNIT AND UNPLUG IT. Do not use the blood/fluid warming unit. Discard the set. Call 3M Patient Warming Customer Service.

Instructions for Use

Preparation and setup of the Ranger infusion fluid warming unit

1. Attach the Ranger infusion fluid warming unit to the I.V. pole. Tighten the pole clamp securely.

CAUTION. To reduce the risks associated with impact and facility medical device damage comply with the following recommendations.

- Clamp the Ranger infusion fluid warming unit to an I.V. pole with a minimum 14" (35.6 cm) radius wheelbase and at a height no higher than 44" (112 cm).

2. Slide the cassette into the slot in the warming unit. The cassette can only fit into the device one way.

3. Prime the warming set. For more information about priming the set, refer to instructions provided with the warming sets.

4. Place the bubble trap in the holder.

5. Plug the power cord into an appropriate outlet. Turn the unit ON. In a few seconds, the alphanumeric display will illuminate. It takes less than two minutes to warm up to the 41°C set point temperature.

6. Begin infusion. When infusion is completed, remove the warming set and discard according to institutional protocol.

Removing the warming set from the Ranger infusion fluid warming unit

1. Close the inlet clamp proximal to the cassette and open all clamps distal to the cassette.
2. Disconnect the warming set from the fluid source, if applicable.
3. Allow fluid to flow to the patient (this may take 2-3 seconds). Close a distal clamp.
4. Remove the cassette from the warming unit and discard according to institutional protocol.
5. Reconnect the patient I.V. line to the fluid source to continue the infusion without warming.

Transferring the warming set from one Ranger warming unit to another

1. Follow steps 1-3 above, then remove the warming set from the first warming unit.
2. During transport, keep the clamps closed and do not infuse fluids while the cassette is outside of the warming unit.
3. Slide the cassette into the second warming unit.
4. Ensure that air is removed from the tubing.
5. Open the clamps and continue the infusion.

Troubleshooting

CONDITION	CAUSE	SOLUTION
-----------	-------	----------

Nothing illuminates on the warming unit panel.	• Unit is not turned on, plugged in, or power cord is not plugged into an appropriate outlet. • Unit failure.	• Turn unit on. Make sure the power cord is plugged into the power entry module of the warming unit. Make sure the warming unit is plugged into a properly grounded outlet. • Call 3M Patient Warming Customer Service.
Alert indicator illuminates and alert sounds, alphanumeric display alternately flashes a temperature of 43°C or higher and the word "HI."	Temporary over-temperature condition because. • An extreme change in flow rates occurred (e.g., from 500 mL/min to stop flow). • Unit was turned on and reached set point temperature before cassette was inserted. • Fluids were prewarmed to above 42°C before being run through the warming unit.	• Open flow to reduce temperature. Alerts will stop when the display reads 41°C. The unit is ready to use. • Alerts will stop when the display reads 41°C. The unit is ready to use. • Turn off unit and unplug it. Discontinue infusion of fluids. Do not warm fluids before infusing them through the Ranger warming unit.
Alert sounds, alphanumeric display and alert indicator light go dark.	Primary controller failure. Unit will no longer operate.	Power to heating plates will shut off if temperature rises to 44°C (warming unit Rev N and newer) or 46°C (warming units Rev A to M). Turn unit off and unplug it. Discontinue use of unit. Discard disposable set. Alert will continue to sound if you do not unplug unit. Call 3M Patient Warming Customer Service.
• Unit alerts soon after plugging it in (unit does not have to be turned on for this condition to occur). • Heater temperature rises to 44°C (warming units Rev N and newer) or 46°C (warming units Rev A to M) and unit shuts down soon after plugging it in (unit does not have to be turned on for this condition to occur).	Test screw on bottom of unit is loose or missing.	Make sure test screw is completely tightened. If it is missing, turn unit off and unplug it. Call 3M Patient Warming Customer Service.
Alert sounds but unit has been turned off.	Independent backup safety system has been activated.	Unplug unit. Call 3M Patient Warming Customer Service.
Cannot remove cassette from unit.	• Cassette is too full, fluids are still being infused, or clamp is open proximal to the cassette. • Warming unit is below patient level, creating	• Make sure fluid is drained from cassette before sliding out the cassette, that fluids are no longer infusing, and that clamp is closed proximal to the warming cassette.

	excessive back pressure.	• Raise unit above patient level.
Alert indicator light illuminates and alert sounds, alphanumeric display alternately flashes a temperature of 33°C or below and the word "LO."	Under temperature condition caused by very high flow using very cold fluid, or defective heater/relay.	Alert should stop when temperature rises above 33°C. If alert continues, turn unit off, unplug unit and discontinue use. Call 3M Patient Warming Customer Service
Alphanumeric display reads "Er 4" or "Open."	Open wire on temperature sensor.	Do not use unit. Call 3M Patient Warming Customer Service.
Alphanumeric display reads "Er 5" or "Open."	Electrical interference.	Remove the unit. Refer to biomedical technician or call 3M Patient Warming Customer Service.

CLEANING

Cleaning the Ranger infusion fluid warming unit

Clean the Ranger warming unit on an as-needed basis.

CAUTION

1. To reduce the risks associated with cross-contamination:

- The cleaning tool provides only superficial cleaning, it does not disinfect or sterilize the interior of the unit.

NOTICE

1. To avoid device damage follow the instructions below.

- Do not immerse the warming unit in cleaning or sterilizing solutions. The unit is not liquid proof.
- Do not clean the warming unit with solvents. Damage to the case, label, and internal components may result.
- Do not insert metallic instruments in the warming unit.
- Do not use abrasive materials or solutions to clean the heater plates.
- Do not allow spills to dry inside the unit, as this may make it more difficult to clean the unit.

TO CLEAN THE EXTERIOR OF THE WARMING UNIT FOLLOW THE INSTRUCTIONS BELOW.


1. Disconnect the Ranger warming unit from the power source.
2. Wipe the outside of the unit with warm, soapy water, nonabrasive cleaning solutions, dilute bleach, or cold sterilants. Do not use abrasive materials.
3. Wipe with a dry, soft cloth.

CLEAN THE HEATING PLATES FOLLOW THE INSTRUCTIONS BELOW.

The Ranger hardware cleaning tool is intended to clean both heating plates of the warming unit. It is not necessary to disassemble the warming unit to use the tool.

METHOD OF CLEANING

1. Unplug the warming unit.
2. Unfold the cleaning tool. Wet the foam pads with a nonabrasive solution, such as Alconox™ brand detergent.
3. Insert the tool from the back of the unit and pull the tool all the way out from the front.
4. Rinse the tool with water and repeat 3 times. Dispose cassette in accordance with the rules of the medical protocol.
5. Wipe off the unit to remove excess fluid.

TO CLEAN RESISTANT, DRIED-ON FLUIDS FOLLOW THE INSTRUCTIONS BELOW.

1. Spray a nonabrasive solution inside the slot of the warming unit and let sit for 15-20 minutes.
2. Clean the unit by using the cleaning tool.

Note: If the cleaning tool cannot be inserted through the slot of the warming unit because of excess dried-on fluids, send the unit to a 3M biomedical technician.

Storage

Store all components in a cool, dry place when not in use.

Servicing

There are no user-serviceable parts in the Ranger infusion fluid warming unit. All service must be performed by 3M Patient Warming or an authorized service technician.

Call 3M Patient Warming Technical Service at 800-228-3957 for service information. Outside of the USA contact your local 3M representative.

SPECIFICATIONS

Recommendations and statement of the equipment manufacturer: electromagnetic radiation		
Device, Model 245, is designed for use in electromagnetic environment described below. The buyer or user of the Device, Model 245, should ensure the device use in this environment.		
Radiation test	Compliance	Electromagnetic environment: recommendations
Radiation test CISPR 11	Group 1	Device, Model 245, uses radio frequency energy only to perform its internal functions. In this regard, the level of radio frequency radiation is extremely low and, most likely, will not affect the operation of nearby electronic equipment.
Radio radiation CISPR 11	Class B	
Emission of harmonic current components IEC 61000-3-2	Class A	
Voltage fluctuations and flicker IEC 61000-3-3	Requirements are met	

Recommendations and statement of the equipment manufacturer: electromagnetic stability

Device, Model 245, is designed for use in electromagnetic environment described below. The buyer or user of the Device, Model 245, should ensure the device use in this environment.			
Electromagnetic stability test	Test level according to IEC 60601	Level of compliance	Electromagnetic environment: recommendations
Electrostatic discharges IEC 61000-4-2	contact ± 6 kV air ± 8 kV	contact ± 6 kV air ± 8 kV	Floors should be wooden, cement or covered with ceramic tiles. If the floors are covered with synthetic material, the relative humidity should be at least 30%.
Fast electrical transients or bursts of pulses IEC 61000-4-4	power lines ± 2 kV	power lines ± 2 kV	The quality of the electricity supply network should meet the standard requirements for commercial or medical institutions.
Dynamic changes in supply voltage IEC 61000-4-5	interphase voltage ± 1 kV phase voltage ± 2 kV	interphase voltage ± 1 kV phase voltage ± 2 kV	The quality of the electricity supply network should meet the standard requirements for commercial or medical institutions.
Voltage dips, short-term faults and voltage fluctuations IEC 61000-4-11	$< 5\% U_T$ (dip $> 95\%$ from U_T) for 0.5 period $40\% U_T$ (dip 60% from U_T) for 6 periods $70\% U_T$ (dip 30% from U_T) for 30 periods $< 5\% U_T$ (dip $> 95\%$ from U_T) for 5 seconds	$< 5\% U_T$ (dip $> 95\%$ from U_T) for 0.5 period $40\% U_T$ (dip 60% from U_T) for 6 periods $70\% U_T$ (dip 30% from U_T) for 30 periods $< 5\% U_T$ (dip $> 95\%$ from U_T) for 5 seconds	The quality of the electricity supply network should meet the standard requirements for commercial or medical institutions. If user of the Device, Model 245, should continuously operate the device during power failure, it is recommended to use an uninterruptible power supply or battery as a power source of the device.
Magnetic field of industrial frequency (50-60 Hz) IEC 61000-4-8	3 A/m	3 A/m	Characteristics of magnetic fields of industrial frequency should be at the level characteristic for the typical room in a typical commercial or medical institution.

NOTE. U_T - AC voltage before test level application.

Recommendations and statement of the equipment manufacturer: electromagnetic stability

Device, Model 245, is designed for use in electromagnetic environment described below. The buyer or user of the Device, Model 245, should ensure the device use in this environment.

Electromagnetic stability test	Test level according to IEC 60601	Level of compliance	Electromagnetic environment: recommendations
Induced radio waves IEC 61000-4-6	3 V (mean square) 150 kHz to 80 MHz	3 V (mean square)	Portable and mobile radio frequency communications equipment should not be located near any part of the device, Model 245, including cables, and the distance between them should be as recommended based on calculations using transmitter wave frequency data. Recommended distance between the device and other equipment $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P}$ within the range from 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ within the range from 800 MHz to 2.5 GHz where P – nominal maximum output power of transmitter in watts (W) according to transmitter manufacturer's data, and d is recommended distance between the device and other equipment in meters (m). The field strength of stationary radio frequency transmitters, determined by examining electromagnetic environment^a, should be below the allowable value of voltage for each frequency bands^b. Interference can be available near equipment marked with the following symbol. 
Radiated radio waves IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

NOTE 1. If the frequency of transmitter is within the range from 80 MHz to 800 MHz, a higher frequency band shall be used.

NOTE 2. These recommendations can be applies not to all cases. The propagation of electromagnetic waves is influenced by the absorption and reflection of waves by structures, objects and people.

^a Intensity of the fields of stationary transmitters, for example, base stations for radiotelephones (mobile and wireless) and terrestrial mobile radio installations, amateur radio equipment, broadcast television sets and radio installations with amplitude or frequency modulation, cannot be calculated theoretically with a high degree of accuracy. To evaluate the electromagnetic environment in the presence of stationary radio-frequency transmitters, it is necessary to take into account the results of electromagnetic environment studies. If the measured field strength in the room where the device, Model 245, is used exceeds the above permissible value, the operation of the device, Model 245, should be monitored to ensure its normal operation. In case of any signs of deviations in the device operation,

it is probably required to take additional measures, for example, to change the orientation or location of the device, Model 245.

- b Within the frequency range from 150 kHz to 80 MHz, the field strength should be less than 3 V/m.

Recommended distances between the device, Model 245, and portable and mobile radio frequency communications equipment

The Device, Model 245, is designed for use in electromagnetic environment with controlled radio frequency interference. The buyer or user of the device, Model 245, can contribute to prevention of electromagnetic interference by maintaining minimum distance between the portable and mobile radio frequency communications equipment (transmitters) and the device, Model 245, in accordance with the recommendations below and maximum output power of the communications equipment.

Nominal maximum output power of transmitter, W	Distance depending on the frequency of transmitter, m		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.30
10	3.69	3.69	7.37
100	11.67	11.67	23.30

If nominal maximum output power of transmitter is not specified above, the recommended distance d in meters (m) can be calculated using equation corresponding to the transmitter frequency, where P is the nominal maximum output power of transmitter in watts (W) according to transmitter manufacturer's data.

- NOTE 1. If transmitter frequency is within the range from 80 MHz to 800 MHz, the recommended distance corresponds to the range of higher frequencies.
- NOTE 2. These recommendations can be applied not to all cases. The propagation of electromagnetic waves is influenced by the absorption and reflection of waves by structures, objects and people.

CERTIFICATIONS

IEC/EN 60601-1; CAN/CSA-C22.2 № 601.1

CLASSIFICATION

MEDICAL—GENERAL MEDICAL EQUIPMENT AS TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH UL 60601-1; CAN/CSA-C22.2, No.601.1; ANSI/AAMI ES60601-1:2005 CSA-C22.2 No. 60601-1:08; Control No.4HZ8

Classified under IEC 60601-1 Guidelines (and other national versions of the Guidelines) as Class I, Type BF, Ordinary equipment, Continuous operation. Classified by Underwriters Laboratories Inc. with respect to electric shock, fire and mechanical hazards only, in accordance with IEC/EN 60601-1 and in accordance with Canadian/CSA C22.2, No. 601.1.

Physical characteristics

Dimensions	4.5 inches in height x 7.5 inches in width x 10 inches in length 11 cm x 19 cm x 25 cm
------------	---

	Accuracy $\pm 15\%$
Weight	7 pounds 7 ounces (3.4 kg) Accuracy $\pm 15\%$
Temperature characteristics	
Set point temperature	41 °C
Set over-temperature (initial)	43 °C
Set under-temperature	33 °C
Set over-temperature (secondary)	44 °C (devices for heating ver. N and more new) 46 °C (devices for heating ver. from A to M)
Electrical characteristics	
Input voltage	100-120 or 220-240 VAC
Operating frequency	50/60 Hz
Maximum heating power	900 W
Fuse	2 x T10A-H (250 V) for operation under voltage conditions 100–120 VAC 2 x T6.3A-H (250 B) for operation under voltage conditions 220–240 VAC
Fuse type	Time delay, high breaking
Leakage current	Meets leakage current requirements in accordance with UL 60601-1 and IEC 601-1.
Environmental Condition	
Operating Temperature Range	15 °C to 40 °C (59 °F to 104 °F)
Storage and Transportation Temperature Range	-20 °C to 45 °C (-4 °F to 113 °F)
Operating Humidity	10 to 85% RH, noncondensing
Atmospheric Pressure Range	50 kPa to 106 kPa

Data on the drug products and materials of animal and(or) human origin contained in the medical product

This medical product does not contain drug products or does not introduce medicinal substances.

Human blood derivatives: Medical product does not contain and does not deliver human blood derivatives.

Materials of animal origin: Materials of animal origin are not intentionally added in the product.

Data on software used

Software	Title	Version	Developer	Date of issue
Controller Model 245; 240V	116E09	7E	West Instruments	17th April 2015

DISPOSAL

Methods of disposal:

Packaging/contents should be disposed in accordance with the local/regional/national/ international regulations.

Medical products should be disposed by the local public law organizations in accordance with SanPiN 2.1.7.2790-10 "Sanitary and Epidemiological Requirements for Medical Waste Treatment" waste Class Г).

Accessories should be disposed by the local public law organizations in accordance with SanPiN 2.1.7.2790-10 "Sanitary and Epidemiological Requirements for Medical Waste Treatment" waste Class B".

WARRANTY

3M ESPE warrants this product will be free from defects in material and manufacture.

Except as expressly provided by the law, 3M shall not be liable for any adverse consequences or damage arising out from the use of this product: direct, indirect, willful, incidental or consequential, regardless of the explanations provided, including warranties, contracts, negligence or objective responsibility.

The warranty period for this Device is 24 months.

The warranty period shall be calculated from the date of equipment installation by the authorized Manufacturer (manufacturer's representative), service organization (if required for equipment operation) or from the date of transfer to the distributor / buyer by 3M.

Manufacturer

3M Company, 3M Health Care, USA

3M Center, 2510 Conway Avenue, Bldg. 275-5W-06, Saint Paul, Minnesota, 55144, USA

Complaints

In case of complaint, contact Authorized Representative of the manufacturer of the medical device on the territory of the Russian Federation:

ZAO "3M Russia"

Russia 121614, Moscow, Krylatskaya str. bldg.17/3

+7 (495) 784 7474 / +7 (495) 784 7475 (fax)

3mra@mmm.com

Data on the last revision of the operational documentation

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APOSTILLE

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1. Country: United States of America

This Public Document: Ranger Irrigation Fluid Warming Unit Model 247 Operator's Manual

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3M

To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M Health Care Business, St. Paul, MN, USA, I certify that the attached Instruction for Use (IFU) for 3M™ Ranger Irrigation Fluid Warming Unit model 247 and disposable warming sets in English is truthful and accurate.

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

**Ranger
Irrigation Fluid Warming Unit
Model 247
Operator's Manual**

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If your RANGER Irrigation Fluid Warming System, Model 247, requires factory servicing, call 3M Customer Support Service. To return a device to 3M for service, first obtain a Return Authorization (RA) number from a customer service representative. Please use the (RA) number on all correspondence when returning a device for service. A shipping carton will be delivered to you at no charge, if needed. Call your local supplier or sales representative to inquire about loaner devices while your device is being serviced.

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INTRODUCTION

RANGER Irrigation Fluid Warming System includes a warming unit, Model 247, and a disposable warming sets.

RANGER Irrigation Fluid Warming System, Model 247, is designed to warm irrigation fluids and deliver these at flow rates under gravity up to 867 mL/min when the bag is hung 100 cm above the scope. At these flow rates, the device maintains fluid output temperatures ranging from 33°C to 41°C (Note: This is for room temperature fluids only). It takes less than 2 minutes to warm up to the 41°C set point temperature.

Alert points on the Model 247 are 48°C and 50°C.

Fluid warming disposable set is sterile, latex-free, single-use-only and is designed to be used with the reusable Model 247 warming unit. For information about using Ranger irrigation fluid warming disposable sets with the Model 247 irrigation fluid warming unit, refer to the "Instructions for Use" included with each disposable component.

The Ranger irrigation fluid warming system should only be used in healthcare facilities by trained medical professionals.

Indications for use

Indication:

RANGER Irrigation Fluid Warming System is intended to warm blood, blood products, and irrigation and infusion solutions.

Indication for use:

- RANGER Irrigation Fluid Warming System is intended to warm irrigation solutions.

Counterindications

Not known

Definition of Symbols

	ON
	OFF
	Fuse
	CAUTION
	Type BF Applied Part
	Voltage, Alternating Current (AC)
	An equipotentiality plug (grounded) conductor other than a protective earth conductor or a neutral conductor, providing a direct connection between electrical equipment and the potential equalization busbar of the electrical installation. Please consult IEC 60601-1; 2005 for requirements.
	This system is subject to European WEEE Directive 2002/96/EC. This product contains electrical and electronic components and must not be disposed of using standard refuse collection. Please consult local directives for disposal of electrical and electronic equipment.
	Protective earth ground
	Consult instructions for use
	Follow instructions for use

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WARNING.

Indicates a hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION.

Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

NOTICE.

Indicates a situation which, if not avoided, could result in property damage only.

WARNING.

1. To reduce the risks associated with hazardous voltage and fire and thermal energy hazards follow the requirements given below.
 - Do not substitute other devices (i.e. Model 245) for the Ranger irrigation fluid warming unit or Ranger irrigation disposable sets.
 - Do not continue use of the unit if the over-temperature alert continues to sound and the temperature does not return to the set point temperature. Immediately stop fluid flow and discard the disposable set. Have the warming unit tested by a biomedical technician or call 3M technical service.
2. To reduce the risks associated with hazardous voltage and fire:
 - Do not modify or service this device, and do not open the warming unit case - there are no user serviceable parts in the Ranger warming unit.
 - Connect power cord to receptacles marked "Hospital Only," "Hospital Grade," or a reliably grounded outlet.

- Use only the power cord specified for this product and certified for the country of use.
 - Do not allow the power cord to get wet.
 - Do not use the Ranger warming system when it appears the warming unit, power cord, or disposable is damaged. Use only 3M specified replacement parts.
 - Keep power cord visible and accessible at all times. The plug on the power cord serves as the disconnect device. The wall socket outlet shall be as close as practical and shall be easily accessible.
3. To reduce the risks associated with air embolism:
- Never use fluids if air bubbles are present in the fluid line.

CAUTION.


1. Not for I.V. use.
2. To reduce the risks associated with cross-contamination:
 - The cleaning tool provides only superficial cleaning, it does not disinfect or sterilize the interior of the unit.
3. To reduce the risks associated with impact and facility medical device damage:
 - Clamp the Ranger irrigation fluid warming unit to an I.V. pole with a minimum 14" (35.6 cm) radius wheelbase and at a height no higher than 44" (112 cm).
4. To reduce the risks associated with environmental contamination:
 - Follow applicable regulations when disposing of this device or any of its electronic components.

NOTICE:

1. Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare professional.
2. To avoid device damage follow indications below:
 - Do not clean the warming unit with solvents. Damage to the case, label, and internal components may result.
 - Do not immerse the warming unit in cleaning or sterilizing solutions. The unit is not liquid proof.
 - Do not insert metallic instruments in the warming unit.
 - Do not use abrasive materials or solutions to clean the heater plates.
 - Do not allow spills to dry inside the unit, as this may make it more difficult to clean the unit.
3. The Ranger irrigation fluid warming unit meets medical electronic interference requirements. If radio frequency interference with other equipment should occur, connect the unit to a different power source.

Product Description

The Ranger irrigation fluid warming system consists of a Model 247 warming unit and a sterile disposable fluid warming set.

Ranger Irrigation Fluid Warming Unit


The Ranger Irrigation Fluid Warming Unit is a compact, lightweight, liquid-resistant device with a clamp located on the side for attachment to an I.V. pole. When mounted to the I.V. pole, the unit fits easily above the 3MTM Bair Hugger 500 or 700 series warming unit. A carrying handle on the top of the unit makes it easy to transport.

Located on the front panel you will find:

- Alphanumeric display that indicates the heater temperature during normal operation. In an over temperature condition, the display alternately flashes a temperature of 48°C or higher and the word "HI." An audible alert also sounds. In an under-temperature condition, the display alternately flashes a temperature of 33°C or lower and the word "LO."
- Alert indicator light that comes on when either an over- or under-temperature condition occurs.

Ranger Irrigation Fluid Warming Disposable Set

Irrigation fluid warming sets include a cassette, tubing, connectors, flow chamber and clamps.

	Irrigation disposable cassettes with air valve	Irrigation disposable cassettes	Irrigation standard patient lines	Irrigation universal patient lines
Model	24760	24770	24775	24795
Flow rate	0-865 ml/min	0-865 ml/min	0-865 ml/min	0-865 ml/min
Priming volume	203 ml	203 ml	195 ml	105 ml
Spike	2 spikes	2 spikes	2 spikes	None
Heat exchanger, dimensions	Length 12.662 inches Width 5.825 mm	Length 12.662 inches Width 5.825 mm	None	None
Length of patient line	None	None	101 inches	86 inches
Inner diameter of patient line	0.282±0.005 inches	0.282±0.005 inches	0.282±0.005 inches	0.282±0.005 inches
Outer diameter of patient line	0.377±0.007 inches	0.377±0.007 inches	0.377 ±0.007 inches	0.377 ±0.007 inches
Tip length	-	-	6±0.025 inches	6±0.025 inches
Inner diameter of tip	-	-	0.182±0.005 inches	0.182±0.005 inches
Outer diameter of tip	-	-	0.37±0.005 inches	0.37±0.005 inches

Drip chamber	Yes	Yes	None	Yes
Bubble trap	None	None	None	None
Filter	None	None	None	None
Sterile	EtO	EtO	EtO	EtO
Sterilization validity	3 years	3 years	3 years	3 years
Single use	No	No	Yes	Yes
Luer	Yes	Yes	Yes	Yes
Clamps	Yes	Yes	Yes	Yes

Ranger Irrigation Bottle Spike Set (Model 24760), Ranger Irrigation Standard Spike Set (Model 24770)

For use with the Ranger patient line with spikes (Model 24775)

Irrigation Bottle Spike Sets (model 24760) and Irrigation Standard Spike Set (model 24770) are designed for use with RANGER irrigation fluid warming device at a flow rate up to 865 ml/min under impact of gravity, with suspension of container 100 cm above the object.

- Irrigation fluid warming set
- Priming volume with patient line: 308 mL
- Not made with natural rubber latex.

STERILE and nonpyrogenic fluid pathway. Do not use if sterile barrier system has been breached. For use with the Ranger irrigation warming unit only.

Caution

- Do not exceed 300 mmHg when using this irrigation set. Follow institutional protocol regarding the use of pressure with irrigation procedures.
- Do not use if the end caps or tubing connectors are not in place. Fluid pathway will not be sterile and/or fluids will leak. Use a new irrigation set.
- Avoid bending and twisting the tubing. A bend or twist in the tubing may stop fluid flow.
- Limit use to 6 hours procedural time, independent of the number of patient lines used. Extended use may cause degradation of materials.
- Must be used with the Ranger patient line with spikes (Model 24775).
- Not for I.V. use.

Notice

Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare professional.

Note: Refer to the Ranger irrigation warming system operator's manual to learn how to operate, troubleshoot, maintain, and store the warming system.

Product Description


The irrigation fluid warming disposable set allows you to warm irrigation fluids to normothermia using the Ranger irrigation warming system.

Illustration of Irrigation Set.


- A. Spikes.
- B. White clamps.
- C. Flow chamber.
- D. Connectors.
- E. Blue clamp.
- F. Fluid warming cassette.

Priming the Irrigation Set.

1. Remove excess air from the irrigation solution container(s).

REF 24760

REF 24770

2. Slide the fluid warming cassette into the slot of the Ranger irrigation warming unit.

Note. Ensure cassette is fully inserted into the fluid warming system. The cassette fits into the warming unit only one way.

3. Close white clamps (B1 and B2).

4. Ensure connectors (D) are tight.

5. Remove end caps from spikes on the irrigation set.

6. Connect the irrigation set to the container(s).

7. Obtain the patient line and twist off end caps from irrigation set and patient line.

8. Connect the patient line to the irrigation cassette.

9. Open white clamps (B1 or B2). Ensure roller clamp is open.

10. Invert the flow chamber until it fills at least 3/4 full, and then turn upright.

Note: Remove residual air bubbles by pinching the tubing.

11. Close the white clamp on the patient line when the line is free of air.

12. Connect the patient line to the scope, Foley catheter, or similar device.

13. Raise I.V. pole.

14. Turn the Ranger irrigation warming unit ON. Allow five minutes for the system to warm up before using.

Removing the Irrigation Set

1. Turn the Ranger irrigation warming unit OFF.

2. Close the blue clamp (E), located below the flow chamber, and open the white clamp between the fluid warming cassette and the end of the patient line.

3. Allow a small amount of fluid to drain from the end of the patient line to reduce pressure in the cassette.

4. Remove the cassette from the warming unit and discard the cassette according to institutional protocol.

Ranger Irrigation Warming System Patient Line with Spikes (Model 24775)

For use with Ranger Irrigation fluid warming Set, Models 24760/24770.

Ranger Irrigation Patient Line Std Spike, Model 24775, is designed for use with Ranger Irrigation Disposable Cassettes, Models 24760/24770, and RANGER device for irrigation fluids warming.

- For flow rates up to 865 mL/min
- For use with a Ranger irrigation warming unit with accessories
- Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare professional
- Not made with natural rubber latex

Sterile and nonpyrogenic fluid pathway. Do not use if sterile barrier system has been breached.


For use with a Ranger irrigation warming Unit.

Warning

This product is designed for single patient use only. To reduce the risk of cross contamination or infection do not reuse.

Cautions

- Do not exceed 300 mmHg when using this irrigation set. Follow institutional protocol regarding the use of pressure with irrigation procedures.
- Do not use if the end caps or tubing connectors are not in place. Fluid pathway will not be sterile and/or fluids will leak. Use a new irrigation set.
- Avoid bending and twisting the tubing. A bend or twist in the tubing may stop fluid flow.
- Limit use to 6 hours procedural time, independent of the number of patient lines used. Extended use may cause degradation of materials.
- Must be used with the Ranger irrigation fluid warming sets (Models 24760/24770/24780).

- Not for I.V. use.

Note: Refer to the Ranger irrigation fluid warming system operator's manual to learn how to operate, troubleshoot, maintain, and store the warming system.

Illustration of Patient Line with Spikes

- A. Scope connector.
- B. Check valve.
- C. Roller clamp.
- D. Patient line.
- E. Pinch clamp.
- F. Warming cassette connector.
- G. Spike connector.

Product Description

The patient line can be attached to a Ranger irrigation fluid warming set or an irrigation solution container with the included spike set.


Instructions

1. Remove tubing cap and attach the patient line to the irrigation fluid warming Cassette.
2. Prime the patient line with irrigation fluid and close one clamp on the patient line.
3. Connect the patient line to the scope, Foley catheter, or similar device.
4. For continued irrigation after warming, remove the patient line from the cassette and attach it to the spike. Attach the spike to the solution container and reprime as needed.

Ranger Irrigation Warming System Patient Line Universal (Model 24795)

For use with 3M Ranger Irrigation Disposable Fluid Warming Set, Models 24760/24770.

Irrigation Patient Line Universal, Model 24795, is designed for use with 3M Ranger Irrigation Disposable Fluid Warming Set, Models 24760/24770, and RANGER device for irrigation fluids warming.

- Priming volume: 105 mL
- Single Use Only
- STERILE and nonpyrogenic fluid pathway.
- Do not use if sterile barrier system has been breached. For use with a Ranger irrigation fluid warming unit.
- Not made with natural rubber latex

Instructions For Use

Cautions

- Do not exceed 300 mmHg when using this patient line for pressure irrigation. Follow institutional protocol regarding the use of pressure with irrigation procedures.
- Do not use if the end caps or tubing connectors are not in place. Fluid pathway will not be sterile and/or fluids will leak. Use a new irrigation set.
- Avoid bending and twisting the tubing. A bend or twist in the tubing may stop fluid flow.
- Limit use to 6 hours procedural time, independent of the number of patient lines used. Extended use may cause degradation of materials.
- Must be used with a 3M Ranger Irrigation Disposable Fluid Warming Set (24760/24770).

- Not for I.V. use.

Note: Refer to the Ranger fluid warming unit Operator's Manual to learn how to operate, troubleshoot, maintain, and store the warming system.

Illustration of Patient Line

- A. Scope connector

- B. Check valve
- C. Roller clamp
- D. Patient line
- E. Pinch clamp
- F. Spike connector

Instructions for Use

1. Remove tubing cap and attach the patient line to the irrigation fluid warming set.
2. Prime the patient line with irrigation fluid and close one clamp on the patient line.
3. Connect the patient line to the scope, Foley catheter, or similar device.

Model 247 Product Safety Features

The following chart describes the safety alarm features of the Ranger irrigation fluid warming unit.

Alert Type	Set-Point	Indication	Description
Primary over-temperature alert system	48°C	LED alternately displays "HI" and temperature, alert sounds, alert light flashes.	Primary Alert: Alert sounds, controller cuts power to the heaters until temperature falls below 48°C. When heater temperature drops below 48°C, audible and visual alerts stop.
Independent back-up over-temperature alert system (cut-off).	50°C	LED displays the temperature or has gone dark, alert sounds.	Primary failure. Temperature will rise to 50°C within a few seconds. The safety backup system is activated at 50°C and power is immediately shut off to the heating plates.
Under-temperature alert	33°C	LED displays "LO" and temperature, alert sounds, alert light flashes.	Heater temperature has fallen to 33°C. Alerts will stop when temperature rises above 33°C. Continue use of unit.

NOTE.

The difference in over-temperature alert system cut-off points does not indicate any difference in safety and is simply a result of tighter controller tolerances.

The Ranger irrigation fluid warming system is designed with very tight controls to regulate its temperature to the 41°C set-point. The microprocessor controller measures the temperature at the fluid path way inside the warming plate four times every second, making it very responsive to changes in flow rate.


Instructions For Use

Preparation and Setup of the Ranger Irrigation Fluid Warming Unit

1. Attach the Ranger irrigation fluid warming unit to the I.V. pole. Tighten the pole clamp securely.


CAUTION:

To prevent tipping, clamp the Ranger irrigation fluid warming unit to an I.V. pole with a minimum 14" (35.6 cm) radius wheelbase and at a height no higher than 44" (112 cm). Failure to do so may result in damage to the product (see Figure 1).

2. Slide the irrigation cassette into the slot in the warming unit. The cassette can only fit into the device one way (see Figure 2).

Figure 2

3. Prime the warming set. For more information about priming the set, refer to instructions provided with the warming sets.
4. Plug the power cord into an appropriate outlet. Turn the unit ON. In a few seconds the alpha-numeric display will illuminate. It takes less than 2 minutes to warm up to the 41°C set point temperature.
5. Begin flow. When flow is completed, remove the warming set and discard according to institutional protocol.

Removing the Irrigation Fluid Warming Set from the Ranger Irrigation Fluid Warming Unit

1. Close the inlet clamp proximal to the fluid warming cassette and open all clamps distal to the cassette.
2. Disconnect the warming set from the fluid source, if applicable.
3. Allow a small amount of fluid to drain from the end of the patient line to reduce pressure in the cassette.
4. Remove the warming set from the warming unit and discard according to institutional protocol.

Troubleshooting

CONDITION	CAUSE	SOLUTION
Nothing illuminates on the warming unit panel.	Unit is not turned on, plugged in, or power cord is not plugged into an appropriate outlet.	Turn unit on. Make sure the power cord is plugged into the power entry module of the warming unit. Make sure the warming unit is plugged into a properly grounded outlet.

	• Unit failure	• Call 3M customer service.
Alert indicator illuminates and alert sounds, alphanumeric display alternately flashes a temperature of 48°C or higher and the word "HI."	Temporary over-temperature condition because. • An extreme change in flow rates occurred (e.g., from 500 mL/min to stop flow). • Unit was turned on and reached set point temperature before cassette was inserted. • Fluids were prewarmed to above 42°C before being run through the warming unit.	• Open flow to reduce temperature. Alerts will stop when the temperature drops below 48°C. The unit is ready to use. • Alerts will stop when the temperature drops below 48°C. The unit is ready to use. • Turn off unit and unplug it. Discontinue infusion of fluids. Do not warm fluids before infusing them through the Ranger irrigation fluid warming unit.
Alert sounds, alpha-numeric display and alert indicator light go dark.	Primary controller failure. Unit will no longer operate.	Power to heating plates will shut off if temperature rises to 50°C. Turn unit off and unplug it. Discontinue use of unit. Discard disposable set. Alert will continue to sound if you do not unplug unit. Call 3M technical support.
Unit alerts soon after plugging it in (unit does not have to be turned on for this condition to occur). Heater temperature rises to 50°C and unit shuts down soon after plugging it in (unit does not have to be turned on for this condition to occur).	Test screw on bottom of unit is loose or missing.	Make sure test screw is completely tightened. If it is missing, turn unit off and unplug it. Call 3M technical support.
Alert sounds but warming unit has been turned off.	Independent backup safety system has been activated.	Unplug unit. Call 3M technical support.
Cannot remove cassette from unit.	• Cassette is too full, fluids are still being infused, or clamp is open proximal to the warming cassette. • Warming unit is below patient level, creating excessive back pressure.	• Make sure fluid is drained from warming cassette before sliding out the cassette, that fluids are no longer infusing, and that clamp is closed proximal to the warming cassette. • Raise unit above patient level
Alerts indicator light illuminates and alert sounds, alpha-numeric	Under temperature condition caused by very high flow using	Alerts should stop when temperature rises above 33°C. If

display alternately flashes a temperature of 33°C or below and the word "LO."	very cold fluid, or defective heater/relay.	alerts continue, turn off unit, unplug unit, and discontinue use. Call 3M technical support.
Alpha-numeric display reads "Er 4" or "Open."	Open wire on temperature sensor.	Do not use unit. Call 3M technical support.
Alpha-numeric display reads "Er 5" or "Open."	Electrical interference.	Remove the unit. Refer to biomedical technician or call 3M technical support.

CLEANING

Cleaning the Ranger Irrigation Fluid Warming Unit

Clean the Ranger irrigation fluid warming unit on an as-needed basis.

CAUTION

- To reduce the risks associated with cross-contamination perform the following:
 - The cleaning tool provides only superficial cleaning, it does not disinfect or sterilize the interior of the unit.

NOTICE



- To avoid device damage follow the instructions below.
 - Do not immerse the warming unit in cleaning or sterilizing solutions. The unit is not liquid proof.
 - Do not clean the warming unit with solvents. Damage to the case, label, and internal components may result.
 - Do not insert metallic instruments in the warming unit.
 - Do not use abrasive materials or solutions to clean the heater plates.
 - Do not allow spills to dry inside the unit, as this may make it more difficult to clean the unit.

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- Disconnect the Ranger irrigation fluid warming unit from the power source.
- Wipe the outside of the unit with warm, soapy water, nonabrasive cleaning solutions, dilute bleach, or cold sterilants. Do not use abrasive materials.
- Wipe with a dry, soft cloth.

TO CLEAN THE HEATING PLATES FOLLOW THE INSTRUCTIONS BELOW.

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METHOD OF CLEANING

- Unplug the warming unit.
- Unfold the cleaning tool. Wet the foam pads with a nonabrasive solution, such as Alconox™ brand detergent.
- Insert the tool from the back of the unit and pull the tool all the way out from the front.

4. Rinse the tool with water and repeat 3 times. Dispose cassette in accordance with the rules of the medical protocol.
5. Wipe off the unit to remove excess fluid.

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NOTE. You may use a nonmetal instrument, such as a cotton swab, to clean the upper channels. If you are unable to adequately clean the unit, call 3M technical service.

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Servicing

There are no user-serviceable parts in the Ranger irrigation fluid warming unit. All service must be performed by 3M Patient Warming or an authorized service technician.

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SPECIFICATIONS

Recommendations and statement of the equipment manufacturer: electromagnetic radiation		
Device, Model 247, is designed for use in electromagnetic environment described below. The buyer or user of the Device, Model 247, should ensure the device use in this environment.		
Radiation test	Compliance	Electromagnetic environment: recommendations
Radiation test CISPR 11	Group 1	Device, Model 247, uses radio frequency energy only to perform its internal functions. In this regard, the level of radio frequency radiation is extremely low and, most likely, will not affect the operation of nearby electronic equipment.
Radio radiation CISPR 11	Class B	
Emission of harmonic current components IEC 61000-3-2	Class A	
Voltage fluctuations and flicker IEC 61000-3-3	Requirements are met	Device, Model 247, is suitable for use in any premises, including residential and directly connected to low-voltage public utility networks, which provide electricity to residential buildings.

Recommendations and statement of the equipment manufacturer: electromagnetic stability			
Device, Model 247, is designed for use in electromagnetic environment described below. The buyer or user of the Device, Model 247, should ensure the device use in this environment.			
Electromagnetic stability test	Test level according to IEC 60601	Level of compliance	Electromagnetic environment: recommendations
Electrostatic discharges IEC 61000-4-2	contact ± 6 kV air ± 8 kV	contact ± 6 kV air ± 8 kV	Floors should be wooden, cement or covered with ceramic tiles. If the floors are covered with synthetic material, the relative humidity should be at least 30%.

Fast electrical transients or bursts of pulses IEC 61000-4-4	power lines ± 2 kV	power lines ± 2 kV	The quality of the electricity supply network should meet the standard requirements for commercial or medical institutions.
Dynamic changes in supply voltage IEC 61000-4-5	interphase voltage ± 1 kV phase voltage ± 2 kV	interphase voltage ± 1 kV phase voltage ± 2 kV	The quality of the electricity supply network should meet the standard requirements for commercial or medical institutions.
Voltage dips, short-term faults and voltage fluctuations IEC 61000-4-11	$< 5 \% U_T$ (dip $> 95 \%$ from U_T) for 0.5 period $40 \% U_T$ (dip 60% from U_T) for 6 periods $70 \% U_T$ (dip 30% from U_T) for 30 periods $< 5 \% U_T$ (dip $> 95 \%$ from U_T) for 5 seconds	$< 5 \% U_T$ (dip $> 95 \%$ from U_T) for 0.5 period $40 \% U_T$ (dip 60% from U_T) for 6 periods $70 \% U_T$ (dip 30% from U_T) for 30 periods $< 5 \% U_T$ (dip $> 95 \%$ from U_T) for 5 seconds	The quality of the electricity supply network should meet the standard requirements for commercial or medical institutions. If user of the Device, Model 245, should continuously operate the device during power failure, it is recommended to use an uninterruptible power supply or battery as a power source of the device.
Magnetic field of industrial frequency (50-60 Hz) IEC 61000-4-8	3 A/m	3 A/m	Characteristics of magnetic fields of industrial frequency should be at the level characteristic for the typical room in a typical commercial or medical institution.

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Recommendations and statement of the equipment manufacturer: electromagnetic stability			
Device, Model 247, is designed for use in electromagnetic environment described below. The buyer or user of the Device, Model 247, should ensure the device use in this environment.			
Electromagnetic stability test	Test level according to IEC 60601	Level of compliance	Electromagnetic environment: recommendations
Induced radio waves IEC 61000-4-6	3 V (mean square) 150 kHz to 80 MHz	3 V (mean square)	Portable and mobile radio frequency communications equipment should not be located near any part of the device, Model 247, including cables, and the distance between them should be as recommended based on calculations using transmitter wave frequency data. Recommended distance between the device and other equipment $d = 1,2 \sqrt{P}$
Radiated radio waves IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

			$d = 1,2 \sqrt{P}$ within the range from 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ within the range from 800 MHz to 2.5 GHz where P – nominal maximum output power of transmitter in watts (W) according to transmitter manufacturer's data, and d is recommended distance between the device and other equipment in meters (m). The field strength of stationary radio frequency transmitters, determined by examining electromagnetic environment ^a , should be below the allowable value of voltage for each frequency bands ^b . Interference can be available near equipment marked with the following symbol. 
NOTE 1.	If the frequency of transmitter is within the range from 80 MHz to 800 MHz, a higher frequency band shall be used.		
NOTE 2.	These recommendations can be applies not to all cases. The propagation of electromagnetic waves is influenced by the absorption and reflection of waves by structures, objects and people.		
^a	Intensity of the fields of stationary transmitters, for example, base stations for radiotelephones (mobile and wireless) and terrestrial mobile radio installations, amateur radio equipment, broadcast television sets and radio installations with amplitude or frequency modulation, cannot be calculated theoretically with a high degree of accuracy. To evaluate the electromagnetic environment in the presence of stationary radio-frequency transmitters, it is necessary to take into account the results of electromagnetic environment studies. If the measured field strength in the room where the device, Model 247, is used exceeds the above permissible value, the operation of the device, Model 247, should be monitored to ensure its normal operation. In case of any signs of deviations in the device operation, it is probably required to take additional measures, for example, to change the orientation or location of the device, Model 247.		
^b	Within the frequency range from 150 kHz to 80 MHz, the field strength should be less than 3 V/m.		

Recommended distances between the device, Model 247, and portable and mobile radio frequency communications equipment

The Device, Model 247, is designed for use in electromagnetic environment with controlled radio frequency interference. The buyer or user of the device, Model 247, can contribute to prevention of electromagnetic interference by maintaining minimum distance between the portable and mobile radio frequency communications equipment (transmitters) and the device, Model 247, in accordance with the **recommendations below and maximum output power of the communications equipment.**

Nominal maximum output power of transmitter, W	Distance depending on the frequency of transmitter, m		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$

		$d = 1,2 \sqrt{P}$	
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.30
10	3.69	3.69	7.37
100	11.67	11.67	23.30

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NOTE 1. If transmitter frequency is within the range from 80 MHz to 800 MHz, the recommended distance corresponds to the range of higher frequencies.

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Physical characteristics

Dimensions	4.5 inches in height x 7.5 inches in width x 10 inches in length 11 cm x 19 cm x 25 cm Accuracy $\pm 15\%$
Weight	7 pounds 7 ounces (3.4 kg) Accuracy $\pm 15\%$
Temperature characteristics	
Set point temperature	41 °C
Set over-temperature (initial)	48 °C
Set under-temperature	33 °C
Set over-temperature (secondary)	50 °C
Electrical characteristics	
Input voltage	100-120 or 220-240 VAC
Operating frequency	50/60 Hz
Maximum heating power	900 W
Fuse	2 x T10A-H (250 V) for operation under voltage conditions 100-120 VAC 2 x T6.3A-H (250 B) for operation under voltage conditions 220-240 VAC
Fuse type	Time delay, high breaking
Leakage current	Meets leakage current requirements in accordance with UL 60601-1 and IEC 601-1.
Environmental Condition	
Operating Temperature Range	15 °C to 40 °C (59 °F to 104 °F)
Storage and Transportation Temperature Range	-20 °C to 45 °C (-4 °F to 113 °F)
Operating Humidity	10 to 85% RH, noncondensing
Atmospheric Pressure Range	50 kPa to 106 kPa

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This medical product does not contain drug products or does not introduce medicinal substances.
Human blood derivatives: Medical product does not contain and does not deliver human blood derivatives.

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Data on software used

Software	Title	Version	Developer	Date of issue
Controller Model 247; 120V/240V	116E09	7E	West Instruments	17th April 2015

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Methods of Accessories disposal: Packaging/contents should be disposed in accordance with the local/regional/national/ international regulations.

Sets should be disposed by the local public law organizations in accordance with SanPiN 2.1.7.2790-10 "Sanitary and Epidemiological Requirements for Medical Waste Treatment" (waste Class B).

WARRANTY

3M Company warrants this product will be free from defects in material and manufacture.

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