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B K B.

Signature of document custodian (affiant)

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State of California

County of Orange


Subscribed and sworn to (or affirmed) before me on this 31st day of August
20 17 by Ashni De Alwis — proved

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Signature

(Seal)

Appendix to the Instruction for use

MindFrame Capture LP revascularization device

1. Name of the medical device	MindFrame Capture LP revascularization device, variants: - MindFrame Capture LP revascularization device, diameter 3 mm, length 20 mm; - MindFrame Capture LP revascularization device, diameter 3 mm, length 30 mm; - MindFrame Capture LP revascularization device, diameter 4 mm, length 20 mm; - MindFrame Capture LP revascularization device, diameter 4 mm, length 30 mm;																												
2.1 Manufacturer of the medical device	Manufacturer: Micro Therapeutics, Inc. DbA Ev3 Neurovascular 9775 Toledo Way, Irvine, CA 92618, USA +1-949-837-3700																												
2.2 Authorized representative of the manufacturer of the medical device	Authorized Representative of Manufacturer: "Medtronic" LLC, 123317, Russia, Moscow, Presnenskaya Naberezhnaya., 10 Tel: (495) 580-73-77, Fax: (495) 580-73-78 www.medtronic.ru																												
3. Indication of the medical device with potential customer	It is intended for temporary use in order to restore blood flow in the cerebral vascular system in patients with acute ischemic stroke.																												
4. Functional characteristics and indication of the medical device	MindFrame Capture LP Revascularization Device have the following main parameters and characteristics: <table border="1"> <thead> <tr> <th></th><th>300011 MindFrame Capture LP 3 x 20</th><th>300010 MindFrame Capture LP 3 x 30</th><th>300013 MindFrame Capture LP 4 x 20</th><th>300012 MindFrame Capture LP 4 x 30</th></tr> </thead> <tbody> <tr> <td>Total length (cm)</td><td>177.0 ± 5.1</td><td>178.0 ± 5.1</td><td>177.0 ± 5.1</td><td>178.0 ± 5.1</td></tr> <tr> <td>Weight (g)</td><td>1.14 ± 0.10</td><td>1.14 ± 0.10</td><td>1.15 ± 0.10</td><td>1.15 ± 0.10</td></tr> <tr> <td>Working diameter (mm)</td><td>3.0 ± 0.5</td><td>3.0 ± 0.5</td><td>4.0 ± 0.5</td><td>4.0 ± 0.5</td></tr> <tr> <td>Working length (cm)</td><td>15.0 ± 1.0</td><td>23.0 ± 2.0</td><td>15.0 ± 1.0</td><td>23.0 ± 2.0</td></tr> </tbody> </table>					300011 MindFrame Capture LP 3 x 20	300010 MindFrame Capture LP 3 x 30	300013 MindFrame Capture LP 4 x 20	300012 MindFrame Capture LP 4 x 30	Total length (cm)	177.0 ± 5.1	178.0 ± 5.1	177.0 ± 5.1	178.0 ± 5.1	Weight (g)	1.14 ± 0.10	1.14 ± 0.10	1.15 ± 0.10	1.15 ± 0.10	Working diameter (mm)	3.0 ± 0.5	3.0 ± 0.5	4.0 ± 0.5	4.0 ± 0.5	Working length (cm)	15.0 ± 1.0	23.0 ± 2.0	15.0 ± 1.0	23.0 ± 2.0
	300011 MindFrame Capture LP 3 x 20	300010 MindFrame Capture LP 3 x 30	300013 MindFrame Capture LP 4 x 20	300012 MindFrame Capture LP 4 x 30																									
Total length (cm)	177.0 ± 5.1	178.0 ± 5.1	177.0 ± 5.1	178.0 ± 5.1																									
Weight (g)	1.14 ± 0.10	1.14 ± 0.10	1.15 ± 0.10	1.15 ± 0.10																									
Working diameter (mm)	3.0 ± 0.5	3.0 ± 0.5	4.0 ± 0.5	4.0 ± 0.5																									
Working length (cm)	15.0 ± 1.0	23.0 ± 2.0	15.0 ± 1.0	23.0 ± 2.0																									

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	Labeled Device Length (mm)	20± 2.0	30 ± 2.0	20 ± 2.0	30 ± 2.0
	Pushwire diameter (mm)	0.33± 5%	0.33± 5%	0.33± 5%	0.33± 5%
	Pushwire length (cm)	175± 1	175± 1	175± 1	175± 1
	Introducer				
	Overall length (cm)	40.0±1,0	40.0±1.0	40.0±1.0	40.0±1.0
	Outer diameter (mm)	< 0.42	< 0,42	< 0,42	< 0.42
	Weight (g)	0.25±0.05	0.25±0.05	0.25±0.05	0.25±0.05
5. Risks of using the medical device, contraindications, expected adverse effects associated with the use of the medical device for intended purpose	In accordance with BS EN ISO 14971:2012 and Covidien's internal Quality procedures, a risk assessment of the MindFrame Capture LP Revascularization Device was performed taking into consideration the possible hazards and associated risks related to the product design and clinical usage of the device. The results of the overall risk analysis demonstrate that the associated risks of using the MindFrame Capture LP Revascularization Device are adequately mitigated through the device design and Instructions for Use.				
6. Technical characteristics of the medical device	MindFrame Capture LP Revascularization Device have the following main parameters and characteristics:				
		300011 MindFrame Capture LP 3 x 20	300010 MindFrame Capture LP 3 x 30	300013 MindFrame Capture LP 4 x 20	300012 MindFrame Capture LP 4 x 30
	Radial force (N/mm)	0.021 – 0.042	0.021 – 0.042	0.024 – 0.029	0.024 – 0.029
	Resheathing force (N)	<2	<2	<2	<2
	Minimum full system tensile Strength (N)	20.334	20.334	20.334	20.334
	Peak force for deployment (N)	<1.1	<1.1	<1.1	<1.1
	Distal Marker Bond Strength (N)	≥ 2.2	≥ 2.2	≥ 2.2	≥ 2.2

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	Expandable Tip Attachment Integrity (N)	≥ 5.2	$\geq 5,2$	$\geq 5,2$	$\geq 5,2$	
	Delivery Wire Tensile Strength (N)	≥ 5.2	≥ 5.2	≥ 5.2	≥ 5.2	
	Coating length (cm)	>127	>127	>127	>127	
	Distal RO markers					
	Length (mm)	0.508 ± 0.127	0.508 ± 0.127	0.508 ± 0.127	0.508 ± 0.127	
	Inner diameter (mm)	0.254 ± 0,025	0.254 ± 0.025	0.254 ± 0.025	0.254 ± 0.025	
	Outer diameter (mm)	0.152 ± 0.025	0.152 ± 0.025	0.152 ± 0.025	0.152 ± 0.025	
	Marker band					
	Length (mm)	0.762 ± 0.025	0.762 ± 0.025	0.762 ± 0.025	0.762 ± 0.025	
	inner diameter (mm)	0.152 ± 0.008	0.152 ± 0.008	0.152± 0,008	0.152 ± 0.008	
	outer diameter (mm)	0.279 ± 0.005	0.279 ± 0.005	0.279 ± 0,005	0.279 ± 0.005	
	Support coil					
	Length (cm)	35.2 ± 0.1	35,2 ± 0.1	35.2 ± 0.1	35.2 ± 0.1	
	inner diameter (mm)	0.18	0.18	0.18	0.18	
	outer diameter (mm)	0.28± 0.03	0.28± 0.03	0.028± 0.03	0.028± 0.03	
	Introducer					
	Inner diameter (mm)	> 0.35	> 0.35	> 0.35	> 0.35	
	Tensile strength (N)	≥ 2.2	≥ 2.2	≥ 2.2	≥ 2.2	
	7. Information on the list of the main characteristics for the operation (application) of the medical device, transportation and storage conditions	Operating conditionals:				
		The device is used in the internal environment of the body				
t°		32 to 42 C°				
Pressure		700-1060 gPa				
Humidity		100 %				
	Storage Conditions					
	Store with other sterile devices in normal operating room temperatures in a dry place and in a manner that protects the integrity of the package and the sterile barrier. Rotate inventory so that					

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device is used before the expiry date on the package.

Air temperature, °C	+10 - +35
Air humidity (relative), non-condensing, %	< 75
Atmosphere pressure, hPa	860 - 1060

Transport Conditions

Air temperature, °C	-35 +40
Air humidity (relative), non-condensing, %	30 - 95
Atmosphere pressure, hPa	700 - 1060

8. List of national standards applied by the manufacturer

International Standards:

№		Standard Title
1	EN ISO 13485:2012 +AC:2012	Medical Devices – Quality Management Systems- Requirements for Regulatory Purposes
2	EN ISO 14971: 2012	Medical Devices – Application of Risk Management to Medical Devices
3	EN 62366:2008	Medical Device – Application of Usability Engineering to Medical Devices
4	ISO 10555-1:2013	Intravascular Catheters – Sterile and Single-Use Catheters – Part 1: General Requirements
5	ISO 11607-1:2006	Packaging for Terminally Sterilized Medical Devices – Part 1: Requirements for Materials, Sterile Barrier Systems and Packaging Systems
6	ISO 10993-1:2009	Biological Evaluation of Medical Devices – Part 1: Evaluation and Testing .
7	ISO 10993-3:2003	Biological Evaluation of Medical Devices – Part 3: Tests for Genotoxicity, Carcinogenicity and Reproductive Toxicity
8	ISO 10993-4:2009	Biological Evaluation of Medical Devices – Part 4: Selection of Tests for Interactions with Blood
9	ISO 10993-5:2009	Biological Evaluation of Medical Devices – Part 5: Tests for In Vitro Cytotoxicity
10	ISO 10993-7:2008	Biological Evaluation of Medical Devices – Part 7: Ethylene Oxide Sterilization Residuals
11	ISO 10993-10:2013	Biological Evaluation of Medical Devices – Part 10: Tests for Irritation and Delayed-Type Hypersensitivity
12	ISO 10993-11:2006	Biological Evaluation of Medical Devices – Part 11: Tests for Systemic Sensitivity
13	ISO 14155: 2011	Clinical Investigation of Medical Devices for Human Subjects – Good Clinical Practices

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	14	EN 556-1: 2001	Sterilization of Medical Devices. Requirements for Medical Devices to be Designated “STERILE.” Requirements for Terminally Sterilized Medical Devices.	
	15	ISO 11737-1:2006	Sterilization of Medical Devices – Microbiological Methods – Part 1: Determination of a Populations of Microorganisms on Products	
	16	ISO 11737-2:2009	Sterilization of Medical Devices – Microbiological Methods – Part 2: Tests of Sterility Performed in the Definition, Validation, and Maintenance of a Sterilization Process	
	17	ISO 11137-1:2006 (2010)	Sterilization of Health Care Products – Radiation – Part 1: Requirements for Development, Validation and Routine Control of a Sterilization Process for Medical Devices	
	18	ISO 11137-2:2012	Sterilization of Health Care Products – Radiation – Part 2: Establishing the Sterilization Dose	
	19	ISO 14644-1:1999	Cleanrooms and Associated Controlled Environments – Part 1: Classification of Air Cleanliness	
	20	ISO 14644-2:2000	Cleanrooms and Associated Controlled Environments – Part 2: Specifications	
	21	EN 1041:2008	Information Supplied by the Manufacturer of Medical Devices	
	22	EN 980:2008	Symbols for Use in the Labeling of Medical Devices	
	23	EN ISO 15233-1:2012 (replaces EN 980)	Medical Devices – Symbols to be Used with Medical Device Labels, Labeling and Information to be Supplied – Part 1: General Requirements	
9. Information on the sterile condition of the medical device, method of its sterilization and the procedure for dealing with a sterile package violation (if the medical device is delivered sterile), or if the medical device is not sterile, indication of the need for sterilization before use	MindFrame Capture LP are sterilized according to the parameters summarized in the table below:			
Sterilization Method		Gamma radiation		
Sterilization Validation Method		per ISO 11137-1: 2006 and ISO 11137-2:2009		
Sterility Assurance Level		25kGy for a 10 ⁻⁶ SAL		
Test	Acceptance Criteria	Result	Pass/Fail	
Pyrogen (LAL)	≤ 0.06	< 0.02 EU/mL	Pass	
Bioburden	< 1000 CFU	< 1000 CFU	Pass	
Do not use if sterile package is damaged.				
10. Information on the procedure for the preparation of the medical device for its reuse, including cleaning,	Not applicable. Single-use only.			

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disinfection, packaging and, if necessary, method of re-sterilization	
11. Information necessary for the identification of medical devices in order to obtain a safe combination, and information on known restrictions on concomitant use of medical devices (for medical devices intended for use with other medical devices and/or supplies)	Capture LP Revascularisation devices are compatible with microcatheters whose ID is 0,042 cm. Microcatheters are not included in the product and are sold separately.
12. Information on precautions taken in case of:	1. The MindFrame Capture LP should only be used by physicians who have received appropriate training in interventional neuroradiology and treatment of acute ischemic stroke. 2. Select a device size and configuration to engage the clot and to maintain sufficient vessel coverage on each side of the embolus along the parent vessel. An incorrectly sized device may result in no blood flow restoration and/or embolus migration. 3. The MindFrame Capture LP family of products is designed for use in vessels $\geq 2\text{mm}$ and $\leq 3.5\text{mm}$ in diameter. Use of the device in vessel diameters outside the recommendation can produce excessive resistance forces on the vessel and device components. If excessive resistance is encountered during the use of the device or any of its components at any time during the procedure, discontinue use. Movement of the device against resistance may result in damage to the vessel or a device component. 4. The device is provided STERILE for single use only. Reusing the device could result in compromised device performance, cross-infection and other safety related hazards. 5. Store in a cool, dry place. 6. Do not re-sterilize. 7. Use the device prior to the 'Use Before' date printed on the package. 8. Carefully inspect the sterile package and device prior to use to verify that neither has been damaged during shipment. Do not use kinked or damaged components. 9. The MindFrame Capture LP has not been shown to be MRI compatible. 10. The device should not be removed or repositioned without recapturing within the introducer Microcatheter unless thrombectomy is being attempted. 11. Exercise caution when crossing the deployed device with adjunctive devices (e.g., micro catheter). 12. Tighten the Rotating Hemostasis Valves sufficiently to create an adequate hemostasis seal without crushing the introducer microcatheter and MindFrame Capture LP shaft. Inadequately tightening the Rotating Hemostasis Valves may lead to premature deployment of the device. 13. After deployment, the distal tip of the device may foreshorten. 14. Do not steam shape or use pre-shaped microcatheters for the 15. MindFrame Capture LP introducermicrocatheter because it could damage the device.
a) Failures of the medical device, malfunction or deviations in the functioning that	Do not use the product in case of damage or sterility impairment!

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may affect the safety of the medical device, including those determined by external signs;	
b) MRI-compatibility.	No information about 8. MRI-compatibility.
13. Information on the circumstances which should prompt the consumer to consult a health professional	A medical device used only by qualified medical personnel who have undergone specialized training and have an appropriate level of surgical experience.
14. Information on the initial release or the latest revision of the operational documentation;	Not applicable
15. Procedure and conditions for the disposal or destruction of a medical device.	<u>After use, this product may be a potential biohazard.</u> Dispose the medical device in accordance with the requirements of the national legislation on environment and health and/or hospital Unused, expired devices are classified as medical waste class «A» Catheters extracted from the patient's body are medical waste Class B, that is epidemiologically hazardous waste. Waste Class B subjects to mandatory disinfection (disinfection)/ disposal. The choice of the method of disinfection/ neutralization is determined by the capabilities an organization engaged in medical and/or pharmaceutical activity and performed in the development of the treatment scheme of medical waste.
16. List and description of medical device materials with direct or indirect contact with the patient's body (human body)	Working part: Nitinol; Distal markers: alloy Pt/W; The proximal markers: alloy Pt/Ir; Flux for nitinol: Indalloy; Acid flux: Kester; Solder: alloy Au/Sn; Support coil: stainless steel; Explorer-follower: nitinol; Coating Polyamidoamine rubber; Coil support: polyethylene terephthalate; The coating sheath: high density polyethylene (HDPE).
17. Shelf Life	2 years

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COPY CERTIFICATION BY DOCUMENT CUSTODIAN

I, Ashni De Alwis, of legal age, being first duly sworn, hereby
Name of document custodian

Swear (or affirm) that the attached reproduction of

Capture IFU
Description of documents

Is a true and exact copy of the correct and complete original document.

[Signature]
Signature of document custodian (affiant)

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.


State of California

County of Orange

Subscribed and sworn to (or affirmed) before me on this 14th day of July,
 20 17 by Ashni De Alwis, proved

to me on the basis of satisfactory evidence to be the person(s) who appeared before me.

[Signature]
 Signature (Seal)

INSTRUCTIONS FOR USE

● MindFrame Capture™ LP

Priya Nankar

70883-001 Rev. 10/14

TABLE OF CONTENTS

MindFrame Capture™ LP

English	3
Français	5
Deutsch	7
Italiano	9
Español	11
Svenska	13
Nederlands	15
Português	18
Suomi	20
Dansk	22
Ελληνικά	24
Česky	26
Magyar	28
Русский	31
Polski	33
Türkçe	36
Norsk	38
Slovenčina	40
Română	42
Български	44
한국어	46
العربية	50
Symbol Glossary	53

English

Instructions for Use

EN

MindFrame Capture™ LP

CAUTIONS

This device is restricted to use by or on the order of a physician.

As with any medical treatment, it is the responsibility of the surgeon/physician to use his or her judgment in utilizing the procedures best suited to the needs of the patient.

Only physicians trained in neuro-interventional procedures should utilize the MindFrame Capture LP.

The MindFrame Capture LP contains (1) MindFrame Capture LP.

INDICATION FOR USE

The MindFrame Capture LP is indicated for temporary use to restore blood flow in the cerebral vasculature of patients suffering from an acute ischemic stroke. The MindFrame Capture LP is positioned across the embolus or blood clot and is used to facilitate the restoration of blood flow and removal of the clot obstruction.

The MindFrame Capture LP is indicated for:

- Endovascular temporary use in patients with acute ischemic stroke
- Endovascular temporary use to restore blood flow in patients who are experiencing symptoms of an acute ischemic stroke caused by an embolus in a cerebral vessel

The MindFrame Capture LP is delivered to and positioned across the embolus and unsheathed creating channels within the clot for blood to flow. Should blood flow stop, the device can be recaptured and redeployed. If the clot does not resolve after three (3) successive captures and device deployments then the device and clot can be slowly withdrawn from the vessel into the guide catheter while aspirating with a syringe.

- Evaluate the risks associated with acute endovascular ischemic stroke therapy (see complications below) and the possible benefits of immediate flow restoration prior to use of the MindFrame Capture LP.

CONTRAINDICATIONS

- Delivery of pharmacological agents not routinely used to treat ischemic stroke.
- Patient presents with nickel allergy.
- Patients with suspected or known allergies to contrast media
- Pregnancy
- Glucose <50mg/dl
- Excessive vessel tortuosity that prevents the placement of the device
- Known hemorrhagic diathesis, coagulation factor deficiency or oral anticoagulant therapy with INR>3.0
- Patient received heparin within 48 hours with a PTT greater than 2 times the lab normal.
- Patient has baseline platelets <30,000
- Evidence of rapidly improving neurological signs of stroke
- Coma
- Pre-existing neurological or psychiatric disease
- Patient has severe sustained hypertension
- CT/MRI scan reveals significant mass effect with midline shift
- Patient's angiogram shows an arterial stenosis >50% proximal to the embolus.

POTENTIAL COMPLICATIONS

Possible complications of the use of the MindFrame Capture LP include but are not limited to:

- Death from perforation of vessel.
- Air embolism.
- Arterial perforation with guidewire.
- Subarachnoid/interventricular hemorrhage due to vessel perforation from guidewire placement or device microcatheter placement.
- Thromboembolic stroke
- Intracerebral hematoma
- Hemiparesis from thromboembolic stroke
- Intraparenchymal bleeding
- Retroperitoneal hematoma
- Confusion

Complications of routine endovascular revascularization include:

- Cerebral ischemia
- Coagulopathy
- Confusion
- Death
- Embolic stroke
- Hematoma, pain, and/or infection at access site
- Intracerebral/intracranial hemorrhage
- Post-procedure bleeding
- Pseudoaneurysm formation
- Renal failure
- Vessel perforation
- Vessel thrombosis

WARNINGS

- The MindFrame Capture LP should only be used by physicians who have received appropriate training in interventional neuroradiology and treatment of acute ischemic stroke.
- Select a device size and configuration to engage the dot and to maintain sufficient vessel coverage on each side of the embolus along the parent vessel. An incorrectly sized device may result in no blood flow restoration and/or embolus migration.
- The MindFrame Capture LP family of products, as noted in the Recommended Sizing Guideline Table on page 3, is designed for use in vessels $\geq 2\text{mm}$ and $\leq 3.5\text{mm}$ in diameter. Use of the device in vessel diameters outside the recommendation can produce excessive resistance forces on the vessel and device components. If excessive resistance is encountered during the use of the device or any of its components at any time during the procedure, discontinue use. Movement of the device against resistance may result in damage to the vessel or a device component.
- The device is provided STERILE for single use only. Reusing the device could result in compromised device performance, cross-infection and other safety related hazards.
- Store in a cool, dry place.
- Do not re-sterilize.
- Use the device prior to the 'Use Before' date printed on the package.
- Carefully inspect the sterile package and device prior to use to verify that neither has been damaged during shipment. Do not use kinked or damaged components.
- The MindFrame Capture LP has not been shown to be MRI compatible.
- The device should not be removed or repositioned without recapturing within the introducer Microcatheter unless thrombectomy is being attempted.
- Exercise caution when crossing the deployed device with adjunctive devices (e.g., micro catheter).
- Tighten the Rotating Hemostasis Valves sufficiently to create an adequate hemostasis seal without crushing the introducer microcatheter and MindFrame Capture LP shaft. Inadequately tightening the Rotating Hemostasis Valves may lead to premature deployment of the device.
- After deployment, the distal tip of the device may foreshorten.
- Do not steam shape or use pre-shaped microcatheters for the MindFrame Capture LP introducer microcatheter because it could damage the device.

Table 1. MindFrame Capture LP
Product Name Numbers and Recommended Sizing Guidelines

Product Number	300010	300011	300012	300013
Product Name	10 MindFrame Capture LP	10 MindFrame Capture LP Short	10 MindFrame Capture LP	10 MindFrame Capture LP Short
Labeled Device Diameter (mm)	3.0	3.0	4.0	4.0
Labeled Device Length (mm)	30	20	30	20
Self Expanded Device Diameter (mm)	3.0	3.0	4.0	4.0
Recommended Vessel Diameter (mm)	≤ 3.0 and ≥ 2.0	≤ 3.0 and ≥ 2.0	≤ 3.5 and ≥ 2.5	≤ 3.5 and ≥ 2.5
Pusher Length	175cm	175cm	175cm	175cm
Introducer Microcatheter Minimum Inner Diameter	.0165"	.0165"	.0165"	.0165"

PROCEDURE

1. The MindFrame Capture LP is delivered endovascularly under fluoroscopic guidance in a manner consistent with other neurovascular catheter-based devices.
2. Antiplatelet and anticoagulation regimen used for interventional intracranial procedures is recommended at the discretion of the treating physician.

MINDFRAME CAPTURE LP PREPARATION AND PROCEDURE

Preparation

1. Administer anti-coagulation and anti-platelet medications per standard institutional guidelines.
2. Aided by angiographic radiography, determine the location and size of the area to be revascularized.
3. Select a MindFrame Capture LP per Table 1.
4. To achieve optimal performance of the MindFrame Capture LP and to reduce the risk of thromboembolic complications, maintain continuous flushing action between a) the femoral arterial sheath and the guide catheter, b) the microcatheter and the guide catheter and c) the microcatheter and the push wire and the MindFrame Capture LP. Check all connections to make sure that during the continuous flush that no air enters the guide catheter or the microcatheter.
5. Position a suitable guide catheter as close to thrombus site as possible employing a standard method. The guide catheter should be appropriately sized to retrieve clot if so desired in subsequent steps. Connect a RHV to the fitting of the guide catheter, and then connect a tube to the continuous flush.
6. With the aid of Table 1, select a microcatheter suitable for advancing the MindFrame Capture LP.
7. Connect a second RHV to the fitting of the microcatheter and then connect a tube to the continuous flush.
8. Set the flush rate per standard institutional guidelines.
9. With the aid of a suitable guide wire, advance the microcatheter until the end of the microcatheter is positioned distal to the thrombus so that the usable length portion of the MindFrame Capture LP will extend past each side of the thrombus in the vessel when fully deployed. Tighten the RHV around the microcatheter.

Delivering the MindFrame Capture LP

10. Insert the distal end of the introducer sheath partially into the RHV connected to the microcatheter. Tighten the RHV and verify that fluid exits the proximal end of the introducer sheath.
11. Loosen the RHV and advance the introducer sheath until it is firmly seated in the hub of the microcatheter. Tighten the RHV around the introducer sheath to prevent back flow of blood, but not so tight as to damage the MindFrame Capture LP during its introduction into the microcatheter. Confirm that there are no air bubbles trapped anywhere in the system.
12. Transfer the MindFrame Capture LP into the microcatheter by advancing the push wire in a smooth, continuous manner. Once the flexible portion of the push wire has entered the microcatheter shaft, loosen the RHV and remove the introducer sheath over the proximal end of the push wire. Once completed, tighten the RHV around the push wire. Leaving the introducer sheath in place will interrupt normal infusion of flushing solution and allow back flow of blood into the microcatheter.
13. Visually verify that the flushing solution is infusing normally. Once confirmed, loosen the RHV to advance the push wire.
14. With the aid of fluoroscopic monitoring, carefully advance the MindFrame Capture LP until its distal markers line up at the end of the microcatheter. The MindFrame Capture LP should be positioned such that the usable length portion of the device will extend past each side of the thrombus in the vessel when the device is fully deployed.

WARNING

- If excessive resistance is encountered during the delivery of the MindFrame Capture LP, discontinue the delivery and identify the cause of the resistance. Advancement of the MindFrame Capture LP against resistance may result in device damage and/or patient injury.

Deploying the MindFrame Capture LP

15. Loosen the RHV around the microcatheter. To deploy the MindFrame Capture LP, fix the push wire to maintain the position of the device while carefully withdrawing the microcatheter in the proximal direction.
16. Retract the microcatheter until it is just proximal to the proximal marker of the MindFrame Capture LP. Tighten the RHV to prevent any movement of the push wire. The usable length of the deployed device should extend past each side of the thrombus.
17. Tighten the RHV around the microcatheter. Angiographically assess the revascularization status of the treated vessel.

- DO NOT rotate the delivery pusher during or after delivery of the device into the parent vessel. Rotating the delivery pusher may result in a kinked device or fracture and detachment of the device from the delivery pusher.

Revascularization Device Recovery

18. If using balloon guide catheter, inflate guide catheter balloon to occlude vessel as specified in Balloon Guide Catheter labeling.
19. To retrieve thrombus, slowly withdraw the microcatheter and MindFrame Capture LP as a unit to the guide catheter tip while applying aspiration to the guide catheter with a 60cc syringe. Never advance the deployed MindFrame Capture LP distally. Note: Ensure microcatheter covers proximal marker
20. Apply vigorous aspiration to the guide catheter using syringe and recover MindFrame Capture LP and microcatheter inside guide catheter. Continue aspirating guide catheter until MindFrame Capture LP and microcatheter are nearly withdrawn from the guide catheter.

- If withdrawal into the guide catheter is difficult, deflate balloon (if using balloon guide catheter) and then simultaneously withdraw guide catheter, microcatheter and MindFrame Capture LP as a unit through the sheath while maintaining aspiration. Remove sheath if necessary.

- If excessive resistance is encountered during recovery of the MindFrame Capture LP, discontinue the recovery and identify the cause of the resistance. Do not perform more than Four recovery attempts in the same vessel using a MindFrame Capture LP.

21. Open the guide catheter RHV to allow the microcatheter and the MindFrame Capture LP to exit without resistance. Use care to avoid interaction with the site of the intervention and to prevent air from entering the system.
22. Aspirate the guide catheter to ensure the guide catheter is clear of any thrombus material.
23. Deflate guide catheter balloon if using balloon guide catheter.
24. If additional flow restoration attempts are desired with:
a new MindFrame Capture LP, then repeat the steps described above starting with the "Preparation" section.
the same MindFrame Capture LP, then:
 - a. Clean the device with saline solution. Do not use solvents or autoclave.
 - b. Carefully inspect the device for damage. If there is any damage, do not use the device and use a new MindFrame Capture LP for subsequent flow restoration attempts following the steps described above starting with the "Preparation" section. Use of a damaged device could result in additional device damage or patient injury.

WARNING

- Do not use each MindFrame Capture LP for more than four flow restoration recoveries.

MindFrame Capture LP Re-Sheathing

If resheathing of the MindFrame Capture LP is necessary (e.g. repositioning), follow these steps:

- Advancing the microcatheter while the device is engaged in dot may lead to embolization of debris. Do not advance the microcatheter against any resistance. Do not reposition more than three times.
1. Loosen the RHV around the microcatheter and around the push wire. With the aid of fluoroscopic monitoring, hold the push wire firmly in its position to prevent the MindFrame Capture LP from moving.
 2. Carefully re-sheath the MindFrame Capture LP by advancing the microcatheter over the MindFrame Capture LP until the distal markers of the MindFrame Capture LP line up at the end of the microcatheter. **If significant resistance is encountered during the re-sheathing process, stop immediately and proceed to the section above entitled "Revascularization Device Recovery".**

HOW SUPPLIED

Each **MindFrame Capture LP** contains one device positioned in an introducer sheath. All are supplied **STERILE (Gamma)** and **FOR SINGLE USE ONLY**. All components should be handled carefully to avoid damaging the device.

Sterile: This device is sterilized with Gamma Radiation.

Contents: One (1) MindFrame Capture™ LP Device.

STORAGE AND HANDLING

Handle with care. Packages should be stored in a manner that protects the integrity of the package, the sterile barrier and prevents exposure to temperatures in excess of 60 degrees C.

WARRANTY DISCLAIMER

ev3 Inc. warrants that reasonable care was used in the design and manufacture of this product. Because ev3 Inc. has no control over the conditions of use, patient selection or handling of the device after it leaves its possession, ev3 Inc. does not warrant either a good effect or against an ill effect following its use. ev3 Inc. shall not be directly or indirectly responsible for any incidental or consequential loss, damage or expenses directly or indirectly arising from the use of this product. ev3 Inc. sole responsibility in the event ev3 Inc. determines the product was defective when shipped by ev3 Inc., shall be the replacement of the product. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including but not limited to any implied warranties of merchantability or fitness for use.

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