
STATE OF MINNESOTA )  
 )ss.  
COUNTY OF HENNEPIN )

I, Deborah L Shaffer, Notary for the State of Minnesota,  
hereby declare this 11<sup>th</sup> day of October, 2016, that the  
attached IFU Fortrex 510331-001a is an Original  
Document.

Signature of Notary Public

# Fortrex™

## 0.035" OTW PTA Balloon Catheter

|    |                                                                                         |
|----|-----------------------------------------------------------------------------------------|
| en | 0.035" OTW PTA Balloon Catheter                                                         |
| fr | Cathéter à Ballonnet de 0,89 mm (0,035") OTW PTA                                        |
| de | 0,89 mm (0,035") OTW PTA Ballonkatheter                                                 |
| it | Catetere a palloncino OTW da 0,89 mm (0,035") per PTA                                   |
| es | Catéter con balón PTA OTW de 0,89 mm (0,035")                                           |
| sv | OTW (Över ledaren)-ballongkateter av PTA, 0,89 mm (0,035 tum)                           |
| nl | 0,89 mm (0,035") OTW PTA ballonkatheter                                                 |
| pt | Cateter de balão OTW PTA de 0,89 mm (0,035")                                            |
| fi | 0,89 mm:n (0,035 tuuman) OTW PTA -pallokatetri                                          |
| da | 0,89 mm (0,035 tommer) OTW PTA-ballonkateter                                            |
| el | Καθετήρας με μπαλόνι OTW PTA 0,89 mm (0,035")                                           |
| cs | Balónkový katétr OTW PTA 0,89 mm (0,035 palců)                                          |
| hu | 0,89 mm (0,035") OTW PTA ballonkatéter                                                  |
| ru | Катетер баллонный 0,89 мм (0,035") OTW PTA                                              |
| pl | Cewnik balonowy 0,89 mm (0,035 cala) OTW do PTA                                         |
| tr | 0,89 mm (0,035") OTW PTA Balon Kateteri                                                 |
| no | 0,89 mm (0,035 tommer) OTW PTA-ballongkateter                                           |
| sk | Balónikový katéter OTW PTA 0,89 mm (0,035 palca)                                        |
| ro | Cateter OTW PTA de 0,89 mm (0,035") cu balon                                            |
| bg | 0,89 mm (0,035") OTW балонен катетър за PTA                                             |
| zh | 0.89 mm ( 0.035" ) OTW PTA 球囊导管                                                         |
| ko | 0.89mm(0.035인치) OTW PTA 풍선 카테터                                                          |
| ar | قسطرة ببالون لعمليات رأب الأوعية الدموية بنظام التوجيه السلكي مقاس 0.89 سم (0.035 بوصة) |

## INSTRUCTIONS FOR USE

### DEVICE DESCRIPTION

Balloon catheters are used to exert radial force to dilate narrow vessel segments. The Fortrex™ 0.035" OTW PTA balloon catheter (Fortrex™ catheter) is an over-the-wire (OTW) 0.035" (0.89 mm) dual lumen catheter (1) with a distally mounted semi-compliant inflatable balloon (2) and an atraumatic tapered tip (3).

The manifold (4) includes a lumen marked "THRU" (5). This is the central lumen of the catheter, which terminates at the distal tip. This lumen is used to pass the catheter over a guidewire with a maximum diameter of 0.035" (0.89 mm). The lumen marked "BALLOON" (6) is the balloon inflation lumen, which is used to inflate and deflate the balloon with a mixture of contrast medium and saline solution. The balloon has two radiopaque markers (7) for positioning the balloon relative to the stenosis. The radiopaque marker bands indicate the dilating or working section of the balloon (8).


Figure 1. Fortrex 0.035" OTW PTA Balloon Catheter

The Fortrex 0.035" OTW PTA balloon catheter is available in multiple balloon sizes. Labeled balloon diameter and length are printed on the manifold (4). Refer to package label for information about catheter length and sheath compatibility.

Consult the balloon compliance chart packaged with the device for the diameters of the balloons at given pressures.

### HOW SUPPLIED

- The Fortrex™ catheter is supplied STERILE for single use only.
- Every Fortrex catheter is packaged with a protective sheath positioned over the balloon.

### STORAGE

The Fortrex catheter should be stored in cool, dry, and dark place. Do not store catheters where they are directly exposed to organic solvents, ionizing radiation, or ultraviolet light. Rotate inventory so that the device is used prior to the "Use by" date on the package label.

### INDICATIONS FOR USE

The Fortrex™ 0.035" OTW PTA balloon catheter is intended to dilate stenoses in the iliac, femoral, ilio-femoral, popliteal, infra-popliteal, and renal arteries, and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. This device is also indicated for stent post-dilatation in the peripheral vasculature.

### CONTRAINDICATIONS

The Fortrex catheter is contraindicated for use in coronary arteries or the neurovasculature, or when unable to cross the target lesion with a guidewire.

### WARNINGS

- This device is supplied STERILE for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilizing could increase the risk of patient infection and risk of compromised device performance.
- This device should only be used by physicians who are experienced and knowledgeable of the clinical and technical aspects of percutaneous transluminal angioplasty.
- To reduce the potential for vessel damage, the inflation diameter of the balloon should approximate the diameter of the lumen at the intended inflation site.
- Precautions to prevent or reduce clotting should be considered. Physician experience and discretion will determine the appropriate anticoagulation therapy for each patient.
- The Fortrex catheter should be used with caution for procedures involving calcified lesions or synthetic vascular grafts.
- Allergic reactions to contrast medium should be identified before PTA.
- When the catheter is exposed to the vascular system, it should be manipulated under high-quality fluoroscopic observation.
- Do not manipulate the balloon in inflated state.
- If resistance is encountered at any time during the insertion procedure, do not force passage. Resistance may cause damage to device or cause vessel damage. Withdraw the catheter carefully.
- Never use air or any gaseous medium to inflate the balloon.
- Do not retract the balloon catheter unless the balloon is free and fully deflated under vacuum, as vessel trauma or product damage may occur.

### PRECAUTIONS

- Do not use contrast media that is contraindicated for intravascular use with this device.
- Carefully inspect the package and catheter prior to use to verify no damage occurred during shipment.
- Do not attempt to pass the Fortrex catheter through a smaller sized sheath introducer than indicated on label (refer to product labeling).

• To avoid banking, advance the PTA catheter slowly, in small increments until the proximal end of the guidewire emerges from the catheter.

- The position of the Fortrex catheter may only be changed with the guidewire in place.
- Balloons should not be inflated in excess of the rated burst pressure.
- After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

### POTENTIAL COMPLICATIONS / ADVERSE EFFECTS

The following complications may result from a balloon dilatation procedure, but may not be limited to:

- |                                                                  |                                                             |
|------------------------------------------------------------------|-------------------------------------------------------------|
| • Arterio-venous fistula                                         | • Hypertension/Hypotension                                  |
| • Allergic reaction to device materials or procedure medications | • Infection                                                 |
| • Aneurysm                                                       | • Inflammation                                              |
| • Arrhythmia                                                     | • Intraluminal thrombus                                     |
| • Bleeding requiring transfusion                                 | • Myocardial infarction                                     |
| • Contrast medium reaction/renal failure                         | • Pseudoaneurysm                                            |
| • Death                                                          | • Sepsis                                                    |
| • Embolism                                                       | • Shock                                                     |
| • Endocarditis                                                   | • Stroke                                                    |
| • Fever                                                          | • Transient Ischemic Attack                                 |
| • Hematoma                                                       | • Venous Thromboembolism                                    |
| • Hemorrhage                                                     | • Vessel injury (e.g., dissection, perforation, or rupture) |

### SELECTION AND PREPARATION OF DEVICE

- The diameter of the expanded balloon should not exceed that of the vessel just distal or proximal to the stenosis.
- Verify that the selected accessories accommodate the Fortrex catheter as labeled.
- Prior to use, carefully inspect the package and catheter to verify no damage occurred during shipment.
- As supplied, the balloon and the guidewire lumen of the Fortrex catheter contain air. This air must be displaced to make certain that only liquid fills the balloon while the catheter is in the bloodstream.

### RECOMMENDED ITEMS

Prepare the following items using sterile technique:

- 10 cc syringe filled with sterile heparinized saline
- Contrast media - the standard inflation medium is a 1:1 mixture of contrast medium and normal saline

**CAUTION: Do not use contrast media that is contraindicated for intravascular use.**

- Appropriately sized exchange guidewire (refer to product labeling)
- Appropriately sized hemostatic introducer sheath (refer to product labeling)
- Inflation device with manometer

### PREPARATION OF PTA CATHETER

#### TO DISPLACE AIR:

1. Carefully remove the catheter from the inner package.

**CAUTION: Carefully inspect the package and catheter prior to use to verify no damage occurred during shipment.**

2. Remove protective tubing of the balloon. DISCARD.
3. Flush out the lumen for the guidewire with heparinized saline solution through the wire lumen marked "THRU".
4. Mix equal volumes of contrast medium and normal saline.
5. Fill a 10 cc syringe with approximately 4 cc of medium solution from step 4.
6. Evacuate air from the balloon segment:
  - a. Attach syringe to the balloon lumen, marked "BALLOON". Hold the syringe pointing down with the plunger up.
  - b. Apply negative pressure and aspirate for 15 seconds. Slowly release the pressure to neutral, allowing contrast to fill the shaft of the Fortrex catheter.
  - c. Disconnect the syringe from the BALLOON port of the dilatation catheter.
  - d. Remove all air from the syringe barrel. Reconnect the syringe to the BALLOON port of the Fortrex catheter. Maintain negative pressure on the balloon until air no longer returns to the device.
  - e. Slowly release the device pressure to neutral.
  - f. On first use, avoid any positive pressure that may unfold and partially inflate the balloon outside the vessel.
7. Repeat steps 6b-6e as necessary to remove all air from the system.
8. Replace the syringe with a manometer controlled dilatation system, without introducing air into the system.

CATHETER INSERTION AND DILATATION

1. The Fortrex catheter can be introduced percutaneously through an appropriate sized introducer sheath.
- WARNING:** When the catheter is exposed to the vascular system, it should be manipulated under high-quality fluoroscopic observation. Do not manipulate the balloon in inflated state.
- CAUTION:** Do not attempt to pass the Fortrex™ catheter through a smaller sized sheath introducer than indicated on the label (refer to product labeling).
- NOTE:** Maintain vacuum during balloon manipulation.
2. Place the prepared catheter over a pre-positioned guidewire, which has been placed through the lesion, and advance the catheter tip to the introduction site. A suitable length guidewire should be used at all times to maintain control and position of the guidewire.
- WARNING:** If resistance is encountered at any time during the insertion procedure, do not force passage. Resistance may cause damage to device or cause vessel trauma. Withdraw the catheter carefully.
- CAUTION:** To avoid kinking, advance the Fortrex catheter slowly, in small increments, until the proximal end of the guidewire emerges from the catheter.
- NOTE:** Use fluoroscopic guidance for all further catheter manipulations during the procedure.
3. Position the catheter with the center of the balloon in the stenotic region. The radiopaque marker bands indicate the working length of the balloon.
- CAUTION:** The position of the Fortrex catheter may only be changed with the guidewire in place.
4. When the marker bands have been appropriately positioned, inflate the balloon to dilate the target area.
- WARNING:** Do not use air or any gaseous substances as a balloon inflation medium.
- CAUTION:** Balloons should not be inflated in excess of the rated burst pressure.
- NOTE:** Repeat inflation of balloon (maximum 10 times) until desired result is achieved.
5. Deflate the balloon by aspirating the inflation device.
6. Maintaining a vacuum in the balloon, withdraw the catheter.
- WARNING:** Do not retract the Fortrex catheter unless the balloon is free and fully deflated under vacuum, as vessel trauma or product damage may occur.
- Multiple balloon inflations and deflations may cause resistance upon device withdrawal. Use a gentle counterclockwise twisting motion to ease withdrawal through the sheath. If the balloon cannot be withdrawn through the sheath, withdraw the catheter and sheath as one unit, while maintaining vessel access with the guidewire.
7. Results should be checked by angiography.
8. Repeat steps 1-7 as needed to achieve desired result.

DISPOSAL

**AUTION:** After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

WARRANTY DISCLAIMER

Although this product has been manufactured under carefully controlled conditions, Covidien has no control over the conditions under which this product is used. Covidien therefore disclaims all warranties, both express and implied, with respect to the product including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Covidien shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind Covidien to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

**CAUTION:** Federal law (USA) restricts this device to sale by or on the order of a physician.

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
STATE OF MINNESOTA )  
 )ss.  
COUNTY OF HENNEPIN )

I, Deborah L Shaffer, Notary for the State of Minnesota, hereby declare this 11<sup>th</sup>  
day of October, 2016, that the attached **Appendix TO IFU for Fortrex OTW**  
**Balloon Catheter**, is a true and complete copy of the original document.

  
Signature of Notary Public

**Appendix to the Instruction for use  
Fortrex OTW PTA Balloon Catheter**

nominal length of the catheter (sm): 40,80,135

nominal diameter of the catheter (mm): 4,5,6,7,8,9,10,12

working length of the catheter (sm): 20,40,80,100

**1.1. With regards to the manufacturer of the medical device - full name and, if available, abbreviated name, including company name, legal form of the legal entity, location address or last name, first name and (if available) patronymic, details of the identity document, address of residence of the individual entrepreneur, as well as phone numbers and (if available) e-mail address of the legal entity or individual entrepreneur;**

**Design facility:**

Covidien LLC  
4600 Nathan Lane North  
Plymouth, MN 55442 USA

**Manufacturer:**

Covidien LLC  
4600 Nathan Lane North  
Plymouth, MN 55442 USA

**Addresses of manufacturing facility of medical device:**

Covidien LLC  
4600 Nathan Lane North  
Plymouth, MN 55442 USA

**1.2. With respect to authorized representative of the manufacturer of the medical device - full and, if available, abbreviated name, including company name, legal form of the legal entity, address (location), telephone numbers and (if available) e-mail address of the legal entity**

**Authorized Representative of Manufacturer:**

"Medtronic" LLC, 123317, Russia, Moscow, Presnenskaya Naberezhnaya, 10  
Tel: (495) 580-73-77, Fax: (495) 580-73-78  
www.medtronic.ru

**Holder of Registration Certificate in Russia**

"Medtronic" LLC, 123317, Russia, Moscow, Presnenskaya Naberezhnaya, 10  
Tel: (495) 580-73-77, Fax: (495) 580-73-78  
www.medtronic.ru

**2. Technical characteristics of the medical device**

| Product Code<br>40 cm | Product Code<br>80 Cm | Product Code<br>135 cm | Ballo<br>on<br>Diam<br>eter<br>(mm) | Ball<br>oon<br>Len<br>gth<br>(mm<br>) | Nominal<br>Pressure     | Rated<br>Burst<br>Pressure | Recomm<br>ended<br>Introduc<br>er<br>Shealth |
|-----------------------|-----------------------|------------------------|-------------------------------------|---------------------------------------|-------------------------|----------------------------|----------------------------------------------|
| A35HPV04<br>020040    | A35HPV04<br>020080    | A35HPV04<br>020135     | 4 mm                                | 20<br>mm                              | 12<br>(atm)/121<br>6kPa | 24(atm)/24<br>32 kPa       | 6 (F)/2.0<br>mm                              |

Appendix to IFU are available upon request

Appendix to the Instruction for use  
Fortrex OTW PTA Balloon Catheter

|                    |                    |                    |      |           |                         |                      |                 |
|--------------------|--------------------|--------------------|------|-----------|-------------------------|----------------------|-----------------|
| A35HPV04<br>040040 | A35HPV04<br>040080 | A35HPV04<br>040135 | 4 mm | 40<br>mm  | 12<br>(atm)/121<br>6kPa | 24(atm)/24<br>32 kPa | 6 (F)/2.0<br>mm |
| A35HPV04<br>080040 | A35HPV04<br>080080 | A35HPV04<br>080135 | 4 mm | 80<br>mm  | 12<br>(atm)/121<br>6kPa | 24(atm)/24<br>32 kPa | 6 (F)/2.0<br>mm |
| A35HPV04<br>100040 | A35HPV04<br>100080 | A35HPV04<br>100135 | 4 mm | 100<br>mm | 12<br>(atm)/121<br>6kPa | 24(atm)/24<br>32 kPa | 6 (F)/2.0<br>mm |
| A35HPV05<br>020040 | A35HPV05<br>020080 | A35HPV05<br>020135 | 5 mm | 20<br>mm  | 12<br>(atm)/121<br>6kPa | 24(atm)/24<br>32 kPa | 6 (F)/2.0<br>mm |
| A35HPV05<br>040040 | A35HPV05<br>040080 | A35HPV05<br>040135 | 5 mm | 40<br>mm  | 12<br>(atm)/121<br>6kPa | 24(atm)/24<br>32 kPa | 6 (F)/2.0<br>mm |
| A35HPV05<br>080040 | A35HPV05<br>080080 | A35HPV05<br>080135 | 5 mm | 80<br>mm  | 12<br>(atm)/121<br>6kPa | 24(atm)/24<br>32 kPa | 6 (F)/2.0<br>mm |
| A35HPV05<br>100040 | A35HPV05<br>100080 | A35HPV05<br>100135 | 5 mm | 100<br>mm | 12<br>(atm)/121<br>6kPa | 24(atm)/24<br>32 kPa | 6 (F)/2.0<br>mm |
| A35HPV06<br>020040 | A35HPV06<br>020080 | A35HPV06<br>020135 | 6 mm | 20<br>mm  | 12<br>(atm)/121<br>6kPa | 24(atm)/24<br>32 kPa | 6 (F)/2.0<br>mm |
| A35HPV06<br>040040 | A35HPV06<br>040080 | A35HPV06<br>040135 | 6 mm | 40<br>mm  | 12<br>(atm)/121<br>6kPa | 24(atm)/24<br>32 kPa | 6 (F)/2.0<br>mm |
| A35HPV06<br>080040 | A35HPV06<br>080080 | A35HPV06<br>080135 | 6 mm | 80<br>mm  | 12<br>(atm)/12<br>16kPa | 23(atm)/23<br>30 kPa | 6 (F)/2.0<br>mm |
| A35HPV06<br>100040 | A35HPV06<br>100080 | A35HPV06<br>100135 | 6 mm | 100<br>mm | 12 (atm)/<br>1216kPa    | 23(atm)/23<br>30 kPa | 6 (F)/2.0<br>mm |
| A35HPV07<br>020040 | A35HPV07<br>020080 | A35HPV07<br>020135 | 7 mm | 20<br>mm  | 9(atm)/91<br>1.9 kPa    | 20(atm)/20<br>26 kPa | 6 (F)/2.0<br>mm |
| A35HPV07<br>040040 | A35HPV07<br>040080 | A35HPV07<br>040135 | 7 mm | 40<br>mm  | 9(atm)/91<br>1.9 kPa    | 20(atm)/20<br>26 kPa | 6 (F)/2.0<br>mm |
| A35HPV07<br>080040 | A35HPV07<br>080080 | A35HPV07<br>080135 | 7 mm | 80<br>mm  | 9(atm)/91<br>1.9 kPa    | 20(atm)/20<br>26 kPa | 6 (F)/2.0<br>mm |
| A35HPV07<br>100040 | A35HPV07<br>100080 | A35HPV07<br>100135 | 7 mm | 100<br>mm | 9(atm)/91<br>1.9 kPa    | 20(atm)/20<br>26 kPa | 6 (F)/2.0<br>mm |
| A35HPV08<br>040040 | A35HPV08<br>040080 | A35HPV08<br>040135 | 8 mm | 40<br>mm  | 9(atm)/91<br>1.9 kPa    | 20(atm)/20<br>26 kPa | 6 (F)/2.0<br>mm |
| A35HPV08<br>080040 | A35HPV08<br>080080 | A35HPV08<br>080135 | 8 mm | 80<br>mm  | 9(atm)/91<br>1.9 kPa    | 19(atm)/19<br>25kPa  | 6 (F)/2.0<br>mm |
| A35HPV08<br>100040 | A35HPV08<br>100080 | A35HPV08<br>100135 | 8 mm | 100<br>mm | 9(atm)/91<br>1.9 kPa    | 18/1824<br>kPa       | 6 (F)/2.0<br>mm |

Appendix to IFU are available upon request

Appendix to the Instruction for use  
Fortrex OTW PTA Balloon Catheter

|                    |                    |                    |          |          |                         |                |                 |
|--------------------|--------------------|--------------------|----------|----------|-------------------------|----------------|-----------------|
| A35HPV09<br>040040 | A35HPV09<br>040080 | A35HPV09<br>040135 | 9 mm     | 40<br>mm | 9(atm)/91<br>1.9 kPa    | 18/1824<br>kPa | 7(F)/ 2.3<br>mm |
| A35HPV09<br>080040 | A35HPV09<br>080080 | A35HPV09<br>080135 | 9 mm     | 80<br>mm | 9(atm)/91<br>1.9 kPa    | 18/1824<br>kPa | 7(F)/ 2.3       |
| A35HPV10<br>040040 | A35HPV10<br>040080 | A35HPV10<br>040135 | 10<br>mm | 40<br>mm | 8(atm)/<br>810.6<br>kPa | 16/1621<br>kPa | mm              |
| A35HPV10<br>080040 | A35HPV10<br>080080 | A35HPV10<br>080135 | 10<br>mm | 80<br>mm | 8(atm)/<br>810.6<br>kPa | 16/1621<br>kPa | 7(F)/ 2.3       |
| A35HPV12<br>040040 | A35HPV12<br>040080 | A35HPV12<br>040135 | 12<br>mm | 40<br>mm | 8(atm)/<br>810.6<br>kPa | 14/1419<br>kPa | mm              |
| A35HPV12<br>080040 | A35HPV12<br>080080 | A35HPV12<br>080135 | 12<br>mm | 80<br>mm | 7(atm)/70<br>9.3 kPa    | 12/1216kP<br>a | 7(F)/ 2.3       |

|                                    |                                                                                                                                                                                                                                                          |
|------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Delivery system                    | OTW                                                                                                                                                                                                                                                      |
| Usable catheter length             | 40, 80, 135(±1,5 mm)                                                                                                                                                                                                                                     |
| Balloon diameter                   | 4 - 12 mm (±0.3 mm)                                                                                                                                                                                                                                      |
| Balloon length                     | 20, 40 MM (±1.5 MM)<br>80,100 MM (±3 MM)                                                                                                                                                                                                                 |
| Coating                            | No coating                                                                                                                                                                                                                                               |
| Guidewire compatibility            | 0.035" (0.89 mm)                                                                                                                                                                                                                                         |
| Sheath compatibility               | 6-7Fr / 2.0-2.3 mm                                                                                                                                                                                                                                       |
| Nominal pressure                   | 7-12 atm /(709.3 – 1216 kPa)                                                                                                                                                                                                                             |
| Rated Burst Pressure               | 12, 14, 16, 18, 19, 20, 23, 24 atm /<br>1216, 1419, 1621, 1824, 1925,<br>2026, 2330, 2432 kPa                                                                                                                                                            |
| Balloon fatigue                    | Each balloon catheter tested must<br>not rupture when inflated to rated<br>burst pressure and deflated, 10<br>times                                                                                                                                      |
| Balloon Bond tensile Strength      | In accordance with ISO 10555-<br>should be >3H <ul style="list-style-type: none"> <li>• Distal balloon bond tensile<br/>force specification &gt;<br/>3.0lbf.</li> <li>• Proximal balloon bond<br/>tensile force specification<br/>&gt;3.0lbf.</li> </ul> |
| Balloon Inflation / Deflation Time | Balloon inflation and deflation<br>time < 60 seconds                                                                                                                                                                                                     |

Appendix to IFU are available upon request

Appendix to the Instruction for use  
Fortrex OTW PTA Balloon Catheter

|                                                  |                                                                                                                                                              |
|--------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The validity of damage to surface irregularities | No damage or negative impact on the product during use                                                                                                       |
| Freedom from leakage                             | No leakage in accordance with ISO 10555-1                                                                                                                    |
| Weight (with packaging)                          | 130±6 g                                                                                                                                                      |
| Weight (without packaging)                       | 5±0,25 g                                                                                                                                                     |
| Surface                                          | On the product surface should be no cracks, scratches, nicks, burrs, foreign inclusions and other defects.                                                   |
| Corrosion resistance                             | Catheter has corrosion resistance during use, transportation, storage in accordance with ISO 10555-1                                                         |
| Radiographic visibility                          | The balloon has two radiopaque markers                                                                                                                       |
| Double Radiopaque Markers                        | Provided on the 150 and 210 mm balloons (1.5 – 4 mm diameters) and 150 and 200 mm balloons (5 – 6 mm diameters) – marking the middle of the longest balloons |
| Diameter Radiopaque Markers                      | 1.3081mm ± 0.00762<br>0,0515''±0.0003                                                                                                                        |
| Tip Configuration                                | The tip of the distal end smooth, rounded, cone-shaped                                                                                                       |
| Measured Catheter Tip Length                     | The measured catheter Tip length is 5.0 +/- 1.0mm                                                                                                            |

**3. Information on the list of the main characteristics for the operation (application) of the medical device, transportation and storage conditions (for example, temperature and humidity, lighting and other characteristics)**

**Operation condition**

These products can be used in standard conditions in operating room. Other additional requirements for installation, installation or maintenance are not available. The product is intended for use by trained medical personnel in the operating room: at room temperature (18 ° C) not higher than 25 ° C, relative humidity of 75%, a pressure of from 86.0 to 106.0 kPa. During the procedure the product is inside of the human body and is operated at a body temperature from 32 to 42 °C

**Storage condition**

The Fortrex catheter should be stored in cool, dry and dark place. Do not store catheters where they are directly exposed to organic solvents, ionizing radiation, or ultraviolet light. Rotate inventory so that the device is used prior to the «Used by» date on the package label

Temperature: from 10 °C till +40 °C

Relative Humidity: from 30% till 75%

Pressure: from 585 till 1060 kPa

**Transportation condition:**

Care should be taken when transporting the devices. Packed by the manufacturer devices f are

Appendix to IFU are available upon request

transported by all modes of transport in closed vehicles in accordance with the rules of carriage of goods, operating in each mode of transport.

Temperature: from -18 °C till +55 °C

Relative Humidity: from 30% till 85%

Pressure: from 585 till 1060 kPa

#### 4. List of national standards applied by the manufacturer

| List of Standards                                     |                                                                                                                                                |
|-------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| Number                                                | Standard Title                                                                                                                                 |
| <b>Quality Systems</b>                                |                                                                                                                                                |
| EN ISO 13485: 2012                                    | Quality Management Systems - Requirements for regulatory purposes                                                                              |
| <b>Risk Management</b>                                |                                                                                                                                                |
| EN ISO 14971:2012                                     | Medical devices – Application of risk management to medical devices                                                                            |
| <b>Device Specific Standard</b>                       |                                                                                                                                                |
| ISO 10555-1:2014                                      | Sterile, single use intravascular catheters—Part 1: General requirements                                                                       |
| ISO 10555-4:2013                                      | Sterile, single use intravascular catheters —Part 4: Balloon dilatation Catheters                                                              |
| ISO 594-1:1986                                        | Conical fittings with a 6% (Luer) taper for syringes, needles and certain other medical equipment—Part 1: General requirements                 |
| ISO 594-2:1998                                        | Conical fittings with a 6% (Luer) taper for syringes, needles and certain other medical equipment—Part 2: Lock fittings                        |
| <b>Packaging, Shipping and Distribution Standards</b> |                                                                                                                                                |
| ISO 11607-1:2006/Amd-1:2014                           | Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems        |
| ISO 11607-2:2006/Amd-1:2014                           | Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes              |
| <b>Labeling Standards</b>                             |                                                                                                                                                |
| EN 1041:2008+A1:2013                                  | Information supplied by the manufacturer with Medical Devices                                                                                  |
| EN 980: 2008                                          | Symbols for use in the labelling of medical device                                                                                             |
| <b>Sterilization Standards</b>                        |                                                                                                                                                |
| ISO 11135-1:2007                                      | Medical devices - Validation and routine control of ethylene oxide sterilization.                                                              |
| ISO 11737-2:2009                                      | Sterilization of medical devices - microbiological methods - Part 2: Tests of sterility performed in the validation of a sterilization process |
| ISO 10993-7:2008/Cor 1:2009                           | Biological evaluation of medical devices – Part 7: Ethylene Oxide sterilization residuals                                                      |

Appendix to IFU are available upon request

Appendix to the Instruction for use  
Fortrex OTW PTA Balloon Catheter

|                                                                                                                                                                                                                                                                             |                                                                                                                                                              |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ISO 11138-1:2006                                                                                                                                                                                                                                                            | Sterilization of health care products --<br>Biological indicators --Part 1: General<br>Requirements                                                          |
| ISO 11138-2:2006                                                                                                                                                                                                                                                            | Sterilization of health care products --<br>Biological indicators --Part 2: Biological<br>indicators for ethylene oxide sterilization<br>processes           |
| <b>Clinical Standards</b>                                                                                                                                                                                                                                                   |                                                                                                                                                              |
| NA                                                                                                                                                                                                                                                                          | NA                                                                                                                                                           |
| <b>Environmental Control Standards</b>                                                                                                                                                                                                                                      |                                                                                                                                                              |
| ISO 14644-1:2015                                                                                                                                                                                                                                                            | Cleanrooms and associated controlled<br>environments. Part 1: Classification of air<br>cleanliness                                                           |
| ISO 14644-2:2015                                                                                                                                                                                                                                                            | Cleanrooms and associated controlled<br>environments. Part 2: Specifications for testing<br>and monitoring to prove continued compliance<br>with ISO 14644-1 |
| <b>Biocompatibility Standards</b>                                                                                                                                                                                                                                           |                                                                                                                                                              |
| ISO 10993-1:2009/C1 2010                                                                                                                                                                                                                                                    | Biological evaluation of medical devices – Part<br>1: Evaluation and testing within a risk<br>management process                                             |
| <b>Other</b>                                                                                                                                                                                                                                                                |                                                                                                                                                              |
| NA                                                                                                                                                                                                                                                                          | NA                                                                                                                                                           |
| <b>5. Warnings and/or precautions taken by the consumer when disposing of the medical device, accessories and consumables used together with it (if any), including information about the infectious, microbial, environmental or physical hazard of the medical device</b> |                                                                                                                                                              |
| <b>Class B</b>                                                                                                                                                                                                                                                              |                                                                                                                                                              |
| <b>6. Procedure and conditions for the disposal or destruction of a medical device.</b>                                                                                                                                                                                     |                                                                                                                                                              |
| After use, this product may be a potential biohazard (hazard class B). Handle and dispose of all such device in accordance with accepted medical practice and applicable local, regional and national laws and regulations                                                  |                                                                                                                                                              |

Name: Djuana Banks

Title: Senior Regulatory Affairs Specialist

Date: 03.21.2017

Djuana Banks

Date

Appendix to IFU are available upon request

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