INSTRUCTION for use

1. PURPOSE: Ophthalmic spectacle frames made of metal, plastic, combined produced by “op Couture Brillen GmbH”, meant for holding the lenses in a correct fixed position over the eyes.

2. MAIN TECHNICAL CHARACTERISTICS AND CONSTRUCTION

Frames, depending on design features and manufacturing technology, are divided into types indicated in Table 1.

Table 1.

Designation of frame type	Name of frame type	Type and designation of temples
FP	Plastic optical frames	Plastic rigid PR; Plastic elastic PE.
FM	Metal optical frames	Metal rigid PR; Metal elastic PE.
FC	Combined optical frames	Metal rigid MR; Metal elastic PE; Plastic rigid PR; Plastic elastic PE

Figure. 1. Metal optical frames (FM)

Figure. 2. Plastic optical frames (FP)

Figure. 3. Combined optical frames (FC)

3.3.2 Ophthalmic spectacle frames consist of a font, temples, nosepads, hinges. Hinges may be plane and spring.


Figure. 4. Structure of the frame

Ophthalmic spectacle frames are characterized by three basic geometric parameters: width (a), width of bridge (b) and length of temples before bending (L).

Applied to the frames of the company "Op Couture Brillen GmbH" dimensions are given in Table 2:

Types of frames for glasses correcting		The distance between the vertical tangents to the beveled facet groove (a). (Mm)	The minimum distance between the vertical tangents to the base of the facet groove of the rims (b). (Mm)	The whole temple length (L) (mm)	Temple length from the hinge axis to the beginning of the bend (L1) (mm)	Temple mass without lenses (g)
Ophthalmic spectacle frames of metal	men	$52 \pm 0,2$	$18 \pm 0,2$	$145 \pm 1,0$	$110 \pm 1,0$	$22 \pm 1,0$
	men	$56 \pm 0,2$	$20 \pm 0,2$	$145 \pm 1,0$	$110 \pm 1,0$	$22 \pm 1,0$
	women	$53 \pm 0,2$	$16 \pm 0,2$	$140 \pm 1,0$	$110 \pm 1,0$	$12 \pm 1,0$
	women	$54 \pm 0,2$	$17 \pm 0,2$	$135 \pm 1,0$	$95 \pm 1,0$	$20 \pm 1,0$
	unisex	$50 \pm 0,2$	$18 \pm 0,2$	$135 \pm 1,0$	$95 \pm 1,0$	$14 \pm 1,0$
	unisex	$56 \pm 0,2$	$20 \pm 0,2$	$150 \pm 1,0$	$110 \pm 1,0$	$18 \pm 1,0$
Ophthalmic spectacle frames of plastic	men	$47 \pm 0,2$	$18 \pm 0,2$	$140 \pm 1,0$	$100 \pm 1,0$	$26 \pm 1,0$
	men	$58 \pm 0,2$	$20 \pm 0,2$	$135 \pm 1,0$	$95 \pm 1,0$	$24 \pm 1,0$
	women	$53 \pm 0,2$	$16 \pm 0,2$	$140 \pm 1,0$	$115 \pm 1,0$	$24 \pm 1,0$
	women	$56 \pm 0,2$	$20 \pm 0,2$	$135 \pm 1,0$	$110 \pm 1,0$	$22 \pm 1,0$
	unisex	$47 \pm 0,2$	$18 \pm 0,2$	$140 \pm 1,0$	$100 \pm 1,0$	$20 \pm 1,0$

	unisex	55 ± 0,2	18 ± 0,2	135 ± 1,0	110 ± 1,0	24 ± 1,0
Ophthalmic spectacle frames of combined	men	58 ± 0,2	20 ± 0,2	140 ± 1,0	110 ± 1,0	24 ± 1,0
	men	56 ± 0,2	22 ± 0,2	135 ± 1,0	100 ± 1,0	20 ± 1,0
	women	54 ± 0,2	14 ± 0,2	140 ± 1,0	110 ± 1,0	22 ± 1,0
	women	53 ± 0,2	15 ± 0,2	130 ± 1,0	90 ± 1,0	20 ± 1,0
	unisex	54 ± 0,2	18 ± 0,2	135 ± 1,0	95 ± 1,0	24 ± 1,0
	unisex	50 ± 0,2	16 ± 0,2	140 ± 1,0	100 ± 1,0	22 ± 1,0

On the surface of the frames there should be no sharp edges and burrs visible to the naked eye of shells, scratches, bubbles, and build-up of solder.

The temples should rotate on the axis smoothly, without jamming.

Decorative finish (if any) of nickel silver (N18 CUNI18ZN20).

The frames of the corrective glasses correspond to the norms: ISO EN 12870, EN 12472, EN 16128.

The frames of the corrective glasses are selected taking into account the patient's anthropometric data.

3. INDICATIONS FOR APPLICATION:

Correction of defects in the optics of the eyes:

- myopia (myopia) - short-sightedness,
- Presbyopia (presbyopia) - hyperopia
- astigmatism (astigmatismus)
- Hypermetropia (hypermetropia), a weak type of clinical refraction
- anisometropia (difference in refractions on the eyes)
- keratoconus

The following may be CONTRAINDICATIONS TO APPLICATION:

- Individual intolerance to vision correction with glasses;
- the inability to correct vision with glasses.

4. MATERIALS:

Rims are made of materials that can not cause irritation, allergic or toxic reactions during their operation in contact with the skin of a person.

Plastic frames are made of cellulose acetate (cellulose acetate polymer (69-75%) and diethyl phthalate (25-31%)), hinges and screws made of nickel silver (N18CUNI18ZN20).

Metal and combined (including hinges and screws) from nickel silver (N18CUNI18ZN20) and cellulose acetate (cellulose acetate (69-75%) and diethyl phthalate

25- 31%) or pure titanium (CAS No. 7740-32- 6; EINECS 241-036-9).

The tips of the temples are made of cellulose acetate (cellulose acetate polymer (69-75%) and diethyl phthalate (25-31%)).

Nose pads are produced of silicone LSR50.

Transportation inserts made of acrylic Polycryl Clear.

Materials comply with the requirements of biological safety.

5. MARKING:

Right temple shall have factory marking with indication of:
model number (1097), technical size (54-16-135), colour or coating code (COL.002).

Figure 5. Frame marking. Right temple.

The left temple indicates the manufacture's name, conformity mark and country of manufacture.


Figure 6. Frame marking. Left temple.

Each frame is packed into an individual consumer plastic packaging.

Figure 7. Individual plastic packaging.

Plastic packaging of each frame or the label placed into plastic packaging should contain:

- nominal value of basic sizes;
- model number;
- model colour;
- bar code
- name of product;
- manufacturer;
- date of manufacture;
- expiration date,
- government standard


Figure 8. Label attached to plastic packaging.

Frames in plastic packaging should be compactly placed into group packaging box. Group packaging box should contain frames of the same type with the same basic sizes, in quantity of 8-10 pcs.

Figure 9. Group packaging box.

Each group packaging box should have a sticker of the manufacturer on it, which should contain following inscriptions:

- model number;
- model colour;
- nominal value of basic sizes;
- quantity of frames

Dimensions of the packaging box, in cm (length x width x height): 23 x 15 x 5,8.

Maximum load on the group packaging box is 10 kg.

Transport packaging is made only by agreement with the consumer.

Frames packed in transport packaging should be resistant to transport shaking with a maximum acceleration of 30 m/s² and frequency of 2 to 3 Hz.

Marking of transport packaging with application of the handling mark "Caution, fragile".

Set of delivery: transport inserts, frame.

6. PREPARING THE ITEM FOR USE: The spectacle frame should not lose its original shape after fitting the lenses. The rims should be symmetrical. No dents, scratches or gaps between spectacle lenses and rims, nor any other assembling defects are allowed.

Spectacle lenses fixed in the frame should not move or fall out of the rims during normal spectacle use.

7. TERMS OF STORAGE AND USE:

- It is not recommended to leave the frame in direct sunlight at high temperatures (for example, into car on the dashboard);

- Need to avoid the fall of the frame on hard surfaces, and mechanical loads that can damage the frame;

- For storage of frames need to use a case, is not permitted the wearing of frame into the pockets of clothing or bags in order to avoid scratches on the cover and mechanical damage;

- To put on and take off frames with both hands is necessary - this prevents their deformation;

- Need to avoid wearing frames during sleep;

- Are not allowed to apply to the frame significant mechanical forces;

- Caution should be exercised when doing sports and the production of works involving sudden movements of the head, the slopes, jumps and bumps.

- Do not apply significant mechanical efforts to the frame;

- Keep the frame in a case or any other package to protect it from mechanical damage;

- Full average term of use is no less than 3 years;

Terms of storage: in heated storages under the environmental temperature from + 5°C to + 50°C, relative air moisture 80%, Atmospheric pressure 700-1060 hPa

- No impact of aggressive substances (acids, alkali) and organic solvents and their fumes is allowed during storage.

- Storage life: plastic and combined frames – 5 years, metal frames – 10 years.

8. TERMS OF USE

Ophtalmic spectacle frames are used in all climate zones at any time of the day.

9. PRECAUTIONS FOR USE

While using ophtalmic spectacle frames it is necessary to exercise caution practising sports and activities involving sudden head movements, bends, jumps and bumps.

10. POSSIBLE SIDE EFFECTS

Allergic reaction in case of individual intolerance of frames material is possible.

11. INFORMATION ON THE POSSIBLE IMPACT OF THE USE ON ABILITY TO DRIVE VEHICLES, MACHINERY

Drivers with vision problems may drive vehicles and machinery using spectacle frames (in the manufacture of frames for corrective glasses).

12. STERILITY

Unsterile.

13. MAINTENANCE AND REPAIRS

Ophthalmic spectacle frames should be repaired in special laboratories licensed for production, maintenance and repairs.

14. WARRANTY OBLIGATIONS

Storage life: plastic and combined frames – 5 years, metal frames – 10 years.

Full storage life – not less than 3 years with temple bends not more than 30 000.

15. TRANSPORTATION

Ophthalmic spectacle frames are transported by all kinds of transport in covered vehicles in accordance with the rules of transportation of goods applicable to this kind of transport.

Temperature from - 45°C to + 50°C

Relative air moisture up to 98%

Atmospheric pressure 700-1060 hPa

16. DISPOSAL

Spectacle frames do not contain harmful substances and components that pose a threat to human health and the environment during and after the end of life and disposal. Ophthalmic spectacle frames should be disposed of as household waste.

Wastes of production and parts with expired year of storage are classified as class A medical waste.

17. CERTIFICATES

The manufacturer uses its own quality system that is subject to verification and updating on an annual basis.

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Legal address: 111024, Moscow, Andronovskoe highway, 26, build. 6, off. 11B


Phone number 8 (495) 660-01-25

Klaus Lippert

General Manager

“op Couture Brillen GmbH”

Signature

Date

10.11. 2017

B 1824/17

Protocol No.

B 1824/17

I, Gottfried Braun, notary public in Passau, hereby certify,
that the above is the true signature, subscribed in my presence, of

Mr. Klaus Lippert,

born on the 10th of March 1974,

address: Mühlbergstr. 16, 94535 Eging am See

who is personally known to me,

acting


on behalf of


op Couture Brillen GmbH

a company with registered office at Passau, Germany.

According to § 21 Section 1 No. 1 BNotO I certify by inspecting the electronic commercial register of the local court of Passau from today, that Mr. Klaus Lippert is sole authorized representative Managing Director of op Couture Brillen GmbH, with company's seat in Passau, registered in the Commercial register of the Local Court of Passau under No. HRB 5423 and he is released from the restrictions of § 181 BGB.

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