Notariat II Tuttlingen

Notar Haller

Bahnhofstraße 103 ♦ 78532 Tuttlingen
Tel.: 07461/98372 ♦ Fax: 07461/98358

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Important information for users of KARL STORZ instruments

1. Important information for users of KARL STORZ instruments

It is recommended to check the suitability of the product for the intended procedure prior to use.

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the cold light fountain power LED 300. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology may occur for this reason. Therefore, no claims may be derived from the information, illustrations, and descriptions in this manual.

© All product pictures, photos and descriptions and texts are the intellectual property of KARL STORZ SE & Co. KG

Utilization and reproduction by third parties require express authorization by KARL STORZ SE & Co. KG.

All rights reserved.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

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Controls, displays, connectors, and their uses

3. Controls, displays, connectors, and their uses

- ① Power switch
- ② TFT touch screen
- ③ Light outlet point (exchangeable)
- ④ SCB connectors*
- ⑤ Potential equalization connector
- ⑥ Power cord socket
- ⑦ Line fuses (fuse holder)

* The KARL STORZ-SCB interface (KARL STORZ Communication Bus), based on the CAN bus, permits remote control of equipment functions, as well as remote display of equipment parameters.

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4. Symbols employed

4.1 Symbols on the device

	Follow instructions for use
	ON
○	OFF
	Potential equalization connector
	Type CF applied part according to IEC 60601-1 (defibrillation-protected)*
	Alternating current
	Electronic information product pollution control (China RoHS)
	This device has been labelled in accordance with the European Directive on waste electrical and electronic equipment (WEEE).
	Manufacturer

* see Section 10.3 for a more detailed explanation

Symbols on label and packaging of the device

For the meanings of the symbols printed on the label or packaging, please refer to the 'Packaging symbols' accompanying instruction leaflet, mat. no. 96216316 DF. can be downloaded from www.karlstorz.com.

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	Menu ①
	Standby ⑦
	Start/Stop (On) ⑦
	Minus ②
	Plus ④
Auto	Automatic light regulation (‘blue’ = off) ⑤
Manual	Manual light regulation (‘orange’ = on) ⑤
	Reject
	Accept
	Remove
	Back
	Continue
	Audio signals off ⑧

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* Depending on the software version of the
IMAGE1 S™, the images may differ slightly from
actual appearance.

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IMAGE1 S™, the images may differ slightly from the actual appearance.

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1. Important information for users of KARL STORZ instruments.....	III
2. Images of the device	IV
3. Controls, displays, connectors, and their uses	V
4. Symbols employed.....	VI
4.1 Symbols on the device.....	VI
4.2 User interface symbols.....	VII
5. General information.....	3
5.1 Device description.....	3
5.2 Property rights.....	3
6. Safety instructions	4
6.1 Explanation of warnings and cautions	4
6.2 Normal use.....	8
6.2.1 Indication	8
6.2.2 Contraindication	8
6.3 Intended use	8
6.4 User qualification.....	8
6.5 Patient profile	9
6.6 User profile.....	9
6.7 Intended conditions of	9
6.7.1 Use	9
6.7.2 Other intended conditions	10
6.8 User position.....	10
6.9 Safety precautions at the site of installation.....	10
6.10 Safety features	11
6.10.1 Self-test.....	11
7. Installation and operating instructions.....	12
7.1 Unpacking the equipment	12
7.2 Basic equipment	12
7.3 Installation	12
7.3.1 Connecting the ground line	13
7.3.2 Connecting the power cord.....	13
7.3.3 KARL STORZ-SCB	13
7.3.4 Connecting light cable.....	14
7.3.5 Connecting the light cable to the endoscope or videoscope.....	14
7.3.6 TL 005 Triple adaptor (optional).....	15
7.4 Commissioning the power LED 300 light source	16
7.4.1 Using the automatic light source control.....	17

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7.4.2 Setting the brightness manually without an SCB connection to the camera control unit	18
7.4.3 Switching the light on and off	18
7.4.4 Decommissioning	18
7.5 Settings	19
7.5.1 Service	19
7.5.2 Language	19
7.5.3 Event protocol	20
7.5.4 Device Info	20
7.5.5 Audio settings	20
8. Maintenance	21
8.1 Fuse replacement	21
8.2 Reprocessing	22
8.2.1 General warnings for medical devices	22
8.2.2 Wipe-down disinfection of device	23
8.3 Maintenance and safety check	24
8.3.1 Maintenance	24
8.3.2 Safety checks/repeat inspections as per IEC 62353	24
8.4 Servicing and repair	25
8.5 Disposal	25
8.6 Repair program	26
8.7 Limitation of liability	27
8.8 Manufacturer's warranty	27
9. Technical description	28
9.1 Information signals	28
9.1.1 Configuration	28
9.1.2 Visual signals	28
9.1.3 Acoustic signals	29
9.1.4 Information signals (general)	29
9.1.5 Availability signal	30
9.2 Troubleshooting	31
9.3 Technical data	33
9.4 Technical documentation	34
10. Spare parts, recommended accessories	35
10.1 List of Spare parts	35
10.2 Recommended accessories	35
10.3 Recommended accessories (for CF as per IEC 60601-1)	35
11. Electromagnetic Compatibility (EMC) Information	37
12. Subsidiaries	54

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5. General information

5.1 Device description

The power LED 300 cold light fountain is designed to supply light for endoscopic diagnostic and surgical procedures.

The device offers a wide range of functions. These include:

- Efficiency: High light intensity at with low power consumption
- Comfortable working environment: Intelligent cooling concept for minimal noise emission
- Durable: No costs or effort required to change the lamp if used properly
- Clear: Direct operation via a color screen with touch-sensitive surface (touchscreen)
- Integration: The SCB function enables the light source to be controlled remotely and the system status information to be retrieved, while providing simple configuration during initialization
- Up-to-date: Software updates can be installed easily and conveniently via an interface

5.2 Property rights

This product is protected in the USA by (at least one of the following) US Patent No(s) 5,788,688; 6,397,286; 6,484,221; 6,824,539.

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6. Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

6.1 Explanation of warnings and cautions

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the product.

NOTE: A Note indicates special information about operating the product, or clarifies important information.

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Safety Instructions

Warnings and cautions

WARNING: The electrical installation of the operating room in which the device is connected and operated must comply with the applicable IEC standards.

WARNING: Keep device out of reach of patients.

WARNING: For safety reasons, do not simultaneously touch the device output sockets and the patient.

WARNING: The instruction manuals and interface specifications for medical devices and/or system components used in combination must be observed precisely.

WARNING: Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended use and interface specifications of the devices used in combination permit this (cf. IEC 60601-1-1 or Chap. 16, IEC 60601-1:2005).

WARNING: Grounding reliability can only be achieved when the unit is connected to a properly installed 'Hospital Only' or 'Hospital Grade' receptacle (i.e. approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

WARNING: To avoid risk of electric shock, this equipment must only be connected to supply mains with protective earth.

WARNING: At the end of its useful operating life, dispose of the unit as electronic scrap.

WARNING: Avoid allowing fluids to penetrate the device. Do not store liquids on or directly above the device.

If liquid penetrates the device, in spite of the precautions, sufficient time should be allowed for evaporation (this also applies for the formation of condensation).

WARNING: Ensure sufficient air circulation for the device.

WARNING: It is not permitted to press more than one touch point simultaneously.

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Safety instructions

Warnings and cautions

WARNING: Only operate the device with the voltage stated on the manufacturer's identification plate.

WARNING: DO NOT use the device in the presence of flammable anesthetics.

WARNING: Servicing may only be performed by specialists authorized by KARL STORZ.

WARNING: Test this device prior to each procedure to ensure that it is functioning correctly. It should not be used if any damage is evident.

WARNING: The device may only be operated with the power cable delivered by KARL STORZ or a similar power cable which has a national inspection seal.

WARNING: Do not use the device if the input controls or the displays are defective.

WARNING: The device must not be used in oxygenated environments.

WARNING: The device may be subject to sudden failures. Keep a spare device at hand or alternate procedures in mind.

WARNING: Applied parts of devices used in combination with this device must be of type BF or CF.

WARNING: Patient leakage currents may be additive when energized endoscopes are used with energized endotherapy devices. This is particularly important if a CF applied part endoscope is used, in which case a type CF applied part energized endotherapy device should be used in order to minimize total patient leakage current.

WARNING: The compatibility of this device with any accessories and/or energized endotherapy devices should be checked in the individual instructions for use.

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Safety instructions

Warnings and cautions

WARNING: The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

WARNING: The signal input and signal output parts of this device are designated by the manufacturer for exclusive connection to equipment which complies with IEC 60601-1.

WARNING: To avoid the risk of electrical shock, do not open the unit. Refer servicing to the manufacturer or to personnel authorized by the manufacturer. Removal of covers by unauthorized personnel will void the equipment's warranty.

WARNING: Always use the lowest possible settings for your intervention.

WARNING: Danger of glare! Never look into the light outlet when the system is switched on.

WARNING: Insufficient ventilation of the unit can cause an internal build-up of heat which leads to a safety shut-off of the light source.

WARNING: The correct configuration of the device must be checked after each software update.

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Safety instructions

Before using the device on the patient, it is imperative that you be acquainted with how the device operates and is controlled.

6.2 Normal use

It is only allowed to use the power LED 300 light source in combination with approved accessories from KARL STORZ or with equipment that is approved as safety related uncritical. Use of the device in fields other than those indicated above is not allowed for safety reasons.

WARNING: No modification of this equipment is allowed.

6.2.1 Indication

The device is intended to create light for endoscopic diagnostic observation and surgical procedures.

6.2.2 Contraindication

The medical device is contraindicated for ophthalmologic interventions. Mandatory protection filters for these kinds of interventions are not implemented for the medical device. For other interventions the responsible physician must decide whether the foreseen application is admissible based on the general condition of the patient.

6.3 Intended use

Cold light sources are designed to supply for endoscopic diagnostic and surgical procedures.

6.4 User qualification

Only qualified doctors and medical assistant staffs are allowed to use the power LED 300 light source after training. The training has to be conducted by skilled persons only. The area agency or the responsible KARL STORZ sales representative is available for trainings, for advices about further training alternatives.

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6.5 Patient profile

The use of the power LED 300 is not limited to a certain patient profile (sex, age, weight, etc.). The attending physician must assess a patient's health for the relevant application and declare it suitable.

6.6 User profile

- Adequate and recognized skills in the relevant application
- Adequate powers of comprehension to rationally assess the surgical situation in hand
- Adequate language skills in at least one of the languages used on the device and in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of contents of the instruction manual
- Have been thoroughly trained at least once in the use of the device
- No physical impairments that could diminish perception of activation alarm signals (visual and acoustic) and that cannot be compensated using auxiliary means

6.7 Intended conditions of use

6.7.1 Use

The device is intended for use in hospitals and doctors' offices. The technical data and ambient conditions are described in the instruction manual. It is the user's responsibility to make sure the device is safe and operates properly before using it. During treatment using the power LED 300 cold light fountain, the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital levels and the anesthetic.

WARNING: Any treatment may only be performed if visual observation of the action of the device is ensured.

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6.8 User position

For operation, the user should always standing in front of the device within a 'cone of vision' with an opening angle of $\pm 45^\circ$ at a distance of approx. 30–70 cm to the front panel.

6.9 Safety precautions at the site of installation

The device may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

It is not intended for use in hazardous zones.

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This also applies to highly combustible and potentially explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants. The device is equipped with a connector for attaching a ground line. It should be connected up in accordance with national regulations.

! WARNING: The device should be placed on a flat surface or a video cart. Do not block air inlets/outlets as overheating may result in fire hazard or device failure during procedure.

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6.10 Safety features

The cold light fountain power LED 300 is equipped with the following safety features:

6.10.1 Self-test

A self-test is performed each time the device is switched on. If any of the tests is negative, a corresponding error message is displayed.

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Functions and features

Once the device has been installed and switched on, the light source emits cold light from its light port. The intensity of the light can be regulated by pressing the + and – buttons. The current intensity is shown on the touchscreen. The device always starts in Standby mode when it is switched on. The most recently used light intensity is always saved.

NOTE: The light intensity can also be controlled by remote control or automatically, provided the technical requirements for this are satisfied.

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7. Installation and operating instructions

7.1 Unpacking the equipment

Carefully unpack the cold light fountain power LED 300 and remove the device and its accessories from their packaging. Check for missing items and evidence of shipping damage. File any complaints with the manufacturer or supplier immediately. If possible, retain the original packing materials for later use; these can come in handy if the device has to be transported

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
7.3 Installation

! NOTE: The cold light fountain power LED 300 including the accessories connected may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

! WARNING: The device is not intended for use in potentially explosive zones. Do not operate the device within the demarcated hazard zone while explosive anesthetic gases are in use.

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Set the device on a horizontal, flat surface (Fig. 1).

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Installation and operating instructions

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7.3.1 Connecting the ground line

The device is equipped with a connector ⑤ for attaching a ground line.
The device's ground line should be installed by a qualified electrician (Fig. 2).

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7.3.2 Connecting the power cord


Connect power cord and insert the power plug fully into the power cord receptacle ⑥ (Fig. 3).

- ⚠ WARNING:** Only operate the device with the voltage stated on the manufacturer's identification plate.
- ⚠ WARNING:** Only connect/disconnect the power plug to/from the power supply outside of areas subject to explosion hazards.
- ⓘ NOTE:** Ensure that access to the power plug is ensured at the site of installation.

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7.3.3 KARL STORZ-SCB

ⓘ NOTE: To prevent the SCB connecting cable from being pulled out accidentally, the SCB connector possesses a protection device. Therefore, to unplug the SCB cable, pull out the connector – do not pull directly on the cable!
Connect the SCB cable 20 0901 70 to an (equivalent) SCB socket on the back of the device (Fig. 4).
Connect the other end of the cable to the KARL STORZ-SCB (KARL STORZ Communication Bus) control device or other SCB devices (see KARL STORZ-SCB control NEO System Instruction Manual).

ⓘ NOTE: The SCB cable 20 0901 70 is an optional accessory (see section 10.2).

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
7.3.4 Connecting light cable

Plug the light cable into the socket ③ as far as it will go (Fig. 5). Only handle the light cable using the handle piece. Never pull on the cable.

NOTE: The outlet ③ is provided with an antiglare flap, which does not allow any direct outlet of light. Light is only emitted when the cable is inserted.

! WARNING: Risk of burning! The light source transmits high energy light which can lead to high surface temperatures (max. 52 °C on the light cable and max. 43 °C at the tip of the endoscope). Tissue damage is possible.

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7.3.5 Connecting the light cable to the endoscope or videoscope

Connect the light cable to the endoscope by a quarter turn of the knurled collar on the screw base (Fig. 6 or 7).

NOTE: The use of original KARL STORZ light cables only is recommended. Light cables from other manufacturers may not be optimally adjusted for light transmission.


NOTE: The degree of protection against electric shock at the applied part as type CF according to IEC 60601-1 (defibrillation-protected) can only be guaranteed by using the KARL STORZ light cables listed in Section 10.3.

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7.3.6 TL 005 Triple adaptor (optional)

The TL 005 triple adaptor (see 'Recommended accessories') can be used for the light cables from the following manufacturers:

- KARL STORZ
- Stryker (position II)
- Olympus (position III)

NOTE: Richard Wolf light cables can be used after lowering the metal sleeve on the connection piece at the KARL STORZ position.

! WARNING: Make sure that the light source is switched off when changing the light cable connection (adaptor).

1. Unscrew the light outlet point ③ (standard adaptor) by turning in a counter-clockwise direction (Fig. 8).
2. Place the triple adaptor onto the screw thread and turn to the left until the adaptor can be felt to click into place (Fig. 9).
3. Then screw in the adaptor hand-tight in a clockwise direction (Fig. 10).
- NOTE:** After being screwed in, it may be the case that the adaptor is positioned asymmetrically inside the device (marking not in 12 o'clock position – see arrow in Fig. 12). This asymmetry has no impact on the use of the adaptor, but should be corrected in order to prevent operating errors.
4. To do so, slightly undo the screw at the edge of the adaptor using a Phillips screwdriver (Fig. 11).
5. Rotate the adaptor (Fig. 12) such that the marking on the surface of the adaptor is in the 12 o'clock position (see arrow).
6. Re-tighten the screw (hand-tight) using the Phillips screwdriver (Fig. 12).
7. Select the required working position (KARL STORZ, II or III) by rotating the inner adaptor component (Fig. 13). This working position should always be in the 12 o'clock position (see note).
8. Insert the light cable into the corresponding socket, and connect (Fig. 14). See also section 7.3.4 and 7.3.5.

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
General applications

The power LED 300 is suitable for use in a wide range of applications. The correct light intensity depends on the particular application, the instrument being used, and the connected light cable. If the camera control unit used is not SCB-compatible, the intensity should always be kept to a minimum.

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**7.4 Commissioning the
power LED 300 light source**

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Once the power LED 300 has been switched on, the first thing that is displayed on the screen is the KARL STORZ logo (Fig. 16).

The device is now running a self test. Once the test has successfully completed, an acoustic availability signal sounds and the light source is ready for use. If the self test fails for critical components, the device goes into 'Malfunction' mode (Safe State) and an acoustic information signal sounds (Fig. 17). Verify that the device is properly connected and that the preparations have been made according to the instructions.


For further information on troubleshooting, refer to the section 'Technical description – Troubleshooting'.

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
The power LED 300 is in manual mode 'Light control' after being switched on, and the light is switched off (Standby mode 5 % light intensity) (Fig. 18).


NOTE: The icon for opening the Settings menu is only shown when the light is switched off (the symbol on the Start/Stop button is 'blue').

Press button 1 to switch to the Settings menu. Further details can be found in the section 'Settings'.

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7.4.1 Using the automatic light source control*

If the power LED 300 light source is connected to the IMAGE1 S™ camera control unit via an SCB connecting cable, automatic light source control can be used. Confirm the information message via the camera control unit. The 'AUTO' button lights up orange (Fig. 19).

*Prerequisite: IMAGE1 S™ with software release version 2.6

Press the Start/Stop button 7 to switch on the light (the symbol on the Start/Stop button 7 is orange) (Fig. 20). The light intensity is then controlled automatically via the IMAGE1 S™, ensuring optimal illumination in every situation.

NOTE: The automatic light regulation is deactivated as soon as the MANUELL 6 button is pressed during operation. The automatic light source control can be switched on again via the camera control unit. Further information can be found in the IMAGE1 S™ instruction manual.

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
* Depending on the software version of the IMAGE1 S™ the images may differ slightly from the actual appearance.

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7.4.2 Setting the brightness manually without an SCB connection to the camera control unit

The power LED 300 is in manual mode after being switched on, and the light is switched off (Standby mode).

The brightness can be modified using the buttons 2 and 4, both when the light is switched on and off (Fig. 21).

Press the button 4 to increase the brightness.

Press the button 2 to reduce the brightness.

The current brightness is shown in % on the display 3.

NOTE: The orange triangles indicate the set setpoint value in relation to the possible setting range.

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7.4.3 Switching the light on and off

The light is switched on and off using the Start/Stop button 7. The current status of the light (on/off) is indicated by the color of the icon on the Start/Stop button (orange = on, blue = off) (Fig. 22).

NOTE: By pressing the + and – buttons, the setpoint value of the light intensity on the device is set. The device switches the preselected setpoint value at the light outlet only after releasing the button. The intensity can be set between 5% and 100%.

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* Depending on the software version of the IMAGE1 S™, the images may differ slightly from the actual appearance.

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7.4.4 Decommissioning

For decommissioning, switch off the unit at the power switch.

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
7.5 Settings


Click on the menu icon ① to call up the Settings menu (Fig. 23).

① **NOTE:** The Settings menu can be called up when the light ⑦ is switched off (icon on the Start/Stop button ⑦ is blue).

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
In the Settings menu, you can select various setting and management fields

- Service (only authorized service technicians)
- Language
- Event protocol
- Device Info
- 'Audio' (on/off)

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7.5.1 Service


The 'Service' window is opened by pressing the 'Service' button in the Settings menu (Fig. 24).

① **NOTE:** This area is password protected (Fig. 25).

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
7.5.2 Language

The 'Language' window is opened by pressing the 'Language' button in the Settings menu (Fig. 24). The following screen languages are available: Danish, English, Finnish, French, German, Italian, Polish, Portuguese, Russian, Korean, Chinese, Swedish, Dutch and Norwegian (displayed in the appropriate language). Select a new language if required, and then confirm.

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7.5.3 Event protocol

The corresponding window is opened by pressing the 'Event protocol' button in the Settings menu (Fig. 24).

The Event protocol contains the saved information and alarms.

Each entry occupies one row and consists of a sequential number, date, time and alarm/ information ID.

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
7.5.4 Device Info


The 'Device Info' window is opened with device-specific data by pressing the 'Device Info' button in the Settings menu (Fig. 24).

In this window, you will find information about your device (e.g., the serial number, software version, and operating hours) (Fig. 28).

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7.5.5 Audio settings

The acoustic signals can be switched on or off (Audio on/Audio off) directly in the Settings menu (Fig. 29).

When the crossed-out bell is orange, the acoustic signals are switched off.

NOTE: When the acoustic signals are switched off, this is indicated by means of an 8 symbol in the top row of the main menu.

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8.2 Reprocessing

8.2.1 General warnings for medical devices

WARNING: Risk of infection: Incorrectly prepared medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments" and the accompanying documentation.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

WARNING: When preparing and using solutions, follow the chemical manufacturer's instructions paying close attention to proper concentration, exposure time, and service life. Incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

WARNING: Country-specific laws and regulations must be observed.

NOTE: The manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments" can be downloaded or requested by visiting www.karlstorz.com.

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CAUTION: Make absolutely sure that no liquid can get into the device housing.

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe ■ surface with a dry low lint cloth.

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8.3 Maintenance and safety check

WARNING: Always unplug the unit and disconnect the endoscope before carrying out any maintenance or cleaning.

8.3.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

8.3.2 Safety checks/repeat inspections as per IEC 62353

WARNING: Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.

Visual inspection

9. Check the device and accessories for any mechanical damage which may impair functionality.
10. Check that inscriptions relevant to safety are legible.

Electric measurements

- Inspect the device safety fuses
- Measure protective ground resistance in accordance with IEC 62353: The limit values can be found in the current standard.
- Measure earth leakage current in accordance with IEC 62353: The limit values can be found in the current standard.
- Touch current measured according to IEC 62353: The limit values can be found in the current standard.
- Patient leakage current measured according to IEC 62353: The limit values can be found in the

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Test for proper functioning

Carry out the functional check as per the instruction manual (see section 7.4).

Documentation

The safety check/repeat inspection and results must be documented.

WARNING: If any defects or shortcomings which could endanger patients, personnel or others are ascertained during these safety checks/repeat inspections, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

NOTE: For detailed information on the range and implementation of the safety check/repeat inspection, please refer to the latest issue of the relevant Service Manual.

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8.4 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts only.

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8.5 Disposal

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

WARNING: At the end of its useful operating life, dispose of the device as electronic scrap.

Please contact either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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8.6 Repair program

For videoscopes and devices, individual repair is necessary. Usually to bridge the repair period, you will receive a device on loan which you then return to KARL STORZ as soon as you receive the repaired device.

In Germany you can refer repairs directly to

KARL STORZ SE & Co. KG
Abt. Reparaturservice
Take-off Gewerbepark 83
78579 NEUHAUSEN

Servicehotline: +49 7461/708 980
E-Mail: technicalsupport@karlstorz.com

In other countries please contact your local KARL STORZ branch or authorized dealer.

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Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff).

We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of instruments or instruments which have been repaired by unauthorized third parties.

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8.7 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the device is connected and operated meet the applicable laws and standards; and
- the device has been used in accordance with its operating instructions at all times.

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8.8 Manufacturer's warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such repair or modification performed during the warranty period shall void all warranty.

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9.1 Information signals

9.1.1 Configuration

All of the information signals created, with the exception of those which display the cause of an inoperable device ('Malfunction', see e.g. page 16 test), are not self-retaining. This means that the signal is only issued whilst the signal conditions are present.

If the signal conditions are only very brief, there are rules as to the minimum visual/acoustic signal duration, with which the described information signals also comply. This is explained in more detail in the following chapters.

Signals are issued depending on priority. Any low priority signal will be overwritten by a signal with a higher priority, or a low priority information signal is suppressed for as long as there is a signal with a higher priority. If numerous signal conditions with the same priority occur, only the condition most recently detected is signaled.

9.1.2 Visual signals

Information signals are displayed inversely (blue font on a white background) compared to normal displays and do not flash.

The visual signal is issued for as long as the signal conditions apply. In the case of signal conditions of a very short duration, it may be difficult for the user to recognize which condition triggered the signal. For this reason, the minimum display duration is set at 5s.

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9.1.3 Acoustic signals


The acoustic signals described below satisfy the requirements in terms of distinctiveness.

The acoustic signal is issued for as long as the signal conditions apply. But at least one complete signal sequence is played.

The acoustic signal can be switched on or off in the Settings menu.

But the acoustic signal can only be switched off temporarily by the user. When the line voltage is switched on, the setting is always 'Audio On'.

9.1.4 Information signals (general)

The pitch is modulated by ± 2 Hz with a modulation frequency of 1.5 Hz. Five different harmonics are generated. In the case of messages which are related to the Safe State, the signal sequence is repeated every 15s, but otherwise it is just played once.

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

9.1.5 Availability signal


On completion of the self-test after switching on, an availability signal sounds. This signal, which is issued following a successful self-test, has the following characteristics:

The pitch is modulated by ± 2 Hz with a modulation frequency of 1.5 Hz. Five different harmonics are generated.

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- Complete failure of the device.

Possible causes:

- Power supply failure.
- Defective power fuse.
- Power cord connector is not properly connected to device socket.

Remedy:

- Check that there is electricity to the wall outlet.
- Change fuses as described in the instruction manual. Use only fuses of the correct rating.
- Push power cord connector firmly into device socket.

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No light is emitted.

Possible causes:

1. Electronics are defective.
2. Overheating of the power LED 300 (LED has switched off due to overheating of the unit caused by covered air vent).

Remedy:

1. Notify Service.
2. Remove objects from in front of the air vents. Switch off the power LED 300, let it cool down (approx. 10 – 15 minutes) and switch it back on. Ensure that it is sufficiently ventilated.

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Technical description

Troubleshooting

- No light output, power switch is lit (on).
- Strong impacts or vibrations during operation. Power supply or LED is defective.
- Turn the unit off and then on again.
- Have the unit repaired.

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- Too little light.
- Ends of light cable and/or Light cable or endoscope is LED service life has been exceeded.
- Clean ends of light cables, and light entrances and exits on the
- Replace light cable or
- Have the unit repaired.

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- Light intensity regulation shows 'errors'.
- Possible causes:
- Internal error in
- Remedy:
- Press the Standby button once, in order to restore 'normal'

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- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part type CF 

According to DIN EN 60825-1:
Class 1 laser

According to IEC 60601-1-2:
Observe the information on electromagnetic compatibility in chapter 11 (p. 36-49).
Directive compliance (for TL 300)
According to Medical Device Directive (MDD):
This medical device belongs to Class I
This medical device bears the CE mark in accordance with MDD 93/42/EEC.

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9.4 Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the device that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the device shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the device.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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Electromagnetic Compatibility (EMC) Information

11. Electromagnetic Compatibility (EMC) Information

WARNING: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this chapter during installation and operation.

The power LED 300 model TL300 corresponds to IEC 60601-1-2 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The power LED 300 model TL300 is a Group 1 device (as CISPR 11). Group 1 contains all the equipment and systems which generate or use RF energy only for their internal functioning.

NOTE: The tables and guidelines that are included in this chapter provide information to the customer or user that is essential in determining the suitability of the device or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the device or system to perform its intended use without causing interferences in other devices or systems or non-medical electrical devices. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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Electromagnetic Compatibility (EMC) Information

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING: Portable HF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the TL 300, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

WARNING: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

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Accessories and cables which have been shown to comply with IEC 60601-1-2:				
Type	Shielded	Length [m]	Ferrite	Used for
PA	No	> 3	No	Potential equalization
Power cord	No	3	No	Power supply

Guidance and manufacturer's declaration – electromagnetic emissions		
power LED 300 model TL300 is intended for use in the electromagnetic environment specified below. The customer or user of the power LED 300 model TL300 should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions acc. to CISPR 11	Group 1	power LED 300 model TL300 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions acc. to CISPR 11	Class B	The power LED 300 model TL300 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. – Professional healthcare facility environment –
Harmonic emissions acc. to IEC 61000-3 2	Class A	
Voltage fluctuations/flicker emissions acc. to IEC 61000 3-3	Complies	

Table 4 – ENCLOSURE PORT (according to IEC 60601-1-2)			
Phenomenon	Basic EMC standard or test method	IMMUNITY TEST LEVELS	
		Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT
ELECTROSTATIC DISCHARGE	IEC 61000-4-2	±8 kV contact ±2 kV, +4 kV, ±8 kV, ±15 kV air	
Radiated RF EM fields	IEC 61000-4-3	3 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz	Not applicable since released only for professional healthcare facility environment
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See table 9	
RATED power frequency magnetic fields	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	

Table 5 – Input a.c. power PORT (according to IEC 60601-1-2)			
Phenomenon	Basic EMC standard	IMMUNITY TEST LEVELS	
		Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT
Electrical fast transients/bursts	IEC 61000 4-4	±2 kV 100 kHz repetition frequency	
Surges Line-to-line	IEC 61000 4-5	±0.5 kV, ±1 kV	
Surges Line to-ground	IEC 61000-4-5	±0.5 kV, ±1 kV, ±2 kV	
Conducted disturbances induced by RF fields	IEC 61000 4-6	3 V 0.15 MHz - 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	Not applicable since released only for professional healthcare facility environment
Voltage dips	IEC 61000-4-11	0 % U_T , 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	
		0 % U_T , 1 cycle and 70 % U_T ; 25/30 cycles Single phase; at 0° Grad	
Voltage interruptions	IEC 61000-4-11	0 % U_T ; 250/300 cycles	

Table 8 – Signal input/output parts PORT (according to IEC 60601-1-2)

Phenomenon	Basic EMC standard	IMMUNITY TEST LEVELS	
		Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT
ELECTROSTATIC DISCHARGE	IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	
Electrical fast transients/bursts	IEC 61000-4-4	±1 kV 100 kHz repetition frequency	
Surges Line-to-ground	IEC 61000-4-5	±2 kV	
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 V 0.15 MHz - 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	Not applicable since released only for professional healthcare facility environment

Table 9 – Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment (according to IEC 60601-1-2)						
Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum Power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380 - 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 - 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
710	704 - 787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1 720	1 700 - 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1 845						
1 970						
2 450	2 400 - 2 570	Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
5 240	5 100 - 5 800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5 500						
5 785						

NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the power LED 300 model TL300 may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

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KARL STORZ Canada Ltd.
7171 Millcreek Drive,
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268 4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421 0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@kse.com

KARL STORZ Veterinary Endoscopy America, Inc.
1 South Los Cerreros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262 8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No.
Col. Granada Del. Miguel Hidalgo
C.P. 11520 Mexico D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing America Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar - Itam Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufnategui 627 6º Piso, B1638 CA
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 158
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsatragrand 14
12739 Skärholmens, Sweden
14108 Kungens Kurva, Sweden
Phone: +46 8 505 648 00
E-Mail: |

KARL STORZ Endoscopy Suomi OY
Tavastia 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania
Phone: +370 68567000
E-Mail: LT-LV info@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlystvej 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopes Nederland B.V.
Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N.V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopes France S.A.S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S.A.
Parque Empresarial San Fernando
Edificio Munch - Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S.r.l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Adria Eos d.o.o.
Capraska 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopja d.o.o.
Verovskova c. 60A
1000 Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

KARL STORZ Endoskop Magyarorszag Kft.
Toberek utca 2. sz. 17/b
HU-1112 Budapest, Hungary
Phone: +36 195 096 31, Fax: +36 195 096 31
E-Mail: info-hu@karlstorz.com

KARL STORZ Endoscopia Romania srl
Str. Prof. Dr. Anton Coloman, nr. 74, Sector 4
041393 Bucuresti, Romania
Phone: +40 (0)31 4250800, Fax: +40 (0)31 4250801
E-Mail: info-ro@karlstorz.com

KARL STORZ Endoskope Greece M.E.P.E.*
Patnariou Gongonou E' 34
54248 Thessaloniki, Greece
Phone: +30 2310 , Fax: +30 2310 304862
E-Mail: info-gr@karlstorz.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG & Service Subsidiary

KARL STORZ Industrial**
Gedik Is Merkezi B Blok
Kat 5, D 38-39, Bagdat Cad. No: 162
Maltepe Istanbul, Turkey
Phone: +90 216 442 9500, Fax: +90 216 442 9030

000 KARL STORZ Endoscopy WOSTOK
Derbenyevskaya nab. 7, building 4
115114 Moscow, Russia
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41
E-Mail: info-ru@karlstorz.com

TOV KARL STORZ Ukraine
Obolonska naberezhna, 15
building 3, office 3
04210 Kiev, Ukraine
Phone: +380 44 42668-14, -15, -19
Fax: +380 44 42668 20
E-Mail: marketing@karlstorz.com.ua

KARL STORZ SE & Co. KG
Representation Office
Sabit Orudschow 1184, apt 23
1025 Baku, Azerbaijan
Phone: +99 450 613 30 60
E-Mail: info-az@karlstorz.com

KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax: +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan
6, Saryarka str. BC "Arman", office 314
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A - Unit 7/A - 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai - United Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy Asia Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, Plaza Saigon
33 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Onon Road Lane Cove NSW 2066
P.O. Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongqiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 610044, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhui District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03 02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03 02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I Park Bundang1, 102dong, 8F
239 Jeongjail Ro, Bundang-Gu
Seongnam, 13556, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8677 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ SE & Co. KG
Representative Office: Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 31745-00
E-Mail: |

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102 0071,
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3J1, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
E-Mail: info@karlstorz.ca

KARL STORZ Endoscopy-America, Inc
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100,
Toll free: 800 421-0837 (USA only)
Fax: +1 424 218-8526
Toll free: 800 321-1304 (USA only)
E-Mail: info@ksca.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@kscla.com

KARL STORZ Endoscopia
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia Mexico S.A de C.V.
Lago Constanza No 326, Col Chapullitero Morales
D.F.C.P. 11520, Mexico, Mexico
Phone: +52 55 1101-1520, Fax: +52 55 55450174

KARL STORZ Marketing America do Sul Ltda
R. Joaquim Floriano, 413, 20º andar -
Itaim Bibi, São Paulo, -SP, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso
B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
P.O. Box 153, Rolf Olsenvei 28
N-2007, (illegible), Norway
Phone: (illegible)

KARL STORZ Endoskop Sverige AB
Storsatragrand 14
12724 Skarholmen, Sweden
Postal address: P.O. Box 8013,
14108 Kungälv Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Valimotie 1B, 4. Kerros, 00380 Helsinki, Finland
Phone: +358 96824774, Fax: +358 968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ Endoskopi Danmark A/S
Skovlystoftvej 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
392 Edinburgh Avenue, Slough
Berkshire, SL1 4UF, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: customerservice@karlstorz-uk.com

KARL STORZ Endoscopia Nederland B.V.
Phone: +31 651 938 738, Fax: +31 135 302 231

KARL STORZ Endoscopy Belgium N.V.
Phone: +32 473 810 451, E-Mail: info@stoplex.be

KARL STORZ Endoscopia France S.A.S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/I/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Iberica S.A.
Parque Empresarial San Fernando
Edificio Munich - Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: marketing@karlstorz.es

KARL STORZ Endoscopia Italia S.r.l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info@karlstorz.it

KARL STORZ Adria Eos d.o.o.
Zadarska 80
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

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KARL STORZ ENDOSKOPE East Mediterranean and Gulf
Building 25 (Al-Hayyan) – Units 205-206
Dubai Health Care City, P.O. Box 118069
Dubai – United Arab Emirates
Further - illegible

KARL STORZ Endoscopy India Private Ltd
D-181, Okhla Industrial Area, Phase-I
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10 – Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
PO Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch Company.
Room 610, China Life Tower No. 6, Chaowai Street,
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 2000300, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F, Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Republic of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch Company.
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 24 6258 9911, Fax: +86 24 6258 9922
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte. Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd
Jeljon Tower 1, Suite 905
16 Neuti-Ro, Bundang-Gu,
Seongnam, 463-847, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: Stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Beisin Rd
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mngwang@karlstorz.tw

KARL STORZ Endoscopy Japan K.K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 (0)3 6380-8633
E-Mail: info@karlstorz.co.jp

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