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
This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document and, where appropriate, the identity of the seal or stamp which the public document bears. This Apostille does not certify the content of the document for which it was issued. To verify the issuance of this Apostille, see www.austlii.edu.au/au/other/dfat/australia/apostille/

COOK MEDICAL EUROPE LTD.
O'Mahon Road,
National Technology Park
Limerick, Ireland.
Company Registration No. 401272
VAT Registration No. 9776234K

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
D – Width of disk stem – 0.080 inch/ 2.032 mm (± 0.010 inch/0.254 mm);

E – Height of stem 0.280 inch/7.112 mm (± 0.010 inch/0.254 mm).

Niall Sheehy
NIALL SHEEHY
Notary Public
6/7 Glentworth St., Limerick.


22-11-2017


Parameters for retention disk for: Percutaneous Nephrostomy Catheter (080208; 080210; 080212):

- A – 6 open holes, diameter – 0.045 inch/1.143 mm (± 0.010 inch/0.254 mm);
- B – Diameter of central opening – 0.200 inch/5.08 mm (± 0.005 inch/ 0.127 mm);
- C – Width of pull tie opening – 0.056 inch/1.422 mm (± 0.010 inch/0.254 mm);
- D – Length of pull tie opening – 0.125 inch/3.175 mm (± 0.016 inch/0.406mm).

Percutaneous Nephrostomy Catheter (080206; 080208; 080210; 080212)

The Percutaneous Nephrostomy Catheter is manufactured with a pigtail curve on the distal end with diameters of 6.0Fr/2mm, 8.3Fr/2.77mm, 10.0Fr/3.3mm, and 12.0Fr/4mm and in a length of 15 or 30 centimeters. The pigtail catheter is characterized by six side ports equally spaced on the inside of the pigtail curve, starting at a location 12 millimeters from the tapered distal tip. The circular pigtail at the distal end of the catheter has a diameter of 16 to 19 millimeters or 20 to 24 millimeters, depending on the catheter French size. Depending on the set, the pigtail catheter may be designed with a set of three ink-marks located near the distal end. The first ink mark is located 5 centimeters from the base of the pigtail curl, with the remaining ink marks spaced 5 centimeters apart. A female Luer-lock adapter and connector cap are attached to the proximal end of the device.


The Catheter may include a retention disc component. The retention disk device includes a black pull tie for securement of the catheter.

Product Catalogue Number	Catheter		Remarks
	Diameter Fr/mm	Length cm	
Tolerance	10%	10%	-
080206	6/2	15	Set consists of pigtail catheter and retention disc.
080208	8.3/2.77	30	Set consists of pigtail catheter and retention disc.
080210	10/3.33	30	Set consists of pigtail catheter and retention disc.
080212	12/4	30	Set consists of pigtail catheter and catheter support

Parameters of FLLA for: Percutaneous Nephrostomy Catheter (080206, 080208, 080210)

Diameters: D – 0.210 inch/5.334mm (± 0.002 inch/0.051mm); E – 0.344 inch/ 8.738mm (± 0.002 inch/0.051mm);

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
Length: D – 0.300 inch/ 7.62mm (± 0.005 inch/0.127mm); E – 0.281 inch/7.137mm (± 0.005 inch/0.127mm); F – 0.055 inch/1.397mm (± 0.005 inch/0.127mm); G – 0.290 inch/7.366 mm (± 0.005 inch/0.127mm); H – 0.725 inch/18.415 mm (± 0.005 inch/0.127mm); I – 0.185 inch/4.699mm (± 0.005 inch/0.127mm).

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Bond strength of components - not less than 15 N;

Tensile strength - not less than 8.1 MPa;

Relative elongation at break point - not less than 150%;

The maximum burst pressure - 81 atm. (1200 psi).

MATERIALS

Device Component	Material Name	Patient Contacting
Catheter	Catheter tubing: Polyurethane Radiopaque Tubing Blue #20393 Blue colorant #20393 Colorants: Titanium dioxide, Bismuth oxychloride CI Pigment Blue 29, CI Pigment Blue 15, Calcium carbonate Carbon black Adapter: Polyamide-6	Direct
Retention Disc only for 080206, 080208, 080210, 080212	Silicone # Q7-4750 Pull tie: Polyamide -6.6 #32131-17-2, Carbon Black # 1333-86-4	Direct

INTENDED USE/INDICATION FOR USE:

Percutaneous Nephrostomy Catheter (080206; 080208; 080210; 080212) - used for percutaneous placement of a pigtail catheter in the renal pelvis for nephrostomy drainage. Intended for one-time use.

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INSTRUCTIONS FOR USE

1. By preliminary plain film, I.V.P., retrograde pyelogram, ultrasound or CAT scan, localize the kidney to be drained.

A prone or slightly oblique position is preferred. Prepare and drape the flank in the usual fashion. Under fluoroscopic control, identify and anesthetize the skin site overlying the collecting system.

2. Introduce the 22 gauge (black hub) needle into the collecting system at a right angle to the skin. Obtain a urine specimen, if desired. Visualization of the pyelocalyceal system may be further enhanced by the injection of contrast medium.

NOTE: Maintain the 22 gauge needle at this location and position throughout the remainder of the procedure. This fixes the kidney in position and allows further injection of the pyelocalyceal system, if desired.

3. Rotate the patient into a more lateral position, with the posterior axillary line of the ipsilateral flank in the superior position. Pass the 18 gauge (pink hub) needle vertically through a small incision in the skin (made with a scalpel blade); insert the needle into the appropriate part of the pyelocalyceal system. Free flow of fluid from the needle after removing the obturator confirms a satisfactory location for the tip of the needle cannula.

4. Pass the flexible end of the Safe-T-J® wire guide through the hub of the 18 gauge needle cannula and advance it into the collecting system. Confirm the position fluoroscopically. Maintain the position of the wire guide as the needle cannula is withdrawn over it.

NOTE: Care must be exercised to avoid withdrawal of the wire guide itself.

5. To facilitate passage of the nephrostomy catheter, dilate the musculofascial tract by progressing from the smallest to the largest dilator provided, advancing the dilator over the wire guide.

6. Pass the pigtail end of the catheter over the external end of the wire guide; gradually advance the pigtail end well into the collecting system. Confirm the position fluoroscopically. Holding the shaft of the pigtail catheter securely in position with one hand, withdraw the wire guide with the other.

7. Tape or suture the retention disc to the skin; secure the catheter to the disc with the plastic tie. A dressing may be applied. Use the connecting tube to connect the pigtail catheter to a drainage bag or leg bag. The one-way stopcock is useful for testing the drainage capability of the collecting system, if desired.

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cutaneous
ephrostomy
Catheter
(080206)


COOK MEDICAL
Percutaneous Nephrostomy Catheter

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For US countries -- in accordance with the US Department of Labor of Occupational Safety and Health Administration (OSHA)

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
Niall Sheehy
Notary Public

Date

NIALL SHEEHY
Notary Public
6/7 Glenthworth St., Limerick.

APOSTILLE (Convention de La Haye du 5 octobre 1961)			
1. Country: Pays/Páis:		IRELAND	
This public document Le présent acte public / El presente documento público			
2. has been signed by a été signé par ha sido firmado por		Niall Sheehy	
3. acting in the capacity of agissant en qualité de quien actúa en calidad de		Notary Public	
4. bears the seal / stamp of est revêtu du sceau / timbre de y está revestido del sello / timbre de		-----	
Certified Attesté / Certificado			
5. at à / en	Cork	6. the le / el día	24/11/2017
7. by par / por	Department of Foreign Affairs and Trade		
8. No sous no bajo el número	9523922017		
9. Seal / stamp: Sceal / taimpe Sello / timbre	10. Signature: Signature: <i>Mary Kelleher</i> Firma:		

This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document and, where appropriate, the identity of the seal or stamp which the public document bears. This Apostille does not certify the content of the document for which it was issued. To verify the issuance of this Apostille, see www.austlii.edu.au/au/other/dfat/ie

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




A Cook-Cope Loop Nephrostomy Catheter includes a catheter, stiffening cannula and stylet. The Cook-Cope Loop Nephrostomy Catheter is manufactured in 8.2Fr/2.73mm, 10.0Fr/3.3mm, 12.0Fr/4mm and 14.0Fr/4.67mm and in a length of 19.5, 30 and 40 centimeters. The catheter is designed with a loop configuration at the distal tip. The distal portion of the loop is characterized by four side ports beginning 7 or 10 millimeters from the tip. A female Luer-lock adapter and connector cap are attached to the proximal end of the Cook-Cope loop catheter. An access hole is placed at the proximal end of the catheter shaft 2 centimeters from the connector cap for thread access for forming the distal curve. The stainless steel stiffening cannula is manufactured in an 18G/1.25mm size and a length of 35, 45.5 and 55.5 centimeters with a 0.040 inch/1.016mm stylet with a length of 37 cm for 085108, 085110, 085112, 085114, 085110-S8. The catheter introduction stiffening cannula is inserted into the catheter during introduction for assistance in catheter placement.

Product Catalogue Number	Catheter		Remarks
	Diameter Fr/mm	Length cm	
Tolerance	10%	10%	-
085108	8.2/2.73	19.5	Set consists of loop catheter and stiffening cannula.
085110	10.0/3.33	19.5	Set consists of loop catheter and stiffening cannula.
085110-S8	10.0/3.33	19.5	catheter only
085112	12.0/4	19.5	Set consists of loop catheter and stiffening cannula.
085114	14.0/ 4.67	19.5	Set consists of loop catheter and stiffening cannula.

Parameters of FLLA for: Cook-Cope Loop Nephrostomy Catheter (085108, 085110, 085110-S8)

Diameter: D – 0.210 inch/5.334mm (± 0.002 inch/0.051mm); E – 0.344 inch/ 8.738mm (± 0.002 inch/0.051mm);

Length: F – 0.172 inch/ 4.369mm (± 0.002 inch/0.051mm); G – 0.281 inch/7.137mm (± 0.002 inch/0.051mm); H – 0.312 inch/7.925mm (± 0.002 inch/0.051mm).

Product Catalogue Number	"A"	"B"	"C"
085108	0.090 inch/2.286mm ± 0.002 inch/0.051mm	0,775 inch/19.685mm ± 0.002 inch/0.051mm	0,062 inch/1.575mm ± 0.002 inch/0.051mm
085110	0.090 inch/2.286mm ± 0.002 inch/0.051mm	0,775 inch/19.685mm ± 0.002 inch/0.051mm	0,062 inch/1.575mm ± 0.002 inch/0.051mm
085110-S8	0.090 inch/2.286mm ± 0.002 inch/0.051mm	0,775 inch/19.685mm ± 0.002 inch/0.051mm	0,062 inch/1.575mm ± 0.002 inch/0.051mm

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Strength characteristics of the stylet:

Bond strength of components - not less than 10 N;

Strength characteristics of the cannula:

Bond strength of components - not less than 10 N;

MATERIALS

Cook-Cope Loop Nephrostomy Catheter (085108; 085110; 085112; 085114; 085110-S8)

Device Component	Material Name	Patient Contacting
Only for 085108		
Catheter	Polyurethane radiopaque White #20123 Colorant: Bismuth oxychloride Suture Holder: Latex MS#147 Tether Spool: Silicone # Q7-4765 Adapter: Polyamide-6	Direct
All except 085108		
Catheter	Polyurethane radiopaque Blue #20393 Colorants: Bismuth oxychloride, Titanium Dioxide, CI Pigment Blue 29, CI Pigment Blue 15, Calcium Carbonate, Carbon Black Suture Holder: Latex Natural Tubing #P0350 Tether Spool: Silicone # Q7-4765 Adapter: Polyamide-6	Direct
Cannula	Cannula: Stainless steel AISI 304 Cannula hub: Polypropylene #3 9003-07-0	Direct
Stylet	Cannula: Stainless steel AISI 304 Cannula hub: Polypropylene #3 9003-07-0	Direct

INTENDED USE/INDICATION FOR USE:

Cook-Cope Loop Nephrostomy Catheter (085108; 085110; 085112; 085114; 085110-S8) -

Used for percutaneous placement of a loop catheter in the renal pelvis for nephrostomy drainage.

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INSTRUCTIONS FOR USE

Catheter Placement

1. By preliminary plain film, I.V.P., retrograde pyelogram, ultrasound, or CT scan, localize the kidney to be drained. A prone or slightly oblique position is preferred. Under fluoroscopic control, identify and anesthetize the skin site overlying the collecting system.
2. Introduce the 21 gage needle below the twelfth rib, about five finger-breadths from the midline. Advance the needle in short steps toward the lower pole until the needle is seen moving with respiration. While the patient holds his or her breath, advance the needle no further than 3 cm. If urine is not aspirated, repeat the process using a slightly different angle until a calyx is successfully punctured.
3. Thread the .018 inch (0.46 mm) diameter wire guide through the needle; and maneuver the wire guide into the pelvis or proximal ureter.
4. Remove the 21 gage needle and advance the 6.3 French introducing a catheter with 20 gage cannula over the wire guide into the calyx. Holding the cannula stationary, further advance the catheter to follow the wire guide into the pelvis and ureter.

5. Remove the wire guide and cannula. Decompress and irrigate the collecting system.
6. Align the curved tip of the Safe-T-J wire guide with the curve of the introducing catheter. Introduce the flexible tip of the wire guide into the catheter and advance it until the wire guide exits the side port of the introducing catheter. If the advancing wire guide encounters resistance, withdraw it slightly and then advance it with a slight rotary motion.
7. When the tip of the Safe-T-J wire guide is within the renal pelvis, carefully withdraw the catheter. **NOTE:** Care must be exercised to avoid withdrawal of the wire guide.
8. To facilitate passage of the Cope nephrostomy catheter, dilate the musculofascial tract by progressing from the smallest to the largest dilator provided.
9. **IMPORTANT: To avoid breaking the catheter tether, unwind several centimeters of tether from the attached spool and gently straighten the catheter loop.** Pass the Cope nephrostomy catheter over the external end of the wire guide; gradually advance the catheter into the collecting system. If it becomes difficult to advance the catheter, insert the catheter introduction stiffening cannula over the wire guide and into the catheter. Some lubricant may be helpful.
10. Confirm catheter position fluoroscopically. While holding the shaft of the catheter securely with one hand, withdraw the wire guide (and stiffening cannula, if used) with the other.
11. While pulling the drawstring, maneuver the catheter back and forth under fluoroscopic monitoring until there is a looping of the distal tip. Wrap the drawstring around the catheter shaft; and cut the loose ends. Cover the knot by unrolling the rubber sleeve.
12. Attach the retention disc to the catheter shaft 1 to 2 cm from the skin surface. To prevent the catheter from being accidentally displaced forward, fix the retention disc in position by gluing or suturing.
13. The connecting tube provides a convenient, secure method of attachment to standard urine drainage pouches. The optional stopcock may be used for intermittent patency check.

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SUPPLY DETAILS

Cook-Cope Loop Nephrostomy Catheter is supplied as follows:

085110-S8: catheter – 1 p. in individual packaging; instructions for use – 1 p.

085108; 085110; 085112; 085114: catheter – 1 p., cannula with stylet – 1 p. in individual packaging; instructions for use – 1 p.

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

- company-manufacturer name and manufacturing address;
- address of the EU representative;
- marking "For specialists only";
- medical device name;
- part number;
- serial number (LOT);
- parameters and general device drawing;
- symbol about CE certification;
- symbol on sterilization;
- symbol on reusability;
- manufacturing date;
- symbol for keep away from sunlight;
- symbol for keep dry;
- bar-code;
- shelf life;
- symbol for "Consult Instruction for Use"
- symbol for caution "Attention, see Instructions for Use"

Cook-Cope
Loop
Nephrostomy
Catheter
(085108)

COOK MEDICAL
Cook-Cope Loop Nephrostomy Catheter

REF 085108
LOT 614494

DISPOSAL RULES

Upon completion of procedure, the devices are to be disposed of as per institutional guidelines for biohazardous medical waste.

For EU countries -in accordance with the guidance documents by World Health Organization guidance document.

For US countries – in accordance with the US Department of Labor of Occupational Safety and Health Administration (OSHA)

For Russian Federation – in accordance with national requirements. Risk class of the medical waste: class B.

Julie Padzorova

Regulatory Affairs Specialist

20 Nov 2017

Date

COOK MEDICAL EUROPE LTD.
O'Halloran Road,
National Technology Park
Limerick, Ireland.
Company Registration No. 492272
VAT Registration No. 9776234K

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
TO ALL TO WHOM THESE PRESENTS SHALL COME, I, Niall Sheehy, Notary Public, duly authorised, admitted, and sworn, practicing at 6/7 Glenthworth St., Limerick, Republic of Ireland, DO HEREBY CERTIFY that I am familiar with certificates as issued from time to time by the **Cook Medical Europe** and that the original certificate hereunto annexed namely, **Nephrostomy Catheters – Instructions for Use** upon which I have signed my name and affixed my Seal of Office is a true and faithful original appears to me to be genuine.

AND I HAVE CERTIFIED as herein at the request of Cook Ireland Limited having it's place of business at O'Halloran Road, National Technology Park, in the city of Limerick, Ireland, and for the benefit of the proper authorities in **Russian Federation** and for all to whom it shall come.

IN FAITH AND TESTIMONY whereof I have hereunto set my hand and affixed my Seal of Office.

Niall Sheehy
Notary Public

Date

6/7 Glenthworth St.
Limerick
6/7 St.

APOSTILLE (Convention de La Haye du 5 octobre 1961)			
1. Country: Pays/Páis:		IRELAND	
This public document Le présent acte public / El presente documento público			
2. has been signed by a été signé par ha sido firmado por		Niall Sheehy	
3. acting in the capacity of agissant en qualité de quien actúa en calidad de		Notary Public	
4. bears the seal / stamp of est revêtu du sceau / timbre de y está revestido del sello / timbre de		-----	
Certified Attesté / Certificado			
5. at à / en	Cork	6. the le / el día	24/11/2017
7. by par / por	Department of Foreign Affairs and Trade		
8. No sous no bajo el número	1434742017		
9. Seal / stamp: 	10. Signature: Signature:  Firma:		
This Apostille certifies the authenticity of the signature and the capacity of the person who has signed the public document. Where appropriate, the identity of the seal or stamp which the public document bears. This Apostille does not certify the content of the document for which it was issued. To verify the issuance of this Apostille, see www.austlii.edu.au/other/dfat			

COOK
MEDICAL

COOK MEDICAL EUROPE LTD.
O'Halloran Road,
National Technology Park
Limerick, Ireland.
Company Registration No. 482272
VAT Registration No. 0776234K

INSTRUCTIONS FOR USE

Nephrostomy Catheters

INTENDED USE

Used to conduct drainage in order to alleviate vesical obstruction, as well as in the range of procedures on the upper urinary tracts, including dilation of the strictures, removal of stones, or stenting of the ureter (Reference **Intended use/Indications for use** in instructions for specific intended uses).


CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Table of conversions

Parameter	Conversion
1 Fr	1/3 mm
1 Inch	25.4 mm

Malecot Catheter (082408; 0082410; 082414; 082416; 082418; 082422; 082424)

Niall Sherry
NIAL SHERRY
Notary Public
6/7 Glentworth St., Limerick.

22-11-2017


The malecot catheter is manufactured in 8.3Fr/2.77mm, 10.0Fr/3.3mm, 14.0Fr/4.67mm, 16.0Fr/5.33mm, 22.0Fr/7.33mm, and 24.0Fr/8mm and in lengths of 30 centimeters. A Luer lock adapter is attached to the proximal end of the catheter. The distal tip of the catheter is tapered 6 to 8 millimeters just proximal to the malecot. The malecot is formed by four slits in the catheter shaft placed 90° apart that facilitate drainage when the malecot is shaped.

Product Catalogue Number	Catheter		Remarks
	Diameter Fr/mm	Length cm	
Tolerance	10%	10%	-
082408	8.3/ 2.77	30	Malecot catheter only
082410	10.0/3.33	30	Malecot catheter only
082414	14.0/ 4.67	30	Malecot catheter only
082416	16.0/5.33	30	Malecot catheter only
082422	22.0/7/33	30	Malecot catheter only
082424	24.0/8	30	Malecot catheter only

Parameters of FLLA for: Malecot Catheter (082408, 082410)

Diameters: D – 0.210 inch/5.334mm (± 0.002 inch/0.051mm); E – 0.344 inch/ 8.738mm (± 0.002 inch/0.051mm);

Length: F – 0.172 inch/ 4.369mm (± 0.002 inch/0.051mm); G – 0.281inch/7.137mm (± 0.002 inch/0.051mm); H – 0.312 inch/7.925mm (± 0.002 inch/0.051mm).


Product Catalogue Number	“A”	“B”	“C”
082408	0.090 inch/2.286mm ± 0.002 inch/0.051mm	0.775 inch/19.685mm ± 0.002 inch/0.051mm	0.062 inch/1.575mm ± 0.002 inch/0.051mm
082410	0.090 inch/2.286mm ± 0.002 inch/0.051mm	0.775 inch/19.685mm ± 0.002 inch/0.051mm	0.062 inch/1.575mm ± 0.002 inch/0.051mm

Parameters of FLLA for Malecot Catheter (082410):

Diameter: F – 0.257 inch/6.528 mm (± 0.005 inch/0.127mm); G – 0.313 inch/7.950 mm (± 0.005 inch/0.127mm);


Length: H – 0.187 inch/4.749 mm (± 0.005 inch/0.127mm); K- 0.457 inch/11.608 mm (± 0.005 inch/0.127mm); L – 0.960 inch/24.384 mm (± 0.010 inch/0.254mm).

Product Catalogue Number	"A"	"B"	"C"	"D"	"E"
082410	0.125inch/3.175mm ± 0.002 inch/0.051mm	0.446 inch/11.328mm ± 0.005 inch/0.127mm	0.096 inch/2.438mm ± 0.005 inch/0.127mm	5°	0.490inch/12.446mm ± 0.005 inch/0.127mm

Parameters of FLLA for Malecot Catheter (082416, 082422, 082424)

Diameters: B – 0.339 inch/8.611 mm (± 0.005 inch/0.127mm); C – 0.300 inch/7.62 mm (± 0.005 inch/0.127mm); D – 0.312 inch/7.925 mm (± 0.005 inch/0.127mm); E – 0.402 inch/10.211 mm (± 0.005 inch/0.127mm); F – 0.550 inch/ 13.97 mm (0.015 inch/0.375 mm); Length: G – 1.010 inch/25.654 mm (± 0.010 inch/0.254mm); H – 0.300 inch/7.62 mm (± 0.005 inch/0.127mm); K - 0.304 inch/7.722 mm (± 0.005 inch/0.127mm); L – 0.406 inch/10.312 mm (± 0.005 inch/0.127mm); M - 0.375 inch/9.525 mm (± 0.005 inch/0.127mm); N – 0.125 inch/3.175 mm (± 0.005 inch/0.127mm).

Product Catalogue Number	"A"
082416	0.150 inch/3.81 mm (± 0.005 inch/0.127mm)
082422	0.250 inch/6.35 mm (± 0.005 inch/0.127mm)
082424	0.250 inch/6.35 mm (± 0.005 inch/0.127mm)

Strength characteristics of the catheter:

Bond strength of components - not less than 15 N;
Tensile strength - not less than 8.1 MPa;
Relative elongation at break point - not less than 150%;
The maximum burst pressure - 81 atm. (1200 psi).

Percutaneous Nephrostomy Malecot Catheter (082004)


The malecot catheter is manufactured in 14.0Fr/4.67mm French sizes and in lengths of 25 centimeters. The distal tip of the catheter is tapered 8 millimeters just proximal to the malecot. The malecot is formed by four slits in the catheter shaft placed 90° apart that facilitate drainage when the malecot is shaped

Product Catalogue Number	Catheter		Remarks
	Diameter Fr/mm	Length cm	
Tolerance	10%	10%	-
082004	14.0/ 4.67	25	Just Malecot

Parameters of FLLA for Percutaneous Nephrostomy Malecot Catheter (082004):

Diameter: F – 0.257 inch/6.528 mm (± 0.005 inch/0.127mm); G – 0.313 inch/7.950 mm (± 0.005 inch/0.127mm);
Length: H – 0.187 inch/4.749 mm (± 0.005 inch/0.127mm); K- 0.457 inch/11.608 mm (± 0.005 inch/0.127mm); L – 0.960 inch/24.384 mm (± 0.010 inch/0.254mm).

Order number	"A"	"B"	"C"	"D"	"E"
082004	0.125inch/3.175mm ± 0.002 inch/0.051mm	0.446inch/11.328mm ± 0.005 inch/0.127mm	0.096inch/2.438mm ± 0.005 inch/0.127mm	5°	0.490inch/12.446mm ± 0.005 inch/0.127mm

Strength characteristics of the catheter:

Bond strength of components - not less than 15 N;
Tensile strength - not less than 8.1 MPa;
Relative elongation at break point - not less than 150%;
The maximum burst pressure - 81 atm. (1200 psi).

MATERIALS

Malecot Catheter (082408; 0082410; 082414; 082416; 082422; 082424)

Device Component	Material Name	Patient Contacting
Only for 082408; 0082410		
Catheter	Catheter tubing: Polyurethane Radiopaque White #20023 Colorants: Bismuth oxychloride, Adapter: Polyamide-6	Direct
Only for 082414; 082416; 082422; 082424		

Device Component	Material Name	Patient Contacting
Only for 082408; 0082410		
Catheter	Catheter tubing: Polyurethane Radiopaque White #20023 Colorants: Bismuth oxychloride, Adapter: Polyamide-6	Direct
Catheter	Catheter tubing: Polyurethane Radiopaque #20393 Blue colorant #20393 Colorants: Titanium dioxide, Bismuth oxychloride CI Pigment Blue 29, CI Pigment Blue 15, Calcium carbonate Carbon black Adapter: Polyamide-6	Direct

Percutaneous Nephrostomy Malecot Catheter (082004)

Device Component	Material Name	Patient Contacting
Catheter	Ethylene-Vinyl Acetate #30213 Colorants: Bismuth oxychloride, CI Pigment Blue 15, Zinc Stearate, Titanium Dioxide, CI Pigment Green 7 Adapter: Polyamide-6	Direct

INTENDED USE/INDICATION FOR USE:

Malecot Catheter (082408; 0082410; 082414; 082416; 082422; 082424) -used for percutaneous placement of a Malecot catheter in the renal pelvis. Intended for one-time use.

Percutaneous Nephrostomy Malecot Catheter (082004) -used for percutaneous placement of a Malecot catheter in the renal pelvis. Intended for one-time use.

CONTRAINDICATIONS:

1. Acute urethritis
2. Acute prostatitis
3. Acute cystitis
4. Urethral abnormalities (diverticula, doubling, etc.) because of the risk of perforation

POTENTIAL ADVERSE EVENTS:

1. Dysuria
2. Fever
3. Postoperative pain or pains in connection with stent placement
4. Incontinence

COMPLICATIONS:

1. Perforation of urinary tracts
2. Infections of the urinary tracts
3. Hematuria
4. Bleeding
5. Pyelonephritis
6. Sepsis

INSTRUCTIONS FOR USE

Malecot Catheter (082408; 0082410; 082414; 082416; 082422; 082424):

Percutaneous Nephrostomy Malecot Catheter (082004):

1. Using preliminary plain film, I.V.P., retrograde pyelogram, ultrasound or CAT scan, localize the kidney to be drained. A prone or slightly oblique position is preferred. Prepare and drape the flank in the usual fashion. Under fluoroscopic control, identify and anesthetize the skin site overlying the collecting system.

2. Introduce the 22 gage (black hub) needle into the collecting system at a right angle to the skin. Obtain a urine sample, if desired. Visualization of the pyelocalyceal system may be further enhanced by the injection of contrast medium.

NOTE: Maintain the 22 gage needle at this location and position throughout the remainder of the procedure. This fixes the kidney in position and allows further injection of the pyelocalyceal system, if desired.

3. Rotate the patient into a more lateral position with the posterior axillary line of the ipsilateral flank in the superior position. Pass the 18 gage (pink hub) needle vertically through

a small incision in the skin (made with a scalpel blade); insert the needle into the appropriate part of the pyelocalyceal system. Free flow of fluid from the needle after removing the obturator confirms a satisfactory location for the tip of the needle cannula.

4. Pass the flexible end of the Safe-T-J® wire guide through the hub of the 18 gage needle cannula; insert the wire guide three more inches into the collecting system. Confirm the position fluoroscopically. Maintain the position of the wire guide as the needle cannula is withdrawn over it. **NOTE:** Care must be exercised to avoid withdrawal of the wire guide itself.

5. To facilitate passage of the nephrostomy catheter, dilate the musculofascial tract by progressing from the smallest to the largest dilator provided.

6. Place the blunt stylet inside the Malecot catheter and secure its position with the Luer lock to straighten the Malecot wings.

7. Pass the Malecot end of the catheter over the external end of the wire guide; gradually advance the catheter end well into the collecting system. Confirm the position fluoroscopically. After releasing the Luer lock in order to open the Malecot wings, carefully remove the blunt stylet.

8. While holding the shaft of the Malecot catheter securely in position with one hand, withdraw the wire guide with the other. When the appropriate position of the Malecot is assured by fluoroscopic visualization, remove the 22 gage needle.

9. Tape or suture the provided retention disc to the skin; secure the catheter to the disc with the plastic pull tie. A dressing may be applied. Use the connecting tube to connect the Malecot catheter to a drainage bag or leg bag. The one-way stopcock is useful for testing the drainage capability of the collecting system, if desired.

NOTE: The individually packaged, flexible catheter stylet should be reserved for straightening the Malecot catheter prior to removal at a later date.

COMPLIANCE WITH RELEVANT REGULATIONS AND TECHNICAL STANDARDS

Cook Incorporated manufactures all medical products in accordance with Medical Devices Directive 93/42/EEC.

Cook Incorporated manufactures all medical products in accordance with the United States Food and Drug Administration (FDA) Quality System Regulations, 21 CFR Part 820, *Quality System Regulations*.

All of Cook facilities are also certified to ISO 13485 (2003/2012) *Quality Management System*.

SHELF LIFE

The “use by” date as denoted by symbols on the labeling. The shelf life is shown by indication the time of manufacture (factory symbol) to the stated expiration date (hourglass symbol). The difference between both dates is the shelf life.

Nephrostomy Catheters have a 3 (three) year shelf-life.

STERILIZATION DETAILS

Nephrostomy catheters are supplied sterile in peel open packages utilizing an EtO (Ethylene Oxide) validated method to a sterility assurance level (SAL) of 10^{-6} . These devices are intended for one-time use, sterile if unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile.

STORAGE RULES

Nephrostomy catheters should be stored in a dark, dry, cool place; avoid extended exposure to light.

TRANSPORTATION RULES

Transportation shall account for the type of transportation employed, the destination, the vendor used for transportation and an estimation of the frequency and harshness of any handling required during transportation.

USAGE CONDITIONS

Nephrostomy catheters are used in the patient care medical institutions.

DELIVERY DETAILS

Nephrostomy Catheters are supplied in individual packaging together with instruction for use. Individual packages are put in the corrugated board boxes. Boxes are being taped with the control tape to cover the sides or along the perimeter.

SUPPLY DETAILS

Nephrostomy catheters are supplied as follows: catheter – 1 p. in individual packaging; instructions for use – 1 p.

TRANSPORTATION LABELLING

Nephrostomy catheters are supplied in the individual packages which are placed in corrugated boxes with sealed with control tape. Transport boxes show the following information: dimensions of the box, unique bar-code, Cook Medical logo and the following symbols:

- symbol for "Fragile"

- symbol for "temperature limitation"

- symbol for "Keep dry"

PERFORMANCE CHARACTERISTICS

Products mentioned in this Technical data summary are stable under vibration loads in the range 10-55 Hz, amplitude displacement of 0.15 mm.

Products mentioned in this Technical data summary, in packaging are stable under vibration loads in the range 10-55 Hz, amplitude displacement of 0.35 mm and shock loads at peak acceleration of 100 m / s² (10g), and duration of shock acceleration of 16 ms.

Products mentioned in this Technical data summary are stable at temperatures in the range of +32 to + 42 ° C.

Products mentioned in this Technical data summary, in packaging are stable at a transportation and storage temperature in the in the range of -50 to + 50 ° C.

Products mentioned in this Technical data summary are stable under sudden changes in temperature.

Products mentioned in this Technical data summary are supplied in sterile packages, do not attempt to clean, disinfect or re-sterilize.

WARRANTY

Manufacturer Cook Medical Europe Ltd., O'Halloran Road, National Technology Park, Limerick, Ireland. Phone: +353 61 334440, fax: +353 61 334441.

Authorized representative on the territory of the Russian Federation -- Company "SCHAG" (ZAO "SCHAG"). Russia, 119002, Moscow, Karmanitzki per. 9. Office 501A

Tel. (495) 956-13-09;
Fax (495) 956-1310,
e-mail: info@schag.ru

Manufacturer's warranty extends to the quality of products.

PACKAGING PARAMETERES

Order number	Height, cm	Width, cm	Weight, g
Tolerance	±10%	±10%	±10%
080206.	42,8	15,2	37,2
080208.	58	15,2	44,4

080210.	58,1	15,2	46
080212.	68,6	15,3	50
082004.	15,4.	43.	30,8.
082408.	15,3	53,3	35.
082410.	53,34.	15,24.	36,4.
082414.	15,4.	58,3.	41.
082416.	15,3	58,4.	48.
082422.	53,34.	15,24.	36,4.
082424.	58,4.	15,3	38.
085108.	63,2.	15,2	59.
085110.	63,1.	15,2	60,8.
085110-S8	58,4.	15,3	47.
085112.	61,5.	15,2	60,8.
085114.	63,3.	15,2	62,8.
023806.	90,7.	12,4.	22.

PACKAGING AND LABELLING

Nephrostomy Catheters are supplied in the individual packaging, which should show no signs of damage or opening. Upon detecting the bends, breaks and other defects the device cannot be used.

The label of the device shows the following device parameters:

- company-manufacturer name and manufacturing address;
- address of the EU representative;
- marking “For specialists only”;
- medical device name;
- part number;
- serial number (LOT);
- parameters and general device drawing;
- symbol about CE certification;
- symbol on sterilization;
- symbol on reusability;
- manufacturing date;
- symbol for keep away from sunlight;
- symbol for keep dry;
- bar-code;
- shelf life;
- symbol for “Consult Instruction for Use”
- symbol for caution “Attention, see Instructions for Use”

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


This Apostille certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears. This Apostille does not certify the content of the document for which it was issued. To verify the issuance of this Apostille, see www.austlii.edu.au/au/other/dfat/ie

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
22-11-2017

Renal Access Cobra Catheter (023806)

The Renal Access Cobra Catheter is 6.0 French/2mm in diameter with a length of 65 centimetres. The catheter is curved towards the distal end in a C-shaped curvature to minimize trauma. It will accept wire guides up to 0.97 millimetres in diameter.

Product Catalogue Number	Catheter		Remarks
	Diameter Fr/mm	Length cm	
Tolerance	10%	10%	-
023806	6/2	65	Accepts a 0.038 inch (0.97 mm) diameter PTFE coated stainless steel wire guide.

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MATERIALS

Renal Access Cobra Catheter (023806)

Device Component	Material Name	Patient Contacting
Catheter	Catheter tubing: Radiopaque Polyamide #40103 Braid: Stainless steel MS#179	Direct
Luer Lock Adapter	Polyamide 6	Direct

INTENDED USE/INDICATION FOR USE:

Renal Access Cobra Catheter (023806) - Used percutaneously to direct a 0.038 inch /0.97 mm diameter wire guide to a specific renal area.

CONTRAINDICATIONS:

1. Acute urethritis

2. Acute prostatitis
3. Acute cystitis
4. Urethral abnormalities (diverticula, doubling, etc.) because of the risk of perforation

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INSTRUCTIONS FOR USE

1. Form the catheter shaft into the required shape.
2. Under fluoroscopic monitoring from the constructed fistula, insert the guide wire into the ureter or the upper ureter loop position in the renal pelvis.
3. Insert the catheter tracking the guide wire and repeat the guide wire insertion in the renal pelvis and confirm that it is in the correct position.

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TRANSPORTATION LABELLING

Nephrostomy catheters are supplied in the individual packages which are placed in corrugated boxes with sealed with control tape. Transport boxes show the following

information: dimensions of the box, unique bar-code, Cook Medical logo and the following symbols:

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Products mentioned in this Technical data summary are stable at temperatures in the range of +32 to +42 ° C.

Products mentioned in this Technical data summary, in packaging are stable at a transportation and storage temperature in the in the range of -50 to +50 ° C.

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20 NOV 2017

Date

COOK MEDICAL EUROPE LTD.
O'Halloran Road,
National Technology Park
Clonsilla, Ireland.
Registration No. 492272
No. 9776234K

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