GE Healthcare

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Chapter 5: China-Specific Information

5.1 China-Specific RoHS Compatibility

This symbol indicates the product contains hazardous materials in excess of the limits established by the Chinese standard SJ/T11363-2006

"Requirements for Concentration Limits for Certain Hazardous Substances in Electronic Information Products".

The number in the symbol is the Environment-friendly Use Period (EFUP), which indicates the period during which the toxic or hazardous substances or elements contained in electronic information products will not leak or mutate under normal operating conditions so that the use of such electronic information products will not result in any severe environmental pollution, any bodily injury or damage to any assets. The unit of the period is **Year**.

In order to maintain the declared EFUP, the product shall be operated normally according to the instructions and environmental conditions as defined in the product manual, and periodic maintenance schedules specified in Product Maintenance Procedures shall be followed strictly.

Consumables or certain parts may have their own label with an EFUP value less than the product. Periodic replacement of those consumable or parts to maintain the declared EFUP shall be done in accordance with the Product Maintenance Procedures.

This product must not be disposed of as unsorted municipal waste, and must be collected separately and handled properly after decommissioning.

Table 5-1: Table of Hazardous Substances’ Name and Concentration

Component name	Hazardous substances’ name					
	(Pb)	(Hg)	(Cd)	(Cr(VI))	(PBB)	(PBDE)
Operator Console	X	O	O	X	O	O
Gantry	X	O	X	X	X	X
LCD Monitor	O	X	O	O	O	O
ECG Cardiac Trigger	X	O	O	X	X	X
Power Distribution Unit	X	O	X	X	X	X
Patient Table	X	O	X	X	X	X
O: Indicates that toxic or hazardous substance contained in all of the homogeneous materials for this part is below the limit requirement in SJ/T11363-2006 X: Indicates that this toxic or hazardous substance contained in at least one of the homogeneous materials used for this part above the limit requirement in SJ/T11363-2006 ■ Data listed in the table represents the best information available at the time of publication ■ Applications of hazardous substances in this medical device are required to achieve its intended clinical uses, and/or to provide better protection to human beings and/or to the environment, due to lack of reasonably (economically or technically) available substances.						

This product consists of devices that may contain mercury, which must be recycled or disposed of in accordance with local, state, or federal laws. (Within this system, the backlight lamps in the monitor display contain mercury).

The Nuclear collimator contains the following potentially hazardous substances:

Lead: Lead salts are toxic and their ingestion may cause serious problems.
The manipulation/handling of lead is subject to regulations.

- Do not discard the collimators among industrial waste or domestic garbage.
- A damaged collimator should not be dispatched through the national postal service.


Technical Publications

Direction 5439352-1EN

Rev. 1

Final - 30-May-2012

NM 600 Series

Nuclear Medicine Imaging Systems

System Description and Safety Manual for Operators

© 2012 by GE Healthcare

© 2012 by GE Healthcare. All Rights Reserved.

Acquisition and/or processing software and the related documentation are confidential and proprietary information of **GE Healthcare**. Only licensees of **GE Healthcare** have a right to use the information contained herein. Only licensees specifically granted copy and/or transfer rights have the right to copy and/or transfer the information. Any unauthorized use, disclosure, assignment, transfer or reproduction of this confidential information will be prosecuted to the full extent of the Law.

Disclaimer

GE Healthcare shall not be liable nor obligated in any manner in respect of bodily injury and/or property damage arising from the use of this system if such use is not in strict compliance with instructions and safety precautions contained in the relevant operating manuals and in all supplements thereto, in all product labels, and according to all terms of warranty and sale of this software, nor if any change not authorized by **GE Healthcare** is made to the software contained herein.

The original manual was written in English.

Warning

User provided programs, scripts or acquisition protocols are NOT validated nor warranted by **GE Healthcare**. The use of data obtained using such user provided programs, scripts or acquisition protocols is the sole responsibility of the party using such programs or protocols.

Users exchanging files and media should beware of the risk of software viruses.

Revision History

Revision	Date	Description of Changes	Chapter/Pages
5439352-1EN r1	June 2012	Updated for O640 release	All O640-related items
5434033-1EN r2	April 2012	■ Generic safety material removed (moved to <i>NM Cameras Safety and Regulatory Manual</i>) ■ B615/D630/D670 IEC Edition 3 conformance ■ System Specifications added	All
5434033-1EN r1	March 2012	■ Manual defined as NM 600 Series ■ Brivo NM 615 and Optima NM/CT 640 systems added, including: ■ System description ■ Safety mechanisms (Safety) ■ Shutdown/Startup procedures	All
5419828-1EN r1	March 2011	■ Discovery NM 630 system added ■ Discovery NM/CT 670 indications added ■ Miscellaneous updates	All
5337418-1EN r2	May 2010	Overall update for D670 release	All
5337418-1EN r1	January 2010	New manual	All

Topics

Before You Start

- Safety Indications in This Document x**
- System Names, Conventions and Images xi**
- System Documentation Set xii**
- Online Access to Operator Manual Set xv**
- Using the Operator Manual Set Effectively xvi**
- Printing this Document xix**
- Conventions in This Document xx**

Chapter 1: System Description

- 1.1 Medical Purpose 1-1**
 - 1.1.1 Brivo NM 615 & Discovery NM 630 Intended Use 1-1
 - 1.1.2 Optima NM/CT 640 Intended Use 1-2
 - 1.1.3 Discovery NM/CT 670 Intended Use 1-4
- 1.2 System Overview 1-6**

1.3	Gantry Description	1-13
1.3.1	NM Detectors	1-14
1.3.2	Collimators	1-17
1.3.2.1	B615 Collimator Carts and Collimator Exchange	1-19
1.3.2.2	Body Contour Apparatus	1-20
1.3.3	System Hand-Held Controller	1-21
1.3.4	O640 CT Rear Panel and Hand-Held Controller	1-26
1.3.5	D670 Rear Control Panel	1-28
1.3.6	Gantry Display	1-29
1.4	Patient Table	1-31
1.5	ECG Trigger Monitor	1-35
1.6	Power Supply Components	1-36
1.6.1	B615,D630 and O640 UPS Support During Power Outages	1-37
1.6.2	D670 Partial UPS Support During Power Outages	1-38
1.7	Operating Consoles	1-39
1.7.1	B615, D630 and O640 Console	1-39
1.7.2	D670 Console	1-40
1.7.3	D670 Working Modes	1-41
1.7.4	D670 Console-dedicated Keyboard	1-43
1.7.5	Application Software	1-45
1.8	System Operating Conditions	1-46

Chapter 2: Standby, Shutdown and Startup Procedures

2.1	Standby	2-3
2.2	B615 & D630 Shutdown and Power OFF	2-4
2.3	O640 Shutdown and Power OFF	2-5
2.4	D670 Shutdown and Power OFF	2-6
2.5	Detectors Temperature When Powering ON from Cold	2-8
2.6	B615 & D630 Powering ON from Cold	2-9
2.7	O640 Powering ON from Cold	2-11
2.8	D670 Powering ON from Cold	2-13

Chapter 3: Emergency/Safety Mechanisms and Procedures

3.1	Emergency Mechanisms for Moving Components	3-1
3.1.1	System Emergency OFF (EMO)	3-2
3.1.2	Emergency Stop Buttons	3-4
3.1.3	Pallet Release	3-7
3.1.4	Detector Manual Release Mechanism	3-8

3.2	Emergency Mechanisms for Moving Components	3-9
3.2.1	Safety Pressure Sensitive Devices (PSD)	3-10
3.2.2	Reset After PSD Activation	3-11
3.2.3	ENABLE Lever on Hand-Held Controller	3-12
3.2.4	Collision Warning Mechanism	3-13
3.3	O640 & D670 X-ray Protection & Emergency Procedures	3-14
3.3.1	O640 Radiation Safety Control Mechanisms	3-14
3.3.2	Radiation Indicators	3-15
3.3.3	Emergency Patient Care During X-ray Exposure	3-17
3.4	Obstruction Free Area and Egress	3-19
Chapter 4: System Specifications		
4.1	Gantry	4-1
4.2	Patient Table	4-4
4.3	Detectors	4-5
4.4	Collimators	4-8
4.5	Acquisition Station	4-12
4.6	Acquisition Modes	4-12
4.6.1	Hybrid SPECT/CT Imaging	4-16
4.7	NM 600 Series Configurations	4-17

4.8	Networking Specifications	4-18
4.8.1	Purpose of Network Connection	4-18
4.8.2	Network Interface Technical Specifications	4-19
4.8.3	Network Information Flow Specifications	4-19
4.9	Environment and Power Requirements	4-26
4.10	System GE Approved Accessories	4-27
4.10.1	UPS Devices	4-28
4.10.2	External ECG Triggers	4-29
4.10.3	Table Accessories	4-29
 Chapter 5: China-Specific Information		
5.1	China-Specific RoHS Compatibility	5-1

List of Figures

1-1	Brivo NM 615 Gantry and Patient Table	1-8
1-2	Discovery NM 630 Gantry and Patient Table	1-9
1-3	Optima NM/CT 640 Gantry and Patient Table.	1-10
1-4	Discovery NM/CT 670 Gantry and Patient Table	1-11
1-5	Other Components in Scan Room	1-12
1-6	NM Detectors	1-14
1-7	Power Supply Components	1-36
1-8	Operator Console	1-40
1-9	D670 Operator Console Screens	1-42
3-1	System Emergency OFF Button	3-2
3-2	Location of Emergency Stop Buttons.	3-5
3-3	Location of Pallet Release Handle	3-7
3-4	Location of Detector Manual Release Mechanism	3-8
3-5	Location of Pressure Sensitive Devices	3-10
3-6	Location of Enable Lever on Hand-Held Controller	3-12
3-7	Location of X-ray ON Indicators	3-16
3-8	Obstruction Free Area	3-20

Before You Start

WARNING

Before any attempt is made to use/service the system, the operator and service personnel must be trained, and must read and be acquainted with all Safety-related Documents, p.xii. This will prepare all users to operate the equipment safely and correctly in order to ensure the well-being of the patient, operator and service personnel.

IMPORTANT

Any safety-related material in this document is classified as a **CAUTION**, unless otherwise stated.

Safety Indications in This Document

This manual uses three safety severity classifications:

DANGER

Danger is used to identify conditions or actions for which a **specific hazard** is known to exist, which **will cause severe or fatal personal injury** or substantial property damage if the instructions are ignored.

WARNING

Warnings are used to identify conditions or actions for which a **specific hazard** is known to exist, which **may cause severe or fatal personal injury** or substantial property damage if the instructions are ignored.

CAUTION

Cautions are used to identify conditions or actions for which a **potential hazard** may exist, which **may cause minor personal injury** or property damage if the instructions are ignored.

System Names, Conventions and Images

NOTE

The following indications are relevant for all documents in the Operator documentation set.

- The images in this manual are for demonstration only. There may be minor differences that do not affect functionality.
- Most of the material provided in this manual is relevant for all NM 600 Series models.
- Some of the described features may be optional, depending on system model/configuration.
- The terms **Hybrid** and **SPECT/CT** are used to designate a combination of SPECT and CT imaging modalities.
- When there are system-specific differences, this is indicated as described below.

Full name	Abbreviation	Markings for system-specific material*	
NM 600 Series or NM600 refers to all of the below systems.			
Brivo NM 615	B615		
Discovery NM 630	D630		
Optima NM/CT 640	O640		
Discovery NM/CT 670	D670		

* When a rectangle indicates more than one system, a combination of the colors for all indicated systems is used, together with the system names.

System Documentation Set

The camera documentation set includes the following documents:

Safety-related Documents

IMPORTANT

To access the electronic safety-related documents, click the **Safety** button on the **How To** (Main Menu) or the **Safety...** bookmark in the documents.

- *NM Cameras Safety and Regulatory Manual*
Contains information about regulatory conformance, safety precautions and emergency procedures for all NM cameras.
- *NM600 System Description and Safety Manual for Operators* (this manual)
Contains an overview description of the NM 600 Series cameras and their components, intended use, safety mechanisms and other system-specific safety information.
- *NM600 Adhesive Labels & Rating Plates & Manufacturing Address User Guides*
Contain all the system safety and regulatory labels.

NOTE

The Rating Plates guides are provided in hard-copy only.

- *Quick Reference Card* (per model; hard-copy and electronic)
Provides emergency procedures and instructions for powering on the system and recovering from emergency situations.

- *Customer Information Sheet (CIS)* (per model; supplied on a separate CD)
Contains open issues and information not covered in other operator documentation (if relevant).
- *Optima NM/CT 640 CT Technical Reference Manual* [REDACTED]
Contains or refers to Optima NM/CT 640 CT-specific safety and regulatory information.

IMPORTANT [REDACTED]

The documentation set does not cover D670 CT-related safety in detail. All D670 CT-related safety procedures, regulatory information and warnings are included in the D670 CT Operator Documents, p.xiv.

Operation-related Documents

- *NM600 NM Clinical Operation Guide*
Contains all operator information required to operate the system for NM-only clinical scans.
- *NM600 Quality Control Guide*
Contains procedures and guidelines for quality control operations performed by the operator.
- *Hybrid Clinical Operation Guide* [REDACTED] [REDACTED] (per model)
Contains all operator information not included in the *NM600 NM Clinical Operation Guide*, that is required to operate the system for hybrid clinical scans.

Additional Reference Documents

- *NM 600 Series User Interface and Acquisition Parameters Reference Guide*
- *Factory Acquisition Protocols Reference Guide* (per model)

NM 600 Series System Description and Safety Manual for Operators
5439352-1EN, Rev. 1, ©2012 GE Healthcare

Patient Accessories Documents

- *NM Cameras Multi-Purpose Arm Support Straps User's Guide*
- *NM Cameras Armband for Arms-Above-Head Scanning User's Guide*
- *NM Cameras Head Holder User's Guide*
- *NM 600 Series Pinhole Collimator User's Guide*

D670 CT Operator Documents

- *BrightSpeed Series Learning and Reference Guide*
Provides all the user information required to operate the CT scanner, as well as related safety information. Contains detailed information as well as step-by-step procedures.
- *BrightSpeed Series Technical Reference Manual*
Contains CT-related safety information, regulatory information, system specifications and technical descriptions.
- *BrightSpeed Series Application Tips and Work-Arounds*

Service Documentation

All service documentation is provided on a separate CD.


Online Access to Operator Manual Set

The operator manual set is accessible online from the **?** in the upper right corner of the **To Do List** screen.

Click any of the links or buttons in the **How To** (Main Menu) to access the relevant information.

Use the **Search** button  or bookmark to search across all documents (see Using Search to Locate Information, p.xvi).

To return to the main menu, use the **How To...** bookmark or **Home** button .

Safety-Related Documents

Click the **Safety** button  on the **How To** (Main Menu) or the **Safety...** bookmark available in all documents to access Safety-related Documents.

Using the Operator Manual Set Effectively

Using Search to Locate Information

- Always use the **Search** function (bookmarks or button) to find information. Search activates a cross-document full-text index that displays a list of the documents and the context in which the word was found.
- Search also provides **Advanced Search Options** where you can query phrases, use boolean operators (AND/OR/NOT, use stemming, “sounds like”, Thesaurus and proximity options, and include bookmarks or comments in the search).
- *Do not* use the Find function (Ctrl+F) as this searches only the current document, and finds the word sequentially without using the full-text index.

NOTE

Documents opened via the Search Index are opened in a separate window, in addition to the currently viewed document.

Commenting

- All documents provide commenting options. Comments are useful to mark and later find frequently-used material, or to send us your comments or requests (you can export comments to a small  file via the **Options** button in the **Comments List**).
- You can add your own comments using the Typewriter  or the Commenting Tools (the most useful tools are **Text Edits** and **Highlight**)

Use the **Show** option to display a list of comments.

- After adding comments, save the file (if you want to retain your comments for future reference inside the document collection, save under the same name and in the same location).

NOTE

After saving a file, you may get an error message stating that the index is not updated. Turn off display of this message.

Final - 30-May-2012


Navigation and Zoom

The navigation bar at the bottom of each page enables you to page back and forth, and to retrace your steps.

■ Navigation:

- The **Home** button  on each page takes you back to the graphic main menu.
- The **Back** button  on each page retraces your browse sequence within and between documents.
- **Zoom:** The default zoom is full-page (fit in window), but you can use the zoom tool  to view the details in any graphic or screenshot. Click  to revert to full page.



Final - 30-May-2012

Printing this Document

This document is created using A5 sheet size. Use the following guidelines when printing:

Print Dialog	Property*	Two pages per sheet (A4/Letter)	One page per sheet (A5 or A4/Letter)	
			A5 is recommended for a compact book and paper saving A4/Letter is recommended when large format and print are needed	
				
Main	Page Scaling	Fit to Printable Area (removes extra white margins)		
Printer Properties (Advanced)	Orientation	Landscape		
	Double-sided	Double-sided or (Print on Both Sides)		
	Page direction / Binding	Open to Side (or Open to Left)	Open to Top or (Flip pages Up)	Open to Side or (Open to Left)
	Number of pages per sheet	2	1	1
	Page Borders	Print Page Borders (adds visual separation between the two pages)	N/A	N/A

* Names of properties and options can differ, depending on your specific printer driver

Conventions in This Document

IMPORTANT

Calls attention to important comments.

NOTE

Contains tips and general comments.

The following conventions are used throughout the manual:

Description	Example
Keys on the operator keyboard, hand-held controller and gantry	<SET>, <Ctrl>
Software interface buttons	[OK], [Apply], [Cancel]
Names of items in the graphical interface including: - Names of dialog boxes, windows, tabs, areas and lists - Menu items - Field and icon labels	Configuration tab; To Do List File menu Gantry icon; Properties field
System messages	–Æ ≠ « ” ± Ω±≤” Σ≤’ a U
System parameters whose actual values must be defined by the user	Type-in the –ø” Σ” ≤” ◇◁
Hyperlinks	Figure 3-1
Paths	Æ±±”, ±∞”, ” øΩØø æ ≥ø≤’ ø¥#
References to other documents	<i>Operator Manual</i>
End of a procedure	◆

Chapter 1: System Description

1.1 Medical Purpose

1.1.1 Brivo NM 615 & Discovery NM 630 Intended Use

NOTE B615 & D630

This section is relevant only for Brivo NM 615 and Discovery NM 630.

The system is a medical tool intended for use by appropriately trained healthcare professionals to aid in detecting, localizing, diagnosing of diseases and in assessment of organ function for the evaluation of diseases, trauma, abnormalities, and disorders such as, but not limited to, cardiovascular disease, neurological disorders and cancer. The system output can also be used by the physician for staging and restaging of tumors, planning, guiding, and monitoring therapy.

The system is a Nuclear Medicine (NM) system, which is intended to yield General Nuclear Medicine imaging procedures for detection of radioisotope tracer uptake in the patient body, using a variety of scanning modes supported by various acquisition types and imaging features designed to enhance image quality.

Scanning modes include Planar mode (Static, Multi-gated, Dynamic and Whole body scanning) and Tomographic mode (SPECT, Gated SPECT, Whole body SPECT).

The acquisition types include single and multi-isotope/multi peak frame/list mode single-photon imaging. The imaging-enhancement features include an assortment of collimators, gating by physiological signals, and real-time automatic body contouring.

The system may include signal analysis and display equipment, patient and equipment supports, components and accessories. The system may be used for patients of all ages.

1.1.2 Optima NM/CT 640 Intended Use

NOTE 0640

This section is relevant only for Optima NM/CT 640.

The Optima NM/CT 640 system is a medical tool intended for use by appropriately trained healthcare professionals to aid in detecting, localizing, diagnosing of diseases and in assessment of organ function for the evaluation of diseases, trauma, abnormalities, and disorders such as, but not limited to, cardiovascular disease, neurological disorders and cancer. The system output can also be used by the physician for staging and restaging of tumors, planning, guiding, and monitoring therapy.

The Optima NM/CT 640 system is a Nuclear Medicine (NM) system, which is intended to yield General Nuclear Medicine imaging procedures for detection of radioisotope tracer uptake in the patient body, using a variety of scanning modes supported by various acquisition types and imaging features designed to enhance image quality. The scanning modes include planar mode (Static, Multi-gated, Dynamic and Whole body scanning) and tomographic mode (SPECT, Gated SPECT, Whole body SPECT). The acquisition types include single and multi-isotope/multi peak frame/list mode single-photon imaging. The imaging-enhancement features include an assortment of collimators, gating by physiological signals, and real-time automatic body contouring.

The Optima NM/CT 640 system includes an integrated CT, which is intended specifically for attenuation correction and anatomical localization purposes only, a processing and review workstation, and may include signal analysis and display equipment, patient and equipment supports, components and accessories.

Utilizing the associated CT portion and post processing, the device can produce attenuation corrected SPECT images using CT-based attenuation maps, as well as functional and anatomical mapping images registered with the SPECT images for purposes of localization and fusion.

The device does not support stand-alone CT operation. The system may be used for patients of all ages.

Final - 30-May-2012

1.1.3 Discovery NM/CT 670 Intended Use

NOTE D670

This section is relevant only for Discovery NM/CT 670.

The Discovery NM/CT 670 system is a medical tool intended for use by appropriately trained healthcare professionals to aid in detecting, localizing, diagnosing of diseases and organ function for the evaluation of diseases, trauma, abnormalities, and disorders such as, but not limited to, cardiovascular disease, neurological disorders and cancer. The system output can also be used by the physician for staging and restaging of tumors, planning, guiding, and monitoring therapy. The Discovery NM/CT 670 system, combining Nuclear Medicine (NM) and Computed Tomography (CT) systems, is intended to produce:

- **NM System:**

General Nuclear Medicine imaging procedures for detection of radioisotope tracer uptake in the patient body, using a variety of scanning modes supported by various acquisition types and imaging features designed to enhance image quality. The scanning modes include planar mode (Static, Multi-gated, Dynamic and Whole body scanning) and tomographic mode (SPECT, Gated SPECT, Whole body SPECT). The acquisition types include single and multi-isotope/multi peak frame/list mode single-photon imaging. The imaging-enhancement features include an assortment of collimators, gating by physiological signals, and real-time automatic body contouring.

- **CT System:**

Cross sectional images of the body by computer reconstruction of x-ray transmission data taken at different angles and planes, including Axial, Cine, Helical, Cardiac, and Gated acquisitions. These images may be obtained with or without contrast. The CT system is indicated for head, whole body, cardiac and vascular x-ray Computed Tomography applications.

- **NM+CT System (Hybrid):**

Combined, hybrid SPECT and CT protocols, for CT-based SPECT attenuation corrected imaging as well as functional and anatomical mapping imaging (localization, registration and fusion).

The GE Discovery NM/CT 670 system may include signal analysis and display equipment, patient and equipment supports, components and accessories. The system may include data and image processing to produce images in a variety of trans-axial and reformatted planes. The images can also be post processed to produce additional images, imaging planes, and analysis results.

Final - 30-May-2012

1.2 System Overview

The NM 600 Series includes the following models:

	Brivo NM 615	Discovery NM 630	Optima NM/CT 640	Discovery NM/CT 670
	Figure 1-1, p.1-8	Figure 1-2, p.1-9	Figure 1-3, p.1-10	Figure 1-4, p.1-11
System type	Single-detector variable-geometry NM-only camera	Dual-detector variable-geometry NM-only camera	Hybrid SPECT-CT system, consisting of: ■ A dual-detector variable-geometry NM camera A low mA fast-rotation-speed CT sub-system	Hybrid SPECT-CT system, consisting of: ■ A dual-detector variable-geometry NM camera ■ A multi-detector CT scanner
Location of operating console	Can be installed in a separate Operator room or inside the Scan room, depending on site preference and subject to local regulations		May require radiation protection. Can be in separate Operator room or inside the Scan room, depending on site preference and subject to local regulations.	Requires radiation protection. Must be in a separate radiation-protected Operator room
			The operator must be at the console to activate the CT part of hybrid scans.	

NOTE

See Ch.4, System Specifications for technical specifications.

The system includes the following components:

Table 1-1: System Components

Scan Room		Operator Area/Room
Main System Components (see Figure 1-4, p.1-11)	Other Components (see Figure 1-5, p.1-12)	
■ Gantry (see page 1-13) ■ Patient table (see page 1-31)	■ Collimators and carts (see page 1-17) ■ Power supply components (see page 1-36) ■ ECG trigger monitor (see page 1-35 and manufacturer documentation) ■ D670 Contrast Injector (see manufacturer documentation) ■ Patient positioning accessories (see Table 1-5, p.1-34)	■ System acquisition computer ■ Operator console (see page 1-39) ■ Processing and reviewing tools: ■ Xeleris 3.x workstation ■ D670 Advantage Windows (AW) workstation

Final - 30-May-2012

Final - 30-May-2012


Figure 1-1: Brivo NM 615 Gantry and Patient Table

NM 600 Series System Description and Safety Manual for Operators
5439352-1EN, Rev. 1, ©2012 GE Healthcare

1-8

Final - 30-May-2012



Figure 1-3: Optima NM/CT 640 Gantry and Patient Table

NM 600 Series System Description and Safety Manual for Operators
5439352-1EN, Rev. 1, ©2012 GE Healthcare

1-10

Final - 30-May-2012

**Figure 1-4: Discovery NM/CT 670 Gantry and Patient Table**

NM 600 Series System Description and Safety Manual for Operators
5439352-1EN, Rev. 1, ©2012 GE Healthcare

1-11

Final - 30-May-2012

One or more collimator carts and Collimators, p.1-17

External ECG Trigger Monitor, p.1-35

Power Supply Components, p.1-36
(specific components depend on system configuration)


Figure 1-5: Other Components in Scan Room

1.3 Gantry Description

The NM detectors are mounted on the front of the gantry, providing various detector modes, that support multiple NM scan types (see B615 NM Detector Modes, p.1-15 and D630, O640, and D670 NM Detector Modes, p.1-16).

The CT sub-system of the O640 and D670 models is located at the rear of the gantry, adding hybrid imaging capabilities and supporting attenuation correction and anatomical localization for SPECT NM studies.

D670 also enables CT-only scan types (see Discovery NM/CT 670 Intended Use, p.1-4).

Final - 30-May-2012

1.3.1 NM Detectors

Each NM detector is covered with a collimator that enables imaging and provides thermal and physical protection for the detector's scintillation crystals (see Collimators, p.1-17).

The D630, O640 and D670 models include two NM detectors each, marked ● (Detector 1) **A** and ● ● (Detector 2) **B**, while B615 includes a single detector.

The NM detectors enable imaging in different modes, depending on the system, as detailed in:

- Table 1-2 "B615 NM Detector Modes", p.15
- Table 1-3 "D630, O640, and D670 NM Detector Modes", p.16

For additional information, see *NM600 NM Clinical Operation Guide*.


Figure 1-6: NM Detectors

Table 1-2: B615 NM Detector Modes

Mode	Description		Mode	Description
Radial 	The detector revolves around the patient		Plantar 	The detector faces upwards
Vertical Inside mode (Hi and Lo) 	Vertical Outside mode (Hi and Lo) 	The detector faces inwards/outwards vertically	Stretcher mode 	The detector faces downwards

Final - 30-May-2012

Table 1-3: D630, O640, and D670 NM Detector Modes

Mode	Description	Mode	Description
L-mode 	Detectors revolve at a fixed 90° angle around the patient	Vertical Inside mode (Hi and Lo) 	Both detectors vertically facing inward
H-mode 	Detectors revolve around the patient, parallel to each other	Stretcher mode 	Both detectors facing downwards
Vertical Outside mode (Hi and Lo) 	Both detectors vertically facing outward	Plantar 	Both detectors facing upwards

1.3.2 Collimators

Each collimator consists of a thick lead plate with holes, designed to allow only those gamma rays traveling in a certain direction to strike the detector, while other gamma rays, with specific energies, are stopped by the lead plate. Besides their role in imaging, the collimators provide thermal and physical protection for the scintillation crystals of the detectors.

 **CAUTION**

The detectors must always be loaded with collimators so that the crystal is not exposed. During procedures that require the collimators to be removed, decoy collimators must be installed to protect the detector.

The collimators available for the system include:

Collimator Type	Examples of isotopes used with collimator
LEHR (Low Energy High Resolution)	Tc ^{99m} , Tl ²⁰¹
ELEGP (Extended Low Energy General Purpose)	Tc ^{99m} , I ¹²³
MEGP (Medium Energy General Purpose)	Ga ⁶⁷ , In ¹¹¹
HEGP (High Energy General Purpose)	I ¹³¹
Pinhole (Pinhole with 3 inserts)	Tl ²⁰¹ , Tc ^{99m} , I ¹²³ , I ¹³¹
Fan-Beam (Low Energy High Resolution Fan-Beam for brain imaging)	Tc ^{99m}
Note: B615 Not available for B615.	

The outer face of each collimator includes:

- Safety Pressure Sensitive Devices (PSD), p.3-10 which prevents inadvertent collision between the patient and the detector during normal imaging procedures.
- The Body Contour Apparatus, p.1-20 housing (optional for **B6**).

When not on the detectors, collimators are stored on collimator carts. For detailed instructions and procedures for replacing collimators and decoy collimators, see *NM600 NM Clinical Operation Guide*.


Final - 30-May-2012

1.3.2.1 B615 Collimator Carts and Collimator Exchange

The B615 collimator cart holds a single collimator, balanced by a heavy mass weight **A** in the cart basis.

Avoid applying additional weight on the balancing mass, as it may interfere with the mechanical structure of the cart.

1.3.2.2 Body Contour Apparatus

The Body Contour apparatus (see Figure 1-6, p.14) consists of a transceiver and receiver pair fitted on each collimator. This option enables the detectors to track the patient's body profile during Tomo, Gated Tomo and Whole Body acquisitions.

NOTE**B615 Only**

Body Contour mechanism is an option. Even if your specific system configuration does not include the mechanism, its housing is still present on the detector.

Final - 30-May-2012


1.3.3 System Hand-Held Controller

The hand-held controller (also known as the RCU) enables you, from alongside the gantry, to do the following:

- Control system movement, including:
 - Table movement controls, p.1-23
 - Detector movement controls, p.1-24
 - Special motion and speed control, p.1-25
- Initialize and stop acquisitions (Acquisition controls, p.1-22)
- Confirm actions and pop-up dialogs displayed on the gantry display (Interaction controls, p.1-22)

All movement and acquisition actions require that the relevant key is pressed simultaneously with pressure applied on the **<ENABLE>** lever **A**.

This ensures that no movement is initiated by unintentionally pressing a key on the controller.

Final - 30-May-2012



Figure 1-2: Discovery NM 630 Gantry and Patient Table

NM 600 Series System Description and Safety Manual for Operators
5439352-1EN, Rev. 1, ©2012 GE Healthcare

Interaction controls

- A Up/Down arrows** Select an action or item from a list
- B <OK>** Confirm the selected action, item or pop-up dialog
- C** (Not in use)
- D <CANCEL>** Close the dialog box or pop-up menu without making any changes

**Acquisition controls****IMPORTANT**

All acquisition controls require simultaneous pressure on the **<Enable>** button **H**.

- I <UNLOAD>** Used to move system for patient removal. At end of acquisition, moves detectors outwards, moves pallet away from the gantry and lowers table.
- F <GO/PAUSE>**
 - During Persistence, starts the acquisition
 - During an NM acquisition with detectors or table motion, pauses the acquisition and opens the **Scan Pause** dialog, allowing selection of commands such as **Continue**, **Store** or **Quit** at the console or via the Controller interaction buttons. Also used to continue a paused NM acquisition involving detectors or table motion.
- G <SET>** Move detectors to the selected mode, prior to acquisition start



Table movement controls

IMPORTANT

All table movement controls require simultaneous pressure on the **<Enable>** button **H**.

- I <OUT>** Moves the pallet away from the gantry
- J Up/Down arrows** Move the patient table upward and downward respectively. If the system detects a potential collision between the table and the gantry, a warning is displayed.
- K <IN>** Moves the pallet towards the gantry

Detector movement controls**IMPORTANT**

All detector movement controls require simultaneous pressure on the **<Enable>** button **ⓘ**.

NOTE B615


All detector movement controls apply to the system's single detector.

NOTE D670 O640

- When the table is in CT position, the detector movement keys are disabled.
- See also O640 CT Rear Panel and Hand-Held Controller, p.1-26 and D670 Rear Control Panel, p.1-28.


- L Up/Down arrows** Move the selected detector/s upward and downward respectively when in Vertical, Stretcher or Plantar modes (select detector/s via **ⓘ**).
- M Rotate or Swivel arrows** Rotate both detectors clockwise or counterclockwise (L-mode or H-mode). In modes Vertical In/Out Lo/Hi, Stretcher and Plantar, when a single detector is selected, swivels that detector.
- N <IN>** Moves selected detector(s) closer to the patient
- P <OUT>** Moves selected detector(s) away from the patient
- O Detector Selector** Click to select detector 1 (●), detector 2 (●●) or both. LEDs indicate the current selection.

Note: B615 Not relevant for B615.

Special motion and speed control

- R Lateral (Right/Left) (pinhole collimator only)
 Moves detectors on lateral axis, when permitted.
 Actual motion direction (right/left or up/down) depends on gantry rotation.
- S **<SLOW/
FAST>** Switch between FAST and SLOW movement of the gantry, detectors and patient table.
 LEDs show whether SLOW or FAST is selected.

Final - 30-May-2012

1.3.4 O640 CT Rear Panel and Hand-Held Controller

The system has two identical rear panels, on either side of the gantry. The rear panels include Emergency Stop Buttons, p.3-4 (red) and reset buttons (green).


The CT hand-held controller is located at the rear of the system, integrated into one of the CT rear panels. The controller is active only when the table telescopic extension is in CT position.

The functionality of the CT hand-held controller is similar to that of the system hand-held controller (see System Hand-Held Controller, p.1-21), as follows:

- Acquisition controls, p.1-22 **A** and Table movement controls, p.1-23 **B** are identical to their use on the gantry display hand-held controller.
- The **Detector Motion C** keys are disabled.

NOTE O640

This section is relevant only for Optima NM/CT 640.


System Description

- In the **Interaction** controls, when the pallet is in CT position:
 - **<OK><MARK POINT>** **E** is used to mark the current table linear position for scout range definition.
 - **<CANCEL>** **D** is reserved for future functionality.
- In the **Special Motion** controls:
 - **F** Turns the sagittal and transaxial planes positioning light on/off.
 - **G** Turns the coronal and transaxial planes positioning light on/off.

1.3 Gantry Description

1.3.4 O640 CT Rear Panel and Hand-Held Controller

1.3.5 D670 Rear Control Panel

The rear control panel provides two sets of the following controls:

Emergency Controls

- **E** Emergency stop (stop all activity)
- **L** Emergency stop reset

Starting CT scans in the room

- **D** Start scan with preparation delay
- **F** Preparation delay timer display

Breathing lights demo

- **G** Start breathing lights demo

Laser Alignment and Landmarks

- **A** Laser alignment light on/off
- **B** Set external landmark
- **C** Set internal landmark

Table Control


- **I**—In **H**—Out **M**—Up **K**—Down
- **J** Table motion with fast speed

NOTE

Controls marked  are not in use.

NOTE **D670**

This section is relevant only for Discovery NM/CT 670.

IMPORTANT

- For detailed usage instructions, refer to D670 CT Operator Documents, p.xiv.
- For safety-related instructions, refer to the *NM Cameras Safety and Regulatory Manual*.


1.3.6 Gantry Display

The gantry display is installed on a boom that can be rotated to the right and left, while the monitor itself can be rotated left and right and tilted up and down.

In addition to its clinical functionality, the monitor can be used for patient entertainment (when the option is available on the system), by pivoting to horizontal position, facing the patient on table.

The gantry display enables the operator to do the following:


- During table and patient setup:
 - View table offset/height and pallet location data **C**
 - View detector orientation **D**
- During idle mode, select gantry modes and activate procedures **A**
-  : Select CT mode **B** **E**

System Description


- During scan persistence, scan settings **F** and persistence images **G**
- During acquisition, view progress information **H** and acquired/acquiring images **I**
- At end of the NM acquisition, view the acquired images **J**

Final - 30-May-2012

1.3 Gantry Description

1.3.6 Gantry Display

NM 600 Series System Description and Safety Manual for Operators
5439352-1EN, Rev. 1, ©2012 GE Healthcare

1-30

1.4 Patient Table


The patient table is designed to support up to 227 kg (500 lbs). Overloading the table may result in injury.

The patient table is used to support and position the patient during scans. The entire table can be raised and lowered, and the pallet can be moved towards and away from the gantry, or through the gantry bore.

Table and pallet motion can be controlled using the System Hand-Held Controller, p.1-21 or automatically by applying one of the pre-defined D630, O640, and D670 NM Detector Modes, p.1-16.

In addition, the table can be swiveled sideways for collimator exchange, and for stretcher and vertical scans.

The table components are described in Table 1-4, p.1-32.

Patient positioning accessories used with the table are described in Table 1-5, p.1-34.

Table 1-4: Patient Table Components

Name	Function
A Pallet	Made of carbon fiber which has reduced attenuation. Slides IN/OUT towards and away from the gantry with a travel range of 200 cm (78"). In an emergency, can be moved manually by pulling or pushing the pallet, using the pallet release handle C. O640 D670: Additional travel to the CT position in the gantry bore is achieved using the table telescopic extension E.
B Mattress	The mattress is placed on top of the pallet. The adhesive strips on the two sides beneath the mattress fix the mattress to the pallet. In addition, four detachable pairs of straps enable more accurate patient positioning and help to keep the patient stationary during scanning procedures.
C Pallet release handle	The pallet release handle is located at the head of the patient table. The handle has to be pressed to enable moving the pallet manually.
D Side ruler	When the pallet is fully retracted, allows reading table longitudinal coordinates (in cm) for the purpose of scan range definition. Optionally, the buttons of the side ruler can be used to interactively define the scan range. Note: B615 The Interactive ruler option is not supported.

Final - 30-May-2012

Table 1-4: Patient Table Components (cont.)

Name	Function
O640 D670 ③ Table Telescopic Extender	Transports the pallet between NM and CT scan positions; a support within the bore ensures consistent table sagging between NM and CT.
④ Base unit	Supports the table, providing up/down motion.
⑥ Foot clutch	The foot clutch is located at the base of the patient table at the end closest to the gantry. To swivel the table, press the foot clutch and push the table aside.

Table 1-5: Patient Positioning Accessories

Name	Function	
■ Note the safety precautions in the <i>MM Cameras Safety and Regulatory Manual</i> ■ Some of the accessories in this list are options ■ Each accessory is supplied with stand-alone documentation (accessible via the How To ? help symbol)		
Arm positioning safety straps	Used to immobilize the patient's arms, for example: ■ Above the head during cardiac SPECT scans ■ Alongside the body during whole body scanning	
Head holder	The head holder is an extender of the pallet, used for head imaging, beyond the limits of the pallet. Maximum load on neck and head support: 34 kg (75 lbs.).	Place the head holder extender under the mattress, ensuring correct placement as detailed in the <i>Head Holder User's Guide</i> , accessible via the Patient Accessories link in the How To page (see Online Access to Operator Manual Set, p.xv).
Detachable pallet extender	Used to support the patient's feet, adding 40 cm to the pallet length. Maximum load: 30 kg (66 lbs.). It is not recommended to include the pallet extender within the scanned Field of View.	To attach the extender, place the extender under the mattress with at least 30 cm (12") overlap over the pallet. Make sure that the fastener strips hold the extender securely in place.
  CAUTION Remember that with the head holder or pallet extender attached, the pallet length is increased. Always verify clearance when the table is fully in so that the extender does not collide with the rear wall of the room.		

1.5 ECG Trigger Monitor

For gated cardiac scanning, the camera operates in conjunction with an ECG trigger monitor.

IMPORTANT

Safety instructions are provided in the *NM Cameras Safety and Regulatory Manual*.

- **630, 630, 640**
These models use an internal ECG trigger monitor, optionally supporting an external one.
The internal ECG trigger monitor supports the following features:
 - Three patient cable leads. When not in use, the cables can be temporarily disconnected via a connector located near the pallet release handle.
 - An R-wave gating trigger signal, used to synchronize the acquisition gated signal.
 - An ECG signal data, displayed alongside the persistence and acquiring data in the **Triggers** area of the console and on the gantry display
- **670** Uses an external ECG trigger monitor, usually mounted on a cart/stand alongside the system.
- For GE-approved external ECG trigger monitors equipment, see External ECG Triggers, p.4-29.

Final - 30-May-2012

1.6 Power Supply Components

The system includes the following power supply components:

System	Component
D670 Only	■ A Main Disconnect & UPS* Control panel (MDC) (A1) which includes the main power disconnect buttons for the entire system ■ B Power Input Distribution Box for the entire system ■ C CT Power Distribution Unit for the CT component ■ D Optional UPS unit (acquisition station only)†
B615 D630	■ E Optional UPS unit (entire system)‡
O640	■ E Optional UPS unit (NM sub-system)‡

* UPS = Uninterruptible Power Supply

† See D670 Partial UPS Support During Power Outages, p.1-38

‡ See B615,D630 and O640 UPS Support During Power Outages, p.1-37


Figure 1-7: Power Supply Components

1.6.1 B615,D630 and O640 UPS Support During Power Outages

NOTE **B615** **D630** **O640**

This section is relevant for Brivo NM 615 and Discovery NM 630 and Optima NM/CT 640 which include the Uninterruptible Power Supply (UPS) option.

In the case of a power outage during a scan, the NM system continues to operate fully. At the UPS, the alarm is heard and the front panel Battery LED lights up.

- Complete the scan and unload the patient
- Power down the system, see:
 - B615 & D630 Shutdown and Power OFF, p.2-4
 - O640 Shutdown and Power OFF, p.2-5
 - D670 Shutdown and Power OFF, p.2-6

The UPS enables 15 to 40 minutes of NM system operation, depending on the scan type and the condition of the UPS battery.

For GE-approved UPS specifications, see UPS Devices, p.4-28.

Final - 30-May-2012

1.6.2 D670 Partial UPS Support During Power Outages

NOTE **D670**

This section is relevant only for Discovery NM/CT 670.

In the case of a power outage during a scan, the gantry stops and error messages are displayed at the acquisition station.

- Unload the patient (see Pallet Release, p.3-7 and Detector Manual Release Mechanism, p.3-8)
- Shut down the acquisition station (see step 1 of D670 Shutdown and Power OFF, p.6).

The UPS provides a partial power supply to enable orderly shutdown of the acquisition station. Data acquired prior to the power outage is saved and can be manually transferred to the remote DICOM station.

For details regarding full-system UPS (NM and CT sub-systems), refer to the GE accessories list.

Final - 30-May-2012

1.7 Operating Consoles

1.7.1 B615, D630 and O640 Console

NOTE **B615** **D630** **O640**

This section is relevant for Brivo NM 615 and Discovery NM 630 and Optima NM/CT 640.

The console can be located in a separate Operator room or in the Scan room. The site is responsible for supplying the console desk.

The B615, D630 and O640 acquisition stations consist of a computer, a single screen and standard keyboard and mouse.

The O640 console also includes a communication device enabling operator-patient dialogue when the operator is seated at the console.

Final - 30-May-2012

1.7.2 D670 Console

NOTE [D670]

This section is relevant only for Discovery NM/CT 670.

The console is located in a dedicated radiation-protected operator room.

The console table (described in detail in the D670 CT Operator Documents, p.xiv) includes the NM and CT acquisition stations, which are operated via the ON/OFF buttons **D**.

On the table are:

- Dual monitors **A** which display either CT or NM interfaces, depending on your D670 Working Modes, p.1-41
- An alphanumeric keyboard **C**
- A D670 Console-dedicated Keyboard, p.1-43 **B**

Figure 1-8: Operator Console

Final - 30-May-2012

1.7.3 D670 Working Modes

NOTE **D670**

This section is relevant only for Discovery NM/CT 670.

The console provides two working modes:

- NM Mode – scheduling, NM scanning interface (see *NM600 NM Clinical Operation Guide* and *Hybrid Clinical Operation Guide*)
- CT Mode – CT scanning interface (see D670 CT Operator Documents, p.xiv)

CT Mode can be activated for clinical scans only if the NM component of the system is fully functional, including:

- NM acquisition station is ON
- Acquisition software is running
- NM Gantry is ON
- NM Startup tests have been passed successfully

Final - 30-May-2012


Figure 1-9: D670 Operator Console Screens

For hybrid scanning, the console automatically switches between NM and CT modes, depending on the current task. Manual switching between NM and CT modes is done by clicking the **CT** icon **A** to move to CT mode or clicking the **SPECT** icon **B** to return to NM mode. This may be required for example to initiate CT warmup, CT FastCal and system setup, or for CT only acquisition or processing. Manual switching can be done regardless of the background color of the icon.


1.7.4 D670 Console-dedicated Keyboard

NOTE **D670**

- This section is relevant only for Discovery NM/CT 670.
- Two keyboard models are shown below; your system may have one or the other. There is no functional difference between the models.

Final - 30-May-2012

The following table describes the dedicated keyboard CT-related controls and indicators:

CT Scan Controls			Patient-Operator Intercom Controls	
A	<Stop Scan> – Aborts the scan immediately	Resume scan by clicking [Resume] in View/Edit screen	J	<Press and Talk> enables operator to communicate with the scan room, using the intercom J and the microphone N
B	<Pause> – pauses scanning after current scan is completed		K	Patient to operator
C	<Start Scan> – flashes green once the tube has reached exposure speed; press to start scan. Flashes for 30 seconds before timing out. When timed out, press to revert to ready state.		L	<Volume Controls> Autovoice to gantry
D	Exposure Indicator – illuminates amber during x-ray exposure		M	Operator to patient
E	<Emergency Stop> – stops all system motion and scanning in progress (see Emergency Stop Buttons, p.3-4).			
F	(not in use)			
G	<Stop Move> – Stops table motion.	Resume scan by clicking [Resume] in View/Edit screen.		
H	<Move to Scan> ■ Flashes green when system is ready to move table to starting position ■ Press to initiate motion ■ Flashes for 180 seconds before timing out			

Final - 30-May-2012

1.7.5 Application Software

The application software, accessed via the operating console, provides functionality for:

- Data acquisition and basic data review
- Data transfer
- System control
- System maintenance

NOTE

Certain capabilities are license-dependent. For more information, contact your GE representative.

While the application software available on the console allows basic data review, image interpretation is usually performed on processing and review workstations connected to the system, like the Xeleris or PACS.

Data transfer between the system and the external workstations is performed according to the DICOM v.3.0 standard. For the DICOM 3.0 Conformance Statement, see:
http://www.gehealthcare.com/user/interoperability/dicom/products/nuclear_dicom.html

Final - 30-May-2012

NOTE

- By default, the NM data acquired by the system automatically includes corrections such as energy, linearity, uniformity (sensitivity), radioisotope physical decay, center-of-rotation (for SPECT).

Although some of these corrections can be cancelled (for example, in order to cancel decay correction for a specific energy session or modify the half-life time of the isotope to zero), this is not advisable because many Xeleris applications consider the corrections as having being applied for all data generated by this system.

- **B615 D630 Only:**

As a stand-alone gamma cameras, Brivo NM 615 and Discovery NM 630 do not use a CT or any other external radioactive source to provide attenuation correction information to be processed in the Xeleris Workstation.

However, the integrated systems with Xeleris enable import of CT images generated on other scanners to the Xeleris workstation and, in specific cases, optionally using them for attenuation correction, after confirmation of proper co-registration between the NM and the CT datasets.

1.8 System Operating Conditions

This system was designed to work under specific operating conditions and environment. For example, maintaining constant temperature and humidity levels is essential in order to ensure system functionality over time. Overheating or underheating or changes in humidity that exceed the requirements specified in Environment and Power Requirements, p.4-26 can cause technical difficulties, system failures or damage to system components. For more details, see *Equipment and Environmental Safety* chapter in the *NM Cameras Safety and Regulatory User Guide*.

Chapter 2: Standby, Shutdown and Startup Procedures

- Standby, p.2-3 is applicable when the system will not be used for a short period of time, for instance overnight or over the weekend. In standby mode you prevent access to the application station and reduce system power usage while keeping the NM detectors warm.
- Shutdown is applicable when power must be disconnected from all subsystems. In shutdown mode the application station, gantry and detectors are powered off and reset. See:
 - B615 & D630 Shutdown and Power OFF, p.2-4
 - O640 Shutdown and Power OFF, p.2-5
 - D670 Shutdown and Power OFF, p.2-6

If any light, monitor or visual indicator stays on after the system is shut down, this may indicate an electrical hazard. In such a case, to prevent possible injury, once the system has been shut down, press the System Emergency OFF (EMO), p.3-2 button and contact your service representative immediately.

Following Shutdown, follow the relevant procedure to reactivate the system:

- B615 & D630 Powering ON from Cold, p.2-9
- O640 Powering ON from Cold, p.2-11
- D670 Powering ON from Cold, p.2-13

WARNING

You may have to wait for the system temperature to stabilize before you can power on the system. For details, see Detectors Temperature When Powering ON from Cold, p.2-8.

IMPORTANT



See *NM600 NM Clinical Operation Guide* for information on Recovering from Errors, including instructions for gantry reset and rebooting the acquisition station and software.

Final - 30-May-2012

2.1 Standby

1. Click [X] **A**.
2. In the **Shut Down Application** dialog, click [Log off].
3. Press the yellow **<Standby>** button **B** on the gantry. ♦



2.2 B615 & D630 Shutdown and Power OFF

NOTE B615 D630

This section is relevant for Brivo NM 615 and Discovery NM 630.

1. Click [X] **A**.
2. In the **Shut Down Application** dialog, click [**Shutdown**].
3. Ignore the logon dialog and wait until acquisition station power shuts down.
4. At the system, push down the main circuit breaker **B** to the **Off** position.
5. If a UPS unit is used, switch off the unit as described in the manufacturer's instructions. ♦



NM 600 Series System Description and Safety Manual for Operators
5439352-1EN, Rev. 1, ©2012 GE Healthcare

2.3 O640 Shutdown and Power OFF

NOTE O640

This section is relevant only for Optima NM/CT 640.

1. Click [X] **A**.
2. In the **Shut Down Application** dialog, click [**Shutdown**].
3. Ignore the logon dialog and wait until acquisition station power shuts down.
4. At the NM system, push down the main circuit breaker **B** to the **Off** position.
5. If a UPS unit is used, switch off the unit as described in the manufacturer's instructions.
6. At the CT gantry, push down the main circuit breaker **C** to the **Off** position. ♦



NM 600 Series System Description and Safety Manual for Operators
5439352-1EN, Rev. 1, ©2012 GE Healthcare

2.4 D670 Shutdown and Power OFF

NOTE **D670**

This section is relevant only for Discovery NM/CT 670.

1. In the Operator room:
 - a. Shut down all currently running CT applications.
 - b. Press the CT shutdown icon **A**.
 - c. Switch to the NM console by sequentially pressing **<Ctrl>-<Ctrl>-1** on the keyboard.
 - d. At the NM console, click **[X]** **B**.
 - e. In the **Shut Down Application** dialog, click **[Shutdown]**.



Final - 30-May-2012

2. In the Scan room:

- a. Push down the **NM gantry main circuit breaker** ③ to the Off position.
- b. If a UPS unit is used, switch off the unit as described in the manufacturer's instructions.
- c. At the **Main Disconnect & UPS Control panel (A1)**, press the red **<Off>** button ④. ◆

2.5 Detectors Temperature When Powering ON from Cold

Before powering on the system, ensure that the environmental conditions specified in Environment and Power Requirements, p.4-26 are reached and have remained stable for at least 10 minutes.

For the NM detectors, a stable operating temperature must be reached to ensure optimal performance.

If the system has been off for:

- 30 minutes wait 2 hours from Power ON till scanning
- 1 hour: wait 4 hours from Power ON till scanning
- 3 or more hours: wait 15 hours from Power ON till scanning

For the CT tube and detectors:

Depending on the system model, you must allow the CT tube and detectors to return to their operating temperature before resuming operation. Wait the following duration after power ON, then perform tube warm-up and resume normal work:

-  If the CT component was powered off for:
 - Less than 10 minutes wait 15 minutes
 - Over 10 minutes wait 60 minutes
-  Wait two hours.


2.6 B615 & D630 Powering ON from Cold

NOTE B615 D630

- This section is relevant for Brivo NM 615 and Discovery NM 630.
- If the system is in **Standby** mode, proceed to step 6.


1. Make sure there is no patient on the patient table.
2. If the **emergency OFF (EMO)** button was pressed, turn the button clockwise to reset it **A**.
3. If a UPS unit is used, switch on the unit as described in the manufacturer's instructions.

Final - 30-May-2012

4. At the gantry, if the **main circuit breaker** **C** was turned off, lift it to the "ON" position.
5. Press the green button **B** to power up the gantry.
6. Switch on the operator console and wait until the login dialog is displayed.
7. Enter the NM login and password, then wait until application loads.
8. Press the green reset button **D** on the hand-held controller holder.
9. Perform the NM start-up test, following the instructions on the gantry display (see *NM600 NM Clinical Operation Guide*). ♦


2.7 O640 Powering ON from Cold

NOTE **O640**

- This section is relevant only for Optima NM/CT 640.
- If the system is in **Standby** mode, proceed to step 6.

1. Make sure there is no patient on the patient table.
2. If the system **emergency OFF (EMO)** button was pressed, turn the button clockwise to reset it **A**.
3. At the CT gantry, lift the main circuit breaker **B** to the ON position.
4. If a UPS unit is used, switch on the unit as described in the manufacturer's instructions.


NM 600 Series System Description and Safety Manual for Operators
5439352-1EN, Rev. 1, ©2012 GE Healthcare

2-11

Final - 30-May-2012

5. At the gantry, if the **main circuit breaker** ① was turned off, lift it to the "ON" position.
6. Press the green button ② to power up the NM gantry.
7. Switch on the operator console and wait until the login dialog is displayed.
8. Enter the login and password, then wait until application loads.
9. Press the green reset button ③ on the hand-held controller holder.
10. Perform the system start-up test, following the instructions on the gantry display (see *NM600 NM Clinical Operation Guide*). ♦


2.8 D670 Powering ON from Cold

NOTE **D670**

- This section is relevant only for Discovery NM/CT 670.
- If the system is in **Standby** mode, proceed to step 6.


1. Make sure there is no patient on the patient table.
2. If the **emergency OFF (EMO)** button was pressed, turn the button clockwise to reset it **(A)**.
3. At the **Main Disconnect & UPS Control panel (A1)**, press the green **<Start>** button **(B)**.
4. If a UPS unit is used, switch on the unit as described in the manufacturer's instructions.

Final - 30-May-2012

5. If the **NM main circuit breaker** **F** was turned off, lift it to the "ON" position.
6. Press the green button **D** to power up the NM gantry.
7. Switch the operator console to "ON" **G** and wait until the login dialog is displayed.
8. Enter the NM login and password, then wait until application loads.
9. Press any reset button (**<Reset>** **H**) on rear control panel, gantry-side green buttons **G**, or green button **E** on the hand-held controller holder).
10. Perform the NM start-up test, following the instructions on the gantry display (see *NM600 NM Clinical Operation Guide*).
11. Warm up the CT tube (see *NM600 Quality Control Guide*). ♦

Chapter 3: Emergency/Safety Mechanisms and Procedures

3.1 Emergency Mechanisms for Moving Components

Every operator must be familiar with the location and operation of all emergency mechanisms. It is recommended that the use of these mechanisms is practised and refreshed periodically so that the operator is prepared to react immediately in a real emergency.

Depending on the type of emergency, you may need to use one or more of the following emergency mechanisms:

Emergency Mechanism	Used in these cases	Functionality
System Emergency OFF (EMO), p.3-2	Fire, flood, earthquake or other major catastrophe	Disconnects the power from the entire system
Emergency Stop Buttons, p.3-4	Patient-related emergencies	Stops all mechanical movement and scan processes
Pallet Release, p.3-7	Emergency egress when there is no power	Releases the table pallet mechanism in order to move the pallet manually
Detector Manual Release Mechanism, p.3-8	Emergency egress when there is no power	Enables manual movement of the detectors

Final - 30-May-2012

3.1.1 System Emergency OFF (EMO)

System emergency OFF buttons are intended for use in the event of a fire, flood, earthquake, or any other emergency.

This button shuts down the power to all system components.

GE Healthcare recommends placing a **System emergency OFF** button near the doorway of every scan and operator room.

However, the facility designer determines the quantity and locations of the **System emergency OFF** buttons.

Final - 30-May-2012

WARNING

Because the power supply to the system is shut down abruptly, the system may not have time to save data. This might corrupt system files or result in loss of patient data.

After using the **System Emergency OFF** button:

Figure 3-1: System Emergency OFF Button

1. Reset the **System emergency OFF** by turning the button **A** clockwise.
2. Follow the system-specific instructions, noting the appropriate detector or CT tube warmup warnings:
 - **B615** **D630** B615 & D630 Powering ON from Cold, p.2-9
 - **O640** O640 Powering ON from Cold, p.2-11
 - **D670** D670 Powering ON from Cold, p.2-13 ♦

Final - 30-May-2012

3.1.2 Emergency Stop Buttons

Several external emergency stop (e-stop) buttons are provided, depending on the system model/type (see Figure 3-2, p.3-5).

Press any of the **Emergency Stop** buttons in the event of a patient-related emergency or if the pallet, table or gantry starts to move unexpectedly.

When pressed, any **Emergency Stop** button halts all table and gantry motion.

NOTE

O640 **D670** In Optima NM/CT 640 and Discovery NM/CT 670, the button also simultaneously stops generation of x-ray radiation and turns off laser alignment lights.

CAUTION

- If the pallet is moving when the **Emergency Stop** button is pressed, the pallet may overrun by up to 10 mm.
- If you press the **Emergency Stop** button during a scan, the system stops data acquisition in progress, and attempts to save all data acquired prior to the emergency stop.

Note: * On some systems, this button may be located on the boom itself and not on the hand-held controller.



Figure 3-2: Location of Emergency Stop Buttons

NM 600 Series System Description and Safety Manual for Operators
5439352-1EN, Rev. 1, ©2012 GE Healthcare

Emergency/Safety Mechanisms and Procedures

To reset the system, press one of the following:

- Green button on the hand-held controller holder **A**
- **O640 Only**: Green reset button on O640 rear panel **B**
- **O640 Only**:
 - Green button on gantry sides **B**
 - Reset key on D670 Rear Control Panel, p.1-28 **C**

3.1 Emergency Mechanisms for Moving Components 3.1.2 Emergency Stop Buttons

3.1.3 Pallet Release

Pallet motion is controlled either automatically when the system moves to a pre-defined position, or by using the hand-held controller. In an emergency:

1. Press the enabler **A** on the pallet release handle.
All automatic system motion is stopped.
2. Manually move the pallet backward or forward. ♦


Figure 3-3: Location of Pallet Release Handle

3.1.4 Detector Manual Release Mechanism

The detectors include a manual mechanism that enables the operator to move the detectors away from the patient in these cases:

- There is no power to the system, for example in case of a power failure or if the operator has pressed System Emergency OFF (EMO), p.3-2
- Failure of the hand-held controller to move the detectors

IMPORTANT

Ensure that the ratchet for the detector manual release mechanism is stored close to the system in an easily accessible and clearly demarcated location.

For each detector separately:

1. At the detector that restricts the patient insert the ratchet **C** into the detector release opening **B**.
2. Wind the handle (Detector1: clockwise; Detector2: counter-clockwise, see arrow **A**) until the patient is no longer restricted by the detector.
3. If necessary, repeat steps 1 and 2 for the other detector. ♦


Figure 3-4: Location of Detector Manual Release Mechanism

3.2 Emergency Mechanisms for Moving Components

The system is fitted with a number of preventive safety devices, including:

- Safety Pressure Sensitive Devices (PSD), p.3-10
- ENABLE Lever on Hand-Held Controller, p.3-12
- Collision Warning Mechanism, p.3-13

Final - 30-May-2012

3.2.1 Safety Pressure Sensitive Devices (PSD)

The face of each collimator is fitted with two **Pressure Sensitive Devices (PSD)** that prevent inadvertent collision between the patient and the face of the detector during normal imaging procedures. On a pinhole collimator, the PSD is inside the cone cover.

WARNING

The PSD does not provide protection against collisions with the front, sides or top of the detector housing. On a pinhole collimator, the lateral surface of the cone is not protected by the PSD.

Always monitor patient and system for potential collisions during positioning and scanning. If necessary, use the Emergency Stop Buttons, p.3-4 to prevent collisions.

Final - 30-May-2012

When pressure is applied to the PSD (for example by the patient's body), the PSD is activated and all movement of patient table, detectors and gantry stops immediately (without disconnecting system power). In this case, all keys on the hand-held controller are disabled, except for the **<OUT>**, **<UP>** and **<ENABLE>** keys on the gantry display hand-held controller, which are required to move the detectors away from the patient.

NOTE

- Each PSD must be tested daily and after every collimator exchange (see *NM 600 NM Clinical Operation Guide*).
- **0640 Only:** PSD activation does not stop CT gantry rotation.




Figure 3-5: Location of Pressure Sensitive Devices

3.2.2 Reset After PSD Activation

1. On the gantry display hand-held controller, press **<OUT>** **B** + **<ENABLE>** **A** or **<UP>** **C** + **<ENABLE>** **A** until there is no longer contact between the obstruction and the PSD.
2. Press **<ENABLE>** **A** and follow the instructions on the gantry display to resume normal operation.

NOTE

If the obstruction that activated the PSD does not enable moving the detectors with the hand-held controller, first move the detectors out manually. ♦

Final - 30-May-2012

3.2.3 ENABLE Lever on Hand-Held Controller

An **<ENABLE>** lever is provided on the right hand side of the hand-held controller.

To prevent movement of the gantry if one of the hand-held controller keys is accidentally pressed, the system responds to movement and acquisition controls only when the **<ENABLE>** lever is depressed (see System Hand-Held Controller, p.1-21).


Figure 3-6: Location of Enable Lever on Hand-Held Controller

3.2.4 Collision Warning Mechanism

The software system controlling gantry motion has information on the position of all the motorized moving parts. It also stores a model of the fixed and moving parts of the system. By correlating these pieces of information, the system is able to detect the eminence of certain types of collisions that might occur between system moving parts and other moving or fixed parts of the system. In addition, the Body Contour Mechanism can be also used in certain circumstances to prevent collision of detectors surface with the patients.

While initially stopping the motion which may cause the collision, the mechanisms described above still leaves to the operator the possibility to continue the motion, despite the warning issued.

In addition, not all the types of collisions can be foreseen by the system.

Therefore, the operator continues to be the sole responsible for collision prevention and the guarantor of patient and system physical integrity.

Final - 30-May-2012

3.3 O640 & D670 X-ray Protection & Emergency Procedures

NOTE O640 D670

This section is relevant for Optima NM/CT 640 and Discovery NM/CT 670.

X-ray equipment, if not properly used, may cause injury. In order to prevent such injury, a number of mechanisms and indications have been incorporated into the system, as follows:

3.3.1 O640 Radiation Safety Control Mechanisms

NOTE O640

- This section is relevant only for Optima NM/CT 640.
- For the D670-related radiation control information, see D670 CT Operator Documents, p.xiv.

- X-rays are only produced when an acquisition protocol specifying x-ray transmission has been started.
- X-rays are generated only by a time-out controlled, dual-step mechanism:
 - To start x-ray:
 - At the hand-held controller, press **<GO>+<ENABLE>** to start a clinical hybrid scan.
 - At the console, click **[Start]** in the **Start X-ray** dialog to enable x-ray generation.
 - If the time-out is reached while system is waiting for **[Start]** to be pressed, **[Start]** is disabled. O640 time-out is 300 seconds
- If any gantry, detector or console problem is detected, x-ray generation is stopped immediately.
- Optional – door interlock mechanism enabling x-ray generation only when the Scan Room door is closed.

3.3.2 Radiation Indicators

- Visual x-ray indicators light up whenever there is danger of x-ray exposure (see Figure 3-7, p.3-16), including:

System	Indicator location
O640 D670	Flashing light on top of gantry*
O640 D670	Indication on gantry display
O640 Only	Flashing x-ray symbol on console screen, together with display of x-ray Tube kV, mA and the on-line updated exposure time
D670 Only	Console dedicated keyboard
D670 Only	Rear control panel

* O640 After system reset, during the automatic testing sequence, the radiation warning light turns on for a few seconds and then turns off. No x-rays are generated during the sequence.

- In accordance with site requirements and local regulations:
 - Leaded walls, windows or partitions to protect the operator from x-ray exposure
 - Safety and x-ray protection devices such as door closed sensors and door lights

Final - 30-May-2012


Figure 3-7: Location of X-ray ON Indicators

NM 600 Series System Description and Safety Manual for Operators
5439352-1EN, Rev. 1, ©2012 GE Healthcare

3.3.3 Emergency Patient Care During X-ray Exposure


WARNING

- While you are in the operator room during x-ray exposure, monitor the patient through the leaded window or using other equipment such as closed-circuit cameras to ensure the patient's safety.
- For detailed x-ray safety instructions, refer to the D670 CT Operator Documents, p.xiv.

- In an emergency, use the Emergency Stop Buttons, p.3-4 **C**.
- **D670 Only:** To stop a specific x-ray scan:
 - To abort the current x-ray scan and stop gantry/table movement, click **<STOP SCAN>** **A**.
 - To pause scanning after the current scan completes, click **<PAUSE>** **B**.

In both cases, click **[Resume]** in the **CT View/ Edit** screen to continue the scan.


Final - 30-May-2012

- **O640 Only:** To stop a specific x-ray scan:
 - At the console, click **[Stop]** .
 - To resume the stopped scan, click **[Resume]** , then click **[OK]**.


The system prompts you to press **<GO>** to move the table to the start position and start the CT acquisition.

3.4 Obstruction Free Area and Egress

When the site is planned, it takes into account all equipment required in the room. Before placing additional items in the room, for example cabinets, ensure that they will not make the room non-compliant with local regulations.

- An unobstructed egress route is required around the system, see red arrows in Figure 3-8, p.3-20.
- An obstruction free area is marked around the system, with yellow lines on the floor, see Figure 3-8, p.3-20.
If objects such as an external ECG trigger monitor are moved into this area for the purpose of an examination, pay special attention that they do not obstruct moving parts of the system. Move any object not being used for the current activity out of the obstruction free area.
- Do not place any objects under the edge of the patient table, as the table up/down movement may crush them and cause damage.

Final - 30-May-2012


Figure 3-8: Obstruction Free Area

NM 600 Series System Description and Safety Manual for Operators
5439352-1EN, Rev. 1, ©2012 GE Healthcare

3-20

NOTE

- An unobstructed egress route is required around the system. The red arrows indicate possible egress routes.
- **D670 Only:** The platform behind the CT gantry can be used to move round the system, but this must be done with care so as not to stumble or fall when crossing the platform.

Chapter 4: System Specifications

This section provides the main system features and specifications for the imaging devices described in this guide, according to the following major system components:

- Gantry
- Patient table
- Detectors
- Collimators
- Acquisition console
- Environment and power requirements

Information that does not refer to a specific system model is valid for all the models described. The values in this section represent typical measurements obtained on the system.

4.1 Gantry

B615 & D630

B615 / D630 Gantry Parameters	Specification
Length	136 cm (53.5")
Width	153 cm (60.2")
Height (including gantry display pole)	211 cm (83.1")
Diameter of gantry bore	70 cm (27.6")
Weight (with LEHR collimators)	B615 1577 kg (3477 lbs) D630 2052 kg (4524 lbs)
Emergency stop buttons	3

O640

O640 Gantry Parameters	Specification
Length	210 cm (82.7")
Width	201 cm (79.1")
Height (including gantry display pole)	211 cm (83.1")
Diameter of gantry bore	70 cm (27.6")
Weight (with LEHR collimators)	2952 kg (6495 lbs)
Emergency stop buttons	4

D670

D670 Gantry Parameters	Specification
Length	220 cm (86.6")
Width	218 cm (85.8")
Height (including gantry display pole)	211 cm (83.1")
Diameter of gantry bore	70 cm (27.6")
Weight (with LEHR collimators)	4345 kg (9581 lbs)
Emergency stop buttons	5

Gantry Motion

NM Gantry Motion Parameters		Specification
Rotational (axial) motion speed		1/60-3.0 rpm automatically, 1.0 rpm or 3.0 rpm manually
Independent radial motion speed		50 or 75 cm/min (19.7 or 29.5"/min)
Rotation		540°
D630, O640, and D670 Detector orientations*	SPECT	■ 90° (L-shape) ■ 180° (opposing)
	Planar	■ 90° (L-shape) ■ 180° (opposing) ■ 0° (horizontal) ■ 180° (vertical)
Swivel Tilt range		0° to 180° around detector axis
Caudal/Cephalic tilt range		Each detector can swivel independently ±45° (available only on hospital bed/gurney)
Detector scan radius (with LEHR collimators)		Minimum: 10 cm (3.9") Maximum: 35 cm (13.8")
Lateral motion speed		50 or 75 cm/min (19.7 or 29.5"/min)

* B615 Not applicable for B615

4.2 Patient Table

Table Parameters	Value (min)	Value (max)
Length	280 cm (110.2")	
Width	40 cm (15.7") (table top)	■ B615, D630 – 53 cm (20.9") ■ O640, D670 – 61 cm (24")
Height	59 cm (23.2")	■ B615, D630 – 98 cm (38.6") ■ O640, D670 – 100 cm (39.4")
Weight	B615 , D630 – 390 kg (860 lbs) O640 , D670 – 565 kg (1246 lbs)	
Patient load		227 kg (500 lbs)
Scan range in Whole-Body (WB) scan		■ 200 cm (78.7") on pallet ■ 230 cm (90.5") with table extender
Attenuation		■ All models: <11% for 140 keV gamma rays ■ O640, D670: Also <14.5% for 120 kVp x-ray.
Vertical travel time (full range)	34 sec (slow)	26 sec (fast)
Horizontal speed (manual, using Hand-held controller)	50 mm/sec (slow)	Up to 100 mm/sec (fast)
O640 , D670 Pallet positioning accuracy per 21 CFR 1020.33 (i)	±1 mm	

4.3 Detectors

D630, O640 and D670 feature two externally mounted large FOV rectangular detectors. B615 has a single externally mounted large FOV rectangular detector.

Each detector includes:

- 59 PMT's per detector: 53 circular with 76 mm (3") diameter and 6 circular with 38 mm (1.5") diameter
- UFOV: 54 cm x 40 cm $\pm$ 0.5 cm (21.25" x 15.75")
- Crystal thickness: 9.5 mm (3/8")
- Energy range: 40-620 keV

Detector Performance

Measured per NEMA NU-1, 2007 using Tc^{99m} where applicable (measurements with $\pm 4\%$ tolerance).

Specifications	Parameter	Range	Data (for 3/8" crystal)
Intrinsic spatial resolution	UFOV FWHM	=<	3.8 mm
	UFOV FWTM	=<	6.9 mm
	CFOV FWHM	=<	3.7 mm
	CFOV FWTM	=<	6.8 mm
Intrinsic energy resolution (Tc^{99m} at 20 kcps)	UFOV FWHM	=<	9.5%

Final - 30-May-2012

Specifications		Parameter	Range	Data (for 3/8" crystal)
Intrinsic flood field uniformity	Integral (max deviation)	UFOV	=<	3.6%
		CFOV	=<	3.0%
	Differential (max deviation)	UFOV	=<	2.3%
		CFOV	=<	2.1%
Intrinsic spatial linearity	Absolute	UFOV	=<	0.40%
		CFOV	=<	0.40%
	Differential	UFOV	=<	0.20%
		CFOV	=<	0.20%
Intrinsic count rate*	Maximum count rate		=>	460 kcps
	Count rate at 20% loss at 20% window		=>	250 kcps
	Maximum count rate at 20% window		=>	400 kcps
Multiple window spatial registration			=<	0.5 mm
SPECT resolution with LEHR collimator (FWHM, mm) with scatter	Central		=<	9.9 mm
	Radial		=<	9.9 mm
	Tangential		=<	7.5 mm

* Not measured per NEMA

Evolution for Bone (optional) SPECT Performance

Specifications	Parameter	Range	Full time	½ time
Evolution for bone SPECT reconstructed spatial resolution with scatter* (LEHR collimators)	Central	=<	6.4 mm	7.0 mm
	Radial	=<	5.7 mm	6.0 mm
	Tangential	=<	5.1 mm	5.4 mm

* Reconstructed without post-filtering

Final - 30-May-2012

4.4 Collimators

Parallel Hole Collimators

- Low Energy High Resolution (LEHR)
- Extended Low Energy General Purpose (ELEGP)
- Medium Energy General Purpose (MEGP)
- (Optional) High Energy General Purpose (HEGP)

Other Collimators

- Pinhole collimator (optional for all models)
- Fan-Beam collimator (optional)
 - B61 – not available

Low energy collimators (LEHR and ELEGP) are the only collimators operational in 90° dual-detector geometry

B615 – configuration geometry not available.

Each collimator set includes a cart on which two collimators are mounted, with the exception of the Pinhole collimator.

B61 – each collimator cart holds a single collimator.

Final - 30-May-2012

Parameter	LEHR	ELEGP	MEGP	HEGP
Field of View (FOV) (mm)	540 x 400	540 x 400	540 x 400	540 x 400
Calculated penetration	0.3% (Tc ^{99m})	0.3% (I ¹²³) 2.3% (Kr ⁸¹)	2% (Ga ⁶⁷)	2% (I ¹³¹)
3/8" System sensitivity (cpm/μCi) @ 100 mm (per detector)*	160 (Tc ^{99m})	320 Tc ^{99m} 224 (I ¹²³)	144 (Ga ⁶⁷)	97 (I ¹³¹)
3/8" System sensitivity (cps/MBq) @ 100 mm (per detector)*	72 (Tc ^{99m})	144 (Tc ^{99m}) 101 (I ¹²³)	65 (Ga ⁶⁷)	43 (I ¹³¹)
3/8" System resolution FWHM (mm) @ 100 mm [†]	7.4	10.3	9.4	12.0
Type of hole	Hex	Hex	Hex	Hex
Hole diameter (mm)	1.5	2.5	3.0	4.0
Septal thickness (mm)	0.2	0.4	1.05	1.8
Hole length (mm)	35	40	58	66
Weight (in kg/lb)	60/132	62/136	103/227	131/289

* Measured with 20% window using relevant isotope for each collimator, with ±10% tolerance.

† Measured at 100 mm distance from collimator face with ±4% tolerance.

Final - 30-May-2012

Pinhole Collimator (Option)

Parameter	Specification		
Field Of View (FOV) (mm)	200 (diameter)*		
Weight (kg/lb)	98/216		
Insert hole diameter	3/8" System sensitivity (cpm/μCi) @ 100 mm (per detector)†	3/8" System sensitivity (cps/MBq) @ 100 mm (per detector)†	3/8" System resolution FWHM (mm) @ 100 mm‡
2 mm	43	19	3.8
4.45 mm	200	90	6.5
8 mm	570	258	11.4

- * Diameter of the circular pinhole cone base (on the detector crystal side).
- † Sensitivity measured with Co⁵⁷ (point source @100 mm from insert center, with PSD cover, 20% window) extrapolated to Tc^{99m}, and with ±10% tolerance.
- ‡ Resolution measured with ±4% tolerance.

Fan-Beam Collimator (Option for D630, O640, D670)

NOTE **D630** **O640** **D670**
This section is relevant for Discovery NM 630 and Discovery NM/CT 670 and Optima NM/CT 640.

Parameter*	Specification
Field Of View (FOV) (mm) [†]	540 x 250
Calculated penetration	0.1% (Tc ^{99m})
Type of hole	Hex
Hole diameter (mm)	1.5
Septal thickness (mm)	0.2
Hole length (mm)	40
Volume Sensitivity 3/8" [‡]	39,960 (cpm/μCi/cm ²) 18,000 (cps/MBq/cm ²)
SPECT resolution 3/8" ^{**}	Central 8.0 mm Peripheral Radial 8.5 mm Peripheral Tangential 6.7 mm
Weight (kg/lb) 1pcs	64/140

* Recommended brain scan radius is 14–16 cm
Focal distance from collimator surface is 350 mm
Effective collimator thickness is 57 mm
† Collimator Field Of View
‡ Sensitivity measured on two detectors, Tc^{99m} with 15% window, with ±10% tolerance.
** SPECT Resolution measured with ±6% tolerance.

4.5 Acquisition Station

- The PC acquisition station is driven by Linux™ Operating System.
- Console screen display matrix 1280 x 1024
- On-the-fly digital corrections for energy, linearity, uniformity, isotope decay, center of rotation (COR) for SPECT.

4.6 Acquisition Modes

Acquisition can be performed after selecting an appropriate gantry mode:

- 180° and 90° detector geometry
 – not available
- On-pallet imaging or a variety of off-pallet scanning orientations (including upright seated or standing patients), appropriate for Static and Dynamic imaging.

Static

- Data acquired by any of the detectors or both in any of the following image matrix sizes: 32x32, 64x64, 128x128, 256x256, and 512x512

Dynamic

- Data acquired by any of the detectors or both in any of the following image matrix sizes: 32x32, 64x64, 128x128, 256x256
- Maximum frame rate: 50 frames/sec

Whole-Body Scan

- Image matrix: 256 x 1024
- Infrared-based real-time Automatic Body Contouring (ABC) system
B615 – optional
- Sequential multi-spot and continuous scanning modes
- Maximum scan range 200 cm (78.7") with additional 30 cm (12") on table extender
- Minimum WB continuous scan speed: 5 cm/min (1.97"/min)
- Maximum WB continuous scan speed: 30 cm/min (11.8"/min)
- Simultaneous anterior/posterior dual view scan
B603 – sequential dual view scan

Multi-gated

- Buffered data for on-line irregular beat rejection
- Equi-time and equi-phase gating modes, with data rebinning into 8, 16, 24, 32 frames per R-R cycle
- User selectable R-R window

SPECT

- Image matrix: 64x64, 128x128, 256x256
- Infrared-based on-line Automatic Body Contouring system
B615 – optional
- Step and Shoot mode, with angular view resolution
B615 (2° to 6°, 9°, 10°, 12°, 15°, 18°, 20°, 30°, 36°, 45°, 60°, 90°)
D630, **O640**, **D670** – including steps of 3°, 5°, 6°, 9°, 10°, 15°, 18°, 20°, 30°, 45°, 60°
- Continuous mode with 3°, 4°, 5°, or 6° angular sampling
- Scan on pallet for 200 cm (78.7") range, or on extender for 30 cm (12") range
- Consecutive multiple SPECT scans combined with table translation enable Volumetrix pasting of up to 6 FOV's (covering most body lengths)
- Either 90° or 180° detector geometries for cardiac SPECT
B615 – not available

Gated SPECT

- Buffered data for on-line irregular beat rejection
- Up to 24 gated bins per cardiac cycle
- User-selectable R-R window
- An extra bin sums all data, both accepted and rejected, securing non-gated study completion in case of exceptionally irregular heart-rate

Dynamic SPECT (Option)

D630, O640, & D670 Only

- Single-SPECT component scan time:
 - Minimum – 1 min/360°
 - Maximum – 15 min/180°
- Up to 40 consecutive single-SPECT scans in H-mode
- Supported collimators: LEHR, ELEGP, MEGP, HEGP and fan-beam
- Continuous mode acquisition, angular sampling: 3°, 4°, 5°, 6°
- Supports pan and zoom
- Supports multi-isotope and multi-peak energy settings

Final - 30-May-2012

4.6.1 Hybrid SPECT/CT Imaging

O640, & **D670** Only

- **D670** CT images are acquired on a standard BrightSpeed CT scanner integrated in the system. Refer to the CT user guides for technical specifications and additional information about the CT sub-system.
- **O640** CT images are acquired on a low-mA fast-rotating-speed CT device dedicated to SPECT attenuation correction and anatomical localization (120 and 140* kV, 10-15 or 10-30* mA, 1* or 2 seconds/rotation, 4x2.5* or 2x5 mm slice thickness).
See *Optima NM/CT 640 CT Technical Reference Manual* for technical specifications and additional information.
Note: * indicates optional items described below
- Hybrid imaging involves parallel-beam collimator SPECT using various scan modes: step and shoot, continuous, gated, dynamic, multiple FOV.
- **D670** D670 CT imaging is performed as:
 - Hybrid acquisition, as described in the *D670 Hybrid Clinical Operation Guide*
 - Standalone CT acquisition, as described in the D670 CT Operator Documents, p.xiv
- **O640** O640 CT imaging is performed only during hybrid scans, as described in the *O640 Hybrid Clinical Operation Guide*

Final - 30-May-2012

4.7 NM 600 Series Configurations

Each of the models in the NM600 family includes various optional components, that are not mentioned explicitly in this manual.

Following are the O640 configurations:

Table 4-1: Optima NM/CT 640 Configurations*

Parameter Name		Supported values per configuration	
		Excel	Custom
Rotation time (seconds)		1, 2	2
Voltage (kV)		120, 140	120
Current (mA) (depends on voltage and rotation time)		■ 10-30 for 120 kV and 1 s ■ 10-15 for 120 kV and 2 s ■ 10-20 for 140 kV and 1 s ■ 10 for 140 kV and 2 s	10-15
Helical pitch		0.75, 1.25, 1.75	1.25
Reconstructed slice thickness (mm)		2.5, 5	5
Hybrid scan length on pallet (cm)		120	80
Clinical hybrid scan flows	Set-up on emission (NM First)	Yes	Yes
	Set-up on scout (CT First)	Yes	No
Active table ruler		Yes	No
Breath-hold support		Yes	No

* The set of factory-recommended hybrid protocols are adapted to the parameters values specific to each configuration.

4.8 Networking Specifications

DICOM 3.0 Compatible Output and Automatic Data Transfer

4.8.1 Purpose of Network Connection

The NM Scanner is intended to be connected to a network in order to support the following functionality:

- DICOM services to push images to other DICOM-compliant machines
- DICOM services to confirm that images have been permanently stored on a DICOM-compliant machine
- DICOM services to obtain DICOM modality worklist information from an information system computer of a remote hospital or radiology department
- DICOM services to verify that the remote DICOM system is connected properly to the NM scanner device
- Remote service diagnostics and configuration tools (InSite)
- Printing to a compatible network printer.

All of the above features are optional on the NM Scanner. This implies that in a nominal environment some or all of the above connections may coexist in any given time frame.

Final - 30-May-2012

4.8.2 Network Interface Technical Specifications

Connection Name: Site Network Port	
Physical network connection type:	IEEE 802.3-1998 1000/100/10 BaseT Ethernet
Speeds and duplex modes supported:	10 Mbps, 100 Mbps, and 1 Gbps half and full duplex Auto-negotiate
Default IP address (from factory):	IP Address: by DHCP Subnet Mask: 255.255.255.0 Gateway: empty
IP addressing:	IPv4 static or DHCP
Qos support:	n/a

4.8.3 Network Information Flow Specifications

Flow Name: DICOM Image Push	
Network Connection on device	Site network port
Usage Type/Function/Purpose	Send a DICOM image or a set of images to a network device
Licensed/optional/required	Optional
Communication Partner Device/ IP Address/Network	Any network device supporting the DICOM application layer protocol(s) listed below
Middle Layer Protocols	TCP/IP
Application Layer Protocol and Encoding	NM Image Storage SC Image Storage
Ports (default)	4006
Traffic characterization and Bandwidth Requirements	On demand, local user initiated. The bandwidth is dependent on the local site.
Latency max	N/A

Final - 30-May-2012

Flow Name: DICOM Storage Commit	
Network Connection on device	Site network port
Usage Type/Function/Purpose	Used to confirm that local DICOM images have been permanently stored on a remote DICOM device
Licensed/optional/required	optional
Communication Partner Device/ IP Address/Network	Any network device supporting the DICOM application layer protocol(s) listed below
Middle Layer Protocols	TCP/IP
Application Layer Protocol and Encoding	Storage Commitment Push Model SOP Class
Ports	4006
Traffic characterization and Bandwidth Requirements	On demand, local user initiated. The bandwidth is dependent on the local site.
Latency max	N/A

Flow Name: DICOM Modality Worklist Information	
Network Connection on device	Site network port
Usage Type/Function/Purpose	Transfer patient information from HIS/RIS system to NM scanner
Licensed/optional/required	Optional
Communication Partner Device/ IP Address/Network	Any network device supporting the DICOM application layer protocol(s) listed below
Middle Layer Protocols	TCP/IP
Application Layer Protocol and Encoding	Basic Modality Worklist Information Model - FIND SOP Class
Ports (default)	4006
Traffic characterization and Bandwidth Requirements	On demand, local user initiated. The bandwidth is dependent on the local site.
Latency max	N/A

Final - 30-May-2012

Flow Name: DICOM Echo	
Network Connection on device	Site network port
Usage Type/Function/Purpose	To verify the remote DICOM system is connected properly to the NM scanner device
Licensed/optional/required	Optional
Communication Partner Device/ IP Address/Network	Any network device supporting the DICOM application layer protocol(s) listed below
Middle Layer Protocols	TCP/IP
Application Layer Protocol and Encoding	Support for C-ECHO is mandatory for all Application Entities which accept associations including NM scanners
Ports (default)	4006
Traffic characterization and Bandwidth Requirements	On demand, local user initiated. The bandwidth is dependent on the local site.
Latency max	n/a

Flow Name: Insite HTTP	
Network Connection on device	Site network port
Usage Type/Function/Purpose	Remote diagnostics and configuration of the scanner
Licensed/optional/required	Optional
Communication Partner Device/ IP Address/Network	Insite OLC
Middle Layer Protocols	TCP/IP
Application Layer Protocol and Encoding	HTTP
Ports (default)	80
Traffic characterization and Bandwidth Requirements	On demand, local user initiated. The bandwidth is dependent on the local site.
Latency max	n/a

Final - 30-May-2012

Flow Name: Insite RFB	
Network Connection on device	Site network port
Usage Type/Function/Purpose	Remote diagnostics and configuration of the scanner, training and support
Licensed/optional/required	Optional
Communication Partner Device/ IP Address/Network	Insite OLC
Middle Layer Protocols	TCP/IP
Application Layer Protocol and Encoding	RFB - Remote FrameBuffer
Ports (default)	5900,5800
Traffic characterization and Bandwidth Requirements	On demand, local user initiated. The bandwidth is dependent on the local site. Bandwidth can range from several KB/s to several MB/s depending on usage.
Latency max	n/a

Flow Name: Insite Telnet	
Network Connection on device	Site network port
Usage Type/Function/Purpose	Remote diagnostics and configuration of the scanner
Licensed/optional/required	Optional
Communication Partner Device/ IP Address/Network	Insite OLC
Middle Layer Protocols	TCP/IP
Application Layer Protocol and Encoding	Telnet
Ports (default)	23
Traffic characterization and Bandwidth Requirements	On demand, local user initiated. The bandwidth is dependent on the local site.
Latency max	n/a

Final - 30-May-2012

Flow Name: Insite SSH	
Network Connection on device	Site network port
Usage Type/Function/Purpose	Remote diagnostics of the scanner
Licensed/optional/required	Optional
Communication Partner Device/ IP Address/Network	ISD
Middle Layer Protocols	TCP/IP
Application Layer Protocol and Encoding	SSH
Ports (default)	22
Traffic characterization and Bandwidth Requirements	On demand, remote user initiated. The bandwidth is dependent on the local site.
Latency max	n/a

Flow Name: Insite FTP	
Network Connection on device	Site network port
Usage Type/Function/Purpose	Remote diagnostics and configuration of the scanner
Licensed/optional/required	Optional
Communication Partner Device/ IP Address/Network	Insite OLC
Middle Layer Protocols	TCP/IP
Application Layer Protocol and Encoding	FTP
Ports (default)	21
Traffic characterization and Bandwidth Requirements	On demand, local user initiated. The bandwidth is dependent on the local site.
Latency max	n/a

Final - 30-May-2012

Flow Name: DNS	
Network Connection on device	Site network port
Usage Type/Function/Purpose	Network services
Licensed/optional/required	Optional
Communication Partner Device/ IP Address/Network	Site DNS Server
Middle Layer Protocols	UDP
Application Layer Protocol and Encoding	DNS
Ports	53
Traffic characterization and Bandwidth Requirements	Initiated by scanner network software sporadically as required to translate a remote host name into an IP address
Latency max	n/a

Flow Name: DHCP	
Network Connection on device	Site network port
Usage Type/Function/Purpose	Network services
Licensed/optional/required	Optional
Communication Partner Device/ IP Address/Network	Site DHCP Server
Middle Layer Protocols	UDP
Application Layer Protocol and Encoding	DHCP
Ports	67, 68
Traffic characterization and Bandwidth Requirements	Initiated by scanner network software sporadically as required to obtain IP address
Latency max	n/a

Final - 30-May-2012

Flow Name: Ping	
Network Connection on device	Site network port
Usage Type/Function/Purpose	Network troubleshooting
Licensed/optional/required	Required
Communication Partner Device/ IP Address/Network	Any network device from which Ping diagnostic may be executed
Middle Layer Protocols	ICMP
Application Layer Protocol and Encoding	n/a
Ports (default)	n/a
Traffic characterization and Bandwidth Requirements	On demand, local initiated bidirectional.
Latency max	n/a

Flow Name: Rexec	
Network Connection on device	Site network port
Usage Type/Function/Purpose	RPC is used for remote administration. Its function is to forward application function call to a remote system over the network.
Licensed/optional/required	optional
Communication Partner Device/ IP Address/Network	Insite OLC
Middle Layer Protocols	TCP/IP
Application Layer Protocol and Encoding	Rexec
Ports	512
Traffic characterization and Bandwidth Requirements	On demand, local/remote user initiated.
Latency max	n/a

4.9 Environment and Power Requirements

Power conditioning is incorporated into the primary power supply of the system.
The system can operate on line voltage based on local conditions and codes.

Power Requirements

B615 D630	208-240 Vac nominal, 1 PH + Neutral 50/60 Hz. Power consumption: 6 kVA maximum, 2.5 kVA average.
O640	380-480 Vac nominal, 3-phase + Neutral, 50/60 Hz. Power consumption: 19 kVA maximum, 10 kVA average
D670	380-480 Vac nominal, 3-phase Delta or Wye, 50/60 Hz. Power consumption: 90 kVA maximum, 22 kVA average.

Environmental Requirements

Parameter	Maximum	Minimum	Recommended	Maximum Change Rate
Temperature	26°C (79°F)	18°C (64°F)	22°C (72°F)	3°C/hr (5°F/hr)
Humidity	60% non-condensing relative	30% non-condensing relative		5%/hr

4.10 System GE Approved Accessories

Use only GE approved equipment with the system.

NOTE

Contact GE Healthcare for an up-to-date list of approved accessories.

With each use, check all accessories for damage and remove them from service if damaged.

**WARNING**

- Do not connect accessories that are not approved as part of the system. Do not use accessories from other modalities.
- Hospital items such as arm boards and catheter bag holders are not secured to the gantry and may interfere with the gantry if not positioned properly.

Final - 30-May-2012

4.10.1 UPS Devices

System	Details	
B615 D630 O640	Manufacturer	Eaton
	Model	EXRTG11000-KL6UMS, 86511
	Rated Voltage	200-240 V
	Power	11 kVA/8 kW
	Standards	FCC class A
D670 (partial UPS)	Manufacturer	Eaton
	Model	9155-10GE
	Rated Voltage	120 Y-208 V
	Power	10 kVA/10 kW
	Standards	FCC class A

Final - 30-May-2012

4.10.2 External ECG Triggers

System	Description	Model	Manufacturer
B615 D630 O640	Cardiac ECG Trigger	IVY 3000 IVY 7600	IVY BioMedical Systems, Inc.
D670	Cardiac ECG Trigger (Nuclear only)	IVY 3000 IVY 7600	
	Cardiac ECG Trigger with recorder and USB data storage (Nuclear and CT)	IVY 3150 IVY 7800	

4.10.3 Table Accessories

Name	Manufacturer/Model
Head Holder	GE H3100NW
Pallet Extender	GE H3100PG
R-made Arm Rest	GE E8500NA
Foam Arm Rest	GE E8500NB
Contoured Leg Rest	GE E8500NC

Chapter 5: China-Specific Information

5.1 China-Specific RoHS Compatibility

This symbol indicates the product contains hazardous materials in excess of the limits established by the Chinese standard SJ/T11363-2006

"Requirements for Concentration Limits for Certain Hazardous Substances in Electronic Information Products".

The number in the symbol is the Environment-friendly Use Period (EFUP), which indicates the period during which the toxic or hazardous substances or elements contained in electronic information products will not leak or mutate under normal operating conditions so that the use of such electronic information products will not result in any severe environmental pollution, any bodily injury or damage to any assets. The unit of the period is **Year**.

In order to maintain the declared EFUP, the product shall be operated normally according to the instructions and environmental conditions as defined in the product manual, and periodic maintenance schedules specified in Product Maintenance Procedures shall be followed strictly.

Consumables or certain parts may have their own label with an EFUP value less than the product. Periodic replacement of those consumable or parts to maintain the declared EFUP shall be done in accordance with the Product Maintenance Procedures.

This product must not be disposed of as unsorted municipal waste, and must be collected separately and handled properly after decommissioning.

Table 5-1: Table of Hazardous Substances’ Name and Concentration

Component name	Hazardous substances’ name					
	(Pb)	(Hg)	(Cd)	(Cr(VI))	(PBB)	(PBDE)
Operator Console	X	O	O	X	O	O
Gantry	X	O	X	X	X	X
LCD Monitor	O	X	O	O	O	O
ECG Cardiac Trigger	X	O	O	X	X	X
Power Distribution Unit	X	O	X	X	X	X
Patient Table	X	O	X	X	X	X
O: Indicates that toxic or hazardous substance contained in all of the homogeneous materials for this part is below the limit requirement in SJ/T11363-2006 X: Indicates that this toxic or hazardous substance contained in at least one of the homogeneous materials used for this part above the limit requirement in SJ/T11363-2006 ■ Data listed in the table represents the best information available at the time of publication ■ Applications of hazardous substances in this medical device are required to achieve its intended clinical uses, and/or to provide better protection to human beings and/or to the environment, due to lack of reasonably (economically or technically) available substances.						

Final - 30-May-2012

This product consists of devices that may contain mercury, which must be recycled or disposed of in accordance with local, state, or federal laws. (Within this system, the backlight lamps in the monitor display contain mercury).

The Nuclear collimator contains the following potentially hazardous substances:

Lead: Lead salts are toxic and their ingestion may cause serious problems.
The manipulation/handling of lead is subject to regulations.

**WARNING**

- Do not discard the collimators among industrial waste or domestic garbage.
- A damaged collimator should not be dispatched through the national postal service.

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