INSTRUCTION for use

1. **PURPOSE:** Ophthalmic spectacle frames made of metal, plastic, combined produced by “Ferdinand Menrad GmbH + Co. KG”, meant for holding the lenses in a correct fixed position over the eyes.

2. **MAIN TECHNICAL CHARACTERISTICS AND CONSTRUCTION**

Frames, depending on design features and manufacturing technology, are divided into types indicated in Table 1.

Table 1.

Designation of frame type	Name of frame type	Type and designation of temples
FP	Plastic optical frames	Plastic rigid PR; Plastic elastic PE.
FM	Metal optical frames	Metal rigid PR; Metal elastic PE.
FC	Combined optical frames	Metal rigid MR; Metal elastic PE; Plastic rigid PR; Plastic elastic PE

Figure. 1. Metal optical frames (FM)


Figure. 2. Plastic optical frames (FP)

Figure. 3. Combined optical frames (FC)

Ophthalmic spectacle frames consist of a front, temples, nose pads, hinges. Hinges may be simple and spring-based.


Figure. 4. Structure of frame

Ophthalmic spectacle frames are characterized by 3 main geometric parameters: width (a), bridge width (b) and temple length to the curve (L1).

The sizes used for spectacle frames by "Ferdinand Menrad GmbH + Co. KG" are given in the table 2:

Ophthalmic spectacle frames		Distance between the vertical tangents to the facet groove of the rim (a). (mm))	Minimum distance between the vertical tangents to the bottom of the facet groove of the rim (b). (mm))	The whole temple length (L) (mm)	Temple length from the hinge axis to the curve (L1) (mm)	Temple mass without lenses (g)
Ophthalmic spectacle frames (metal)	men	$54 \pm 0,2$	$16 \pm 0,2$	$135 \pm 1,0$	$95 \pm 1,0$	$22 \pm 1,0$
	men	$58 \pm 0,2$	$18 \pm 0,2$	$135 \pm 1,0$	$95 \pm 1,0$	$21 \pm 1,0$
	women	$56 \pm 0,2$	$17 \pm 0,2$	$135 \pm 1,0$	$95 \pm 1,0$	$22 \pm 1,0$
	women	$55 \pm 0,2$	$17 \pm 0,2$	$135 \pm 1,0$	$95 \pm 1,0$	$21 \pm 1,0$
	unisex	$55 \pm 0,2$	$16 \pm 0,2$	$135 \pm 1,0$	$95 \pm 1,0$	$22 \pm 1,0$
	unisex	$54 \pm 0,2$	$18 \pm 0,2$	$140 \pm 1,0$	$100 \pm 1,0$	$22 \pm 1,0$
Ophthalmic spectacle frames (plastic)	men	$57 \pm 0,2$	$15 \pm 0,2$	$140 \pm 1,0$	$100 \pm 1,0$	$21 \pm 1,0$
	men	$59 \pm 0,2$	$17 \pm 0,2$	$145 \pm 1,0$	$105 \pm 1,0$	$20 \pm 1,0$
	women	$55 \pm 0,2$	$15 \pm 0,2$	$155 \pm 1,0$	$115 \pm 1,0$	$21 \pm 1,0$
	women	$59 \pm 0,2$	$15 \pm 0,2$	$155 \pm 1,0$	$115 \pm 1,0$	$20 \pm 1,0$
	unisex	$59 \pm 0,2$	$17 \pm 0,2$	$140 \pm 1,0$	$100 \pm 1,0$	$21 \pm 1,0$
	unisex	$59 \pm 0,2$	$16 \pm 0,2$	$140 \pm 1,0$	$100 \pm 1,0$	$21 \pm 1,0$
Ophthalmic spectacle frames (combined)	men	$59 \pm 0,2$	$14 \pm 0,2$	$155 \pm 1,0$	$115 \pm 1,0$	$22 \pm 1,0$
	men	$59 \pm 0,2$	$15 \pm 0,2$	$150 \pm 1,0$	$110 \pm 1,0$	$21 \pm 1,0$
	women	$55 \pm 0,2$	$13 \pm 0,2$	$135 \pm 1,0$	$95 \pm 1,0$	$22 \pm 1,0$
	women	$59 \pm 0,2$	$19 \pm 0,2$	$135 \pm 1,0$	$95 \pm 1,0$	$22 \pm 1,0$
	unisex	$59 \pm 0,2$	$18 \pm 0,2$	$140 \pm 1,0$	$100 \pm 1,0$	$21 \pm 1,0$
	unisex	$59 \pm 0,2$	$15 \pm 0,2$	$155 \pm 1,0$	$115 \pm 1,0$	$20 \pm 1,0$

On the surface of the frames there may be no sharp edges or burrs, shells, scratches, bubbles, sagging solder visible to the naked eye.

Temples should rotate on the axis smoothly, without jamming.

Decorative trim (if any) is from nickel silver (N18 CuNi18Zn20).

Ophthalmic spectacle frames correspond to the European norms (standards): DIN ISO EN 12870, DIN EN 12472, DIN EN 16128.

Ophthalmic spectacle frames are chosen according to the anthropometric parameters of the patient. It is necessary to consult a medical worker to decide whether to apply this method of vision correction.

3. INDICATIONS FOR APPLICATION:

Correction of defects of the optics of eyes:

- myopia
- presbyopia
- astigmatismus
- hypermetropia (small type of clinical refraction)
- anisometropia (difference of refraction on eyes)
- keratoconus and others.

4. The following may be considered CONTRAINDICATIONS:
- individual intolerance to sight correction by means of spectacles;
 - impossibility of sight correction by means of spectacles.

5. MATERIALS:

Frames are made of materials that do not cause irritation, allergic or toxic reactions during their exploitation upon contact with human skin.

Plastic frames are produced of cellulose acetate (cellulose acetate polymer (69-75 %) and diethyl phthalate (25-31 %)), hinges and screws made of nickel silver (N18 CuNi18Zn20).

Metal and combined frames are produced of nickelsilver (N18 CuNi18Zn20), stainless steel X5CrNi18-10 (chromium 10-30 %, nickel 12-14 %, molybdenum 2-3 %), hinges and screws made of nickel silver (N18 CuNi18Zn20).

Tips of temples are produced of cellulose acetate (cellulose acetate polymer (69-75 %) and diethyl phthalate (25-31 %)). Nose pads are produced of PVC (phthalate free).

The transport insert made of acrylic Polycryl Clear.

Materials comply with the requirements of biological safety.

6. MARKING:

The right temple indicates the manufacture's name, conformity mark and country of manufacture.


Figure 5. Frame marking. Right temple.

Left temple shall have factory marking with indication of:

model number (16043), technical size (54-18-140), colour or coating code (6100).


Figure 6. Frame marking. Left temple.

Each frame is packed into an individual consumer plastic packaging. Each consumer packaging (plastic bag) is marked with a conformity mark.


Figure 7. Individual plastic packaging.

Plastic packaging of each frame or the label placed into plastic packaging should contain:

- nominal value of basic sizes;
- model number;
- model colour;
- bar code
- name of product;
- manufacturer;
- date of manufacture;
- expiration date.


Figure 8. Label attached to plastic packaging.

Frames in plastic packaging should be compactly placed into group packaging box. Group packaging box should contain frames of the same type with the same basic sizes, in quantity of 8 - 10 pcs.


Figure 9. Group packaging box.

Each group packaging box should have a sticker of the manufacturer on it, which should contain following inscriptions:

- model number;
- model colour;
- nominal value of basic sizes;
- quantity of frames

Dimensions of the packaging box, in cm (length x width x height): 23 x 15 x 5,8.

Maximum load on the group packaging box is 10 kg.

Transport packaging is made only by agreement with the consumer.

Frames packed in transport packaging should be resistant to transport shaking with a maximum acceleration of 30 m/s² and frequency of 2 to 3 Hz.

Marking of transport packaging with application of the handling mark "Caution, fragile".

Set of delivery: transport inserts, frame.

7. PREPARING THE ITEM FOR USE: The spectacle frame should not lose its original shape after fitting the lenses. The rims should be symmetrical. No dents, scratches or gaps between spectacle lenses and rims, nor any other assembling defects are allowed.

Spectacle lenses fixed in the frame should not move or fall out of the rims during normal spectacle use.

8. TERMS OF STORAGE AND USE:

- It is not recommended to leave the rim in direct sunlight at high temperatures (eg, car control panel);
- Need to avoid the fall of the rim on hard surfaces, and mechanical loads that can damage the rim;
- Storage frames need to use a case or pouch, is not permitted the wearing of eyeglasses in the pockets of clothing or bags in order to avoid scratches on the cover and mechanical damage;

- To put on and take off glasses with both hands is necessary - this prevents their deformation;

- Need to avoid wearing eyeglasses during sleep;

- Are not allowed to apply to the rim significant mechanical forces;

- Caution should be exercised when doing sports and the production of works involving sudden movements of the head, the slopes, jumps and bumps.

- Do not apply significant mechanical efforts to the frame;

- Keep the frame in a case or any other package to protect it from mechanical damage;

- Full average term of use is no less than 3 years..

- Terms of storage: in heated storages under the environmental temperature from + 5°C to + 50°C, relative air moisture 80%, atmospheric pressure 700-1060 hPa

- No impact of aggressive substances (acids, alkali) and organic solvents and their fumes is allowed during storage.

- Storage life: plastic and combined frames – 5 years, metal frames – 10 years.

9. TERMS OF USE

Ophthalmic spectacle frames are used in all climate zones at any time of the day.

10. PRECAUTIONS FOR USE

While using ophthalmic spectacle frames it is necessary to exercise caution practising sports and activities involving sudden head movements, bends, jumps and bumps.

11. POSSIBLE SIDE EFFECTS

Allergic reaction in case of individual intolerance of frames material is possible.

12. INFORMATION ON THE POSSIBLE IMPACT OF THE USE ON ABILITY TO DRIVE VEHICLES, MACHINERY

Drivers with vision problems may drive vehicles and machinery using spectacle frames (in the manufacture of frames for corrective glasses).

13. STERILITY

Unsterile.

14. MAINTENANCE AND REPAIRS

Ophthalmic spectacle frames should be repaired in special laboratories licensed for production, maintenance and repairs.

15. WARRANTY OBLIGATIONS

Storage life: plastic and combined frames – 5 years, metal frames – 10 years.

Full storage life – not less than 3 years with temple bends not more than 30 000.

16. TRANSPORTATION

Ophthalmic spectacle frames are transported by all kinds of transport in covered vehicles in accordance with the rules of transportation of goods applicable to this kind of transport.

Temperature from - 45°C to + 50°C

Relative air moisture up to 98%

Atmospheric pressure 700-1060 hPa

17. DISPOSAL

Spectacle frames do not contain harmful substances and components that pose a threat to human health and the environment during and after the end of life and disposal. Ophthalmic spectacle frames should be disposed of as household waste.

Wastes of production and parts with expired year of storage are classified as class A medical waste.

18. CERTIFICATES

The manufacturer uses its own quality system that is subject to verification and updating on an annual basis.

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Legal address: 111024, Moscow, Andronovskoe highway, 26, build. 6, off. 11B
Phone number 8 (495) 660-01-25

Hermann Müller-Menrad
Managing Director
Hermann Müller-Menrad GmbH + Co. KG

Signature

Date

8-12-2017

File No. 12657 / 2017

I, Dr. Joachim Schervier, hereby certify that this document was signed in my presence by

Mr. Hermann Müller-Menrad ,
born on 19.03.1961,
80804 München, Mottlstraße 29, Germany,
identified by his identity card.

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