
OPERATIONAL DOCUMENTATION FOR MEDICAL DEVICE

INSTRUCTIONS FOR USE


IN VITRO REAGENTS KIT FOR RAPID DETECTION OF PROCALCITONIN
(RAMP PROCALCITONIN)

Name: Uma Gaur

Position: Associate Director, QA/RA

Signed before me on this 3 day of
October, 2022 by Uma Gaur

Signature: 

Nikolay N. Eskin, Notary Public in and
for the Province of British Columbia

My Commission is Permanent

Nikolay N. Eskin, Notary Public

201-223 West Broadway

Vancouver, BC V5Y 1P2

Tel: 604-416-0154

My Commission is Permanent Response Biomedical Corporation

Last version of the IFU: 09.2022

TEST NAME

In vitro reagents kit for rapid detection of procalcitonin (RAMP Procalcitonin) (hereafter test kit, RAMP Procalcitonin)

One kit includes materials to run 25 tests.

WARNING!

For in vitro diagnostic use only

Failure to follow RAMP Procalcitonin test procedure may result in invalid and/or erroneous results. Read the entire Instructions For Use prior to performing test.

INTENDED USE

RAMP Procalcitonin assay is a qualitative immunochromatographic test indicated for use as in vitro diagnostic product used with a RAMP Reader or RAMP 200 instrument to measure levels of prohormone procalcitonin (PCT) in human EDTA anti-coagulated whole blood. RAMP Procalcitonin assay is intended for professional use only and is to be used in conjunction with other laboratory findings and clinical assessments as an aid in the risk assessment of progression to severe sepsis and septic shock.

POTENTIAL USERS

Physicians and diagnostic laboratory personnel who have been trained and admitted to work with RAMP Procalcitonin test in accordance with the established procedure.

INDICATIONS AND CONTRA-INDICATIONS FOR USE, SIDE EFFECTS

Indications for use: RAMP Procalcitonin assay is intended to be used in conjunction with other laboratory findings and clinical assessments as an aid in the risk assessment of progression to severe sepsis and septic shock.

Contra-indications for use, side effects: not found.

SUMMARY

In the absence of infection, the prohormone procalcitonin (PCT) is expressed in thyroid C-cells and subsequently undergoes post-translational cleavages to yield mature calcitonin [1]. The

circulating plasma concentration of PCT in healthy individuals is very low [2,3]. Conversely, microbial infection triggers a substantial increase in the expression of PCT in non-neuroendocrine cells throughout the body which are unable to perform the necessary post-translational cleavages resulting in increased serum PCT concentrations [1].

Serum concentrations of PCT increase proportionally to the severity of infection, ranging from slightly elevated in infections with minor systemic inflammatory response to very high values in cases of severe sepsis and septic shock [2-4]. Sepsis is a common and frequently fatal condition that is the result of the body's excessive inflammatory response to an infection. Despite advances in critical care medicine, sepsis and septic shock are leading causes of death in intensive care units [4-6].

The speed at which appropriate therapy is provided in sepsis cases impacts patient outcomes [6,7]. The traditional method of sepsis diagnosis however, relies on culturing blood, urine or cerebrospinal fluid which takes upwards of 24 hours, lacks sensitivity and specificity and frequently does not correlate to the manifestation of clinical symptoms [2,3]. Alternatively over 20 years of research suggests PCT is an excellent biochemical marker for sepsis as it is highly specific for bacterial infections and correlates strongly with their extent and severity [4,8-12].

Research suggests that a PCT concentration <0.5 ng/mL can be considered normal while concentrations greater than 2 ng/mL are highly suggestive of severe sepsis [3,13,14]. The RAMP Procalcitonin assay should not be used in the diagnosis or prognosis of conditions other than sepsis such as non-specific complaints, dementia, malaria, tuberculosis, syncope, urinary tract infections and post-cardiac transplant graft failure [2,15].

TEST EXPLANATION END PRINCIPLE

The RAMP Procalcitonin assay is a quantitative immunochromatographic test for the determination of procalcitonin levels in EDTA whole blood. Mixed EDTA whole blood is applied into the sample well of the test cartridge. The red blood cells are retained in the sample pad, and the separated plasma migrates along the strip. Fluorescent-dyed particles coated with anti-procalcitonin antibodies bind to procalcitonin, if present in the sample. As the sample migrates along the strip, procalcitonin-bound particles are captured at the detection zone and excess fluorescent-dyed particles are captured at the control zone.

The RAMP instrument then measures the amount of fluorescence emitted by the complexes bound at the detection zone and at the control zone. Using a ratio between the two fluorescence values, a quantitative reading is calculated. For further information on the use of the instrument, refer to the RAMP Reader and RAMP 200 User Manuals.

STABILITY, STORAGE, TRANSPORTATION

Store RAMP Procalcitonin test kit at +2 °C – +8 °C. Unopened kit is stable until the indicated expiry date. The storage at +15 °C – +25 °C for 14 days is permissible.

The opening of the kit box does not influence on the expiry date. Once the pouch with test cartridge and assay tip opened, they must be used or discarded within 60 minutes.

Do not freeze. Do not use the kit beyond the stated expiry date or in the case of storage conditions violation.

The RAMP Procalcitonin test kit is to be transported by all kinds of covered vehicles in accordance with corresponding rules under storage conditions.

If transportation temperature conditions were violated the kits cannot be used.

RAMP Procalcitonin test kit reagents

The RAMP test kit contains all the reagents necessary for the quantification of procalcitonin in EDTA whole blood.

MATERIALS PROVIDED

Each RAMP Procalcitonin kit contains enough materials to run 25 tests, including:

- Test cartridge and test tip in individual pouch - 25 pcs.
- Sample buffer; 0,15 ml/vial - 25 vials.
- Transfer device.
- Lot card.
- IFU.

Content	Quantity
1 Test Cartridge and 1 Assay Tip in foil pouch, with desiccant. The Test Cartridge contains monoclonal antibodies to procalcitonin (active component) immobilized on the nitrocellulose membrane The Assay Tip contains latex fluorescent-dyed particles coated with anti-procalcitonin antibodies and sodium azide (conservative)	25
Sample Diluent, 0,15 ml in a vial Contains animal protein, surfactant and ProClin300/ProClin950 (preservative)	25
Transfer Device* for sample intake, dilution, and transfer into the sample well	1

Lot Card in anti-static pack; contains lot specific information to be inserted into the instrument	1
Instructions for use**	1

*- the volume is indicated on the device

** - the Instructions for Use (hard-copy form) may be presented separately

MATERIALS REQUIRED (BUT NOT PROVIDED)

- **RAMP Reader System** instrument, with accessories, **Response Biomedical Corporation**, registered in Russian Federation, Registration Certificate No P3H 2014/2188, or
- **RAMP® 200 System** including **RAMP® 200 Control Module** and **RAMP® 200 Test Module**, with accessories, **Response Biomedical Corporation**, registered in Russian Federation, Registration Certificate No ΦC3 2010/06545.
- Specimen collection tubes: EDTA (Venous Whole Blood)

LOT CARD CALIBRATION

Each RAMP test kit includes a lot card that is individually packaged in an anti-static pouch. The lot card provides information specific to the kit test cartridge lot, including lot number, expiration date, and standard curve information. For further details on loading lot-specific information, see the RAMP instrument User Manual. No additional calibration beyond insertion of the lot card is necessary. This operation is required only once per test kit lot.

For each new lot, remove the lot card from its pouch and insert it into the lot card slot on the instrument. Once the lot card has been uploaded, return to its pouch and do not discard. Avoid touching the contacts at the end of the lot card.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use.
- For professional use only.
- For use by qualified personnel per local, state, or Federal regulations or accrediting agency requirements.
- A medical device is intended for single use only. Do not reuse.
- Read the entire Instructions for Use (IFU) prior to use. Directions should be read and followed carefully, or invalid or erroneous results may occur.
- Diagnostic and prognostic uses of the **RAMP Procalcitonin** assay in non-septic diseases (e.g. non-specific complaints, dementia, malaria, tuberculosis, syncope, urinary tract infections and

post-cardiac transplant graft failure) or non-infectious conditions (e.g. burn trauma, extensive surgery, heatstroke and pancreatitis) have not been tested.

- Do not interchange or mix components of different RAMP tests, lots or components from other manufacturers.
- Do not use the kit or any kit component beyond the stated expiry date.
- Do not use any visibly damaged components.
- Do not insert a cartridge on which blood or any other fluid is spilled, into the instrument.
- Exercise standard precautions required for handling all laboratory reagents and patient samples.

Use individual protection equipment

- The device contains material of animal origin and should be handled as a potential biohazard.
- The sample buffer provided contains ProClin®, a potential skin sensitizer. Avoid spilling or splashing reagents containing ProClin on skin or clothing. In case of contact, thoroughly flush with water.

SYSTEM QUALITY CONTROL

The RAMP instrument has error checking and self-diagnostic functions (Internal Quality Control (IQC)) that assure system integrity. These include algorithms and measurements used to confirm acceptable operator technique, sample handling, and test performance. Frequency of IQC may be programmed at desired intervals.

Valid results are displayed only after all performance requirements have been met.

Each RAMP assay has built-in controls. Test cartridges have a control zone that is scanned as part of the test protocol to ensure proper sample flow.

Control limits for each lot of test cartridges are established during the manufacturing process and are incorporated in the test-specific lot parameters. If a control result does not meet specifications, the sample result is not reported and a message is displayed.

TEST RUN MESSAGES

When the RAMP instrument is unable to continue a specific task it will emit an audio alarm and display a message. Refer to the RAMP User Manual 'Troubleshooting Guide' section for a full description of all messages. If repeated tests give unexpected or inconsistent results, contact the authorized representative of Response Biomedical Corporation in Russian Federation.

SAMPLE COLLECTION AND PREPARATION

- Use ONLY EDTA Whole Blood (K₂EDTA tubes are recommended).
- Other sample types and anticoagulants have not been evaluated. Do not use samples that have been frozen.

- Avoid short draws and blood samples that show gross hemolysis as these may interfere with the test and cause erroneous results. If this occurs, another blood sample should be obtained and tested.
- Testing must be completed within 2 hours of phlebotomy. However, if this is not possible, the EDTA whole blood can be stored for up to 2 days at +2 to +8°C. If stored, allow blood samples to equilibrate at +18 to +25°C for at least 15 minutes prior to use.

PROCEDURE

Prior to sample preparation allow all components to come to room temperature (at +18 to +25°C) for at least 15 minutes.

- Keep the test cartridge and test tip in the sealed foil pouch until ready for use. Once opened, test cartridges and test tips must be used or discarded within 60 minutes.
 - The test cartridge, test tip, and buffer vial should be discarded after a single-use. Do not reuse.
1. Prepare RAMP instrument for test cartridge. Refer to the appropriate RAMP User Manual for detailed instructions on Starting a Test.
 2. Ensure that the EDTA whole blood or sample is well mixed by gentle inversion.
 3. Uncap the buffer vial and place upright on a clean, dry level surface, or in a holder.
 4. Open a test pouch and remove the test cartridge and tip. Place the test cartridge on a clean, level surface. Firmly attach the test tip to the supplied transfer device.
 5. Before inserting the test tip into the sample, fully depress the transfer device plunger.
 6. Insert the tip into sample and fully release plunger. The test tip should fill with 75 µL of blood.
 7. Immediately transfer the filled test tip into the buffer vial close to, but not touching, the bottom.
 8. Mix sample slowly by fully pressing and releasing the plunger 10 times; while keeping the tip submerged in the buffer for optimal mixing and to minimize air bubbles.
 9. Once mixing is complete, draw 75 µL of sample into the test tip by releasing the plunger one final time and immediately dispense liquid into the sample well of the test cartridge. Small droplets may remain in the tip; this is expected.
 10. Immediately insert the test cartridge fully into the instrument and press until firm resistance is felt.
 11. The instrument will draw the cartridge in and test development will begin.
 12. The instrument will analyze the cartridge and report the result in approximately 15 minutes.
 13. Record the result, if required. For additional information on printing and/or uploading results, please refer to the User Manual.

14. Remove the used test cartridge and discard all used test components according to local biohazard procedures (see "Utilisation"). DO NOT reuse.

For additional information on the general operation and troubleshooting of the instrument, please refer to the RAMP User Manual.

LIMITATIONS

- Factors such as technical or procedural errors or the presence of substances in blood specimens other than those that have been evaluated (see Interference section of this IFU), may interfere with the RAMP test and cause erroneous results.
- As with any immunoassay, patient specimens may contain heterophilic antibodies that may result in either falsely elevated or depressed results. Presence of these antibodies may be due to elevated levels of rheumatoid factor, treatment with mouse monoclonal antibodies for diagnostic or therapeutic purposes, or other undetermined factors. The RAMP test has been formulated to reduce the effects of heterophilic antibodies, but complete elimination of heterophilic interference from all samples cannot be guaranteed.

EXPECTED VALUES

The RAMP Procalcitonin Expected Values (reference range) study was conducted at one clinical site and included 125 apparently healthy individuals (59 women, 66 men). The normal range for the RAMP Procalcitonin assay (determined as the 95th percentile of these results) includes values ≤ 0.36 ng/mL. No statistically significant bias due to gender, age or ethnicity was observed in this study. Each laboratory should investigate the transferability of the expected values to its own patient population and, if necessary, determine its own reference range.

MEASUREMENT RANGE

0.20 to 200.0 ng/mL.

PCT levels in excess of 200 ng/mL are reported as greater than ($>$) 200.0 ng/mL, values below 0.20 ng/mL are reported as less than ($<$) 0.20 ng/mL.

DETECTION LIMIT

To determine the Limit of Blank (LoB), 60 replicates of a blank material were tested. The LoB, estimated as the non-parametric determination of the 95th percentile of all 60 replicates, was determined to be 0.20 ng/mL.

To determine the Limit of Detection (LoD), 15 replicates each of 4 low level EDTA whole blood samples were tested over 4 runs. The LoD was calculated to be 0.36 ng/mL.

TRUENESS

Has not been determined because of absence of International standard sample.

PRECISION

Repeatability and precision estimates were obtained using three levels of frozen control materials (spiked pooled plasma samples) and three lots of RAMP Procalcitonin assay, tested in duplicate over 20 days, two runs per day for a total of 80 replicates per level. The precision of the RAMP Procalcitonin test is shown in the following table.

Mean Value (ng/mL)	Repeatability	Within-Laboratory Precision
	%CV	
0.58	16.0 – 18.7	15.8 – 19.3
2.72	8.0 – 9.0	7.7 – 9.3
92.7	6.6 – 8.5	7.2 – 9.4

Precision estimates for whole blood samples were obtained using spiked EDTA whole blood and one lot of RAMP Procalcitonin assays. Whole blood samples were evaluated at three levels across the reportable range of the RAMP Procalcitonin assay and were each tested in replicates of 10 in a single run. The precision performance of the RAMP Procalcitonin assay was found to be 14.1, 6.8 and 8.2 %CV at 0.75, 3.00 and 103.2 ng/mL respectively.

ANALYTICAL SPECIFICITY

No cross-reactivity was observed for calcitonin (up to 5 ng/mL), katacalcin (up to 10 ng/mL) or α -CGRP and β -CGRP (up to 30 ng/mL).

HOOK EFFECT

No high dose hook effect was observed for the RAMP Procalcitonin assay up to the highest level tested (~2000 ng/mL).

INTERFERENCE

None of the following factors has been found to influence the RAMP Procalcitonin assay

- Hemolysis (hemoglobin to 500 mg/dL)
- Bilirubinemia (conjugate bilirubin to 40 mg/dL; unconjugated bilirubin to 40 mg/dL)
- Hyperlipidemia (triglycerides to 3260 mg/dL)

- Human serum albumin (to 2.4 g/dL)

However it is recommended not to use samples that appear to be clearly hemolyzed, lipemic or icteric and to collect a new sample if possible.

11 commonly used pharmaceutical compounds were evaluated for potential interference in the **RAMP Procalcitonin** assay. Compounds were evaluated by spiking different concentrations of each individual potential interferent into EDTA whole blood with procalcitonin added to provide levels of ~0.50 and ~10 ng/mL. The pharmaceutical compounds listed in the following table were tested at the indicated concentration and found to have no interference with the **RAMP Procalcitonin** assay.

Tested Compound	Tested Concentration
Acetaminophen	180 mg/dL
Acetylsalicylic Acid	240 mg/dL
Ibuprofen	144 mg/dL
Cefotaxime	200 mg/dL
Dobutamine	25 µg/mL
Dopamine	30 mg/dL
Heparin Sodium Salt	28.8 U/mL
Noradrenalin	5 µg/mL
Vancomycin	1.8 mg/mL
Imipenem	100 mg/dL
Furosemide	40 mg/dL

TECHNICAL SERVICE AND REPAIR

RAMP Procalcitonin test kit is not the subject of technical service and repair

WASTE UTILIZATION

Remove the used test cartridge from the reader. All the wastes (reagents, sample material, package) as well as the non-used expired reagents must be utilized in accordance with SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for the maintenance of the territories of urban and rural settlements, for water bodies, drinking water and drinking water supply, atmospheric air, soils, residential premises, operation of industrial and public premises, organization and conduct of sanitary and anti-epidemic (preventive) measures". Unused reagents and the wastes that did not contact with sample material belong to class A of medical wastes.

The samples used as well as every waste contacted or contaminated with biological materials belong to class B of medical wastes.

EXPIRY TIME

18 months.

The kits or its components should be used until the expiry date on the kit label. Do not use beyond it.

GUARANTEE

The Manufacturer's guarantee is valid if the test is used in accordance with the IFU. All warning and precautions should be taken into consideration. The guarantee is valid until the expiry date indicated on the kit label.

In the case of kit package or components found damaged, contact with Analytica Ltd. company, the authorized representative of Response Biomedical Corporation in Russian Federation.

MANUFACTURER AND DEVELOPER

Response Biomedical Corporation

1781 – 75th Avenue West, Vancouver,

British Columbia V6P 6P2 Canada,

Tel. +1 604 456 6010, fax +1 604 456 6066

AUTHORIZED REPRESENTATIVE IN RUSSIAN FEDERATION

ANALYTICA Ltd.

Serebryakova Proezd 2, Bld. 1,

Moscow, 129343 Russian Federation

For correspondence: 129343, Moscow, post box 93

Tel (495) 737-03-63, (800) 200-19-89

Fax: (495) 737-03-65

www.analytica.ru


info@analytica.ru

SYMBOLS USED


Symbol	Meaning
Box label	
	Batch Code
	Catalogue number
	In vitro diagnostic Medical Device
	Harmful, irritant
	CE Mark
	Manufacturer
	Temperature Limit
	Use-by Date
	Authorized Representative in European Community
	Consult Instructions for Use
	Contains Sufficient for 25 Tests
	Bar Code
Test Cartridge/Assay Tip Pouch Label (Foiled Pouch with Desiccant)	
	Do Not Reuse
	Batch Code
	Consult Instructions for Use
Test Cartridge Marking	
	Do Not Reuse
RAMP Sample Diluent Vial Label (Vial, 0,15 ml)	
	Harmful, irritant
	Batch Code
Lot Card Label (Anti-Static Pack)	
	Attention – Observe Precautions For Handling Electrostatic Sensitive Devices
Lot Card Symbols	
	Batch Code
	Use-by Date
	Caution

SHORT PROCEDURE

1. Prepare whole venous blood samples to be tested. Prepare your RAMP instrument for test cartridge.
2. Uncap the buffer vial and place upright on a clean, dry level surface, or in a holder.
3. Open a test pouch and remove the test cartridge and tip. Firmly attach the test tip to the supplied transfer device.
4. Fully depress the transfer device plunger, insert the tip into sample and fully release plunger.

5. Transfer the filled test tip into the buffer vial. Mix sample slowly by fully pressing and releasing the plunger 10 times.
6. Dispense 75 μ L of diluted sample into the sample well of the test cartridge.
7. Immediately insert the test cartridge fully into the instrument. It will analyze the cartridge and report the result in approximately 15 minutes.
8. Treat and utilize the wastes.

REFERENCES

- [1]. Muller et al. 2001. J Clin Endocrinol Metab. 86: 396-404
- [2]. Jin et al. 2010. Lab Med. 41(3):173-7
- [3]. McGee et al. 2009. Clin Lab News. 35(7)
- [4]. Muller et al. 2000. Crit Care Med. 28(4): 977-983
- [5]. Angus et al. 2001. Crit Care Med. 29(7): 1303-10
- [6]. Clec'h et al. 2004. Crit Care Med. 32(5) : 1166-1169
- [7]. Dellinger et al. 2013. Crit Care Med. 41(2) : 580-637
- [8]. Assicot et al. 1993. Lancet 341: 515-518
- [9]. Schuetz et al. 2011. BMC Medicine. 9:107-15
- [10]. Christ-Crain et al. 2004. Lancet 363: 600-607
- [11]. Uzzan et al. 2006. Crit Care Med. 34: 1996-2003
- [12]. Kutz et al. 2015. Critical Care 19: 74
- [13]. Becker et al. 2008. Crit Care Med. 36(3): 941-952
- [14]. Meisner 2014. Ann Lab Med. 34(4): 263-273
- [15]. Schneider et al. 2007. Pathology 39(4): 383-390

**ADDENDUM TO
INSTRUCTIONS FOR USE OF
IN VITRO REAGENTS KIT FOR RAPID DETECTION OF PROCALCITONIN
(RAMP PROCALCITONIN)**

This addendum contains the information concerning specific pathology, population/demographic aspects of the test, the list of regulation documents the test corresponds to, as well as the diagnostic sensitivity and specificity data obtained during the clinical evaluation of the test within the Russian registration procedure.

The diseases, pathologic conditions or risk factors associated with the test: sepsis, septic shock

Population and demography aspects associated with the test: The test has not any population or demography limitations.

Diagnostic sensitivity and specificity

The investigation was performed during the Russian test registration procedure and included the measurement of PCT concentration of 22 patients with sepsis (10 men, 12 women) and 125 healthy persons (blood donors). Cut-off level of 0.5 ng/ml was used.

Both sensitivity and specificity were 100%.

The list of International Standards and other regulation documents the test corresponds to:

Document	Name
Directive 98/79/EC	Directive 98/79/EC of the European Parliament and of the Council of 27 October 1998 on in vitro diagnostic medical devices.
ISO 13485	Medical devices - Quality management systems - Requirements for regulatory purposes.
ISO 14971	Medical devices – Application of risk management to medical devices.
ISO 15223-1	Medical devices - Symbols to be used with medical device labels, labeling and information to be supplied - Part 1: General requirements.
ISO 18113-1	In vitro diagnostic medical devices — Information supplied by the manufacturer (labelling) - Part 1: Terms, definitions and general requirements

Document	Name
ISO 18113-2	In vitro diagnostic medical devices — Information supplied by the manufacturer (labelling) - Part 2: In vitro diagnostic reagents for professional use.
ISO 23640	In vitro diagnostic medical devices - Evaluation of stability of in vitro diagnostic reagents.
SOR/98-282	Canadian Medical Devices Regulations.
CLSI C28-A3c	Clinical Laboratory Standards Institute. Defining, establishing, and verifying reference intervals in the clinical laboratory; Approved guideline - Third Edition.
CLSI EP5-A3	Clinical Laboratory Standards Institute. Evaluation of precision of quantitative measurement procedures; Approved guideline - Third Edition.
CLSI EP7-A2	Clinical Laboratory Standards Institute. Interference testing in clinical chemistry; Approved guideline – Second Edition.
CLSI EP17-A	Clinical Laboratory Standards Institute. Protocols for determination of limits of detection and limits of quantitation; Approved guideline.
CLSI EP25-A	Clinical Laboratory Standards Institute. Evaluation of stability of in vitro diagnostic reagents; Approved guideline.

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