State of Minnesota

SECRETARY OF STATE

APOSTILLE

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1. Country: United States of America

This Public Document: 3M Document

2. has been signed by Alison Leigh Chulla

3. acting in the capacity of Notary Public, State of Minnesota

4. bears the seal/stamp of Irma Vargas, Notary Public, State of Minnesota

CERTIFIED

5. at: St. Paul, Minnesota

6. Date: September 8, 2017

7. by: Secretary of State, State of Minnesota

8. File No: 9640422-6

9. Seal/Stamp:

10. Signature:

Steve Simon
Secretary of State

3M Health Care Business


3M Center
2510 Conway Ave., Bldg. 275-5W-06
St. Paul, MN 55144 U.S.A.
651 733 1110

To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M Health Care Business, St. Paul, MN, USA, I certify that the attached Instruction for Use (IFU) for 3M™ Tegaderm™ CHG Chlorhexidine Gluconate I.V. Securement Dressing and 3M™ Tegaderm™ CHG Chlorhexidine Gluconate I.V. Securement Dressing with PICC/CVC Securement Device in English is truthful and accurate.


Regina Feferman
International Regulatory Affairs Associate
3M Health Care


Date

State of Minnesota)
County of Dakota) ss

Subscribed and sworn to before me this 6th day of September, 2017.

Alison Leigh Chulla, Notary Signature
My Commission Expires: January 31, 2021

Instruction for Use

3M™ Tegaderm™ CHG Chlorhexidine Gluconate I.V. Securement Dressing

1. 3M™ Tegaderm™ CHG Chlorhexidine Gluconate I.V. Securement Dressing in the following sizes 7 cm x 8,5 cm, 8,5 cm x 11,5 cm, 10 cm x 12 cm, 10 cm x 15,5 cm, 12 cm x 12 cm.

2. 3M™ Tegaderm™ CHG Chlorhexidine Gluconate I.V. Securement Dressing with PICC/CVC Securement Device in the following sizes 8,5 cm x 11,5 cm, 10 cm x 15,5 cm.

REF	1657R
REF	1658R
REF	1659R
REF	1660R
REF	1877R-2100
REF	1879R-2100

Description:

3M Tegaderm™ CHG Chlorhexidine Gluconate I.V. Securement Dressing is used to cover and protect catheter sites and to secure devices to skin. It is available in a variety of shapes and sizes.

3M Tegaderm™ CHG dressing consists of a transparent adhesive dressing and an integrated gel pad containing 2% w/w chlorhexidine gluconate (CHG), a well-known antiseptic agent with broad spectrum antimicrobial and antifungal activity. The gel pad absorbs fluid. The transparent film provides an effective barrier against external contamination including fluids (waterproof), bacteria, viruses* and yeast, and protects the I.V. site.

In vitro testing (time kill and zone of inhibition) demonstrates that the Tegaderm CHG gel pad in the dressing has an antimicrobial effect against a variety of gram-positive and gram-negative bacteria, and yeast.

3M Tegaderm™ CHG dressing is transparent, allowing continual site observation, and is breathable, allowing good moisture vapor exchange.

*In vitro testing shows that the transparent film of Tegaderm™ CHG dressing provides a viral barrier from viruses 27 nm in diameter or larger while the dressing remains intact without leakage. The barrier to viruses is due to the physical properties of the dressing rather than ancillary properties of CHG

3M™ PICC/CVC Securement Device + Tegaderm™ I. V. Advanced Securement Dressing ("Securement System") is used to secure Peripherally Inserted Central Venous Catheters (PICCs) and short term Central Venous Catheters (CVCs) to the skin and cover and protect catheter insertion sites. The Securement System consists of a device and a dressing. The molded plastic device is integrated onto a breathable base with a gentle silicone adhesive. The soft cloth bordered transparent film dressing is made of a thin film backing with a non-latex adhesive. A large, notched, film-covered soft cloth tape strip is included in the system.

The Tegaderm™ I. V. Advanced dressing is breathable, allowing good moisture vapor exchange and has a transparent window, allowing continuous site observation. The transparent film provides an effective barrier against external contamination including fluids (waterproof), bacteria, viruses* and protects the I.V. site.

*In vitro testing shows that the transparent film of the Tegaderm CHG dressing provides a viral barrier from viruses 27 nm in diameter or larger while the dressing remains intact without leakage.

Indications:

3M Tegaderm™ CHG Chlorhexidine Gluconate I.V. Securement Dressing can be used to cover and protect catheter sites and to secure devices to skin. Common applications include central venous and arterial catheters, other intravascular catheters and percutaneous devices. Tegaderm™ CHG Dressing is intended to reduce skin colonization and catheter colonization and to suppress regrowth of microorganisms commonly related to bloodstream infections. Tegaderm CHG is intended to reduce catheter-related bloodstream infections (CRBSI) in patients with central venous or arterial catheters.

3M™ PICC/CVC Securement Device + Tegaderm™ I. V. Advanced Securement Dressing ("Securement System") is used to secure Peripherally Inserted Central Venous Catheters (PICCs) and Short Term Central Venous Catheters (CVCs) to the skin and to cover and protect catheter insertion sites. The device will accommodate the majority of single, double and triple lumen CVC catheters up to and including 12 French.

Warnings: CHG Chlorhexidine Gluconate I.V. Securement Dressing

- DO NOT USE TEGADERM CHG DRESSINGS ON PREMATURE INFANTS OR INFANTS YOUNGER THAN 2 MONTHS OF AGE. USE OF THIS PRODUCT ON PREMATURE INFANTS MAY RESULT IN HYPERSENSITIVITY REACTIONS OR NECROSIS OF THE SKIN.
- THE SAFETY AND EFFECTIVENESS OF TEGADERM™ CHG DRESSINGS HAS NOT BEEN EVALUATED IN CHILDREN UNDER 18 YEARS OF AGE. FOR EXTERNAL USE ONLY. DO NOT ALLOW THIS PRODUCT TO CONTACT EARS, EYES, MOUTH OR MUCOUS MEMBRANES.
- DO NOT USE THIS PRODUCT ON PATIENTS WITH KNOWN HYPERSENSITIVITY TO CHLORHEXIDINE GLUCONATE.
- THE USE OF CHLORHEXIDINE GLUCONATE CONTAINING PRODUCTS HAS BEEN REPORTED TO CAUSE IRRITATIONS, SENSITIZATION, AND GENERALIZED ALLERGIC REACTIONS. IF ALLERGIC REACTIONS OCCUR, DISCONTINUE USE IMMEDIATELY, AND IF SEVERE, CONTACT A PHYSICIAN.

Hypersensitivity reactions associated with topical use of chlorhexidine gluconate have been reported in several countries. The most serious reactions (including anaphylaxis) have occurred in patients treated with lubricants containing chlorhexidine gluconate, which were used during urinary tract procedures. Caution should be used when using chlorhexidine gluconate containing preparations, and the patient should be observed for the possibility of hypersensitivity reactions.

Warnings: CHG Chlorhexidine Gluconate I.V. Securement Dressing with PICC/CVC Securement Device

- DO NOT USE 3M™ PICC/CVC SECUREMENT DEVICE + TEGADERM™ I. V. ADVANCED SECUREMENT DRESSING ("SECUREMENT SYSTEM") WHERE LOSS OF ADHESION COULD OCCUR, SUCH AS WITH A CONFUSED PATIENT, DIAPHORETIC OR NON-ADHERENT SKIN, OR WHEN THE ACCESS DEVICE IS NOT MONITORED DAILY.
- FAILURE TO FOLLOW THE MANUFACTURER'S INSTRUCTIONS FOR USE MAY RESULT IN COMPLICATIONS INCLUDING INADEQUATE SECUREMENT AND SKIN

IRRITATION.

- THIS IS A SINGLE USE SECUREMENT SYSTEM. REUSE AND/OR REPACKAGING MAY CREATE A RISK OF PATIENT OR USER INFECTION, COMPROMISE THE STRUCTURAL INTEGRITY AND/OR ESSENTIAL MATERIAL AND DESIGN CHARACTERISTICS OF THE SECUREMENT SYSTEM, WHICH MAY LEAD TO DEVICE FAILURE.

Precautions:

- 3M Tegaderm CHG dressing should not be placed over infected wounds. It is not intended to be used as a treatment of percutaneous device-related infections.
- In the case of clinical wound infection, systemic antibacterials should be used if indicated. Any active bleeding at the insertion site should be stabilized before applying the dressing.
- Do not stretch the dressing during application. Mechanical skin trauma may result if the dressing is applied with tension.
- The skin should be clean, dry and free of detergent residue. Allow all preps and protectants to dry completely before applying the dressing to prevent skin irritation and to ensure good adhesion.
- Do not reuse. Reuse may result in compromising product integrity and lead to device failure.
- The 3M™ PICC/CVC Securement Device + Tegaderm™ I. V. Advanced Securement Dressing ("Securement System") may be used on an infected site, only when under the care of a health care professional.
- Observe universal blood and body fluid precautions and infection control procedures, during application and removal of the Securement System.
- Any active bleeding at the insertion site should be stabilized before applying the Securement System.
- The skin should be dry and free of detergent residue to prevent skin irritation and to ensure good adhesion. Allow all preps and protectants to dry completely before applying the Securement System to prevent skin irritation and to ensure good adhesion.
- Do not use on patients with known tape or adhesive allergies.
- Avoid Securement System contact with alcohol or acetone; both can weaken bonding of components and pad adherence.
- Minimize catheter manipulation during application and removal of the Securement System.
- Orient the device so arrows point toward the catheter insertion site.
- Securement System adherence and catheter position should be routinely inspected.
- Securement System should be replaced at least every 7 days.
- Do not stretch the Tegaderm™ I. V. Advanced dressing during application. Mechanical skin trauma may result if the dressing is applied with tension.
- Antimicrobial ointments containing polyethylene glycols may compromise the strength of the Tegaderm™ I. V. Advanced dressing.
- • Do not re-sterilize

Clinical Trial Results for Tegaderm CHG Chlorhexidine Gluconate I.V. Securement Dressing:

A randomized, controlled clinical trial consisting of 1879 subjects with 4163 central venous and arterial catheter insertion sites was conducted at 11 hospitals (1). Results showed that the

use of Tegaderm™ CHG resulted in a statistically significant 60% reduction in the incidence of catheter - related bloodstream infections (P=0.02). Study results also demonstrate a statistically significant reduction in skin colonization (P<0.001) and catheter colonization (P<0.0001) in the chlorhexidine vs. non-chlorhexidine group.

Variable	Non-chlorhexidine vs. chlorhexidine dressings (941 patients/2055 catheters vs. 938 patients/2108 catheters)
Catheter-related bloodstream infection	
Incidence densities (n per 1000 catheter-days)	1.3 vs. 0.5
Hazard ratio	0.402 [0.186 to 0.868], P=0.02
Catheter colonization	
Incidence densities (n per 1000 catheter-days)	10.9 vs. 4.3
Hazard ratio	0.412 [0.306 to 0.556], P<0.0001

(1) Timsit JF et al Randomized Controlled Trial of Chlorhexidine Dressing and Highly Adhesive Dressing for Preventing Catheter-Related Infections in Critically Ill Adults Am. J. Respir. Crit. Care Med. 2012; 186:1272-1278

Instructions for Use: Tegaderm CHG Chlorhexidine Gluconate I.V. Securement Dressing:

Failure to follow the manufacturer's instructions for use may result in complications including skin irritation and/or maceration.

Dressing Selection: Choose a dressing large enough to provide at least one inch margin of adherence on dry, healthy skin around the catheter site.

Site Preparation: Prepare the site according to institution protocol. Clipping of hair at the site may improve dressing adhesion. Shaving is not recommended. The skin should be clean, dry and free of detergent residue. Allow all preps and protectants to dry completely before applying the dressing to prevent skin irritation and to ensure good adhesion.

Any active bleeding at the insertion site should be stabilized before applying the dressing.

Application:

- 1.Open package and remove sterile dressing.
- 2.Peel liner from dressing, exposing adhesive surface.
- 3.Avoid stretching the dressing during application to reduce the risk of mechanical skin trauma.
- 4.Center the gel pad over the catheter insertion site. Apply firm pressure to entire dressing starting in the center to the outer frame edges to enhance adhesion.
- 5.Slowly remove frame while smoothing down transparent film dressing edges.
- 6.Smooth the transparent film dressing from the center towards the dressing edges, using firm pressure to enhance adhesion.
- 7.After dressing has been applied, apply the sterile tape strip(s) to further secure I.V. tubing or to stabilize catheter. Refer to figures on packaging.
- 8.Document dressing change information on label according to facility protocol. Remove label from frame and place on dressing.

Site Care:

1. The site should be observed daily for signs of infection or other complications. If infection is suspected, remove the dressing, inspect the site directly, and determine appropriate medical intervention. Infection may be signaled by fever, pain, redness, swelling, or unusual odor or discharge.
 2. Inspect the dressing daily and change the dressing as necessary, in accordance with facility protocol; dressing changes should occur at a minimum of every 7 days, per current Centers for Disease Control and Prevention (CDC) recommendations. Dressing changes may be needed more frequently with highly exudative sites. The Tegaderm™ CHG dressing should be changed as necessary:
 - If the dressing becomes loose, soiled or compromised in any way
 - If the site is obscured or no longer visible
 - If there is visible drainage outside the gel pad
 - If the dressing appears to be saturated or overly swollen
 - To test if the dressing is fully saturated, lightly press down on a corner of the gel pad with your finger. If the gel pad remains displaced once your finger is removed, the dressing should be changed.
- Note: Tegaderm CHG dressing is not designed to absorb large quantities of blood or fluid.

Removal:**Stabilize catheter during removal of the 3M Tegaderm CHG Dressing**

1. Remove documentation label and securement tape strip(s) from top of dressing.
2. Using a low and slow removal technique, start removing the dressing from where the catheter or tubing exits the dressing toward the catheter insertion site. Avoid skin trauma by peeling the dressing back, rather than pulling it up from the skin.
3. When the CHG gel pad is exposed, grasp a corner of the gel pad and the transparent film dressing between thumb and finger.
4. Apply sterile alcohol swabs or wipes, or sterile solutions (i.e., sterile water or normal saline) between gel pad and skin to facilitate removal of the gel pad dressing. If needed, a medical adhesive solvent can be used to help remove the dressing border.
5. Continue the low and slow removal method until the dressing is completely removed.

3M™ PICC/CVC Securement Device + Tegaderm™ I. V. Advanced Securement Dressing ("Securement System") Selection:


The device will accommodate the majority of single, double and triple lumen CVC catheters up to and including 12 French. Select the appropriate size Securement System where the Tegaderm™ I. V. Advanced dressing is large enough to provide at least one-inch margin of adherence on dry, healthy skin around the catheter site.

Site Preparation: Prepare the site according to institution protocol. Clipping of hair where the Securement System will be placed may improve adhesion of the Securement System. Shaving is not recommended. The skin should be clean, dry and free of detergent residue. Allow all preps and skin protectants to dry completely before applying the Securement System to prevent skin irritation and to ensure good adhesion. Any active bleeding at the insertion site should be stabilized before applying the Securement System.

Refer to Figure 1 for device application instructions.

Device Application:

1. Before placing device on the skin, orient the device with arrows pointing toward the insertion site as shown in **Figure 1a**.
2. Place single, double or triple lumen CVC catheter into device as shown in **Figure 1b** and weave lumen(s) under the single plastic arm.
3. Remove liner of attached tape strip and secure lumen(s) to the device base as shown in **Figure 1c**.
4. Position device on skin at the desired location. Pull and remove liner from one side of base to expose adhesive while holding the device in place. Pull and remove liner from the other side of base.
5. Apply pressure to the device base to establish good adhesion to the skin,
6. Apply Tegaderm™ I. V. Advanced dressing following the instructions below.

Tegaderm™ I. V. Advanced Securement Dressing Application:

1. Peel the liner from the dressing, exposing the adhesive surface.
2. Place the dressing so the transparent film covers the insertion site and the border of the dressing covers the single plastic arm on the device. **Figure 2**
Do not stretch dressing during application. Mechanical skin trauma may result if the dressing is applied with tension.
3. Apply firm pressure on dressing including edges to enhance adhesion to the skin.

4. Slowly remove the frame while smoothing down the dressing edges.
5. Smooth the dressing from the center towards the edges, using firm pressure to enhance adhesion.


Figure 2

Tape Strip Application

1. Remove liner from sterile tape strip.
2. Grasp non-adhesive tab of tape strip, and bend slightly with thumb. **Figure 3**
3. Lift catheter lumen(s) and apply the notch end of tape strip under catheter lumen(s) and over the dressing edge. Push tape strip notch forward, abutting up against catheter lumen(s). **Figure 4**
4. Apply pressure on tape strip to enhance adhesion.
5. Slowly remove the frame from the tape strip while smoothing down the edges.
6. On the label, document dressing change information according to facility protocol. Place the label on top of dressing over the catheter lumens. **Figure 5**

Figure 3

Figure 4

Figure 5

Site Care:

1. The site should be observed daily for signs of infection or other complications. If infection is suspected, remove the Tegaderm™ I. V. Advanced dressing, inspect the site directly, and determine appropriate medical intervention. Infection may be signaled by fever, pain, redness, swelling, or unusual odor or discharge.
2. Inspect the Securement System daily and change the system as necessary, in accordance with facility protocol. Securement System changes should occur at least every 7 days and may be needed more frequently with highly exudative sites or if integrity of the dressing is compromised.

Removal of Securement System

Minimize catheter manipulation during removal of the 3M™ PICC/CVC Securement Device +Tegaderm™ I. V. Advanced Securement Dressing ("Securement System").

1. Remove documentation label from top of dressing.
2. Remove the dressing tape strip by gently peeling the notched ends away from the insertion site.
3. Lift the catheter lumen(s) with one hand and place gloved index finger on base of device. Utilizing the low and slow removal technique, slowly start removing the dressing towards the insertion site.
4. When the catheter hub is exposed, move your gloved finger to secure the catheter hub and continue to remove the dressing until device is uncovered. Leave remainder of dressing in place over catheter insertion site. **Figure 6**
Avoid skin trauma by peeling the dressing back, rather than pulling it up from the skin.
5. Remove device tape strip from catheter lumen(s). **Figure 7**
6. Use gloved finger to stabilize the catheter hub and gently remove the catheter lumen(s) out from under the plastic arm of device.
7. Secure the catheter with one hand and use your other hand to remove the device from the patient's skin. **Figure 8**
8. Stabilize the catheter with one hand and carefully remove the remainder of the dressing over the insertion site using the low and slow removal technique.


Figure 6


Figure 7


Figure 8

Shelf Life and Storage Information:

For best results, store in a cool, dry place. 3M Tegaderm CHG Dressings are considered validated to a two (2) year shelf life. For shelf life, refer to the expiration date on the package. Sterility of the dressing is guaranteed unless individual package is damaged or open.

Information about repair and maintenance

This medical device isn't subject to maintenance and repair.

Utilization

Dispose of the device in accordance with applicable local/regional/national/international regulations. Medical devices in Russian Federation must be disposed of in the local public entities in accordance with SanPiN 2.1.7.2790-10 "Sanitary-epidemiological requirements for the handling of medical waste" to the waste class B

Transport conditions

Shelf life of the products is guaranteed as labeled, if stored according to the recommended storage conditions. During transportation, the products may experience short term temperature variations without affecting product safety and performance.

Transportation of medical device is carried out by all type of transport.

For further information outside the United States, contact your local 3M representative or contact us at www.3M.com and select your country.

Catalog #	Dressing Size	Average Amount of CHG per Dressing (mg based on gel pad size)	Quantity per box (primary)
1657R	8.5 cm x 11.5 cm (3- 1/2 X 4- 1/2 in)	45	25
1658R	10 cm X 12 cm (4 X 4-3/4 in)	45	25
1659R	10 cm x 15.5 cm (4 X 6- 1/8 in)	78	25
1660R	7 cm x 8.5 cm (2-3/4 x 3-3/8 in)	15	25
1877R-2100	8.5 cm x 11.5 cm (3- 1/2 x 4- 1/2 in)	45	20
1879R-2100	10 cm x 15.5 cm (4 x 6-1/8 in)	78	20

Explanation of Symbols

-  Not Made With Natural Rubber Latex
-  Caution, see instructions for use
-  Do not use if package is damaged or open
-  Do not reuse
-  Use by date
-  Batch code
-  Manufacturer
-  Date of manufacture
-  Sterilized using ethylene oxide
-  Do not resterilize

Authorised representative of the
 manufacturer in Russia:
 ZAO "3M Russia", Russia 121614,
 Moscow, Krylatskaya str. bldg.17/3
 T.+74957847474

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To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M Health Care Business, St. Paul, MN, USA, I certify that the attached Instruction for Use (IFU) for 3M™ Tegaderm™ I.V. Port Dressing in English is truthful and accurate.

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
3M Tegaderm CHG I.V. Port Dressing

Description:

3M Tegaderm CHG I.V. Port Dressing is used to cover and protect vascular access sites and to secure devices to the skin. The dressing consists of a bordered transparent film with an adhesive-free window. The border has multiple notch configurations and laminated tape strips to provide a better seal around catheters, tubing and other devices. The dressing is transparent, allowing continual site observation, and is breathable, allowing good moisture vapor exchange.

The Tegaderm dressing provides an effective barrier against external contamination, including fluids (waterproof), bacteria, viruses* and yeast.

**In vitro* testing shows that the transparent film provides a viral barrier from viruses 27 nm in diameter or larger while the dressing remains intact without leakage.

Indications:

3M Tegaderm CHG I.V. Port Dressing can be used to cover and protect vascular access sites and to secure devices to skin.

Warnings:

DO NOT PUNCTURE CHG GEL WITH NONCORING NEEDLE

DO NOT USE TEGADERM™ CHG DRESSINGS ON PREMATURE INFANTS OR INFANTS YOUNGER THAN 2 MONTHS OF AGE. USE OF THIS PRODUCT ON PREMATURE INFANTS MAY RESULT IN HYPERSENSITIVITY REACTIONS OR NECROSIS OF THE SKIN.

THE SAFETY AND EFFECTIVENESS OF TEGADERM™ CHG DRESSINGS HAS NOT BEEN EVALUATED IN CHILDREN UNDER 18 YEARS OF AGE. FOR EXTERNAL USE ONLY. DO NOT ALLOW THIS PRODUCT TO CONTACT EARS, EYES, MOUTH OR MUCOUS MEMBRANES.

DO NOT USE THIS PRODUCT ON PATIENTS WITH KNOWN HYPERSENSITIVITY TO CHLORHEXIDINE GLUCONATE. THE USE OF CHLORHEXIDINE GLUCONATE CONTAINING PRODUCTS HAS BEEN REPORTED TO CAUSE IRRITATIONS, SENSITIZATION, AND GENERALIZED ALLERGIC REACTIONS. IF ALLERGIC REACTIONS OCCUR, DISCONTINUE USE IMMEDIATELY, AND IF SEVERE, CONTACT A PHYSICIAN.

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
Precautions:

1. Any active bleeding at the insertion site should be stabilized before applying the dressing.
2. Do not stretch the dressing during application. Mechanical skin trauma may result if the dressing is applied with tension.
3. The skin should be dry and free of detergent residue to prevent skin irritation and to ensure good adhesion.

4. Antimicrobial ointments containing polyethylene glycols may compromise the strength of the cover dressing film.
5. Use care and cover the port dressing to protect from water when showering.

Instructions for Use: Refer to figures.

Site Preparation: Prepare the skin and insert noncoring needle according to facility protocol. Clipping of the hair may improve dressing adhesion. Shaving is not recommended. The skin should be clean, dry and free of detergent residue. Allow all antiseptic preps and skin protectants to dry completely before application of the dressing to prevent skin irritation and to ensure good adhesion.

Application of the 3M Tegaderm CHG I.V. Port Dressing

1. Open package and remove sterile dressing.
2. Peel the liner from the dressing exposing the border adhesive surface.
3. Center the dressing over the accessed noncoring needle. Do not stretch the dressing during application, as mechanical skin trauma may result.
4. Apply dressing to the skin with firm yet gentle pressure, beginning in the dressing region opposite the selected tubing notch.
5. While removing the paper frame, pinch the dressing notch around the tubing and apply firm yet gentle pressure to dressing border to enhance adhesion to skin.
6. The laminated sterile tape strips and documentation label may be used to further secure the I.V. tubing.
7. Document the dressing change information on the label according to the facility protocol.

Site Care

1. The insertion site should be visually observed daily for signs and symptoms of infection or other complications in accordance with facility protocol.
2. The dressing changes should occur at least every 7 days and may be needed more frequently with highly exudative sites or if integrity of the dressing is compromised.
3. Use care and cover the port dressing to protect from water when showering.

Removal

1. Remove tape strips.
2. Always secure I.V. tubing during dressing changes with either tape or gloved finger.
3. In a low and slow manner, gently peel the dressing back upon itself starting at the tubing notch, until the dressing has been removed completely.

Shelf Life and Storage Information

For best results, store in a cool, dry place. For shelf life, refer to the expiration

date on the package. 3M Tegaderm I.V. CHG Advanced Securement Dressing are considered validated to a two (2) year shelf life. Supplied in boxes of individually packaged sterile dressings – 25 per box. Sterility of the dressing is guaranteed unless individual package is damaged or open.

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- This medical device isn't subject to maintenance and repair.

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Transport conditions

Shelf life of the products is guaranteed as labeled, if stored according to the recommended storage conditions. During transportation the products may experience short term temperature variations without affecting product safety and performance.

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For further information outside of the United States, contact your local 3M representative or contact us at 3M.com and select your country.

Explanation of Symbols

Not Made With Natural Rubber Latex

Caution, see instructions for use

Do not use if package is damaged or open

Do not reuse

Use by date

Batch code

Manufacturer

Date of manufacture

Sterilized using ethylene oxide

Do not resterilize

CE 0086

Made in U.S.A. with
globally sourced materials by
3M Health Care
2510 Conway Ave.
St. Paul, MN 55144

Authorised representative of the
manufacturer in Russia:
ZAO "3M Russia", Russia 121614,
Moscow, Krylatskaya str. bldg.17/3
T.+74957847474

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3M and Tegaderm are trademarks of 3M.
© 2016, 3M. All rights reserved.