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This Apostille is not valid for use anywhere within the United States of America, its territories or possessions.

This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears.

This Apostille does not certify the content of the document for which it was issued.

Apostille

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1. District of Columbia, United States of America
AHMAD MISKI
2. This public document has been signed by _____
NOTARY PUBLIC IN AND FOR THE DISTRICT OF COLUMBIA
3. acting in the capacity of _____
AHMAD MISKI, NOTARY PUBLIC IN AND FOR THE

DISTRICT OF COLUMBIA
4. bears the seal/stamp of _____
CERTIFIED
5. at Washington, D.C.
6. the 22 day of JANUARY 2018
7. by Secretary of the District of Columbia
8. No. 466903
9. Seal/Stamp
10. Signature:

LAUREN C. VAUGHAN

Lauren C. Vaughan
Secretary of the District of Columbia

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January 17, 2018


To Whom it May Concern:

This is to certify that the attached document is the Instruction for use for:

Product: Xprecia Systems PT/INR Test Strips

For and on behalf of

Siemens Healthcare Diagnostics Inc.

Bobby Zinek

Regulatory Specialist 3

Siemens Healthcare Diagnostics Inc.


511 Benedict Avenue
Tarrytown, NY, 10591, USA

914 631-8000
www.siemens.com/diagnostics


Based upon available information

US Legalization

I, Yanadi Ahami, Vice President, hereby
certify that is an original document.

Authorized Signature

District of Columbia: SS
Subscribed and sworn to
before me,
this 19 day of Jan, 2018


Ahmad Miski, Notary Public,
District of Columbia

His Commission expires Jan 31st, 2019

Commonwealth of Massachusetts
County of Norfolk

On this **18th** day of **January**, 2018, before me, the undersigned notary public,
personally appeared **Bobby Zinke**, proved to me through satisfactory evidence of
identification, which was **Company Identification Card** to be the person whose name
is signed on the attached document in my presence.

Jody A. Soucie, Notary Public

Expiration Date: January 20, 2023

SIEMENS

Xprecia™ System

PT/INR TEST STRIPS

Xprecia™ Systems PT/INR Strips are designed for use only with the Siemens Xprecia Stride™ Coagulation Analyzer. Other test strips won't work with Stride.

INTENDED USE

The Xprecia Stride Coagulation System is intended for multiple patient use by professional healthcare providers for quantitative prothrombin time (PT) testing for the monitoring of oral anticoagulation therapy with warfarin, a vitamin K antagonist. The system uses fresh capillary (finger-stick) whole blood. It's intended for professional *in vitro* diagnostic use at the point-of-care.

CLEANING AND DISINFECTING STRIDE

You should clean and disinfect the exterior of Stride, and the test strip port protective cap, after patient and quality control tests using a Siemens recommended germicidal wipe. For detailed cleaning instructions, see the User Guide.

PRINCIPLES OF THE PROCEDURE

A PT test is a blood test that measures the time it takes for blood to clot, reporting results using either the International Normalized Ratio (INR) or seconds. The INR was developed to standardize the PT results.

Stride measures PT values in whole blood. To begin testing, the contact end of the test strip is inserted into the test strip port on the analyzer. Then, a blood sample is applied to the test strip target area. The blood sample is automatically drawn by capillary action into the reaction chamber of the strip where the blood mixes with reagents and activates the coagulation cascade. When Stride senses that the blood has clotted, the testing stops. A PT result is calculated and appears on the analyzer as either INR or seconds.

For complete instructions, see the user guide that came with your Xprecia Stride System. If you have questions, contact your local technical support provider.

REAGENTS

The test strips are packaged in 1 vial with 25 test strips. Each test strip contains the reagent Dade® Innovin®, which is a preparation of purified recombinant human tissue factor combined with synthetic phospholipids, calcium, and stabilizers.

WARNINGS AND CAUTIONS

- Always follow the safety procedures and precautions listed here and throughout the User Guide when using Stride.
- All parts of Stride are potentially infectious and capable of transmitting blood-borne pathogens between patients and healthcare professionals. For more information, refer to "Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Healthcare Settings 2007," found at <http://www.cdc.gov/hicpac/2007IP/2007isolationPrecautions.html>.
- You should disinfect the device after each patient use. You can only use Stride for testing multiple patients when all standard precautions and the recommended cleaning and disinfection procedures are followed.

- Use only whole blood from a finger-stick for testing.
- Only auto-disabling, single use lancing devices may be used with Stride.
- Don't consume test strips or drink control solutions.
- Always store the test strips in the original test strip vial with the cap closed.
- Use the test strip within 10 minutes of when you remove it from the test strip vial.
- Apply the blood sample to the test strip within 120 seconds after the Apply Sample screen displays.
- Always wear protective gloves and follow your facility's policies and procedures when performing tests involving biological samples and control solutions.
- Only apply blood to a test strip after it's inserted in, and blood is requested by, Stride.
- Within 15 seconds of sticking the fingertip, apply the drop of blood to the test strip target area.
- Never use bent, scratched, or damaged test strips.
- Don't scan a test strip barcode and then use a different test strip from another vial.
- Each test strip is single-use only; never perform a second test using the same test strip.
- Don't touch or move the test strip after you apply the drop of blood. Don't move the test strip during the test.
- Never add more blood to the test strip after the test has begun.
- A best practice is to make sure there is adequate lighting to conduct a test.
- Refer to the following general safety reference materials:
 - "Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline-Third Edition" Clinical and Laboratory Standards Institute (CLSI) M29-A3.
 - "FDA Public Health Notification: Use of Fingertick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens: Initial Communication" (2010), found at <http://www.fda.gov/medicaldevices/safety/alertsandnotices/ucm224025.html>.
 - "CDC Clinical Reminder: Use of Fingertick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens" (2010), found at <http://www.cdc.gov/injectionsafety/Fingertick-DevicesBGM.html>.

TEST STRIP STORAGE AND STABILITY

Always store the test strips as packaged and use within the expiration date printed on the test strip vial: don't refrigerate them.

- Unopened Storage: 4–30°C (39.2–86°F) @ 0-75% RH
- Opened Storage: 4–30°C (39.2–86°F), use before the 24 month lot expiration date after first opening vial

OPTIONAL LIQUID QUALITY CONTROL

For quality control information, see the instructions for use that came in your control solution package.

ELECTRONIC ONBOARD QUALITY CONTROL

When each test strip is inserted, Stride automatically conducts 2 on-strip quality control checks designed to help ensure test strip integrity. The first control checks the presence of adequate sample reagent on the test strip, and the second control detects test strip degradation due to exposure to environmental conditions.

HOLDING STRIDE

Always handle Stride with care and don't mishandle it. Rough treatment or impact with hard objects, such as dropping, may damage parts and lead to incorrect operational results.

You can use Stride in the following positions:

- Flat on a tabletop or other hard surface (display side up).
- In your hand in a level position.
- In your hand within a 45 degree angle up or down.

- Don't run a patient test or quality control (QC) test when holding Stride at extreme angles.

MATERIALS PROVIDED

- Xprecia Systems PT/INR Strips

MATERIALS REQUIRED BUT NOT PROVIDED

- Xprecia analyzer
- Alcohol wipe, cotton ball or tissue
- Single-use lancing device to collect capillary blood samples obtained from a finger-stick.
Prepare the lancing device according to manufacturer instructions.

PREPARING THE PATIENT SAMPLES

Prepare the patient to collect a blood sample by following your facility's standard procedure. For example, clean the finger with an alcohol wipe, or have the patient wash his or her hands in warm, soapy water. Dry the fingertip before taking a blood sample. Use only whole blood from a finger-stick for testing.

IMPORTANT

- Within 15 seconds of sticking the fingertip, apply the drop of blood to the test strip target area.
- Don't add more blood to the test strip once the test has begun.
- Don't touch the test strip while the test is in progress.
- Discard expired test strip vials.

 **WARNING** Clinical studies have shown that Xprecia Stride PT/INR results from a venous sample were falsely low when compared to a capillary sample from the same patient. Testing of venous samples should not be performed on Xprecia Stride.

PERFORMING A PATIENT TEST

INSTRUCTIONS

1. To turn Stride on, press the On/Off button.
2. If you are prompted, enter an Operator ID (OID).

3. On the Home screen, tap Patient Test.
4. If you are prompted, enter a Patient ID (PID).
5. Open the test strip vial and remove only 1 test strip.
6. Immediately close the vial. Make sure the vial cap seals tightly.
7. To scan the test strip barcode, aim the barcode reader at the barcode on the test strip, then tap the Scan button.
Tip An audible tone sounds, a check mark displays, and the screen changes when the barcode is accepted.
8. Gently, but firmly, insert the test strip with the printed side up, into the test strip port until it stops.
9. (Optional) If Stride requests that you scan the test strip vial, aim the barcode reader at the barcode on the test strip vial, then tap the Scan button.
10. The Wait screen displays as Stride warms the test strip until it reaches operating temperature, in approximately 20 seconds.
11. After the Apply Sample screen displays, prepare and apply the patient blood sample. See the following:
 - Finger-stick sample collection method below.
 The test will begin when the sample is drawn into the test strip by capillary action.
Requirement Be sure to fill the entire target area with blood. Don't overfill the target area.
12. After the Test in Progress screen displays, don't touch the test strip or add more blood. The INR value, time, and date display when the test is complete.
Tip You can set Stride down to tend to the patient while the test is in progress.
13. Press the Test Strip Eject button to discard the test strip according to your facility's biohazard control policies.
Requirement When ejecting a used test strip, always point Stride down facing your biohazard container before you press the Test Strip Eject button.
14. Discard the lancet according to your facility's biohazard control policies.
15. Clean and disinfect the entire exterior surface of Stride, and the test strip port protective cap if necessary, with a Siemens recommended germicidal wipe.
Requirement Siemens recommends you remove, clean, disinfect, and dry the cap when you observe blood or other contaminants inside the opening in the cap.
16. Remove your gloves, thoroughly wash and dry your hands, and put on a new pair of gloves before performing another patient test.
 For complete instructions, see the user guide that came with your Xprecia System. If you have questions, contact your local technical support provider.

COLLECTING A FINGER-STICK BLOOD SAMPLE INSTRUCTIONS

1. To stick the finger, firmly place the lancet against the finger and press the lancet trigger.
2. Gently squeeze from the base of the finger to form a round drop of blood. If the blood smears or runs, wipe it off with a tissue and gently squeeze another round drop of blood.
Requirement The drop of blood should be about the same size as the test strip target area (a minimum of 6 μ L in volume). Low sample volume will cause an error message.
3. Return to *Performing a Patient Test* to follow steps 12–16.

IMPORTANT

- Apply the drop of blood to the test strip target area within 60 seconds of lancing the fingertip.

RESULTS

Results display either as the International Normalized Ratio (INR) or as seconds. Desired INR values may vary depending upon the clinical practice and test methodologies, and the optimum therapeutic range for this method should be established by each user.¹ Normal values can vary and each site should establish its own reference interval. Based on test samples from 120 non anticoagulated subjects, the normal INR range for the Xprecia Stride PT/INR system was 0.9 – 1.1.

Each lot of Xprecia Systems PT/INR test strips is calibrated to a reference lot of human recombinant thromboplastin traceable to the World Health Organization International Reference Preparation.²

Stride also displays PT results in seconds. The reported time is derived from the INR result and the equation below. The calculation is performed with an International Sensitivity Index (ISI) of 1.0 and a typical Mean Normal Plasma Prothrombin Time of 12.0 seconds.

$$\text{INR} = \left(\frac{[\text{Patient Prothrombin Time (sec)}]}{[\text{Mean Normal Prothrombin Time (sec)}]} \right)^{\text{ISI}}$$

LIMITATIONS

The PT/INR assay measures results in patients on warfarin therapy. Don't use this assay to analyze patients on heparin therapy. *In vitro* studies showed Stride tests results are not affected by heparin levels up to 3 U/mL, low molecular weight heparin up to 2 U anti factor Xa activity/mL, ascorbic acid levels up to 1.5 mg/dL, triglycerides up to 3270 mg/dL, hemoglobin up to 200 mg/dL, acetaminophen up to 20 mg/dL, uric acid up to 24 mg/dL, unconjugated bilirubin up to 20 mg/dL, Oritavancin up to 5 mg/L, Fondaparinux up to 2.5 mg/L, Clopidogrel concentration up to 40 mg/L, Daptomycin concentration up to 300 mg/L, and conjugated bilirubin up to 29 mg/dL.

Hematocrit values between 25–55% don't significantly affect test results.

Note Patients with long clotting times (>8 INR, >96 seconds) are outside the reporting range. Confirm results by using an alternative test method. Contact the patient's physician.

EXPECTED VALUES

During the initiation phase of oral anticoagulant therapy, it's prudent to monitor patient status daily or at least 4 to 5 times per week, until some degree of consistency of the INR response to a stable dose is noted.³

The frequency of testing in stabilized anticoagulated patients should be determined on an individual patient basis. In general, don't exceed intervals of 4 weeks when determining INR.

UNUSUAL RESULTS

If an unexpected result is reported or if you are concerned that the result doesn't match the clinical symptoms or patient history, repeat the patient test and then perform a QC test using the Xprecia Systems PT Controls (as described in the Xprecia Systems PT Controls Instructions for


use) to confirm the device is working correctly. If the results are confirmed, more in-depth testing may be necessary. INR results inconsistent with the patient's clinical presentation may be an indication of improper test strip storage, interference with certain drugs, or changes in the patient's diet. Follow up any unexpected results with additional coagulation studies and testing. If the patient's PT value seems unusually low or high, run a QC test using the Xprecia Systems PT Controls as described in the instructions for use that came in your control solution package.

PERFORMANCE CHARACTERISTICS


Reporting Range: 0.8–8.0 INR. Results outside this range will not be displayed and a Results out of Range error message appears.

Sensitivity: The Factor Sensitivity was determined using nine (9) Xprecia Stride™ Coagulation Analyzers, three (3) lots of Xprecia™ System PT/INR test strips, and normal human red blood cells mixed with various concentrations of normal human plasma. Each level of Factor Sensitivity was carried out twice on each analyzer. The study verifies Factor sensitivity at the following levels (% of normal factor level; in vitro testing): Factor II < 36%; Factor V < 58%; Factor VII < 52%; and Factor X < 68%.

Xprecia Stride Factor Sensitivity

Accuracy: Accuracy was evaluated using 412 subjects across 4 sites using 3 lots of test strips per site. The INR of capillary samples measured on the analyzer were compared against the INR of venous plasma samples measured on Siemens BCS XP laboratory analyzer using Dade® Innovin® reagent. Results are shown below:

Repeatability: Whole blood precision was determined for finger-stick blood from sample duplicates over 8 weeks at multiple point-of-care sites using 3 lots of test strips run by typical point-of-care operators. Typical CVs are shown below.

INR Range	n	SD	CV%
0.8–1.9	20	0.06	5.9
2.0–3.0	53	0.10	4.1
3.1–4.5	18	0.14	4.2
4.6–8.0	9	0.20	3.6

Precision: Two levels of Liquid Quality Control were analyzed in duplicate twice a day for 20 days at 3 point-of-care sites run by typical point-of-care operators. Typical CVs were less than 7%.

TECHNICAL ASSISTANCE

For technical assistance, contact your local technical support provider.

SIEMENS

Xprecia™ System

PT/INR STRIPS

English

SIEMENS

REFERENCES

1. Hirsh J, Dalen J E, Anderson DR, Poller L, Bussey H J, Ansell J, Deykin D: Oral anticoagulants: mechanism of action, clinical effectiveness, and optimal therapeutic range. Chest 119: 85-215, 2001
2. WHO Expert Committee on Biological Standardization: Sixty-second Report. Geneva, World Health Organization, 2011 (WHO Technical Report Series, No. 979)
3. Fairweather R B, Ansell J, van den Besselaar A, Brandt J T, Bussey H J, Poller L, Triplett D A, White R H: College of American pathologists conference XXXI on laboratory monitoring of anticoagulant therapy: Laboratory monitoring of oral anticoagulant therapy 122: 768-781, 1998.

MADE IN AU

SIEMENS HEALTHCARE DIAGNOSTICS INC.
511 BENEDICT AVENUE
TARRYTOWN, NY 10591-5097 USA

SIEMENS HEALTHCARE DIAGNOSTICS LTD.
SIR WILLIAM SIEMENS SQ.
FRIMLEY, CAMBERLEY, UK GU16 8QD

10714628 REV. D

DATE OF ISSUANCE NOVEMBER 2017

10714623

WWW.SIEMENS.COM/POC

© 2014–2017 SIEMENS HEALTHCARE DIAGNOSTICS | XPRECIA AND XPRECIA STRIDE ARE TRADEMARKS OF SIEMENS HEALTHCARE DIAGNOSTICS
ALL RIGHTS RESERVED | ALL OTHER TRADEMARKS ARE THE PROPERTY OF THEIR RESPECTIVE OWNERS

Addendum to Instruction for use
Xprecia Systems PT/INR Test Strips

1. MD Name

Xprecia Systems PT/INR Test Strips:

- Xprecia Systems Test Strips (25 test strips in a plastic vial) - 4 vials
- Package insert

2. MD intended Use

It's intended for professional in vitro diagnostic use at the point-of-care. Test-strips are intended for use with the Xprecia Stride coagulation analyzers for quantitative prothrombin time / International Normalized Ratio (PT/INR) testing for the monitoring of oral anticoagulation therapy with warfarin, a vitamin K antagonist.

3. Field of Application

Clinical laboratory diagnostics

4. Contraindications

This test is contraindicated for patients receiving heparin.

5. Possible Adverse Effects of Medical Device Use

Not applicable to this MD.

6 Information on Potential MD Consumers

The medical device is held for use in treatment and prevention facilities and the research centers for carrying out researches only by laboratory and medical personnel.

7. Indications

PT/INR test strips for Xprecia systems are designed for use with Xprecia Stride coagulation analyzers (hereinafter, "analyzer" and "Xprecia analyzer" terms may occur). The analyzer is intended for quantification of prothrombin time (PT) during the monitoring of oral anticoagulation therapy with warfarin, a vitamin K antagonist.

Warfarin is prescribed to prevent blood clots from forming or growing larger in blood or blood vessels. When using anticoagulation medication, its blood level have to stay within a specific therapeutic range, as determined by a doctor. The doctor needs to monitor the warfarin activity

to ensure that the medication dosage is correct. To monitor the activity, the doctor orders a PT test.

8. Operating Conditions

Temperature range	15-35 °C
Ambient relative humidity	20-80% non-condensing.
Pressure	atmospheric

9. Analyzing Sample Type

Whole blood sample.

10. Target Analyte Description

Analyte: protrombin;

Assay type: qualitative

11. Manufacturer Information

MD Developer

Name: Siemens Healthcare Diagnostics Inc.

Address: 511 Benedict Avenue Tarrytown, NY 10591-5097, USA

Phone: +1 302-631-8000

Website: www.siemens.com/diagnostics

MD Manufacturer

Name: Siemens Healthcare Diagnostics Inc.

Address: 511 Benedict Avenue Tarrytown, NY 10591-5097, USA

Phone: +1 302-631-8000

Website: www.siemens.com/diagnostics

MD Manufacturing Site Address

Universal Biosensor Pty Ltd.

1 Corporate Avenue, Rowville, Victoria, 3178, Australia

12. Shipping and Storage condition

Store in accordance with local regulations. Store in original container protected from direct sunlight in a dry, cool and well-ventilated place, away from incompatible materials, foods and

drinks. The container should be tightly sealed before use. Opened containers should be resealed and kept upright to avoid possibility of leakages. Do not store in containers with no labeling. Use appropriate measures to avoid environmental contamination.

Shipping Conditions


Temperature: -20 to +40°C

Humidity: 0% to 85% RH non-condensing @40°C

Atmospheric pressure

13. Accuracy.

Accuracy was evaluated using 412 subjects across 4 sites using 3 lots of test strips per site. The INR of capillary samples measured on the analyzer were compared against the INR of venous plasma samples measured on Siemens BCS XP laboratory analyzer using Dade Innovin reagent. Results are shown below:

The following statistics will be derived to compare the method vs. **laboratory measurement procedure** using whole capillary blood (and whole venous blood target value only).

Over measurement range

- Slope: confidence interval 95% within 0.80–1.20
- Intercept: from +0.3 to -0.3
- Determination coefficient (r^2) ≥ 0.82 ***; target value only: $r^2 \geq 0.90$

For range over 6.0 INR *

INR range	Values			
	Median Bias		Permissible difference (90% bias)	
	Target value	Mandatory level	Target value	Mandatory level
0.8–1.9	± 0.2	± 0.3	± 0.5	± 0.5
2.0–4.5	± 0.2	± 0.3	± 0.5	± 30%
4.6–6.0	± 0.4	± 0.7	N/A	± 30%

* If measuring range reaches 8.0, then for range 2.0–8.0: 97% bias values should be within ± 40%; applicable for combined data sets from multiple clinical sites only.

** Applicable for combined data sets from multiple clinical sites or for data sets with $n \geq 10$ within this INR analyzed range only.

*** For combined data sets from multiple clinical sites: $r^2 \geq 0.85$

For submitting in US regulation authority only, the following statistics will be derived to compare the method vs. **reference device** (capillary blood samples).

Over measurement range

- Slope: 0.80–1.20
- Intercept: from +0.5 to –0.5
- Detemonation coefficient (r^2) > 0.80

Target values only

The following statistics will be derived to compare the method during measurement of capillary (reference) and venous blood with Xprecia Stride.

Over measurement range

- Slope: within 0.90–1.10
- Intercept: from +0.3 to –0.3
- Detemonation coefficient (r^2) > 0.90, target value: 0.95

14. Precision

Two levels of liquid quality control were tested in duplicate bid for 20 days in 3 therapeutic sites; the testing was performed by ordinary medical stuff. Exemplary CVs were less than 7%.

Precision — Repeatability	The value is less than or equal to 10% CV for samples falling into the following INR ranges: • INR less than 2.0 • INR 2.0–3.0 • INR 3.1–4.5 • INR 4.6–8.0
Precision — Intermediate	The value is less than or equal to 10% CV for samples falling into the following ranges. • PT LQC 1 • PT LQC 2

15. Information on MD Disinfection, Cleaning and Sterilization

Not applicable to this MD

16. Environmental Protection

16.1 Environmental Requirements when Using the MD

This device does not harm the environment if properly operated, shipped and stored.

16.2 Disposal and Destruction Procedure

Follow the next instructions for the disposal of hazardous waste:

- Dispose the used test strips in an approved hazardous waste container.
- Follow an institutional policy for the disposal of hazardous waste.
- When ejecting a used test strip, always point the analyzer down facing your biohazard container before you press the Test Strip Eject button.

Waste formation should be avoided or minimized if possible. This product, solutions and product derivatives should be disposed in accordance with the requirements of environmental and waste disposal legislation as well as any requirement of local authorities. The excessive amounts of product and non-recyclable products should be disposed by contractor having license for waste disposal. Do not stream unprocessed wastes in sewers, excluding the cases of full conformance with all authorized authorities requirements. The package is subject to recycling. Incineration or landfill should only be performed when recycling is not feasible.

Utilize all materials by a safe and acceptable method in according with the operating normative regulations. Used medical products have to be utilized in accordance with the local legislation as epidemiologically dangerous wastes.

17. Warranty Obligations

This device has been developed, manufactured, tested, and packaged in compliance with relevant safety procedures. Siemens Healthcare Diagnostics Inc provides warranty for the product for the period of shelf life on condition the package is not damaged.

The manufacturer guarantees the safety of medical device, the absence of unacceptable risks of hazards to human life, health and the environment when using test strips as intended in the conditions provided for by the manufacturer.

In the case of issues regarding Manufacturer's warranty liabilities please contact a legal representative of Manufacturer in Russia - Siemens Healthcare LLC:

18. Complaints

On the product quality issues, please contact:

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
1. Hirsh J, Dalen J E, Anderson DR, Poller L, Bussey H I, Ansell J, Deykin D: Oral anticoagulants: mechanism of action, clinical effectiveness, and optimal therapeutic range. Chest 119: 85-215, 2001
2. WHO Expert Committee on Biological Standardization. Sixty-second Report Geneva, World Health Organization, 2011 (WHO Technical Report Series, No. 979).
3. Fairweather R B, Ansell J, van den Besselaar A, Brandt J T, Bussey H I, Poller L, Triplett D A, White R H: College of American pathologists conference XXXI on laboratory monitoring of anticoagulant therapy. Laboratory monitoring of oral anticoagulant therapy 122: 768-781, 199

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
Apostille

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1. District of Columbia, United States of America
2. This public document has been signed by AHMAD MISKI
3. acting in the capacity of NOTARY PUBLIC IN AND FOR THE DISTRICT OF COLUMBIA
AHMAD MISKI, NOTARY PUBLIC IN AND FOR THE
DISTRICT OF COLUMBIA
4. bears the seal/stamp of _____
5. at Washington, D.C.
6. the 16, day of AUGUST 2017
7. by Secretary of the District of Columbia
8. No. 452625
9. Seal/Stamp
10. Signature:

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August 15, 2017

To Whom it May Concern:

This is to certify that the attached document is the Instruction for Use for:

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Jody A. Soucie
Regulatory Technical Specialist

Commonwealth of Massachusetts
On this 15th day of August, 20 17
I certify that the IFU
document is a true, exact, complete and unaltered copy
of the original.
 Sandra L. Giampa
SANDRA L. GIAMPA, Notary Public
My Commission Expires November 5, 2021

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SIEMENS

REFERENCES

1. Quick A., Stanly-Brown M. and Bancroft FW: A study of the coagulation defect in hemophilia and in jaundice. Am. J. Med. Sci 190: 501-511, 1935.
2. Hirsh J, Daen J.E, Anderson DR, Peller L, Bussey - I, Ansel J, Daykin C: Oral anticoagulants: median sir of action, clinical effectiveness, and optimal therapeutic range. Chest 115: 85-2'S 2001.
3. Hirsh J, Daen J.E, Anderson DR, Peller L, Bussey - I, Ansel J, Daykin C: Oral anticoagulants: median sir of action, clinical effectiveness, and optimal therapeutic range. Chest 115: 85-2'S 2001.
4. WHO Expert Committee on Biological Standardization. Thirty-third Report. Geneva, World Health Organization, 1983 (WHO Technical Report Series, No. 687).
5. Fairweather R.B, Ansel J, Baselaer A, Brandt J.T, Bussey - I, Peller L, Tiplatt D.A, White R.H: College of American pathologists conference XXXI on laboratory monitoring of anticoagulant therapy. Laboratory monitoring of oral anticoagulant therapy 22: 76B-781, 1999.

CE MADE IN AU

SIEMENS HEALTHCARE DIAGNOSTICS INC.
511 BENEDICT AVENUE
TARRYTOWN, NY 10591-5097 USA

SIEMENS HEALTHCARE DIAGNOSTICS LTD.
SIR WILLIAM SIEMENS SQ.
FRIMLEY, CAMBERLEY, UK GU16 8QD

11065099 REV. A



DATE OF ISSUANCE JULY 2015

11064957

WWW.SIEMENS.COM/POC

© 2015 SIEMENS HEALTHCARE DIAGNOSTICS | X-PEC A AND X-PEC A STRIDE ARE TRADEMARKS OF SIEMENS HEALTHCARE DIAGNOSTICS
ALL RIGHTS RESERVED. ALL OTHER TRADEMARKS ARE THE PROPERTY OF THEIR RESPECTIVE OWNERS

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