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
Kit for endovascular diagnostic and therapeutic pro-
cedures, including angiography, angioplasty, consist-
ing of:

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Michael Keaveny,
Notary Public,
Headsfort Place, Kells,
County Meath, Ireland.

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"АркРойал"(ArcRoyal), Ирландия
Virginia Road, Kells, Co. Meath, Ireland

Manufactured by

ArcRoyal, Ireland.
Virginia Road, Kells, Co. Meath, Ireland

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Kit for endovascular diagnostic and therapeutic procedures (PTCA Accessory Kit) consists of Y-connector, Torque device and guidewire insertion tool

Y- Connectors are intended for application during endovascular procedures. Y-connector could be used with diagnostic catheter, guiding catheter or sheath introducer. Y-connector is manufactured of transparent plastic to control formation of air bubbles, provided with one or two adjustable haemostatic valves allowing to introduce instruments with diameter from 3F to 9F, preventing blood reflux and air aspiration. Y-connector is provided with a sideport, which allows washing out the instrument which is placed the lumen of catheter/introducer, and could be used as additional infusion line for injection of contrast or other medicinal solutions.

The Torque Device is intended for *short-term non-invasive use* (not intended in whole or in part, to penetrate inside the body, either through a body orifice or through the surface of the body), used to control a guidewire during cardiovascular procedures to introduce guidewires into the vascular system, allowing them to be positioned and locked in place.

Guidewire insertion tool is used to facilitate the introduction of a guide wire by providing sup-

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Common side effects, known at endovascular angiographic and another intravascular treatment procedures with application of guiding catheter procedures, include but are not limited to:

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- Vasospasm
- Thrombus
- Arterial dissection
- Vessel perforation
- Death

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- Keep the Y-connector filled with an appropriate flushing solution when it is in use
- When withdrawing fluid from the sideport, ensure the adjustable Haemostatic valve is/are closed and the Y-connector is monitored to prevent air from entering the Y-connector

Procedure:

- Prior to use, attach an infusion line to the female luer on the Y-connectors sideport
- Using normal saline, flush the Y-connector to remove any air from all channels. To flush the valve section, open the adjustable valve by rotating the thumb wheel counter clockwise and continue to fill the assembly
- Connect the rotating male luer adapter of the Y-connector to the guiding catheter. Ensure the connection is airless.
- Open the adjustable Haemostatic valve and gently insert the Guidewire and/or the catheter. Use the Guidewire insertion tool to facilitate Guidewire entry, removing the tool once entry is accomplished.
- Carefully tighten the adjustable thumb wheel by rotating clockwise around the Guidewire and/or dilation catheter to create a seal
- Advance the Guidewire and/or catheter into the guiding catheter and then the desired vessel
- Perform the dilation by following standard angioplasty procedures
- Following angioplasty, deflate the catheter and loosen the adjustable Haemostatic valve/s of the Y-connector Before retracting the balloon and/or Guidewire
- Disconnect the Y-connector from the guiding catheter and discard

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package.

2. Loosen the cap to the cylinder, by turning several times counterclockwise (about 720 degrees).
3. Introduce the proximal end of the guidewire into hole of the cap; move the torque device along the guidewire for a necessary length (15-45 cm) and fix the device by turning the cap a few times in a clockwise direction.

Guidewire insertion tool

Use the guide wire insertion tool to facilitate the introduction of a guide wire during general intravascular procedures. This is achieved by providing support during insertion and passage through hemostasis valves and protecting the guide wire tips.

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damaged or opened, remove from the operative area to prevent unintended use, and return to the manufacturer.

- This device is for single patient use only. It is not intended for disinfection and subsequent re-use, which may result in device failure or create the risk of contamination.
- This device has not been evaluated for reprocessing or re-sterilization. Reprocessing or re-sterilization may damage the device, rendering it unusable or may lead to device failure, which could result in patient illness, injury or death.
- Do not expose them to organic solvents or ionizing radiation. Do not autoclave.
- Do not use after the expiration date stated on the packaging labels.
- Consider the following to prevent or reduce the risk of Clot formation:
 - Use systemic anticoagulants
 - Finish or rinse the Y-connector with an appropriate flushing solution prior to use to remove all air

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The manufacturer guarantees conformity of medical devices to the requirements, indicated in this instruction for use provided that usage, operation, transportation and storage terms set forth herein are met.

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The medical device is transported using all available modes of transportation in the covered vehicles according to the transportation rules forced on each transportation mode, in cool, dry, dark place, at temperature from 0°C to +40°C and air relative humidity from 20% to 80% (without condensation).

Medical device is stored at a room temperature from +15°C to +25°C, at a place protected from direct sunlight.

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After use all used components and packing materials may be potentially biologically dangerous. Their treatment and disposal are performed according to the accepted medical practice, and also according to local, state and federal laws and provisions.

The products which had a contact with blood and / or body fluids are subject to disposal and / or destruction according to SanPiN 2.1.7.2790-10 (B class).

Unused products (which had no contact with blood and / or body fluids), including those with expired shelf-life, are subject to disposal and / or destruction, according to SanPiN 2.1.7.2790-10 (Class A).

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These medical devices are for single use only and are not subject to maintenance and repair.

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Authorized Representative of the manufacturer
in the Russian Federation:

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Johnson & Johnson LLC

Russian Federation, Krylatskaya street -
17, building 2, Moscow, 121614

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