(EN) Instructions for use

bredent
Weissehorne Strasse 2
49250 Senden - Germany
GmbH & Co. KG
Gerald Micko // 25.10.2018
General Manager

009486EX-20181025

CE
0483

The mass of dental thermoplastic Polyan IC

Before using the product, please read these instructions for use carefully.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Dental thermoplastic mass Polyan IC, options:

1. **Polyan IC transparent** (540PI024 in cartridges of 24 g each - 20 pcs., 540PI030 in cartridges of 30 g-20 pcs., 540PI005 in a bank of 500 g)
2. **Polyan IC pink 1** (540PI124 in cartridges of 24 g each - 20 pcs., 540PI130- in cartridges of 30 g-20 pcs., 540PI105 in a bank of 500 g)
3. **Polyan IC pink 2** (540PI224 in cartridges of 24 g each - 20 pcs., 540PI230 in cartridges of 30 g-20 pcs., 540PI205 in a bank of 500 g)
4. **Polyan IC pink 3** (540PI324 in cartridges of 24 g each - 20 pcs., 540PI330 in cartridges of 30 g-20 pcs., 540PI305 in a bank of 500 g)
5. **Polyan IC pink 4, with streaks** (540PI424 in cartridges of 24 g each - 20 pcs., 540PI430 in cartridges of 30 g-20 pcs., 540PI405 in a bank of 500 g)
6. **Polyan IC pink 5** (540PI524 in cartridges of 24 g each - 20 pcs., 540PI530 in cartridges of 30 g-20 pcs., 540PI505 in a bank of 500 g)

Description of the product

The product Polyan IC is a thermoplastic resins which can be process with the resin injection moulding system „thermopress 400“ of bredent GmbH & Co. KG. The processing parameters to ensure reproducible quality are preset in the thermopress 400 unit.

These so-called Polyan IC is available in different colors. All thermoplastics are biocompatible resins, which can be used to fabricate metal-free, removable or fixed dental restorations depending on the type of resin that is used. For the approved dental indications, please refer to the table listed under item 2. Polyan IC is prefilled in aluminium cartridges as a granular material and is available in different filling quantities or in a bulk quantity of 500 g. No medicine is included into the medical device.

2. Dental-technical indications


Patient target group

Adults with partial or complete edentulous upper and / or lower jaw

Mechanical retentions:

To achieve reliable bonding of denture tooth and denture base, the contact surfaces of all prefabricated teeth need to be sandblasted with 110 µm aluminium oxide. A circumferential groove (1) needs to be milled in the cervical area and mechanical retention (undercut) (2) needs to be included in the basal area (see drawing).

High-quality extensions or repairs are carried out by rebasing in the injection technique (retentive). When using Polyan IC, it is also possible to perform repairs using conventional PMMA-based cold-curing resins, such as uni.lign and top.lign professional for bre.dentan HP.

Product	Removable/Removable with restrictions	Fixed	Additional information
Polyan IC (PMMA) hard-rigid	Denture bases for full and partial dentures and implant-supported restorations, splints, resin saddles, simple orthodontic appliances (such as Bionator)	not recommended	Colors: crystal-clear, pink 1, pink 2, pink 3, pink 4 veined and pink 5. To ensure reliable bonding of denture tooth and Polyan IC base material, the teeth must have mechanical retentions (grooves) and need to be conditioned with our special bonding agent poly.link IC [REF polylnk5] and used in accordance with the instructions for use of this product. After completion of the injection moulding process, allow the flask to cool down to room temperature for at least 30 min. and then devest carefully.

General information on proper preparation of acrylic teeth to be anchored in thermoplastic materials

Measures	** Polyan IC	
	poly.link IC	visio.link
*Chemical bond between denture tooth and base material	✓	-
Increase of microretentive abonding to the base materal	-	✓
mechanical retentions required	✓	✓

* No chemical bond of base material to base material (rebasing)!

** To achieve red-white esthetics with crea.lign, it is sufficient to sandblast with 110 µm and condition with visio.link (visio.lign system) on the Polyan IC base material.

3. Dental–technical contraindications

Polyan IC	Retention elements such as retainers, attachments, telescopes, crowns and bridges
------------------	---

The products must not be used for patients with a suspected allergy until allergy tests have been completed to confirm that the patient is not allergic to the materials contained in the product.

Material samples for allergy sufferers

Material sample	REF	Color	Qty	Type
Polyan IC	99250PIC	crystal-clear	1 piece	PMMA
Polyan IC	99250PI1	pink 1	1 piece	PMMA
Polyan IC	99250PI2	pink 2	1 piece	PMMA
Polyan IC	99250PI3	pink 3	1 piece	PMMA
Polyan IC	99250PI4	pink 4	1 piece	PMMA
Polyan IC	99250PI5	pink 5	1 piece	PMMA

Thermopress materials – potential allergenic substances			
The following information is not binding and based only on professional assessment. It is not based on results of analyses or information provided by our raw material manufacturers or suppliers. The respective material should not be used for patients with a suspected allergy until allergy tests have been completed to confirm that the patient is not allergic to the materials contained in the product.			
Traces of polymer type can be contained:	Polyan IC		
	PMMA (Polymethylmethacrylate)		
	yes	no	n.s.
Methyl methacrylat	x		
Ethylene glycol dimethacrylate (EDGMA)		x	
Triethylene glycol dimethacrylate (TEGDMA)		x	
Urethane diol		x	
BIS-GMA (Bisphenol A-dimethacrylate) 2,2-Bis-(2-hydroxy-ethacryloxypropoxy)-phenyl)propane		x	
Butane diol dimethacrylate (BDDMA)		x	
2-hydroxyethyl metacrylate (HEMA)		x	
Camphor quinone		x	
Benzoyl peroxide (BPO)			x
Dimethyl - p -toluidin		x	
N,N-Bis 2-hydroxyethyl)-p-toluidine		x	
N,N-D4T: N,N-Dimethyl-4 toluidine		x	
Hydroxy methoxy benzophenone		x	
Butyl-hydroxy-toluene (BHT, DBPC)			x
Formaldehyde	x		
Colophony		x	
Cobalt chloride			x
Phthalates			x

Please consult a doctor if you feel any paraesthesia (inflammation/irritation) or if the denture breaks.

4. Hazard and safety information

Do not inhale vapors during the melting phase. Use caution in case of spilled granular material! Risk of accident by slipping! For further safety information about the product, please refer to the safety datasheet.

The lot numbers of all components should be recorded to ensure traceability and to handle complaints.

Information – 500 g bulk quantity:

Please ensure dry storage to prevent the material from absorbing any moisture. Water vapor resulting from residual moisture during the heating phase of the cartridge in the unit may result in bursting of the cartridge!

Prior to the use, please predry the non-mixed granular material, the integration of foreign substances or contaminated particles must be avoided when filling the cartridges. Different thermoplastic granular materials must not be mixed.

5. Storage, transport

Open package must be protected from moisture. The granular material needs to be stored under dry and clean conditions of transport, see samples.

Sterile

Not sterile, not sterilisable

6. Processing / use

Installation
No installation necessary

Single use
Based on the usage of the product Polyan IC is a single use device
Polyan IC is processed by thermopress 400 injection molding unit.

Model fabrication:

When processing Polyan IC, a formaldehyde-free class IV stone is used for the fabrication of plaster models. Exakto-Rock S (REF 5700SE50 ivory or 5700SB50 white) is particularly suitable; alternatively, Thixo Rock can be used.

Mixing ratios:

Utilization should be carried out in according to norms and rules approved by the Russian legastatic.

Use of formaldehyde-free class IV plaster

Material	Model	Exakto-Rock S (or Thixo Rock)	Distilled water	Hardening time
Polyan IC	per	100 g	20 ml	1 h


Please invest the finished model into the flask once the total expansion time for the plaster has elapsed.

Wax model:

The dimensions used for wax models should be very generous. In other word, clasps and other delicate elements should be made thicker. A modelling wax with a low melting point, e. g. Biotec modelling wax (REF 510 0061 O), is ideal for the boiling out phase with boiling water. As a result, any color residues on the injection-moulded thermoplastic part are avoided.

Film sprue:

A film sprue must be attached to the wax model in order to benefit from the outstanding mechanical properties of thermoplastics. Apply this film sprue using a 1.5 mm thick wax plate. Using a 10 mm thick injection moulding channel (REF 430 0172 O), attach the molten thermoplastic directly to the model via the 1.5 mm wax plate.

Processing parameters:

Investment technique	Material	Set temp. in °C	Heating time minutes	Pressure time seconds	Injection Speed	Injection Force	Flask/ring temperature in °C	Program No. in thermo-press 400 version 2.4 / *2.57
Flask	*Polyan IC	265 °C	15 min.	60 sec.	8 M/SEC	145 KN	40 °C	1

Please adjust the following parameters for compressed air operated injection moulding devices:

Material	Parameter – pressure (force) in bar
Polyan IC	9.5

To facilitate opening the halves of the flask and avoid damaging the moulded object, the press-out device and punch (REF 14000904) should be used. Any larger pieces of plaster must be carefully removed from the moulded object using a pneumatic devesting chisel. Fine cross-cut silicone burs (REF S263QG60) have proven to be perfectly suitable for finishing and reducing the moulded object.

The special bur set (REF 33000830) is recommended for easier finishing of the thermoplastics.

Polishing

First of all, smooth the structure down using sandpaper with a 280 grain, before evenly smoothing it down with the handpiece using a 400. Then the structure is prepolished with pumice powder and a goat-hair brush (REF 35000610) using the polishing motor without applying any pressure. Complete high gloss polishing with the handpiece using a goat-hair brush (REF 35000550) and the high-gloss polishing paste Abraso-Star K50 (REF 52000161).

7. Cleaning

When using heated cleaning devices (ultrasonic), make sure not to exceed a temperature of 50 °C for the cleaning liquid. Do not use any cleaning agents that contain alcohol or acids.

8. Symbols

8.1 Symbols according to DIN EN ISO 15223-1

REF	Reference number
LOT	Lot number
Qty	Quantity
	Dry storage
	Attention
	Manufacturer

8.2 Other symbols

	Preheating temperature
	Preheating time
	Injection speed
	Injection force

9. Additional information

Shelf life / storage

If the pellets are properly stored in the original packaging, the shelf life is not limited.

Service life (period of serviceability)

The service life of a prosthesis made of Polyan IC material, manufactured technically correctly, can reach 10 years.

Recycling

On the basis of the necessary technical regulation and after consulting with disposal agents and competent authorities, it can be disposed of with household waste or incinerated with household waste. Small amount can be disposed of together with municipal waste.

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used by dentists, dental technicians and adequately trained dental staff and is intended for the indication described under item 2. The processing parameters described under item 6 must be strictly adhered to in order to avoid material degradation and the possibility of chipping or fracture of the completed restorations. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products. If components of other systems are used, all claims for warranty or compensation are excluded!

Medical devices „Polyan IC“ shall be conveyed using all types of covered vehicles in accordance with the requirements and rules of cargo transportation applicable to the corresponding mode of transport. Type of shipping – by parcel, motor vehicles, planes, containers.

Transportation conditions: ambient temperature; -10°C to 45°C, relative humidity not more than 80%.

Product boxes should not be exposed to excessive shocks and atmospheric precipitation during handling and transportation, in accordance with handling marks bearing the respective signs „Fragile, handle with care“, „Top side“, „Protect from moisture“.

The manufacturer guarantees the conformity Polyan IC to all regulatory requirements subject to the conditions of operation, transportation and storage.

The Manufacturer's warranty relates to the unopened, undamaged and correctly stored product.

Name of manufacturer:

bredent GmbH & Co. KG, Weißenhorner Straße 2, 89250 Senden, Germany,
+49 7309 872-443, e-mail: info@bredent.com

Name of product representative company in Russia:

Limited Liability Company "ALADENTT" (LLC <<ALADENTT>>), 123103, Moscow, Marshal Zhukov Avenue, 74, building 1, office 49, +7 (495) 585-51-10, e-mail: office@aladent.ru

bredent
Weissenhorner Strasse 2
89250 Senden Germany
bredent GmbH & Co. KG
Gerald Micko // 25.10.2018
General Manager

File-No. 3568/2018

I hereby certify that the document affixed and sealed hereto was signed on Page 1 and 6 in my presence by

Mr. Gerald Micko, born on 04. May 1966,
88471 Laupheim, Zwischen den Wegen 4,
personally known to me

acting on behalf of

bredent Beteiligungs- und Verwaltungs-GmbH in Senden (89250 Senden, Weißenhorner Str. 2)
as this company's director,
the latter acting as general partner of

bredent GmbH & Co KG in Senden (89250 Senden, Weißenhorner Str. 2)

Upon my inspection of the commercial register of 24th October 2018 I further hereby certify that bredent Beteiligungs- und Verwaltungs-GmbH in Senden and bredent GmbH & Co. KG in Senden are registered with the Municipal Court of Memmingen - Commercial Register - under nos. HRB 12912 and HRA 11358, respectively, and that Mr. Peter Brehm is entitled to act individually as legal representative of bredent Beteiligungs- und Verwaltungs-GmbH in Senden as well as the latter is entitled to act individually as legal representative of bredent GmbH & Co. KG in Senden.

And I do hereby further certify, that

- the aforesaid Corporation is duly incorporated
- it is in good standing under the laws of the Federal Republic of Germany
- it has a legal corporate existence;
- its registered corporate domicile is Senden

so far as the commercial register shows.

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