
Instruments to the TCM 3000 BL – Morcellator for Surgical Procedures

Operating Manual

CE 1275

Nouvag AG – St.Gallerstrasse, 23-25, CH 9403 Goldach, Switzerland

Tel. +41 (0)71 846 66 00, Fax +41 (0)71 846 66 70

info@nouvag.com www.nouvag.com

Date: 20.09.2017

Leszek Kazana, Area Sales Manager


Table of contents

	page
INTRODUCTION	4
1. PURPOSE	4
2. INDICATIONS	4
3. CONTRAINDICATIONS	4
4. WARNINGS AND PRECAUTIONS	5
5. SIDE EFFECTS	5
6. OPERATION CONDITIONS	5
7. General information about composition of the TCM 3000 BL Morcellator for Surgical Procedures with Appliances	8
8. INFORMATION ABOUT MP	9
9. CONSTRUCTION AND OPERATION	36
1. Principles of Operation	36
2 Application	36
3. Method of application	36
4. MP overview	37
5. Operation Plan for Instruments for the TCM 300 BL – Morcellator for Surgical Procedures	39
6 Before use	41
6.1 Preparation of MP	41
6.1.1 Assembly of the morcellator set with protective tube (option 1)	41
6.1.2 Assembly of the Morcellator set with trocar cannula (option 2)	45
6.1.3 Assembly of the trocar cannula with myoma drill and mounting gripper (option 3)	49
7 Use	49
7.1 Using the control panel	50
7.2 Operation	51
7.2.1 Option 1: morcellation with protective tube	53
7.2.2 Option 2: morcellation with trocar cannula	55
7.2.3 Option 3: myomectomy	56
8. After use	57
10. SCOPE OF DELIVERY	60
11. CLEANING OF MP COMPONENTS	63
General information about the disinfection, presterilization cleaning and Sterilization of MP	64
General recommendations on preparation of MP	65

	Instruments preparation	67
11.1	Preparation of the trocar cannula	68
11.2	Preparation of cutting tubes, protective tubes, obturators, grippers and myoma drill	69
12.	MAINTENANCE AND REPAIR	71
12.1	General maintenance	71
13.	SAFETY INFORMATION	71
13.1	Modification and misuse	72
13.2	Use precautions	72
13.3	List of operational hazards, risk analysis	72
13.4	Safety requirements, operation conditions and rules	72
13.5	Safety check	73
14.	MEDICAL PRODUCT CLASSIFICATION	73
15.	COMPLIANCE WITH INTERNATIONAL STANDARDS	73
16.	PACKING AND MARKING	74
17.	TRANSPORTATION AND STORAGE CONDITIONS	84
18.	SERVICE LIFE	84
19.	INFORMATION ON DISPOSAL	84
20.	WARRANTY COVERAGE	85
21.	CLAIMS	85
	Appendix 1. Parts of the MP	86
	Appendix 2. Table of components of the MP	87
	Appendix 3. Certificate of Analysis	94

General

INTRODUCTION

Nowadays newly developed medical products (MP) make hysterectomy possible by a laparoscopic intervention. The advantage for the patient of the laparoscopic method is less pain and a much faster recovery. Patients appreciate much the minimal invasive character of this operation because it leaves almost no scars on the patient's skin, reduces the healing of the wounds and saves money and time for all involved parties. A morcellator being a special device for morcellation and further extraction of mass organs and tissues directly in the body cavity without additional incision of tissues to remove the organ is used for such procedures. In order to satisfy with requirements for each particular procedure Nouvag AG offers instruments intended for use with TCM 3000 BL Morcellator for Surgical Procedures.

1. PURPOSE

Instruments to the TCM 3000 BL - Morcellator for Surgical Procedures are used for morcellation and extraction of tissue and for removing myomas of the uterus during laparoscopic procedures in gynecology..

2. INDICATIONS.

Instruments to the TCM 3000 BL - Morcellator for Surgical Procedures are used used for surgical laparoscopic intervention during general surgical and gynecological procedures for morcellation and extraction of tissues in the following cases:

- myoma of the uterus;
- fibroma of the uterus;
- uterus neoplasms requiring partial hysterectomy;
- benign ovarian neoplasms.

3. CONTRAINDICATIONS

Relative or absolute contraindications may result from the general medical findings or in special cases in which the patient's risk for motor-driven tools is significantly increased. Ovaries, fallopian tubes, myomata and other structures must be devascularized and dissected before morcellation.

The laparoscopic deployment of morcellators is contraindicated at:

- Treatment of malignant tumors
- Treatment of vascularized tissue

- Preparation of tissue
- In gynecologic surgery in which tissue to be morcellated is known or suspected to contain malignancy
- For removal of uterine tissue containing fibroids in patients who are...
 - Pre or postmenopausal, or
 - Candidates for the bloc tissue removal, for example through the vagina or via a mini-laparotomy incision

4. WARNINGS AND PRECAUTIONS

- Inadvertent morcellation in the peritoneal cavity can result in severe damage of organs and/or abdominal wall
- The sharp edges of the MP bear a risk of tissue damage.
- During the procedure, tissue can spread in the peritoneal cavity.
- The morcellator produces heat when continuously used for a longer period.
- It is not allowed to perform a morcellation without protective tube or trocar cannula.

5. SIDE EFFECTS

Tissue morcellation with electric-powered surgical cutting MP is a well-established procedure. There have never been reported any unacceptable side effects provided the MP is properly used as intended. Acceptable side effects include sensation of pain and anesthesia reactions.

6. OPERATION CONDITIONS

To ensure proper operation of the Instruments to the TCM 3000 BL - Morcellator for Surgical Procedures one must meet the following conditions:

Ambient Conditions:

	Operation:
Relative humidity:	30% - 80%.
Temperature:	10°C - 30°C
Atmospheric pressure:	800 hPa - 1060 hPa

	Before use read this manual and the manuals for all other equipment used in the procedure thoroughly.		The products described in this manual must be used by trained and qualified medical personnel only in accordance with the requirements and with reference to the manual.
	For protection of the patient, the user and third parties the products used must be prepared as specified in the manual.		Improper use or repair of the MP, or failure to observe these instructions, will relieve us from any obligation arising from warranty provisions or other claims.
	The use of third-party products is the responsibility of the operator. Functionality and patient safety cannot be guaranteed with third-party accessories.		Prior to using the MP, before startup, and before operation, the user must always ensure that the Instruments are in good working order and are clean, sterile and operational.
	The MP must be protected against direct sunrays, heat supplied by centrally batteries, excessive dust, vibration and mechanical shocks.		The MP is re-usable. This MP is not repairable.
	Uterine tissue may contain unsuspected cancer. The use of laparoscopic power morcellators during fibroid surgery may spread cancer, and decrease the long-term survival of patients. This information should be shared with patients when considering surgery with the use of these MPs.		Instruments to the TCM 3000 BL Morcellator for Surgical Procedures should be operated under constant surveillance of qualified medical professionals.
	The MP is not sterile on delivery. All sterilizable parts must be sterilized before use (see Section Instrument Preparation).		Never operate the clamping mechanisms of the cutting tubes while the system is running. This could result in instrument damage.
	Operate the MP only outside the danger zone of explosive, flammable mixtures or gases.		Keep hands clear of the cutting tubes when operating the MP. Danger of serious injury.
	Disregarding the guidelines on intermittent operation (1 min on, 1 min off for 8 cycles) may cause burns at contact with the gear unit.		Before starting use and operation of the MP obligatory read the Instructions for Use for the medical device TCM 3000 BL Morcellator for Surgical Procedures with Accessories.

Modification and misuse

- Replacement of MP parts with external ones or modification of Instruments for to the TCM 3000 BL - Morcellator for Surgical Procedures is prohibited. The manufacturer is not liable for any damages resulting from unauthorized modifications or manipulations. The warranty will be cancelled.
- Use of the Instruments for to the TCM 3000 BL - Morcellator for Surgical Procedures outside the indications described in Section **PURPOSE** is prohibited. The user or operator is solely responsible for any adverse consequences resulting from such use.

Requirements to the operator:

Education:

- a university-educated physician

Health:

- good physical and psychophysical condition

Language proficiency:

- English or
- the language of the Operation Manual.

Medical application:

Instruments for to the TCM 3000 BL - Morcellator for Surgical Procedures are designed solely for use in special conditions of surgical emergency units in medical facilities by medically qualified persons. The MP is stable during operation at the temperature of 10°C to 30°C and relative humidity of 30% to 80%. The MP is resistant to corrosion during operation, transportation and storage. The MP is resistant to biological liquids and body tissues it comes in contact with during operation. The MP is resistant to pre-sterilization cleaning, disinfection and sterilization. **The operator must wear medical gloves and a surgical mask while using this medical product.**

7. General information about composition of the TCM 3000 BL Morcellator for Surgical Procedures with Appliances


Fig. 1. The TCM 3000 BL - Morcellator for Surgical Procedures with Accessories

Fig. 1 shows general view of the TCM 3000 BL - Morcellator for Surgical Procedures with Accessories, for which this MP Instruments for the TCM 3000 BL - Morcellator for Surgical Procedures are intended. It consists of a compact TCM 3000 BL - Morcellator control unit (1) with a simple and intuitive panel to control the speed of the attached instruments, a foot switch with precise speed control (3), a gear unit (4), an ergonomic handle, attached to the gear unit (5) and an E-Motor 21 electric motor (6).

Fig. 2 shows how to assemble the TCM 3000 BL - Morcellator for Surgical Procedures before starting its use.


Fig. 2. How one can assemble the Morcellator.

8. Information about MP

To suit the needs of every procedure the model range of the Instruments for the TCM 3000 BL - Morcellator for Surgical Procedures consist of instruments with three different diameters. They are available in Ø12, Ø15 and Ø20 mm versions. By carefully composed and aligned instruments the system as a whole delivers excellent results at morcellation and guarantees for maximum safety at operation. Parts of Instruments to the TCM 3000 BL - Morcellator for Surgical Procedures are described in Appendix 2.

Table 1. Intended use of the parts of Instruments to the TCM 3000 BL - Morcellator for Surgical Procedures.

№	Part	Intended use
1.	Protective tube, Ø 12 mm Ø 15 mm Ø 20 mm	To insert the cutting tube and protect tissues during the procedure
2.	Obturator, Ø 12 mm Ø 15 mm Ø 20 mm	To forward the cutting tube to the tissues to be removed
3.	Cutting tube, Ø 12 mm Ø 15 mm Ø 20 mm	To resect the tissues to be removed
4.	Trocar cannula, Ø 12 mm	To insert the cutting tube and protect tissues during the

	Ø 15 mm	procedure
5.	Myoma drill, Ø 10 mm	To conduct myomectomy (myoma removal)
6.	Grasping forceps, Claw gripper, Ø 10 mm	To grasp the resected tissue for its further removal
7.	Grasping forceps, Tenaculum gripper, Ø 10 mm Ø 5 mm Ø 13,5/15,5 mm	To remove small tissue pieces left

Obturator

5053nou Obturator, Ø 20 mm

5152nou Obturator, Ø 15 mm

5151nou Obturator, Ø 12 mm

Table 2. Technical specifications: obturator

Nose shape	cone
Taper angle	max 30°
Applied part blunting radius	max 0.3 mm

How to use obturator

Obturator has a conical nose part. When the obturator is inserted through a trocar cannula or protective tube, it gradually moves tissues apart.

Fig. 3. General view of the MP

Fig. 4. Labeling on the MP

Labeling is made by laser on the non-applied part of the MP. Labeling contains the following information:

- MP ref. number 5152
- Lot number Lot 1937C1602R
- CE certification

Fig. 5. Obturator case

Fig. 6. Labeling on the obturator case

NOUVAG AG
St. Gallerstrasse 23-25
CH-9403 Goldach

Manufacturer's name and address

Typ : Obturator 15 mm für Morzellator
Type : Obturator 15 mm for morcellator

MP description

REF 5152nou MP Ref. number

LOT 1937C1602R

Lot number

CE mark

Non-sterile

Autoclaving at 135°C max

Suitable for thermal disinfection

Indication of the country of manufacture – Switzerland

Qty.: 1

Quantity of MPs in a pack: pc.

2D data matrix (QR-code)

Cutting tube

5156nou Cutting tube, Ø 20 mm

5155nou Cutting tube, Ø 15 mm

5154nou Cutting tube, Ø 12 mm

Table 3. Technical specifications: cutting tube

Blade thickness, max	0.4 mm
Width of the cutting edge, max	0.2 mm
Fitting diameter to fix in a gear unit, max	20 mm

Cutting tube (Blade)

High-tech cutting tubes manufactured by Nouvag AG, fulfill all requirements to equipment intended for surgery procedures in Laparoscopic Hysterectomy as well as Myomenuklation. High class refinement of the blade of the cutting tube enables multiple use with consistent cutting effect, without having to sharpen the blade.

Fig. 7. How the tube works

2 pilot pins (locking pieces) on the cutting tube are designed to fit it in handle of the gear unit. They should clamp in the gear unit slots to insert the instrument into its working position.

Fig. 8. General view of the MP

Fig. 9. Labeling on the MP

Labeling is electro-eroded directly onto the non-applied surface of the MP. Labeling contains the following information:


- Information about manufacturer: NOUVAG AG
- MP Ref. number: 51555

Fig. 10. Cutting tube pack

Fig. 11. Labeling on the cutting tube pack

	Manufacturer's name and address
Typ : Klinge 15 mm für Morzellator Type : Blade 15 mm for morcellator	MP description
REF 5155nou	MP ref. number
LOT 7406C1511R	MP lot number
	CE mark and notified body
	Non-sterile
	Autoclaving at 135°C max
	Suitable for thermal disinfection
	Indication of the country of manufacture – Switzerland
Qty.: 1	Quantity of MPs in a pack: pc.

Stainless steel

Material of the MP: stainless steel

2D data matrix (QR-code)

Protecting tube

5139nou Protecting tube, Ø 20 mm

5138nou Protecting tube, Ø 15 mm

5137nou Protecting tube, Ø 12 mm

Table 4. Technical specifications: protective tube

Positioner, protrusion, max	5 mm
Clamping nut, length x diameter, max	50 mm x 30 mm


Fig. 12 Tissue positioner

Positioner at the tube's distal end performs 2 functions: protection of underlying tissues during morcellation, and tissue "supply", which allows continuous performance of morcellation. Also prevents tissue pieces from rotating at morcellation.

The tissue positioner acts also as a guiding element at morcellation

During insertion the cutting edge may be retracted inside of the protective tube to prevent complications (CUT-NO-CUT function)

Fig. 13. General view of the MP.

Fig. 14. General view of the MP.

Fig. 15. Labeling on the MP.

Labeling on the MP is made by electric discharge machining. Labeling contains the following information:

- Information about manufacturer: NOUVAG AG
- MP diameter 15 mm

icon "Autoclaving at 135°C max"

icon "Suitable for thermal disinfection"

Fig. 16. Protective tube casing



Fig. 17. Labeling on the protective tube casing.

Manufacturer's name and address

Indication of the country of manufacture – Switzerland

Typ : Schutzrohr 15 mm
Type : Protection sleeve 15 mm

MP description: Protective tube, diameter 15 mm

REF 5138nou

MP ref. number

LOT 1936C1602R

MP lot number

CE mark and notified body

Non-sterile

Autoclaving at 135°C max

Suitable for thermal disinfection

2D data matrix (QR-code)

Qty.: 1

Quantity of MPs in a pack: pc.

Stainless steel

Material of the MP: stainless steel

Trocar cannula

5142nou Trocar cannula for accessories, Ø 15 mm

5141nou Trocar cannula for accessories, Ø 12 mm

Table 5. Technical specifications: trocar cannula

Force (moment) for opening/closing	0.2-0.4 N, with two fingers
Holder diameter x length, max	35 mm x 50 mm
Connecting thread	M24

How to use trocar cannula

Optimal spatial position of troacar cannulas is important for succesful performance of laparoscopy procedures. If more than one trocar cannula is used, they should be positioned 7-10 cm from each other (or palm width). Trocar cannulas should be positioned in a triangular form, aimed to the target organ with its vertex. This allows a surgeon use both hands along the laparoscope axis. Avoid such position of trocar sleeves when the surgeon and its assistants will have to work contrary to visual axis of the laparoscope.

Proximal end of the trocar cannula has a valve ensuring minimum gas loss during insertion and removal of the instrument. Manual regulation of the valve allows removing any objects from the surgical area or introduce necessary instruments to the abdominal cavity.

Fig. 18. General view of the MP.

During insertion the cutting edge may be retracted inside of the protective tube to prevent complications (CUT-NO-CUT function)

Fig. 19. Valve to open the trocar cannula.

Fig. 20. How the trocar cannula valve works.

Fig. 21. How to disassemble/assemble the trocar cannula

Fig. 22. Labeling on the MP

Labeling is electro-eroded directly onto the non-applied surface of the MP. Labeling contains the following information:

- Information about manufacturer: NOUVAG AG
- MP Ref. number: 5142
- Lot number 1939C1602R
- MP diameter 15 mm

icon "Autoclaving at 135°C max"

icon "Suitable for thermal disinfection"

- CE certification mark.

Fig. 23. Trocar tube pack.


Fig. 24. Labeling on the trocar tube pack.

NOUVAG AG
St. Gallerstrasse 23-25
CH-9403 Goldach

Manufacturer's name and address

Typ : Trokarhulse 15 mm
Type : Trocar sleeve 15 mm

MP description

REF 5142nou

MP ref. number

LOT 1939C1602R

MP lot number

CE mark and notified body

Non-sterile

Autoclaving at 135°C max

Suitable for thermal disinfection

Indication of the country of manufacture – Switzerland

Qty.: 1

Quantity of MPs in a pack: pc.

Stainless steel

Material of the MP: stainless steel

2D data matrix (QR-code)

5193 Mioma drill, Ø10 mm

Fig. 25. General view of the MP.

Table 6. Technical specifications: myoma drill

Drill	Right-hand, operating length NMT 30 mm
Handle	6 edges, length 100 mm, diameter 35 mm
Dimensions, max	Length 330 mm, diameter 10 mm

Myoma drill is inserted into surgery area through the trocar cannula during tissue removal (myomectomy) using myoma drill through the trocar cannula.

Fig. 26. Labeling on the MP.

Labeling is electro-eroded directly onto the non-applied surface of the MP. Labeling contains the following information:

- Information about manufacturer: NOUVAG AG
- MP Ref. number: 5193
- Material of the MP Germany Stainless
- CE certification mark.


Fig. 27. Myoma drill pack.

Fig. 28. Labeling on the myoma drill pack

5193	MP Ref. number
Myomborher D=10mm, NL=330mm	MP description
Anzahl: 1	Number of MPs in a pack: 1 pc.
Datum: 26.02.2016	Date of manufacture
Nouvag AG St. Gallenstrasse 22-23 CH-8400 Grenchen Tel.: +41 71 642 45 50	Manufacturer's name and address
Nouvag GmbH St. Gallenstrasse 22-23 CH-8400 Grenchen Tel.: +41 71 642 45 50	Name and address of distributor in the EU.

Gripper

No. 5195

Claw gripper grasping forceps, Ø10 mm/420 mm

No. 5196

Tenaculum gripper grasping forceps, Ø10 mm/420 mm

No. 5194


Tenaculum gripper grasping forceps, Ø5 mm/420 mm

No. 5192

Tenaculum gripper grasping forceps, Ø13.5 mm/420 mm

Table 7. Technical specifications: gripper

Force (moment) for smooth opening/closing	0.2-0.68 N, with two fingers
Jaws free travel angle	does not exceed 25-35°
Head protrusion above surface	max 0.5 mm, riveted-over
Blunting radius	max 0.1 mm

Cutting tube cuts tissue for removal. Claw gripper grasping forceps, diameter 10 mm, are introduced through the tube. Cut tissues are grasped by claw gripper forceps, and removed with them. Ensure atraumatic, solid grasping of tissues.

How to use the gripper (example)

Fig. 29. General view of the MP

Fig. 30. Labeling on the MP

Labeling is electro-eroded directly onto the non-applied surface of the MP. Labeling contains the following information:

- Information about manufacturer: NOUVAG AG
- MP Ref. number: 5192
- Material of the MP Germany Stainless
- CE certification mark.

Fig. 31. Gripper pack.

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


Non-sterile

Date of manufacture

All surgical instruments contact with a patient. **Operator must wear medical gloves and a surgical mask while using this medical product (not included in the scope of delivery).** Materials used for MPs are provided in Table of the Appendix 2. Mechanical characteristics of materials are provided in the Table 8.

Table 8. Mechanical characteristics of materials.

Material	Strength limit	Yield limit	Elongation	Hardness	Roughness (µm), max	Density g/cm ³	Water absorption at 23°C, %	Permissible air operation temperature, °C
Stainless steel AISI 304	≥515 MPa	≥205 MPa	≥40 %	≤201 HB	Ra 0.8			
Stainless steel 1.4305	≥505 MPa	≥215 MPa	≥40 %	240 HB	Ra 0.8			
Polyester-amide PEI MT blue		105 MPa			Ra 0.8	1.27	1.35	-50...170

9. CONSTRUCTION AND OPERATION

1. Principles of Operation

Instruments to the TCM 3000 BL – Morcellator for Surgical Procedures consists of an electric motor, which drives a cutting tube (blade). You can adjust speed and direction (to the right or to the left) of cutting tube rotation. Cutting tubes with blades at the distal end are inserted through an abdominal wall that has been enlarged by an Insufflator (user-supplied), by the means of a protective tube with an obturator or through a trocar cannula. The tissue that is to be removed is guided to the cutting tube, cut off and removed through the tube. The cutting tube is installed in the gear unit attached to the electric motor which is located in an ergonomic handle. The cable connects the electric motor with the control panel displaying various morcellator operation modes for a physician. The physician adjusts morcellator speed using the control panel or a foot switch.

The set of morcellation instruments includes cutting tubes of various diameters, protective tubes, trocar cannula, laparoscopic instruments (gripper, a myoma drill).

2. Application

Instruments to the TCM 3000 BL – Morcellator for Surgical Procedures are used in surgical gynecology.

3. Method of application

Surgical procedure of morcellation is conducted under local or general anesthesia. The patient body shall be placed in Trendelenburg position.

The cutting tube of the morcellator is inserted through an abdominal wall with an obturator placed into the protective tube or through a trocar cannula. By rotating rapidly (rotation speed is adjusted by the foot switch) the cutting tube cuts off the tissue that is to be removed. The cut off tissue is grasped with the claw gripper. The tip of the cutting tube is always located inside the trocar cannula or the protective tube in the lower abdominal wall to avoid deflection or wandering of the cutting tube inside the abdominal cavity. The cut off tissues fall onto the cutting edge of the cutting tube without sucking inside. Large segments of the cut off tissue (10 mm long and larger) are removed with the gripper through an outlet of the protective tube or trocar cannula and placed into a sterile tray.

After large tissue segments are removed with Tenaculum gripper small pieces left are to be removed. The patient shall be moved out of Trendelenburg position so that intra-abdominal liquid starts coming into the lower pelvis. There intra-abdominal liquid is aspirated and replaced with heparinized Ringer's normal saline solution.

Due to effective cutting off the structure of removed tissues is just slightly damaged ensuring high reliability of histological analysis.

This manual describes three different options for removing tissue by a laparoscopic procedure using the MP **Instruments to the TCM 3000 BL – Morcellator for Surgical Procedures**.

Option 1: Tissue removal by morcellation with a protective tube

Instruments (cutting tubes, gripper) are inserted through a protective tube sealed by a valve system. The valve can be covered by a thumb for removing the gripper from the sealed protective tube to prevent gas loss.

Option 2: Tissue removal by morcellation with a trocar cannula


Instruments (cutting tubes, gripper) are inserted through the trocar cannula. Trocar cannula contains a valve that prevents gas loss when withdrawing the instruments.

Option 3: Tissue removal (myomectomy) with a myoma drill through the trocar cannula

Myoma drill is inserted through the trocar cannula. Trocar cannula can be used with or without a cutting tube. No gas is lost after removal of instruments with the trocar cannula.

4. MP overview


Protective tube with obturator

- 1. Protective tube
- 2. Lock nut
- 3. Nose for tissue positioning (tissue positioner)


- 4. Obturator

Cutting tube


- 1. Cutting edge
- 2. Locking pins

Trocar cannula

- | | |
|--------------------|---|
| 1. Protective tube | 4. Valve |
| 2. Seal cap holder | 5. Main body |
| 3. Valve actuator | 6. Nose for tissue positioning
(tissue positioner) |

Gripper

1. Mouth insert
2. Gripper mouth
3. Sleeve shaft
4. Fixed section
5. Flexible leg
6. Locking device
7. Spring component

} **Handle**

Fig. 33. Configuration of MP accessories

5. Operation plan for Instruments to the TCM 3000 BL – Morcellator for Surgical Procedures Option 1


Fig.34. Configuration of the instruments used.

Option 1. Morcellation with a protective tube

Section Option 1 (Fig.34) represents the configuration with protective tube. At the hysterectomy in laparoscopic surgery the morcellator is inserted into the abdominal cavity with the obturator inserted in the protective tube. After removing the obturator, the cutting tube can be inserted. To expose the cutting tube to the tissue that has to be removed, the protective tube has to be unlocked by holding the gear unit with one hand and performing a counter movement of the protective tube with the other hand.

Option 2. Morcellation with a trocar cannula

When morcellating with a trocar cannula (Fig. 34, option 2) the abdominal wall is kept in safe distance from the operation area by inflating the abdominal cavity. This gives also safe clearance for the surgical operation. The trocar cannula stays flexible while operating, which means it can be rotated against the organ or tissue to be protected. Only then the obturator can be removed and the cutting tube can be inserted.

Configuration of instruments needed for this procedure is shown on the chart (Fig.35).


Fig 35. Image-guided morcellation (endoscope).

Option 3


Fig. 36. Instruments used for option 3.

Option 3. Tissue removal (myomectomy) with a myoma drill through a trocar cannula

Myoma drill is inserted through the trocar cannula. Large segments of the cut off tissue are grasped and removed with the gripper through an outlet of trocar cannula and placed into a sterile tray. The Ø10 mm claw gripper are inserted through the trocar cannula to grasp a myomatous nodule and remove it.

6. Before use

TCM 3000 BL - Morcellator control unit and all required accessories and instruments must be placed on a level nonslip surface. For reliable operation of the device, the latter requires enough space to avoid undesirable interference so remove all the limitations and everything affecting manipulations. The display, keypad and indicator lights must be fully visible at all times. Instruments to the TCM 3000 BL Morcellator for Surgical Procedures should be placed within reach, allowing their comfortable use during surgical manipulations. It is very important to ensure that no objects can fall on the foot switch. The power plug at the rear of the device must be accessible at all times.

6.1. Preparation of MP

The MP is not sterile on delivery. All sterilizable parts must be sterilized before use (see Section **Instrument Preparation**).

If MPs have already been sterilized: when removing them from the sterile packaging, ensure that the sterile packaging is not damaged and that sterility indicator confirms sterility and shelf life of the sterile item has not expired.

The Morcellator (as per option 1 to option 3) must be assembled under sterile conditions (**wear gloves and a mask, place Morcellator and Instruments on a sterile surface**).

6.1.1 Assembly of the morcellator set with protective tube (option 1)

There is a risk of injury to tissue by the sharp edge at the distal end of the cutting tube.

- Handle the cutting tube carefully.
- If the cutting tube is not used, rotate the protective tube to the NO CUT position.
- If the cutting tube is not used, it must be placed on a flat, sterile surface so it cannot fall and cause injury.


- Check cutting tubes.
- Make sure that the cutting edge of the blade is sharp and undamaged (e.g. no cracks, deformation or burrs).
- Use only cutting tubes in perfect condition.

1. Inserting the cutting tube through the gear unit

Remove the seal unit from the gear unit.


Insert the cutting tube through the gear unit. The two guide pins on the cutting tube must fit into the slot on the flange of the gear unit.

2. Fixing the cutting tube to the gear unit

Hold the cutting tube with one hand.

Connect the sealing unit to the gear unit with the bayonet closure with the other hand.


3. Mounting protective tube

Select a suitable protective tube.

Push the protective tube over the distal end of the cutting tube. The guide pin of the protective tube must be in the distal slot of the gear unit.

Lock the protective tube in position with a click.


4. Protecting the cutting tube

The distal end of the cutting tube is very sharp.

To prevent injury, the cutting tube must be covered at all times when it is not in use.

To cover the cutting tube, rotate the protective tube until the words "NO CUT" are visible in the distal slot of the gear unit.


The cutting tube is completely covered by the protective tube.

5. Replacing seal holder

Remove the seal holder from the gear unit by rotating it slightly.


Unscrew the membrane screw from the seal unit.

Check the seal holder before use. Do not use damaged seals (frayed, discolored, yellowed, porous, non-elastic).

Position the edge of the cross seal (1) to point to the inside stop (2) of the seal holder (3) and press the cross seal in.


Position the edge of the membrane seal (4) to point to the slot (5) of the membrane screw (6) and press the membrane seal into the slot.

Screw the membrane screw (6) into the thread of the seal holder (3) to the stop and connect the seal unit to the gear unit.


6. Inserting obturator

Insert the obturator through the seal holder, gear unit and protective tube to the stop.

7. Attaching ergonomic handle

Lift and hold the release lever (1) of the ergonomic handle (2). Push the handle (2) over the motor connection of the gear unit. Release the handle release lever. The handle release lever must lock in position.

- Do not open the handle release lever to remove the ergonomic handle during the procedure.

8. Connecting electric motor (optional)


Insert the electric motor coupling into the bottom end of the handle.

Continue to insert the electric motor until it locks into the gear unit's motor connection.

Align the red dot on the device connection plug of the electric motor with the red dot on the electric motor connection of the control unit.

Plug the electric motor device connection plug into the electric motor connection of the control unit. The plug must lock in position.


9. Connecting foot switch (optional)

Align the red dot on the device connection plug of the foot switch with the red dot on the foot switch connection of the control unit.

Plug the foot switch device connection plug into the foot switch connection of the control unit.

The plug must lock in position.

Functional check

- Check Morcellator for loose parts
- Check that Morcellator is correctly assembled
- Are all labels easily readable?
- Check cutting tube: cutting tube evenly ground, no broken edges or bends?
- Switch on control unit
- Switch on electromotor (press foot switch).

- Run electric motor with Morcellator for about 20 seconds; it must reach the set speed.
- If it does not run evenly, switch off electric motor with Morcellator.
- Switch off control unit.

 **The control unit will only operate if the foot switch has been connected to the control unit.**

- Check that Morcellator is correctly assembled.
- Check that the gear unit runs smoothly; if necessary treat it with Nou-Clean care spray (shall be user-supplied, as it is not included into the scope of delivery).
- Check that the electric motor is clean and runs smoothly; if necessary clean it and treat with Nou-Clean care spray (shall be user-supplied, as it is not included into the scope of delivery).
- Repeat the functional check.

6.1.2 Assembly of the Morcellator set with trocar cannula (option 2)

Using the trocar cannula

 The trocar cannula is used instead of the protective tube. The integrated valve flap prevents the escape of gas when instruments are moved in and out.

The nose of the trocar cannula partially extends beyond the cutting tube to protect organs in the abdominal cavity from the sharp blade. This requires the Morcellator to be rotated under visual supervision (endoscopy) so the nose of the trocar cannula covers the cutting tube blade to protect the organ.

1. Visual check

- Check the trocar cannula for sealing, corrosion and sterility.
- Ensure that it is free from rust, dents and scratches.
- The distal end of the obturator must not be damaged.
- The sealing flap and the seal must not be damaged.
-


Risk of injury

Select blade. Make sure that the cutting edge of the blade is sharp and undamaged (e.g. no cracks or burrs).


2. Inserting the cutting tube through the gear unit

Remove the seal unit from the gear unit.
Insert the cutting tube through the gear unit. The two guide pins on the cutting tube must fit into the slot on the flange of the gear unit..

3. Fixing the cutting tube to the gear unit

Hold the cutting tube with one hand.
Connect the sealing unit to the gear unit with the bayonet closure with the other hand.


4. Mounting trocar cannula

Select a suitable trocar cannula. Open the trocar cannula valve by actuating the lever to prevent damage to the cutting tube by pushing the trocar cannula too far.

Push the trocar cannula over the distal end of the cutting tube. The trocar cannula must enter and be in the distal slot of the gear unit.


Lock the trocar cannula in position with a click.

5. Replacing seal holder

Remove the seal holder from the gear unit by rotating it slightly.


Unscrew the membrane screw from the seal unit. Select a suitable seal holder.


Keep a cross seal and a membrane seal ready.

6. Position the edge of the cross seal (1) to point to the inside stop (2) of the seal holder (3) and press the cross seal in.

Position the edge of the membrane seal (4) to point to the slot (5) of the membrane screw (6) and press the membrane seal into the slot.


7. Screw the membrane screw (6) into the thread of the seal holder (3) to the stop and connect the seal unit to the gear unit.

8. Inserting obturator


Insert the obturator through the seal holder, gear unit and the trocar cannula.

9. Attaching handle

Lift and hold the release lever (1) of the ergonomic handle (2). Push the handle (2) over the motor connection of the gear unit.

Release the handle release lever. The handle release lever must lock in position.

10. Connecting electric motor (optional)

Insert the electric motor coupling into the bottom end of the handle.

Continue to insert the electric motor until it locks into the gear unit's motor connection.

Align the red dot on the device connection plug of the electric motor with the red dot on the electric motor connection of the control unit.


Plug the electric motor device connection plug into the electric motor connection of the control unit. The plug must lock in position.

11. Connecting foot switch (optional)

Align the red dot on the device connection plug of the foot switch with the red dot on the foot switch connection of the control unit.

Plug the foot switch device connection plug into the foot switch connection of the control unit. The plug must lock in position.

Functional check

- Check Morcellator for loose parts
- Check that Morcellator is correctly assembled
- Are all labels easily readable?
- Check cutting tube: cutting tube evenly ground, no broken edges or bends?
- Switch on control unit
- Switch on electromotor (press foot switch).
- Run electric motor with Morcellator for about 20 seconds; it must reach the set speed.
- If it does not run evenly, switch off electric motor with Morcellator.
- Switch off control unit.

The control unit will only operate if the foot switch has been connected to the control unit.

- Check that Morcellator is correctly assembled.
- Check that the gear unit runs smoothly; if necessary treat it with Nou-Clean care spray (shall be user-supplied, as it is not included into the scope of delivery).
- Check that the electric motor is clean and runs smoothly; if necessary clean it and treat with Nou-Clean care spray (shall be user-supplied, as it is not included into the scope of delivery).
- Repeat the functional check.

Risk of injury by electric motor failure

Problem: the electric motor becomes very hot.

Cause: the electric motor was damaged during preparation.

Remedy: do not continue to use the electric motor. **Always have a spare electric motor available.**

6.1.3 Assembly of the trocar cannula with myoma drill and mounting gripper (option 3)

Assembly of the trocar cannula with myoma drill is similar to that of the Morcellator set with the protective tube (option 1) and with the trocar cannula (option 2), except that myoma drill and gripper are used instead of cutting tubes.

7. Use

Before use

Make sure that the product has been correctly prepared and inspected.

Risk of injury

TCM 3000 BL - Morcellator control panel has been designed for operation solely with the instruments manufactured by Nouvag AG. The user is solely responsible if different products are used with the set. Incompatibility may risk injury to patients.

Switching on and off

Switch on the control panel. The display lights up. A short signal tone sounds. The control panel displays the value "200 rpm".

Operating mode (electric motor)

The motor is designed for intermittent operation (1 minute on, 1 minute off for 8 cycles).

7.1 Using the control panel

Setting the speed

The speed can be set manually (within a range of 50 to 100 rpm).

Press ▲ to increase the speed.

Press ▼ to reduce the speed.

The set speed is shown on the display.

Starting electric motor

Press foot switch.

The blade rotates at the set speed.

The electric motor LED lights up.

To stop the electric motor remove foot from foot switch.

Risk of injury and damage due to heat!

The Morcellator may generate heat if operated for a long period.

Remedy:

Use Morcellator only in intermittent mode of "1 minute on, 1 minute off for 8 cycles".

Stop the Morcellator immediately if one or more components of the Morcellator set become very warm.

Precaution when changing parameters. Unusual behavior of instruments while operating may provoke false reactions and jeopardize the patient.

Every new setting and the new behavior of the instrument has to be verified and tested.

Available rotary modes

-  (100 – 400 rpm)
-  (300-700 rpm)
-  (500 – 1000 rpm)

Activating mode setting

Press  to activate mode setting.

The Vario LED lights up.

The Low LED lights up (Default setting)

The display shows «VARIO LOW, 100-400 rpm».

Selecting mode

Press ,  or  to select the corresponding mode. The display shows the mode and the corresponding speed range that can be selected with the foot switch. The corresponding LED lights up.

Starting electric motor

Press foot switch. The electric motor LED lights up. The blade rotates at a speed within the range set by how far down the foot switch is pressed.

To stop the electric motor remove foot from foot switch.

Exit mode setting

Press .

Changing direction of rotation

To change direction of rotation, press .

The display shows the mode and the corresponding speed range that can be selected with the foot switch. The corresponding LED lights up.

7.2 Operation

Risk of death or serious injury!

Unsupervised cutting may injure the abdominal wall and organs and cause the death of the patient.

Remedy:

- Before every procedure make sure that an endoscopic procedure is preferable to a conventional procedure.
- Use only the products described in this manual.
- During the procedure use the protective tube or trocar cannula as a precautionary measure. Follow the instructions for "Option 1: Morcellation with protective tube" and "Option 2: Morcellation with trocar cannula" in Section **Preparation of Components**.
- Use Morcellator under visual control (endoscope) only.
- Use Morcellator only for morcellation of tissue that has been fully prepared and is clearly visible. Nouvag AG accepts no liability for any other applications.
- Nouvag also recommends the use of a second gripper or a similar retaining MP to prevent

uncontrolled movement of large morcellated pieces of tissue. This will require an additional percutaneous access.

- Fix the Morcellator in a ventrolateral position to prevent movements inside the uterus.
- Prevent movement of the lateral blood vessels, kidneys and retroperitoneum.
- Pull the tissue towards the rotating blade for morcellation with the gripper.

Risk of injury due to gas loss!

Gas may be lost during the procedure. This may result in the abdominal wall dropping.

- Select an Insufflator (not included into the scope of delivery) with adequate flow to prevent the abdominal wall from dropping.
- Use the Morcellator with protective tube as in option 1: when removing or replacing instruments, immediately close the proximal opening of the seal unit with the thumbs.
- Use the Morcellator with trocar cannula as in option 2: pull the Morcellator with cutting tube out of the trocar cannula. The trocar cannula valve closes to prevent gas loss.
- After removal of the morcellated tissue, insert the Morcellator into the trocar cannula and repeat the process.

Risk of injury!

The Morcellator system is designed for removing tissue with the greatest possible safety. For this purpose the protective tube and the trocar cannula are both fitted with a valve flap.

Morcellation without protective tube or trocar cannula is prohibited.

Tissue escapes the gripper.

Tissue may escape the gripper during the procedure and fall into the abdominal cavity.

Causes:

- The tissue is not gripped tightly enough.

Remedy:

- Grip the tissue at a different position.
- Press the jaws of the gripper together as tightly as possible.

Risk of injury due to heat generated by friction!

The user must ensure that as little friction heat as possible is generated during a procedure on a patient.

High speeds and high application pressure may cause thermal necrosis of tissue.

Risk of injury!

The distal end of the cutting tube is extremely sharp.

Precautions:

- Handle the cutting tube particularly cautiously and carefully.
- Where possible, have the Morcellator with protective tube in protected position (NO CUT), particularly if the cutting tube is not used for any period.
- Insert and remove the protective tube and cutting tube under visual control (endoscopy).

Risk of injury!

When inserting the obturator through the percutaneous access, there is risk of injury to blood vessels, intestinal loops or the bladder. This can result in significant complications, which will require a Laparotomy.

Risk of damage to the blade of the cutting tube!

If the jaws of the gripper are not completely closed when it is pulled from the protective tube, the blade will be damaged (see fig. 15).


Fig. 15. Use of gripper.

Risk of injury!

Damaged blades and cutting tubes must not be used.

Cutting tubes must not be sharpened.

7.2.1 Option 1: morcellation with protective tube.

1. Hysterectomy

Select the procedure for the hysterectomy in accordance with current medical standards.

2. Inserting cutting tube


Make sure that the obturator is inside the cutting tube.

Make sure that the protective tube covers the blade on the cutting tube (set position to NO CUT; "PROTECTED" is also shown on the black ring between the gear unit and protective tube).

Insert Morcellator (with obturator) into an existing percutaneous access under visual control (endoscopy).


3. Using cutting tube

To use the cutting tube, rotate the protective tube to uncover the blade of the cutting tube (unprotected CUT position is displayed). To cover the blade again, rotate the protective tube (protected NO CUT position is displayed; "PROTECTED" is also shown on the black ring between the gear unit and protective tube).

4. Inserting gripper

Pull the obturator out of the protective tube. Close the proximal opening of the gear unit with the thumbs immediately. Insert gripper carefully through the proximal opening of the gear unit.

10. Removing tissue


Grip the tissue with the gripper under good visual control (endoscopy).

The handle lock is set.

Start the electric motor of the Morcellator set. Pull the tissue towards the rotating blade for morcellation with the gripper.

As the cut progresses pull the gripper with the tissue further out of the Morcellator.

Carefully pull the complete gripper with the separated tissue out of the Morcellator and the patient.

The retracted mouth must not touch the blade of the cutting tube. Repeat the process until the tissue has been completely removed.

Make sure that the proximal opening of the seal unit is closed immediately with the thumb when there is no instrument in the Morcellator to prevent large quantities of gas from escaping.

7.2.2 Option 2: morcellation with trocar cannula

1. Hysterectomy

Select the procedure for the hysterectomy in accordance with current medical standards.

2. Inserting trocar cannula

Insert trocar cannula with obturator into an existing access under visual control (endoscopy).

Remove obturator.

3. Insert the cutting tube through the trocar cannula.

Actuate the valve flap lever of the trocar cannula and hold in this position.



Carefully insert prepared cutting tube with gear unit and seal unit into the proximal opening of the trocar cannula.

Cutting tube must be fully inserted. Release the valve flap lever of the trocar cannula.

The guide pin of the trocar cannula must be in the distal slot of the gear unit.


Lock the trocar cannula in position with a click.

Make sure that the cutting tube remains in the NO CUT position.

4. Inserting gripper

Insert gripper carefully through the proximal opening of the gear unit and trocar cannula.


5. Preparing cutting tube

Rotate the trocar cannula so the CUT position is visible from above.

If the cutting tube is completely inserted, the cutting edge is now only partially covered by the trocar cannula.

6. Protecting tissue

The nose on the trocar cannula, which partially covers the sharp edge of the blade, must protect tissue that must not be injured. Rotate the Morcellator to the desired position for this.

7. Removing tissue


Grip the required tissue with the gripper. The handle lock is set.

Start the Morcellator set motor.

Pull the tissue towards the rotating blade with the gripper.

Pull the gripper out of the patient through the gear unit. The retracted mouth of the gripper must not touch the blade.

Repeat the process until the morcellated tissue has been completely removed.

7.2.3 Option 3: myomectomy

The myoma drill has been designed exclusively for removing tissue from a myoma.

1. Inserting Morcellator

Insert Morcellator with protective tube through the percutaneous access.


Keep gripper ready to remove the tissue from the myoma drill.

2. Inserting myoma drill

Insert trocar cannula through an additional percutaneous access. Open valve flap of trocar cannula and insert myoma drill.

3. Uncovering myoma

Remove tissue from the myoma with the myoma drill.

4. Inserting gripper

Insert gripper through the second access and remove separated tissue from the myoma drill. Remove myoma drill.

5. Removing separated tissue

If the piece of tissue is larger than the diameter of the cutting tube, start the Morcellator set electric motor. Pull the tissue towards the rotating blade with the gripper.

(Smaller pieces of tissue can be removed without morcellation).

Carefully pull the complete gripper out of the Morcellator and the patient. The retracted mouth must not touch the blade of the cutting tube.

6. Removing additional tissue

If necessary, re-insert the myoma drill and repeat the above process.

8. After use

1. Switching off

Switch off the power switch on the rear panel of the device ("O" position).

2. Removing protective tube (when using Morcellator with protective tube)

Rotate protective tube to NO CUT position.

Close the proximal opening of the gear unit with the thumb to prevent gas loss. Hold the protective tube and carefully disconnect gear unit with cutting tube and remove from patient. The protective tube can also be removed from the percutaneous access.

If applicable, remove any other instruments.

3. Removing cutting tube (when using Morcellator with trocar cannula)

Close the proximal opening of the gear unit with the thumb to prevent gas loss.

Hold the trocar cannula and carefully disconnect gear unit with cutting tube from the trocar cannula and remove from patient. The valve flap closes automatically as soon as the cutting tube is pulled out.

The trocar cannula can remain as access throughout the remaining procedure.


4. Disconnecting gear unit connections

Pull the electric motor coupling from the gear unit's motor connection.

Pull the electric motor device connection plug from the electric motor connection of the control unit.


Pull the electric motor out of the handle. Lift and hold the release lever of the handle.

Pull the handle from the gear unit's electric motor connection.

5. Disconnecting foot switch connections

Pull the foot switch device connection plug from the foot switch connection of the control unit.

6. Preparing for preparation at the location of use

Take the product to the preparation area immediately. Prepare the product as specified in this manual.


7. Disassembling trocar cannula

Disassembling trocar cannula (four steps)

Disassemble trocar cannula before decontamination.

Disassemble trocar cannula and rinse in deionized water during disassembly.

1. Unscrew sealing cap holder from the proximal end of the trocar cannula.
2. Remove the internal seal from the sealing cap holder.
3. Unscrew the valve body from the main body.
4. Remove the red O-ring from the valve body.


8. Disassembly of electric motor (optional)

Disassembling electric motor (three steps).

Disassemble electric motor before contamination.

1. Unscrew handpiece holder from electric motor
2. Unscrew motor cap from electric motor
3. Pull motor plug from electric motor.


9. Disassembly of gear unit and Morcellator with protective tube

Disassembling seal unit

Remove membrane screw from seal unit and remove membrane seal.

Remove seal holder from gear unit and remove cross seal.

Risk of injury for the user!

The distal end of the cutting tube is very sharp, handle the cutting tube very carefully.


Removing protective tube and cutting tube

Make sure that the lock nut of the protective tube covers the distal end of the cutting tube (protected NO CUT position is visible).

Hold the gear unit with one hand so the distal end of the cutting tube does not point to the user or any other person.

With the other hand carefully disconnect the protective tube from the gear unit and pull it off.

Hold the cutting tube with the same hand and press it lightly against the gear unit until it releases itself from the handle and slides out smoothly. It can now be pulled out of the proximal opening of the gear unit.

10. SCOPE OF DELIVERY

Instruments to the TCM 3000 BL – Morcellator for Surgical Procedures

1. 5154nou Cutting tube, Ø 12 mm
2. 5155nou Cutting tube, Ø 15 mm

3. 5156nou Cutting tube, Ø 20 mm
4. 5141nou Trocar cannula for instruments Ø12 mm
5. 5142nou Trocar cannula for instruments Ø15 mm
6. 5151nou Obturator, Ø12 mm
7. 5152nou Obturator, Ø15 mm
8. 5153nou Obturator, Ø20 mm
9. 5137nou Protective tube, Ø12 mm
10. 5138nou Protective tube, Ø15 mm
11. 5139nou Protective tube, Ø20 mm
12. 5193 Myoma drill, Ø10 mm
13. 5194 Tenaculum gripper, Ø 5mm
14. 5195 Claw gripper, Ø10 mm
15. 5196 Tenaculum gripper, Ø 10 mm
16. 5192 Tenaculum gripper, Ø 13.5/15 mm


On MP receipt the following actions shall be completed:

- Ensure that the shipping package contains all parts listed in Section Scope of Delivery.
- If any parts are damaged or missing, contact the Nouvag agent from which the MP was ordered or Nouvag AG.

To carry out immediate inspection of the TCM 3000 BL – Morcellator for Surgical Procedures with Instruments compare the received MP with the below pictures of the set components:

Cutting tube

Protective tube

Obturator

Trocac cannula

Myoma drill


Gripper

Fig. 37. Visual check of the scope of delivery

11. CLEANING OF MP COMPONENTS

Schematic flow chart for reesterialization as per EN ISO 17664

General information on disinfection, pre-sterilization cleaning and sterilization of MP

- Reusable MP are subject to re-sterilization and subsequent storage in conditions excluding secondary microbial contamination.
- Sterilization of such MP should be performed immediately before the first use of a MP, and after each patient to prevent cross-contamination.
- Disassembly, washing and rinsing of MP which contacted patients' blood and serum may should be performed in rubber gloves, goggles, medical mask, medical gown and waterproof apron. Personal hygiene should be strictly observed.
- Containers with working solutions of detergents and disinfectants should have caps, clear inscriptions indicating name of a product, its concentration, purpose, date of preparation, end of shelf life. For products ready for use name of the product and its purpose should be indicated, if multiple use of a product is permitted, date when its use began should be additionally provided.
- During disinfection and pre-sterilization cleaning in chemical solutions MP are submerged into a product solution to fill all channels and cavities. Dismountable MP are submerged disassembled, instruments with locking parts should be soaked open, with making several operational movements with such MPs.
- Disinfection by wiping may be used for MDs that does not immediately contact a patient, or where their design does not permit using the submersion method.
- After disinfection reusable MDs should be washed from residues of disinfectants according to recommendation provided in instructions for use of the appropriate disinfectant.
- Degree of cleaning is confirmed by a test: amidopyrine test for blood residues and phenolphthalein test for detergent residues. Quality control should be performed every day. Control should be performed on 1% of simultaneously treated instruments of each position, but not less than three. Check both external surfaces, and internal surfaces of tools! Document proving fact of such inspection is a pre-sterilization quality log.
- This MP is sterilized by physical method (autoclaving).
- Sterilization should be performed in a functionally separated premise – "clean room".
- MP should be transported into the clean room in closed sterilizable containers.
- Personnel performing pre-sterilization cleaning and disinfection of the MD should wear special garments, and use mask, goggles and rubber gloves, waterproof apron over the gown.
- Sterilized MPs should be kept under conditions excluding any secondary microbial contamination in a special cabinet, protected from dust, humidity and temperature fluctuations.
- Shelf life should be indicated on package and determined by type of a packing materials, according to instructions for its use.
- It is not allowed to use sterilized MPs with expired term of storage after sterilization.
- Sterilization of MDs should be recorded in a log.

General recommendations on preparation of MP

- Preparation procedures include cleaning, washing disinfection, washing and drying.
- Instruments should be cleaned only in a cleaning solution (e.g. Schulke perform or 3% DESOMEDAN ID) with observance of personal security measures. Strictly adhere to the manufacturer's instructions for use when using solutions!
- Use only demineralized water (VE-water) for rinsing and washing.
- Use only daily prepared working solutions. For strong contaminations it is recommended to replace solution frequently.
- It is recommended to use for cleaning soft lint-free cloths, plastic brushes or washing guns. Use of metal items is forbidden. Treatment should be performed in plastic or enamel-lined (without damaged enamel lining) containers, closing with covers.
- Manual external cleaning is performed by a disposable cloth, disinfected brushes. Brush cleaning is performed over a solution surface to prevent from splashing contaminated fluid.
- For effective cleaning open hinged instruments (gripper) to minimize overlapping surfaces. Dismountable instruments should be disassembled.
- Hollow instrument is cleaned by a rinsing gun (rinsing pressure 3 bar- 4 bar).
- MP disinfection is performed jointly with cleaning.
- Machine cleaning and disinfection are performed in a disinfector wash (e.g., thermal disinfector Miele G7835 CD) with performance of the thermal disinfection procedure.
- For the avoidance of damages MP parts should be properly fixed in the disinfector wash. Loading MPs should ensure thorough washing, including for internal cavities of instruments through the appropriate valves.
- Instrument drying is a part of machine cleaning and disinfection.
- Upon completion instruments should be removed manually, using disposable sterile gloves treated with an antiseptic.
- Inspection is performed visually. Critical parts, such as handles, hinges and cutting surfaces should be controlled with special attention. All instruments with internal channels should be checked for passability. If any instruments do not satisfy the inspection, they should be replaced.
- Insufficiently cleaned instruments should be re-cleaned and then washed thoroughly.
- Instruments with micro cracks near hinges and/or damaged, bent or otherwise worn instruments should be replaced, as they cannot or can improperly perform their functions.
- If instruments require treatment actions, the appropriate treatment actions are performed before function test, as a rule.
- Before function test lubricate ball and threaded instruments. Proper operation of instrument should be ensured by the check. For this purpose any disassembled instruments should be assembled. If necessary, after the check instruments may be disassembled again for sterilization. Perform disassembly/assembly according to manufacture's instructions.
- It is recommended to pack MPs into a sterilization pack. Follow manufacturer's instructions for use and time of sterilization validity!
- Package should ensure the possibility for labeling with the following information:
 - sterilization date,

- packer,
- shelf life,
- contents.
- If sterilization is performed without craft packs, storage time for sterilized MPs is: 6 hours if sterilized in baskets, and 20 days, if sterilized in boxes with filters; craft packs are stored up to three days in sterile conditions. Any items sterilized unpacked should be used immediately after sterilization.
- Sterilize in a vacuum autoclave (class B or S according to EN 13060) at 135°C. Sterilization parameters

Sterilization mode					
Steam pressure in sterilization chamber, MPa (kgf/cm ²)		Working temperature in sterilization chamber, °C		Time of sterilization hold, min	
rated value	limit value	rated value	limit value	rated value	limit value
0.20	±0.02	135	±2	5-25	±2

- The manufacturer recommends using indicator tapes to control sterilization (e.g. indicator tapes Comply™ by 3M).
- The manufacturer confirms validation of the re-sterilization procedure by a Certificate (see Appendix 3).
- Returned instruments are MP, which – regardless of whether they have been used – are returned to the manufacturer.
- Possible causes of return may be detected MP fault, failure during the warranty period or the need for maintenance. If a MP is returned to the manufacturer or its authorized representative, they may be sent back to the user only if:
 - passed cleaning and disinfection according to manufacturer's instructions, are dry and declared "sanitarily clean", and/or
 - clearly stated as not contaminated, and properly packed.
- Decontamination of returned products – as in common circulation – should be performed without delay to avoid damaging of an instrument. Decontamination should not be performed if it may change or destruct the product, falsify results of a test or may a test impossible. In case of doubts please consult with the manufacturer.
- The following information should be provided in accompanying documents for each individual MP returned to the manufacturer:
 - intended use of the MP;
 - decontamination method;
 - date of treatment;
 - last name of the person who performed the treatment;
 - reason for return.

Cleaning of components of the MD “Instruments to the TCM 3000 BL – Morcellator for Surgical Procedures” is performed as follows:

Instruments preparation

Risk of infection!

Prepare the product before its first use and before every subsequent use as specified in this manual.

Inadequate and/or incomplete preparation of the product may cause an infection of the patient. Observe the following specific instructions:

- At unpacked steam sterilization of single items they should be fixed in the tray to prevent them from rolling around, causing damage.
- Allow Accessories to the 3000 BL – Morcellator for Surgical Procedures to cool down after steam sterilization!

Immersing instruments in cleaning solution

- Never exceed the maximum concentration and exposure time specified by the manufacturer of the cleaning and disinfection solution.
- Air bubbles must not adhere to the instruments.
- All components of the instruments must be completely immersed in the cleaning solution.
- All lumina of the instruments must be completely filled with cleaning solution and free from bubbles.

- Do not use cleaning agents that contain solvents.
- If sterilized MP components are not used immediately, a sterilization indicator must be placed on the packaging and it must be labeled with the sterilization date.

- Perform cleaning, disinfection and sterilization after every treatment.
- To autoclave single items, such as cutting tubes, Nouvag AG recommends the use of sterilization packages to prevent of damage.

- Make sure that sterilization packaging is no more than 80 % full.
- Autoclave material for at least five minutes at max. 135 °C.
- If sterilized MP components are not used immediately, the sterilization packaging must be labeled with the sterilization date.
- Nouvag AG recommends including a sterility indicator.

11.1 Preparation of the trocar cannula

Essential requirements

Avoid long wait times until preparation due to risk of drying and corrosion. Maximum time between use and preparation: 2 hours. The trocar cannula is designed for 500 sterilization cycles.

Disassembly

Before decontamination disassemble the trocar cannula as described in section "**Disassembling trocar cannula**". Remove seal cap, internal seal and O-ring.

Manual cleaning

- Rinse and clean all surfaces and lumina thoroughly with cleaning solution and suitable brushes.
- Place trocar cannula in an ultrasonic bath with cleaning solution, if necessary.
- Rinse all components of the trocar cannula thoroughly with cleaning gun and deionized water.

Machine cleaning and disinfection

- Arrange all components of the trocar cannula in the cleaning machine.
- Use a fine-mesh basket for small parts, such as seals and O-rings.
- Run a standard machine cleaning and thermal disinfection cycle.

Autoclaving

Autoclave in vacuum autoclave (class B or S as per EN 13060) at 135°C for at least five minutes*. Do not exceed the maximum load of the sterilizer when sterilizing multiple instruments in one sterilization cycle. Autoclaving without vacuum must be followed by a drying phase. Allow the trocar cannula to dry in a bag for at least one hour at room temperature with the paper side up.

** Temperature holding times conform to the national guidelines and standards with the maximum holding time 25 minutes.*

Assembly

- Make sure that all parts of the trocar cannula are completely dry after preparation.
- The flap in the valve body must be open when screwing the spacer with internal seal to the valve body.
- Assemble the trocar cannula in reverse order of disassembly (see Section "**Disassembling trocar cannula**").

11.2. Preparation of cutting tubes, protective tubes, obturators, grippers and myoma drill

Clean all products with suitable cleaning brushes.

Reprocessing restrictions

The end of the product service life is normally determined by wear and tear and damage through use. The cutting tubes are designed for a service life of 20 cycles.

At the location of use

Remove surface soiling with a napkin/paper towel.

Storage and transport

No special requirements - avoid long periods before preparation to prevent the danger of drying-out and corrosion, maximum time between use and preparation is 2 hours. If sterilized instruments are not used immediately, the sterilization packaging must be labeled with the sterilization date. Nouvag AG recommends using sterility indicators.

Manual cleaning

- Use disinfectant solution, neutral cleaning agent, suitable cleaning brushes and demineralized water. Wipe all surfaces thoroughly with disinfectant (exposure time as specified by the disinfectant manufacturer).
- Immerse all parts in the neutral cleaning agent (exposure time as specified by the cleaning agent manufacturer). Use suitable brushes that can reach all parts of the lumen for cavities.
- Rinse all parts thoroughly under running water.

Automatic cleaning and disinfection

Equipment: Cleaning/disinfection MP with a special load carrier which facilitates the connection of the instruments to the cleaning/disinfection MP and the channels. Use a fine-mesh basket for small parts. Use neutral cleaning agent at the concentration recommended by the manufacturer.

- Place the cutting tubes, protective tubes, myoma drills, grippers and ergonomic handle in the load carrier (ensure that the channels can be rinsed).
- Place seal holder, membrane screw, membrane seal, cross seal, obturator in the fine-mesh basket and close the basket.
- Set a cleaning cycle that offers sufficient cleaning and rinsing. Perform the final rinse with deionized water.
- Perform a 10-minute rinse cycle at 93°C to facilitate thermal disinfection.
- When removing the products, check whether visible contamination remains in the gaps and the grooves. If necessary, repeat the cycle or clean manually.

Drying

If a drying program is not provided by the washer-disinfector, all parts must be dried in a drying cabinet at 60°C.

Inspection and maintenance

- Perform a visual inspection to check for damage, corrosion and wear and tear.
- Check cutting edge for damage and rotary components for bending.
- Check that rotary components move smoothly.
- Check seals and seal holders for elasticity, cracking, deformation and porosity.

Packaging

- Storage of single unit: Pack the parts in individual packaging for sterile items. The bags must be large enough to ensure that the seal is not subject to strain. When using sterilization container, observe the instructions of the manufacturer.
- Storage: Sort the instruments on trays intended for this purpose or place them on all-purpose sterilization trays.

Risk of damage

Blade of the cutting tube may damage the sterile packaging.

Sterilization of gripper, cutting tubes, protective tubes, obturator, myoma drills

Autoclave in the vacuum autoclave (Class B or S as per EN 13060) at 135°C for at least 5 minutes*. When sterilizing several instruments during one sterilization cycle, do not exceed the maximum sterilizer load. A drying cycle must be added in the case of autoclaves without a post-vacuum function. Allow the instruments to dry in the bag for at least 1 hour at room temperature with the paper side facing upwards.

** Temperature exposure times are based on country-specific guidelines and standards. Maximum exposure time: 25 minutes.*

The MP is not sterile on delivery. All sterilizable parts must be sterilized before use (see Section Instrument Preparation).

If MPs have already been sterilized: when removing them from the sterile packaging, ensure that the sterile packaging is not damaged and that sterility indicator confirms sterility and shelf life of the sterile item has not expired.

The Morcellator (as per option 1 to option 3) must be assembled under sterile conditions (wear gloves and a mask, place Morcellator and Instruments on a sterile surface).

12. Maintenance and repair

This medical product is non-repairable. If MP components have any defects, they should be replaced with new ones.

12.1 General maintenance

Before every use perform a visual inspection of the MP to check for damage and wear and tear. Check cutting edges and rotary components for bending. Check rotary components for gaps.

Users can perform the following maintenance procedures themselves:

Maintenance, protective tube

- Replace damaged protective tubes with new ones.
Protective tube, Ø 12 mm (Ref. 5137nou)
Protective tube, Ø 15 mm (Ref. 5138nou)
Protective tube, Ø 20 mm (Ref. 5139nou)

Maintenance, obturator

- Replace damaged obturators.
Obturator, Ø 12 mm (Ref. 5151nou)
Obturator, Ø 15 mm (Ref. 5152nou)
Obturator, Ø 20 mm (Ref. 5153nou)

Maintenance, trocar cannula

- Replace damaged trocar cannulas
Trocar cannula Ø 12 mm (Ref. 51484nou), 1 item per package
Trocar cannula Ø 15 mm (Ref. 51502nou), 1 item per package.

Maintenance, trocar seal holder

- Replace damaged seals with new ones.

13. Safety information

Your safety, the safety of your team, and of course that of your patients is very important to us. It is therefore essential to bear the following information in mind

13.1 Modification and misuse

- Replacement of MP parts with external ones or modification of **Instruments to the TCM 3000 BL – Morcellator for Surgical Procedures** is prohibited. The manufacturer is not liable for any damages resulting from unauthorized modifications or manipulations. The warranty will be cancelled.
- Use of the MP **Instruments to the TCM 3000 BL – Morcellator for Surgical Procedures** outside the indications described in Section **PURPOSE** is prohibited. The user or operator is solely responsible for any adverse consequences resulting from such use.

13.2 Use precautions

In order to avoid risks and complications you must follow the below precautions:

- The cutting tube of the morcellator must be moved against the abdomen in one direction only;
- Cut off tissues to be removed must be extracted through a hollow tube;
- Uncontrolled movement of morcellator tubes and/or instruments towards lateral vessels, bowels or retroperitoneum must be excluded;
- The morcellator must be ventrolaterally positioned during the whole procedure in order to avoid bowels or vessels injuries;
- Sharp edge of the morcellator cutting tube must be visible during the whole procedure;
- The MP must be used by qualified and experienced surgeons only;
- The MP must be operated as per the Operation Manual and be functional at all times;
- Rotating tube must be set into motion under continuous visual control only.

13.3 List of operational hazards, risk analysis

The manufacturer runs the system of risk analysis and monitoring. Risk analysis is performed during post marketing monitoring of results of MP use. Risk analysis is a part of the device technical documents.

Residual risks of use of **Instruments to the TCM 3000 BL – Morcellator for Surgical Procedures** include injuries of vessels and other anatomical sites.

13.4 Safety requirements, operation conditions and rules

This product is designed solely for use by medically qualified, experienced and trained persons.

Surgical procedure using **Instruments to the TCM 3000 BL – Morcellator for Surgical Procedures** must be completed under continuous supervision of medically qualified staff able to react to any possible adverse events and complications immediately.

To avoid injuries no MP manipulations must be done while the electric motor is switched on.

When the MP is not used, the protective tube must be pulled out to completely cover the blade of the cutting tube to avoid injuries.

Ambient conditions

	Operation:
Relative humidity:	30% - 80%.
Temperature:	10°C - 30°C
Atmospheric pressure:	800 hPa - 1060 hPa

The medical product is always safe for the environment while being used.

13.5. Safety check

Safety check for the MP **Instruments to the TCM 3000 BL – Morcellator for Surgical Procedures** is obligatory before each use. The MP is non-repairable and should be disposed after expiration of its service life.

If you have any questions regarding MP maintenance and repair as well as after-sales servicing please contact our authorized representative in the Russian Federation: **ООО MedTestConsulting, Russia, 123181, Moscow, 1 Kulakova Str., bldg. 1, office 193, phone/fax +7 (495) 750-29-01, +7 (495) 729-77-14, E-mail: medtestconsult@mail.ru.**

14. MEDICAL PRODUCT CLASSIFICATION

Purpose	Cutting instruments
Frequency of application	Multiple application
Operating mode	device for short-time modes of operation
Scope of use	Product for professional use in healthcare institutions
Possibility of repairs	non-repairable
Type of perceived mechanical impact	portable, for use in healthcare institutions
Type of human body contact	Short-term contact (less than 24 h) with milieu interior

15. COMPLIANCE WITH INTERNATIONAL STANDARDS

Instruments to the TCM 3000 BL – Morcellator for Surgical Procedures comply with Council Directive 93/42/EEC on Medical Products:

EN ISO 10993-1:2009	Biological evaluation of medical devices. Part 1: Evaluation and testing within a risk management
ISO 3964	Dental handpieces. Coupling dimensions.

EN ISO 17664-2004

Sterilization of medical devices. Information provided by manufacturer for treatment of re-sterilizable medical devices.

ISO 13485:2012


Design and development, production, distribution and servicing of medical devices

CE certificate

Compliance with Council Directive 93/42/EEC on Medical Devices

16. PACKING AND LABELING

Table 9. Explanation of symbols

	Important information		Autoclavable at max. 135° C
	Warning		Suitable for thermal disinfection
	Manufacturer		Keep dry
	Keep away from direct sunlight		Non-sterile
	Certified by the Canadian Standards Association (CSA).		Electrical and electronic devices that have reached the end of their service life comprise hazardous waste and may not be disposed of together with household waste. Local disposal regulations apply.
	CE mark with Notified Body		Symbol indicating the order number
	Intermittent operation, 1 min on, 1 min off for 8 cycles		Symbol indicating the serial number and the date of manufacture (year/month)
	2D data matrix (QR-code)		Biohazard
	Non-sterile		Indication of country of manufacture: Switzerland

Individual package of the MP

Fig. 38. Individual package of the trocar cannula

Fixpac box size 255x155x100 mm

Internal PE-LD bag size 200x500x0.075 mm


Fig. 39. Individual package of the obturator

DP 40 250 PP transparent packing tubes, size Ø45x255 mm.

Foam inserts size Ø45x10 mm.

Fig. 40. Individual package of the protective tube


DP 30 200 PP transparent packing tubes, size Ø35x210 mm.
Foam inserts size Ø35x10 mm.

Fig. 41. Labeling on individual package.


Manufacturer's name and address


Information on the country of origin: Switzerland

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Type : Protection sleeve 15 mm Name of the MP: Protection sleeve 15 mm

Reference number of the MP.

Number of MP lot

CE mark with notified body

Non-sterile


Autoclaving at 135°C max

Suitable for thermal disinfection

2D datamatrix (QR-code).

Number of MPs in a package: 1 item.

Material of MP: stainless steel

Fig. 42. Labeling of Instruments for the MP.

Labeling is electro-eroded directly onto the non-applied surface of the MP. The marking includes the following information:

- information on the manufacturer: NOUVAGAG
- Reference number of the MP: 5193
- material, the MP is made of: Germany Stainless (stainless steel, produced in Germany)
- CE certification mark.

3. Manufacturer's carton is marked with the following data:

A shipping package (picture 43) protects the MP against mechanical stress associated with transportation. The MP packed for transportation is resistant to vibration and shock. The shipping package is a 395x395x245 mm (WxDxH) carton marked with the below information (fig. 43, 44, 45, 46).

Fig. 43. Shipping package

Fig. 44. Information on the shipping package.



Fig.45. Marking on the shipping package

Symbols:

Name and address of the manufacturer

Type : TCM 3000 BL Morcellator Name of the MP

REF 3281 Reference number

SN Serial number

Power : 100V~/115V~/230V~ 50-60Hz 110VA Symbol indicating power consumption

Fuse : T1AL 250VAC Fuses used

CE 1275 c Symbol indicating that CE Certificate has been obtained and that the MP is certified by the Canadian Standards Association (CSA)

 2D data matrix (QR-code)

 Indicator of hazardous waste which may not be disposed of together with household waste

 Protection class II

 Information on the country of manufacture


Fig. 46. Informational block on the shipping package.

Informational block on the shipping package includes:

1. Corporate logo with manufacturer's address and phone numbers.
2. The following notice: «The operation manual is not available in printed form with this product. The operation manual as PDF is on the CD-ROM that was included in the delivery. A CD-ROM drive and the Adobe Acrobat software is needed to read or print this operation manual. To order the operation manual in printed form please check on the type plate for the address or the website (www.nouvag.com). If you send us the article number (31946 RUS) with the language code and the note "Printed copy", a printed manual will be sent to you free of charge».

The shipping package is marked with the following handling symbols.

Fig.47 Handling symbols on the shipping package.

Table 10. Handling symbols on the package.

	Keep dry The device must not be located near moisture sources or in a humid facility
	Handle with care. Fragile goods. Goods shall be handled with care. Such a symbol is used to mark goods which may be damaged in the result of tossing or shock.

Transportation conditions

The symbol is used to mark the transportation requirements (humidity and atmospheric pressure).

Transportation conditions

The symbol is used to mark the transportation requirements (temperature).

The MP is packed as follows inside the shipping carton:

Fig.48. Shipping package, inside view.

The MP is packed into a 370x220x130 mm polyamide box, its accessories are packed into a 375x155x130 mm cardboard box. The box also contains a CD-ROM with the Operation Manual and a warranty card and both boxes are packed into the shipping package.

17. TRANSPORTATION AND STORAGE CONDITIONS.

There are no special storage requirements. The MP may be stored at the below conditions. The MP packed in the manufacturer's package should be stored in the warehouse and put on the racks. Avoid putting any other equipment or goods on the MP.

There are no special transportation requirements as well. The MP doesn't require disassembling for transportation. The MP may be shipped by any covered vehicles or in multi-purpose containers as per the relevant shipping regulations at the below ambient conditions. Arrangement and attachment of shipping packages must ensure their stable position preventing packages from moving and bumping into each other.

Transportation and storage:

Relative humidity:	Max. 90 %
Temperature:	0° C – 60°
Atmospheric pressure:	700 gPa – 1060 gPa

18. SERVICE LIFE.

Shelf life of the MP is 5 years. End of service life of products, generally, is determined by their wear and tear, and damages due to their use. Cutting tubes are intended for 20 working cycles (use-sterilization), trocar cannulas are designed for 500 working cycles (use-sterilization). After the service life is over this MP must be disposed of. After exploitation resource is over the MP must be disposed of as per national regulations of the importing country as well as recommendations of the manufacturer.

19. INFORMATION ON DISPOSAL

The MP shall be disposed of and destroyed as per the classification and rules of assembling, operation, decontamination, handling, storage, transportation, recording and disposal enforced by the competent authority of the importing country (Safety class of medical wastes as per Sanitary Regulations SanPin 2.1.7.2790-10 according to Russian Classification).

Do not dispose the MP with household waste!

When disposing of MP the requirements specified in legislation must be followed.

To ensure environmental protection, old MP can be returned to the dealer or manufacturer.

Sterilize instruments before disposal. If they are not sterilized before disposal, please observe the national and local regulations for disposal of infectious waste.

User is liable for disposal of the MP to specially dedicated collection facilities for recycling of medical waste.

Separate collection and treatment of reworked equipment during disposal helps saving natural resources. Besides, use of the latest methods in waste recycling promotes protection of health of service personnel and the environment. If you have any questions on disposal of the MP after expiration of their service life, please contact the supplier.

20. WARRANTY COVERAGE

The manufacturer Nuvag AG provides warranty coverage of the MP. Manufacturer's address is St.Gallerstrasse, 23-25, CH 9403 Goldach, Switzerland. Guaranteed term of use is 1 year from the date of sale. If you return the warranty card for registration within four weeks of the date of purchase, warranty coverage will be extended for a further 6 months. Improper use or repair, or failure to observe these instructions, will relieve the manufacturer from any obligation arising from warranty provisions or other claims.

Please send any claims or questions regarding the circulation of this MP in the Russian Federation to:

MedTestConsulting Ltd., Russia, 123181, Moscow, 1 Kulakova Str., bldg. 1, office 193, phone/fax +7 (495) 750-29-01, +7(495) 729-77-17, E-mail: medtestconsult@mail.ru.

21. CLAIMS



Please send claims to quality of the MP to the manufacturer Nuvag AG, St.Gallerstrasse, 23-25, CH 9403 Goldach, Switzerland, and to its representative in the Russian Federation: MedTestConsulting Ltd., Russia, 123181, Moscow, 1 Kulakova Str., bldg. 1, office 193, phone/fax +7 (495) 750-29-01, +7(495) 729-77-17, E-mail: medtestconsult@mail.ru.

Appendix 1. Parts of Instruments to the TCM 3000 BL – Morcellator for Surgical Procedures:


1. 5154nou Cutting tube, Ø 12 mm
2. 5155nou Cutting tube, Ø 15 mm
3. 5156nou Cutting tube, Ø 20 mm
4. 5141nou Trocar cannula for Ø 12 mm instruments
5. 5142nou Trocar cannula for Ø 15 mm instruments
6. 5151nou Obturator, Ø 12 mm
7. 5152nou Obturator, Ø 15 mm
8. 5153nou Obturator, Ø 20 mm
9. 5137nou Protective tube, Ø 12 mm
10. 5138nou Protective tube, Ø 15 mm
11. 5139nou Protective tube, Ø 20 mm
12. 5193 Myoma drill, Ø 10 mm
13. 5194 Tenaculum gripper, Ø 5 mm
14. 5195 Claw gripper, Ø 10 mm
15. 5196 Tenaculum gripper, Ø 10 mm
16. 5192 Claw gripper, Ø 13.5/15 mm


Appendix 2. Table of components of the MP


Table 1. Dimensions and weight of MP components


No.	Component name	Appearance	Weight (gr), max	Size (mm), max	Material grade
1.	5139nou Protecting tube, Ø20 mm		130	Length 185 Diameter 20	Medical stainless steel 1.4305
2.	5138nou Protecting tube, Ø15 mm		100	Length 185 Diameter 15	Medical stainless steel 1.4305
3.	5137nou Protecting tube, Ø12 mm		75	Length 185 Diameter 12	Medical stainless steel 1.4305

4.	5153nou Obturator, Ø20 mm		83	Length 340 Diameter 20	Polyesterimide PEI MT
5.	5152nou Obturator, Ø15 mm		82	Length 340, Diameter 15	Polyesterimide PEI MT
6.	5151nou Obturator, Ø12 mm		81	Length 340, Diameter 12	Polyesterimide PEI MT
7.	5156nou Cutting tube, Ø20 mm		630	Length 250 Diameter 20 Thickness left 0.4	Medical stainless steel AISI 304
8.	5155nou Cutting tube, Ø15 mm		610	Length 251 Diameter 15 Thickness left 0.4	Medical stainless steel AISI 304

9.	5154nou Cutting tube, Ø12 mm		610	Length 251 Diameter 12 Thickness left 0.4	Medical stainless steel AISI 304
10.	5142nou Trocacannula for Ø15 mm instruments		320	Length 185 Diameter 15	Medical stainless steel 1.4305
11.	5141nou Trocacannula for Ø12 mm instruments		312	Length 185 Diameter 12	Medical stainless steel 1.4305
12.	5193 Myoma drill, Ø10mm		100	Full length 330, Diameter 10, handle diameter 100, Drill working	Medical stainless steel AISI 304

				length 30	
				Right-hand drill	
13.	5195 Claw gripper, Ø10mm		215	Full length 525 Working length 420 Jaw length 35 Diameter 10	Medical stainless steel AISI 304

14.	5196 Tenaculum gripper, Ø10mm		217	Full length 525 Working length 420 Jaw length 35 Diameter 10	Medical stainless steel AISI 304
-----	--	--	-----	---	-------------------------------------

15.	5194 Tenaculum gripper, Ø5mm		198	Full length 525 Working length 420 Jaw length 35 Diameter 5	Medical stainless steel AISI 304
-----	------------------------------------	--	-----	---	----------------------------------

16.	5192 Tenaculum gripper, Ø13.5/15.00 mm		232	Full length 525 Working length 420 Jaw length 35 Diameter large 15 Diameter small 13,5	Medical stainless steel AISI 304
-----	--	--	-----	---	-------------------------------------

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[Signature]

Zwisler
smart lab solutions

Study site complies GLP and GMP
Accredited according to ISO 17025
AKS P 20804-EU / ZIG-P 436 05 01

CERTIFICATE OF ANALYSIS

MORCELLATOR

REVISION:2

STERILITY VALIDATION ACCORD. ISO 17664:2004

Product Codes: Morcellator
Manufacturer: NOUVAG AG

According to AAMI TIR12:2004 /ISO 17665-1:2006/ 11737-2:2009 Morcellators were tested for efficacy of the sponsors sterilization process (31932 41/10) for reuse (Zwisler report 1010.0717). The test specimen were analysed reassembled. For the sterilisation process #1 and #2, two test specimen each were contaminated with spore stripes of *Geobacillus stearothermophilus* (with worst case by placement of 10 biological indicators on each device). For the sterilization process #3, three test specimen were contaminated with a challenge suspension of *Geobacillus stearothermophilus* spores and Soybean-casein digest medium (with worst case by placement of 10 biological indicators on each device). This microbial challenge was let dry for approx. one hour. Thereafter, the test samples were sterilized in a partial cycle (2.5 min, 134°C) of the recommended sterilizing instructions, whereas reference test samples were kept unclean and the bacterial load was determined by bioburden.

Recommended Sterilization Process

- reassembled test specimen were double packed and sealed in foil for sterilization
- sterilization parameters (moisted heat, type B autoclave): two pre-vacuum steps, holding time 5 minutes at 134°C and drying step for ten minutes (validated after: partial cycle, holding time 2.5 minutes at 134°C and drying step for ten minutes)

After the entire procedure the test specimen were tested for sterility by direct immersion of the samples for 14 days in Soybean-casein digest medium. The test was replicated three times concerning sterilization efficiency (product code: Morcellator). In all cases the sterilization process applied resulted in a reduction of more than 10^6 bacteria per device.

The sterilizing procedure applied is
adequate for the tested medical device.

Date: 20.09.2017

Leszek Kazana, Area Sales Manager

[Signature: L. Kazana]

[Signature: Christine Sauter]

Christine Sauter, Quality Assurance Personnel

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