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RPR

**KITS FOR THE DETECTION
OF NON TREPONEMAL SYPHILIS-ASSOCIATED, REAGIN, ANTIBODIES
IN HUMAN SERUM AND/OR PLASMA BY AGGLUTINATION**

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 500

REF 72516

KITS FOR THE QUALITATIVE AND SEMI-
QUANTITATIVE DETECTION OF NON TREPONEMAL
SYPHILIS-ASSOCIATED, REAGIN, ANTIBODIES IN
HUMAN SERUM OR PLASMA BY MACROSCOPIC
AGGLUTINATION ON DISPOSABLE TEST CARDS

IVD **CE**

883684 - 2014/12

BIO-RAD

TABLE OF CONTENT

1.	INTENDED USE.	5
2.	SUMMARY AND EXPLANATION OF THE TEST.	5
3.	PRINCIPLES OF THE PROCEDURE.	5
4.	REAGENTS	6
5.	WARNING AND PRECAUTIONS.	6
6.	SPECIMENS	8
7.	PROCEDURE.	8
8.	TEST LIMITATION	10
9.	PERFORMANCES CHARACTERISTICS	11
10.	BIBLIOGRAPHY REFERENCES..	12

1. INTENDED USE

RPR kits are intended for use for the qualitative and semi-quantitative detection of non treponemal syphilis-associated antibodies (reagin) in human serum and plasma as an aid in the diagnosis in syphilis infection.

2. SUMMARY AND EXPLANATION OF THE TEST

Syphilis is a chronic infection that progresses through distinct stages of infection: primary, secondary, tertiary, and quaternary. These stages produce diverse clinical symptoms, typically producing initial chancres then syphilitic rash followed by long periods of dormancy. Untreated infection may eventually result in cardiovascular problems and neurosyphilis.

The infection is caused by the spirochaete *Treponema pallidum*, and is usually acquired by sexual contact, although the disease may be transmitted by transfusion of infected blood. Intra-uterine infection also occurs. The organism has proved virtually impossible to culture in artificial media, and diagnosis of the infection usually depends on the demonstration of antibodies in the blood, which appear soon after initial infection, and may persist for many years.

Tests for Syphilis fall into four categories: direct microscopic examination, treponemal antibody tests, non- treponemal antibody tests and direct antigen tests.

The RPR (Rapid Plasma Reagin) is a “non-treponemal” test, in that the antibodies detected are not specific for *T pallidum*, although their presence in patient’s serum or plasma is strongly associated with infection by the organism. This type of test measures antibody (IgG and IgM) produced in response to lipoidal material released from damaged host cells as well as to lipoprotein-like material released from the spirochaetes. These antibodies tend to disappear after successful cure of the infection.

3. PRINCIPLES OF THE PROCEDURE

RPR kits use carbon particles coated with a mixture of lipid antigens, which will combine with reagin antibody present in patient’s serum or plasma. The particles are suspended in a medium containing components to eliminate non-specific reactions.

Positive reactions are shown by macroscopic agglutination (clumps) of the particles.

Although the kit is intended for use primarily as a qualitative test, antibody levels may be titrated by doubling dilution.

Agglutination patterns are interpreted by eye.

4. REAGENTS

4.1. Description

Identification on label	Description	Presentation	
		72515 100 tests	72516 500 tests
R1 RPR Antigen	RPR Antigen Carbon particles coated with cardiolipin, lecithin and cholesterol antigens in phosphate buffer	1 vial 2 ml	1 vial 10 ml
R2 Positive control	Positive control Human serum containing antibodies associated to <i>T. pallidum</i> , negative for HBs Antigen, anti-HIV1/2, and anti-HCV antibodies diluted in phosphate buffer	1 vial 1 ml	1 vial 2 ml
R3 Negative control	Negative control Rabbit serum in phosphate buffer	1 vial 1 ml	1 vial 2 ml
	Dispensing bottle (re-usable)	1	1
	Dispensing needle (re-usable)	1	1
	Test cards (10 circles)	10	50

4.2. Storage and handling requirements

This kit should be stored at +2-8°C. Store bottles up-right. Do not freeze.

Each item of the kit preserved at +2-8°C can be used up to the expiry date mentioned on the package. After opening and in the absence of contamination, the R1, R2 & R3 reagents preserved at +2-8°C can be used up to the expiry date shown on the label.

5. WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. For healthcare professional use.

5.1. Health and Safety precautions:

- This test kit should be handled only by qualified personnel trained in laboratory procedures and familiar with their potential hazards. Wear appropriate protective clothing, gloves and eye/face protection and handle appropriately with the requisite Good Laboratory Practices.

- The control materials supplied are derived from human serum. They have been tested at donor level and found negative for HBs Antigen, anti-HIV1/2, and anti-HCV antibodies. No known test method can offer complete assurance that infectious agents are absent. Therefore, all human blood derivatives, reagents and human specimens should be handled as if capable of transmitting infectious disease, following recommended Universal Precautions for blood borne pathogens as defined by local, regional and national regulations.

- Biological spills: Human source material spills should be treated as potentially infectious.

Spills not containing acid should be immediately decontaminated, including the spill area, materials and any contaminated surfaces or equipment, with an appropriate chemical disinfectant that is effective for the potential biohazards relative to the samples involved (commonly a 1:10 dilution of household bleach, 70-80% Ethanol or Isopropanol, an iodophor [such as 0.5% Wescodyne™ Plus, etc.), and wiped dry.

Spills containing acid should be appropriately absorbed (wiped up) or neutralized, the area flushed with water and wiped dry; materials used to absorb the spill may require biohazardous waste disposal. Then the area should be decontaminated with one of the chemical disinfectants.

Note: Do not place solutions containing bleach into the autoclave!

- Dispose of all specimens and material used to perform the test as though they contain an infectious agent. Laboratory, chemical or biohazardous wastes must be handled and discarded in accordance with all local, regional and national regulations.

- The Safety Data Sheet is available on www.bio-rad.com.

5.2. Precautions related to the procedure

5.2.1. Preparing

The reliability of the results depends on correct implementation of the following Good Laboratory Practices:

- Do not mix or associate reagents from different lots within a test run.
- Do not use expired reagents.
- Before use wait for 30 minutes for the reagents to stabilize at room temperature (18-30°C).

5.2.2. Processing

- Do not change the assay procedure.
- Use a new distribution tip for each sample.
- Do not touch the reaction surface of the agglutination cards.

6. SPECIMENS

Serum or plasma (EDTA, Sodium Citrate, Sodium Heparin and ACD) specimens should be free of blood cells. They may be stored at +2-8°C for up to 7 days before testing. Specimens needing longer storage should be frozen at -20°C or lower. Frozen specimens should be thawed and well mixed before testing. Do not repeat more than 5 freeze/thaw cycles.

Heated samples at 56°C during 1 hour do not impact the results. Specimens containing up to 120 g/l of albumin, 200 mg/l of bilirubin, 33 g/l of triolein and samples containing up to 2 g/l of hemoglobin do not affect the results. However, it is not recommended to use hyperlipemic and hyperhemolysed samples.

If the specimens are to be shipped, they must be packaged in accordance with the regulations in force regarding the transport of etiological agents and preferably transport frozen.

7. PROCEDURE

7.1. Materials required

7.1.1. Material provided

Properly calibrated disposable dropping bottle and dropping needle to dispense carbon particles.

7.1.2. Materials required but not provided

Rotator to rotate test cards at 100 rpm, in a circle of diameter approximately 1 cm.

Properly calibrated and maintained pipette and tips to deliver a volume of 50 µl.

0,9% saline solution (for semi-quantitative procedure)

7.2. Assay Procedure

The kit positive (R2) and negative (R3) controls must be run with each run of tests.

7.2.1. Qualitative test

1. Place 50 µl of specimen or control into a circle on the test card.
2. Spread the specimen and controls evenly over the test circle area.
3. Shake the vial of RPR antigen to ensure thorough mixing, **just before use** to avoid sedimentation.

4. Attach the dropping needle to the dropping bottle and take up the RPR antigen by suction.
5. Invert the dropper and gently squeeze to expel air from the needle.
6. Holding the dropping bottle vertically over the test specimen, dispense a single drop of RPR antigen.
7. Place test card on a card rotator and rotate at 100 rpm for 8 minutes.
8. Read **immediately** and interpret results visually in good light. (Refer to 7.4)
Note: If it is not possible to read immediately, card must be maintained on the card rotator up to 15 minutes.
9. Return unused antigen from dropping bottle to glass vial.
10. Clean out dropping bottle and needle with distilled water and allow drying before reusing.

7.2.2. Semi-quantitative test

1. Make doubling dilutions from undiluted to 1:16 in 0,9% saline solution.
2. Place 50 µl of each dilution and controls in to a separate circle on the test card.
3. Spread each dilution evenly over the test circle.
4. Continue as from Qualitative Test section 3.

The titre of the specimen is expressed as the reciprocal of the highest dilution showing agglutination of the carbon particles.

If the highest dilution tested (1:16) is reactive, proceed with a further dilution series by preparing doubling dilutions of the sample from 1:32 to 1:512 using 0,9% saline solution.

Mix well and continue as from step 2 in the semi-quantitative test.

7.3. Quality Control

Use Positive Control (R2) and the Negative Control (R3) in each run to validate the assay.

7.4. Interpretation of the results and Test Validation criteria


Strong Reactive (SR):

Large clumps of carbon particles with a clear background.


Reactive (R):

Large clumps of carbon particles somewhat more dispersed than in Strong Reactive.



Weak Reactive (WR):

Small clumps of carbon particles with light grey background.

Trace Reactive (TR):

Slight clumping of carbon particles typically seen as a button of aggregates in the centre of the test circle or dispersed around the edge of the test circle.

Non-Reactive (NR):

Typically a smooth grey pattern or a button of non-aggregated carbon particles in the centre of the test circle or a large circle of carbon particles without any clumping inside.

Reactive samples should be recorded as antibody positive and must (in view of the nonspecific nature of the antibodies detected) be subjected to further tests to determine the presence or absence of specific anti-treponemal antibodies.

For the assay to be valid, the Positive Control should give a Strong Reactive pattern and the Negative Control should give a clearly Non-Reactive result.

8. TEST LIMITATION

No single test or definitive reference standard is available for every stage of the disease. Thus, Syphilis diagnosis relies predominantly on serological testing, requiring results from both non-treponemal and treponemal methods.

False positive results can be observed when reading does not occur immediately after rotating.

RPR carbon test is not specific for Syphilis. All reactive samples should be retested with treponemal methods such as TPHA to confirm the results.

9. PERFORMANCES CHARACTERISTICS

9.1. Precision Study

A sample panel constituted of one negative sample, 1 low positive sample (titer 1:4) and 1 positive sample (titer 1:16) were tested for repeatability in 8 replicates during the same run. For intermediate precision and inter lot reproducibility studies, the samples were tested in 2 replicates per day during 5 days (reading by two different operators) and on two different lots.

9.1.1. Repeatability

All the negative sample replicates gave negative results and all the positive samples replicates gave positive results.

9.1.2. Intermediate precision / Inter Lot reproducibility

All the negative sample replicates gave negative results and all the positive samples replicates gave positive results whatever the conditions.

9.2. Clinical performance

9.2.1. Specificity

The specificity study was a retrospective study performed on 102 frozen serum samples from the routine of the laboratory of a Sexually Transmitted Disease center, France.

The results from Syphilis RPR assay were compared to a CE marked RPR/VDRL assay.

Total number specimens	Initial Reactive (IR)	Repeat Reactive (RR)	RR Specificity (%)	CI 95%
102	2 (equivocal)	0	100% 102/102	[96.4% - 100.0%]

On 102 samples, 2 samples were found equivocal. After re-test, all samples were negative.

The diagnostic specificity on the retrospective samples was equal to $102/102 = 100\%$ with a confidence interval at 95 % of [96.4% – 100.0%].

9.2.2. Sensitivity

The sensitivity study was a retrospective study performed on 101 frozen serum samples from the routine of the laboratory of a Sexually Transmitted Disease center, France.

The results from Syphilis RPR assay were compared to a CE marked RPR/VDRL assay.

Total number specimens	Initial Reactive (IR)	Repeat Reactive (RR)	RR Specificity (%)	CI 95%
101	100 (1 equivocal)	101	100% 101/101	[96.4% - 100.0%]

On 101 samples, 1 sample was found equivocal. After re-test, all samples were positive.

The diagnostic sensitivity on the retrospective samples was equal to 100% (101/101) with a confidence interval at 95 % of [96.4% - 100.0%].

9.3. Analytical Specificity - Cross Reactivity Study

The analytical specificity of the Syphilis RPR assay was evaluated on a total of 65 cross-reacting samples chosen to be from other infectious diseases and/or other conditions that can give a false-positive result due to non-specificity (Rheumatoid factor, Lyme disease, EBV, Rubella, leptospirosis, SLE (Lupus), Hepatitis B, Hepatitis C, HIV 1/2, multipara women, pregnant women).

The assay specificity obtained is equal to 100.0% (65/65) with a 95% confidence interval of [94.5% - 100.0%].

9.4. Prozone Effect

Three (3) very high Ab Syphilis positive samples (>1:128) were tested undiluted and diluted to verify the absence of prozone effect.

All non-diluted positive samples are positive.

No Prozone effect is seen with titers up to 1:256.

10. BIBLIOGRAPHY REFERENCES

1. Larsen SA., Pettit , et coll., EDTA -treated plasma in the Rapid Plasma Reagin card test and the toluidine red unheated serum test for serodiagnosis of syphilis. J Clin Microbiology 1983;17;431-5
2. Larsen S.A., Pope V., et coll., A manual of Tests for Syphilis. 9th Edition ;1998; 193 – 207
3. Portnoy J. Modifications of the rapid plasma reagin (RPR) card test for syphilis, for use in large-scale testing. Am J Clin Pathol;1963;40;473-9
4. Singh AE and Romanowski. Syphilis: Review with emphasis on clinical, epidemiological, and some biologic features. Clin Microbiol Rev. 1999 Apr; 12(2):187–209.
5. Stability of selected serum proteins after long-term storage in the Janus Serum Bank. Clin Chem Lab Med. 2009. 47:596-606.

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1. Larsen SA., Pettit, et coll., EDTA –treated plasma in the Rapid Plasma Reagin card test and the toluidine red unheated serum test for serodiagnosis of syphilis. J Clin Microbiology 1983;17;431-5
2. Larsen S.A., Pope V., et coll., A manual of Tests for Syphilis. 9th Edition ;1998; 193 – 207
3. Portnoy J. Modifications of the rapid plasma reagin (RPR) card test for syphilis, for use in large-scale testing. Am J Clin Pathol;1963;40;473-9
4. Singh AE and Romanowski. Syphilis: Review with emphasis on clinical, epidemiological, and some biologic features. ClinMicrobiol Rev. 1999 Apr; 12(2):187–209.
5. Stability of selected serum proteins after long-term storage in the Janus Serum Bank. Clin Chem Lab Med. 2009. 47:596-606.

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