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DECLARATION

I, Masashige Kuwayama, CEO of REPROLIFE Inc. (2-5-3 9F, Shinjuku, Shinjuku-ku, 160-0022 Tokyo, Japan), the undersigned, do hereby solemnly and sincerely declare:


That the attached documents:

Instructions for use

are true copies prepared by REPROLIFE Inc.

And I make this solemn declaration conscientiously believing the same to be true and correct.

Dated this 19th day of September, 2017


Masashige Kuwayama
CEO

For Oocytes and Embryos

VITRIFICATION KIT (101) : For 1 Patient(3 times Uses)

Contents:

- 1 vial of 1.0ml Equilibration Solution (ES)
- 2 vials of 1.0ml Vitrification Solution (VS)
- 4 Cryotops:
- 1 Cryotop can contain more than one oocyte/embryo (recommend 3-4 maximum for oocytes and 4-cells embryos, and one for blastocyst for single embryo transfer).
- 3 Vial Plates with 3 wells for each.

Instructions:

Preparation

- The whole process should be performed under the room temperature (25-27°C).
- Fill a nitrogen container.
- Compare the thickness of zona pellucida with the post-thawing space, and take note for oocytes.

Attention: Use a right diameter size of Pasteur pipette for oocytes and embryos (140-160 µm), and blastocysts (160-220 µm).

Equilibration of oocytes and embryos

1. Fill in the 1st well of a Vial Plate with 300 µl of ES, and 300 µl of VS in both of the 2nd and 3rd well.
2. Put an oocyte/embryo on the surface of ES in the 1st well.
3. The oocyte/embryo will sink and begin to shrink, and gradually return to the original size (maximum time is 15 min for oocytes and blastocysts, and 12 min for other stages of embryos).

Vitrification

Attention: The following steps must be taken in more than 25 sec and less than 60 sec.

4. Transfer the oocyte/embryo from the 1st well to the half depth of the 2nd well with VS. (Not with the minimum volume of ES at the first step) The oocyte/embryo immediately floats to the surface of VS while being

washed.

5. After washing the inside wall of the pipette with fresh VS media, take only the oocyte/embryo and transfer it to the bottom of the well. Wait until the oocyte/embryo stops floating in VS.

6. Transfer the oocyte/embryo to the half depth of the 3rd well with VS, and mix the media by pipette around for 5 times.

7. Take the oocyte/embryo only at the tip of the pipette, and put it on the end of the Cryotop seal with a minimum volume of VS.

8. Immediately submerge the Cryotop into liquid nitrogen.

9. Place the cap, and store the Cryotop in a nitrogen tank.

Please leave the Cryotop in liquid nitrogen at all times.

Quality Control Tests:

- Solutions
Lot is labeled on vial.
- Vial Plates
Lot is labeled on pack.
- Cryotops
Lot is labeled on pack.

Successfully passed the following controls.

- Sterility : sterility test.
- Endotoxin by KT methodology (Each component).
- Efficiency : survival of 50/50 mouse embryos and porcine oocytes.

Storage and stability:

Solutions and kits can be transported under the room temperature, and then must be kept in the fridge at 2-8°C until the expiration date.

Use before the expiration date indicated on each of the labels.

Composition:

- Modified HEPES Buffered MEM
- Hydroxy Propyl Cellulose
- Ethylene Glycol
- Dimethyl Sulfoxide
- Endotoxin free Trisulfoxide

References:

- Kamegama M. Highly efficient vitrification for cryopreservation of human oocytes and embryos: The Cryotop method. *Theriogenology* 67, 73-80, 2007.
- Cote A, Kamegama M. Comparison of conventional outcomes achieved with fresh and cryopreserved donor oocytes vitrified by the Cryotop method. *Fertil Steril*, 2009; 92: 1057-64, 2007.
- Ando M, Ueda S, Oda G, Cornejo F, Vizard G, Ando S. Cryotop vitrification of human oocytes results in high survival rate and healthy offspring. *Reproductive Biomedicine Online* 14, 5-057, 2007.
- Velez G, Kamegama M. Improving cryopreservation systems. *Theriogenology* 68(1), 234-44, 2008.
- Kamegama M. Highly efficient vitrification method for cryopreservation of human oocytes. *Reproductive Biomedicine Online* 11:288-294, 2005.
- Ueda S, Kamegama M. High survival rate of bovine oocytes matured in vivo following vitrification. *J Reprod Dev*. 50:395-400, 2004.
- Patel V, Kamegama M. Effect of oocyte type and donor's sexual maturity on vitrification of mice whole oocytes at germinal vesicle stage. *Zygote* 12, 333-335, 2004.
- Hoshi S, Kamegama M. Improved Survival of Vitrified in vivo-derived porcine embryos. *J. Reprod. Develop.* 50, 491-500, 2004.
- Hoshi R, Kamegama M. Cryopreservation of porcine embryos derived from in vitro-matured oocytes. *Biological Reproduction* 71, 498-507, 2004.

Product only for in vitro use.


2010/11/01 Ver. 2.1

Cryotech Japan 2-3-3-8F Shinjuku, Shinjuku-ku, Tokyo 160-0022 JAPAN contact@cryotech-japan.jp <http://cryotech-japan.jp/>

For Oocytes and Embryos

WARMING KIT (102) : For 1 time Use

Contents:

- 1 vial of 1.8 ml Warming Solution (TS)
- 1 vial of 0.5ml Diluent Solution (DS)
- 1 vial of 1ml Washing Solution (WS)
- 1 Warm Plate with 4 wells

Instructions:

Preparation

- The whole process should be performed under the room temperature (25-27°C).

Attention: Use a right diameter size of Pasteur pipette for oocytes and embryos (140-150 µm), and blastocysts (190-220 µm).

- Place a Warm plate and TS vial (with rd) in an incubator at 37°C 3 hours before the use (overnight storage is preferable).

- Expose DS and WS vials to the room temperature air at least 1 hour before the use.

- Take a patient's cane out of a liquid nitrogen tank, and take off a cover cap of Cryotecs and prop up a Cryotec to the inside wall of a cooling rack with liquid nitrogen tearing it against the wall.

Warming and dilution of CPAs

1. Take the Warm Plate out of the incubator, and fill the second well with 300µl of DS.
2. Take the TS vial out of the incubator, and expel TS to the 1st square well.
3. Quickly (within 1 sec) put the Cryotec into the 1st square well with TS, and wait for 1 min.
4. Aspirate the oocyte/embryo and 3mm long of TS after it into the pipette, and expel them most slowly to the bottom of the 2nd well (DS), and wait for 3 min.
5. While waiting, fill the 3rd(WS1) and 4th well(WS2) for each with 300µl of WS.

6. Aspirate the oocyte/embryo and 3mm long of DS after it into the pipette, and expel them most slowly to the bottom of the 3rd well (WS1), and wait for 5 min.

7. Give a survival judgment at the end of this step if the shrunk oocyte/embryo is recovered or not.

8. Put the oocyte/embryo on the surface of the 4th well (WS2). When they sink and reach to the bottom, put them again on the surface of the same WS2 to wash for 2 times in total.

9. Put the oocyte/embryo into the droplet of the culture media for the recovery for ICSI or ET.

Note: 2 to 4 hours of culture for oocytes, and 3 hours for embryos.

Quality Control Tests:

- Solutions
Lot is labeled on Vial.
- Warm Plates
Lot is labeled on pack.

Successfully passed the following controls.

- Sterility : sterility test.
- Endotoxin by KT methodology (Each component).
- Efficiency : survival of 50/50 mouse embryos and porcine oocytes.

Storage and stability:

Solutions and kits can be transported under the room temperature, and then must be kept in the fridge at 2-6°C until the expiration date. Use before the expiration date indicated on each of the labels.

Composition:

- Modified HEPES Buffered MEM
- Hydroxy Propyl Cellulose
- Endotoxin free TrisHCl

References:


- Kawayama M. Highly efficient vitellization for cryopreservation of human oocytes and embryos: The CryoTec method. *Theriogenology* 67, 73-80, 2007.
- Cebal A, Kawayama M. Comparison of conventional oocyte vitellization with fresh and cryopreserved donor oocytes vitellized by the CryoTec method. *Fertil Steril*, 88(4): 1057-64, 2007.
- Aubert M, Luetjens J, Durr G, Camus F, Vimeau C, Arnaud B. Cryotop vitellization of human oocytes results in high survival rate and healthy daughters. *Reproductive Biomedicine Online* 14, 5-607, 2007.
- Vello G, Kawayama M. Improving cryopreservation systems. *Theriogenology* 68(1), 226-44, 2008.
- Kawayama M. Highly efficient vitellization method for cryopreservation of human oocytes. *Reproductive Biomedicine Online* 11:300-303, 2005.
- Uchida J, Kawayama M. High survival rate of bovine oocytes matured in vitro following vitellization. *J Reprod Dev*. 50:288-95, 2004.
- Riedl Y, Kawayama M. Effect of oocyte type and donor's sexual maturity on vitellization of murine whole oocytes at germinal vesicle stage. *Zygote* 12, 333-338, 2004.
- Hoshi S, Kawayama M. Improved Survival of Vitellized in vivo-derived porcine embryos. *J. Reprod. Develop.* 50, 481-485, 2004.
- Basil R, Kawayama M. Cryopreservation of porcine embryos derived from in vitro-matured oocytes. *Biology of Reproduction*. 71, 432-437, 2004.

Product only for in vitro use.

Cryotech Japan 2-6-3-8F Shinjuku, Shinjuku-ku, Tokyo 168-0022 JAPAN contact@cryotech-japan.jp <http://cryotech-japan.jp/>

2010/11/01 Ver. 2.1


Cryotech®

For Oocytes and Embryos

VITRIFICATION SOLUTION SET (110) : For 10 times Uses

Contents:

- 2 vials of 1.8 ml Equilibration Solution (ES)
- 4 vials of 1.8 ml Vitrification Solution (VS)

Instructions:

Preparation

- The whole process should be performed under the room temperature (25-27°C).
- Fill a nitrogen container.
- Compare the thickness of zona pellucida with the perivitelline space, and take note for oocytes.

Attention: Use a right diameter size of Pasteur pipette for oocytes and embryos (140-150 µm), and blastocysts (160-220 µm).

Equilibration of oocytes and embryos

1. Fill in the 1st well of a Vitrification Plate with 300 µl of ES, and 300 µl of VS in both of the 2nd and 3rd well.
2. Put an oocyte/embryo on the surface of ES in the 1st well.
3. The oocyte/embryo will sink and begin to shrink, and gradually return to the original size (maximum time is 15 min for oocytes and blastocysts, and 12 min for other stages of embryos).³

Vitrification

Attention: The following steps must be taken in more than 25 sec and less than 90 sec.

4. Transfer the oocyte/embryo from the 1st well to the half depth of the 2nd well with VS. (Not with the minimum volume of ES at the first step) The oocyte/embryo immediately floats to the surface of VS while being washed.
5. After washing the inside wall of the pipette with fresh VS media, take only the oocyte/embryo and transfer it to the bottom of the well. Wait until the oocyte/embryo stops

floating in VS.

6. Transfer the oocyte/embryo to the half depth of the 3rd well with VS, and mix the media by pipette around for 5 times.
7. Take the oocyte/embryo only at the tip of the pipette, and put it on the end of the Cryotec seat with a minimum volume of VS.
8. Immediately submerge the Cryotec into liquid nitrogen.
9. Place the cap, and store the Cryotec in a nitrogen tank.

Please leave the Cryotec in liquid nitrogen at all times.

Quality Control Tests:

-Solutions

Lot is labeled on vial.

Successfully passed the following controls.

- Sterility : sterility test .
- Endotoxin by KT methodology (Each component).
- Efficiency: survival of 50/50 mouse embryos and porcine oocytes.

Storage and stability:

Solutions and kits can be transported under the room temperature, and then must be kept in the fridge at 2-8°C until the expiration date. Use before the expiration date indicated the each of the labels.

Composition:


- Modified HEPES Buffered MEM
- Hydroxy Propyl Cellulose
- Ethylene Glycol
- Dimethyl Sulfoxide
- Endotoxin free Trehalose

References:

- Kuwayama M. Highly efficient vitrification for cryopreservation of human oocytes and embryos: The CryoTop method. *Theriogenology* 67, 73-80, 2007.
- Cobo A, Kuwayama M. Comparison of concomitant outcome achieved with fresh and cryopreserved donor oocytes vitrified by the Cryotop method. *Fertil Steril*. J89(8): 1657-64, 2007.
- Antinori M, Licata E, Danti G, Cerusio F, Versed C, Antinori S. Cryotop vitrification of human oocytes results in high survival rate and healthy deliveries. *Reproductive Biomedicine Online* 14, 5-667, 2007.
- Vajta G, Kuwayama M. Improving cryopreservation systems. *Theriogenology* 65(1), 236-44, 2006.
- Kuwayama M. Highly efficient vitrification method for cryopreservation of human oocytes. *Reproductive Biomedicine Online* 11:300-308, 2005.
- Ushijima J, Kuwayama M. High survival rate of bovine oocytes matured in vitro following vitrification. *J Reprod Dev*. 50:885-96, 2004.
- Fukui Y, Kuwayama M. Effect of cryodevice type and donor's sexual maturity on vitrification of minke whale oocytes at germinal vesicle stage. *Zygote* 12, 333-338, 2004.
- Hochi S, Kuwayama M. Improved Survival of Vitrified in vivo-derived porcine embryos. *J. Reprod. Develop.* 50, 481-486, 2004.
- Esaki R, Kuwayama M. Cryopreservation of porcine embryos derived from in vitro-matured oocytes. *Biology of Reproduction*. 71, 432-437, 2004.


Product only for in vitro use.

For Oocytes and Embryos

WARMING SOLUTION SET (205) : For 5 times Uses

Contents:

- 5 vials of 1.8 ml Warming Solution (TS)
- 1 vial of 1.8 ml Diluent Solution (DS)
- 2 vials of 1.8 ml Washing Solution (WS)

Instructions:

Preparation

-The whole process should be performed under the room temperature (25-27°C).

Attention: Use a right diameter size of Pasteur pipette for oocytes and embryos (140-150 µm), and blastocysts (160-220 µm).

-Place a Warm plate and TS vial (with rid) in an incubator at 37°C 3 hours before the use (overnight storage is preferable).

-Expose DS and WS vials to the room temperature air at least 1 hour before the use.

-Take a patient's cane out of a liquid nitrogen tank, and take off a cover cap of Cryotec and prop up a Cryotec to the inside of a cooling rack with liquid nitrogen learning it against the wall.

Warming and dilution of CPAs

1. Take the Warm Plate out of the incubator, and fill the second well with 300µl of DS.
2. Take the TS vial out of the incubator, and expel TS to the 1st square well.
3. Quickly (within 1 sec) put the Cryotec into the 1st square well with TS, and wait for 1 min.
4. Aspirate the oocyte/embryo and 3mm long of TS after it into the pipette, and expel them most slowly to the bottom of the 2nd well (DS), and wait for 3 min.
5. While waiting, fill the 3rd(WS1) and 4th wells(WS2) for each with 300µl of WS.

6. Aspirate the oocyte/embryo and 3mm long of DS after it into the pipette, and expel them most slowly to the bottom of the 3rd well (WS1), and wait for 5 min.

7. Give a survival judgment at the end of this step if the shrunk oocyte/embryo is recovered or not.

8. Put the oocyte/embryo on the surface of the 4th well (WS2). When they sink and reach to the bottom, put them again on the surface of the same WS2 to wash for 2 times in total.

9. Put the oocyte/embryo into the droplet of the culture media for the recovery for ICSI or ET.

Note: 2 to 4 hours of culture for oocytes, and 3 hours for embryos.

Quality Control Tests:

-Solutions

Lot is labeled on vial.

Successfully passed the following controls.

- Sterility : sterility test .
- Endotoxin by KT methodology (Each component).
- Efficiency: survival of 50/50 mouse embryos and porcine oocytes.

Storage and stability:

Solutions and kits can be transported under the room temperature, and then must be kept in the fridge at 2-8°C until the expiration date. Use before the expiration date indicated on each of the labels.

Composition:

- Modified HEPES Buffered MEM
- Hydroxy Propyl Cellulose
- Endotoxin free Trehalose

References:


- Kuwayama M. Highly efficient vitrification for cryopreservation of human oocytes and embryos: The CryoTop method. *Theriogenology* 67, 73-80, 2007.
- Cobo A, Kuwayama M. Comparison of concomitant outcome achieved with fresh and cryopreserved donor oocytes vitrified by the Cryotop method. *Fertil Steril*. J89(6): 1657-64, 2007.
- Antinori M, Licata E, Dani G, Cerusico F, Versaci C, Antinori S. Cryotop vitrification of human oocytes results in high survival rate and healthy deliveries. *Reproductive BioMedicine Online* 14, 5-667, 2007.
- Valta G, Kuwayama M. Improving cryopreservation systems. *Theriogenology* 65(1), 236-44, 2006.
- Kuwayama M. Highly efficient vitrification method for cryopreservation of human oocytes. *Reproductive BioMedicine Online* 11:300-308, 2005.
- Ushijima J, Kuwayama M. High survival rate of bovine oocytes matured in vitro following vitrification. *J Reprod Dev*. 50:685-86, 2004.
- Fukui Y, Kuwayama M. Effect of cryodevice type and donor's sexual maturity on vitrification of minke whale oocytes at germinal vesicle stage. *Zygote* 12, 333-338, 2004.
- Hochi S, Kuwayama M. Improved Survival of Vitrified in vitro-derived porcine embryos. *J. Reprod. Develop.* 50, 481-488, 2004.
- Esaki R, Kuwayama M. Cryopreservation of porcine embryos derived from in vitro- matured oocytes. *Biology of Reproduction*. 71, 432-437, 2004.

Product only for in vitro use.

2016/11/01 Ver. 2.1

Cryotech Japan 2-5-3-9F Shinjuku, Shinjuku-ku, Tokyo 160-0022 JAPAN contact@cryotech-japan.jp <http://cryotech-japan.jp/>

Drafts of the Instructions For Use for Russian market

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- Kuwayama M. Highly efficient vitrification for cryopreservation of human oocytes and embryos: The CryoTop method. *Theriogenology* 67, 732007.
- Cobo A, Kuwayama M. Comparison of concomitant outcome achieved with fresh and cryopreserved donor oocytes vitrified by the Cryotop method. *Fertil Steril.* J89(6): 1657-64, 2007.
- Antinori M, Licata E, Dani G, Cerusico F, Versaci C, Antinori S. vitrification of human oocytes results in high survival rate and healthy deliveries. *Reproductive BioMedicine Online* 14, 5
- Vajta G, Kuwayama M. Improving cryopreservation systems. 65(1), 236-44, 2006.
- Kuwayama M. Highly efficient vitrification method for cryopreservation of human oocytes. *Reproductive BioMedicine Online* 11:300
- Ushijima J, Kuwayama M. High survival rate of bovine oocytes matured invitro following vitrification. *J Reprod Dev.* 50:685
- Fukui Y, Kuwayama M. Effect of cryodevice type and donor's sexual maturity on vitrification of minke whale oocytes at germinal vesicle stage. *Zygote* 12,333-338, 2004.
- Hochi S, Kuwayama M. Improved Survival of Vitrified in vivo embryos. *J. Reprod. Develop.* 50, 481-486, 2004.

- Esaki R, Kuwayama M. Cryopreservation of porcine embryos derived from invitro-matured oocytes. *Biology of Reproduction*. 71, 432

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- Kuwayama M. Highly efficient vitrification for cryopreservation of human oocytes and embryos: The CryoTop method. *Theriogenology* 67, 732007.
- Cobo A, Kuwayama M. Comparison of concomitant outcome achieved with fresh and cryopreserved donor oocytes vitrified by the Cryotop method. *Fertil Steril*. J89(6): 1657-64, 2007.
- Antinori M, Licata E, Dani G, Cerusico F, Versaci C, Antinori S. vitrification of human oocytes results in high survival rate and healthy deliveries. *Reproductive BioMedicine Online* 14, 5
- Vajta G, Kuwayama M. Improving cryopreservation systems. 65(1), 236-44, 2006.

- Kuwayama M. Highly efficient vitrification method for cryopreservation of human oocytes. *Reproductive BioMedicine Online* 11:300
- Ushijima J, Kuwayama M. High survival rate of bovine oocytes matured invitro following vitrification. *J Reprod Dev.* 50:685
- Fukui Y, Kuwayama M. Effect of cryodevice type and donor's sexual maturity on vitrification of minke whale oocytes at germinal vesicle stage. *Zygote* 12,333-338, 2004.
- Hochi S, Kuwayama M. Improved Survival of Vitrified in vivo embryos. *J. Reprod. Develop.* 50, 481-486, 2004.
- Esaki R, Kuwayama M. Cryopreservation of porcine embryos derived from invitro- matured oocytes. *Biology of Reproduction*. 71, 432

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- Kuwayama M. Highly efficient vitrification for cryopreservation of human oocytes and embryos: The CryoTop method. Theriogenology 67, 732-737, 2007.
- Cobo A, Kuwayama M. Comparison of concomitant outcome achieved with fresh and cryopreserved donor oocytes vitrified by the Cryotop method. Fertil Steril. 89(6): 1657-64, 2007.
- Antinori M, Licata E, Dani G, Cerusico F, Versaci C, Antinori S. vitrification of human oocytes results in high survival rate and healthy deliveries. Reproductive BioMedicine Online 14, 5
- Vajta G, Kuwayama M. Improving cryopreservation systems. 65(1), 236-44, 2006.
- Kuwayama M. Highly efficient vitrification method for cryopreservation of human oocytes. Reproductive BioMedicine Online 11:300
- Ushijima J, Kuwayama M. High survival rate of bovine oocytes matured invitro following vitrification. J Reprod Dev. 50:685

- Fukui Y, Kuwayama M. Effect of cryodevice type and donor's sexual maturity on vitrification of minke whale oocytes at germinal vesicle stage. Zygote 12,333-338, 2004.
- Hochi S, Kuwayama M. Improved Survival of Vitrified in vivo embryos. J. Reprod. Develop. 50, 481-486, 2004.
- Esaki R, Kuwayama M. Cryopreservation of porcine embryos derived from invitro-matured oocytes. Biology of Reproduction. 71, 432

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- Kuwayama M. Highly efficient vitrification for cryopreservation of human oocytes and embryos: The CryoTop method. *Theriogenology* 67, 732007.
- Cobo A, Kuwayama M. Comparison of concomitant outcome achieved with fresh and cryopreserved donor oocytes vitrified by the Cryotop method. *Fertil Steril*. J89(6): 1657-64, 2007.
- Antinori M, Licata E, Dani G, Cerusico F, Versaci C, Antinori S. vitrification of human oocytes results in high survival rate and healthy deliveries. *Reproductive BioMedicine Online* 14, 5
- Vajta G, Kuwayama M. Improving cryopreservation systems. 65(1), 236-44, 2006.
- Kuwayama M. Highly efficient vitrification method for cryopreservation of human oocytes. *Reproductive BioMedicine Online* 11:300
- Ushijima J, Kuwayama M. High survival rate of bovine oocytes matured invitro following vitrification. *J Reprod Dev*. 50:685

- Fukui Y, Kuwayama M. Effect of cryodevice type and donor's sexual maturity on vitrification of minke whale oocytes at germinal vesicle stage. *Zygote* 12,333-338, 2004.
- Hochi S, Kuwayama M. Improved Survival of Vitrified in vivo embryos. *J. Reprod. Develop.* 50, 481-486, 2004.
- Esaki R, Kuwayama M. Cryopreservation of porcine embryos derived from invitro-matured oocytes. *Biology of Reproduction*. 71, 432

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- Kuwayama M. Highly efficient vitrification for cryopreservation of human oocytes and embryos: The CryoTop method. *Theriogenology* 67, 732-737, 2007.
- Cobo A, Kuwayama M. Comparison of concomitant outcome achieved with fresh and cryopreserved donor oocytes vitrified by the Cryotop method. *Fertil Steril*. 89(6): 1657-64, 2007.
- Antinori M, Licata E, Dani G, Cerusico F, Versaci C, Antinori S. vitrification of human oocytes results in high survival rate and healthy deliveries. *Reproductive BioMedicine Online* 14, 565-569, 2007.
- Vajta G, Kuwayama M. Improving cryopreservation systems. *Reproductive BioMedicine Online* 11:300-304, 2006.
- Kuwayama M. Highly efficient vitrification method for cryopreservation of human oocytes. *Reproductive BioMedicine Online* 11:300-304, 2006.
- Ushijima J, Kuwayama M. High survival rate of bovine oocytes matured in vitro following vitrification. *J Reprod Dev*. 50:685-688, 2004.
- Fukui Y, Kuwayama M. Effect of cryodevice type and donor's sexual maturity on vitrification of minke whale oocytes at germinal vesicle stage. *Zygote* 12, 333-338, 2004.
- Hochi S, Kuwayama M. Improved Survival of Vitrified in vivo embryos. *J. Reprod. Develop.* 50, 481-486, 2004.
- Esaki R, Kuwayama M. Cryopreservation of porcine embryos derived from invitro-matured oocytes. *Biology of Reproduction*. 71, 432-437, 2004.

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- Kuwayama M. Highly efficient vitrification for cryopreservation of human oocytes and embryos: The CryoTop method. *Theriogenology* 67, 732-737, 2007.
- Cobo A, Kuwayama M. Comparison of concomitant outcome achieved with fresh and cryopreserved donor oocytes vitrified by the Cryotop method. *Fertil Steril*. J89(6): 1657-64, 2007.
- Antinori M, Licata E, Dani G, Cerusico F, Versaci C, Antinori S. Vitrification of human oocytes results in high survival rate and healthy deliveries. *Reproductive BioMedicine Online* 14, 5, 2007.
- Vajta G, Kuwayama M. Improving cryopreservation systems. *Reproductive BioMedicine Online* 11:300-304, 2006.
- Kuwayama M. Highly efficient vitrification method for cryopreservation of human oocytes. *Reproductive BioMedicine Online* 11:300-304, 2006.
- Ushijima J, Kuwayama M. High survival rate of bovine oocytes matured in vitro following vitrification. *J Reprod Dev*. 50:685-688, 2004.
- Fukui Y, Kuwayama M. Effect of cryodevice type and donor's sexual maturity on vitrification of minke whale oocytes at germinal vesicle stage. *Zygote* 12, 333-338, 2004.
- Hochi S, Kuwayama M. Improved Survival of Vitrified in vivo embryos. *J. Reprod. Develop.* 50, 481-486, 2004.
- Esaki R, Kuwayama M. Cryopreservation of porcine embryos derived from in vitro-matured oocytes. *Biology of Reproduction*. 71, 432-437, 2004.

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- Kuwayama M. Highly efficient vitrification for cryopreservation of human oocytes and embryos: The CryoTop method. *Theriogenology* 67, 732-737, 2007.
- Cobo A, Kuwayama M. Comparison of concomitant outcome achieved with fresh and cryopreserved donor oocytes vitrified by the Cryotop method. *Fertil Steril*. J89(6): 1657-64, 2007.
- Antinori M, Licata E, Dani G, Cerusico F, Versaci C, Antinori S. vitrification of human oocytes results in high survival rate and healthy deliveries. *Reproductive BioMedicine Online* 14, 565-569, 2007.
- Vajta G, Kuwayama M. Improving cryopreservation systems. *Reproductive BioMedicine Online* 11:300-304, 2006.
- Kuwayama M. Highly efficient vitrification method for cryopreservation of human oocytes. *Reproductive BioMedicine Online* 11:300-304, 2006.
- Ushijima J, Kuwayama M. High survival rate of bovine oocytes matured in vitro following vitrification. *J Reprod Dev*. 50:685-688, 2004.
- Fukui Y, Kuwayama M. Effect of cryodevice type and donor's sexual maturity on vitrification of minke whale oocytes at germinal vesicle stage. *Zygote* 12, 333-338, 2004.
- Hochi S, Kuwayama M. Improved Survival of Vitrified in vivo embryos. *J. Reprod. Develop.* 50, 481-486, 2004.
- Esaki R, Kuwayama M. Cryopreservation of porcine embryos derived from in vitro-matured oocytes. *Biology of Reproduction*. 71, 432-437, 2004.

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- Kuwayama M. Highly efficient vitrification for cryopreservation of human oocytes and embryos: The CryoTop method. *Theriogenology* 67, 732-737, 2007.
- Cobo A, Kuwayama M. Comparison of concomitant outcome achieved with fresh and cryopreserved donor oocytes vitrified by the Cryotop method. *Fertil Steril*. J89(6): 1657-64, 2007.
- Antinori M, Licata E, Dani G, Cerusico F, Versaci C, Antinori S. vitrification of human oocytes results in high survival rate and healthy deliveries. *Reproductive BioMedicine Online* 14, 5
- Vajta G, Kuwayama M. Improving cryopreservation systems. 65(1), 236-44, 2006.
- Kuwayama M. Highly efficient vitrification method for cryopreservation of human oocytes. *Reproductive BioMedicine Online* 11:300
- Ushijima J, Kuwayama M. High survival rate of bovine oocytes matured invitro following vitrification. *J Reprod Dev.* 50:685
- Fukui Y, Kuwayama M. Effect of cryodevice type and donor's sexual maturity on vitrification of minke whale oocytes at germinal vesicle stage. *Zygote* 12, 333-338, 2004.
- Hochi S, Kuwayama M. Improved Survival of Vitrified in vivo embryos. *J. Reprod. Develop.* 50, 481-486, 2004.
- Esaki R, Kuwayama M. Cryopreservation of porcine embryos derived from invitro-matured oocytes. *Biology of Reproduction*. 71, 432

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