STATE OF NORTH CAROLINA

Department of The
Secretary of State

APOSTILLE

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1. COUNTRY: UNITED STATES OF AMERICA
2. This Public Document
has been signed by TIFFANNY A THOMAS
3. acting in the capacity of NOTARY PUBLIC
COMMISSION EXPIRES 9/5/2022
4. bears the seal/stamp of COUNTY OF GUILFORD, NORTH CAROLINA

CERTIFIED

5. at Raleigh, North Carolina
6. the 1ST DAY OF NOVEMBER, 2017
7. by Secretary of State or Deputy Secretary of State, State of North Carolina
8. No. 17711
9. Seal/Stamp
10. Signature

Scan to verify online.

by:

Hailey Haynes
Hailey Haynes
Deputy Secretary of State

Claine J. Marshall

Certification# 5706249623

Verify the issuance of this Apostille at <http://www.sosnc.gov/apostille>

Secretary of State

content of the document for which it was issued. This Apostille is not valid for use anywhere within the United States of America, its territories or possessions. [Verify the issuance of this Apostille at <http://www.sosnc.gov/apostille>]

INSTRUCTION FOR USE

Enteral Tube Feeding Sets – Sterile, with accessories

USAGE Used for enteral feeding and/or administering drugs and also aspiration of the contents of the gastrointestinal tract by minimally invasive procedures gastrostomy of the gastroenteroanastomosis.

DEVICE DESCRIPTION

Parameter	Conversion
1Fr/3	1 mm
1 inch	25,4 mm

Percutaneous Endoscopic Gastrostomy Set

(PEG-24-PUSH-I-S, PEG-24-PULL-I-S, FLOW-20-PUSH-I-S, FLOW-20-PULL-I-S)


Figure 2. FLOW-20-PULL-I-S & PEG-24-PULL-I-S

These devices are designed for use in endoscopy applications. They are invasive, supplied sterile and disposable. Percutaneous Endoscopic Gastrostomy (PEG) is the modality of choice for long-term enteral access and was first described in 1980 by Ponsky and Gauderer. PEG tubes

are an established supportive treatment for a variety of medical conditions. The PEG system is intended for the delivery of enteral nutrition to patients requiring nutritional support. In patients who are unable to maintain sufficient oral intake and have a functional G.I. tract, PEG tubes provide long-term enteral access. Percutaneous endoscopic gastrostomy (PEG) devices - initial placement tubes are placed percutaneously in conjunction with upper gastrointestinal endoscopy. These devices are available in Pull or Push methods. In the pull technique, the PEG tube is pulled through the oropharynx with the aid of an attached loop that pulls it through the stomach and abdominal wall. Using the push technique, the tube is percutaneously inserted directly into the stomach over a wire guide.

The set includes feeding tube, wire guide, bolster kit (the composition described below). In addition preparation tray (the composition described below), bolster kit includes bolster, tubing clamp.

Set includes:

- Bolster: diameter 20 Fr/ 6.67 mm ($\pm 1,5/0,5$), length 15.85mm or 24 Fr/ 8 mm ($\pm 1,5/0,5$), length 15.88mm.
- Twist lock: length 1.2 inches/ 30.48 mm ($\pm 0,02/0,5$); diameter 0.450 inch/11.43mm;
- Cable tie: 4.0 inches/ 101.6 mm +/- 0.250 inches/6.35 mm ($\pm 0,250$ дюйма/6.35 мм),
- Scissors: 5 ½ inches/ 139.7 mm ($\pm 0,02/0,5$),
- Universal feeding adapter (UFA-24-S; UFA-20-S) -please see further below in this TDS
- Bolus feeding adapter (BFA-20-S) -please see further below in this TDS
- Adapter for G/L tube: body 1.240 inch/ 31.5 mm +/- 0.063 inches/1.6 mm in length (+/- 0.063 дюйма/ 1.6 мм),
- Haemostat: 5 inch/ 127 mm ($\pm 0,02/0,5$)
- Tubing clamp: fits 0.268 inch Outer Diameter Tubing (6.8 mm)

Preparation tray is included which has needle cannula, scalpel, gauze pads, drape and cold snare.

- Scalpel: length of handle and length of cutting surface: Scalpel has a 7.5 cm ($\pm 0,5$ мм) length and its sharp part is 1 cm. Roughness of the sharp surface of the cutting edges, blades for scalpels: from 0.05 to 0.63 mkm.
- Gauze pads: 101.6 mm x 101.6 mm ($\pm 0,5$ мм),
- Drape: 24 inch/ 609.6 mm x 24 inch/ 609.6 mm ($\pm 0,02/0,5$),
- Cold snare: Snare size: 1.5cm x 4.3cm; length 150cm
- Needle cannula:

PEG/FLOW - PUSH

Needle – diameter 19 G / 1,07 mm (±0,03mm), length 2.86 inch/ 72.6 mm(± 0,02/0,5)

Cannula - diameter 18 G / 1,25 (±0,05mm), length 2.75 inch/ 69.9 mm(± 0,02/0,5)

PEG/FLOW - PULL

Needle – diameter 14 G / 2,05mm (±0,1 mm), length 2.86 inch/ 72.6 mm(± 0,02/0,5)

Cannula - diameter 14 G / 2,05mm (±0,1 mm), length 3.125 inch/ 79.4 mm(± 0,02/0,5)

Technical Specifications

Product Code/Model	Tube		Wire Guide	
	Diameter Fr/mm	Length, cm	Diameter Inch/ mm	Length cm
Tolerance	± 0.9/0.3mm	+2.54/- 5.08cm	± 0.9/0.3 mm	± 2cm
FLOW-20-PULL-I-S	20/6.67	76.2	0.035/0.89	260
FLOW-20-PUSH-I-S	20/6.67	76.2	0.035/0.89	260
PEG-24-PULL-I-S	24/8	76.2	0.035/0.89	260
PEG-24-PUSH-I-S	24/8	76.2	0.035/0.89	260

**Bolus Feeding Adapter (BFA-20-S)
Universal Feeding Adapter UFA-20-S, UFA-24-S)**

Adapters are non-invasive, supplied non-sterile and disposable. A universal and bolus adapter is included with the PEG feeding tubes to allow physician to choose adapter of choice for delivering nutrition. These adapters are used with PEG/FLOW-PULL/PUSH-I-S devices.


Figure 3. BFA – Bolster Adapter

Figure 4.UFA–Bolster Adapter

Technical Specifications

Product Code/Model	Diameter, Fr/ mm	Length, Inch / mm	Adapter
Tolerance	± 0.9/0.3mm	± 0.063 /1.6mm	—
BFA-20-S	20 / 6.67	1.240 /31.5	Bolus Adapter
UFA-20-S	20/ 6.67	1.240 / 31.5	Universal Adapter
UFA-24-S	24/ 8	1.240 / 31.5	Universal Adapter

INTENDED USE/INDICATION FOR USE

Percutaneous Endoscopic Gastrostomy Sets (PEG-24-PUSH-I-S, PEG-24-PULL-I-S, FLOW-20-PUSH-I-S, FLOW-20-PULL-I-S) are used for percutaneous endoscopic gastrostomy placement to provide enteral nutrition to patients requiring nutritional support.

CONTRAINDICATIONS

Percutaneous Endoscopic Gastrostomy Sets (PEG-24-PUSH-I-S, PEG-24-PULL-I-S, FLOW-20-PUSH-I-S, FLOW-20-PULL-I-S):

- sepsis,
- severe gastroesophageal reflux,
- large ascites or diffuse inflammatory,
- infection,
- neoplastic disease involving the walls of the abdomen or anterior stomach,
- gastrointestinal obstruction
- proximal small bowel fistulae
- disorders of coagulation
- potential infiltration in the abdominal wall or stomach wall
- lack of stomach function.

POTENTIAL ADVERSE EVENTS

Percutaneous Endoscopic Gastrostomy Sets (PEG-24-PUSH-I-S, PEG-24-PULL-I-S, FLOW-20-PUSH-I-S, FLOW-20-PULL-I-S).

- The failure of gastrostomy, depressurization (for gastrostomy)
- Erosive - ulcerative changes in the mucosa of the stomach (gastrostomy for) / duodenum (for jejunostomy)
- Local bleeding
- Peritonitis
- Bumper - syndrome (gastrostomy tube migration into the wall of the stomach , followed by epithelialization gastrostomy hole)
- Infection
- Perforation
- Damage to the esophagus
- Aspiration / pneumonia

- Hyperemia
- Puffiness
- Pneumomediastinum
- Pneumoperitoneum
- Pain

Additional complications associated with endoscopic procedures include:

- Aspiration
- Respiratory distress or respiratory arrest
- Cardiac arrhythmia and cardiac arrest

POTENTIAL COMPLICATIONS

Percutaneous Endoscopic Gastrostomy Sets (PEG-24-PUSH-I-S, FLOW-20-PUSH-I-S, PEG-24-PULL-I-S, FLOW-20-PULL-I-S):

- bronchopulmonary aspiration and pneumonia,
- respiratory distress or airway obstruction,
- peritonitis or septic shock,
- colocutaneous, gastrocolocutaneous or small bowel fistula,
- gastric dilatation,
- sigmoid intra-abdominal herniation and volvulus,
- persistent fistula following tube removal, .
- esophageal injury,
- necrotizing fasciitis,
- candida cellulitis,
- improper placement or inability to place tube,
- tube dislodgment or migration,
- hemorrhage, and tumor metastasis.
- pneumoperitoneum,
- peristomal wound infection and purulent drainage,
- stomal leakage,
- bowel obstruction,
- gastroesophageal reflux (GERD)
- blockage or deterioration of tube.
- death

Note: Patients with cirrhosis have an increased risk of developing ascites which is a contraindication to PEG placement.

INSTRUCTIONS FOR USE

Percutaneous Endoscopic Gastrostomy Set (PEG-24-PUSH-I-S, FLOW-20-PUSH-I-S)

Tube Placement

1. Upon removing the device and its components from the package, visually inspect with particular attention to kinks, bends or breaks in the feeding tube assembly. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.
2. After introducing the gastroscope, insufflate the stomach and examine the mucosa. Determine if the mucosa is free of ulcerations or bleeding before proceeding.
3. Position the tip of the gastroscope and illuminate the left anterior wall of the stomach. Note: It may be necessary to dim the lights in the procedure room to better visualize the light through the abdominal wall. Important: Ideally, placement of the PEG should be low in the body of the stomach or high in the antrum towards the pylorus to allow for easier insertion of the PEG-J (jejunal) feeding tube into the PEG, if indicated as an option.
4. Manipulate the gastroscope tip until the light trans-illuminates the desired

PEG position. (See fig. 1)

5. Lightly depress the illuminated area with a finger while viewing the site with the gastroscope. The depression of this area should be clearly visualized with the scope. After determining that the mucosa is healthy, proceed with this procedure.
6. Prepare the site following surgical guidelines as determined by your institution. Drape the area using the enclosed surgical drape. Inject local anesthetic into the PEG site.
7. Using the enclosed scalpel, make a 1 cm long incision through the skin, subcutaneous tissue. Caution: A smaller incision may contribute to extreme resistance of the gastrostomy feeding tube when exiting the fascia.
8. Endoscopically observe the site.
9. While maintaining stomach insufflation, insert the needle and cannula unit through the skin incision and into the stomach. Leave the cannula in place to maintain access to the stomach while removing the inner needle.
10. Place the floppy tip of the wire guide through the needle cannula and into the stomach.
11. Maintain stomach insufflation to obtain close proximity of the stomach and abdominal walls. Place a snare or non-spiked biopsy forceps through the channel of the gastroscope and grasp the end of the wire guide. Caution: Do not tighten the snare around the needle cannula after removal of the inner stylet as this may interfere with the passage of the wire guide. Note: For wire guide removal grasp the wire with a snare or non-spiked biopsy forceps at least 5 cm from the tip of the wire guide. (See fig. 2) Pull the wire guide to the tip of the endoscope, not into the endoscope channel. (See fig. 3)

12. While maintaining the snare or non-spiked biopsy forceps securely around the wire guide, remove the gastroscope and wire from the patient's mouth. The wire guide will now be protruding from both the patient's mouth and incision site.
13. Using water-soluble lubricant (supplied separately) and gauze, thoroughly lubricate the dilator and entire external length of the tube including the internal bumper.
14. Advance the dilator tip over the wire guide and through the patient's mouth. (See fig. 4) Note: Firm tension on both ends of the wire guide will ease passage through the oropharynx.
15. When the end of the dilator meets the cannula in the stomach, push the dilator through the abdominal wall. Note: Maintaining wire guide tension will aid this process.
16. After the tapered end of the dilator passes through the abdominal wall, release the wire guide and gently pull the tube and wire guide through the abdominal wall.
17. When the internal bumper of the PEG tube enters the mouth, reintroduce the gastroscope and view the tip as it advances through the esophagus and into the stomach. Monitor the patient for respiratory distress as you advance the internal bumper through the esophagus. (See fig. 5)
18. While observing centimeter increments, slowly pull the dilating portion of the tube through the abdominal incision.
19. Bring the internal bumper in contact with the abdominal wall, carefully avoiding excess tension.
20. Apply gentle pressure to the exiting portion of the feeding tube and remove the wire guide. Caution: Blanching of the site indicates excessive pressure on the mucosa and should be avoided.
21. Slide the bolster onto the tube and toward the skin surface. (See fig. 6) The hemostats provided with the kit may be used to facilitate placement of the bolster over the PEG tube. Note: Antiseptic ointment (supplied separately) may be applied to the surrounding tissue prior to sliding the bolster into position. Warning: The bolster should sit close to the skin but not tight against the skin. Excessive traction on the tube may cause premature removal, fatigue or failure of the device.
22. Cut the tube at the X mark.
23. Secure the twist lock or cable tie around the bolster collar, being careful not to crimp it. Important: Use the twist lock or cable tie to secure the bolster to the tube. This will help prevent future migration of the tube and reduce the need to constantly reposition or pull on the tube.
24. Cut off the excess length of the cable tie (if applicable) using the scissors provided.
25. If desired, slide the tubing clamp onto the tube leaving a gap between bolster and tubing clamp.
26. Plug in the adapter of choice (universal or bolus) and close the caps. (See fig. 7) Note: The adapter may be secured to the tube with a cable tie.

27. Note the centimeter marking on the tube that is closest to the bolster and record it on the patient's chart and on the patient information sheet in the Patient Care Manual.

Note: The Patient Care Manual enclosed in the kit is intended as a reference for caregivers of the patient. It is essential for the Patient Care Manual to accompany the patient and be explained to all people responsible for the care of the patient.

28. Dispose of residual kit materials per institutional guidelines for biohazardous medical waste.

29. Record physician's instructions for feeding and administration of prescribed medications in the Patient Care Manual.

30. The patient should remain NPO for 24 hours unless otherwise directed by a physician.

TUBE REMOVAL


The PEG tube has been designed for removal using the external/traction method shown below. If this method of removal is not possible, another method such as endoscopic or surgical should be utilized.

Caution: If a replacement tube is desired, it must be placed immediately after removal. External tube removal may result in minimal trauma or bleeding that may require treatment.

External/Traction Method

The feeding tube may be removed without the need for endoscopy in the following manner:

1. Grasp the gastric feeding tube near the stoma site. While slowly rotating the tube, gently push 1-2 cm of the tube into the stomach to separate the tube from the stoma tract. Warning: If the tube does not rotate freely within the tract, do not attempt to use traction as a method of removal.
2. Hold the gastric tube near the stoma site and apply counter pressure by placing the fingers of the other hand around the base of the tube.
3. Loosely cover the stoma site with either a towel or a drape.
4. Holding the gastric feeding tube straight, apply steady traction to the tube until the internal dome emerges through the abdominal wall. Caution: The tube must be pulled straight out of the stoma tract.
5. Dispose of the feeding tube per institutional guidelines for biohazardous medical waste. The stoma tract should heal and close within 24 hours.

Percutaneous Endoscopic Gastrostomy Set (PEG-24-PULL-I-S, FLOW-20-PULL-I-S)

Tube Placement

1. Upon removing the device and its components from the package, visually inspect with particular attention to kinks, bends or breaks in the feeding tube assembly. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.
2. After introducing the gastroscope, insufflate the stomach and examine the mucosa. Determine if the mucosa is free of ulcerations or bleeding before proceeding.
3. Position the tip of the gastroscope and illuminate the left anterior wall of the stomach. Note: It may be necessary to dim the lights in the procedure room to better visualize the light through the abdominal wall. Important: Ideally, placement of the PEG should be low in the body of the stomach or high in the antrum towards the pylorus to allow for easier insertion of the PEG-J (jejunal) feeding tube into the PEG, if indicated as an option.
4. Manipulate the gastroscope tip until the light trans-illuminates the desired PEG position. (See fig. 1)

5. Lightly depress the illuminated area with a finger while viewing the site with the gastroscope. The depression of this area should be clearly visualized with the scope. After determining that the mucosa is healthy, proceed with this procedure.
6. Prepare the site following surgical guidelines as determined by your institution. Drape the area using the enclosed surgical drape. Inject local anesthetic into the PEG site.
7. Using the enclosed scalpel, make a 1 cm long incision through the skin, subcutaneous tissue. Caution: A smaller incision may contribute to extreme resistance of the gastrostomy feeding tube when exiting the fascia.
8. Endoscopically observe the site.
9. While maintaining stomach insufflation, insert the needle and cannula unit through the skin incision and into the stomach. Leave the cannula in place to maintain access to the stomach while removing the inner needle.
10. Place the looped insertion wire through the needle cannula and into the stomach. (See fig. 2)
11. Maintain stomach insufflation to obtain close proximity of the stomach and abdominal walls. Place a snare or non-spiked biopsy forceps through the channel of the gastroscope and grasp the looped end of the wire. Caution: Do not tighten the snare around the needle cannula after removal of the inner stylet as this may interfere with passage of wire loop.
12. While maintaining the snare or non-spiked biopsy forceps securely around the looped insertion wire, remove the gastroscope and the wire from the patient's mouth. (See fig. 3) The insertion wire will now be protruding from both the patient's mouth and incision site.
13. Pass the insertion wire loop through the looped wire on the dilator end of the feeding tube. (See fig. 4)
14. Place the internal bumper of the feeding tube through the extended insertion wire loop. (See fig. 4)
15. Pull the feeding tube gently through the insertion wire loop and form a knotless connection by applying simultaneous gentle traction on both wire loops. (See fig. 4)
16. Using water-soluble lubricant (supplied separately) and gauze, thoroughly lubricate the dilator and entire external length of the tube including the internal bumper.
17. Advance the dilator tip through the patient's mouth by pulling on the wire exiting the abdominal incision. Depression of the patient's tongue may facilitate the initial introduction of the feeding tube. Gentle pressure on either side of the incision may prevent excessive pulling on the incision site.
18. Continue to pull the insertion wire until the dilator protrudes through the abdominal wall. Note: The tube is extra-long to permit control of the insertion set at all times. Caution: Do not push the tube from the oral portion as this may result in the tube curling in the patient's stomach.

19. When the internal bumper of the PEG tube enters the mouth, reintroduce the gastroscope and view the tip as it advances through the esophagus and into the stomach. (See fig. 5) Monitor the patient for respiratory distress as you advance the internal bumper through the esophagus.
20. While observing centimeter increments, slowly pull the introduction dilator and tube through the abdominal incision. Bring the internal bumper in contact with the stomach wall, carefully avoiding excess tension.
21. Apply gentle pressure to the exiting portion of the feeding tube. Caution: Blanching of the site indicates excessive pressure on the mucosa and should be avoided.
22. Slide the bolster over the dilation catheter loop and onto the tube, past the X mark. The hemostats provided with the kit may be used to facilitate placement of the bolster over the PEG tube. Cut the tube at the X mark. Note: Antiseptic ointment (supplied separately) may be applied to the surrounding tissue prior to sliding the bolster into position. Warning: The bolster should sit close to the skin but not tight against the skin. Excessive traction on the tube may cause premature removal, fatigue, or failure of the device. (See fig. 6)
23. Secure the twist lock or cable tie around the bolster collar, being careful not to crimp it. Important: Use the twist lock or cable tie to secure the bolster to the tube. This will help to prevent future migration of the tube and reduce the need to constantly reposition or pull on the tube.
24. Cut off the excess length of the cable tie (if applicable) using the scissors provided.
25. If desired, slide the tubing clamp onto the tube leaving a gap between bolster and tubing clamp
26. Plug in the adapter of choice (universal or bolus) and close the caps. (See fig. 7) Note: The adapter may be secured to the tube with a cable tie.
27. Note the centimeter marking on the tube that is closest to the bolster and record it on the patient's chart and on the patient information sheet in the Patient Care Manual. Note: The Patient Care Manual enclosed in the kit is intended as a reference for the caregivers of the patient. It is essential for the Patient Care Manual to accompany the patient and be explained to all people responsible for the care of the patient.
28. Dispose of residual kit materials per institutional guidelines for biohazardous medical waste.
29. Record physician's instructions for feeding and administration of prescribed medications in the Patient Care Manual.
30. The patient should remain NPO for 24 hours unless otherwise directed by a physician.

Tube Removal


The PEG tube has been designed for removal using the external/traction method shown below. If this method of removal is not possible, another method such as endoscopic or surgical should be utilized.


Caution: If a replacement tube is desired, it must be placed immediately after removal. External tube removal may result in minimal trauma or bleeding that may require treatment.

External/Traction Method

The feeding tube may be removed without the need for endoscopy in the following manner:

1. Grasp the gastric feeding tube near the stoma site. While slowly rotating the tube, gently push 1-2 cm of the tube into the stomach to separate the tube from the stoma tract. Warning: If the tube does not rotate freely within the tract, do not attempt to use traction as a method of removal.
2. Hold the gastric tube near the stoma site and apply counter pressure by placing the fingers of the other hand around the base of the tube.
3. Loosely cover the stoma site with either a towel or a drape.
4. Holding the gastric feeding tube straight, apply steady traction to the tube until the internal dome emerges through the abdominal wall. Caution: The tube must be pulled straight out of the stoma tract.
5. Dispose of the feeding tube per institutional guidelines for biohazardous medical waste. The stoma tract should heal and close within 24 hours.

Universal Feeding Adapters (UFA-20-S, UFA-24-S) and Bolus Feeding Adapter (BFA-20-S) - these are accessories and no IFU applies to these products.

MATERIALS

Percutaneous Endoscopic Gastrostomy Set (PEG-24-PUSH-I-S, PEG-24-PULL-I-S, FLOW-20-PUSH-I-S, FLOW-20-PULL-I-S)

Component	Material	Patient Contacting
Feeding Tube	Silicone Dow Corning 94-599; Black Ink – Marabu TampaStar TPR	Yes
Wire guide	AISI 304 Stainless Steel PTFE MS#171	Yes
Scalpel	AISI 304 Stainless Steel	Yes
Obturator	AISI 304 Stainless Steel	Yes
Cold Snare	AISI 304 Stainless Steel PTFE MS#171 Polyethylene # 276-73	Yes

Haemostat	AISI 304 Stainless Steel	No
Gauze	Cotton	Yes
Snare	AISI 304 Stainless Steel PTFE MS#171	Yes
Scissors	AISI 304 Stainless Steel	No

Bolus Feeding Adapter (BFA-20-S, BFA-24-S) and Universal Feeding Adapter (UFA-20-S, UFA-24-S)

Component	Material	Patient Contacting
Adapters	PVC	No

BIOCOMPATIBILITY TESTING

Product Name/Description	Biocompatibility Test
PEG/FLOW-PULL/PUSH-I - Feeding Tubes w/ ink markings, Bolsters	Cytotoxicity Study; Sensitization Study; Intracutaneous Study, and/or Biocompatibility Testing and Risk Assessment; and/or Chemical Characterization
PEG/FLOW-PULL/PUSH-I, BFA, UFA - Introduction Components and Adapters	Cytotoxicity Study; Sensitization Study; Intracutaneous Study; and/or Biocompatibility Testing and Risk Assessment

Conclusion: Results of the biological testing and our product history pertaining to the products have shown no indications of a toxic reaction. We therefore anticipate no adverse reactions against any of the components comprising these products. Thus we consider the recommendation as laid out in EN ISO 10993-1 as met.

PERFORMANCE CHARACTERISTICS

Sets for enteral feeding sterile, with accessories sustainable:

- vibratory loads in the range of 10-55 Hz,
- at acceleration 100 m/S2 (10g) and duration of impact acceleration 16 MS,
- at temperatures from +32 to +42°C,

Sets for enteral feeding sterile, with accessories in the transport packaging sustainable:

- vibratory loads in the range 10-55 Hz, displacement amplitude 0,35 mm shock peak acceleration

100 m/S² (10g) and duration of impact acceleration 16 MS,
-at a temperature of transportation and storage -50 to +50°C
Sets for enteral feeding sterile, with accessories are not susceptible to sudden changes in temperature.

Gauze pads soak - in not less than 5 sec

Drape

Surface density 42 g/m² (+- 10%)

Water resistance Not less than 10 cm H₂O

Tensile strength in crosswise direction in dry condition Not less than 40 kPa

Tensile strength in lengthwise direction in dry condition Not less than 40 kPa

Tensile strength in crosswise direction in wet condition Not less than 20 kPa

Tensile strength in lengthwise direction in wet condition Not less than 20 kPa

Catheter measuring device

Length of scale line 4.4 cm

Graduation

Marking from 1,2 to 4,4 cm:

from 1,2 to 1,7 cm – 0,5 cm

from 1,7 to 2,4 cm – 0,7 cm

from 2,4 to 4,4 cm – 1,0 cm

Hardness/firmness of the component made from stainless steel 40-42 HRC

The strength of the connection tubing connectors – not less than 40 N

Color of adapter:

input – red

output – blue

Syringe 60 ml

The minimum length of the plunger protruding from the cylinder -- 5 mm

Diameter of the tip -- not less than 5 mm

Dead space – not more than 0.07ml

Parameters of the scale with tolerances:

Less than half the nominal capacity: $\pm$ (1.5% of the nominal volume + from 2% of the merged volume) is equal to or greater than the nominal capacity: $\pm$ 5% of the nominal volume

COMPLIANCE WITH RELEVANT REGULATIONS AND TECHNICAL STANDARDS

Wilson Cook Medical, Inc., Cook Incorporated manufacture all medical products in accordance with the United States Food and Drug Administration (FDA) Quality System Regulations, 21 CFR Part 820, *Quality System Regulations*.

All Cook facilities also certified to *ISO 13485 (2003/2012) Quality Management System*.

WARRANTY

An authorized representative of the manufacturer in the territory of the Russian Federation is the ZAO "SHAG".

Russia, 119002, Moscow, Karmanitskiy per., 9.office 501a.

phone: (495) 956-13-09;

Fax: (495) 956-1310,

e-mail: info@schag.ru

Manufacturer's warranty binding on quality products.

INFORMATION ABOUT DEVICE PACKAGING

Percutaneous Endoscopic Gastrostomy Set (PEG-24-PUSH-I-S, PEG-24-PULL-I-S, FLOW-20-PUSH-I-S, and FLOW-20-PULL-I-S) are supplied in sterile packages. FLOW-I and PEG-I devices are individually packaged in a heat sealed Tyvek and PETG tray. The Wire Guides included with the FLOW-PUSH-I and PEG-PUSH-I devices are inserted into a wire guide holder prior to the PETG tray package. The Kit components are packaged in an Ethylene Vinyl Acetate Bag prior to the PETG tray package.

Universal Feeding Adapters (UFA-20-S, UFA-24-S) and Bolus Feeding Adapter (BFA-20-S) are packaged in a heat sealed Tyvek and Mylar pouch.

SHELF LIFE

The “use by” date as denoted by symbols on the labelling. The shelf life is shown by indication the time of manufacture (factory symbol) to the stated expiration date (hourglass symbol). The difference between both dates is the shelf life.

Percutaneous Endoscopic Gastrostomy Set (PEG-24-PUSH-I-S, PEG-24-PULL-I-S, FLOW-20-PUSH-I-S, FLOW-20-PULL-I-S) have a three (3) year shelf life.

Universal Feeding Adapters (UFA-20-S, UFA-24-S) and Bolus Feeding Adapter (BFA-20-S) have a three (3) year shelf life.

STERILIZATION DETAILS

All devices falling under the scope of this Technical Data Summary with the exception of accessories -Universal Feeding Adapters (UFA-20-S, UFA-24-S) and Bolus Feeding Adapter (BFA-20-S) - are supplied sterile in peel open packages utilizing an EtO (Ethylene Oxide) validated method to a sterility assurance level (SAL) of 10^{-6} . Do not use the sterile product if there is doubt as to whether the product is sterile.

Universal Feeding Adapters (UFA-20-S, UFA-24-S) and Bolus Feeding Adapter (BFA-20-S) are non-sterile and disposable.

STORAGE RULES

These devices should be stored in a dry location, away from temperature extremes.

PARAMETERS OF PACKAGED DEVICE – tolerance + - 10%

RPN	Height, cm	Width, cm	Weight, Gr
PEG-24-PULL-I-S	35.5	20.3	427.4
PEG-24-PUSH-I-S	35.5	20.3	430.4
FLOW-20-PULL-I-S	35.5	20.4	412.5
FLOW-20-PUSH-I-S	35.5	20.2	430.1
UFA-24-S	22.7	12.7	14.2
UFA-20-S	22.7	12.6	14.1


BFA-20-S	22.7	12.6	10.2
----------	------	------	------

LABELING


Enteral Tube Feeding Sets – Sterile are supplied in the individual packaging, which should show no signs of damage or opening. Upon detecting the bends, breaks and other defects the device cannot be used.

The label of the device above shows the following device parameters:


- company-manufacturer and manufacturing address;
- device name (trade mark);
- serial number;
- marking about CE certification;
- marking about sterilization;
- address of the EU representative;
- trade mark of the manufacturer;
- manufacturing date;
- sizes;
- schematic picture;
- lot number;
- marking for “Single use only”
- marking for “Attention. See Instructions for Use” (except Universal Feeding Adapters and Bolus Feeding Adapters)
- marking for consult Instructions for Use (only for Entuit™ Start Initial Placement Gastrostomy Sets)
- shelf life


FLOW-20-PULL-I-S
PERCUTANEOUS ENDOSCOPIC
GASTROSTOMY SET
G22761

7/2/14

PEG-24-PUSH-I-S
PERCUTANEOUS ENDOSCOPIC
GASTROSTOMY SET
G22758

These sets are single use only, do not attempt to repair, carry out maintenance or reuse.

DISPOSAL RULES

Upon completion of procedure, the devices are to be disposed of as per institutional guidelines for biohazardous medical waste.

For EU countries -in accordance with Council Directive 93/42/EEC

For Russian Federation – in accordance with national requirements, regulated by sanitary rules and norms N2.1.7.2790-10 of December 12, 2010 "Sanitary-epidemiological requirements for the treatment of medical waste". Risk class of medical waste: class B.

Doris A. Hawks

Doris A. Hawks
Global Regulatory Affairs Specialist
Wilson-Cook Medical, Inc.
State of North Carolina, Forsyth COUNTY

Oct 27, 2017
Date

I, Tiffany A. Thomas a Notary Public certify that,

Doris A. Hawks personally came before me and acknowledged that she is Regulatory Affairs Specialist of Wilson-Cook Medical Inc. and acknowledged on behalf of the corporation, the due execution of the foregoing instrument.

Witness my hand and notarial seal this 27 day of October 20 17.

Tiffany A. Thomas

Notary Public

My Commission expires: 9-5-2022

Tiffany A Thomas
Notary Public, North Carolina
County of Guilford
My Commission Expires September 5, 2022

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